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

100467389 | Rev:5
Released: 13 Oct 2020
CO: 100563705
Release Level: 4. Production

2

—

100467389 | Rev:5
Released: 13 Oct 2020
CO: 100563705

[illegible][illegible][illegible][illegible][illegible][illegible]


Introducer (Figure 4) / Gaine d'introduction (figure 4) / Einführsystem (Abbildung 4) / Introduuttore (Figura 4) / Introducitor (Figura 4) / Introduktor (afbeelding 4) / Indfører (figur 4) / Sisäänviejä (kuva 4) / Introducer (figur 4) / Innfører (figur 4) / Εισαγωγικός (Εικόνα 4) / Zaváděč (Obrázek 4) / Bevezető (4. ábra) / Introduktor (Rycina 4) / Zavádzač (obrázok 4) / Интродьюсер (Фигура 4) / Dispositif de introduction (Figura 4) / Sisestusseade (Joonis 4) / Levaditaja (4. attēls) / Stömikilis (4 pav.) / Introduuser (Sekil 4) / Интродьюсер (рисунок 4) / Интродьюсер (малюнок 4) / Интродьюсер (4-сурет) / 導入器 (図 4) / 引導套組 (圖 4) / 삽입기 (그림 4) / 導入器 (圖 4) / Introducer (Gambar 4) / Uvodač (slika 4.) / Uvodnik (slika 4) / Проводник (слика 4) / Uvajalo (Slika 4) / 4. (شك 4) / 4. (شك 4)

Vessel Dilator (Figure 5) / Dilatateur vasculaire (figure 5) / Gefäßdilator (Abbildung 5) / Dilatore per vasi sanguigni (Figura 5) / Dilatador vascular (Figura 5) / Dilatador de vaso (Figura 5) / Bloedvatdilatator (afbeelding 5) / Kardilator (figur 5) / Suonen laajennin (kuva 5) / Kärldilator (figur 5) / Kardilator (figur 5) / Διαστολικός αγγείο (Εικόνα 5) / Cévné dilatátor (Obrázek 5) / Értágító (5. ábra) / Rozszerzacz naczyniowy (Rycina 5) / Cievny dilatátor (obrázok 5) / Судовий дилататор (Фигура 5) / Dilatator vascular (Figura 5) / Veresoonne laiendaja (Joonis 5) / Asinsvadu dilatators (5. attēls) / Kraujagysdės plėtiklis (5 pav.) / Damar Dilatoru (Sekil 5) / Сосудистый дилататор (рисунок 5) / Судинний дилататор (малюнок 5) / Тамыр кеңейткіші (5-сурет) / 血管扩张器 (图 5) / 血管擴張器 (圖 5) / 혈관 확장기 (그림 5) / 血管擴張器 (圖 5) / Dilator Pembuluh (Gambar 5) / Dilator žila (slika 5.) / Dilator krvnog suda (slika 5) / Дилататорот на крвни садови (слика 5) / Dilator žile (Slika 5) / 5. (شك 5) / 5. (شك 5)

Mini-Guidewire (Figure 6) / Mini-guide (figure 6) / Mini-Führungsdraht (Abbildung 6) / Mini filo guida (Figura 6) / Miniguia (Figura 6) / Mini fio-guida (Figura 6) / Mini-voerdraad (afbeelding 6) / Mini-guidewire (figur 6) / Miniohjainvaimeri (kuva 6) / Minileidare (figur 6) / Mini-leidesonde (figur 6) / Μικροεγώγιμος οδηγός (Εικόνα 6) / Mini vodíci drát (Obrázek 6) / Mini vezetődrót (6. ábra) / Miniprowadnik (Rycina 6) / Vodiac minidrót (obrázok 6) / Минипроводчик (Фигура 6) / Mini-fir de ghidare (Figura 6) / Minijuhtedraht (Joonis 6) / Mini vadstiga (6. attēls) / Mažoji kreipiamoji viela (6 pav.) / Mini Kılavuz tel (Sekil 6) / Минипроводчик (рисунок 6) / Mini-provodičnik (малюнок 6) / Шабарын бағыттағыш сым (6-сурет) / 微导线 (图 6) / 微導絲 (圖 6) / 소형 가이드 와이어 (그림 6) / 小形ガイドワイヤ (圖 6) / Kawat Pemandu Mini (Gambar 6) / Mini žica vodilica (slika 6.) / Mini vodíci-žica (slika 6) / Мини жица-водилка (слика 6) / Mini vodilna žica (Slika 6) / 6. (شك 6) / 6. (شك 6)

Needle (Figure 7) / Aiguille (figure 7) / Nadel (Abbildung 7) / Ago (Figura 7) / Aguja (Figura 7) / Agulha (Figura 7) / Naald (afbeelding 7) / Nål (figur 7) / Neula (kuva 7) / Nål (figur 7) / Nål (figur 7) / Βελόνα (Εικόνα 7) / Jehla (Obrázek 7) / Tű (7. ábra) / Iгла (Rycina 7) / Ihla (obrázok 7) / Игла (Фигура 7) / Ac (Figura 7) / Nēl (Joonis 7) / Adats (7. attēls) / Adats (7 pav.) / Iğne (Sekil 7) / Игла (рисунок 7) / Игла (малюнок 7) / Игла (7-сурет) / 針 (图 7) / 針頭 (圖 7) / 바늘 (그림 7) / 針 (圖 7) / Jarum (Gambar 7) / Iгла (slika 7.) / Iгла (slika 7) / Игла (слика 7) / Iгла (Slika 7) / 7. (شك 7) / 7. (شك 7)

English

STERILE. Sterilized with ethylene oxide gas. Nonpyrogenic. For single use only. Do not re-sterilize.
Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

I. System Components/Description

Catheter Sheath Introducer (CSI, Included)

A catheter sheath introducer (CSI) facilitates percutaneous entry of an intravascular device. The sheath surface has a hydrophilic coating to enhance entry and withdrawal during vessel access. Color coding indicates the French size of the largest intravascular device that will pass through the CSI. Refer to Figure 4.

French	Color
4	Red
5	Gray
6	Green
7	Orange

Vessel Dilator (Included)

A vessel dilator facilitates the percutaneous entry of the CSI by forming an atraumatic transition from the skin through the subcutaneous tissue to the vessel. Refer to Figure 5.

Mini-Guidewire (Included)

A 0.021" diameter, mini-guidewire (Refer to Figure 6) is provided for establishing vessel access. The mini-guidewire is bare (stainless steel), nitinol, polymer, or hydrophilic coated depending on the kit configuration. Note that the hydrophilic coated, polymer, and nitinol wires contain nickel, and should not be used in patients in which this may cause an allergic response.

IV Catheter Needle (Included)

An IV Catheter access needle (Refer to Figure 7) is provided with the hydrophilic wire and polymer wire kit options to facilitate initial entry into the artery. Removal of the needle leaves the IV cannula in place to facilitate insertion of the mini-guidewire into the target vessel. Flush the needle per the Preparation section. Note that this needle contains chrome-nickel steel, and should not be used in patients in which this may cause an allergic response.

Bare Needle (Included)

A bare access needle (Refer to Figure 7) is provided for the bare wire and nitinol wire kit options. It also facilitates initial entry into the artery. After initial entry of the bare needle into the target vessel the mini-guidewire can be placed through the hub of the needle. Flush the needle per the Preparation section. Note: For catalog codes where the 0.021" mini-guidewire is hydrophilic or polymer, the bare needle is not included.

II. Indications

RAIN Sheath™ Transradial Introducer is indicated to facilitate placing a catheter through the skin into a radial artery.

III. Contraindications

- None Known

IV. Warnings

- Use of alcohol, antiseptic solutions, or other solvents should be avoided, as they may adversely affect the device.
- For the IV Catheter needle, do not reinsert the needle into the IV catheter at any time. The needle could damage the IV catheter, resulting in an IV catheter embolus.
- Do not leave the CSI in place for extended periods of time without a catheter in place.
- If using a hydrophilic or polymer wire, do not use with a bare needle, as this may damage the integrity of the coating or jacket.
- Manipulate the mini-guidewire slowly and carefully to avoid damage to the vessel wall, while monitoring the tip position and movement using standard catheterization technique.
- Once the vessel dilator is removed, manipulate the sheath introducer slowly and carefully to minimize the chances of kinking.
- Persons with allergic reactions to nickel may suffer an allergic response to components of this device.
- During the procedure, provide a proper anticoagulant or antiplatelet therapy to the patient.
- Do not use power injector for contrast media injection from the side port.
- Prior to radial access procedures, it is recommended to verify adequate collateral flow through the ulnar artery, such as with an Allen test. If collateral blood supply to the hand is considered inadequate, an alternate access site should be considered.
- Do not manually re-shape the tip of the mini-guidewire by applying external force intended to bend or affect the shape of mini-guidewire.

V. Precautions

- For single use only. Do not re-process or re-sterilize. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.
- Store in a dry, dark, cool place.
- Do not use if package is open or damaged.
- Use prior to the "use-by" date.
- Prior to use, confirm that the sheath size is appropriate for the access vessel to be used.
- The entire procedure, from skin stick to product removal, must be performed aseptically.
- Remove the product from the packaging tray with care to avoid damage to the sheath.
- If resistance is felt upon insertion or withdrawal, investigate the cause before continuing.
- Do not pull the side port with excessive force.
- This product is intended to be used by a trained professional using a standard catheterization technique.
- Use the product immediately after opening the package.
- Do not flush liquid such as contrast media or heparinized saline solution from side tube while a device such as a dilator or a catheter is inside of the sheath.
- Do not inject drugs including oil components such as lipid emulsion, castor oil, interfacial active agent or solubilization agent such as alcohol, through the side tube. It may cause cracks on the sheath.
- Be careful not to cut the side tube when holding it with forceps, or not to cut with scissors and knives.
- Do not pinch the plastic cannula and/or the side tube with forceps. It may cause scratches on it. Attention should be paid not to damage the plastic cannula with forceps or sharp edged tools.
- Do not scratch the sheath with needle point, cutting tool, or other edged tools.
- Do not incline a guide wire and / or a catheter while inserting them via the valve of the sheath.

VI. Complications

Possible complications include, but are not limited to:

- Abrupt vessel closure
- Additional intervention
- Allergic reaction (device, contrast medium and medications)
- Air embolism
- Infection
- Intimal tear
- Hematoma at puncture
- Hemorrhage
- Inflammation / Infection / Sepsis
- Ischemia
- Perforation of vessel wall
- Thrombus formation
- Peripheral nerve injury
- Pain
- Vascular complications (e.g. intimal tear, dissection, pseudo aneurysm, perforation, rupture, spasm, occlusion)

VII. Instructions for Use

a. Preparation

1. Select an appropriately sized CSI.
2. Carefully remove CSI from package using sterile technique.
CAUTION: To prevent damage while removing sheath from package, gently pull up the side port (at the sheath hub) to remove from the package.
3. Inspect the system contents for any signs of damage; do not use if any damage is present.
4. Soak the CSI in sterile heparinized saline or similar isotonic solution to activate the hydrophilic coating. If using a hydrophilic access wire, soak the wire as well to activate the coating.
5. Prior to use, remove air from CSI through the side port extensions by inverting and flushing with heparinized saline or suitable isotonic solution.
CAUTION: Failure to flush the device prior to use may result in air embolism.
6. Insert vessel dilator through the CSI, hemostasis valve, snapping it into place at the hub. Flush the vessel dilator with heparinized saline or another isotonic solution.
7. Use a syringe filled with heparinized saline or similar isotonic solution to flush the access needle (IV Catheter needle or bare needle) prior to use.

b. Recommended Procedure

1. Introduce the needle into the vessel using an aseptic technique. Holding the needle in place, insert the flexible end of the mini-guidewire through the needle and into the vessel. Gently advance the guidewire to the desired depth.
2. Do not withdraw the guidewire back through the metal cannula of the needle after it has been inserted as it may damage the guidewire.
3. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the CSI is inserted into the vasculature.
4. Thread the CSI with vessel dilator over the guidewire, grasping the sheath close to the skin to prevent buckling. Advance the assembly through the tissue and into the vessel.
CAUTION: To prevent damage to the CSI tip or kinking of the CSI body, do not withdraw the vessel dilator while advancing and positioning the CSI in the vessel.
5. If injection is required at this point in the procedure, remove the mini-guidewire and use the hub of the vessel dilator for injection.
6. Detach the vessel dilator from the CSI by releasing the snap fit ring at the hub. Slowly withdraw the mini-guidewire and vessel dilator.
7. Consider establishing a heparin drip via the side port extension.
8. Introduce the selected guide wire through the CSI.
9. Introduce the selected catheter through the CSI using one of the following methods:
 - a. Straighten the catheter tip by hand
 - b. Insert a guidewire through the catheter until the tip is straight.**NOTE: Hold the sheath in place when inserting, positioning, or removing the catheter.**
CAUTION: If there is a reduction in the valve hemostasis, insert the tip of the vessel dilator into the valve and withdraw it slowly.
10. To exchange catheters, slowly withdraw the catheter from the vessel and repeat steps 7 to 9.
11. Once all catheters and wires are removed, remove the sheath when clinically indicated by placing compression on the vessel above the puncture site, and slowly withdraw the CSI. Discard the sheath appropriately.
Important: Upon removing any intravascular device, aspirate via the side port extension to collect any fibrin that may have been deposited within or at the tip of the device.

VIII. DISCLAIMER OF WARRANTY AND LIMITATIONS OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE CORDS PRODUCT(S) DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL CORDS BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND CORDS TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

Descriptions or specifications in Cords printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Cord Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

18

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

44

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

4

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

■

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

**THIS PAGE HAS BEEN
INTENTIONALLY LEFT
BLANK**

PPE Specification
Labeling Specification
Spec for Merlin IFU

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Cardinal Corporation
64200 North West 62nd Avenue
Miami Lakes, Florida 33074 USA

Cardinal Castrol
Castrol Brand, Castrol
Castrol, Ireland

© 2020 Cardinal Health. All Rights Reserved. CDMHE, the Castrol LOGO, CARDINAL HEALTH and BIAH are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Patent www.cardinalhealth.com/patents

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6. RAIN Sheath Transradial Hydrophilic Thin-walled Introducer, diameter 2 mm, length 16 cm, consists of:

6.1. Transradial thin-walled introducer with hydrophilic coating RAIN Sheath Transradial, diameter and length variants: 2 mm (6F) x 16 cm - 1 pc.

6.2. Mini-Guidewire 0.021" (0.53 mm) - 1 pc., options:

6.2.1. Bare (Stainless Steel) Mini-Guidewire.

6.2.2. Nitinol Mini-Guidewires.

6.3. Vessel dilator, diameter and length variants: 2.13 mm x 21.60 cm - 1 pc.

6.4. Stainless steel Vygon Bare Needle 21G (08x40 mm) - 1 pc.

7. RAIN Sheath Transradial Hydrophilic Thin-walled Introducer, diameter 2.3 mm, length 10 cm, consists of:

7.1. Transradial thin-walled introducer with hydrophilic coating RAIN Sheath Transradial, diameter and length variants: 2.3 mm (7F) x 10 cm - 1 pc.

7.2. Mini-Guidewire 0.021" (0.53 mm) - 1 pc., variants:

7.2.1. Bare (Stainless Steel) Mini-Guidewire.

7.2.2. Nitinol Mini-Guidewires.

7.3. Vessel dilator, diameter and length variants: 2.44 mm x 15.80 cm - 1 pc.

7.4. Stainless steel Vygon Bare Needle 21G (08x40 mm) - 1 pc.

8. RAIN Sheath Transradial Hydrophilic Thin-walled Introducer, diameter 2.3 mm, length 16 cm, consists of:

8.1. Transradial thin-walled introducer with hydrophilic coating RAIN Sheath Transradial, diameter and length variants: 2.3 mm (7F) x 16 cm - 1 pc.

8.2. Mini-Guidewire 0.021" (0.53 mm) - 1 pc., variants:

8.2.1. Bare (Stainless Steel) Mini-Guidewire.

8.2.2. Nitinol Mini-Guidewires.

8.3. Vessel dilator, diameter and length variants: 2.44 mm x 21.80 cm - 1 pc.

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outweigh the risks associated in the appropriate patient population and as per its intended use, by a physician. Overall, the Cordis RAIN Sheath™ Transradial Introducer remains a safe and effective diagnostic and interventional tool for endovascular procedures.

9. GENERAL DESCRIPTION

Each RAIN Sheath™ device is supplied sterile and consists of a sheath introducer, a vessel dilator (0.021" guidewire compatible), an IV catheter needle or a bare access needle, and a 0.021" mini guidewire (either bare stainless steel or bare nitinol).

The RAIN Sheath™ facilitates placing a catheter through the skin into a radial artery.

Catheter Sheath Introducer

The catheter sheath introducer, also referred to as sheath cannula or sheath, has a hydrophilic coating, a smooth transition to the vessel dilator, and a low-profile diameter (OD) to facilitate entry to the vessel, minimize trauma and prevent spasm and occlusion of the radial artery. The hydrophilic coating facilitates handling of the product with reasonable force when attempting to insert and withdraw from the vasculature. The hydrophilic coating must be activated by soaking the sheath introducer in saline prior to use. The sheath has a hemostasis valve to minimize or prevent blood leakage during use.

The hub is over molded to the sheath and is color coded by French size.

French Size	Color
4F	Red
5F	Gray
6F	Green
7F	Orange

RAIN Sheath Transradial Hydrophilic Thin-walled Introducer

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Vessel Dilator

The RAIN Sheath™ vessel dilator facilitates the percutaneous entry of the catheter sheath introducer by forming a smooth atraumatic transition between the mini guidewire and the sheath introducer. The vessel dilator maintains the integrity of the sheath and dilates the blood vessel during insertion. The distal section of the RAIN Sheath™ dilator is lubricated with MDX, a lubricant, to enhance insertion of the dilator into the vasculature. The distal section of the dilator is tapered to create a smooth shoulder transition between the vessel dilator tip and the mini guide wire. The OD of the dilator is optimized to fit with the ID of the sheath cannula. The hub is over molded to the dilator which includes the locking mechanism between the hub of the vessel dilator and the hub of the cannula to prevent the vessel dilator from backing out of the sheath cannula while the device is inserted.

Stainless steel Vygon Bare Needle

The bare needle is a 21-gauge (21G) stainless steel access needle. After entry of the bare needle into the target vessel, the mini-guidewire can be placed through the hub of the needle.

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Mini-Guidewire

A 0.021" diameter mini-guidewire is provided in each RAIN Sheath™ device for maintaining vessel access. The mini-guidewire is either bare stainless steel (hereafter referred to as the bare wire) or bare nitinol (hereafter referred to as the nitinol wire), depending on the device configuration. The bare stainless steel wires are 45 cm long, and the bare nitinol wire is 43 cm long. These devices are used in the RAIN Sheath™ in a manner consistent with their intended use.

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10. TECHNICAL DESCRIPTION OF THE MEDICAL DEVICE

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The RAIN Sheath™ device has a traditional sheath with a bleed-back valve and a stopcock. The sheath has a hemostasis valve that minimizes or prevents blood leakage during use. Additionally, the introducer sheath has a hydrophilic coating which imparts lubricity and facilitates ease of insertion of the device. During radial access, the dilator taper and close fit of the dilator with the sheath, combined with the hydrophilic coating on the sheath, act to expand the arteriotomy and facilitate atraumatic entry into the radial artery.

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The product configurations of RAIN Sheath™ Transradial device are based on sheath French size, sheath length and the specific wire and needle provided with the device. RAIN Sheath™ is offered in French sizes of 4F, 5F, 6F and 7F, with lengths of 10 cm and 16 cm. The sheath hub and dilator are color coded by French size to align with accepted conventions for the guiding catheters to which the sheaths are compatible. The RAIN Sheath™ Transradial device is available in eight (8) product configurations with French sizes between 4F – 7F.

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Accessories

The following accessories are packaged with the RAIN Sheath™ device: a 0.021" mini-guidewire (either bare (stainless steel) or bare Nitinol) and a needle (Stainless steel Vygon Bare Needle 21G (08x40 mm)).

Scope of delivery:

1. Transradial thin-walled introducer with hydrophilic coating RAIN Sheath Transradial, one of the diameters and lengths in the set - 5 pc.
2. Instructions for use – 1 pc.

Weight and size characteristics of the MD

All the specified in the section values are nominal, the permissible error may be $\pm 10\%$, unless otherwise specified.

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11. MEDICINAL SUBSTANCE STATEMENT

The Cordis RAIN Sheath™ Transradial Introducer is not composed of substances or of combinations of substances that are absorbed by or locally dispersed in the human body.

12. WARNINGS

- Use of alcohol, antiseptic solutions, or other solvents should be avoided, as they may adversely affect the device.
- For the IV Catheter needle, do not reinsert the needle into the IV catheter at any time. The needle could damage the IV catheter, resulting in an IV catheter embolus.
- Do not leave the CSI in place for extended periods of time without a catheter in place.
- If using a hydrophilic or polymer wire, do not use with a bare needle, as this may damage the integrity of the coating or jacket.
- Manipulate the mini-guidewire slowly and carefully to avoid damage to the vessel wall, while monitoring the tip position and movement using standard catheterization technique.
- Once the vessel dilator is removed, manipulate the sheath introducer slowly and carefully to minimize the chances of kinking.

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- Persons with allergic reactions to nickel may suffer an allergic response to components of this device.
- During the procedure, provide a proper anticoagulant or antiplatelet therapy to the patient.
- Do not use power injector for contrast media injection from the side port.
- Prior to radial access procedures, it is recommended to verify adequate collateral flow through the ulnar artery, such as with an Allen test. If collateral blood supply to the hand is considered inadequate, an alternate access site should be considered.
- Do not manually re-shape the tip of the mini-guidewire by applying external force intended to bend or affect the shape of mini-guidewire.
- This product is designed and intended for single use. It is not designed to undergo reprocessing and resterilization after initial use. Reuse of this product, including after reprocessing and/or resterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- Reuse of this product, including after reprocessing and/or resterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- This device contains the following substance defined as CMR 1A and/or CMR 1B and/or endocrine disrupting substances in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0.
- Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

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13. PRECAUTIONS

- For single use only. Do not re-process or re-sterilize. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.
- Store in a dry, dark, cool place.
- Do not use if package is open or damaged.
- Use prior to the "use by" date.
- Prior to use, confirm that the sheath size is appropriate for the access vessel to be used.
- The entire procedure, from skin stick to product removal, must be performed aseptically.
- Remove the product from the packaging tray with care to avoid damage to the sheath.
- If resistance is felt upon insertion or withdrawal, investigate the cause before continuing.
- Do not pull the side port with excessive force.
- This product is intended to be used by a trained professional using a standard catheterization technique.
- Use the product immediately after opening the package.
- Do not flush liquid such as contrast media or heparinized saline solution from side tube while a device such as a dilator or a catheter is inside of the sheath.
- Do not inject drugs including oil components such as lipid emulsion, castor oil, interfacial active agent or solubilization agent such as alcohol, through the side tube. It may cause cracks on the stopcock.

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CAUTION: Failure to flush the device prior to use may result in air embolism.

6. Insert vessel dilator through the CSI's hemostasis valve, snapping it into place at the hub. Flush the vessel dilator with heparinized saline or another isotonic solution.

7. Use a syringe filled with heparinized saline or similar isotonic solution to flush the access needle (IV Catheter needle or bare needle) prior to use.

b. Recommended Procedure

1. Introduce the needle into the vessel using an aseptic technique. Holding the needle in place, insert the flexible end of the mini guidewire through the needle and into the vessel. Gently advance the guidewire to the desired depth.

2. Do not withdraw the guidewire back through the metal cannula of the needle after it has been inserted as it may damage the guidewire.

3. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the CSI is inserted into the vasculature.

4. Thread the CSI with vessel dilator over the guidewire, grasping the sheath close to the skin to prevent buckling. Advance the assembly through the tissue and into the vessel.

CAUTION: To prevent damage to the CSI tip or kinking of the CSI body, do not withdraw the vessel dilator while advancing and positioning the CSI in the vessel.

5. If injection is required at this point in the procedure, remove the mini-guidewire and use the hub of the vessel dilator for injection.

6. Detach the vessel dilator from the CSI by releasing the snap fit ring at the hub. Slowly withdraw the mini-guidewire and vessel dilator.

7. Consider establishing a heparin drip via the side port extension.

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PETG tray with sr n features to hold product components in place. Figure below provides a representation of the components assembled in the RAIN Sheath™ tray.

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Each tray is placed into a labeled Tyvek and Nylon/HDPE film chevron-style pouch. Five (5) pouches and one (1) IFU are inserted into a labeled carton

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and sealed with a closure label. The RAIN Sheath™ device is provided sterile. The RAIN Sheath™ Transradial packaging is represented as shown in Picture 5.

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Transport packaging is a corrugated cardboard box in which any required number of secondary (consumer) packaging is placed. The dimensions and weight of the transport packaging may vary depending on the number of products supplied.

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23. LABELING

The following marking symbols are applied to packaging:

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