Serial No.4791-03-22

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Hereby certify that I was presented with:

- ☒ the original document
☐ a true copy of the original document
☐ a copy that is not the original and is not a true copy of the original

☐ computerized information (please specify)
 from the database of

Drawn up in the language English
 in which I am fluent.

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I confirm, with my signature and seal, that the
 attached document marked with the letter/
 number A is an exact copy of the
 document or information presented to me.

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
A
Q

INSTRUCTION FOR USE		AlphaBio
Approved Name: <u>Asher Pizak</u> Position: <u>QA TRA Manager</u>		Signature, Stamp
	 <u>Shani Amich</u> Product Manager	

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INSTRUCTION FOR USE

Dental Implant Spiral with cover screw

Dental Implant Spiral with cover screw, in variations:

- | | | |
|-----|------|---------------------------|
| 1. | 1300 | Spiral IH D3.3mm L10.0mm |
| 2. | 1301 | Spiral IH D3.3mm L11.5mm |
| 3. | 1303 | Spiral IH D3.3mm L13.0mm |
| 4. | 1306 | Spiral IH D3.3mm L16.0mm |
| 5. | 1308 | Spiral IH D3.3mm L8.0mm |
| 6. | 1330 | Spiral IH D4.2mm L10.0mm |
| 7. | 1331 | Spiral IH D4.2mm L11.5mm |
| 8. | 1333 | Spiral IH D4.2mm L13.0mm |
| 9. | 1336 | Spiral IH D4.2mm L16.0mm |
| 10. | 1338 | Spiral IH D4.2mm L8.0mm |
| 11. | 1340 | Spiral IH D5.0mm L10.0mm |
| 12. | 1341 | Spiral IH D5.0mm L11.5mm |
| 13. | 1343 | Spiral IH D5.0mm L13.0mm |
| 14. | 1346 | Spiral IH D5.0mm L16.0mm |
| 15. | 1348 | Spiral IH D5.0mm L8.0mm |
| 16. | 1350 | Spiral IH D3.75mm L10.0mm |
| 17. | 1351 | Spiral IH D3.75mm L11.5mm |
| 18. | 1353 | Spiral IH D3.75mm L13.0mm |
| 19. | 1356 | Spiral IH D3.75mm L16.0mm |
| 20. | 1358 | Spiral IH D3.75mm L8.0mm |
| 21. | 1360 | Spiral IH D6.0mm L10.0mm |
| 22. | 1361 | Spiral IH D6.0mm L11.5mm |
| 23. | 1363 | Spiral IH D6.0mm L13.0mm |
| 24. | 1368 | Spiral IH D6.0mm L8.0mm |

1. DESCRIPTION

The Alpha-Bio Tec Spiral implant in variation include different types and sizes of dental implants made of Titanium alloy Ti 6Al 4V ELI and undergo a unique surface treatment. The implants are provided sterile by gamma irradiation.

The implantation is a surgical procedure which includes opening a flap or punching through the soft tissue, drilling and inserting the implant into the jawbone.

Restoration on the implants is done using abutments that enable fixed or removable restoration in various options: a single crown, a number of connected crowns, partial or full dentures.

Cover screw is supplied with each and every implant. The cover screw is used to seal the implant concealed under the soft tissue for the osteointegration period. Cover screw is blue colored for Spiral implant. Made of titanium alloy ELI 4V-6Al-Ti.

2. DESIGNATION

The product is intended for surgical implantation in the alveolar and / or basal bone of the lower jaw or upper jaw to install and fix denture retainer

3. APPLICATION AREA

Dentistry, maxillofacial surgery

4. INDICATIONS

Alpha-Bio Tec implants are intended for rehabilitating fully or partially edentulous patients. Dental implants are intended to be used in a manner that they will integrate the bone (osseointegration).

5. RECOMMENDATIONS

An essential condition for placing dental implants is the existence of adequate quantity and quality of bone. This data is gathered during the implantation planning stage by making appropriate radiographs (panoramic and computer tomography) of the implantation site. Anatomic areas near the implantation site such as: blood vessels, nerves, maxillary sinus and nasal cavity must be identified in order to prevent accidental damage. The performance of surgical procedures is subject to the patient's systemic condition.

After opening the package and removing the implant Alpha-Bio Tec it is recommended to place Alpha-Bio Tec implants using torque up to 50 Ncm.

In case of immediate loading it is recommended to place the implants using torque of at least 35Ncm.

6. CONTRAINDICATIONS

Patients who take steroid-based anticonvulsant and anticoagulant drugs. Patients receiving radiotherapy, chemotherapy or any other immunosuppressive treatment or who have been administered radiotherapy in the last 5 years. Metabolic bone disorders, uncontrolled bleeding disorders such as: hemophilia, thrombocytopenia, Degenerative diseases, osteoradionecrosis, renal failure, organ transplant recipients, HIV positive, malignant diseases, diseases that compromise the immune system, unbalanced diabetes mellitus, psychotic diseases, hypersensitivity to one of the components of the implant in general and titanium in particular, pregnancy, inability of the patient to maintain reasonable oral hygiene, lack of patient cooperation, use of alcohol, narcotics and uncontrolled endocrine diseases. Any systemic condition that is unbalanced and therefore precludes surgical procedures. Lack of appropriate experience or training of the attending dentist. Particular caution should be used while placing narrow (physical diameter less than 3.7mm) implants and internal connection narrow platform implants in the posterior region. These implants should not be overloaded as single restoration in posterior region.

Contraindications to implantation - Absolute contraindications.

1) General - severe general somatic diseases: diseases of the cardiovascular system in the stage of decompensation; diseases of the blood and blood-forming organs (lymphogranulomatosis, leukemia, hemolytic anemia); mental and behavioral disorders; immunopathological diseases and conditions (defects in the complement system with a pronounced decrease in the body's resistance, phagocytic disorders, humoral insufficiency syndromes, hypoplasia of the thymus and parathyroid glands); some diseases of the musculoskeletal system and connective tissue (rheumatic and rheumatoid processes, congenital osteopathies, bone dysplasias, conditions after radiation and drug therapy); diseases of the skeletal system and other pathological conditions that cause a violation of trophism and a weakening of the regenerative capacity of bone tissue (congenital osteopathies, bone dysplasias, conditions after radiation and drug therapy); diseases of the endocrine system, nutritional disorders and metabolic disorders (type I diabetes mellitus, dysfunction of the thyroid and parathyroid glands, diseases of the pituitary gland and adrenal glands); malignant neoplasms; tuberculosis; AIDS; sexually transmitted diseases; some skin diseases (dermatoses and scleroderma), regular use of narcotic drugs in the anamnesis; bisphosphonate treatment; as well as a number of diseases, provided that implantation is not authorized by an appropriate specialist: congenital defects and prosthetics of heart valves.

2) Local: some diseases of the oral mucosa (chronic recurrent aphthous stomatitis, lupus erythematosus, pemphigus, Sjogren's syndrome, Behcet's syndrome); generalized periodontitis, severe.

7. RELATIVE CONTRAINDICATIONS

Previously irradiated bone, treatment with anticoagulant drugs or bisphosphonates, bruxism, parafunctional habits, untreated and/or uncontrolled periodontal disease, temporomandibular joint disease, various pathologies of the oral mucosa.

Relative contraindications to implantation.

1) General: osteoporosis; low estrogen in women - for example, after ovariectomy; benign neoplasms; chronic infectious diseases; bad habits (alcohol and smoking abuse, drug addiction); age under 18; pregnancy and lactation.

2) Local: poor oral hygiene.

8. POST-OPERATION INSTRUCTIONS FOR PATIENTS

Avoid actions that subject the body to physical stress after implantation. Eat only cold, soft food. Take antibiotics, oral rinses and analgesics as instructed by the attending dentist. Do not spit or brush your teeth at the day of implantation. It is advisable to place an ice pack intermittently on the external area of the implantation on the day of the procedure.

9. POSSIBLE SHORT-TERM ADVERSE EFFECTS

Pain, swelling, infection, bleeding, hemorrhages.

10. POSSIBLE SHORT-TERM COMPLICATIONS AND RISKS

Nerve damage causing temporary or permanent loss of sensation of the lower lip and/or tongue, damage to blood vessels (may cause life threatening bleeding), perforation to nasal cavity and/or maxillary sinus and local and systemic infection.

11. POSSIBLE LONG-TERM COMPLICATIONS AND RISKS

Chronic pain, lack of sensation in the lip or tongue, local or systemic infection, perforation to the sinus, rupture of nasal cavity, damage to adjacent teeth, fracture of implant, fracture of restoration, fracture of bone, esthetic disorders and infectious endocarditis. Implants can fail (failure to integration).

12. CAUTION

Patients should be instructed to contact the attending dentist immediately in case of any complications or post-operation adverse effects.

The law restricts the sale of this device to, or on the order of, a license physician or dentist.

13. WARNINGS!

An error in patient evaluation, in pre-treatment diagnosis and in the treatment plan resulting from inexperience or lack of inappropriate training of the attending dentist may endanger the health of the patient and cause severe, irreversible and even life threatening damage.

Therefore, implantation must not be performed without appropriate training of the implanting dentist at a recognized institute!

- For Internal Hex Implants: The torque applied must not exceed 30N/cm when using the FID (Fast Implant Delivery) holder to insert the implant into the bone.

- The implant is designed for a single use only. Do not attempt to re sterilize or re-use due to the risk of cross infection.

- The implant is sterile unless its package is open or damaged!

- Performing implantations may require prophylactic antibiotics!

- Alpha-Bio Tec. devices are intended for trained practitioners only.

14. STERILITY

Implants and cover screws have been sterilized by gamma irradiation and are ready for use.

Open the implant and cover screw case onto a sterile surface.

All implants are for single use only and are for use prior to the labeled expiration date.

Do not use implants if the packaging has been damaged or previously opened.

The product is non-toxic, pyrogen-free.

15. OPERATING AND STORAGE CONDITIONS. SHELF LIFE. PATIENT CONTACT.

The Implant system parts should be kept in their original package, in a dry area at room temperature.

Storage

Temperature	20°±5°C
Atmospheric pressure	NA
Relative humidity	NA

Operation

Temperature	20 $\pm$ 5°C
Atmospheric pressure	NA
Relative humidity	NA

The shelf life (sterility) of implants is 5 years.

Permanent contact with the internal environment of the body at a temperature: +32°C - +42°C

The implantation period is unlimited, with adequate implantation and hygiene of the oral cavity, implants are installed for life.

16. WARRANTY OBLIGATION

In the case of failure of an implant, Alpha-Bio Tec will replace/provide another implant in exchange free of charge, in accordance with the following conditions:

- Filling in a report form supplied by Alpha-Bio Tec and attaching a radiogram taken before and after implantation.

- Submitting the report no later than 6 months from the beginning of the event, with the failed implant.

This is the full scope of the warranty for implantation provided by Alpha-Bio Tec that lists the sole remedies related to implantation.


For questions about the quality of a medical product, manufactured by Alpha-Bio Tec Ltd. please contact:


LLC «Nobel Biocare Russia»,

Russian Federation, 109004 Moscow, Stanislavskogo str., 21, bld. 2,

Telephone: +7 (495) 974-77-55 fax +7 (495) 974-77-86

17. SYMBOLS ON THE PACKAGE DESCRIPTION

SYMBOL	DESCRIPTION
	Sterile. Sterilized using irradiation
	Do not reuse
	Caution! refer to the instructions for use
	Consult instructions for use
	Do not use if package is damaged
	Use by
	Batch code
	Catalogue number

	Manufacturing date
	Manufacturer
	CE Mark
	The law restricts the sale of this device to, or on the order of, a li-censed physician or dentist.
	The symbol of Rostest

18. GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE

Number	REF	Name, Connection, Diameter, Length
1.	1300	Spiral IH D3.3mm L10.0mm
2.	1301	Spiral IH D3.3mm L11.5mm
3.	1303	Spiral IH D3.3mm L13.0mm
4.	1306	Spiral IH D3.3mm L16.0mm
5.	1308	Spiral IH D3.3mm L8.0mm
6.	1330	Spiral IH D4.2mm L10.0mm
7.	1331	Spiral IH D4.2mm L11.5mm
8.	1333	Spiral IH D4.2mm L13.0mm
9.	1336	Spiral IH D4.2mm L16.0mm
10.	1338	Spiral IH D4.2mm L8.0mm
11.	1340	Spiral IH D5.0mm L10.0mm
12.	1341	Spiral IH D5.0mm L11.5mm
13.	1343	Spiral IH D5.0mm L13.0mm
14.	1346	Spiral IH D5.0mm L16.0mm
15.	1348	Spiral IH D5.0mm L8.0mm
16.	1350	Spiral IH D3.75mm L10.0mm
17.	1351	Spiral IH D3.75mm L11.5mm
18.	1353	Spiral IH D3.75mm L13.0mm
19.	1356	Spiral IH D3.75mm L16.0mm
20.	1358	Spiral IH D3.75mm L8.0mm
21.	1360	Spiral IH D6.0mm L10.0mm
22.	1361	Spiral IH D6.0mm L11.5mm
23.	1363	Spiral IH D6.0mm L13.0mm
24.	1368	Spiral IH D6.0mm L8.0mm

Cover screw diameter 3.7 mm, length 4.9 mm

19. CONDITIONS FOR TRANSPORTATION

Transportation of products should be carried out by all types of covered transport in compliance with the conditions and requirements established for this type of transport.

Temperature	20°±5°C
Atmospheric pressure	NA
Relative humidity	NA

20. POTENTIAL USERS

Medical devices are intended for use by a dentist with experience in this procedure (dentists-restorers, surgeons).

Procedure conditions:

The product is intended for use in a medical institution (in dental clinics)

Patient Category:

Used for men and women over 18 years old.

21. PACKAGE. SCOPE OF DELIVERY

The implant in the package together with the cover screw is supplied in 10, 20 or 50 pcs. in a group cardboard box along with the instruction for use.

22. REQUIREMENTS FOR THE MAINTENANCE AND REPAIR OF MEDICAL DEVICES

Not applicable.

23. REQUIREMENTS FOR THE DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES

Medical devices after used are potentially hazardous medical waste and are class B medical waste — epidemiologically hazardous waste.

Disinfection and disposal of used implant should be carried out in accordance with the requirements of sanitary regulations of SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, living quarters, the operation of industrial, public premises, the organization and implementation of sanitary and anti-epidemic (preventive) measures"

Dispose in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

24. THE NAME AND LEGAL ADDRESS OF THE MANUFACTURER OF MEDICAL DEVICES AND AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RF:

Alpha Bio Tec Ltd., Israel

7 Hatnufa St., P.O.B 3936 Kiryat Arye, Petach Tikva 4951025 Israel

Manufacturing site

1. Alpha Bio Tec Ltd., 21 HaReches Blvd., POB 264, Modiin, 7171102 Israel

25. AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:

LLC «Nobel Biocare Russia»,

Russian Federation, 109004 Moscow, Stanislavskogo str., 21, bld. 2.

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