02/2021

Barcelona, March 7, 2023

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I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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
Laura López
Regulatory Affairs Manager
BIOSYSTEMS S.A.

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User Manual BioSystems BTS

Dear Customer,

Thank you for purchasing our semi-automatic analyser BTS. This semi-automatic analyser is one of the most technically advanced and easiest to operate in the market. We are sure that its features will make it a valued instrument for your laboratory. Although it can be operated in a logical and simple way through its menu options, we recommend you to read this manual carefully. It will help you for performing the installation and the maintenance correctly, and will allow you to get the maximum benefit from its many possibilities.

User Manual BioSystems BTS

1. INTRODUCTION

2. INSTALLATION

3. DESCRIPTION

4. GENERAL OPERATING

5. MEASUREMENT PROCEDURES

6. MAINTENANCE

7. TROUBLESHOOTING GUIDE

8. ACCESSORIES AND SPARE PARTS

9. TECHNICAL SPECIFICATIONS

10. END OF PRODUCT LIFE

11. MANUAL VERSION

1. INTRODUCTION

1. INTRODUCTION

1.1 About Semi-Automatic Analyzer

The semi-automatic analyzer is an in vitro diagnostics device designed to perform clinical biochemical and turbidimetric tests.

The semi-automatic analyser offers high reliability of its entire measurement system thanks to its optical system based on LED technology (340 to 670 nm) and great versatility in reading devices such as **flow cell, macro, micro, and semi-cuvettes and tubes**. All these Reading devices can be thermostatted quickly and accurately.

1. INTRODUCTION

1.1 About Semi-Automatic Analyzer

The functionality of this semi-automatic analyzer is completed by a series of automatic calculation modes that cover most of the possibilities of photometric analysis: **end point, differential mode, fixed time and kinetic, as well as non-linear techniques**. The calculations also allow the use of multi-point calibrators and calibration factor, the measurement of bichromatic test and QC management.

BTS' software allows the programming of up to 39 different parameters, which define a specific technique. The analyzer incorporates all the official techniques of BioSystems reagents - **more than 75 pre-programmed tests** - both for biochemistry and turbidimetry tests. In addition, the user can define new techniques unlimitedly.

1. INTRODUCTION

1.2. Intended Use

The semi-automatic analyzer is used to determine analyte concentrations using in vitro biochemistry, turbidimetry measurements on **human samples of serum, urine, plasma, cerebrospinal fluid, seminal plasma or whole blood**. The operation is based on the measurement of the concentration of the analyte present in the sample by measuring its absorbance after the use of one or more reagents.

The semi-automatic analyzer **is for professional use only**, for users with adequate training and capacity for its use. In conjunction with the installation of the instrument, it is recommended to read the user manual carefully and follow all the indications, warnings and precautions detailed therein.

2. INSTALLATION


2. INSTALLATION

2.1. Unpacking and Setting Up the Instrument

Unpack the instrument and place it together with its accessories on a flat work surface.

- a) Connect the instrument through its power adapter to the mains.
- b) Connect the waste bottle to the instrument through the fluid outlet fitting of the pump. To do this, use the supplied tube.
- c) **Install the flow cell in the instrument, as indicated in the image.**
- d) **Pass the aspiration tube (c) through the inside of the device until it exits in the aspiration zone.**
- e) **Then connect the waste outlet tube (d) from the cuvette to the plastic connector. You can start the analyzer and use it.**

NOTE: You can also refer to the Quick Guide sheet found in your packaging.

2. INSTALLATION

2.2. Packaging and Storage

If the semi-automatic analyzer must be stored for some time, it is advisable to clean the fluid circuit thoroughly. First, with washing solution and then with distilled water. The cannal must be drained of liquid once the cleaning is finished. **The semi-automatic analyzer must be protected from environmental aggressions, such as dust, excessive humidity, or direct sunlight.**

If it is necessary to re-send the semi-automatic analyzer or to move it, it is essential to use the original packaging.


3. DESCRIPTION

3. DESCRIPTION

3.1. Analyzer Components - Front

1. Display Color: Color Touch Screen to interact with user application.

2. Aspiration Tube: Teflon tube that aspirates the sample to flow cell.


3. Printer

4. Flow Cell Cover

3. DESCRIPTION

3.1. Analyzer Components - Back

5. ON-OFF Push Button:
Press 1 second to turn on
Press 5 seconds to turn off

6. Peristaltic Pump

7. Power-in Adapter

8. Fan Filter

9. USB and Ethernet Ports


10. External Waste Tube from fitting to pump

3. DESCRIPTION

3.1. Analyzer Components - Aspiration System

The aspiration system has the function of moving the liquids when we are using a flow cuvette.

3. DESCRIPTION

3.2. Analyzer Components - Tests with standard cuvettes or tubes

Standard disposable cuvettes can be used for absorbance readings. The cuvettes can be macro, semi-micro or micro, glass, quartz or plastic. The optical system of this semiautomatic analyzer is designed to be able to do readings with test tubes with a diameter of 12 mm. Its height should not exceed 75mm. To read with test tubes, remove the flow cuvette, place it in the housing intended for this purpose and insert the tube adapter into the cuvette holder.

The semi-automatic analyzer is set for 10 mm square cuvettes so that the absolute absorbance values read using test tubes will not match the readings of the same sample made with 10 mm cuvettes. For this reason, only tubes for absorbance measurements should be used in techniques programmed with CALIBRATOR.


4. GENERAL OPERATION

4. GENERAL OPERATION

4.1 Starting Up

When turning on, the semi-automatic analyser shows the main screen:

01. CONCENTRATION TEST: To measure concentrations in accordance with a pre-defined test, that has been programmed previously in the semi-automatic analyser.

02. ABSORBANCE TEST: To measure absorbance according to user preferences.

03. QUALITY CONTROL: To manage quality control and control chart.


04. TEST CONFIGURATION: To edit tests or add new tests into the analyser.

05. HISTORICAL DATA: To review the test results of samples previously performed.

4. GENERAL OPERATION

4.2.1 General Options

HOME: Main screen direct access from any screen and tests for all systems components

UTILITIES: To perform adjustments and testing system components.

SETTINGS: Instrument configuration options.

USER MANUAL: to consult the user manual

SHUTDOWN: This process turn-off analyser. The analyser drives a user to do end day maintenance.

QUICK RESTART: Fast option to restart BTS software

ALARMS: The pending alarms icon shows a list of warnings or advice to the user.

4. GENERAL OPERATION

4.2.2 General Buttons: Utilities & Settings

UTILITIES

To perform adjustments and testing system components.

- **Show IP:** Shows the instrument IP to facilitate communication through digital devices (if it is connected to an internet network).
- **Software Update:** Automatic procedure to update the BTS Software.
- **Pump Adjustment:** The screen allows users to check and adjust the movement of the pump to placing the sample in flow cell.
- **Temperature Adjustment:** The screen allows users to check and adjust the sensor temperature gap.
- **Light Reference Values:** The screen shows the current LED light values from the algorithm to adjust the light intensity and integration time for sample readings.

- **Functionality Test:** test to check the mechanical functionality of the different equipment systems (printer, fan, pumps).

- **Recover Default Values:** Procedure to recover the original factory values.

SETTINGS

Allows users to adapt some of the semi-automatic analyser aspects:

- **Setup clock/date:** Adjustment time and date.
- **Wash volume:** (microL) Volume aspirated when the user presses the wash button.
- **Patient code:** Enable or disable patient code. When enabled, the semi-automatic analyser requests the ID patient before each sample test.

4. GENERAL OPERATION

4.2.3 General Buttons: Action Buttons

Go to the next menu or screen.

Return to the previous menu or screen.

WASH BUTTON:

By pressing this button, the analyser aspirates a fixed volume of liquid into the aspiration system. It can be used to perform washing, cleaning and priming procedures.


To print tests results when printer is available.

4. GENERAL OPERATION

4.3 Main Menu Description

01. CONCENTRATION TEST

This option allows the calculation of concentrations by using any of the TEST previously programmed and stored in the application. On selecting this option, a screen will pop up with a list of the programmed tests.


When a test is selected, the instrument will request the input of a baseline to adjust the ZERO. **Always use distilled water to do the baseline.**

When the adjustment to ZERO is complete, the display shows information about reagent blank and standard values saved for the selected test.

4. GENERAL OPERATION

4.3 Main Menu Description

01. CONCENTRATION TEST

At this point, it is possible to do a new reagent blank, calibrations and quality control for the selected test.


Next, the program guides the user through the test procedure by means of different messages displayed on the screen.

When a reading is done, the display shows the value. If there is more than one replica, the display shows the average value. Abnormal values can be removed in case they appear.

4. GENERAL OPERATION

4.3 Main Menu Description

02. ABSORBANCE TEST

Option to perform direct absorbance readings.

4. GENERAL OPERATION

4.3 Main Menu Description

Absorbance		We 09/10 19:27	
Reading mode	Monochromatic	Stabilization time	1 s
Wavelength	340	Temperature	25 °C
Sample volume	300 uL		

READ BASELINE

Waiting for baseline

02. ABSORBANCE TEST


The work section is divided into two parts: On the top side, the user can select reading conditions. On the bottom, the user can see the absorbance readings.

Once the Reading parameters have been determined, the analyser requests a baseline, and subsequently, the user can take absorbance readings of the samples.

4. GENERAL OPERATION

4.3 Main Menu Description

03. QUALITY CONTROL

This menu option gives access to the quality control management program, based on the Levey-Jennings charts.

4. GENERAL OPERATION

4.3 Main Menu Description


03. QUALITY CONTROL

This option shows a list with the names of all programmed techniques. The user can select any test and check the evolution of all quality controls levels in control charts.

4. GENERAL OPERATION

4.3 Main Menu Description

04. TEST CONFIGURATION

This option allows the modification of any of the TESTS previously programmed and stored in the application. On selecting this option, a screen will pop up with a list of the programmed tests. When a test is selected, user can modify any parameter of pre-programmed test or create a new test.

IT IS RECOMMENDED NOT TO MODIFY ANY OF THE PRE-PROGRAMMED TESTS IN ORDER TO GET ACURATE AND RELIABLE RESULTS.

4. GENERAL OPERATION

4.3 Main Menu Description

≡ CALCIUM CPC - Concentration test Th 03/10 12:14

GENERAL CONTROLS CALIBRATION

Name
CALCIUM CPC

Technique Sample Units	Measure type	Reaction type
mg/dL	End point	Increasing

Reading mode
Monochromatic

Main wavelength
505

Test conditions

Aspiration volume Temperature Decimal positions

← SAVE

04. TEST CONFIGURATION

PARAMETERS: This tab shows the main characteristics of each test.

To obtain the most accurate and reliable results, BioSystems recommends to not modify the parameters specified in the official applications.

4. GENERAL OPERATION

4.3 Main Menu Description

04. TEST CONFIGURATION - GENERAL

NAME

Test name

TECHNIQUE SAMPLE UNITS

Units of test results

MEASURE TYPE

Users can choose between different types (see section 5 for more information about measure types):

- End point
- Kinetics
- Differential
- Fixed Time

REACTION TYPE

The software displays a window with the following two options: INCREASING / DECREASING. Indicates whether the absorbance of test increases or decreases with time.

READING MODE

The software displays a window with the following two options:

MONOCHROMATIC: A single measurement is made of the sample at the selected wavelength.

BICHROMATIC: Two readings are taken of the sample. One at one wavelength and the other at a reference wavelength. The absorbance value of the sample is the difference between both absorbances:

$$ABS = ABS \lambda_{main} - ABS \lambda_{reference}$$

MAIN/REFERENCE WAVELENGTH

Wavelengths in nanometres at which test to do the readings.

4. GENERAL OPERATION

4.3 Main Menu Description

04. TEST CONFIGURATION - GENERAL

ASPIRATION VOLUME: The volume in μL of sample, blank and/or standard to be aspirated for the measurements. The permitted margin is between 100 and 5000 μL .

Less the volume aspirated, more the carryover contamination between two consecutive samples will occur.

TEMPERATURE: Thermostat showing temperature of cuvette. Programmable between 25 and 40°C with a resolution of one degree.

INCUBATION TIME: Warming time before the first measurement.

STABILIZATION TIME

Time elapsed between the end of the aspiration cycle and the measurement in ENDPOINT and DIFFERENTIAL modes. It is usually five second, though other values can be programmed, allowing the measurement to be delayed until the liquid homogenises perfectly in the cuvette, especially when working with high viscosity.

DILUTION FACTOR

Enter a dilution factor for the sample. The result of the concentration will be multiplied by the factor input. The user has to dilute the sample manually by the factor entered in the programming.

INITIAL ABSORBANCE

For KINETIC and FIXED TIME modes only. For an increasing reaction, it is the maximum admissible absorbance value for the first reading. For a decreasing reaction, it is the minimum admissible absorbance value for the first reading.

4. GENERAL OPERATION

4.3 Main Menu Description

04. TEST CONFIGURATION - GENERAL

DETECTION LIMIT

Minimum detection limit of the test. If sample exceed this limit, program gives a warning message.

LINEARITY LIMIT

Maximum linearity limit. If it exceeds that limit, the program gives a warning message.

MINIMUM NORMALITY LIMIT

The minimum limit for normal values. That value is printed with the results and is shown on screen.

MAXIMUM NORMALITY LIMIT

The maximum limit for normal values. That value is printed with the results and is shown on screen.

BLANK LIMIT

The absorbance limit value that a blank can take. If it exceeds that limit, the program gives a warning message. It is possible to enable or disable to do blank test. It is mandatory for END POINT and DIFFERENTIAL measure types.

FACTOR (A)

The multiplicative value applied to the concentration results.

CONSTANT (B)

The additive value applied to the concentration results.

4. GENERAL OPERATION

4.3 Main Menu Description

☰ CALCIUM CPC - Concentration test Th 03/10 19:56

GENERAL CONTROLS CALIBRATION

Type
STANDARD

Available calibrators

Calibrator	Status	Concentration	LOT
Calibrator 1	ENABLED	300 mg/dL	a1222
Calibrator 2	ENABLED	600 mg/dL	uw822
Calibrator 3	ENABLED	900 mg/dL	21289s

← SAVE

04. TEST CONFIGURATION - CONTROLS

CONTROLS: Users can enable or disable QC levels (up to 3 levels) and edit QC values (QC range and lot n°).

4. GENERAL OPERATION

4.3 Main Menu Description

CALCIUM CPC - Concentration test Th 03/10 12:19

GENERAL CONTROLS CALIBRATION

Quality control concentrations

Level 1
ENABLED
LOT 23dj892
Lower limit 73,6 mg/dL
Upper limit 99,6 mg/dL

Level 2
ENABLED
LOT 23dj892
Lower limit 218 mg/dL
Upper limit 296 mg/dL

Level 3

SAVE

04. TEST CONFIGURATION - CONTROLS

CALIBRATION: Users can select the calibration type; by factor or standards.


1) Calibrating by a factor, users must enter a calibration factor.

2) Calibrating by standard, users can edit the number of calibration points (up to five points), lot number and concentration values.

4. GENERAL OPERATION

4.3 Main Menu Description

05. HISTORICAL DATA

It shows a history of all operations performed with the semi-automatic analyser. It allows visualising the stored values of the results of the concentration of the patients.

The information shown can be sorted by:

- Date
- Sample ID
- Test results
- Measure type (Blank, Calibration, QC or Samples)

5. MEASUREMENT PROCEDURES

5. MEASUREMENT PROCEDURES & CALCULATIONS

This chapter reviews the different measuring procedures based on the analysis mode and the calculations that are made to obtain the concentration values. The different formulas used are shown. Measurement procedures are the same for blanks, standards, controls and samples. Controls are treated like samples for every calculation.

A Blank	Blank absorbance
A Standard	Standard absorbance
A Sample	Absorbance of the sample
[...] λ_{main}	Absorbance at main wavelength
[...] $\lambda_{reference}$	Absorbance at reference wavelength
[...] R1+R2	Absorbance of the blank, standard or sample with the two reagents (global reaction)
[...] R1	Absorbance of the blank, standard or sample with the first reagent (reaction of the sample blank)
[...] T1	Absorbance of the blank, standard or sample measured after incubation 1
[...] T2	Absorbance of the blank, standard or sample measured after incubation 2
A FRAC	Absorbance of the fraction
A TOTAL	Total absorbance
$\Delta A/min$	Velocity of change of absorbance of the standard
C Sampe	Concentration of the sample
C Standard	Programmed of oncentration of the standard
n	Number of replicates
f	Standard factor
TR	Fixed factor, depending on the type of reaction. Its value is (+1) for increasing reactions and (-1) for decreasing reactions.

5. MEASUREMENT PROCEDURES & CALCULATIONS

5.1 End Point – Definition

ABSORBANCE

The absorbance of the reaction is measured just once against the distilled water baseline. In this procedure, one or two reagents can be used and the absorbance can be measured at one or two wavelengths. Calibration is based on the use of standards (one or several) or a programmed factor. A blank is prepared for each test with reagent alone or the reagent with water. The absorbance of the blank is also measured against the distilled water baseline.

$$A_{\text{Sample}} = [A_{\text{Sample}}] - [A_{\text{Blank}}]$$

$$A_{\text{Standard}} = [A_{\text{Standard}}] - [A_{\text{Blank}}]$$

$$A_{\text{QC}} = [A_{\text{QC}}] - [A_{\text{Blank}}]$$

MONOCHROMATIC / BICHROMATIC

The absorbance can be measured at one or two wavelengths. In the case of bichromatic readings, the absorbance value is taken as the difference between the absorbance at the main wavelength and the absorbance at the reference wavelength.

$$A_{\text{Sample}} = [A_{\text{Sample}}]_{\lambda_{\text{principal}}} - [A_{\text{Sample}}]_{\lambda_{\text{reference}}}$$

$$A_{\text{Standard}} = [A_{\text{Standard}}]_{\lambda_{\text{principal}}} - [A_{\text{Standard}}]_{\lambda_{\text{reference}}}$$

$$A_{\text{Blank}} = [A_{\text{Blank}}]_{\lambda_{\text{principal}}} - [A_{\text{Blank}}]_{\lambda_{\text{reference}}}$$

5. MEASUREMENT PROCEDURES & CALCULATIONS

5.1 End Point – Concentration

Based on the value of the absorbance obtained, the concentration of the analyte in the sample can be calculated:

FACTOR

The concentration is calculated using a programmed factor with the formula:

$$C_{\text{Sample}} = TR \cdot F \cdot (A_{\text{Sample}} - A_{\text{Blank}})$$

SINGLE STANDARD

If a standard is used at a single point, the concentration is calculated with the formula:

$$C_{\text{Sample}} = \frac{[A_{\text{Sample}}] - [A_{\text{Blank}}]}{[A_{\text{Standard}}] - [A_{\text{Blank}}]} \cdot C_{\text{Standard}}$$

MULTI-STANDARD

If a multi-point standard is used, the concentration is calculated using a standard curve or function. This curve is obtained from the programmed concentration values of the standards and the absorbance values measured for each one with respect to the baseline, with the formula:

$$C_{\text{Sample}} = \text{Func} [A_{\text{Sample}}]$$

5. MEASUREMENT PROCEDURES & CALCULATIONS

5.2 Differential bi-reagent – Definition

ABSORBANCE

For each test, a blank is prepared using distilled water instead of the sample or just reagents. The absorbances of this blank, with the first reagent (R1) and with both reagents (R2), are also measured against the distilled water baseline. Calibration is based on the use of standards (one or several) or a programmed factor. The absorbance value is the difference between the absorbance measured with both reagents and the absorbance measured with just the first reagent:

$$A_{\text{Sample}} = [A_{\text{Sample}}]_{R1+R2} - [A_{\text{Sample}}]_{R1}$$

$$A_{\text{Standard}} = [A_{\text{Standard}}]_{R1+R2} - [A_{\text{Standard}}]_{R1}$$

$$A_{\text{Blank}} = [A_{\text{Blank}}]_{R1+R2} - [A_{\text{Blank}}]_{R1}$$

MONOCHROMATIC / BICHROMATIC

The absorbance can be measured at one or two wavelengths. In the case of bichromatic readings, the absorbance value is taken as the difference between the absorbance at the main wavelength and the absorbance at the reference wavelength.

$$A_{\text{Sample}} = [A_{\text{Sample}}]_{\lambda_{\text{principal}}} - [A_{\text{Sample}}]_{\lambda_{\text{reference}}}$$

$$A_{\text{Standard}} = [A_{\text{Standard}}]_{\lambda_{\text{principal}}} - [A_{\text{Standard}}]_{\lambda_{\text{reference}}}$$

$$A_{\text{Blank}} = [A_{\text{Blank}}]_{\lambda_{\text{principal}}} - [A_{\text{Blank}}]_{\lambda_{\text{reference}}}$$

5. MEASUREMENT PROCEDURES & CALCULATIONS

5.2 Differential bi-reagent – Concentration

Based on the value of the absorbance obtained, the concentration of the analyte in the sample can be calculated.

FACTOR

The concentration is calculated using a programmed factor with the formula:

$$C_{\text{Sample}} = TR \cdot F \cdot (A_{\text{Sample}} - A_{\text{Blank}})$$

SINGLE STANDARD

If a standard is used at a single point, the concentration is calculated with the formula:

$$C_{\text{Sample}} = \frac{[A_{\text{Sample}}] - [A_{\text{Blank}}]}{[A_{\text{Standard}}] - [A_{\text{Blank}}]} \cdot C_{\text{Standard}}$$

MULTI-STANDARD

If a multi-point standard is used, the concentration is calculated using a standard curve or function. This curve is obtained from the programmed concentration values of the standards and the absorbance values measured for each one with respect to the baseline, with the formula:

$$C_{\text{Sample}} = \text{Func} [A_{\text{Sample}}]$$

5. MEASUREMENT PROCEDURES & CALCULATIONS

5.3 Fixed Time – Definition

ABSORBANCE

The absorbance of the reaction is read in two specific times with respect to the distilled water baseline at just one wavelength. Calibration is based on the use of standards (one or several) or a programmed factor. The absorbance value is the difference between the absorbance measured at time interval T1 and the absorbance measured at time interval T2.

$$A_{\text{Sample}} = [A_{\text{Sample}}] T2 - [A_{\text{Sample}}] T1$$

$$A_{\text{Standard}} = [A_{\text{Standard}}] T2 - [A_{\text{Standard}}] T1$$

$$A_{\text{Blank}} = [A_{\text{Blank}}] T2 - [A_{\text{Blank}}] T1$$

5. MEASUREMENT PROCEDURES & CALCULATIONS

5.3 Fixed Time – Concentration

Based on the value of the absorbance obtained, the concentration of the analyte in the sample can be calculated.

FACTOR

The concentration is calculated using a programmed factor with the formula:

$$C_{\text{Sample}} = TR \cdot F \cdot (A_{\text{Sample}} - A_{\text{Blank}})$$

SINGLE STANDARD

If a standard is used at a single point, the concentration is calculated with the formula:

$$C_{\text{Sample}} = \frac{[A_{\text{Sample}}] - [A_{\text{Blank}}]}{[A_{\text{Standard}}] - [A_{\text{Blank}}]} \cdot C_{\text{Standard}}$$

MULTI-STANDARD

If a multi-point standard is used, the concentration is calculated using a standard curve or function. This curve is obtained from the programmed concentration values of the standards and the absorbance values measured for each one with respect to the baseline, with the formula:

$$C_{\text{Sample}} = \text{Func} [A_{\text{Sample}}]$$

5. MEASUREMENT PROCEDURES & CALCULATIONS

5.4 Kinetics – Definition

Variation in the **ABSORBANCE** by time unit

The kinetic mode is used to measure the catalytic activity of an enzyme. The absorbance of the reaction against the distilled water baseline is measured regularly during several cycles, between time intervals **T_i** and **T_f** as programmed in the test. The readings are taken at only one wavelength. Based on these absorbance readings, the analyser calculates the variation in the absorbance of the reaction by time unit. Calibration is based on the use of standards (one or several) or a programmed factor.

LINEARITY CHECK

The catalytic activity is measured by the velocity of the reaction, which is proportional to the slope of the absorbance curve as opposed to the time. This slope is calculated using the linear regression method for all the absorbances measured between time intervals **T_i** and **T_f**. The most common units for measuring the slope are $\Delta A/\text{min}$. The semiautomatic analyser reads at regular intervals and automatically calculates the best straight line regression using the least squares method.

5. MEASUREMENT PROCEDURES & CALCULATIONS

5.4 Kinetics – Concentration

FACTOR

The concentration is calculated using a programmed factor with the formula:

$$C_{\text{sample}} = \text{TR. F.} \left(\frac{\Delta A}{\Delta t} \Big|_{\text{sample}} - \frac{\Delta A}{\Delta t} \Big|_{\text{blank}} \right)$$

SINGLE STANDARD

If a standard is used at a single point, the concentration is calculated with the formula:

$$C_{\text{sample}} = \text{TR. F.} \frac{\frac{\Delta A}{\Delta t} \Big|_{\text{sample}}}{\frac{\Delta A}{\Delta t} \Big|_{\text{standard}} - \frac{\Delta A}{\Delta t} \Big|_{\text{blank}}} \cdot C_{\text{standard}}$$

6. MAINTENANCE

6. MAINTENANCE

6.1 General Maintenance Rules

The semi-automatic analyser is designed for trouble free operation and requires minimum maintenance. To ensure reliable capability to prolong the instrument's life, perform periodical maintenance for the instrument according to the user manual.

Do not perform any maintenance other than what is specified in this chapter; otherwise, human injury and instrument damage may occur. The unauthorized maintenance for instrument may cause system damage and human injury, and invalidate the terms of maintenance service stated in contract. After maintenance, make sure the systems are in working order.

CAUTION

6. MAINTENANCE

6.1 General Maintenance Rules

- Never use detergents or abrasive products to clean the outside of the semi-automatic analyser. Use only a damp cloth dampened in water and neutral soap.
- If a reagent or corrosive product spills onto the semi-automatic analyser, wipe it immediately with a damp cloth in water.
- **Always leave the flow cell filled with distilled water or air. Never leave any reaction mixture in it.**
- The cuvette holder tray is equipped with watertight joints in order to prevent penetration of liquid into the inside of the instrument. If liquid is spilt onto the tray, wipe it with a damp cloth.

BIOLOGICAL DANGEROUS

Be sure to put on protective clothes, gloves, and even goggles when necessary before maintenance.

6. MAINTENANCE

6.2 Maintenance Schedule

FREQUENCY	ACTION	WHAT TO DO
At the beginning of working day	Wash	Aspirate washing solution
	Rinse	Aspirate distilled water
Between methods	Rinsing	Aspirate distilled water only
At the end of working day	Wash	Aspirate washing solution
	Rinse	Aspirate distilled water + air
	Clean	Clean outside the analyzer
Every week	Wash	Aspirate washing solution
	Rinse	Aspirate abundant distilled water + air
Every month	Verify	Pump calibration
	Verify	Flow cell external cleanness, and clean if needed
	Clean	Filter fan

6. MAINTENANCE

6.3 Daily Maintenance

It is recommended to do the initial start up and shut-down steps to optimize the functionality of the analyzer.

At the beginning of the day, once the analyzer has been switched on, it is recommended to aspirate 5 ml of washing solution and then plenty of distilled water to prepare the aspiration system.

On completing a series of measurements or after a long stand-by, it is recommended to rinse aspiration system with washing solution and abundant distilled water in order to clean and remove any bubble.

At the end of the working day, wash thoroughly with a washing solution. Lastly, rinse with distilled water and empty the aspiration systems by performing wash cycles with air only.

6. MAINTENANCE

6.4 Preventive Maintenance

It is recommended to do a weekly preventive maintenance to optimize the functionality of the analyzer.

CLEANING ASPIRATION SYSTEM

It is necessary to clean the aspiration system properly after each series of measurements and at the end of the day. Wash the aspiration system thoroughly with washing solution. Keep washing solution inside aspiration system for 5 minutes. Lastly, rinse with abundant distilled water and empty the circuit by performing wash cycles with air only.

If the outer end of the aspiration tube is damaged, replace it. Always use original spares.

PUMP calibration

Verify if peristaltic pump is calibrated by using UTILITY function CALIBRATE PUMP. In case peristaltic pump is not well calibrated, proceed to adjust by using UTILITY function CALIBRATE PUMP.

ADVICE: Change the tube of the peristaltic pump once a year

6. MAINTENANCE

6.4 Preventive Maintenance

It is recommended to do a weekly preventive maintenance to optimize the functionality of the analyzer.

CLEANING FLOW CELL

Check for cleanness of flow cell. If flow cell is dirty or has splashes on it, proceed to clean as described.

- a) To clean the inside, proceed as described in CLEANING ASPIRATION SYSTEM.
- b) To clean the outside, first use a paper towel dampened with distilled water and then use a paper towel dampened with alcohol (70%) to remove any dust or splashes and let it air dry.

ADVICE: It is recommended not to touch flow cell optical sides to avoid dirtying the optical parts of the cell. It could affect optical readings.

CLEANING FAN FILTER

- 1) The filter is on the back of the semi-automatic analyzer. Remove the plastic cover by pulling lightly to access the filter.
- 2) Change or clean the dust filter of the fan. To clean the filter, wash it with water and let it dry before replacing it. Replace the filter in its housing.
- 3) Put the plastic cover back in position.

6. MAINTENANCE

6.5 Disposable Waste

The semi-automatic analyzer has a waste bottle where the remains are stored. To carry out the waste disposal correctly, should follow good laboratory practices in accordance with the local or national legislation as per the country.

BIOLOGICAL DANGER

The waste contained in the waste bottle carries biological risk for the operator. Wear gloves, gown, and lab goggles to handle the waste bottle.

7. TROUBLESHOOTING GUIDE

7. TROUBLESHOOTING GUIDE

Troubleshooting directions are designed to help the user search and eliminate the malfunctions in the operation, and also offer ways to obtain timely and convenient technical help from BioSystems. For correct and quick troubleshooting, read through the user manual carefully and be familiar with the operation and maintenance measures.

GENERAL PROBLEM	PROBLEM DESCRIPTION	LIKELY CAUSE or CORRECTIVE ACTION
ASPIRATION SYSTEM	Liquid is not aspirated when the Dynamic key <WASH> switch is pressed	• Check if peristaltic pump tubing is connected properly, if not reconnects correctly. • Perform pump calibration again. • Replace tubings • Call Technical Service
	Peristaltic Pump operates but liquid is not aspirated.	• Ensure that tubings are properly connected. • Check tubing for blockage, crimps or leaks. • Perform pump calibration again. • Replace tubings. • Call Service Engineer.

7. TROUBLESHOOTING GUIDE

GENERAL PROBLEM	PROBLEM DESCRIPTION	LIKELY CAUSE or CORRECTIVE ACTION
CONTAMINATION PROCESS	Contamination between samples	It is not necessary to wash the instrument between different samples for the same method. However: a) In case of high risk of cross-over, you could pre-rinse the flow cell by aspirating distilled water, followed by air (to avoid dilution) between the samples. b) In case of a strongly colored reagent, it might be necessary to pre-rinse the flow cell with (working) reagent before aspirating the first measurement. Otherwise the first measurement could be influenced by the distilled water used to zero the instrument.
	Contamination between methods	Always rinse the flow cell with distilled water to avoid contamination problems.

7. TROUBLESHOOTING GUIDE

GENERAL PROBLEM	PROBLEM DESCRIPTION	LIKELY CAUSE or CORRECTIVE ACTION
TEST RESULTS	Poor reproducibility and erratic readings.	• Check if the aspiration volume is programmed and whether the amount of liquid aspirated is as programmed. • If the aspiration is not proper, calibrate pump. • Check expiry date of the reagent and control material used. • Check if all the necessary precautions have been taken. Always use clean glassware. • Check flow cell for air bubbles. • Call service engineer.
	Control values not within range	• Check if the flowcell is placed in the light path of the cuvette. • Wash flowcell with washing solution and rinse with distilled water. • Check QC lot n° and limit ranges. • Prepare new reagent and fresh calibrator and QC • Perform a new reagent BLANK and CALIBRATION and test QC again • Tested again with new QC and reagent. • Call technical service.

7. TROUBLESHOOTING GUIDE

GENERAL PROBLEM	PROBLEM DESCRIPTION	LIKELY CAUSE or CORRECTIVE ACTION
TEST RESULTS	The result reproducibility is low	• There are bubbles inside the flow-cell, and it should be washed. Check the aspiration condition of flow-cell. • The usage amount of reagent is too little, and the dosage should be increased. Clean with washing solution and distilled water
	Analyzer working but suddenly LCD display becomes faint and alphanumeric characters are not seen.	• Check for any heavy fluctuation in mains supply. Should not be more than +/-10%. • Check loose connection of main cable or DC jack pin.
SEMI-AUTOMATIC ANALYZER	External keyboard is not working	• Check the cable connected to the rear panel of the analyzer. • In the event that USB-C ports fail, contact Technical Service • Connect the analyzer through the internet and work through your PC.

7. TROUBLESHOOTING GUIDE

GENERAL PROBLEM	PROBLEM DESCRIPTION	LIKELY CAUSE or CORRECTIVE ACTION
REACTION TEMPERATURE	The temperature attained in the flow cell is not within $\pm 0,5^{\circ}\text{C}$ of the selected temperature.	• Ensure that the flow cell is completely inserted in the slot and is in its proper position. • Refrigerated reagents should be allowed to attain Room Temperature before use. • Adjust Temperature by using utility function menu for Temperature adjustment • Call Technical Service.

7. TROUBLESHOOTING GUIDE

GENERAL PROBLEM	PROBLEM DESCRIPTION	LIKELY CAUSE or CORRECTIVE ACTION
THERMAL PRINTER	Internal printer is not working	• Make sure that the paper roll is installed pro-berly with the heat sensitive surface up. (The analyzer will not print if the paper is fed reverse). • Check if paper is thermal paper as per the specifications. • Check printer lever locking system. It should be in closed position. • Call Technical Service.




8. ACCESSORIES AND SPARE PARTS

8. ACCESSORIES AND SPARE PARTS

8.1 End User



AC17433 Paper Roll 	AC17435 Flow Cell 	AC17436 Flow Cell's Aspiration Tube 	AC17437 Flow Cell's Waste Tube 	AC17438 Peristaltic Pump Tube 	AC17440 Waste External Tube 
AC17441 Bottle Waste 	AC17442 Fan Filter Kit 	AC10453 Test Tube Adaptor 	AC17443 Washing Solution Bottle (4x100mL) 		

If any of the accessories are lost or damaged, or if consumable material needs to be re-stocked, always use original material.
In this case, contact the usual dealer and order the spare parts using their corresponding codes.

8. ACCESSORIES AND SPARE PARTS

8.1 End User - electrical components








AC17439	AC14694	CA10456	CA10455
Power-in Adapter 12V3.5A	Mains cable UK	Mains cable America	Mains cable Europe
			

If any of the accessories are lost or damaged, or if consumable material needs to be re-stocked, always use original material.
In this case, contact the usual dealer and order the spare parts using their corresponding codes.

8. ACCESSORIES AND SPARE PARTS

8.2 Technical Service

AC17444 Printer 	AC17445 Peristaltic Pump 	AC17446 CPU 	AC17448 Photometer PCB 	AC17452 Optical Bench 	AC17450 General FAN 
AC17451 Fittings KIT 	AC17447 Display 	AC17449 Flow Cell Cover 			

If any of the accessories are lost or damaged, or if consumable material needs to be re-stocked, always use original material.
In this case, contact the usual dealer and order the spare parts using their corresponding codes.

8. ACCESSORIES AND SPARE PARTS

8.3 Warranty Limits

Improper use (falls, negligence, electrical mains outside tolerances, inappropriate environmental or location conditions) and the internal manipulation of the instrument by persons not authorised by the manufacturer or without using original consumables and spares (tubes, etc.) will render the warranty invalid.

Consumer components such as the peristaltic pump tubes, Teflon tubes are not included in the warranty.

ATENCION

9. TECHNICAL SPECIFICATIONS

9. TECHNICAL SPECIFICATIONS

Compliance with applicable directives and legislation: 98/79/CE Directive on in vitro diagnostic medical devices.

Optical System

Nominal range	-0.2 to 3.5 A
N° of wavelengths	8
Wavelengths installed	340, 405, 505, 535, 560, 600, 635, 670
Band width of each wavelength	10 ± 2 nm
Wavelength error (1)	± 2 nm
Digital resolution	0.0001 A for 2.0 A
Screen resolution	0.001 A
Baseline stability	better than 0.001A in 60 minutes
Accuracy of the measurement (1)	CV ≤ 1% at 0.1 A CV ≤ 0.1% at 2 A
Trueness of the measurement (1)	±5% at 0.1 A ± 2% at 1.0 A ± 2% at 2.0 A ± 5% at 3.0 A

Thermostatic System

Thermostatisation range	25° C to 37° C
Accuracy of the temperature	± 0.5°C
Temperature stability	± 0.2° C in 30 minutes

(1) para 340, 405 y 505 nm

9. TECHNICAL SPECIFICATIONS

Compliance with applicable directives and legislation: 98/79/CE Directive on in vitro diagnostic medical devices.

Aspiration System

Aspiration tube (external length)	approx. 80 mm
-----------------------------------	---------------

Flow Cuvette

Length of light passage	10 mm
Diameter of light passage	1.5 mm
Volume	18 µL
Quality of the glass	Quartz
Standard cuvettes	Macro, semi-micro, micro
Test tubes	Round base with a diameter of 12 mm x 75 mm in length

Peristaltic Pump

Type of operation	Stepper motor
Nominal flow	10 mL/min
Aspiration volume	100 µL - 5000 µL
Cross contamination	< 5% at 200 µL
	< 3% at 300 µL
	< 0,7% at 400 µL
Waste bottle capacity	1L

9. TECHNICAL SPECIFICATIONS

Compliance with applicable directives and legislation: 98/79/CE Directive on in vitro diagnostic medical devices.

Processing Unit

CPU	H6 Quad-core 64-bit 1.8GHZ ARM Cortex™-A53
Memory+Onboard Storage	1GB LPDDR3 (shared with GPU) +8GB EMMC Flash
USB3.0 Ports	2 x USB 3.0 Host
Onboard Network	10/100M/1000M, Ethernet RJ45
Database capacity (results)	>1.000.000 results

Dimensions and Weight

Size (H x W x D)	180mm x 245mm x 438mm (7.1in x 9.6in x 17.3in)
Weight	4290 gr. (9.43 lb)

Environmental Conditions

Use	Indoor
Height	< 2000 m
Room temperature	10° C – 35° C

External Power Supply

Voltage	100 V to 240 V AC
Frequency	50/60 Hz
Output voltage	15V CC
Output power	30W

Instrument Power Consumption

Average	5 W performing measurements
Minimum	2 W on standby

10. END OF PRODUCT LIFE

10. END OF PRODUCT LIFE

Compliance with applicable directives and legislation: 98/79/CE Directive on in vitro diagnostic medical devices.

At the end of the product life of the analyser, disposal of the unit must be done in accordance with the environmental legislation enforced in each country. If that country is an EU member state, the WEEE directive on electrical and electronic appliances will apply. In other words, when the appliance's useful life has ended, it is converted into waste and must be separated from household waste for correct recycling. For this purpose, contact your regular distributor for them to carry out the recycling process.

Use gloves to handle any of the analyzer that has been in contact with biological material: tips, containers, tubes, etc.

Dispose of this material in accordance with current legislation for the disposal of hazardous biological waste from your national or local government.

11. MANUAL VERSION

11. MANUAL VERSION

Manual Version	Software Version	Date Revision	Modification
1.0	1.0	January 2021	Initial version
1.1	1.0	May 2021	Changes in the illustrations - Position of the Flow Cell - pages 9 and 14

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