
I, PHILLIP JONES, of Windsor House, Victoria Street, Windsor,
Berkshire, England, Notary Public duly admitted and sworn,

DO HEREBY CERTIFY AND ATTEST:

That I have no reason to doubt the authenticity of the document
hereunto annexed namely PICO Single Use Negative Pressure Wound
Therapy System, dated 03rd December 2019 and signed for and on
behalf of Smith & Nephew Medical Limited by Helen May, Snr
Regulatory Affairs Manager.

IN TESTIMONY WHEREOF I have hereunto set my hand and
affixed my Seal of Office, in the Town of Windsor aforesaid, this day
Friday, 6th December 2019.

PHILLIP JONES
Notary Public

+44 (0) 1753 851 591

My Commission expires on death

Smith & Nephew Medical Limited
101 Hessle Road
Hull
HU3 2BN
England

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101 Hessle Road, Hull, HU3 2BN, United Kingdom

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PICO Single Use Negative Pressure Wound Therapy System

INSTRUCTION FOR USE

System variants

Low to moderately exuding wounds - 7-day system

(2 dressings, 1 pump, secondary fixation strips):

Multisite 15x20	66801356
Multisite 20x25	66801357
10cm x 20cm	66801358
10cm x 30cm	66801359
10cm x 40cm	66801360
15cm x 15cm	66801361
15cm x 20cm	66801362
15cm x 30cm	66801363
20cm x 20cm	66801364
25cm x 25cm	66801365

This user guide is not intended as a guarantee or warranty. It is intended only as a guide. For medical questions please consult a physician.

The product must be used in accordance with this user guide and all applicable labelling.

1. Description

The PICO Single Use Negative Pressure Wound Therapy System consists of a pump and sterile dressing(s). The PICO pump maintains negative pressure wound therapy (NPWT) at 80 mmHg (nominal) to the wound surface. Exudate is managed by the dressing through a combination of absorption and evaporation of moisture through the outer film.

PICO is intended for use in wound sizes (surface area x depth) up to 400 cm³ which are considered to be low to moderately exuding.

2. Indications for use

PICO is indicated for patients who would benefit from a suction device (negative pressure wound therapy) as it may promote wound healing via removal of low to moderate levels of exudate and infectious materials.

Examples of appropriate wound types include:

- Chronic
- Acute
- Traumatic
- Subacute and dehiscent wounds
- Partial-thickness burns

- Ulcers (such as diabetic or pressure)
- Flaps and grafts
- Surgically closed incision sites

PICO Single Use Negative Pressure Wound Therapy System is suitable for use both in a hospital and homecare setting.

3. Contraindications

The use of PICO is contraindicated in the presence of:

- Patients with malignancy in the wound bed or margins of the wound (except in palliative care to enhance quality of life).
- Previously confirmed and untreated osteomyelitis.
- Non-enteric and unexplored fistulas.
- Necrotic tissue with eschar present.
- Exposed arteries, veins, nerves or organs.
- Exposed Anastomotic sites.
- Emergency airway aspiration.
- Pleural, mediastinal or chest tubing drainage.
- Surgical suction.

4. Warnings

1. Certain patients are at high risk of bleeding complications which, if uncontrolled, could potentially be fatal.

Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately disconnect pump, leave dressing in place, take appropriate measures to stop bleeding and seek immediate medical assistance.

2. The use of anticoagulants does not deem a patient inappropriate for treatment with PICO; however, haemostasis must be achieved

before applying the dressing. Patients suffering from difficult haemostasis or who are receiving anticoagulant therapy have an increased risk of bleeding. During therapy, avoid using haemostatic products that may increase the risk of bleeding, if disrupted. Frequent assessment must be maintained throughout the therapy.

3. At all times care should be taken to ensure that the pump and tubing and connectors do not:

- Lie in a position where it could cause pressure damage to the patient.
- Trail across the floor where it could present a trip hazard or become contaminated.
- Present a risk of strangulation or a tourniquet to patients.
- Rest on or pass over a source of heat.
- Become twisted or trapped under clothing or bandages so that the negative pressure is blocked.

4. Sharp edges or bone fragments in a wound must be covered or removed prior to using PICO due to risk of puncturing organs or blood vessels while under negative pressure.

5. In the event that defibrillation is required, disconnect the pump from the dressing prior to defibrillation. Remove the dressing if it is positioned in a location that will interfere with defibrillation.

6. MR Unsafe. PICO is not MRI compatible. Do not take PICO into the MRI suite.

7. PICO has not been studied on paediatric patients. Patient size and weight should be considered when prescribing this therapy.

8. PICO is unsuitable for use in areas where there is danger of explosion (e.g. hyperbaric oxygen unit).

9. The battery compartment door could represent a choking hazard for

young children. Keep out of the reach of children.

10. PICO is not suitable for use in the presence of flammable anesthetic mixture with oxygen or nitrous oxide.

5. Precautions

1. Precautions should be taken in the following types of patients who are at high risk of bleeding complications:

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2. PICO dressings should only be applied by a healthcare professional. Dressings are not to be removed or changed by the patient.

3. Where PICO is used on infected wounds, more frequent dressing changes may be required. Regular monitoring of the wound should be maintained to check for signs of infection.

4. If deemed clinically appropriate, care should be taken that the application of a circumferential dressing or the use of negative pressure on ischaemic limbs does not compromise circulation.

5. PICO does not contain audible alerts. The pump should be carried so that it is accessible and the patient/ healthcare professional can check the status routinely.

6. Although PICO can be used under clothing/bedding, it is important that occlusive materials e.g. film dressings, are not applied over the pad area of the dressing as this will impair the intended evaporation of moisture through its outer layer.

7. The PICO dressing should not be covered by rigid immobilisation devices or casts which might apply excessive pressure and cause tissue injury at the wound site, especially where the tubing enters the dressing (i.e. the device's soft port).

8. Prolonged placement of rigid or opaque materials over the PICO dressing may prevent the regular inspection and assessment of the wound, and disrupt scheduled or required dressing changes.

9. Where PICO is used on patients with fragile skin, a skin protectant such as SECURA™ should be used on areas of skin where fixation strips are to be applied. Inappropriate use or repeated application of fixation strips may otherwise result in skin stripping.

10. If reddening or sensitisation occurs, discontinue use and contact the treating healthcare professional.

11. Do not use PICO with oil-based products such as petrolatum as it may compromise establishing an effective seal.

12. The use of negative pressure presents a risk of tissue ingrowth into foam when this is used as wound filler. When using foam filler with PICO, tissue ingrowth may be reduced by using a non-adherent wound contact layer or by increasing the frequency of dressing changes.

13. PICO may be used in conjunction with surgical drains provided the dressing is not placed over tubing where it exits the skin. Any surgical drain should be routed under the skin away from the edge of the dressing and function independently of the PICO Single

Use Negative Pressure Wound Therapy System.

14. Pump must be protected from sources of fluid e.g. from incontinence or spillages. Discontinue device use if fluid ingress is observed.

15. When showering, the PICO pump should be disconnected from the dressing. Ensure the end of the tubing attached to the dressing is facing down so that water does not enter the tubing.
 16. Do not take the pump apart.
 17. The dressing should not be used with any other suction pump.
 18. Do not alter or cut tubing configuration or pull on the tubing or soft port.
 19. Do not cut the dressing pad as this may lead to loss of NPWT application. For conformability in awkward areas the film border may be cut up to 1cm away from the edge of the pad.
 20. Always ensure that the dressing is positioned centrally over the wound. The soft port should be positioned uppermost on intact skin and not extend over the wound so that the risk of fluid collecting around the soft port and potentially blocking the negative pressure is minimised.
 21. CT scans and x-ray have the potential to interfere with some electronic medical devices. Where possible, move the pump out of the x-ray or scanner range. If the pump has been taken into the CT scan or x-ray range, check that the system is functioning correctly following the procedure.
 22. The system is single use only. Use of any part of this system on more than one patient may result in cross- contamination that may lead to infection.
 23. High temperatures and humidity may reduce wear time of dressings.
 24. During transport, there is a potential for radio frequency interference that could affect PICO performance. If the pump malfunctions, replace batteries. If not corrected, contact your healthcare professional to replace the system.
- PICO is not intended for use aboard aircraft, the batteries should be removed during air travel.
25. The potential for electromagnetic interference in all environments cannot be eliminated. Use caution if PICO is near electronic equipment such as RFID (Radio Frequency Identification) readers, anti-theft equipment or metal detectors.
 26. Where PICO is used to bolster skin grafts, it is important to visually inspect the system regularly, especially in the first week of treatment to ensure that negative pressure is continually applied and a seal is maintained.

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7. Instructions for use

7.1. Guidance on wound suitability for management with PICO

PICO should be used on wounds which fit comfortably within the area of the pad, observing precautions on port positioning (on intact skin and not extending over the wound).

As a guide:

- Depth - Wounds greater than 0.5cm in depth are likely to require a foam or gauze NPWT filler to ensure adequate treatment of all the wound surfaces. Wounds treated with the larger sizes of the PICO system should generally be no more than 2cm in depth.
- Exudate – When used over a moderately exuding wound, the size of the wound should generally be no more than 25% of the dressing pad area.

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2. Using a clean technique, peel off the central release handle and place the dressing centrally over the wound to reduce the chance of wound fluid coming into contact with the soft port. Ensure the dressing lies flat to the wound and surrounding skin.

The soft port should be uppermost from the wound (depending on the patient's primary position), placed on intact skin and not extending over the wound to prevent fluid pooling around the soft port and blocking the negative pressure. Remove the other two handles and smooth the dressing around the wound to prevent creasing. Reposition if required to ensure border is not creased.

3. Once the dressing is in place, remove the pump and the batteries from the tray. The direction in which the batteries should be placed is indicated inside the battery compartment. Insert the batteries. Replace the cover.

Following this all three indicators should flash once (Refer to section (i) of Table 1).

4. Join the pump to the dressing by twisting together the tubing connectors. Press the orange button to start the application of negative pressure. The green indicator will start to flash (indicates pump working to establish NPWT). Depending on the size of the wound, the pump should take up to 30 seconds to establish negative pressure wound therapy. If NPWT has not been established the pump will wait for 15 seconds and then take a further 30 seconds to establish NPWT (the green indicator will continue to flash during this time). If unsuccessful the yellow air leak indicator will illuminate. To troubleshoot refer to section (ii) of Table 1.

5. If using SECURA prior to application of the fixation strips (see section 5. Precautions) wipe the area surrounding the dressing and allow skin to dry.

6. Apply the fixation strips to each of the four sides of the dressing. Remove top carrier on the strip after each one has been applied. These strips maintain the seal over the wear time of the dressing.

In awkward areas, it may be useful to apply the strips to help achieve a seal prior to switching on the pump.

Place each strip so that it overlaps the dressing border by approximately 1cm.

Ensure tubing is not twisted or trapped between clothing.

Please note that if at any time the fixation strips are removed, the dressing should also be replaced.

7.3. Dressing change

1. Dressings should be changed in line with standard wound management guidelines, typically every 3-4 days. More frequent dressing changes may be required depending on the level of exudate, condition of the dressing or other patient considerations; e.g. when PICO is used on infected wounds. At the healthcare professional's discretion a PICO dressing may be left in place for up to 7 days.

2. Inspect the dressing regularly. If the dressing appears ready for changing; press the orange button and disconnect the dressing from

the pump. The fixation strips should be stretched away from the skin and the dressing lifted at one corner and peeled back until it has been fully removed. Apply another dressing as per section 7.2. Application, connect to the pump and press the orange button to reinitiate the therapy.

3. Based on dressing change frequency, a new PICO Single Use Negative Pressure Wound Therapy System kit will be required dependent on whichever of the following occurs first - either when the supplied dressing(s) have been used or when the pump life ends and it automatically stops functioning (all the indicators will turn off at this point).

4. The dressing should be disposed of as clinical waste. The batteries should be removed from the pump; and

both batteries and pump disposed of according to local regulations.

5. For additional information on disposal requirements see: www.mypico.com or speak to your Smith & Nephew representative.

7.4. Use with fillers and wound contact layers

PICO is compatible **with standard gauze and foam fillers** used in traditional NPWT where this is clinically appropriate – for example on a defect wound. When filler is used, the filler and the PICO dressing should be changed 2 to 3 times a week, according to local clinical protocol and manufacturer's instructions. Gauze should loosely fill to the surface of the wound.

Avoid over packing.

PICO may be used over the top of a non- adherent layer if required, for example over a skin graft. Gel strips may be used around the border to help achieve a seal.

On infected wounds or wounds at risk of infection, ACTICOAT™ Flex antimicrobial dressings may be used under PICO.

PICO is compatible with **compression therapy** if clinically required; however, the wound should be first accurately assessed to ensure the exudate is of low

to moderate volume. Application of active NPWT to a wound may cause production of additional exudate; so following the first application of PICO under a compression system, check the PICO dressing after 2-3 days to see if fluid handling is acceptable or if a dressing change is required and that it is appropriate to continue with compression therapy.

When applying compression therapy over PICO, ensure tubing connecting the PICO pump and dressing runs over the top of the first layer and beneath subsequent layers to avoid it being held against

the skin. If using compression hosiery, run tubing outside of stockings. For guidance on the correct application of a compression therapy, refer to the relevant manufacturer's instructions.

8. General Use

8.1. Showering and bathing

Light showering is permissible; however, the pump should be disconnected (see section 5. Precautions) and placed in a safe location where it will not get wet.

The dressing should not be exposed to a direct spray or submerged in water. Ensure the end of the tubing attached to the dressing is facing down so that water does not enter the tubing.

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
9. Faults and technical assistance

If your pump develops a fault or there are signs of damage, refer to section (ii) of Table 1.


Table 1 – Pump status indication and faults


PICO has visual indicators to let the user know when there is an issue. PICO does not contain audible alerts. The pump should be carried so that it is accessible and the patient/ healthcare professional can check the status routinely.

Section (i) – Normal function

Display status	Indicator status	Possible cause	Comments/trouble shooting
	All indicators off.	The pump is OFF. The pump has reached the end of its life. The batteries are no longer functional.	The therapy has been paused. Pressing the orange button will restart the therapy and the green indicator will flash. When the pump reaches the end of its life, it will automatically cease functioning, in this case all the indicators will turn off. Pressing the orange button will not provide a green flashing indicator. If the pump has had less than the specified days usage, the batteries may not be functional and should be replaced.
	All indicators flash once.	This reflects the pump self-test once batteries have been inserted and the cover has been replaced.	This is expected.
	Green 'OK' indicator flashes	Dressing applied, and full system is functioning properly. No issues.	The pump may be heard running occasionally as it maintains the negative pressure. This is normal. If this occurs frequently (several times an hour) smooth down the dressing to remove any creases that may be allowing air into the system. NPWT is still being applied in this situation.

Section (ii) – Alerts and faults

Display status	Indicator status	Possible cause	Comments/trouble shooting
	Yellow 'leak' indicator flashes	Air leak detected possibly due to a	Smooth down the dressing and the strips to remove any creases that are allowing air into the system.

		creased dressing/ border/strip. Pump is in auto pause. NPWT is not being applied to the wound. The pump will auto pause for 1 hour and then will automatically try to re-establish therapy if no remedial action is taken.	Press the orange button to restart the therapy. The green "OK" indicator will flash as the pump tries to establish therapy. Depending on the size of the wound, the pump should take up to 30 seconds to establish negative pressure wound therapy. If NPWT has not been established the pump will wait for 15 seconds and then take a further 30 seconds to establish NPWT (the green indicator will continue to flash during this time). If unsuccessful the yellow air leak indicator will illuminate. If this happens, repeat smoothing actions and press the orange button. If the leak is resolved the green light will continue to flash.
	Green 'OK' indicator flashes. Yellow 'battery low' indicator flashes	System on and functioning properly. Battery power low	Change of batteries required in <24 hours. Pause therapy by pressing the orange button Push open the battery cover at the top of the pump and remove the old batteries. Insert 2 new AA lithium (L91) batteries. Replace the cover. Press the orange button, the therapy will re-start and the green indicator will flash.
	Yellow 'leak' indicator flashes. Yellow 'battery low' indicator flashes.	Creased dressing/ border/strip. Battery power low.	Address air leak as above. Change of batteries or device required in <24 hours as above..
	All indicators solidly	Pump failed. System is not usable.	Contact S&N representative. Apply new pump and dressing.

	illumina ted.		
--	------------------	--	--

10. Specifications

Maximum Dimensions	Height (mm) 63.96 ± 0.5 Front side width (mm) $12.5 +0.4/-0.1$ Back side width (mm) 9.75 ± 0.5 Length (mm) 71.52 ± 0.5
Weight	<120 g
Maximum Vacuum	13,34 kPa (100 mmHg)
Mode of Operation	Continuous
Storage/Transport	5 – 25°C, 10 – 75% RH 700 to 1060 mbar atmospheric pressure
Operating Environment	5 – 35°C, 10 – 95% RH 700 to 1060 mbar atmospheric pressure

11. Safety and electromagnetic compatibility

When used in accordance with the manufacturer's instructions, PICO complies with the general requirements for safety of electrical medical equipment IEC 60601-1 and the electromagnetic safety requirements of electrical medical equipment IEC 60601-1-2.

Electromagnetic compatibility

This equipment has been tested and found to comply with the limits for medical devices to IEC 60601-1-2. These limits are designed to provide reasonable protection against electromagnetic interference typical of medical installations and home use environment. Portable and mobile RF communication equipment can affect Medical Electrical Equipment. This equipment generates, uses, and can radiate radio frequency energy and if not used in accordance with the instructions, may cause harmful interference to other devices in the vicinity.

For additional information on electromagnetic immunity and electromagnetic emissions see: www.mypico.com or ask your Smith & Nephew representative for a hardcopy.

PICO is packed in the UK with individual components made in the following countries:

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Trade Marks of Smith & Nephew www.smith-nephew.com

©Smith & Nephew

Date of issue 06/2017

UK Croxley Park, Building 5, Hatters Lane, Watford, Hertfordshire WD18 8YE Tel: 08000 157 573.

AUSTRALIA Smith & Nephew Pty Ltd.,

85 Waterloo Road, North Ryde NSW 2113. Tel +61 2 9857 3999.

AUSTRIA Smith & Nephew GmbH,

Concorde Business Park C3, A-2320 Schwechat. Tel: +43 (0)1 707 91 02.

BELGIUM Smith & Nephew N.V.-S.A., Hector Henneaulaan 366 - 1930 Zaventem.

Tel: +32 2 702 29 11.

CANADA Smith & Nephew Inc.,

2280 Argentia Road, Mississauga, ON L5N 6H8. Tel: 1 800 463 7439.

CHINA Smith & Nephew Ltd,

12F Headquarters Building, 168 Xizang Road (M), Huangpu District, Shanghai 200001.

Tel: +86 21 2330 3000.

DENMARK Smith & Nephew A/S,

Slotsmarken 14, 2970 Hørsholm, Tel: +45 45 80 61 00.

EASTERN EUROPE Smith & Nephew GmbH, Concorde Business Park C3, A-2320 Schwechat. Tel: +43 (0)1 707 91 02.

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26th km, Paeonias - Markopoulou Ave, PO Box 182, 19400 Koropi. Tel: 662 6201-7.

HONG KONG Smith & Nephew Ltd.,

Unit 813-818, 8/F, Delta House, 3 On Yiu Street, Shatin, N.T. Tel: (852) 2648 7700.

INDIA Smith and Nephew HealthCare Pvt. Ltd.

B-501-509, Dynasty Business Park, Andheri-Kurla Road, Andheri (E), Mumbai-400059. Tel: 022 4005 5090.

IRELAND Smith & Nephew Limited, Carraig Court, George's Avenue, Blackrock, Co. Dublin, Tel: 01-217 0444.

ITALY Smith & Nephew S.r.l. - via De Capitani 2A,

20864 Agrate Brianza (MB). Tel: +39 039 60941.

JAPAN Smith & Nephew KK, 2-4-1 Shiba-Koen, Minato-ku, Tokyo 105-0011 Japan Tel: 03-5403-8930

KOREA Smith & Nephew Ltd,

13F, ASEM Tower 517, Yeongdong-daero, Gangnam-gu,

Seoul, Korea (135-798). Tel: 2 6001 7575.

MALAYSIA Smith & Nephew Healthcare Sdn. Bhd. Suite 11.01B & 11.02, Level 11, Menara AmFIRST, No 1 Jalan 19/3, 46300 Petaling Jaya, Selangor. Tel: +603-7958 7103.

MALTA Cherubino Ltd,

Delf Building, Sliema Road, Gzira. Tel: 21 343270.

MEXICO Smith & Nephew,

SA de CV, Patricio Sanz No. 442, Colonia del Valle, C.P. 03100, México D.F., Tel. 53 40 22 60.

MIDDLE EAST Smith & Nephew FZE,

4th Floor National Bank of Abu Dhabi Building,

Al Jazari Building 47, Dubai Healthcare City - District 8, PO Box 9715, Dubai, U.A.E. Tel: 971 4 4299111.

NETHERLANDS Smith & Nephew B.V.,

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2132 NB Hoofddorp, Tel: (+31) (0) 20 6543999.

NEW ZEALAND Smith & Nephew Limited, Unit A, 36 Hillside Road, Wairau Valley, Auckland 0627 Tel: +64 (0) 9 820 2840.

NORWAY Smith & Nephew A/S, PO Box 224,
Nye Vakas vei 64, 1379 Nesbru, Tel: +47 66842020.

PAKISTAN BSN Medical (PVT) Ltd,
A/69 SITE, Manghopir Road, PO Box 3659, Karachi 75700. Tel: + 92 21 32561176-77.

POLAND Smith & Nephew sp. z o.o. Ul. Osmańska 12 02-823 Warszawa
Tel: +48 22 360 41 20.

PORTUGAL Smith & Nephew Lda.,
Estrada Nacional n° 10, Km 131, Parque Tejo - Bloco C, 2625-445 Forte da Casa. Tel: +351 214 460 650.

SINGAPORE Smith & Nephew Pte Ltd.,
1A International Business Park, #09-03 Tolaram 609933. Tel: (+65) 6270 0552.

SOUTH AFRICA Smith & Nephew (Pty) Limited, 30 The Boulevard, West End Office Park, Westville, 3629. Tel: 031 2428111.

SPAIN Smith & Nephew SAU, Fructuós Gelabert, 2-4, 08970 Sant Joan Despí, Barcelona. Tel: 93 373 73 01.

SWEDEN Smith & Nephew AB, Box 143, Kråketorpsgratan
20, 43122 Mölndal. Tel: 031-746 58 00.

SWITZERLAND Smith & Nephew Schweiz AG
Advanced Wound Management, Oberneuhofstrasse 10d, 6340 Baar Tel. +41 41 766 22 44.

TAIWAN (R.O.C.) Smith & Nephew (Overseas) Ltd.
Taiwan Branch, 9F-2, No. 50, Sec.1, Xincheng South Road Taipei City. Tel: 886 2 3322 5020.

THAILAND Smith & Nephew Ltd.
16th Floor Building A, The 9th Towers Grand Rama 9, 33/4 Rama 9 Road, Huaykwang Sub-District, Huaykwang District, Bangkok 10310. Tel: 66(0)2-116-5581.

TURKEY Smith ve Nephew Medikal Cihazlar Tic. LTD.
ti. Bahçeköy Mah. Orkide Sokak No:8/E 84473 Sarıyer/ ISTANBUL Tel: +90 212 226 22 55.

UNITED STATES Smith & Nephew, Inc.
5600 Clearfork Main Street, Suite 600, Fort Worth, TX 76109.

Customer Care Center: 1-817-900-4000.

VENEZUELA Ace Sales Corp, 8085 NW 68th Street, Miami, Florida 33166. Tel: 305.835.0310.

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**PICO Single Use Negative Pressure Wound Therapy
System**

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1. PICO Single Use Negative Pressure Wound Therapy System contents:
 - Pump -1pcs.
 - Multisite dressings 15cm x 20cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
2. PICO Single Use Negative Pressure Wound Therapy System contents:
 - Pump-1pcs.
 - Multisite dressings 20cm x 25cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
3. PICO Single Use Negative Pressure Wound Therapy System contents:
 - Pump-1pcs.
 - dressings 10cm x 20cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
4. PICO Single Use Negative Pressure Wound Therapy System contents:
 - Pump-1pcs.
 - dressings 10cm x 30cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
5. PICO Single Use Negative Pressure Wound Therapy System contents:
 - Pump-1pcs.
 - dressings 10cm x 40cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
6. PICO Single Use Negative Pressure Wound Therapy System contents:
 - Pump -1pcs.
 - dressings 15cm x 15cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
7. PICO Single Use Negative Pressure Wound Therapy System contents:

- Pump-1pcs.
- dressings 15cm x 20cm - 2 pcs.
- Secondary fixation strips - 12 pcs.
- extension tubing - 1 pc.
- AA batteries - 2 pcs.
- instruction for use - 8 pcs.

8. PICO Single Use Negative Pressure Wound Therapy System contents:

- Pump-1pcs.
- dressings 15cm x 30cm - 2 pcs.
- Secondary fixation strips - 12 pcs.
- extension tubing - 1 pc.
- AA batteries - 2 pcs.
- instruction for use - 8 pcs.

9. PICO Single Use Negative Pressure Wound Therapy System contents:

- Pump -1pcs.
- dressings 20cm x 20cm - 2 pcs.
- Secondary fixation strips - 12 pcs.
- extension tubing - 1 pc.
- AA batteries - 2 pcs.
- instruction for use - 8 pcs.

10. PICO Single Use Negative Pressure Wound Therapy System contents:

- Pump-1pcs.
- dressings 25cm x 25cm - 2 pcs.
- Secondary fixation strips - 12 pcs.
- extension tubing - 1 pc.
- AA batteries - 2 pcs.
- instruction for use - 8 pcs.

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Fixation strips: Polyethylene woodless bleached paper with solventless silicone coating, orange print ink, A8 acrylic adhesive, EU45 polyurethane film, anti-blocking agent, polyester film with silicone and matte polyethylene coating, orange print ink, polypropylene tape

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Article/catalogue number:

Name: PICO Single Use Negative Pressure Wound Therapy System

Size*: 10x20 /10x30 /10x40 /15x15 /15x20 /15x30/ 20x20/ 25x25/ Multisite 15x20/ Multisite 20 x25

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Precautions: -

Symbol of conformity to Russian standards (PCT mark);




Instruction for Use: <http://www.smith-nephew.com/about-us/where-we-operate/europe/russia/ifus/>

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3.2 Data on labeling

Symbol/text	Meaning/translation
PICO ®	PICO
	First remove protective liner from the fixation strip
	Then remove top carrier on the fixation strip
	"Alarm: air leak" indicator
	"Battery low" indicator
OK	"OK" indicator - the pump is functioning properly
	Manufacturer
	Date of manufacture
 	Packaging contains 2 AA batteries E91 (PICO 7), L 91 (PICO, PICO 7Y)
IP24	Protection degree of the housing of electrical equipment: the device is protected from finger penetration and cannot be damaged by splashes of water falling in any direction.
3V	Device voltage
STERILE EO	The product is EtO sterilized

	Single use. Do not reuse
	Warning!
	Incompatible with magnetic resonance - keep away from magnetic resonance equipment
	International classification
	EC: Do not dispose of as household waste
	Classification of equipment: BF type, isolated
	EU Mark of Compliance
	Follow the instructions for use
	Use by medical prescription only
	Start / Pause / Resume treatment
	Follow the instructions for use
	Keep away from direct sunlight
	Protect from moisture
	Batch number
	Storage temperature
	Do not use if package is opened or damaged
	1 pump included

	The service life of the device is up to 7 days
	Dressing size
	Shelf life
	Unique product identifier
	International labeling and accounting code
	Barcode
	Catalog number
	Recyclable packaging material
	European authorized representative
	Packaging contains 6 fixation strips
	Open by folding back the paper edge
	Schematic instructions for use
	PICO dressing included
	Schematic instructions for use
4. MD DISPOSAL OR DESTRUCTION PROCEDURE AND CONDITIONS	

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8. LIST OF STANDARDS

Smith & Nephew Medical Ltd. declares that **PICO Single Use Negative Pressure Wound Therapy System, variations** meets the requirements of Directive 93/42/EEC, a full list of international requirements is available on request.

9. MD MAINTENANCE AND REPAIR REQUIREMENTS

Not applicable.

10. MD PERFORMANCE AND SAFETY CHECKING

Not applicable.

11. WARRANTY AND CLAIMS

Claimee: Smith & Nephew LLC, 105120, Russia, Moscow, 2nd Syromyatnichesky Lane, 1, 9th floor, office 1, room 1. tel/fax: + 7 495 755 55 03

WARRANTY: The manufacturer guarantees compliance with the declared characteristics, the period of preservation of sterility, the safety and integrity of the packaging subject to proper transportation, operation, and storage of the medical device for 2 years from the date of manufacture/sterilization. **Shelf life** - 24 months from the date of manufacture / sterilization indicated on the package, provided that it is stored in the original packaging, in a dry, well-ventilated place, protected from direct sunlight, at a temperature not exceeding 25°C.

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MANUFACTURING SITE:

1. Smith & Nephew Medical Limited, 101 Hessle Road, Hull, HU3 2BN, United Kingdom

2. Flextronics Technology (Shenzhen) Co., Ltd. 3-6/F, Building 1 and 6/F, Building 4, Yusheng Industrial Park, 467 Xixiang Section, National Highway 107, Xixiang, Baoan District, 518126 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA.

3. Avery Dennison Belgie BVBA, Tieblokkenlaan 1, Turnhout B-2300, Belgium.

4. Smith & Nephew Medical (Suzhou) Limited, No. 12, Wuxiang Road Comprehensive Free Zone, West Zone Suzhou Industrial Park Suzhou Jiangsu 215021 China.

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Smith & Nephew LLC, 105120, Russia, Moscow, 2nd Syromyatnichesky Lane, 1, 9th floor, office 1, room 1. tel/fax: + 7 495 755 55 03

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КИТАЙ Smith & Nephew Ltd, 12F Headquarters Building, 168 Xizang Road (M), Huangpu District, Shanghai 200001.

Tel: +86 21 2330 3000.

ДАНИЯ Smith & Nephew A/S, Slotsmarken 14, 2970 Hørsholm,
Tel: +45 45 80 61 00.

ВОСТОЧНАЯ ЕВРОПА Smith & Nephew GmbH, Concorde Business Park C3, A-2320 Schwechat.
Tel: +43 (0)1 707 91 02

ФРАНЦИЯ Smith & Nephew SAS, 25 bd Alexandre Oyon 72019 LE MANS Cedex 2.
Tel: 02 43 83 23 23.

ФИНЛЯНДИЯ Smith & Nephew Oy, Äyritie 12 C, 01510 Vantaa.
+358 207 866 300.

ГЕРМАНИЯ Smith & Nephew GmbH, Friesenweg 4/ Haus 21, 22763 Hamburg,
Tel: 040 879744-0.

ГРЕЦИЯ Costas A. Papaellinas (Hellas) SA, 26th km, Paeonias – Markopoulou Ave, PO Box 182, 19400 Koropi.
Tel: 662 6201-7.

ГОНКОНГ Smith & Nephew Ltd., Unit 813-818, 8/F, Delta House, 3 On Yiu Street, Shatin, N.T.
Tel: (852) 2648 7700.

ИНДИЯ Smith and Nephew HealthCare Pvt. Ltd. B-501-509, Dynasty Business Park, Andheri-Kurla Road, Andheri (E), Mumbai-400059.
Tel: 022 4005 5090.

ИРЛАНДИЯ Smith & Nephew Limited, Carraig Court, George's Avenue, Blackrock, Co. Dublin,
Tel: 01-217 0444.

ИТАЛИЯ Smith & Nephew S.r.l. – via De Capitani 2A, 20864 Agrate Brianza (MB).
Tel: +39 039 60941.

ЯПОНИЯ Smith & Nephew KK, 2-4-1 Shiba-Koen, Minato-ku, Tokyo 105-0011 Japan
Tel: 03-5403-8930.

КОРЕЯ Smith & Nephew Ltd, 13F, ASEM Tower 517, Yeongdong-daero, Gangnam-gu, Seoul, Korea (135-798).
Tel: 2 6001 7575.

МАЛАЙЗИЯ Smith & Nephew Healthcare Sdn. Bhd. Suite 11.01B & 11.02, Level 11, Menara AmFIRST, No 1 Jalan 19/3, 46300 Petaling Jaya, Selangor.
Tel: +603-79587103.

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Smith & Nephew FZE, 4th Floor National Bank of Abu Dhabi Building, Al Jazari Building 47, Dubai Healthcare City – District 8, PO Box 9715, Dubai, U.A.E.
Tel: 971 4 4299111.

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ШВЕЙЦАРИЯ Smith & Nephew Schweiz AG Oberneuhofstrasse 10d, CH – 6340 Baar.

Tel. +41 41 766 22 44.

ТАЙВАНЬ Smith & Nephew (Overseas) Ltd. Taiwan Branch, 9F-2, No. 50, Sec.1, Xinsheng South Road Taipei City.

Tel: 886 2 3322 5020.

ТАЙЛАНД Smith & Nephew Ltd. 16th Floor Building A, The 9th Towers Grand Rama 9, 33/4 Rama 9 Road, Huaykwang Sub-District, Huaykwang District, Bangkok 10310.

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5600 Clearfork Main Street, Suite 600, Fort Worth, TX 76109.

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1. Smith & Nephew Medical Ltd., 101 Hessle Road, Hull, HU3 2BN, United Kingdom
2. Flextronics Technology (Shenzhen) Co., Ltd. 3-6/F, Building 1 and 6/F, Building 4, Yusheng Industrial Park, 467 Xixiang Section, National Highway 107, Xixiang, Baoan District, 518126 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA.
3. Avery Dennison Belgie BVBA, Tieblokkenlaan 1, Turnhout B-2300, Belgium.
4. Smith & Nephew Medical (Suzhou) Limited, No. 12, Wuxiang Road Comprehensive Free Zone, West Zone Suzhou Industrial Park Suzhou Jiangsu 215021 China.

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[Handwritten signature of V.V. Ivanova]

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DO HEREBY CERTIFY AND ATTEST:

That I have no reason to doubt the authenticity of the document
hereunto annexed namely PICO 7 Single Use Negative Pressure Wound
Therapy System, dated 03rd December 2019 and signed for and on
behalf of Smith & Nephew Medical Limited by Helen May, Snr
Regulatory Affairs Manager.

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+44 (0) 1753 851 591

My Commission expires on death

Regulated by the Faculty Office of the Archbishop of Canterbury

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PICO 7 Single Use Negative Pressure Wound Therapy System

INSTRUCTION FOR USE

Hand this Leaflet to Patient after Dressing Application

PICO™ 7 is intended for use by or on the direction of a trained and licensed physician in accordance with these instructions for use.

1. Description

PICO 7 consists of a pump and sterile dressing(s). The PICO 7 pump maintains negative pressure wound therapy (NPWT) at 80 mmHg (nominal) to the wound surface. Exudate is managed by the dressing through a combination of absorption and evaporation of moisture through the outer film.

The PICO 7 dual dressing kit is intended to be used for up to 7 days on low to moderately exuding wounds. Low exuding wounds are considered to be up to 0.6g of liquid exudate/cm² of wound area/24 hours. Moderate exuding wounds are considered to be up to 1.1g of liquid exudate/cm² of wound area/24 hours. 1g of exudate is approximately equal to 1ml of exudate.

The frequency of dressing changes can be affected by multiple factors such as wound type, wound size, rate or volume of exudate, orientation or environmental conditions.

2. Indications for use

PICO 7 is indicated for patients who would benefit from a suction device (NPWT) as it may promote wound healing via removal of low to moderate levels of exudate and infectious materials.

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PICO 7 single use negative pressure systems are suitable for use both in a hospital and homecare setting.

3. Contraindications

PICO 7 is contraindicated for:

- Patients with malignancy in the wound bed or margins of the wound (except in palliative care to enhance quality of life).
- Previously confirmed and untreated osteomyelitis.
- Non-enteric and unexplored fistulas.
- Necrotic tissue with eschar present.
- Exposed arteries, veins, nerves or organs.
- Exposed anastomotic sites

PICO 7 should not be used for the purpose of:

- Emergency airway aspiration.
- Pleural, mediastinal or chest tube drainage.
- Surgical suction.

4. Important Information

Important information for Patients and Healthcare Professionals

Pump contains a MAGNET. You must keep the pump at least 10cm away from **implantable devices** such as **pacemakers, defibrillators and shunt valves**.

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5. Warnings

1. Certain patients are at high risk of bleeding complications which, if uncontrolled, could potentially be fatal. Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately disconnect pump, leave dressing in place, take appropriate measures to stop bleeding and seek immediate medical assistance.
2. Haemostasis must be achieved before applying the dressing, although the use of anticoagulants does not deem a patient inappropriate for treatment with PICO 7. Patients suffering from difficult haemostasis or who are receiving anticoagulant therapy have an increased risk of bleeding. During therapy, avoid using haemostatic products that may increase the risk of bleeding, if disrupted. Frequent assessment must be maintained throughout the therapy.
3. At all times care should be taken to ensure that the pump and tubing and connectors do not:
 - Lie in a position where it could cause pressure damage to the patient.
 - Trail across the floor where it could present a trip hazard or become contaminated.
 - Present a risk of strangulation or a tourniquet to patients.
 - Rest on or pass over a source of heat.
 - Become twisted or trapped under clothing or bandages so that the NPWT is blocked.
4. Sharp edges or bone fragments in a wound must be covered or removed prior to using PICO 7 due to risk of puncturing organs or blood vessels while under NPWT.
5. In the event that defibrillation is required, remove the dressing if it is positioned in a location that will interfere with defibrillation.
6. MR Unsafe. The PICO 7 pump is not MRI compatible. Remove the PICO 7 pump from the dressing before entering the MRI suite.
7. PICO 7 has not been studied on paediatric patients. Patient size and weight should be considered when prescribing this therapy.
8. PICO 7 is unsuitable for use in areas where there is danger of explosion (e.g. hyperbaric oxygen unit).
9. The system contains small parts which could represent a choking hazard for young children. Keep out of the reach of children.
10. PICO 7 is not suitable for use in the presence of flammable anaesthetic mixture with oxygen or nitrous oxide.
11. PICO dressings should only be applied and changed by a healthcare professional. The patient must not remove their dressing unless the healthcare professional believes it to be clinically appropriate.
12. Each PICO dressing (including Multisite) must be used to dress one wound only.

6. Precautions

1. PICO 7 pump contains a magnet. Care should be taken to keep pump at least 10cm from

implantable devices, such as pacemakers, defibrillators and shunt valves.

2. Precautions should be taken in the following types of patients who are at high risk of bleeding complications:

- Receiving anticoagulant therapy or platelet aggregation inhibitors or actively bleeding.
- Having weakened or friable blood vessels or organs in or around the wound as a result of, but not limited to; anastomoses, infection, trauma or radiation.
- Suffering from difficult wound haemostasis.
- Untreated for malnutrition.
- Non-compliant or combative.
- Suffering from wounds in close proximity to blood vessels or delicate fascia.

3. If pain, reddening, odour, sensitisation or a sudden change in the volume or colour of wound fluid occurs during use, contact your healthcare professional

4. Where PICO 7 is used to bolster skin grafts, it is important to visually inspect the system regularly, especially in the first week of treatment to ensure that NPWT is continually applied and a seal is maintained.

5. Where PICO dressings are used on infected wounds, more frequent dressing changes may be required. Regular monitoring of the wound should be maintained to check for signs of infection.

6. If deemed clinically appropriate, care should be taken that the application of a circumferential dressing or the use of NPWT on ischaemic limbs does not compromise circulation.

7. PICO 7 does not contain audible alerts. The pump should be carried so that it is accessible and the patient/ healthcare professional can check the status routinely.

8. Although PICO dressings can be used under clothing/bedding, it is important that occlusive materials e.g. film dressings, are not applied over the pad area of the dressing as this will impair the intended evaporation of moisture through its outer layer.

9. The PICO dressing should not be covered by rigid immobilisation devices or casts which might apply excessive pressure and cause tissue injury at the wound site, especially where the tubing enters the dressing.

10. Prolonged placement of rigid or opaque materials over the PICO dressing may prevent the regular inspection and assessment of the wound, and disrupt scheduled or required dressing changes.

11. Where PICO dressings are used on patients with fragile skin, a skin protectant such as SECURA™ No-Sting Barrier Film should be used on areas of skin where fixation strips are to be applied. Inappropriate use or repeated application of fixation strips may otherwise result in skin stripping.

12. Do not use PICO dressings with oil- based products such as petrolatum as it may compromise establishing an effective seal.

13. The use of NPWT presents a risk of tissue ingrowth into foam when this is used as wound filler. When using foam filler with PICO 7, tissue ingrowth may be reduced by using a non-adherent wound contact layer or by increasing the frequency of dressing changes.

14. PICO 7 may be used in conjunction with surgical drains provided the dressing is not placed over tubing where it exits the skin. Any surgical drain should be routed under the skin away from the edge of the dressing and function independently of the PICO 7 System.

15. The pump must be protected from sources of fluid e.g. from incontinence or spillages. Discontinue PICO 7 use if fluid ingress is observed.

16. When showering, the PICO 7 pump should be disconnected from the dressing. Whilst disconnected, ensure the end of the tubing attached to the dressing is facing down so that water does not enter the tube.

17. Do not take the pump apart.

18. The PICO dressing should only be used with PICO pumps.

19. Do not alter or cut tubing configuration or pull on the tubing or soft port.

20. Do not cut the PICO dressing pad as this may lead to loss of NPWT application.

21. Always ensure that the PICO dressing is positioned centrally over the wound. The soft port should be positioned uppermost on intact skin and not extend over the wound so that the risk of fluid collecting around the soft port and potentially blocking the NPWT is minimised.
22. CT scans and x-ray have the potential to interfere with some electronic medical devices. Where possible, move the pump out of the x-ray or scanner range. If the pump has been taken into the CT scan or x-ray range, check that the system is functioning correctly following the procedure.
23. The PICO 7 system is single use only. Use of any part of this system on more than one patient may result in cross contamination that may lead to infection.
24. High temperatures and humidity may reduce wear time of PICO dressings.
25. The PICO 7 system is intended for use in a hospital and homecare setting. The system can also be used in aircraft, car, train, and boat transportation. During transport, there is a potential for radio frequency interference that could affect PICO 7 performance. If the pump malfunctions, replace batteries. If not corrected, contact your healthcare professional to replace the system.
26. When applying dressings next to one another, ensure the dressing borders do not overlap.

7. Adverse reactions

Excessive bleeding is a serious risk associated with the application of suction to wounds which may result in death or serious injury. Careful patient selection, in view of the above stated contraindications, warnings and precautions is essential. Carefully monitor the wound and dressing for any evidence of a change in the blood loss status of the patient. Notify the healthcare professional of any sudden or abrupt changes in the volume or the colour of exudate.

8. Instructions for use

8.1. Guidance on wound suitability

PICO dressings should be used on wounds which fit comfortably within the area of the pad, observing precautions on soft port positioning (on intact skin and not extending over the wound).

PICO Multisite dressings are designed to enhance conformability when dressing awkward anatomical areas. PICO dressings (including Multisite) must be used to dress one wound only. As a guide:

- **Depth** - Wounds greater than 0.5cm in depth are likely to require a foam or gauze NPWT filler to ensure adequate treatment of all the wound surfaces. Wounds treated with the larger sizes of the PICO 7 system should generally be no more than 2cm in depth.
- **Exudate** – PICO 7 is intended for use on wounds where the level of exudate is low (up to 0.6g of liquid exudate/cm² of wound area/24 hours) to moderate (up to 1.1g of liquid exudate/cm² of wound area/24 hours). 1g of exudate is approximately equal to 1ml of exudate. When used over a moderately exuding wound, the size of the wound should generally be no more than 25% of the dressing pad area.

8.2. Application

The dressing should only be applied by a healthcare professional.

1. Remove any excess hair to ensure close approximation of the dressing to the wound. If necessary, irrigate the wound with sterile saline and pat the wound dry.

2. Using a clean technique, peel off the first release handle and place the dressing centrally over the wound to reduce the chance of wound fluid coming into contact with the soft port. Ensure the dressing lies flat to the wound and the surrounding skin. The soft port should be uppermost from the wound (depending on the patient's primary position), placed on intact

skin and not extending over the wound to prevent fluid pooling around the soft port and blocking the NPWT. Remove the other remaining handle(s) and smooth the dressing around the wound to prevent creasing. Reposition if required to ensure border is not creased.

3. Once the dressing is in place, remove the pump and the batteries from the tray.

Caution: Pump contains a MAGNET.


You must keep the pump at least 10cm away from implantable devices such as pacemakers, defibrillators and shunt valves. The direction in which the batteries should be placed is indicated inside the battery compartment. Insert the batteries.

Replace the cover. Following this all four indicators should illuminate for 3 seconds.


4. Join the pump and dressing tubes by twisting together the connectors. The soft tube can be directly connected to the pump if long tubing is not required. Press the orange button to start the application of NPWT. The green OK indicator and orange air leak indicator will start to flash together (indicates pump working to establish NPWT). Depending on the size of the wound, the


pump should take up to 65 seconds to establish NPWT. If after 65 seconds the system has not established NPWT, just the orange air leak indicator will flash. To troubleshoot refer to [section 13](#).

5. If using SECURA No-Sting Barrier Film prior to application of the fixation strips (see Precaution 11) wipe the area surrounding the dressing and allow skin to dry.

6. Apply the fixation strips all the way around the dressing border. Remove top carrier on the strip after each one has been applied. These strips maintain the seal over the wear time of the dressing. In awkward areas, it may be useful to apply the strips to help achieve a seal prior to switching on the pump. Place each strip so that it overlaps the dressing border by approximately 1cm.

Ensure tubing is not twisted or trapped between clothing. Please note that if at any time the fixation strips are removed, the dressing should also be replaced.


If desired, gel patches may be applied in addition to the fixation strips to help achieve or maintain a seal.

8.3. Dressing change

Dressings should only be changed by a healthcare professional.

1. Dressings should only be changed in line with standard wound management guidelines, typically every 3-4 days. At the healthcare professional's discretion a PICO dressing may be left in place for up to 7 days. The orange dressing full indicator on the pump will flash if it detects a full dressing or blocked filter. More frequent dressing changes may be required depending on the level of exudate, condition of the dressing or other patient considerations; e.g. when PICO 7 is used on infected wounds. Additional dressings for PICO 7 are available for purchase separately.

2. Inspect the PICO dressing regularly. If the dressing appears ready for changing (see diagrams A-C), press the orange button and disconnect the dressing from the pump. The fixation strips should be stretched away from the skin and the dressing lifted at one corner and peeled back until it has been fully removed. Apply another dressing as per [section 8.2](#). Application, connect to the pump and press the orange button to reinitiate the therapy.

- (A) Dressing properly positioned and is acceptable to be left in place
 - (B) Dressing requires change – Port may block with fluid
 - (C) Dressing requires change – Absorbent area is full
3. The PICO dressing should be disposed of as clinical waste.
 4. The pump life ends and it automatically stops functioning at 7 days (all the indicators will turn off at this point). The batteries should be removed from the pump; and both batteries and pump disposed of according to local regulations.
 5. For additional information on disposal requirements see: www.mypico.com or speak to your Smith & Nephew representative.

8.4. Use of PICO dressings with fillers

PICO dressings are compatible with standard gauze and foam fillers used in traditional NPWT where this is clinically appropriate – for example on a defect wound. When filler is used, the filler and the PICO dressing should be changed 2 to 3 times a week, according to local clinical protocol and manufacturer's instructions. Gauze should loosely fill to the surface of the wound. Avoid over-packing.

8.5. Use of PICO dressings with non-adherent layers

PICO dressings may be used over the top of a non-adherent layer if required, for example over a skin graft. On infected wounds or wounds at risk of infection, ACTICOAT™ Flex may be used under PICO dressings.

8.6. Use of PICO 7 with compression

PICO 7 is compatible with compression therapy if clinically required; however, the wound should first be accurately assessed to ensure the exudate is of low to moderate volume. Application of active NPWT to a wound may cause production of additional exudate; so following the first application of PICO dressing under a compression system, check the dressing after 2-3 days to see if fluid handling is acceptable or if a dressing change is required and that it is appropriate to continue with compression therapy.

When applying compression therapy over PICO dressings, ensure tubing connecting the PICO 7 pump and dressing runs over the top of the first layer and beneath subsequent layers to avoid it being held against the skin. If using compression hosiery, run tubing outside of stockings. Do not store the pump within layers of compression. For guidance on the correct application of a compression therapy, refer to the relevant manufacturer's instructions.

8.7. Dressing removal

Where a healthcare professional believes it to be clinically appropriate, the patient may be instructed to remove the PICO dressing at the end of therapy. (See [Section 17](#) – Patient Dressing Removal Instructions.)

9. General Use

9.1. Showering and bathing

Light showering is permissible; however, the PICO 7 pump should be disconnected and placed in a safe location where it will not get wet. The PICO dressing should not be exposed to a direct spray or submerged in water. Whilst

disconnected, ensure the end of the tubing attached to the dressing is facing down so that water does not enter the tube.

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11. Safety of PICO 7

When used in accordance with the manufacturer’s instructions, PICO 7 complies with the General Requirements for Safety of Electrical Medical Equipment (IEC 60601-1). PICO 7 is intended for uncontrolled environments e.g. home use (IEC60601-1-11).

12. Electromagnetic compatibility of PICO 7




PICO 7 has been tested and found to comply with the limits for medical devices to IEC 60601-1-2-2014 4th edition. These limits and test levels are intended to provide reasonable safety with regard to electromagnetic disturbances when PICO 7 is used in a typical medical installation and uncontrolled environment like home use.

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
13. Troubleshooting

PICO 7 has visual indicators to let the user know when there is an issue. PICO 7 does not contain audible alerts. The pump should be carried so that it is accessible and the patient/ healthcare professional can check the status routinely in case there is a fault or in case of damage.

Display/Indicator status	Possible cause	Comments/troubleshooting
--------------------------	----------------	--------------------------

 All indicators off	The pump is in standby. The pump has completed its course of NPWT. The batteries have depleted.	NPWT is paused. Press the orange button to restart NPWT. Pressing the orange button will not restart NPWT. Healthcare professional to apply new pump and dressing if further NPWT is required. If the pump has not yet completed its course of NPWT, replace the batteries.
 Green 'OK' and orange 'leak' indicators flash	The pump is working to achieve NPWT but has not reached the intended pressure.	Wait up to 65 seconds. Assess whether NPWT has been established.
 Green 'OK' indicator flashes	System is functioning properly. No issues.	The pump may be heard running occasionally as it maintains the negative pressure. This is normal.
 Green 'OK' and orange 'battery low' indicators flash	System is functioning properly but the batteries are low.	Replace the batteries and press the orange button to restart therapy.
 Orange 'leak' indicator flashes	A high air leak has been detected. NPWT is not being applied. (Note: the pump will automatically try to restart NPWT after 1 hour).	Smooth down the dressing and strips to remove any creases. Press the orange button to restart NPWT. If the air leak remains, the orange 'leak' indicator will flash again after approximately 60 seconds. Ensure that the tube connectors have been twisted together securely.
 Orange 'leak' and orange 'battery low' indicators flash	A high air leak has been detected and the batteries are low. NPWT is not being applied. (Note: the pump will automatically try to restart NPWT after 1 hour).	Resolve the air leak according to instructions above. Also replace the batteries and press the orange button to restart therapy.
 Orange 'dressing full' indicator flashes	Dressing is saturated or filter is blocked. NPWT is not being applied. (Note: the pump will automatically try to restart NPWT after 1 hour).	Healthcare professional to replace the dressing with a new one and press the orange button to restart NPWT.

 Orange 'dressing full' and orange 'battery low' indicators flash	Dressing is saturated or filter is blocked and the batteries are low. NPWT is not being applied. (Note: the pump will automatically try to restart NPWT after 1 hour).	Healthcare professional to replace the dressing with a new one. Also replace the batteries and press the orange button to restart therapy.
 All indicators solidly illuminated	A pump error has been detected. The pump can no longer apply NPWT.	Healthcare professional to apply a new pump and dressing.

14. Cautions

This user manual is not intended as a guarantee or warranty. It is intended only as a guide. For medical questions please consult a physician. The product must be used in accordance with this user manual and all a PICO 7 is packed in the UK with individual components being made in the following countries:

Dressing - UK
 Fixation strips - Origin as marked
 Pump - Origin as marked
 Pump clip – China
 Batteries – Origin as marked

UK Smith & Nephew, Croxley Park, Building 5, Hatters Lane, Watford, Hertfordshire WD18 8YE. Tel: 08000 157 573.

AUSTRALIA Smith & Nephew Pty Ltd.,
 85 Waterloo Road, North Ryde NSW 2113. Tel +61 2 9857 3999.

AUSTRIA Smith & Nephew GmbH, Concorde Business Park C3,
 A-2320 Schwechat. Tel: +43 (0)1 707 91 02.

INDONESIA PT Smith & Nephew Healthcare, Total Building, 5th Floor, J1. Let. Jend. S. Parman Kav
 106A, Jakarta 11440. Tel: 21 5680521.

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Smith & Nephew Medical Limited, 101 Hessle Road, Hull, HU3 2BN, England.

Contact details for other countries can be found at: www.smith-nephew.com

™Trade Marks of Smith & Nephew

©Smith & Nephew Date of issue 12/2018

This product may be covered by one or more US Patents. See www.smith-nephew.com/patents

15. System variants

PICO 7 Single Use Negative Pressure Wound Therapy system

Low to moderately exuding wounds – up to 7-day system

2 dressings, 1 pump & pump clip, secondary fixation strips

Multisite 66802000

Multisite 66802001

10cm x 20cm	66802002
10cm x 30cm	66802003
10cm x 40cm	66802004
15cm x 15cm	66802005
15cm x 20cm	66802006
15cm x 30cm	66802007
20cm x 20cm	66802008
25cm x 25cm	66802009

16. Patient Dressing Removal Instructions

If you have been directed by your healthcare professional to remove the dressing at the end of therapy follow these instructions:

1. Stop the PICO 7 pump by pressing the orange button - All indicators will turn off.
2. Remove the pump from the dressing by untwisting the connectors.
3. Remove the dressing by stretching the fixation strips away from the skin. Lift the dressing at one corner and peel back until it has been fully removed.
4. If pain, reddening, odour or sensitisation occurs, or if wound fluid is present contact your healthcare professional.
5. The PICO dressing and fixation strips should be disposed of as advised by your healthcare professional. The batteries should be removed from the pump; and both batteries and pump disposed of according to local regulations.

**ADDENDUM TO THE INSTRUCTION FOR USE
OF MEDICAL DEVICE**

**PICO 7 Single Use Negative Pressure Wound Therapy
System**

The information contained in this Addendum to the instruction for use is a priority when the product is supplied to the Russian Federation

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1. PICO 7 Single Use Negative Pressure Wound Therapy System contents:
 - Pump – 1pc
 - Multisite dressings 15cm x 20cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tube - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
2. PICO 7 Single Use Negative Pressure Wound Therapy System includes:
 - Pump– 1pc
 - Multisite dressings 20cm x 25cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
3. PICO 7 Single Use Negative Pressure Wound Therapy System contents:
 - Pump– 1pc
 - dressings 10cm x 20cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
4. PICO 7 Single Use Negative Pressure Wound Therapy System contents:
 - Pump– 1pc
 - dressings 10cm x 30cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
5. PICO 7 Single Use Negative Pressure Wound Therapy System contents:
 - Pump– 1pc
 - dressings 10cm x 40cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.
6. PICO 7 Single Use Negative Pressure Wound Therapy System contents:
 - Pump – 1pc
 - dressings 15cm x 15cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - extension tubing - 1 pc.
 - AA batteries - 2 pcs.
 - instruction for use - 8 pcs.

7. PICO 7 Single Use Negative Pressure Wound Therapy System contents:

- Pump – 1pc
- dressings 15cm x 20cm - 2 pcs.
- Secondary fixation strips - 12 pcs.
- extension tubing - 1 pc.
- AA batteries - 2 pcs.
- instruction for use - 8 pcs.

8. PICO 7 Single Use Negative Pressure Wound Therapy System contents:

- Pump – 1pc
- dressings 15cm x 30cm - 2 pcs.
- Secondary fixation strips - 12 pcs.
- extension tubing - 1 pc.
- AA batteries - 2 pcs.
- instruction for use - 8 pcs.

9. PICO 7 Single Use Negative Pressure Wound Therapy System contents:

- Pump – 1pc
- dressings 20cm x 20cm - 2 pcs.
- Secondary fixation strips - 12 pcs.
- extension tubing - 1 pc.
- AA batteries - 2 pcs.
- instruction for use - 8 pcs.

10. PICO 7 Single Use Negative Pressure Wound Therapy System contents:

- Pump – 1pc
- dressings 20cm x 25cm - 2 pcs.
- Secondary fixation strips - 12 pcs.
- extension tubing - 1 pc.
- AA batteries - 2 pcs.
- instruction for use - 8 pcs.

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Dressing: Translucent, one-side siliconized high-density polyethylene film, wound contacting adhesive, perforated EU30 polyurethane, anti-blocking concentrate, A8 acrylic adhesive, polyester, super absorbent, K5 acrylic adhesive, EU33 polyurethane film, anti-blocking agent

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Article/catalogue number:

Name: PICO 7 Single Use Negative Pressure Wound Therapy System

Size*: 10x20 /10x30 /10x40 /15x15 /15x20 /15x30/ 20x20/ 25x25/ Multisite 15x20/ Multisite 20 x25

Pieces per package, pcs.: 1 pc.

Date of manufacture, best before, lot/batch number: See the original packaging.

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Authorized representative in the Russian Federation: Smith & Nephew LLC, 105120, Russia, Moscow, 2nd Syromyatnichesky Lane, 1, 9th floor, office 1, room 1. tel/fax: + 7 495 755 55 03

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






*Only one size corresponding to original article/catalogue number will be specified

3.2 Data on labeling

Symbol/text	Meaning/translation
	PICO 7
	First remove protective liner from the fixation strip
	Then remove top carrier on the fixation strip
OK	"OK" indicator - the pump is functioning properly
	"Air leak" indicator
	"Dressing full" indicator
	"Battery low" indicator
	Manufacturer
	Date of manufacture
STERILE EO	The product is EtO sterilized
	Single use. Do not reuse

	Warning!
	Incompatible with magnetic resonance - keep away from magnetic resonance equipment
	International classification
	EC: Do not dispose of as household waste
	EU Mark of Compliance
	Follow the instructions for use
	The service life of the device is up to 7 days
	Warning! Magnetic field! The pump contains a magnet: Keep the pump at least 10cm/4" away from the implanted medical devices such as pacemakers, defibrillators, and shunt valves
	Protection degree of the housing of electrical equipment: the device is protected from finger penetration and cannot be damaged vertically or almost vertically by dripping water
	Use by medical prescription only
	Classification of equipment: applied part type BF, defibrillation proof
	Start / Pause / Resume treatment
	1 pump included
	The service life of the device is up to 7 days
	Batch number

	Dressing size Length 25cm (9.8 inch) - 20cm (7.9 inch) Width 20cm (7.9 inch) - 15cm (5.9 inch) (in accordance with the design variant))
	Shelf life
	Catalog number
	Unique product identifier
	International labeling and accounting code
	Barcode
	Recyclable packaging material
	European authorized representative
	Do not use if package is opened or damaged
	Keep away from direct sunlight
	Storage temperature
	Protect from moisture
	Follow the instructions for use
	Packaging contains 6 fixation strips
	Open by folding back the paper edge

 x1	Packaging contains 1 clip
 x1	Packaging contains 1 extension tubing
 x 2	Packaging contains AA E91 battery
	Schematic instructions for use
	PICO dressing included
	Schematic instructions for use

4. MD DISPOSAL OR DESTRUCTION PROCEDURE AND CONDITIONS

Disposal of packages / contents must be carried out in a safe and authorized manner in compliance with the requirements of regional / national / international regulations regarding the clinical waste handling procedure.

The used medical device disposal or destruction procedure - disposal is carried out in accordance with the applicable legislation of the Russian Federation. Waste class according to SanPiN 2.1.7.2790-10 "B"

The new medical device disposal or destruction procedure

- disposal is carried out in accordance with the applicable legislation of the Russian Federation. Waste class according to SanPiN 2.1.7.2790-10 "A"

5. DESCRIPTION OF STERILIZATION METHOD

Sterile. EtO sterilized.

6. LIST OF MATERIALS OF ANIMAL AND/OR HUMAN ORIGIN

Not applicable.

7. CONTENT OF MEDICAL SUBSTANCES

The section does not apply to the considered products since the Smith & Nephew Medical Limited products do not contain medical substances, and no medical substance is used in the manufacturing process.

8. LIST OF STANDARDS

Smith & Nephew Medical Ltd. declares that PICO Single Use Negative Pressure Wound Therapy System, variations meets the requirements of Directive 93/42/EEC, a full list of international requirements is available on request.
9. MD MAINTENANCE AND REPAIR REQUIREMENTS
Not applicable.
10. MD PERFORMANCE AND SAFETY CHECKING
Not applicable.
11. WARRANTY AND CLAIMS
Claimee: Smith & Nephew LLC, 105120, Russia, Moscow, 2nd Syromyatnichesky Lane, 1, 9th floor, office 1, room 1. tel/fax: + 7 (495) 755 55 03 WARRANTY: The manufacturer guarantees compliance with the declared characteristics, the period of preservation of sterility, the safety and integrity of the packaging subject to proper transportation, operation, and storage of the product for 2 years from the date of manufacture/sterilization of the product. Shelf life - 24 months from the date of manufacture / sterilization indicated on the package, provided that it is stored in the original packaging, in a dry, well-ventilated place, protected from direct sunlight, at a temperature not exceeding 25°C. Warranty for PICO Single Use Negative Pressure Wound Therapy System, variations applies to materials, functioning and lack of production defects, provided that the device is used for only one patient. Reuse is prohibited.
12. NAME AND LEGAL ADDRESS OF THE MANUFACTURER, AUTHORIZED REPRESENTATIVE OF MD
MANUFACTURER: Smith & Nephew Medical Limited, 101 Hessle Road, Hull, HU3 2BN, United Kingdom PLACE OF MANUFACTURE: 1. Smith & Nephew Medical Limited, 101 Hessle Road, Hull, HU3 2BN, United Kingdom 2. Flextronics Technology (Shenzhen) Co., Ltd. 3-6/F, Building 1 and 6/F, Building 4, Yusheng Industrial Park, 467 Xixiang Section, National Highway 107, Xixiang, Baoan District, 518126 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA. 3. Avery Dennison Belgie BVBA, Tieblokkenlaan 1, Turnhout B-2300, Belgium. 4. Smith & Nephew Medical (Suzhou) Limited, No. 12, Wuxiang Road Comprehensive Free Zone, West Zone Suzhou Industrial Park Suzhou Jiangsu 215021 China. AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA: Smith & Nephew LLC, 105120, Russia, Moscow, 2nd Syromyatnichesky Lane, 1, 9th floor, office 1, room 1. tel/fax: + 7 (495) 755 55 03

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PICO 7Y Single Use Negative Pressure Wound Therapy System

Healthcare professional user manual

PICO™ 7Y is intended for use by or on the direction of a trained and licensed physician in accordance with these instructions for use.


This user manual contains information specific for use by a healthcare professional and is not appropriate for use by patients and caregivers.

Dressings should only be applied, changed and removed by a healthcare professional.

1. Description

The PICO 7Y pump maintains negative pressure wound therapy at 80 mmHg (nominal) to two wound surfaces simultaneously. Exudate is managed by the dressings through a combination of absorption and evaporation of moisture through the outer film.

The PICO 7Y system is intended to be used for up to 7 days on low exuding wounds. For moderate exuding wounds the system is intended to be used for up to 4 days. For 7 days' use on moderate exuding wounds, additional dressings will be required (available for purchase separately).

A	1x PICO 7Y Pump
B	Y-connector
C	OK Indicator
D	Leak Indicator
E	Check Dressing Indicator
F	Battery Low Indicator
G	Pump Clip
H	2x AA Lithium Batteries
I	2x Large Multisite Dressings
J	Soft Port
K	Fixation Strips
L	2x Extension Tubing Pump

M	Pump & Tubing Connectors

2. Indications for use

PICO 7Y is indicated for patients who would benefit from a suction device (Negative Pressure Wound Therapy) as it may promote wound healing via removal of low to moderate levels of exudate and infectious materials. Appropriate wound types include:

- Chronic
- Acute
- Traumatic
- Subacute and dehiscent wounds
- Partial-thickness burns
- Ulcers (such as diabetic or pressure)
- Flaps and grafts
- Surgically closed incision sites

PICO 7Y systems are suitable for use both in a hospital and homecare setting

3. Contraindications

PICO 7Y is contraindicated for:

- Patients with malignancy in the wound bed or margins of the wound (except in palliative care to enhance quality of life).
- Previously confirmed and untreated osteomyelitis.
- Non-enteric and unexplored fistulas.
- Necrotic tissue with eschar present.
- Exposed arteries, veins, nerves or organs.
- Exposed anastomotic sites.

4. Important information

Pump Placement Warning

The PICO 7Y pump contains a **MAGNET**. Keep the PICO 7Y pump at least 10cm away from other medical devices at all times. Failure to do so can cause the other medical device to fail which can result in serious harm including death.

See Warning number 1. Magnet Warning

For more information on electromagnetic immunity and electromagnetic emissions see: www.mypico.com or ask your Smith & Nephew representative for a hardcopy.

5. Warnings

1. The PICO 7Y pump contains a **MAGNET** that can cause other medical devices in close proximity to fail, leading to serious harm including death. The PICO 7Y pump

must be positioned at least 10cm away from other medical devices that could be affected by magnetic interference. These include but are not limited to:

- Implantable Cardioverter - Defibrillator (ICD)
- Pacemakers
- Insulin Pumps
- Shunt Valves
- Neurostimulators
- Cochlear Implants

THIS WARNING APPLIES AT ALL TIMES TO ALL USERS.

This applies to both patients and caregivers.

You must keep the PICO 7Y pump at least 10cm away from other devices:

- If you have an electronic medical device and are helping take care of somebody else using the PICO 7Y system.
 - If the patient is wearing the PICO 7Y pump in a public area where they may come in close contact with someone else who has an electronic medical device.
2. Certain patients are at high risk of bleeding complications which, if uncontrolled, could be fatal. Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately disconnect pump, leave dressings in place, take appropriate measures to stop bleeding and seek immediate medical assistance.
 3. Haemostasis must be achieved before applying the dressings, although the use of anticoagulants does not deem a patient inappropriate for treatment with PICO 7Y. Patients suffering from difficult haemostasis or who are receiving anticoagulant therapy have an increased risk of bleeding. During therapy, avoid using haemostatic products that may increase the risk of bleeding, if disrupted. Frequent assessment must be maintained throughout the therapy.
 4. ssPinICgOs™ should only be applied, changed and removed by a healthcare professional.
 5. PICO 7Y is unsuitable for use in areas where there is danger of explosion (e.g. hyperbaric oxygen unit).
 6. PICO 7Y is unsuitable for use in the presence of flammable anaesthetic mixture with oxygen or nitrous oxide.
 7. Sharp edges or bone fragments in a wound must be covered or removed prior to using PICO 7Y due to risk of puncturing organs or blood vessels while under Negative Pressure Wound Therapy.
 8. Each PICO dressing (including Multisite) must be used to dress one wound only.
 9. At all times care should be taken to ensure that the pump and tubing and connectors do not:
 - Lie in a position where they could cause pressure damage to the patient.
 - Trail across the floor where they could present a trip hazard or become contaminated.
 - Present a risk of strangulation or a tourniquet to patients.
 - Rest on or pass over a source of heat.
 - Become twisted or trapped under clothing or bandages so that the air-path delivering Negative Pressure is blocked.
 10. MR unsafe. The PICO 7Y pump is not MRI compatible. Remove the PICO 7Y pump from the dressings before entering the MRI suite.
 11. PICO 7Y has not been studied on paediatric patients. Patient size and weight should be considered when prescribing this therapy.

12. The system contains small parts which could represent a choking hazard. Keep out of the reach of children.
13. Keep PICO 7Y away from pets, pests and other animals that could damage the PICO 7Y device.

6. Precautions

1. The PICO 7Y pump contains a magnet. See Warning number 1. Magnet Warning.
2. **Precautions should be taken in the following types of patients who are at high risk of bleeding complications:**
 - Receiving anticoagulant therapy or platelet aggregation inhibitors or actively bleeding.
 - Having weakened or friable blood vessels or organs in or around the wounds as a result of, but not limited to; anastomoses, infection, trauma or radiation.
 - Suffering from difficult wound haemostasis.
 - Untreated for malnutrition.
 - Non-compliant or combative.
 - Suffering from wounds in close proximity to blood vessels or delicate fascia.
3. Monitor for pain, reddening, odour, sensitisation or a sudden change in the volume or colour of wound fluid during use.
4. Where PICO 7Y is used to bolster skin grafts, it is important to visually inspect the system regularly, especially in the first week of treatment to ensure that Negative Pressure Wound Therapy is continually applied and a seal is maintained.
5. Where PICO dressings are used on infected wounds, more frequent dressing changes may be required. Regular monitoring of the wounds should be maintained to check for signs of infection.
6. If deemed clinically appropriate, care should be taken that the application of a circumferential dressing or the use of Negative Pressure Wound Therapy on ischaemic limbs does not compromise circulation.
7. PICO 7Y does not contain audible alerts. The pump should be carried so that it is accessible and the patient/healthcare professional can check the status routinely.
8. Although PICO dressings can be used under clothing/bedding, it is important that occlusive materials e.g. film dressings, are not applied over the pad area of the dressings as this will impair the intended evaporation of moisture through their outer layers.
9. The PICO dressings should not be covered by rigid immobilisation devices or casts which might apply excessive pressure and cause tissue injury at the wound sites, especially where the tubing enters the dressings.
10. Prolonged placement of rigid or opaque materials over the PICO dressings may prevent the regular inspection and assessment of the wounds, and disrupt scheduled or required dressing changes.
11. Where PICO dressings are used on patients with fragile skin, a skin protectant such as SECURA™ No-Sting Barrier Film should be used on areas of skin where fixation strips are to be applied. Inappropriate use or repeated application of fixation strips may otherwise result in skin stripping.
12. Do not use PICO dressings with oil based products such as petrolatum as it may compromise establishing an effective seal.
13. The use of Negative Pressure Wound Therapy presents a risk of tissue ingrowth into foam when this is used as wound filler. When using foam filler with PICO 7Y, tissue ingrowth may be reduced by using a non-adherent wound contact layer or by increasing the frequency of dressing changes.
14. PICO 7Y may be used in conjunction with surgical drains provided the dressings are not placed over tubing where it exits the skin. Any surgical drains should be routed under the skin away from the edge of the dressings and function independently of the PICO 7Y system.
15. The pump must be protected from sources of fluid e.g. from incontinence or

spillages. Discontinue PICO 7Y use if fluid ingress is observed.

16. When showering, the PICO 7Y pump should be disconnected from the dressings. Whilst disconnected, ensure the end of the tubing attached to the dressings is facing down so that water does not enter the tubes.
17. Do not take the pump apart.
18. The dressings supplied should only be used with PICO 7Y pumps.
19. Do not alter or cut tubing configuration or pull on the tubing or soft ports.
20. Do not cut the PICO dressing pad as this may lead to loss of Negative Pressure Wound Therapy application.
21. When applying dressings next to one another, ensure the dressing borders do not overlap.
22. Always ensure that the PICO dressing is positioned centrally over the wound. The soft port should be positioned uppermost on intact skin and not extend over the wound so that the risk of fluid collecting around the soft port and potentially blocking the air-path is minimised.
23. CT scans and x-ray have the potential to interfere with some electronic medical devices. Where possible, move the pump out of the x-ray or scanner range. If the pump has been taken into the CT scan or x-ray range, check that the system is functioning correctly following the procedure.
24. The potential for electromagnetic interference in all environments cannot be eliminated. Use caution if PICO is near electronic equipment such as RFID (Radio Frequency Identification) readers, anti-theft equipment or metal detectors.
25. The PICO 7Y system is single use only. Use of any part of this system on more than one patient may result in cross contamination that may lead to infection.
26. High temperatures and humidity may reduce wear time of PICO dressings.
27. The PICO 7Y system is intended for use both in a hospital and homecare setting. The system can also be used in aircraft, car, train, and boat transportation. Special care must be taken regarding pump positioning when in close proximity to other people (see magnet warning).
28. During transport there is a potential for radio frequency interference that could affect PICO 7Y performance. If the PICO 7Y pump malfunctions, replace the system.

7. Adverse reaction

Excessive bleeding is a serious risk associated with the application of suction to wounds which may result in death or serious injury. Careful patient selection, in view of the above stated contraindications, warnings and precautions is essential. Carefully monitor the wounds and dressings for any evidence of a change in the blood loss status of the patient or any sudden or abrupt changes in the volume or the colour of exudate.

8. Instructions for use (healthcare professional only)

8.1. Guidance on wound suitability

PICO dressings should be used on wounds which fit comfortably within the area of the pad (generally no more than 25% of the dressing pad area), observing precautions on soft port positioning (on intact skin and not extending over the wound). PICO Multisite dressings are designed to enhance conformability when dressing awkward anatomical areas. Each PICO Multisite dressing must be used to dress one wound only. As a guide:

- Depth - Wounds greater than 0.5cm in depth are likely to require a foam or gauze Negative Pressure Wound Therapy filler to ensure adequate treatment of all the wound surfaces. Wounds treated with the PICO 7Y system should generally be no more than 2cm in depth.
- Exudate – PICO 7Y is intended for use on wounds where the level of exudate is low (up to 0.6g of liquid exudate/cm² of wound area/24 hours) to moderate (up to 1.1g of liquid exudate/cm² of wound area/24 hours). 1g of exudate is approximately equal to 1ml of

exudate.

8.2. Use of PICO dressings with fillers

PICO dressings are compatible with standard gauze and foam fillers used in traditional Negative Pressure Wound Therapy where this is clinically appropriate – for example on a defect wound. When filler is used, the filler and the PICO dressing should be changed 2 to 3 times a week, according to local clinical protocol and manufacturer's instructions. Gauze should loosely fill to the surface of the wound. Avoid over-packing.

The use of Negative Pressure Wound Therapy presents a risk of tissue ingrowth into foam when this is used as wound filler. When using foam filler with PICO 7Y, tissue ingrowth may be reduced by using a non-adherent wound contact layer or by increasing the frequency of dressing changes.

8.3. Use of PICO dressings with non-adherent layers

PICO dressings may be used over the top of a non-adherent layer if required. On infected wounds or wounds at risk of infection, ACTICOAT™ Flex 3 and ACTICOAT™ Flex 7 may be used under PICO dressings.

8.4. Use of PICO dressings with compression

PICO 7Y is compatible with compression therapy if clinically required; however, the wounds should first be accurately assessed to ensure the exudate is of low to moderate volume. Application of active Negative Pressure Wound Therapy to a wound may cause production of additional exudate; so following the first application of PICO dressings under a compression system, check the dressings after 2-3 days to see if fluid handling is acceptable or if a dressing change is required and that it is appropriate to continue the use of PICO with compression therapy. When applying compression therapy over PICO dressings, ensure tubing connecting the PICO 7Y pump and dressings runs over the top of the first layer and beneath subsequent layers to avoid it being held against the skin. If using compression hosiery, run tubing outside of stockings. Do not store the pump within layers of compression. For guidance on the correct application of a compression therapy, refer to the relevant manufacturer's instructions.


8.5. Application

1. Dressings should only be applied, changed or removed by healthcare professionals.

If part of your clinical guidelines remove any excess hair to ensure close approximation of the dressing to the wound. If necessary, irrigate the wound with sterile saline and pat the wound dry with sterile gauze.

2. Using a clean technique, peel off the first piece of backing plastic. If possible ensure that the soft port is positioned higher than the wound (depending on the patient's primary position). Ensure that the soft port is placed over intact skin and not extending over the wound. These steps will prevent fluid blocking the air-path delivering negative pressure. Ensure the dressing lies flat to the wound and the surrounding skin. If the soft port is not placed higher than the wound, the system will function as normal but the wear time of the dressing may be reduced.

3. Remove the remaining backing plastic and smooth the dressing around the wound to prevent creasing. Reposition if required to ensure border is not creased.


4. Repeat steps 2 - 3 for the second dressing. Ensure the dressings do not overlap.

The PICO 7Y pump contains a MAGNET. Keep the PICO 7Y pump at least 10cm away from other medical devices at all times. Failure to do so can cause the other medical device to fail which can result in serious harm including death.


See Warning number 1. Magnet Warning.


5. Once the dressings are in place, remove the pump and the batteries from the tray. The direction in which the batteries should be placed is indicated inside the battery compartment. Insert the batteries. Following this all four indicators should illuminate for 3 seconds. Replace the cover
6. Join the pump and dressing by twisting together the connectors. Extension tubes are included or the soft port tubes can be directly connected to the pump. If extension tubes are not used for initial application, they should be retained for dressing change.

7. Press the Orange Button to start the application of the therapy. The green OK indicator and orange leak indicator will start to flash together (indicates pump working to establish therapy). Depending on the size of the wounds, the pump can take up to 100 seconds to establish therapy. Once therapy is established just the OK indicator will flash. If after 100 seconds the system has not established therapy, just the orange air leak indicator will flash.

8. If using SECURA No-Sting Barrier Film prior to application of the fixation strips (see Precaution 11) wipe the area surrounding the dressings and allow skin to dry.
9. Apply the fixation strips all the way around the dressing border. Remove the printed top layer on the strip after each one has been applied. These strips maintain the seal over the wear time of the dressings. In awkward areas, it may be useful 100 sec to apply the strips to help achieve a seal prior to switching on the pump. Place each strip so that it overlaps the dressing borders by approximately 1cm.


10. Ensure tubing is not twisted or trapped between clothing. Please note that if at any time the fixation strips are removed, the dressing(s) should also be replaced. If desired, gel patches may be applied in addition to the fixation strips to help achieve or maintain a seal.
11. If extension tubes have not been used for initial application, remember to retain for dressing change. A pump clip is supplied to attach the pump to clothing if required.

9. Check dressing indicator

When this indicator flashes, the condition of the dressings needs to be checked. The indicator will flash routinely every 24 hours as a reminder to check the dressings.

(See **Section 11. Dressing change**)


If no dressing change is required, refresh the device by pressing the Orange Button TWICE.

10. Dressing change

Dressings should only be applied, changed or removed by healthcare professionals.

Dressings should only be changed in line with standard wound management guidelines, typically every 3-4 days. At the healthcare professional's discretion a PICO™ dressing may be left in place for up to 7 days. More frequent dressing changes may be required depending on the level of exudate, condition of the wound(s) or other patient considerations, e.g. when PICO 7Y is used on infected wounds. Additional dressings for use with the PICO 7Y pump are available for purchase separately. Inspect the PICO dressings regularly. If either dressing appears ready for changing, remove the dressing as per Section 13 and apply another dressing as per Section 9.5.

A – OK, Acceptable to be left in place

B – Change, Soft port may block with fluid

C – Change, Absorbent area is full

Based on dressing change frequency, further dressings may be required, see Section 20. System variants. PICO Multipack does not contain extension tubes. Extension tubes must be retained for dressing change.

11. Replacing the batteries

If the Battery Low indicator flashes, the batteries in PICO will require changing. To change the batteries:

1. Remove the back cover from PICO to access the battery compartment.
2. Remove the used batteries and dispose of in accordance with local regulations.
3. Insert new batteries in the orientation indicated inside the battery compartment.
4. Replace the back cover, all four indicators should illuminate for 3 seconds.

5. Press the Orange Button to resume therapy.

12. Dressing removal

Dressings should only be applied, changed or removed by healthcare professionals.

1. Stop the PICO 7Y pump by pressing the Orange Button. All indicators will turn off.
2. Remove the pump from the dressings by untwisting the connectors.
3. Remove the PICO dressings by stretching the fixation strips away from the skin. Lift the dressings at one corner and peel back until it has been fully removed.
4. The PICO dressings and fixation strips should be disposed of as clinical waste in accordance with local protocol. The batteries should be removed from the pump, and both batteries and pump disposed of according to local regulations.
5. For additional information on disposal requirements speak to your Smith & Nephew representative.

13. Cleaning, showering and bathing




Adherence to clinical directives concerning hygiene is of prime importance. The pump may be wiped clean with a damp cloth using soapy water or a weak disinfectant solution.



13.1. Showering and bathing



Light showering is permissible; however, the PICO 7Y pump should be disconnected and placed in a safe location where it will not get wet. The PICO dressings should not be exposed to a direct spray or submerged in water. Whilst disconnected, ensure the ends of the dressing tubes are facing down so that water does not enter the tubes.

14. Troubleshooting

PICO 7Y has visual indicators to let the user know when there is an issue. PICO 7Y does not contain audible alerts. The pump should be carried so that it is accessible and the patient/healthcare professional can check the status routinely in case there is a fault or in case of damage.

Indicator status	Possible cause	Comments/troubleshooting
 All indicators off	The pump is in standby. The pump has completed its course of therapy. The batteries have depleted.	Therapy has been paused. Press the Orange Button to resume therapy. Pressing the Orange Button will not restart therapy. Healthcare professional to apply new pump and dressings if further therapy is required. If the pump has not yet completed its course of therapy, replace the batteries. See section 12. Replacing the batteries.
 Green 'OK' and orange 'Leak' indicators flash	The pump is working to achieve therapy but has not reached the intended pressure.	 Depending on the size of the wound the pump

		may take up to 100 seconds to establish therapy. Pump will then flash green 'OK' indicator OR orange 'Leak' indicator.
 Green 'OK' indicator flashes	Pump is functioning. Dressing may be full. Refer to section 11. Dressing change. Pump may be heard running occasionally as it maintains negative pressure.	Dressing may be full.  • A OK- Acceptable to be left in place • B Change required • C Absorbent area is full Dressings should only be applied, changed or removed by healthcare professionals. Dressings should be regularly monitored.
 Green 'OK' and orange 'Battery Low' indicators flash	System is functioning properly but the batteries are low.	Replace the batteries and press the Orange Button to restart therapy. See section 12. Replacing the batteries.
 Orange 'Leak' indicator flashes	A high air leak has been detected. Therapy is not being applied. (Note: the pump will automatically try to restart therapy after 1 hour).	1. Smooth down the dressings and strips to remove any creases. 2. Ensure that the tube connectors have been twisted together securely. 3. Press the Orange Button to restart the pump. 4. If the air leak remains the orange 'Leak' indicator will flash again after approximately 100 seconds (depending on the size of the wound). One of the dressings may have a leak. To identify which of the dressings has a leak, follow these steps: 1. Untwist one of the tubes from the Y-connector. 2. Cover the open branch of the Y-connector with a finger and press

		Orange Button to resume therapy.  3. If the green 'OK' indicator flashes, connect the dressing and repeat steps 1-3 for remaining dressing. However, if the leak remains, the orange 'Leak' indicator will flash again after approximately 100 seconds. If the leak persists the dressing may need changing. Refer to Section 11 Dressing change.
 Orange 'Check Dressing' indicator flashes	Pump is functioning properly however check dressings. See section 10. Check dressing indicator.	Check the dressings. See section 10. Check dressing indicator.
 Orange 'Leak' and orange 'Battery LoW' indicators flash	A high air leak has been detected and the batteries are low. Therapy is not being applied.	Replace batteries. See section 12. Replacing the batteries. Then troubleshoot air leak as above..
 All indicators solidly illuminated	A pump error has been detected. The pump can no longer apply therapy	Apply new pump and dressings.

15. Cautions

This user manual is not intended as a guarantee or warranty. It is intended only as a guide. For medical questions please consult a physician. The product must be used in accordance with this user manual and all applicable labelling.

16. Specifications

Pump Dimensions	Height (mm) 71,58 ± 0,17 Front side width (mm) 11,8 + 0,15/-0,05 Back side width (mm) 9,1 ± 0,11 Length (mm) 65.58 ± 0.15
Weight	<100g

Battery Type	2 x AA 1.5V (FR6)
Ingress Protection	IP22
Maximum Vacuum	13.34 kPa (100 mmHg)
Mode of Operation	Continuous
Storage/Transport	5°C to 25°C (-25°C to 5°C allowable for up to 7 days), 10 – 75% relative humidity, 700 to 1060 mbar atmospheric pressure
Operating Environment	5°C – 40°C, 10 – 95% relative humidity, 700 to 1060 mbar atmospheric pressure

17. Safety

When used in accordance with the manufacturer's instructions, PICO™ 7Y complies with the General Requirements for Safety of Electrical Medical Equipment (IEC 60601-1). PICO 7Y is intended for uncontrolled environments e.g. home use (IEC 60601-1-11).

18. Electromagnetic compatibility of PICO 7Y

PICO 7Y has been tested and found to comply with the requirements for medical devices to IEC 60601-1-2. These limits and test levels are intended to provide reasonable safety with regard to electromagnetic disturbances when PICO 7Y is used in a typical medical installation and uncontrolled environment like home use.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. For additional information on electromagnetic immunity and electromagnetic emissions see: www.mypico.com or ask your Smith & Nephew representative for a hardcopy.

19. System variants

PICO 7Y Single Use Negative Pressure Wound Therapy system Low to moderately exuding wounds – up to 7 day system. 2 x large Multisite dressings, 1 pump & pump clip, fixation strips, 2 extension tubes.

PICO 7Y 66802031

PICO 7Y is packed in the UK with individual components being made in the following countries:

Dressing – UK,
Fixation strips – Origin as marked,
Pump – China,
Pump clip – China,
Batteries – Origin as marked,
Tubing – UK

Smith & Nephew Medical Limited, 101 Hessle Road, Hull, HU3 2BN, United Kingdom.

Contact details for other countries can be found at www.smith-nephew.com

Trade Marks of www.smith-nephew.com

©Smith & Nephew Date of issue 12/2018

This product may be covered by one or more US patents.

See: www.smith-nephew.com/patents

**ADDENDUM TO THE USER MANUAL
OF MEDICAL DEVICE**

**PICO 7Y Single Use Negative Pressure Wound Therapy
System**

The information contained in this Addendum to the instruction for use is a priority when the product is delivered to the Russian Federation

1. MD NAME, MD MODELS OR VARIATIONS

PICO Single Use Negative Pressure Wound Therapy System, variations:

1. PICO 7Y Single Use Negative Pressure Wound Therapy System contents:
 - Pump- 1pcs.
 - Pump clip – 1 pcs.
 - Multisite dressings 20cm x 25cm - 2 pcs.
 - Secondary fixation strips - 12 pcs.
 - Extension tubing - 2 pcs.
 - AA batteries - 2 pcs.
 - Instruction for use - 3 pcs.

2. LIST AND DESCRIPTION OF MD MATERIALS

Pump: Bisphenol A-based polycarbonate, colourant, black print ink, grey print ink, polyester, silicone rubber, polypropylene, matte polyester protective film, stainless steel

Extension tubing: PVC, polypropylene

Dressing: Translucent, one-side siliconized high-density polyethylene film, wound contact adhesive, perforated EU30 polyurethane, anti-blocking concentrate, A8 acrylic adhesive, polyester, super absorbent, K5 acrylic adhesive, EU33 polyurethane film, anti-blocking agent

Soft port: EU94 high-strength polyurethane film, anti-blocking agent, UV curable adhesive, nylon acrylic copolymer, polyester, curable adhesive, PVC, aromatic polyurethane thermoplastic

Fixation strips: Polyethylene woodless bleached paper with insoluble silicone coating, orange print ink, A8 acrylic adhesive, EU45 polyurethane film, anti-blocking agent, polyester film with silicone and matte polyethylene coating, orange print ink, polypropylene tape

3. PACKAGING AND LABELING

3.1. Draft of labeling in Russian

The medical device package labeling in the Russian Federation contains the following information in Russian:

Article/catalogue number:

Name: PICO 7Y Single Use Negative Pressure Wound Therapy System

Size: Multisite 20 x25

Pieces per package, pcs.: 1 pc.

Date of manufacture, best before, lot/ batch number: See the original packaging.

Name and address of manufacturer: Smith & Nephew Medical Ltd, 101 Hessle Road Hull, HU3 2BN United Kingdom

Authorized representative in the Russian Federation: Smith & Nephew LLC, 105120, Russia, Moscow, 2nd

Syromyatnichesky Lane, 1, 9th floor, office 1, room 1. tel/fax: + 7 495 755 55 03

Sterilization method: Ethylene oxide sterilization. Sterile.

Registration certificate No.:

Intended use, method of application, contraindications: See Instructions for Use

Storage conditions: store in a dry, dark place, at a temperature of +5°C to + 25°C




Precautions:


Symbol of conformity to Russian standards (PCT mark)

Instruction for Use: <http://www.smith-nephew.com/about-us/where-we-operate/europe/russia/ifus/>


3.2 Data on labeling

Symbol/text	Meaning/translation
	PICO 7Y
	First remove protective liner from the fixation strip
	Then remove top carrier on the fixation strip
OK	"OK" indicator - the pump is functioning properly
	"Air leak" indicator
	"Dressing full" indicator
	"Battery low" indicator
	Manufacturer
	Date of manufacture
STERILE EO	The product is EtO sterilized
	Single use. Do not reuse
	Warning!
	Incompatible with magnetic resonance - keep away from magnetic resonance equipment
	International classification
	EC: Do not dispose of as household waste
CE 0086	EU Mark of Compliance

	Follow the instructions for use
	The service life of the device is up to 7 days
	Warning! The PICO 7Y pump contains a MAGNET. Always keep the PICO 7Y pump at least 10cm/4" away from the other medical devices. Failure to do so can cause the other medical device to fail, which can result in serious harm including death.
	Protection degree of the housing of electrical equipment: the device is protected from finger penetration and cannot be damaged vertically or almost vertically by dripping water
	"Medical worker" symbol
	Relative humidity limits
	Atmospheric pressure limits
	Use by medical prescription only
	Classification of equipment: applied part type BF, defibrillation proof
	Start / Pause / Resume treatment
	1 pump included
	The service life of the device is up to 7 days
	Batch number
	Catalog number

	Dressing size Length 25cm (9.8 inch) - 20cm (7.9 inch) Width 20cm (7.9 inch) - 15cm (5.9 inch)
	Shelf life
UDI	Unique product identifier
GTIN	International labeling and accounting code
	Barcode
	Recyclable packaging material
EC REP	European authorized representative
	Instruction for use included
	Do not use if package is opened or damaged
	Keep away from direct sunlight
	Storage temperature
	Protect from moisture
	Follow the instructions for use
	Packaging contains 6 fixation strips
	Open by folding back the paper edge

 x1	Packaging contains 1 clip
 x2	Packaging contains 2 extension tubes
"AA"/L91  x 2	Packaging contains AA L91 battery
	Schematic instructions for use
	Low noise pump (below 35 dB)
	High air leakage resistance
	"Dressing full" indicator
	Separation layer
	Evaporates 80% of the exudate
	Intended for low to moderate wound exudate
	Silicone adhesive

	Water-proof dressing
	Schematic instructions for use
4. MD DISPOSAL OR DESTRUCTION PROCEDURE AND CONDITIONS	
Disposal of packages / contents must be carried out in a safe and authorized manner in compliance with the requirements of regional / national / international regulations regarding the clinical waste handling procedure.	
The used medical device disposal or destruction procedure - disposal is carried out in accordance with the applicable legislation of the Russian Federation. Waste class according to SanPiN 2.1.7.2790-10 "B"	
The new medical device disposal or destruction procedure - disposal is carried out in accordance with the applicable legislation of the Russian Federation. Waste class according to SanPiN 2.1.7.2790-10 "A"	
5. DESCRIPTION OF STERILIZATION METHOD	
Sterile. EtO sterilized.	
6. LIST OF MATERIALS OF ANIMAL AND/OR HUMAN ORIGIN	
Not applicable.	
7. CONTENT OF MEDICAL SUBSTANCES	
The section does not apply to the considered products since the Smith & Nephew Medical Limited products do not contain medical substances, and no medical substance is used in the manufacturing process.	
8. LIST OF STANDARDS	
Smith & Nephew Medical Ltd. declares that PICO Single Use Negative Pressure Wound Therapy System, variations meets the requirements of Directive 93/42/EEC, a full list of international requirements is available on request.	
9. MD MAINTENANCE AND REPAIR REQUIREMENTS	
Not applicable.	
10. MD PERFORMANCE AND SAFETY CHECKING	
Not applicable.	
11. WARRANTY AND CLAIMS	
Claimee: Smith & Nephew LLC, 105120, Russia, Moscow, 2nd Syromyatnichesky Lane, 1, 9th floor, office 1, room 1. tel/fax: + 7 (495) 755 55 03 WARRANTY: The manufacturer ensures compliance with the declared characteristics, the period of preservation	

of sterility, the safety and integrity of the packaging subject to proper transportation, operation, and storage of the product for 2 years from the date of manufacture/sterilization of the product. **Shelf life** - 24 months from the date of manufacture / sterilization indicated on the package, provided that it is stored in the original packaging, in a dry, well-ventilated place, protected from direct sunlight, at a temperature not exceeding 25°C.

Warranty for **PICO Single Use Negative Pressure Wound Therapy System, variations** applies to materials, functioning and lack of production defects, provided that the device is used for only one patient. Reuse is prohibited.

12. NAME AND LEGAL ADDRESS OF THE MANUFACTURER, AUTHORIZED REPRESENTATIVE OF MD

MANUFACTURER:

Smith & Nephew Medical Limited, 101 Hessle Road, Hull, HU3 2BN, United Kingdom

PLACE OF MANUFACTURE:

1. Smith & Nephew Medical Limited, 101 Hessle Road, Hull, HU3 2BN, United Kingdom
2. Flextronics Technology (Shenzhen) Co., Ltd. 3-6/F, Building 1 and 6/F, Building 4, Yusheng Industrial Park, 467 Xixiang Section, National Highway 107, Xixiang, Baoan District, 518126 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA.
3. Avery Dennison Belgie BVBA, Tieblokkenlaan 1, Turnhout B-2300, Belgium.
4. Smith & Nephew Medical (Suzhou) Limited, No. 12, Wuxiang Road Comprehensive Free Zone, West Zone Suzhou Industrial Park Suzhou Jiangsu 215021 China.

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA:

Smith & Nephew LLC, 105120, Russia, Moscow, 2nd Syromyatnichesky Lane, 1, 9th floor, office 1, room 1. tel/fax: + 7 (495) 755 55 03

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