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Covidien llc., USA

Regulatory Affairs Operations Specialist

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"Operating documents" for submission in the Russian Federation:

Patient Return Electrode Valleylab REM Polyhesive with accessories

Manufactured by:

Covidien llc., USA. 15 Hampshire Street, Mansfield MA 02048, USA

Manufacturing sites:

Covidien llc, 5920 Longbow Drive, Boulder, Colorado, 80301, USA;

New Deantronics Taiwan LTD, 12f. No. 51, Sec. 4, Chong Yang Rd., Tu Cheng Dist., New Taipei City,
Taiwan, R.O.C.


State of Colorado

County of Boulder

} §

Subscribed and sworn to before me this 17th day of August 2020. Proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same for purposes therein stated.

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COVIDIEN™

Valleylab™

Cordless Patient Return Electrodes

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REF E7508

**REM Polyhesis™
Adult Cordless
Patient Return
Electrode**

REF E7509

**REM Polyhesis™
Adult Cordless
Patient Return
Electrode**

REF E7509B

**REM Polyhesis™
Adult Cordless
Patient Return
Electrode**

A cordless, single-use patient return electrode with conductive adhesive for use with E0560 REM™ Patient Return Electrode Cord. The patient return electrodes are for use with Covidien REM-equipped generators.

Do not use if package is opened or damaged.

Intended Use

All Covidien patient return electrodes are designed for use in traditional monopolar surgical procedures with duty cycles of 25% or less (for example: 10 seconds on followed by 30 seconds off).

These patient return electrodes are intended for use with electrosurgical generators equipped with a contact quality monitoring system such as the Covidien REM system.

REM System Considerations

The REM system is designed to minimize the risk of burns at the patient return electrode site due to a reduction in patient contact area during monopolar electrosurgery.

During a surgical procedure, the amount of current delivered in a given time determines the amount of heat that occurs under the patient return electrode. It is not possible to foresee what combination of current and duty cycle may be safely used in every situation.

When using a patient return electrode, the user should consult the generator and accessory manufacturers' recommendations and technical specifications regarding recommended maximum duty cycles.

Notice

Ensure the patient return electrode selected for the procedure is appropriate. Use the following chart as a guide to patient return electrode selection.

Patient Weight & Recommended Patient Return Electrode Model	
Patient Weight	Electrode Model
1 lb – 6 lb 0.45 kg – 2.7 kg	E7512
6 lb – 30 lb 2.7 kg – 13.6 kg	E7510 – 25 E7510 – 25DB
>30 lb >13.6 kg	E7506 E7507 E7507DB E7508 E7509 E7509B

Warning

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

Confirm proper electrosurgical generator settings before proceeding with surgery. Use the lowest power settings that will achieve the desired surgical effect.

Do not use the patient return electrode if it is damaged or modified. Product performance may be compromised. Replace before proceeding.

Non-traditional procedures that utilize high current, long duty cycles, or both (for example: tissue lesioning, tissue ablation, tissue vaporization, and procedures in which conductive fluids such as saline or lactated Ringer's solution are introduced into the surgical site for distention or to conduct the RF current) increase the risk of excessive heating under a fully applied patient return electrode to the point of injuring the patient. Use of more than one patient return electrode may help mitigate the increased risk.

Caution

Do not use patient return electrode if package seal is broken or conductive adhesive is dry.

Caution

Do not attempt to sterilize patient return electrodes. Sterilization of this product may compromise its performance and safety.

Do not attempt to relocate the patient return electrode after initial application.

Electrode gel should not be used with this patient return electrode. Gel is incompatible with the Polyheseive surface and will compromise patient return electrode performance.

Position surgical electrode cables to avoid contact with the patient or other leads.

Do not use the patient return electrode beyond its expiration date.

Do not reuse the patient return electrode. It is single use only.

E7509/E7509B - To avoid the possibility of a dry patient return electrode, use the package contents within four weeks after opening. Reseal after each use.

When applying the patient return electrode, do not allow adjoining edges to touch or overlap.

Refer to the generator user's guide for additional cautions and warnings.

Select and Prepare the Application Site

1. Select a well vascularized, convex area in close proximity to the surgical site for patient return electrode application. Avoid scar tissue, bony prominences, excessive adipose tissue, and areas where fluid may pool.

**Caution**

Special clinical circumstances may require alternate application sites. If alternate sites are necessary, ensure maximum patient-to-patient return electrode contact.

2. The patient return electrode site should be free of excessive hair. Remove hair from the selected application site in accordance with the policies and procedures of your facility.


Caution

Excessive hair at the application site may result in high impedance and poor electrical conduction. High impedance can result in a REM alarm.

3. Clean and dry the application site as needed. The application site must be free of oils, lotions, and other topically applied products to ensure proper adhesion between the patient's skin and the patient return electrode.

Apply the Patient Return Electrode

1. Check the patient return electrode's expiration date. If expired, discard the patient return electrode.
2. **E7508** - Remove the patient return electrode from the pouch.
E7509/E7509B - Open the patient return electrode pouch at the zip lock. Remove a patient return electrode and reseal the pouch.
3. Inspect the cord prior to each use to make sure that the clamp, connector, and cable are not broken or damaged. Discard damaged cords.
4. Open the clamp by lifting the latch. Do not force the latch beyond 90° when in the lifted position.
5. Fold back the liner from the patient return electrode tab to expose the foil contacts.
6. Insert the tab completely into the clamp with the foil contacts facing up so that the entire foil-contact area is covered by the clamp.

7. Lower the latch until it is level with the body of the clamp.

8. Remove the liner from the patient return electrode and lightly touch the surface of the conductive adhesive to ensure that it is moist. Do not use dry patient return electrodes.

9. Apply the patient return electrode to the patient and lightly massage the entire surface of the patient return electrode to ensure secure contact with the patient's skin.

10. Prior to connecting the patient return electrode to the REM-equipped generator, turn the generator on. If operating properly, a brief self-test will be followed by an audible tone and the REM-alarm indicator displays red. This indicates that no applied patient return electrode is detected, resulting in a REM alarm.
11. Insert the patient return electrode connector into the generator to correct the alarm condition. If corrected, the warning tone ceases and the REM-alarm indicator changes to green. If the alarm condition is not corrected, refer to the section Correcting REM-Alarm Conditions.


Caution

Whenever the generator is on and not in use, keep active accessories away from the patient.

Remove the Patient Return Electrode

1. After the procedure, disconnect the patient return electrode connector from the electrosurgical generator.

2. Slowly remove the patient return electrode with one hand while supporting the skin with the other to avoid skin trauma.

3. Remove the cord from the patient return electrode by lifting the latch.
4. Discard the patient return electrode.

Cord Care

If the cord requires cleaning, wipe with a mild cleaning solution or blood-dissolving detergent. Rinse the cord thoroughly by wiping with a clean, damp cloth. *Do not* immerse the cord. Wipe the cord completely dry before reuse.

Do not autoclave as the process will damage the cord.

Important

Discard and replace any reusable cord that cannot be adequately cleansed of contaminants.

Correcting REM-Alarm Conditions

A REM alarm condition occurs when the REM Contact Quality Monitoring System detects compromised contact between the patient and the patient return electrode. During a REM alarm, the generator's REM-fault indicator illuminates red, an audible alarm sounds, and the generator does not produce monopolar output power. The REM-alarm condition clears when the REM Contact Quality Monitoring System detects that the patient return electrode and the patient contact resistance are within the acceptable range.

A REM alarm can result from:

- Failure to properly insert the REM patient return electrode connector into the generator
- Failure to properly attach the clamp to the patient return electrode
- A damaged cord
- Poor contact quality (high impedance) between the patient and the patient return electrode
- Extremely low impedance between the patient and the patient return electrode

High impedance (greater than 135 ohms) between the patient and the patient return electrode may result from the patient return electrode not being placed firmly on the skin or being placed on high-impedance tissue—such as adipose or scar tissue, bony areas, excessive hair, or an emaciated patient.

Low impedance (less than 5 ohms) between the patient and the patient return electrode may result from placement on edematous tissue.

REM-Alarm Troubleshooting Instructions

Use these troubleshooting instructions to correct the problem causing the REM alarm. It may be necessary to choose an alternate site to achieve adequate contact between the patient and patient return electrode.

1. Ensure that the patient return electrode is applied to a well vascularized area close to the surgical site. Avoid scar tissue, adipose tissue, bony prominences, and areas where fluid may pool.
2. Turn on the generator. (Press (I) on the power switch.)

REM alarm not silenced

1. Disconnect the patient return electrode cord from the generator.
2. Inspect the center pin of the patient return electrode connector to ensure that it is not bent or missing. (If damaged, replace the cord before proceeding.)
3. Insert the patient return electrode connector into the generator's receptacle. Make sure the center pin goes into the corresponding hole and the connector is fully inserted.

REM alarm not silenced

1. Apply firm pressure over the entire surface area of the patient return electrode, particularly the center.

REM alarm not silenced

1. Disconnect the patient return electrode cord connector from the generator. (Do not remove the cord from the patient return electrode or remove the patient return electrode from the patient.)
2. Select an alternative site on the patient's thigh, upper biceps, calf, or lower back for a second patient return electrode.
3. Attach a second patient return electrode to another cord, and apply to the patient.

4. Insert the connector of the second patient return electrode cord into the generator.

REM alarm not silenced

1. Disconnect the patient return electrode cord connector from the generator. (Do not disconnect the cord from the patient return electrode or remove attached patient return electrodes from the patient.)
2. Connect a third patient return electrode to another cord.
3. Apply the patient return electrode to the patient's thigh, upper biceps, calf, or lower back.


4. Insert the connector of the third patient return electrode into the generator.

REM alarm not silenced

1. Disconnect the patient return electrode cord connector from the generator. (Do not disconnect the cord from the patient return electrode or remove attached patient return electrodes from the patient.)
2. Insert a Valleylab Multiple Return/S Cord Adapter (E0507B) into the generator's patient return electrode receptacle.
3. Select the cords from the two patient return electrodes that are situated on the most well vascularized, convex areas, and are in closest proximity to the surgical site.
4. Insert the connectors of two of the patient return electrode cords into the adapter.

REM alarm not silenced

1. Use a backup electrostimulation generator and repeat these troubleshooting steps.

REM alarm not silenced

1. Remove all patient return electrodes from the patient, and use bipolar output and appropriate accessories for the surgery.

After Correcting the REM Alarm Condition

Once the REM alarm condition is corrected, leave the generator on during draping to ensure that any attached patient return electrodes are not disturbed. Remove any patient return electrodes from the patient that are not in use.

When ever the generator is on and not in use, keep all cables and accessories away from the patient.

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Part No. 1012058

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG.

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© 2011 Covidien.

Covidien Inc, 15 Hampshire Street,
Mansfield, MA 02048 USA.

Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made in USA. Printed in USA.

www.covidien.com

REV 02/2011

Rx
ONLY

0086

NEAR

[REF] E7507 REM Polyhesis™ Adult Patient Return Electrode
9' (2.7 m) Cord

[REF] E7507DB REM Polyhesis™ Adult Patient Return Electrode
15' (4.6 m) Cord

A single-use patient return electrode with conductive adhesive and integrated cable for use with Covidien REM equipped generators.

Intended Use

All patient return electrodes are designed for use in traditional monopolar surgical procedures with duty cycles of 25% or less (for example: 10 seconds on followed by 30 seconds off).

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REM System Considerations

The REM system is designed to minimize the risk of burns at the return-electrode site due to a reduction in the patient contact area during monopolar electrosurgery.

During a surgical procedure, the amount of current delivered in a given time determines the amount of heat that occurs under the electrode. It is not possible to foresee what combination of current and duty cycle may be safely used in every situation.

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Do not use the return electrode if it is damaged or modified. Product performance may be compromised. Replace before proceeding.

Non-traditional procedures that utilize high current, long duty cycles, or both (for example: tissue lesioning, tissue ablation, tissue vaporization, and procedures in which conductive fluids such as saline or lactated Ringer's solution are introduced into the surgical site for distention or to conduct the RF current) increase the risk of excessive heating under a fully applied return electrode to the point of injuring the patient. Use of more than one patient return electrode may help mitigate the increased risk.

Caution

Do not use the electrode if package seal is broken or conductive adhesive is dry.

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Do not attempt to relocate the return electrode after initial application.

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Caution

Electrode gel should not be used with this return electrode. Gel is incompatible with the Polyhesiva surface and will compromise return electrode performance.

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Do not use a return electrode beyond its expiration date.

Do not reuse the return electrode. It is single use only.

When applying the return electrode, do not allow adjoining edges to touch or overlap.

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1. Select a well vascularized, convex area in close proximity to the surgical site for electrode application. Avoid scar tissue, bony prominences, excessive adipose tissue, and areas where fluid may pool.

Caution

Special clinical circumstances may require alternate application sites. If alternate sites are necessary, ensure maximum patient-to-return-electrode contact.

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
3. Clean and dry the application site as needed. The application site must be free of oils, lotions, and other topically applied products to ensure secure contact between the patient's skin and the return electrode.

Apply the Return Electrode

1. Check the return electrode's expiration date. If expired, discard the return electrode.
2. Remove the return electrode from the pouch.
3. Remove the liner from the return electrode, and lightly touch the surface of the conductive adhesive to ensure that it is moist. Do not use dry return electrodes.

4. Apply the patient return electrode to the patient and lightly massage the entire surface of the return electrode to ensure secure contact with the patient's skin.

5. Prior to connecting the return electrode to the REM-equipped generator, turn the generator on. If operating properly, a brief self test will be followed by an audible tone and the REM-alarm indicator displays red. This indicates that no patient-applied return electrode is detected, resulting in a REM alarm.
6. Insert the return-electrode connector into the generator to correct the alarm condition. If corrected, the warning tone ceases and the REM-alarm indicator changes to green. If the alarm condition is not corrected, refer to the following section, Correcting REM-Alarm Conditions.

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Remove the Return Electrode

1. After the procedure, disconnect the patient electrode connector from the electrosurgical generator.
2. Slowly remove the return electrode with one hand while supporting the skin with the other to avoid skin trauma.

3. Dispose of the used patient return electrode.

Correcting REM-Alarm Conditions

A REM alarm condition occurs when the REM Contact Quality Monitoring System detects compromised contact between the patient and the return electrode. During a REM alarm, the generator's REM-fault indicator illuminates red, an audible alarm sounds, and the generator does not produce monopolar output power. The REM alarm condition clears when the REM Contact Quality Monitoring System detects that the patient return electrode and the patient contact resistance are within the acceptable range.

A REM alarm can result from:

- Failure to properly insert the REM return-electrode connector into the generator
- A damaged cable
- Poor contact quality (high impedance) between the patient and the return electrode
- Extremely low impedance between the patient and the return electrode

High impedance (greater than 135 ohms) between the patient and the return electrode may result from the electrode not being placed firmly on the skin or being placed on high-impedance tissue—such as adipose or scar tissue, bony areas, excessive hair, or an emaciated patient.

Low impedance (less than 5 ohms) between the patient and the return electrode may result from placement on edematous tissue.

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Use these troubleshooting instructions to correct the problem causing the REM alarm. It may be necessary to choose an alternate site to achieve adequate contact between the patient and return electrode.

1. Ensure that the return electrode is applied to a well vascularized area close to the surgical site. Avoid scar tissue, adipose tissue, bony prominences, and areas where fluid may pool.
2. Turn on the generator (Press **ON** on the power switch.)

REM alarm not silenced

1. Disconnect the return-electrode cable from the generator.
2. Inspect the center pin of the return-electrode connector to ensure that it is not bent or missing. (If damaged, replace the return electrode before proceeding.)
3. Insert the return-electrode connector into the generator. Make sure the center pin goes into the corresponding hole and the connector is fully inserted.

REM alarm not silenced

Apply firm pressure over the entire surface area of the return electrode, particularly the center.

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REM alarm not silenced

1. Disconnect the electrode-cable connector from the generator. (Do not remove the return electrode from the patient.)
2. Select an alternative site on the patient's thigh, upper biceps, calf, or lower back for a second return electrode.
3. Apply a second return electrode to the patient.


4. Insert the connector of the second return electrode into the generator.

REM alarm not silenced

1. Disconnect the electrode-cable connector from the generator. (Do not remove attached return electrodes from the patient.)
2. Apply a third return electrode to the patient's thigh, upper biceps, calf, or lower back.


3. Insert the connector of the third return electrode into the generator.

REM alarm not silenced

1. Disconnect the electrode cable connector from the generator. (Do not remove attached return electrodes from the patient.)
2. Insert an E0507B Multiple Return/S Cord Adapter into the generator's return-electrode receptacle.
3. Select the cables from the two return electrodes that are situated on the most well vascularized, convex areas, and are in closest proximity to the surgical site.
4. Insert the connectors of two of the return-electrode cables into the adapter.

REM alarm not silenced

Use a backup electrosurgical generator and repeat these troubleshooting steps.

REM alarm not silenced

Remove the return electrode(s) from the patient, and use bipolar output and appropriate accessories for the surgery.

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Once the REM-alarm condition is corrected, leave the generator on during draping to ensure that any attached return electrodes are not disturbed. Remove any return electrodes from the patient that are not in use.

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REF E7510-25DB

**REM Polyhesive™
Infant Patient
Return Electrode**
15" (4.6 m) Cord

A single-use patient return electrode with conductive adhesive for use with Covidien REM™-equipped generators.

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Type BF applied part

Indications for Use

The REM Polyhesive Infant Patient Return Electrode is a single use, non-sterile dispersive electrode with a preattached cord. The electrode adheres to the patient over its entire surface. Its purpose is to complete the electrosurgical circuit between the generator, the active electrode, and the patient.

The product is used for general monopolar electrosurgery on newborns, infants, and children weighing between approximately 6 and 30 lb.

REM System Considerations

This split-plate patient return electrode is intended for use only with electrosurgical generators equipped with a contact quality monitoring system such as the Covidien REM system. See the compatibility table on the front cover of this manual for a list of generators known to be compatible with the Polyhesive Infant Return Electrode.

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Performance Considerations

All Covidien patient return electrodes are designed for use in traditional monopolar electrosurgical procedures with duty cycles of 25% or

less (for example: 10 seconds on followed by 30 seconds off).

During a surgical procedure, the amount of current delivered in a given time determines the amount of heat that occurs under the patient return electrode. It is not possible to foresee what combination of current and duty cycle may be safely used in every situation. Always use the lowest power levels possible to achieve the intended tissue effect and strictly observe the 25% duty cycle established for the patient return electrode.

When using a patient return electrode, the user should consult the generator and active accessory labeling and technical specifications regarding recommended maximum duty cycles for the generator and active accessories.

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Product Limitations This patient return electrode is only intended for use on infants weighing 6 lb. to 30 lb. (2.7 kg to 13.6 kg).

Warning

Do not use in high-current (> 500 mA) and/or long-duty cycle (activation time greater than 60 seconds) procedures (for example: tissue lesioning, tissue ablation, tissue vaporization, and procedures in which conductive fluids such as saline or lactated Ringer's solution are introduced into the surgical site for distention or to conduct the RF current). These conditions increase the risk of excessive heating under a fully applied patient return electrode to the point of seriously injuring the patient.

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Precaution

Surgeons should be familiar with the literature on electrosurgical effects on small patients and may wish to consider bipolar electrosurgery, which does not require a patient return electrode.

Use of this product with the adhesive border is recommended. When using this product without the adhesive border, the patient return electrode should be consistently monitored for fluid invasion and should not be placed in areas where fluids are present.

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Do not attempt to sterilize Covidien patient return electrodes. Sterilization of this product may compromise its performance and safety.

Do not attempt to relocate the patient return electrode after the initial application.

Electrode gel should not be used with this patient return electrode. Gel is incompatible with the Polyhesive surface and will compromise patient return electrode performance.

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1. Whenever patient size permits, use an adult-size patient return electrode to reduce the possibility of patient burns. The patient return electrode must adhere completely to the application site without overlapping.
2. Select a well vascularized, convex area in close proximity to the surgical site for patient return electrode application. Avoid scar tissue, bony prominences, excessive adipose tissue, and areas where fluid may pool. The preferred application site for an infant is the back or torso.

Precaution

Special clinical circumstances may require alternate application sites. If alternate sites are necessary, ensure maximum patient contact between the patient and the patient return electrode.

3. Clean and dry the application site as needed. The application site must be free of oils, lotions, and other topically applied products to ensure secure contact between the patient's skin and the patient return electrode.

Apply the Patient Return Electrode

1. Check the patient return electrode's expiration date. If expired, discard the patient return electrode.
2. Remove the patient return electrode from the pouch.
3. Remove the liner, and lightly touch the surface of the patient return electrode to ensure the conductive adhesive is moist. Do not use dry patient return electrodes.
4. Use the patient return electrode with or without the adhesive border around the exterior edge. To use with the adhesive border, remove the separate border liner.
5. Place the patient return electrode on the patient. For maximum patient safety, the entire conductive surface of the patient return electrode should be in contact with the patient.
6. Lightly massage the entire surface of the patient return electrode to ensure secure contact with the patient's skin.

7. Prior to connecting the patient return electrode to the REM-equipped generator, turn the generator on. If operating properly, a brief self test will be followed by an audible tone and the REM-alarm indicator displays red. This indicates that no applied patient return electrode is detected, resulting in a REM alarm.
8. Insert the patient return electrode connector into the generator to correct the alarm condition. If corrected, the warning tone ceases and the REM-alarm indicator changes to green. If the alarm condition is not corrected, refer to the following section, Correcting REM-Alarm Conditions.

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Remove the Patient Return Electrode

1. After the procedure, disconnect the patient return electrode connector from the electrosurgical generator.
2. Slowly remove the patient return electrode with one hand while supporting the skin with the other to avoid skin trauma.
3. Dispose of the used patient return electrode in compliance with local regulations for the disposal of medical waste.

Correcting REM-Alarm Conditions

A REM-alarm condition occurs when the REM Contact Quality Monitoring System detects compromised contact between the patient and the patient return electrode. During a REM alarm, the generator's REM-fault indicator illuminates red, an audible alarm sounds, and the generator does not produce monopolar output power. The REM-alarm condition clears when the REM Contact Quality Monitoring System detects that the patient return electrode and the patient contact resistance are within the acceptable range.

A REM alarm can result from:

- Failure to properly insert the REM patient return electrode connector into the generator
- A damaged cord
- Poor contact quality (high impedance) between the patient and the patient return electrode
- Extremely low impedance between the patient and the patient return electrode

High impedance (greater than 135 ohms) between the patient and the patient return electrode may result from the patient return

being placed on high-impedance tissue—such as adipose or scar tissue, bony areas, excessive hair, or an emaciated patient.

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REM alarm not silenced

1. Disconnect the patient return electrode cable from the generator.
2. Inspect the center pin of the patient return electrode connector to ensure that it is not bent or missing. (If damaged, replace the patient return electrode before proceeding.)
3. Insert the patient return electrode connector into the generator. Make sure the center pin goes into the corresponding hole and the connector is fully inserted.

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REM alarm not silenced

1. Disconnect the patient return electrode from the generator.
2. Remove the patient return electrode from the patient and replace it with another E7510-25 patient return electrode.
3. Insert the patient return electrode connector into the generator.

REM alarm not silenced

1. Use a backup electrosurgical generator and repeat these troubleshooting steps.

REM alarm not silenced

1. Remove the patient return electrode from the patient, and use bipolar output and appropriate accessories for the surgery.

After Correcting the REM-Alarm Condition

Once the REM-alarm condition is corrected, leave the generator on during draping to ensure that any attached patient return electrodes are not disturbed. Remove any patient return electrodes from the patient that are not in use.

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Part No. PT00033894

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and Internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company,  brands are trademarks of their respective owner.

© 2011 Covidien,
 Covidien Inc,
15 Hampshire Street, Mansfield, MA 02048 USA.
 Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made In USA. Printed in USA.

www.covidien.com

REV 04/2016

Single Use

Rx
ONLY

Consult
Instructions
for use

0086

REF E7512

Valleylab™

REM Polyhesive™ Neonatal Patient Return Electrode

- en** Instructions for Use
- fr** Mode d'emploi
- pt** Instruções de uso
- de** Gebrauchsanleitung
- es** Instrucciones de uso
- it** Istruzioni per l'uso
- nl** Gebruiksaanwijzing
- sv** Bruksanvisning
- ru** Инструкция по применению
- zh** 使用说明

Part No. 1012119

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

Other brands are trademarks of a Covidien company.

© 2011 Covidien.
Covidien Inc, 15 Hampshire Street,
Mansfield, MA 02048 USA.
Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made In USA. Printed in USA.

www.covidien.com

REV 02/2011

Rx
ONLY

en

Valleylab™

REF E7512

REM Polyheseive™ Neonatal Patient Return Electrode 9' (2.7 m) Cord

A reduced-size, single-use patient return electrode with conductive adhesive for patients weighing 1 lb. to 6 lb. (0.45 kg to 2.7 kg). For use with Covidien REM™-equipped generators.

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This split-plate patient return electrode is intended for use with electrosurgical generators equipped with a contact quality monitoring system such as the Covidien REM system.

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Patient Weight & Recommended Patient Return Electrode Model	
Patient Weight	Electrode Model
1 lb. – 6 lb. 0.45 kg – 2.7 kg	E7512

Patient Weight & Recommended Patient Return Electrode Model

6 lb. – 30 lb. 2.7 kg – 13.6 kg	E7510 – 25 E7510 – 25DB
>30 lb. >13.6 kg	E7506 E7507 E7507DB E7508 E7509 E7509B

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Product Limitations—This patient return electrode is only intended for use on infants weighing 1 lb. to 6 lb. (0.45 kg to 2.7 kg). For infants larger than the recommended weight, use Covidien's E7510-25 Infant REM Polyheseive patient return electrode.

Do not use in high-current and/or long-duty cycle procedures (for example: tissue lesioning, tissue ablation, tissue vaporization, and procedures in which conductive fluids such as saline or lactated Ringer's solution are introduced into the surgical site for distention or to conduct the RF current). These conditions increase the risk of excessive heating under a fully applied patient return electrode to the point of seriously injuring the patient.

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Applied current should not exceed 300 milliamperes nor be applied for longer than 30 seconds continually.

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Caution

Surgeons should be familiar with the literature on electrosurgical effects on small patients, and may wish to consider bipolar electrosurgery, which does not require a patient return electrode.

Caution

Use of this product with the adhesive border is recommended. When using this product without the adhesive border, the patient return electrode should be consistently monitored for fluid invasion, and should not be placed in areas where fluids are present.

Do not use the patient return electrode if the package seal is broken or conductive adhesive is dry.

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Procedure to Determine Maximum Power Settings

Prior to initial pad use, the electrosurgical generator should be calibrated by a qualified biomedical engineer or technician to determine the control setting that represents the maximum current level (300 milliamps) recommended for this patient return electrode. A standard electrosurgical analyzer that uses a 125 ohm load should be used for this calibration. A reading of 300 milliamps should be identified in each cut, coag, and blend mode. Instructions should be mounted on the electrosurgical generator that state the unit must not exceed these settings when using the E7512 Neonatal REM Polyhesive Patient Return Electrode.

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1. The recommended application site on premature infants 1 lb. to 6 lb. (0.45 kg to 2.7 kg) is on the back, inferior to the shoulder blades and superior to the sacrum.
2. Select a well vascularized, convex area in close proximity to the surgical site for patient return electrode application. Avoid scar tissue, bony prominences, excessive adipose tissue, and areas where fluid may pool.

Caution

Special clinical circumstances may require alternate application sites. If alternate sites are necessary, ensure maximum patient-to-patient return electrode contact and minimize current levels.

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Apply the Patient Return Electrode

1. Check the patient return electrode's expiration date. If expired, discard the patient return electrode.
2. Remove the patient return electrode from the pouch.
3. Remove the liner, and lightly touch the surface of the patient return electrode to ensure the conductive adhesive is moist. Do not use dry patient return electrodes.
4. Use the patient return electrode with or without the adhesive border around the exterior edge. To use with the adhesive border, remove the separate border liner.
5. Place the patient return electrode on the patient. For maximum patient safety, the entire conductive surface of the patient return electrode should be in contact with the patient.
6. Lightly massage the entire surface of the patient return electrode to ensure secure contact with the patient's skin.
7. Prior to connecting the patient return electrode to the REM-equipped generator, turn the generator on. If operating properly, a brief self test will be followed by an audible tone and the REM-alarm indicator displays red. This indicates that no patient-applied patient return electrode is detected, resulting in a REM alarm.
8. Insert the patient return electrode connector into the generator to correct the alarm condition. If corrected, the warning

tone ceases and the REM-alarm indicator changes to green. If the alarm condition is not corrected, refer to the following section, Correcting REM-Alarm Conditions.

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Remove the Patient Return Electrode

1. After the procedure, disconnect the patient return electrode connector from the electrosurgical generator.
2. Slowly remove the patient return electrode with one hand while supporting the skin with the other to avoid skin trauma.
3. Dispose of the used patient return electrode.

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A REM alarm can result from:

- Failure to properly insert the REM patient return electrode connector into the generator
- A damaged cable
- Poor contact quality (high impedance) between the patient and the patient return electrode
- Extremely low impedance between the patient and the patient return electrode

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Use these troubleshooting instructions to correct the problem causing the REM alarm. It

may be necessary to choose an alternate site to achieve adequate contact between the patient and patient return electrode.

1. Ensure that the patient return electrode is applied to a well vascularized area close to the surgical site. Avoid scar tissue, adipose tissue, bony prominences, and areas where fluid may pool.
2. Turn on the generator. (Press [I] on the power switch.)

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REM alarm not silenced

1. Disconnect the patient return electrode cable connector from the generator.
2. Remove the patient return electrode from the patient and replace with another E7512 patient return electrode.
3. Insert the patient return electrode connector into the generator.

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Valleylab™

REM Patient Return Electrode Cord and Clamp

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Part No. 1016717

COVIDIEN, COVIDIEN with logo and the Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG.

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© 2011 Covidien.

Manufactured for Covidien LLC, 15 Hampshire Street, Mansfield, MA 02048 USA.
EC/REP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.

Made in Taiwan. Printed in USA.

www.covidien.com

REV 02/2011

Rx
ONLY

REF E0560 REM Patient Return Electrode Cord, Reusable
15' (4.6 m)

REF E0560E Cord Clamp, ERBE, Reusable
15' (4.6 m)


For use with E7509 REM Polyhesive™ Adult Cordless Patient Return Electrodes

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Before Surgery

1. Inspect the cord prior to each use to make sure that the clamp, plug, and cable are not broken or damaged. Discard damaged cords.
2. Open the clamp by lifting the latch. Do not force the latch beyond 90° when in the lifted position.
3. Remove the release liner from the return electrode tab area, exposing the foil contacts.

4. Insert the tab completely into the clamp with the foil contacts facing up so that the entire foil contact area is covered by the clamp.
5. Lower the latch until it is level with the body of the clamp.

6. Apply the patient return electrode according to the manufacturer's instructions.
7. Connect the cord plug to the generator.

For Covidien generators: prior to inserting the cord plug into the REM™ generator, turn the generator to the ready (on) mode to ensure that the alarm system is operating properly. Refer to the generator instruction manual for additional information.

Connect the cord plug to the generator to correct the alarm condition (audible alarm silenced, alarm light off). Return the generator to the standby (off) mode until ready for use. If the alarm condition is not corrected, refer to "Correcting REM Alarm Conditions" in the REM Polyhesive Patient Return Electrode instructions for use.

For non-Covidien generators: refer to the generator instruction manual.

After Surgery

1. To avoid skin trauma when removing the return electrode, peel off slowly with one hand while supporting the underlying tissue with the other hand.
2. To release the clamp, lift the latch.

Cord Care

If cleaning is necessary, wipe the cord with a mild cleaning solution or blood dissolving detergent. Rinse the cord thoroughly by wiping with a clean, damp cloth. Do not immerse. Wipe the cord completely dry before reuse.

Do not autoclave. This process will damage the cord.

Detailed Cord Inspection Instructions

Follow your institution's procedure for periodic inspection of the cord by the hospital biomedical department.

1. Inspect the cable for nicks or cuts. Bend the cable back and forth at the entrance to the clamp. Look for broken strands or corrosion on the wire. Repeat this inspection at the plug end. If any of these problems are found, discard the cord.
2. Inspect the clamp and plug for cracks. Check to make sure that the pin on the plug is in place and firmly attached. Discard the cord if there are any cracks or a missing or loose pin.
3. Using a milliohm meter, attach one lead to the clamp end of the cord and the other lead to the plug end. Measure the resistance of the cord while flexing the cable where it enters the clamp and the plug. The resistance must be less than 900 milliohms. Discard the cord if the resistance is higher than 900 milliohms.

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For Valleylab generators: prior to inserting the cord plug into the REM generator, turn the generator to the ready (on) mode to ensure that the alarm system is operating properly. Refer to the generator instruction manual for additional information.

Connect the cord plug to the generator to correct the alarm condition (audible alarm silenced, alarm light off). Return the generator to the standby (off) mode until ready for use. If the alarm condition is not corrected, refer to "Correcting REM Alarm Conditions" in the PolyHesive II Patient Return Electrode instructions for use.

For non-Valleylab generators: refer to the instruction manual.

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The cord can be sterilized with ethylene oxide gas up to 5 times if necessary. Do not autoclave. This process will damage the cord.

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1. Inspect the cable for nicks or cuts. Bend the cable back and forth at the entrance to the clamp. Look for broken strands or corrosion on the wire. Repeat this inspection at the plug end. If any of these problems are found, discard the cord.
2. Inspect the clamp and plug for cracks. Check to make sure that the pin on the plug is in place and firmly attached. Discard the cord if there are any cracks or a missing or loose pin.
3. Using a milliohmmeter, attach one lead to the clamp end of the cord and the other lead to the plug end. Measure the resistance of the cord while flexing the cable where it enters the clamp and the plug. The resistance must be less than 900 milliohms. Discard the cord if the resistance is higher than 900 milliohms.

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
Сurgical cable must be kept away from people and other electrical equipment.

Сurgical before

1. Each time before use, check the cable. Confirm the clip, plug and cable are not cracked or damaged. Throw away the damaged cable.
2. Turn the clip and open the clip. Do not force the clip to rotate more than 90°.
3. Remove the protective cover of the clip area. Expose the contact surface.

4. Let the contact surface face up. Insert the clip completely into the clip. Make the clip stay in the contact surface.
5. Put the clip. Make it flat with the clip.

6. Follow the manufacturer's instructions to attach the patient's return electrode.
7. Connect the clip to the electrode.

For Valleylab electrode: Insert the clip into the REM electrode. Turn the electrode to the "REM" mode. Confirm the system is working. For other electrodes, refer to the electrode user manual.

Connect the clip to the electrode to correct the alarm condition (alarm sound, alarm light). Turn the electrode to the "REM" mode. Confirm the system is working. If the alarm condition cannot be corrected, refer to the PolyHesive II patient return electrode user manual for the "Correct REM alarm condition".

For non-Valleylab electrode: Refer to the electrode user manual.

Сurgical after

1. To avoid skin injury when removing the return electrode, use one hand to hold the return electrode while the other hand holds the tissue below it.
2. Turn the clip and open the clip.

Сurgical maintenance


If necessary, clean the cable with a mild detergent or alcohol solution. Clean the cable with a clean cloth. Do not immerse the cable in liquid. After use, dry the cable thoroughly.

If necessary, use ethylene oxide gas to sterilize the cable. Sterilize the cable for 5 cycles. Do not use high pressure steam sterilization. This method will damage the cable.

Сurgical inspection details

Follow the hospital's medical department's regular inspection of the cable.

1. Check the cable for cracks or cuts. Bend the cable at the connection point. Check for cracks or corrosion. Check the connection point. If any problem is found, throw away the cable.
2. Check the clip and the plug for cracks. Check the clip on the needle. If the clip is not fixed in place, if there are cracks or the needle is loose or missing, throw away the cable.
3. Use a multimeter. Connect one lead to the clip and the other lead to the plug. Measure the resistance of the cable. The resistance should be less than 900 ohms. If the resistance is higher than 900 ohms, throw away the cable.

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