
Phillip H Jones – Notary Public
Two Chimneys, North Street
Winkfield Windsor SL4 4TE England
Telephone +44 (0) 7860 398485
Email: info@notarypublicassistance.co.uk

I, PHILLIP JONES, Notary Public duly admitted and sworn, carrying on practice at Two Chimneys North Street Winkfield Windsor SL4 4TE England.

CERTIFY THAT the document hereto attached entitled "Instruction for Use of the Medical Device", bearing date 07th July 2022, the original of which purports to be signed by Ken Fergusson, Senior Regulatory Affairs Manager of Smith & Nephew Medical Limited is to the best of my knowledge and belief a true copy of the original.

In faith and testimony whereof I, the said Notary, have subscribed my name and set and affixed my seal of office at Windsor aforesaid this Tuesday, 12th July 2022.

Notary Public

My Commission expires on death

Protocol no. 9229 /22

Regulated by the Faculty Office of the Archbishop of Canterbury

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101 Hessle Road, Hull, HU3 2BN, United Kingdom

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Senior Regulatory Affairs Manager

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
smith&nephew

PICO[◊]
Dressings
User Manual

Dressings are designed for use with the PICO single use negative pressure wound therapy system. Any information on the PICO systems refer to the respective PICO system user manual. PICO is intended for use by or on the direction of a trained and licensed physician in accordance with the PICO system user manual.

1. Description

PICO dressings are packaged sterile and are for single-use only.

The PICO dressing manages exudate through a combination of absorption and evaporation of moisture through the outer film. The PICO dressing is intended to be used for up to 7 days on low exuding wounds.

For moderate exuding wounds dressing changes will be required in line with standard wound management guidelines.

Low exuding wounds are considered to be up to 0.6g of liquid exudate/cm² of wound area/24 hours. Moderate exuding wounds are considered to be up to 1.1g of liquid exudate/cm² of wound area/24 hours. 1g of exudate is approximately equal to 1ml of exudate.

The frequency of dressing changes can be affected by multiple factors such as wound type, wound size, rate or volume of exudate, orientation or environmental conditions.

2. Indications for use

The PICO single use negative pressure systems are indicated for patients who would benefit from a suction device (NPWT) as it may promote wound healing via removal of low to moderate levels of exudate and infectious materials. Appropriate wound types include:

- Chronic
- Acute
- Traumatic
- Subacute and dehiscent wounds
- Partial-thickness burns
- Ulcers (such as diabetic or pressure)
- Flaps and grafts
- Surgically closed incision sites

PICO single use negative pressure systems are suitable for use both in a hospital and homecare setting.

3. Contraindications

PICO is contraindicated for:

- Patients with malignancy in the wound bed or margins of the wound (except in palliative care to enhance quality of life).
- Previously confirmed and untreated osteomyelitis.
- Non-enteric and unexplored fistulas.
- Necrotic tissue with eschar present.
- Exposed arteries, veins, nerves or organs.
- Exposed anastomotic sites.

PICO should not be used for the purpose of:

- Emergency airway aspiration.
- Pleural, mediastinal or chest tube drainage.
- Surgical suction.

Important Information

Important information for Patients and Healthcare Professionals

The pump contains a MAGNET.

You must keep the pump at least 10cm away from implantable

devices such as pacemakers, defibrillators, shunt valves, neurostimulators, cochlear

implants, or other life sustaining equipment such as insulin pumps. This warning applies

to ALL users, including patients, caregivers, healthcare professionals, and all people who may come within

10cm of the device.

For more information on electromagnetic immunity and electromagnetic emissions see: www.mypico.com or ask your Smith & Nephew representative for a hardcopy.

5. Warnings

1. Magnet Warning

The PICO pump contains a magnet. Interaction between your device and other medical devices can result in device failure and lead to serious harm, including death. For this reason the device must be kept at least 10cm away from all implantable or life sustaining medical devices. Extension tubes are included to enable the pump to be positioned more than 10cm away from any devices that could be affected. This warning applies to ALL users including patients, caregivers, healthcare professionals and all persons close to the device.

2. Certain patients are at high risk of bleeding complications which, if uncontrolled, could potentially be fatal. Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately disconnect pump, leave dressing in place, take appropriate measures to stop bleeding and seek immediate medical assistance.

3. Haemostasis must be achieved before applying the dressing, although the use of anticoagulants does not deem a patient inappropriate for treatment with PICO. Patients suffering from difficult haemostasis or who are receiving anticoagulant therapy have an increased risk of bleeding. During therapy, avoid using haemostatic products that may increase the risk of bleeding, if disrupted. Frequent assessment must be maintained throughout the therapy.

4. At all times care should be taken to ensure that the pump, tubing and connectors do not:

Lie in a position where they could cause pressure damage to the patient.

Trail across the floor where they could present a trip hazard or become contaminated.

Present a risk of strangulation or a tourniquet to patients.

Rest on or pass over a source of heat.

Become twisted or trapped under clothing or bandages so that the therapy is blocked.

5. Sharp edges or bone fragments in a wound must be covered or removed prior to using PICO due to risk of puncturing organs or blood vessels while under negative pressure wound therapy.

6. In the event that defibrillation is required, remove the dressing if it is positioned in a location that will interfere with defibrillation.

7. PICO has not been studied on paediatric patients. Patient size and weight should be considered when prescribing this therapy.

8. The system contains small parts which could represent a choking hazard for young children. Keep out of the reach of children.

9. PICO is unsuitable for use in areas where there is danger of explosion (e.g. hyperbaric oxygen unit)

10. PICO dressings should only be applied and changed by a healthcare professional. The patient must not remove their dressing unless the healthcare professional believes it to be clinically appropriate.

11. Each PICO dressing (including multisite) must be used to dress one wound only.

Precautions

Precautions should be taken in the following types of patients who are at high risk of bleeding situations:

Receiving anticoagulant therapy or platelet aggregation inhibitors or actively bleeding.

Having weakened or friable blood vessels or organs in or around the wound as a result of, but not limited to; anastomoses, infection, trauma or radiation.

Suffering from difficult wound hemostasis.

Untreated for malnutrition.

Non-compliant or combative.

Suffering from wounds in close proximity to blood vessels or delicate fascia.

2. If pain, reddening, odour, sensitization or a sudden change in the volume or colour of wound fluid occurs during use, contact your healthcare professional right away.

3. **Where PICO is used to bolster skin grafts, it is important to visually inspect the system regularly, especially in the first week of treatment to ensure that NPWT is continually applied and a seal is maintained.**

4. Where PICO dressings are used on infected wounds, more frequent dressing changes may be required. Regular monitoring of the wound should be maintained to check for signs of infection.

5. If deemed clinically appropriate, care should be taken that the application of a circumferential dressing or the use of negative pressure wound therapy on ischaemic limbs does not compromise circulation.

6. Although PICO dressings can be used under clothing/bedding, it is important that occlusive materials e.g. film dressings, are not applied over the pad area of the dressing as this will impair the intended evaporation of moisture through its outer layer.

7. The PICO dressing should not be covered by rigid immobilization devices or casts which might apply excessive pressure and cause tissue injury at the wound site, especially where the tubing enters the dressing.

8. Prolonged placement of rigid or opaque materials over the PICO dressing may prevent the regular inspection and assessment of the wound, and disrupt scheduled or required dressing changes.

9. Where PICO dressings are used on patients with fragile skin, a skin protectant such as SECURA No-Sting Barrier Film should be used on areas of skin where fixation strips are to be applied. Inappropriate use or repeated application of fixation strips may otherwise result in skin stripping.

10. Do not use PICO dressings with oil-based products such as petrolatum as it may compromise establishing

11. The use of negative pressure wound therapy presents a risk of tissue ingrowth into foam when this is used as a wound filler. When using foam filler with PICO, tissue ingrowth may be reduced by using a non-adherent wound contact layer or by increasing the frequency of dressing changes.

12. PICO may be used in conjunction with surgical drains provided the dressing is not placed over tubing where it exits the skin. Any surgical drain should be routed under the skin away from the edge of the dressing and function independently of the PICO System.

13. When showering, the PICO pump should be stopped and disconnected from the dressing. Whilst disconnected, ensure the end of the tubing attached to the dressing is facing down so that water does not enter the tube.

14. The PICO dressing should only be used with PICO pumps.

15. Do not alter or cut tubing configuration or pull on the tubing or soft port.

16. Do not cut the PICO dressing pad as this may lead to loss of negative pressure wound therapy application.

17. Always ensure that the PICO dressing is positioned centrally over the wound. The soft port should be positioned uppermost on intact skin and not extend over the wound so that the risk of fluid collecting around the soft port and potentially blocking the therapy is minimized.

18. The PICO system is single use only. Use of any part of this system on more than one patient may result in cross contamination that may lead to infection.

19. High temperatures and humidity may reduce wear time of PICO dressings.

When applying dressings next to another, ensure the dressing borders do not overlap.

Adverse reactions

Excessive bleeding is a serious risk associated with the application of suction to wounds which may result in death or serious injury. Careful patient selection, in view of the above stated contraindications, warnings and precautions is essential. Carefully monitor the wound and dressing for any evidence of a change in the blood loss status of the patient. Notify the healthcare professional of any sudden or abrupt changes in the volume or the colour of exudate.

8. Instructions for use

8.1. Guidance on wound suitability

PICO dressings should be used on wounds which fit comfortably within the area of the pad, observing precautions on soft port positioning (on intact skin and not extending over the wound). PICO Multisite dressings are designed to enhance conformability when dressing awkward anatomical areas. PICO dressings (including Multisite) must be used to dress one wound only.

Depth - Wounds greater than 0.5cm (1/4in) in depth are likely to require a foam or gauze negative pressure wound therapy filler to ensure adequate treatment of all the wound surfaces. Wounds greater than 2cm (3/4in) in depth must be treated with the use of a filler along with the PICO 14 system to ensure adequate contact with the wound. **Exudate** - PICO is intended for use on wounds where the level of exudate is low (up to 0.6g of liquid exudate/cm² of wound area/24 hours) to moderate (up to 1.1g of liquid exudate/cm² of wound area/24 hours). 1g of exudate is approximately equal to 1ml of exudate. When used on a moderately exuding wound, the size of the wound should generally be no more than 25% of the dressing pad area.

8.2. Application

The dressing should only be applied by a healthcare professional.

1. Remove any excess hair to ensure close approximation of the dressing to the wound. If necessary, irrigate the wound with sterile saline and pat the wound dry.
2. Using a clean technique, peel off the first release handle and place the dressing centrally over the wound to reduce the chance of wound fluid coming into contact with the soft port. Ensure the dressing lies flat to the wound and the surrounding skin. The dressing should be applied with the soft port positioned higher than the wound (depending on the patient's primary position), placed on intact skin and not extending over the wound to prevent fluid pooling around the soft port and blocking the therapy. Remove the other remaining handle(s) and smooth the dressing around the wound to prevent creasing. Reposition if required to ensure border is not creased.

Warning: The PICO pump contains a MAGNET.

Keep the PICO pump at least 10cm (4 inches) away from other medical devices at all times.

(See Section 5.1 Magnet Warning)

3. Once the dressing is in place, remove the pump and the batteries from the tray. The direction in which the batteries should be placed is inside the battery compartment. Insert the batteries. Replace the cover. Following this all four indicators should illuminate for 3 seconds.
4. Join the pump and dressing tubes by twisting together the connectors. The soft tube can be directly connected to the pump if long tubing is not required.


Press the orange button to start the application of NPWT. For any information or troubleshooting refer to

pective PICO system user manual.

If using SECURA No-Sting Barrier Film prior to application of the fixation strips (see Precaution 9) wipe the area surrounding the dressing and allow skin to dry.

6. Apply the fixation strips all the way around the dressing border. Remove top carrier on the strip after each one has been applied. These strips maintain the seal over the wear time of the dressing. In awkward areas, it may be useful to apply the strips to help achieve a seal prior to switching on the pump. Place each strip so that it overlaps the dressing border by approximately 1cm.


Ensure tubing is not twisted or trapped between clothing. Please note that if at any time the fixation strips are removed, the dressing should also be replaced. If desired, gel patches may be applied in addition to the fixation strips to help achieve or maintain a seal.

8.3. Dressing change

Dressings should only be changed by a healthcare professional

1. Dressings should only be changed in line with standard wound management guidelines, typically every 3-4 days. At the healthcare professional's discretion a PICO dressing may be left in place for up to 7 days. The orange dressing full indicator on the pump will flash if it detects a full dressing or blocked filter. More frequent dressing changes may be required depending on the level of exudate, condition of the dressing or other patient considerations; e.g. when PICO is used on infected wounds.

2. Inspect the PICO dressing regularly. If the dressing appears ready for changing (see diagrams A-C), press the orange button and disconnect the dressing from the pump. The fixation strips should be stretched away from the skin and the dressing lifted at one corner and peeled back until it has been fully removed. Apply another dressing as per section 8.2. Application, connect to the pump and press the orange button to reinitiate the therapy.

(A) Dressing properly positioned and is acceptable to be left in place

(B) Dressing requires change – Port may block with fluid

(C) Dressing requires change – Absorbent area is full

3. The PICO dressing should be disposed of as clinical waste.

8.4. Use of PICO dressings with fillers

PICO dressings are compatible with standard gauze and foam fillers used in traditional NPWT where this is clinically appropriate – for example on a defect wound. When filler is used, the filler and the PICO dressing should be changed 2 to 3 times a week, according to local clinical protocol and manufacturer's instructions. The filler should loosely fill to the surface of the wound. Avoid over-packing.

8.5. Use of PICO dressings with non-adherent layers

dressings may be used over the top of a non-adherent layer if required, for example over a skin on infected wounds or wounds at risk of infection, ACTICOAT™ Flex may be used under PICO dressings.

Use of PICO with compression

ICO is compatible with compression therapy if clinically required; however, the wound should first be accurately assessed to ensure the exudate is of low to moderate volume. Following the first application of ICO dressing under a compression system, check the dressing after 2-3 days to see if fluid handling is acceptable or if a dressing change is required and that it is appropriate to continue with compression therapy.

When applying compression therapy over PICO dressings, ensure tubing connecting the PICO pump and dressing runs over the top of the first layer and beneath subsequent layers to avoid it being held against the skin. If using compression hosiery, run tubing outside of stockings. Do not store the pump within layers of compression. For guidance on the correct application of compression therapy, refer to the relevant manufacturer's instructions.

8.7. Dressing removal

Where a healthcare professional believes it to be clinically appropriate, the patient may be instructed to remove the PICO dressing at the end of therapy. Refer to PICO system user manual for patient dressing removal instructions.

Follow these instructions for dressing removal:

1. Stop the PICO pump by pressing the orange button - All indicators will turn off.
2. Remove the pump from the dressing by untwisting the connectors.
3. Remove the PICO dressing by stretching the fixation strips away from the skin. Lift the dressing at one corner and peel back until it has been fully removed.
4. Dispose of the dressing and fixation strips as clinical waste.

1. Showering and bathing

ight showering is permissible; however, the PICO pump should be stopped and disconnected (see section i. Precautions) and placed in a safe location where it will not get wet. The PICO dressing should not be exposed to a direct spray or submerged in water. While disconnected, ensure the end of the tubing attached to the dressing is facing down so that water does not enter the tube. After showering and bathing the PICO pump should be reconnected to the dressing and restarted by pressing the orange button.

.0. General use

During use of your PICO dressing, for any information or troubleshooting on your pump refer to the respective PICO system user manual.

1. Cautions

This user manual is not intended as guarantee or warranty. It is only intended as a guide. For medical questions please consult a physician. The product must be used in accordance with this user manual and all applicable labelling.

ICO is packed in the UK with individual components being made in the following countries:

Dressing – UK

ixation strips – Origin as marked

K Smith & Nephew, Croxley Park, Building 5, Hatters Lane, Watford, Hertfordshire WD18 8YE. Tel: 08000 57 573.

USTRALIA Smith & Nephew Pty Ltd., 85 Waterloo Road North Ryde NSW 2113. Tel +61 2 9857 3999.

SMITH & NEPHEW GmbH, Concorde Business Park C3, A-2320 Schwechat. Tel: +43 (0)1 707 91 02.

INDONESIA PT Smith & Nephew Healthcare, Total Building, 5th Floor, J1. Let. Jend. S. Parman Kav 106A, Jakarta 11440. Tel: 21 5680521.

NETHERLANDS Smith & Nephew Nederland CV, Bloemlaan 2, 2132 NP Hoofddorp, The Netherlands Tel: (+31) (0) 20 6543999.

SINGAPORE Smith & Nephew Pte Ltd., 29 Media Circle, #06-05 Alice @ Mediapolis, Singapore 138565 Tel: (+65) 6270 0552

SWITZERLAND Smith & Nephew Schweiz AG, Advanced Wound Management, Oberneuhofstrasse 10d, 6340 Baar. Tel: +41 41 766 22 44.


Smith & Nephew Medical Limited,
101 Hessle Road, Hull, HU3 2BN England
Contact details for other countries can
be found at: www.smith-nephew.com

Trade Marks of Smith & Nephew
©Smith & Nephew
Date of issue 10/2020
This product may be covered by one or more
patents. See www.smith-nephew.com/patents

- 12. System variants**
- PICO Multipack of Single Dressings for use with PICO single use Negative Pressure Wound Therapy Systems**
- | | |
|---|----------|
| 5 x individually packaged dressings, secondary fixation strips | |
| 10cm x 20cm | 66802022 |
| 10cm x 30cm | 66802023 |
| 10cm x 40cm | 66802024 |
| 15cm x 15cm | 66802025 |
| 15cm x 20cm | 66802026 |
| 15cm x 30cm | 66802027 |
| 20cm x 20cm | 66802028 |
| 25cm x 25cm | 66802029 |
| Multisite 15cm x 20cm | 66802020 |
| Multisite 20cm x 25cm | 66802021 |

13. Glossary of symbols

	Warning magnetic field. Pump contains a MAGNET. You must keep the pump at least 10cm away from implantable devices such as pacemakers, defibrillators, shunt valves, neurostimulators, cochlear implants, or other life sustaining equipment such as insulin pumps
	Caution: Federal (USA) law restricts this device to sale by or on order of a physician
	Single Use. Do not reuse
	EU: Not for general waste

STERILE EO

Product is sterilized by Ethylene Oxide

Caution

Do not use if the package is damaged

Storage temperature

Manufacturer

Date of manufacture

Consult instructions for use

Keep product out of sunlight

Keep dry

LOT

Lot number

CE Mark

ADDENDUM TO MEDICAL DEVICE INSTRUCTIONS FOR USE

"PICO dressings for single use negative pressure wound therapy system, variations"

delivering the device to the Russian Federation, the information contained in this addendum to the instructions takes priority. The addendum to medical device instructions for use is an integral part of the instructions for

MD NAME, MODELS OR VARIATIONS OF MD

The following MD names used in the text of related documents are identical to the name of the medical device "PICO Dressings for Single Use Negative Pressure Wound Therapy System, variations":

- PICO Dressing Packs/ PICO Dressing Kits
- PICO Dressings
- Dressings
- Device (s)

The names of medical device can be mentioned in the text in abbreviated form for the convenience of perceiving information.

PICO dressings for single use negative pressure wound therapy system, variations:

1. PICO Dressing 10cm x 20cm, contents:
 - Dressing – 5 pcs.
 - Fixation strips – 30 pcs.
 - IFU – 8 pcs.
 - Addendum to IFU – 8 pcs.
2. PICO Dressing 10cm x 30 cm, contents:
 - Dressing – 5 pcs.
 - Fixation strips – 30 pcs.
 - IFU – 8 pcs.
 - Addendum to IFU – 8 pcs.
3. PICO 10cm x 40xm, contents:
 - Dressing – 5 pcs.
 - Fixation strips – 30 pcs.
 - IFU – 8 pcs.
 - Addendum to IFU – 8 pcs.
4. PICO 15cm x 15xm, contents:
 - Dressing – 5 pcs.
 - Fixation strips – 30 pcs.
 - IFU – 8 pcs.
 - Addendum to IFU – 8 pcs.
5. PICO 15cm x 20cm, contents:
 - Dressing – 5 pcs.
 - Fixation strips – 30 pcs.
 - IFU – 8 pcs.
 - Addendum to IFU – 8 pcs.
6. PICO 15cm x 30xm, contents:
 - Dressing – 5 pcs.
 - Fixation strips – 30 pcs.
 - IFU – 8 pcs.
 - Addendum to IFU – 8 pcs.
7. PICO 20cm x 20cm, contents:

Dressing – 5 pcs.
 Fixation strips – 30 pcs.
 IFU – 8 pcs.
 Addendum to IFU – 8 pcs.

PICO 25cm x 25cm, contents:

Dressing – 5 pcs.
 Fixation strips – 30 pcs.
 IFU – 8 pcs.
 Addendum to IFU – 8 pcs.

9. PICO Multisite 15cm x 20cm, contents:

Dressing – 5 pcs.
 Fixation strips – 30 pcs.
 IFU – 8 pcs.
 Addendum to IFU – 8 pcs.

10. PICO Multisite 20cm x 25cm, contents:

Dressing – 5 pcs.
 Fixation strips – 30 pcs.
 IFU – 8 pcs.
 Addendum to IFU – 8 pcs.

MD Classification

Local Russian Nomenclature classification type in acc. with MoH 4n Order – 287470

Potential risk class in acc. with Directive EEC of MD 93/42 – IIb

Frequency of use – The product is single use only.

Product classification by duration of contact - The product has a prolonged (<30 days) contact

Sterility – sterile.

Invasiveness – Non-invasive medical device

2 INTENDED USE

Intended for use with PICO single use negative pressure wound therapy (NPWT) system

3 TECHNICAL SPECIFICATIONS

See section 12 of the Instructions for Use.

Name	Description/ Component characteristics (Dimensions tolerance $\pm 5\%$, unless otherwise specified)	
Dressing weight	Dressing dimensions WxL (cm)	Dressing weight (g)
	10 x 20	15 $\pm$ 10%
	10 x 30	21 $\pm$ 10%
	10 x 40	27 $\pm$ 10%
	15 x 15	17 $\pm$ 10%
	15 x 20	22 $\pm$ 10%
	15 x 30	33 $\pm$ 10%
	20 x 20	29 $\pm$ 10%
	25 x 25	49 $\pm$ 10%
	Multisite 15 x 20	21 $\pm$ 10%
	Multisite 20 x 25	35 $\pm$ 10%
Dressing	Designed to evacuate and evaporate exudate from the wound site in NPWT conditions	
	Dressing dimensions WxL (cm)	Absorbent layer – maximum dimensions WxL (cm)
	10x20	5.6x15
	10x30	5.6x25
	10x40	5.6x35
	15x15	10x10
	15x20	10x15
	15x30	10x20

	20x20	15x15
	25x25	20x20
Multisite Dressing	Designed to conform to curved areas of the body	
	Dressing dimensions WxL (cm)	Absorbent layer – maximum dimensions WxL (cm)
	Multisite 15x20	10x15
	Multisite 20x25	15x20
Dressing release handle	Designed to protect the adhesive layer of the dressing during storage and transportation	
	Dressing dimensions WxL (cm)	Dressing release handle dimensions WxL (cm)
	10x20	10x20
	10x30	10x30
	10x40	10x40
	15x15	15x15
	15x20	15x20
	15x30	15x30
	20x20	20x20
	25x25	25x25
	Multisite 15x20	15x20
	Multisite 20x25	20x25
Secondary fixation strip	Designed to maintain the seal over the wear time of the dressing	
	Dressing dimensions WxL (cm)	Secondary fixation strip dimensions LxW (cm)
	10x20	25x5
	10x30	
	10x40	
	15x15	
	15x20	
	15x30	
	20x20	
	25x25	
	Multisite 15x20	
	Multisite 20x25 /	
Secondary fixation strip release handle	To be removed before applying the secondary fixation strip	
	Dressing dimensions WxL (cm)	Release handle dimensions LxW (cm)
	10x20	14x5 (central part) 10,2 x 5 (side parts)
	10x30	
	10x40	
	15x15	
	15x20	
	15x30	
	20x20	
	25x25	
	Multisite 15x20	
	Multisite 20 x25	
Secondary fixation strip fixing film	Fixing the dressing directly	
	Dressing dimensions WxL (cm)	Fixing film dimensions LxW (cm)
	10x20	25x5
	10x30	
	10x40	
	15x15	
	15x20	
	15x30	
	20x20	
	25x25	
	Multisite 15x20	
	Multisite 20 x25	
Secondary fixation strip top carrier	To be removed after applying the secondary fixation strip	
	Dressing dimensions WxL (cm)	Top carrier dimensions LxW (cm)
	10x20	26x5
	10x30	
	10x40	

	15x15	
	15x20	
	15x30	
	20x20	
	25x25	
	Multisite 15x20	
	Multisite 20 x25	
Secondary fixation strip top carrier handle bars	Designed to remove the secondary fixation strip top carrier	
	Dressing dimensions WxL (cm)	Secondary fixation strip top carrier handle bar dimensions WxL (cm)
	10x20	1x5
	10x30	
	10x40	
	15x15	
	15x20	
	15x30	
	20x20	
	25x25	
	Multisite 15x20	
	Multisite 20x25	
Soft port	Designed to maintain negative pressure in the wound by connecting the dressing and the pump (either directly or using extension tubing)	
	Dressing dimensions WxL (cm)	Soft port length (cm)
	10x20	33.3 +0.9/-1.1
	10x30	
	10x40	
	15x15	
	15x20	
	15x30	
	20x20	
	25x25	
	Multisite 15x20	
	Multisite 20 x25	
Dressing tubing Male connector	Designed to connect dressing soft port and either extension tubing or pump connector	
	Dressing tubing length (mm): 37 Dressing tubing inner diameter (mm): 1,5 Dressing tubing outer diameter (mm): 3 Connection dimensions of the male connector (mm): - Connector tip diameter (mm): 4,5 - External thread diameter (mm): 9	

Dressing	
Appearance	An adhesive, absorbent dressing with a baking film is free from foreign bodies. No colour leaching.
Absorbency	50 g/100 cm ² minimum
Absorbency (compressed)	50 g/100 cm ² minimum
Fluid Handling Capacity (Absorbtion)	40 g/100 cm ² minimum
Fluid Handling Capacity (Permeability)	80 g/100 cm ² minimum
Extensibility (Machine Direction)	1.4 kgf/cm maximum (1,373 N/mm)
Extensibility (Transverse Direction)	1.4 kgf/cm maximum (1,373 N/mm)
Waterproofness	Complies (inc. filter/top film seal)*
Moisture Vapour Permeability	100 – 3500 g/m ² /24hrs
Perforation size	0.4 – 1.8 mm
Perforation Pattern	Holes to be distributed. 8 holes/cm ² minimum
Port-tube Tensile Strength (reverse pull)	Holds 15 N minimum (for 15 sec)
Tube connector Tensile Strength	7N minimum disconnection force

*The method is based upon a pharmacopoeia method, contained in BP 1993 Appendix XX method K, and as it is a pass/fail test based on a simple physical characteristic carried out under controlled conditions formal validation is not considered necessary.

Fixation strips	
Appearance	The dressing shall be free from foreign bodies. A film pre-lift must be present on either side of the carrier to facilitate the removal of the carrier. No excessive delamination of the carrier must be visible following removal of central handle
Weight of Base Film	40 – 50 g/m ²
Weight of Adhesive	27 – 43 u/m ²
MVP	600 g/m ² /24hrs minimum
Extensibility	0.2 kgf/cm maximum (0.196 N/mm)
Permanent Set	5 % maximum
Tensile Strength	0.5 kgf/cm minimum (0.490 N/mm)
Adhesion to Steel	150-800 gf/cm (0.147-0.785 N/mm)
Fim/Carrier Stripping load	1.0 – 6.0 gf/cm (0.000981-0.005884 N/mm)
Adhesive/Liner Stripping Load	1.0gf/cm minimum (0.981 N/mm)
Shear Adhesion	2.5 mm maximum
Waterproofness	Complies*
Strip width	Width tolerance: 50 mm ±1 mm

*The method is based upon a pharmacopoeia method, contained in BP 1993 Appendix XX method K, and as it is a pass/fail test based on a simple physical characteristic carried out under controlled conditions formal validation is not considered necessary.

4 LIST AND DESCRIPTION OF MD MATERIALS

Device materials:

Component		Component material	Brand/CAS No/material characteristics	Type of contact with the patient/person
Dressing	Release handle	Semi-transparent , one-side siliconised, hogh density polyethylene film	Mondi HDPE CT(70tHDut) U00B6	Prolonged contact (no more then 7 days) with breached or compromised surfaces
	Wound contacting adhesive		Silpuran 2111 A/B Part A: ▪ Polydimethylsiloxane, vinyl terminated CAS-No: 68083-19-2 ▪ MQ resin CAS No: 104199-38-4 ▪ Platinum catalyst (<50ppm) CAS No: 68478-92-2 ▪ <0.02% toluene (impurity) Part B: ▪ Polydimethylsiloxane, hydrogen terminated CAS No: 70900-21-9. ▪ Polydimethylsiloxane, vinyl terminated CAS-No: 68083-19-2 ▪ Dimethylsiloxane-methylhydrogensiloxanecopolymer, trimethylsilyl-terminated. CAS No: 68037-59-2	
	Wound contact layer	Perforated polyurethane	Elastollan SP806	
		Antiblock masterbatch	Konzentrat 978	
	Constructional Adhesive	A8 acrylic adhesive	▪ 14.661% (2-EHA) 2-ethylhexyl acrylate ▪ 1.892% glacial acrylic acid ▪ 14.661% (n-BA) n-butyl acrylate ▪ 68.624% acetone ▪ 0.162% Di (4-tert-butylcyclohexyl) peroxydicarbonate (known as BCHPC)	
		K5 acrylic adhesive	▪ 22.15% n-butyl acrylate	

			▪ 24.72% 2-ethylhexyl acrylate ▪ 0.50% 2-hydroxyethyl methacrylate ▪ 2.50% n-butyl methacrylate ▪ 49.90% deionised water ▪ 0.10% Ginasul 46P Surfactant ▪ 0.13% Ammonium persulphate	
	Spacer layer	Polyester	Baltex, thickness 2 mm	
	Absorbent layer	Superabsorbent	Airlaid Chem-Posite 11C-450, thickness 3-5 mm	
	Top film	EU33 Polyurethane film	SP9109	
		Antiblocking agent	Konzentrat 978	
Dressing soft port	Top film	EU94 High strength polyurethane film	Elastollan 685A 10N	Prolonged contact (no more than 7 days) with breached or compromised surfaces
		Antiblocking agent	Konz 927 EBS free	
	Constructional Adhesive	UV curable adhesive	Loctite 3341	
		Curable adhesive	Loctite 4011-SG	
	Filter membrane	Acrylic copolymer on a nylon support	Versapor R	
	Spacer layer	Polyester	Baltex, thickness 2 mm	
	Dressing tubing	PVC	Unichem 7088G-015 MED	
	Soft Port Dressing Connector	Aromatic thermoplastic polyurethane	Desmopan 9370A	
Secondary fixation strips	Release handle	Polyethylene wood-free bleached paper with a solventless silicon coating	Loparex 111/60	Limited (up to 24 h) contact with intact skin
		Orange print	PMS 151	
	A8 acrylic adhesive		▪ 14.661% (2-EHA) 2-ethylhexyl acrylate ▪ 1.892% glacial acrylic acid ▪ 14.661% (n-BA) n-butyl acrylate ▪ 68.624% acetone ▪ 0.162% Di (4-tert-butylcyclohexyl) peroxydicarbonate (known as BCHPC)	Prolonged contact (no more than 7 days) with intact skin
	Top film	Polyurethane film	SP9109 resin	
		Antiblocking agent	Konzentrat 978	
	Top carrier	Polyester film with silicone and matt polyethylene	Loparex 111/23	Limited (up to 24 h) contact with intact skin
		Orange print	PMS 151	

Packing materials:

Component	Component material	Brand or characteristics of material
Individual packaging (pouch) materials :		
Top web	Paper	KG9111
Base Web	Nylon	NLS90
Secondary packaging (Carton Box) materials :		
Carton Box	Carton	650 micron

5 USE, TRANSPORTATION AND STORAGE CONDITIONS OF MD

Use conditions

Operational temperature	+5-+40°C
Relative humidity	10-95%

Transportation and storage conditions

Operational temperature

not exceeding 25°C

Relative humidity

not exceeding 60%

LABELLING INFORMATION**1.1 Draft of the labeling**






Draft of labeling medical device contains following information in Russian:

Article/catalogue number:**Name:** PICO dressings for single use negative pressure wound therapy system, variations**Size*:** 10x20 /10x30 /10x40 /15x15 /15x20 /15x30/ 20x20/ 25x25/ Multisite 15x20/ Multisite 20x25**Pieces per package:** 5 pcs.**Date of manufacturing, use by, lot/batch number:** see on the package**Name and address of manufacturer:** Smith & Nephew Medical Ltd., 101 Hessle Road, Hull, HU3 2BN, United Kingdom**Information about authorized representative in the Russian Federation:** Smith & Nephew LLC, 1, 2nd Syromyatnichesky Lane, 9 floor, prem.1, room 1, Moscow, 105120, Russia. Tel: +7 (495) 755-55-03**Sterilization method:** Sterilized using ethylene oxide. Sterile**Registration Certificate Number:****Indications, mode of application, contraindications:** see IFU**Storage conditions:** in a dry, dark place at a temperature not exceeding 25 ° C.**Precautions:** -**Symbol of conformity to Russian standards (if necessary).****Instruction for use:** <http://www.smith-nephew.com/about-us/where-we-operate/europe/russia/ifus/>

* There will only be one size listed for this article/catalogue number

6.2 Labelling data

Symbol/ Text	Meaning
	Peel off the release handle of the secondary fixation strip first
	Then peel off the top carrier of secondary fixation strip
PICO° Dressings	PICON dressings
	Follow instructions for use
REF	REF number
	Use by
	Turn away the paper edge to open
Multisite 20cm x 25cm 7.9in. x 9.8in.	Dressing dimensions: Multisite 20cmx25cm; 7,9 In. x 9,8 In. (according to variation)

	Dressing dimensions: Length 25 cm(9,8 in.) – 20 cm (7,9 in.) Width 20 cm(7,9 In.) - 15 cm (5,9 In.) (according to variation)
	Kit contains 6 fixation strips
	Device schematic IFU (PICO pouch)
	Smith & Nephew company logo
PICO[®]	PICO
Packaged in UK with individual components made in the following countries Dressing – UK Fixation strips - UK	-
Smith & Nephew Medical Limited, 101 Hessle Road, Hull HU3 2BN England Trade Mark of / Marque de commerce de Smith & Nephew This Product may be covered by one or more US Patents See: www.smith-nephew.com/patents	-
	European authorised representative
	Not intended for reuse
	QR-code
	Barcode
	The packaging material is recyclable
	Kit contains 5 PICO dressing

7 MD DISPOSAL RULES AND CONDITIONS

The disposal of the packaging/contents shall be carried out in a safe and permitted manner in compliance with the requirements of regional/national/international regulations concerning the treatment of medical waste.

Procedure of disposal or destruction of used medical product-disposal is carried out in accordance with the legislation in force in the Russian Federation. Class of waste on SanPiN 2.1.3684-21 "Б"

Procedure of disposal or destruction of unused medical product-disposal is carried out in accordance with the legislation in force in the Russian Federation. Class of waste on SanPiN 2.1.3684-21 "А"

INFORMATION ABOUT THE STERILE CONDITION OF THE MEDICAL DEVICE, THE METHOD OF ITS STERILIZATION AND THE PROCEDURE FOR ACTION IN CASE OF VIOLATION OF THE STERILE PACKAGING

Medical device is sterilized with ethylene oxide.

Do not use if the packaging is damaged.

9 DEVICES INCORPORATING NON-VIABLE MATERIALS OF ANIMAL ORIGIN

Not applicable.

10 DEVICES INCORPORATING A MEDICINAL SUBSTANCE

Smith & Nephew Medical Ltd. products do not incorporate a medicinal substance, nor is it present in the manufacturing process or present in the final device, therefore this section does not apply

11 REFERENCES TO PREVIOUS MODIFICATIONS OF THE MEDICAL PRODUCT OR SIMILAR MODIFICATIONS OF MEDICAL DEVICES

Not applicable

12 DESCRIPTION OF ACCESSORIES, MEDICAL DEVICES OR NON-MEDICAL DEVICES INTENDED FOR USE IN COMBINATION WITH THE CLAIMED MEDICAL DEVICE

See section 8 of the Instructions for use.

Standard foam fillers:

Length 125±5 mm

Width 100±5 mm

Thickness 15±5 mm

Pore size: 40-50 PPI

Density: 1.30-1.60 lbs/ft³ (20.82-25.63 kg/m³)

Tensile strength: min 10 psi (0.069 N/mm²)

Elongation: min 110 %

Tensile strength: min 2.5 lb/in (0.438 N/mm)

Standard gauze fillers:

Surface density not less than 36 g/m

pH aqueous extract - 6-8

Breaking force:

- warp – not less than 78,0 N (8,0 kgs),

- weft - not less than 34,0 N (3,5 kgs),

Wetting speed – not more than 10 sec

Foam and gauze fillers are not supplied with the PICO dressings.

Compatible PICO devices:

- PICO single use Negative Pressure Wound Therapy System (registration certificate № RZN 2013/1110),

- PICO 7 single use Negative Pressure Wound Therapy System, (registration certificate № RZN 2013/1110),

- PICO 7Y single use Negative Pressure Wound Therapy System, variations (registration certificate № RZN 2013/1110),

- PICO14 single use Negative Pressure Wound Therapy System, variations (the product is in the process of registration)

13 LIST OF NATIONAL STANDARDS APPLIED BY THE MANUFACTURER OF A MEDICAL DEVICE

Standard	Description
Packaging	
BS EN ISO 11607-1: 2009 + A1:2014	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
BS EN ISO 11607-2: 2017	Packaging for terminally sterilized medical devices - Part 2:

	Validation requirements for forming, sealing and assembly processes
BS EN ISO 780:2015	Packaging. Distribution packaging. Graphical symbols for handling and storage of packages
ISTA 2A:2011	Partial Simulation Performance Tests - ISTA 2A tests for packaged-products 150 lb. (68 kg) or less
Device applicable standards	
BS EN 13726-1:2002/AC:2003	Test methods for primary wound dressings. Aspects of absorbency
BS EN 13726-2: 2002	Test methods for primary wound dressings. Moisture vapour transmission rate of permeable film dressings
BS EN 13726-3: 2003	Test methods for primary wound dressings. Waterproofness
BS EN 13726-4: 2003	Test methods for primary wound dressings. Conformability
Biocompatibility	
BS EN ISO 10993-1:2018	Biological evaluation of medical devices. Part 1. Evaluation and testing within a risk management process
BS EN ISO 10993-2:2006	Biological Evaluation of Medical Devices - Part 2: Animal Welfare Requirements
BS EN ISO 10993-3:2014	Biological evaluation of medical devices. Part 3. Tests for genotoxicity, carcinogenicity and reproductive toxicity
BS EN ISO 10993-5:2009	Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity
BS EN ISO 10993-7:2008	Biological evaluation of medical devices. Part 2. Ethylene oxide sterilization residuals
BS EN ISO 10993-10:2013	Biological evaluation of medical devices -- Part 10: Tests for irritation and skin sensitization
BS EN ISO 10993-11:2009	Biological evaluation of medical devices. Part 11: Tests for systemic toxicity
BS EN ISO 10993-12:2012	Biological evaluation of medical devices. Part 12. Sample preparation and reference materials
BS EN ISO 10993-18:2009	Medical devices. Biological evaluation of medical devices. Part 18. Chemical characterization of materials
Restriction of Hazardous Substances (RoHS) Directive - 2011/65/EU	Directive 2011/65/EU. Restriction of hazardous substances directive
ISO 17075-1:2017	Leather — Chemical determination of chromium(VI) content in leather — Part 1: Colorimetric method
Sterility/Microbiology	
EN 556-1: 2001 + AC:2006	Sterilization of medical devices - Requirements for medical devices to be designated 'STERILE' - Part 1: Requirements for terminally sterilized medical devices
BS EN ISO 11737-1 2018	Sterilization of medical devices - Microbiological methods - Part 1: Determination of a population of microorganisms on products
BS EN ISO 11737-2:2009	Sterilization of medical devices. Microbiological methods. Part 2. Tests of sterility performed in the definition, validation and maintenance of a sterilization process
BS EN ISO 11135:2014	Sterilization of health-care products. Ethylene oxide. Part 1. Requirements for the development, validation and routine control of a sterilization process for medical devices
BS EN ISO 13485:2016	Medical devices. Quality management systems. Requirements for regulatory purposes
BS EN ISO 9001:2015	Quality management systems. Requirements
BS EN ISO 14971:2019	Medical devices. Application of risk management to medical devices
BS EN ISO 14644-1:2015	Cleanrooms and associated controlled environments. Part 1. Classification of air cleanliness by particle concentration
Labeling	
BS EN ISO 1041: 2008+A1:2013	Symbols to be used with medical device labels, labelling and information to be supplied-Information supplied by the manufacturer of medical devices.
BS EN ISO 15223-1:2016	Symbols to be used with medical device labels, labelling and information to be supplied Part 1: General requirements
14 REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD	
Не применимо.	Not applicable.
15 CHECKING THE FUNCTIONALITY AND SAFETY OF MD	

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Not applicable.

WARRANTY AND CLAIM

Importing organization: Smith & Nephew LLC; 1, 2nd Syromyatnichesky Lane, 9 floor, prem.1, room 1, Moscow, 105120, Russia; Tel: +7 (495) 755-55-03

WARRANTY: The manufacturer guarantees compliance with the stated characteristics, terms preservation of sterility, safety and integrity of packaging with proper conditions of transportation, operation and storage of the medical device for 2 years since manufacturing/sterilization. Shelf life - 24 months from the date of manufacture /sterilization indicated on the package, provided that it is stored in its original packaging, in a dry, well-ventilated place, protected from direct sunlight, at a temperature not exceeding 25 ° C.

The guarantee for the medical device **"PICO dressings for single use negative pressure wound therapy system, variations"** should cover the materials, performance, and absence of manufacturing faults subject to use for one patient only. Do not reuse.

17 NAME AND LEGAL ADDRESS OF MANUFACTURER

MANUFACTURER:


Smith & Nephew Medical Ltd., 101 Hessle Road, Hull, HU3 2BN, United Kingdom

MANUFACTURING SITE:

1. Smith & Nephew Medical Ltd., 101 Hessle Road, Hull, HU3 2BN, United Kingdom
2. Smith & Nephew Medical (Suzhou) Limited, No. 12, Wuxiang Road Comprehensive Free Zone, West Zone Suzhou Industrial Park Suzhou Jiangsu 215021 China.

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA:

Smith & Nephew LLC, 1, 2nd Syromyatnichesky Lane, 9 floor, prem.1, room 1, Moscow, 105120, Russia
Tel: +7 (495) 755-55-03

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