**HAEMONETICS®**
THE Blood Management Company
Haemonetics Corporation
400 Wood Road
Braintree, MA 02184-9114
Tel: 781.848.7100
www.haemonetics.com

HAEMONETICS CORPORATION
400 WOOD ROAD
BRAINTREE, MA 02184 USA
1-800-225-5242

"APPROVED"
Julia PERLMAN
Julia PERLMAN
INTERNATIONAL Regulatory
Affairs specialist
6-OCTOBER-2014

Postscript to the **MCS®+** - Operation Manual

Yulia PERIMAN
~~Yulia PERIMAN~~
INTERNATIONAL Regulatory Affairs Specialist
6-OCTOBER-2014

HAEMONETICS CORPORATION
400 WOOD ROAD
BRAINTREE, MA 02184 USA
1-800-225-5242

Postscript to the MCS[®]+ - Operation Manual -

Printed in France/USA
Haemonetics Corporation
400 Wood Road
Braintree, Massachusetts 02184, USA

CE 0123
EC REP HAEMONETICS U.K. LTD.
5 Ashley Drive, Bothwell
Scotland G71 8BS
Great Britain

P/N 85270-30, Manual revision: B
June 2006

©2002, 2006, Haemonetics Corporation. All rights reserved.

P/N 85270-30, Manual revision: B

CONSUMER INFORMATION

Proprietary rights

The contents of this guide are property of the Haemonetics Corporation. Haemonetics® and eLynx™ are trademarks or registered trademarks of the Haemonetics Corporation in the United States and/or other countries.

Any information or descriptions contained in this guide may not be reproduced and released to any of the general public, or used in conjunction with any professional instruction without written consent of Haemonetics Corporation, USA. Please direct any written inquiries to the appropriate address.

International Headquarters

Haemonetics SA
Signy Centre, rue des Fléchères
P.O. Box 262, 1274 Signy 2, Switzerland
Tel. [+41 22] 363 90 11
Fax [+41 22] 363 90 54

Corporate Headquarters

Haemonetics Corporation
400 Wood Road
Braintree, MA 02184, USA
Tel. [+1 781] 848 7100
Fax [+1 781] 848 5106

What is the purpose of this guide?

This postscript is intended to supplement information contained in the device and procedure operation manuals as supplied by the Haemonetics Corporation.

Use this postscript in conjunction with instruction and training as supplied by qualified Haemonetics personnel. Haemonetics guarantees its merchandise when correctly used by a properly trained operator. Haemonetics accepts no responsibility for problems resulting from failure to comply with prescriptions as outlined by the company. In addition, the operator must fully understand and implement the local standard operating procedures concerning blood-contaminated material as well as blood products.

It remains the full responsibility of the customer to assess and ensure the safety of any products obtained from Haemonetics prescribed procedures, prior to further application or use. Haemonetics declines any responsibility for choices made by the customer concerning the utilization of these products and by-products.

Table of Contents

Section 1 Listing NOTICES and Relevant HELP Messages

MCS+ SAFETY SYSTEM: SOFTWARE DETECTION S/1-2

MCS+ SAFETY SYSTEM: HARDWARE DETECTION..... S/1-21

Section 2 Using MCS+ Data Acquisition Features

DESCRIBING THE MCS+ COMMUNICATION OPTIONS...: S/2-2

Printer S/2-2

Wireless S/2-2

Network..... S/2-2

Disk S/2-3

None S/2-3

ENTERING PROCEDURE DATA WITH A BAR-CODE READER..... S/2-4

Disposable set list and lot numbers..... S/2-4

Anticoagulant and solution lot numbers S/2-5

Donation, donor and operator code numbers..... S/2-6

ACCESSING STORED PROCEDURE DATA..... S/2-7

Data management menu S/2-7

DOWNLOADING, CLEARING AND PRINTING PROCEDURE DATA. . . S/2-9

Downloading data S/2-9

Clearing data S/2-9

Printing data..... S/2-10

PRESENTING THE WIRELESS FEATURE AND COMPONENTS S/2-11

What does wireless mean..... S/2-11

Describing the wireless device components S/2-11

How does the Haemonetics device interface with the eLynx system? . S/2-12

PRESENTING THE ELYNX SYSTEM S/2-13

Section 1

Listing NOTICES and Relevant HELP Messages

MCS+ SAFETY SYSTEM: SOFTWARE DETECTION	1-2
MCS+ SAFETY SYSTEM: HARDWARE DETECTION.....	1-21

MCS+ SAFETY SYSTEM: SOFTWARE DETECTION

The MCS+ is designed with sophisticated software capable of detecting irregularities in the function of the MCS+ components.

The following table provides a list of the numerical NOTICE messages provided by this software detection system which an operator could encounter while working with the MCS+ device.

Note: Protocol-specific NOTICE messages exist. The operator should consult the relevant MCS+ protocol manual for further information.

#	NOTICE message	HELP message
0	SYSTEM PROBLEM THE COMPUTER HAS ENCOUNTERED AN ERROR.	A CPU error caused the MCS+ to reset. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
5	AIR DETECTED BY THE BLOOD LINE AIR DETECTOR (BLAD).	The Air Detector detects air. 1. Check tubing installation in detector. 2. Remove air bubbles. 3. Resume operation.
6	AIR DETECTED BY THE DONOR LINE AIR DETECTOR 1 (DLAD 1).	The Air Detector detects air. 1. Check tubing installation in detector. 2. Remove air bubbles. 3. Resume operation.
7	AIR DETECTED BY THE DONOR LINE AIR DETECTOR 2 (DLAD 2).	The Air Detector detects air. 1. Check tubing installation in detector. 2. Remove air bubbles. 3. Resume operation.
15 to 24	SYSTEM PROBLEM THE COMPUTER HAS ENCOUNTERED AN ERROR.	A CPU error caused the MCS+ to reset. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
25	CENTRIFUGE PROBLEM THE CENTRIFUGE IS NOT TURNING.	1. Check that bowl is installed properly. 2. Check that Centrifuge rotates freely. 3. Check header of bowl. 4. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
26	CENTRIFUGE PROBLEM THE CENTRIFUGE WAS TURNING TOO FAST.	1. Wait until Centrifuge has stopped. 2. Check Centrifuge for obstructions. 3. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.

#	NOTICE message	HELP message
27	CENTRIFUGE PROBLEM THE CENTRIFUGE VACUUM IS NOT DETECTED.	The vacuum is too low to hold the bowl. 1. Remove bowl from the Centrifuge. 2. Re-seat bowl. 3. Ensure that bowl is properly seated. 4. Resume operation.
28	SYSTEM PROBLEM FLUID HAS BEEN DETECTED IN THE CENTRIFUGE.	A leak may have occurred. This procedure MUST be discontinued! 1. Discontinue the procedure. 2. Note the condition and the lot number of the disposable before disposing of it. 3. Contact Haemonetics Field Service.
29 30	SYSTEM PROBLEM THE COMPUTER HAS ENCOUNTERED AN ERROR.	A CPU error caused the MCS+ to reset. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
31	CENTRIFUGE PROBLEM THE CENTRIFUGE SPEED WAS TOO SLOW.	1. Wait until Centrifuge has stopped. 2. Check that bowl is installed properly. 3. Check Centrifuge for obstructions. 4. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
32 to 36	SYSTEM PROBLEM THE COMPUTER HAS ENCOUNTERED AN ERROR.	A CPU error caused the MCS+ to reset. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
37	AIR DETECTED BY AIR DETECTED BY SYSTEM MONITOR.	The Air Detector detects air. 1. Check tubing installation in detector. 2. Remove air bubbles. 3. Resume operation.
38	SYSTEM PROBLEM THE COMPUTER HAS ENCOUNTERED AN ERROR.	A CPU error caused the MCS+ to reset. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
39	CENTRIFUGE PROBLEM THE CENTRIFUGE SPEED WAS TOO SLOW.	1. Wait until Centrifuge has stopped. 2. Check that bowl is installed properly. 3. Check Centrifuge for obstructions. 4. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
40	PUMP PROBLEM AC PUMP WAS TURNING TOO SLOW.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.

#	NOTICE message	HELP message
41	PUMP PROBLEM BLOOD PUMP WAS TURNING TOO SLOW.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
42	PUMP PROBLEM TRANSFER PUMP WAS TURNING TOO SLOW.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
44	PUMP PROBLEM AC PUMP WAS TURNING TOO FAST.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
45	PUMP PROBLEM BLOOD PUMP WAS TURNING TOO FAST.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
46	PUMP PROBLEM TRANSFER PUMP WAS TURNING TOO FAST.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
48	PUMP PROBLEM AC PUMP WAS TURNING IN THE WRONG DIRECTION.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
49	PUMP PROBLEM BLOOD PUMP WAS TURNING IN THE WRONG DIRECTION.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
50	PUMP PROBLEM TRANSFER PUMP WAS TURNING IN THE WRONG DIRECTION.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
52	PUMP PROBLEM AC PUMP POSITION IS NOT CORRECT.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
53	PUMP PROBLEM BLOOD PUMP POSITION IS NOT CORRECT.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
54	PUMP PROBLEM TRANSFER PUMP POSITION IS NOT CORRECT.	A restriction may have occurred. 1. Check pump for restrictions. 2. Check pump for kinked tubing. 3. Resume operation.
56 57	PUMP PROBLEM THE AC/BLOOD PUMP RATIO IS INCORRECT.	An occlusion may have occurred at the Blood pump or AC pump. 1. Check both pumps for occlusions. 2. Resume operation.

#	NOTICE message	HELP message
58	VALVE PROBLEM THE RED VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
61	SYSTEM PROBLEM THE COMPUTER HAS ENCOUNTERED AN ERROR.	A CPU error caused the MCS+ to reset. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
62	ECV LIMIT THE DONOR ECV LIMIT IS REACHED.	The allowed Extra Corporeal Volume (ECV) limit was reached. A complete Return cycle must be performed. 1. Press RETURN to continue.
63	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
64	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	A Watchdog error caused the MCS+ to reset. 1. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
65	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	A DPM/SPM connector caused the MCS+ to reset. 1. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
66	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	A pressure error caused the MCS+ to reset. 1. Check installation of DPM and SPM. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
67	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	An Air Detector caused the MCS+ to reset. 1. Check tubing installation in detectors. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
68	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	A State Fault caused the MCS+ to reset. 1. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
69	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	An ECV limit caused the MCS+ to reset. 1. Check tubing installation at Blood Pump. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.

#	NOTICE message	HELP message
70	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	A pump ratio caused the MCS+ to reset. 1. Check tubing installation at Blood Pump. 2. Check tubing installation at AC Pump. 3. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
71	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	A pump motor caused the MCS+ to reset. 1. Check tubing installation at all pumps. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
72	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	A prime error caused the MCS+ to reset. 1. Check tubing installation at all pumps. 2. Check proper installation of disposable. 3. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
73	SYSTEM PROBLEM THE SAFETY SYSTEM HAS ENCOUNTERED AN ERROR.	A power interruption caused the MCS+ to reset. 1. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
74 to 94	SYSTEM PROBLEM THE COMPUTER HAS ENCOUNTERED AN ERROR.	A CPU error caused the MCS+ to reset. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
100	CENTRIFUGE PROBLEM CENTRIFUGE SAFETY SWITCH NOT OPERATIONAL.	The centrifuge safety switch must be operational to operate the MCS+. 1. Power off, and try again. NOTE: If problem persists, call Haemonetics Field Service.
101	PROCESS PROBLEM RETURN CYCLE TAKING TOO LONG.	Emptying bowl took longer than expected. An occlusion may have occurred. 1. Check tubing at pumps and bowl. 2. Check bowl. 3. Resume operation.
102	PROCESS PROBLEM BLAD DETECTED EARLIER THAN EXPECTED.	Emptying bowl took less time than expected. 1. Check the tubing at the BLAD. 2. Check the bowl for occlusions. 3. Check tubing for sources of air. 4. Clear air from BLAD. 5. Resume operation.

#	NOTICE message	HELP message
103	AC FLOW ERROR LOW AC FLOW DETECTED BY DRIP MONITOR.	1. Check that the AC bag contains fluid. 2. Check AC tubing for occlusions. 3. Check AC tubing for clamps. 4. Check Drip Chamber installation. 5. Check Drip Chamber fluid level. 6. Resume operation.
104	AC FLOW ERROR HIGH AC FLOW DETECTED BY DRIP MONITOR.	1. Check tubing at AC Pump. 2. Check Drip Chamber installation. 3. Check Drip Chamber fluid level. 4. Resume operation.
105	AC FLOW ERROR NO AC FLOW DETECTED BY DRIP MONITOR.	1. Check that the AC bag contains fluid. 2. Check donor tubing for occlusions. 3. Check donor tubing for closed clamps. 4. Check AC tubing for occlusions. 5. Check AC tubing for closed clamps. 6. Resume operation.
106	SYSTEM PROBLEM NO CUFF PRESSURE DETECTED.	1. Check that Cuff is attached. 2. Check Cuff for leaks. 3. Resume operation.
107	PROCESS PROBLEM NO AIR / PLASMA INTERFACE DETECTED.	1. Check pump tubing installation. 2. Check tubing for occlusions. 3. Check venipuncture. 4. Press RETURN to continue.
108	DATA CANNOT BE SENT CURRENTLY TO HAEMONET PLEASE CHECK THE CONNECTION: RS-232 CABLE CONNECTED TO THE MCS+. HAEMONET DATA ACQUISITION SERVER STARTED ON PC.	<i>No Help Screen is provided.</i>
109	PRIMING PROBLEM ANTICOAGULANT PUMP IS NOT LOADED.	Please verify the following: 1. AC pump is correctly loaded. 2. DPM tubing is not clamped.
110	INSTALLATION PROBLEM UNEXPECTED FLUID DETECTED IN DISPOSABLE SET	Please verify the following: 1. Before Autoload: AC/ Saline bags not spiked. After Autoload: Pumps are loaded correctly. 2. Discard the set and power off if a significant amount of fluid is in the disposable set.
115	SYSTEM PROBLEM NO PRESSURE CHANGE DETECTED AT THE SPM.	1. Check the SPM tubing for occlusions. 2. Check the SPM tubing for clamps. 3. Ensure the SPM filter is dry. 4. Resume operation. NOTE: The product shelf life becomes 24 hours if the SPM filter is wet.

#	NOTICE message	HELP message
116	SYSTEM PROBLEM NO PRESSURE CHANGE DETECTED AT DPM.	The DPM cannot detect venous pressure change. 1. Check the DPM tubing for occlusions. 2. Check the DPM tubing for clamps. 3. Ensure the DPM filter is dry. 4. Resume operation.
162	STOP	
175	AIR DETECTED BY THE ANTICOAGULANT AIR DETECTOR (ACAD).	The Air Detector detects air. 1. Check tubing installation in detector. 2. Remove air bubbles. 3. Resume operation.
223	CENTRIFUGE PROBLEM THE CENTRIFUGE LID IS NOT LOCKED.	1. Close centrifuge lid. 2. Tighten cover knob completely. 3. Resume operation.
224	VALVE PROBLEM THE PURPLE VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
225	VALVE PROBLEM THE ORANGE VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
226	VALVE PROBLEM THE BLUE VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
227	VALVE PROBLEM THE YELLOW VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
228	VALVE PROBLEM THE WHITE VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
229	VALVE PROBLEM THE GREEN VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
230	VALVE PROBLEM THE CLEAR VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.

#	NOTICE message	HELP message
252	SYSTEM PROBLEM THE COMPUTER HAS ENCOUNTERED AN ERROR.	A CPU error caused the MCS+ to reset. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
254	VALVE PROBLEM THE RED VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
255	CONNECTOR THE DPM FILTER HAS BEEN DISCONNECTED.	Please install the filter. 1. Attach filter. 2. Ensure that the filter is not wet. 3. Resume operation.
256	CONNECTOR THE SPM FILTER HAS BEEN DISCONNECTED.	Please install the filter. 1. Attach filter. 2. Ensure that the filter is not wet. 3. Resume operation.
257	PRESSURE PROBLEM HIGH PRESSURE ERROR AT THE DPM.	The DPM high pressure limit was exceeded. 1. Check tubing for occlusions. 2. Check tubing for closed clamps. 3. Check tubing at valves. 4. Check tubing at pumps. 5. Resume operation.
258	PRESSURE PROBLEM LOW PRESSURE ERROR AT THE DPM.	The DPM low pressure limit was exceeded. 1. Check tubing for occlusions. 2. Check tubing for closed clamps. 3. Check tubing at valves. 4. Check tubing at pumps. 5. Resume operation.
259	PRESSURE PROBLEM SPM HIGH PRESSURE ERROR.	The SPM high pressure limit was exceeded. 1. Check tubing for occlusions. 2. Check tubing for closed clamps. 3. Check tubing at valves. 4. Relieve pressure by opening the effluent valve. NOTE: The product shelf life becomes 24 hours if the SPM filter is wet.
260	PRESSURE PROBLEM THE SPM PRESSURE WAS TOO LOW.	The SPM low pressure limit was exceeded. 1. Check tubing for occlusions. 2. Check tubing for closed clamps. 3. Check tubing at valves. 4. Check that product bags have not collapsed.
261	PROCESS PROBLEM DRAW NOT ALLOWED, PRESS RETURN TO CONTINUE.	Cycle was interrupted with a full bowl. The MCS+ cannot achieve re-separation. Press RETURN to continue. NOTE: Check venipuncture before pressing RETURN.

#	NOTICE message	HELP message
262	DOOR OPEN THE PROTOCOL CARD DOOR IS OPEN.	1. Close the Protocol Card door. 2. Resume operation. NOTE: DO NOT remove/insert cards during operation of the MCS+.
263	CENTRIFUGE PROBLEM THE CENTRIFUGE LID IS NOW LOCKED.	
264	DOOR CLOSED THE PROTOCOL CARD DOOR IS CLOSED.	1. Ensure the door stays closed. 2. Continue operation. NOTE: DO NOT remove/insert cards during operation of the MCS+.
265	SYSTEM PROBLEM PUMP MANIFOLD NOT INSTALLED.	The white pump manifold is not snapped in place. 1. Snap manifold in place. 2. Check pump tubing installation. 3. Resume operation.
266	VALVE PROBLEM THE BLUE VALVE WAS DETECTED OUT OF POSITION.	1. Toggle the valve manually. 2. Resume operation. NOTE: If problem persists, call Haemonetics Field Service.
280	POLE WEIGHER PROBLEM THE WEIGHT OF THE WEIGHER BAGS CHANGED.	An unexpected weight change occurred. Do NOT touch the weigher during operation. 1. Check that the bags hang free. 2. Check tuning at valves. 3. Resume operation.
281	POLE WEIGHER PROBLEM THE WEIGHT OF THE WEIGHER BAGS CHANGED. PRESS YES TO TARE THE WEIGHER.	An unexpected weight change occurred. Do NOT touch the weigher during operation. 1. Check that the bags hang free. 2. Check tuning at valves. 3. Press YES to tare the weigher and resume operation.
282	POLE WEIGHER PROBLEM PLEASE REMOVE ALL WEIGHT FROM POLE WEIGHER.	1. Remove all weight from the pole weigher. 2. Resume operation. 3. Replace bag when requested.
283	POLE WEIGHER PROBLEM EXPECTED FLOW TO POLE WEIGHER NOT DETECTED	No change in bag weight was detected. 1. Check that bags are properly installed. 2. Check tubing for occlusions. 3. Check tubing at valves. 4. Resume operation.

#	NOTICE message	HELP message
284	POLE WEIGHER PROBLEM EXPECTED FLOW FROM POLE WEIGHER NOT DETECTED	No change in bag weight was detected. 1. Check that bags are properly installed. 2. Check tubing for occlusions. 3. Check tubing at valves. 4. Resume operation.
285	POLE WEIGHER PROBLEM UNEXPECTED WEIGHT VARIATION	An unexpected weight change occurred. Do NOT touch the weigher during operation. 1. Check that the bag hangs free. 2. Check tubing at valves. 3. Resume operation.
286	SELF TEST PROBLEM THE [TYPE] WEIGHER IS OUT OF TOLERANCE.	The weigher signal is out of the allowed range. 1. Remove all weight from weigher. 2. Extend weigher arm. 3. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
301	STOP	
331	AIR DETECTED BY THE BLOOD LINE AIR DETECTOR (BLAD).	The Air Detector detects air. 1. Check tubing installation in detector. 2. Remove air bubbles. 3. Resume operation.
338	NO FLOW PRIME FLOW CONTROL PROBLEM.	1. Check that tubing between the AC bag and bowl inlet is correctly installed and has no occlusions. 2. Check that AC bag seal is open. 3. Turn Blood pump to bring DPM pressure to the third of its range.
339	NO FLOW RETURN FLOW CONTROL PROBLEM.	Caused by high pressure or vacuum at DPM. 1. Check donor tubing for occlusions. 2. Check donor tubing for closed clamps. 3. Check venipuncture. 4. Relieve pressure. 5. Resume operation.
340	AIR REMOVED AIR REMOVED FROM THE AIR DETECTOR.	The Air Detector detects fluid. 1. Check tubing for residual air. 2. Resume operation.
341	NO FLOW NO FLOW FROM THE DONOR IS DETECTED.	Caused by high pressure or vacuum at DPM. 1. Check donor tubing for occlusions. 2. Check donor tubing for closed clamps. 3. Check venipuncture. 4. Relieve pressure. 5. Resume operation.
345	PRESSURE PROBLEM THE DPM PRESSURE IS OUT OF RANGE.	1. Close needle clamp. 2. Open RED valve. 3. Turn Blood pump CLOCKWISE. 4. Resume operation.

#	NOTICE message	HELP message
366	PROCESS PROBLEM TERMINATE THE PROCEDURE NOW? PRESS YES TO END PROCEDURE. PRESS NO TO RETURN TO THE READY MODE.	No Help Screen is provided.
367	TERMINATING PROCEDURE PRESS YES TO RETURN BOWL CONTENT AND END. PRESS NO TO END THE PROCEDURE IMMEDIATELY.	No Help Screen is provided.
368	PROCESS PROBLEM INITIATE RETURN NOW? PRESS RETURN TO INITIATE RETURN CYCLE. PRESS SAVE TO RETURN BOWL CONTENT AND END. PRESS DRAW TO RESUME DRAW CYCLE.	No Help Screen is provided.
370	SYSTEM PROBLEM THE EFFLUENT TUBING IS NOT INSTALLED.	The effluent tubing may not be installed. 1. Check that each effluent line is loaded in a valve of the same color. 2. Check slide clamp at SPM. 3. Resume operation.
371	SYSTEM PROBLEM THE PLASMA BAG IS CLAMPED.	The tubing to the product bag is clamped. 1. Check tubing at Yellow valve. 2. Check tubing for closed clamps. 3. Check tubing for kinks. 4. Resume operation.
372	SYSTEM PROBLEM THE PROTOCOL CARD DOES NOT MATCH DISPOSABLE.	The installed Protocol Card and disposable do not match. 1. Check that the correct disposable is installed. 2. Check that the correct Protocol Card is inserted. 3. Clean the disposable ID window.
373	PRESSURE PROBLEM NO DONOR PRESSURE CHANGE DETECTED AT THE DPM.	The DPM cannot detect venous pressure change. 1. Check the DPM tubing for occlusions. 2. Check the DPM tubing for clamps. 3. Ensure the DPM filter is dry. 4. Resume operation.
374	PRESSURE PROBLEM NO SYSTEM PRESSURE CHANGE DETECTED AT THE SPM.	1. Check the SPM tubing for occlusions. 2. Check the SPM tubing for clamps. 3. Ensure the SPM filter is dry. 4. Resume operation. NOTE: The product shelf life becomes 24 hours if the SPM filter is wet.

#	NOTICE message	HELP message
375	PROCEDURE RECOVERY THE POWER HAS BEEN OFF FOR A SHORT TIME.	A recovery sequence can be performed to CONTINUE the last procedure if desired, or a NEW procedure may be initiated.
376	PROCESS PROBLEM RETURN CYCLE WAS NOT COMPLETED.	The bowl is not empty. 1. Press RETURN to complete return cycle. 2. Press STOP to stop draw cycle. 3. Press DRAW to start draw cycle.
377	FRONT WEIGHER PROBLEM EXPECTED FLOW TO FRONT WEIGHER NOT DETECTED.	No change in bag weight was detected. 1. Check that bags are properly installed. 2. Check tubing for occlusions. 3. Check tubing at valves. 4. Resume operation.
378	FRONT WEIGHER PROBLEM EXPECTED FLOW FROM FRONT WEIGHER NOT DETECTED.	No change in bag weight was detected. 1. Check that bags are properly installed. 2. Check tubing for occlusions. 3. Check tubing at valves. 4. Resume operation.
380	LINE SENSOR PROBLEM EXPECTED FLUID FLOW FROM BOWL NOT DETECTED.	Line sensor detects air. Please verify the following: Line sensor tubing is correctly installed.
382	WEIGHER PROBLEM THE WEIGHT OF THE WEIGHER BAGS CHANGED.	An unexpected weight change occurred. Do NOT touch the Weigher during operation. 1. Check that the bag hangs free. 2. Check tubing at valves. 3. Press YES to tare the weigher and resume operation.
383	BOWL TYPE MISMATCH	A different bowl size was detected. Confirm whether the MCS+ detected the correct bowl: 1. Press YES to accept the new bowl size and continue. 2. Press NO to keep the existing bowl size and continue.
384	TRANSFER PROBLEM NO WEIGHER INCREASE DURING TRANSFER.	The transfer was not completed correctly. Please check that: 1. The RBC bags are on the weigher. 2. The RBC bags are open. 3. The additive bag is open. 4. The tubing is correctly loaded.
385	PROCESS PROBLEM THE SOLUTION LINES MAY BE IN WRONG VALVES.	The solution lines may be in wrong valves. 1. Check line and valve colors for mismatch. 2. If mismatch, clamp lines and place them in correct valves. 3. Prime compensation (purple) line manually to "Y" connector. 4. Press PRIME to continue.
394	NOTICE THE AC BAG MAY BE ALMOST EMPTY.	

#	NOTICE message	HELP message
395	NOTICE THE SALINE BAG MAY BE ALMOST EMPTY.	
396	PROCESS PROBLEM EXCEED TOTAL PRODUCT VOLUME LIMIT? PRESS RETURN TO REDUCE THE PRODUCT TARGET VOLUMES AND REMAIN WITHIN THE LIMIT. PRESS YES TO COLLECT THE TARGETED VOLUMES AND EXCEED THE PRODUCT VOLUME LIMIT.	The FDA collection volume limit will be reached. 1. Press RETURN to stay within volume limit. Volume of plasma will be reduced. If no plasma product, collection will be ended. 2. Press YES to exceed the max volume limit. Collection will continue and exceed limit.
399	PROCESS PROBLEM RETURN NOT ALLOWED, PRESS DRAW TO CONTINUE.	No Help Screen is provided.
400	SELF TEST PROBLEM THE COMPUTER CPU TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
402	SELF TEST PROBLEM THE COMPUTER BIOS CRC TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
404	SELF TEST PROBLEM THE PROTOCOL CARD CRC TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
405	SELF TEST PROBLEM THE CALIBRATION CRC TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
406	SELF TEST PROBLEM THE CONFIGURATION CRC TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
407	NOTICE THE PROTOCOL CARD DOES NOT MATCH DISPOSABLE.	The installed Protocol Card and disposable do not match. 1. Check that the correct disposable is installed. 2. Check that the correct Protocol Card is inserted. 3. Clean the disposable ID window.
408	SELF TEST PROBLEM ANTICOAGULANT AIR DETECTOR (ACAD) TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
409	SELF TEST PROBLEM BLOOD LINE AIR DETECTOR (BLAD) TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.

#	NOTICE message	HELP message
410	SELF TEST PROBLEM DONOR LINE AIR DETECTOR 1 (DLAD1) TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
411	SELF TEST PROBLEM DONOR LINE AIR DETECTOR 2 (DLAD2) TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
413	VALVE PROBLEM PURPLE VALVE SELF TEST FAILED.	1. Toggle the valve manually. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
414	VALVE PROBLEM ORANGE VALVE SELF TEST FAILED.	1. Toggle the valve manually. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
415	VALVE PROBLEM RED VALVE SELF TEST FAILED.	1. Toggle the valve manually. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
416	VALVE PROBLEM BLUE VALVE SELF TEST FAILED.	1. Toggle the valve manually. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
417	VALVE PROBLEM YELLOW VALVE SELF TEST FAILED.	1. Toggle the valve manually. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
418	VALVE PROBLEM WHITE VALVE SELF TEST FAILED.	1. Toggle the valve manually. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
419	VALVE PROBLEM GREEN VALVE SELF TEST FAILED.	1. Toggle the valve manually. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
420	SELF TEST PROBLEM A/D CONVERTER 20% TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.

#	NOTICE message	HELP message
421	SELF TEST PROBLEM A/D CONVERTER 50% TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
422	SELF TEST PROBLEM A/D CONVERTER 80% TEST FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
423	SELF TEST PROBLEM PLEASE CLAMP THE DPM LINE AND REMOVE FILTER.	During self test, the DPM and SPM filters can not be installed. 1. Clamp tubing to the appropriate filter. 2. Remove filter from connector. 3. Resume operation.
424	SELF TEST PROBLEM PLEASE CLAMP THE SPM LINE AND REMOVE FILTER.	During self test, the DPM and SPM filters can not be installed. 1. Clamp tubing to the appropriate filter. 2. Remove filter from connector. 3. Resume operation.
425	SELF TEST PROBLEM PLEASE INSTALL THE DPM FILTER.	Please install the filter. 1. Attach filter. 2. Ensure that the filter is not wet. 3. Resume operation.
426	SELF TEST PROBLEM PLEASE INSTALL THE SPM FILTER.	Please install the filter. 1. Attach filter. 2. Ensure that the filter is not wet. 3. Resume operation.
427	SELF TEST PROBLEM THE DPM CONNECTOR TEST FAILED.	The MCS+ could not detect the pressure monitor filter was removed. 1. Remove pressure filter from connector. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
428	SELF TEST PROBLEM THE SPM CONNECTOR TEST FAILED.	The MCS+ could not detect the pressure monitor filter was removed. 1. Remove pressure filter from connector. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
429	SELF TEST PROBLEM THE DPM IS OUT OF TOLERANCE.	The Pressure Monitor signal is out of the allowed range. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.

#	NOTICE message	HELP message
430	SELF TEST PROBLEM THE SPM IS OUT OF TOLERANCE.	The Pressure Monitor signal is out of the allowed range. 1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
431	SELF TEST PROBLEM THE WEIGHER IS OUT OF TOLERANCE.	The Weigher signal is out of the allowed range. 1. Remove all weight from Weigher. 2. Extend Weigher arm. 3. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
432	SELF TEST PROBLEM THE CUFF IS OUT OF TOLERANCE.	The Cuff pressure is out of the allowed range. 1. Check Cuff tubing for occlusions. 2. Check the Cuff connection. 3. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
433 434	SELF TEST PROBLEM THE SAFETY SYSTEM FAILED THE SELF TEST.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
435	SELF TEST PROBLEM THE LINE SENSOR OUTPUT IS TOO LOW.	1. Remove tubing from Line Sensor. 2. Check Line Sensor for contamination. 3. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
436	SELF TEST PROBLEM THE LINE SENSOR OUTPUT IS TOO HIGH.	1. Remove tubing from Line Sensor. 2. Check Line Sensor for contamination. 3. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
437	SELF TEST PROBLEM THE SAFETY SYSTEM FAILED THE SELF TEST.	
438	SELF TEST PROBLEM THE WATCHDOG FAILED THE SELF TEST.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
439	SELF TEST PROBLEM LOW PRESSURE ERROR AT THE DPM.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.

#	NOTICE message	HELP message
441	SELF TEST PROBLEM HIGH PRESSURE ERROR AT THE DPM.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
443	SELF TEST PROBLEM STATE TEST ON CHANNEL 1 FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
444	SELF TEST PROBLEM STATE TEST ON CHANNEL 2 FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
447	SELF TEST PROBLEM TIME BASE TEST ON CHANNEL 1 FAILED.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
448	SELF TEST PROBLEM TIME BASE TEST ON CHANNEL 2 FAILED.	1. Turn power OFF, then ON and try again NOTE: If problem persists, call Haemonetics Field Service.
449 452 453	SELF TEST PROBLEM THE SAFETY SYSTEM FAILED THE SELF TEST.	1. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
454	VALVE PROBLEM PURPLE VALVE MAY NOT BE LOADED.	The disposable test indicates that the valve may not be loaded. 1. Check tubing installation at valve. 2. Resume operation.
455	VALVE PROBLEM ORANGE VALVE MAY NOT BE LOADED.	The disposable test indicates that the valve may not be loaded. 1. Check tubing installation at valve. 2. Resume operation.
456	VALVE PROBLEM RED VALVE MAY NOT BE LOADED.	The disposable test indicates that the valve may not be loaded. 1. Check tubing installation at valve. 2. Resume operation.
457	VALVE PROBLEM BLUE VALVE MAY NOT BE LOADED.	The disposable test indicates that the valve may not be loaded. 1. Check tubing installation at valve. 2. Resume operation.
458	VALVE PROBLEM WHITE VALVE MAY NOT BE LOADED.	The disposable test indicates that the valve may not be loaded. 1. Check tubing installation at valve. 2. Resume operation.
459	VALVE PROBLEM YELLOW VALVE MAY NOT BE LOADED.	The disposable test indicates that the valve may not be loaded. 1. Check tubing installation at valve. 2. Resume operation.

#	NOTICE message	HELP message
460	VALVE PROBLEM GREEN VALVE MAY NOT BE LOADED.	The disposable test indicates that the valve may not be loaded. 1. Check tubing installation at valve. 2. Resume operation.
461	SELF TEST PROBLEM THE CENTRIFUGE VACUUM SYSTEM IS OBSTRUCTED.	Centrifuge vacuum could not be released. 1. Remove bowl. 2. Check that white filter in Centrifuge Chuck base is dry. 3. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
462	SELF TEST PROBLEM DISPOSABLE ID READER FAILED SELF TEST.	1. Remove white pump manifold. 2. Clean disposable ID reader window. 3. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
466	WEIGHER PROBLEM PLEASE REMOVE ALL WEIGHT FROM WEIGHER.	1. Remove any weight on Weigher. 2. Resume operation. 3. Replace bag AFTER self test.
467	SELF TEST PROBLEM PLEASE DISCONNECT THE CUFF FROM REAR PANEL.	Excessive Cuff pressure was measured. 1. Disconnect cuff from rear panel. 2. Deflate cuff completely. 3. Connect cuff to rear panel. 4. Resume operation.
468	SELF TEST PROBLEM IS LINE SENSOR CLEAN AND TUBING REMOVED?	1. Remove tubing from Line Sensor. 2. Check Line Sensor for contamination. 3. Press YES to re-test the Line Sensor. NOTE: If problem persists, call Haemonetics Field Service.
469	SELF TEST PLEASE UNLOCK THE CENTRIFUGE LID.	The centrifuge knob must be tested for safe operation. 1. OPEN and CLOSE the centrifuge knob as requested.
470	SELF TEST PLEASE LOCK THE CENTRIFUGE LID.	The centrifuge knob must be tested for safe operation. 1. OPEN and CLOSE the centrifuge knob as requested.
471	NOTICE THE PROTOCOL PARAMETERS HAVE BEEN RESET.	The saved protocol parameters have been reset. Before starting the next procedure: 1. Re-program the Haemo Calculator. 2. Re-program Modify.
472	VALVE PROBLEM CLEAR VALVE SELF TEST FAILED.	1. Toggle the valve manually. 2. Turn power OFF, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.

#	NOTICE message	HELP message
475	SELF TEST PROBLEM THE CALIBRATION REVISION IS INCOMPATIBLE.	1. Turn power off, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
476	SELF TEST PROBLEM THE CONFIGURATION REVISION IS INCOMPATIBLE.	1. Turn power off, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
477	NOTICE THE PROTOCOL PARAMETERS HAVE BEEN RESET.	The saved protocol parameters have been reset. Before starting the next procedure: 1. Re-program the Haemo Calculator. 2. Re-program Modify.
478	SELF TEST PROBLEM THE CENTRIFUGE RELAY TEST FAILED.	.
482	SELF TEST PROBLEM DUAL OPTICS ERROR.	1. Turn power off, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.

MCS+ SAFETY SYSTEM: HARDWARE DETECTION

The MCS+ device is designed with a special electronic safety circuit board which will oversee all device functions, including the MCS+ software programming. The *WATCHDOG* system will completely shut down all device function if a problem is detected at this level, guaranteeing security for the donor/patient.

The following information could appear to the operator if this type of situation should occur.

Note: These messages will appear on the MCS+ display screen in English.

	Notification	Instructions
–	SOFTWARE WATCHDOG INVALID ID	The computer has sensed an invalid Software Watchdog ID code. Turn power off, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
–	SOFTWARE WATCHDOG INVALID ID SEQUENCE	The computer has sensed an incorrect Software Watchdog ID sequence. Turn power off, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
–	SOFTWARE WATCHDOG ID TIME OUT	The computer has sensed a Software Watchdog ID time out. Turn power off, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.
–	SOFTWARE WATCHDOG NOT ENABLED	The computer has sensed that the Software Watchdog was not enabled. Turn power off, then ON and try again. NOTE: If problem persists, call Haemonetics Field Service.

S/1-22

Listing NOTICES and Relevant HELP Messages

P/N 85270-30, Manual revision: B

Section 2

Using MCS+ Data Acquisition Features

DESCRIBING THE MCS+ COMMUNICATION OPTIONS	S/2-2
Printer	S/2-2
Wireless	S/2-2
Network	S/2-2
Disk	S/2-3
None	S/2-3
ENTERING PROCEDURE DATA WITH A BAR-CODE READER	S/2-4
Disposable set list and lot numbers	S/2-4
Anticoagulant and solution lot numbers	S/2-5
Donation, donor and operator code numbers	S/2-6
ACCESSING STORED PROCEDURE DATA	S/2-7
Data management menu	S/2-7
DOWNLOADING, CLEARING AND PRINTING PROCEDURE DATA ..	S/2-9
Downloading data	S/2-9
Clearing data	S/2-9
Printing data	S/2-10
PRESENTING THE WIRELESS FEATURE AND COMPONENTS	S/2-11
What does wireless mean	S/2-11
Describing the wireless device components	S/2-11
How does the Haemonetics device interface with the eLynx system? ..	S/2-12
PRESENTING THE ELYNX SYSTEM	S/2-13

DESCRIBING THE MCS+ COMMUNICATION OPTIONS

The MCS+ device can be equipped with a wireless antenna, an external communication box or internal data card and a bar-code reader. These elements allow the operator to introduce data directly into the MCS+ device and subsequently transfer procedure data to an external non-medical device such as a printer, or link the device to a central monitoring computer by using HaemoNet or eLynx. The MCS+ protocols are designed to request the data inputs at specific moments, then transfer the stored information at the appropriate time.

Warning: Only Haemonetics-approved data acquisition devices and cables should be connected to the MCS+ device. To comply with the requirements of IEC 60601-1-1 Standard for Medical Electrical Equipment, any non-Haemonetics approved device *not* powered from the MCS+ device must be placed at least 1.5 meters from the donor and any cable must not exceed a length of 15 meters.

For cleaning procedures of any MCS+ data acquisition accessories not described in the relevant accompanying manual, Haemonetics recommends using a disinfectant solution appropriate for cleaning the MCS+ device. Further information is provided in the chapter of the MCS+ device operation manual, "Maintaining the MCS+ Equipment".

The MCS+ device will be configured with one of the following options at the time of installation. The apheresis center should contact the local Haemonetics representative if a configuration change is desired.

Note: The operator is responsible for making sure that the protective earthing and the leakage current of the final configuration in normal and single fault conditions comply with IEC 60601-1-1 Standard, Medical electrical equipment.

All ancillary equipment and cables should be set-up in safe manner to prevent any mechanical damage.

Printer

This configuration will allow the MCS+ device to send a record of the complete procedure statistics to a printer at the end of each collection procedure.

Wireless

The wireless feature allows the transmission of procedure data to a wireless access point connected to a handheld device, internal PC or network, the eLynx communication system, or to an external computer network over the Internet. Further information is available in "Presenting the wireless feature and components" on page S/2-11.

Network

HaemoNet is the name of the Haemonetics communication network. This configuration will connect the MCS+ device to a central monitoring computer. The

MCS+ device will receive prompts from the monitoring computer and will respond by transmitting requested data. This transmitted data can be stored in a database and/or viewed directly on the central monitor computer screen.

The MCS+ device will register the prompts sent by the monitoring computer and respond by transferring the requested data. The data transmitted will be stored in a database or displayed on the computer screen for central monitoring purposes.

Disk

This configuration is used by Haemonetics technicians to record certain procedure data. It is designed to log information for technical purposes only.

None

The MCS+ device will not transfer any data to or from the communication box or data card.

ENTERING PROCEDURE DATA WITH A BAR-CODE READER

The bar-code reader is used to scan information contained on disposable packaging and solution bags, as well as the donation number, donor number and operator number for each collection procedure. The MCS+ bar-code reader can scan any type of bar-code. This information can be grouped as:

- Disposable set lot and list numbers.
- Anticoagulant and other solution codes.
- Donation number, donor number and operator codes.

The MCS+ device can be configured to "request" either all or only certain of these data groups. There is a specific sequence in which the information will appear on the MCS+ display screen. If the device is not configured to request certain information, the operator will not receive the related screen message, as depicted in the following examples.

Disposable set list and lot numbers

The first series of screen messages will appear once the disposable set elements have been installed, but prior to the MCS+ pump autoload sequence. The following type of screen will be displayed for "unbundled" disposable sets:

DISPOSABLE	
LIST #	LOT #
HARNESS IN	
BOWL	
HARNESS OUT	
NEEDLE	
Press DRAW to confirm	
Press MODIFY to switch between the fields. Press STOP to clear and start again.	

Figure 2-1, Example of the unbundled disposable set message

The following type of screen will be displayed for "closed" disposable sets:

DISPOSABLE	
LIST #	LOT #
DISPOSABLE	
Press DRAW to confirm	
Press MODIFY to move between the fields. Press STOP to clear and start again.	

Figure 2-2, Example of the closed disposable set message

When either of these screens appears, the operator should:

- Scan the appropriate information using the bar-code reader. Once the data has been registered, a “beep” will be heard.
- Press the Draw key to access the PRIME mode screen. Two “beeps” will be heard. The operator can proceed to prime the disposable set.

Anticoagulant and solution lot numbers

Once the operator has primed the disposable set, the following screen will be displayed for MCS+ protocol options not using additional solutions.

ANTICOAGULANT	
AC #	LOT # [REDACTED]
AC #2	
Press DRAW to confirm	
Press MODIFY to move between the fields. Press STOP to clear and start again.	

Figure 2-3, Example of the anticoagulant message

The following screen will be displayed for MCS+ protocols using additional solutions (example: saline, SAG-M). The operator can enter one or two solution in addition to the AC solution code, for each collection procedure.

AC & SOLUTIONS			
AC #	LOT # [REDACTED]	Solution #1	LOT #
AC #2		Solution #2	
Press DRAW to confirm			
Press MODIFY to move between the fields. Press STOP to clear and start again.			

Figure 2-4, Example of the AC & solutions message

When either of these screens appears, the operator should:

- Scan the appropriate information using the bar-code reader. Once the data has been registered, a “beep” will be heard.
- Press the Draw key to access the next data input screen. Two “beeps” will be heard.

Donation, donor and operator code numbers

Prior to entering the READY mode, at which time the collection procedure can be initiated, the following screen will be displayed. At this time the operator can enter the donation number, donor number and operator code.

Note: The donation number is the link to all other data entered for each procedure. The donation number identifies each individual procedure and serves as a reference for the entire data acquisition system. If the donation number is not entered at this point, the MCS+ device will prompt the operator to enter it at the end of the collection procedure.

DONATION & OPERATOR	
DONATION #	
DONOR #	
OPERATOR #	
Press DRAW to confirm	
Press MODIFY to move between the fields. Press STOP to clear and start again.	

Figure 2-5, Example of the donation & operator message

When this screen appears, the operator should:

- Scan the appropriate information using the bar-code reader. Once the data has been registered, a "beep" will be heard.
- Press the Draw key to access the READY screen. Two "beeps" will be heard. The operator can initiate blood component collection.

Note: The STOP key can be pressed to clear the any of the screens when entering data with the bar-code reader, if the Draw key has not been pressed to confirm registration of the data.

ACCESSING STORED PROCEDURE DATA

The MCS+ device will retain all data from the last ten collection procedures performed. If the operator wishes to produce a printed record, the device must be connected to a printer. The operator will need to access the DATA MANAGEMENT menu and proceed as follows:

Data management menu

- Power on the MCS+ device.
- Wait for the following screen message to be displayed:


Figure 2-6, Example of the power on message

- Press the STOP key immediately to enter the CONFIGURATION menu which will be displayed as follows:


Figure 2-7, Example of the configuration message

- Press the Modify key and scroll through the options until reaching DATA MANAGEMENT.
- Press the Draw key to access the following screen:

DATA MANAGEMENT
DOWNLOAD DATA
CLEAR DATA
Press STOP to return to the previous screen. Press MODIFY to change selection, DRAW to execute.

Figure 2-8, Example of the data management message

- Press the Modify key until DOWNLOAD DATA is highlighted.
- Press the Draw key to receive the list of stored procedure data.

Note: The STOP key can be pressed to return to the previous screen when selecting an option from the DATA MANAGEMENT menu.

DOWNLOADING, CLEARING AND PRINTING PROCEDURE DATA

Downloading data

The upper portion of this screen will indicate to the operator the configuration of the MCS+ device as either *printer* or *network*.

Note: This selection can be modified only by an authorized Haemonetics technician.

DOWNLOAD DATA			
Procedure:	1	Device:	Printer
1	431821446	13/03/1996	10:44
2	831451185	13/03/1996	17:13
3	1A31BCD438	15/03/1996	23:05
4		16/03/1996	9:19
5	891000191	17/03/1996	7:32
6	23435 885	19/03/1996	21:01
7	0188258137	20/03/1996	8:14
8		21/03/1996	10:49
9			
10			

Press +/- to select the procedures.
Press DRAW to download DATA to the selected device.

Figure 2-9, Example of the download data message

Once the operator has accessed the list of stored procedure data, the selection can be made to download this information as follows:

- Press the "+" or "-" key to highlight a procedure on the display screen.
- Press the Draw key to download the selected procedure to a printer or to HaemoNet, the Haemonetics network.

Clearing data

After downloading stored procedure data, the operator can clear the data from the MCS+ device as follows:

- Press the STOP key to return to the DATA MANAGEMENT menu.
- Press the Modify key until CLEAR DATA is highlighted.

The following screen will appear:

CLEAR DATA
READY TO CLEAR DATA
Press DRAW to clear all recorded data. Press STOP to abort.

Figure 2-10, Example of the clear data message

- Press the Draw key to remove the procedure data stored in the device.

Printing data

When the MCS+ device is configured to provide a printed record of the collection procedure, the information will be presented, as depicted in the following example of an MCS+ procedure print-out.

PROCEDURE RECORD			
DONATION #	77021220124	DONOR #	
Date:	31-08-1999		
MACHINE TYPE:	MCS+	Protocol:	LDP
SERIAL NB.:	94A125	VERSION:	MCS+2-LDP-X.X-LL
DISPOSABLE			
	LIST #	LOT #	*
ANTICOAGULANT			
	LOT #		
	96tKLJH3455		
SOLUTIONS			
	LOT #		
PROCEDURE DATA			
Number of Cycles	8		
Elapsed Time	81	min	
Volume Processed	4100	ml	
AC Volume Used	483	ml	
Plasma Weight	10	g	
Platelet Weight	360	g	
NaCl Volume Used	0	ml	
HAEMOCALCULATOR DATA			
Sex	M	Target Plasma Wgt	0 g
Height	170 cm	End Criteria	2
Weight	65 kg	Target Yield	5.0 10e11
Blood Volume	4800 ml	Process Volume	4000 ml
HCT	40 %	Target Cycles	8
Plt Pre-Count	250 10e3	Time (estimated)	80 min.
COMMENTS			
PROCEDURE RESULTS			
PLATELETS	10e11	WBC COUNT	10E6
OPERATOR #	0426		
PHYSICIAN _____ (if required by local rules)			
Date / Signature			

PRESENTING THE WIRELESS FEATURE AND COMPONENTS

What does wireless mean

This section describes the addition of wireless capability to the Haemonetics devices and advantages of this feature.

Wireless technology provides the following benefits:

- Enables the Haemonetics devices to transmit data electronically to a hand-held device, PC or a central location over the Internet.
- Allows electronic data collection and management with an appropriate system, such as eLynx.
- Eliminates the need for cabling which allows easier movement of devices in the center.

Note: The wireless function is independent of the operator and of the selected protocol. No additional operator training is required to use the wireless feature.

Note: A short overview of the eLynx system is provided in the following section Presenting the eLynx system.

Describing the wireless device components

The wireless hardware components are:

- An antenna located on the Haemonetics device rear panel, designed for connection to wireless devices.
- A wireless module communication component, designed to provide connectivity in several type of applications (not visible, the module is inside the device).

Figure 11, Example of the wireless antenna located on the device rear panel

Note: The wireless components of the Haemonetics device comply with all relevant IEC and ISO standards for medical electrical equipment.

The operator should refer to the specific Haemonetics operation manuals for all applicable electrical and procedural information and warnings.

How does the Haemonetics device interface with the eLynx system?

The antenna, located on the rear panel and protected by a hard plastic cover, allows the transmission of procedure data to a wireless access point connected to the communication system, eLynx.

The eLynx system collects and consolidates this information, along with donor information entered on the eLynx handheld device.

The eLynx system can then print or transmit the data to a host database management system and/or a network.

Note: The network installation and maintenance is the customer's responsibility. If the network uses wireless technology, the customer should consider and mitigate hazards associated with electromagnetic interference (EMI) and compatibility (EMC).

PRESENTING THE ELYNX SYSTEM

eLynx is a computerized system to electronically record and consolidate procedure information that would normally be collected by hand.

eLynx is a wireless system that enables the operator:

- ! To manage donors, devices and units in an automated environment.
- ! To enter procedure data into the system through a handheld PDA or a standard PC.
- ! To get procedure information in specified XML formats from external systems or devices.
- ! To review online and print the collected procedure information.

eLynx is a temporary data storage system intended to facilitate the collection and consolidation of procedure activities, and its subsequent download or transfer to the system of record. Possible permanent storage systems include donor management systems, simple electronic file archives, and optical or paper based archives.

The following figure depicts a typical network architecture in a donor blood center, using the eLynx system.

1. Wireless access point
2. Wireless enabled devices
3. eLynx handheld device (pocket PC based PDA)
4. Hub/switch
5. Terminal server
6. Serial devices
7. eLynx console (standard PC or laptop)
8. eLynx server
9. External computer system(s)

Type of connections

- A. Wireless connection
- B. Wired ethernet connection
- C. Wired serial connection
- D. Network connection (eg. LAN, WAN, VPN)


Figure 12, The eLynx™ system interaction and connection

HAEMONETICS CORPORATION
400 WOOD ROAD
BRAINTREE, MA 02184 USA
1-800-225-5242

Yulia Perelman
Yulia Perelman
INTERNATIONAL Regulatory Affairs
Specialist
6-OCTOBER-2014

P/N 85270-30, Manual revision: B

This Document contains 40 pages
that are bound and stamped pages.

Yulia PerIman

Yulia PERIMAN

INTERNATIONAL REGULATORY
AFFAIRS SPECIALIST

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