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INSTRUCTIONS FOR USE FOR A MEDICAL DEVICE

«Ophthalmic blades»

Productions of Microtrack Surgicals (Mikrotrek Serdzhikals)

A-38: Adarsh-2, Industrial Estate, Behind Ashish Cinema, Odhav, Ahmedabad-
382 415,

Gujarat, India.

Phone No: - 079-22-970972


Email ID: - microtrack10@yahoo.com

S.R. No. 2392-2024

ILABEN R. PATEL
NOTARY
GOVT. OF GUJARAT
(INDIA)

11 JUL 2024

Date of last approval June 2024

2 IV

1. Purpose of MI, potential consumers and scope of application.
The type and duration of contact with the human body.

Purpose.

The medical product "Ophthalmic blade" is intended for use in surgical operations on the organs of vision, as well as for punctures, incisions and delamination of eye tissues during ophthalmic surgical interventions with minimal damage to the tissues of the eyeball.

Potential consumers The product is intended for use in medical institutions, in operating rooms, by qualified medical personnel (ophthalmic surgeon).

Scope - ophthalmic surgery.

The type and duration of contact of the medical device with the human body.

Short-term contact with the mucous membrane of the eye, short-term contact with intact skin.

2. Product description, design, technical and functional characteristics.

2.1. Execution options.

Ophthalmic blade, consisting of:

1. Ophthalmic blade - 1 pack. (1pack - 10 pcs), execution options:

1. Lancet-type ophthalmic blade, design options:


- 1.1. Ophthalmic blade Lancet type 15°
- 1.2. Ophthalmic blade Lancet type 22.5°
- 1.3. Ophthalmic blade Lancet type 30°
- 1.4. Ophthalmic blade Lancet type 45°

2. Ophthalmic sickle blade, design options:

- 2.1. Ophthalmic sickle blade straight with bevel up, size 20 G
- 2.2. Ophthalmic sickle blade straight with bevel down, size 20 G
- 2.3. Ophthalmic sickle blade curved with bevel up, size 20 G
- 2.4. Ophthalmic sickle blade curved with bevel down, size 20 G
- 2.5. Ophthalmic sickle blade straight with bevel up, size 19 G
- 2.6. Ophthalmic sickle blade straight with bevel down, size 19 G
- 2.7. Ophthalmic sickle blade curved with bevel up, size 19 G
- 2.8. Ophthalmic sickle blade curved with bevel down, size 19 G

3. Ophthalmic keratome blade for cutting, design options:

- 3.1. Ophthalmic keratome blade for cutting, blade width 1.2 mm
- 3.2. Ophthalmic keratome blade for cutting, blade width 1.3 mm
- 3.3. Ophthalmic keratome blade for cutting, blade width 1.4 mm
- 3.4. Ophthalmic keratome blade for cutting, blade width 1.5 mm
- 3.5. Ophthalmic keratome blade for cutting, blade width 1.6 mm
- 3.6. Ophthalmic keratome blade for cutting, blade width 1.7 mm
- 3.7. Ophthalmic keratome blade for cutting, blade width 1.8 mm
- 3.8. Ophthalmic keratome blade for cutting, blade width 1.9 mm
- 3.9. Ophthalmic keratome blade for cutting, blade width 2.0 mm
- 3.10. Ophthalmic keratome blade for cutting, blade width 2.1 mm
- 3.11. Ophthalmic keratome blade for cutting, blade width 2.2 mm
- 3.12. Ophthalmic keratome blade for cutting, blade width 2.3 mm
- 3.13. Ophthalmic keratome blade for cutting, blade width 2.4 mm

- 3.14. Ophthalmic keratome blade for cutting, blade width 2.5 mm
- 3.15. Ophthalmic keratome blade for cutting, blade width 2.6 mm
- 3.16. Ophthalmic keratome blade for cutting, blade width 2.65 mm
- 3.17. Ophthalmic keratome blade for cutting, blade width 2.75 mm
- 3.18. Ophthalmic keratome blade for cutting, blade width 2.8 mm
- 3.19. Ophthalmic keratome blade for cutting, blade width 3.0 mm
- 3.20. Ophthalmic keratome blade for cutting, blade width 3.2 mm
- 3.21. Ophthalmic keratome blade for cutting, blade width 3.5 mm
- 3.22. Ophthalmic keratome blade for cutting the transparent cornea, blade width 2.65 mm
- 3.23. Ophthalmic keratome blade for cutting the transparent cornea, blade width 2.8 mm
- 3.24. Ophthalmic keratome blade for cutting the transparent cornea, blade width 3.2 mm

4. Ophthalmic keratome extended blade, design options:

- 4.1. Ophthalmic keratome extended blade, blade width 3.5 mm
- 4.2. Ophthalmic keratome extended blade, blade width 4.1 mm
- 4.3. Ophthalmic keratome extended blade, blade width 5.1 mm
- 4.4. Ophthalmic keratome extended blade, blade width 5.2 mm
- 4.5. Ophthalmic keratome extended blade, blade width 5.5 mm
- 4.6. Ophthalmic keratome extended blade, blade width 6.0 mm

5. Ophthalmic microvitrectomy blade, design options:

- 5.1. Ophthalmic microvitrectomy blade straight, size 19 G
- 5.2. Ophthalmic microvitrectomy blade curved, size 19 G
- 5.3. Ophthalmic microvitrectomy blade straight, size 20 G
- 5.4. Ophthalmic microvitrectomy blade curved, size 20 G
- 5.5. Ophthalmic microvitrectomy blade straight, size 23G
- 5.6. Ophthalmic microvitrectomy blade curved, size 23G
- 5.7. Ophthalmic microvitrectomy blade straight, size 24 G
- 5.8. Ophthalmic microvitrectomy blade curved, size 24 G

6. Ophthalmic blade with controlled depth, design options:

- 6.1. Ophthalmic blade with controlled depth, 300 microns
- 6.2. Ophthalmic blade with controlled depth, 350 microns
- 6.3. Ophthalmic blade with controlled depth, 400 microns
- 6.4. Ophthalmic blade with controlled depth, 450 microns
- 6.5. Ophthalmic blade with controlled depth, 500 microns
- 6.6. Ophthalmic blade with controlled depth, 550 microns
- 6.7. Ophthalmic blade with controlled depth, 600 microns

II. Instructions for use - 1 pc.

2.2. Description of the medical device.

The medical device consists of a blade, a handle, a plastic case or a protective cap.

The blades have different sizes and shapes required for certain surgical procedures.

A plastic case or protective cap protects medical personnel from accidental cuts during disposal. The case body has a cutout for easier removal of the blade from the case.

The protective cap covers the blade:


before use, the cap slides along the handle / after use, the cap slides back onto the blade.

The color coding of the handle corresponds to the blade version.

To identify the blade with a cap, the letter "C" is added to the model number in the catalog.

For example: M922201C.


The version of the "Ophthalmic blade with controlled depth" is available only in a case.

The version "Ophthalmological sickle blade" has a certain blade sharpening. The diagram shows the direction of sharpening the knife blade.

In the versions "Ophthalmic sickle blade", "Ophthalmic microvitreoretinal blade" the working part of the blade can be straight or curved. The figure shows the types of straight and curved working part of the blade using the example of the version "Ophthalmic microvitreoretinal blade".


	
straight working part of the blade	curved working part of the blade

1. Lancet-type ophthalmic blade


The blade is used for cutting the sclera, for ultra-precise cuts and punctures.

2. Ophthalmic sickle blade


The blade is used to create a tunnel (channel).
The blades have a straight or curved working part.

3. Ophthalmic keratome blade for cutting


The blade serves to create a corneal incision.
The blades have a curved working part.

4. Ophthalmic keratome extended blade


A trapezoidal blade is used to widen the incision during lens implantation. The blades have a curved working part.

5. Ophthalmic microvitreoretinal blade

The blade is used for microvitreoretinal procedures for precise, uniform and sealable punctures.
The blades have the option of a straight and curved type of the working part.


6. Ophthalmic blade with controlled depth

The blade is used for gentle incision of the sclera.

2.3 Medical device design.

1. Blade. 2. Handle. 3. Plastic case. 4. Plastic cap.

2.4 Technical characteristics of the medical device


General scheme of measurements of technical characteristics.

1. Blade length
2. Blade width.
3. The angle of the working part of the blade.
4. Outer diameter of the blade tube.
5. Handle diameter.
6. Handle length.
7. Cap length
8. Cap diameter

General scheme for measuring technical characteristics.

1. Lancet-type ophthalmic blade	Blade			Handle (yellow color)	
	Width, mm, ± 0,18	Cutting angle, ± 3°	Length, mm, ± 2,5	Diameter r, mm, ±0,4	Length, mm, ± 5
1.1. Ophthalmic blade Lancet type 15°	1,75	15	11,2	6	120
1.2. Ophthalmic blade Lancet type 22.5 °	1,75	22,5	11,2	6	120
1.3. Ophthalmic blade Lancet type 30°	1,75	30	11,2	6	120
1.4. Ophthalmic blade Lancet type 45°	1,75	45	11,2	6	120
Weight, g, 3 ± 0,3					
2.Ophthalmic sickle blade	Blade		Handle (pink color)		
	Width, mm, ± 0,18	Length, mm, ±2,5	Diameter, mm, ±0.4	Length, mm, ± 5	
2.1. Ophthalmic sickle blade straight with bevel up, size 20 G	2,5	13	6	120	
2.2. Ophthalmic sickle blade straight with bevel down, size 20 G	2,5	13	6	120	
2.3. Ophthalmic sickle blade curved with bevel up, size 20 G	2,5	13	6	120	
2.4. Ophthalmic sickle blade curved with bevel down, size 20 G	2,5	13	6	120	
2.5. Ophthalmic sickle blade straight with bevel up, size 19 G	2,7	13	6	120	
2.6. Ophthalmic sickle blade straight with bevel down, size 19 G	2,7	13	6	120	
2.7. Ophthalmic sickle blade curved with bevel up, size 19 G	2,7	13	6	120	
2.8. Ophthalmic sickle blade curved with bevel down, size 19 G	2,7	13	6	120	
Angle of bending of the working part of the blade 45°, ± 3°					
Weight, g, 3 ± 0,3					
3. Ophthalmic keratome blade for cutting	Blade		Handle (blue color)		
	Width, mm, ± 0,18	Length, mm, ±2,5	Diameter, mm, ±0.4	Length, mm, ± 5	
					

3.1. Ophthalmic keratome blade for cutting, blade width 1.2 mm	1,2	13	6	120
3.2. Ophthalmic keratome blade for cutting, blade width 1.3 mm	1,3	13	6	120
3.3. Ophthalmic keratome blade for cutting, blade width 1.4 mm	1,4	13	6	120
3.4. Ophthalmic keratome blade for cutting, blade width 1.5 mm	1,5	13	6	120
3.5. Ophthalmic keratome blade for cutting, blade width 1.6 mm	1,6	13	6	120
3.6. Ophthalmic keratome blade for cutting, blade width 1.7 mm	1,7	13	6	120
3.7. Ophthalmic keratome blade for cutting, blade width 1.8 mm	1,8	13	6	120
3.8. Ophthalmic keratome blade for cutting, blade width 1.9 mm	1,9	13	6	120
3.9. Ophthalmic keratome blade for cutting, blade width 2.0 mm	2,0	13	6	120
3.10. Ophthalmic keratome blade for cutting, blade width 2.1 mm	2,1	13	6	120
3.11. Ophthalmic keratome blade for cutting, blade width 2.2 mm	2,2	13	6	120
3.12. Ophthalmic keratome blade for cutting, blade width 2.3 mm	2,3	13	6	120
3.13. Ophthalmic keratome blade for cutting, blade width 2.4 mm	2,4	13	6	120
3.14. Ophthalmic keratome blade for cutting, blade width 2.5 mm	2,5	13	6	120
3.15. Ophthalmic keratome blade for cutting, blade width 2.6 mm	2,6	13	6	120
3.16. Ophthalmic keratome blade for cutting, blade width 2.65 mm	2,65	13	6	120
3.17. Ophthalmic keratome blade for cutting, blade width 2.75 mm	2,75	13	6	120
3.18. Ophthalmic keratome blade for cutting, blade width 2.8 mm	2,8	13	6	120
3.19. Ophthalmic keratome blade for cutting, blade width 3.0 mm	3,0	13	6	120
3.20. Ophthalmic keratome blade for cutting, blade width 3.2 mm	3,2	13	6	120
3.21. Ophthalmic keratome blade for cutting, blade width 3.5 mm	3,5	13	6	120
3.22. Ophthalmic keratome blade for cutting the transparent cornea, blade width 2.65 mm	2,65	13	6	120
3.23. Ophthalmic keratome blade for cutting the transparent cornea, blade width 2.8 mm	2,8	13	6	120
3.24. Ophthalmic keratome blade for cutting the transparent cornea, blade width 3.2 mm	3,2	13	6	120

Angle of bending of the working part of the blade $45^\circ, \pm 3^\circ$

Blade point angle $50^\circ, \pm 3^\circ$

Weight, g, $3 \pm 0,3$

Ophthalmic keratome extended

	Blade		Handle (green color)	
	Width, mm, $\pm 0,18$	Length, mm, $\pm 2,5$	Diameter, mm, $\pm 0,4$	Length, mm, ± 5
4.1. Ophthalmic keratome extended blade, blade width 3.5 mm	3,5	13	6	120
4.2. Ophthalmic keratome extended blade, blade width 4.1 mm	4,1	13	6	120
4.3. Ophthalmic keratome extended blade, blade width 5.1 mm	5,1	13	6	120
4.4. Ophthalmic keratome extended blade, blade width 5.2 mm	5,2	13	6	120
4.5. Ophthalmic keratome extended blade, blade width 5.5 mm	5,5	13	6	120
4.6. Ophthalmic keratome extended blade, blade width 6.0 mm	6,0	13	6	120

Angle of bending of the working part of the blade $45^\circ, \pm 3^\circ$

Weight, g, $3,6 \pm 0,3$

5. Ophthalmic microvitrectomy blade

	Blade		Handle (pink color)	
	Width, mm, $\pm 0,18$	Length, mm, $\pm 2,5$	Diameter, mm, $\pm 0,4$	Length, mm, ± 5
5.1. Ophthalmic microvitrectomy blade straight, size 19 G	1,2	22	6	120
5.2. Ophthalmic microvitrectomy blade curved, size 19 G	1,2	22	6	120
5.3. Ophthalmic microvitrectomy blade straight, size 20 G	1,1	22	6	120
5.4. Ophthalmic microvitrectomy blade curved, size 20 G	1,1	22	6	120
5.5. Ophthalmic microvitrectomy blade straight, size 23G	1,0	22	6	120
5.6. Ophthalmic microvitrectomy blade curved, size 23G	1,0	22	6	120
5.7. Ophthalmic microvitrectomy blade straight, size 24 G	0,95	22	6	120
5.8. Ophthalmic microvitrectomy blade curved, size 24 G	0,95	22	6	120

Angle of bending of the working part of the blade $45^\circ, \pm 3^\circ$

Weight, g, $3,25 \pm 0,3$

6. Ophthalmic blade with controlled depth	Blade	Handle (pink color)	
	Punching accuracy, micron	Diameter, mm, ± 0.4	Length, mm, ± 5
6.1. Ophthalmic blade with controlled depth, 300 microns	300	6	97
6.2. Ophthalmic blade with controlled depth, 350 microns	350	6	97
6.3. Ophthalmic blade with controlled depth, 400 microns	400	6	97
6.4. Ophthalmic blade with controlled depth, 450 microns	450	6	97
6.5. Ophthalmic blade with controlled depth, 500 microns	500	6	97
6.6. Ophthalmic blade with controlled depth, 550 microns	550	6	97
6.7. Ophthalmic blade with controlled depth, 600 microns	600	6	97
Weight, g, $3,2 \pm 0,3$			
Cap			
Length, mm - 62 ± 2			
diameter mm - 9 ± 2			

The surface roughness of the cutting edges R_a is not more than 0.63 microns.
The roughness of the remaining surfaces R_a is not more than 1.25 μm .

Blade hardness according to Vickers - 52-54 HRC.

The width of the cutting edge of the products should be no more than 3 microns.

3. Indications for use:

Surgical operations on the organs of vision, as well as for punctures, incisions and delamination of eye tissues during surgical ophthalmic interventions with minimal damage to the tissues of the eyeball.

Contraindications for use.

There are no contraindications to the use of the medical device.

Precautionary measures:

- When using a medical device, precautions inherent in the procedure for performing an ophthalmic operation must be observed.
- For single use only.
- Do not use if package is opened or damaged.

4. How to work with the product.

1. Check the integrity of the blister pack.
2. Remove the medical device from the blister pack.
3. Remove the blade from the plastic case / slide the plastic cap along the handle.
4. Inspect the blade for defects and integrity violations.
5. Perform surgical manipulations in accordance with generally accepted medical practice.
6. Insert the blade back into the plastic case / slide the plastic cap back onto the blade.
7. After use, the medical device must be disposed of.

5. Terms of use, storage and transportation..

- Application conditions.

The product is operated at temperatures from 15 °C to 35 °C and humidity 50-60%.

- Storage conditions.

Store the medical product in its original packaging at a temperature of 0°C to 25°C in a dry room protected from direct sunlight at a humidity not exceeding 70%.

- Conditions of transportation.

The medical product is transported by all modes of transport in covered vehicles in accordance with the rules for the carriage of goods; applicable to every mode of transport.

Ambient temperature during transportation 0 from °C to 45°C, humidity not more than 80%. Avoid damage to corrugated cardboard boxes during transportation

6. Sterilization.

Medical device sterilized with ethylene oxide.

7. Packing.

- The blade without a protective cap is placed in a plastic case and packed in an individual package - a blister of 1 pc.



The blade with a protective cap is packed in an individual package - a blister of 1 pc.

- Individual (blister) packaging is placed in consumer packaging - a cardboard box of 10 pcs. Inserted instructions for use - 1 pc. Each box is sealed with a laminating film.

- Consumer packaging is placed in transport packaging - corrugated cardboard boxes - 50 pcs.

8. Marking.

Symbols used :

	Manufacturer		Date of manufacture		Do not re-sterilize
	Batch code		Use up to		Do not re-apply
	Marking CE approved		Sterilization ethylene oxide		Catalog number
	Non-toxic		Storage temperatur e range		Do not use if packaging is damaged
	Refer to instructions for use		Top		Keep away from moisture

	Do not expose to direct sunlight		Carefully! Refer to supporting documentation
H Stainless steel			

9. Standards applied in the manufacture of a medical device.

EN ISO 14971:2012; ISO 15223-1:2020; EN 1041:2008; EN ISO 14155: 2011; EN ISO 10993-1:2009; EN ISO 10993-5:2009; EN ISO 10993-7: 2008; ISO 10993-10: 2010; ISO 14644-1: 2015; EN ISO 11607-1:2009; EN ISO 11607-2:2006; EN ISO 15004-1:2009; EN ISO 11135-1: 2007; EN ISO 11138-2: 2009 EN ISO 11140-1: 2009; EN ISO 11737-1: 2006; EN ISO 11737-2: 2009; ISO 13485: 2003.

10. Disposal.

When handled on the territory of the Russian Federation after use, the medical device is disposed of as medical waste of class B in accordance with SanPIN 2.1.3684-21. A medical device with broken packaging and expired shelf life is disposed of as medical waste of Class A in accordance with SanPIN 2.1.3684-21.

11. Requirements for maintenance and repair.

Maintenance and repair requirements do not apply to this medical device. The medical product is disposable. After use, the medical device must be disposed of.

12. Manufacturer's Warranty

The manufacturer guarantees the sterility of the product, as well as its compliance with the requirements of international and national standards throughout the entire shelf life of the product, subject to the rules of transportation, storage and conditions of use in accordance with the instructions for use established by the manufacturer.

The shelf life of the product is 5 years from the date of sterilization.

Manufacturer's name and address

Microtrack Surgicals

A - 38, Adarsh - 2, Industrial Estate, B/h. Ashish Cinema, Odhav, Ahmedabad - 382 415, Gujarat, India. Phone No: - 079-22970972 Email ID: microtrack10@yahoo.com


Authorized representative on the territory of the Russian Federation:

Kom Span LLC.

Russia, 115035, Moscow, Kosmodarnianskaya nab., building 40-42, str.3, fl. 1, premises XVII, room. 1, of. 2.
tel./fax: (495) 434-45-06.

Date of last approval June 2024

FOR, MICROTRACK SURGICALS
[Signature]
PARTNER


EXECUTANT PARTNER
SIGNED BEFORE ME
[Signature]
Uday R. Patel
NOTARY
GOVT. OF GUJARAT
INDIA

ILABEN R. PATEL
NOTARY PUBLIC
GUJARAT
REGD. No. NTR/507/02
2428 3347 0311
मेरा आधाार, मेरी पहचान

ILABEN R. PATEL
NOTARY PUBLIC
GUJARAT
REGD. No. NTR/507/02
2428 3347 0311
मेरा आधाार, मेरी पहचान

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