"Approved"
Julia PERIMAN
INTERNATIONAL RA
Specialist
Julia PERIMAN
11-August-2015

Haemonetics® ACP®215 Cell Wash Procedure – Operation Manual –

HAEMONETICS CORPORATION
400 WOOD ROAD
BRAINTREE, MA 02184 USA
1-800-225-5242

CE 0123

Authorized European Representative
HAEMONETICS U.K. LTD.
5 Ashley Drive, Bothwell
Scotland G71 8DA
Great Britain

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Contents

ACP 215 Cell Wash Procedure

PRESENTING THE ACP 215 CELL WASH PROCEDURE.	1
Intended use	1
Listing procedure material.	1
Understanding procedure guidelines.	1
DESCRIBING THE CELL WASH DISPOSABLE SET	2
Solution lines	2
Bacteria barrier filter	2
Pump tubing with pump stops.	3
Blow molded bowl	3
DPM filter	3
SPM filter	3
Heat sealed tubing	3
Final product bag	3
Product waste bag.	3
PREPARING THE ACP 215 DEVICE FOR A PROCEDURE	4
System testing	4
Selecting a procedure	4
INSTALLING THE DISPOSABLE SET	6
Inspecting the disposable material	6
Time limits for disposable sets and solutions	6
Installing the final product bag and tubing	7
Installing the product waste bag and tubing	7
Installing the disposable bowl	7
Installing the pump tubing	8
Installing the solution lines and line sensor tubing	9
Installing the pressure monitor filters.	9
PERFORMING THE BOWL SEATING TEST	11
CONNECTING THE SOLUTIONS AND RBC PRODUCT.	14
VERIFYING THE DISPOSABLE SET INSTALLATION.	15
Initiating the disposable set checklist.	15
Color-coded valve tubing	16
Line sensor tubing.	16
Wash solution bag	17
Additive solution bag	17
SLAD tubing	17

P/N 85280-30, Manual revision: DA

Solution pump tubing	18
Blood pump tubing	18
BLAD tubing	18
Centrifuge bowl position	19
Additive solution seal	19
RBC bag tubing	19
ENTERING OPERATOR INITIALS	20
EXPLAINING MODIFIABLE PROCEDURE PARAMETERS	21
DESCRIBING THE CELL WASH PROCEDURE.	24
Verifying the sterile connection	24
Receiving the printed procedure summary	27
Depicting the procedure summary printout	29
Removing the final product	30
APPENDIX A: DEFINING THE MODIFIABLE PARAMETERS	32
APPENDIX B: TROUBLESHOOTING DURING A PROCEDURE	34
Listing notice messages	34

PRESENTING THE ACP 215 CELL WASH PROCEDURE

Intended use

The Haemonetics ACP 215 device is designed to wash red blood cells for removing plasma proteins. According to the validated conditions for CE marking, red blood cells suitable for washing with the ACP 215 device are CPD/SAG-M RBC units derived from whole blood or collected by apheresis and stored for up to 14 days. Red cells units after washing and stored in SAG-M at $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$ may be labeled for extended storage for 14 days.

Listing procedure material

The following material is needed to perform an ACP 215 cell wash procedure:

- The ACP 215 device and the printer.
- A sterile connection device.
- The Haemonetics disposable set LN235.
- Wash solution: saline glucose (0.9% saline 0.2% glucose) or saline (0.9% saline).

Optional:

- SAG-M solution (when extended storage is requested).

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The following guidelines should be respected concerning the RBC product to be washed:

- The initial unit of RBC product should be leukoreduced.
- The hematocrit of the initial RBC product should be between 50% and 80%.

DESCRIBING THE CELL WASH DISPOSABLE SET

1. Wash solution line
2. Additive solution line
3. Bacteria barrier filter
4. Solution pump tubing
5. Blood pump tubing
6. Blow molded bowl
7. DPM filter
8. SPM filter
9. Heat sealed tubing
10. Final product bag
11. Product waste bag


Figure 1, LN235 Cell wash disposable set

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Additive solution line

The orange-striped solution line contains a male luer connector to attach the selected additive solution bag to the disposable set.

Note: The orange-striped line is not used during certain cell wash procedures and the blue-striped line is not used at all for any cell wash procedure. The slide clamp on any unused solution line should remain closed at all times.

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Blow molded bowl

The disposable set contains a 275 ml blow molded bowl (BMB). A diverter is located in the bottom of the bowl, which introduces all fluids into the separation chamber at the outer diameter. The head of the bowl contains an *inlet port* connected to the blood line and an *outlet port* connected to the effluent line.

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DPM filter

A 0.2 μ m filter located on the blood line is used as a bacteria barrier between the draw pressure monitor and the disposable set. It is coupled with the receptor on the DPM to monitor pressure within the blood line from the blood pump to the RBC bag.

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SPM filter

A 0.2 μ m filter located on the effluent line is used as a bacteria barrier between the system pressure monitor and the disposable set. It is coupled with the receptor on the SPM to monitor pressure within the bowl and effluent tubing.

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ACP 215 SYSTEMS TESTS			
Rom	P	Air Sensors	P
Calibration	P	Analog	P
Config	P	Valves	P
Protocols	P	Watch Dog	P

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- ➔ Press the MODIFY PROGRAM key to scroll the list until the arrow indicates CELL WASH.
- ➔ Press the YES key to validate the selection.

Once the cell wash procedure has been selected, the following information will be displayed.

CELL WASH PROTOCOL	
PLEASE SELECT PROTOCOL OPTION	
--->	WITH ADDITIVE
	WITHOUT ADDITIVE

Figure 4, Cell wash procedure option menu

The operator should select if the washed RBC product will be resuspended with an additive solution or not.

- ➔ Press the MODIFY PROGRAM key to scroll the list until the arrow indicates the selection.
- ➔ Press the YES key to validate the selection.

Warning: If the option "WITHOUT ADDITIVE" is selected, it is the responsibility of the operator to ensure that, the quality of the final RBC product fulfills the local requirements in accordance with the local SOP.

Once the operator has pressed the YES key from the option menu, all valves will open and the following message will be displayed, along with target values based on the procedure programming in the row **TARG**.

CELL WASH PROTOCOL				
PLEASE LOAD DISPOSABLE SET PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 5, Disposable set installation message

The actual statistics will be updated and displayed during the procedure in the row **REAL**:

- PUMP = pump speed.
- TIME = procedure duration in minutes.
- WASH = wash solution used in ml.
- ADDI = additive solution used in ml.

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Note: This message will remain displayed during the disposable set installation sequence. The operator should install the disposable set completely, prior to pressing the YES key. Pressing the YES key at this point will display the bowl seating test, which should be initiated after the disposable set has been installed.

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The operator can proceed to install the disposable set as follows:

Figure 6, Position of the final product bag and tubing

- ➔ Place the open disposable set tub against the control panel with the head of the bowl facing down.
- ➔ Remove the final product bag from the tub and hang it from the small bag pins located on the left side of the front panel.
- ➔ Install the red-striped RBC bag tubing in the red-coded valve.
- ➔ Install the final product printed tubing into the white-coded valve.
- ➔ Remove the coiled solution lines and place them temporarily over the machine cover, to the left of the tub.

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- ➔ Remove the product waste bag and green-striped tubing from the tub.
- ➔ Hang the product waste bag on the large bag pins located on the right side of the front panel with the port facing out.

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- ➔ Loop the pump tubing segment between the blood pump stops around the groove of the blood pump rotor.
- ➔ Loop the pump tubing segment between the solution pump stops around the groove of the solution pump rotor with the clear pump stop towards the front of the device.
- ➔ Secure the tubing in each tubing guide so that each pump stop is located to the **right** of the tubing guide.
- ➔ Rotate each pump clockwise until the notch in the pump rotor catches the tubing and draws it completely into the pump groove.
- ➔ Rotate each pump clockwise another quarter of a turn to ensure proper placement of the tubing in the pump groove.
The pump should be able to rotate freely.

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
Figure 8, Pump stop and bacteria barrier filter orientation

Note: The bacteria barrier filter must be placed between the front solution pump tubing guide and the Y-connector leading to the bowl as depicted.

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- ➔ Install all three of the color-coded solution lines into the respective color-coded valves.
 - Yellow for the wash solution line.
 - Orange for the additive solution line.
 - Blue for the non-used solution line.

Note: The blue-striped line will not be connected to a solution bag during any cell wash procedure and the orange-striped line will not be connected to a solution bag during certain procedures. Nonetheless, these tubing sections must be installed in the appropriate valves for all of the cell wash procedures.

- ➔ Verify that the solution lines from the Y-connector are not restricted or kinked.
- ➔ Verify that there is some slack in the tubing on either side of the valves to ensure solution flow.
- ➔ Install the clear effluent line tubing into the line sensor, ensuring that the tubing is fully inserted into the groove.
- ➔ Close the line sensor cover.
- ➔ Install the green-striped tubing in the green-coded valve.

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Once the disposable set has been completely installed, the operator can proceed to the bowl seating test sequence.

P/N 85280-30, Manual revision: DA

PERFORMING THE BOWL SEATING TEST

The ACP 215 device will perform a test to ensure that the bowl is properly seated in the centrifuge chuck. The operator should initiate this test once the disposable set has been completely installed, as follows:

CELL WASH PROTOCOL				
PLEASE LOAD DISPOSABLE SET PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 9, Disposable set installation message

➔ Press the YES key.

All of the valves will close and the following message will be displayed.

BOWL SEATING TEST				
IS CENTRIFUGE COVER LOCKED PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 10, Bowl seating test message (1)

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BOWL SEATING TEST				
TESTING CENTRIFUGE				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 11, Bowl seating test message (2)

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BOWL SEATING TEST				
DID BOWL PRODUCE UNUSUAL NOISE PRESS YES OR NO TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 12, Bowl seating test message (3)

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➔ Press the NO key to continue.

The ACP 215 device will prompt the operator to initiate the disposable set checklist. At this point, the operator should proceed to the section *Connecting the solutions and RBC product*, **prior to** initiating the checklist.

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BOWL SEATING TEST				
CHANGE DISPOSABLE PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 13, Bowl seating test message (4)

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- ➔ Unlock and open the centrifuge cover.
- ➔ Remove the set completely.
- ➔ Install a new disposable set.
- ➔ Press the YES key again to receive the following message.

BOWL SEATING TEST				
PRESS DOWN ON BOWL PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 14, Bowl seating test message (5)

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BOWL SEATING TEST				
IS CENTRIFUGE COVER LOCKED PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 15, Bowl seating test message (1)

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- ➔ Press the NO key to continue.

The ACP 215 device will prompt the operator to initiate the disposable set checklist. At this point, the operator should proceed to the section *Connecting the solutions and RBC product*, **prior to** initiating the checklist.

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Using aseptic technique, the operator should:

- ➔ Spike the wash solution bag onto the yellow-striped line.
- ➔ Attach the luer lock connector of the additive solution to the orange-striped tubing connector (when applicable) and crack the breakable seal.
- ➔ Verify that the breakable seal is fully separated internally by bending the component back and forth in multiple axes until clearance is observed between the two broken halves.

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- ➔ Verify that the sterile connection weld is not occluded by squeezing the weld area between the fingers and rolling, to fully open the weld.
- ➔ Hang the RBC bag (port facing down) next to the final product bag on the front bags pins.
- ➔ Verify that each color-coded tubing segment is properly installed in the appropriate color-coded valve.
- ➔ Verify again that the sterile connection weld is fully open and that no kinks are evident in the tubing.

At this point, the operator can verify the installation of the disposable set according to the messages displayed during the disposable set checklist.

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Figure 16, Disposable set installation message

➔ Press the YES key to initiate the checklist.

A Full or Minimum checklist will be displayed, depending on the technical configuration of ACP 215 device. In addition, certain messages will be displayed only if the operator has selected to use an additive solution. The following table indicates with the letter **X** which messages will be displayed for each disposable set checklist configuration.

Table 1, Disposable set checklist configurations: Cell Wash procedure

Message	Full checklist		Minimum checklist	
	with additive solution	without additive solution	with additive solution	without additive solution
Color-coded valve tubing (6)	X	X		
Line sensor tubing	X	X	X	X
Wash solution bag	X	X		
Additive solution bag	X			
SLAD tubing	X	X		
Solution pump tubing	X	X		
Blood pump tubing	X	X		
BLAD tubing	X	X		
Centrifuge bowl position	X	X	X	X
Additive solution seal	X		X	
RBC bag tubing	X	X	X	X

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- Check, Press Yes to Continue				
BLUE striped tubing in BLUE valve				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 17, Example of a valve tubing message

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Line sensor tubing

- Check, Press Yes to Continue				
Tubing in line sensor Line sensor cover closed				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 18, Line sensor tubing message

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Wash solution bag

- Check, Press Yes to Continue				
YELLOW striped tubing connected to WASH solution bag				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 19, Wash solution connection message

- ➔ Verify that the yellow-striped solution tubing segment is connected to the wash solution bag.
- ➔ Correct any connection if incorrect.
- ➔ Press the YES key to continue.

Additive solution bag

- Check, Press Yes to Continue				
ORANGE striped tubing connected to additive solution bag				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 20, Additive solution connection message

- ➔ Verify that the orange-striped solution tubing is correctly connected to the additive solution bag.
- ➔ Correct the connection if incorrect.
- ➔ Press the YES key to continue.

SLAD tubing

- Check, Press Yes to Continue				
Solution tubing laced in solution air detector				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 21, SLAD tubing message

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Solution pump tubing

- Check, Press Yes to Continue				
SOLUTION pump loaded correctly filter downstream pump (CCW)				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 22, Solution pump tubing message

- ➔ Verify that the solution pump tubing is correctly installed.
- ➔ Reposition the tubing if incorrectly installed.
- ➔ Press the YES key to continue.

Blood pump tubing

- Check, Press Yes to Continue				
BLOOD pump loaded correctly filter upstream pump (CCW)				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 23, Blood pump tubing message

- ➔ Verify that the blood pump tubing is correctly installed.
- ➔ Reposition the tubing if incorrectly installed.
- ➔ Press the YES key to continue.

BLAD tubing

- Check, Press Yes to Continue				
BLOOD tubing laced in BLOOD air detector				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 24, BLAD tubing message

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Centrifuge bowl position

- Check, Press Yes to Continue				
Bowl is seated correctly in the centrifuge.				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 25, Centrifuge bowl position message

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Additive solution seal

- Check, Press Yes to Continue				
Additive solution seal is broken				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 26, Additive solution seal message

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RBC bag tubing

- Check, Press Yes to Continue				
Sterile connect RBC bag to RED Striped tubing.				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 27, RCB bag sterile connection message

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ENTER OPERATOR INITIALS A A				
PRESS YES TO SCROLL PRESS SAVE PROGRAM TO SELECT				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 28, Operator initials message (1)

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- Press the YES key to scroll the alphabet, letter by letter from A to Z.
- Press the SAVE PROGRAM key the when the appropriate letter is displayed.

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OPERATOR INITIALS J D				
PRESS YES TO CONTINUE PRESS NO TO RE-ENTER				
	PUMP	TIME	WASH	ADDI
REAL:	0	0	0	0
TARG:	0	32	1840	240

Figure 29, Operator initials message (2)

At this point the operator can change the entry as follows:

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Warning: Any modification to the Haemonetics default values remains the responsibility of the local site and should be in accordance with the local SOP. The device will retain in memory any settings which the operator decides to save, and will not automatically revert to the Haemonetics default values. However, when using SAG-M additive solution, the operator should change the "Final Preservative Vol" parameter to 80 ml and save it as the default value.

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- ➔ Press the START key.

The following table lists the full set of modifiable parameters. The parameters in **bold type** will be displayed when the device is configured for the **minimum set**. An asterisk (*) denotes the parameters displayed only during a procedure using additive solution.

A list of parameter definitions is provided in *Appendix A*.

Note: The default value provided for each wash solution volume is based on a four-step wash sequence.

Table 2, ACP 215 Modifiable parameters for the Cell Wash procedure

Parameter	Unit of measure	Default settings	Range of settings	Increment of change
RBC Unit Weight	g	550	200-600	5
RBC Unit Hematocrit	%	55	30-90	5
Max Bowl Fill Rate	ml/min	150	50-250	10
RBC Bowl Volume	ml	140	50-250	10
Centrifuge Speed	rpm	8000	3000-8000	50
Num Wash Steps	–	4	1-5	1
Wash Soln Vol (#1)	ml	230	10-1000	10
Wash Soln Vol (#2)	ml	230	10-1000	10
Wash Soln Vol (#3)	ml	230	10-1000	10
Wash Soln Vol (#4)	ml	230	10-1000	10
Wash Soln Vol (#5)	ml	0	0-1000	10
Wash Flow Rate (#1)	ml/min	200	10-250	10
Wash Flow Rate (#2)	ml/min	200	10-250	10
Wash Flow Rate (#3)	ml/min	200	10-250	10
Wash Flow Rate (#4)	ml/min	200	10-250	10
Wash Flow Rate (#5)	ml/min	200	10-250	10
Equilibration Time	s	30	5-600	5
Re-Centrifugation Time	s	60	15-600	5
* Preservative Soln Vol (#1)	ml	140	10-500	10
* Preservative Soln Vol (#2)	ml	100	10-500	10
* Preservative Soln Rate	ml/min	200	10-250	10

Table 2, ACP 215 Modifiable parameters for the Cell Wash procedure

Parameter	Unit of measure	Default settings	Range of settings	Increment of change
* Final Preservative Vol	ml	100	10-500	10
Preservative Ratio	%	75	5-95	5
Bowl Conc. Rate	ml/min	60	10-250	10
Line Sensor Trip%	%	70	1-99	1
Bowl Empty Rate	ml/min	200	10-250	10
Min Empty Bowl Vol	ml	130	25-130	5
Conc LS Delay Vol	ml	80	30-100	5

DESCRIBING THE CELL WASH PROCEDURE

Once the operator has pressed the START key, the ACP 215 device will automatically perform the cell wash procedure according to the programmed parameters.

First, air will be purged from the final product bag and the disposable set will be primed with wash solution. Then, blood will be drawn from the RBC bag to fill the bowl up to a maximum RBC volume.

The wash solution will then be sent to the bowl over a series of steps to remove the plasma protein from the packed RBC layer and flush the effluent into the product waste bag. If the operator has selected to use additive solution, this will also be sent to the bowl according to the procedure programming.

In the case of a larger RBC unit, the bowl will be partially emptied to permit another complete wash of the remaining unwashed cells. Once all of the RBCs have been washed, a concentration cycle will be performed to optimize the hematocrit of the final product.

Verifying the sterile connection

Before the final wash cycle, once the RBC bag is empty, the ACP 215 device will evaluate the integrity of the sterile connection weld by holding and monitoring a negative pressure for a period of time. During this test, the following information will be displayed.

CELL WASH PROTOCOL				
STERILE CONNECTION CHECK IN PROGRESS. PLEASE WAIT				
	PUMP	TIME	WASH	ADDI
REAL:	0	20	920	140
TARG:	0	32	1840	240

Figure 31, Sterile connection test message (1)

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Occasionally, a naturally occurring pocket of air in the product bag may move as a result of the negative pressure in the tubing and cause a fluctuation of pressure. This pressure change might also indicate a leak at the connection site, so the integrity of the sterile connection will be controlled again.

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CELL WASH PROTOCOL				
STERILE CONNECTION CHECK IN PROGRESS. PLEASE WAIT				
	PUMP	TIME	WASH	ADDI
REAL:	0	20	920	140
TARG:	0	32	1840	240

Figure 32, Sterile connection test message (1)

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STERILE CONNECTION FAILURE				
PLEASE CHECK BLOOD BAG FOR AIR PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	20	920	140
TARG:	0	32	1840	240

Figure 33, Sterile connection test message (2)

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- Reposition the bag to ensure that the air is displaced to the top of the bag.

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CELL WASH PROTOCOL				
STERILE CONNECTION CHECK IN PROGRESS. PLEASE WAIT				
	PUMP	TIME	WASH	ADDI
REAL:	0	20	920	140
TARG:	0	32	1840	240

Figure 34, Sterile connection test message (1)

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CELL WASH PROTOCOL				
STERILE CONNECTION FAILURE PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	20	920	140
TARG:	0	32	1840	240

Figure 35, Sterile connection test message (3)

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→ Press the YES key to receive the following message.

STERILE CONNECTION FAILURE				
PRODUCT LIMITED TO 24 HOURS PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	20	920	140
TARG:	0	32	1840	240

Figure 36, Sterile connection failure message

Once this message is displayed:

→ Press the YES key to receive the following message.

PRODUCT LIMITED TO 24 HOURS				
PRESS NO TO ABORT PROCEDURE PRESS START TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	20	920	140
TARG:	0	32	1840	240

Figure 37, Product limitation warning message

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CELL WASH PROTOCOL				
PROCEDURE ABORTED PLEASE UNLOAD DISPOSABLE SET				
	PUMP	TIME	WASH	ADDI
REAL:	0	20	920	140
TARG:	0	32	1840	240

Figure 38, Procedure termination message

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- ➔ Press the START key.

A statement will appear on the procedure summary printout informing the operator that the storage of the blood product is limited to 24 hours as follows:

WARNING *****				
BLOOD PRODUCT SHELF LIFE				
LIMITED TO 24 HOURS.				

Figure 39, Product limitation warning on the procedure printout

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Once the concentrated, washed RBC product has been completely transferred to the final product bag, along with any final additive solution, all valves will close and the pumps will stop. The ACP 215 device will beep and the following message will be displayed:

CELL WASH PROTOCOL				
IS PRINTER ONLINE PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	32	1840	240
TARG:	0	32	1840	240

Figure 40, Printer verification message

At this point, the operator should:

- ➔ Verify that the printer is attached to the rear panel of the device.
- ➔ Verify that the printer is online.
- ➔ Press the Yes to receive the printed procedure summary.

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IS THE PRINTOUT ACCEPTABLE?				
PRESS YES TO CONTINUE PRESS NO TO RETRY				
	PUMP	TIME	WASH	ADDI
REAL:	0	32	1840	240
TARG:	0	32	1840	240

Figure 41, Printout confirmation message

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CELL WASH PROTOCOL				
IS PRINTER ONLINE PRESS YES TO CONTINUE				
	PUMP	TIME	WASH	ADDI
REAL:	0	32	1840	240
TARG:	0	32	1840	240

Figure 42, Printer verification message

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The following example depicts the information found on the cell wash procedure summary printout. For certain items, the values are generated automatically by the ACP 215 device. For other items, a blank space is provided for the operator to write the information. The name of the manufacturer should be included when documenting the lot number of each item used for the procedure.

HAEMONETICS ACP 215	
SERIAL NUMBER: 95B072	
WASH PROCEDURE SUMMARY	
Date & Time:	_____
Unit Number:	_____
RBC Unit Weight:	250 g
RBC Unit Hematocrit:	65%
RBC Unit Col. Date:	_____
RBC Unit Exp Date:	_____
Estimated Wash Volume:	1840 ml
Estimated Additive Solution:	240 ml
Final Product Weight:	_____ g
Final Product Hct:	_____ %
Procedure Time:	32 Min.
Operator Initials	A A
Disposable Lot Num:	_____
Exp. Date:	_____
0.9% NaCl Lot Num:	_____
Exp. Date:	_____
Additive Solution:	_____
Lot Num:	_____
Exp. Date:	_____
Wafer Lot Num:	_____
Operator Signature:	_____

Figure 43, Example of a procedure summary printout

Removing the final product

Once the procedure summary has been printed, the operator will receive the PROCEDURE COMPLETE message and should remove the final product as follows.

- ➔ Heat seal the tubing 10 cm [4 inches] from the bag.
- ➔ Detach the final product bag from the set and package it in accordance with the local SOP.
- ➔ Store the product at 2° C - 6° C in accordance with the local SOP.
- ➔ Remove the used disposable set and dispose of in accordance with the local SOP.

Under certain circumstances, sterile air may not be completely displaced through the disposable set, therefore some blood may remain in the bowl at the end of the procedure. In this case, the operator can manually transfer the blood in the bowl to the final product bag as follows.

- ➔ Remove the green-striped tubing from the green valve.
- ➔ Manipulate the waste bag as follows.
 - ➔ *Grasp the waste bag on each side near the bottom corners.*
 - ➔ *Gently squeeze the bag to push any air pockets above the fluid towards the inlet port of the bag, ensuring that no waste fluid reaches the port and exits the bag.*
 - ➔ *Release the bag once the air is pocket is observed at the port.*
- ➔ Remove the final product bag tubing from the white valve.
- ➔ Manually rotate the blood pump counterclockwise until the bowl is empty.

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
Figure 44, ACP 215 procedure selection menu

APPENDIX A: DEFINING THE MODIFIABLE PARAMETERS

RBC UNIT WEIGHT	The weight in grams of the red blood cells to be processed. Determined by weighing the RBC bag and its contents prior to the wash procedure, then taring the weight of the empty bag. Any samples removed prior to the wash procedure wash should be excluded from this weight.
RBC UNIT HEMATOCRIT	The measured hematocrit of the well-mixed unit of RBCs prior to the wash procedure.
MAX BOWL FILL RATE	The rate in ml/ min at which the blood is pumped from the unwashed RBC bag into the centrifuge bowl.
RBC BOWL VOLUME	The absolute RBC volume in ml that is pumped into the bowl during the first wash cycle, calculated from the RBC Unit Weight and Hct.
CENTRIFUGE SPEED	The number of revolutions per minute at which the centrifuge spins when operating at full speed.
NUMBER WASH STEPS	The number of times a batch of cells is washed before being pumped into the final washed product bag.
WASH SOLN VOL (#1-5)	The volume in ml of wash solution that is used to wash cells between one cycle of starting and stopping the centrifuge.
WASH FLOW RATE (#1-5)	The rate in ml/min of wash solution that is pumped through the bowl during each wash step.
EQUILIBRATION TIME	The time in seconds that the bowl and centrifuge are stopped in between wash steps.
RE-CENTRIFUGATION TIME	The time in seconds that the bowl spins before wash solution is pumped through it.
PRESERV. SOLN VOL (#1)	The volume in ml of additive solution that is pumped through the bowl before the first batch of washed blood is pumped into final product bag.
PRESERV. SOLN VOL (#2)	In a two cycle wash procedure, the volume in ml of additive solution that is pumped through the bowl before the second batch of washed blood is pumped into the final product bag.
PRESERVATIVE SOLN RATE	The rate in ml/min of additive solution that is pumped through the bowl at the end of each series of wash steps.
FINAL PRESERVATIVE VOL	The volume of additive solution that is pumped into the final washed product bag after all the blood has been pumped into the product bag.

PRESERVATIVE RATIO	The ratio of additive solution to total volume of supernate, indicating what percent of the supernate is additive solution after 100 ml have been added.
BOWL CONCENTRATION RATE	The rate in ml/min at which the first batch of RBCs is pumped back into the bowl.
LINE SENSOR TRIP %	The % transmittance of a threshold value that triggers a software response to detection of RBCs passing through the line sensor tubing.
BOWL EMPTY RATE	The rate in ml/min at which the washed product is pumped from the bowl to the final product bag.
MIN EMPTY BOWL VOL	The volume at which the bowl contents will be emptied to the final product bag, to allow the remaining volume of unwashed RBCs to be pumped into the bowl and processed.
CONC LS DELAY VOL	A programmed volume that is not monitored by the line sensor during the concentration process.

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Listing notice messages

During the cell wash procedure, the ACP 215 device may interrupt operation and display particular messages, to indicate a situation that can easily be corrected by the operator. The following table list these notice messages along with a brief explanation of the situation and the recommended operator actions.

INSTALLATION NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- CHECK CORRECT DISPOSABLE LOADED - INSTALL DPM FILTER PRESS YES TO CONTINUE	The DPM filter has not been detected on the donor pressure monitor.	➔ Press the filter firmly onto the DPM while providing a quarter turn clockwise to ensure a tight seal. ➔ Ensure DPM sleeve is properly attached. ➔ Press the YES key to resume operation.
- CHECK CORRECT DISPOSABLE LOADED - INSTALL SPM FILTER PRESS YES TO CONTINUE	The SPM filter has not been detected on the system pressure monitor.	➔ Press the filter firmly onto the SPM while providing a quarter turn clockwise to ensure a tight seal. ➔ Ensure SPM sleeve is properly attached. ➔ Press the YES key to resume operation.
- CHECK CORRECT DISPOSABLE LOADED - INSTALL TUBING IN LINE SENSOR PRESS YES TO CONTINUE	The effluent tubing has not been detected as the line sensor.	➔ Insert the tubing completely in the groove of the line sensor and close the cover.
- BOWL SEATING TEST - BOWL NOT SEATED PROPERLY PRESS YES TO CONTINUE	The centrifuge has spun the bowl and the test has not passed.	➔ Verify that the bowl is completely seated in the centrifuge chuck. ➔ Reposition if necessary. ➔ Press the YES key to resume operation.

INSTALLATION NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- IS FINAL PRODUCT BAG EMPTY? - Press YES to start prime Press NO to continue purge	Residual air is in the final product bag during the pre-Prime state.	The operator has two options: ➔ If the final product bag is empty, press the YES key to proceed to prime the disposable set. or ➔ If the final product bag still contains air, press the NO key to continue purging air from the bag.

PRIME STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- AIR AT SLAD - Check tubing installation & flow Press START to continue	The SLAD did not detect fluid during the Prime state.	➔ Inspect the solution line for occlusions and ensure that the tubing path to the SLAD is clear. ➔ Ensure that the wash solution bag is spiked. ➔ Press the START key to resume operation.
- AIR AT BLAD - Check tubing installation & flow Press START to continue	The BLAD detected air during the Prime state.	➔ Inspect the blood line for occlusions and ensure that the tubing path is clear. ➔ Ensure the BLAD is loaded properly. ➔ Press the START key to resume operation.
<i>Initial message</i> - INLET PRESSURE OUTSIDE LIMITS - Check for occlusion in line ***** <i>Corrected state</i> - PRESSURE RELIEVED - Press START to continue	The DPM detected pressure outside of limits during the Prime state.	➔ Inspect the inlet harness lines for occlusions. ➔ Ensure the BLAD is loaded properly. ➔ If there are no occlusions, manually open the red valve briefly. Once the pressure is relieved: ➔ Press the START key to resume operation.

PRIME STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- AIR AT BLAD - Check flow in solution line Press START to continue	The BLAD detected air before starting to fill the bowl.	➔ Inspect the solution line for occlusions and ensure that the tubing path is clear. ➔ Press the START key to resume operation. If the message is repeated: ➔ Inspect the disposable set for a leak. If a leak is detected: ➔ Discontinue the procedure and handle the unprocessed cells according to the local SOP.

FILL STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- AIR AT BLAD - Press START to continue Press NO to skip to next step	The BLAD detected air before a programmed bowl blood volume was reached.	The operator has two options: ➔ Press the START key to continue filling the bowl. or ➔ Press the NO key to advance to the next step of the procedure.
- FILL LONGER THAN EXPECTED - Press START to continue Press NO to skip to next step	Blood remains in the RBC bag after a programmed bowl blood volume was reached.	The operator has two options: ➔ If the bag is not empty, press the START key to continue filling the bowl. or ➔ If the bag is empty, press the NO key to advance to the next step of the procedure.
<i>Initial message</i> - BOWL PRESSURE OUTSIDE LIMITS - Check for effluent line occlusion <i>Corrected state</i> - PRESSURE RELIEVED - Press START to continue	The SPM detected pressure outside of acceptable limits as the bowl was filling.	➔ Inspect the effluent harness lines for occlusions. Once the pressure is relieved: ➔ Press the START key to resume operation.

FILL STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
<i>Initial message</i> -INLET PRESSURE OUTSIDE LIMITS- Check for occlusion in line ***** <i>Corrected state</i> - PRESSURE RELIEVED - Press START to continue	The DPM detected pressure outside of acceptable limits as the bowl was filling.	➔ Inspect the inlet harness line for occlusions. ➔ Check the blood line tubing for occlusions. Once the pressure is relieved: ➔ Press the START key to resume operation.

WASH STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- AIR AT SLAD - Check flow in solution line Press START to continue	The SLAD detected air in the wash solution line.	➔ Inspect the solution line for occlusions and ensure that the tubing path is clear. ➔ Ensure that the wash solution bag is spiked correctly and has fluid present. ➔ Press the START key to resume operation.
- AIR AT SLAD - Check flow in preservative line Press START to continue	The SLAD detected air in the additive solution line.	➔ Inspect the additive solution line for occlusions and ensure that the tubing path is clear. ➔ Ensure that the additive solution bag is correctly connected. ➔ Press the START key to resume operation.

EMPTY STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
<i>Initial message</i> -INLET PRESSURE OUTSIDE LIMITS- Check for occlusion in line ***** <i>Corrected state</i> - PRESSURE RELIEVED - Press START to continue	The DPM detected pressure outside of acceptable limits as the bowl was emptying.	→ Inspect the inlet harness lines for occlusions. → If there are no occlusions, manually open the red valve briefly. Once the pressure is relieved: → Press the START key to resume operation.
- EMPTY LONGER THAN EXPECTED - Press START to continue Press NO to skip to next step	Blood remains in the bowl when it should be completely emptied.	The operator has two options: → If the bowl is not empty, press the START key to continue emptying the bowl. or → If the bowl is empty, press the NO key to advance to the next step of the procedure.
- AIR DETECTED EARLY - Press Start to continue Press NO to skip to next step	BLAD detected air early during the Empty state.	The operator has two options: → Press the START key to continue operation. or → Press the NO key to advance to the next step of the procedure.
- RBCs IN EFFLUENT LINE? - Press YES to end concentration Press NO to resume filling	The line sensor has detected foam or RBCs during the concentration process.	→ Inspect the effluent tubing to determine the presence of RBCs or foam. The operator has two options: → If the line is opaque, press the YES key to end the concentration process. or → If no RBCs are observed in the line, press the NO key to continue the concentration process.

Yulia PERIMAN
INTERNATIONAL RA Specialist
Yulia PERIMAN
11- August -2015

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HAEMONETICS 215 SERIAL NUMBERS: 95B072
WASH PROCEDURE SUMMARY

Date & Time:
Unit Number:
RBC Unit Weight: 250 g
RBC Unit Hematocrit: 65%
RBC Unit Col. Date:
RBC Unit Exp Date:
Estimated Wash Volume: 1840 ml
Estimated Additive Solution: 240 ml
Final Product Weight: _____ g
Final Product Hct: _____ %
Procedure Time: 32 Min.
Operator Initials A A
Disposable Lot Num:
Exp. Date:
0.9% NaCl Lot Num:
Exp. Date:
Additive Solution:
Lot Num:
Exp. Date:
Wafer Lot Num:
Operator Signature:

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CE 0123

EC REP

HAEMONETICS U.K. LTD.
5 Ashley Drive, Bothwell
Scotland G71 8BS
Great Britain

Printed in France

Haemonetics Corporation

400 Wood Road, Braintree, MA 02184, USA

©2001, 2004, 2007 Haemonetics Corporation. All rights reserved.

P/N 85279-30, Manual revision: DA
July 2007

Table of Contents

ACP 215 Deglycerolization Procedure

PRESENTING THE ACP® 215 DEGLYCEROLIZATION PROCEDURE	1
Intended use	1
Listing procedure material.	1
Understanding procedure guidelines.	1
DESCRIBING THE DEGLYCEROLIZATION DISPOSABLE SET	4
Solution lines	4
Bacteria barrier filter	4
Pump tubing with pump stops.	5
Blow molded bowl	5
DPM filter	5
SPM filter	5
Heat sealed tubing	5
Final product bag	5
Product waste bag.	5
PREPARING THE ACP 215 DEVICE FOR A PROCEDURE	6
System testing	6
Selecting a procedure	6
INSTALLING THE DISPOSABLE SET	8
Inspecting the disposable material	8
Time limits for disposable sets and solutions	8
Installing the final product bag and tubing	9
Installing the product waste bag and tubing	9
Using the chuck adaptor (LN235 only).	10
Installing the disposable bowl	11
Installing the pump tubing	11
Installing the solution lines and line sensor tubing	13
Installing the pressure monitor filters.	13
PERFORMING THE LINE SENSOR CALIBRATION TEST	14
Explaining the line sensor test sequence	14
PERFORMING THE BOWL SEATING TEST	18
CONNECTING THE SOLUTIONS AND RBC PRODUCT.	21
VERIFYING THE DISPOSABLE SET INSTALLATION.	22
Initiating the disposable set checklist.	22
Color-coded valve tubing	23

P/N 85279-30, Manual revision: DA

Line sensor tubing	23
Hypertonic solution bag	24
Wash solution bag	24
Additive solution bag	24
SLAD tubing	25
Solution pump tubing	25
Blood pump tubing	25
BLAD tubing	26
Centrifuge bowl position	26
Additive solution seal	26
RBC bag on the shaker	27
ENTERING OPERATOR INITIALS	28
EXPLAINING MODIFIABLE PROCEDURE PARAMETERS	29
DESCRIBING THE DEGLYCEROLIZATION PROCEDURE	32
Verifying the sterile connection	32
Receiving the printed procedure summary	35
Depicting the procedure summary printout	36
Depicting the estimated free hemoglobin chart	38
Removing the final product	39
APPENDIX A: DEFINING THE MODIFIABLE PARAMETERS	41
APPENDIX B: TROUBLESHOOTING DURING A PROCEDURE	43
Listing notice messages	43

PRESENTING THE ACP[®] 215 DEGLYCEROLIZATION PROCEDURE

Intended use

The ACP 215 is intended to glycerolize and deglycerolize red blood cells derived from whole blood, or collected by apheresis, that have been stored in an approved anticoagulant/additive solution for up to 7 days at 1–6°C prior to glycerolization.

The conditions that have been used for the CE marking validation are as follows: red cells derived from whole blood units and collected in CPD, glycerolized by the ACP 215, and frozen at –65° C or colder, deglycerolized by the ACP 215 and stored in SAG-M at 4° C ± 2° C may be labeled for extended storage for 72 hours.

Listing procedure material

The following material is needed to perform an ACP 215 deglycerolization procedure:

- The ACP 215 device, the shaker and the printer.
- A sterile connection device.
- The Haemonetics[®] disposable set LN235 (275 ml) or LN236 (325 ml).
- An 800-1000 ml RBC storage bag (qualified for the storage of frozen cells) for the transfer of a frozen RBC product stored in a 600 ml bag.
- Wash solution (0.9% saline/0.2% dextrose).
- Hypertonic Solution (12% NaCl).
- SAG-M solution.

Understanding procedure guidelines

- The hematocrit of the thawed RBC unit should be approximately 60% (range 44-70%).
- The initial glycerolized RBC unit must have been stored at –60° C to –80° C and must be properly thawed according to local regulations and the local SOP.

Note: Haemonetics recommends placing a secondary, watertight cover around the glycerolized RBC unit when warmed in a water bath, to prevent the ports from becoming contaminated.

- If the glycerolized RBC unit has been frozen in a 600 ml storage bag, the unit must be transferred to a qualified 800-1000 ml storage bag after thawing, prior to the deglycerolization procedure.

Warning: Improper thawing of the frozen RBC product may damage the red blood cells thereby producing high supernatant hemoglobin content. All welds between the thawed RBC bag, the transfer bag and the deglycerolization disposable set must be made using a sterile connection device.

Guidelines for choosing Deglycerolization Disposable/Bowl Size

The choice of bowl size depends on the red cell mass of the collected product.

In order to conserve Red Cells and avoid spillage into waste, it is recommended that only units at or below an RBC mass of 180ml be washed using the LN235 (275 ml bowl) deglycerolization disposable. The LN236 (325 ml bowl) deglycerolization disposable should be used for units with a larger RBC mass.

Table 1 below summarizes when the 325 ml versus the 275 ml bowl should be used and also details the expected final hematocrit of the product. The shaded area indicates where the 325 ml bowl (LN236) should be used for deglycerolization.

Table 1, Expected final product Hct

RBC mass (mL) in	Final HCT	
	325 ml Bowl	275 ml Bowl
150	37	42
160	39	44
170	41	47
180	44	50
190	46	Use 325 ml bowl
200	49	Use 325 ml bowl
210	51	Use 325 ml bowl
220	51	Use 325 ml bowl
230	51	Use 325 ml bowl
240	51	Use 325 ml bowl
250	51	Use 325 ml bowl
260	51	Use 325 ml bowl
270	51	Use 325 ml bowl
280	51	Use 325 ml bowl

Note: Depending on recovery of the red cells the hematocrit of the final product may be slightly lower.

Red Cell Mass estimate when the donor hematocrit is known.

When the hematocrit of the donor is known, the table below can be used to determine whether the 275 ml bowl (LN235) or the 325 ml bowl (LN236) should

be used. The red cell mass is based on the weight of the whole blood drawn and the donor hematocrit. To determine the red cell mass of the unit, locate the donor's hematocrit in the first column of the table.

Move across the row until you locate the appropriate column that denotes the weight of the whole blood unit drawn. **If the RBC mass >180 mL: use the 325 ml bowl.** The shaded yellow area indicates where the 325 ml bowl (LN236) should be used for deglycerolization.

Table 2, Red cell mass (donor hematocrit known)

		Weight of whole blood (gm)															
HCT		400	410	420	430	440	450	460	470	480	490	500	510	520	530	540	550
38		143	147	151	154	158	161	165	168	172	176	179	183	186	190	194	197
40		151	155	158	162	166	170	174	177	181	185	189	192	196	200	204	208
42		158	162	166	170	174	178	182	186	190	194	198	202	206	210	214	218
44		166	170	174	178	183	187	191	195	199	203	208	212	216	220	224	228
46		174	178	182	187	191	195	200	204	208	213	217	221	226	230	234	239
48		181	186	190	195	199	204	208	213	217	222	226	231	235	240	245	249
50		189	193	198	203	208	212	217	222	226	231	236	241	245	250	255	259
52		196	201	206	211	216	221	226	231	235	240	245	250	255	260	265	270

If the red cell mass is not known and the donor's hematocrit is not known.

Haemonetics recommends that the site determine the post-thaw RBC volume and post-thaw RBC hematocrit of the unit prior to deglycerolization. The red cell mass can be estimated by using the formula provided below.

$$\text{Weight} / 1.09 = \text{Volume}$$

$$\text{Volume} \times \text{Hematocrit} = \text{Red Cell Mass}$$

If the red cell mass of the unit is greater than 180mLs, the 325 ml bowl (LN 236) should be utilized for the deglycerolization process.

DESCRIBING THE DEGLYCEROLIZATION DISPOSABLE SET

1. Hypertonic solution line
2. Wash solution line
3. Additive solution line
4. Bacteria barrier filter
5. Solution pump tubing
6. Blood pump tubing
7. Blow molded bowl (LN235, 275 ml bowl & LN236, 325 ml bowl)
8. DPM filter
9. SPM filter
10. Heat sealed tubing
11. Final product bag
12. Product waste bag


Figure 1, LN235 Deglycerolization disposable set

Solution lines

Wash solution line

The yellow-striped solution line contains a spike which is used to connect the wash solution bag to the disposable set.

Additive solution line

The orange-coded solution line contains a male luer connector to attach the selected additive solution bag to the disposable set.

Hypertonic solution line

The blue-striped solution line contains a spike which is used to connect the hypertonic solution bag to the disposable set.

Bacteria barrier filter

A 0.2 μm is integrated on the set, downstream of the solution pump. This filter acts as a barrier to prevent potential contaminants from passing into the system after spiking or connecting any of the solution bags.

Pump tubing with pump stops

The pump tubing sections contain stops which will secure the pump tubing during pump loading and operation.

- The **blood** pump tubing section contains **two blue** pump stops.
- The **solution** pump tubing contains **one clear** stop and **one blue** stop.

Blow molded bowl

The disposable set contains either a 275 ml blow molded bowl (BMB) or a 325 ml blow molded bowl. A diverter is located in the bottom of the bowl, which introduces all fluids into the separation chamber at the outer diameter. The head of the bowl contains an *inlet port* connected to the blood line and an *outlet port* connected to the effluent line.

When installing a bowl, the inlet port should be positioned facing the left side of the centrifuge. The outlet port, located below the inlet port should be positioned facing the right side of the centrifuge.

DPM filter

A 0.2 µm filter located on the blood line is used as a bacteria barrier between the draw pressure monitor and the disposable set. It is coupled with the receptor on the DPM to monitor pressure within the blood line from the blood pump to the RBC bag.

In case of any occlusion in the line or any obstruction generated by the sterile connection process, the pressure monitor will stop the pump and alert the operator.

SPM filter

A 0.2 µm filter located on the effluent line is used as a bacteria barrier between the system pressure monitor and the disposable set. It is coupled with the receptor on the SPM to monitor pressure within the bowl and effluent tubing.

In case of high pressure, the blood or solution flow into the bowl will be stopped, preserving the integrity of the bowl seal.

Heat sealed tubing

The end of the blood line tubing is heat sealed to achieve a closed system and is connected to the RBC bag using a sterile connection device.

Final product bag

This bag has a 600 ml capacity and is used to accommodate the final product at the end of the procedure.

Product waste bag

This bag has a 5 liter capacity and is used to collect the effluent that exits the bowl. This effluent is the waste of the different solutions used to wash the RBCs.

Warning: The operator must remain attentive to the capacity of the product waste bag when modifying any of the Haemonetics default values for the procedure parameters. The total volume of solutions used for the complete procedure must not exceed the 5 liter capacity.

PREPARING THE ACP 215 DEVICE FOR A PROCEDURE

Prior to initiating a procedure, the operator should:

- ➔ Verify that the power cord is attached to the device.
- ➔ Verify that the biohazard waste bag is attached and hanging freely.
- ➔ Power-on the device.

System testing

Immediately after the ACP 215 device has been powered-on, a series of internal self-diagnostic tests will be conducted. The letter **P** will be displayed next to each test after it has passed.

ACP 215 SYSTEMS TESTS				
Rom	P		Air Sensors	P
Calibration	P		Analog	P
Config	P		Valves	P
Protocols	P		Watch Dog	P

Figure 2, Example of the ACP 215 system test display

Once the self-tests have been completed, the operator will need to verify for the safety system that the centrifuge cover can be properly locked. A series of messages will instruct the operator to lock, unlock then re-lock the cover until this safety check is complete.

When the operator has completed the test, the ACP 215 software revision level will be displayed briefly, followed by the procedure selection menu.

Selecting a procedure

Once the procedure selection menu is displayed, the operator should:

HAEMONETICS ACP 215	
PLEASE SELECT PROTOCOL	
==>	CELL WASH
	GLYCEROLIZATION
	DE-GLYCEROLIZATION

Figure 3, ACP 215 procedure selection menu

- ➔ Press the MODIFY PROGRAM key to scroll the list until the arrow indicates DEGLYCEROLIZATION as follows.

HAEMONETICS ACP 215	
PLEASE SELECT PROTOCOL	
	CELL WASH
	GLYCEROLIZATION
→	DE-GLYCEROLIZATION

Figure 4, ACP 215 procedure selection menu

➔ Press the YES key to validate the selection.

Once the deglycerolization procedure has been selected, the following message will be displayed, along with target values based on the procedure programming in the row **TARG**.

DE-GLYCEROLIZATION PROTOCOL					
PLEASE LOAD DISPOSABLE SET PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 5, Disposable set installation message

The actual statistics will be updated and displayed during the procedure in the row **REAL**.

- PUMP = pump speed.
- TIME = procedure duration in minutes.
- HYPE = hypertonic solution used in ml.
- NACL = wash solution used in ml.
- ADDI = additive solution used in ml.

At this point, the operator can proceed to install the disposable set.

Note: This message will remain displayed during the disposable set installation sequence. The operator should install the disposable set completely, prior to pressing the YES key. Pressing the YES key at this point will display the line sensor calibration test, which should be initiated after the disposable set has been installed.

INSTALLING THE DISPOSABLE SET

Inspecting the disposable material

The operator should inspect the disposable material prior to, as well as during installation, using the following guidelines:

- ➔ Verify that the disposable set corresponds with the selected procedure.
- ➔ Verify that neither the packaging tub nor Tyvek® cover has been damaged.
- ➔ Remove the Tyvek cover from the disposable set tub and inspect all tubing sections.
- ➔ Ensure that no kinks or occlusions are present which could obstruct the flow through the disposable set.
- ➔ Verify that there are no visible defects or particulate within the set.

If a disposable set is damaged, it should be returned to Haemonetics Corporation as specified in the section, *Returned Goods Authorization system*, of the *ACP 215 Automated Cell Processing System Operation Manual*.

Warning: A disposable set should not be used if the packaging has been damaged during shipping or storage, since this could compromise the sterility of the collected blood products.

Time limits for disposable sets and solutions

Disposable sets may be removed from the protective packaging, installed on the device, and primed for up to 8 hours before a procedure is initiated. If solutions have not been attached to the disposable set, the opened, dry set should be used within 24 hours.

The protective outer wrap on the solution bags controls water evaporation. In order to assure that the solution remains within its specification limits, the solutions should be used and the procedure initiated within 8 hours after the outer wrap has been removed.

Installing the final product bag and tubing

The operator should secure the thawed RBC unit on the shaker with the shaker magnets as depicted and proceed to install the disposable set as follows:


Figure 6, Position of the thawed RBC unit, final product bag and tubing

- ➔ Place the open disposable set tub against the control panel with the head of the bowl facing down.
- ➔ Remove the final product bag from the tub and hang it from the small bag pins located on the left side of the front panel.
- ➔ Install the red-striped RBC bag tubing in the red-coded valve.
- ➔ Install the final product printed tubing into the white-coded valve.
- ➔ Remove the coiled solution lines and place them temporarily to the left of the tub.

Installing the product waste bag and tubing

- ➔ Remove the product waste bag and green-striped tubing from the tub.
- ➔ Hang the product waste bag from the large bag pins located on the right side of the front panel.
- ➔ Install the green-striped tubing in the green-coded valve.

Warning: The bowl will not be emptied correctly if there is an occlusion in the green-striped tubing disturbing the free path of sterile air between the product waste bag and the bowl, or if the product waste bag is removed prior to the completion of the procedure.

Using the chuck adaptor (LN235 only)

The LN235 disposable set (275 ml bowl), requires a chuck adaptor. The chuck adaptor is a white plastic ring that snaps into the centrifuge chuck. Install the chuck adaptor before loading the disposable set. Removing the chuck after a procedure requires a special removal tool.

The chuck adaptor (P/N 51348-00) and removal tool (ETC-6136), are available separately.

1. Chuck adaptor
2. Removal tool
3. Arms


Figure 7: Chuck adaptor and removal tool

Use the following steps to install the chuck adaptor:

- Place the chuck adaptor inside the centrifuge chuck.
- Press down on the adaptor until it clicks into place.

1. Chuck adaptor
2. Removal tool


Figure 8: Chuck adaptor and removal tool

Use the following steps to remove the chuck adaptor:

- Compress the arms of the removal tool until it fits into the chuck adaptor.
- Push down on the tool until it clicks into place.
- Using two hands, pull upwards to remove the chuck adaptor.
- Push the tool through the bottom of the chuck adaptor to separate.

Installing the disposable bowl

At this point, the operator can easily remove the bowl and install it as follows:

Figure 9, Installing the bowl

- ➔ Open the centrifuge cover.
- ➔ Insert the bowl in the centrifuge chuck, ensuring that the effluent port is directed to the **right**.
- ➔ Press firmly down on the shoulder of the bowl until the bowl is completely seated in the chuck.
- ➔ Close and secure the centrifuge cover.
- ➔ Verify that the centrifuge knob is in the **locked** position.

Installing the pump tubing

At this point, the operator can install the pump tubing.

- ➔ Loop the pump tubing segment between the blood pump stops around the groove of the blood pump rotor.
- ➔ Loop the pump tubing segment between the solution pump stops around the groove of the solution pump rotor.
- ➔ Secure the tubing in each tubing guide so that each pump stop is located to the **right** of the tubing guide.
- ➔ Rotate each pump clockwise until the notch in the pump rotor catches the tubing and draws it completely into the pump groove.
- ➔ Rotate each pump clockwise another quarter of a turn to ensure proper placement of the tubing in the pump groove.
The pump should be able to rotate freely.

1. SLAD
2. Solution tubing pump stops
3. Solution pump
4. Pump tubing guides
5. Bacteria barrier filter
6. Blood pump
7. Blood tubing pump stops
8. BLAD


Figure 10, Pump stop and bacteria barrier filter orientation

Note: The bacteria barrier filter must be placed between the solution pump tubing guide and the Y-connector leading to the bowl as depicted.

- ➔ Thread the clear tubing between the blood pump and the Y-connector attached to the DPM filter through the blood line air detector (BLAD).
- ➔ Thread the clear tubing between the solution pump and the Y-connector of the solution lines through the solution line air detector (SLAD).
- ➔ Verify that the pump stops are situated against the tubing guides once the tubing has been installed in the air detectors.

Warning: The operator should respect all standard precautions applicable to rotating machinery when installing the pump tubing.

Installing the solution lines and line sensor tubing

The operator should install the remaining tubing segments as follows:

- ➔ Install all three of the color-coded solution lines into the respective color-coded valves.
 - *Yellow for the wash solution line.*
 - *Orange for the additive solution line.*
 - *Blue for the hypertonic solution line.*
- ➔ Verify that the solution lines from the Y-connector are not restricted or kinked.
- ➔ Verify that there is some slack in the tubing on either side of the valves to ensure solution flow.
- ➔ Install the clear effluent line tubing into the line sensor, ensuring that the tubing is fully inserted into the groove.
- ➔ Close the line sensor cover.

Installing the pressure monitor filters

- ➔ Connect the DPM filter to the draw pressure monitor.
- ➔ Connect the SPM filter to the system pressure monitor.
- ➔ Press each filter firmly onto the DPM and SPM while providing a quarter turn clockwise to ensure a tight seal.

Warning: If for any reason, the DPM filter becomes wet during the procedure, the device may not be able to correctly detect pressure when emptying the RBC bag, resulting in a message **FILL LONGER THAN EXPECTED**.

If, at any time during a procedure, the SPM filter becomes wet, the device may not be able to monitor the pressure in the bowl. Due to a potential inability to inadequately monitor pressure in the bowl, the sterile pathway should be considered compromised and the product life should be reduced to 24 hours.

After completion of the procedure, Haemonetics recommends that the operator evaluate and record on the procedure record whether the SPM filter has become wet.

Once the disposable set has been completely installed, the operator can proceed to the line sensor calibration test sequence.

PERFORMING THE LINE SENSOR CALIBRATION TEST

The optical line sensor monitors the effluent exiting the bowl and evaluates the free hemoglobin during the deglycerolization procedure. The line sensor has the capability to categorize the free hemoglobin level to be in the following ranges:

- Less than 50 mg%.
- Between 50 and 100 mg%.
- Between 100 and 150 mg%.
- Greater than 150 mg%.

The purpose of the line sensor calibration test is to ensure that the optical sensor provides readings within specified limits.

*Note: The Line Sensor is intended to assist the operator in **estimating** the level of free hemoglobin in the deglycerolized product. Although the ACP 215 system provides an estimate of free hemoglobin level on the procedure summary printout, the final decision that the free hemoglobin level is acceptable is the responsibility of the operator, not the ACP 215 device. The operator should follow the local SOP when determining free hemoglobin levels and document the decision on the deglycerolization summary printout.*

Explaining the line sensor test sequence

Once the disposable set has been completely installed, the operator should:

DE-GLYCEROLIZATION PROTOCOL					
PLEASE LOAD DISPOSABLE SET PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 11, Disposable set installation message

➔ Press the YES key to receive the following message.

CALIBRATION OF LINE SENSOR					
PRESS YES TO CALIBRATE Press NO TO BYPASS CALIBRATION					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 12, Line sensor test message (1)

At this point, all valves will close and the ACP 215 device will prompt the operator to decide whether or not to perform the line sensor calibration test at this point.

Haemonetics recommends that the line sensor calibration test be performed at least once each day of use, for example prior to starting the first procedure of the day. If the line sensor passes the calibration check, then the operator can bypass the test for any subsequent procedures performed during that day as follows.

➔ Press the NO key.

The line sensor calibration test will be bypassed and the device will proceed directly to the bowl seating test, described in the following section of the manual.

If the operator selects to perform the test at this point:

➔ Press the YES key to initiate the test sequence.

The following message will be displayed and the operator should:

CALIBRATION OF LINE SENSOR					
OPEN COVER AND REMOVE TUBING PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 13, Line sensor test message (2)

➔ Open the line sensor cover and remove the tubing.

➔ Press the YES key while the cover is open.

The line sensor will automatically conduct a self-test using the standard atmosphere as a filter. The following message will be displayed and the operator should continue as follows:

CALIBRATION OF LINE SENSOR					
FILTER IS IN RANGE PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 14, Line sensor test message (3)

➔ Press the YES key to receive the following message.

CALIBRATION OF LINE SENSOR					
INSERT FILTER 2 INTO LINE SENSOR PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 15, Line sensor test message (4)

- ➔ Place the special Haemonetics filter #2 into the line sensor groove.
- ➔ Press the YES key, ensuring that the filter blocks the line sensor optical beam when the key is pressed.

The line sensor will automatically conduct a second self-test and the following message will be displayed:

CALIBRATION OF LINE SENSOR					
FILTER 2 IS IN RANGE PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 16, Line sensor test message (5)

- ➔ Press the YES key to receive the following message.

CALIBRATION OF LINE SENSOR					
INSERT FILTER 3 INTO LINE SENSOR PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 17, Line sensor test message (6)

- ➔ Place the special Haemonetics filter #3 into the line sensor groove.
- ➔ Press the YES key, ensuring that the filter blocks the line sensor optical beam when the key is pressed.

The line sensor will automatically conduct a third self-test and the following message will be displayed.

CALIBRATION OF LINE SENSOR					
LINE SENSOR TEST COMPLETE YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 18, Line sensor test message (7)

➔ Press the YES key to receive the following message.

CALIBRATION OF LINE SENSOR					
INSTALL TUBING IN LINE SENSOR PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 19, Line sensor test message (8)

When this message is displayed, the operator should:

- ➔ Re-install the line sensor tubing completely in the line sensor groove.
- ➔ Close the line sensor cover and press the YES key to continue.

At this point, the ACP 215 device will proceed to the bowl seating test.

However, if during any of the three self-tests, the line sensor is detected to be out of range, the device will stop operation and display the following message:

CALIBRATION OF LINE SENSOR					
FILTER IS NOT IN RANGE CALL HAEMONETICS FIELD SERVICE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 20, Line sensor out of range message

If this message is displayed, the operator should contact Haemonetics Field Service to calibrate the line sensor.

Note: When ordering replacement line sensor calibration test filters from Haemonetics, the following part numbers should be specified:

- Filter #2 as P/N 50062-00.
- Filter #3 as P/N 50339-00.

PERFORMING THE BOWL SEATING TEST

The ACP 215 device will perform a test to ensure that the bowl is properly seated in the centrifuge chuck. Once the line sensor calibration test has been completed, or bypassed, the following message will be displayed.

BOWL SEATING TEST					
IS CENTRIFUGE COVER LOCKED PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 21, Bowl seating test message (1)

At this point, the operator may verify the installation of the bowl prior to pressing the YES key to initiate the test, as follows:

- ➔ Unlock and open the centrifuge cover.
- ➔ Press firmly on the shoulder of the bowl to confirm that the bowl is completely seated in the chuck.
- ➔ Close and lock the centrifuge cover.
- ➔ Press the YES key to perform the bowl seating test.

Warning: The operator must NOT open the centrifuge cover while the centrifuge is operating.

The centrifuge will spin the bowl at 8000 rpm for approximately 30 seconds as the following message is displayed.

BOWL SEATING TEST					
TESTING CENTRIFUGE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 22, Bowl seating test message (2)

After completing the bowl seating test, the centrifuge will stop automatically. The device will prompt the operator to continue with the following message.

BOWL SEATING TEST					
DID BOWL PRODUCE UNUSUAL NOISE PRESS YES OR NO TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 23, Bowl seating test message (3)

If the operator determines that the bowl seating test ran smoothly without any unusual noise:

- ➔ Press the NO key to continue.

The ACP 215 device will prompt the operator to initiate the disposable set checklist. At this point, the operator should proceed to the following section, *Connecting the solutions and RBC product*, **prior to** initiating the checklist.

However, if the operator determines that the bowl was producing a loud or high-pitched noise during the bowl seating test:

- ➔ Press the YES key to receive the following message.

BOWL SEATING TEST					
CHANGE DISPOSABLE PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 24, Bowl seating test message (4)

At this point, the operator should change the disposable set as follows:

- ➔ Unlock and open the centrifuge cover.
- ➔ Remove the set completely.
- ➔ Install a new disposable set.
- ➔ Press the YES key to receive the following message.

BOWL SEATING TEST					
PRESS DOWN ON BOWL PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 25, Bowl seating test message (5)

To ensure proper seating of the bowl inside of the centrifuge chuck:

- ➔ Press firmly on the shoulder of the bowl.
- ➔ Press the YES key to continue.

Warning: An incorrectly loaded bowl may generate particulate.

The following message will be re-displayed.

BOWL SEATING TEST					
IS CENTRIFUGE COVER LOCKED PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 26, Bowl seating test message (1)

The operator should:

- ➔ Verify that the centrifuge cover is closed and locked.
- ➔ Press the YES key to continue.

The ACP 215 device will repeat the bowl seating test with the new disposable set. If the operator determines that the bowl seating test ran smoothly without any unusual noise:

- ➔ Press the NO key to continue.

The ACP 215 device will prompt the operator to initiate the disposable set checklist. At this point, the operator should proceed to the following section, *Connecting the solutions and RBC product*, **prior to** initiating the checklist.

CONNECTING THE SOLUTIONS AND RBC PRODUCT

Using aseptic technique, the operator should:

- ➔ Spike the wash solution bag onto the yellow-striped line.
- ➔ Spike the hypertonic solution bag onto the blue-striped line.
- ➔ Attach the luer lock connector of the additive solution to the orange-striped tubing connector and crack the breakable seal.
- ➔ Verify that the breakable seal is fully separated internally by bending the component back and forth in multiple axes until clearance is observed between the two broken halves.

Warning: Any color-coded tubing that is connected to an improper solution bag may result in product damage.

Any incorrect connection or restriction in the solution flow may result in product damage.

The operator must verify the expiration date of the solution prior to use to ensure that it is sterile

- ➔ Perform a sterile connection weld between the red-striped line and the tubing of the RBC bag.

When using a sterile connecting device, the operator must ensure that the diameters of the tubing being connected meet the diameter specifications as listed in the operating instructions of the sterile connection device.

Note: If the RBC bag has two tubing segments, the operator should select the segment that does not contain a breakable seal. If this segment is not long enough to perform a proper weld, the segment with a breakable seal can be used, however, the operator should verify that the seal is fully separated prior to performing a weld.

- ➔ Verify that the sterile connection weld is not occluded by squeezing the weld area between the fingers and rolling, to fully open the weld.
- ➔ Verify that each color-coded tubing segment is properly installed in the appropriate color-coded valve.
- ➔ Verify again that the sterile connection weld is fully open and that no kinks are evident in the tubing.

At this point, the operator can verify the installation of the disposable set according to the messages displayed on the disposable set checklist.

Warning: An improperly installed disposable set may result in product damage.

VERIFYING THE DISPOSABLE SET INSTALLATION

Initiating the disposable set checklist

The ACP 215 device will prompt the operator to verify the installation of the set through a sequence of messages referred to as the *disposable set checklist*. Once the bowl seating test has been successfully completed, the operator should:

DE-GLYCEROLIZATION PROTOCOL					
START DISPOSABLE CHECK LIST PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 27, Disposable set installation message

➔ Press the YES key to initiate the checklist.

A Full or Minimum checklist will be displayed, depending on the technical configuration of ACP 215 device. The following table indicates with the letter **X** which messages will be displayed for each disposable set checklist configuration.

Table 3, Disposable set checklist configurations: Deglycerolization procedure

Message	Full checklist	Minimum checklist
Color-coded valve tubing (6)	X	
Line sensor tubing	X	X
Hypertonic solution bag	X	
Wash solution bag	X	
Additive solution bag	X	
SLAD tubing	X	
Solution pump tubing	X	
Blood pump tubing	X	
BLAD tubing	X	
Centrifuge bowl position	X	X
Additive solution seal	X	X
RBC bag tubing	X	X
RBC bag/shaker	X	X

Color-coded valve tubing

The ACP 215 device will display a sequence of messages prompting the operator to verify the correct installation of each color-coded tubing segment in the appropriate color-coded valve, as depicted in the following example:

- Check, Press Yes to Continue					
BLUE striped tubing in BLUE valve					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 28, Example of a valve tubing message

For each message, the operator should:

- ➔ Verify that the color-coded tubing is installed in the appropriate valve.
- ➔ Reposition the tubing if incorrectly installed.
- ➔ Press the YES key to proceed.

The list will proceed in the following order.

- Blue-striped tubing in the blue valve.
- Yellow-striped tubing in the yellow valve.
- Orange-striped tubing in the orange valve.
- Red-striped tubing in the red valve.
- Clear tubing in the white valve.
- Green-striped tubing in the green valve.

Line sensor tubing

- Check, Press Yes to Continue					
Tubing in line sensor. Line sensor cover closed.					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 29, Line sensor tubing message

- ➔ Verify that the line sensor tubing is fully inserted in the line sensor groove.
- ➔ Reposition the tubing if incorrectly installed.
- ➔ Verify that the line sensor cover is closed.
- ➔ Press the YES key to continue.

Hypertonic solution bag

- Check, Press Yes to Continue					
BLUE Striped tubing connected to Hypertonic Solution bag.					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		57	50	1470	240

Figure 30, Hypertonic solution connection message

- ➔ Verify that the blue-striped solution tubing segment is connected to a hypertonic solution bag.
- ➔ Correct any connection if incorrect.
- ➔ Press the YES key to continue.

Wash solution bag

- Check, Press Yes to Continue					
YELLOW striped tubing connected to WASH solution bag					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 31, Wash solution connection message

- ➔ Verify that the yellow-striped solution tubing segment is connected to the wash solution bag.
- ➔ Correct the connection if incorrect.
- ➔ Press the YES key to continue.

Additive solution bag

- Check, Press Yes to Continue					
ORANGE striped tubing connected to additive solution bag					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 32, Additive solution connection message

- ➔ Verify that the orange-striped solution tubing is connected to the additive solution bag.
- ➔ Correct the connection if incorrect.
- ➔ Press the YES key to continue.

SLAD tubing

- Check, Press Yes to Continue					
Solution tubing laced in solution air detector					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 33, SLAD tubing message

- ➔ Verify that the tubing is correctly installed in the SLAD.
- ➔ Reposition the tubing if incorrectly installed.
- ➔ Press the YES key to continue.

Solution pump tubing

- Check, Press Yes to Continue					
SOLUTION pump loaded correctly. filter downstream pump (CCW).					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 34, Solution pump tubing message

- ➔ Verify that the solution tubing is correctly installed in the solution pump.
- ➔ Reposition the tubing if incorrectly installed.
- ➔ Press the YES key to continue.

Blood pump tubing

- Check, Press Yes to Continue					
BLOOD pump loaded correctly. filter upstream pump (CCW).					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 35, Blood pump tubing message

- ➔ Verify that the blood tubing is correctly installed in the blood pump.
- ➔ Reposition the tubing if incorrectly installed.
- ➔ Press the YES key to continue.

BLAD tubing

- Check, Press Yes to Continue					
BLOOD tubing laced in BLOOD air detector					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 36, BLAD tubing message

- ➔ Verify that the tubing is correctly installed in the blood line air detector.
- ➔ Reposition the tubing if incorrectly installed.
- ➔ Press the YES key to continue.

Centrifuge bowl position

- Check, Press Yes to Continue					
Bowl is seated correctly in the centrifuge.					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 37, Centrifuge bowl position message

- ➔ Visually verify that the bowl is correctly installed.
- ➔ Press the YES key to continue.

Additive solution seal

- Check, Press Yes to Continue					
Additive solution seal is broken					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 38, Additive solution seal message

- ➔ Verify that the seal is broken at the port of the additive solution bag.
- ➔ Break the seal if it is unbroken.
- ➔ Press the YES key to continue.

RBC bag tubing

- Check, Press Yes to Continue					
Sterile connect RBC bag to RED Striped tubing.					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 39, RCB bag sterile connection message

- ➔ Verify that the RBC product bag is sterile connected to the red-striped tubing segment of the disposable set.
- ➔ Perform a sterile connection if not present.
- ➔ Ensure that the sterile connection weld is fully open.
- ➔ Press the YES key to continue.

RBC bag on the shaker

- Check, Press YES To Continue					
Secure RBC Bag to the shaker.					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		57	50	1470	240

Figure 40, RBC bag/shaker message

- ➔ Verify that the thawed glycerolized RBC unit is properly secured to the shaker.
- ➔ Correct the placement of the bag and secure with the magnets if necessary.
- ➔ Press the YES key to continue.

ENTERING OPERATOR INITIALS

Following the disposable set checklist, the device will immediately prompt the operator to enter his/her initials in the space provided on the top portion of the message (AA in the following example).

ENTER OPERATOR INITIALS A A					
PRESS YES TO SCROLL PRESS SAVE PROGRAM TO SELECT					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 41, Operator initials message (1)

The letter which can be modified will be blinking and the operator should:

- ➔ Press the YES key to scroll the alphabet, letter by letter from A to Z.
- ➔ Press the SAVE PROGRAM key the when the appropriate letter is displayed.

Once the first letter has been saved, the cursor will move to the second letter. The operator should repeat the procedure for the second letter.

After the operator initials has selected two letters, the device will prompt the operator to confirm the initials with the following message.

OPERATOR INITIALS J D					
PRESS YES TO CONTINUE PRESS NO TO RE-ENTER					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 42, Operator initials message (2)

At this point the operator can change the entry as follows.

- ➔ Press the NO key.

The device will repeat the procedure and the cursor will blink next to the first initial. The operator can repeat the sequence with the corrected entry.

If the initials are acceptable:

- ➔ Press the YES key to continue.

EXPLAINING MODIFIABLE PROCEDURE PARAMETERS

Once the operator has confirmed the initials for the current procedure, the following message will be displayed.

DE-GLYCEROLIZATION PROTOCOL					
PRESS MODIFY TO SET PARAMETERS PRESS START TO BEGIN					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	0	0	0	0
TARG:		55	50	1470	240

Figure 43, Modifiable parameters message

At this point, certain procedure parameters are accessible to the operator. Each ACP 215 device is delivered with pre-installed settings tested by Haemonetics, called the *default* values. The operator can either use the default values or adjust them for the present procedure. Depending on the technical configuration of the device, either a **minimum set** or a **full set** of parameters will be available.

Warning: Any modification to the Haemonetics default values remains the responsibility of the local site and should be in accordance with the local SOP. The device will retain in memory any settings which the operator decides to save, and will not automatically revert to the Haemonetics default values.

To initiate the current procedure without consulting the modifiable parameters:

- ➔ Press the START key.

To access and verify or modify the available parameters for the current procedure:

- ➔ Press the MODIFY PROGRAM key to view the current settings.
- ➔ Use the Arrow keys to change the value as necessary.
- ➔ Press the MODIFY PROGRAM key to access the next parameter available.

To retain any modified values for the subsequent procedure(s):

- ➔ Press the SAVE PROGRAM key.

To initiate the procedure at this point:

- ➔ Press the START key.

The following table lists the full set of modifiable parameters. The parameters in **bold type** will be displayed when the device is configured for the **minimum set**. A list of parameter definitions is provided in *Appendix A*.

Note: The default value provided for each wash solution volume is based on a five-step wash sequence.

Table 4, ACP 215 Modifiable parameters for the Deglycerolization procedure

Parameter	Unit of measure	Default setting	Range of settings	Increment of change
Thawed Unit Weight	g	340	115-570	5
Thawed Unit HCT	%	60	40-100	10
Min Flow Rate	ml/min	50	10-250	10
Bowl Volume	ml	275	125-500	5
Centrifuge Speed	rpm	8000	3000-8000	50
Thawed Unit Osmolality	mOsm/kgH ₂ O	5100	100-10000	100
Hyper Soln Volume	ml	50	0-500	10
Hyper Soln Flow Rate	ml/min	150	10-250	10
Hyper Equalize Delay	s	150	10-1000	10
Dilution to Wash Delay	s	60	5-1000	5
Osmolality Rate	mOsm/kgH ₂ O/min	500	50-1000	50
Num Saline Dilution	–	2	0-5	1
NaCl Dilute Vol (#1)	ml	340	10-1000	10
NaCl Dilute Vol (#2)	ml	400	10-1000	10
Constant Dilution Rate	–	NO	YES-NO	–
Dilution Rate	ml/min	60	10-250	10
NaCl Rinse Vol (#1)	ml	50	10-500	10
Max Bowl Fill Rate	ml/min	150	100-250	10
Min Bowl Fill Rate	ml/min	50	10-100	10
Num Wash Cycles	–	5	0-9	1
Wash Soln Vol (#1)	ml	60	10-1000	10
Wash Soln Vol (#2)	ml	70	10-1000	10
Wash Soln Vol (#3)	ml	100	10-1000	10
Wash Soln Vol (#4)	ml	150	10-1000	10
Wash Soln Vol (#5)	ml	300	10-1000	10
Wash Flow Rate (#1)	ml/min	50	10-250	5
Wash Flow Rate (#2)	ml/min	50	10-250	5

Table 4, ACP 215 Modifiable parameters for the Deglycerolization procedure

Parameter	Unit of measure	Default setting	Range of settings	Increment of change
Wash Flow Rate (#3)	ml/min	75	10-250	5
Wash Flow Rate (#4)	ml/min	75	10-250	5
Wash Flow Rate (#5)	ml/min	100	10-250	5
Equilibration Time	s	45	5-1000	5
Re-Centrifugation Time	s	60	5-1000	5
Num Resuspensions	–	1	0-1	1
Preserv. Soln Vol	ml	240	10-500	10
Preservative Soln Rate	ml/min	120	10-250	10
Bowl Empty Rate	ml/min	200	10-250	10

DESCRIBING THE DEGLYCEROLIZATION PROCEDURE

Once the operator has pressed the START key, the ACP 215 device will automatically perform the deglycerolization procedure according to the programmed parameters.

First, air will be purged from the final product bag and the disposable set will be primed with wash solution. Then, the glycerolized RBC unit on the shaker will be diluted with hypertonic and wash solutions, before being transferred to the spinning bowl for centrifugation.

The RBCs in the bowl will be rinsed with wash solution as the effluent is sent to the product waste bag. The packed RBC layer in the bowl will then be washed over a series of steps before being suspended in additive solution and sent to the final product bag.

Verifying the sterile connection

At the end of the first fill cycle, the ACP 215 device will evaluate the integrity of the sterile connection weld by holding and monitoring a negative pressure for a period of time. During this test, the following information will be displayed.

DE-GLYCEROLIZATION PROTOCOL					
STERILE CONNECTION CHECK IN PROGRESS. PLEASE WAIT					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	17	50	340	0
TARG:		57	50	1470	240

Figure 44, Sterile connection test message (1)

A steady pressure level indicates a robust weld and operation will continue normally.

Occasionally, a naturally occurring pocket of air in the product bag may move as a result of the negative pressure in the tubing and cause a fluctuation of pressure. This pressure change might also indicate a leak at the connection site, so the integrity of the sterile connection will be checked again.

The blood pump will rotate at a rate to produce the required negative pressure level and the following information will be displayed.

DE-GLYCEROLIZATION PROTOCOL					
STERILE CONNECTION CHECK IN PROGRESS. PLEASE WAIT					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	17	50	340	0
TARG:		57	50	1470	240

Figure 45, Sterile connection test message (1)

If the negative pressure is maintained as required, the procedure will continue normally.

However, if the pressure level is again **not** maintained according to the specific criteria, the following information will be displayed.

STERILE CONNECTION FAILURE					
PLEASE CHECK BLOOD BAG FOR AIR PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	17	50	340	0
TARG:		57	50	1470	240

Figure 46, Sterile connection test message (2)

At this point, the operator can assess the blood bag for a valid leak or excess air. If excess air is observed:

- ➔ Reposition the bag to ensure that the air is displaced to the top of the bag.

To continue the procedure:

- ➔ Press the YES key.

The ACP 215 device will repeat the sterile connection test as the following information is displayed.

DE-GLYCEROLIZATION PROTOCOL					
STERILE CONNECTION CHECK IN PROGRESS. PLEASE WAIT					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	17	50	340	0
TARG:		57	50	1470	240

Figure 47, Sterile connection test message (1)

If the test is successfully completed, the procedure will continue normally.

However, if the test is again **unsuccessful**, the following information will be displayed to notify the operator of a leak in the system.

DE-GLYCEROLIZATION PROTOCOL					
STERILE CONNECTION FAILURE PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	17	50	340	0
TARG:		57	50	1470	240

Figure 48, Sterile connection test message (3)

At this point the operator should proceed as follows:

- ➔ Press the YES key to receive the following message.

STERILE CONNECTION FAILURE					
PRODUCT LIMITED TO 24 HOURS PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	17	50	340	0
TARG:		57	50	1470	240

Figure 49, Sterile connection failure message

Once this message is displayed:

- ➔ Press the YES key to receive the following message.

PRODUCT LIMITED TO 24 HOURS					
PRESS NO TO ABORT PROCEDURE PRESS START TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	17	50	340	0
TARG:		57	50	1470	240

Figure 50, Product limitation warning message

When this message is displayed, the operator can either choose to terminate the procedure directly, or choose to continue the procedure knowing that the storage life of the product will be reduced to 24 hours.

To terminate the procedure:

- ➔ Press the NO key to receive the following message.

DE-GLYCEROLIZATION PROTOCOL					
PROCEDURE ABORTED PLEASE UNLOAD DISPOSABLE SET					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	17	50	340	0
TARG:		57	50	1470	240

Figure 51, Procedure termination message

To continue with the procedure normally, the operator should:

- ➔ Press the START key.

A statement will appear on the procedure summary printout informing the operator that the storage of the blood product is limited to 24 hours as follows.

WARNING *****
BLOOD PRODUCT SHELF LIFE
LIMITED TO 24 HOURS.

Figure 52, Product limitation warning on the procedure printout

Receiving the printed procedure summary

Once the deglycerolized RBC product has been completely transferred to the final product bag, along with any final additive solution, all valves will close and the pumps will stop. The ACP 215 device will beep and the following message will be displayed:

DE-GLYCEROLIZATION PROTOCOL					
IS PRINTER ONLINE PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	57	50	1470	240
TARG:		57	50	1470	240

Figure 53, Printer verification message

At this point, the operator should:

- ➔ Verify that the printer is attached to the rear panel of the device.
- ➔ Verify that the printer is online.
- ➔ Press the YES key to receive the printed procedure summary.

At the end of the printing operation, the ACP 215 device will display the following message, prompting the operator to respond about the quality of the printout.

IS THE PRINTOUT ACCEPTABLE?					
PRESS YES TO CONTINUE PRESS NO TO RETRY					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	57	50	1470	240
TARG:		57	50	1470	240

Figure 54, Printout confirmation message

At this point the operator will have the following options:

- ➔ Press the YES key if the printout is acceptable and terminate the procedure.
- ➔ Press the NO key if the printout is not acceptable.

The following message will be re-displayed.

DE-GLYCEROLIZATION PROTOCOL					
IS PRINTER ONLINE PRESS YES TO CONTINUE					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	57	50	1470	240
TARG:		57	50	1470	240

Figure 55, Printer verification message

- ➔ Press the YES key to print another copy of the procedure summary.

Depicting the procedure summary printout

The following example depicts the information found on the deglycerolization procedure summary printout. For certain items, the values are generated automatically by the ACP 215 device. For other items, a blank space is provided for the operator to write the information. The name of the manufacturer should be included when documenting the lot number of each item used for the procedure.

HAEMONETICS ACP 215	
SERIAL NUMBER: 95B072	
THAWED BLOOD DEGLYCEROLIZATION	
PROCEDURE SUMMARY	
Date and Time:	_____
Unit Number:	_____
Hypertonic Volume:	50 ml
Shaker:	ON
First Dilution Volume:	340 ml
Shaker:	ON
Bowl Fill:	800 ml
RBC Spillage:	NO
Second Dilution Volume:	400 ml
Shaker:	ON
Bowl Fill:	900 ml
RBC Spillage:	NO
1st Wash Volume:	60 ml
2nd Wash Volume:	70 ml
3rd Wash Volume:	100 ml
4th Wash Volume:	150 ml
5th Wash Volume:	300 ml
Additive Solution:	240 ml
Final Product: Estimated	
Free Hemoglobin Level < :	50 mg%
Procedure Time:	56 Min.
Operator Initials	A A
Thawed unit temperature:	_____
12% NaCl Lot Num:	_____
Exp. Date:	_____
0.9% NaCl Lot Num:	_____
Exp. Date:	_____
Additive Solution:	_____
Lot Num:	_____
Exp. Date:	_____
Wafer Lot Num:	_____
Supt. Hemoglobin level in tubing is:	
	ACCEPTABLE
Operator Signature:	_____

Figure 56, Example of a procedure summary printout

Depicting the estimated free hemoglobin chart

The following example depicts the information found on the printed estimated free hemoglobin chart, which represents the levels of free hemoglobin released into the product waste bag during the dilution, wash and additive solution steps.


Figure 57, Example of an estimated free hemoglobin printout

Removing the final product

Once the procedure summary has been printed, the ACP 215 device will notify the operator with the following message that the the final product bag tubing is being prepared for sampling segments.

DE-GLYCEROLIZATION PROTOCOL					
MAKING SEGMENTS...					
	PUMP	TIME	HYPE	NACL	ADDI
REAL:	0	57	50	1470	240
TARG:		57	50	1470	240

Figure 58, Sampling segments message

After this message, the operator will receive the PROCEDURE COMPLETE message and should remove the final product as follows:

- ➔ Heat seal the tubing 10 cm [4 inches] from the final product bag.
- ➔ Detach the final product bag from the set and package it in accordance with the local SOP.
- ➔ Store the product at 2° C - 6° C in accordance with the local SOP.
- ➔ Heat seal the tubing at each marked segment for sampling segments.
- ➔ Remove the used disposable set and dispose of in accordance with the local SOP.

Under certain circumstances, sterile air may not be completely displaced through the disposable set, therefore some blood may remain in the bowl at the end of the procedure. In this case, the operator can manually transfer the blood in the bowl to the final product bag as follows:

- ➔ Remove the green-striped tubing from the green valve.
- ➔ Manipulate the waste bag as follows:
 - ➔ Grasp the waste bag on each side near the bottom corners.
 - ➔ Gently squeeze the bag to push any air pockets above the fluid towards the inlet port of the bag, ensuring that no waste fluid reaches the port and exits the bag.
 - ➔ Release the bag once the air is pocket is observed at the port.
- ➔ Remove the final product bag tubing from the white valve.
- ➔ Manually rotate the blood pump counterclockwise until the bowl is empty.

Once the blood has been transferred to the final product bag:

- ➔ Detach the final product bag from the set, package it and store at 2° C - 6° C in accordance with the local SOP.

Warning: Upon completion of any procedure, the used disposable set should be considered as biologically contaminated. It must be disposed of according to local standard operating procedure for the removal of such material and should not be mixed with non-biologically contaminated waste.

Once the DPM/SPM filters have been removed from the device, the device will automatically display the procedure selection menu and is prepared to initiate a new procedure.


Figure 59, ACP 215 procedure selection menu

APPENDIX A: DEFINING THE MODIFIABLE PARAMETERS

THAWED UNIT WEIGHT	The weight in grams of the thawed 40% W/V glycerolized RBC unit to be processed. Determined by weighing the freezer bag and its contents prior to the wash procedure, then taring the weight of the empty freezer bag. Any samples removed prior to the deglycerolization procedure should be excluded from this weight.
THAWED UNIT HEMATOCRIT	The measured hematocrit of the well-mixed unit of thawed 40% W/V glycerolized blood prior to the deglycerolization procedure should be constant at a $60\% \pm 5\%$ Hct based on the validation by the local facility.
MIN FLOW RATE	The minimum rate in ml/min in which the solution pump rotates during the procedure.
BOWL VOLUME	The volume in ml that is pumped into the bowl during the first dilution, calculated from the thawed RBC Unit Weight and Hct.
CENTRIFUGE SPEED	The number of revolutions per minute at which the centrifuge spins when operating at full speed.
THAWED UNIT OSMOLALITY	The osmolality of the 40% W/V glycerolized RBC unit to be processed.
HYPER SOLN VOLUME	The volume of the hypertonic (12% saline) solution that will be added to the thawed deglycerolized RBC unit.
HYPER SOLN FLOW RATE	The rate in ml/min in which hypertonic (12% saline) solution is pumped to the bowl.
HYPER EQUALIZE DELAY	The amount of time set for the volume of hypertonic (12% saline) solution to equilibrate with the thawed glycerolized RBC unit.
DILUTION TO WASH DELAY	The amount of time set for the volume of the wash (0.9% saline/0.2% dextrose) solution to equilibrate with the diluted thawed RBC unit.
OSMOLALITY RATE	The rate of change to decrease the osmolality of the thawed unit and achieve an osmolality of blood with the minimum amount of glycerol present.
NUM SALINE DILUTION	The number of dilution cycles performed.
NACL DILUTE VOL (#1)	The volume in ml of wash (0.9% saline/0.2% dextrose) solution used during the first dilution cycle.
NACL DILUTE VOL (#2)	The volume in ml of wash (0.9% saline/0.2% dextrose) solution used during the second dilution cycle.

CONSTANT DILUTION RATE	A parameter, that when set to Yes, will send solution at a constant rate during each dilution cycle, according to the default setting of the Dilution Rate parameter.
DILUTION RATE	The rate in ml/min at which solution is sent during a dilution cycle.
NACL RINSE VOLUME (#1)	The volume in ml of wash (0.9% saline/0.2% dextrose) solution used during the rinse cycle.
MAX BOWL FILL RATE	The maximum rate in ml/min at which the blood pump rotates during the procedure.
MIN BOWL FILL RATE	The minimum rate in ml/min at which the blood pump rotates during the procedure.
NUM WASH CYCLES	The number of wash cycles required to remove the residual glycerol from the blood being processed.
WASH SOLN VOL (#1-5)	The volume in ml of wash (0.9% saline/0.2% dextrose) solution that is used to wash cells between one cycle of starting and stopping the centrifuge.
WASH FLOW RATE (#1-5)	The rate in ml/min of wash (0.9% saline/0.2% dextrose) solution that is pumped through the bowl during each wash step.
EQUILIBRATION TIME	The time in seconds that the bowl and centrifuge are stopped in between wash steps.
RE-CENTRIFUGATION TIME	The time in seconds that the bowl spins before wash (0.9% saline/0.2% dextrose) solution is pumped through it.
NUM RESUSPENSIONS	The number of times that the bowl must stop to facilitate red cell suspension in the bowl.
PRESERV. SOLN VOL	The volume in ml of additive solution that is pumped through the bowl before the final batch of washed blood is pumped into the final product bag.
PRESERVATIVE SOLN RATE	The rate in ml/min of additive solution that is pumped through the bowl to the final washed product bag.
BOWL EMPTY RATE	The rate in ml/min at which the washed product is pumped from the bowl to the final product bag.

APPENDIX B: TROUBLESHOOTING DURING A PROCEDURE

Listing notice messages

During the deglycerolization procedure, the ACP® 215 device may interrupt operation and display particular messages, to indicate a situation that can easily be corrected by the operator. The following table list these notice messages along with a brief explanation of the situation and the recommended operator actions.

INSTALLATION NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- CHECK CORRECT DISPOSABLE LOADED - INSTALL DPM FILTER PRESS YES TO CONTINUE	The DPM filter has not been detected on the donor pressure monitor.	➔ Press the filter firmly onto the DPM while providing a quarter turn clockwise to ensure a tight seal. ➔ Press the YES key to resume operation.
- CHECK CORRECT DISPOSABLE LOADED - INSTALL SPM FILTER PRESS YES TO CONTINUE	The SPM filter has not been detected on the system pressure monitor.	➔ Press the filter firmly onto the SPM while providing a quarter turn clockwise to ensure a tight seal. ➔ Press the YES key to resume operation.
- CHECK CORRECT DISPOSABLE LOADED - INSTALL TUBING IN LINE SENSOR PRESS YES TO CONTINUE	The effluent tubing has not been detected as the line sensor.	➔ Insert the tubing completely in the groove of the line sensor and close the cover.
- BOWL SEATING TEST - BOWL NOT SEATED PROPERLY PRESS YES TO CONTINUE	The centrifuge has spun the bowl and the test has not passed.	➔ Verify that the bowl is completely seated in the centrifuge chuck. ➔ Reposition if necessary. ➔ Press the YES key to resume operation.
- IS FINAL PRODUCT BAG EMPTY? - Press YES to start prime Press NO to continue purge	Residual air is in the final product bag during the pre-Prime state.	The operator has two options: ➔ If the final product bag is empty, press the YES key to proceed to prime the disposable set. or ➔ If the final product bag still contains air, press the NO key to continue purging air from the bag.

PRIME and DILUTE STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- AIR AT SLAD - Check flow in hypertonic line Press START to continue	The SLAD did not detect fluid during the Prime state.	➔ Inspect the solution line for occlusions and ensure that the tubing path to the SLAD is clear. ➔ Ensure that the hypertonic solution bag is spiked. ➔ Press the START key to resume operation.
- AIR AT SLAD - Check flow in NaCl-Glucose line Press START to continue	The SLAD did not detect fluid during the Prime state.	➔ Inspect the solution line for occlusions and ensure that the tubing path to the SLAD is clear. ➔ Ensure that the wash solution bag is spiked. ➔ Press the START key to resume operation.
- AIR AT BLAD - Check flow in solution line Press START to continue	The BLAD detected air during the Prime state.	➔ Inspect the blood line for occlusions and ensure that the tubing path is clear. ➔ Press the START key to resume operation.
<i>Initial message</i> - INLET PRESSURE OUTSIDE LIMITS - Check for occlusion in line leading to thawed unit bag ***** <i>Corrected state</i> - PRESSURE RELIEVED - Press START to continue	The DPM detected pressure outside of limits during the Prime state.	➔ Inspect the inlet harness lines for occlusions. ➔ If there are no occlusions, manually open the red valve briefly. Once the pressure is relieved: ➔ Press the START key to resume operation.

PRIME and DILUTE STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- AIR AT BLAD - Check flow in solution line Press START to continue	The BLAD detected air before starting to fill the bowl.	➔ Inspect the solution line for occlusions and ensure that the tubing path is clear. ➔ Press the START key to resume operation. If the message is repeated: ➔ Inspect the disposable set for a leak. If a leak is detected: ➔ Discontinue the procedure and handle the unprocessed cells according to the local SOP.

FILL STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- AIR AT BLAD - Press START to continue Press NO to skip to next step	The BLAD detected air before a programmed bowl blood volume was reached.	The operator has two options: ➔ Press the START key to continue filling the bowl. or ➔ Press the NO key to advance to the next step of the procedure.
- FILL LONGER THAN EXPECTED - Press START to continue Press NO to skip to next step	Blood remains in the RBC bag after a programmed bowl blood volume was reached.	The operator has two options: ➔ If the bag is not empty, press the START key to continue filling the bowl. or ➔ If the bag is empty, press the NO key to advance to the next step of the procedure.

FILL STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
<i>Initial message</i> - BOWL PRESSURE OUTSIDE LIMITS - Check for effluent line occlusion ***** <i>Corrected state</i> - PRESSURE RELIEVED - Press START to continue	The SPM detected pressure outside of acceptable limits as the bowl was filling.	➔ Inspect the effluent harness lines for occlusions. Once the pressure is relieved: ➔ Press the START key to resume operation.
<i>Initial message</i> -INLET PRESSURE OUTSIDE LIMITS- Check for occlusion in line leading to thawed unit bag ***** <i>Corrected state</i> - PRESSURE RELIEVED - Press START to continue	The DPM detected pressure outside of acceptable limits as the bowl was filling.	➔ Inspect the inlet harness line for occlusions. ➔ Check the blood line tubing for occlusions. Once the pressure is relieved: ➔ Press the START key to resume operation.

DEGLYCEROLIZATION STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
- AIR AT SLAD - Check flow in solution line Press START to continue	The SLAD detected air in the wash solution line	➔ Inspect the solution line for occlusions and ensure that the tubing path is clear. ➔ Ensure that the wash solution bag is spiked correctly. ➔ Press the START key to resume operation.
- AIR AT SLAD - Check flow in preservative line Press START to continue	The SLAD detected air in the additive solution line.	➔ Inspect the additive solution line for occlusions and ensure that the tubing path is clear. ➔ Ensure that the additive solution bag is correctly connected. ➔ Press the START key to resume operation.

EMPTY STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
<i>Initial message</i> -INLET PRESSURE OUTSIDE LIMITS- Check for occlusion in line leading to thawed unit bag ***** <i>Corrected state</i> - PRESSURE RELIEVED - Press START to continue	The DPM detected pressure outside of acceptable limits as the bowl was emptying.	→ Inspect the inlet harness lines for occlusions. → If there are no occlusions, manually open the red valve briefly. Once the pressure is relieved: → Press the START key to resume operation.
- EMPTY LONGER THAN EXPECTED - Press START to continue Press NO to skip to next step	Blood remains in the bowl when it should be completely emptied.	The operator has two options: → If the bowl is not empty, press the START key to continue emptying the bowl. or → If the bowl is empty, press the NO key to advance to the next step of the procedure.
- AIR DETECTED EARLY - Press Start to continue Press NO to skip to next step	BLAD detected air early during the Empty state.	The operator has two options: → Press the START key to continue operation. or → Press the NO key to advance to the next step of the procedure.

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HAEMONETICS ACP 215	
SERIAL NUMBER: 95B072	
THAWED BLOOD DEGLYCEROLIZATION	
PROCEDURE SUMMARY	
Date and Time:	
Unit Number:	
Hypertonic Volume:	50 ml
Shaker:	ON
First Dilution Volume:	
Shaker:	ON
Bowl Fill:	800 ml
RBC Spillage:	NO
Second Dilution Volume:	400 ml
Shaker:	ON
Bowl Fill:	900 ml
RBC Spillage:	NO
1st Wash Volume:	
2nd Wash Volume:	70 ml
3rd Wash Volume:	100 ml
4th Wash Volume:	150 ml

5th Wash Volume: 300 ml

Additive Solution: 240 ml

Final Product: Estimated

Free Hemoglobin Level < : 50 mg%

Procedure Time: 56 Min.

Operator Initials A A

Thawed unit temperature:

12% NaCl Lot Num:

Exp. Date:

0.9% NaCl Lot Num:

Exp. Date:

Additive Solution:

Lot Num:

Exp. Date:

Wafer Lot Num:

Supt. Hemoglobin level in tubing is :
ACCEPTABLE

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"Approved"
Yulia PERMAN
INTERNATIONAL RA
Specialist
11. August - 2015

Haemonetics® ACP®215 Glycerolization Procedure – Operation Manual –

**HAEMONETICS CORPORATION
400 WOOD ROAD
BRAINTREE, MA 02184 USA
1-800-225-5242**

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Table of Contents

ACP 215 Glycerolization Procedure

PRESENTING THE ACP® 215 GLYCEROLIZATION PROCEDURE	1
Intended use	1
Listing procedure material	1
Understanding procedure guidelines	1
DESCRIBING THE GLYCEROLIZATION DISPOSABLE SET	3
Glycerol solution line	3
Bacteria barrier filter	3
Pump tubing	3
DPM filter	3
Heat sealed tubing	3
PREPARING THE ACP 215 DEVICE FOR A PROCEDURE	4
System testing	4
Selecting a procedure	4
INSTALLING THE DISPOSABLE SET	6
Inspecting the disposable material	6
Time limits for disposable sets and solutions	6
Installing the pump tubing	6
Installing the BLAD tubing and DPM filter	7
Connecting the solutions and RBC product	8
VERIFYING THE DISPOSABLE SET INSTALLATION	9
Initiating the disposable set checklist	9
Bacteria barrier filter	10
RBC bag tubing	10
Drip chamber vent open	10
Drip chamber spiked to bottle and primed	11
RBC bag on shaker	11
ENTERING OPERATOR INITIALS	12
EXPLAINING MODIFIABLE PROCEDURE PARAMETERS	13
DESCRIBING THE GLYCEROLIZATION PROCEDURE	15
Receiving the printed procedure summary	15
Depicting the procedure summary printout	16
Removing the final product	17
APPENDIX A: DEFINING THE MODIFIABLE PARAMETERS	18
APPENDIX B: TROUBLESHOOTING DURING A PROCEDURE	19
Listing notice messages	19

P/N 85278-30, Manual revision: DA

PRESENTING THE ACP® 215 GLYCEROLIZATION PROCEDURE

Intended use

The ACP 215 is intended to glycerolize and deglycerolize red blood cells derived from whole blood or collected by apheresis that have been stored in an approved anticoagulant/additive solution for up to 7 days at 1–6°C prior to glycerolization.

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Listing procedure material

The following material is needed to perform an ACP 215 glycerolization procedure:

- The ACP 215 device, the shaker and the printer.
- A sterile connection device.
- The Haemonetics® disposable set LN225.
- An 800-1000 ml RBC transfer bag qualified for the storage of frozen cells.
- Glycerol solution, 57%W/ V (weight to volume).

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The following guidelines should be respected concerning the RBC product to be glycerolized:

- The hematocrit of the initial RBC product to be glycerolized should be 75% ± 5%.

Note: To achieve the appropriate hematocrit, the cells should be centrifuged at 1615 x g for 4 minutes.

- If the site intends to use the LN235 disposable set (275 ml), for deglycerolization, the total red cell mass of the product should not exceed 180 mls. The LN236 disposable (325 ml bowl) should be used when the red cell mass is greater than 180 mls.

The red cell mass is based on the donor's hematocrit and the weight of the whole blood drawn from the donor. For example, if the donor's hematocrit is 44% the site should not collect more than 430 mls of whole blood. Refer to Table 2, "Red cell mass (donor hematocrit known)" in the Deglycerolization Procedure Manual for further information on red cell mass.

- The initial unit of RBC product should be glycerolized and frozen within six to seven days of collection, depending on the local regulations and guidelines.
- The temperature of the glycerol and RBC unit to be glycerolized should be brought to 20° C - 30° C prior to the procedure.

Note: Haemonetics recommends placing a secondary, watertight cover around the RBC unit when warmed in a water bath, to prevent the ports from becoming contaminated.

- The RBC product must be transferred to a qualified 800-1000 ml storage bag, prior to the glycerolization procedure. The quantity of glycerol solution to be added could not be contained in the RBC collection bag.

DESCRIBING THE GLYCEROLIZATION DISPOSABLE SET

1. Glycerol solution spike
2. Ratchet clamp
3. Blood pump tubing
4. Bacteria barrier filter
5. DPM filter
6. Slide clamp
7. Heat sealed tubing


Figure 1, LN225 Glycerolization disposable set

Glycerol solution line

The glycerol solution line contains a spike which is used to connect the glycerol solution bottle to the disposable set and a ratchet clamp.

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Pump tubing

The blood pump tubing sections contain **two blue** pump stops which will secure the pump tubing during pump loading and operation.

DPM filter

A 0.2 μm filter located on the blood line is used as a bacteria barrier between the draw pressure monitor and the disposable set. It is coupled with the receptor on the DPM to monitor pressure within the blood line from the blood pump to the RBC bag. In case of any occlusion in the line or any obstruction generated by the sterile connection process, the pressure monitor will stop the pump and alert the operator.

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ACP 215 SYSTEMS TESTS			
Rom	P		Air Sensors P
Calibration	P		Analog P
Config	P		Valves P
Protocols	P		Watch Dog P

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HAEMONETICS ACP 215	
PLEASE SELECT PROTOCOL	
➔➔➔➔	CELL WASH
	GLYCEROLIZATION
	DE-GLYCEROLIZATION

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- ➔ Press the MODIFY PROGRAM key to scroll the list until the arrow indicates GLYCEROLIZATION.

HAEMONETICS ACP 215	
PLEASE SELECT PROTOCOL	
	CELL WASH
•==>	GLYCEROLIZATION
	DE-GLYCEROLIZATION

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Once the glycerolization procedure has been selected, the following message will be displayed, along with target values based on the procedure programming in the row **TARG**.

GLYCEROLIZATION PROTOCOL				
PLEASE LOAD DISPOSABLE SET PRESS YES TO CONTINUE				
	PUMP	TIME	GLYSOL	BLOOD
REAL:	0	00:00	0	0
TARG:	70	10:40	442	680

Figure 5, Disposable set installation message

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- PUMP = pump speed.
- TIME = procedure duration in minutes.
- GLYSOL = glycerol solution used in ml.
- BLOOD = blood volume processed in ml.

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Note: This message will remain displayed during the disposable set installation sequence. The operator should install the disposable set completely, prior to pressing the YES key. Pressing the YES key at this point will display the disposable set checklist.

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- Verify that the disposable set corresponds with the selected procedure.
- Verify that the Tyvek® pouch has not been damaged.
- Open the Tyvek pouch and inspect all tubing sections.
- Ensure that no kinks or occlusions are present which could obstruct the flow through the disposable set.
- Verify that there are no visible defects or particulate within the set.

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- Remove the harness from the Tyvek pouch.
- Close the ratchet clamp on the glycerol solution line.
- Loop the pump tubing segment between the blood pump stops around the groove of the blood pump rotor.
- Secure the tubing in each tubing guide so that each pump stop is located to the **right** of the tubing guide.
- Rotate the pump clockwise until the notch in the pump rotor catches the tubing and draws it completely into the pump groove.

- ➔ Rotate the pump clockwise another quarter of a turn to ensure proper placement of the tubing in the pump groove.
The pump should be able to rotate freely.

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1. Blood pump
2. Pump tubing guides
3. Blood tubing pump stops
4. BLAD
5. Bacteria barrier filter
6. DPM filter


Figure 6, Tubing installation and bacteria barrier filter position

Note: The bacteria barrier filter must be positioned as depicted.

Installing the BLAD tubing and DPM filter

- ➔ Thread the clear tubing between the blood pump and the Y-connector attached to the DPM filter through the blood line air detector (BLAD).
- ➔ Verify that the pump stops are situated against the tubing guides once the tubing has been installed in the BLAD.
- ➔ Connect the DPM filter to the draw pressure monitor.
- ➔ Press the filter firmly onto the DPM while providing a quarter turn clockwise to ensure a tight seal.

Warning: If for any reason, the DPM filter becomes wet during the procedure, the device may not be able to correctly detect pressure in the blood line.

Once the disposable set has been completely installed, the operator can proceed to the connect the glycerol solution and the blood product.

Connecting the solutions and RBC product

Using **aseptic technique**, the operator should:

- ➔ **Spike the glycerol solution bottle onto the solution line and suspend the bottle from the solution pole.**
- ➔ Ensure that the air vent on the spike is open, to allow air to flow into the bottle as it empties.
- ➔ Perform a sterile connection weld between the blood line tubing and the tubing of the RBC bag.

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- ➔ Verify that the sterile connection weld is not occluded by squeezing the weld area between the fingers and rolling, to fully open the weld.
- ➔ Secure the RBC bag on the shaker using the magnets, as depicted.

Figure 7, Securing the RBC bag on the shaker

- ➔ Verify again that the sterile connection weld is fully open and that no kinks are evident in the tubing.

Warning: An improperly installed disposable set may result in product damage. Any leak in or potential opening in the set reduces the product life to 24 hours.

At this point, the operator can proceed to the disposable set checklist and verify the installation of the disposable set.

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Figure 8, Disposable set installation message

- ➔ Press the YES key to initiate the checklist.

A Full or Minimum checklist will be displayed, depending on the technical configuration of ACP 215 device.

The following table indicates with the letter **X** which messages will be displayed for each disposable set checklist configuration.

Table 1, Disposable set checklist configurations: Glycerolization procedure

Message	Full checklist	Minimum checklist
Bacteria barrier filter	X	
RBC bag tubing	X	
Drip chamber vent open	X	
Drip chamber spiked to bottle/ primed	X	X
RBC bag on shaker	X	X

Bacteria barrier filter

- Check, Press YES To Continue				
Bacteria filter between BLOOD <FRONT> pump and RBC Bag.				
	PUMP	TIME	GLYSOL	BLOOD
REAL:	0	00:00	0	240
TARG:	70	10:40	442	680

Figure 9, Bacteria filter position message

- ➔ Verify that the bacteria barrier filter is correctly positioned.
- ➔ Correct the position if necessary.
- ➔ Press the YES key to continue.

RBC bag tubing

- Check, Press YES To Continue				
Sterile dock RBC Bag.				
	PUMP	TIME	GLYSOL	BLOOD
REAL:	0	00:00	0	240
TARG:	70	10:40	442	680

Figure 10, Sterile connection message

- ➔ Verify that the RBC product bag is sterile connected to the disposable set.
- ➔ Perform a sterile connection if not present.
- ➔ Ensure that the sterile connection weld is fully open.
- ➔ Press the YES key to continue.

Drip chamber vent open

- Check, Press YES To Continue				
RED TAB for vent above drip chamber is open.				
	PUMP	TIME	GLYSOL	BLOOD
REAL:	0	00:00	0	240
TARG:	70	10:40	442	680

Figure 11, Vent tab message

- ➔ Verify that the vent tab on the glycerol solution spike drip chamber is open.
- ➔ Open it if it is closed.
- ➔ Press the YES key to continue.

Drip chamber spiked to bottle and primed

- Check, Press YES To Continue				
Spike connected and prime drip chamber with Glycerol.				
	PUMP	TIME	GLYSOL	BLOOD
REAL:	0	00:00	0	240
TARG:	70	10:40	442	680

Figure 12, Drip chamber spiked to bottle/prime message

- Verify that the glycerol solution bottle is spiked and the drip chamber is primed.
- Spike the bottle and prime the drip chamber if necessary.
- Press the YES key to continue.

Note: The operator should also verify that the target amount of glycerol solution displayed (i.e. the amount required to complete the procedure) is less than the amount of solution available in the glycerol solution bottle.

RBC bag on shaker

- Check, Press YES To Continue				
Secure RBC Bag to the shaker.				
	PUMP	TIME	GLYSOL	BLOOD
REAL:	0	00:00	0	240
TARG:	70	10:40	442	680

Figure 13, RBC bag/shaker message

- Verify that the RBC bag is correctly installed on the shaker.
- Reposition if incorrectly installed.
- Press the YES key to continue.

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- ➔ Press the YES key to scroll the alphabet, letter by letter from A to Z.
- ➔ Press the SAVE PROGRAM key when the appropriate letter is displayed.

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- ➔ Press the YES key to continue.

EXPLAINING MODIFIABLE PROCEDURE PARAMETERS

Each ACP 215 device is delivered with pre-installed procedure parameter settings called the *default* settings. Certain parameters, referred to as the *modifiable parameters* are accessible to the operator. The operator can either use the default settings, or adjust them for the present procedure.

Warning: Any modification to the Haemonetics default values remains the responsibility of the local site and should be in accordance with the local SOP. The device will retain in memory any settings which the operator decides to save, and will not automatically revert to the Haemonetics® default values.

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Figure 16, Modifiable parameters message

The ACP 215 device retains in memory any modified settings which were saved from the previously performed procedures. If the operator decides to use these settings for the current procedure:

- Press the START key to initiate the glycerolization procedure.

However, the operator may decide to access and/or modify certain parameters according to the needs of the current procedure. Depending on the technical configuration of the ACP 215 device, either a **minimum set** or a **full set** of modifiable parameters will be available to the operator.

To access the modifiable parameters for the glycerolization procedure:

- Press the MODIFY PROGRAM key.
- Verify the current parameter settings for the procedure.
- Use the Arrow keys to change the value as necessary.
- Press the MODIFY PROGRAM key to access the next parameter available.

To retain the modified settings for the subsequent procedure(s):

- Press the SAVE PROGRAM key.

To initiate the glycerolization procedure at this point:

- Press the START key.

The following table lists the full set of modifiable parameters for the glycerolization procedure. The parameter in **bold type** will be displayed when the device is configured for the **minimum set**.

A list of definitions for each parameter is provided in *Appendix A*.

Table 2, ACP 215 modifiable parameters for the Glycerolization procedure

Parameter	Unit of measure	Default setting	Range of settings	Increment of change
Unit Weight	g	225	135-555	5
RBC Unit Hematocrit	%	75	20-100	5
RBC Unit Osmolality	mOsm/kgH ₂ O	320	10-600	10
Glycerol Concentration	g%	57	20-60	1
Glycerol Conc. Target	g%	40	10-50	1
Glycerol Osmolality * *minimum value = (RBC Unit Osmolality + 100)	mOsm/kgH ₂ O	7230	200-9000	10
Target Osmolality Rate	mOsm/kgH ₂ O/min	500	50-900	50
Max Glycerol Flow Rate	ml/min	70	10-70	10
Constant Flow	—	NO	YES-NO	—
Glycerol Flow Rate	ml/min	60	10-70	10

DESCRIBING THE GLYCEROLIZATION PROCEDURE

Once the operator has pressed the START key, the ACP 215 device will automatically perform the glycerolization procedure according to the programmed parameters.

First, the disposable set will be primed with glycerol solution. Once the set has been primed, the volume of glycerol used for the procedure will be counted and displayed to the operator with the procedure statistics. The glycerol solution flow rate will progressively increase according to the procedure parameters until the procedure is complete.

Note: The operator should verify the volume in the glycerol solution bottle during the procedure and compare the consumed volume to the displayed volume.

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Once the RBC product has been completely glycerolized, the pumps will stop, the device will beep and the following message will be displayed. The shaker will **stop after** thirty more seconds.

GLYCEROLIZATION PROTOCOL				
IS PRINTER ONLINE PRESS YES TO CONTINUE				
	PUMP	TIME	GLYSOL	BLOOD
REAL:	0	10:40	442	680
TARG:	70	10:40	442	680

Figure 17, Printer verification message

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At the end of the printing operation, the device will display the following message, prompting the operator to respond about the quality of the printout.

IS THE PRINTOUT ACCEPTABLE?				
PRESS YES TO CONTINUE PRESS NO TO RETRY				
	PUMP	TIME	GLYSOL	BLOOD
REAL:	0	10:40	442	680
TARG:	70	10:40	442	680

Figure 18, Printout confirmation message

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Figure 19, Printer verification message

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The following example depicts the information found on the glycerolization procedure summary printout. For certain items, the values are generated automatically by the ACP 215 device. For other items, a blank space is provided for the operator to write the information. The name of the manufacturer should be included when documenting the lot number of each item used for the procedure.

HAEMONETICS ACP 215	
SERIAL NUMBER: 95B072	
GLYCEROLIZATION PROTOCOL	
PROCEDURE SUMMARY	
Date & Time:	_____
Blood Unit Number:	_____
Glycerol Volume:	442 ml
Shaker:	ON
Procedure Time:	12 Min.
Operator Initials:	A A
Disposable Lot Num:	_____
Glycerol temperature:	_____
Lot Num:	_____
Exp. Date:	_____
Wafer Lot Num:	_____
Initial Unit Net Weight:	_____
Operator Signature:	_____

Figure 20, Example of a procedure summary printout

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- ➔ Heat seal the tubing 10 cm [4 inches] from the bag.
- ➔ Detach the final product bag from the set and package it in accordance with the local SOP.
- ➔ Remove the supernatant glycerol from the product.
- ➔ Store and freeze the product accordance with the local SOP.

Warning: After the glycerolization process has been completed, the supernatant glycerol must be removed from the product prior to freezing, using a centrifugation technique in accordance with the local SOP, to result in a $60\% \pm 5\%$ hematocrit. This is accomplished by centrifuging at $1248 \times g$ for 10 minutes.

Improper RBC glycerolization, freezing, or storage may damage the red blood cells and produce high free hemoglobin content.

- ➔ Remove the used disposable set and dispose of in accordance with the local SOP.

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Once the DPM filter has been removed from the device, the device will automatically display the procedure selection menu and is prepared to initiate a new procedure.

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Figure 21, ACP 215 procedure selection menu

APPENDIX A: DEFINING THE MODIFIABLE PARAMETERS

UNIT WEIGHT	The weight in grams of the RBC unit to be processed. Determined by weighing the storage bag and its contents prior to the glycerolization procedure, then taring the weight of the empty storage bag. Any samples removed prior to the procedure should be excluded from this weight.
RBC UNIT HEMATOCRIT	The measured hematocrit of the concentrated RBC product prior to the glycerolization procedure.
RBC UNIT OSMOLALITY	The osmolality of the RBC unit to be processed.
GLYCEROL CONCENTRATION	The concentration of glycerol used to glycerolized the RBC unit, expressed in weight to volume percentage.
GLYCEROL CONC. TARGET	The target of the final glycerol concentration, expressed in weight to volume percentage.
GLYCEROL OSMOLALITY* *minimum value = (RBC UNIT OSMILARITY + 100)	The osmolality of the glycerol solution, expressed in mOsm/Kg of water.
TARGET OSMOLALITY RATE	The rate at which the glycerol solution is added to the RBC unit, expressed in mOsm/Kg of water per minute.
MAX GLYCEROL FLOW RATE	The maximum rate in ml/min. at which the blood pump rotates during the procedure.
CONSTANT FLOW	A parameter, that when set to YES, will send glycerol solution at a constant rate according to the default setting of the Glycerol Flow Rate parameter.
GLYCEROL FLOW RATE	The rate in ml/min. at which solution is sent during a cycle.

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Listing notice messages

During the glycerolization procedure, the ACP® 215 device may interrupt operation and display particular messages, to indicate a situation that can easily be corrected by the operator. The following table list these notice messages along with a brief explanation of the situation and the recommended operator actions.

INSTALLATION NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
CHECK CORRECT DISPOSABLE LOADED INSTALL DPM FILTER PRESS YES TO CONTINUE	The DPM filter has not been detected on the donor pressure monitor.	➔ Press the filter firmly onto the DPM while providing a quarter turn clockwise to ensure a tight seal. ➔ Press the YES key to resume operation.
CHECK CORRECT DISPOSABLE LOADED REMOVE SPM FILTER PRESS YES TO CONTINUE	An SPM filter has been detected on the system pressure monitor.	➔ Verify the disposable set which has been installed. ➔ Remove the SPM filter and the rest of the disposable set, if it is not a glycerolization set. ➔ Install a glycerolization set. ➔ Press the YES key once an appropriate set has been installed, to receive the disposable set checklist.
CHECK CORRECT DISPOSABLE LOADED REMOVE TUBING FROM LINE SENSOR PRESS YES TO CONTINUE	Tubing has been detected in the line sensor.	➔ Verify the disposable set which has been installed. ➔ Remove the tubing and the rest of the disposable set, if it is not a glycerolization set. ➔ Install a glycerolization set. ➔ Press the YES key once an appropriate set has been installed, to receive the disposable set checklist.
PROCEDURE TIME > 30 MIN PLEASE MODIFY GLY PARAMETERS PRESS MODIFY TO CONTINUE	A parameter value is incorrectly set and the procedure time will require 30 or more minutes to complete.	➔ Press the MODIFY key to access and verify the settings. ➔ Correct any values as necessary and press the MODIFY key to continue.

INSTALLATION NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
GLYCEROL VOLUME > 1000 ML PLEASE MODIFY GLY PARAMETERS PRESS MODIFY TO CONTINUE	A parameter value is incorrectly set and the glycerol volume required for the procedure will be greater than or equal to 1000 ml.	➔ Press the MODIFY key to access and verify the settings. ➔ Correct any values as necessary and press the MODIFY key to continue.

PRIME STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
GLYCEROLIZATION PROTOCOL CHECK FOR LINE OCCLUSION PRESS START TO CONTINUE	The BLAD has not detected fluid during the Prime state.	➔ Inspect the solution line for occlusions and ensure that the tubing path to the BLAD is clear. ➔ Press the START key to resume operation.
<i>Initial message</i> GLYCEROLIZATION PROTOCOL HIGH PRESSURE DETECTED CHECK FOR LINE OCCLUSION ***** <i>Corrected state</i> - PRESSURE RELIEVED - PRESS START TO CONTINUE	The DPM detected pressure outside of limits during the Prime state	➔ Inspect the harness for occlusions and ensure that the tubing path is clear. Once the pressure is relieved: ➔ Press the START key to resume operation.

GLYCEROLIZATION STATE NOTICE MESSAGE	SITUATION	CORRECTIVE ACTION
<i>Initial message</i> GLYCEROLIZATION PROTOCOL HIGH PRESSURE IN LINE CHECK LINE TO RELIEVE PRESSURE ***** <i>Corrected state</i> PRESSURE RELIEVED PRESS START TO CONTINUE	The DPM detected pressure outside of limits during the Glycerolization state.	The operator has two options: → Relieve the pressure and continue the procedure. or → Press the STOP key to end the procedure. Once the pressure is relieved: → Press the START key to resume operation.

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CONSUMER INFORMATION

Proprietary rights

The contents of this manual are property of the Haemonetics Corporation. Haemonetics® and ACP® are trademarks or registered trademarks of Haemonetics Corporation in the United States, other countries, or both.

The following are trademarks or registered trademarks of their respective owners in the United States, other countries, or both: Seiko — Seiko Instruments USA Inc.; Tyvek — E.I. du Pont de Nemours and Company; Terumo — Terumo Medical Corporation.

Any information or descriptions contained in this manual may not be reproduced and released to any of the general public, or used in conjunction with any professional instruction without written consent of Haemonetics Corporation, USA. Please direct any written inquiries to the appropriate address:

International Headquarters	Corporate Headquarters
Haemonetics SA	Haemonetics Corporation
Signy Centre, rue des Fléchères 6	400 Wood Road
P.O. Box 262, 1274 Signy-Centre	Braintree, MA 02184
Switzerland	USA
Tel. [+41 22] 363 90 11	Tel. [+1 781] 848 7100
Fax [+41 22] 363 90 54	Fax [+1 781] 848 5106

Legal Disclaimer

This manual is intended for use as a guide for material supplied by the Haemonetics Corporation. It provides the operator with necessary information to carry out specific procedures and maintain Haemonetics produced equipment. Use this manual in conjunction with instruction and training supplied by qualified Haemonetics personnel.

Haemonetics provides a limited product warranty when its products are used in accordance with published instructions by a properly trained operator. Any failure to follow the procedures described in this manual could result in impaired equipment function, injury to the operator or patient, or the voiding of applicable product warranties. Haemonetics accepts no responsibility for problems resulting from failure to comply with prescriptions as outlined by the company.

Any modifications estimated as necessary by the customer should be evaluated by a Haemonetics clinical specialist (or the Haemonetics Customer Care Center).

Safe utilization of Haemonetics material and equipment requires the operator to correctly handle and dispose of blood-contaminated material. The operator of any Haemonetics equipment must fully understand and implement the local prevailing policies and procedures of each facility in which Haemonetics products are used concerning blood-contaminated material as well as blood products.

It remains solely the responsibility of the customer to fully assess and ensure the safety of any blood products obtained from Haemonetics prescribed procedures, prior to further application or use. Haemonetics declines any responsibility for choices made by the consumer concerning the utilization of these products and by-products.

P/N 85271-30, Manual revision: EB

Haemonetics worldwide contact information

Corporate headquarters

Haemonetics Corporation
400 Wood Road
Braintree, MA 02184
U.S.A.
Tel.: +1 781 848 7100
Fax: +1 781 848 5106

Haemonetics Austria

Handelsges.m.b.H.
Mariahilfer Strasse 123
3rd Floor
1060 Wien, Austria
Tel.: 0800 29 2777
Fax: 0800 29 2820

Haemonetics Belgium-NV

Braine-L'Alleud Parc de L'Alliance 9
Boulevard de France, bât A
1420 Braine L'Alleud
Belgium
Tel.: (FR): 0800 754 80
Fax: 0800 75512
Tel. (NL): 0800 754 82
Fax: 0800 755 12

Haemonetics BV

Claudius Prinsen Hof
Verlengde Poolseweg 34-46
4818 CL Breda
The Netherlands
Tel.: 0800 0222 707
Fax: 0800 0223 066

Haemonetics CZ, spol. s.r.o.

Ptašínského 8
60200 Brno, Czech Republic
Tel.: 800 143 243
Fax: 800 143 250

Haemonetics (Hong Kong) Ltd.

Rooms 2202, 22/F
Harbour Centre, 25 Harbour Road
Wanchai, Hong Kong
Tel.: +852 2868 9218
Fax: +852 2801 4380

International headquarters

Haemonetics S.A.
Signy Centre
Rue des Fléchères 6
P.O. Box 262
1274 Signy Centre, Switzerland
Tel: +41 22 363 9011
Fax: +41 22 363 9054

Haemonetics France SARL

46 bis, rue Pierre Curie
Z.I. Les Gâtines
78370, Plaisir, France
Tel.: 0800 90 11 58
Fax: 0800 91 48 76

Haemonetics GmbH

Wolfratshauser Straße 84
81379 Munich, Germany
Tel.: 0800 180 8890
Fax: 0800 182 8064

Haemonetics Hospitaler LTDA

Avenida Bernardino de Campos, 98
Sobrelaja, Paraisio, CEP 04004-040
City of São Paulo
State of São Paulo, Brazil

Haemonetics India

JMD Regent Square
Regus - Level 6
Mehrauli Gurgaon Road
Gurgaon - 122 001
Haryana (Delhi-NCR) India
Tel.: +91 124 471 1819
Fax: +91 981 027 3237

Haemonetics Italia SA

Via Senigalia 18/2 - Torre A
Milan 20161
Italy
Tel.: 800 870 200
Fax: 800 870 375

Haemonetics Japan GK

Kyodo Building, 16
Ichiban-cho, Chiyoda-ku
Tokyo 102-0082, Japan
Tel.: +81 3 3237 7260
Fax: +81 3 3237 7330

Haemonetics Korea

30th Floor ASEM Tower
159-1 Samsung-dong, Kangnam-ku
Seoul 135-798, Korea
Tel.: +82 2 6001 3280
Fax: +82 2 6001 3281

**Haemonetics Medical Devices
(Shanghai) International Trading
Company**

Room 1103-06 Evergo Mansion
No. 1325 Middle Huaihai Road
200031 Shanghai, China
Tel.: +86 21 3406 0700
Fax: +86 21 5466 8852

Haemonetics Scandinavia AB

Ideon
223 70 Lund, Sweden
Tel.: 020 797 150
Fax.: 020 794 670

Haemonetics Asia

Incorporated Taiwan Branch

26F-1, No. 102, Roosevelt Road Sec. 2
Taipei, Taiwan
Tel.: +886 2 2369 0722
Fax: +886 2 2364 3698

Haemonetics Ltd.

5 Ashley Drive
Bothwell
Glasgow G71 8BS
United Kingdom
Tel.: +44 113 273 7711
Fax: +44 113 273 4055

Lebanon Representative Office

Haemonetics SA
c/o The Regus Group
Azarieh Bldg., Block 3,
5th floor
Regus, Beirut, Lebanon
PO Box 11-503

Russia Representative Office

Haemonetics SA
Smolenskaya Passage
Smolenskaya Square 3
121099 Moscow Russia
Tel.: +7 495 937 8239

Table of Contents

Chapter 1, Explaining General Information

PROVIDING AN OVERVIEW	1-2
What is the ACP® 215 system?	1-2
What is the purpose of this manual?	1-2
What is required to perform a procedure?	1-2
UNDERSTANDING THE USE OF SYMBOLS	1-3
Symbols found in this document	1-3
Symbols found on the device	1-3
LISTING ACP 215 DEVICE SPECIFICATIONS	1-5
Dimensions and weight	1-5
Environmental conditions	1-5
Power requirements	1-6

Chapter 2, Presenting the ACP® 215 Equipment

DESCRIBING THE ACP 215 CONTROL PANEL	2-3
Display screen	2-3
Keypad	2-3
DESCRIBING THE ACP 215 CABINET COMPONENTS	2-4
Centrifuge system	2-4
Pumps	2-5
Pressure monitors	2-5
Valves	2-6
Air detectors	2-6
Optical line sensor	2-6
Solution bag pole	2-7
Power entry module	2-7
Power cord	2-7
Biohazard waste bag	2-7
DESCRIBING THE SHAKER AND THE PRINTER	2-8
Shaker	2-8
Printer	2-9

Chapter 3, Maintaining the ACP® 215 Equipment

CLEANING PROCEDURES	3-2
Listing tools and cleaning supplies	3-2
Exterior surfaces, valves, shaker and printer	3-3
Pressure monitors	3-3

P/N 85271-30, Manual revision EB

Air detectors. 3-3

Optical line sensor. 3-4

Pumps 3-4

Fluid detector. 3-5

Centrifuge components 3-5

Filter screens 3-6

CUSTOMER SERVICE 3-7

 Clinical training 3-7

 Field service. 3-7

 Product Return Guidelines 3-7

Chapter 4, Ensuring Safety and Quality for an ACP®215 Procedure

HANDLING THE ACP 215 EQUIPMENT 4-2

 Storing the device and material 4-2

 Inspecting the material. 4-2

 Installation 4-2

PREVENTING PROBLEMS DURING AN ACP 215 PROCEDURE 4-3

 Understanding the risk of hemolysis. 4-3

 Avoiding the consequences of flow restriction 4-3

 Avoiding bowl misalignment 4-4

 Avoiding overheating due to mechanical situations 4-4

WARNINGS FOR THE OPERATOR 4-5

 Electrical shock hazards. 4-5

 Leakage current control 4-5

 Mechanical hazards/rotating parts 4-5

 Power outlet connection 4-5

 Communicable disease precautions 4-6

Chapter 5, Troubleshooting for the ACP® 215 Device

LISTING SYSTEM ERROR MESSAGES 5-2

LISTING INTERRUPT SAFETY SYSTEM MESSAGES 5-5

Chapter 6, Reference information

APPENDIX A: IEC 60601-1-2:2001 STANDARD REQUIREMENTS 6-2

 Operation precautions 6-2

 Electromagnetic immunity 6-2

Chapter 1

Explaining General Information

Providing an overview	1-2
What is the ACP® 215 system?	1-2
What is the purpose of this manual?	1-2
What is required to perform a procedure?	1-2
Understanding the use of symbols	1-3
Symbols found in this document	1-3
Symbols found on the device	1-3
Symbols found on disposable packaging	1-4
Listing ACP 215 device specifications	1-7
Dimensions and weight	1-7
Environmental conditions	1-7
Power requirements	1-7

P/N 85271-30, Manual revision EB

PROVIDING AN OVERVIEW

What is the ACP[®] 215 system?

In the tradition of its innovative technology, Haemonetics has produced the ACP 215 - a compact, lightweight, automated cell processing system which is simple to operate.

The ACP 215 system can perform the following types of procedures:

- Prepare collected red blood cell products with glycerol for freezing.
- Deglycerolize thawed red blood cell products, removing extra cellular components, and resuspend the red blood cells in an additive solution.
- Wash red blood cell units stored in commonly used AC and additive solutions to remove plasma and extra cellular components.

The ability to process-control the addition and removal of cryoprotectant solution as well as wash banked blood in a closed system significantly improves the safety and availability of blood products.

What is the purpose of this manual?

This manual is intended to supply anyone involved in using the ACP 215 equipment with information to attain:

- An awareness of the purpose of the device and the implications of its procedures.
- An understanding of how to safely operate the system, correctly install the appropriate disposable material, and troubleshoot any difficulties.
- An ability to consistently apply the principles behind safe operation, proper maintenance and correct handling, to ensure optimal results.

What is required to perform a procedure?

The following material is required to perform an ACP 215 procedure:

- The ACP 215 device, to process the collected red blood cell products.
- A dedicated shaker, to properly mix the blood with procedure-specific solutions.
- A disposable set designed for the selected procedure.
- Appropriate solutions.
- A printer, to provide the operator with printed procedure data.
- A sterile connection device, to ensure proper tubing welds.

UNDERSTANDING THE USE OF SYMBOLS

Symbols found in this document

The terms *Note*, *Caution* and *Warning* are used in this manual with the following symbols to emphasize certain details for the operator.

Note: provides useful information regarding a procedure or operating technique when using Haemonetics material.

Caution: advises the operator against initiating an action or creating a situation which could result in damage to equipment, or impair the quality of the by-products; personal injury is unlikely.

Warning: advises the operator against initiating an action or creating a situation which could result in serious personal injury to either the operator or the recipient of the product.

- Text preceded by this bullet indicates an item on a list of information for the operator.
- ➔ Text preceded by this bullet indicates an action for the operator.

Symbols found on the device

The following symbols are found on the ACP 215 device.

- IEC Standard 60601-1, Medical Electrical Equipment, Part 1: General requirements for safety.
- IEC Standard 60417-1, Graphical symbols for use on equipment, Part 1: Overview and application.

Type B

Protection against electric shock, particularly regarding: allowable leakage current and reliability of the protective earth connection.

Protective earth (ground)

Used to identify any terminal intended for connection to an external conductor, for protection against electrical shock in case of a fault.

Alternating current

Used to indicate on the rating plate that the device is suitable for alternating current only.

Fuse symbol
Used to identify fuse boxes or the location of a fuse box.

Power OFF
Position of the main power switch indicating disconnection from the mains.

Power ON
Position of the main power switch indicating connection to the mains.

IPX1

Protection against ingress of liquid
Indicates that the enclosure of the device is designed to provide a specified degree of protection against harmful ingress of water or liquid into the equipment (under applicable conditions).

Attention (Consult accompanying documents)

Pollution control mark
Pollution control mark for products containing any of the six referenced substances (Lead, Mercury, Cadmium, etc...) according to new Chinese regulations.

Electrical and Electronic Equipment Waste
Dispose of the device using a separate collection method (according to EU and local regulation for waste electrical and electronic equipment).

Fragile, handle with care

Storage conditions, KEEP DRY

Storage conditions, temperature level

LISTING ACP 215 DEVICE SPECIFICATIONS

Dimensions and weight

The approximate weight and dimensions of the ACP 215 device are as follows:

Characteristic	Value (approximate)	
	Cabinet cover open	Cabinet cover closed
Height	67.5 cm [26.5 inches]	43 cm [17 inches]
Width	55 cm [21.5 inches]	
Depth	55 cm [21.5 inches]	30.5 cm [12 inches]
Weight	25 kg [56 lb]	

Environmental conditions

The following environmental conditions should be respected pertaining to operation and storage of the ACP 215 device.

Condition	Value
Operating Temperature	18° C to 35° C
Humidity	95% maximum, non-condensing
Storage temperature range	-20° C to 50° C<85% RH

Warning: Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

Note: Store disposables in a dry place away from solvent vapors and extreme temperatures. Storing disposable sets at extreme temperatures or humidity may cause defects. Inspect all disposable sets and disposable set packaging carefully before use.

Power require-
ments

The power is supplied to the device from an external AC power source. The power source must be properly grounded. Haemonetics will regulate the proper voltage setting upon installation. The power is supplied to the machine from an external AC source.

Electrical input power			
Rated voltage	Rated current	Fuse	Frequency
100–120 V~	2.6 A	F5 A, 250 V	50-60 Hz
230 V~±10%	1.9 A	F2.5 A, 250 V	50-60 Hz

Warning: For continued protection against fire hazard, replace fuse with the same type and rating.

Enclosure/chassis leakage current specifications*			
Condition	Polarity	Ground	Max value
Normal	Normal	Normal	100 µA
	Reverse	Normal	100 µA
Single fault	Reverse	Open	500 µA
	Normal	Open	500 µA

Note: The power source used must be properly grounded.

Caution: The ACP 215 device must be operated in an environment compatible to the requirements of the IEC 60601-1-2:2001 Standard, Electromagnetic compatibility.

Mobile RF communication equipment not approved by Haemonetics and portable communication equipment can affect the ACP 215 device. Any accessories and cables not approved by Haemonetics used in conjunction with the device may increase hazards and influence compatibility with EMC requirements. Therefore, non-approved accessories and cables must not be used.

In addition, the ACP 215 device and accessories must not be placed directly adjacent to, or top of other equipment, unless specifically approved by Haemonetics.

Chapter 2

Presenting the ACP® 215 Equipment

DESCRIBING THE ACP 215 CONTROL PANEL	2-3
Display screen	2-3
Keypad	2-3
DESCRIBING THE ACP 215 CABINET COMPONENTS	2-4
Centrifuge system	2-4
Pumps	2-5
Pressure monitors	2-5
Valves	2-6
Air detectors	2-6
Optical line sensor	2-6
Solution bag pole	2-6
Power entry module	2-7
Power cord	2-7
Biohazard waste bag	2-7
DESCRIBING THE SHAKER AND THE PRINTER	2-8
Shaker	2-8
Printer	2-9

P/N 85271-30, Manual revision EB


Figure 2-1, The ACP 215 device with the shaker

P/N 85271-30, Manual revision EB

DESCRIBING THE ACP 215 CONTROL PANEL

The control panel is located on the inside of the ACP 215 cabinet cover. It consists of a **display screen** and a **keypad**, containing groups of keys used to operate the device. There is a protective plastic coating on the keypad, which allows for efficient cleaning and disinfecting.

Display screen

The ACP 215 display screen provides procedural information and posts the volume of the processed solutions, pump speed, and procedure time. It is used to guide the operator through the procedure and display warning messages.

- 1. Display screen
- 2. Start/Stop keys
- 3. Pump control keys
- 4. Programming keys
- 5. Shaker key


Figure 2-2, ACP 215 control panel

Keypad

START and STOP keys

These two keys are active throughout the entire procedure and will provide the operator with the basic control of starting or interrupting the procedure.

PUMP control keys

These keys allow the operator to START, to STOP, or to increase (arrow up) or decrease (arrow down) the speed of the pumps during a procedure. The pump control keys are active when the pump is actually rotating.

PROGRAMMING keys

The MODIFY PROGRAM, YES, NO, SAVE PROGRAM keys will allow the operator to adjust a limited set of parameters, which are presented in each relevant procedure operation manual.

SHAKER key

This key is used only to start and stop the shaker; speed or amplitude cannot be modified.

DESCRIBING THE ACP 215 CABINET COMPONENTS

- 1. Centrifuge system
- 2. Solution pump
- 3. Blood pump
- 4. DPM
- 5. SPM
- 6. Red valve
- 7. White valve
- 8. Green valve
- 9. Blue valve
- 10. Yellow valve
- 11. Orange valve
- 12. SLAD
- 13. BLAD
- 14. Line Sensor
- 15. Solution bag pole


Figure 2-3, ACP 215 top deck view

Centrifuge system

The **centrifuge chuck**, seated in the **centrifuge well**, is designed to hold the disposable blow molded bowl (BMB) in place during operation and spin the bowl up to 8000 revolutions per minute (rpm).

The **centrifuge lid**, referred to as the *cover*, contains a locking knob; it will secure the disposable bowl in the chuck while isolating the spinning bowl from the operator. The **centrifuge drain tube**, located at the rear of the device, must be connected at all times to the biohazard waste bag supplied by Haemonetics®.

When installing a bowl, the operator should exert a downward pressure on the head of the bowl and ensure that the bowl is completely seated. The bowl will be completely secured once the operator has locked the centrifuge lid. To remove a bowl at the end of a procedure, the operator should open the centrifuge lid and pull upward on the bowl.

Caution: If any abnormality or noise occurs during the bowl seating test related to the spinning bowl, the operator should interrupt the procedure and remove and dispose of the disposable set. Initiate a new procedure using a new disposable set.

Warning: Under no circumstances, should the operator attempt to open the centrifuge lid if the bowl is still spinning. The bowl must come to a complete stop before attempting to open the centrifuge for any reason.

Do not use any bowl which cannot be properly seated in the centrifuge chuck.

Fluid detector

The centrifuge well is equipped with an electronic fluid detection system designed to detect the presence of liquid. The detector is mounted on the wall of the centrifuge well. The ACP 215 safety system will automatically stop the centrifuge and the pumps if there is contact between liquid of any sort and the fluid detector.

Pumps

Located on the left side of the ACP 215 top deck are two pumps which use peristaltic movements to displace fluids through the disposable set tubing.

Figure 2-4, ACP 215 pump rotor

Blood Pump

The blood pump will rotate during certain procedures simultaneously with the solution pump to direct solution to the RBC bag. As the bowl is filling, the blood pump directs the RBC product to the bowl.

At the end of certain procedures, the blood pump directs the final product from the bowl to the final product bag.

Solution Pump

The solution pump will rotate at certain points during the procedures to transfer solution to the RBC bag.

As the bowl is filling, or when transferring the final product from the bowl back to the RBC bag or the final product bag, the solution pump is stopped.

Pressure monitors

There are two electronically controlled pressure monitors on the ACP 215 device, which function with the correlating filter on the disposable set to measure pressure in the disposable tubing. The pressure monitors provide feedback to the system about the flow of blood components.

The draw pressure monitor (**DPM**), located on the left side of the top deck, measures pressure in the blood line tubing.

The system pressure monitor (**SPM**), located on the right side of the top deck, measures pressure in the effluent tubing. This measurement verifies that the sterile seal, between the head and the body of the centrifuge bowl, remains intact.

If a pressure monitor detects that pressure is outside of acceptable limits, the ACP 215 safety system will interrupt the procedure and provide an explanatory message.

Valves

There are six pneumatic pinch valves located on the ACP 215 top deck which automatically control the flow of fluids through the disposable set tubing. These valves are color-coded to facilitate accurate tubing installation.

The ACP 215 safety system will control the valves during the self-diagnostic tests. Once the operator has selected a procedure, the appropriate valves will automatically open, in preparation for installing the disposable set tubing. In case of power failure, the valves revert to a closed status.

Air detectors

The ACP 215 device is equipped with two ultrasonic sensors designed to detect the presence of air, bubbles or foam in the fluids flowing through the disposable set tubing. The solution line air detector (**SLAD**) is located to the right of the solution pump. The blood line air detector (**BLAD**) is located to the right of the blood pump.

In the case of an occlusion in the tubing, or an interruption in the fluid flow, the pump(s) will stop and a message will be displayed with information about the condition.


Figure 2-5, Example of an air detector

Optical line sensor

Located on the right side of the ACP 215 top deck is the optical line sensor which monitors the blood components passing through the effluent tubing.

Caution: The line sensor will not provide accurate readings if the optical lens it is obstructed in any way; thus the lens must be cleared of any extraneous substances to ensure proper functioning of the system.

Further information about specific functions of the line sensor is provided in each procedure operation manual.

Solution bag pole

Located on right side of the ACP 215 cabinet is a height-adjustable pole, used to hang solution bags during the procedure. There is a knob located on the base of the pole that permits the operator to adjust the height of the pole by pulling the knob outward to extend or retract the pole.

Power entry module

The power entry module (**PEM**) is located on the left side panel of the device. Externally, the module consists of an ON/OFF power switch and a power-input receptacle for the power cord. In the case of an emergency, the power switch can be used to stop all device function.

Internally, the PEM contains the fuse panel. It will interrupt power supply to the system in the event of an electrical current surcharge. The design of the PEM also provides a filter-effect for the device against the effects of a power surge.

- 1. ON/OFF power switch
- 2. Power input receptacle


Figure 2-6, Power entry module (PEM)

Power cord

The power cord provided by Haemonetics is designed to connect the ACP 215 device with an external power source via the power-input receptacle, located on the power entry module on the side panel.

Biohazard waste bag

The biohazard waste bag is designed to collect any biologically contaminated material from the centrifuge well in the rare case of a spill or leak. Two biohazard waste bags are supplied with the delivery of each ACP 215 device.

A bag must be attached at all times to the centrifuge drain tube, located at the rear of the device. The bag must hang freely, with the clamp open, visible to the operator.

Warning: The biohazard waste bags are not to be used to collect or store blood products. When a bag contains evacuated waste products, it must be clamped, removed and properly disposed of, according to the local policies concerning biologically contaminated material. A new bag must be placed before resuming operation.

DESCRIBING THE SHAKER AND THE PRINTER

- 1. Shaker connector port
- 2. Printer connector port


Figure 2-7, ACP 215 rear panel

Shaker

During certain procedures, a shaker device manufactured by Haemonetics properly mixes the blood in the RBC bag with the added solutions. The ACP 215 device provides power to the shaker and controls the shaker operation. The RBC bag is secured on the shaker platform with magnets.

The operator should connect the shaker to the rear panel of the ACP 215 device as follows.

- ➔ Place the shaker device to the left of the ACP 215 device.
- ➔ Attach the shaker power cord to the connector port labeled SHAKER on the rear panel of the ACP 215 device.
- ➔ Secure the connection using the two thumb screws.


Note: When ordering any additional shaker magnets from Haemonetics®, the following part number should be specified.

- ACP 215 shaker magnets, P/N 09437-00.

Printer

The printer is a satellite unit which provides a procedure summary printout automatically at the end of the procedure. Examples of a printout are provided in each procedure operation manual. The ACP 215 device controls the printer, however the printer has an independent power source via a wall unit adapter that converts AC into DC.

The printer is connected to the rear panel of the device through an RS232 serial port, using a 9 pin connector. The function of each pin is defined as follows.

Pin Number	Definition
1	RS-232C Data Carrier Detect
2	RS-232C Receive Data
3	RS-232C Transmit Data
4	RS-232C Data Terminal Ready
5	RS-232C RS-232C Ground
6	RS-232C Data Set Ready
7	RS-232C Request to Send
8	RS-232C Clear to Send
9	RS-232C Ring Indicator

Connecting the printer

The operator should connect the printer to the rear panel of the device as follows.

- ➔ Attach one end of the printer serial cable to the serial port labeled PRINTER on the rear panel of the ACP 215 device.
- ➔ Secure the printer cable to the rear panel connector using the two thumb screws.
- ➔ Attach the other end of the printer serial cable to the serial port located on the rear panel of the printer.
- ➔ Secure the cable in to the printer using the two thumb-screws.
- ➔ Insert the DC plug of the printer AC adaptor into the power supply jack located at the rear of the printer.
- ➔ Plug the printer AC adaptor power cord into an electric outlet.
- ➔ Place the printer to the right of the ACP® 215 device.
- ➔ Power-on the printer.
- ➔ Ensure that there is enough paper remaining in the printer.

To remove the printer AC adaptor, the operator should:

- ➔ Power off the printer first.
- ➔ Unplug the AC adaptor and the DC plug.

The operator should not touch the pins of the DC plugs.

Note: Further information about replacing the paper or cleaning the printer is provided in the Operation Manual for the DPU-414 Thermal Printer, Seiko® Instruments Inc.

Chapter 3

Maintaining the ACP® 215 Equipment

CLEANING PROCEDURES.	3-2
Listing tools and cleaning supplies	3-2
Exterior surfaces, valves, shaker and printer	3-3
Pressure monitors	3-3
Air detectors	3-3
Optical line sensor	3-3
Pumps.	3-4
Fluid detector	3-4
Centrifuge components.	3-4
Filter screens.	3-5
CUSTOMER SERVICE.	3-6
Clinical training	3-6
Field service	3-6
Returned Goods Authorization system	3-6
Haemonetics® Cleaning and maintenance record	3-8

CLEANING PROCEDURES

To ensure consistent performance of the device, the operator needs to primarily perform routine cleaning of the components. A record of routine cleaning should be kept along with any routine or preventive service maintenance performed by a Haemonetics® representative. An example of a cleaning and routine maintenance form is provided at the end of this chapter. Annual preventive maintenance should be performed by a Haemonetics service representative.

The frequency and type of cleaning each ACP 215 device will depend on the number of procedures performed. Special cleaning needs may arise and should be dealt with promptly. Haemonetics recommends the following routine cleaning schedule for each device, **based on an average of three procedures per day, or approximately sixty per month**. Adjustments should be made to the cleaning schedule according to the actual number of procedures performed.

- **Daily:** Clean the exterior surfaces.
- **Weekly:** Clean the pressure monitors, air detectors, optical sensors (line sensor and optical bowl sensor), the fluid detector and the inside of the centrifuge well.
- **Monthly:** Clean the pump rotors and the pump wells.
- **Quarterly:** Clean the filter screens.

Warning: To eliminate the potential danger of electrical shock, the operator must clean the ACP 215 device only when it is disconnected from an external power source.

Listing tools and cleaning supplies

A tool kit has been included with each device containing the following articles.

- Phillips-head screwdriver, to remove the retainer plate from the air vent.
- Allen (hexagonal head) wrench, to remove the pump heads from the pump assembly.
- A brush, to remove loose debris within the roller assembly.
- Silicon lubricant, used for the O-ring gasket of the centrifuge chuck.

The following list describes the basic material required for routine cleaning.

- Disinfectant cleaning solution, specific for blood borne pathogens.
- Clear warm water.
- 70% Isopropyl alcohol.
- Lint-free gauze or cloth for cleaning and drying.
- Cotton swabs.
- Protective gloves.
- Spray bottle, recommended for cleaning the centrifuge well.

Caution: Electric shock hazard, do not remove back panel. Refer to servicing to qualified service personnel.

Caution: Do not use solvents, strong alcohol solutions, or abrasive cleaning agents. The following cleaning solutions or active ingredients have been tested for use on the ACP 215 device:

- 12.5% Benzalkonium Chloride solution
- 70% Isopropyl Alcohol
- 8% Glutaraldehyde solution
- Coverage Plus solution
- 10% bleach/90% water solution (cannot be used on air detector)

10% bleach/90% water solution may fade the color of the device over time with continued use.

If there is no established institutional policy for decontamination, Haemonetics recommends that blood spills be cleaned with a cleaning solution followed by a wipe down with distilled/sterile water. Use lint-free gauze or cloth to apply the cleaning solution and water.

Exterior surfaces, valves, shaker and printer

The exterior cabinet, keypad, display screen, valves, shaker and printer should be wiped daily, as well as following any spill, using an appropriate cleaning solution.

Pressure monitors

The pressure monitors (DPM/SPM) should be cleaned in the following manner.

- Depress and hold the white ring as if installing the disposable filter.
- Wipe the silver rod thoroughly in a circular motion, using isopropyl alcohol or warm water.
- Dry the rod and release the pressure on the ring.

Caution: It is very important to completely dry the pressure monitor rod before installing the disposable set. Alcohol or cleaning solution residue could cause a reaction with the plastic material of the corresponding disposable set filter and affect the function of the filter.

Air detectors

The air detectors (SLAD/DLAD) are designed with a groove to hold the disposable tubing. The contents of the tubing are monitored by the sensors which are located internally on either side of this groove. If a procedure is interrupted due to an air detector alarm, the operator should remove the tubing and clean the groove before continuing the procedure.

The operator should use **only** warm water and lint-free gauze to clean and dry inside of the tubing groove. The groove should be kept free of any particles, such as powder residue from disposable gloves, since this could lead to an erroneous detection of air.

Caution: *Alcohol may cause the plastic housing to degrade and must not be used to clean the air detectors.*

Optical line sensor

The line sensor contains two very small lenses which are centered on either side of the disposable tubing groove. The lenses must be kept completely free of particles or debris which could produce inaccurate readings and influence the device performance.

The operator should carefully pass lint-free gauze moistened with warm water **only** through the groove to clean, then dry the lenses.

Caution: *If cleaning solution should come into contact with the optical sensor lenses, the lenses should be cleaned immediately with lint-free gauze and warm water, then thoroughly dried. Dried cleaning solution can leave an opaque film on the lens.*

Pumps

Remove debris from the rotors and the pump wells on a routine basis, as well as after any spills, to ensure efficient operation.

Use the following steps for routine cleaning:

- Loosen the pump screw and remove the pump rotor from the housing.
- Clean the pump well using warm water and lint-free cloth or gauze.
- Dry with lint-free cloth (or compressed air, if available).
- Spray the rotor with a center approved cleaning solution.
- Hold the rotor under running water, rolling the guides and rollers back and forth to allow the water to flush any accumulated debris.

Caution: *Do not allow the rotors to sit immersed in fluid.*

- Gently shake the excess water out of the rotor.
- Dry the rotor with a lint-free cloth and allow to air dry at least an hour.
- Ensure that all of the rollers spin freely.
- Place the rotor onto the pump shaft, rotating until the grooves on the bottom of the rotor drop onto the drive pin on the pump shaft.
- Replace and tighten the pump screw.

Note: The pump screw will not tighten if the pump head is not properly aligned with the pump shaft.

In the case of a fluid spill, use the same cleaning method but use disinfectant cleaning solution, followed by a clear water rinse. Do not immerse the pump rotor in water.

Fluid detector

The fluid detector is located inside of the centrifuge well. The surface of the detector should be cleaned using a cotton swab moistened with 70% alcohol.

Centrifuge components

The centrifuge well, chuck, hinged lid and locking knob can be wiped routinely using the cleaning solution and a lint-free cloth.

After a major cleaning, the operator should apply a small amount of the silicon lubricant to the O-ring gasket to prevent it from cracking. It is not necessary to remove the gasket when applying the lubricant.

Cleaning after a blood spill

In the case of a blood spill, the operator should:

- Power off the device and disconnect it from the external power source before cleaning.
- Ensure that the equipment waste bag is attached to the drain tube and hangs freely.
- Wipe the centrifuge lid with cleaning solution.
- Clean the centrifuge chuck and well using the disinfectant solution and a lint-free cloth until all traces of blood components are removed.
- Lubricate the O-ring gasket with a small amount of the silicon lubricant.

Haemonetics recommends that the operator wear protective gloves to avoid direct contact with the cleaning solution and/or any spilled blood which may be present.

In the case of a larger spill, fluid and/or blood may be evacuated into the biohazard waste bag. The operator should complete the following additional steps and contact the local Haemonetics representative for further instructions before using the device.

- Wipe away as much blood as possible with gauze pads.
- Irrigate the centrifuge drain holes with cleaning solution, until the drain tube is rinsed clear of the spilled material.
- Remove the bag and replace it with a new bag.
- Dispose of the used waste bag according to local established policies concerning the disposal of biohazard waste.

Note: A 50 ml syringe attached to a 20 cm section of disposable tubing placed in the drain holes can be used for irrigation. The biohazard waste bag should be monitored to avoid overfilling.

Warning: An authorized Haemonetics technician should perform a leakage current control after any major fluid spill involving the ACP 215 device. Leakage current represents a primary indication of electrical shock hazard and should be checked according to guidelines as described in local standard operating procedures.

Filter screens

The ACP 215 device is equipped with filter screens on the bottom of the cabinet, which eliminate dust from incoming cool air. The filters should be cleaned routinely, especially if dust becomes visible on the screens.

To clean the filters, the operator should:

- Remove the retainer plates using a Phillips-head screwdriver.
- Remove the filter screens from the panel.
- Rinse the screens under running water - DO NOT use any cleaning agents.
- Gently squeeze the screens to remove excess water.
- Place the screens on a clean, dry cloth and allow to dry completely.
- Reinsert the screens into the panel, ensuring that all openings are completely covered by the filter.
- Replace the retainer plates and tighten the screws.

Warning: To avoid electrical shock, the filter screen should be completely dry before it is reinstalled on the ACP 215 cabinet.

CUSTOMER SERVICE

Clinical training

Haemonetics employs a staff of Clinical Specialists to provide training in the use of the ACP 215 equipment. The local Haemonetics representative will schedule staff training upon delivery of the ACP 215 equipment and should be contacted to organize further instruction when needed.

Field service

Haemonetics maintains a worldwide network of company-trained service representatives responsible for responding to technical needs concerning equipment. These technical specialists are available to diagnose and repair any malfunctions, as well as provide routine annual or semi-annual maintenance of the equipment, including leakage current tests. If service beyond the routine maintenance and cleaning described in this manual is required, the local Haemonetics representative should be contacted to provide specific instruction.

Product Return Guidelines

If for any reason merchandise must be returned to the company, the customer should contact the local Haemonetics representative to arrange for repairs or returns using procedures to ensure proper handling and subsequent analysis. No returns will be accepted without advanced authorization.

Warning: Haemonetics products must be properly cleaned and packaged prior to their return. It remains an important responsibility of the customer to reduce potential health hazards by being aware of the risks involved in the shipping, handling and testing of this material. Units returned to Haemonetics for repair are subject to biohazard charges if any component is contaminated with blood or blood products.

Haemonetics® Cleaning and maintenance record

Device serial number: Year:

Action	January	Febr.	March	April	May	June	July	August	Sept.	Oct.	Nov.	Dec.
Clean cabinet, control panel, valves, shaker and printer												
Clean air detectors												
Clean pump rotors												
Clean centrifuge cover and well												
Clean optical line sensor												
Clean air filters												
Inspect O-ring: apply silicon lubricant												
Clean DPM and SPM												
Verify biohazard waste bag												
Maintenance performed by: (date and initials)												
Reviewed by: (date and initials)												

Annual preventive maintenance should be scheduled by a supervisor when appropriate and performed by a Haemonetics® service representative or a qualified engineer.

Annual maintenance performed by: Date performed:

Reviewed by: Date/initials:

Chapter 4

Ensuring Safety and Quality for an ACP® 215 Procedure

HANDLING THE ACP 215 EQUIPMENT	4-2
Storing the device and material	4-2
Inspecting the material	4-2
Installation	4-2
PREVENTING PROBLEMS DURING AN ACP 215 PROCEDURE	4-3
Understanding the risk of hemolysis	4-3
Avoiding the consequences of flow restriction.	4-3
Avoiding bowl misalignment.	4-4
Avoiding overheating due to mechanical situations.	4-4
WARNINGS FOR THE OPERATOR	4-5
Electrical shock hazards	4-5
Leakage current control	4-5
Mechanical hazards/rotating parts.	4-5
Power outlet connection	4-5
Communicable disease precautions	4-6

P/N 85271-30, Manual revision EB

HANDLING THE ACP 215 EQUIPMENT

Safe and successful ACP 215 operation will depend in part on the proper routine handling of the ACP 215 equipment. The operator should be aware of the problems which could result if the device or disposable material is stored, installed or used incorrectly.

Storing the device and material

The ACP 215 device must not be operated or stored in an area where flammable gases or vapors are present. The ACP 215 disposable set material should be kept in a dry, well-ventilated area and isolated from any chemical vapors. The operator should handle the disposable set elements with clean, dry hands or gloves.

Warning: *Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.*

Inspecting the material

Prior to installation, the operator should complete a visual inspection of the disposable set elements and control for twisted or flattened sections.

After installing the disposable set, the operator should verify the correct placement of the individual elements, prior to initiating a procedure. It is important that the tubing remain free of any twists or occlusions which could cause a flow obstruction.

Warning: *Upon completion of any procedure, the used disposable should be considered as biologically contaminated. It must be disposed of according to local standard operating procedure for the removal of such material and should not be mixed with non-biologically contaminated waste.*

Installation

Installation of the disposable set is explained in detail in each procedure operation manual.

The following disposables may be used with the ACP 215:

List No.	Protocol
LN 225	Glycerolization
LN 235	Deglycerolization and Cell Wash
LN 236	Deglycerolization

PREVENTING PROBLEMS DURING AN ACP 215 PROCEDURE

Understanding the risk of hemolysis

Hemolysis involves the destruction of red blood cell membranes, with the release of free hemoglobin into the plasma portion of the blood. Free hemoglobin does not have the capacity to transport oxygen and can produce serious problems. The remnants of the red cell can stimulate clot formation and damage the vascular nature of the lungs and the kidneys. This could lead to respiratory complications and/or renal failure.

Hemolysis of red cells can occur during a procedure in the rare event of a mechanically induced situation, such as overheating or excessive pressure.

Warning: Forcing a pump to work against a severe flow restriction can lead to hemolysis, and thus, consequently high levels of free hemoglobin. It is important that the operator remain attentive to this fact in the case of any “high pressure” alarms during operation.

If there is any suspicion that hemolysis has occurred, the operator should follow the local SOP applicable in this type of situation. The operator should also contact the local Haemonetics® representative if there is any suspicion of a mechanically induced problem.

Avoiding the consequences of flow restriction

As the bowl is filling, a flow restriction in the effluent tubing can create pressure on the outlet port of the disposable bowl. This unrelieved pressure can deform the rotary seal of the disposable bowl. If the functional characteristics of the rotary seal are altered, the increased friction and excessive heat can make the contents of the bowl unsuitable for further use.

As the bowl is emptying, a flow restriction in the effluent tubing can cause the pressure in the centrifuge bowl to drop severely. This sudden drop in pressure could potentially produce hemolysis.

To eliminate these potential problems, the operator should not clamp the effluent tubing.

P/N 85271-30, Manual revision EB

Avoiding bowl misalignment

An improperly installed disposable bowl can become misaligned as it spins. This can create excessive friction, and consequently overheat the bowl contents. The operator should follow the bowl seating procedure to ensure proper placement of the bowl in the centrifuge.

Warning: The operator must not use any bowl which cannot be properly seated in the centrifuge chuck. Overheating can occur, subsequently lead to hemolysis and make any blood being processed unsafe for reinfusion. During operation the operator should interrupt the ACP 215 procedure if an abnormality or noise occurs, related to the spinning bowl.

Avoiding overheating due to mechanical situations

Overheating could also result from a mechanical or maintenance-related problem, such as a defective bearing or seal within the centrifuge well. In this case, the operator should contact the local Haemonetics representative and discontinue use of the device until it is serviced.

Warning: If any component of the ACP 215 device has become overheated during a procedure, and thus overheats the blood being processed, the blood components cannot be considered safe for re-infusion.

WARNINGS FOR THE OPERATOR

Electrical shock hazards

The operator should always use the ACP 215 device with clean, dry hands, or gloves. The internal parts of the device contain various electrical components. Contact with any of these components, when the device is connected to an external powered source, could result in an electrical shock to the operator.

The operator should never remove any of the ACP 215 cabinet panels. Maintenance requiring access to the inner cabinet remains the responsibility of a Haemonetics-trained technician.

Leakage current control

Each ACP 215 device receives a careful inspection for leakage current prior to leaving the factory. Haemonetics recommends that a control be performed for leakage current by an authorized representative as part of the annual preventative maintenance. It remains the responsibility of the facility to ensure that this control is performed.

In the event of any major spill in which fluid may enter the cabinet, the operator is responsible to ensure that a leakage current test is performed before re-using the device. The test is necessary to avoid the risk of electrical shock and should be conducted by a Haemonetics trained technician.

Mechanical hazards/rotating parts

As with any equipment containing rapidly rotating parts, the potential for severe injury exists if contact is made, or if clothing becomes entangled with the moving parts. The ACP 215 device contains a safety feature, designed to prevent the centrifuge from spinning if the system has not been properly secured. However, the operator should respect the usual precautions taken when working with equipment containing rotating mechanical parts.

Power outlet connection

The power cord supplied by Haemonetics is a hospital grade power cord. Do not replace the power cord with a substitute. If necessary, contact the Haemonetics Customer Care Center for a replacement. Always ensure the power cord is connected to an appropriately grounded power source. To comply with the IEC 60601-1-1:2000 Standard for Medical Electrical Equipment, general requirements for safety, an additional Multiple Socket-Outlet or extension cord shall not be connected to the ACP 215 device.

Communicable disease precautions

Despite testing and screening to detect communicable diseases such as hepatitis, syphilis or HIV, the risk remains that the blood being processed may be infected. The operator must take the appropriate precautions when handling blood products and disposing of blood-contaminated material, to ensure personal safety as well as the safety of others who may come in contact with the material.

Proper handling of blood contaminated material

If a leak or blood spill should occur, it should be cleaned immediately. The operator should follow the local standard operating procedure outlining the steps to follow and product(s) to be used for the disinfection of blood-contaminated material.

If any blood-contaminated material must be returned to Haemonetics® for further inspection, the operator should consult the product return guidelines as described in the chapter *Maintaining the ACP® 215 Equipment*.

Proper disposal of biologically contaminated materials

Any disposable material used during an ACP 215 procedure is considered as biologically contaminated. It must be disposed of according to local standard operating procedure for the removal of such material and should not be mixed with non-biologically contaminated waste.

Chapter 5

Troubleshooting for the ACP® 215 Device

LISTING SYSTEM ERROR MESSAGES	5-2
LISTING INTERRUPT SAFETY SYSTEM MESSAGES	5-5

P/N 85271-30, Manual revision EB

LISTING SYSTEM ERROR MESSAGES

When the ACP 215 device is powered-on, several internal tests are performed. Any test failure will interrupt the test sequence and a system error message will be displayed. If a message of this type is displayed, the operator should record the message, power off the device and contact Haemonetics® Field Service.

If for any reason the system has been powered off, or an unexpected power failure has occurred during a procedure, there is no automatic recovery procedure. The operator should follow the local SOP applicable in this type of situation.

The following table lists the system error messages with a corresponding code number and relevant operator actions.

System Error Message	Corrective Action
CPU OPERATION ERROR CODE = 4 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
ROM TEST FAILURE CODE = 5 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
CALIBRATION CRC FAILURE CODE = 7 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
DPM FAILED 20% ANALOG TEST CODE = 10 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
DPM FAILED 50% ANALOG TEST CODE = 10 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
DPM FAILED 80% ANALOG TEST CODE = 10 Please POWER OFF.	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
SPM FAILED 20% ANALOG TEST CODE = 11 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
SPM FAILED 50% ANALOG TEST CODE = 11 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.

System Error Message	Corrective Action
SPM FAILED 80% ANALOG TEST CODE = 11 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
XXX* VALVE POSITION FAULT CODE = 12 PLEASE POWER OFF * XXX DENOTES THE COLOR OF THE VALVE	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
WATCHDOG ERROR CODE = 14 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
WATCHDOG ERROR CODE = 15 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
LINE SENSOR VOLTAGE TOO LOW CODE = 50 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
LINE SENSOR VOLTAGE TOO HIGH CODE = 51 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
CONFIGURATION CRC FAILURE CODE = 52 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
SPM FILTER INTERLOCK FAILURE CODE = 502 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
DPM FILTER INTERLOCK FAILURE SPM FILTER INTERLOCK FAILURE CODE = 504 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
DPM FILTER INTERLOCK FAILURE CODE = 507 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
AIR DETECTOR FAULT: [SLAD] CODE = 509 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.

System Error Message	Corrective Action
AIR DETECTOR FAULT: [BLAD] CODE = 510 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
DPM OUT OF TOLERANCE CODE = 514 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
SPM OUT OF TOLERANCE CODE = 515 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
CPU RELAY CLOSING FAILURE CODE = 701 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
CPU RELAY OPENING FAILURE CODE = 702 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.

LISTING INTERRUPT SAFETY SYSTEM MESSAGES

Specific inputs are checked at 50 ms intervals during an ACP 215 procedure. If an error is detected, all peripherals are shut down, and an error message is displayed.

For most these of situations, the procedure can be automatically recovered. However, in some of these situations, the operator will be required to power off the device and the procedure cannot be recovered.

The following table list these messages which may be displayed during a procedure and the relevant operator actions.

Interrupt Safety System Message	Corrective Action
CENTRIFUGE OVERSPEED PRESS YES TO CONTINUE	1. Press the YES key to continue. 2. If the message persists, stop the procedure and contact Haemonetics Field Service.
CENTRIFUGE UNDERSPEED PRESS YES TO CONTINUE	1. Press the YES key to continue. 2. If the message persists, stop the procedure and contact Haemonetics Field Service.
CENTRIFUGE COVER UNLOCKED PRESS YES TO CONTINUE	1. Verify the position of the cover lock and ensure that it is in the locked position. 2. Press the YES key to resume the procedure. 3. If the message persists, stop the procedure and contact Haemonetics Field Service
DPM FILTER NOT INSTALLED PRESS YES TO CONTINUE	1. Ensure that the DPM filter is correctly installed on the draw pressure monitor and reinstall if necessary. 2. Press the YES key to resume the procedure. 3. If the message persists, stop the procedure and contact Haemonetics Field Service.
SPM FILTER NOT INSTALLED PRESS YES TO CONTINUE	1. Ensure that the SPM filter is correctly installed on the system pressure monitor and reinstall if necessary. 2. Press the YES key to resume the procedure. 3. If the message persists, stop the procedure and contact Haemonetics Field Service
XXX* VALVE POSITION FAULT PRESS YES TO CONTINUE * XXX DENOTES THE COLOR OF THE VALVE	1. Press the YES key to continue. 2. If the message persists, stop the procedure and contact Haemonetics Field Service.
BLOOD PUMP POSITION FAULT PRESS YES TO CONTINUE	1. Ensure that the tubing is properly installed in the blood pump and reinstall if necessary. 2. Press the YES key to resume the procedure. 3. If the message persists, stop the procedure and contact Haemonetics Field Service

Interrupt Safety System Message	Corrective Action
SOLUTION PUMP POSITION FAULT PRESS YES TO CONTINUE	1. Ensure that the tubing is properly installed in the solution pump and reinstall if necessary. 2. Press the YES key to resume the procedure. 3. If the message persists, stop the procedure and contact Haemonetics Field Service
BLOOD PUMP WRONG DIRECTION PRESS YES TO CONTINUE	1. Press the YES key to continue. 2. If the message persists, stop the procedure and contact Haemonetics Field Service.
SOLUTION PUMP WRONG DIRECTION PRESS YES TO CONTINUE	1. Press the YES key to continue. 2. If the message persists, stop the procedure and contact Haemonetics Field Service.
CHECK SHAKER CONNECTION PRESS YES TO CONTINUE	1. Ensure that the shaker connector is attached to the rear panel of the ACP 215 device. 2. Press the YES key to resume the procedure 3. If the message persists, stop the procedure and contact Haemonetics Field Service.
FLUID IN CENTRIFUGE DETECTED Please POWER OFF	If there is a spill in the centrifuge, the centrifuge and pump(s) will stop. 1. Power off the ACP 215 device. 2. Refer to the appropriate section of the manual for guidance on handling a blood spill. 3. Contact Haemonetics Field Service.
WATCHDOG ERROR Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
SOFTWARE WATCHDOG FAULT INVALID STATE CODE = 802 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
SOFTWARE WATCHDOG FAULT CURRENT MODULE ID = XXX* EXPECTED MODULE ID = XXX CURRENT STATE = XXX CODE = 800 Please POWER OFF	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
* XXX DENOTES THE VALUE DISPLAYED	

Interrupt Safety System Message	Corrective Action
===== INTERNAL SOFTWARE ERROR MODULE: XXX* LINE: XXX STATE: XXX PARAM: XXX ===== * XXX DENOTES THE VALUE DISPLAYED	1. Power off the ACP 215 device. 2. Contact Haemonetics Field Service.
INTERNAL SOFTWARE ERROR CALL HAEMONETICS FIELD SERVICE	1. Power off the ACP® 212 device. 2. Contact Haemonetics® Field Service.

Chapter 6

Reference information

APPENDIX A: IEC 60601-1-2:2001 STANDARD REQUIREMENTS. 6-2

 Operation precautions 6-2

 Electromagnetic immunity 6-2

APPENDIX A: IEC 60601-1-2:2001 STANDARD REQUIREMENTS

Caution: The ACP 215® device must be operated in an environment compatible to the requirements of the IEC 60601-1-2 Standard, Electromagnetic compatibility (EMC).

Operation precautions

A power cord is supplied with the device. Do not replace the power cord with a substitute. If necessary, contact the Haemonetics® Customer Care Center for a replacement. Always ensure the power cord is connected to an appropriately grounded power source.

Mobile radio frequency (RF) communication equipment not approved by Haemonetics and portable communication equipment can affect the system. Any accessories and cables not approved by Haemonetics used in conjunction with the device may increase hazards and influence compatibility with EMC requirements. Therefore, non-approved accessories and cables must not be used.

In addition, the Cymbal device and accessories must not be placed directly adjacent to, or on top of other equipment, unless specifically approved by Haemonetics.

Caution: Grounding reliability can only be achieved when the equipment is connected to a properly grounded outlet.

Electromagnetic immunity

Note: The Cymbal device is intended for use in the electromagnetic environment specified below. The customer or operator of the Cymbal device should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The Cymbal system uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Note: The Cymbal system is intended for use in the electromagnetic environment specified below. The customer or operator of the Cymbal system should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8 kV air	±6kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, voltage variations on power supply input lines. EC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec.	<5% U_T (>95% dip in U_T) for 0,5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the operator of the Cymbal system requires continued operation during power mains interruptions, it is recommended that the Cymbal system be powered from an uninterruptable power supply or battery or Haemonet-ics approved battery.
Power frequency (50/60 Hz) magnetic field EC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the AC mains voltage prior to application of the test level.			

Note: The Cymbal system is intended for use in the electromagnetic environment specified below. The customer or operator of the Cymbal system should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environmental guidance
			Portable and mobile RF communication equipment should be used no closer to any part of the Cymbal system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio AM and FM broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Cymbal system is used exceeds the applicable RF compliance level above, the Cymbal system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Cymbal system. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environmental guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3Vrms	$d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio AM and FM broadcast and TV broadcast cannot be predicted theoretically with accuracy. To asses the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Cymbal system is used exceeds the applicable RF compliance level above, the Cymbal system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Cymbal system.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Note: The Cymbal system is intended for use in the electromagnetic environment specified below. The customer or operator of the Cymbal system should ensure that it is used in such an environment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz, the separation for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Julia Periman
INTERNATIONAL
RASPECTINGST
MAY 19, 2016

HAEMONETICS CORPORATION
400 WOOD ROAD
BRAINTREE, MA 02184 USA
1-800-225-5242

This Document contains 50 pages
that are signed, stamped and
bound pages.

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HAEMONETICS U.K. LTD
5 Ashley Drive, Bothwell
Scotland G71 8BS
Great Britain

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