VALIDATED

15 / 10 2014

**Additional information to the maintenance documentation
(instruction for use)**

Medical devices: Dianatal® obstetric gel (Dianatal® obstetric set)

Manufacturer of the medical device (Private Label Manufacturer): HCB Women's Health GmbH,
Surenweg 7, CH-6318 Walchwil, Switzerland
tel.: +41 41 544 22 20, fax: +41 41 544 81 76

e-mail: andreas.f.schaub@bluewin.ch

Manufacturer of the medical device (Original Equipment Manufacturer): Medena AG, Aemmenmattstrasse
2, CH-3123 Belp, Switzerland

tel.: +41 31 810 46 00, fax: +41 31 810 46 50

e-mail: michael.grassberger@medena.ch

General descriptions of the product:

Dianatal® obstetric gel is an odourless, almost colourless, highly viscous gel with a slightly acidic pH value of 5.5 – 6.5 and an electrical conductivity of minimal 5 mS·cm⁻¹ to facilitate vaginal delivery. When applied it forms a low-friction bioadhesive film on the vaginal birth canal, reducing the friction forces between the child and the vaginal birth canal.

Dianatal® obstetric gel contains following ingredients: propylene glycol, carbomer, hydroxyethylcellulose and purified water.

Dianatal® obstetric gel is filled in sterile single-use tubes with a long, thin insertion nozzle and a breakable seal and packed in single-use steribag. The seal must be removed before use. Each tube contains 5 ml of Dianatal® obstetric gel.

Dianatal® is non-allergenic, bioadhesive, electroconductive, slightly acidic, free of preservatives, latex, and phthalates and does not cause irritation to mucous membranes and or eyes and it does not contain any pharmaceutically active ingredients.

The containment consists of a white polyethylene, single-use tube (5 ml).

The packaging materials are sterilized by a validated-irradiation process.

No other accessory equipment, adapters or devices are necessary for the use.

Packaging:

Each tube is sterile, sealed in a steribag and packed in a cardboard-box with six units, including an instruction leaflet.

Storage and transport recommendations:

Storage condition recommendation:

Room temperature: 5°C to 30°C

Tolerated temperature excess (transport, etc.):

Excess to 40°C: max. 1 week

Excess not below - 10°C: max. 1 week

Disposal considerations:

Waste from residues/unused products: Dispose of in accordance with local regulations.

Contaminated packaging: Offer rinsed packaging material to local recycling facilities.

Unused gel in tubes, empty tubes as well as unused gel remained after the procedure must be disposed in accordance with applicable health rules and regulations of the Russian Federation (SanPIN) for the handling of medical waste and the Guidelines for disinfection, and sterilization of medical devices (MU).

Warranty and liability limitations:

HCW WOMEN'S HEALTH GMBH warrants that these products are in compliance with the legal requirements, in particular those set forth in the EU Directive 93/42/EEC concerning medical devices. Therefore, in case the product exhibits non-conformities or defects that can be directly attributed to material or processing errors that occurred during its manufacture or packaging ("non-conformities or defects"), and said non-conformities or defects appeared prior to the expiration date (UBD), HCB WOMEN'S HEALTH GMBH will, at its discretion, replace the product:

- with one of equal construction, or
- with one exhibiting the identical function.

All the following conditions must have been met in order to make claims within the context of this warranty:

- the product must not have been used after its UBD.
- the non-conformity or defect did not initially occur after the UBD.
- if the product in question is a disposable product or is a product, which HCB WOMEN'S HEALTH GMBH designates as being for single use, the fault must be detected at the first use.
- the non-conformity or defect was immediately reported to HCB WOMEN'S HEALTH GMBH or its representative.
- the non-compliant or defective product was forwarded to HCB WOMEN'S HEALTH GMBH or its representative immediately after the non-conformity or defect was reported.

The above warranty in no way limits and applies in conjunction with any and all regulations governing warranties on the part of the seller and the product liability laws.

For the remainder, this warranty applies only to the extent defined above. Additional claims arising from this warranty are excluded. This warranty, as well as the cited legal claims applies in place of, and to the exclusion of, all other warranties, either verbal or written, or express or implicit in the part of HCB WOMEN'S HEALTH GMBH. None of our representatives or employees is empowered to reach differing agreements.

Due to numerous factors such as shipping, storage, and handling on the part of the user, all of which lie outside the control of HCB WOMEN'S HEALTH GMBH, HCB WOMEN'S HEALTH GMBH is not liable, in particular, for damages of any type resulting directly or indirectly from the product or its employment unless said damages can be shown to be the result of a failure on the part of the product.

In particular, HCB WOMEN'S HEALTH GMBH assumes no liability for damages resulting from sterilization or reuse of products approved for single use, or which HCB WOMEN'S HEALTH GMBH designates as being for single use.

Instructions for use must be accurately studied before using of the product.

Only highly qualified medical staff with special training to carry out the relevant procedures should operate medical devices manufactured by HCB WOMEN'S HEALTH GMBH.

All the pending disputes should be resolved in the courts at the location of HCB WOMEN'S HEALTH GMBH, in accordance with Swiss law.

Zug.

15.10.2014

Dr. Andreas F. Schaub
CEO HCB Women's Health GmbH

Official Certification

Seen for authentication of the above signature
affixed by

Andreas Schaub

born

20.02.1965, Pass X 040 9977

The Notary Public assumes no responsibility
regarding the content or authenticity of the
undersigned document.

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