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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of BLUE SAIL MEDICAL CO., LTD. on the annexed INSTRUCTION FOR USE ON MEDICAL DEVICE is genuine.

China Council for the Promotion
of International Trade

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Authorized

Signature:

Sun Danfeng

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
(Date: May. 10, 2023)

证明网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTION FOR USE ON MEDICAL DEVICE:
Medical diagnostic (examination) disposable non-sterile gloves
“BENOVY” nitrile finger textured powder free
chlorinated (beaded cuff)

Sales Manager

(signature, stamp)

 **Mr. Bin Li /**

Initial issue of operational documentation
(instruction for use) for medical device: May 2020

Ver. 02

2023

TABLE OF CONTENTS

Data on the medical device	3
1.1 Name of the medical device	3
1.2 Manufacturer (producer) of the medical device	3
1.3 Purpose of the medical device	3
1.4 Scope of application of the medical device	3
1.5 Intended users of the medical device	3
1.6 Operation (application) conditions of the medical device	3
1.7 Indications for use of the medical device	3
1.8 Contradictions of the medical device	4
1.9 Shelf life of the medical device	4
2 Instruction for use of the medical device	5
3 Warning and precautions	6
4 Transportation, storage and disposal conditions of the medical device	7
4.1 Transportation conditions of the medical device	7
4.2 Storage conditions of the medical device	7
4.3 Disposal conditions of the medical device	7
5 Methods and means of disinfection and cleaning of the medical device	7
6 Specification of the medical device	8
6.1 Description of the medical device	8
6.2 Physical specifications of the medical device	9
6.2.1 Physical and mechanical specifications of the medical device	9
6.2.2 Inspection levels and acceptable quality level of the medical device	9
7 Circumstances, in which the consumer must consult healthcare professional	10
8 Manufacturer's warranties	10
9 Data on the authorized representative of the manufacturer	10
10 Scope of delivery of the medical device	10
11 List of national standards applied by the manufacturer (producer) of the medical device	11

1. Data on the medical device

1.1 Name of the medical device

Medical diagnostic (examination) disposable non-sterile gloves "BENOVY" nitrile finger textured powder free chlorinated (beaded cuff), color: blue, lilac-blue, black, white, pink, green, fuchsia, violet size: XS, S, M, L, XL.

1.2 Manufacturer (producer) of the medical device

Manufacturer (producer) of the medical device: Blue Sail Medical Co., Ltd. (No.21 Qingtian Road, Qilu Chemical Industrial Park, Zibo City, Shandong Province, China).

Address of the medical device manufacture: Blue Sail Medical Co., Ltd. No.21 (Qingtian Road, Qilu Chemical Industrial Area, 255414, Zibo, Shandong Province, P.R. China

1.3 Purpose of the medical device

The gloves are intended for single use as protective barrier in the process of wearing on healthcare professional's hands during examination/treatment of a patient or for other sanitary goals.

1.4 Scope of application of the medical device

Gloves are used during medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces, providing hand hygiene, as also protection of patient and healthcare professional from cross-contamination.

1.5 Intended users of the medical device

Medical personnel, who carry out medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces and other activities, which require compliance with the requirements for protection of patient/healthcare professional from cross-infection.

1.6 Operation (application) conditions of the medical device

Medical examination disposable gloves are used in medical or healthcare facilities indoors in a dry place, away from direct sunlight, at temperature from 0°C to 25°C and relative humidity from 65% up to 85% in accordance with the instruction for use on medical device.

Gloves for external use only (the device is non-invasive).

1.7 Indications for use of the medical device

The indications for use of the medical device are the necessity of conduction of medical researches, diagnostic, therapeutic, dental and treatment procedures, sanitation and laboratory researches, as well as handling contaminated medical materials/devices/surfaces.

1.8 Contradictions of the medical device

The medical device has no contradictions.

1.9 Shelf life of the medical device

The established shelf life of examination nitrile gloves is 3 years.

2 Instruction for use of the medical device

Preliminary hygienic treatment of your hands is required before use of the gloves (treat hands with skin antiseptic or wash with soap and water). Gloves are put on dry hands only.

How to put the gloves on: Take a glove from the box. It is recommended to touch the bounded surface of the glove corresponding to the wrist (upper edge of the cuff). Pull on the first glove. Take the second glove with your bare hand, touching the bounded surface of the glove corresponding to the wrist. To avoid touching the skin of the forearm with a gloved hand, rotate the outer surface of the glove to be put on with folded fingers in the glove so that you could pull the glove over the other hand. Having put the gloves on, it is required to touch only those objects provided for by the indications for use and conditions of the use of the gloves. To avoid contamination, it is required to change gloves before proceeding with the new procedure.

The used gloves should be put off carefully, avoiding splashing from the gloves' surfaces where microbial contamination of hands and environment is possible.

How to take the gloves off: To take a glove off, pick it up at the level of the wrist without touching the skin of the forearm and pull it off your hand, turning it inside out. Holding the removed glove with the gloved hand, it is necessary to place the fingers of your bare hand between the glove and wrist. Remove the second glove by rolling it down the hand and put into the first glove. Discard the gloves and wash your hands thoroughly.

3 Warning and precautions

The device is intended for single use;

- For individual use only;
- Every time it is necessary to perform hygienic treatment of your hands before putting on or changing gloves, as well as after removing the gloves;
- It is necessary to choose the gloves of the right size - they should be easily put on and should not impede movement or squeeze the hands. Do not use in case of size mismatch;
- Touching is recommended to the surface of the glove corresponding to the wrist (upper edge of the cuff);
- Depending on the type of the performed manipulation, to keep an optimal level of protection, it is recommended to change gloves every 60-120 minutes;
- It is necessary to change the gloves in case of any defect, puncture, damage, contamination, fluid contact with skin under the glove, as well as stickiness/ If one glove is damaged, one should change both;
- It is necessary to put the gloves on dry hands only;
- Do not wet and keep away from moisture;
- Avoid exposure to direct sunlight;
- Gloves contain residual process chemicals that may cause an allergic reaction;
- Gloves do not guarantee 100% hand protection against punctures, cuts and damage caused by exposure to chemicals;
- It is recommended to consult a healthcare professional in the event of dermal hypersensitivity reactions when using the device.

4 Transportation, storage and disposal conditions of the medical device

4.1 Transportation conditions of the medical device

The transportation of the device may be made by all means of transport (sea, land, railway or air transport) according to the rules of transportation, which are applied to such means of transport. During transportation, the storage conditions should not be violated.

4.2 Storage conditions of the medical device

Devices should be stored in a dry place protected from direct sunlight, heat, ignition and ultraviolet sources, as well as ozone, at temperature from 0°C to 25°C and relative humidity - from 65% up to 85%.

Products should be handled with caution when they are stored at temperatures below 15°C, as during careless handling, they can be stiff and sensitive to deformation. After storage at a low temperature, gloves should be warmed up to 30°C before use throughout all product volume.

4.3 Disposal conditions of the medical device

The disposal of the medical device in the Russian Federation shall be carried out by medical organizations as class B waste in accordance with SanPiN 2.1.3684-21 "Sanitary-epidemiological requirements for medical waste management" and approved guideline of the organization.

5 Methods and means of disinfection and cleaning of the medical device

Examination gloves are a disposable medical device. Reuse of the medical device is a gross violation of application methods.

6 Specification of the medical device

6.1 Description of the medical device

Form. Straight fingers with a thumb located on the same plane with other fingers, ambidextrous (non-anatomic). The glove's cuff has a standard length. Rolled cuff.

Surface. Medical gloves have a powder-free surface. Powder-free gloves are treated with chlorine (chlorinated).

Texture of the external surface. Gloves have a textured surface on fingers.

Sizes. Depending on the size, table 1 gives the values of width, length and thickness of the medical device.

Table 1 – Sizes of the medical device:

Value name	Size				
	XS	S	M	L	XL
Width W (figure 1), mm	≤80	80±10	95±10	110±10	≥110
Length l (figure 1), mm, min	220		230		
Thickness (in points, shown on figure 2), mm, min: - smooth surface - textured surface	0.08 0.11				
Thickness of rolled cuff, mm, max	2.50				


Figure 1 - Points for width and length measurement of the glove


Figure 2 – Points for wall thickness measurement of the glove

6.2 Physical specifications of the medical device

6.2.1 Physical and mechanical specifications of the medical device

Physical and mechanical parameters of the medical device are represented in table 2.

Table 2 - Values of physical and mechanical parameters of the medical device:

DESCRIPTION	SPECIFICATION
<u>1. Before accelerated ageing:</u>	
- force at break, N, min	7,0
- elongation at break, %, min	500
<u>2. After accelerated ageing:</u>	
- force at break, N, min	6,0
- elongation at break, %, min	400

6.2.2 Inspection levels and acceptable quality level of the medical device

For the purposes of quality control, medical gloves are selected and examined in accordance with standard ISO 2859-1.

Hermeticity. Medical gloves are hermetic. Acceptable quality level (AQL) is ≤ 1.5 .

7 Circumstances, in which the consumer must consult healthcare professional

The patient is recommended to consult healthcare professional in the event of dermal hypersensitivity reactions, when using the device.

8 Manufacturer's warranties

The manufacturer (producer) guarantees preservation of all the parameters and characteristics of the medical device taking into account acceptable quality level (AQL) ≤ 1.5 , as well as the quality of the medical gloves within all shelf life (3 years) from the date of manufacture upon the condition of package integrity, observance of transportation, storage and operation conditions in accordance with the instruction.

9 Data on the authorized representative of the manufacturer

Please contact the authorized representative of the manufacturer for any and all purposes, related to circulation of the medical device on the territory of the Russian Federation:

"BENOVY" Company Limited ("BENOVY" Co. Ltd.)

Russia, 194292, St. Petersburg, 1, 5-y Verkhniy pereulok, litera J, office 205

Tel.: +7 (812) 309-16-32

www.benovy-company.ru

e-mail: benovy@benovy-company.ru

10 Scope of delivery of the medical device

Gloves may be delivered in the following scopes of supply:

- bulk packaging - 2000 gloves, 1000 pairs in a consumer package;
- 400 gloves, i.e. 200 pairs in a consumer package;
- 200 gloves, i.e. 100 pairs in a consumer package;
- 100 gloves, i.e. 50 pairs in a consumer package;
- 50 gloves, i.e. 25 pairs in a consumer package;
- 40 gloves, i.e. 20 pairs in a consumer package;
- 20 gloves, i.e. 10 pairs in a consumer package;
- 10 gloves, i.e. 5 pairs in a consumer package;
- 4 gloves, i.e. 2 pairs in a consumer package;
- 2 gloves, i.e. 1 pair in a consumer package.

List of national standards applied by the manufacturer (producer) of the medical device

Table 3 – List of standards, applied by the manufacturer:

Designation	Name
ISO 13485	Medical devices. Quality management systems. Requirements for regulatory purposes.
ISO 14971	Medical devices. Application of risk management to medical devices.
ISO 15223-1	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1. General requirements.
ISO 2859-1	Statistical methods. Sampling procedures for inspection by attributes. Part 1. Sampling schemes indexed by acceptance quality limit for lot-by- lot inspection.
ASTM D6319	Standard specification for nitrile examination gloves for medical application.
EN 455-1	Medical gloves for single use. Part 1. Method for determination of tightness.
EN 455-2	Medical gloves for single use. Part 2. Methods for determination of physico-mechanical properties.
EN 455-4	Medical gloves for single use. Part 4. Requirements and testing for shelf life determination.

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