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Twintree Tower (B) #101,6 Yulgok-ro

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Twintree Tower(B) #101,6 Yulgok-ro,
Jongno-gu, Seoul, Republic of Korea

Cellprep

The next generation of Liquid-based Cytology

Cellprep Urine, CSF
Solution for Liquid-Based
Cytology Processor

Instructions for Use

www.biodyne.asia

Bio dyne

Kyungeun Jang
Chief Quality Officer

Approved _____

(Signature)

Biodyne

24-4, Achesan-ro 5-gil
Seongdong-gu, Seoul, Republic of Korea

Table of contents	
1. General information.....	5
2. Description of the Cellprep Urine, CSF Solution.....	6
2.1. Intended use, indications and contra-indications	6
2.2. Description of the solution	7
2.3. Principles of operation	7
2.4. Delivery set	7
3. Type of the sample analyzed.....	7
4. Information on conditions necessary for collection, processing and preparation of samples, data on stability of the samples analyzed, including storage conditions and duration, transportation conditions, restrictions on freezing (thawing) cycles	7
5. Information on consumers and potential users	8
6. Installation, preparation	8
7. Application of the Cellprep Urine, CSF Solution	9
7.1. Description of the medical or non-medical devices intended for use in combination with Cellprep Urine, CSF Solution.....	9
7.1.1. <i>Equipment for slide preparation for thin layer cytology or liquid-based cytology</i>	9
7.1.2. <i>Slides</i>	9
7.1.3. <i>Vortex</i>	9
7.1.4. <i>Centrifuge</i>	9
7.2. Application of the Cellprep Urine, CSF Solution	10
7.2.1. <i>Urinary System: Kidney, Ureter, Urinary bladder, Urethra</i>	10
7.2.2. <i>CSF</i>	10
7.3. Insertion of the vial with the Cellprep Urine, CSF Solution into equipment for slide preparation for thin layer cytology/liquid-based cytology	11
7.4. Interfering substances or sample limitations that may interfere with study result	19
7.5. Identification of the medical devices in order to obtain a safe combination and (or) information on known restrictions on the combined use of medical devices	19
8. Specifications of the Cellprep Urine, CSF Solution	20
9. Recommendations for staining the slides	21
10. A. B. C. Suction Numerical Standard Values for LBC Processors.....	22
11. Cautions and warnings	22
12. Transportation and storage conditions.....	22
13. Packaging	23

14. Symbol identification	24
15. Manufacturer's responsibility.....	25
16. Repair and Maintenance	25
17. Cleaning and Disinfection	25
18. Warranty and shelf life	25
19. Disposal	26
20. Information on the initial release and revision of the Instructions for Use.....	26
21. List of the standards applied.....	26
Appendix No. 1.....	29
Certificate of Analysis - Cellprep Urine, CSF Solution (sample)	29
Appendix No. 2.....	30
Material Safety Data Sheet.....	30

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Cellprep

1. General information

Name of the medical device - **Cellprep Urine, CSF Solution for Liquid-Based Cytology Processor** (hereafter - "Cellprep Urine, CSF Solution", "Urine, CSF Solution", "Urine, CSF Solution", "Cellprep Solution(s)", "medical device").

Use:

For detecting oncological and inflammatory diseases of the urinary system (kidneys, bladder, ureter, urethra, prostate) and brain.

Collected cells:


Urine (From daily urine, Cystoscopic)

CSF (lumbar puncture, suboccipital puncture)

The specially designed vial, in which the solution is kept, is optimized for Cellprep LBC devices.

The vial cap consists of two parts: middle cap and top cap. This design provides that after collecting the sample from the patient, the physician can open only the middle cap and place the sample in a vial with solution. The top cap of the vial can be opened only by the Cellprep LBC Processor during slide preparation. The middle cap is equipped with an integrated nano-polyester filter that helps to remove impurities and obtain a monolayer preparation.

The medical device is intended for clinical laboratory diagnostics (in vitro diagnostics).

The device is for single use only.

Class of potential risk of solution application – **2a**

Type of the medical device in accordance with the Russian Nomenclature Classification – **214660**.

Solution name	Catalog No.	Volume
Cellprep Urine, CSF Solution for Liquid-Based Cytology Processor	UC-1000	10 ml

The medical device is intended for professional use only.

The Cellprep Urine, CSF Solution for Liquid-Based Cytology Processor relates to the medical devices which have no direct or indirect contact with the human body.

The medical device is not sterile, and it is not a subject to sterilization.

2. Description of the Cellprep Urine, CSF Solution

2.1. Intended use, indications and contra-indications

Intended use and indications:

The Cellprep Urine, CSF Solution is an auxiliary device for in vitro diagnostics and is intended for preservation of urothelial urine and CSF cells for subsequent preparation of slides by liquid-based cytology using the Cell Prep Plus LBC processor or the Cellprep AUTO LBC system.

Monolayered slides prepared are stained in accordance with approved protocols, followed with the results interpretation.

Contra-indications:

Not intended for preservation of histological samples (biopsy).

Side effects:

- not applicable.

Description of target analyte	Not applicable
Functional use	Auxiliary device for in vitro diagnostics
Specific pathology, conditions or risk factor which detection, determination or differentiation the IVD medical device is intended for	Suspected neoplasms or precancerous condition of the kidneys, bladder, ureter, urethra, prostate, brain.
Population, demographic aspects of the medical device use	The Cellprep Urine, CSF Solution is intended for optimal preservation of urothelial urine and CSF cells, collected from patients, regardless of their gender, race and age differences.
Information on analytical sensitivity (limit of detection), analytical specificity, diagnostic sensitivity and diagnostic specificity.	Not applicable

2.2. Description of solution

Colorless transparent liquid without foreign impurities

Solution composition

Components:

CHEMICAL NAME	CAS No.	% (ppm)
Ethanol	64-17-5	42%
Deionised water	7732-18-5	51%
Glycerine	56-81-5	6%
Phosphorous acid	7664-38-2	0.45%
Ammonium acetate	631-61-8	0.45%
Sodium Citrate Tribasic Dihydrate	6132-04-3	0.10%

The Cellprep Urine, CSF Solution does not contain medicinal products for medical use, materials of animal and human origin.

2.3. Functional principle

- The solution preserves and fixes the cells collected.
- The urine fixed in the vial, while keeping the temperature regime from +10 °C to +30 °C, can be used for liquid-based cytology studies for up to 3 months.
- The CSF fixed in the vial shall be kept at temperature from +10 °C to +30 °C no more than 5 hours.

2.4. Delivery set

Cellprep Urine, CSF Solution for Liquid-Based Cytology Processor

Including:

1. Cellprep Urine, CSF Solution for Liquid-Based Cytology Processor, vial, 10 ml - 25 pcs/pack.
2. Instructions for Use in paper - 1 pc. for a shipping box (200 vials)

3. Type of the sample analyzed

Urine (From daily urine, Cystoscopic)

CSF (lumbar puncture, suboccipital puncture)

4. Information on conditions necessary for collection, processing and preparation of samples, data on stability of the samples analyzed, including storage conditions and duration, transportation conditions, restrictions on freezing (thawing) cycles

Recommendations for organization of preliminary (pre-laboratory) stage of work with urine and CSF samples.

Sample collection	Urine. For cytological study of urine sediment, it is better to use the first portion of first void urine, which contains the largest number of cellular elements. Additionally, urine obtained at other times of the day can be used. Pour collected sample (not less than 30 ml urine) into 50 ml tube and centrifuge at 2000 rpm for 5 minutes. Remove supernatant, then pour the sediment into Urine, CSF Solution Vial.
-------------------	---

	CSF. In most cases, a lumbar puncture is used to obtain CSF. A lumbar puncture is performed by a physician in a medical institution. The patient lies on his/her side with knees drawn up to the chest, then a needle is inserted into the space between 4 and 5 vertebrae in the subarachnoid cavity. The first five CSF drops are removed, as they contain bloody tap from blood vessels damaged during the procedure. The fluid is collected in the sterile tubes. Safe volume to be collected: adults - 8-10 ml, children - 5-7 ml, infants - 2-3 ml. Prohibited to shake the specimen collected, expose it to temperature differences, as this changes its parameters. The collected fluid is transferred from the tube into the Urine, CSF Solution vial. Specimen should be prepared for 2 hours.
Preservation	The Cellprep Urine, CSF Solution preserves and fixes the cells collected.
Label	Label the vial with sample in the designated area on the vial label (hand lettering according to referral or barcode label).
Storage	The urine fixed in the vial, while keeping the temperature regime from +10 °C to +30 °C, can be used for liquid-based cytology studies for up to 3 months. The CSF fixed in the vial shall be kept at temperature from +10 °C to +30 °C no more than 5 hours. Freeze/Thaw Prohibited. Cooling, especially freezing, of both non-fixed and fixed material is strictly prohibited.
Transportation	Transportation temperature regime: from +10 °C to +30 °C.

5. Information on consumers and potential users

For professional use only (e.g., pathologist, laboratory assistant).

Before using the Cellprep solutions, personnel should carefully read the Instructions for Use and complete the training.

6. Installation, preparation

Not applicable.

The Cellprep Urine, CSF Solution is completely ready to use and does not require any preparation.

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7. Application of the Cellprep Urine, CSF Solution

7.1. Description of the medical or non-medical devices intended for use in combination with Cellprep Urine, CSF Solution

7.1.1. Equipment for slide preparation for thin layer cytology/liquid-based cytology

1) Cell Prep Plus LBC Processor with accessories by Biotdyne Co., Ltd., Korea (registration certificate No. ΦC3 2010/07117) is intended for fixation and transfer of cells of various body fluids and tissue samples from a vial to a slide creating a cell monolayer so that the samples are suitable for further staining and microscopy analysis. Hereafter referred to as the Cell Prep Plus.

2) Cellprep AUTO LBC System by Biotdyne Co., Ltd., Korea (registration certificate No. ΦC3 2020/11462) is an auxiliary equipment intended for preparation of a monolayered cytological sample using LBC method for its further staining and examination with the microscope for the purpose of early diagnosis of precancerous conditions, detection of oncopathology and the degree of its malignance. Hereafter referred to as the Cellprep AUTO.

The vial containing the Cellprep Urine, CSF Solution and the sample inside must be placed in one of the above devices. Through the process, the device smears the cells on a slide producing a mono-layered circle of cells for further examination. It minimizes misdiagnosis by reducing the amount of overlapping cells on a slide.

7.1.2. Slides

In addition to the slides included in the delivery set of Cell Prep Plus LBC Processor, it is allowed to use the slides that have similar parameters and characteristics, such as "Microscopic Slides and Coverslips" by Leica Biosystems Richmond, Inc., USA, registration certificate No. P3H 2013/1157.

7.1.3. Vortex

To shake and mix the vial with the Cellprep Urine, CSF Solution and the sample, it is recommended to use:

- Vortex mixer V-3, SIA "ELMI", Latvia, registration certificate No. P3H 2017/5466;

or other medical devices of vortex type with similar technical characteristics, approved for circulation within the Russian territory.

7.1.4. Centrifuge

For centrifugation of samples, it is recommended to use the centrifuge model "Z 206A" of the medical device "Laboratory centrifuges Z 100M, Z 160M, Z 200A, Z 206A, Z 233M-2, Z 216MK, Z 300, Z 300K, Z 323, Z323K, Z 383, Z 383K, Z 400, Z 400K, Z 513, Z 513K, Z 366, Z 36HK, SIEVA-2 with/without accessories", HERMLE Labortechnik GmbH, Germany, registration certificate No. P3H 2016/4723;

or other centrifuges with similar technical characteristics, approved for circulation within the Russian territory.

7.2. Application of the Cellprep Urine, CSF Solution

7.2.1. Urinary System: Kidney, Ureter, Urinary bladder, Urethra

Fig. 1 Collecting sample	Fig. 2 Centrifuge	Fig. 3 Pouring to vial, vortex	Fig. 4 Insert vial to Cellprep Processor
			
	Tube (50ml)/2000rpm, for 5 min	Collect sediment (Cellprep vial)	Pre-fix, then insert

- 1) Pour collected sample (not less than 30 ml) into 50 ml tube and centrifuge at 2000 rpm for 5 minutes (Fig.2).
- 2) Remove supernatant, then pour the sediment into Cellprep Urine, CSF Solution Vial, then vortex for 5 secs and leave to fix for more than 20 minutes.(Fig 3)

Only one patient sample is placed per vial.

After this stage, put the vial with urine sample in LBC Processor.

Depending on the equipment,

1. Cellprep AUTO: place the vial in a tray in order. As processor has auto vortexing system, there is no need to use vortex outside of processor.
2. Cell Prep PLUS: before placing the vial in a processor, vortex it for 10-20 seconds or shake 5-10 times.

7.2.2. CSF

- 1) Since the usual amount of collected CSF specimen is a very small, insert all collected specimen into the Cellprep Urine, CSF Solution Vial without any pre-preparation as there is a possibility of cell damage or loss during this step. However, if the specimen is too bloody, centrifuge it then pour sediment into vial.
- 2) Thoroughly mix the specimen, then leave to pre-fix for more than 20 minutes.

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- Pour the specimen in Cellprep Urine, CSF Vial, prepare a slide within 2 hours. Storage of CSF sample in the vial should not exceed 5 hours.

Only one patient sample is placed per vial.

When to use Vortex mixer:

1. Cellprep AUTO: Due to characteristic of CSF specimen, user needs to mix the vial instead of using vortex function of device.
2. Cell Prep PLUS: Due to characteristic of CSF specimen, user needs to mix the vial.

7.3. Insertion of the vial with Cellprep Urine, CSF Solution into equipment for slide preparation for thin layer cytology/liquid-based cytology


The innovative patented air transfer technology (Blowing Technology™) used in the Cellprep processors allows to get a high quality monolayer slide and reduce the loss of cells during cytological examination to zero.

Depending on the equipment,

Cellprep AUTO: after turning on the instrument, select the operation mode on the touch screen (Auto Turbidity Smear Mode or Manual Turbidity Mode).

Detailed description of how to operate the Cellprep AUTO LBC System is provided in the corresponding user manual.


Auto Turbidity Smear Mode

Please click on the **SETUP** button for 3 seconds or longer in the main screen.


The screen will switch to User Setup Screen.


Select Turbidity ON.

Using the auto- turbidity selection.

Place the slide in the rack with the barcode facing towards the slide holder. The barcode should face the barcode scanner, and make sure the slide holder is in the correct position.


Prepare vials and load them in the tray.


Up to 40 vials can be loaded at once. The slide, completed a smear, is come out according to the order of Vial's work.

Since the equipment has an integrated vortex, there is no need to specially mix up the vials.


Vortex should not be used for CSF specimens.

Put the vials from the front so no vacancy remains.


Push the tray by aligning the groove on the location of the next section when mounting the vial tray.

Select the types of specimens


Set the mode according to the type of specimen by pressing the buttons on the touch screen.

Press the green "U" button for urine and CSF samples.

Specify the quantity of vials for operating.

Number press menu comes out when press red number on **Working Count**.

Enter the desired quantity of vials in the field.

Cellprep AUTO System

Cell Status (Manual Mode)

A B C 88888

Type of Specimen

Gyne

G U B E R

Filter Mgt Forward

Status

Message

START STOP RESET SETUP MANUAL INITIAL

Total Vial Processing Result

Left Vial Tray

Right Vial Tray

Main Slide Result

Slide Glass Processing Result-1

Number	Gyne/NonGyne	Type	Barcode Data
Slide-01	A	A	ABDEFGHIJKLMNOP
Slide-02	A	A	ABDEFGHIJKLMNOP
Slide-03	A	A	ABDEFGHIJKLMNOP
Slide-04	A	A	ABDEFGHIJKLMNOP
Slide-05	A	A	ABDEFGHIJKLMNOP
Slide-06	A	A	ABDEFGHIJKLMNOP
Slide-07	A	A	ABDEFGHIJKLMNOP
Slide-08	A	A	ABDEFGHIJKLMNOP
Slide-09	A	A	ABDEFGHIJKLMNOP
Slide-10	A	A	ABDEFGHIJKLMNOP

Initial & Data Clear History Solution Result

Press **START** button. Cap opening, suction, blowing and slide transfer is automatically operated.

If you want to stop the operation, please push '**STOP**' button. Push the **START** button to continue with the operation. If you want to start from the beginning after stopping the operation, press **STOP** and then **HOME** to reset the position.


Check the result of smear and fix the slide.

When all smearing completes, **TOTAL VIAL PROCESSING RESULT** displays in the touch screen. If the work completed in normal, the vial displays with Green, otherwise it displays with Black.

Press the arrow on **TOTAL VIAL PROCESSING RESULT** or below of **SLIDE RESULT**, it moves to **GLASS PRECESSING RESULT** and the **SLIDE** work can be checked.


Take out the Smearing completed Slide that hold in the cell fixing alcohol container, and fix for 30 minutes in separate ethanol container then dye.

Manual Turbidity Mode


Please click on the **SETUP** button for 3 seconds or longer in the main screen.

The screen will switch to User Setup Screen.


Select Turbidity OFF.

Disable auto-turbidity selection function.


Vortex ON/OFF Selection Function.

Vortex should not be used for CSF specimens.


Prepare vials and load them in the tray.


Up to 40 vials can be loaded at once. The slide, completed a smear, is come out according to the order of Vial's work.

Put the vials from the front so no vacancy remains.

Push the tray by aligning the groove on the location of the next section when mounting the vial tray.

Select **A**, **B**, **C** on the touch screen according to the density and condition of prepared specimen.

Mode A:


- Used when amount of collected cells in solution vial are normal

Mode B:

- Used when too less amount of collected cells in solution vial

Mode C:


- Used when too concentrated amount of collected cells in solution vial

Select the types of specimens.


Set the mode according to the type of specimen by pressing the buttons on the touch screen.

Press the green "U" button for urinary and CSF samples.

Specify the quantity of vials for operating.


Number press menu comes out when press red number on **Working Count**.

Press **START** button.


Cap opening, suction, blowing and slide transfer is automatically operated.

If you want to stop the operation, please push 'STOP' button. Push the **START** button to continue with the operation. If you want to start from the beginning after stopping the operation, press **STOP** and then **HOME** to re-set the position.

Check the result of smear and fix the slide.


When all smearing completes, **TOTAL VIAL PROCESSING RESULT** displays in the touch screen. If the work completed in normal, the vial displays with Green, otherwise it displays with Black.

Press the arrow on **TOTAL VIAL PROCESSING RESULT** or below of **SLIDE RESULT**, it moves to **GLASS PRECESSING RESULT** and the SLIDE work can be checked.

Take out the Smearing completed Slide that hold in the cell fixing alcohol container, and fix for 30 minutes in separate ethanol container then dye.


Cell Prep PLUS: after turning on the instrument, select the sample type on the touch screen. Press "URINE" button for urine and CSF samples.

Before inserting the vial containing the sample into the processor, shake it well 5-10 times and place it upside down in the instrument with the vial cap facing down.

Place the vial in the holder. It is necessary to shake the vial several times before inserting it into the instrument to ensure that the cells are evenly distributed in the vial. Label and place the slide in the slide holder unit of the processor. The smeared side of the slide should be facing down.

Select A, B, C mode on the touch screen according to the amount of cells and conditions.

Cellprep

YYYY/MM/DD HH:MM:SS

STATUS		WORK 12345	
SUCTION		BLOWER	
00.1 kPa		1.23 kg/cm	
1 LINE IN	GYNE	A	12 CYCLE

Bio-dyne

GYNE URINE BODY

FNA RPT CUSTOM

A

B

C

MESSAGE : AUTO RUNNING (PLEASE PUSH START BUTTON)

MANUAL

SETUP

RELEASE

HOME

START

STOP

Mode A: Used when amount of collected cells in solution vial are normal

Mode B: Used when too less amount of collected cells in solution vial

Mode C: Used when too concentrated amount of collected cells in solution vial

Push 'START' button on the touch screen. Vial cap opening, suction, blowing and slide out is automatically operated.

Detailed description of how to operate the Cell Prep Plus LBC Processor with accessories is provided in the corresponding user manual.

After smearing insert the slide into 95% ethanol immediately and let it fix for 30 mins.

Screen out the lesion from the smeared slide after it is stained and mounted.

Interpretation of the results is carried out according to the clinical guidelines.

7.4. Interfering substances or sample limitations that may interfere with study result

Samples with a high content of blood and mucous, and samples with an insufficient content of material in all fields of view of the slide are considered inadequate, since they may cause subsequent difficulties in making a diagnosis. If abnormal cells are present in the sample, then any smear is considered adequate.


When collecting samples of the material, it is necessary to follow established clinical guidelines and sample preparation protocols that are appropriate for the type and cellularity of the sample. Failure to comply with this requirement may result in incorrect results.

7.5. Identification of the medical devices in order to obtain a safe combination and (or) information on known restrictions on the combined use of medical devices

The medical and/or non-medical devices listed in Section 7.1 may be used with Cellprep Urine, CSF Solution to ensure a safe combination.

It is NOT ALLOWED to use the Cellprep Solution with other products than those listed in Section 7.1.

8. Specifications of the Cellprep Urine, CSF Solution

Name of characteristic / parameter	Reference value
Appearance	The vial should not contain holes, cracks, deformations, scratches, and other defects that affect the quality of the device use. Appearance of solution - colorless transparent liquid without foreign impurities
Foreign impurities	
Fungi	No
Bacteria	No
Parasites	No
Others	1 unit / 10 ml
Reference volume	10 ml $\pm$ 10%
Torque effect of a vial middle cap	19 kgf·cm $\pm$ 5 kgf·cm
Torque effect of a vial top cap	14 kgf·cm $\pm$ 5 kgf·cm
pH of solution	pH=5.1 $\pm$ 0.25
Dimensions of vial:	
Total vial height:	89.0 mm $\pm$ 1.0 mm
Diameter of a vial middle cap:	42.7 mm $\pm$ 1.0 mm
Diameter of a vial middle part:	36.1 mm $\pm$ 1.0 mm
Diameter of a vial bottom part:	38.1 mm $\pm$ 1.0 mm
Weight of vial with solution	35.0 g $\pm$ 10%
Marking	There should be no soiling, poorly typed text. Lot number, expiration date, product name should correspond to the set parameters. The vial label must contain the following symbols and signs:  The shipping package with a cardboard box must contain the following symbols and signs: 

9. Recommendations for staining the slides

Step No.	Name of chemical for staining	Duration	
		Gyn	Non-Gyn
Fixation	95% Ethanol	30 minutes	30 minutes
1	80% Ethanol	10 dips	10 dips
2	70% Ethanol	10 dips	10 dips
3	Water	2 minutes	2 minutes
4	Hematoxylin	4 minutes	2 minutes
5	Water	1 minute 30 seconds	1 minute 30 seconds
6	0.5% HCl	3-5 dips	3-5 dips
7	Water	5 minutes	5 minutes
8	95% Ethanol	10 dips	10 dips
9	95% Ethanol	10 dips	10 dips
10	OG-6	3 minutes	1 minute 30 seconds
11	95% Ethanol	10 dips	10 dips
12	95% Ethanol	10 dips	10 dips
13	EA-50	3 minutes	1 minute 30 seconds
14	95% Ethanol	10 dips	10 dips
15	95% Ethanol	10 dips	10 dips
16	100% Ethanol	10 dips	10 dips
17	100% Ethanol	10 dips	10 dips
18	Xylene	10 dips	10 dips
19	Xylene	1 minute (10 dips)	1 minute (10 dips)

Desired reagent	Expected manufacture concentration	95% Ethanol (ml)	DW (ml)
50% Ethanol	50%	526	474
70% Ethanol	70%	737	263
80% Ethanol	80%	840	160
Desired reagent	Expected manufacture concentration	70% Ethanol (ml)	HCl (ml)
0.5% HCl-alcohol	0.5%	995	5

10. A. B. C. Suction Numerical Standard Values for LBC Processors

Pressure & Time

Specimen	Mode	Suction			Blower		
		Pressure		Time (sec)	Pressure		Time (sec)
URINE, CSF	A	20.0	40.0	1.0	0.05	0.05	0.5
	B	27.0	40.0	1.0	0.05	0.05	0.5
	C	15.0	35.0	1.0	0.05	0.05	0.5

11. Cautions and warnings

<Cautions>

1. Not for internal use.
2. Do not drink or inhale vapours.
3. Avoid media contact with eyes or skin.
4. When contacting with eyes, rinse eyes immediately and thoroughly with clean water and seek medical attention.
5. Colorless volatile liquid. In case of prolonged exposure, it affects the liver and central nervous system.
6. Use chemical resistant gloves, goggles and a lab coat/apron for your safety.
7. After use, the Cellprep solutions must be handled as a biohazard material.
8. Do not use the Cellprep Solution after the expiration date.

The Cellprep solutions should not be used in damaged or cracked vials, if the color of solution changes, or if foreign impurities are visible.

12. Transportation and storage conditions

Transportation by all means of transport is allowed.

Transportation and storage temperature: from +10°C to +30°C

Store in a closed container in an area that prevents direct contact with sunlight.

Store and use only in a well-ventilated area.

13. Packaging

The Cellprep solutions are delivered in polypropylene vials. Vials (25 pcs) are placed between two plastic vial trays and covered with polyethylene film. Eight packages (25 vials each) covered with film (200 vials in total), together with the Instructions for Use, are placed in a cardboard box. Then a cardboard box with several sealing non-combustible inserts is additionally placed in a cardboard package for hazardous materials and put on a pallet.

It is allowed to provide the Instructions for Use in electronic form by sending by e-mail at the consumer's request.

Basic information about the Cellprep Urine, CSF Solution packaging

Packaging	Material	Dimensions / quantity
Cellprep Urine, CSF Solution Vial The vial is a primary packaging for Cellprep solutions.	Polypropylene	Total vial height: 89.0 mm $\pm$ 1 mm Weight of vial with solution: 35.0 g $\pm$ 10%
Plastic vial trays covered with polyethylene film. Secondary consumer package	Plastic trays (polypropylene) Polyethylene film	Dimensions: 25.5 x 25.5 x 8.5 cm $\pm$ 10% Weight of packaging with vials (gross) - 1.06 kg $\pm$ 10% Quantity: 25 vials
Shipping package of Cellprep Urine, CSF Solution The shipping package includes an inner cardboard box and an outer fiberboard box for hazardous materials.	Cardboard box, inner	Dimensions: 525 x 270 x 370 mm $\pm$ 10% Weight of box with vials (gross): 9.40 kg $\pm$ 10% Quantity: 8 packages covered with film. 200 vials in total.
Package for hazardous materials. Shipping package	Outer fibreboard box for hazardous materials Several inserts from polymer non-combustible material or foam	Box dimensions: 530 x 320 x 450 mm $\pm$ 10% or 590 x 395 x 470 mm $\pm$ 10% Weight of box with vials (gross): 12.10 kg $\pm$ 10% Quantity: 1 inner cardboard box with 200 vials.
Pallet	Plastic or wooden pallet	Dimensions: 110 x 110 x 132 cm $\pm$ 10%

14. Symbol identification

Symbol	Description
--------	-------------

Causes serious eye irritation

Carcinogenic substances

Flammable liquid and vapours.

European Community
CE logo: This device complies with the EC Directives

Temperature range
(temp. +10 °C ~ +30 °C)

Biohazard

In vitro diagnostic medical device

Manufacturer's lot number.

Symbol	Description
--------	-------------

Best before...

Keep dry

Do not use hooks

Fragile, handle with care

Top

Manufacturer

Authorized representative
in the European
Community

Refer to the Instructions
for Use

15. Manufacturer's responsibility

Biodyne Co., Ltd. is a legal manufacturer responsible for the quality, efficiency and safety of the registered medical device for use, normal operation and possible negative results of using them, except for cases of negligence, incorrect use through the fault of user, or force majeure.

In case of any adverse reactions or events while using the product, please contact customer support service or your local distributor immediately.

16. Repair and Maintenance

The medical device does not require maintenance or repair.

For additional services or support, please contact Support Services or your local representative.

Support Service

Tel. +82-1522-4811

E-mail: info@biodyne.asia

Authorized representative in the Russian Federation:

BioLine, LLC

Legal address: Russia, 197022, Saint Petersburg, municipal district Aptekarsky Ostrov, 23 Professora Popova Str., Let. E, structure 3H/5H, room 208

Postal address: Russia, 197022, Saint Petersburg, 23 Professora Popova Str., Let. E

E-mail: main@bioline.ru

Tel.: +7(812) 320-49-49.

Hot line: 8-800-555-49-40.

17. Cleaning and Disinfection

Not applicable.

The Cellprep Urine, CSF Solution does not require cleaning and disinfection.

18. Warranty and shelf life

Biodyne Co., Ltd., Korea, warrants that the device has passed a comprehensive quality control procedure based on our internal testing standards and that the device is free from defects and conforms to all specifications and/or agreed specifications.

The manufacturer warrants the high quality of the device throughout the shelf life, subject to the recommended storage conditions.

Cellprep Urine, CSF Solution shelf life - 2 years (while observing recommended storage conditions).

The expiry date is indicated on the label.

Use of the Cellprep Solutions after the expiry date is prohibited.

The urine fixed in the vial, while keeping the temperature regime from +10 °C to +30 °C, can be used for liquid-based cytology studies for up to 3 months.°

The CSF fixed in the vial shall be kept at temperature from +10 °C to +30 °C no more than 5 hours.

19. Disposal

After use, the Cellprep solutions delivered by Biotdyne Co., Ltd., Korea, must be handled as biohazard waste.

Waste material should be disposed of in an approved incinerator in compliance with all federal, municipal and local government regulations.

For the Russian Federation, it is recommended to dispose of the device in accordance with the rules implemented by SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste referred to class B (epidemically hazardous waste).

20. Information on the initial release and revision of the Instructions for Use

Initial release. Version 01 from 09-2023.

21. List of the standards applied

a) Quality Management System applied by Biotdyne Co., Ltd., Republic of Korea, is certified to the requirements of EN ISO 13485:2016 by the TUV Rheinland LGA Products GmbH certification body in relation to: "Distribution and Service of Liquid-based cytology system. Design and Development, Manufacture and Distribution of Cell preservation solution and Membrane filters".

b) The manufacturer (Biotdyne Co., Ltd., Korea) complies with the principles of Good Manufacturing Practice (GMP) for medical devices. GMP Certificate is issued by Korea Conformity Laboratories.

The Cellprep Solutions comply with the following standards:

- Directive 98/79/EC on in vitro diagnostic medical devices
- EN 980:2008 Symbols for use in the labelling of medical devices
- EN ISO 14971:2012 Medical devices. Application of risk management to medical devices

Cellprep

28

Biodyne Co., Ltd.

24-4, Ahasan-ro 5-gil, Seongdong-gu, Seoul, 04793, Republic of Korea

Tel. +82-1522-4811 **Fax**+82-2-468-0910

E-mail: overseas@biodyne.asia

www.biodyne.asia

For the safe use of the device, please read the instructions for use carefully.
Photos in this manual may differ from actual device parts.

Appendix No. 1.

Certificate of Analysis - Cellprep Urine, CSF Solution (sample)

Cellprep Urine, CSF Solution for Liquid-based Cytology Processor						
Manufacturer: Biodyne Co., Ltd. 24-4, Ahasan-ro 5-gil, Seongdong-gu, Seoul, 04793, Republic of Korea Tel. +82-1522-4811 E-mail: overseas@biodyne.asia						
LOT No.		XXXX		Exp. date:		XXXXX
Certificate of Analysis No. XXXX dd. XXXXX						
No.	Parameter	Equipment	Standard	Compliance	Interval	Specification
1	External surface	Identification table Eyes	1. There should be no foreign impurities inside the vial.	Complies	At once after production	Observation from a distance of 30 mm
			2. The vial, middle cap and top cap should match.	Complies		
			3. Lot number must be indicated correctly.	Complies		
			4. The solution must not leak away.	Complies		
			5. There should be no cracks.	Complies		
			6. Solution and label should match.	Complies		
			7. There should be no foreign impurities on the vial exterior.	Complies		
			8. Check that the language on label and type are indicated correctly.	Complies		
2	Color	Eyes	1. Content should be transparent.	Complies		Visual inspection
3	Foreign impurities	Eyes	1. There should be no foreign impurities.	Complies		Visual assessment
4	Volume	Graduated cylinder	1. Volume of the solution should be determined by the volume for each type. (□ 10ml)	Complies		Measurement
5	pH of solution	pH-meter	5.1±0.25	Complies		Measurement
Transportation and storage temperature: from +10°C to +30°C. Store in a closed container in an area that prevents direct contact with sunlight. Store and use only in a well-ventilated area. Shelf life - 2 years (while observing recommended storage conditions).						

XXXX – variable data.

Cellprep LBC Solution

| 1. | CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

- A. PRODUCT NAME: Cellprep LBC Solutions (Urine CSF Solution)
- B. APPLICATION: Ethanol based, Medical preservation liquid for cell cultures.
- C. COMPANY IDENTIFICATION BIODYNE CO.,LTD.
24-4, Ahasan-ro 5-gil, Seongdong-gu,
Seoul, 04793, Republic of Korea
Phone: +82-1522-4811 / Fax: +82-70-4903-0901

| 2. | HAZARDS IDENTIFICATION

A. NATIONAL FIRE PROTECTION ASSOCIATION (NFPA) CLASSIFICATION

HEALTH HAZARDS = 0~1

FLAMMABILITY = 0

REACTIVITY = 0

SPECIFIC HAZARD = N/A

Rating: 0 = Insignificant 1 = Slight 2 = Moderate 3 = High 4 = Extreme

PRECAUTIONARY STATEMENTS

PREVENTION

- Obtain special instructions before use (P201).
Do not handle until all safety precautions have been read and understood (P202).
Keep container tightly closed (P233).
Ground/bond container and receiving equipment (P240).
Take precautionary measures against static discharge (P243).
Do not breathe dust/fume/gas/mist/vapors/spray (P260).
Wash eye thoroughly after handling (P264).
Do not eat, drink or smoke when using this product (P270).
Wear protective gloves/protective clothing/eye protection/face protection (P280).
Use personal protective equipment as required (P281).

RESPONSE

IF ON SKIN (or hair): Rinse skin with water/shower.

In case of fire: Use manufacturer/supplier or the competent authority to specify appropriate media for extinction. IF

IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

If eye irritation persists: Get medical advice/attention.

IF exposed or concerned: Get medical advice/attention.

IF exposed: Call a doctor/physician.

Get medical advice/attention if you feel unwell.

STORAGE

Store in cool/well-ventilated place.

Store locked up.

DISPOSAL

Dispose of contents/container in accordance with local/regional/national/international regulation.

| 3. | COMPOSITION / INFORMATION ON INGREDIENTS

Cellprep LBC Solution

CHEMICAL NAME	CAS No.	ENCS Number*	%(w/w)
Ethanol	64-17-5	2-202	42%
Di-water	7732-18-5		51%
Glycerine	56-81-5		6%
Phosphorous acid	7664-38-2		0.45%
Ammonium acetate	631-61-8		0.45%
Sodium Citrate Tribasic Dihydrate	6132-04-3		0.10%

* ENCS Number: Japanese Existing and New Chemical Substances.

4. FIRST AID MEASURES

A. EYE CONTACT

Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, get medical advice/attention.

B. SKIN CONTACT

Rinse skin with water/shower or wash with plenty of soap and water. If skin irritation or rash occurs, seek medical advice/attention. Wash contaminated clothing before reuse.

C. INHALATION

Not a likely route of exposure. If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing. If experiencing respiratory symptoms or feel unwell, call a doctor/physician.

D. INGESTION

Not a likely route of exposure. Do not induce vomiting without medical advice. There may be irritation to the gastrointestinal tract with nausea and vomiting.

E. MAJOR TOXIC SYMPTOMS AND EFFECTS

Inhalation of high concentration of vapor may cause mild irritation of the respiratory tract. Frequent or prolonged contact with product may defect and dry the skin, leading to discomfort and dermatitis.

F. NOTE TO PHYSICIAN

Prolonged or frequently repeated skin contact may aggravate an existing dermatitis condition in some individuals. Based on the individual reactions of the patient, the physician's judgment should be used to control symptoms and clinical condition.

5. FIRE FIGHTING MEASURES

A. SUITABLE EXTINGUISHING MEDIA

General extinguishing media is usable.

B. SPECIFIC HAZARDS ARISING FROM THE CHEMICAL

No cause of harmful hazards arose from the chemical.

C. SPECIAL PROTECTIVE EQUIPMENT AND PRECAUTIONS FOR FIRE FIGHTERS

In case of fire, wear a full face positive-pressure self contained breathing apparatus and protective suit. Keep containers cool with water spray.

6. ACCIDENTAL RELEASE MEASURES

A. PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES

Use personal protective equipment recommended in Section 8 (Exposure Controls/Personal Protection). Stop or reduce any leaks if it is safe to do so. Have emergency equipment (for fires, spills, leaks, etc.) readily available.

B. ENVIRONMENTAL PRECAUTIONS

Prevent material from entering sewers or waterways.

C. METHODS AND MATERIALS FOR CONTAINMENT AND CLEANING UP

Clean with towel or paper towel on its spill.

7. HANDLING AND STORAGE

A. PRECAUTIONS FOR SAFE HANDLING

Do not get in eyes, on skin, on clothing. Do not take internally. Use with adequate ventilation. Try not to breathe vapors/gases/dust. Keep the containers closed when not in use.

B. CONDITIONS FOR SAFE STORAGE, INCLUDING ANY INCOMPATIBILITIES

Store in suitable labeled containers. Store the containers tightly closed. Store away from heat. Keep containers placed in cool, well-ventilated areas. Have appropriate fire extinguishers available in and near the storage area.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

A. OCCUPATIONAL EXPOSURE LIMITS

KOREA

Korea. OELs (ISHL Article 42; MOL Public Notice No. 1986-45, as amended through MOL Public Notice No. 2007-25, June 8, 2007)

CAS RN: 64-17-5

Name: Ethanol

The 8-hour TWA: 1,000 ppm

The 8-hour TWA: 1,900 mg/m³

The 15-minute STEL: - ppm

The 15-minute STEL: - mg/m³

JAPAN

Japan. OELs - JSOH (Japan Society of Occupational Health: Recommendation of Occupational Exposure Limits, 2007)

CAS RN: 64-17-5

Name: Ethanol

Not established.

OSHA

OSHA Table Z-1 Limits for Air Contaminants (June 30, 1993) (29 CFR 1910.1000) (1971 Permissible Exposure Limits (PELs))

CAS RN: 64-17-5

Name: Ethanol

The 8-hour PEL: 1,000 ppm

The 8-hour PEL: 1,900 mg/m³

NIOSH

NIOSH. Pocket Guide to Chemical Hazards, 2005

CAS RN: 64-17-5

Name: Ethanol

NIOSH Recommended exposure limit (REL) [for up to a 10-hour workday during a 40-hour workweek]: 1,000 ppm.

NIOSH Recommended exposure limit (REL) [for up to a 10-hour workday during a 40-hour workweek]: 1,900 mg/m³.

NIOSH Immediately dangerous to life or health (IDLH) concentration: 3,300 ppm

B. ENGINEERING MEASURES

General ventilation is recommended. Use local exhaust ventilation if necessary to control airborne mist and vapor. Provide mechanical ventilation of confined spaces.

C. INDIVIDUAL PROTECTION MEASURES

GENERAL ADVICE

The use and choice of personal protection equipment is recommended when use of the product, the workplace and the way the product is handled. In general, we recommend as a minimum precaution as wearing gloves.

Cellprep LBC Solution

RESPIRATORY PROTECTION

If significant mists, vapors or aerosols are generated an approved respirator is recommended. A suitable filter material depends on the amount and type of chemicals being handled. If respiratory protection is required, institute a complete respiratory protection program including selection, fit testing, training, maintenance and inspection.

EYE PROTECTION

When handling this product, the use of splash chemical goggles can be used when necessary.

HAND PROTECTION

Wearing appropriate chemical resistant gloves are recommended.
Regular glove replacement is recommended.

SKIN PROTECTION

No specific harm on skin. Wearing appropriate chemical resistant cloths are recommended.

| 9. | PHYSICAL AND CHEMICAL PROPERTIES

A. PHYSICAL APPEARANCE:	Clear, colorless liquid.
B. ODOR:	Mild Characteristic odor..
C. ODOR THRESHOLD:	No data available.
D. pH:	No data available.
E. MELTING/FREEZING POINT:	No data available.
F. INITIAL BOILING POINT/BOILING POINT RANGE:	No data available.
G. FLASH POINT:	About 30°C *
H. EVAPORATION RATE:	No data available.
I. FLAMMABILITY (GAS, SOLID):	Not applicable.
J. UPPER AND LOWER EXPLOSION LIMITS:	No data available
K. VAPOUR PRESSURE:	No data available.
L. SOLUBILITY:	Soluble in water.
M. RELATIVE VAPOUR DENSITY:	No data available.
N. SPECIFIC GRAVITY:	No data available.
O. n-OCTANOL/WATER PARTITION COEFFICIENT:	No data available.
P. AUTOIGNITION TEMPERATURE:	No data available.
Q. DECOMPOSITION TEMPERATURE:	No data available.
R. VISCOSITY:	No data available.

* Sample was not flashed until boiled (about 30°C)

| 10. | STABILITY AND REACTIVITY

A. CHEMICAL STABILITY

Stable under normal conditions.

B. POSSIBILITY OF HAZARDOUS REACTIONS Hazardous polymerization will not occur.

C. CONDITIONS TO AVOID

No specific conditions to avoid.

D. INCOMPATIBLE MATERIALS

No materials to be incompatible. _____

| 11. | TOXICOLOGICAL INFORMATION

A. INFORMATION ON THE LIKELY ROUTES OF EXPOSURE

EYE CONTACT: Yes.

SKIN CONTACT: Yes.

INGESTION: No.

INHALATION: No.

B. SYMPTOMS RELATED TO THE PHYSICAL, CHEMICAL AND TOXICOLOGICAL CHARACTERISTICS

No data available.

ACUTE TOXICITY:

- Oral: None known. No data available.
- Skin: None known. No data available.
- Inhalation: None known. No data available.

SKIN CORROSION/IRRITATION:

Can cause mild irritation to skin. Frequent or prolonged contact with product may defect and dry the skin, leading to discomfort and dermatitis.

SERIOUS EYE DAMAGE/IRRITATION:

Can cause mild irritation to eye.

RESPIRATORY OR SKIN SENSITIZATION:

This product is not expected to be a sensitizer.

CARCINOGENICITY:

ACGIH Carcinogen Category: A4 (not classifiable as to carcinogenicity in humans).

None of the substances in this product are listed as carcinogens in humans by the International Agency for Research on Cancer (IARC), American Conference of Governmental Industrial Hygienists (ACGIH).

ASPIRATION HAZARD:

None known. No data available.

| 12 | **ECOLOGICAL INFORMATION**

* No toxicity studies have been conducted on this product. The following results are for the hazardous components.

A. AQUATIC AND ECOLOGIC TOXICITY

Occurrence of ecologic toxicity is possible when exposed or deposited in soil for long period of time.

B. PERSISTENCE AND DEGRADABILITY

The organic portion of the substances is expected to be biodegrade to a moderate extent.

C. BIOACCUMULATIVE POTENTIAL

The product not expected to significantly bioaccumulation due to low bioaccumulation factor (BCF) and octanol/water partition coefficient.

D. MOBILITY IN SOIL

None known. No data available.

E. OTHER ADVERSE EFFECTS

None known. No data available.

| 13. | **DISPOSAL CONSIDERATIONS**

A. DISPOSAL METHODS

Cellprep LBC Solution

Hazardous wastes must be transported by a licensed hazardous waste transporter and disposed of or treated in a properly licensed hazardous waste treatment, storage, disposal or recycling facility.
Consult local, state, and federal regulations for specific requirements.

B. PRECAUTIONS (ON DISPOSAL OF COTAMINATED CONTAINERS AND PACKAGES)

Do not dispose of wastes in local sewer or with normal garbage. Combustion products are carbon monoxide, carbon dioxide and water.

| 14. | TRANSPORT INFORMATION

- A. UN No.: 1170.
- B. Proper Shipping Name: ETHANOL SOLUTION
- C. CLASS: 3
- D. PACKING GROUP. III
- E. Not restricted under IATA regulations

| 15. | REGULATORY INFORMATION**INTERNATIONAL CHEMICAL CONTROL LAWS****EUROPE**

The disclosed substances in this product are included in or exempted from the EINECS or ELINCS inventories.
U.S.

The disclosed substances in this product are included on or exempted from the TSCA 8(b) Inventory (40 CFR 710).

JAPAN

The disclosed substances in this product are comply with the Law Regulating the Manufacture and Importation of Chemical Substances and are listed on the Ministry of Economy, Trade and Industry (METI).

KOREA

Observational Chemicals List - None of the components of this product are regulated under TCCL.

| 16. | OTHER INFORMATION

This product material safety data sheet provides health and safety information. The above information is believed to be correct but does not propose to be all inclusive and shall be used only as a guide. BIODYNE CO.,LTD. shall not be held liable for any damage resulting from handling or from contact with the above product. Each individual should make a determination as to the suitability of the information for their particular purpose(s). Various government agencies may have specific regulations regarding the transportation, handling, storage, use, or disposal of this product which may not be covered by this MSDS. The user is responsible for full compliance.

A. REFERENCES

Regulated chemicals info. produced by the National Institute of Technology and Evaluation, Japan.

Chemicals Evaluation and Research Institute, Japan (CERI)

WHO/IPCS : International Chemical Safety Cards (ICSC).

EU European Chemicals Bureau (ECB): International Uniform Chemical Information Database (IUCLID). Ariel

WebInsight DB (3E Company).

Korea Testing & Research Institute

B. OTHER INFORMATION

- Date issued: April 01, 2024
- Version number: 6.3
- Technical Support and Revised from Attached Research Center of Biodyne Co., Ltd.

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As a result of checking at my office, I have
found that the attached INSTRUCTIONS
FOR USE copy exactly corresponds with
the original.

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This is hereby attested on this 13th day of
May, 2024 at this office.

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Name of the office
HANMI LAW AND NOTARY OFFICE

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
Belong to
Seoul Central District Prosecutors' Service

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
Address of the office
Twintree Tower(B) #101,6 Yulgok-ro,
Jongno-gu, Seoul, Republic of Korea

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Notary Public HAN MYUNG SUB

(Signature of Notary Public)

This office has been authorized by the
Minister of Justice, the Republic of Korea,
to act as Notary Public since March 14,
2007 under Law No. 7428.

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