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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of LIFETECH SCIENTIFIC
(SHENZHEN) CO., LTD. on the annexed STATEMENT is genuine.

China Council for the Promotion
of International Trade

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17 January, 2019

We

Lifetech Scientific (Shenzhen) Co., Ltd
Floor 1-5, Cybio Electronic Building, Langshan 2nd Street,
North Area of High-tech Park, Nanshan District,
Shenzhen, 518057
P.R. China

Here state that:

The document only used for product registration in Russia and not for other use. All the contents are true and lawful.

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
CeraFlex™ PDA Closure System, in different versions

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1 Device Description

The CeraFlex™ Patent Ductus Arteriosus (PDA) Closure Systems including CeraFlex™ PDA Occluder and SteerEase™ Introducer. The CeraFlex™ PDA Occluder is a percutaneous transcatheter occluding device for the non-surgical closure of patent ductus arteriosus defects. The device consists of a knitted nitinol wire mesh, which has been heated and set to expand into a T-shape. The disc on the aorta side is designed to hold the plug at the orifice of the ductus arteriosus, while the waist expands to fill and occlude it. A PTFE membrane is sewn inside the disc with nylon thread on the aorta side, and two PTFE membranes are sewn inside the waist to help block the abnormal blood flow from the aorta into the pulmonary artery. All metallic structures are plated with a titanium nitride (TiN) bioceramic coating (applied by vacuum plasma deposition technology), to effectively reduce the release of nickel ions into blood and surrounding tissues in order to mitigate thrombogenicity of the device and to speed up endothelial cell growth on the surface of device.


Figure 1. The CeraFlex™ PDA Occluder

The CeraFlex™ PDA Closure System consists of PDA occluder, delivery cable, loader, Haemostatic valve, sheath and dilator. The system (Refer to Figure 1 and Figure 2) is comprised of:

- **PDA Occluder:** A device for closing PDA. The CeraFlex™ PDA Occluder is connected with delivery cable by a loop connection through the holes. The CeraFlex™ PDA Occluder is preassembled with delivery cable.
- **Sheath:** The sheath is used to advance the device to its appropriate position.
- **Dilator:** The dilator is used to ease penetration of tissues and the vessel wall.
- **Haemostatic valve:** The haemostatic valve at the proximal end of the sheath minimizes bleeding. The sideport with the flexible extension tube and stopcock are used to flush the system.
- **Loader:** The loader is used to introduce the occluder with the attached delivery cable into the sheath.
- **Delivery cable with Plastic vise/Handle (proximal):** The delivery cable is used to advance (push) the occluder through the sheath. The handle and delivery cable are also used to retrieve and/or reposition the occluder, if the size, position or the expansion of the occluder is deemed to be unsatisfactory. The handle is used to facilitate direction control and serve as "release control handle" for disconnecting (releasing) the occluder from the delivery cable.


Figure 2. CeraFlex™ PDA Occluder and the SteerEase™ Introducer

2 Indications and Usage

The CeraFlex™ PDA Occluder is a percutaneous, transcatheter closure device that is intended for the non-surgical closure of patent ductus arteriosus (PDA).

Indications:

- Patients have PDA or PDA accompanies with other mild cardiac disease;
- Patients' weight more than 6 kgs, and 6 months old and above,
- The narrowest portion of the PDA is 2 mm or larger.

Contraindications:

- Patients have PDA which is needed for patients' survival due to other cardiac anomalies;
- Patients' weight less than 6 kgs, or less than 6 months old,
- Presence of thrombus at the intended site of implant, or documented evidence of venous thrombus in the vessels through which access to the defect is gained,
- Active endocarditic or other infections producing bacteremia.

English 1

- Patients vasculature, through which access to the defect, is inadequate to accommodate the appropriate sheath size;
- Patients with pulmonary hypertension higher than 0.4 R_p/R_s or with pulmonary vascular resistance more than 8 woods units

4 Warnings:

- The device should be removed if more than half of the left pulmonary artery lumen is obstructed by the device, or if more than 3 mm length extends into the pulmonary artery
- Patients allergic to nickel may suffer an allergic reaction to the device.
- The PDA Closure System must only be used in hospitals, which have experience with invasive treatment of congenital heart disease, and used by physicians who are trained in transcatheter defect closure techniques.
- Physicians must be prepared to deal with emergency cases which require removal of embolized devices that result in a critical hemodynamic compromise. This includes the availability of an on-site surgeon
- Embolized devices must be withdrawn via the delivery track within a sheath.
- The PDA Closure System is for single use only. Do not reuse or restenilize. Structural integrity and/or function may be impaired or be lost through cleaning, restenilization, or reuse and may cause adverse patient reactions. Lifetech will not be responsible for any direct or consequential damages or expenses resulting from reuse of any of the components in the PDA Closure System
- The PDA Occluder should not be released from the delivery cable if the device does not match with the defect or if the device position is unstable. Withdraw the device into the introduce sheath and redeploy it if positionings still unsatisfactory, withdraw the device completely and replace with a new one.
- Ensure the connection safety before loading the device into the loader
- Retrieve the device when it is hard to release the device from the delivery cable.

5 Precautions:

- Both the PDA Occluder and the SteerEase™ Introducer are for single use only, and the valid period of sterilization is indicated on the labels. Do not reuse or restenilize
- Post-procedure
 - Patients should take appropriate endocarditis prophylaxis for 6 months following the device implantation. The decision to continue endocarditis prophylaxis beyond 6 months can be made by physicians
 - Lung perfusion scan should be done if flow speeds is greater than 3m/s, or if the Z-score is -2 for the left pulmonary artery diameter
- MR conditional
A patient with an implanted PDA Occluders can be scanned safely immediately after placement of the device under the following conditions:
 - Static magnetic field of 3 T or less
 - Spatial gradient magnetic field of 720 G/cm or less
 - Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning

Note: MRI image quality may be compromised if the scan area is in the same area or relatively close to the position of the device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device is necessary

- Use in specific populations
 - Pregnancy – Care must be taken to minimize the radiation exposure to the fetus and the mother.
 - Nursing mothers – Although appropriate biocompatibility testing has been conducted for this implant device, there has been no quantitative assessment of the presence of leachables in breast milk.

6 Potential Adverse Events

Placement of the CeraFlex™ PDA Occluders involves using standard interventional cardiac catheterization techniques. The following adverse events might be expected from interventional cardiac catheterization techniques.

- | | |
|--------------------------------------|-----------------------|
| ◊ Allergic reaction | ◊ Arrhythmias |
| ◊ Cardiac tamponade | ◊ Thrombus formation |
| ◊ Septal tearing | ◊ Device migration |
| ◊ Injury to the nerve or vessel wall | ◊ Cardiac perforation |
| ◊ Incomplete sealing of the defect | ◊ Infection |
| ◊ Pulmonary artery dissection | ◊ Fever |

2English

7 Product Features

The CeraFlex™ PDA Closure system is available with specifications in the following table (Table 1).

Table 1. List of the available PDA Occluders and recommended introducers

Spec.	ΦD1mm	ΦD2mm	ΦD1mm	L/mm	Recommended Introducers
CeraFlex™ PDA Occluders					SteerEase™ Introducers
LT-PDAI-0406	10	6	4	4	SFP6F-I
LT-PDAI-0608	12	8	6	6	SFP7F-I
LT-PDAI-0810	14	10	8	7	SFP7F-I
LT-PDAI-1012	16	12	10	7	SFP8F-I
LT-PDAI-1214	20	14	12	7	SFP9F-I
LT-PDAI-1416	22	16	14	8	SFP9F-I
LT-PDAI-1618	24	18	16	8	SFP10F-I
LT-PDAI-1820	26	20	18	9	SFP12F-I
LT-PDAI-2022	28	22	20	9	SFP12F-I
LT-PDAI-2224	30	24	22	10	SFP14F-I


Figure 3. CeraFlex™ PDA Occluder

8 Directions for Use

- Perform a right heart catheterization in a routine way.
- Use angiography to size the heart defect. There are two options for angiography performed of the patent ductus arteriosus. The first is to introduce an exchange guidewire through the ductus and pass a pigtail catheter with side holes in the communication. Perform a biplane angiogram to opacify the PDA. The second option is to pass a pigtail catheter into the proximal descending aorta via the femoral artery and perform the biplane angiogram to opacify the PDA. Under the angiography, measure ductal size both on the pulmonary and aortic ends of the duct.
- Select the proper size occluder. Select a CeraFlex™ PDA Occluder based on the smallest diameter measured in the PDA. It is recommended to select a device of which the smaller end is at least 2 mm larger than the narrowest portion of the PDA.
- Introduce a 0.038-inch J-tipped exchange guide wire.
- Flush the dilator, sheath with saline solution, and then insert the dilator into the sheath, tighten the clamp nut to ensure the dilator fasten with sheath securely. Advance the introducer sheath with dilator over the exchange wire via the pulmonary artery into the aorta and position the sheath in the descending aorta while withdrawing the dilator.
- Clockwise screw the loader into the Haemostatic valve until it stops, and then check connection is reliable.
- Infuse the saline solution from the stopcock to flush the haemostatic valve and the loader, while keeping the tip of the loader immersed under the saline solution.
- Carefully check the connection between the occluder and the delivery cable (Refer to Figure 4). Immerse the device and the loader in saline solution and pull the CeraFlex™ PDA Occluder into the loader with a quick jerking motion (Refer to Figure 5).
- Twist the seal cover of the haemostatic valve to suitable degree (Refer to Figure 6), and then repeatedly push and pull to make sure the tension is appropriate.

Warning:

If the seal cover is too tight, then the cable will be hard to push. However, if the seal cover is too loose, there will be a risk of aeroembolism.

- Inject sterile heparinized saline from stopcock of the loader to eliminate air from the loader.

English3


Figure 4. Illustration of CeraFlex™ PDA Occluder System

Figure 5. Illustration of loading the device into a loader


Figure 6. Illustration of clockwise screwing the seal cover of Haemostatic valve

- Keep the sterile heparinized saline outflowing from the loader to eliminate air of the loader. Simultaneously, insert the loader into the Haemostatic Valve of the sheath gently until stop, and then tighten the clamp nut of loader. Advance the occluder into the descending aorta by delivery cable without any rotation. Do not remove the loader from the sheath when advancing the device.

Warning:

If there is air in the loader and sheath when advancing the device, it will result in an air embolism. If the loader is removed from the sheath when advancing the device, there will also be risk of an air embolism.

- Deploy the retention skirt only and pull firmly against the orifice of the PDA. This can be observed under fluoroscopy, or it can be clearly felt as a tugging sensation in synchrony with the aortic pulsation. The position of the device is confirmed with repeated angiograms in the aorta using the pigtail catheter. The device can be adjusted until the retention skirt is well seated at the orifice. Retract the delivery sheath and deploy the cylindrical portion of the device securely in the patent ductus arteriosus while applying slight tension.
- Perform an aortography to verify the correct position of the device. If the position is unsatisfactory, stabilize the sheath and pull the delivery cable until the device is completely within the sheath. It can be repositioned and deployed, or removed completely from the patient.

Warning:

When withdrawing the occluder back into the sheath, pull back on the delivery cable, and do not push the sheath. Otherwise, the tip of the occluder may be stuck outside the sheath, and then the occluder cannot be resheathed.

- Release the device when the position is satisfactory. To release the device, press the blue button and move back, while maintaining the orange button press state.
- Hold the sheath to touch the hub of the occluder, and then pull back the cable gently until it separates from the occluder. Finally remove the delivery cable and sheath from the patient. (Refer to Figure 7)
- Repeat aortography.


Figure 7. Illustration of releasing the occluder

- 1 Pressing the lock button as the arrow of 1 (the color of the button is orange);
- 2 Moving the control button as the arrow of 2 (the color of the button is blue)

Warning:

4English

Remove the device if over 3 mm length extends into the pulmonary artery, or if more than half of the left pulmonary artery lumen is obstructed by the device. In questionable cases, perform TTE before releasing the device with Doppler measurement of left pulmonary artery flow velocity. The device should be removed if the left pulmonary artery flow is greater than 3.0 m/s (or greater than 75% of the LPA velocity before cardiac catheterization).

Note

The device should not be released from the delivery cable if the device does not match the defect or if device position is unstable. Withdraw the device into the sheath and redeploy. If still unsatisfactory, withdraw the device completely and replace it with a new one.

If the CeraFlex™ occluder becomes detached from the delivery cable during loading, it cannot be used and must be discarded.

9 Post-procedure

- All patients should be kept overnight for observation
- Antibiotic therapy
- Re-examination by TTE at 24 hours prior to discharge

10 Expiry Date

Both the CeraFlex™ PDA Occluder and the SteerEase™ introducer are sterilized with ethylene oxide. The shelf life is indicated on the label. The manufacturing date and use by date are marked on the label. Do not use an expired product.

11 Package and Label

The CeraFlex™ PDA Closure Systems are supplied sterile. The CeraFlex™ PDA Occluder is pre-connected with delivery cable by loop connection and pre-assembled with its Loader and Haemostatic valve. The sheath and dilator are packaged separately. The package content is shown in Table 2.

Table 2. The packaging contents of the PDA Closure Systems

Package	CeraFlex™ PDA Closure Systems	
	CeraFlex™ PDA Occluders	SteerEase™ Introducers
Contents	a PDA Occluder, a Delivery Cable, a Loader, and a Haemostatic Valve	a Sheath and a Dilator

The PDA device is fixed and protected on a PETG tray, then sealed into inner and outer Tyvek/linpouches, where a label and a sterilization indicator are attached. The product is sterilized by ethylene oxide gas and put in a box with the Instructions for Use, a patient card, a customer feedback form and a certificate of compliance. An outer label is attached to the box.

The sheath and dilator of the SteerEase™ Introducer are fixed and protected on a scale board, then sealed in two dialyzing pouches, where a label and a sterilization indicator are attached. The product is sterilized and put in a box with the Instructions for Use, a customer feedback form, and a certificate of compliance. A label is attached to the box.

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[illegible]

4. **Figure 1** shows the following:
- A-axis:** the length of the longest projection of the polygonal skeleton (TD) measured over the vertices of the polygonal skeleton.
 - B-axis:** the length of the longest projection of the skeleton (TD) measured over the midpoints of the edges of the polygonal skeleton.
 - Perimeter:** the length of the polygonal skeleton.
 - Area:** the area of the polygonal skeleton.

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of LIFETECH SCIENTIFIC
(SHENZHEN) CO., LTD. on the annexed DOCUMENT is genuine.

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证书验证网址: Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Lifetech Scientific (Shenzhen) Co., Ltd.,
Floor 1-5, Cybio Electronic Building,
Langshan 2nd Street, North Area
of High-tech Park, Nanshan District,
Shenzhen, 518057, P.R. China

Ling Li/Regulatory Affairs Supervisor

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Scientific Corporation)
Scientific India Private Limited)