AEGIS LIFESCIENCES

PROTECTION THROUGH PERFECTION

Authorized representative of the manufacturer in
Russian Federation

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General Manager of MRPC Avis Life Science LLC/

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INSTRUCTIONS FOR USE OF MEDICAL DEVICE

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INSTRUCTIONS FOR USE
Absorbable Haemostat (Oxidized Regenerated Cellulose)
Surgi-ORC

MEDICAL DEVICE NAME

Surgi-ORC surgical hemostatic absorbable material (oxidized regenerated cellulose):

I. Surgi-ORC surgical hemostatic absorbable material (oxidized regenerated cellulose), in the required design variants – 12 pcs.

Design variants:

1) Surgi-ORC Original/Standard: SOO-502; SOO-66; SOO-0101; SOO-0102; SOO-0203; SOO-0205; SOO-0214; SOO-0304; SOO-0408;

2) Surgi-ORC Knit: SOK-0101; SOK-0102; SOK-0103; SOK-0135; SOK-0203; SOK-0204; SOK-0304; SOK-0404; SOK-5608; SOK-0609;

3) Surgi-ORC Fibril: SOF-0101; SOF-0102; SOF-0202; SOF-0203; SOF-0204; SOF-0304; SOF-0404;

4) Surgi-ORC Non-Woven/SNOW: SON-0102; SON-0202; SON-0203; SON-0204; SON-0304; SON 0404.

II. Instructions for use – 1 pc.

The potential users of the medical device are surgeons, nursing staff.

DESCRIPTION

Surgi-ORC material is the absorbable (biodegradable) tissue obtained by controlled chemical oxidation of regenerated cellulose. The product is sterile. The fabric is white pale yellow cast and has a faint, caramel-like aroma. The product retains its properties for the specified shelf life at room temperature. However, there may be mild discolouration over time, but this does not affect the efficacy.

Surgi-ORC material is available in ORIGINAL/STANDARD, KNIT, FIBRIL and NON-WOVEN/SNOW variants (models), made of the same oxidized regenerated cellulose.

Surgi-ORC material meets all the international requirements for oxidized regenerated (recovered) cellulose.

Parameter	Acceptance criterion
Identification	a precipitate in the form of white flakes should appear
Loss on drying	NMT 15.0 %
PH index	1.0...3.0
Residue on ignition	NMT 0.15 %
Limit of Nitrogen	NMT 0,5 %
Limit of Formaldehyde	NMT 0,5 %
Cellulose quantification	18.0% — 24.0 % of Carboxyl group (COOH)
Sterility	negative microbial growth

Parameter	Surgi-ORC Original/Standard	Surgi-ORC Knit	Surgi-ORC Fibril	Surgi-ORC Non-Woven/SNOW
Colour	Light yellow	Light yellow	Light yellow	Light yellow
Odour	Faint, caramel-like aroma	Faint, caramel-like aroma	Faint, caramel-like aroma	Faint, caramel-like aroma
Appearance	Pale yellow gauze napkin (knitted mesh).	Pale yellow gauze napkin (knitted mesh).	Light yellow cotton drape (fibrous material)	Light yellow cotton drape (fibrous material)
Breaking Strength	(50...70) N	(65...85) N	N/A	N/A
Elongation	(1,5...2,5) cm	(1,0...2,0) cm	N/A	N/A

INTENDED USE

Surgi-ORC material is used adjunctively in surgical procedures to assist in the control of capillary, venous and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective. Surgi-ORC material can be cut to size for use in endoscopic procedures.

OPERATING PRINCIPLE

Surgi-ORC absorbs the human blood and forms a brown or black gelatinous mass. This gelatinous mass acts as a physical matrix for the platelets adherence. Platelets aggregation leads to the formation of a platelet-fibrin plug, the swelling material applies pressure on the surrounding tissues, promoting haemostasis.

Anti-Bacterial effect of the hemostatic agents made of Oxidized Regenerated Cellulose (ORC) are accounted for by the characteristics of the oxidized cellulose itself. Low pH of oxidized cellulose, which can inhibit the growth of microorganisms, leads to the bactericidal effect.

In addition to its local haemostatic properties, Surgi-ORC material exhibits in vitro bactericidal properties against a wide range of Gram positive and Gram negative organisms including aerobes and anaerobes, such as:

- methicillin-resistant *Staphylococcus aureus* (MRSA)
- penicillin-resistant *Streptococcus pneumoniae* (PRSP)
- vancomycin-resistant *Enterococcus* (VRE)
- methicillin-resistant *Staphylococcus epidermidis* (MRSE)
- Staphylococcus aureus* *Bacillus subtilis* *Staphylococcus epidermidis* *Proteus vulgaris*
- Micrococcus luteus* *Corynebacterium xerosis*
- Streptococcus pyogenes* Group A *Mycobacterium phlei*
- Streptococcus pyogenes* Group B *Clostridium tetani*
- Streptococcus salivarius* *Clostridium perfringens*
- Branhamella catarrhalis* *Bacteroides fragilis*
- Escherichia coli* *Enterococcus*
- Klebsiella aerogenes* *Enterobacter cloacae*
- Lactobacillus* sp. *Pseudomonas aeruginosa*
- Salmonella enteritidis* *Pseudomonas stutzeri*
- Shigella dysenteriae* *Proteus mirabilis*
- Serratia marcescens*

However, Surgi-ORC hemostat is not intended as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent post-operative infections.

No adjuvants or additives are added to the product. Surgi-ORC material consists entirely of oxidized regenerated cellulose, and thus it is a medical device.

For risk management purposes, the instructions for use warn against the use of Surgi-ORC as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents.

Thus, the haemostatic effect of Surgi-ORC is based on its swelling into a gelatinous mass when saturated with blood, resulting in the formation of a stable clot. The action mechanism of Surgi-ORC is independent from ® the blood coagulation mechanism of the body. The Surgi-ORC material does not achieve its intended action by pharmacological, immunologic or metabolic means.

INDICATIONS FOR USE

All four Surgi-ORC material options — ORIGINAL/STANDARD, KNIT, FIBRIL and NON-WOVEN/SNOW — consist of oxidized reduced cellulose (ORC) and are intended to be used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial hemorrhage when ligation or other conventional methods of control are impractical or ineffective. Options are designed with different weave patterns to control bleeding depending on the type of surgery and the surface of the bleeding site.

Surgi-ORC Original/Standard	Composition — 100% oxidized regenerated cellulose (ORC) Characteristics — Fine weave of ORC, high flexibility. Indications: Surgi-ORC Original/Standard is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective. Haemostat can be cut to size for use in endoscopic procedures. Surgi-ORC Original/Standard can be used in many in many areas of surgery. e.g., cardiovascular surgery, haemorrhoidectomy, implantation of vascular prostheses, biopsies, lung operations, surgery to the face and jaw, gastric resection, operations on the throat or nose, liver and gallbladder operations, gynecological operations, thoracic and abdominal sympathectomies, neurosurgery, especially cerebral operations, thyroid operations, skin transplantations, and treatment of superficial injuries. Surgi-ORC Original/Standart is indicated also for adjunctive use in dental application to assist in the control of bleeding in exodontia and oral surgery. It may also be used to help achieve haemostasis after single or multiple teeth extraction, alveoloplasty, gingival haemorrhage, impaction, biopsies, and other procedures in the oral cavity. Dimensions: REF SOO-502 – 0.5 x 2.0 inches (1.3 x 5.1 cm); REF SOO-66 – 0.6 x 0.6 inches (1.5 x 1.5 cm); REF SOO-0101 – 1 x 1 inch (2.5 x 2.5 cm); REF SOO 0102 – 1 x 2 inches (2.5 x 5.1 cm); REF SOO-0203 – 2 x 3 inches (5.1 x 7.6 cm); REF SOO-0205 – 2 x 5 inches (5.1 x 12.7 cm); REF SOO-0214 – 2 x 14 inches (5.1 x 35.6 cm); REF SOO-0304 – 3 x 4 inches (7.6 x 10.1 cm); REF SOO-0408 – 4 x 8 inches (10.2 x 20.3 cm).
Surgi-ORC Knit	Composition — 100% oxidized regenerated cellulose (ORC)

	Characteristics – Denser weave/densely woven knit of ORC (compared to the ORIGINAL/STANDARD), Soft; pliable weave structure designed to hold a suture and for delicate tissue; provides matrix for platelet adhesion and aggregation for heavier bleeding. Indications: Surgi-ORC Knit is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective. Haemostat can be cut to size for use in endoscopic procedures. Surgi-ORC Knit can be used in many in many areas of surgery. e.g., cardiovascular surgery, haemorrhoidectomy, implantation of vascular prostheses, biopsies, lung operations, surgery to the face and jaw, gastric resection, operations on the throat or nose, liver and gallbladder operations, gynecological operations, thoracic and abdominal sympathectomies, neurosurgery, especially cerebral operations, thyroid operations, skin transplantations, and treatment of superficial injuries. Surgi-ORC Knit is indicated also for adjunctive use in dental application to assist in the control of bleeding in exodontia and oral surgery. It may also be used to help achieve haemostasis after single or multiple teeth extraction, alveoloplasty, gingival haemorrhage, impaction, biopsies, and other procedures in the oral cavity. Dimensions: REF SOK-0101 – 1 x 1 inch (2.5 x 2.5 cm); REF SOK-0102 – 1 x 2 inches (2.5 x 5.1 cm); REF SOK-0103 – 1 x 3 inches (2.5 x 7.6 cm); REF SOK-0135 – 1 x 3.5 inches (2.5 x 8.8 cm); REF SOK-0203 – 2 x 3 inches (5.1 x 7.6 cm); REF SOK-0204 – 2 x 4 inches (5.1 x 10.2 cm); REF SOK-0304 – 3 x 4 inches (7.6 x 10.2 cm); REF SOK-0404 – 4 x 4 inches (10.2 x 10.2 cm); REF SOK-5608 – 5.6 x 8 inches (14.2 x 20.3 cm); REF SOK-0609 – 6 x 9 inches (15.2 x 22.9 cm)
Surgi-ORC Fibril	Composition — 100% oxidized regenerated cellulose (ORC) Characteristics — Lightweight, tufted, soft and layered structure. 7 layers of fibrous texture. Can be easily peeled off in any thickness and quantity required. It provides high flexibility and conformability to surgical sites. It is suitable for irregular shaped wounds and hard-to-reach areas. Intended for use directly for hemostasis. Indications: Material is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective. Surgi-ORC Fibril can be cut to size for use in endoscopic procedure.

		Surgi-ORC Fibril can be used in many areas of surgery. e.g., cardiovascular surgery, haemorrhoidectomy, implantation of vascular prostheses, biopsies, lung operations, surgery to the face and jaw, gastric resection, operations on the throat or nose, liver and gallbladder operations, gynecological operations, thoracic and abdominal sympathectomies, neurosurgery (especially cerebral operations), thyroid operations, skin transplantations, and treatment of superficial injuries. Dimensions: REF SOF-0101 – 1 x 1 inch (2.5 x 2.5 cm); REF SOF-0102 – 1 x 2 inches (2.5 x 5.1 cm); REF SOF-0202 – 2 x 2 inches (5.1 x 5.1 cm); REF SOF-0203 – 2 x 3 inches (5.1 x 7.6 cm); REF SOF-0204 – 2 x 4 inches (5.1 x 10.2 cm); REF SOF-0304 – 3 x 4 inches (7.6 x 10.2 cm); REF SOF-0404 – 4 x 4 inches (10.2 x 10.2 cm)
Surgi-ORC Woven/SNOW	Non-	Composition — 100% oxidized regenerated cellulose (ORC) Characteristics – Structured non-woven ORC matrix, needle punched with interlocking fibres and feasible in both open and minimally invasive procedures. Increased surface contact to the bleeding site. Indications: Material is intended to be used as an adjunctive haemostat in surgical procedures to assist in the control of capillary, venous, and small arterial hemorrhage when ligation or other conventional methods of control are impractical or ineffective. Surgi-ORC Non-Woven/SNOW can be cut to size for use in endoscopic procedure. Surgi-ORC Non-Woven/SNOW material can be used in many areas of surgery. e.g., cardiovascular surgery, haemorrhoidectomy, implantation of vascular prostheses, biopsies, lung operations, surgery to the face and jaw, gastric resection, operations on the throat or nose, liver and gallbladder operations, gynecological operations, thoracic and abdominal sympathectomies, neurosurgery (especially cerebral operations), thyroid operations, skin transplantations, and treatment of superficial injuries. Dimensions: REF SON-0102 – 1 x 2 inches (2.5 x 5.1 cm); REF SON-0202 – 2 x 2 inches (5.1 x 5.1 cm); REF SON-0203 – 2 x 3 inches (5.1 x 7.6 cm); REF SON-0204 – 2 x 4 inches (5.1 x 10.2 cm); REF SON-0304 – 3 x 4 inches (7.6 x 10.2 cm); REF SON-0404 – 4 x 4 inches (10.2 x 10.2 cm)

CONTRAINDICATIONS

Surgi-ORC material can be used as a form of tamponade. When used in this manner, it is to be removed after hemostasis is achieved.

Surgi-ORC material should not be used for implantation in bone defects, such as fractures, since there is a possibility of interference with callus formation and a theoretical chance of cyst formation.

When Surgi-ORC material is used to help achieve hemostasis in, around, or in proximity to bone holes, areas of bone junction, during the spinal cord, or the optic nerve and chiasm interventions, it must always be removed after hemostasis is achieved since it will swell and could exert unwanted pressure.

Surgi-ORC material should not be used to control hemorrhage from large arteries.

Surgi-ORC material should not be used on non-hemorrhagic serous oozing surfaces, since body fluids other than whole blood, such as serum, do not react with Surgi-ORC material to produce satisfactory hemostatic effect.

Surgi-ORC material is an absorbable hemostat, and should not be used as an adhesion prevention product.

Surgi-ORC material should not be used in patients with known hypersensitivity to cellulose.

POSSIBLE SIDE EFFECTS

Cases of "Encapsulation" of fluid and foreign body reactions have been reported.

There have been reports of stenotic effect when ORC material has been applied as a wrap during vascular surgery.

Although it has not been well established that stenosis was directly related to the use of ORC haemostatic material, it is important to exercise caution and avoid tight application of the material as a wrap.

Paralysis and nerve damage have been reported when ORC hemostatic material were used around, in, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm.

Blindness has been reported in connection with surgical repair of a lacerated left frontal lobe when ORC hemostat was placed in the anterior cranial fossa.

Possible prolongation of drainage in cholecystectomies and difficulty passing urine per urethra after prostatectomy have been reported.

Burning, stinging sensations, and sneezing may occur when ORC haemostatic material is used as a tamponade for nasal bleeding, which are believed to be due to the low pH of the product.

WARNINGS

Surgi-ORC material is supplied sterile and must not be re-sterilized.

Surgi-ORC material is intended for single use and cannot be reused. Reuse of the Surgi-ORC material may result in cross-infection.

Closing with Surgi-ORC material in a contaminated wound without drainage may lead to complications and should be avoided.

Surgi-ORC material should not be moistened with water or saline prior to application, as it may influence the hemostatic effect.

Surgi-ORC material should not be impregnated with anti-infective agents or with other materials such as buffering or hemostatic substances.

It is advisable to remove Surgi-ORC material once hemostasis is achieved. It must always be removed from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm, regardless of the type of surgical procedure because Surgi-ORC material, by swelling, may exert pressure resulting in paralysis and/or nerve damage.

Dislodgement of Surgi-ORC material could possibly occur by means such as repacking, further intraoperative manipulation, lavage, exaggerated respiration, etc.

There have been reports that in procedures such as lobectomy, laminectomy and repair of a frontal skull fracture and lacerated lobe that Surgi-ORC Absorbable Hemostat migrated from the site of application.

Special care must be taken by physicians, regardless of the type of surgical procedure, to consider the advisability of removing Surgi-ORC Absorbable Hemostat after hemostasis is achieved.

Although Surgi-ORC Absorbable Hemostat is bactericidal against a wide range of pathogenic microorganisms, it is not intended as a substitute for systemically administered therapeutic or prophylactic antimicrobial agents to control or prevent post-operative infections.

Surgi-ORC Absorbable Hemostat is not intended for the primary treatment of coagulation disorders.

The safety and efficacy of Surgi-ORC Absorbable Haemostat in pregnant women have not been established.

PRECAUTIONS

Use only as much Surgi-ORC Absorbable Hemostat as is necessary for hemostasis, holding it firmly in place until bleeding stops. Remove any excess before surgical closure in order to facilitate absorption and minimize the possibility of foreign body reaction.

In urological procedures, minimal amounts of Surgi-ORC Absorbable Hemostat should be used and care must be exercised to prevent plugging of the urethra, ureter, or a catheter by dislodged portions of the product.

Surgi-ORC Absorbable Hemostat could be prevented in chemically cauterized areas.

If Surgi-ORC Absorbable Hemostat is used temporarily to line the cavity of large open wounds, it should be placed so as not to overlap the skin edges.

Surgi-ORC Absorbable Hemostat should also be removed from open wounds by forceps or by irrigation with sterile water or saline solution after bleeding has stopped.

Precautions should be taken in otorhinolaryngologic surgery to ensure that none of the material is aspirated by the patient.

Care should be taken not to apply Surgi-ORC Absorbable Hemostat too tightly when it is used as a wrap during vascular surgery.

Packing should be avoided, especially in rigid cavities where swelling may interfere with normal function or cause necrosis.

Endoscopic procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques.

DOSAGE AND ADMINISTRATION

Sterile technique should be followed when removing Surgi-ORC Absorbable Haemostat from the sterile container. Minimal amounts of appropriately sized Surgi-ORC Absorbable Haemostat is applied on the bleeding site or pressed tightly against the tissue until haemostasis is achieved. The necessary amount of Surgi-ORC Absorbable Haemostat depends on the nature and intensity of the bleeding to be stopped. It is not recommended to moisten the material with water or saline solution before use. Opened but unused Surgi-ORC Absorbable Haemostat should be disposed of as it cannot be re-sterilized.

The medical device is intended for use (operation) under standard operating conditions in a hospital operating room.

The medical device is not subject to maintenance and repair. The product is disposable.

SET COMPLETENESS

The packaging of the medical device contains:

- Surgi-ORC material of one design variant — 12 pcs;
- IFU — 1 pc.

STORAGE

The medical device is transported by all means of transport at the temperatures from 0 to 30 °C. The relative humidity should not exceed 85 %.

Store the material in its original packaging in a clean and dry room at a temperature not exceeding 30 °C. The relative humidity should not exceed 85 %.

The outer packaging of the product should be closed immediately after use.

Do not refrigerate or freeze the product.

The shelf life of the medical device is 36 months from the manufacturing date.

The shelf life is indicated on the packaging of the medical device; the product must not be used after the expiry date.

DISPOSAL

Products coming into contact with blood and/or biological fluids are class B medical devices under medical wastes classification and may be recycled and/or disposed in accordance with applicable national regulations.

Unused products (not having contact with blood and/or biological fluids), including those with expired shelf life, belong to class A medical device under medical wastes classification and may be recycled and/or disposed in accordance with applicable national regulations (SanPiN 2.1.3684-21).

MEDICAL DEVICE MANUFACTURER

Aegis Lifesciences Private Limited

215/216, Mahagujarat Industrial Estate, Moraiya, Ahmedabad-382213, Gujarat
India

Tel.: +91 9726903131, +91 9727103131

E-mail: info@aegis-lifesciences.com, admin@aegis-lifesciences.com

Web-site: www.aegis-lifesciences.com

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER OF A MEDICAL DEVICE IN THE TERRITORY OF THE RUSSIAN FEDERATION

Medical Research and Production Company Avis Life Science Limited Liability Company
(MRPC Avis Life Science LLC)

1D, office 315, Leningradskaya St., Khimki, urban district Khimki, Moscow Region,
Russia

Tel.: +7 (965) 120 20 85

E-mail: mnpkavis@bk.ru

SYMBOLS USED IN THE LABELLING

Prohibition on reuse

Warning!

Manufacturer

Do not re-sterilize

Do not use when
the package is
damaged

Date of
manufacture

Protect from
moisture

Keep away from
sunlight

Radiation
sterilization

Refer to
the Instructions for
Use

Expiry date


Catalogue
number

Upper temperature
range

Batch code

Authorized
representative
in the
European
Community

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