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#501 WonPoong B/D, 301 Hwagok-ro, Gangseo-gu, Seoul Korea

Statement

PROTIA Inc. declares that we have issued User manuals as attached document:

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User manuals

In vitro diagnostic assay kit for qualitative detection of
autoantigen-specific IgG antibodies in human serum or plasma

PROTIA ANA Profile 18 Autoimmune panel

Signature

Kook Jin Lim
President / CEO

Instructions for use of a medical device

In vitro diagnostic assay kit for qualitative detection of autoantigen-specific IgG antibodies in human serum or plasma "PROTIA ANA Profile 18 Autoimmune Panel"

The name of medical device (full)

In vitro diagnostic assay kit for qualitative detection of autoantigen-specific IgG antibodies in human serum or plasma "PROTIA ANA Profile 18 Autoimmune Panel"

Name of the medical device (short)

"PROTIA ANA Profile 18 Autoimmune Panel"

Intended use

In vitro diagnostic test kit for qualitative detection of autoantigen-specific IgG antibodies in human serum or plasma "PROTIA ANA Profile 18 Autoimmune Panel" (further in the text as "PROTIA ANA Profile 18 Autoimmune panel") is intended for the qualitative determination of autoantigen-specific IgG in human serum or plasma (anticoagulant: lithium heparin, sodium citrate, EDTA-K2) using enzyme-linked immunosorbent assay (ELISA).

The kit is used as an aid in the diagnosis of autoimmune diseases.

The use of the product has no population or demographic aspects.

Functional purpose

The kit is used as an aid in the diagnosis of autoimmune diseases.

Indications for use

Mixed connective tissue disease (MCTD=SLE+SS+SSD+myositis), systemic lupus erythematosus (SLE), Sjogren's syndrome (SS), myositis (poly-/dermato-) (PM/DM), systemic sclerosis (systemic scleroderma) (SSc), primary biliary cirrhosis (PBC).

Contraindications

None.

Possible side effects

No side effects were identified.

Potential consumers of a medical product

Health care facilities of any profile.

Information on population and demographic aspects of the use of a medical device

The use of the product has no population or demographic aspects.

Professional level of potential users

Clinical laboratory diagnostics doctor, medical laboratory technician.

Application area

Clinical laboratory diagnostics in vitro.

List and description of medical device materials that come into direct or indirect contact with the patient's body

The product contains no materials that come into direct or indirect contact with the body of the patient and personnel using the product, subject to the requirements of the operational documentation (instructions for use).

List of materials of animal and (or) human origin, medicines and pharmaceutical substances
The product contains no materials of human origin, medicines and pharmaceutical substances.

Methods for sterilization and disinfection of products

The product is non-sterile and does not require sterilization or disinfection.

Frequency of use of a medical device

Single use product.

Product operation software

No software available.

Product maintenance and repair

The product does not require maintenance and cannot be repaired.

Brief information

An autoimmune disease is an abnormal immune response that occurs when the immune system's self-tolerance does not develop, or when tolerance is broken and T cells are activated to attack themselves in response to the body's own cells or tissues, resulting in the production of autoantibodies that attack their own cells or tissues. Specifically, systemic autoimmune disease activates immune cells to attack self-antigens, causing extensive damage to inflamed areas and various organs. Because pathogenic autoantigens are produced differently in different tissues or organs of the human body, inflammatory diseases progress differently in each patient. Typical systemic autoimmune diseases include systemic lupus erythematosus (SLE), Sjögren's syndrome, systemic sclerosis (SSc), myositis, and mixed connective tissue disease (MCTD), and each disease produces autoantibodies against a specific autoantigen. Therefore, the diagnosis and treatment of patients with suspected systemic autoimmune disease is aided by the measurement of various autoantigen-specific IgG antibodies, including antinuclear antibodies.

Test principle

PROTIA ANA Profile 18 Autoimmune Panel, based on the principle of enzyme immunoassay, includes nitrocellulose membranes containing various autoantigens, which are adsorbed on strips located at regular intervals, allowing testing for 18 autoantigen-specific IgG antibodies in a single test. If autoantigen-specific IgG antibodies bind to autoantigens, they are immobilized on the membrane after the washing step. The immobilized IgG antibodies bind to biotin anti-human IgG antibodies, and the biotin is captured by streptavidin conjugated to alkaline phosphatase. The color develops after the addition of the substrate solution in the last step of enzyme incubation, and the color intensity is analyzed using a color measuring device.

Composition of the kit

In vitro diagnostic assay kit for the qualitative detection of autoantigen-specific IgG antibodies in human serum or plasma "PROTIA ANA Profile 18 Autoimmune Panel":

Name	Content	Quantity
Autoantigen panel	Autoantigen panel	20 panels (2 packs of 10 panels)
Sample diluent	Sample diluent	15 ml x 2 E.A.
Antibody solution	Antibody solution	20 ml x 1 E.A.
Enzyme solution	Enzyme solution	20 ml x 1 E.A.
Substrate solution	Substrate solution	20 ml x 1 E.A.
Wash solution 20x	Washing solution 20x	45 ml x 1 E.A.
Instructions	Instructions	1 PC.

An in vitro diagnostic test kit for the qualitative detection of autoantigen-specific IgG antibodies in human serum or plasma "PROTIA ANA Profile 18 Autoimmune Panel" is designed for the study of 20 samples and is intended for use in healthcare institutions.

Assaying is possible manually, semi-automatically and automatically using the following equipment:

1. Manual method: Immunological photometric analyzer "Q-SMART" for in vitro diagnostics on "PROTIA" panels and a rocking shaker.
2. Semi-automatic method: Automated device "Q-PROCESSOR" for processing "PROTIA" panels for in vitro diagnostics and Immunological photometric analyzer "Q-SMART" for in vitro diagnostics on "PROTIA" panels.
3. Automatic immunological photometric analyzer "Q-STATION ELITE" for in vitro diagnostics on "PROTIA" panels.

Description of components of the kit for in vitro diagnostic analysis for the qualitative detection of autoantigen-specific IgG antibodies in human serum or plasma "PROTIA ANA Profile 18 Autoimmune Panel":


1. Autoantigen panel: Sky blue plastic panel with three white membranes attached to it
2. Sample diluent: Blue-violet liquid
3. Antibody solution: Yellow liquid
4. Enzyme solution: Colorless liquid
5. Substrate solution: Light yellow liquid
6. Washing solution 20x: Colorless liquid

Design and main functional elements.

Composition of components and their quantity in the PROTIA ANA Profile 18 Autoimmune panel

No.	Name	Description
1	Autoantigen panel	A plastic cassette with a plastic frame containing membranes with adsorbed autoantigens. A plastic panel of sky blue color, to which three white membranes are attached (S-1, S-2, S-3 - standard lines; from 1 to 18 - autoantigen lines) 20 panels (2 packs of 10 panels). Each membrane has a different colored top strip and a white bottom strip.
2	Sample diluent	Two white 20 ml bottles with a blue screw cap each contain 15 ml of blue-violet liquid
3	Antibody solution	White 20 ml bottle with yellow screw cap containing 20 ml of yellow liquid
4	Enzyme solution	White 20 ml bottle with red screw cap containing 20 ml of colorless liquid
5	Substrate solution	Black 20 ml bottle with black screw cap containing 20 ml of light yellow liquid
6	Washing solution 20x	White 50 ml bottle with white screw cap containing 45 ml of colorless liquid
7	Instructions	1 PC

Autoantigen panel

Description of autoantigen panel

No.(1)	Autoantigens	Code
1	Histidyl-tRNA synthetase	Jo-1
2	PM-Scl 100	PM-Scl
3	Centromere protein B	CENP-B
4	Scl-70	Scl-70
5	Ro/SS-A (60 kDa)	SSA 60
6	Ro/SS-A (52 kDa)	SSA 52
7	La/SS-B	SSB
8	M2	M2
9	dsDNA	dsDNA
10	Nucleosome	Nuc.
11	Histone	Histone
12	Proliferating Cell Nuclear Antigen	PCNA
13	Ribosomal Phosphoprotein P0	RPP
14	Sm	Sm
15	U1-snRNP 68/70 kDa	RNP 70
16	U1-snRNP A	RNP A
17	U1-snRNP C	RNP C
18	RNP/Sm	RNP/Sm

(1) The strip numbers are shown in the figure.

Analyzed samples

Human serum or plasma (anticoagulant: lithium heparin, sodium citrate, EDTA-K2)

Test procedure.**1. Preparation before the test**

All reagents should be brought to room temperature 30 minutes before use and mixed well immediately before use. After bringing to room temperature, open the package with the autoantigen panel.

NOTE. Autoantigen panels not used for testing should be immediately placed in an aluminum bag, sealed and refrigerated.

2. Preparation of reagents and samples**(1) Preparation of wash solution (1x)**

Before testing, dilute the wash solution 20 times with deionized water.

Example: 19 ml deionized water + 1 ml 20x rinse solution = 20 ml 1x rinse solution

(2) Sample preparation

The test uses human blood serum or plasma (anticoagulant: lithium heparin, sodium citrate, EDTA-K2). Before testing, remove blood cells or any solids by centrifugation. Hemolyzed or contaminated samples may result in incorrect results. Store serum/plasma samples at 2~8°C for short-term use (within 2 weeks) and -20°C or lower for longer-term use. Repeated freezing and thawing of serum/plasma samples should be avoided. Sample volume required for testing: 9 µl.

**Test Methods****Manual method**

1) The optimal conditions for the laboratory are to conduct the experiment at + 25°C ± 2 °C and at a humidity of 50% ± 20%.

2) All reagents should be brought to room temperature 30 minutes before use and mixed thoroughly immediately before use. Open the packaging of the allergen panel after its temperature has equalized to room temperature.

NOTE: Autoantigen panels not used for testing should be immediately placed in an aluminum bag, sealed and refrigerated.

3) Wet the test membrane completely by adding 600 µL of diluted 1x wash solution to the allergen

panel, rock the panel for 5 minutes, and then remove any remaining wash solution (rocking frequency recommended: up to 100 oscillations/min).

4) Add in to autoantigen panel 891 μL of sample diluent.

5) Dispense 9 μL of patient sample into the panel and incubate while rocking the panel at room temperature for 45 minutes.

6) Remove any remaining solutions from the panel and rinse the membrane twice with diluted washing solution 1x. At each washing step, add 900 μL of diluted 1x wash solution, incubate while rocking the panel for 5 minutes and then remove 1x wash solution from the panel, the solution should not remain in the panel.

7) Add 500 μL of antibody solution to the panel and incubate while shaking the panel at room temperature for 30 minutes.

8) Remove the antibody solution from the panel and rinse the membrane twice with diluted washing solution 1x. At each washing step, add 600 μL of diluted 1x washing solution, incubate while rocking the panel for 5 minutes and then remove 1x washing solution from the panel, the solution should not remain in the panel.

9) Add 500 μL of enzyme solution to the panel and incubate while shaking the panel for 30 minutes.

10) Remove the enzyme solution from the panel and rinse the membrane as indicated in step (8).

11) Add 500 μL of substrate solution and incubate while shaking the panel at room temperature in a dark place for 20 minutes.

12) Remove the substrate solution from the panel and rinse the membrane once with 1000 μL of deionized water.

13) Leave the panel on the table until it dries naturally.

NOTE: Make sure the membrane is completely dry before moving on to the next step.

14) Insert the panels into the "Q-SMART" immunological photometric analyzer for in vitro diagnostics on the "PROTIA" panels and evaluate the results. For more details, see the Manual for Immunological photometric analyzer "Q-SMART" for in vitro diagnostics on "PROTIA" panels

* for incubation and shaking of the panel, use any shaking shaker with a swing speed of up to 100 oscillations / min. There are no restrictions on the choice of device model.

Semi-automatic method

1) Automated device "Q-PROCESSOR" for processing "PROTIA" panels for in vitro diagnostics.

2) Refer to the instruction manual on Automated device "Q-PROCESSOR" for processing "PROTIA" panels for in vitro diagnostics

3) Insert the panels into the Q-SMART analyzer and evaluate the results. For more details, see the manual for the Immunological Photometric Analyzer "Q-SMART" for in vitro diagnostics on "PROTIA" panels.

Automatic method

1) Automatic immunological photometric analyzer "Q-STATION ELITE" for in vitro diagnostics on "PROTIA" panels.

2) Refer to the instruction manual for Automatic immunological photometric analyzer "Q-STATION ELITE" for in vitro diagnostics on "PROTIA" panels.

3. Assessment and interpretation

The level of autoantigen-specific IgG antibodies is qualitatively analyzed using measuring devices in AU/mL.

The test results are defined as negative or positive as shown in the table below.

AU/mL	Level of autoantigen-specific IgG antibodies	Results
< 11.00	Not found or weak manifestation	Negative
≥ 11.00	Strong manifestation	Positive

Quality control

After completion of the test, the PC strip (S1) in the middle of the membrane should be clearly visible. If the color intensity of the reference line is weak, it will not be read correctly by optical measuring instruments. In this case, it is recommended to retest.

Analytical and diagnostic characteristics:

1) Threshold value: 11.00 AU/mL*

*AU/mL was arbitrarily set as the signal level value obtained on the measuring device for the company's product.

2) Analytical specificity: No noticeable interference or cross-reactivity is observed at the following concentrations.

Substances that cause interference and cross-reactivity	Concentration
Lithium heparin	330 units/dL
Sodium citrate	3.2%
EDTA-K2	0.1 mg/dL
Hemoglobin	1000 mg/dL
Bilirubin (unconjugated)	40 mg/dL
Triglyceride	1500 mg/dL
Glucose	1000 mg/dL
Human serum albumin	6 g/dL
Rheumatoid factor	60 IU/mL
Biotin	0.351 mg/dL
Human IgE	300 IU/mL
Human IgM	4 mg/mL
Human IgA	5 mg/mL
Human IgD	0.09 mg/mL

3) Accuracy: Repeatability (within batches, within a day, between test dates) and reproducibility (between batches, test sites and testers, instrumentation, testing equipment) of test results meet the acceptance criteria (CV within 20%).

4) Diagnostic sensitivity 95% and more.

5) Diagnostic specificity 95% and more.

Precautions for use

- 1) For diagnostic purposes only *in vitro* and professional use.
- 2) The "PROTIA ANA Profile 18 Autoimmune panel" kit may be used as a diagnostic aid, and the final clinical diagnosis or dosing regimen for immunotherapy should be determined by the physician after evaluation of all clinical and laboratory results.
- 3) Between results of "PROTIA ANA Profile 18 Autoimmune panel" and the results of *in vivo* and/or other *in vitro* tests may be variations due to the lack of a national or international standard, which may cause the levels of autoantigens to differ between tests.
- 4) False-positive test results may be obtained due to cross-reactivity of the test autoantigens with other autoantigens.
- 5) Smoking, eating and drinking are prohibited in areas where samples or test reagents are used.
- 6) All samples may potentially contain unknown infectious materials. When handling samples, you must wear disposable gloves and wash your hands after testing.
- 7) For safety reasons, it is prohibited to place needles, knives or any other objects that could cause injury near areas where human blood samples and reagents are handled.
- 8) All patient blood samples and used kit components must be considered as biohazards. They must be disposed in according to appropriate guidelines.
- 9) Do not use kits after the expiration date.
- 10) Autoantigen panels are supplied with a desiccant and must be properly sealed after each use. Each strip is for single use.
- 11) Some reagents included in the kit contain sodium azide as a preservative. Sodium azide is known to react with lead or copper to form a potentially explosive metal azide when drained. Therefore, it is always necessary to rinse with plenty of water to avoid the formation of metal azide.
- 12) The substrate solution may contain a black precipitate. This is not an abnormal phenomenon and does not affect the test results in any way.
- 13) Be careful not to form bubbles. Especially when using automatic devices, bubbles should be removed before starting the test as they can affect the volume of distribution.
- 14) If the membranes are not completely dried after the final reaction, the test results analyzed by the measuring instruments may be incorrect.
- 15) It is prohibited to mix or use solutions from different batches.
- 16) Products must be stored in the refrigerator at a temperature of 2~8°C and avoid freezing them.
- 17) If any serious incident occurs in connection with this medical device, the manufacturer should be contacted.

Disposal and destruction

The product must be disposed of in accordance with local, regional and national regulations*.

*Products during use, storage and transportation do not emit harmful substances that pollute the environment. Waste generated as a result of using a set of reagents for the intended purpose established by the manufacturer, as well as unused products (expired expiration date, damaged consumer packaging/labeling, damaged packaging/labeling of reagents, etc.) are classified as medical waste of class "A", "B" in accordance with Art. 49 of the Federal Law of November 21, 2011 N 323-FZ "On the fundamentals of protecting the health of citizens in the Russian Federation."

Medical waste is subject to collection, neutralization, placement, storage, transportation, accounting and disposal in accordance with current rules and regulations (SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of territories of urban and rural settlements, for water bodies, drinking water and drinking water supply to the population, atmospheric air, soil, residential premises, operation of industrial and public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" in the current version, "Guidelines for disinfection, pre-sterilization cleaning and sterilization of medical devices" (approved by the Ministry of Health of Russia on December 30, 1998 N MU-287-113).

There are no environmental or physical risks if disposal methods are followed.

Packaging unit: 1 kit (for 20 tests).

Storage conditions: Store at 2~8°C.

Best before date: 18 months from date of manufacture (each panel is for single use). After opening, it can be stored at a temperature of 2 to 8°C for no more than 2 months.

Transportation allowed at temperatures from -15°C to +30°C for 5 days; buy transport of all types in covered vehicles in accordance with the transportation rules in force for transport of this type.

Symbols used on packaging and labels

Symbol	Description
	Manufacturer
	Carefully!
	Temperature Range
	Prohibition on re-use
	Date of manufacture

Symbol	Description
	Medical device for in vitro diagnostics
	Refer to instructions for use
	Enough content to run <n> tests
	Biological risk

Symbol	Description
	Use before
	Batch code
	Catalog number
	Authorised representative

[Manufacturer]

PROTIA Inc.

A-702, 401, Yangcheon-ro, Gangseo-gu, Seoul, 07528, Republic of Korea

Tel.: +82-2-6968-5278, fax: +82-2-6442-8345

www.protia.co.kr, Email: protia@protia.co.kr

{Authorised representative}

On the territory of the Russian Federation, all questions regarding the quality and circulation of a medical device should be addressed to: 109387, Moscow, st. Tikhaya, 28, office. 216; Limited Liability Company "BIO FOCUS" (BIO FOCUS LLC), info@biofocus.ru. Tel. (495) 352-8777.

PROTIA Inc.

PRI-ANP01-en002 (20240517)

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Registered No. 2024 - 00964

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Notarial Certificate

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HUANG YUXUAN

attorney-in-fact of

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PROTIA INC.

KOOK JIN LIM / President

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appeared before me and admitted said
principal's subscription to the
STATEMENT

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This is hereby attested on this 28 th
day of May., 2024 at this office.

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HANLIM LAW AND NOTARY OFFICE

BELONG TO SEOUL SOUTHERN
DISTRICT PROSECUTOR'S OFFICE

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HOON OH

Attorney at law acting as Notary Public

This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public since
February7, 2020 under Law No. 15150.

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