Registered No. 2024-8868

## NOTARIAL CERTIFICATE

HANMI LAW AND NOTARY OFFICE

Twintree Tower(B) #101,6 Yulgok-ro,  
Jongno-gu, Seoul, Republic of Korea

# Cellprep

The next generation of Liquid-based Cytology

**Cellprep Respiratory Tract  
Solution for Liquid-Based  
Cytology Processor**

## **Instructions for Use**

[www.biodyne.asia](http://www.biodyne.asia)

Kyungeun Jang  
Chief Quality Officer

Approved

  
**Biodyne** (Signature)

24-4, Achesan-ro 5-gil  
Seongdong-gu, Seoul, Republic of Korea

## Table of contents

|                                                                                                                                                                                                                                                                |    |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----|
| 1. General information.....                                                                                                                                                                                                                                    | 5  |
| 2. Description of the Cellprep Respiratory Tract Solution .....                                                                                                                                                                                                | 6  |
| 2.1. <i>Intended use, indications and contra-indications</i> .....                                                                                                                                                                                             | 6  |
| 2.2. <i>Description of the solution</i> .....                                                                                                                                                                                                                  | 6  |
| 2.3. <i>Principles of operation</i> .....                                                                                                                                                                                                                      | 7  |
| 2.4. <i>Delivery set</i> .....                                                                                                                                                                                                                                 | 7  |
| 3. Type of the sample analyzed .....                                                                                                                                                                                                                           | 7  |
| 4. Information on conditions necessary for collection, processing and preparation of samples, data on stability of the samples analyzed, including storage conditions and duration, transportation conditions, restrictions on freezing (thawing) cycles ..... | 7  |
| 5. Information on consumers and potential users .....                                                                                                                                                                                                          | 8  |
| 6. Installation .....                                                                                                                                                                                                                                          | 8  |
| 7. Application of the Cellprep Respiratory Tract Solution .....                                                                                                                                                                                                | 8  |
| 7.1. <i>Description of the medical or non-medical devices intended for use in combination with the Cellprep Respiratory Tract Solution</i> .....                                                                                                               | 8  |
| 7.1.1. <i>Equipment for slide preparation for thin layer cytology or liquid-based cytology</i> .....                                                                                                                                                           | 8  |
| 7.1.2. <i>Slides</i> .....                                                                                                                                                                                                                                     | 9  |
| 7.1.3. <i>Vortex</i> .....                                                                                                                                                                                                                                     | 9  |
| 7.2. <i>Application of the Cellprep Respiratory Tract Solution</i> .....                                                                                                                                                                                       | 9  |
| 7.2.1. <i>Sputum</i> .....                                                                                                                                                                                                                                     | 9  |
| 7.2.1. <i>Bronchial washing, bronchial brushing</i> .....                                                                                                                                                                                                      | 9  |
| 7.3. <i>Insertion of the vial with the Cellprep Respiratory Tract Solution into equipment for slide preparation for thin layer cytology/liquid-based cytology</i> .....                                                                                        | 9  |
| 7.4. <i>Interfering substances or sample limitations that may interfere with study result</i> .....                                                                                                                                                            | 17 |
| 7.5. <i>Identification of the medical devices in order to obtain a safe combination and (or) information on known restrictions on the combined use of medical devices</i> .....                                                                                | 17 |
| 8. Specifications of the Cellprep Respiratory Tract Solution.....                                                                                                                                                                                              | 18 |
| 9. Recommendations for staining the slides .....                                                                                                                                                                                                               | 19 |
| 10. A. B. C. Suction Numerical Standard Values for LBC Processors .....                                                                                                                                                                                        | 19 |
| 11. Cautions and warnings .....                                                                                                                                                                                                                                | 20 |
| 12. Transportation and storage conditions.....                                                                                                                                                                                                                 | 20 |
| 13. Packaging .....                                                                                                                                                                                                                                            | 20 |
| 14. Symbol identification .....                                                                                                                                                                                                                                | 22 |
| 15. Manufacturer's responsibility.....                                                                                                                                                                                                                         | 23 |

|                                                                                      |    |
|--------------------------------------------------------------------------------------|----|
| 16. Repair and Maintenance .....                                                     | 23 |
| 17. Cleaning and Disinfection .....                                                  | 23 |
| 18. Warranty and shelf life .....                                                    | 23 |
| 19. Disposal .....                                                                   | 24 |
| 20. Information on the initial release and revision of the Instructions for Use..... | 24 |
| 21. List of the standards applied .....                                              | 24 |
| Appendix No. 1.....                                                                  | 27 |
| Certificate of Analysis - Cellprep Respiratory Tract Solution (sample).....          | 27 |
| Appendix No. 2 Material Safety Data Sheet .....                                      | 30 |

# *Cellprep*

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## 1. General Information

Name of the medical device - **Cellprep Respiratory Tract Solution for Liquid-Based Cytology Processor** (hereafter - "Cellprep Respiratory Tract Solution", "Respiratory Tract Solution", "Respiratory Tract Solution", "Cellprep Solution(s)", "medical device").


The Cellprep Respiratory Tract Solution ensures optimal preservation of cells of sputum, bronchial washings, bronchial brushing, bronchoalveolar lavage, preserving their morphology, without deforming the nuclei and without causing hyperchromasia.

### Use:

For studying the lungs, pharynx, larynx, respiratory tract for cancer

The specially designed vial, in which the solution is kept, is optimized for Cellprep LBC devices.

The vial cap consists of two parts: middle cap and top cap. This design provides that after collecting the sample from the patient, the physician can open only the middle cap and place the sample in a vial with solution. The top cap of the vial can be opened only by the Cellprep LBC Processor during slide preparation. The middle cap is equipped with an integrated nano-polyester filter that helps to remove impurities and obtain a monolayer preparation.

The medical device is intended for clinical laboratory diagnostics (in vitro diagnostics).

The device is for single use only.

Class of potential risk of solution application – 2a

Type of the medical device in accordance with the Russian Nomenclature Classification – 214660.

| Solution name                       | Catalog No. | Volume |
|-------------------------------------|-------------|--------|
| Cellprep Respiratory Tract Solution | RPT-1000    | 20 ml  |

The medical device is intended for professional use only.

The Cellprep Respiratory Tract Solution for Liquid-Based Cytology Processor relates to the medical devices which have no direct or indirect contact with the human body.

The medical device is not sterile, and it is not a subject to sterilization.

## **2. Description of the Cellprep Respiratory Tract Solution**

### **2.1. Intended use, indications and contra-indications**

#### **Intended use and indications:**

The Cellprep Respiratory Tract Solution is an auxiliary device for in vitro diagnostics and is intended for preservation of exfoliated cells of sputum, bronchial washing, bronchial brushing, bronchoalveolar lavage for subsequent preparation of slides by liquid-based cytology using the Cell Prep Plus LBC processor or the Cellprep AUTO LBC system.

Monolayered slides prepared are stained in accordance with approved protocols, followed with the results interpretation.

#### **Contra-indications:**

Not intended for preservation of histological samples (biopsy).

#### **Side effects:**

- not applicable.

|                                                                                                                                        |                                                                                                                                                                                                                                                           |
|----------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description of target analyte                                                                                                          | Not applicable                                                                                                                                                                                                                                            |
| Functional use                                                                                                                         | Auxiliary device for in vitro diagnostics                                                                                                                                                                                                                 |
| Specific pathology, conditions or risk factor which detection, determination or differentiation the IVD medical device is intended for | Suspected neoplastic tumor of the respiratory system.                                                                                                                                                                                                     |
| Population, demographic aspects of the medical device use                                                                              | The Cellprep Respiratory Tract Solution is intended for optimal preservation of exfoliated cells of sputum, bronchial washing, bronchial brushing, bronchoalveolar lavage, collected from patients, regardless of their gender, race and age differences. |
| Information on analytical sensitivity (limit of detection), analytical specificity, diagnostic sensitivity and diagnostic specificity. | Not applicable                                                                                                                                                                                                                                            |

### **2.2. Description of the solution**

Colorless transparent liquid without foreign impurities

### **Solution composition**

#### **Components:**

| CHEMICAL NAME                     | CAS No.   | % (ppm) |
|-----------------------------------|-----------|---------|
| Ethanol                           | 64-17-5   | 34%     |
| Deionised water                   | 7732-18-5 | 49%     |
| Glycerine                         | 56-81-5   | 16%     |
| Phosphorous acid                  | 7664-38-2 | 0.48%   |
| Ammonium acetate                  | 631-61-8  | 0.48%   |
| Sodium Citrate Tribasic Dihydrate | 6132-04-3 | 0.04%   |

The Cellprep Respiratory Tract Solution does not contain medicinal products for medical use, materials of animal and human origin.

### **2.3. Principle of operation**

- The solution preserves and fixes the cells collected.
- The material (sample) fixed in the vial, while keeping the temperature regime from +10 °C to +30 °C, can be used for liquid-based cytology studies for up to 3 months.

### **2.4. Delivery set**

Cellprep Respiratory Tract Solution for Liquid-based Cytology Processor  
Including:

1. Cellprep Respiratory Tract Solution for Liquid-based Cytology Processor, 20 ml - 25 pcs/pack.
2. Instructions for Use in paper - 1 pc. for a shipping box (200 vials)

## **3. Type of the sample analyzed**

Exfoliated cells of sputum, bronchial washing, bronchial brushing, bronchoalveolar lavage.

## **4. Information on conditions necessary for collection, processing and preparation of samples, data on stability of the samples analyzed, including storage conditions and duration, transportation conditions, restrictions on freezing (thawing) cycles**

Recommendations for organization of preliminary (pre-laboratory) stage of work with samples.

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sample collection | <b>1. Sputum</b><br>Recently coughed up sputum (about 3 cc) is collected and placed in a Cellprep Respiratory Tract solution vial. To collect an adequate material, the cough must be "deep-seated", it is necessary to collect sputum before eating, after brushing the teeth.<br><br><b>2. Bronchial washing, bronchial brushing</b><br>Material for cytological examination: material obtained by bronchoscopy (tumor puncture sample, biopsy print, brushing material, bronchioloalveolar lavage (BAL), transbronchial puncture (under visual or ultrasound control). If a neoplastic process is suspected and the morphological study of this material fails, a transthoracic puncture is performed. In addition, it is possible to study the material collected during mediastinoscopy, thoracoscopy, surgery by the cytological method. |
| Preservation      | The Cellprep Respiratory Tract Solution preserves and fixes the cells collected.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

|                |                                                                                                                                                                                                                                                                                                        |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Label          | Label the vial with sample in the designated area on the vial label (hand lettering according to referral or barcode label).                                                                                                                                                                           |
| Storage        | The material (sample) fixed in the vial, while keeping the temperature regime from +10 °C to +30 °C, can be used for liquid-based cytology studies for up to 3 months.<br><b>Freeze/Thaw Prohibited.</b><br>Cooling, especially freezing, of both non-fixed and fixed material is strictly prohibited. |
| Transportation | Transportation temperature regime: from +10 °C to +30 °C.                                                                                                                                                                                                                                              |

## **5. Information on consumers and potential users**

For professional use only (e.g., pathologist, laboratory assistant).

Before using the Cellprep solutions, personnel should carefully read the Instructions for Use and complete the training.

## **6. Installation**

Not applicable.

The Cellprep Respiratory Tract Solution is completely ready to use and does not require any preparation.

## **7. Application of the Cellprep Respiratory Tract Solution**

### **7.1. Description of the medical or non-medical devices intended for use in combination with Cellprep Respiratory Tract Solution**

#### **7.1.1. Equipment for slide preparation for thin layer cytology or liquid-based cytology**

1) Cell Prep Plus LBC Processor with accessories by Biotdyne Co., Ltd., Korea (registration certificate No. ΦC3 2010/07117) is intended for fixation and transfer of cells of various body fluids and tissue samples from a vial to a slide creating a cell monolayer so that the samples are suitable for further staining and microscopy analysis. Hereafter referred to as the Cell Prep Plus.

2) Cellprep AUTO LBC System by Biotdyne Co., Ltd., Korea (registration certificate No. ΦC3 2020/11462) is an auxiliary equipment intended for preparation of a monolayered cytological sample using LBC method for its further staining and examination with the microscope for the purpose of early diagnosis of precancerous conditions, detection of oncopathology and the degree of its malignance. Hereafter referred to as the Cellprep AUTO.

The vial containing the Cellprep Respiratory Tract Solution and the sample inside must be placed in one of the above devices. Through the process, the device smears the cells on a slide producing a monolayered circle of cells for further examination. It minimizes misdiagnosis by reducing the amount of overlapping cells on a slide.

### **7.1.2. Slides**

In addition to the slides included in the delivery set of Cell Prep Plus LBC Processor, it is allowed to use the slides that have similar parameters and characteristics, such as "Microscopic Slides and Coverslips" by Leica Biosystems Richmond, Inc., USA, registration certificate No. P3H 2013/1157.

### **7.1.3. Vortex**

To shake and mix the vial with the Cellprep Respiratory Tract Solution and the sample, it is recommended to use:

1) Vortex mixer V-3, SIA "ELMI", Latvia, registration certificate No. P3H 2017/5466;

2) Personal vortex V-1 plus, SIA "Biosan", Latvia, registration certificate No. ΦC3 2011/09797;

or other devices with similar technical characteristics, approved for circulation within the Russian territory.

## **7.2. Application of the Cellprep Respiratory Tract Solution**

### **7.2.1. Sputum**

Insert sputum sample (about 3 cc) with the part of blood and discoloration (yellow) into Cellprep Respiratory Tract Solution vial, then vortex for over 5 secs and leave to fix for 20 mins.

**Only one patient sample is placed per vial.**

### **7.2.1. Bronchial washing, bronchial brushing**

1) Pour collected sample into 50ml tube and centrifuge at 2000rpm for 5 minutes.

2) Remove supernatant, then pour the sediment into Cellprep Respiratory Tract Vial and vortex for 5 sec and leave to pre-fix for more than 20 minutes.

▷ About 3cc of the sample should be put in the vial.

**Only one patient sample is placed per vial.**

## **7.3. Insertion of the vial with the Cellprep Respiratory Tract Solution into equipment for slide preparation for thin layer cytology/liquid-based cytology**

The innovative patented air transfer technology (Blowing Technology™) used in the Cellprep processors allows to get a high-quality monolayered slide and reduce the loss of cells during cytological examination to zero.


Depending on the equipment,

1. **Cellprep AUTO**: place the vial in a tray in order. As processor has auto vortexing system, there is no need to use vortex outside of processor.
2. **Cell Prep PLUS**: before placing the vial in a processor, vortex it for 10-20 seconds or shake 5-10 times.

**Cellprep AUTO**: after turning on the instrument, select the operation mode on the touch screen (Auto Turbidity Smear Mode or Manual Turbidity Mode).


Detailed description of how to operate the Cellprep AUTO LBC System is provided in the corresponding user manual.

## Auto Turbidity Smear Mode


Please click on the **SETUP** button for 3 seconds or longer in the main screen.

The screen will switch to User Setup Screen.


Select **Turbidity ON**.

Using the auto- turbidity selection.


Place the slide in the rack with the barcode facing towards the slide holder. The barcode should face the barcode scanner, and make sure the slide holder is in the correct position.


Prepare vials and load them in the tray.

Up to 40 vials can be loaded at once. The slide, completed a smear, is come out according to the order of Vial's work.


Since the equipment has an integrated vortex, there is no need to specially mix up the vials.

Put the vials from the front so no vacancy remains.


Push the tray by aligning the groove on the location of the next section when mounting the vial tray.

Select the types of specimens


Set the mode according to the type of specimen by pressing the buttons on the touch screen.

Press the purple "R" button for respiratory samples.


Specify the quantity of vials for operating.

Number press menu comes out when press red number on **Working Count**.

Enter the desired quantity of vials in the field.

Press **START** button. Cap opening, suction, blowing and slide transfer is automatically operated.

If you want to stop the operation, please push '**STOP**' button. Push the **START** button to continue with the operation. If you want to start from the beginning after stopping the operation, press **STOP** and then **HOME** to reset the position.


Check the result of smear and fix the slide.

When all smearing completes, **TOTAL VIAL PROCESSING RESULT** displays in the touch screen. If the work completed in normal, the vial displays with Green, otherwise it displays with Black.

Press the arrow on **TOTAL VIAL PROCESSING RESULT** or below of **SLIDE RESULT**, it moves to **GLASS PRECESSING RESULT** and the **SLIDE** work can be checked.


Take out the Smearing completed Slide that hold in the cell fixing alcohol container, and fix for 30 minutes in separate ethanol container then dye.

## Manual Turbidity Mode


Please click on the **SETUP** button for 3 seconds or longer in the main screen.

The screen will switch to User Setup Screen.


Select **Turbidity OFF**.

Disable auto-turbidity selection function.


Vortex ON/OFF Selection Function.


Prepare vials and load them in the tray.


Up to 40 vials can be loaded at once. The slide, completed a smear, is come out according to the order of Vial's work.

Put the vials from the front so no vacancy remains.

Push the tray by aligning the groove on the location of the next section when mounting the vial tray.

Select A, B, C on the touch screen according to the density and condition of prepared specimen.

#### Mode A:

- Used when amount of collected cells in solution vial are normal


#### Mode B:

- Used when too less amount of collected cells in solution vial

#### Mode C:


- Used when too concentrated amount of collected cells in solution vial

Select the types of specimens.


Set the mode according to the type of specimen by pressing the buttons on the touch screen.

Press the purple "R" button for respiratory samples.

Specify the quantity of vials for operating.


Number press menu comes out when press red number on **Working Count**.

Press **START** button.


Cap opening, suction, blowing and slide transfer is automatically operated.

If you want to stop the operation, please push 'STOP' button. Push the **START** button to continue with the operation. If you want to start from the beginning after stopping the operation, press **STOP** and then **HOME** to reset the position.

Check the result of smear and fix the slide.


When all smearing completes, **TOTAL VIAL PROCESSING RESULT** displays in the touch screen. If the work completed in normal, the vial displays with Green, otherwise it displays with Black.

Press the arrow on **TOTAL VIAL PROCESSING RESULT** or below of **SLIDE RESULT**, it moves to **GLASS PRECESSING RESULT** and the SLIDE work can be checked.

Take out the Smearing completed Slide that hold in the cell fixing alcohol container, and fix for 30 minutes in separate ethanol container then dye.

**Cell Prep PLUS:** after turning on the instrument, select the sample type on the touch screen. Press the "RPT" button for respiratory samples.

Before inserting the vial containing the sample into the processor, shake it well 5-10 times and place it upside down in the instrument with the vial cap facing down.

Place the vial in the holder. It is necessary to shake the vial several times before inserting it into the instrument to ensure that the cells are evenly distributed in the vial.

Label and place the slide in the slide holder unit of the processor. The smeared side of the slide should be facing down.

Select A, B, C mode on the touch screen according to the amount of cells and conditions.

# Cellprep

YYYY/MM/DD HH:MM:SS

|                                                   |  |            |  |         |  |          |  |         |  |
|---------------------------------------------------|--|------------|--|---------|--|----------|--|---------|--|
| STATUS                                            |  | WORK 12345 |  | GYNE    |  | URINE    |  | BODY    |  |
| SUCTION                                           |  | BLOWER     |  | FNA     |  | RPT      |  | CUSTOM  |  |
| 00.1 kPa                                          |  | 1.23 kg/cm |  | A       |  | B        |  | C       |  |
| 1 LINE IN                                         |  | GYNE       |  | A       |  | 12 CYCLE |  | Blodyne |  |
| MESSAGE : AUTO RUNNING (PLEASE PUSH START BUTTON) |  |            |  |         |  |          |  |         |  |
| MANUAL                                            |  | SETUP      |  | RELEASE |  | HOME     |  | START   |  |
|                                                   |  |            |  |         |  | STOP     |  |         |  |

**Mode A:** Used when amount of collected cells in solution vial are normal

**Mode B:** Used when too less amount of collected cells in solution vial

**Mode C:** Used when too concentrated amount of collected cells in solution vial

Push 'START' button on the touch screen. Vial cap opening, suction, blowing and slide out is automatically operated.

Detailed description of how to operate the Cell Prep Plus LBC Processor with accessories is provided in the corresponding user manual.

After smearing insert the slide into 95% ethanol immediately and let it fix for 30 mins.

Screen out the lesion from the smeared slide after it is stained and mounted.

Interpretation of the results is carried out according to the clinical guidelines.

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## **7.4. Interfering substances or sample limitations that may interfere with study result**

Samples with a high content of blood and mucus, and samples with an insufficient content of material in all fields of view of the slide are considered inadequate, since they may cause subsequent difficulties in making a diagnosis. If abnormal cells are present in the sample, then any smear is considered adequate.

When collecting samples of the material, it is necessary to follow established clinical guidelines and sample preparation protocols that are appropriate for the type and cellularity of the sample. Failure to comply with this requirement may result in incorrect results.


## **7.5. Identification of the medical devices in order to obtain a safe combination and (or) information on known restrictions on the combined use of medical devices**

The medical and/or non-medical devices listed in Section 7.1 may be used with Cellprep Respiratory Tract Solution to ensure a safe combination.

It is NOT ALLOWED to use the Cellprep Solution with other products than those listed in Section 7.1.

18

## 8. Specifications of the Cellprep Respiratory Tract Solution

| Name of characteristic / parameter | Reference value                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Appearance                         | The vial should not contain holes, cracks, deformations, scratches, and other defects that affect the quality of the device use.<br><br>Appearance of the solution - colorless transparent liquid without foreign impurities                                                                                                                                                                                                                                                          |
| Foreign impurities                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Fungi                              | No                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Bacteria                           | No                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Parasites                          | No                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Others                             | 1 unit / 10 ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Reference volume                   | 20 ml $\pm$ 10%                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Torque effect of a vial middle cap | 19 kgf·cm $\pm$ 5 kgf·cm                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Torque effect of a vial top cap    | 14 kgf·cm $\pm$ 5 kgf·cm                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| pH of solution                     | pH=5.2 $\pm$ 0.25                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Dimensions of vial:                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Total vial height:                 | 89.0 mm $\pm$ 1.0 mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Diameter of a vial middle cap:     | 42.7 mm $\pm$ 1.0 mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Diameter of a vial middle part:    | 36.1 mm $\pm$ 1.0 mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Diameter of a vial bottom part:    | 38.1 mm $\pm$ 1.0 mm                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Weight of vial with solution       | 45.0 g $\pm$ 10%                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Marking                            | <p>There should be no soiling, poorly typed text. Lot number, expiration date, product name should correspond to the set parameters.</p> <p>The vial label must contain the following symbols and signs:</p>  <p>The shipping package with a cardboard box must contain the following symbols and signs:</p>  |

## 9. Recommendations for staining the slides

| Step No  | Name of chemical for staining | Duration            |                     |
|----------|-------------------------------|---------------------|---------------------|
|          |                               | Gyn                 | Non Gyn             |
| Fixation | 95% Ethanol                   | 30 minutes          | 20 minutes          |
| 1        | 80% Ethanol                   | 10 dips             | 10 dips             |
| 2        | 70% Ethanol                   | 10 dips             | 10 dips             |
| 3        | Water                         | 2 minutes           | 2 minutes           |
| 4        | Water                         | 1 minute 30 seconds | 1 minute 30 seconds |
| 5        | 0.5% HCl                      | 3-5 dips            | 3-5 dips            |
| 7        | Water                         | 5 minutes           | 5 minutes           |
| 8        | 95% Ethanol                   | 10 dips             | 10 dips             |
| 9        | 95% Ethanol                   | 10 dips             | 10 dips             |
| 11       | 95% Ethanol                   | 10 dips             | 10 dips             |
| 12       | 95% Ethanol                   | 10 dips             | 10 dips             |
| 14       | 95% Ethanol                   | 10 dips             | 10 dips             |
| 15       | 95% Ethanol                   | 10 dips             | 10 dips             |
| 16       | 100% Ethanol                  | 10 dips             | 10 dips             |
| 17       | 100% Ethanol                  | 10 dips             | 10 dips             |
| 18       | Xylene                        | 10 dips             | 10 dips             |
| 19       | Xylene                        | 1 minute (10 dips)  | 1 minute (10 dips)  |

| Desired reagent  | Expected manufacture concentration | 95% Ethanol (ml) | DW (ml)  |
|------------------|------------------------------------|------------------|----------|
| 50% Ethanol      | 50%                                | 526              | 474      |
| 70% Ethanol      | 70%                                | 737              | 263      |
| 80% Ethanol      | 80%                                | 840              | 160      |
| Desired reagent  | Expected manufacture concentration | 70% Ethanol (ml) | HCl (ml) |
| 0.5% HCl-alcohol | 0.5%                               | 995              | 5        |

## 10. A. B. C. Suction Numerical Standard Values for LBC Processors

### Pressure & Time

| Suction                    | Mouth | Suction  |            | Blower   |            |
|----------------------------|-------|----------|------------|----------|------------|
|                            |       | Pressure | Time (sec) | Pressure | Time (sec) |
| RPT<br>(respiratory tract) | A     | 20.0     | 35.0       | 1.0      | 0.05       |
|                            | B     | 37.0     | 35.0       | 1.0      | 0.05       |
|                            | C     | 10.0     | 35.0       | 1.0      | 0.05       |

## **11. Cautions and warnings**

### **<Cautions>**

1. Not for internal use.
2. Do not drink or inhale vapours.
3. Avoid media contact with eyes or skin.
4. When contacting with eyes, rinse eyes immediately and thoroughly with clean water and seek medical attention.
5. Colorless volatile liquid. In case of prolonged exposure, it affects the liver and central nervous system.
6. Use chemical resistant gloves, goggles and a lab coat/apron for your safety.
7. After use, the Cellprep solutions must be handled as a biohazard material.
8. Do not use the Cellprep Solution after the expiration date.

**The Cellprep solutions should not be used in damaged or cracked vials, if the color of solution changes, or if foreign impurities are visible.**

## **12. Transportation and storage conditions**

Transportation by all means of transport is allowed.

Transportation and storage temperature: from +10°C to +30°C

Store in a closed container in an area that prevents direct contact with sunlight.

Store and use only in a well-ventilated area.

## **13. Packaging**

The Cellprep solutions are delivered in polypropylene vials. Vials (25 pcs) are placed between two plastic vial trays and covered with polyethylene film. Eight packages (25 vials each) covered with film (200 vials in total), together with the Instructions for Use, are placed in a cardboard box. Then a cardboard box with several sealing non-combustible inserts is additionally placed in a cardboard package for hazardous material and put on a pallet.

It is allowed to provide the Instructions for Use in electronic form by sending by e-mail at the consumer's request.

### **Basic information about the Cellprep Respiratory Tract Solution packaging**

| <b>Packaging</b>                                                                                        | <b>Material</b> | <b>Dimensions / quantity</b>                                                                  |
|---------------------------------------------------------------------------------------------------------|-----------------|-----------------------------------------------------------------------------------------------|
| Cellprep Respiratory Tract Solution vial<br><br>The vial is a primary packaging for Cellprep solutions. | Polypropylene   | Total vial height: 89.0 mm $\pm$ 1.0 mm<br><br>Weight of vial with solution: 45.0 g $\pm$ 10% |

|                                                                                                                                                                        |                                                                                                               |                                                                                                                                                                                                 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Plastic vial trays covered with polyethylene film<br>Secondary consumer package                                                                                        | Plastic trays (polypropylene)<br>Polyethylene film                                                            | Dimensions:<br>25.5 x 25.5 x 8.5 cm $\pm 10\%$<br>Weight of packaging with vials (gross)<br>- 1.30 kg $\pm 10\%$<br>Quantity: 25 vials                                                          |
| Shipping packaging of Cellprep Respiratory Tract Solution<br>The shipping package includes an inner cardboard box and an outer fiberboard box for hazardous materials. | Cardboard box, inner                                                                                          | Dimensions:<br>525 x 270 x 370 mm $\pm 10\%$<br>Weight of box with vials (gross):<br>11.35 kg $\pm 10\%$ (for volume of 20 ml)<br>Quantity: 8 packages covered with film<br>200 vials in total. |
| Packaging for hazardous materials (outer packaging).<br>Shipping package                                                                                               | Outer fibreboard box for hazardous materials<br>Several inserts from polymer non-combustible material or foam | Box dimensions: 530 x 320 x 450 mm $\pm 10\%$ or 590 x 395 x 470 mm $\pm 10\%$<br>Weight of box with vials (gross):<br>14.10 kg $\pm 10\%$<br>Quantity: 1 inner cardboard box with 200 vials.   |
| Pallet                                                                                                                                                                 | Plastic or wooden pallet                                                                                      | Dimensions: 110 x 110 x 132 cm $\pm 10\%$                                                                                                                                                       |

11032

## 14. Symbol Identification

| Symbol                                                                              | Description                                                                | Symbol                                                                              | Description                                         |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------|
|    | Causes serious eye irritation                                              |    | Best before...                                      |
|    | Carcinogenic substances                                                    |    | Keep dry                                            |
|    | Flammable liquid and vapours.                                              |    | Do not use hooks                                    |
|    | European Community<br>CE logo: This device complies with the EC Directives |    | Fragile, handle with care                           |
|  | Temperature range<br>(temp. +10 °C ~ +30 °C)                               |  | Top                                                 |
|  | Biohazard                                                                  |  | Manufacturer                                        |
|  | In vitro diagnostic medical device                                         |  | Authorized representative in the European Community |
|  | Manufacturer's lot number.                                                 |  | Refer to the Instructions for Use                   |

## **15. Manufacturer's responsibility**

Biodyne Co., Ltd. is a legal manufacturer responsible for the quality, efficiency and safety of the registered medical device for use, normal operation and possible negative results of using them, except for cases of negligence, incorrect use through the fault of user, or force majeure.

**In case of any adverse reactions or events while using the product, please contact customer support service or your local distributor immediately.**

## **16. Repair and Maintenance**

The medical device does not require maintenance or repair.

For additional services or support, please contact Support Services or your local representative.

### **Support Service**

Tel. +82-1522-4811

E-mail: [info@biodyne.asia](mailto:info@biodyne.asia)

### **Authorized representative in the Russian Federation:**

#### **BioLine, LLC**

**Legal address:** Russia, 197022, Saint Petersburg, municipal district Aptekarsky Ostrov, 23 Professora Popova Str., Let. E, structure 3H/5H, room 208

**Postal address:** Russia, 197022, Saint Petersburg, 23 Professora Popova Str., Let. E

E-mail: [main@bioline.ru](mailto:main@bioline.ru)

Tel.: +7(812) 320-49-49.

Hot line: 8-800-555-49-40.

## **17. Cleaning and Disinfection**

Not applicable.

The Cellprep Respiratory Tract Solution does not require cleaning and disinfection.

## **18. Warranty and shelf life**

Biodyne Co., Ltd., Korea, warrants that the device has passed a comprehensive quality control procedure based on our internal testing standards and that the device is free from defects and conforms to all specifications and/or agreed specifications.

The manufacturer warrants the high quality of the device throughout the shelf life, subject to the recommended storage conditions.

Cellprep Respiratory Tract Solution shelf life - 2 years (while observing the recommended storage conditions).

The expiry date is indicated on the label.

Use of the Cellprep Solutions after the expiry date is prohibited.

The material (sample) fixed in the vial, while keeping the temperature regime from +10 °C to +30 °C, can be used for liquid-based cytology studies for up to 3 months.°

## **19. Disposal**

**After use, the Cellprep solutions delivered by Biodyne Co., Ltd., Korea, must be handled as biohazard waste.**

Waste material should be disposed of in an approved incinerator in compliance with all federal, municipal and local government regulations.

For the Russian Federation, it is recommended to dispose of the device in accordance with the rules implemented by SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste referred to class B (epidemically hazardous waste).

## **20. Information on the initial release and revision of the Instructions for Use**

Initial release. Version 01 from 09-2023.

## **21. List of the standards applied**

a) Quality Management System applied by Biodyne Co., Ltd., Republic of Korea, is certified to the requirements of EN ISO 13485:2016 by the TUV Rheinland LGA Products GmbH certification body in relation to: "Distribution and Service of Liquid-based cytology system. Design and Development, Manufacture and Distribution of Cell preservation solution and Membrane filters".

b) The manufacturer (Biodyne Co., Ltd., Korea) complies with the principles of Good Manufacturing Practice (GMP) for medical devices. GMP Certificate is issued by Korea Conformity Laboratories.

The Cellprep Solutions comply with the following international standards:

- Directive 98/79/EC on in vitro diagnostic medical devices
- EN 980:2008 Symbols for use in the labelling of medical devices
- EN ISO 14971:2012 Medical devices. Application of risk management to medical devices

The Cellprep Solutions comply with the following Russian standards:

- GOST R 51088-2013 "In vitro diagnostic medical devices. Reagents, kits, the test - systems, control materials, culture media. Requirements to devices, and to supporting documentation";
- GOST R ISO 18113-1-2015 "In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 1. Terms, definitions and general requirements";
- GOST R ISO 18113-2-2015 "In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 2. In vitro diagnostic reagents for professional use";
- GOST R ISO 23640-2015 "In vitro diagnostic medical devices. Evaluation of stability of in vitro diagnostic reagents".
- GOST R ISO 15223-1-2023 "Medical devices. Symbols applied in medical device labelling, marking and in supporting documentation. Part 1. General requirements".

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**Biodyne Co., Ltd.**

24-4, Ahasan-ro 5-gil, Seongdong-gu, Seoul, 04793, Republic of Korea

**Tel. +82-1522-4811 Fax+82-2-468-0910**

**E-mail: [overseas@biodyne.asia](mailto:overseas@biodyne.asia)**

**[www.biodyne.asia](http://www.biodyne.asia)**

For the safe use of the device, please read the instructions for use carefully.  
Photos in this manual may differ from actual device parts.

## Appendix No. 1.

### Certificate of Analysis - Cellprep Respiratory Tract Solution (sample)

| Cellprep Respiratory Tract Solution for Liquid-Based Cytology Processor                                                                                        |                  |                                  |                                                                       |            |                          |                                      |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|----------------------------------|-----------------------------------------------------------------------|------------|--------------------------|--------------------------------------|
| Manufacturer: Biodyne Co., Ltd.<br>24-4, Ahasan-ro 5-gil, Seongdong-gu, Seoul, 04793, Republic of Korea<br>Tel. +82-1522-4811<br>E-mail: overseas@biodyne.asia |                  |                                  |                                                                       |            |                          |                                      |
| LOT No. XXXX                                                                                                                                                   |                  |                                  | Exp. date: XXXXX                                                      |            |                          |                                      |
| Certificate of Analysis No. XXXX dd. XXXXX                                                                                                                     |                  |                                  |                                                                       |            |                          |                                      |
| No                                                                                                                                                             | Parameter        | Equipment                        | Standard                                                              | Compliance | Interval                 | Specification                        |
| 1                                                                                                                                                              | External surface | Identification table<br><br>Eyes | 1. There should be no foreign impurities inside the vial.             | Complies   | At once after production | Observation from a distance of 30 mm |
|                                                                                                                                                                |                  |                                  | 2. The vial, middle cap and top cap should match.                     | Complies   |                          |                                      |
|                                                                                                                                                                |                  |                                  | 3. Lot number must be indicated correctly.                            | Complies   |                          |                                      |
|                                                                                                                                                                |                  |                                  | 4. The solution must not leak away.                                   | Complies   |                          |                                      |
|                                                                                                                                                                |                  |                                  | 5. There should be no cracks.                                         | Complies   |                          |                                      |
|                                                                                                                                                                |                  |                                  | 6. Solution and label should match.                                   | Complies   |                          |                                      |
|                                                                                                                                                                |                  |                                  | 7. There should be no foreign impurities on the vial exterior.        | Complies   |                          |                                      |
|                                                                                                                                                                |                  |                                  | 8. Check that the language on label and type are indicated correctly. | Complies   |                          |                                      |

|                                                                                                                                                                                                                                                                         |                    |                    |                                                                                      |          |  |                   |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|--------------------|--------------------------------------------------------------------------------------|----------|--|-------------------|
| 2                                                                                                                                                                                                                                                                       | Color              | Eyes               | 1. Content should be transparent.                                                    | Complies |  | Visual inspection |
| 3                                                                                                                                                                                                                                                                       | Foreign impurities | Eyes               | 1. There should be no foreign impurities.                                            | Complies |  | Visual assessment |
| 4                                                                                                                                                                                                                                                                       | Volume             | Graduated cylinder | 1. Volume of the solution should be determined by the volume for each type. (□ 20ml) | Complies |  | Measurement       |
| 5                                                                                                                                                                                                                                                                       | pH of solution     | pH-meter           | 5.2±0.25                                                                             | Complies |  | Measurement       |
| Transportation and storage temperature: from +10°C to +30°C. Store in a closed container in an area that prevents direct contact with sunlight. Store and use only in a well-ventilated area.<br>Shelf life - 2 years (while observing recommended storage conditions). |                    |                    |                                                                                      |          |  |                   |

XXXX – variable data.

100

## MSDS (MATERIAL SAFETY DATA SHEET)

|                       |         |
|-----------------------|---------|
| PRODUCT               | Rev 6.3 |
| Cellprep LBC Solution |         |

## | 1. | CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

- A. PRODUCT NAME: Cellprep LBC Solutions (Cellprep Respiratory Tract Solution)
- B. APPLICATION: Ethanol based, Medical preservation liquid for cell cultures.
- C. COMPANY IDENTIFICATION      BIODYNE CO.,LTD.  
 24-4, Ahasan-ro 5-gil, Seongdong-gu,  
 Seoul, 04793, Republic of Korea  
 Phone: +82-1522-4811 / Fax: +82-70-4903-0901

## | 2. | HAZARDS IDENTIFICATION

## A. NATIONAL FIRE PROTECTION ASSOCIATION (NFPA) CLASSIFICATION

HEATH HAZARDS = 0~1

FLAMMABILITY = 0

REACTIVITY = 0

SPECIFIC HAZARD = N/A

Rating: 0 = Insignificant 1 = Slight 2 = Moderate 3 = High 4 = Extreme

## PRECAUTIONARY STATEMENTS

PREVENTION

Obtain special instructions before use (P201).

Do not handle until all safety precautions have been read and understood (P202).

Keep container tightly closed (P233).

Ground/bond container and receiving equipment (P240).

Take precautionary measures against static discharge (P243).

Do not breathe dust/fume/gas/mist/vapors/spray (P260).

Wash eye thoroughly after handling (P264).

Do not eat, drink or smoke when using this product (P270).

Wear protective gloves/protective clothing/eye protection/face protection (P280).

Use personal protective equipment as required (P281).

RESPONSE

IF ON SKIN (or hair): Rinse skin with water/shower.

In case of fire: Use manufacturer/supplier or the competent authority to specify appropriate media for extinction. IF

IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

If eye irritation persists: Get medical advice/attention.

IF exposed or concerned: Get medical advice/attention.

IF exposed: Call a doctor/physician.

Get medical advice/attention if you feel unwell.

STORAGE

Store in cool/well-ventilated place.

Store locked up.

DISPOSAL

Dispose of contents/container in accordance with local/regional/national/international regulation.

## | 3. | COMPOSITION / INFORMATION ON INGREDIENTS

**MSDS (MATERIAL SAFETY DATA SHEET)**

|                       |         |
|-----------------------|---------|
| PRODUCT               | Rev 6.3 |
| Cellprep LBC Solution |         |

| CHEMICAL NAME                     | CAS No.   | ENCS Number* | %(w/w) |
|-----------------------------------|-----------|--------------|--------|
| Ethanol                           | 64-17-5   | 2-202        | 34%    |
| Di-water                          | 7732-18-5 | -            | 49%    |
| Glycerine                         | 56-81-5   | -            | 16%    |
| Phosphorous acid                  | 7664-38-2 | -            | 0.48%  |
| Ammonium acetate                  | 631-61-8  | -            | 0.48%  |
| Sodium Citrate Tribasic Dihydrate | 6132-04-3 | -            | 0.04%  |

\* ENCS Number: Japanese Existing and New Chemical Substances.

**4. FIRST AID MEASURES****A. EYE CONTACT**

Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, get medical advice/attention.

**B. SKIN CONTACT**

Rinse skin with water/shower or wash with plenty of soap and water. If skin irritation or rash occurs, seek medical advice/attention. Wash contaminated clothing before reuse.

**C. INHALATION**

Not a likely route of exposure. If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing. If experiencing respiratory symptoms or feel unwell, call a doctor/physician.

**D. INGESTION**

Not a likely route of exposure. Do not induce vomiting without medical advice. There may be irritation to the gastrointestinal tract with nausea and vomiting.

**E. MAJOR TOXIC SYMPTOMS AND EFFECTS**

Inhalation of high concentration of vapor may cause mild irritation of the respiratory tract. Frequent or prolonged contact with product may defect and dry the skin, leading to discomfort and dermatitis.

**F. NOTE TO PHYSICIAN**

Prolonged or frequently repeated skin contact may aggravate an existing dermatitis condition in some individuals. Based on the individual reactions of the patient, the physician's judgment should be used to control symptoms and clinical condition.

**5. FIRE FIGHTING MEASURES****A. SUITABLE EXTINGUISHING MEDIA**

General extinguishing media is usable.

**B. SPECIFIC HAZARDS ARISING FROM THE CHEMICAL**

No cause of harmful hazards arose from the chemical.

**C. SPECIAL PROTECTIVE EQUIPMENT AND PRECAUTIONS FOR FIRE FIGHTERS**

In case of fire, wear a full face positive-pressure self-contained breathing apparatus and protective suit. Keep containers cool with water spray.

**6. ACCIDENTAL RELEASE MEASURES****A. PERSONAL PRECAUTIONS, PROTECTIVE EQUIPMENT AND EMERGENCY PROCEDURES**

Use personal protective equipment recommended in Section 8 (Exposure Controls/Personal Protection). Stop or reduce any leaks if it is safe to do so. Have emergency equipment (for fires, spills, leaks, etc.) readily available.

**B. ENVIRONMENTAL PRECAUTIONS**

Prevent material from entering sewers or waterways.

**C. METHODS AND MATERIALS FOR CONTAINMENT AND CLEANING UP**

**MSDS (MATERIAL SAFETY DATA SHEET)**

|                              |         |
|------------------------------|---------|
| PRODUCT                      | Rev 6.3 |
| <b>Cellprep LBC Solution</b> |         |

Clean with towel or paper towel on its spill.

**| 7. | HANDLING AND STORAGE****A. PRECAUTIONS FOR SAFE HANDLING**

Do not get in eyes, on skin, on clothing. Do not take internally. Use with adequate ventilation. Try not to breathe vapors/gases/dust. Keep the containers closed when not in use.

**B. CONDITIONS FOR SAFE STORAGE, INCLUDING ANY INCOMPATIBILITIES**

Store in suitable labeled containers. Store the containers tightly closed. Store away from heat. Keep containers placed in cool, well-ventilated areas. Have appropriate fire extinguishers available in and near the storage area.

**| 8. | EXPOSURE CONTROLS/PERSONAL PROTECTION****A. OCCUPATIONAL EXPOSURE LIMITS****KOREA**

Korea. OELs (ISHL Article 42; MOL Public Notice No. 1986-45, as amended through MOL Public Notice No. 2007-25, June 8, 2007)

CAS RN: 64-17-5

Name: Ethanol

The 8-hour TWA: 1,000 ppm

The 8-hour TWA: 1,900 mg/m<sup>3</sup>

The 15-minute STEL: - ppm

The 15-minute STEL: - mg/m<sup>3</sup>

**JAPAN**

Japan. OELs - JSOH (Japan Society of Occupational Health: Recommendation of Occupational Exposure Limits, 2007)

CAS RN: 64-17-5

Name: Ethanol

Not established.

**OSHA**

OSHA Table Z-1 Limits for Air Contaminants (June 30, 1993) (29 CFR 1910.1000) (1971 Permissible Exposure Limits (PELs))

CAS RN: 64-17-5

Name: Ethanol

The 8-hour PEL: 1,000 ppm

The 8-hour PEL: 1,900 mg/m<sup>3</sup>

**NIOSH**

NIOSH. Pocket Guide to Chemical Hazards, 2005

CAS RN: 64-17-5

Name: Ethanol

NIOSH Recommended exposure limit (REL) [for up to a 10-hour workday during a 40-hour workweek]: 1,000 ppm.

NIOSH Recommended exposure limit (REL) [for up to a 10-hour workday during a 40-hour workweek]: 1,900 mg/m<sup>3</sup>.

NIOSH Immediately dangerous to life or health (IDLH) concentration: 3,300 ppm

**B. ENGINEERING MEASURES**

General ventilation is recommended. Use local exhaust ventilation if necessary to control airborne mist and vapor. Provide mechanical ventilation of confined spaces.

**C. INDIVIDUAL PROTECTION MEASURES****GENERAL ADVICE**

The use and choice of personal protection equipment is recommended when use of the product, the workplace and the way the product is handled. In general, we recommend as a minimum precaution as wearing gloves.

## MSDS (MATERIAL SAFETY DATA SHEET)

|                       |         |
|-----------------------|---------|
| PRODUCT               | Rev 6.3 |
| Cellprep LBC Solution |         |

### RESPIRATORY PROTECTION

If significant mists, vapors or aerosols are generated an approved respirator is recommended. A suitable filter material depends on the amount and type of chemicals being handled. If respiratory protection is required, institute a complete respiratory protection program including selection, fit testing, training, maintenance and inspection.

### EYE PROTECTION

When handling this product, the use of splash chemical goggles can be used when necessary.

### HAND PROTECTION

Wearing appropriate chemical resistant gloves are recommended.  
Regular glove replacement is recommended.

### SKIN PROTECTION

No specific harm on skin. Wearing appropriate chemical resistant cloths are recommended.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

|                                               |                            |
|-----------------------------------------------|----------------------------|
| A. PHYSICAL APPEARANCE:                       | Clear, colorless liquid.   |
| B. ODOR:                                      | Mild Characteristic odor.. |
| C. ODOR THRESHOLD:                            | No data available.         |
| D. pH:                                        | No data available.         |
| E. MELTING/FREEZING POINT:                    | No data available.         |
| F. INITIAL BOILING POINT/BOILING POINT RANGE: | No data available.         |
| G. FLASH POINT:                               | About 30°C *               |
| H. EVAPORATION RATE:                          | No data available.         |
| I. FLAMMABILITY (GAS, SOLID):                 | Not applicable.            |
| J. UPPER AND LOWER EXPLOSION LIMITS:          | No data available          |
| K. VAPOUR PRESSURE:                           | No data available.         |
| L. SOLUBILITY:                                | Soluble in water.          |
| M. RELATIVE VAPOUR DENSITY:                   | No data available.         |
| N. SPECIFIC GRAVITY:                          | No data available.         |
| O. n-OCTANOL/WATER PARTITION COEFFICIENT:     | No data available.         |
| P. AUTOIGNITION TEMPERATURE:                  | No data available.         |
| Q. DECOMPOSITION TEMPERATURE:                 | No data available.         |
| R. VISCOSITY:                                 | No data available.         |

\* Sample was not flashed until boiled (about 30°C)

## 10. STABILITY AND REACTIVITY

### A. CHEMICAL STABILITY

Stable under normal conditions.

B. POSSIBILITY OF HAZARDOUS REACTIONS Hazardous polymerization will not occur.

### C. CONDITIONS TO AVOID

No specific conditions to avoid.

### D. INCOMPATIBLE MATERIALS

No materials to be incompatible. \_\_\_\_\_

## 11. TOXICOLOGICAL INFORMATION

## MSDS (MATERIAL SAFETY DATA SHEET)

|                       |         |
|-----------------------|---------|
| PRODUCT               | Rev 6.3 |
| Cellprep LBC Solution |         |

### A. INFORMATION ON THE LIKELY ROUTES OF EXPOSURE

EYE CONTACT: Yes.  
SKIN CONTACT: Yes.  
INGESTION: No.  
INHALATION: No.

### B. SYMPTOMS RELATED TO THE PHYSICAL, CHEMICAL AND TOXICOLOGICAL CHARACTERISTICS

No data available.

#### ACUTE TOXICITY:

- Oral: None known. No data available.
- Skin: None known. No data available.
- Inhalation: None known. No data available.

#### SKIN CORROSION/IRRITATION:

Can cause mild irritation to skin. Frequent or prolonged contact with product may defect and dry the skin, leading to discomfort and dermatitis.

#### SERIOUS EYE DAMAGE/IRRITATION:

Can cause mild irritation to eye.

#### RESPIRATORY OR SKIN SENSITIZATION:

This product is not expected to be a sensitizer.

#### CARCINOGENICITY:

ACGIH Carcinogen Category: A4 (not classifiable as to carcinogenicity in humans).  
None of the substances in this product are listed as carcinogens in humans by the International Agency for Research on Cancer (IARC), American Conference of Governmental Industrial Hygienists (ACGIH).

#### ASPIRATION HAZARD:

None known. No data available.

## | 12 | ECOLOGICAL INFORMATION

\* No toxicity studies have been conducted on this product. The following results are for the hazardous components.

### A. AQUATIC AND ECOLOGIC TOXICITY

Occurrence of ecologic toxicity is possible when exposed or deposited in soil for long period of time.

### B. PERSISTENCE AND DEGRADABILITY

The organic portion of the substances is expected to be biodegrade to a moderate extent.

### C. BIOACCUMULATIVE POTENTIAL

The product not expected to significantly bioaccumulation due to low bioaccumulation factor (BCF) and octanol/water partition coefficient.

### D. MOBILITY IN SOIL

None known. No data available.

### E. OTHER ADVERSE EFFECTS

None known. No data available.

## | 13. | DISPOSAL CONSIDERATIONS

### A. DISPOSAL METHODS

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## MSDS (MATERIAL SAFETY DATA SHEET)

|                       |         |
|-----------------------|---------|
| PRODUCT               | Rev 6.3 |
| Cellprep LBC Solution |         |

Hazardous wastes must be transported by a licensed hazardous waste transporter and disposed of or treated in a properly licensed hazardous waste treatment, storage, disposal or recycling facility. Consult local, state, and federal regulations for specific requirements.

### B. PRECAUTIONS (ON DISPOSAL OF COTAMINATED CONTAINERS AND PACKAGES)

Do not dispose of wastes in local sewer or with normal garbage. Combustion products are carbon monoxide, carbon dioxide and water.

## | 14. | TRANSPORT INFORMATION

- A. UN No.: 1170.
- B. Proper Shipping Name: ETHANOL SOLUTION
- C. CLASS: 3
- D. PACKING GROUP. III
- E. Not restricted under IATA regulations

## | 15. | REGULATORY INFORMATION

### INTERNATIONAL CHEMICAL CONTROL LAWS

#### EUROPE

The disclosed substances in this product are included in or exempted from the EINECS or ELINCS inventories. U.S.

The disclosed substances in this product are included on or exempted from the TSCA 8(b) Inventory (40 CFR 710).

#### JAPAN

The disclosed substances in this product are comply with the Law Regulating the Manufacture and Importation of Chemical Substances and are listed on the Ministry of Economy, Trade and Industry (METI).

#### KOREA

Observational Chemicals List - None of the components of this product are regulated under TCCL.

## | 16. | OTHER INFORMATION

This product material safety data sheet provides health and safety information. The above information is believed to be correct but does not propose to be all inclusive and shall be used only as a guide. BIODYNE CO.,LTD. shall not be held liable for any damage resulting from handling or from contact with the above product. Each individual should make a determination as to the suitability of the information for their particular purpose(s). Various government agencies may have specific regulations regarding the transportation, handling, storage, use, or disposal of this product which may not be covered by this MSDS. The user is responsible for full compliance.

### A. REFERENCES

Regulated chemicals info. produced by the National Institute of Technology and Evaluation, Japan.

Chemicals Evaluation and Research Institute, Japan (CERI)

WHO/IPCS : International Chemical Safety Cards (ICSC).

EU European Chemicals Bureau (ECB): International Uniform Chemical Information Database (IUCLID). Ariel

WebInsight DB (3E Company).

Korea Testing & Research Institute

### B. OTHER INFORMATION

- Date issued: April 01, 2024
- Version number: 6.3
- Technical Support and Revised from Attached Research Center of Biodyne Co., Ltd.

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As a result of checking at my office, I have  
found that the attached INSTRUCTIONS  
FOR USE copy exactly corresponds with  
the original.

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This is hereby attested on this 29th day of  
May, 2024 at this office.

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Name of the office  
HANMI LAW AND NOTARY OFFICE

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
Belong to  
Seoul Central District Prosecutors' Service

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Address of the office  
Twintree Tower(B) #101,6 Yulgok-ro,  
Jongno-gu, Seoul, Republic of Korea

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Notary Public JUNG HYUNG KEUN

    
(Signature of Notary Public)

This office has been authorized by the  
Minister of Justice, the Republic of Korea,  
to act as Notary Public since March 14,  
2007 under Law No. 7428.

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