12/2019

TIMBRE
DEL ESTADO
BioSystems

Barcelona, July, 20, 2020

FE 3575198

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I, Laura López, Regulatory Affairs Manager at BioSystems S.A., hereby certify, that attached documents are true and accurate copies of User Manual of the product and translation of this document into Russian language in accordance with the documents.

Also, I hereby certify that the attached documents are for purposes of product registration in Russia and in the countries of Eurasian Economic Union.

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Laura López
Regulatory Affairs Manager
BIOSYSTEMS S.A.

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BIOSYSTEMS, S.A.
Costa Brava, 30,
08030 BARCELONA
Tel. +34 933 11 00 00

Regulatory Affairs Department
BIOSYSTEMS S.A.
llopez@biosystems.es
Barcelona, July 20th, 2020

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This reage

CLINICAL

Urea is sy
in the urin
Elevated
protein ca
or treatme
Post-rena
prostatic
variability
Clinical di
both clinic

PRINCIPL

Urea in the
can be me

CONTENTS

A. Reag
B. Reag

COMPOSITI

A. Reagent:
dehydroge

B. Reagent: M

WARNING

P301+P311;
unwell. P33

STORAGE AN

Store at 2-8°C.

Components ar

well closed and

On board stabili

stable 30 days.

Indications of

Parameters".

WARNING AND

Exercise the nor

available for prof

with local guidel

reported to BioS

ADDITIONAL M

Biochemistry Cal

cod. 18044).

REAGENT PREF

Reagents are pro

SAMPLES

Serum, plasma or

Urea in serum

anticoagulant⁵

Urea in urine is st

CALIBRATION

A reagent blank s

lot change or as r

QUALITY CONTR

It is recommende

cod. 18042) and II

18054 and cod. 18

Each laboratory s

corrective action if

UREA/BUN-UV

COD 21516 8 x 60 mL + 8 x 15 mL	COD 23516 4 x 60 mL + 4 x 15 mL
Only for <i>in vitro</i> use in the clinical laboratory	

BioSystems

UREA/BUN-UV
UREASE/GLUTAMATE DEHYDROGENASE

INTENDED USE

Reagent for the measurement of urea concentration in human serum, plasma or urine. The obtained values are useful as an aid in the diagnosis and treatment of chronic renal failure, and to evaluate the function of renal glomerulus.

This reagent is for use in the BioSystems BA analyzers.

CLINICAL SIGNIFICANCE

Urea is synthesized in the liver as a by-product of the deamination of amino acids. Its elimination in the urine represents the major route for nitrogen excretion.


Elevated urea concentration in plasma is found as a result of a high-protein diet, increased protein catabolism, after a gastrointestinal hemorrhage, mild dehydration, shock and heart failure or treatment with glucocorticoids (pre-renal uremia)^{1,2}.

Post-renal uremia is caused by conditions that obstruct urine outflow: nephrolithiasis, tumor or prostatic hypertrophy. The usefulness of urea as an indicator of renal function is limited by the variability of its plasma concentration as a result of nonrenal factors^{1,2}.

Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE OF THE METHOD

Urea in the sample consumes, by means of the coupled reactions described below, NADH that can be measured by spectrophotometry^{3,4}.

CONTENTS

	COD 21516	COD 23516
A. Reagent	8 x 60 mL	4 x 60 mL
B. Reagent	8 x 15 mL	4 x 15 mL

COMPOSITION

A. Reagent: Tris 100 mmol/L, 2-oxoglutarate 5.6 mmol/L, urease > 140 U/mL, glutamate dehydrogenase > 140 U/mL, ethyleneglycol 220 g/L, sodium azide 0.95, pH 8.0.

B. Reagent: NADH 1.5 mmol/L, sodium azide 9.5 g/L.

WARNING: H302 Harmful if swallowed. EUH031: Contact with acids liberates toxic gas. P301+P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. P330 Rinse mouth.

STORAGE AND STABILITY

Store at 2-8°C.

Components are stable once opened until the expiry date marked in the label if they are stored well closed and care is taken to prevent contamination during their use.

On board stability: Reagents open and kept in the refrigerated compartment of the analyzer are stable 30 days.

Indications of deterioration: Absorbance of the blank below the limit indicated in "Test Parameters".

WARNING AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

ADDITIONAL MATERIALS REQUIRED (NOT PROVIDED)

Biochemistry Calibrator (BioSystems cod. 18011) or Biochemistry Calibrator Human (BioSystems cod. 18044).

REAGENT PREPARATION

Reagents are provided ready to use.

SAMPLES

Serum, plasma or urine collected by standard procedures.

Urea in serum or plasma is stable for 7 days at 2-8°C. Heparin is recommended as anticoagulant⁵.

Urea in urine is stable for 2 days at room temperature if microbial growth is prevented⁵.

CALIBRATION

A reagent blank should be done every day and a calibration at least every 30 days, after reagent lot change or as required by quality control procedures.

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum level I (cod. 18005, cod. 18009 and cod. 18042) and II (cod. 18007, cod. 18010 and 18043) and the Biochemistry Control Urine (cod. 18054 and cod. 18066) to verify the accuracy of the measurement procedure.

Each laboratory should establish its own internal Quality Control scheme and procedures for corrective action if control results are not within the acceptable limits.

REFERENCE VALUES

Serum and plasma¹: 12.8 - 42.8 mg/dL urea = 6 - 20 mg/dL BUN = 2.14 - 7.14 mmol/L urea. Concentrations in the neonatal period are lower, and in adults over 60 years of age are higher than in adults. Concentrations also tend to be slightly higher in males than in females.

Urine¹: 26 - 43 g/24-h urea = 12 - 20 g/24 h BUN = 428 - 714 mmol/24-h urea.

These ranges are given for orientation only, each laboratory should establish its own reference ranges.

METROLOGICAL CHARACTERISTICS

The metrological characteristics described below have been obtained using a BA400 analyzer and following the guidelines of the Clinical & Laboratory Standards Institute (CLSI).

- Detection limit: 3.69 mg/dL urea = 1.72 mg/dL BUN = 0.614 mmol/L urea.
- Linearity limit: 300 mg/dL urea = 140 mg/dL BUN = 50 mmol/L urea. For samples with higher values, dilute manually or refer to the Test Parameterization for Automatic dilution (note that all these samples will be diluted with the same dilution ratio).
- Precision:

Serum. Mean urea concentration	Repeatability (CV)	Within-laboratory (CV)
26.8 mg/dL = 4.47 mmol/L	3.5 %	5.0 %
137 mg/dL = 22.9 mmol/L	1.1 %	1.7 %

Urine. Mean urea concentration	Repeatability (CV)	Within-laboratory (CV)
1291 mg/dL = 215 mmol/L	3.1 %	4.3 %
1771 mg/dL = 295 mmol/L	2.9 %	3.1 %

- Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

LIMITATIONS OF THE PROCEDURE

- Interferences: bilirubin (up to 30 mg/dL), hemolysis (hemoglobin up to 500 mg/dL) and lipemia (triglycerides up to 1625 mg/dL) do not interfere. Other drugs and substances may interfere⁶.

BIBLIOGRAPHY

- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2012.
- Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AAC Press, 2001.
- Talke H and Schubert GE. Enzymatische Hamstoffbestimmung in blut und serum im optischen test nach Warburg. Klinische Wochenschrift 1965, 43: 174-175.
- Gutmann I, Bergmeyer HU. Methods of enzymatic Analysis, ed Bergmeyer HU. Academic Press, NY, 1974, 4: 1794-1798.
- World Health Organization (WHO). Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1, Rev 2, 2002.
- Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AAC Press, 2000.

TEST PARAMETERS

R1: use Reagent A

R2: use Reagent B

	BA200	BA400
GENERAL		
Name	UREA-BUN-UV	UREA-BUN-UV
Short name	UREA	UREA
Sample type	serum / plasma / urine	serum / plasma / urine
Analysis mode	fixed-time bireagent	fixed-time bireagent
Unit	mg/dL	mg/dL
Decimals	0	0
Reaction type	decreasing	decreasing
PROCEDURE		
Reading mode	monochromatic	monochromatic
Main filter	340	340
Reference filter	-	-
Sample	2	2
Vol. R1	240	240
Vol. R2	60	60
Reading 1 (cycle)	20	40
Reading 2 (cycle)	23	45
Dilution factor	- / - / 50	- / - / 50
Reduced factor	10	10
CALIBRATION AND BLANK		
Blank type	distilled water	distilled water
Calibration mode	experimental calibrator	experimental calibrator
Number of calibrators	1	1
Calibration curve	-	-
OPTIONS		
Blank absorbance limit	1.100	1.100
Kinetic blank limit	-	-
Linearity limit	300 / 300 / 15000	300 / 300 / 15000
Substrate depletion	-	-

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DISPOSAL AND DESTRUCTION OF MEDICAL DEVICES

National standards must be followed when transferring and disposing of the product, its components and packaging. Chemicals and finished and used reagents must be disposed of in accordance with regional regulations.

LIABILITY AND GUARANTEES

Any change to the set and/or exchange or replacement of any components of different batches from one set to another may negatively affect the expected results and the validity of the analysis as a whole. In the event of such changes and/or substitutions, any requirement to replace the set is invalid. Claims submitted due to incorrect interpretation of the analysis results by the client are also invalid.

However, in case of any claims, the manufacturer's liability should not exceed the cost of the kit. Damage caused to the kit during transportation is not the responsibility of the manufacturer. The manufacturer's warranties of BioSystems S. A. apply to the fact that this product meets the specifications stated on the label, but only if the product is used by qualified personnel in accordance with the instructions set out on the label and the materials are used before the expiration date indicated on the label.

MAINTENANCE

No requirement

MAINTENANCE SERVICE

Not required

POTENTIAL CONSUMER

Professional user: a laboratory specialist who has the appropriate work permit

INFORMATION ABOUT THE PRESENCE OF A MEDICINAL PRODUCT IN A MEDICAL DEVICE

Does not contain

INDICATIONS FOR USE

Not applicable

CONTRAINDICATIONS TO USE

Not applicable

STERILITY

It is delivered in non-sterile form. Sterility is not required

INFORMATION ABOUT THE POSSIBILITY OF REUSE

Not re-usable

SYMBOLS and LABELS

Symbol	Name
	See manual
	Catalogue number
	Medical device <i>in vitro</i>
	Lot
	European Community symbol
	Temperature
	Expire before
	Manufacturer

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AUTHORISED REPRESENTATIVE of the MANUFACTURER

"MD-Consulting and Development" LLC, 129337, Russian Federation, Moscow, Yaroslavskoye Shosse, 6/2, office 106

MANUFACTURER AND ADDRESS OF PRODUCTION

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