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APPROVED

Deputy Executive Officer,
Aurena Laboratories AB

 Anders Bared

January 10.2023


INSTRUCTION FOR USE OF THE MEDICAL DEVICE

Product for washing and irrigation of the nasal cavity for children and adults
Aqualor® Active Forte

Aurena Laboratories AB, Sweden

1. Name of the medical device
2. General properties
3. Purpose
4. Therapeutic indications
5. Contraindications
6. Side effects
7. Precautions
8. Methods and means of cleaning
9. Procedure of use
10. Conditions of use, storage, and transportation
11. Disposal procedure
12. Manufacturer / Manufacturing site
13. Authorized representative of the manufacturer in the Russian Federation
14. Warranty
15. Explanation of symbols on the package

SIGNED BEFORE MF

Christian Nilsson
Notary Public, City of Karlstad

Identification: ☐ driver's license ☐ passport ☐ identity card

Number: n/a ☒ Bank ID

2023-10-20

1. Name of the medical device

Product for washing and irrigation of the nasal cavity for children and adults Aqualor® Active Forte

Dosage forms:

1. Product for washing and irrigation of the nasal cavity for children and adults Aqualor® Active Forte Mini 50 mL;
2. Product for washing and irrigation of the nasal cavity for children and adults Aqualor® Active Forte 150 mL

How supplied:

Option 1

- Product for irrigation and rinsing of the nasal cavity for children and adults Aqualor® Active Forte Mini 50 mL in an individual package (can) – 1 pc.;
- Instruction for Use – 1 pc.;
- Carton – 1 pc.

Option 2

- Product for irrigation and rinsing of the nasal cavity for children and adults Aqualor® Active Forte 150 mL in an individual package (balloon) – 1 pc.;
- Instruction for Use – 1 pc.;
- Carton – 1 pc.

2. General properties

Name of the parameter	Description
Device description	The product for washing and irrigation of the nasal cavity for children and adults Aqualor® Active Forte Mini 50 mL, Aqualor® Active Forte 150 mL (hereinafter, product, medical device) consists of a valve-package system containing a solution for irrigation and washing of the nasal cavity, in an aerosol can with an anatomic spray nozzle for the nose. The can with a bag inside is under pressure. The can contains propellant gas around the bag, thus separated from the liquid to be sprayed. The gas always remains inside the can. Anatomic spray nozzle tightly fixed on the can makes it possible to conveniently position the fingers on it to adjust the pressing force during application. In addition, the anatomic spray nozzle has the soft rounded form to minimize the effect of the nozzle surface contact with the nasal cavity walls. The can is equipped with a valve which allows regulating the delivery of the product. This is not a medicinal product.
Filling volume	50 mL, 150 mL
Description of the device content	Transparent, colorless, odorless liquid with a slightly salty taste, which is a seawater solution (64–80% of natural seawater) with hypertonic NaCl concentration.
Content composition	Hypertonic seawater solution (2.1% NaCl concentration), containing natural seawater, purified water, carbon dioxide (0.4% CO₂ concentration). pH 3.5–5.5. Osmolality of 640–720 mOsm/L. Preservative free.
Technical parameters	Type of spraying: “shower” (continuous fine-dispersed); pressure of 5.8 bar (0.58 MPa).
Sterilization	Sterilization by radiation. Remains sterile for the entire period of use.

Package	The medical device in an individual package (balloon) is placed in consumer packaging (carton) with the instruction for use.
----------------	---

3. Purpose

The product is designed for prevention and complex treatment of ENT diseases in children aged 2 years and adults by irrigation and washing of the nasal cavity.

The product is designed to remove excess mucus, reduce secretions, soften and remove crusts, moisturize the nasal cavity mucosa. The product has a positive effect on the nasal cavity and nasopharyngeal mucosa, helps reduce nasal congestion and restore nasal breathing. Provides the most thorough washing of the nasal cavity, facilitates removal of any pathogens (viral, bacterial, fungal), as well as allergens and dust particles from the nasal cavity and nasopharyngeal mucosa due to the "shower" spraying and adding carbon dioxide (CO₂) bubbles to the seawater solution, without disturbing the biocenosis (normal microflora) of the nasal cavity mucosa. CO₂ bubbles also contribute to the destruction of bacterial biofilms capable to protect themselves from antibiotics (causing antibiotic resistance), antiseptics, human immune system factors, and cause acute and chronic inflammatory diseases. The salts contained in seawater contribute to the liquefaction of mucus, normalization of its production in mucosal goblet cells, hamper formation and maturation of new biofilms, which reduces the bacterial load on the nasal cavity mucosa and increases the efficacy of antibacterial drugs used. Calcium and magnesium ions contained in the salt significantly improve the functioning of the nasal mucosa ciliated epithelium cells, normalize the rheological properties of mucus, increase the resistance against the introduction of viruses and bacteria. Iodine and sodium chloride have antiseptic properties. Zinc and selenium ions create conditions for the activation of the nasal mucosa and paranasal sinuses local immunity. Hypertonic seawater solution reduces nasal mucosa edema due to the removal of excess fluid from the intercellular space due to the difference in osmotic pressures.

After washing the nasal cavity with the product, the therapeutic efficacy of the drugs applied to the nasal mucosa increases, the duration of respiratory diseases is reduced, the risk of infection spreading into the paranasal sinuses and the ear cavity is decreased.

The product is intended for individual use in one patient, in order to avoid infectious contamination. Conditions for use: at home and in the medical institutions.

Area of application: otorhinolaryngology.

4. Therapeutic indications

The product is used for:

- prevention and complex treatment of acute and chronic inflammatory diseases of the nasal cavity, paranasal sinuses and rhinopharynx (infectious, allergic, atrophic):
 - acute and chronic rhinitis (running nose);
 - acute and chronic sinusitis (including maxillary sinusitis);
 - acute and chronic adenoiditis;
 - allergic and vasomotor rhinitis;
- prevention and complex treatment of respiratory infections (influenza, SARS, etc.);
- prevention and complex treatment after surgical interventions in the nasal cavity and paranasal sinuses;
- preparation of the nasal cavity mucosa for the use of nasal drugs.

The product is indicated for children over 2 years of age and adults.

5. Contraindications:

Individual intolerance of the can contents components.

Children under 2 years of age.

Special warnings:

Pregnancy and breastfeeding are not contraindications.

6. Side effects

No side effects revealed.

7. Precautions

Inspect the product and package before use. Do not use the medical device in case of damage to the integrity of the package.

To ensure the product good quality and dosing system operation, do not disassemble the medical device (e.g., do not remove the anatomic nozzle from the can).

The can is under pressure: may burst if heated.

- Do not puncture or heat, even after use.
- Keep away from heat, hot surfaces, sparks, open flames, and other sources of ignition. Do not smoke during application.

- Do not expose to the temperature above 50 °C.

- Do not expose to sunlight.

- Keep out of the reach of children.

8. Methods and means of cleaning

The product is intended for individual use. After each application, clean the anatomic nozzle with soap and warm water, then rinse.

9. Procedure of use

Wash the nasal cavity as needed, 2–4 times a day on average, and more often if necessary. Duration and frequency of treatment are not limited.

Important! During the procedure, the person subjected to washing should not breathe in, talk, or shout.

For children from 2 years

Perform nasal washing in children in an upright position of the body, preferably sitting on the arms of an adult. In this case, the child's head should be placed straight and be slightly tilted forward (Fig. 1). Ensure that the child does not throw back his head during the wash.

1. Remove the protective cap before use. At the first application, it may be necessary to press the anatomical nozzle several times before spraying the product. Insert the anatomical nozzle into the right nasal passage (Fig. 2).

2. Press the nozzle and hold it for several seconds, thus washing the right nasal passage (Fig. 3).

3. Remove residual flushing fluid by blowing it out.

4. Repeat the procedure with the left nasal passage.

5. if necessary, repeat the procedure.

Fig. 1.

Fig. 2.

Fig. 3.

For adults

Perform nasal washing in adults with the body upright (sitting or standing). In this case, the head should be located straight and be slightly tilted forward. Make sure that during washing, the head is not thrown back (Fig. 4).

1. Remove the protective cap before use. At the first application, it may be necessary to press the anatomical nozzle several times before spraying the product. Insert the anatomical attachment into the right half of the nose (Fig. 5).
2. Press the nozzle and hold it for a few seconds, thus washing the right nasal passage (Fig. 6).
3. Remove residual flushing fluid by blowing it out (Fig. 7).
4. Repeat the procedure with the left nasal passage.
5. if necessary, repeat the procedure.

Fig. 4.

Fig. 5.


Fig. 6.

Fig. 7.

10. Conditions of use, storage, and transportation

When stored properly, the shelf life of the medical device is 3 years. Do not use the product after expiry date.

Store and use at a temperature from 15 °C to 30 °C and humidity not above 70 %. Keep out of reach of children. Protect from heating.

The medical device is transported by all types of transport in covered vehicles in accordance with the regulations for freight transportation applicable for each means of transport. Ambient temperature during transportation should be from -20 °C to 40 °C, humidity of not more than 80 %. When transporting, avoid damaging the cardboard boxes.

11. Disposal procedure

The medical device should be disposed of in accordance with the applicable legislation at the time of disposal in the country where the medical device is being used. Class A of medical waste hazard. –

12. Manufacturer / Manufacturing site / Developer

Aurena Laboratories AB, Sweden
Fjärrviksvägen 22, SE-653 50 Karlstad, Sweden
www.aurena.se, www.bagonvalve.com
Tel.: +46 54 44 44 800; e-mail: contact@aurena.se

13. Authorized representative of the manufacturer in the Russian Federation

Nizhpharm JSC, Russia
7 Salganskaya St., Nizhny Novgorod 603105
Tel.: +7 (831) 278-80-88; Fax: +7 (831) 430-72-28
E-mail: med@stada.ru, website: www.aqualor.ru

14. Warranty

The manufacturer guarantees the sterility of the medical device and its compliance with the requirements of international and national standards during the medical device entire shelf life, subject to the compliance with the rules of transportation, storage, and operating conditions established by the manufacturer and when used in accordance with the instruction for use.

15. Explanation of symbols on the package

	Manufacturer		Sterilization by radiation
	Best before...		Batch code
	Date of manufacture		


Deputy Executive Officer,
Aurena Laboratories AB

Anders Bared

SIGNED BEFORE ME

Christian Nilsson
Notary Public, City of Karlstad

Identification: ☐ driver's license ☐ passport ☐ identity card
Number: n/a ☒ Bank ID
2023-10-20

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