
11.06, 2019

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Production sites

Tokuyama Dental Corporation Kashima Factory
26, Sunayama, Kamisu-shi, Ibaraki 314-0255, Japan

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«PROTECO», JSC

Legal address 192102, Saint-Petersburg, Salova street, 57, building1, letter "I" and 196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"

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Dental composite, light cured filling material ESTELITE Σ QUICK- radiopaque composite material with submicron filler. Contains 82% by weight of silicon-zirconium and composite filler. A high filler load offers decreased polymerization shrinkage. The inorganic filler contained in the material is a spherical submicron filler (average particle size of 0.2 μm ; fractional composition of particles from 0.1 to 0.3 μm), providing excellent gloss stability and abrasion resistance. The monomeric base contains Bisphenol A di (2-hydroxy propoxy) dimethacrylate (the same compound as Bis-GMA (bisphenol A-

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glycidyl methacrylate)) and triethylene glycol dimethacrylate. The product used radical-enhanced polymerization activation system (Radical-Amplified Photopolymerization initiator technology - RAP technology). When using a photopolymerizer with a wavelength in the range of 400-500 nm (peak - 470 nm), the curing time can be reduced to one third compared with traditional composite materials. The radical-enhanced polymerization activation system provides an extended working time of product. It is recommended to be used for direct restoration of anterior and posterior teeth, including the occlusal surface of teeth and carious cavities of all classes (according to Black's classification), restoration of teeth with direct composite veneers, diastema closure, and repairing ceramic and composite restorations.

ESTELITE Σ QUICK 3 SYRINGE KIT includes the filling material in 3 shades: A2, A3, OA3; ESTELITE Σ QUICK 3 SYRINGE KIT II – includes shades OA1, OA2, OA3; ESTELITE Σ QUICK 3 SYRINGE KIT III – shades A2, OA2, OA3. ESTELITE Σ QUICK SYRINGE PROMO KIT includes the filling material in 6 shades: A2, A3, OA1, OA2, OA3, OPA2.

6. Precautionary measures.

- 1) Use the devices only in strict conformance with the below listed recommendations.
- 2) The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
- 3) Do not use the devices for any purposes other than those listed in the Instruction for use.
- 4) Do not use the devices if the safety seals are broken or appear to have been tampered with.
- 5) If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.
- 6) Use examination gloves (plastic, vinyl, or latex) at all times when handling the devices to avoid the possibility of allergic reactions from methacrylic monomers. If the materials penetrate the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 7) Do not mix the materials with other brands of resin composite or other shades of the same material to avoid incomplete cure or entrapment of air bubbles.
- 8) Avoid contact of the devices with eyes, mucosal membrane, skin and clothing. - If the material comes in contact with eyes, thoroughly flush eyes with water and immediately refer to a medical care. If the material comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. If the device comes into contact with the skin or clothing, immediately

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saturate the area with an alcohol soaked cotton swab or gauze. Instruct the patient to rinse his mouth immediately after treatment.

- 9) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 10) The device should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 11) To avoid the unintentional ingestion of the device, do not leave it unsupervised within the reach of patients and children.
- 12) Clean the placement instruments and brushes with alcohol after use.
- 13) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

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Cavity Preparation.

- Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in erasing demarcations between the cavity margins and

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restoration, thereby enhancing both esthetics and retention.

- Scalloped (undulant) bevel could be desirable depending on the esthetics and retention.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

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Bonding system.

- Apply a bonding system according to manufacturer's instructions.
- Do not use self-cured bonding systems. We advise that light-cured bonding systems should be used with ESTELITE Σ QUICK since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

Dispensing

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately.

Do not apply unnecessary force to the syringe immediately after removing from the refrigerator.

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8. Performance Specifications

Table 2. Performance specifications of Dental composite, light cured filling material ESTELITE Σ QUICK

Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 $\pm$ 5 seconds. sample material remains physically uniform

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Polymerization depth ¹⁾	At least 1.0 mm (Opaque shades) At least 1.5 mm (Shades, except Opaque)
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa
Water absorption ¹⁾	Maximum: 40 µg / mm ³
Solubility ¹⁾	Maximum: 7.5 µg /mm ³
Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ¹⁾	After 10 hours of light irradiation only slight changes in color are allowed. ¹⁾
Radiocontrast ¹⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

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ESTELITE Σ QUICK is functioning during operation at nominal temperatures from 32 to 42 ° C and is resistant to the environmental effects of the oral cavity.

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*The Instruction for use inserted in the package of each version of the medical devices contain the shortened information, in accordance with the requirements of the Order 11N of the Ministry of health. of the Russian Federation dd. 19.01.2017 , Section II, point 8.

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-Temperature +18 °C-+30 °C

15. Shelf life.

Shelf life: 3 years from a date of production.

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Dental composite, light cured filling material ESTELITE Σ QUICK in syringes

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depending on the degree of their epidemiological, toxicological and radiation hazard, are classified as epidemiologically safe waste, close in composition to solid household waste (MSW). To dispose the medical devices residues safely, an unused material from the syringe should be removed and light-cured before disposing according to the local regulations.

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

APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by TORAO TSUKUMA
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of TORAO TSUKUMA, Notary
- Certified
5. at Tokyo
6. DEC. 06, 2019
7. by the Ministry of Foreign Affairs
8. 19-Nº 069956
9. Seal/stamp:
10. Signature

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11.06, 2019

A handwritten signature in black ink, appearing to read "M. Sakamoto".

Makoto Sakamoto
President
Tokuyama Dental Corporation

INSTRUCTION FOR USE

For the medical device

Dental composite, light cured filling material ESTELITE in kits

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Revision 1

1. Name of Medical Device.

Dental composite, light cured filling material ESTELITE in kits

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II:

- Dental composite, light cured filling material ESTELITE Σ QUICK in syringes, 2,0 ml (3,8 g), shades: A1, A2, A3, A3.5, OA3, OPA2 - 6 pcs.;
- Dental composite, light cured filling material ESTELITE FLOW QUICK in syringes 1,0 ml (1,8 g), shades: A2, A3, OA3 - 3 pcs.;
- Dental bonding material TOKUYAMA BOND FORCE II (5 ml) - 1 pc.; (registration certificate № RZN 2016/4080 dd. 11.05.2016)
- Dispensing palette for bond Dispensing Well - 1 pc. (registration certificate № RZN 2016/4080 dd. 11.05.2016);
- ESTELITE Σ QUICK Shade Guide - 1 pc.;
- Disposable Applicator (Superfine) - 50 pcs. (registration certificate № RZN 2016/4080 dd. 11.05.2016);
- Plastic tips for syringes ESTELITE FLOW QUICK Syringe Tip - 27 pcs.

2. Details of the developer, manufacturer and authorized representative

Developer

Tokuyama Dental Corporation

(38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan)

Tel: +81-3-3835-7201

URL: <http://www.tokuyama-dental.com>

Manufacturer

Tokuyama Dental Corporation

(38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan)

Tel: +81-3-3835-7201

URL: <http://www.tokuyama-dental.com>

Production sites

Tokuyama Dental Corporation Kashima Factory

26, Sunayama, Kamisu-shi, Ibaraki 314-0255, Japan

Authorized representative of the manufacturer in the Russian Federation

«PROTECO», JSC

Legal address 192102, Saint-Petersburg, Salova street, 57, building1, letter "I" and

196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"

Phone +7 (812) 635-88-90; Email: info@protecodent.ru

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

3. Intended use of the product with an indication of the potential consumer, field of application.

Intended use: the product is intended to fill the cavity of a tooth (restore damaged hard tooth tissues).

Potential consumers of a medical device: dentists.

Field of application: dentistry.

4. Indications, contraindications and possible side effects associated with the intended use of the medical device.

Indication:

For ESTELITE Σ QUICK

- Direct anterior and posterior restorations including the occlusal surface and all carious classes
- Direct bonded composite veneer
- Diastema closure
- Repair of porcelain/composite

For ESTELITE FLOW QUICK,

- Direct anterior and posterior restorations
- As the first layer of direct restorations
- Blocking out cavity undercuts before fabricating indirect restorations
- Repair of porcelain/composite

Contraindications: the product contains methacrylic monomers. Do not use to treat patients who are allergic or hypersensitive to methacrylic and similar to them monomers.

Possible side effects: allergic reactions to methacrylic monomers.

5. Description and operating principal of the medical device.

The kit consists of Dental composite, light cured filling material ESTELITE Σ QUICK in syringes (see subpoint 6.2.) and Dental composite, light cured filling material ESTELITE FLOW QUICK in syringes (see subpoint 6.1.), TOKUYAMA BOND FORCE II (registration certificate № RZN 2016/4080 dd. 11.05.2016), bond dispensing pallets Dispensing Well (registration certificate № RZN 2016/4080 dd. 11.05.2016), ESTELITE Σ QUICK Shade Guide, Disposable Applicator (Superfine) (registration certificate № RZN 2016/4080 dd. 11.05.2016), plastic tips for syringes ESTELITE FLOW QUICK Syringe Tip. Dental composite, light cured filling material ESTELITE FLOW QUICK is a low viscosity, medium flow, light cured, radiopaque composite resin. Putting ESTELITE FLOW QUICK directly into the prepared cavity makes working with the material easier.

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Photopolymerization acceleration technology (Radical - Amplified Photopolymerization initiator technology - RAP technology) is used in the material. The photopolymerization time can be reduced to one third compared to traditional light-cured materials containing camphorquinone (CQ) as a photoinitiator when using a photopolymerizer with a wavelength of 400-500 nm (peak - 470 nm). Contains 71% by weight of specific submicron silicon-zirconium and silicon-titanium filler. The filler particles have a spherical shape, the average particle size is 0.3 micron (from 0.04 to 0.6 micron), which ensures a durable gloss of the material and its high abrasion resistance. The monomeric base contains Bisphenol A polyethoxymethacrylate (Bis-MPEPP), triethylene glycol dimethacrylate (TEGDMA), and 1,6-bis (methacryl ethyloxycarbonylamino) trimethyl hexane (UDMA). It is recommended to use for direct restoration of anterior and posterior teeth (in particular, for small cavities, surface defects, and tunneling technique of restoration), as a first layer when performing direct restorations, to eliminate undercuts of the cavity before performing indirect restorations, to restore defects in ceramic and composite restorations.

Dental composite, light cured filling material ESTELITE Σ QUICK- radiopaque composite material with submicron filler. Contains 82% by weight of silicon-zirconium and composite filler. A high filler load offers decreased polymerization shrinkage. The inorganic filler contained in the material is a spherical submicron filler (average particle size of 0.2 μm ; fractional composition of particles from 0.1 to 0.3 μm), providing excellent gloss stability and abrasion resistance. The monomeric base contains Bisphenol A di (2-hydroxy propoxy) dimethacrylate (the same compound as Bis-GMA (bisphenol A-glycidyl methacrylate)) and triethylene glycol dimethacrylate. The product used radical-enhanced polymerization activation system (Radical-Amplified Photopolymerization initiator technology - RAP technology). When using a photopolymerizer with a wavelength in the range of 400-500 nm (peak - 470 nm), the curing time can be reduced to one third compared with traditional composite materials. The radical-enhanced polymerization activation system provides an extended working time of product. It is recommended to be used for direct restoration of anterior and posterior teeth, including the occlusal surface of teeth and carious cavities of all classes (according to Black's classification), restoration of teeth with direct composite veneers, diastema closure, and repairing ceramic and composite restorations.

ESTELITE Σ QUICK Shade Guide allows you to compare the color of the teeth and choose the appropriate shade of the relevant product. ESTELITE FLOW QUICK Syringe Tip is designed to optimize the dispensing of a flowable composite material and apply it on various surfaces.

6. Precautionary measures.

- 1) Use the devices only in strict conformance with the below listed recommendations.
- 2) The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

- 3) Do not use the devices for any purposes other than those listed in the Instruction for use.
- 4) Do not use the devices if the safety seals are broken or appear to have been tampered with.
- 5) If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.
- 6) Use examination gloves (plastic, vinyl, or latex) at all times when handling the devices to avoid the possibility of allergic reactions from methacrylic monomers. If the materials penetrate the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 7) Do not mix the materials with other brands of resin composite or other shades of the same material to avoid incomplete cure or entrapment of air bubbles.
- 8) Avoid contact of the devices with eyes, mucosal membrane, skin and clothing. - If the material comes in contact with eyes, thoroughly flush eyes with water and immediately refer to a medical care. If the material comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. If the device comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze. Instruct the patient to rinse his mouth immediately after treatment.
- 9) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 10) The device should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 11) To avoid the unintentional ingestion of the device, do not leave it unsupervised within the reach of patients and children.
- 12) Clean the placement instruments and brushes with alcohol after use.
- 13) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

7. Order of application.

Dental composite, light cured filling material ESTELITE Σ QUICK.

Cleaning. Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

Shade selection.

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Select the appropriate shade of ESTELITE Σ QUICK prior to placing the rubber dam using the shade guide provided for ESTELITE Σ QUICK.

- Complete the shade selection within 5 minutes; desiccated teeth are lighter than wet or moist teeth.
- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade.
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

Isolation.

A rubber dam is the preferred method to isolate the tooth from contamination.

Cavity Preparation.

- Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in erasing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.
- Scalloped (undulant) bevel could be desirable depending on the esthetics and retention.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

Pulp protection.

Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. Do not use eugenol based materials to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE Σ QUICK.

Bonding system.

- Apply a bonding system according to manufacturer's instructions.
- Using the disposable applicator, apply TOKUYAMA BOND FORCE II on the surface to be bonded and wait 10 seconds or more. Using an oil-free air/water syringe, apply mild air continuously to the TOKUYAMA BOND FORCE II surface until the runny TOKUYAMA BOND FORCE II stays in the same position without any movement (usually for 5 seconds).

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Light-cure the surface for 10 seconds or more. Confirm that the light-curing unit has sufficient intensity ($>300\text{mW/cm}^2$) before using.

- Do not use self-cured bonding systems. We advise that light-cured bonding systems should be used with ESTELITE Σ QUICK since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

Dispensing

- Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately.
- Do not apply unnecessary force to the syringe immediately after removing from the refrigerator.

Filling and contouring

Fill the cavity incrementally. Increments should not exceed the indicated depth (please refer to the table 1).

Do not mix ESTELITE Σ QUICK with other brands of resin composite or other shades of ESTELITE Σ QUICK to avoid incomplete cure or voids.

Curing

Light-cure each increment for at least the indicated time (please see table 1), keeping the curing light tip within a distance of 2 mm from the increment. If other brands of composite resins are layered over the cured composite, follow their instructions.

Finishing

Shape and polish the restoration. For finishing, use fine finishing diamond points. To erase the demarcation between the composite and the natural tooth rotate a 12 fluted carbide bur slowly along the margins without water. Use finishing metal strips or vinyl strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool.

Dental composite, light cured filling material ESTELITE FLOW QUICK

Cleaning

Teeth should be cleaned with pumice and water to remove surface stains and plaque.

Shade Selection

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Select the appropriate shade of ESTELITE FLOW QUICK.

Isolation

A rubber dam is the preferred method of isolation.

Cavity Preparation

Prepare the cavity. In case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.

Pulp Protection


Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. **DO NOT USE EUGENOL BASED MATERIALS** to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE FLOW QUICK.

Bonding system

- Apply a bonding system according to its manufacturer's instructions.
- Using the disposable applicator, apply TOKUYAMA BOND FORCE II on the surface to be bonded and wait 10 seconds or more. Using an oil-free air/water syringe, apply mild air continuously to the TOKUYAMA BOND FORCE II surface until the runny TOKUYAMA BOND FORCE II stays in the same position without any movement (usually for 5 seconds). Light-cure the surface for 10 seconds or more. Confirm that the light-curing unit has sufficient intensity ($>300\text{mW/cm}^2$) before using.

Attaching the syringe tip

- 1) Hold the nozzle and turn the cap counterclockwise for removal. [1]
- 2) Hold the nozzle and place the Syringe Tip onto it. Turn the Syringe Tip clockwise until it is securely locked. [2]

For Dental composite, light cured filling material ESTELITE FLOW QUICK (1 ml)

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Filling and contouring

Cavities can be directly filled from the syringe, or indirectly filled with an instrument after extruding the paste on the mixing pad. Increments should not exceed the indicated depth at a time (please refer to the table 2).

- Do not apply unnecessary force to the syringes, immediately after removal from the refrigerator.
- Do not mix ESTELITE FLOW QUICK with other brands of composite resin to avoid incomplete cure or air bubble contamination. In order to avoid air bubble contamination, do not mix ESTELITE FLOW QUICK with another shade of the paste.
- After extruding the paste, remove and discard the Syringe Tip, wipe the nozzle with gauze and attach the cap immediately.

Curing

Light-cure each increment for at least the indicated time in the table 2.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

Finishing

Shape and polish the restoration. For finishing, use fine finishing diamond points. To erase the demarcation between the composite and the natural tooth rotate a 12 fluted carbide bur slowly along the margins without water. Use finishing metal strips or vinyl strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool.

Indication for filling and curing.

Be sure to light-cure the material extraorally and check the time needed for complete hardening of it with your light curing unit before performing the bonding procedure. The following table illustrates the correlation between curing time and increment depth.

Dental composite, light cured filling material ESTELITE Σ QUICK

Table 1. Relationship between curing time and increment depth.

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Lamp Type	Lamp power (mW /cm ²)	Cure time (seconds)	Layer thickness (mm) ¹⁾	
			A1, A2, A3, A3.5, B1, B2, B3, C1, C2, OA1, OA2, CE(Inc.), BW, WE	A4, A5, B4, C3, OA3, OPA2
Halogen	400	10	1.9	1.6
		20	2.2	1.9
	600	10	2.0	1.7
		20	2.5	2.0
	800	5	1.8	1.6
		10	2.2	1.9
LED	260	10	1.7	1.4
		20	2.0	1.7
	900	5	1.7	1.5
		10	2.0	1.8
Plasma arc	950	3	1.9	1.5

Dental composite, light cured filling material ESTELITE FLOW QUICK

Table 2. Relationship between curing time and increment depth:

Lamp Type	Lamp power (mW /cm ²)	Cure time (seconds)	Layer thickness (mm) ¹⁾		
			A1, A2, B1, B2, C1, OA1, OA2, BW, CE (Inc.)	A3, A 3.5, B3, C2, OA3	A4 (Cerv.) B4, C3
Halogen	400	10	2.1	2.0	1.8
		20	2.5	2.3	1.9
	600	10	2.3	2.2	1.9
		20	2.5	2.4	2.1
	800	5	2.0	-	1.4
		10	2.3	2.2	1.9
20		2.6	2.5	2.2	
LED	120	30	2.1	-	1.6
		40	2.2	2.1	1.6
		60	2.4	2.2	1.9
	260	10	2.0	-	1.4
		20	2.3	2.1	1.8
		30	2.4	2.2	2.0
	900	5	2.1	-	1.4
		10	2.3	2.2	1.8
Plasma arc	950	3	2.2	2.0	1.8

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Performance Specifications

**Table 3. Performance specifications of Dental composite, light cured filling material
ESTELITE Σ QUICK, ESTELITE FLOW QUICK**

Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 ± 5 seconds. sample material remains physically uniform
Polymerization depth ¹⁾	At least 1.0 mm (Opaque shades) At least 1.5 mm (Shades, except Opaque)
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa
Water absorption ¹⁾	Maximum: 40 $\mu\text{g} / \text{mm}^3$
Solubility ¹⁾	Maximum: 7.5 $\mu\text{g} / \text{mm}^3$
Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ¹⁾	After 10 hours of light irradiation only slight changes in color are allowed. ¹⁾
Radiocontrast ¹⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

1) These tests are performed with reference to the ISO 4049 standard.

2) This test is performed according to our internal standard.

Table 4. Performance specifications of Dental bonding material

Test item	Specification
Appearance	Homologous mixed solution
Light curing depth	minimum 0.1 mm
Bond strength ¹⁾	minimum 12 MPa to enamel, minimum 10 MPa to dentine
Coating property	Form the uniform thickness film

TOKUYAMA BOND FORCE II.

1) This test is performed by referring to ISO/TS 11405: 2003.

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

ESTELITE Σ QUICK and ESTELITE FLOW QUICK are functioning during operation at nominal temperatures from 32 to 42 ° C and is resistant to the environmental effects of the oral cavity.

9. The standards of compliance of the medical device

ISO 13485:2016, ISO 14971:2012, ISO 15223-1:2016, ISO 7405:2008, EN ISO 10993-1-2009, EN ISO10993-3-2014, EN ISO 10993-5-2009, ISO 10993-10:2010, EN ISO 10993-12-2012, ISO 4049:2009, EN ISO 10993-11-2009, EN ISO 10993-18-2009, EN 1041-2008, EN 1641-2009

10. Scope of supply

Dental composite, light cured filling material ESTELITE in kits


REF NO.	Components
13201	ESTELITE Σ QUICK SYRINGE SYSTEM KIT II: - Dental composite, light cured filling material ESTELITE Σ QUICK in syringes, 2,0 ml (3,8 g), shades: A1, A2, A3, A3.5, OA3, OPA2 - 6 pcs.; - Dental composite, light cured filling material ESTELITE FLOW QUICK in syringes 1,0 ml (1,8 g), shades: A2, A3, OA3 - 3 pcs.; - Dental bonding material TOKUYAMA BOND FORCE II (5 ml) - 1 pcs.; (registration certificate № RZN 2016/4080 dd. 11.05.2016) - Dispensing palette for bond Dispensing Well - 1 pc. (registration certificate № RZN 2016/4080 dd. 11.05.2016); - ESTELITE Σ QUICK Shade Guide - 1 pc.; - Disposable Applicator (Superfine) - 50 pcs. (registration certificate № RZN 2016/4080 dd. 11.05.2016); - Plastic tips for syringes ESTELITE FLOW QUICK Syringe Tip - 27 pcs. Instructions for Use*

*The Instruction for use inserted in the package of each version of the medical devices contain the shortened information, in accordance with the requirements of the Order 11N of the Ministry of health of the Russian Federation dd. 19.01.2017 , Section II, point 8.

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

11. Interpretation of symbols used in labeling of consumer and transport packaging:

Symbol	Description	Symbol	Description
	Manufacturer		LOT number
	Catalogue number		"Use before"
	Warning symbol to show the possibility for the occurrence of skin allergic reaction		«Caution! Refer to the instructions for use»
	Highly flammable liquid and vapor.		Authorized Representative in the EU
	Temperature range		conformity to GOST R
	Compliance with European standards		

12. Requirements for maintenance and repair of the medical device.

Medical devices do not required maintenance and repair.

13. Sterilization, cleaning and disinfection.

The medical device is supplied non-sterile and is not intended to sterilize the product before using it.

Cleaning methods of ESTELITE Σ QUICK SHADE GUIDE is presented in the table below:

	Clean method
Correct	- washing in soapy water -ultrasonic cleaning when immersed in water -immersion in a 3% solution of glutarald - wiping with a cotton swab dipped in alcohol solution (80%)
Wrong	- Steam sterilizer ¹⁾ -Sterilization by dry heat ¹⁾

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

	- Immersion in alcohol solution (80%) ²⁾ - Immersion in sodium hypochlorite solution ³⁾ Immersion in a solution of hydrogen peroxide ³⁾
--	--

- 1) Shade samples and base plates melt
- 2) Base plate will turn white and crack
- 3) Shade samples may discolor

Dispensing Well cleaning: after use the dispensing pallette should be cleaned by alcohol-containing disinfecting wipes, intended for disinfection of devices made of polymer materials.

14. Conditions of transportation, storage and operation.

Transportation of products packaged in shipping containers must be carried out by all types of closed vehicles in compliance with the precautions specified on the shipping container and in accordance with the rules for the carriage of goods in force on this type of transport.

Transportation terms and conditions:

- temperature - from 0 to + 25°C
- relative humidity of air - 10- 100%
- Atmosphere pressure: 800 ~ 1,060 hPa

Storage conditions:

- temperature - from 0 to + 25°C
- relative humidity of air - 10-70 %
- Atmosphere pressure: 900 ~ 1,060 hPa

Operating Conditions

- Operation temperature from +32 ° C to +42 ° C

Room conditions during operation of the medical devices:

-Temperature +18 °C-+30 °C

15. Shelf life.

Shelf life: 2 years from a date of production

16. The procedure for the disposal and destruction of medical devices

Medical devices: Dental composite, light cured filling material ESTELITE in kits

Instruction for use for the medical device:

ESTELITE Σ QUICK SYRINGE SYSTEM KIT II

Dental composite, light cured filling material ESTELITE Σ QUICK and ESTELITE FLOW QUICK in syringes,
Dental bonding material TOKUYAMA BOND FORCE II (registration certificate № RZN 2016/4080 dd. 11.05.2016)
Plastic tips for syringes ESTELITE FLOW QUICK Syringe Tip,
ESTELITE Σ QUICK Shade Guide,
Dispensing palette for bond Dispensing Well (registration certificate № RZN 2016/4080 dd. 11.05.2016)

depending on the degree of their epidemiological, toxicological and radiation hazard, are classified as epidemiologically safe waste, close in composition to solid household waste (MSW). To dispose the medical devices residues safely, an unused material from the syringe should be removed and light-cured before disposing according to the local regulations, Dental bonding material TOKUYAMA BOND FORCE II should be absorbed by a neutral absorbing material such as gauze or cotton wool.

Disposable Applicator (Superfine) depending on the degree of their epidemiological, toxicological and radiation hazard are classified as epidemiologically hazardous medical wastes: materials and instruments, things, contaminated by blood and/or other biological liquids. After use applicators should be placed in specific pouches for medical wastes, different by color, depending on the hazard class of wastes (in this case yellow color pouches are used). After a pouch is full for $\frac{3}{4}$ of its volume, pouches should be tied by screws and transported to a place of temporary storage.

The products must be collected and transferred to specialized licensed organizations involved in the disposal of medical waste. Transportation of garbage from institutions to the place of disposal is carried out by specialized companies.

Medical devices must be disposed in accordance with medical procedures in medical institutions and local, state or federal regulations. The choice of methods for safe disposal and destruction of products depends on the capacity and profile of the medical organization, the availability of facilities for decontamination / neutralization of waste, method of disposal / destruction of waste, adopted in the administrative territory (incineration, landfill, disposal).

Instruction for use for the medical device:

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
APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by TORAO TSUKUMA
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of TORAO TSUKUMA, Notary

Certified

5. at Tokyo
6. DEC. 06. 2019
7. by the Ministry of Foreign Affairs
8. 19-№ 069949
9. Seal/stamp:
10. Signature

T. TANAKA

Toshie TANAKA

For the Minister for Foreign Affairs

Registered No. 866 of 2019

Certificate of Acknowledgment of Notary

On this 6th day of DEC , 2019 , before me , Torao Tsukuma , a notary in and for Tokyo Legal Affairs Bureau , personally appeared Tomoko Mifune, with satisfactory evidence of her identification and of her being the agent of Makoto Sakamoto, President of Tokuyama Dental Corporation, said Makoto Sakamoto , had acknowledged to have affixed his signature and seal to the attached document.

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Indication:

For PALFIQUE ESTELITE PASTE

- Direct anterior and posterior restorations including the occlusal surface and all carious classes;

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Dental composite, light cured filling material PALFIQUE ESTELITE PASTE is a light-cured, radiopaque, submicron filled composite resin. Contains 82% by weight of silicon-zirconium with spherical particles of 0.2 µm and composite filler. The monomeric base contains Bisphenol A di (2-hydroxy propoxy) dimethacrylate (the same compound as Bis-GMA (bisphenol A-glycidyl methacrylate)) and triethylene glycol dimethacrylate. It is recommended to be used for direct restoration of the front and posterior teeth, including the occlusal surface of the teeth and carious cavities of all classes (according to Black's classification).

The kit includes the filling material in 5 shades: A1, A2, A3, OA2 (2 syringes), OA3.

6. Precautionary measures.

- 1) Use the devices only in strict conformance with the below listed recommendations.
- 2) The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
- 3) Do not use the devices for any purposes other than those listed in the Instruction for use.
- 4) Do not use the devices if the safety seals are broken or appear to have been tampered with.
- 5) If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.

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PALFIQUE ESTELITE PASTE SYRINGE INTRO KIT II

- 6) Use examination gloves (plastic, vinyl, or latex) at all times when handling the devices to avoid the possibility of allergic reactions from methacrylic monomers. If the materials penetrate the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 7) Do not mix the materials with other brands of resin composite or other shades of the same material to avoid incomplete cure or entrapment of air bubbles.
- 8) Avoid contact of the devices with eyes, mucosal membrane, skin and clothing. - If the material comes in contact with eyes, thoroughly flush eyes with water and immediately refer to a medical care. If the material comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. If the device comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze. Instruct the patient to rinse his mouth immediately after treatment.
- 9) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 10) The device should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 11) To avoid the unintentional ingestion of the device, do not leave it unsupervised within the reach of patients and children.
- 12) Clean the placement instruments and brushes with alcohol after use.
- 13) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

7. Order of application.

1. Prepare the cavity, thoroughly rinse with water and dry. If the pulp protection is required, use a lining material. Do not use eugenol based materials, as these materials will inhibit curing of bonding agent and/or composite resin.
2. Apply light-cured bonding system according to manufacturer's instructions. Using the disposable applicator, apply TOKUYAMA BOND FORCE II on the surface to be bonded and wait 10 seconds or more. Using an oil-free air/water syringe, apply mild air continuously to the TOKUYAMA BOND FORCE II surface until the runny TOKUYAMA BOND FORCE II stays in the same position without any movement (usually for 5 seconds). Light-cure the surface for 10 seconds or more. Confirm that the light-curing unit has sufficient intensity ($>300\text{mW/cm}^2$) before using.
3. Place an adequate amount of the material on the mixing pad and fill the cavity. Light irradiation should be for a minimum of 30 seconds for each layer ($\leq 2\text{ mm}$).

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Dental composite, light cured filling material PALFIQUE ESTELITE PASTE Light irradiation should be for a minimum of 30 seconds for each layer (≤ 2 mm).

8. Performance Specifications

Table 1. Performance specifications of Dental composite, light cured filling material PALFIQUE ESTELITE PASTE

Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 ± 5 seconds. sample material remains physically uniform
Polymerization depth ¹⁾	At least 1.5 mm
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa
Water absorption ¹⁾	Maximum: $40 \mu\text{g} / \text{mm}^3$
Solubility ¹⁾	Maximum: $7.5 \mu\text{g} / \text{mm}^3$
Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ¹⁾	After 10 hours of light irradiation only slight changes in color are allowed. ¹⁾
Radiocontrast ¹⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

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PALFIQUE ESTELITE PASTE SYRINGE INTRO KIT II

**Table 2. Performance specifications of Dental bonding material
TOKUYAMA BOND FORCE II.**

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PALFIQUE ESTELITE PASTE is functioning during operation at nominal temperatures from 32 to 42 ° C and is resistant to the environmental effects of the oral cavity.

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
*The Instruction for use inserted in the package of each version of the medical devices contain the shortened information, in accordance with the requirements of the Order 11N of the Ministry of health of the Russian Federation dd. 19.01.2017, Section II, point 8.

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Dispensing Well cleaning: after use the dispensing pallet should be cleaned by alcohol-containing disinfecting wipes, intended for disinfection of devices made of polymer materials.

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- temperature - from 0 to + 23°C
- relative humidity of air - 10- 100%
- Atmosphere pressure: 800 ~ 1,060 hPa

Storage conditions:

- temperature - from 0 to + 23°C
- relative humidity of air - 10-70 %
- Atmosphere pressure: 900 ~ 1,060 hPa

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
APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by TORAO TSUKUMA
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of TORAO TSUKUMA, Notary
5. at Tokyo
6. DEC. 06. 2019
7. by the Ministry of Foreign Affairs
8. 19- № 069957
9. Seal/stamp:
10. Signature

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The kit consists of Dental composite, light cured filling material ESTELITE FLOW QUICK in 3 shades: A2, A3, A3.5 and ESTELITE FLOW QUICK Metal Syringe Tips.

Dental composite, light cured filling material ESTELITE FLOW QUICK is a low viscosity, medium flow, light cured, radiopaque composite resin. Putting ESTELITE FLOW QUICK directly into the prepared cavity makes working with the material easier. Photopolymerization acceleration technology (Radical - Amplified Photopolymerization initiator technology - RAP technology) is used in the material. The photopolymerization time can be reduced to one third compared to traditional light-cured materials containing camphorquinone (CQ) as a photoinitiator when using a photopolymerizer with a wavelength of 400-500 nm (peak - 470 nm). Contains 71% by weight of specific submicron silicon-zirconium and silicon-titanium filler. The filler particles have a spherical shape, the average particle size is 0.3 micron (from 0.04 to 0.6 micron), which ensures a durable gloss of the material and its high abrasion resistance. The monomeric base contains Bisphenol A polyethoxymethacrylate (Bis-MPEPP), triethylene glycol dimethacrylate (TEGDMA), and 1,6-bis (methacryl ethyloxycarbonylamino) trimethyl hexane (UDMA). It is recommended to use for direct restoration of anterior and posterior teeth (in particular, for small cavities, surface defects, and tunneling technique of restoration), as a first layer when performing direct restorations, to eliminate undercuts of the cavity before performing indirect restorations, to restore defects in ceramic and composite restorations.

ESTELITE FLOW QUICK Metal Syringe Tips are designed to optimize the dispensing of a flowable composite material and apply it on various surfaces.

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Pulp Protection


Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. **DO NOT USE EUGENOL BASED MATERIALS** to protect the pulp as these materials will inhibit curing of bonding agent and/or **ESTELITE FLOW QUICK**.

Bonding system

Apply a bonding system according to its manufacturer's instructions.

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Filling and contouring

Cavities can be directly filled from the syringe, or indirectly filled with an instrument after extruding the paste on the mixing pad. Increments should not exceed the indicated depth at a time (please refer to the table 1).

- Do not apply unnecessary force to the syringes, immediately after removal from the refrigerator.
- Do not mix **ESTELITE FLOW QUICK** with other brands of composite resin to avoid incomplete cure or air bubble contamination. In order to avoid air bubble contamination, do not mix **ESTELITE FLOW QUICK** with another shade of the paste.
- After extruding the paste, remove and discard the Syringe Tip, wipe the nozzle with gauze and attach the cap immediately.

Curing

Light-cure each increment for at least the indicated time in the table 1.

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Table 1. Relationship between curing time and increment depth:

Lamp Type	Lamp power (mW /cm ²)	Cure time (seconds)	Layer thickness (mm) ¹⁾		
			A1, A2, B1, B2, C1, OA1, OA2, BW, CE (Inc.)	A3, A 3.5, B3, C2, OA3	A4 (Cerv.) B4, C3
Halogen	400	10	2.1	2.0	1.8
		20	2.5	2.3	1.9
	600	10	2.3	2.2	1.9
		20	2.5	2.4	2.1
	800	5	2.0	-	1.4
		10	2.3	2.2	1.9
LED	120	20	2.6	2.5	2.2
		30	2.1	-	1.6
		40	2.2	2.1	1.6
	260	60	2.4	2.2	1.9
		10	2.0	-	1.4
		20	2.3	2.1	1.8
		30	2.4	2.2	2.0
	900	5	2.1	-	1.4
		10	2.3	2.2	1.8
Plasma arc	950	3	2.2	2.0	1.8

8. Performance Specifications

Table 2. Performance specifications of Dental composite, light cured filling material ESTELITE FLOW QUICK

Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 ± 5 seconds. sample material remains physically uniform
Polymerization depth ¹⁾	At least 1.5 mm
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa
Water absorption ¹⁾	Maximum: 40 µg / mm ³
Solubility ¹⁾	Maximum: 7.5 µg /mm ³

Instruction for use for the medical device

ESTELITE FLOW QUICK 3 Syringe Kit

Page 6 of 9

Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ¹⁾	After 10 hours of light irradiation only slight changes in color are allowed. ¹⁾
Radiocontrast ¹⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

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ESTELITE FLOW QUICK is functioning during operation at nominal temperatures from 32 to 42 ° C and is resistant to the environmental effects of the oral cavity.

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10. Scope of supply

Dental composite, light cured filling material **ESTELITE** in kits

REF NO.	Components
12100	ESTELITE FLOW QUICK 3 Syringe Kit : - Dental composite, light cured filling material ESTELITE FLOW QUICK in syringes 2,0 ml (3,6 g), shades: A2, A3, A3.5 - 3 pcs.; - ESTELITE FLOW QUICK Metal Syringe Tip - 45 pcs. Instructions for Use*.

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Instruction for use for the medical device

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Instruction for use for the medical device

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This public document

2. has been signed by TORAO TSUKUMA

3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau

4. bears the seal/stamp of TORAO TSUKUMA, Notary

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