
15.06, 2019
«Approved by»

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«PROTECO», JSC

Legal address 192102, Saint-Petersburg, Salova street, 57, building1, letter "I" and 196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"
Phone +7 (812) 635-88-90; Email: info@protecodent.ru

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Indications:

For PALFIQUE ESTELITE PASTE, PALFIQUE ESTELITE PASTE MINI

- Direct anterior and posterior restorations including the occlusal surface and all carious classes;

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Dental composite, light cured filling material PALFIQUE ESTELITE PASTE, PALFIQUE ESTELITE PASTE MINI - is a light-cured, radiopaque, submicron filled composite resin. Contains 82% by weight of silicon-zirconium with spherical particles of 0.2 µm and composite filler. The monomeric base contains Bisphenol A di (2-hydroxy propoxy) dimethacrylate (the same compound as Bis-GMA (bisphenol A-glycidyl methacrylate)) and triethylene glycol dimethacrylate. It is recommended to be used for direct restoration of the front and posterior teeth, including the occlusal surface of the teeth and carious cavities of all classes (according to Black's classification).

6. Precautionary measures.

- 1) Use the devices only in strict conformance with the below listed recommendations.
- 2) The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
- 3) Do not use the devices for any purposes other than those listed in the Instruction for use.
- 4) Do not use the devices if the safety seals are broken or appear to have been tampered with.
- 5) If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.
- 6) Use examination gloves (plastic, vinyl, or latex) at all times when handling the devices to avoid the possibility of allergic reactions from methacrylic monomers. If the materials penetrate the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 7) Do not mix the materials with other brands of resin composite or other shades of the same material to avoid incomplete cure or entrapment of air bubbles.
- 8) Avoid contact of the devices with eyes, mucosal membrane, skin and clothing. - If the material comes in contact with eyes, thoroughly flush eyes with water and immediately refer to a medical care. If the material comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. If the device comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze. Instruct the patient to rinse his mouth immediately after treatment.

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Page 3 of 7

- 9) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 10) The device should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 11) To avoid the unintentional ingestion of the device, do not leave it unsupervised within the reach of patients and children.
- 12) Clean the placement instruments and brushes with alcohol after use.
- 13) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

7. Order of application.

- 1) Prepare the cavity, thoroughly rinse with water and dry. If the pulp protection is required, use a lining material. Do not use eugenol based materials, as these materials will inhibit curing of bonding agent and/or composite resin.
- 2) Apply light-cured bonding system according to manufacturer's instructions.
- 3) Place an adequate amount of the material on the mixing pad and fill the cavity. Light irradiation should be for a minimum of 30 seconds for each layer (≤ 2 mm).
- 4) Shape and polish the restoration.

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Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 ± 5 seconds. sample material remains physically uniform
Polymerization depth ¹⁾	At least 1.0 mm (Opaque shades) At least 1.5 mm (Shades, except Opaque)
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa

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Page 4 of 7

Water absorption ¹⁾	Maximum: 40 µg / mm ³
Solubility ¹⁾	Maximum: 7.5 µg /mm ³
Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ^{1),2)}	After 10 hours of light irradiation only slight changes in color are allowed. ¹
Radiocontrast ²⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

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
Page 5 of 7

	Compliance with European standards		mark of conformity GOST R
---	---------------------------------------	---	------------------------------

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- temperature - from 0 to + 23°C;
- relative humidity of air - 10- 100%
- Atmosphere pressure: 800 ~ 1,060 hPa

Storage conditions:

- temperature - from 0 to + 23°C;
- relative humidity of air - 10-70 %
- Atmosphere pressure: 900 ~ 1,060 hPa

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Medical devices must be disposed in accordance with medical procedures in medical institutions and local, state or federal regulations. The choice of methods for safe disposal and destruction of products depends on the capacity and profile of the medical organization, the

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Page 6 of 7

availability of facilities for decontamination / neutralization of waste, method of disposal / destruction of waste, adopted in the administrative territory (incineration, landfill, disposal).

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

This public document

2. has been signed by TORAO TSUKUMA

3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau

4. bears the seal/stamp of TORAO TSUKUMA, Notary

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
Tokuyama Dental Corporation

Instruction For Use

Dental composite, light cured filling material ESTELITE in individual packages (syringes), with accessories

Versions

Dental composite, light cured filling material ESTELITE Σ QUICK in a syringe, 2.0 ml (3.8 g), shades: A1; A2; A3; A 3,5; A4; A5; B1; B2; B3; B4; C1; C2; C3; OA1; OA2; OA3; OPA2; CE (Inc); BW; WE.

Dental composite, light cured filling material ESTELITE Σ QUICK MINI in a syringe, 1.0 g, shades: A2; A3.

Revision 1

1. Name of medical device.

Dental composite, light cured filling material ESTELITE in individual packages (syringes), with accessories (hereinafter referred to as ESTELITE material, material, device)

Versions:

1. Dental composite, light cured filling material ESTELITE Σ QUICK in a syringe, 2.0 ml (3.8 g), shades: A1; A2; A3; A 3,5; A4; A5; B1; B2; B3; B4; C1; C2; C3; OA1; OA2; OA3; OPA2; CE (Inc); BW; WE.
2. Dental composite, light cured filling material ESTELITE Σ QUICK MINI in a syringe, 1.0 g, shades: A2; A3.

Accessories: ESTELITE Σ QUICK Shade Guide

2. Details of the developer, manufacturer and authorized representative

Developer

Tokuyama Dental Corporation
(Tokuyama Dental Corporation)
(38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan)
Tel: +81-3-3835-7201
URL: <http://www.tokuyama-dental.com>

Manufacturer

Tokuyama Dental Corporation
(38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan)
Tel: +81-3-3835-7201
URL: <http://www.tokuyama-dental.com>

Production sites

Tokuyama Dental Corporation Kashima Factory
26, Sunayama, Kamisu-shi, Ibaraki 314-0255, Japan

Authorized representative of the manufacturer in the Russian Federation

«PROTECO», JSC
Legal address 192102, Saint-Petersburg, Salova street, 57, building1, letter "I" and 196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"
Phone +7 (812) 635-88-90; Email: info@protecodent.ru

3. Intended use of the product with an indication of the potential consumer, field of application

Intended use: the product is intended to fill the cavity of a tooth (restore damaged hard tooth tissues).

Potential consumers of a medical device: dentists.

Field of application: dentistry.

4. Indications, contraindications and possible side effects associated with the intended use of the medical device.

Indications:

For ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI

- Direct anterior and posterior restorations including the occlusal surface and all carious classes
- Direct bonded composite veneer
- Diastema closure
- Repair of porcelain/composite

Contraindications: the product contains methacrylic monomers. Do not use to treat patients who are allergic or hypersensitive to methacrylic and similar to them monomers.

Possible side effects: allergic reactions to methacrylic monomers.

5. Description and operating principal of the medical device.

Dental composite, light cured filling material ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI - radiopaque composite material with submicron filler. Contains 82% by weight of silicon-zirconium and composite filler. A high filler load offers decreased polymerization shrinkage. The inorganic filler contained in the material is a spherical submicron filler (average particle size of 0.2 μm ; fractional composition of particles from 0.1 to 0.3 μm), providing excellent gloss stability and abrasion resistance. The monomeric base contains Bisphenol A di (2-hydroxy propoxy) dimethacrylate (the same compound as Bis-GMA (bisphenol A-glycidyl methacrylate)) and triethylene glycol dimethacrylate. The product used radical-enhanced polymerization activation system (Radical-Amplified Photopolymerization initiator technology - RAP technology). When using a photopolymerizer with a wavelength in the range of 400-500 nm (peak - 470 nm), the curing time can be reduced to one third compared with traditional composite materials. The radical-enhanced polymerization activation system provides an extended working time of product. It is recommended to be used for direct restoration of anterior and posterior teeth, including the occlusal surface of teeth and carious cavities of all classes (according to Black's classification), restoration of teeth with direct composite veneers, diastema closure, and repairing ceramic and composite restorations.

ESTELITE Σ QUICK Shade Guide allows you to compare the color of the teeth and choose the appropriate shade of the relevant product.

6. Precautionary measures.

- 1) Use the devices only in strict conformance with the below listed recommendations.
- 2) The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
- 3) Do not use the devices for any purposes other than those listed in the Instruction for use.
- 4) Do not use the devices if the safety seals are broken or appear to have been tampered with.

Instruction For Use for the medical device:

Dental composite, light cured filling material ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI in syringe

Page 3 of 9

- 5) If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.
- 6) Use examination gloves (plastic, vinyl, or latex) at all times when handling the devices to avoid the possibility of allergic reactions from methacrylic monomers. If the materials penetrate the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 7) Do not mix the materials with other brands of resin composite or other shades of the same material to avoid incomplete cure or entrapment of air bubbles.
- 8) Avoid contact of the devices with eyes, mucosal membrane, skin and clothing. - If the material comes in contact with eyes, thoroughly flush eyes with water and immediately refer to a medical care. If the material comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. If the device comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze. Instruct the patient to rinse his mouth immediately after treatment.
- 9) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 10) The device should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 11) To avoid the unintentional ingestion of the device, do not leave it unsupervised within the reach of patients and children.
- 12) Clean the placement instruments and brushes with alcohol after use.
- 13) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

7. Order of application.

Cleaning. Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

Shade selection.

Select the appropriate shade of ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI prior to placing the rubber dam using the shade guide provided for ESTELITE Σ QUICK.

- Complete the shade selection within 5 minutes; desiccated teeth are lighter than wet or moist teeth.
- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade.
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

Isolation.

Instruction For Use for the medical device:

Dental composite, light cured filling material ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI in syringe

Page 4 of 9

A rubber dam is the preferred method to isolate the tooth from contamination.

Cavity Preparation.

- Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in erasing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.
- Scalloped (undulant) bevel could be desirable depending on the esthetics and retention.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

Pulp protection.

Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. Do not use eugenol based materials to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI.

Bonding system.

Apply a bonding system according to manufacturer's instructions.

Do not use self-cured bonding systems. We advise that light-cured bonding systems should be used with ESTELITE Σ QUICK since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

Dispensing

Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately.

Do not apply unnecessary force to the syringe immediately after removing from the refrigerator.

Filling and contouring

Fill the cavity incrementally. Increments should not exceed the indicated depth (please refer to the table 1).

Do not mix ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI with other brands of resin composite or other shades of ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI to avoid incomplete cure or voids.

Curing

Light-cure each increment for at least the indicated time (please see table 1), keeping the curing light tip within a distance of 2 mm from the increment. If other brands of composite resins are layered over the cured composite, follow their instructions.

Finishing

Instruction For Use for the medical device:

Dental composite, light cured filling material ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI in syringe

Page 5 of 9

Shape and polish the restoration. For finishing, use fine finishing diamond points. To erase the demarcation between the composite and the natural tooth rotate a 12 fluted carbide bur slowly along the margins without water. Use finishing metal strips or vinyl strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool.

Indication for filling and curing.

Be sure to light-cure the material extraorally and check the time needed for complete hardening of it with your light curing unit before performing the bonding procedure. The following table illustrates the correlation between curing time and increment depth.

Dental composite, light cured filling material ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI.

Table 1. Relationship between curing time and increment depth:

Lamp Type	Lamp power (mW /cm ²)	Cure time (seconds)	Layer thickness (mm) ¹⁾	
			A1, A2, A3, A3.5, B1, B2, B3, C1, C2, OA1, OA2, CE(Inc.), BW, WE	A4, A5, B4, C3, OA3, OPA2
Halogen	400	10	1.9	1.6
		20	2.2	1.9
	600	10	2.0	1.7
		20	2.5	2.0
	800	5	1.8	1.6
		10	2.2	1.9
LED	260	10	1.7	1.4
		20	2.0	1.7
	900	5	1.7	1.5
		10	2.0	1.8
Plasma arc	950	3	1.9	1.5

8. Performance Specifications

Table 2. Performance Specifications: Dental composite, light cured filling material ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI.

Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 ± 5 seconds. sample material remains physically uniform
Polymerization depth ¹⁾	At least 1.0 mm (Opaque shades) At least 1.5 mm (Shades, except Opaque)
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa
Water absorption ¹⁾	Maximum: 40 µg / mm ³
Solubility ¹⁾	Maximum: 7.5 µg /mm ³
Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions

Instruction For Use for the medical device:

Dental composite, light cured filling material ESTELITE Σ QUICK, ESTELITE Σ QUICK MINI in syringe

Page 6 of 9

Color stability ^{1),2)}	After 10 hours of light irradiation only slight changes in color are allowed. ¹
Radiocontrast ³⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

1) These tests are performed with reference to the ISO 4049 standard.

2) This test is performed according to our internal standard.

Dental composite, light cured filling material ESTELITE is functioning during operation at nominal temperatures from 32 to 42 ° C and is resistant to the environmental effects of the oral cavity.

9. The standards of compliance of the medical device

ISO 13485:2016, ISO 14971:2012, ISO 15223-1:2016, ISO 7405:2008, EN ISO 10993-1-2009, EN ISO 10993-3-2014, EN ISO 10993-5-2009, ISO 10993-10:2010, EN ISO 10993-12-2012, ISO 4049:2009, EN ISO 10993-11-2009, EN ISO 10993-18-2009, EN 1041-2008, EN 1641-2009

10. Interpretation of symbols used in labeling of consumer and transport packaging:

Symbol	Description	Symbol	Description
	Manufacturer		Lot number
	Catalogue number		"Use by date..."
	Warning symbol to show the possibility for the occurrence of skin allergic reaction		«Caution! Refer to the instructions for use»
	Temperature range		Authorized Representative to the EU
	Compliance with European standards		mark of conformity GOST R

11. Sterilization, cleaning and disinfection.

The medical device is supplied non-sterile and is not intended to sterilize the product before using it.

Cleaning methods of ESTELITE Σ QUICK SHADE GUIDE is presented in the table below:

	Clean method
Correct	- washing in soapy water -ultrasonic cleaning when immersed in water -immersion in a 3% solution of glutarald - wiping with a cotton swab dipped in alcohol solution (80%)
Wrong	- Steam sterilizer ¹⁾ -Sterilization by dry heat ¹⁾ - Immersion in alcohol solution (80%) ²⁾ - Immersion in sodium hypochlorite solution ³⁾ Immersion in a solution of hydrogen peroxide ³⁾

- 1) Shade samples and base plates melt
- 2) Base plate will turn white and crack
- 3) Shade samples may discolor

12. Conditions of transportation, storage and operation.

Transportation of products packaged in shipping containers must be carried out by all types of closed vehicles in compliance with the precautions specified on the shipping container and in accordance with the rules for the carriage of goods in force on this type of transport.

Transportation terms and conditions:

- temperature - from 0 to + 25°C
- relative humidity of air - 10- 100%
- Atmosphere pressure: 800 ~ 1,060 hPa

Storage conditions:

- temperature - from 0 to + 25°C
- relative humidity of air - 10-70 %
- Atmosphere pressure: 900 ~ 1,060 hPa

Operating Conditions

- Operation temperature from +32 ° C to +42 ° C

Room conditions during operation of the medical devices:

- Temperature +18 °C-+30 °C

13. Shelf life

Shelf life: 3 years from a date of production.


14. The procedure for the disposal and destruction of medical devices

Medical devices Dental composite, light cured filling material ESTELITE in individual packages (syringes), with accessories”, depending on the degree of their epidemiological, toxicological and radiation hazard, are classified as epidemiologically safe waste, close in composition to solid

household waste (MSW). To dispose the medical devices residues safely, an unused material from the syringe should be removed and light-cured before disposing.

The products must be collected and transferred to specialized licensed organizations involved in the disposal of medical waste. Transportation of garbage from institutions to the place of disposal is carried out by specialized companies.

Medical devices must be disposed in accordance with medical procedures in medical institutions and local, state or federal regulations. The choice of methods for safe disposal and destruction of products depends on the capacity and profile of the medical organization, the availability of facilities for decontamination / neutralization of waste, method of disposal / destruction of waste, adopted in the administrative territory (incineration, landfill, disposal).

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2. Details of the developer, manufacturer and authorized representative.

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«PROTECO», JSC

Legal address 192102, Saint-Petersburg, Salova street, 57, building 1, letter "I" and 196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"
Phone +7 (812) 635-88-90; Email: info@protecodent.ru

3. Intended use of the product with an indication of the potential consumer, field of application.

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Indications:

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- Direct posterior restorations including the occlusal surface

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Dental composite, light cured filling material ESTELITE POSTERIOR - radiopaque composite material. Contains 84% by weight of silicon-zirconium filler (average particle size: 2 microns, fractional composition of particles: from 0.1 to 10 microns). A high filler load offers decreased polymerization shrinkage. The monomeric base contains Bisphenol A di (2-hydroxy propoxy) dimethacrylate (the same compound as Bis-GMA (bisphenol A-glycidyl methacrylate)) triethylene glycol dimethacrylate (TEGDMA) and bisphenol A polyethoxy methacrylate (Bis-MPEPP). The product used radical-enhanced polymerization activation system (Radical-Amplified Photopolymerization initiator technology - RAP technology). When using a photopolymerizer with a wavelength in the range of 400-500 nm (peak - 470 nm), the curing time can be reduced to one third compared with traditional composite materials. The radical-enhanced polymerization activation system provides an extended working time of product. It is recommended to be used for direct restoration of posterior teeth, including the occlusal surface of the teeth, repairing ceramic and composite restorations.

6. Precautionary measures.

1. Use the devices only in strict conformance with the below listed recommendations.
2. The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
3. Do not use the devices for any purposes other than those listed in the Instruction for use.
4. Do not use the devices if the safety seals are broken or appear to have been tampered with.
5. If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.
6. Use examination gloves (plastic, vinyl, or latex) at all times when handling the devices to avoid the possibility of allergic reactions from methacrylic monomers. If the materials penetrate the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
7. Do not mix the materials with other brands of resin composite or other shades of the same material to avoid incomplete cure or entrapment of air bubbles.
8. Avoid contact of the devices with eyes, mucosal membrane, skin and clothing. - If the 9. material comes in contact with eyes, thoroughly flush eyes with water and immediately refer to a medical care. If the material comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. If the device comes into contact with the skin or clothing, immediately

saturate the area with an alcohol soaked cotton swab or gauze. Instruct the patient to rinse his mouth immediately after treatment.

9. When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
10. The device should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
11. To avoid the unintentional ingestion of the device, do not leave it unsupervised within the reach of patients and children.
12. Clean the placement instruments and brushes with alcohol after use.
13. This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

7. Order of application.

Cleaning

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

Shade Selection

Select the appropriate shade of ESTELITE POSTERIOR prior to placing the rubber dam.

- Complete the shade selection within 5 minutes; desiccated teeth are lighter than wet or moist teeth.
- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade –
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

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Cavity Preparation

Prepare the cavity and rinse with water. Add chamfers to the margins of posterior preparations (class I, II) because chamfers assist in diminishing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.

- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

Pulp Protection

Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. Do not use the eugenol based materials to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE POSTERIOR.

Bonding system.

- Apply a bonding system according to manufacturer's instructions.
- Do not use self-cured bonding systems. We advise that light-cured bonding systems should be used with ESTELITE POSTERIOR since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

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Filling and contouring

Fill the cavity incrementally using a suitable placing instrument. Increments should not exceed the indicated depth (please refer to the table 2).

Do not mix ESTELITE POSTERIOR with other brands of resin composite or other shades of ESTELITE POSTERIOR to avoid incomplete cure or entrapment of air bubbles.

Curing

Light-cure each increment for at least the indicated time (please see table 1), keeping the curing light tip within a distance of 2 mm from the increment.

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Finishing

Shape and polish the restoration. For finishing, use fine finishing diamond points. To diminish the demarcation between the composite and the natural tooth, rotate a 12 fluted carbide bur slowly along the margins without water. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool.

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Dental composite, light cured filling material ESTELITE POSTERIOR

Table 1. Relationship between curing time and increment depth.

Lamp Type	Lamp power (mW /cm ²)	Cure time (seconds)	Layer thickness (mm) ¹⁾
Halogen	400	10	1.9
		20	2.3
	600	10	2.0
		20	2.4
	800	5	1.8
		10	2.4
LED	260	10	1.9
		20	2.0
	900	5	1.8
		10	2.0
Plasma arc	950	3	2.1

8. Performance Specifications

Table 2. Performance Specifications : Dental composite, light cured filling material ESTELITE POSTERIOR in syringe

Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 ± 5 seconds. sample material remains physically uniform
Polymerization depth ¹⁾	At least 1.5 mm
Ultimate strength of flexure ¹⁾	Minimum: 120 MPa
Water absorption ¹⁾	Maximum: 40 µg / mm ³
Solubility ¹⁾	Maximum: 7.5 µg /mm ³
Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ^{1),2)}	After 10 hours of light irradiation only slight changes in color are allowed. ¹⁾
Radiocontrast ¹⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

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11. Sterilization, cleaning and disinfection

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-Temperature +18 °C-+30 °C

13. Shelf life

Shelf life: 3 years from a date of production

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Medical devices Dental composite, light cured filling material ESTELITE in individual packages (syringes), with accessories", depending on the degree of their epidemiological, toxicological and radiation hazard, are classified as epidemiologically safe waste, close in composition to solid household waste (MSW). To dispose the medical devices residues safely, an unused material from the syringe should be removed and light-cured before disposing.

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
APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by TORAO TSUKUMA
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of TORAO TSUKUMA, Notary
5. at Tokyo
6. SEP. 17. 2019
7. by the Ministry of Foreign Affairs
8. 19- N^o 117264
9. Seal/stamp:
10. Signature

T. TANAKA

Toshie TANAKA

For the Minister for Foreign Affairs

Registered No.669 of 2019

Certificate of Acknowledgment of Notary

On this 17th day of SEP , 2019 , before me , Torao Tsukuma , a notary in and for Tokyo Legal Affairs Bureau , personally appeared Tomoko Mifune, with satisfactory evidence of her identification and of her being the agent of Makoto Sakamoto, President of Tokuyama Dental Corporation, said Makoto Sakamoto , had acknowledged to have affixed his signature and seal to the attached document.

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15.06, 2019
«Approved by»

Makoto Sakamoto
President
Tokuyama Dental Corporation

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Version

Dental composite, light cured filling material ESTELITE FLOW QUICK

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Version:

Dental composite, light cured filling material ESTELITE FLOW QUICK consisting of:

- Dental composite, light cured filling material ESTELITE FLOW QUICK in a syringe, 2.0 ml (3.6 g): shades A1; A2; A3; A3.5; A4 (Cerv.); B1; B2; B3; B4; C1; C2; C3; OA1; OA2; OA3; CE (Inc.); BW.
- ESTELITE FLOW QUICK Metal Syringe Tip - 15 pcs.

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«PROTECO», JSC

Legal address 192102, Saint-Petersburg, Salova street, 57, building1, letter "I" and 196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"

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Dental composite, light cured filling material ESTELITE FLOW QUICK in syringe

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Dental composite, light cured filling material ESTELITE FLOW QUICK is a low viscosity, medium flow, light cured, radiopaque composite resin. Putting ESTELITE FLOW QUICK directly into the prepared cavity makes working with the material easier. Photopolymerization acceleration technology (Radical - Amplified Photopolimerization initiator technology - RAP technology) is used in the material. The photopolymerization time can be reduced to one third compared to traditional light-cured materials containing camphorquinone (CQ) as a photoinitiator when using a photopolymerizer with a wavelength of 400-500 nm (peak - 470 nm). Contains 71% by weight of specific submicron silicon-zirconium and silicon-titanium filler. The filler particles have a spherical shape, the average particle size is 0.3 micron (from 0.04 to 0.6 micron), which ensures a durable gloss of the material and its high abrasion resistance. The monomeric base contains Bisphenol A polyethoxymethacrylate (Bis-MPEPP), triethylene glycol dimethacrylate (TEGDMA), and 1, 6-bis (methacryl ethyloxycarbonylamino) trimethyl hexane (UDMA). It is recommended to use for direct restoration of anterior and posterior teeth (in particular, for small cavities, surface defects, and tunneling technique of restoration), as a first layer when performing direct restorations, to eliminate undercuts of the cavity before performing indirect restorations, to restore defects in ceramic and composite restorations.

ESTELITE FLOW QUICK Metal Syringe Tip is designed to optimize the dispensing of a flowable composite material and apply it on various surfaces.

6. Precautionary measures.

- 1) Use the devices only in strict conformance with the below listed recommendations.

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Dental composite, light cured filling material ESTELITE FLOW QUICK in syringe

- 2) The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
- 3) Do not use the devices for any purposes other than those listed in the Instruction for use.
- 4) Do not use the devices if the safety seals are broken or appear to have been tampered with.
- 5) If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.
- 6) Use examination gloves (plastic, vinyl, or latex) at all times when handling the devices to avoid the possibility of allergic reactions from methacrylic monomers. If the materials penetrate the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 7) Do not mix the materials with other brands of resin composite or other shades of the same material to avoid incomplete cure or entrapment of air bubbles.
- 8) Avoid contact of the devices with eyes, mucosal membrane, skin and clothing. - If the material comes in contact with eyes, thoroughly flush eyes with water and immediately refer to a medical care. If the material comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. If the device comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze. Instruct the patient to rinse his mouth immediately after treatment.
- 9) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 10) The device should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 11) To avoid the unintentional ingestion of the device, do not leave it unsupervised within the reach of patients and children.
- 12) Clean the placement instruments and brushes with alcohol after use.
- 13) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

7. Order of application.

Cleaning

Teeth should be cleaned with pumice and water to remove surface stains and plaque.

Shade Selection

Select the appropriate shade of ESTELITE FLOW QUICK.

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Instruction For Use for the medical device:
Dental composite, light cured filling material ESTELITE FLOW QUICK in syringe

Pulp Protection

Bonding system

Attaching the syringe tip (Metal Syringe Tip)

-

Do not apply unnecessary force to the syringes, immediately after removal from the refrigerator.

After extruding the paste, remove and discard the Syringe Tip, wipe the nozzle with gauze and attach the cap immediately.

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Finishing

Shape and polish the restoration. For finishing, use fine finishing diamond points. To erase the

Dental composite, light cured filling material ESTELITE FLOW QUICK in syringe

demarcation between the composite and the natural tooth rotate a 12 fluted carbide bur slowly along the margins without water. Use finishing metal strips or vinyl strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool.

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Dental composite, light cured filling material ESTELITE FLOW QUICK

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8. Performance Specifications.

Table 2. Performance Specifications: Dental composite, light cured filling material ESTELITE FLOW QUICK

Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 ± 5 seconds. sample material remains physically uniform
Polymerization depth ¹⁾	At least 1.0 mm (Opaque shades) At least 1.5 mm (Shades, except Opaque)
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa
Water absorption ¹⁾	Maximum: 40 µg / mm ³
Solubility ¹⁾	Maximum: 7.5 µg /mm ³

Instruction For Use for the medical device:
Dental composite, light cured filling material ESTELITE FLOW QUICK in syringe

Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ^{1),2)}	After 10 hours of light irradiation only slight changes in color are allowed. ¹⁾
Radiocontrast ²⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

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
9. The standards of compliance of the medical device.

ISO 13485:2016, ISO 14971:2012, ISO 15223-1:2016, ISO 7405:2008, EN ISO 10993-1-2009, EN ISO10993-3-2014, EN ISO 10993-5-2009, ISO 10993-10:2010, EN ISO 10993-12-2012, ISO 4049:2009, EN ISO 10993-11-2009, EN ISO 10993-18-2009, EN 1041-2008, EN 1641-2009

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Symbol	Description	Symbol	Description
	Manufacturer		Lot number
	Catalogue number		"Use by date..."
	Warning symbol to show the possibility for the occurrence of skin allergic reaction		«Caution! Refer to the instructions for use»
	Temperature range		Authorized Representative to the EU

Instruction For Use for the medical device:
Dental composite, light cured filling material ESTELITE FLOW QUICK in syringe

Symbol	Description	Symbol	Description
	Compliance with European standards		mark of conformity GOST R

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
Tokuyama Dental Corporation

Instruction For Use

**Dental composite, light cured filling material ESTELITE in individual packages
(syringes), with accessories**

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«PROTECO», JSC
Legal address 192102, Saint-Petersburg, Salova street, 57, building 1, letter "I" and 196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"
Phone +7 (812) 635-88-90; Email: info@protecodent.ru

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Dental composite, light cured filling material ESTELITE BULK FILL Flow in syringe

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Dental composite, light cured filling material ESTELITE BULK FILL Flow - low-viscosity radiopaque composite material. This low polymerization stress flowable material can be placed 4 mm increments. The material incorporates Radical-amplified polymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time. Contains 70% by weight of spherical silicon-zirconium filler and composite filler (SR). The spherical shape of the composite filler provides low shrinkage stress during polymerization. The spherical filler contained in the material (average particle size 0.2 μm ; fractional composition of particles from 0.1 to 0.3 μm) provides excellent gloss stability and wear resistance. The monomeric matrix contains Bisphenol A di (2-hydroxy-propoxy) dimethacrylate (the same compound as Bis-GMA (bisphenol A-glycidyl methacrylate)), triethylene glycol dimethacrylate (TEGDMA) and bisphenol A polyethoxy methacrylate (Bis-MPEPP), mequinol, dibutyl hydroxy toluene and UV absorber. It is recommended to be used for direct restoration of the front and posterior teeth, to remove undercuts of the cavity before performing indirect restorations, filling cavities, and for repairing defects of ceramic and composite restorations.

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6. Precautionary measures.

- 1) Use the devices only in strict conformance with the below listed recommendations.
- 2) The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
- 3) Do not use the devices for any purposes other than those listed in the Instruction for use.
- 4) Do not use the devices if the safety seals are broken or appear to have been tampered with.
- 5) If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.

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Dispensing

- 1) Hold the nozzle and turn the cap counterclockwise for removal. [1]
- 2) Hold the nozzle and place the Dispensing Tip onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

Filling and contouring

Cavity can be directly filled from the Dispensing Tip. Increments should not exceed 4mm at a time. Do not apply excessive force to extrude the material from the syringe immediately after removal from the refrigerator. Do not mix ESTELITE BULK FILL Flow with other brands of resin composite or other shades of ESTELITE BULK FILL Flow to avoid incomplete cure or entrapment of air bubbles.

- After extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.

Curing

Light-cure each increment for at least the indicated time (please see the table 1). Light-cure each layer (up to 4mm increment) using Halogen light-curing or LED light-curing unit with a wavelength of 400-500 nm. Be sure the light-curing unit with an intensity of 800mW/cm² or more. Hold the light guide tip as close to the restorative as possible during light exposure.

Finishing

Shape and polish the restoration. For finishing, use finishing discs and/ or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

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Dental composite, light cured filling material ESTELITE BULK FILL Flow in syringe

Light-cure each layer (up to 4mm increment) using Halogen light-curing or LED light-curing unit with a wavelength of 400-500 nm. Hold the light guide tip as close to the restorative as possible during light exposure.

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Test parameter	Specification
Sensitivity to ambient light ¹⁾	After light exposure for 60 ± 5 seconds sample material remains physically uniform.
Polymerization depth ¹⁾	Universal shade: At least 3.5 mm Shades, except universal: At least 3.0 mm
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa
Water absorption ¹⁾	Maximum: 40 µg /mm3
Solubility ¹⁾	Maximum: 7.5 µg /mm3
Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide.
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ¹⁾	After 24 hours of exposure to Xenon light, only slight changes in color are allowed.
Radiocontrast ¹⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness.

- 1) These tests are performed with reference to the ISO4049 standard.
- 2) This test is performed according to our internal standard.

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Dental composite, light cured filling material ESTELITE BULK FILL Flow in syringe

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

APOSTILLE

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- Country: JAPAN
This public document
 - has been signed by TORAO TSUKUMA
 - acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
 - bears the seal/stamp of TORAO TSUKUMA, Notary
- Certified
- at Tokyo
 - SEP. 17. 2019
 - by the Ministry of Foreign Affairs
 - 19- NO 117266
 - Seal/stamp:
 - Signature

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On this 17th day of SEP, 2019, before me, Torao Tsukuma, a notary in and for Tokyo Legal Affairs Bureau, personally appeared Tomoko Mifune, with satisfactory evidence of her identification and of her being the agent of Makoto Sakamoto, President of Tokuyama Dental Corporation, said Makoto Sakamoto, had acknowledged to have affixed his signature and seal to the attached document.

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Versions:

1. Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW High L consisting of:
 - Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW High L in a syringe, 1.8 ml (3.0 g), shades A1; A2; A3; A3.5; OPA2; OPA3; OPA4.
 - Tokuyama Dispensing Tip - 10 pcs.
2. Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW Medium L consisting of:
 - Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW Medium L in a syringe, 1.8 ml (3.0 g), shades A1; A2; A3; A3.5; A4; A5; B3; OPA2; OPA3; OPA4; CE; BW.
 - Tokuyama Dispensing Tip - 10 pcs.
3. Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW SuperLow L consisting of:
 - Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW SuperLow L in a syringe, 1.8 ml (3.0 g), shades A1; A2; A3; A3.5; A4; A5.
 - Tokuyama Dispensing Tip - 10 pcs.

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«PROTECO», JSC

Instruction For Use for the medical devices:

Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L in syringe

Page 2 of 9

Legal address 192102, Saint-Petersburg, Salova street, 57, building 1, letter "I" and
196128, Saint-Petersburg, Varshavskaya street, 5/2, office 401, business center "Varshavsky"
Phone +7 (812) 635-88-90; Email: info@protecodent.ru

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- Direct anterior and posterior restorations
- As the first layer of direct restorations
- Blocking out cavity undercuts before fabricating indirect restorations
- Repair of porcelain/composite

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Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L – low-viscosity radiopaque composite material. Putting of product directly into the prepared cavity simplifies work with the product. Includes radical-enhanced polymerization activator technology (RAP technology). RAP technology provides shorter curing times and longer working hours. It has three degrees of fluidity: High, Medium and SuperLow. Each type contains a spherical silica-zirconium filler and a composite filler (SR). Due to its spherical shape, the composite filler reduces shrinkage during polymerization. The spherical filler contained in the product (average particle size: 200 nm, particle size range: from 100 to 300 nm) provides excellent gloss retention and wear resistance. The SR filler is the powder which is thermal polymerized mixture of its ingredients. The monomer matrix comprises Bisphenol A di (2-hydroxy-propoxy) dimethacrylate (the same compound as Bis-GMA (bisphenol A-glycidyl methacrylate)), bisphenol A polyethoxy methacrylate (Bis-MPEPP), triethylene glycol dimethacrylate (TEGDMA), 1, 6 bis- (methacryl-ethyloxycarbonylamino)trimethyl hexane (UDMA), Mequinol, dibutylhydroxytoluene and UV absorber. It is recommended to be used for direct restoration of the anterior and posterior teeth, filling cavities, for eliminating cracks and preparing the stump before performing indirect restorations, for repairing defects of ceramic and composite restorations.

Instruction For Use for the medical devices:

Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L in syringe

Page 3 of 9

ESTELITE UNIVERSAL FLOW Shade Guide allows you to compare the color of the teeth and choose the appropriate shade of the relevant product.

TOKUYAMA Dispensing Tip is designed to optimize the dispensing of a flowable composite material and apply it on various surfaces.

6. Precautionary measures.

- 1) Use the devices only in strict conformance with the below listed recommendations.
- 2) The devices are designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
- 3) Do not use the devices for any purposes other than those listed in the Instruction for use.
- 4) Do not use the devices if the safety seals are broken or appear to have been tampered with.
- 5) If the devices cause an allergic reaction or hypersensitivity, discontinue its use immediately.
- 6) Use examination gloves (plastic, vinyl, or latex) at all times when handling the devices to avoid the possibility of allergic reactions from methacrylic monomers. If the materials penetrate the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- 7) Do not mix the materials with other brands of resin composite or other shades of the same material to avoid incomplete cure or entrapment of air bubbles.
- 8) Avoid contact of the devices with eyes, mucosal membrane, skin and clothing. - If the material comes in contact with eyes, thoroughly flush eyes with water and immediately refer to a medical care. If the material comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. If the device comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze. Instruct the patient to rinse his mouth immediately after treatment.
- 9) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- 10) The device should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- 11) To avoid the unintentional ingestion of the device, do not leave it unsupervised within the reach of patients and children.
- 12) Clean the placement instruments and brushes with alcohol after use.
- 13) This product is designed to be used at room temperature (18 - 30°C). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

7. Order of application.

Cleaning.

Thoroughly clean the tooth surface with a rubber cup and free-fluoride paste and rinse with water.

Shade Selection

Prior to isolation of tooth, select the appropriate shade of ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L using shade guide provided for ESTELITE UNIVERSAL FLOW.

Isolation

A rubber dam is the preferred method of isolation.

Cavity Preparation

Prepare the cavity. In case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.

Pulp Protection

Calcium hydroxide should be applied if the cavity is in close proximity to the pulp. Do not use eugenol-based materials to protect the pulp as these materials will inhibit curing ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L. ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L are not indicated for direct pulp capping.

Bonding System

Apply a bonding system according to its manufacturer's instructions.

Dispensing

- Hold the nozzle and turn the cap counterclockwise for removal. [1]
- Hold the nozzle and place the Dispensing Tip (metal tip) onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

Filling and contouring

Cavities can be directly filled from the Dispensing Tip, or indirectly filled with an instrument after extruding the paste on the mixing pad. Increments should not exceed the indicated curing depth at a time (please refer to the table 1).

- Do not apply excessive force to extrude the material from the syringe immediately after removal from the refrigerator.

- Do not mix ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L with other brands of resin composite or other shades of ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L to avoid incomplete cure or entrapment of air bubbles.

- In case of syringe, after extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.

Curing

Light-cure each increment for at least the indicated time in the table 1.

If other brands of composite resins are layered over the cured composite, follow their instructions.

Finishing

Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

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Table 1. Relationship between curing time and increment depth

Lamp Type	Lamp power (mW /cm ²)	Cure time (seconds)	Layer thickness (mm) ¹⁾		
			CE, BW, A1, A2, A3, A3.5, A4, B3	A5	OPA2, OPA3, OPA4
Halogen	600	10	1.8	1.6	1.3
		20	2.3	2.0	1.6
	800	10	2.0	1.7	1.5
		20	2.4	2.1	1.7
	1200	5	1.8	1.6	1.3
		10	2.2	2.0	1.6
LED	600	10	1.7	1.5	1.2
		20	2.1	1.8	1.4
	800	10	1.9	1.7	1.3
		20	2.2	2.0	1.6
	1200	10	2.0	1.8	1.5
	2300	3	1.7	1.5	1.2
Plasma arc	1200	3	1.7	1.5	1.1

¹⁾ The thickness of the material layer was set based on the results of tests performed in accordance with clause 7.10 "Light curing depth" of ISO 4049.

8. Performance Specifications

Table 2. Performance Specifications

Dental composite, light cured filling material ESTELITE UNIVERSAL FLOW High L, ESTELITE UNIVERSAL FLOW Medium L, ESTELITE UNIVERSAL FLOW SuperLow L

Test parameter	Specification
Sensitivity to ambient light ¹⁾	After exposure to light for 60 ± 5 seconds. sample material remains physically uniform
Polymerization depth ¹⁾	At least 1.0 mm (Opaque shades) At least 1.5 mm (Shades, except Opaque)
Ultimate strength of flexure ¹⁾	Minimum: 80 MPa
Water absorption ¹⁾	Maximum: 40 $\mu\text{g} / \text{mm}^3$
Solubility ¹⁾	Maximum: 7.5 $\mu\text{g} / \text{mm}^3$
Shade ¹⁾	Shade of a given product, immersed in water for 7 days, should correspond as closely as possible to the shade prescribed by the manufacturer's guide
Appearance ²⁾	The given product must be homologous, without bubbles and foreign inclusions
Color stability ^{1),2)}	After 24 hours of light irradiation, only slight changes in color are allowed. ²⁾
Radiocontrast ¹⁾	Radiocontrast of a given product must meet or exceed this indicator for aluminum of the same thickness

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
2) This test is performed according to our internal standard

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Cleaning methods of ESTELITE UNIVERSAL FLOW SHADE GUIDE is presented in the table below:

	Clean method
Correct	- washing in soapy water -ultrasonic cleaning when immersed in water -immersion in a 3% solution of glutarald - wiping with a cotton swab dipped in alcohol solution (80%)
Wrong	- Steam sterilizer¹⁾ -Sterilization by dry heat¹⁾ - Immersion in alcohol solution (80%)²⁾ - Immersion in sodium hypochlorite solution³⁾ - Immersion in a solution of hydrogen peroxide³⁾

1) Shade samples and base plates melt

2) Base plate will turn white and crack

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Medical devices Dental composite, light cured filling material ESTELITE in individual packages (syringes), with accessories”, depending on the degree of their epidemiological, toxicological and radiation hazard, are classified as epidemiologically safe waste, close in composition to solid household waste (MSW). To dispose the medical devices residues safely, an unused material from the syringe should be removed and light-cured before disposing.

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