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Registered No. 2024-2047


NOTARIAL CERTIFICATE

Seowon Law & Notary Office Inc.


84, Wonheong-ro, Seowon-gu, Cheongju, Korea

Date: 06.05.2024.

«Approved»

  (Signature)
You Jae Hoon
CEO

DiaDent Group International

DiaDent Group International

16, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu,
Cheongju-si, Chungcheongbuk-do, 28161, Korea
Tel : 82-43-266-2315 Fax : 82-43-262-6053

Version 2
INSTRUCTION FOR USE

FOR THE MEDICAL DEVICE

Root canal permanent filling dental material Healtex® in kits

We, the company “DiaDent Group International”, Republic of Korea, legal address: 16, Osongsaengmyeong 4-ro Osong-eup Heungdeok-gu, Cheongju-si Chungcheongbuk-do, 28161, Republic of Korea, being the manufacturer of the medical device “Root canal permanent filling dental material Healtex® in kits” refers to To the Russian Federal Service on Surveillance in Healthcare (Roszdravnadzor) with the request to conduct the state registration of the abovementioned medical device and issue the Registration Certificate in order to release safe, efficient and high quality medical devices in turnover on the Russian market.

The attached document is provided to DiaDent Group International.

This document certifies that the contents of attached true and is confirmed by DiaDent Group International.

Contents

1. **Medical device name**
2. **Information about the manufacturer, developer and authorized representative of the manufacturer**
3. **Intended use of the medical device**
4. **Potential consumers of the medical device**
5. **Indications, contraindications, and potential adverse reactions**
6. **Device description**
7. **Precautions**
8. **Application procedure**
9. **Functional properties.**
10. **Main materials of the medical device.**
11. **List of international standards that the medical device complies with**
12. **Description of symbols used in labeling of consumer and transport packaging.**
13. **Maintenance and repair**
14. **Sterilization, disinfection**
15. **Transport and storage conditions**
16. **Shelf life**
17. **Disposal and/or destruction**

1. Medical device name

Root canal permanent filling dental material Healtex® in kits

1. HealSealer Regular Kit Type A consisting of:
 - dental material Healtex® HealSealer in syringes 16 g - 1pc. (dual-chamber syringe: base 6.5 g, catalyst 9.5 g)
 - mixing pad - 1 pc.
 - spatula - 1 pc.
 - Instruction for use – 1 pc.
2. HealSealer Regular Kit Type B consisting of:
 - dental material Healtex® HealSealer in syringes 4 g - 1pc. (dual-chamber syringe: base 1.625 g, catalyst 2.375 g)
 - mixing pad - 1 pc.
 - spatula - 1 pc.
 - Instruction for use – 1 pc.
3. HealRoot BioSealer Regular Kit consisting of:
 - dental material Healtex® HealRoot BioSealer in syringes 2 g - 1 pc.
 - disposable tips for syringes HP002 - 20 pcs.
 - silicone cap - 1 pc.
 - Instruction for use – 1 pc.
4. HealRoot BioSealer Intro Kit consisting of:
 - dental material Healtex® HealRoot BioSealer in syringes 0.5 g - 1 pc.
 - disposable tips for syringes HP002 - 4 pcs.
 - silicone cap - 1 pc.
 - Instruction for use – 1 pc.
5. Healtex® HealRoot BioMTA Regular Kit consisting of:
 - dental material Healtex HealRoot BioMTA in Eppendorf tube 0,5 g. - 1 pc.
 - mixing pad - 1 pc.
 - spatula - 1 pc.
 - Instruction for use – 1 pc.

Hereinafter – Root canal filling dental material Healtex® in kits, dental material, kits, medical device, device.

2. Information about the manufacturer, developer and authorized representative of the manufacturer

Manufacturer

DiaDent Group International

16, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, Republic of Korea

Tel. +82-43-266-2315

E-mail: diadent@diadent.co.kr

Developer

DiaDent Group International

16, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, Republic of Korea

Tel. +82-43-266-2315

E-mail: diadent@diadent.co.kr

Authorized representative of manufacturer (producer)

PROTECO, Joint Stock Company

PROTECO JSC

57 Salova st., build. 1, lit. I 192102, St. Petersburg,

Tel./fax: +7 (812) 779-30-90;

E-mail: info@protecodent.ru

3. Intended use of the medical device

It is used for permanent obturation of the root canals of teeth.

4. Potential consumers of the medical device

Dentists.

5. Indications, contraindications, and potential adverse reactions

Indications

Root canal filling dental material Healtex® HealSealer

Root canal filling with joint use of gutta-percha points.

Root canal filling dental material Healtex® HealRoot BioSealer

It is used for permanent root canal filling, in combination with or without root canal filling points.
It is used for filling the root canal from the crown side.

Root canal filling dental material Healtex® HealRoot BioMTA

- Treatment of perforations and bifurcation of the root canal caused by carious lesions
- Treatment of root canal perforations due to internal resorption
- Apexogenesis (induction of root development in a viable tooth with an inflamed crown pulp)
- Apexiation (induction of the formation of a mineralized barrier at the end of the root of a young permanent tooth with incomplete root development and necrotic pulp)
- Surgical treatment of root canal perforations due to internal resorption
- Pulp capping
- Removal of the pulp of the affected part of the crown while maintaining the viability and function of the remaining radicular part

Contraindications to the use of a medical device:

Intolerance to the components of raw materials of a medical device by patients.

Potential adverse reactions:

Allergic reaction to the components of the devices.

The material may cause hypersensitivity reactions very rarely. In such cases, discontinue use and consult a doctor.

6. Device description

HealSealer Regular Kit Type A, HealSealer Regular Kit Type B

The kits consist of one dual-chamber syringe (with a protective cap) with dental material Healtext® HealSealer (the weight of material in the syringe depends on the version of the kit), a mixing pallet and a spatula.

Dental material Healtext® HealSealer is a radiopaque hardening composite based on composite resin for permanent filling of the root canal of a tooth. It is a yellow base and a white catalyst placed in the corresponding chambers of the syringe. Base and catalyst are mixed in a 1:1 ratio during application and are used with points for permanent root canal filling of molars. The material has excellent sealing ability with low shrinkage, antimicrobial activity, high biocompatibility, resistance to degradation over time with long-term dimensional stability.

The mixing pallet and spatula are designed to mix the base and catalyst in order to obtain a homogeneous final material.

HealRoot BioSealer Regular Kit, HealRoot BioSealer Intro Kit

The kits consist of one syringe (with a protective cap) with dental material Healtext® HealRoot BioSealer (the weight of the material in the syringe depends on the version of the kit), disposable syringe tips HP002 (the number of tips depends on the version of the kit) and a silicone cap.

Dental material Healtext® HealRoot BioSealer is a white radiopaque hardening material for complete permanent root canal filling based on MTA (mineral trioxide aggregate). It can be used with or without root canal obturation materials. The material seals the root canal by curing when it reacts with the moist environment in the oral cavity. The lack of shrinkage and high adhesion make the material an ideal sealant for permanent obturation and is suitable for all gutta-percha obturation methods.

Disposable syringe tips HP002 are designed to optimize the dosing of the material and its application on various surfaces.

The silicone cap is designed to prevent contact of the material with oxygen and moisture in the air.

Healtext® HealRoot BioMTA

The kit consists of an Eppendorf tube with dental material Healtext® HealRoot BioMTA (the mass of the material in the tube is 0.5 g), a spatula and a mixing pallet.

Dental material Healtext® HealRoot BioMTA is a white hydrophilic biocompatible powder that promotes cementogenesis and has excellent sealing ability. The material is used to seal the root canal.

The mixing pallet and spatula are designed to mix the material with water to form a paste.

7. Precautions

General Precautions:

1. This product shall only be used by a qualified dentist with careful observance of the instructions for use and precautions.
2. Check the product for abnormalities, foreign matter or impurities, expiration date, conformity with the description before using the device. Do not use material in case of discoloration or hardening.
3. Do not use to treat patients allergic to any component of the product and patients with current inflammation.
4. The device is not intended for ingestion.
5. Be sure to install a rubber dam before the procedure to protect the respiratory tract and to avoid swallowing.
6. Do not use the material for purposes other than the intended use.

7. Do not use an expired device.
8. Do not mix with other products.
9. Avoid contact of the product with skin and eyes. In case of contact, rinse immediately with plenty of water and seek medical attention; if dental material **Healtex® HealSealer** comes in contact with the skin, remove immediately the material with cotton and alcohol, then wash thoroughly the affected area with soap and water.
10. Improper storage may result in poor performance and reduced shelf life.

Dental material Healtex® HealSealer

1. Overfilling the canal with filling material can lead to damage to the periapical tissues.
2. Ingestion is prohibited. If accidentally swallowed, give the patient plenty of water and seek medical attention.
3. It is highly recommended to close the cap tightly after each use.
4. Do not use heat when mixing the base material and catalyst.

Dental material Healtex® HealRoot BioSealer

1. Do not use the material for purposes other than the intended use.
2. Do not store or use the contents in a container other than that provided by the manufacturer.
3. Keep the product closed until use to avoid premature reaction with atmospheric moisture.
4. Use with caution as pH increases to 10~12 when reacting with moisture.

Dental material Healtex® HealRoot BioMTA

1. The material may cause hypersensitivity in very rare cases. In this case, discontinue use and consult a doctor.
2. Avoid contact with water other than for the intended use as the product hardens on contact with moisture.
3. Increased humidity leads to a deterioration in the effectiveness of Healtex® HealRoot BioMTA. The product shall be stored in an airtight container prior to immediate use to prevent premature reaction with the humidity in the air.
4. This product shall be used immediately after mixing.
5. Remove excess material from the end of the spatula to clean and reuse the spatula using a disposable cotton pad or a disposable cotton swab moistened with alcohol. Wipe the spatula after removing the material three or more times with a disposable cotton pad/swab moistened with alcohol.

8. Application procedure

Healtex® HealSealer

1. Isolate the area with a rubber dam. Prepare and shape the root canal
2. Remove foreign matter from the root canal. Disinfect and clean thoroughly the root canal with sodium hypochlorite or EDTA and dry completely.
3. Open the protective cap and squeeze out the base (yellow) and catalyst (white) in a 1:1 ratio from the double syringe onto the mixing pallet (they can be distributed at different speeds due to the different viscosities of the base and catalyst).
4. Mix base (yellow) and catalyst (white) for 10-20 seconds using a spatula or metal spatula until completely mixed and acquiring ivory color.
5. Measure the length of the root canal using an electronic apex locator or file. Choose the appropriate gutta-percha point. Apply and distribute evenly a small amount of the mixture on the gutta-percha post and insert slowly the post in the direction of the apex. Dentists can apply alternatively the mixture with a paper point, reamer, or spiral filler (Lentulo spiral). Insert the coated point very slowly in a clockwise direction during obturation and remove it in a counterclockwise direction to prevent the formation and entrapment of air bubbles and overfilling

Page 5 of 11

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE
Root canal permanent filling dental material Healtex® in kits

of the channel. If a Lentulo coil is used, advance and withdraw it slowly towards the apex at extremely low speed.

6. Take an x-ray to check the filling after completing the procedure.

Note.

Clean thoroughly the root canal before applying the material. If hydrogen peroxide was used for irrigation, be sure to check that there are no traces of flushing solution in the channel. This will avoid reaction with the material, which can lead to the formation of bubbles. If sodium hypochlorite and hydrogen peroxide are used alternately to irrigate the canal, hypochlorite shall be used last.

Healtex® HealRoot BioSealer

1. Prepare and flush the root canal with standard endodontic techniques.
2. Put a disposable tip on the syringe and insert into the root canal.
3. Fill the root canal by slowly pressing the plunger. Do not apply too much pressure as this may cause bubbles and/or overflow.
4. When using a gutta-percha point, it must be inserted into the root canal.
5. Use a dental instrument to cut off the top. Excess material can be removed with a cotton pad or cotton swab moistened with water.
6. Remove and dispose of the disposable tip at the end of the procedure and insert a new disposable tip.
7. Close the new disposable tip tightly with the supplied silicone cap.

Healtex® HealRoot BioMTA

Direct pulp capping and perforation treatment

- 1) Isolate the tooth with a rubber dam and form an indentation using a burr at high speed, without touching the tissues close to the pulp.
- 2) Remove dead tissue near the pulp carefully using a new sterilized diamond bur at high speed.
- 3) Rinse the prepared cavity, open pulp, perforation site with 5.25% NaOCl solution.
- 4) If bleeding is observed in the patient, control it with a cotton ball moistened with NaOCl and applied to the hemorrhagic area.
- 5) Mix Healtex® HealRoot BioMTA with distilled water by rubbing in the indicated ratio* for 30 seconds. Add water depending on the ambient temperature and humidity if the mixture hardens or dries out. Fill the prepared cavity after thorough mixing.

***Mixing ratio = 0.5 g : 0.225 g (0.225 cm³)**

(Healtex® HealRoot BioMTA: distilled water)

- 6) If the patient is bleeding, apply a thin layer of Healtex® HealRoot BioMTA and press down with a cotton ball. Place the remaining material in the cavity after that.
- 7) Aim a powerful jet of water at the hardened material with a blaster to check for flushing after placing the material in the cavity. Remove the material completely if at least a part of the material has separated from the cavity, and refill again.
- 8) Using suitable instruments, create a cavity using Healtex® HealRoot BioMTA that will not disturb the linear angle of the tooth cavity. During this procedure, all dentinal tubules located close to the pulp shall be covered with a material of at least 3 mm. It is necessary to perform the restoration of the crown part of the tooth if this condition is difficult to fulfill.
- 9) Take an x-ray picture to make sure that the material has not entered the pulp area. If this happens, prescribe painkillers to reduce pain. If there are visible pockets in the image that are not filled with material, it is necessary to completely remove all material and repeat the whole procedure again;
- 10) Use polymer resin or glass ionomer cement to seal the Healtex® HealRoot BioMTA layer and prevent saliva contamination.

Root apex filling

- 1) Use a surgical burr, after gaining access to the apical lesion with a flap, to excise the root by 3 mm.
- 2) Prepare the cavity to a depth of 3 mm using an ultrasonic handpiece (low power level);
- 3) If the patient is bleeding, control it by applying pressure to the area of bleeding (for example, using a cotton ball with a hemostatic agent such as NaOCl).
- 4) Stain the area to check for cracks or additional canals in the apical area using methylene blue and examine under a microscope.
- 5) Remove excess moisture from the formed cavity.
- 6) Mix Healtex® HealRoot BioMTA with distilled water by rubbing in the indicated ratio* for 30 seconds. Add water depending on the ambient temperature and humidity if the mixture hardens or dries out.

***Mixing ratio = 0.5 g : 0.225 g (0.225 cm³)**

(Healtex® HealRoot BioMTA: distilled water)

- 7) Place the required amount of material into the prepared cavity using a special tool.
- 8) Compact the material gently with a fine plugger. Add small portions of material and compact. Repeat until you have to use an instrument smaller than the size of the cavity;
- 9) Remove excess material and clean the root surface with a damp gauze or cotton ball.
- 10) Rinse the cavity thoroughly after curing with saline and check for washout. The material shall not be washed out.
- 11) Take an x-ray to make sure the material is properly placed. Remove all material if there is not enough material and repeat the filling again.
- 12) Close the surgical opening after confirming the correct filling of the root apex.

Partial pulpectomy of anterior teeth

- 1) Isolate the tooth with a rubber dam and form an indentation using a burr at high speed, without touching the tissues close to the pulp.
- 2) Remove dead tissue near the pulp using a new sterilized diamond bur at high speed.
- 3) Open the pulp chamber, rinse and remove the pulp to the required depth with a new sterilized diamond bur at high speed.
- 4) Flush the prepared cavity with copious amounts of 5.25% NaOCl solution until bleeding stops.
- 5) Mix Healtex® HealRoot BioMTA with distilled water by rubbing in the indicated ratio* for 30 seconds. Add water if the mixture hardens or dries out, depending on the ambient temperature and humidity. Fill the prepared cavity after thorough mixing.

***Mixing ratio = 0.5 g : 0.225 g (0.225 cm³)**

(Healtex® HealRoot BioMTA: distilled water)

- 6) Remove excess water with a sterile cotton ball, and tamp the material into the pulp chamber.
- 7) Fill in the remaining amount of material into the cavity without allowing Healtex® HealRoot BioMTA to harden
- 8) Aim a powerful jet of water at the hardened Healtex® HealRoot BioMTA material with the emptier after inserting the material into the cavity for washout test. Remove the material completely if at least a part of the material has separated from the cavity, and proceed with the root canal treatment, because re-insertion can lead to infection of the pulp.
- 9) Remove the top of the Healtex® HealRoot BioMTA material after curing and place the polymer resin material or glass ionomer cement over it as instructed.

Partial pulpectomy of the posterior teeth

- 1) Isolate the tooth with a rubber dam and form an indentation using a burr at high speed, without touching the tissues close to the pulp.
- 2) Remove dead tissue near the pulp using a new sterilized diamond bur at high speed.
- 3) Open all pulp horns of the tooth and remove the pulp with a new sterilized diamond bur at high

speed. Remove carefully the remaining pulp from the axial wall of the cavity.

4) Flush the prepared cavity with copious amounts of 5.25% NaOCl solution until bleeding stops.

5) Mix Healtex® HealRoot BioMTA with distilled water with rubbing movements in the indicated ratio for 30 seconds. Add water if mixture hardens or dries out, depending on an ambient temperature and humidity. Fill the prepared cavity after thorough mixing.

***Mixing ratio = 0.5 g : 0.225 g (0.225 cm³)**

(Healtex® HealRoot BioMTA: distilled water)

6) Remove excess water with a sterile cotton ball, and tamp the material into the pulp chamber.

7) Fill in the remaining amount of material into the cavity without allowing Healtex® HealRoot BioMTA to harden.

8) Aim a powerful jet of water at the hardened Healtex® HealRoot BioMTA material with the emptier after inserting the material into the cavity for washout test. Remove the material completely, if at least a part of the material has separated from the cavity, and proceed with the root canal treatment, because re-insertion can lead to infection of the pulp;

9) Remove the outer part of the Healtex® HealRoot BioMTA material after curing and place the polymer resin material or glass ionomer cement on it according to the instructions. The minimum thickness of the remaining Healtex® HealRoot BioMTA material shall not be less than 3 mm to prevent secondary infection of the pulp.

9. Functional properties.

Table 1.

Parameter	Value
Healtex® HealSealer	
Weight in syringe	
HealSealer Regular Kit Type A	16 g ± 5 %
HealSealer Regular Kit Type B	4 g ± 5 %
Appearance	Homogeneous, without foreign matter
Working time	1 hour ± 10 %
Flow property	Not less than 17 mm
Material curing time	Within 3-4 hours
Film thickness	Less than 50 µm
Radiopacity	Equivalent to a layer of aluminum with a thickness of more than 3 mm
Solubility	Less than 3.0 %
Healtex® HealRoot BioSealer	
Weight	
HealRoot BioSealer Regular Kit	2 g ± 5 %
HealRoot BioSealer Intro Kit	0.5 g ± 5 %
Appearance	Homogeneous, without foreign matter
Working time	20 min ± 10 %
Flow property	Not less than 17 mm
Material curing time	6 hours 30 minutes ± 10 %
Film thickness	Less than 50 µm
Radiopacity	Equivalent to a layer of aluminum with a thickness of more than 3 mm
Solubility	Less than 3.0 %
Healtex® HealRoot BioMTA	
Weight in tube	0.5 g ± 5 %

Appearance	White powder, without foreign matter
Working time	20 min $\pm$ 10 %
Material curing time	40 $\pm$ 10 minutes
Radiopacity	Equivalent to a layer of aluminum with a thickness of more than 3 mm
Solubility	Less than 3.0 %
Change in size after curing	The shrinkage rate after curing shall be 1.0% and the expansion rate shall be 0.1% or less.

10. Main materials of the medical device.

Dental material Healtex® HealSealer

1. Base
 - Diglycidyl ether of bisphenol A
 - Diglycidyl ether of bisphenol F
 - Zirconium dioxide
2. Catalyst
 - 1,3-benzenedimethanamine
 - Poly(1,4-butanediol)bis(4-aminobenzoate)
 - Polyoxypropylene diamine
 - Zirconium dioxide
 - Calcium tungstate

Dental material Healtex® HealRoot BioSealer

- Calcium cement powder
- Calcium Aluminate Cement Powder
- Ytterbium Fluoride Powder
- Zirconium oxide

Dental material Healtex® HealRoot BioMTA

- Calcium cement powder
- Zirconium oxide

11. List of international standards that the medical device complies with

HealSealer Regular Kit Type A, HealSealer Regular Kit Type B.

EN ISO 13485; EN ISO 14971; ISO 6876; ISO 7405; EN ISO 10993-1; EN ISO 10993-3; EN ISO 10993-5; EN ISO 10993-10; EN ISO 10993-11; IEC 62366-1; EN ISO 15223-1

HealRoot BioSealer Regular Kit, HealRoot BioSealer Intro Kit

EN ISO 13485; EN ISO 14971; EN ISO 10993-1; EN ISO 10993-3; EN ISO 10993-5; EN ISO 10993-10; EN ISO 10993-11; IEC 62366-1; ISO 6876; ISO 7405; EN ISO 15223-1

Healtex® HealRoot BioMTA Regular Kit

EN ISO 13485; EN ISO 14971; EN ISO 10993-1; EN ISO 10993-3; EN ISO 10993-5; EN ISO 10993-10; EN ISO 10993-11; IEC 62366-1; EN ISO 15223-1; ISO 6876; ISO 7405; EN 1041; EN 1641

Reference standards are used in the latest edition: the organization maintains and updates the list of reference standards as a record of its own management system in accordance with ISO 13485.

12. Description of symbols used in labeling of consumer and transport packaging.

	Manufacturer		Temperature limit
	Date of manufacture		Refer to instructions for use or to instruction for use in electronic form
	Use before		Keep away from moisture
	Lot number		Caution!
	Catalog number		Keep away from sunlight
	Top		Limit on the number of tiers in a stack
	Fragile, handle with care		Do not use hooks/Do not pierce

13. Maintenance and repair

The device is disposable. Maintenance and repair are not required.

14. Sterilization, disinfection

Medical devices "Root canal permanent filling dental material Healtex® in kits" are non-sterile, designed for single use. Sterilization and disinfection are not provided.

15. Transport and storage conditions

Storage. Keep away from direct sunlight, excessive heat and high humidity. Store in a dry, cool place tightly closed.

Transportation. Transportation of devices packed in a shipping container shall be carried out by all types of covered vehicles, observing the precautions indicated on the shipping container and in accordance with the rules for the carriage of goods in force for this type of transport.

Storage, transportation and operation conditions indoors:

- Temperature: from 15 to 25 °C for kits with dental material Healtex® HealRoot BioSealer; from 2 to 27 °C for other devices
- Relative humidity: not more than 60% for kits with dental materials Healtex® HealRoot BioMTA, Healtex® HealRoot BioSealer; there are no special requirements for other devices;
- Atmospheric pressure: no special requirements.

16. Shelf life

Shelf life of the devices "HealSealer Regular Kit Type A", "HealSealer Regular Kit Type B" is 2 years from the date of manufacture, subject to storage conditions.

Shelf life of the devices "HealRoot BioSealer Regular Kit", "HealRoot BioSealer Intro Kit" is 2 years from the date of manufacture, subject to storage conditions.

Shelf life of the device "Healtex® HealRoot BioMTA Regular Kit" is 3 years from the date of manufacture, subject to storage conditions.

17. Disposal and/or destruction

Medical devices “Root canal permanent filling dental material Healtex® in kits” after their intended use are classified as epidemiologically hazardous waste (devices contaminated with blood and/or other biological fluids of patients).

In this case, the medical organization only collects and temporarily stores the above waste. Collection is carried out in single-use bags or containers with color coding accepted for hazardous medical waste in a room or in places for temporary storage of medical waste until further transportation to a disposal site.

Transportation of medical waste from medical institutions, their disposal are carried out by specialized organizations licensed to handle epidemiologically hazardous medical waste.

Consumer packaging, defective devices and unused devices with an expired shelf life are epidemiologically safe waste (non-hazardous waste that does not have contact with the biological fluids of patients), close in composition to municipal solid waste (MSW). For safe disposal of curable device residues, unused material shall be removed from the syringe or Eppendorf tube and cured before being disposed of in accordance with local regulations.

Non-hazardous wastes are collected in reusable containers or disposable bags (except for hazardous waste bags). The procedure for collection and storage of non-hazardous waste is determined in accordance with the waste management scheme in each specific organization.

Transportation and disposal are carried out taking into account the requirements adopted for a given territory, in accordance with the sanitary legislation of a particular country for the maintenance of populated areas and the handling of manufacture and consumption waste.

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HA, JUNG MIN, attorney-in-fact of DIADENT GROUP INTERNATIONAL CEO JAE HOON YOU, appeared before me and admitted said principal's subscription to the attached VERSION2 INSTRUCTION FOR USE FOR THE MEDICAL DEVICE.

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This is hereby attested on this 19th day of June, 2024 at this office.

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Name of the office
Seowon Law & Notary Office Inc.
Belong to
Cheongju District Prosecutors' Office
Address of the office
84, Wonheong-ro, Seowon-gu, Cheongju, Korea

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PARK JUNG MIN

(Signature of Notary Public)

This office has been authorized by the Minister of Justice, the Republic of Korea, to act as Notary Public since FEB, 7. 2020 under Law No. 173.

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