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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of ANQING BAOJIE
PACKAGING CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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验证网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

APPROVED BY:

General Director of Anqing Baojie Packaging Co., Ltd.

(No.123 XingYe Road, Development District, Anqing, Anhui Province, China)

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INSTRUCTION FOR USE

on medical device:

Packing material disposable for sterilization DENTAL FORMULA

Initial issue of operational documentation (instructions for use) for medical device:

Version: 1.0 dated 25.06.2024

1 Data on the medical device

1.1 Name of the medical device

Packing material disposable for sterilization DENTAL FORMULA, comprising:

I. Design versions:

1. Roll combined flat disposable DENTAL FORMULA RKP for steam and gas sterilization, dimensions: width from 50 to 500 mm with a step of 5 mm, length from 100 to 200 m with a step of 50 m – from 1 to 12 pcs/pack.
2. Roll combined gusseted disposable DENTAL FORMULA RKS for steam and gas sterilization, dimensions: width from 50 to 500 mm with a step of 5 mm, fold from 20 to 100 mm with a step of 5 mm, length from 100 to 200 m with a step of 50 m – from 1 to 12 pcs/pack.
3. Pouch combined flat self-sealing disposable DENTAL FORMULA PKPS for steam and gas sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm – from 5 to 3000 pcs/pack.
4. Pouch combined flat disposable DENTAL FORMULA PKP for steam and gas sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm – from 5 to 3000 pcs/pack.
5. Kraft paper pouch flat self-sealing disposable DENTAL FORMULA PBPS for heat and steam sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm – from 5 to 3000 pcs/pack.
6. Bleached kraft paper pouch flat self-sealing disposable DENTAL FORMULA PBBPS for heat and steam sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm – from 5 to 3000 pcs/pack.
7. Kraft paper pouch combined flat self-sealing disposable DENTAL FORMULA PBKPS for heat and steam sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm – from 5 to 3000 pcs/pack.
8. Bleached kraft paper pouch combined flat self-sealing disposable DENTAL FORMULA PBBKPS for heat and steam sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm – from 5 to 3000 pcs/pack.

II. Shortened consumer instruction for use – 1 pc.

1.2 Manufacturer (producer) of the medical device

Manufacturer (producer) of the medical device: Anqing Baojie Packaging Co., Ltd., China, No.123 XingYe Road, Development District, Anqing City, Anhui Province, China.

Address of the manufacturing facility of the medical device: Anqing Baojie Packaging Co., Ltd., China, No.123 XingYe Road, Development District, 246003 Anqing City, Anhui Province, People's Republic of China.

1.3 Purpose of the medical device

Packing material disposable for sterilization DENTAL FORMULA (hereinafter referred to as "packing material") is intended for disposable use for finishing sterilization of medical devices, maintaining sterility to the place of their use and observing aseptic properties.

1.4 Scope of the medical device

The field of application of packing material is sterilization of medical devices by the personnel of medical and preventive institutions and medical organizations operating sterilization equipment.

1.5 Indications for use of the medical device

Indications for the use of medical devices are the need to sterilise the products, to maintain sterility to the place of use and to comply with aseptic properties. The use of packaging material in medical organisations is a mandatory requirement.

1.6 Contraindications for use of the medical device

There are no contraindications and side effects of the medical device.

1.7 Intended users of the medical device

The packing material can be used by various specialists involved in sterilization of medical devices in medical organizations. The use of the packing material does not require special qualifications or special education.

1.8 Operating principle of the medical device

The principle of action of a medical device is to create a hermetic physical barrier between the sterilized article and the external environment. The packaging material is permeable to sterilizing agents, allowing them to affect the packaged articles, completely freeing them from pathogenic and non-pathogenic microorganisms, including their spore forms, and impermeable to viruses, bacteria and fungi of the external environment. Such a combination of properties of the packing material allows to sterilize the packaged products and preserve it until the moment of their direct use. The packing material is provided with chemical indicators of sterilization processes, reacting to one or more critical variables of the sterilization process and allowing to distinguish sterilized packaged articles from unsterilized ones. Thus, the medical device provides prevention of nosocomial infections in patients and personnel of medical and preventive institutions.

1.9 Application conditions of the medical device

The packing material is used in therapeutic or therapeutic-preventive medical institutions in a closed room in a dry place, away from direct sunlight, at temperatures from 10°C to 30°C and relative humidity not more than 80%, in accordance with the instructions for the medical device.

1.10 Classification of the medical device

Potential risk of use of the packaging material in accordance with Regulation (EU) 2017/745 - Class I. In accordance with the Order of the Ministry of Health of the Russian Federation No. 4n dated 06.06.2012. "On approval of the nomenclature classification of medical devices" – class 1.

1.11 Type of contact with human body

Short-term contact with intact human skin.

1.12 Shelf life of the medical device

The established shelf life of the medical device is 5 years.

2 Instruction for use of the medical device

1) Before packaging medical devices, the packing material is inspected to check its integrity and assess its suitability with regard to the expiration date (final date of use). It is forbidden to use the packaging if it has holes, tears in the material, foreign inclusions, violations of seam integrity, signs of wetting (warping, streaks).

2) Select a pouch of such size that the medical device to be sterilised can be placed freely and occupies approximately $\frac{3}{4}$ of the internal volume of the pouch. This is necessary to prevent damage to the pouch when packing the article, for ease of sealing and to compensate for the expansion of the package during the vacuum stage (in case of steam sterilisation).

2a) From the rolls, pouches are prepared in advance in such a size that the medical device to be sterilised fits freely and occupies approximately $\frac{3}{4}$ of the internal volume of the pouch. Special cutting devices are used to cut sections corresponding to the length of the articles to be packed, taking into account the allowance for the free positioning of the article and the width of the seams.

Attention! Since all inscriptions and other markings, including information on the date of sterilization, may be indicated on the package only outside the area where the sterilized items are placed, it is recommended to leave a free margin of 1.5-2 cm behind the seam for applying such information.

One side of the roll section shall be heat-sealed using approved heat-welding machines designed for sealing medical sterilisation packaging materials. The width of the heat seam must be at least 6 mm. A properly made heat seam must be at least as strong as the finished side seams and the paper and plastic film must not be burnt.

3) The clean, dry, pre-sterilized, washed and dried medical devices to be sterilized shall be placed in optimally sized pouches as described above. The pouches should be filled to no more than $\frac{3}{4}$ of their volume to avoid tearing during sterilization. Medical devices that have a working and a non-working part (e.g. surgical instruments) are placed in the pouches with the non-working part facing towards the opening of the pouch. To minimize the risk of damage to the package by the edges of stabbing and cutting instruments when packing them, gauze, cloth or paper napkins are used to wrap sharp parts of the product.

4) Before closing the pouches, remove as much air as possible from the pouches by ironing it by hand in the direction from the closed end to the open end.

The unsealed side of pouches prepared from rolls as well as heat-sealing pouches packed with sterilisable articles shall be heat sealed using approved heat-welding machines designed for sealing medical sterilisation packaging materials. The width of the heat seam must be at least 6 mm. A properly executed heat seam should be at least as strong as the finished side seams and the paper and plastic film should not be burnt.

To facilitate sealing of the packages, additional space should be provided between the product and the sealed edge on the filling side, ensuring a distance between them of at least 30 mm. To avoid creasing during heat-welding of double packs, an additional 20-30 mm length allowance should be left on the outer package.

4a) Before closing self-sealing pouches, as much air as possible shall be removed from the pouch by ironing it by hand from the closed end to the open end. Then, pressing the smoothed edges of the pouch with one hand against the table, remove the protective paper strip from the adhesive strip

of the pouch. Fold the protruding part along the marked bend line and press the adhesive layer to the package surface, starting from the centre of the strip and smoothing to the edges. Iron the seam until all air bubbles are removed.

5) Pouches with sterilisable articles, especially heavy ones, are packed into sterilisation baskets, cassettes, pallets and placed vertically ("sideways"). The rule "paper to paper, film to film" shall be observed when combined pouches are placed side by side. Sterilisation of articles shall be carried out in accordance with the established procedure in accordance with the current regulatory documentation and according to the operating instructions. Recommended conditions:

- Autoclave steam sterilisation: 121°C, 20 min
- Ethylene oxide sterilisation: 55°C, 60 min
- Formaldehyde sterilisation: 55°C, 120 min
- Dry air sterilisation: 160°C, 150 min.

6) During sterilization processing there is a visually discernible change in the color of the chemical indicator for the corresponding method of sterilization, applied to the packaging material, which indicates the fact of sterilization by this method and allows to distinguish the sterilized product from non-sterilized. Full or partial preservation of the original color of the indicator, indicates non-compliance with the required conditions of effective sterilization. Medical devices with indicators that have not reached the final state are considered non-sterile and are not subject to repeated (additional) sterilization.

Packages on or in which there are moisture residues after sterilization by steam method, immediately after sterilization shall be dried, without opening, in a desiccator or air sterilizer at a temperature not exceeding 85°C.

7) After sterilization and before direct use sterilized products are stored in sealed bags in closed cabinets in clean dry rooms, at a temperature of +5 ° C to +40 ° C and relative humidity of 30% to 50%, avoiding direct sunlight and moisture.

Packages are opened by cutting off the edge of the package or by delamination. To do this, the loose ends of the paper and film are clamped with fingers and the adhesive seam is torn.

3 Warning and precautions

- The product is intended for single use only.
- Do not use a damaged product.
- Do not use a wet or damp product.
- Use of the product does not require special qualifications or special training.
- Do not use the product after its expiration date.
- Sterilization of packaged medical devices is carried out in sterilizers approved in the established order for use, in accordance with the current regulatory documentation for the sterilization method used and in accordance with the operating instructions for a particular type of sterilizer.
- Packages, on or in which there are moisture residues after sterilization by steam method, immediately after sterilization dried, without opening, in a desiccator or air sterilizer at a temperature not exceeding 85 ° C.

4 Specifications of the medical device

4.1 Description of the medical device

Rolls combined flat are flat sleeves made of medical paper (opaque side) and transparent polymer film (transparent side), joined together by thermal and/or adhesive seams on two sides (side seams). The paper side or the outer side of the sleeve film is appropriately labeled, as well as chemical indicators of class 1. In rolls combined with a fold, the increased surface area of the film is collected in folds on the edges of the roll, which allows to place inside the bulk sterilizable products.

Combined flat pouches are rectangular flat envelopes made of medical paper (opaque side) and transparent polymer film (transparent side), joined together by thermal and/or adhesive seam on three sides (side and bottom seams).

Flat kraft paper pouches are rectangular envelopes made of kraft paper, joined together by thermal and/or adhesive seam on three sides (side and bottom seams).

Combined flat kraft paper pouches are rectangular envelopes made of kraft paper (opaque side) and transparent polymer film (transparent side), joined together by thermal and/or adhesive seam on three sides (side and bottom seams).

On the inner surface of the protruding paper part of the self-sealing pouches there shall be a strip of heat-resistant adhesive, covered with a protective paper tape.

Sterility preservation period of medical devices after sterilisation - up to 5 years, but not more than the shelf life of the enclosed medical device and not more than the remaining shelf life of the packaging material.

Dimensions of the medical device are shown in Table 1.

Permissible dimensional deviations:

- for width and (or) length up to 150 mm, permissible dimensional deviation $\pm 15\%$;
- for width and (or) length from 150 mm to 500 mm, the permissible dimensional deviation is $\pm 10\%$;

Allowable dimensional deviation in length for all roll widths $\pm 5\%$.

Table 1 – Dimensions of the medical device:

Design version	Dimensions								
	Width, W, mm			Fold height, H, mm			Length, L, m		
	from	to	step	from	to	step	from	to	step
1. Roll combined flat disposable DENTAL FORMULA RKP	50	500	5	-	-	-	100	200	50
2. Roll combined gusseted disposable DENTAL FORMULA RKS	50	500	5	20	100	5	100	200	50
Design version	Dimensions								
	Width, W, mm			Fold height, H, mm			Length, L, mm		
	from	to	step	from	to	step	from	to	step
3. Pouch combined flat self-sealing disposable DENTAL FORMULA PKPS	30	1000	5	-	-	-	30	1000	5
4. Pouch combined flat disposable DENTAL FORMULA PKP	30	1000	5	-	-	-	30	1000	5

5. Kraft paper pouch flat self-sealing disposable DENTAL FORMULA PBPS	30	1000	5	-	-	-	30	1000	5
6. Bleached kraft paper pouch flat self-sealing disposable DENTAL FORMULA PBBPS	30	1000	5	-	-	-	30	1000	5
7. Kraft paper pouch combined flat self-sealing disposable DENTAL FORMULA PBKPS	30	1000	5	-	-	-	30	1000	5
8. Bleached kraft paper pouch combined flat self-sealing disposable DENTAL FORMULA PBBKPS	30	1000	5	-	-	-	30	1000	5

Along one edge of the rolls and pouches on the paper side and on the front of the kraft paper pouches are chemical indicators of class 1 according to ISO 11140-1 for the respective sterilization methods, which change their color after sterilization:

- for the dry heat method of sterilization - pink, with the explanation that it turns green after treatment with this agent;
- for the steam method of sterilization - blue, with the explanation that it turns grey-black after treatment with steam;
- for the gas sterilization method (using ethylene oxide) - pink, with the explanation that after treatment with this agent it becomes yellow;
- for the gas sterilization method (using formaldehyde) - crimson, with the explanation that after treatment with this agent it becomes green (the color shade can be from green to grey-green).

Complete or partial preservation of the original colour of the indicator indicates non-compliance with the required conditions of effective sterilisation. Medical devices with indicators that have not reached the final state are considered non-sterile and are not subject to repeated (additional) sterilisation.

5 Transportation, storage and disposal conditions of the medical device

5.1 Transportation conditions

The transportation of the device may be made by all means of transport (sea, land, railway or air transport) according to the rules of transportation, which are applied to such means of transport. During transportation, the storage conditions should not be violated.

5.2 Storage conditions

Packing material should be stored in the packaging of the manufacturer in a dry place, at a distance of at least 1 m from heating devices at a temperature from +5 ° C to +40 ° C and relative humidity not more than 85% during the entire shelf life. In case of partial consumption the remains of packaging materials should be stored in clean rooms under similar conditions.

5.3 Disposal conditions

Products that have fallen into disrepair, including due to expiration of their shelf life, as well as in case of package integrity failure, shall be disposed of as Class A waste in accordance with the requirements of SanPiN 2.1.3684-21.

Used products are classified as Class A waste and disposed of in accordance with the requirements of SanPiN 2.1.3684-21.

6 Information on completeness and packaging of the medical device

Packaging provides protection from mechanical and climatic factors during transportation and storage.

Packing material should be packed in bags made of polymeric materials, or wrapped in plastic film (including stretch film). Packaging material in group packing should be packed in a transport box made of corrugated cardboard, the transport packaging is glued with adhesive tape.

7 Scope of delivery of the medical device

The delivery set of packing material includes group packs of products in transportation containers and shortened consumer instruction for use, which defines the purpose and conditions of use of materials.

In a group package the packing material is included in the quantity from 1 to 12 rolls and from 5 to 3000 pouches.

8 Symbols used in medical device labelling

Refer to the instructions for use

Use only once

Do not use if the packaging is damaged

Manufacturer, name and address

Date of manufacture

Batch code


Use before

Keep away from moisture

Temperature range from +5°C to +40°C

Humidity range from 0% to 85%

Packaging material - high density polyethylene

Open here

9 List of national standards, applied by manufacturer (producer) of the medical device

Table 2 – List of standards applied by the Manufacturer.

Index	Name
ISO 9001:2015	Quality Management Systems. Requirements.
ISO 13485:2016	Medical Devices. Quality Management Systems. Requirements for Regulatory Purposes.
EN ISO 14971:2019	Medical Devices - Application of Risk Management to Medical Devices.
ISO 15223-1:2021	Medical Device-Symbols to be used with information to be supplied by the manufacturer - Part 1: General Requirements
ISO 11140-1:2014	Sterilization of medical products. Chemical indicators. Part 1. General requirements
ISO 11607-1-2019	Packaging for medical devices subject to final sterilization. Part 1. Requirements for materials, sterilization barrier systems and packaging systems

10 Manufacturer's warranty

The manufacturer (producer) guarantees the conformity of the quality of the packaging material to the requirements of this technical documentation in compliance with the rules of use, transportation and storage.

Safety and quality of the product are guaranteed during the entire shelf life.

The manufacturer is responsible for defects of the product, except for defects resulting from violation of the rules of use, transportation and storage conditions, or the action of third parties, or force majeure.

11 Data on Authorized representative of the manufacturer

Please contact the Authorized representative of the manufacturer for any and all purposes, related to circulation of the medical device on the territory of the Russian Federation:

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网站 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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ANQING BAOJIE PACKAGING CO., LTD

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SHORTENED CONSUMER INSTRUCTION FOR USE

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Version: 1.0 dated 25.06.2024

Packing material disposable for sterilization DENTAL FORMULA. Shortened consumer instruction for use.

Roll combined flat disposable DENTAL FORMULA RKP for steam and gas sterilization, dimensions: width from 50 to 500 mm with a step of 5 mm, length from 100 to 200 m with a step of 50 m

Kraft paper pouch flat self-sealing disposable DENTAL FORMULA PBPS for heat and steam sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm

PURPOSE. Packing material in rolls disposable for sterilization of medical devices DENTAL FORMULA is intended for disposable use for finishing sterilization of medical devices, maintaining sterility to the place of their use and observing aseptic properties.

RISK OF USE. The potential risk of use is class 1. There are no contraindications of the medical device and no side effects.

INSTRUCTION FOR USE. 1) Before packaging medical devices, the packing material is inspected to check its integrity and assess its suitability with regard to the expiration date (final date of use). It is forbidden to use the packaging if it has holes, tears in the material, foreign inclusions, violations of seam integrity, signs of wetting (warping, streaks).

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3) The clean, dry, pre-sterilized, washed and dried medical devices to be sterilized shall be placed in optimally sized pouches as described above. The pouches should be filled to no more than $\frac{3}{4}$ of their volume to avoid tearing during sterilization. Medical devices that have a working and a non-working part (e.g. surgical instruments) are placed in the

Roll combined gusseted disposable DENTAL FORMULA RKS for steam and gas sterilization, dimensions: width from 50 to 500 mm with a step of 5 mm, fold from 20 to 100 mm with a step of 5 mm, length from 100 to 200 m with a step of 50 m

Bleached kraft paper pouch flat self-sealing disposable DENTAL FORMULA PBBPS for heat and steam sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm

Pouch combined flat self-sealing disposable DENTAL FORMULA PKPS for steam and gas sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm

Kraft paper pouch combined flat self-sealing disposable DENTAL FORMULA PBKPS for heat and steam sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm

Pouch combined flat disposable DENTAL FORMULA PKP for steam and gas sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm

Bleached kraft paper pouch combined flat self-sealing disposable DENTAL FORMULA PBBKPS for heat and steam sterilization, dimensions: width from 30 to 1000 mm with a step of 5 mm, length from 30 to 1000 mm with a step of 5 mm

pouches with the non-working part facing towards the opening of the pouch. To minimize the risk of damage to the package by the edges of stabbing and cutting instruments when packing them, gauze, cloth or paper napkins are used to wrap sharp parts of the product.

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4a) Before closing self-sealing pouches, as much air as possible shall be removed from the pouch by ironing it by hand from the closed end to the open end. Then, pressing the smoothed edges of the pouch with one hand against the table, remove the protective paper strip from the adhesive strip of the pouch. Fold the protruding part along the marked bend line and press the adhesive layer to the package surface, starting from the centre of the strip and smoothing to the edges. Iron the seam until all air bubbles are removed.

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7) After sterilization and before direct use sterilized products are stored in sealed bags in closed cabinets in clean dry rooms, at a temperature of +5 °C to +40 °C and relative humidity of 30% to 50%, avoiding direct sunlight and moisture.

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APPLICATION CONDITIONS. The packing material is used in therapeutic or therapeutic-preventive medical institutions in a closed room in a dry place, away from direct sunlight, at temperatures from 10°C to 30°C and relative humidity not more than 80%, in accordance with the instructions for the medical device.

Sterility preservation period of medical devices after sterilization - up to 5 years, but not more than the shelf life of the enclosed medical device and not more than the remaining shelf life of the packaging material.

TRANSPORTATION AND STORAGE CONDITIONS. The transportation of the device may be made by all means of transport (sea, land, railway or air transport) according to the rules of transportation, which are applied to such means of transport. During transportation, the storage conditions should not be violated.

Packing material should be stored in the packaging of the manufacturer in a dry place, at a distance of at least 1 m from heating devices at a temperature from +5 °C to +40 °C and relative humidity not more than 85% during the entire shelf life. In case of partial consumption the remains of packaging materials should be stored in clean rooms under similar conditions. Do not use the products after their expiration date.

DISPOSAL. Products that have fallen into disrepair, including due to expiration of their shelf life, as well as in case of package integrity failure, shall be disposed of as Class A waste in accordance with the requirements of SanPiN 2.1.3684-21. Used products are classified as Class A waste and disposed of in accordance with the requirements of SanPiN 2.1.3684-21.

WARRANTY. The manufacturer (producer) guarantees the conformity of the quality of the packaging material to the requirements of this technical documentation in compliance with the rules of use, transportation and storage.

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Registration certificate P3H XXXX/XXXX of XX.XX.XXXX.

In the event of an adverse event with signs of an adverse event (incident), it is necessary to notify the authorized representative of the manufacturer in the territory of the Russian Federation about it: "BE PLUS" Co., Ltd., Russia, 194292, St. Petersburg, ext. ter. c. municipal district Sergievskoe, 12, 6-y Verkhniy pereulok, litera A, part of office 5.3.1.. Tel.: +7 (812) 309-36-62; e-mail: beplus@beplus.spb.ru

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
Open here

Instructions for use of the medical device "Packing material disposable for sterilization DENTAL FORMULA" is issued for the first time.

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