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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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
THIS IS TO CERTIFY THAT: the seal of HUQIU IMAGING
(SUZHOU) CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Released for the first time

Table of Contents

Introduction	4
1. Medical device name.....	5
2. Information about the manufacturer of the medical device and its authorized representative, including the name of the legal entity, postal address, telephone, fax numbers, and E-mail address	7
3. Purpose of the medical device indicating the user	7
4. Product description and operating principle	8
5. Product functional characteristics	11
6. Generalized results of clinical trials conducted for the purpose of registration of a medical device, or a link to a source where such information is available to the user	12
7. Technical characteristics necessary for the user to use the product for the intended purpose determined by the manufacturer.....	13
8. Information on the availability of a drug, biological material and nanomaterial.....	14
9. Installation and commissioning information.....	14
10. Premises requirements	14
11. Product Operation	14
12. Operational Maintenance.....	23
13. Information on cleaning, disinfecting and sterilizing the product	26
14. Precautions	27
15. Transporting the product.....	28
16. Product storage	29
17. Product disposal	29
18. Warranty.....	29
19. Service Manual.....	30
20. List of international standards to which a medical device meets	34
21. Scope of delivery.....	34
ACCEPTANCE AND SALE CERTIFICATE.....	36
WARRANTY CERTIFICATE.....	37
EMC INFORMATION.....	38

Introduction

Instructions for use the product are provided to the user in Russian. The instructions for use compiled by the manufacturer include an operating manual, instructions for installation, maintenance, repair, transportation, storage, and disposal of the product.

Instructions for use the product are provided to the user in electronic form on a CD. Instructions for use the product are also provided to the user on paper.

The method of providing instructions for use is suitable and accessible to all users. With the Instructions for use provided on the CD, the user has the opportunity to:

1) View the Instructions for use on the user's computer screen, which is used in conjunction with the product. The minimum technical requirements for the user's computer to view the operational document correspond to the minimum technical requirements for the user's computer working in conjunction with the product;

2) Receive the current version of the Instructions for use, identical to the Instructions for use on paper;

3) Obtain an additional paper version of the Instructions for use using office printing tools.

Symbols used in the document:

Symbol	Description
	Attention!
	Note
	EC Declaration of Conformity Marking
	Manufacture date
	Manufacturer
	Serial number
	"Protect from moisture".
	"Top"
	"Stacking is prohibited"
	"Don't turn over"
	«ON»
	«OFF»

1. Medical device name**Dry Imager Printer****Versions:**

- 1. HQ-430DY Dry Imager Printer consisting of:**
 - 1.1.HQ-430DY Dry Imager Printer – 1 pc.
 - 1.2.Power cable – 1 pc.
 - 1.3.Network cable – 1 pc.
 - 1.4.Plates – 2 pcs.
 - 1.5.Rubber rollers – 4 pcs.
 - 1.6.Cloth for the printer's thermal head – 1 pc.
 - 1.7.CD with instructions – 1 pc.
 - 1.8.Operating manual – 1 pc.
 - 1.9.Data sheet – 1 pc.
- 2. HQ-450DY Dry Imager Printer consisting of:**
 - 2.1.HQ-450DY Dry Imager Printer – 1 pc.
 - 2.2.Power cable – 1 pc.
 - 2.3.Network cable – 1 pc.
 - 2.4.Plates – 2 pcs.
 - 2.5.Rubber rollers – 4 pcs.
 - 2.6.Cloth for the printer's thermal head – 1 pc.
 - 2.7.CD with instructions – 1 pc.
 - 2.8.Operating manual – 1 pc.
 - 2.9.Data sheet – 1 pc.
- 3. HQ-460DY Dry Imager Printer consisting of:**
 - 3.1.HQ-460DY Dry Imager Printer – 1 pc.
 - 3.2.Power cable – 1 pc.
 - 3.3.Network cable – 1 pc.
 - 3.4.Plates – 2 pcs.
 - 3.5.Rubber rollers – 4 pcs.
 - 3.6.Cloth for the printer's thermal head – 1 pc.
 - 3.7.CD with instructions – 1 pc.
 - 3.8.Operating manual – 1 pc.
 - 3.9.Data sheet – 1 pc.
- 4. HQ-762DY Dry Imager Printer consisting of:**
 - 4.1.HQ-762DY Dry Imager Printer – 1 pc.
 - 4.2.Power cable – 1 pc.
 - 4.3.Network cable – 1 pc.
 - 4.4.Plates – 2 pcs.
 - 4.5.Rubber rollers – 4 pcs.
 - 4.6.Cloth for the printer's thermal head – 1 pc.
 - 4.7.CD with instructions – 1 pc.
 - 4.8.Operating manual – 1 pc.
 - 4.9.Data sheet – 1 pc.
- 5. HQ-760DY Dry Imager Printer consisting of:**
 - 5.1.HQ-760DY Dry Imager Printer – 1 pc.
 - 5.2.Power cable – 1 pc.
 - 5.3.Network cable – 1 pc.
 - 5.4.Plates – 2 pcs.
 - 5.5.Rubber rollers – 4 pcs.
 - 5.6.Cloth for the printer's thermal head – 1 pc.

- 5.7.CD with instructions – 1 pc.
- 5.8.Operating manual – 1 pc.
- 5.9.Data sheet – 1 pc.

hereinafter referred to as: printer, product

For more information about the characteristics and differences between product versions, see 5.1.


Figure 1 – Dry Imager Printer, general view of the HQ-430DY, HQ-460DY, HQ-760DY versions


Figure 2 – Dry Imager Printer, version HQ-450DY, general view

Figure 3 – Dry Imager Printer, version HQ-762DY, general view

2. Information about the manufacturer of the medical device and its authorized representative, including the name of the legal entity, postal address, telephone, fax numbers, and E-mail address

Developer and Manufacturer	Huqiu Imaging (Suzhou) Co., Ltd (Хуцю Имэджинг (Сучжоу) Ко., Лтд.) No. 19 Tianmushan Road, SND, Suzhou, 215163, Jiangsu, People's Republic of China
Production site	Huqiu Imaging (Suzhou) Co., Ltd (Хуцю Имэджинг (Сучжоу) Ко., Лтд.) No. 19 Tianmushan Road, SND, Suzhou, 215163, Jiangsu, People's Republic of China
Authorised representative	Limited Liability Company "Salute-Service" (Salute-Service LLC) 432071, Ulyanovsk region, Ulyanovsk, Federation st., 63 Tel.: 8 (8422) 58-16-00; 8 800 511 65 38 E-mail: salute73@inbox.ru

3. Purpose of the medical device indicating the user

This product is designed to capture and display digital images generated by diagnostic imaging systems such as computed tomography (CT), magnetic resonance imaging (MRI), positron emission tomography (PET), gamma camera, or ultrasound, as well as for their reproduction in various image formats on X-ray thermographic film.

3.1. Scope and potential users:

Above printers may only be used for general radiography examinations. The HQ-762DY and HQ-760DY printers can also be used for mammography examinations.

The printers are designed for high-performance operation and can be used as central network printers for printing medical images transmitted over the hospital network in DICOM format.

Image sources include images from CT, MRI, digital and computed radiography systems, fluoroscopy, PACS, nuclear medicine, mobile radiography, dental and other medical imaging systems.

The product is intended for use in medical institutions, clinics, clinical dispensaries and hospitals.

Potential user: medical personnel of appropriate qualifications.

3.2. Indications for use

The need to print DICOM medical images on films of various formats when images are created using diagnostic imaging systems, for example: X-ray, computed tomography (CT), magnetic resonance imaging (MRI), positron emission tomography (PET), gamma camera, or ultrasound.

3.3. Contraindications

Not found

3.4. Side effects

Not found

3.5. Residual risks

Residual risks are insignificant and acceptable in accordance with the requirements adopted by regulatory documents.

4. Product description and operating principle

4.1. Operating principle of the product

The product implements the principle of direct thermal printing to reproduce high-quality digital diagnostic images obtained on various digital diagnostic devices with digital image recording, including radiographic systems, CT and MRI systems, via a network. DICOM files are used to print medical digital diagnostic images. For printing, dry type thermographic film must be used, registered in the register of medical devices. The printer allows you to print high quality black and white digital medical images on thermographic film with high resolution.

4.2. Explanation of the properties and characteristics of the medical device

Using the principle of printing images on thermographic film allows you to effectively perform work in natural light, without the use of additional reagents, because printing method is dry one.

The product complies with IEC 61010-1(2010) protection class IP 20 against moisture and dust.

An overcurrent protection circuit is built into the device's power supply. When the device experiences overcurrent, the fuse (type BOJACK F6.3AL250V, size 5x20mm, 6.3A, 250V) opens and stops supplying power.

Note:

1) Consumer independently ensures the availability of a fuse with appropriate characteristics for replacement;

2) If necessary, the fuse is replaced by the consumer (see cl. 19.4).

4.3. Description of components

Along with the product, the manufacturer supplies component parts intended for use in combination with the printer:

- Power cable – designed to connect the product to the power supply;

Figure 4 – Power cable

- Network cable – designed for connecting to a local/Internet network;

Figure 5 – Network cable

- Plates – used when it is necessary to adjust the height of the tray: the plates are placed under the film cartridge in the printer tray.

Figure 6 – Plates

- Rubber rollers – designed to press the film for better printing. When erased, the rollers installed in the printer must be replaced.

Figure 7 – Rubber rollers

- Cloth for printer's thermal printhead – designed to clean the thermal head of the printer from accumulated dust and dirt that remain after heat treatment and thermal transfer, accumulate and lead to premature wear of the thermal elements;

- CD with instructions – CD with electronic instructions for user convenience.

Figure 9 – CD with instructions

- Instructions for use – instructions on paper.
- Data sheet – passport with technical data

5. Product functional characteristics

5.1. Product characteristics

Unspecified permissible parameter deviations are $\pm 10\%$

No.	Parameter name	Product modifications				
		HQ-430DY	HQ-450DY	HQ-460DY	HQ-762DY	HQ-760DY
1	Printing method	dry, film installation in normal lighting				
2	Print resolution, dpi (pixel/mm)	320 (12.6)			508 (20)	
3	Print speed, sheets per hour	up to 70 pcs. (14"x17")	up to 80 pcs. (14"x17")	up to 70 pcs. (14"x17")	up to 60 pcs. (14"x17")	up to 60 pcs. (14"x17")
		up to 80 pcs. (11"x14")	up to 85 pcs. (11"x14")	up to 80 pcs. (11"x14")	up to 65 pcs. (11"x14")	up to 65 pcs. (11"x14")
		up to 85 pcs. (10"x12")	up to 90 pcs. (10"x12")	up to 85 pcs. (10"x12")	up to 70 pcs. (10"x12")	up to 70 pcs. (10"x12")
		up to 90 pcs. (8"x10")	up to 100 pcs. (8"x10")	up to 90 pcs. (8"x10")	up to 80 pcs. (8"x10")	up to 80 pcs. (8"x10")
4	Operational readiness time after switching on, minutes, no more	3				
5	Grayscale contract resolution, bit	14				
6	Film feeding system	sheet-by-sheet				
7	Number of feed trays	1 tray	2 trays	4 trays	2 trays	
8	Film parameters	medical thermographic film (blue) conforming to ISO 4090				
9	Capacity, sheets: feed cassette/tray	50... 100 / 100				
10	Standard brand and size of film, inch	14"x17", 11"x14", 10"x12", 8"x10"				
11	Standard film sizes, cm	35x43, 28x35, 25x30, 20x25				
12	Interfaces	10/100/1000 Base-T Ethernet (RJ-45)				
13	Image file format/network protocols	DICOM/DICOM 3.0, network printing				
14	Adjusting the image quality	automatic calibration using built-in densitometer				
15	Provided image control	by contrast, halftone, polarity, rotation, saturation, scale				
16	Control panel	LCD touch screen (4"), online display, modes: Ready, Printing, Notification, Error				
17	Processor brand	Intel atom				
18	Amount of RAM/independent memory (hard drive), Gb	2 / 500				
19	Built in software: name, version, release date	HQPrint version: V2 and higher, release date: 18.05.2020				
20	Input power during printing / in standby mode, W	600±75				
21	Rated supply voltage, V	100 – 240	100 – 240	100 – 240	220 – 240	100 – 240
22	Supply frequency, Hz	50±60				
23	Maximum current consumption, mA	500				
24	Mass, kg	50	50	50	75	50
25	Weight in transport packaging, kg, no more	63	63	63	100	63
26	Overall dimensions (length×width×height), mm	740×566×445	730×566×415	740×566×445	730×525×1006	740×566×445
27	Overall dimensions of transport packaging (length×width×height), mm	870×700×495	870×700×495	870×700×495	810×620×1170	870×700×495
28	Continuous operation time, h, no more	8				
29	Fuse type	type BOJACK F6.3AL250V, size 5x20 mm, 6.3A, 250V				

Noise level: During operation - Max. 64 dBA; Standby - Max. 54 dBA

5.2. Product components characteristics

No.	Name	Parameters
1	Power cable, 1 pc.	
1.1	Type of input/output connections	IEC 320 C14 / CEE 7/7
1.2	Enclosure protection class according to IEC 60529	IEC 53(RVV), CCC
1.3	Rated voltage, V	300±500
1.4	Length, m	2
1.5	Weight, g	230
2	Network cable, 1 pc.	
2.1	Twisted pair type	UTP
2.2	Category	5e
2.3	Connector type	RJ-45
2.4	Length, m	5
2.5	Weight, g	150
3	Plates, 2 pcs.	
3.1	Overall dimensions (length×width×height), mm	(140×100×03)
3.2	Weight, g	61
4	Rubber rollers, 4 pcs.	
4.1	Mark	AVD-25
4.2	Overall dimensions (length×diameter), mm	(18×20)
4.3	Weight, g	14
5	Cloth for printer's thermal head, 1 pc.	
5.1	Material	Non-woven lint-free
5.2	Overall dimensions (length×width×height), mm	(225×225×01)±2 mm
5.3	Weight, g	13
6	CD with instructions, 1 pc.	
6.1	Size	CD in cardboard packaging
6.2	Overall dimensions (length×width×height), mm	(125×125×01)
6.3	Weight, g	20

6. Generalized results of clinical trials conducted for the purpose of registration of a medical device, or a link to a source where such information is available to the user

The Dry Imager Printer has passed clinical trials.

During clinical trials, the following was confirmed:

- compliance of the medical device with regulatory documentation, technical and operational documentation;
- compliance of the documentation with the established purpose and indications for use;
- completeness and reliability of the characteristics of the medical device established by the regulatory documentation, technical and operational documentation of the manufacturer;
- quality of the medical device, effectiveness and safety of its use.

7. Technical characteristics necessary for the user to use the product for the intended purpose determined by the manufacturer

To print medical digital diagnostic images, the product must be connected to medical equipment with digital image registration via a network, the user's computer or the network to which the user's (operator's) work computer is connected. Via the Ethernet interface (RJ-45), the printer receives data for printing medical digital diagnostic images. The printer prints on thermographic film.

7.1. Minimum technical specifications of the user's computer used with the product

Element	Parameter
Monitor	19 inches
Computer	CPU: Intel Pentium D HDD: 500 Gb RAM: 2 Gb
Operating system	Windows 7 Professional, service pack 1, Mac OS X 10.9; Ubuntu 14.04 (64-bit);
Data input/control devices	Mouse, keyboard
Type of input and output connections	Ethernet (RJ-45)
Additional software	Program for viewing images in DICOM format: View Form (Philips), Inobitek DICOM viewer, E-Film, Asteris, Horos, etc.

7.2. Medical digital diagnostic images' characteristics for printing in DICOM format

Tags ID	Name	Comments
(0018.5101)	View Position	Correct values are CC/MLO
(0028.1351)	Partial View Description	
(0008, 103E)	Series Description	
(0054.0220)	View Code Sequence	
(0008.0100)	Code Meaning	
(0020.0020)	Patient Orientation	Correct values – R/L
(0020.0062)	Image Laterality	
(0020.0060)	Laterality	
(0020.0020)	Patient Orientation	Correct values – MG
(0008.0060)	Modality	
(0028.0004)	Photometric Interpretation	Correct values – MONOCHROME1/MONOCHROME2
(0020.0062)	Image Laterality	Correct values – U
(0018.5101)	View Position	Correct values – AP
(0008.0060)	Modality	Correct values – DX
(0028.0004)	Photometric Interpretation	Correct values – MONOCHROME1/MONOCHROME2

7.3. Characteristics of the consumables used with the product – thermographic film*

No.	Name Characteristics	Options for thermographic medical film sizes			
		8"×10" (20×25 cm)	11"×14" (28×35 cm)	14"×17" (35×43 cm)	10"×12" (25×30 cm)
1	Size error, cm	±1			
2	Thickness of one sheet, microns	(197±1)			
3	Density, g/m ²	260±20			
4	Maximum print resolution (px size), microns	20			
5	Transmission of halftones, bits	14			
6	Grayscale display	4096			
7	Range of reproducible optical density, B:	≥2.9			
8	Minimum optical density	≤0.08			
9	Average gradient, G	1.5~2.1			
10	Sensitivity, S	25~60			
11	Sensitivity to light	not sensitive			
12	Colour	Blue			
13	Regulatory Compliance	ISO 4090			

*not included in the main product, supplied separately

8. Information on the availability of a drug, biological material and nanomaterial

The product and its components do not contain drugs, biological materials or nanomaterials.

9. Installation and commissioning information

Installation of the printer in winter should be done only after keeping it unpacked in a closed heated room for at least 4 hours.

The following requirements apply to printer placement:

- Do not install it in an inclined position;
- The power socket used to connect the power plug must be in an accessible location so that you can quickly disconnect the printer from the network;
- Do not turn on the printer without protective grounding;
- It is recommended to use an uninterruptible power supply device;
- It is not recommended to place or throw heavy objects on the power cord or damage it in any way.

Only dry thermographic film that meets specifications 5.1 should be used for printing. Using film other than the recommended one may cause it to jam as it passes through the thermal printhead, which will negatively affect print quality and cause the printer's thermalhead to fail.

10. Premises requirements

The product is intended for indoor use only and is typically used in a laboratory/laboratory environment.

In the rooms intended for operating the printer, there must be no aggressive environments, the mass concentration of dust in the air should not exceed 0.75 mg/m², the electrical component of the electromagnetic field should not exceed 3 V/m in the frequency range of 150 kHz to 300 MHz.

Permissible ambient temperatures during operation: +5...+40 °C, humidity: ≤80% Atmosphere Pressure Range: 86 – 106 kPa

It is not recommended to place the printer in the following places:

- with mechanical vibration;
- with high humidity;
- with excessive exposure to direct sunlight;
- with very low or high temperature.

It is not recommended to use the printer in close proximity to flammable materials.

Moving the printer to another location, if necessary, must be carried out by at least two people.

11. Product Operation

11.1. Preparation for use

When preparing the printer for use, in order to ensure safety precautions, it is prohibited:

- start work without studying these Instructions for use;
- use elements, fixtures, tools and devices not provided for in the design.

Before starting operation, you should make sure by external inspection that the printer is in good condition and the following components are included in the package:

- Printer – 1 pc.
- Power cable – 1 pc.
- Network cable – 1 pc.
- Plates – 2 pcs.
- Rubber rollers – 4 pcs.
- Cloth for the printer's thermal head – 1 pc.
- CD with instructions – 1 pc.
- Instructions for use – 1 pc.
- Data sheet – 1 pc.

The locations of components and controls are shown in the Figures 10-14.


Figure 10 – Main internal components of the HQ-450DY, HQ-460DY, HQ-760DY versions


Figure 11 – Main internal components of the HQ-430DY version


Figure 12 – Main internal components of the HQ-762DY version


Figure 13.1 – Main external components of the product (HQ-450DY, HQ-460DY, HQ-760DY, HQ-430DY)


Figure 13.2 – Main external components of the product (HQ-762DY)

LCD screen	Displays the size and number of remaining film sheets in the trays. Equipped with a touch panel that navigates through the printer menu. Error messages are displayed when an error condition occurs.
Film exit hole (tray)	Printed sheets of film (prints) are ejected from the upper and lower film cassettes into this tray.
POWER button	Press the switch to turn on the device's power. When the power is turned on, the indicator turns green. To turn off the device, press the switch again. The indicator will flash and then go off. Pressing the switch while the indicator is flashing turns the device's power back on.
PAUSE button	Used to pause printing images
Upper Film Cassette Compartment	Designed for use with KX410 A/F film sheets
Lower Film Cassette Compartment	Designed for use with KX410 A/F film sheets
Printer top cover	Opens to clear jamming the film loaded in the film cassette
Printer back cover	Opens to clear jamming the film loaded in the film cassette


Figure 14.1 – Connectors located on the product (HQ-450DY, HQ-460DY, HQ-760DY, HQ-430DY)


Figure 14.2 – Connectors located on the product (HQ-762DY)

USB interface	Connector for connecting a keyboard and mouse to the printer for making changes to settings, USB interface type is USB 2.0
Network interface	Compliant with 100BASE-TX standard. Used to connect to a LAN network (RJ-45 modular connector)
Power connector	AC connector – AC IN Use an appropriate local power cord (supplied).
Power switch	This switch is the main power switch for the device (turning the device on or off). When the switch is set to on, power is supplied to the device when you press the POWER button on the front panel.

11.2. Turn on

Connect the device to a LAN network of 100BASE-TX LAN standard, having previously specified the required IP address.

Connect the power cable with the appropriate connector to the printer, and the other connector to the power supply.

Switch the power switch “I/O” on the rear panel of the printer to the “I” position.

The device does not turn on even if the power switch “I/O” on the rear panel of the printer is turned to the “I” position.

Press the POWER button on the printer's front control panel to turn on. The appearance and location of the control panel are shown in Figures 15 and 16. Functionally, the control panels of all versions are the same.


Figure 15 – Front control panel of the Dry Imager Printer, versions HQ-430DY, HQ-460DY, HQ-760DY

Figure 16 – Front control panel of the Dry Imager Printer, versions HQ-450DY, HQ-762DY

The printer's LCD touch screen will light up, the printer is preparing for operation and warming up (see Figure 17)


Figure 17 – Warm-up mode on the printer's LCD touch screen

BoxUp	Upper feed cassette, size
BoxDn	Bottom feed cassette, size
Num	Film quantity
Mission	Current queue of jobs waiting to be printed
Wait	Expect

The printer then goes into ready-to-use mode (Figure 18).

Figure 18 – The printer is ready for use, the upper and lower trays are shown as full.

BoxUp	Upper feed cassette, size
BoxDn	Bottom feed cassette, size
Num	Film quantity
Mission	Tasks
Ready	Ready for work

When using the printer for the first time, there will be zero film remaining in the upper and lower feed cassettes because the film sheets have not yet been loaded into the printer.

When reusing the product, the indicator shows the number of sheets of film loaded in the upper and lower cassettes. After this, the printer will be ready to print.

11.3. Turn off

To turn off the printer, press the POWER button to enter the Close Image menu (Turning off the printer).

Press the "OK" button to turn off or "Cancel" to return to the previous menu. The printer will turn off.


Figure 19 – Close Image Submenu (Turning off the printer)

You can also turn off the printer from the Main Menu (Figure 22) by clicking the Close Image button (Turning off the printer).

11.4. Product Application

Inserting film into a cassette:

- pull out the printer cassette compartment (Figure 20).


Figure 20 – Removing the printer cassette compartment

- Remove the film cassette from the cardboard packaging, remove the protective packaging and place it in the printer cassette compartment (Figure 21).


Figure 21 – Placing the film cassette in the printer cassette compartment

Hold the film using a moisture-proof bag or gloves. Do not touch the film itself. Dust, fingerprints and sweat can cause poor print quality.

- Reinstall the printer cassette compartment.

11.5. Print

Before printing, make sure the product is properly connected to the LAN.

Make sure the film sheets are loaded correctly into the printer cassette.

Press the POWER button on the printer's front control panel to turn on. After a few seconds, a message will appear on the LCD touch screen. If the dimensions and number of remaining sheets of film (Figure 16) are displayed, the printer is ready to print.

Transfer the image from the digital image registration device or the user's computer.

When the printer status changes, the message on the LCD touch screen will change to a new one.

Message	Printer status
Processing	Image processing in progress
Printing	Printing in progress
Reserved	Expectation

When printing is complete, the printer will return to standby mode. Information about the size and amount of remaining film will be displayed on the LCD touch screen (Figure 16).

It is not recommended to remove the film while the printer is in the printing state.

If you remove the film by force, the surface of the film may become dirty or damaged.

It is recommended not to open the printer covers while the printer is in the printing state.

If the outer covers are open while the printer is printing, they will close automatically if damaged. At this time, the printer emits a characteristic sound signal, and a corresponding message appears on the LCD touch screen indicating that the printer covers are open.

11.6. Consumption of film sheets during printing.

If the film cassette runs out of film while the printer is running, "SET FILM" will appear on the LCD touch screen.

When the current image has finished printing, the printer will stop and you can add film to the cassette. After inserting the film cassette into the printer, printing will resume.

Do not expose printed images to direct sunlight or bright room lighting. Store images according to the film manufacturer's recommendations. Unless otherwise stated, it is recommended to store images in a polypropylene bag in a closet at a temperature not exceeding +25°C.

For other messages on the printer's LCD touch screen, see cl. 20.3.

11.7. Working with the printer menu

1) Main menu

When you touch the icon  on the printer's LCD touch screen (Figure 18), the Main Menu (Figure 22) opens. To select the appropriate submenu, press the corresponding on-screen button.


Figure 22 – Main Menu

Mission Manage	Managing received tasks
Sample Print	Test print (of sample)
Service Mode	Service section
Close Imager	Turning off the printer

To return to the printer ready screen (Figure 18), press .

2) Managing received tasks

From the Main Menu screen (Figure 22), click the Mission Manage button. This submenu is designed to receive information about current print jobs received by the printer from connected medical devices/computers and control their printing (Figure 23).


Figure 23 – Mission Manage submenu

To return to the printer ready screen (Figure 18), press .

To return to the previous screen, press the sign .

3) Test print (of sample)

Clicking the Sample Print button from the Main Menu (Figure 22) opens the corresponding submenu (Figure 24).


Figure 24 – Sample print submenu

Print Size	Film size for printing
------------	------------------------

Select the film size (in inches, see cl. 8.3 for size matching) for printing from data presented. The printer will make a test print.

To return to the printer ready screen (Figure 18), press .

To return to the previous screen, press the sign .

4) Service Mode

When you click on the Service Mode button from the Main Menu (Figure 22), the corresponding submenu (Figure 25) opens.


Figure 25 – Service Mode submenu

Service	Service
Sensor Status1	Status of the sensor 1
Sensor Status2	Status of the sensor 2
Motor Control	Feed control
Error Log	Error log

This submenu is INFORMATIONAL for the user. To change settings/configure/reconfigure the printer, contact the manufacturer (Authorized Representative) or an authorized service center.

To return to the printer ready screen (Figure 18), press .

To return to the previous screen, press the sign .

12. Operational Maintenance

12.1. Contents and interval of maintenance, including cleaning and disinfection of the medical device

Cleaning and disinfecting the printer are carried out according to section 13.

- Before maintaining electrical equipment, make sure that it is not energized.
- Do not remove identification labels from the printer: they are marked with a serial number and other useful information.
- Changes and modifications to the printer by the user that are not approved by the manufacturer (Authorized Representative) are prohibited.

1) Maintenance of the printer as a whole consists of regularly inspecting all working components for serviceability and integrity; special attention should be paid to searching for irregularities in electrical connections, damage to protective coatings, integrity of the housing, insulation of electrical wiring, and grounding devices.

2) Maintenance of the printer is carried out only after it has been de-energized.

3) When using the printer, you should regularly, but at least once a week, clean the external surfaces of dust and dirt (see cl. 13.2).

4) During daily inspection, you should ensure that the printer is being used correctly. Inspection before using the printer should include the following checks:

- the printer starts in normal mode;
- network connectors are connected correctly;
- status of the power switch (checked upon completion of operation).

5) Routine maintenance:

It is recommended to perform the following maintenance tasks every 1000 pages printed:

- Cleaning the thermal printhead;

- Cleaning the rollers.
- 6) Unscheduled cleaning of the thermal printhead is carried out in the following cases:
 - if contamination is visually detected;
 - vertical stripes of white or black color are visible on the film after printing.
- 7) To clean the thermal printhead, proceed as follows:
 - Turn off the printer power and allow it to cool completely before cleaning the thermal printhead;
 - Open the back cover and loosen the fastening bolt (Figure 26);


Figure 26 – Removing the back cover

- Tilt the thermal printhead against the mounting axis (Figure 27);


Figure 27 – Removing the thermal printhead

- Clean the thermal printhead. To clean the thermal printhead, use a cloth (not linen) completely soaked in surgical spirit (Figure 28).


Figure 28 – Cleaning the thermal printhead

- Set the thermal printhead to its original state;
 - Tighten the mounting bolt;
 - Close the back cover.
- Cleaning the rollers is carried out as follows:
- Turn off the printer power and allow it to cool completely before cleaning the rollers;
 - Open the back cover and loosen the fastening bolt (Figure 29).


Figure 26 – Removing the back cover

- Tilt the thermal printhead against the mounting axis (Figure 27);


Figure 27 – Removing the thermal printhead

- Clean the rollers (Figure 31).

To clean the rollers, use a cloth (not linen) from the supply kit or bologna.


Figure 31 – Roller placement

- Set the thermal printhead to its original state;
- Tighten the mounting bolt;
- Close the back cover of the printer.

The manufacturer (Authorized Representative) does not take responsibility for unprofessional repair attempts and the consequences caused by this!

The manufacturer supplies rubber rollers with the printer. When erased, the rollers installed in the printer must be replaced.

12.2. The need for calibration to ensure proper and safe operation of a medical device during its service life. The estimated service life of the product is 15 years.

There is no need for calibration.

Image quality is controlled by automatic calibration using the printer's built-in densitometer.

12.3. Methods to reduce risks associated with the installation, calibration, or maintaining of a medical device.

Residual risks are insignificant and acceptable in accordance with the requirements adopted by regulatory documents. Only personnel who have studied these instructions for use and have experience in operating office equipment are allowed to service the printer.

13. Information on cleaning, disinfecting and sterilizing the product

13.1. Precautions for cleaning and disinfection

Risk of electric shock if liquid enters the device.

Before cleaning or disinfecting, turn off the product and unplug it from the electrical outlet. Do not allow liquids to enter the housing. Do not use spray for cleaning/ disinfect the housing. Connect the printer to the power supply only after it is completely dry inside and out.

Risk of product damage due to aggressive substances.

Do not use aggressive substances such as strong or weak alkalis, strong acids, acetone, formaldehyde, halogenated hydrocarbons, or phenol to clean the printer or its components. If the device becomes contaminated with aggressive substances, immediately clean it with a mild cleaning agent.

Corrosion of product parts due to the use of aggressive cleaning and disinfection agents.

Do not use caustic cleaners, aggressive solvents or abrasives for polishing. Do not expose printer components to prolonged exposure to aggressive cleaning and disinfection agents.

Printer damage from ultraviolet or other high-energy radiation.

Do not disinfect it using UV, beta or gamma rays or other high energy radiation. Do not store the printer in areas exposed to UV radiation.

13.2. Cleaning and disinfecting the printer

If you have questions about cleaning and disinfection, as well as disinfection and the cleaning products used, please contact the Technical Support Center of Huqiu Imaging (Suzhou) Co., Ltd., or the Authorized Representative of the manufacturer in Russia, Salute-Service LLC, 432071, Russian Federation, Ulyanovsk region, Ulyanovsk, st. Federation, 63. Telephone: 8 (8422) 58-16-00; 8 800 511 65 38. E-mail: salute73@inbox.ru.

General recommendations for cleaning the device:

1. Turn off the device using the power switch.
2. Disconnect the power plug from the power supply.
3. Clean and disinfect all accessible surfaces of the printer, including the power cord, using a clean microfiber or chamois cloth and 96.4-97% surgical spirit.
4. Connect the device to the power supply only after it is completely dry.

External surface

The outer surface and housing of the printer are treated with a clean microfiber or chamois cloth and 96.4-97% surgical spirit.

Do not sterilize the printer with ethylene oxide gas or vapor.

13.3. Criteria for product unsuitability

The printer is unusable if, after work carried out in accordance with the instructions in section 19, the following occurs:

- the printer does not turn on;
- the control panel is inactive (LCD touch screen is dark);
- the printer does not pick up film from the tray;
- printed images have defects or streaks.

13.4. Product sterilization

The product is non-sterile and does not require sterilization during use.

14. Precautions

14.1. Precautions for use

1) The printer meets safety requirements when used for the purposes specified in this Operational Document. When preparing for use and using the printer, safety precautions must be observed.

2) Safety precautions:

When operating the printer, the following must be ensured:

- periodic monitoring of it;
- timely implementation of preventive maintenance and repair regulations;
- compliance with the requirements for safe placement.

3) Repairs to printer equipment must be carried out by the manufacturer (Authorized Representative) or an organization authorized by it. Independent troubleshooting (except for those regulated in section 19) and performing repair and adjustment work is not allowed.

In case of violation of the rules and regulations of operation, the requirements of safety measures established in this Operational Document, the manufacturer (Authorized Representative) does not bear any responsibility for the quality and technical condition of the product, as well as for any consequences that occurred during preparation for operation and/or during the operation of the printer, including those resulting in damage to human health and life, damage to the environment and human habitat.

Changes and modifications to the printer by the user that are not approved by the manufacturer (Authorized Representative) are prohibited.

4) In order to use the printer correctly and prevent danger or damage, be sure to execute the following instructions:

- Do not connect several devices to one multisectional socket at the same time. This may cause overheating and fire.
- If the power plug is wet or dirty, wipe it off before plugging it into the outlet.
- If the socket is worn out or the plug does not fit into the socket, do not insert it. Use only standard multisectional sockets.
- If you notice any strange smoke, odor, or noise coming from the printer, turn off the power. Then: Unplug the printer. After waiting for the smoke to stop, seek help for warranty repairs.
- When removing the plug from the socket, do not pull on the cord. This could damage it, which could result in a fire or printer failure.
- Do not connect or pull out the plug with wet hands. There may be a risk of electric shock.
- Do not bend or place heavy loads on the soft power cord as this may cause a fire.

Failure to follow these instructions may result in serious injury or death.

- Do not disassemble or repair the printer yourself.
- Do not touch the thermal printhead with your hand. This may cause burns or affect print quality.
- Use only recommended accessories, otherwise it may cause damage.
- Install the product in a stable place. A fall may cause damage and cause injury.
- Do not use the printer when it is broken, as this may result in fire or electric shock. Immediately turn off the power and contact your place of purchase or the manufacturer (Authorized Representative).

14.2. Warning about the impact of external factors on the functioning of a medical device.

Semiconductor devices can fail due to static electricity. Therefore, to protect it from discharges, first plug the plug into the socket, and only then turn on the printer's power switch, and not vice versa.

When turning off, perform these steps in reverse order. If the printer fails due to static electricity, turn off the product and contact the manufacturer (Authorized Representative).

Condensation: When the temperature changes significantly, for example when moving from a cold room to a heated room, the latter will create a large amount of moist air, and this may cause the printer to fog up inside. Under these conditions, the printer cannot operate normally, and if used for a long time in such conditions, it may fail. If the temperature changes significantly, the thermographic film may fog up inside the device cassette. This may ultimately damage the printer. If the printer is still used under these conditions, white spots may appear on the output film.

If fogging occurs, you must turn off the power to the printer. After the fogging goes away, you need to restart the printer.

- Do not allow water or foreign objects to enter the printer. If this happens, turn off the printer, dry it, and only then turn it on again.

14.3. Warning or precautions to be taken by the user when disposing of a medical device

The product does not pose an infectious or microbial hazard.

Printer materials do not have irritating, toxic or sensitizing effects.

When disposing of the product, apply the requirements of section 17.

14.4. Physical hazard of a medical device

The printer does not have accessible metal parts galvanically connected to the power supply. The design of the product components eliminates the possibility of cuts and injuries during their use and repair.

15. Transporting the product

1) During transportation, loading and unloading, packaging with products must be protected from any mechanical damage, vibrations, shocks, falls and direct moisture.

2) Products must be transported in the manufacturer's packaging, in an upright position, in no more than one tier. The placement and fastening of the printer in vehicles should ensure their stable position, exclude the possibility of the products moving and hitting each other, the walls of the transport container and vehicles; for this it is recommended to use inserts made of foam rubber materials.

3) Products transportation can be carried out by all types of transport in closed vehicles, over any distance, in accordance with the rules for the transportation of goods in force for this type of transport.

4) Temperature: -10 – 50 °C.

5) Relative Humidity Range: 10 – 80 %.

6) Atmosphere Pressure Range: 50 – 106 kPa

16. Product storage

- 1) Protect the product from vibrations, any mechanical impact, shocks and falls.
- 2) The product should be stored on stable surfaces, away from sources of moisture and/or heat. The distance from the packaged printer placed on a rack to any object (including the walls and floor of the storage) must be at least 0.1 m. The distance from the packaged product to the nearest heating device must be at least 1.0 m.
- 3) Stacking of products or placing extraneous cargo on them is prohibited.
- 4) Keep the device away from direct sunlight.
- 5) Prepare the storage room in such a way that it is free of dust, vapors of acids and alkalis, aggressive gases and other harmful impurities that cause corrosion.
- 6) Temperature: -10 – 50 °C.
- 7) Relative Humidity Range: 10 – 80 %.
- 8) Atmosphere Pressure Range: 50 – 106 kPa

17. Product disposal

- 1) The packaging is disposed of in household waste collection areas.
- 2) Before disposing of the device, render it inoperative so that it cannot become a source of danger. To do this, disconnect the device from the AC power supply.

3) The device or its parts must not be disposed of with other medical waste.

Note - in order to prevent environmental pollution, all waste generated during the disposal of products is subject to mandatory collection and subsequent disposal, in the prescribed manner and in accordance with the current requirements and standards of local regulatory documentation, as epidemiologically safe waste.

Disposal of products used in the infectious diseases department of a medical facility, as well as parts of the product that have contact with infected patients, is carried out in the prescribed manner and in accordance with the current requirements and standards of local regulatory documentation, as epidemiologically hazardous waste. Disposal of products should be carried out as disposal of class A waste according to SanPiN 2.1.3684-21 (epidemiologically safe waste, similar in composition to solid waste).

4) For more information about centralized recycling dates and locations, please contact your local government or waste disposal company.

5) Since regulations for disposal may differ in different countries, if necessary, contact the authorized representative of the manufacturer in the Russian Federation Limited Liability Company "Salute-Service" (LLC "Salute-Service"), 432071, Ulyanovsk region, Ulyanovsk, Federation st., 63. E-mail salute73@inbox.ru (Tel.: +7 (8422) 58-16-00). Salute-Service LLC carries out, if necessary, separate collection of the old printer, aimed at preventing potential damage that may be caused to the environment and health, for reuse and recycling of used products.

18. Warranty

18.1. The manufacturer (Authorized Representative) guarantees the normal operation of the printer for 1 year from the date of sale through a retail network or receipt by the user (customer) subject to the conditions of operation, transportation and storage.

Guaranteed shelf life is 1 year from the date of manufacture.

Service life – 15 years

18.2. During the warranty period, the manufacturer (Authorized Representative) carries out warranty repairs of the printer or failed equipment.

After repair and/or warranty service, the printer will be returned to the owner in its original condition at the time of delivery.

18.3. The components used in the design have a residual service life at the time of installation in the printer that is not less than its warranty period.

The warranty period for components is set equal to the warranty period for the printer and expires simultaneously with the expiration of its service life.

18.4. Complaints are made during the warranty period, provided that the consumer complies with the requirements of the printer's Operational Document and draws up a complaint report.

18.5. Warranty obligations do not apply to:

- for scratches, normal wear and tear, intentional damage or damage due to careless handling or failure to follow the instruction manual, improper installation or use of inappropriate voltage;
- for damage due to chemical or electrochemical reaction, rust, corrosion or exposure to water; for damage due to exposure to an unusual (inappropriate) environment;
- for accidental damage caused by foreign objects or substances, as well as due to cleaning, for repairs performed by unauthorized service organizations or unauthorized service partners; or repair using non-original spare parts;
- on cables, fuses, clamps and electrical connectors;
- when attempting self-repair without permission from the manufacturer (Authorized Representative);
- in case of damage caused due to force majeure (typhoon, fire, flood, military action, etc.);
- if the printer is used by unauthorized persons or in conflict with these instructions.

18.6. Submitting requests to the manufacturer's Service Center (Authorized Representative) to eliminate defects in the printer is possible only with the simultaneous presentation of a correctly completed warranty card. At the same time, it must clearly indicate information about the printer, the defects of which are subject to elimination (name, serial number), as well as information about the sale of the printer (date of transfer to the buyer, name and address of the seller), certified by the signature and seal (stamp) of the seller, as well as the buyer's signature.

19. Service Manual

19.1. Before contacting the service center, it is recommended to study this section to determine the possible cause of the problem. A few simple checks and adjustments can eliminate it and return the printer to normal operation.

19.2. The contents and interval of maintenance, including cleaning and disinfection of the medical device, are given in cl. 12.1.

19.3. Possible malfunctions and methods for eliminating them are indicated in the Table 1.

Failure to do so may result in personal injury, damage to the printer, and void the warranty!

Performing any operation that is expressly prohibited in the Operational Document, as well as any settings or assembly actions not recommended or prohibited in this Operational Document, will void the warranty!

Table 1 – Possible malfunctions and solutions

Error message on the printer's LCD touch screen	Possible reason	Solution
SET FILM	Out of film	Insert film into cassette
SET TRAY 1	The upper cassette compartment is missing or not inserted all the way.	Press the upper cassette compartment until it stops.
		The detector sensor of the upper cassette compartment is not working correctly.
SET TRAY 2	The lower cassette compartment is missing or not inserted all the way.	Press the lower cassette compartment until it stops.
		The detector sensor of the lower cassette compartment is not working correctly.
CLOSE COVER	The back and/or top cover is not closed.	Close the back and/or top cover.
		The back and/or top cover detector sensor is not working correctly.
REMOVE FILM FRONT	The film on the input gets pinched before it reaches the thermal printhead.	Remove the cassette compartment and remove jammed films.
REMOVE FILM REAR	The film on the input gets pinched after passing through the thermal printhead.	Open the back cover and remove any jammed films.
REMOVE FILM SOME PLACES	The film is jammed somewhere else in the printer.	Open the top and back covers and remove any jammed films.
CT TROUBLE	There may be a system software error	Turn off the power and turn on the printer again. The error may be resolved. If the error still occurs, immediately turn off the power and contact the service center.
MECHA TROUBLE	There may be a system software error	Turn off the power and turn on the printer again. The error may be resolved. If the error still occurs, please immediately turn off the power and contact your dealer.
HEAD IN	The thermal printhead is	Wait until the temperature of the thermal printhead is normal.

Error message on the printer's LCD touch screen	Possible reason	Solution
COOLING PLEASE WAIT	overheating. The printer will automatically stop printing.	The error message will disappear. You can continue printing.
No picture on LCD screen none	The POWER button on the front panel of the printer is not pressed.	Press the switch to turn on the device.
The device does not turn on when you press the POWER button on the front of the printer.	The main power switch I/O on the back of the printer is not set to the "I" position.	Move the main power switch I/O on the back of the printer to the "I" position.
The product does not turn on	Violation of recommended power supply parameters	Replace the fuse in the device's power supply.

19.4. To replace the fuse, remove the fuse holder from the slot (see Figure 33) in the device's power supply (see Figure 32).


Figure 32 – Location of the device's power supply on the back side


Figure 33 – Procedure for replacing the fuse in the device's power supply

Note – the device's power supply contains 2 fuses.

19.5. If the above methods for solving problems that have arisen do not help, contact the manufacturer's service department (authorized representative)!

20. List of international standards to which a medical device meets

Designation	Name
1	2
EN 1041:1998	Information provided by the manufacturer along with medical devices
IEC 62304-2022	Medical devices. Software. Life cycle processes
IEC 61326-1:2012	Electrical equipment for measurement, control and laboratory use. Electromagnetic compatibility requirements. Part 1. General requirements
IEC 61010-1-2014	Safety of electrical instrumentation and laboratory equipment. Part 1. General requirements
ISO 13485-2017	Medical devices. Quality management systems. Requirements for regulatory purposes
EN 55011:2009+A1:2010 (CISPR 11)	Industrial, scientific and medical (ISM) radio frequency equipment. Electromagnetic interference characteristics, limits and measurement methods
ISO 15223-1:2016	Medical devices. Symbols used for marking on medical devices, labels and accompanying documentation. Part 1. Basic Requirements
IEC 62366-1:2020	Medical devices. Part 1. Designing medical devices with respect to usability
ISO 14971-2021	Medical devices. Application of risk management to medical devices

21. Scope of delivery

Versions:	Name of components	Quantity, pcs
Dry Imager Printer HQ-430DY	Dry Imager Printer HQ-430DY	1
	Power cable	1
	Network cable	1
	Plates	2
	Rubber rollers	4
	Cloth for printer's thermal printhead	1
	CD with instructions	1
	Instructions for use	1
	Data sheet	1
Dry Imager Printer HQ-450DY	Dry Imager Printer HQ-450DY	1
	Power cable	1
	Network cable	1
	Plates	2
	Rubber rollers	4
	Cloth for printer's thermal printhead	1
	CD with instructions	1
	Instructions for use	1
	Data sheet	1
Dry Imager Printer HQ-460DY	Dry Imager Printer HQ-460DY	1
	Power cable	1
	Network cable	1

	Plates	2
	Rubber rollers	4
	Cloth for printer's thermal printhead	1
	CD with instructions	1
	Instructions for use	1
	Data sheet	1
Dry Imager Printer HQ-762DY	Dry Imager Printer HQ-762DY	1
	Power cable	1
	Network cable	1
	Plates	2
	Rubber rollers	4
	Cloth for printer's thermal printhead	1
	CD with instructions	1
	Instructions for use	1
	Data sheet	1
Dry Imager Printer HQ-760DY	Dry Imager Printer HQ-760DY	1
	Power cable	1
	Network cable	1
	Plates	2
	Rubber rollers	4
	Cloth for printer's thermal printhead	1
	CD with instructions	1
	Instructions for use	1
	Data sheet	1

ACCEPTANCE AND SALE CERTIFICATE

Printer _____

_____ modification

Manufacturing date _____

Quality Control Department stamp (acceptance stamp)

Price

Sold _____ Sale date _____
name of trade enterprise

WARRANTY CERTIFICATE

WARRANTY CERTIFICATE STUB NO. _____ Seized « ____ » _____ Executor _____ Full Name	Warranty Certificate No. _____ name of the manufacturer's authorized representative and his address WARRANTY CERTIFICATE No. _____ Printer _____ version and serial number Sold by store _____ Store name and number and its address Sale date _____ Store stamp _____ _____ Seller's personal signature Completed work _____
	Executor _____ Owner _____ surname, first name, patronymic signature name of the company that performed the repairs, Stamp . and its address position and signature of the repair company manager

EMC INFORMATION

Emissions Guide for EMC

Electromagnetic compatibility test	Conformity	Electromagnetic Environment - Guidance
Radio interference according to CISPR 11 (EN 55011:2009 +A1:2010)	Group 1	The product uses RF energy only for its internal operation. Thus, the RF emissions are very low and are not likely to interfere with nearby electronic devices.
Radio interference according to CISPR 11 (EN 55011:2009 +A1:2010)	Class B	Suitable for use in residential areas and premises directly connected to low-voltage electrical distribution networks used for domestic purposes.
Harmonic current components according to IEC 61000-3-2	Class A	Suitable for use in buildings of any public purpose, with the exception of domestic houses and buildings directly connected to the public low-voltage power distribution network supplying residential buildings.
Voltage fluctuations and flicker according to IEC 61000-3-3	Corresponds	

Electromagnetic immunity

Immunity test	Test level according to IEC 61326	Compliance level	Electromagnetic Environment - Leadership
Electrostatic discharge (ESD) according to IEC 61000-4-2	±4 kV – contact discharge ±8 kV – air discharge	±4 kV – contact discharge ±8 kV – air discharge	The floor in the room is made of wood, concrete or ceramic tiles. For floors covered with synthetic material, the relative air humidity is at least 30%.
Nanosecond burst noise – as required by IEC 61000-4-4	± 0.5 kV for I/O lines ± 1 kV for power lines	± 0.5 kV for I/O lines ± 1 kV for power lines	Mains power quality is consistent with typical conditions in a commercial or hospital environment
High energy microsecond burst noise – as required by IEC 61000-4-5	± 0.5 kV when applying interference according to the wire-to-wire circuit (DM) interference ± 1 kV when applying interference according to the wire-to-ground circuit (CM)	± 0.5 kV when applying interference according to the wire-to-wire circuit (DM) interference ± 1 kV when applying interference according to the wire-to-ground circuit (CM)	Mains power quality should be consistent with that of a typical commercial or hospital environment.
Voltage drops, short interruptions and dynamic changes in power supply voltage are as required by IEC 61000-4-11	0% UT during half cycle 0% UT during 1 cycle 70% UT during 25/30 ¹ cycles 0% UT during 250/300 ² cycles	0% during half cycle 0% during 1 cycle 70% during 25/30 ¹ cycles 0% during 250/300 ² cycles	Mains power quality is consistent with typical conditions in a commercial or hospital environment. If the user of the product needs to ensure continuous operation in conditions of possible interruptions in the mains voltage, it is recommended that the product be powered from an uninterruptible power supply or battery.

Immunity to power frequency magnetic fields according to IEC 61000-4-8	3 A/m, 50 and 60 Hz	3 A/m, 50 and 60 Hz	Power frequency magnetic field levels should be consistent with typical commercial or hospital environments
Note: Note: UT is the voltage level of the AC electrical network before the test stimulus is applied. ¹ 25/30 cycles means 25 cycles for testing at 50Hz and 30 cycles for testing at 60Hz. ² 250/300 cycles means 250 cycles for testing at 50Hz and 300 cycles for testing at 60Hz.			

Electromagnetic immunity

Immunity test	Test level according to IEC 61326	Compliance level	Electromagnetic Environment - Leadership
			Portable and mobile equipment should be used further away from machine parts, including cables, taking into account the recommended distance calculated from the corresponding transmission frequency equation.
Conducted interference induced by radio frequency electromagnetic fields is as required by IEC 61000-4-6	3V rms voltage in the band of 150 kHz to 80 MHz	3V	$d = 1.2\sqrt{P}$
RF Electromagnetic Field is according to the requirements of IEC 61000-4-3	3 V/m (80 MHz to 1 GHz band) 3 V/m (1.4 GHz to 2 GHz band) 1 V/m (2.0 GHz to 2.7 GHz band)	3 V/m (80 MHz to 1 GHz band) 3 V/m (1.4 GHz to 2 GHz band)	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz
			$d = 2.3\sqrt{P}$ 800 MHz - 2.7 GHz where P is the rated maximum output power of the transmitter in watts (W) as specified by the manufacturer; d is the recommended spatial separation in meters (m). The field strength of radio waves from fixed radio transmitters based on observations of the electromagnetic environment ¹ should be lower than the compliance level in each frequency band ² Interference may occur in the vicinity of equipment marked with: 

Note. At frequencies of 80 MHz and 800 MHz, a higher field strength is used.

Note. This guidance may not apply in all cases. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.

¹ The field strength of radio waves from stationary radio transmitters, such as base stations of radiotelephone networks (cellular or wireless) and land mobile radio stations, amateur radio stations, AM and FM radio broadcast transmitters, television transmitters, cannot be determined by calculation with sufficient accuracy. To assess the electromagnetic environment associated with fixed RF transmitters, field strength measurements must be made. If the measured values at the product location exceed the

applicable compliance levels above, the product should be monitored to verify that it is functioning properly. If a deviation from normal operation is detected, additional measures may need to be taken, such as reorienting or repositioning the product.

² Outside the band of 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

Electromagnetic immunity

Recommended spatial separation distances between portable and mobile RF communications equipment and the printer

Rated maximum transmitter output power (W)	Spatial separation depending on transmitter frequency, m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.7	3.8	7.3
100	11.7	12	23

When determining the recommended values of spatial separation d in meters (m) for transmitters with a rated maximum output power not listed in the table, substitute the rated maximum output power P in watts (W) specified in the transmitter manufacturer's documentation into the given expression corresponding to the frequency of the transmitter.

Note. At frequencies of 80 MHz and 800 MHz, the spatial separation for the larger field strength is applied.

Note. This guidance may not apply in all cases. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.

FOR NOTES

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

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