Date: 25th of December, 2023

M.D.L. SRL Approved by
Marcello Dell'Oca Via Tavani, 1/A (signature)
23014 DELEBIO (SO) Italia CEO
C.F./P.IVA 00656810140 M.D.L. S.r.l.

**INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE**

Biopsy Gun PLURI-GUN

1. Name of medical product

Biopsy Gun PLURI-GUN, comprising:

1. Biopsy Gun PLURI-GUN, 1 pce
2. Storage and carrying case, 1 pce
3. Instructions for use, 1 pce

2. Information about the developer, manufacturer of the medical device and an authorized representative of the manufacturer

Developer

M.D.L S.r.l.

Registered address: Via Tavani, 1/a – IT 23014 DELEBIO (SO), Italy

Tel.: +39 0342/682130

Fax: +39 0342/691316

E-mail: info@mdlsrl.com

Website: www.mdlsrl.com

Manufacturer

M.D.L S.r.l.

Registered address: Via Tavani, 1/a – IT 23014 DELEBIO (SO), Italy

Tel.: +39 0342/682130

Fax: +39 0342/691316

E-mail: info@mdlsrl.com

Website: www.mdlsrl.com

Manufacturing site

M.D.L S.r.l.

Via Tavani, 1/a – IT 23014 DELEBIO (SO), Italy

Authorized representative of the manufacturer:

Aktimed Plus Limited Liability Company (Aktimed Plus LLC)

Registered address: Lunacharskogo Prospect, 43, Block 1, Bld. A, Rooms 21-32, St. Petersburg, 194291, Russia

Tel.: +7 (812) 703-50-71

E-mail: info@aktimed-plus.ru

3. Intended use and intended users of the medical device

Intended use: The product is intended for the automated cutting biopsy of soft tissues from the anatomical structures of the patient's body for the purpose of pathological examination/testing of the biopsy material.

Intended users of the medical device: Dedicated and certified medical personnel who meet the qualification requirements applied to medical staff with higher professional (medical) education (medical doctors) in the relevant area of training (specialty), and who have professional knowledge, skills and abilities to perform cutting biopsy of soft tissues.

4. Indications, contraindications for use, side effects, precautions

Indications:

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE

Biopsy Gun PLURI-GUN

Performing the cutting biopsy of soft tissues from anatomical structures of the body, such as the liver, kidneys, prostate, breast, spleen, lymph nodes and various soft tissue tumours, for the purpose of tissue sample examination (tissue sample obtained through a biopsy procedure) in case of a suspected malignant tumour in any of the above tissues and/or suspected inflammatory process of unknown origin.

Caution!

- Soft tissue biopsy should be performed under ultrasound, X-ray or CT guidance to achieve the most precise placement of the biopsy needle.
- A MRI-guided biopsy must not be performed.

Contraindications to use:

Absolute contraindications (manipulations by means of the product are strictly contraindicated):

- Biopsy of bone, bone marrow, central nervous system and other hard tissues; and
- Purulent inflammatory process at the biopsy site (a tissue sample obtained via a biopsy procedure).

Relative contraindications (manipulations by means of the product are possible if its diagnostic benefit outweighs the risk, but complications are highly probable):

- Uncontrolled coagulation (severe bleeding disorders);
- Severe general condition of the patient (impaired consciousness to the point of a deep stupor; severe disorders of the respiratory or cardiovascular systems);
- Heart failure in the decompensation stage;
- Diabetes mellitus in the decompensation stage; and
- Acute infectious diseases (influenza, ARVI, COVID-19.)

Caution! The need to perform the procedure, its risks and contraindications are to be considered individually for each patient, and depend on his/her condition and presumptive diagnosis.

Note:

- *If the patient is taking anticoagulants (drugs to reduce blood clotting) or non-steroidal anti-inflammatory drugs, such as aspirin, a biopsy by means of the product can only be performed at least five days after their withdrawal.*
- *Typically, a cutting biopsy of soft tissues is performed under local anaesthesia. However, general anaesthesia must be used in some cases, for example, if the patient is unable to control his/her movements or experiences severe anxiety. In such cases, no eating or drinking is allowed within 8–12 hours before the procedure.*

Side effects:

- Pain in the biopsy sampling site;
- Minor swelling, bleeding, hematoma;
- Haemorrhage, infection, injury to adjacent tissues.

Caution! The risk of damaging the viscera during the procedure is minimal; however, it may increase because of patient's motions and very deep breaths.

After taking a biopsy sample (tissue sample obtained via a biopsy procedure), rest and immobilization are recommended to the patient. Until the complete healing of the tissues, the patient should avoid drinking alcoholic beverages and undertaking any physical activity; also, visiting a bathhouse, sauna, swimming pool, or heating the biopsy site are not recommended. Furthermore, the patient should not remove the bandages or wet the biopsy needle insertion site without the doctor's permission within the first 1-2 days after the procedure.

Caution! If the patient experiences a severe lightning or throbbing pain, severe swelling, redness of the biopsy site, bleeding, or deterioration of the general health condition, he/she should consult a doctor.

Precautions

- Use the product only in accordance with these Instructions for Use.
- Use the product only after completing proper cleaning, disinfection and sterilization procedures.
- Users must employ the cleaning, disinfection and sterilization methods recommended by the manufacturer. Follow the manufacturer's recommended disinfection and sterilization parameters for the product; otherwise, a failure to comply with the recommendations may result in the reduced service life, damage and/or failure of the product.
- Before each use, always check the product for (internal and external) mechanical damage and proper operation.
- Do not use the product in case of damage to any of its parts or its malfunction.
- If the product is faulty, see Section 11 "Troubleshooting".
- Follow the warning signs and take precautions when handling the product.
- Never try to disassemble or repair the product yourself, as this may lead to its breakdown and mechanical injury to the user.
- When you bring the product from a cold place to a warm room, keep it at the room temperature for at least 4 hours; otherwise, moisture settling on the product and its parts may damage it.
- Use the product only for its intended purpose, that is, to perform an automated soft tissue biopsy procedure on the patient's anatomical structure for tissue sample examination.
- Do not make any modifications in the product. Any modifications to the design may result in the product failure, malfunction, or injury of the user/patient.
- Protect the product from falls and hard impacts.
- The product can be used only by dedicated, certified medical personnel meeting the qualification requirements that apply to the medical staff with higher professional (medical) education (medical doctors) in the relevant area of training (specialty), who have professional knowledge, skills and abilities to perform a cutting biopsy of soft tissues and have studied these instructions for use in detail. The doctor is responsible for the correctness of the biopsy procedure and the techniques used with the product.
- It is recommended to use the product in conjunction with a custom-designed soft tissue biopsy needle for a reusable automatic system, manufactured by M.D.L. S.r.l. (Marketing Authorisation No. FSZ 2012/13006 dated 16 December 2016), and soft tissue biopsy needles for a reusable automatic system manufactured by M.D.L., having the following characteristics: diameter of 12G to 20G and length of 80 mm to 300 mm. The biopsy needle is not included in the product delivery package.
- Medical personnel should handle the product only while using personal protective equipment for hands (medical gloves).
- A soft tissue biopsy should be performed under ultrasound, X-ray or CT guidance to achieve the most accurate position of the biopsy needle.

Caution! MRI-guided biopsy is prohibited.

- The biopsy needle must be installed in the product with a protective cap (sheath) on, so as to avoid injury to the user. The protective cap (sheath) should be removed only immediately before insertion of the biopsy needle to take a biopsy sample (tissue sample obtained through a biopsy procedure).
- In order to prevent damage, never get a proof shot or test the product with a biopsy needle installed, as this may result in an injury and damage to the biopsy needle. Always perform empty shots (simulating the biopsy process in the air outside of the tissue) without a biopsy needle.

- When the safety of the product is released, it is ready for use. Take care not to press the release button inadvertently and trigger the device.
- Do not remove the biopsy needle from the product in order to extract a biopsy sample (tissue sample obtained through a biopsy procedure).
- To prevent the potential cross-contamination, do not use the same biopsy needle on multiple patients. A new biopsy needle should be used for each patient. Using the same biopsy needle on multiple patients can lead to infection.
- Re-cleaning, disinfection and sterilization of a biopsy needle after its use is strictly forbidden. The waste needle must be disposed of in accordance with the manufacturer's recommendations set out in the accompanying documentation and the regulations in force in the country of use at the time of disposal.
- Never perform a biopsy using the device if the patient's condition disqualifies him/her from the procedure.
- The risk of damaging the viscera during the procedure may increase because of patient's motions and very deep breaths.
- Make sure that the biopsy device with a new biopsy needle has been disinfected and sterilized for each patient.

5. Operating principle

The Biopsy Gun PLURI-GUN is a spring-loaded, reusable device for performance of a cutting biopsy of soft tissues and viscera, such as the liver, kidneys, mammary glands, lymph nodes, spleen and prostate glands. It is recommended to perform this type of biopsy under X-ray, ultrasound or CT guidance in order to achieve the most precise placement possible of the biopsy needle. The automatic biopsy device has a special-purpose spring mechanism, which, after triggering, introduces a biopsy needle into the abnormal mass and cuts off a tissue. The penetration depth can be selected: 15 mm or 22 mm.

6. Product description

The Biopsy Gun PLURI-GUN is a spring-loaded reusable device for performance of a cutting biopsy of soft tissues and viscera, such as the liver, kidneys, mammary glands, lymph nodes, spleen and prostate gland, which can be used under X-ray, ultrasound or CT guidance.

The surface of the PLURI-GUN device body has the following marking (Figure 1):

- Manufacturer's logo
- Trade name of the product
- "Open " ("Open here");
- Biopsy needle penetration (cutting) depth: "22mm/15mm"
- Version "4.0"

Figure 1 – Biopsy Gun PLURI-GUN

The Biopsy Gun PLURI-GUN is recommended to be used in conjunction with a custom-designed soft tissue biopsy needle for a reusable automatic system manufactured by M.D.L. S.r.l. (Marketing Authorisation No. FSZ 2012/13006 dated 16 December 2016), and soft tissue biopsy needles for a reusable automatic system manufactured by M.D.L., having the following characteristics: diameter of 12G to 20G and length of 80 mm to 300 mm. The biopsy needle is not included in the product delivery package. The biopsy needle diameter and length should be chosen depending on requirements of the procedure in each specific case.


The Biopsy Gun PLURI-GUN is provided with indicators of the charge (cocking) of the biopsy needle cannula and the stylet (see Table 1)

Table 1. Indicator of charge (cocking) of the cannula and stylet

	The device is idle, both indicators are colourless.
	- Charging (cocking) the biopsy needle cannula via the first movement (pulling back) of the plunger — the charge (cocking) indicator of the biopsy needle cannula turns white.
	- Charging (cocking) the biopsy needle stylet via the second movement (pulling back) of the plunger — the charge (cocking) indicator of the biopsy needle stylet turns white. - The product is ready for operation — BOTH charge (cocking) indicators of the biopsy needle cannula and stylet turn white.

The PLURI-GUN device provides adjustment of the shot (penetration) depth of the biopsy needle using a penetration (cutting) depth switch (Table 2):

Table 2 – Setting the biopsy needle shot depth

22 mm $\pm$ 0.5 mm		
15 mm $\pm$ 0.5 mm		

A manually switched lever safety protects the product against inadvertent activation (Table 3). The product cannot be activated, unless the safety is released.

Table 4 – Releasing the PLURI-GUN safety

	The Biopsy Gun PLURI-GUN is at safety. The product cannot be activated unless the safety is released.
---	---

Before activation, release the safety by moving the safety latch to the appropriate position.

After the shot, the safety returns to its original position.

7. Procedure of use


Figure 2 – How to use the Biopsy Gun PLURI-GUN

A. How to prepare the PLURI-GUN device

1. Before using the Biopsy Gun PLURI-GUN, make sure that there is no damage to the storage and carrying cases and no visible damage to the product's body.
2. Before the first use and after each intended use, the biopsy device must be cleaned, disinfected and sterilized by autoclaving (for cleaning, disinfection and sterilization of the PLURI-GUN device, see Section 12 of this document).
3. After completing the cleaning, disinfection and sterilization steps, sanitize the skin of your hands, then put on medical gloves and remove the *PLURI-GUN* device from its sterile packaging.
2. Before installing a biopsy needle into the device, checking its functionality (correct operation) is recommended:
 - Hold the device in the palm of your hand with your fingers on the white recess, the plunger. (Figure 2).
 - Twice firmly pull the charging (cocking) plunger of the cannula and stylet toward the palm of your hand, until both indicators on the product body turn white. The PLURI-GUN device has two indicators: at the first movement of the plunger, the colour of the distant indicator changes, which indicates the charge (cocking) of the cannula; at the second movement, the near indicator turns white, indicating that the stylet is charged (cocked). (Fig.2)
 - When both indicators turn white, the Biopsy Gun PLURI-GUN is charged and ready for use.
 - Then, check the functioning of the penetration (cutting) depth selection in the cocked position by turning the special switch to the appropriate value: 15 mm or 22 mm.
 - Check the functionality of the product by switching the lever safety from a horizontal to a vertical position (hold the Biopsy Gun PLURI-GUN with the white recess down).
 - Press the release button to activate the Biopsy Gun PLURI-GUN. Once released, the safety lever will automatically return to the safe position.

B. Preparing the biopsy needle

After completing the functional testing of the Biopsy Gun PLURI-GUN, prepare the biopsy needles following the rules of antiseptics.

- It is recommended to use the product in conjunction with a custom-designed soft tissue biopsy needle for a reusable automatic system, manufactured by M.D.L. S.r.l. (Marketing Authorisation No. FSZ 2012/13006 dated 16 December 2016), and soft tissue biopsy needles for a reusable automatic system manufactured by M.D.L., having the following characteristics: diameter of 12G to 20G and length of 80 mm to 300 mm. The biopsy needle is not included in the product delivery package. The biopsy needle diameter and length should be custom-selected depending on requirements of the procedure.
- The needle has the depth mark and is packaged together with a spacer. The doctor should select a needle depending on the biopsy procedure parameters. During the entire procedure, the doctor is to hold the needle only by the special fixing bar (spacer).
- Before use, visually inspect the biopsy needle for signs of damage and check that the disposable biopsy needle package is not opened and/or damaged. Remove the needle from the package while wearing medical gloves. Inspect the needle components for damage. Each biopsy needle is placed in a protective sheath.

Caution! DO NOT USE if needle components are damaged.

C. Biopsy procedure

1. A biopsy needle can be used in two ways:

- I. Insert the selected biopsy needle into the fully charged PLURI-GUN device as follows:
Open the cover of the cocked device. Insert the biopsy needle into the corresponding "seats" of the device so that the tip of the biopsy needle extends through the slot in the front wall of the device (needle outlet). Pull the spacer to remove it and close the cover. Introduce the biopsy needle through an incision in the skin under ultrasound, x-ray or CT guidance, until the tip of the needle is as close to the biopsy site as possible.
 - II. Introduce the biopsy needle with the spacer through an incision in the skin under ultrasound, X-ray or CT guidance, until the tip of the needle is close to the biopsy site. Once the position has been ascertained, insert the needle into the cocked PLURI-GUN device. Insert the needle into the corresponding "seats" of the device so that the tip of the biopsy needle extends into the slot in the front wall of the device (needle outlet). Partially close the cover, while keeping the needle in place. Remove the spacer from the needle, holding the product firmly. Close the cover.
2. The depth of the biopsy needle shot (penetration, cutting) is selected and set by turning the special switch to the appropriate value: 15 mm or 22 mm (Table 2). The introduction depth of the biopsy needle can be set only when the device is fully cocked.
 3. While maintaining product position and needle orientation, move the lever safety (Table 3) in order to prepare the device for precision biopsy, and press the release button to release the cannula and stylet.
 4. Remove the needle from the patient's body.
 5. A biopsy sample (tissue sample obtained through a biopsy procedure) can be extracted from the biopsy needle by moving (pulling back) the plunger of the device. The cannula will open the "window" of the biopsy needle stylet with the obtained biopsy sample. The biopsy device design allows extracting the obtained tissue sample without removing the biopsy needle from the device.
 6. If another biopsy sample (tissue sample obtained through a biopsy procedure) needs to be taken from the patient after extraction of the tissue sample, the device can be prepared for use only after re-cocking by a single pull of the plunger back. This movement activates the product for the second biopsy.
Note: After one biopsy, if additional samples need to be taken, always assess the integrity of the needle and inspect it for damage before starting the next biopsy.
- Caution!** Do not use the needle if any defects are noticed.
7. If additional biopsies are not needed, put on the protective cap (sheath), and remove the needle by opening the product cover and reinserting the spacer.

8. After removing the needle, dispose of waste needles in an appropriate biohazard container according to special hospital procedures.

Note: Clean, disinfect and sterilize the PLURI-GUN device before each use as instructed in Section 12 of this document.

Table 4 – Specifications of the Biopsy Gun PLURI-GUN

Specification		Value
Dimensions	Length, mm ($\pm 10\%$)	122
	Width, mm ($\pm 10\%$)	35
	Height, mm ($\pm 10\%$)	35
Weight, kg		≤ 0.3
Shot depth, mm (± 0.5 mm)		15
		22

8. Scope of delivery

Biopsy Gun PLURI-GUN, comprising:


1. Biopsy Gun PLURI-GUN, 1 pce;
2. Storage case, 1 pce;
3. Instructions for Use*, 1 pce.

*Available in hard copy or in soft copy documentation posted online.

9. Packaging

Packaging material and design do not adversely affect the content and ensure integrity of devices subject to their proper storage, the minimum risk of contamination of the content when opening the packaging, and protection of the content under normal transportation and storage conditions.

The Biopsy Gun PLURI-GUN is packed in a labelled plastic case. Instructions for Use for the device are packed in an individual clear LDPE bag. The case containing the device and the Instructions for Use are packaged in an additional clear LDPE bag (Figure 3).

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE
Biopsy Gun PLURI-GUN

Figure 3 – Packaging of the Biopsy Gun PLURI-GUN

Biopsy Gun PLURI-GUNS are packed in shipping cartons in quantities requested by the buyer. The size of a shipping carton is as ordered by the buyer.

Explanation of symbols of package labelling:

Table 5 – Designations and symbols

Symbol	Description	Symbol	Description
	Manufacturer		Manufacture date
	Part number		Keep dry
	Serial number		Do not use if the package is damaged and read <i>Instructions for Use</i>
	Handle with care		Read <i>Instructions for Use</i> or to electronic <i>Instructions for Use</i>
	Not MRI-compatible		
	Upper temperature limit		Protect against sunlight

10. Storage, shipping and operation conditions

Storage, shipping and operation conditions

- temperature: below 50 °C
- relative humidity: no requirements
- atmospheric pressure: no requirements

11. Troubleshooting

This section provides a list of possible malfunctions of the Biopsy Gun PLURI-GUN (Table 6) during its preparation for use and intended use, and possible causes of their occurrence and recommendations for actions in the event of malfunctions. If all troubleshooting steps are taken but the trouble remains, please contact the authorized representative in the Russian Federation.

Table — 6 List of possible malfunctions of the Biopsy Gun PLURI-GUN during preparation for use and intended use

Malfunction	Recommended corrective actions
The device does not charge (does not cock) or charges (cocks) with difficulty	- Properly position the biopsy device in your hand. - Pull the plunger back to the end of its stroke for charging (cocking) of the biopsy needle cannula and stylet. - Make sure that you do not press the release button while charging (cocking) the device. - Clean, disinfect and sterilize the Biopsy Gun PLURI-GUN (see Section 19 of this document).

A biopsy needle cannot be installed in the device	- Check if the PLURI-GUN device is fully charged (cocked) or completely uncharged (uncocked). - Check that the "seats" (grooves) for installation of the biopsy needle are not distorted. - Make sure that a fully compatible biopsy needle is inserted into the device.
The device cover cannot be closed	- Check if the cover can be closed without the biopsy needle installed. - Make sure that the biopsy needle is installed correctly. - Clean, disinfect and sterilize the PLURI-GUN device (see Section 12 of this document). - Make sure that a fully compatible biopsy needle is inserted into the device.
Biopsy procedure (shot of the biopsy needle) is constrained	- Make sure the device is fully charged (cocked). - Make sure that the safety of the biopsy device is released. - Make sure that you press the release button firmly. - Clean, disinfect and sterilize the device (see Section 12 of this document).
The penetration (cutting) depth cannot be switched	- Make sure the device is fully charged (cocked). - Clean, disinfect and sterilize the biopsy device (see Section 12 of this document).

12. Cleaning, disinfection and sterilization of the device

The Biopsy Gun PLURI-GUN is supplied non-sterile. Before first use and after each intended use, clean, disinfect and sterilize the device by autoclaving.

1. Preparation of the device: Prepare the device for cleaning by opening the top cover and turning the penetration depth switch to 22 mm for optimal access to the internal parts of the device. The introduction depth of the biopsy needle can be set only when the device is fully cocked.
2. Disinfection and pre-sterilization cleaning: Prepare aqueous solution of 0.3% enzyme concentrate. Wash the device with the above solution and leave to soak for 10 minutes. Clean all visible areas of the device with a soft brush so that no visible traces of contaminants remain.
Caution!
 - DO NOT USE 2 % glutaraldehyde or other antibacterial solutions;
 - DO NOT USE highly alkaline or highly acidic cleaners;
 - DO NOT USE concentrated chlorhexidine gluconate or any alcohol-containing agent to clean the device.
3. Washing and drying: Take the device from the solution and rinse with running water. Place the device on a dry, sterile cloth. After drying, perform a visual inspection to ensure that cleaning is effective and that all visually detectable contaminants have been removed. If visual inspection results are unsatisfactory and contaminants remain on the device, repeat the cleaning. Do not use the device if traces of contaminants cannot be cleaned due to deterioration of the device, such as corrosion, discoloration, metal fatigue, peeling of the outer coating, and cracking.
4. Sterilization: Prepare the device as specified in paragraph 1 of this section. Pack the biopsy device in a sterilization package, then place the biopsy device in an autoclave. Parameters for sterilization by autoclaving are provided in Table 7. The sterilization process takes 10 minutes. Allow the device to cool. After every sterilization procedure, the biopsy device should be inspected for integrity before its further use.

Table 7 – Sterilization parameters of the Biopsy Gun PLURI-GUN

Parameters	Values
Time	10 min
Temperature	134 °C ± 3 °C

Pressure	2 bar
----------	-------

Caution!

- DO NOT STERILIZE with radiation;
- USE personal protection means
- The maximum possible number of reprocessing cycles is not more than 600.

13. Disposal

Medical devices "Biopsy Gun PLURI-GUN" are classified as epidemiologically hazardous waste.

In this case, a medical organization only performs the collection, disinfection and temporary storage of the above waste. Disinfection must be carried out in compliance with sanitary rules and other regulations of a particular country in all medical institutions, regardless of their profile. Waste is collected in single-use bags or containers with colour markings, adopted for hazardous medical waste, in a room or in places suitable for temporary storage of medical waste before their subsequent transportation to a destruction or disposal site.

Transportation of medical waste (used implants or their components) from medical facilities and its disposal are carried out by specialized organizations licensed to treat epidemiologically hazardous medical waste.

Consumer packaging is classified as epidemiologically safe waste (non-hazardous waste not contacting body fluids of patients or infectious patients), which is close by its composition to municipal solid wastes (MSW). However, it is recommended to dispose of them separately as industrial waste through specialized organizations licensed for this type of activity.

Epidemiologically safe waste is collected in reusable containers or disposable bags (excluding bags of the colour adopted for hazardous waste). The procedure for washing and disinfection of reusable containers is determined in accordance with the waste management scheme in each specific organization. The procedure for collection of epidemiologically safe waste is determined in accordance with the waste management scheme in each specific organization.

Transportation and disposal are performed considering requirements adopted for a given region, in accordance with the sanitary legislation of a particular country for the maintenance of populated areas and treatment of production and consumption waste.

14. Warranty period of operation and storage

Table 12 – Warranty period of operation and storage

Parameters	Biopsy Gun PLURI-GUN
Possibility of the device reprocessing	≤600 cycles
Number of activations of charging and releasing mechanisms	≤12000
Guaranteed shelf life	>1 year
Warranty period	1 year

15. List of international standards that the medical device complies with

DIRECTIVE 93/42/EEC, MDR 2017/745, EN ISO 13485/A11, EN ISO 14971, EN ISO 15223-1, ASTM F2503-20, EN 62366-1 + A2, ISO/TR 20416, EN ISO 11737-1, EN ISO

INSTRUCTION FOR USE FOR THE MEDICAL DEVICE
Biopsy Gun PLURI-GUN

11737-2, AAMI TIR 12, AAMI TIR 30, EN ISO 15883-1, ISO TS 15883-5, EN ISO 17664, MEDDEV 2.7/1, Rev. 4, MEDDEV 2.12/1, Rev. 8.

Note: The standards are used in the latest edition: the organization maintains and updates the list of reference standards as a record of its own management system in accordance with ISO 13485.

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