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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of HUQIU IMAGING
(SUZHOU) CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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
“Thermographic medical-grade film “Dry Film HQ-KX410 A/F” in design versions

Instruction for use

02.04.2018

1 Name and purpose

1.1 Purpose of medical device

«Thermographic medical-grade film «Dry Film HQ-KX410 A/F» » in design versions (hereinafter- device, thermographic film, or film) is intended for use in medical imaging systems in printers using the principle of imaging with the help of temperature change.

Thermographic medical-grade film is intended for use in medical imaging systems in printers using the principle of imaging with the help of temperature change.

Consumable single-use material.

«Thermographic medical-grade film «DRY FILM HQ-KX410 A/F» in design versions:

- 1 Thermographic medical-grade film «Dry Film HQ-KX410 A/F», size 8x10 inches (20x25 cm) – 1 cartridge (100 sheets);
- 2 Thermographic medical-grade film «Dry Film HQ-KX410 A/F», size 10x12 inches (25x30 cm) – 1 cartridge (100 sheets);
- 3 Thermographic medical-grade film «Dry Film HQ-KX410 A/F», size 11x14 inches (28x35 cm) – 1 cartridge (100 sheets);
- 4 Thermographic medical-grade film «Dry Film HQ-KX410 A/F», size 14x17 inches (35x43 cm) – 1 cartridge (100 sheets).

After opening of an individual package (foil-coated package) of cartridge with film, it is recommended to load it at once in the imagesetter (printer)

Attention!

Safety measures and additional information on use of this device are given in operating manual for the printer.

1.2 Indications for use.

For use in medical imaging systems in printers using the principle of image acquisition by temperature change

1.3 Contraindications to use

There are no contraindications

1.4 Side effects

No side effects

1.5 Risks of using

Risks of using	Predictable the occurrence of the risk	Problem solution
The error in the image (hall)	Deformation (bending/crumpling) of the film when loaded into the image output device	Study the instructions for use before use
	Touching the printed side of the film	
	Placement of the print between the objects containing the plasticizer, contact of the films with the eraser	
	Erroneous use of films	
Error in the image with the ability to diagnose	Violation of storage conditions	Study the instructions for use before use
	Violation of the order of application of films	
	Contact of films with alcohol-containing liquids or other volatile organic solvents	
	Violation of the order of application of films	
	Erroneous use of films	Two-stage quality control during production: current control and quality control of manufactured products
	Flawed research	
The inability to use	Violation of the order of application of films	Study the instructions for use before use
	Incompatibility with the printer	

2 Manufacturer and authorized person

Manufacturer:	Huqiu Imaging (Suzhou) Co., Ltd, China
Manufacture address:	22 Building, No. 369 Lushan Road, Suzhou New District, 215129 China
Authorized representative in the RF:	Limited liability company "Salute-Service" ("Salute-Service" Ltd)
Legal address:	63, Federacii str., Ulyanovsk, 432071
Telephone:	8 (8422) 58-16-00; 8 (8422) 58-04-00
E-mail:	salute73@inbox.ru

3 Description of the device

Thermographic film possesses such characteristics as insensitivity to light, low level of grey colour ("fog"), high contrast, high resolution, high optical density. The film makes it possible to obtain high quality images due to optimized half-tones (displays to 4096 tone colours of grey) and has a special technology of increase of image resolution (resolution- 508 dpi). Smoothness of transmission of film half-tones provides for excellent diagnostic sharpness of reproduced image.

During printing on printers (or other medical imagesetters), oxide of activator is formed in the thermo emulsion layer under the influence of temperature, and silver is reduced creating the image.

The device can be used directly in daylight (natural illumination), the process of development is absent, there's no need in a dark room and/or in chemical reagents. All these characteristics correspond to high level of convenience for use.

The device is a consumable material for printing and is intended for use with printers and other medical imagesetters using similar consumable material.


Figure 1 – Thermographic medical-grade film «Dry Film HQ-KX410 A/F» in individual package

Figure 2 – Thermographic medical-grade film «Dry Film HQ-KX410 A/F»


Figure 3 – Appearance of the film itself, without cartridge

4 Description of materials

Thermographic film is made on 191-micron support of polyethylene terephthalate of blue colour, coated with thermal emulsion coating, has a protective layer and rounded angles (see figure 3).

Thermo emulsion layer is not sensitive to light, but reacts to high temperatures providing for delivery of the image on the film with method of direct thermographic printing with low “fog” level, high optical density and contrast of diagnostic imaging, along with its excellent stability.

Protective layer gives to the film stability to scratches, action of moisture and other adverse external effects.

Thermographic film is made of polyethylene terephthalate (CAS № 25038-59-9), manufactured by the company «Shanghai Jieli Insulation Materials Co., Ltd.» (China), coloured with blue colouring grade 16900-vms, manufactured by the company «Ampacet Corporation» (USA), coated with mixture, manufactured by the company «Huqiu Imaging (Suzhou) Co., Ltd» (China) of the following composition: silver bromide (CAS № 7785-23-1) and sodium thiosulphate (activator) (CAS № 7772-98-7) with secondary coating with protective mixture, manufactured by the company «Huqiu Imaging (Suzhou) Co., Ltd» (China) of the following composition: trifluoro-vinyl chloride copolymer with vinyl-butyl ether and monovinyl ether of ethylene glycol, dimethylamine hydrochloride (CAS № 506-59-2), 2,4-toluilen diisocyanate (CAS № 584-84-9), 2,6-toluilen diisocyanate (CAS № 91-08-7), hydroxyphenyl benzotriazole grade Tinuvin®479-DW, polymethyl siloxane (CAS № 8050-81-5).

The device is not sterile.

Any materials of human and/or animal origin, medicinal products are not used in device manufacturing.

5 Technical characteristics

Name of the characteristic	Dimensions of thermographic medical-grade film «Dry Film HQ-KX410 A/F»			
	8×10 inches (20×25 cm)	11×14 inches (28×35 cm)	14×17 inches (35×43 cm)	10×12 inches (25×30 cm)
Thickness of one sheet, µm	197	197	197	197
Density, g/m ²	260	264	264	264
Maximum print resolution (pixels value), µm	20	20	20	20
Transmission of half-tones, bit	14	14	14	14
Display of grey tone colours	4096	4096	4096	4096
Range of reproduced optical density, B:	≥2,9	≥2,9	≥2,9	≥2,9
Sensitivity to light	Not sensitive			
Number of sheets in a cartridge	100	100	100	100

6 The order of work with the film

ATTENTION! Never open the film with a knife (for example, a paper cutter). It can damage the film and it cannot be used for printing.

Before installation in the imagesetter take away the white sheet of paper (see figure 4), and then divide the sheets of film quickly thumbing them through.

Figure 4 – Cartridge with films and white sheet of paper which it is necessary to remove prior to printing

ATTENTION! Try to not to touch the printing side of the film. Be accurate to not to bend or crease the films.

ATTENTION! Preliminary tests (one-time test printing) should be carried out in an attempt to use this film in imagesetters for the purpose to achieve the desired results.

ATTENTION! After opening a cartridge with film it is recommended to load it into the imagesetter at once. During film loading its bending, deformation is not allowed.

Partly used films should be removed from the printer, closed up in the original individual package (foil-coated package) and stored in a cool dark place (in accordance with storage conditions). To avoid condensation, before unpacking the devices again and using them after storage, wait till they are warmed up to operating temperature.

For storage of partly used film it is recommended to use the original individual package (foil-coated package).

Never glue adhesive tape on the films. Also, do not place a plastic rubber on top of a print, and do not put a print between subjects containing plasticizer (table mat, etc.). Take care to ensure that any alcohol containing liquids or other volatile organic solvents could not get on the prints.

7 Operation, transportation and storage conditions of thermographic film.

Shelf life

Operation conditions of the device:

- temperature- from plus 18 to 24 °C;
- humidity: from 50 to 65%.

Transportation and storage conditions of the device:

- temperature- from plus 10 to 23 °C;
- relative humidity: from 30 to 65%.

The film shall be stored in cool, dry, non-dusty place, away from direct sunrays, acids and alkaline gases, in vertical position, to avoid adverse effect of external pressure.

ATTENTION! Store the film away from sources of heat, as well as acid and alkaline gases, such as hydrogen sulphide, ammonia, sulphur dioxide, formaldehyde, etc.

Terms of archivation correspond to the requirements of ANSI IT 9.11 and IT 9.19

Shelf life of the device: 24 months from the day of manufacture, with observance of storage conditions;

Print storage period: not less that 10 years, with observance of storage conditions.

8 List of EU harmonized standards

No.	Standard	Name of standard
1	EN ISO 13485:2016	Medical devices- Quality management systems- Requirements for regulatory purposes
2	EN ISO 14971:2012	Medical devices- Risk management for medical devices (ISO 14971:2007, Revised version 2007-10-01)
3	ISO 9236-1:2004	Sensitometry of screen/film systems for medical radiography- Part 1: Determination of sensitometric curve shape, speed and average gradient
4	ISO 4090:2001	Medical radiographic cassettes/ screens/films and hard-copy imaging films. Dimensions and specifications

9 The order of sterilization and disinfection

Thermographic film is a non-invasive consumable material for medical imagesetters (printing), it is not subject to sterilization and disinfection.

10 Maintenance and care

Thermographic film is not subject to maintenance and repair.

11 The order of disposal

Disposal must be carried out in accordance with the norms accepted in medical practice, as well as corresponding local, state and federal regulations, and instructions according to SanPiN 2.1.7.2790-10, class A.

12 Scope of supply

Scope of supply: 100 sheets of film in a cartridge (vacuum-sealed in a foil-coated package);

Number of packages in transport package: 5 cartridges 100 sheets each.

Weight and overall dimensions of transport package (length x width x height)

- 1) Thermographic medical-grade film «Dry Film HQ-KX410 A/F», size 8x10 inches (20x25 cm):
weight of transport package with 5 cartridges– 10,2 kg, transport package dimensions 470 × 240 × 510 mm;
- 2) Thermographic medical-grade film «Dry Film HQ-KX410 A/F», size 10x12 inches (25x30 cm):
weight of transport package with 5 cartridges – 13,9 kg, transport package dimensions 470 × 240 × 510 mm;
- 3) Thermographic medical-grade film «Dry Film HQ-KX410 A/F», size 11x14 inches (28x35 cm):
weight of transport package with 5 cartridges – 17,3 kg, transport package dimensions 470 × 240 × 510 mm;
- 4) Thermographic medical-grade film «Dry Film HQ-KX410 A/F», size 14x17 inches (35x43 cm):
weight of transport package with 5 cartridges – 24,9 kg, transport package dimensions 470 × 240 × 570 mm.

13 Warranty obligations

The manufacturer is not responsible for damages or accidents in the result of misuse, user's mistake or external emergency situations. On the matters of product quality please contact the authorized representative of the manufacturer.

14 Reclamation

In case if defective products are detected, the developer and manufacturer Huqiu Imaging (Suzhou) Co., Ltd, China, address: 22 Building, No. 369 Lushan Road, Suzhou New District, 215129 China guarantees the quality of products. If any defective products are detected, please contact the authorized representative in the Russian Federation.

The authorized representative of the developer and manufacturer in the Russian Federation- Limited liability company "Salute-Service" ("Salute-Service" Ltd), 63, Federacii str., Ulyanovsk, 432071, telephone: 8 (8422) 58-16-00; 8 (8422) 58-04-00, E-mail: salute73@inbox.ru – carries out reception of reclamations.

Damages of the device which occurred during operation and mishandling shall not be deemed defects.

Annex 1. List of national standards used by the manufacturer of medical devices

EN ISO 13485:2016 Medical devices – Quality management systems – Requirements for regulatory purposes

EN ISO 14971:2012 Medical devices – Application of risk management to medical devices

ISO 9236-1:2004 Photography – Sensitometry of screen/film systems for medical radiography – Part 1:
Determination of sensitometric curve shape, speed and average gradient

ISO 4090:2001 Photography. Medical radiographic cassettes/screen/films and hard-copy imaging films. Dimensions
and specifications

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