Instruction Manual

V 1.3

Tyromotion

tyromotion

Manufacturer and Distributor:

Tyromotion GmbH
Bahnhofgürtel 59
8020 Graz - AUSTRIA
Tel.: +43 316 908909
Fax.: +43 316 908909 – 90
E-Mail: office@Tyromotion GmbH.com
Internet: www.Tyromotion GmbH.com

Table of Content

1 Introduction.....	7
1.1 Safety	8
1.1.1 Intended Use.....	8
1.1.2 Indications.....	8
1.1.3 Contraindications	9
1.1.4 Side Effects.....	11
1.1.5 Environmental Conditions.....	11
1.1.6 Proper Use and Precautions.....	11
1.1.7 Responsibility of the Owner	13
1.1.8 Errors and Incompleteness.....	14
1.1.9 Property of Tyromotion GmbH.....	14
1.1.10 Guarantee & Exclusion of Liability.....	14
1.1.11 Service Life of medical device.....	16
1.1.12 Environmental Protection.....	16
1.1.13 Utilization	16
2 Technical Specifications	17
2.1 Tymo Therapy Plate.....	17
2.2 Markings on the Tymo Therapy Plate	19
3 Preparation/Basic Information	20
3.1 Operation, transportation, storage	20
3.2 Installation of the Software	24
3.3 Installation of the Bluetooth Connector	25

3.3.1 Installation on Windows 8	27
3.3.2 Installation on Windows 7	30
3.3.3 Installation on Windows XP	31
3.4 Preparation of the Tymo Therapy Plate.....	35
3.5 Explanation of the LED Display on the Therapy Plate.....	36
3.6 Battery and Energy-Saving Mode	37

4 Software41

4.1 Software start	41
4.2 Work program.....	42
4.3 Therapy Chapter	43
4.4 The main components	45
4.5 The patient account	46
4.6 The patient details	49
4.7 Assumptions and documentation	50
4.8 Assessment tools	52
4.8.1 Weight distribution.....	53
4.8.2 Force	53
4.8.3 Position holding test	54
4.8.4 Active weight transfer	55
4.8.5 Pedometer	56
4.9 Therapy programs	57
4.9.1 Control Mode settings	58
4.9.1.1 Power Mode 1	59
4.9.1.2 Power Mode 2 or dynamic mode for therapy 1D and 2D.....	59
4.9.2 During therapy	60
4.9.3 After therapy	61
4.9.4 Balloon.....	62
4.9.5 Chicken and worm	62
4.9.6 Applehunter	63

4.9.7	Firefighters	63
4.9.8	Cars	64
4.9.9	Labyrinth	64
4.9.10	Dinning.....	65
4.9.11	Trampoline	65
4.9.12	Symbols	66
4.9.13	Math.....	66
4.9.14	Words	67
4.10	Therapy report	68
4.11	Settings.....	69
4.12	Hardware configuration	70
4.12.1	Calibrate sensors	71
4.12.2	Upgrade module.....	72
4.12.3	Hardware testing	72
5	Software	73
5.1	Rolling element 1D and 2D	73
5.1.1	Description.....	73
5.1.2	Technical specifications.....	73
5.1.3	Images.....	74
5.2	Tyrostation electric height adjustable desk with motor and Tyrostation height adjustable stool	75
5.2.1	Description.....	75
5.2.2	Technical specifications.....	75
5.2.3	Images.....	76
5.3	All-in-One PC, wireless mouse and wireless keyboard	77
5.3.1	Description.....	77
5.3.2	Technical specifications.....	77
5.3.3	Images.....	79
5.3.4	Notes	80

tyromotion

5.4 Tymo X-Pad	81
5.4.1 Description	81
5.4.2 Technical specifications	81
5.4.3 Images	81
6 Delivery set	82
7 Tymo Therapy Plate Maintenance	84
7.1 Cleaning the Tymo Therapy Plate	84
7.2 Repair Tymo Therapy Plate and its elements	85
8 FAQ - Frequently Asked Questions	86

1 Introduction

Dear Customer!

Thank you very much for deciding to purchase the Tymo Therapy Plate.

Due to the use of special sensors and the high-quality processing of the product, this system makes it possible to record and process movements, process information of weight and the targeted application of strength by patients, for therapeutic and training purposes. It can thus be used in the guidance of specially developed therapy programs and for the determining of progress values during the rehabilitation process.

In these instructions, you will find all the directions you will need for the operation and handling of the Tymo Therapy Plate. For further information on our products, visit our website or contact us personally. www.Tyromotion GmbH.com

1.1 Safety

1.1.1 Intended Use

Tymo Therapy Plate for testing and rehabilitation is a therapeutic device that is above all to be used with patients with restricted physical and/or mental abilities. The device is used to test and rehabilitate balance and equilibrium functions.

1.1.2 Indications

1. Stroke.
2. Traumatic brain injury (TBI).
3. Cerebral hemorrhage.
4. Spinal cord injury.
5. Brain tumor.
6. Multiple sclerosis.
7. Cerebral palsy.
8. Parkinson's disease.
9. Muscular dystrophies.
10. Herniated vertebral discs of the spine.

11. As part of neurorehabilitation, device can also be used to train the cardiovascular system, as well as the muscles of the upper extremities.
12. The device can be used for orthopedic rehabilitation or physiotherapy.

1.1.3 Contraindications

As with every other therapy, the doctor in charge of your treatment is always responsible for the indication of therapy. The same indications and contraindications as for manually conducted physiotherapy treatment always apply to therapy and/or training with the Tymo Therapy Plate. Typical contraindications are for example:

1. Substantial pain symptoms despite conventional pain therapy in the area of the affected members.
2. In case of severe cognitive disorders, aggressive behavior, psychotic illnesses or neurotic disturbances the possibility to train on the Tymo Therapy Plate should be carefully considered in each individual case.
3. Epilepsy: the device cannot be used.
4. Severe ataxia: the device cannot be used.
5. Severe osteoporosis: there exists the risk of fractures. The device cannot be used.
6. Severe spasticity (Ashworth scale 4-5).
7. Unstable fractures: do not carry out any training if you have unstable or still inadequately consolidated fractures.
8. Arthritis of wrist and/or finger joints.
9. Swelling/infection/fracture.

tyromotion

10. Pacemaker: patients with pacemaker may use Tymo Therapy Plate only if distance between them is at least 15 cm. Modern pacemakers are resistant to this type of influence.
11. Infections.
12. Joint problems: the strain on joints during training can cause pain and irritation in states of reduced load-bearing strength.
13. Orthostatic circulatory problems.
14. Skin problems: before and after every training session, check for injuries and pressure marks that existed before or have been caused by the training.
15. Apraxia, symptoms of neglect.

Use of the device in acute diseases and unstable fractures is strictly forbidden!

Please be aware that your patients may have other contraindications that are not listed here. In the case of questions or feedback, please contact Tyromotion GmbH (for contact data see start of the document).

In order to ensure risk-free operation, the following measures and instructions are to be applied.

Is the system suitable for the intended therapy:

1. Continual checking of the therapy success of the patient.
2. Observation of current studies on the subject area.
3. Ceasing of the treatment if there is a sign of disproportional worsening of the patient's condition.

1.1.4 Side Effects

No Side Effects.

1.1.5 Environmental Conditions

The Tymo Therapy Plate is only intended for use in therapeutic rehabilitation centers and private households with moderate environmental conditions. Training outside the premises intended for the purpose is always to be avoided.

1.1.6 Proper Use and Precautions

The Tymo Therapy Plate is a therapeutic device that is above all to be used with patients with restricted physical and/or mental abilities. In order to prevent problems with the use of the device, and in addition to the contraindications described in the Chapter 1.1.3 the following instructions for proper use are also to be adhered to:

- Only use the Tymo Therapy Plate on a firm surface.
- With dynamic measurements, only use the Tymo Therapy Plate with the Tymo Therapy Plate -rolling elements specially developed for this. Placing another object underneath the Tymo Therapy Plate is expressly not permitted by the manufacturer.

tyromotion

- Only use the Tymo Therapy Plate in closed rooms intended for therapy.
- Make sure that the patients who train standing on the Tymo Therapy Plate are wearing clean shoes.
- While you assist a patient in the therapy on the Tymo Therapy Plate, wear shoes that prevent your feet from getting underneath the Therapy Plate during dynamic measuring.
- Switch off the Tymo Therapy Plate when the Tymo Therapy Plate is not in use - the battery life will increase to more than a week

No claim is made for the completeness of this list; it is instead intended to serve as a basis for correct use.

1.1.7 Responsibility of the Owner

The owner is responsible for ensuring that all persons operating the system have read and understood this manual. No guarantee is made, however, that a person who has read this manual is qualified to operate, inspect, test, calibrate, repair, modify or remedy system faults in the Tymo Therapy Plate. The owner is to ensure that the installation, maintenance, calibration and repair of the system and the remedying of faults are only carried out by properly trained, fully qualified employees.

If the Tymo Therapy Plate does not function correctly or does not react to the commands described in this manual, the operator is to contact Tyromotion GmbH.

Despite its robust and high-quality construction, the Tymo Therapy Plate only has a limited ability to withstand shocks. If the device falls from table onto the floor, defects can occur in the sensors or the battery. These defects may possibly not become recognizable until sometime after a fall. Therefore, do not use the Tymo Therapy Plate after a fall; instead, immediately contact Tyromotion GmbH or the dealer from whom you bought the Tymo Therapy Plate.

Address for reclamations:

ZAO "SP Beka-Hospitec", Russia, 124489, Moscow, Zelenograd, Alleya Sosnovaya, 6A, bld.1.

Phone: +7 (495) 742-44-30,

Fax: +7 (495) 742-44-35,

E-mail: info@beka.ru.

1.1.8 Errors and Incompleteness

If you detect errors in this manual or find that something is missing from it, please contact Tyromotion GmbH.

1.1.9 Property of Tyromotion GmbH

The content of this user manual, including all data and illustrations is copyrighted information of Tyromotion GmbH and is provided exclusively for operational, maintenance and repair purposes. Distribution for other purposes and copying without prior written consent from Tyromotion GmbH is forbidden.

1.1.10 Guarantee & Exclusion of Liability

Tyromotion GmbH guarantees the purchaser of the system that the system for a period of 12 months from the date of installation at the premises of the owner, in the case of normal use, will be free from defects in materials and quality of manufacture, and that it will correspond to the mechanical and electrical specifications published by Tyromotion GmbH (unless the warranty period is extended by an optional service contract). This guarantee is granted on the condition that the system is correctly installed, operated and maintained in accordance with the instructions for use.

The customer is to direct all warranty claims to Tyromotion GmbH in writing within 60 days after the occurrence of a problem and before the elapsing of the warranty period. The only obligation of Tyromotion GmbH in accordance with the warranty consists of repairing, replacing or correcting, in its sole discretion, parts that are defective or non-compliant. After the repair or replacing of the defective or non-compliant parts, Tyromotion GmbH has no further obligations towards the owner as regards such defective or non-compliant parts. All repairs and maintenance work in accordance with

this guarantee are to be carried out by an authorized Tyromotion GmbH service representative.

If repairs, maintenance or other work is carried out by third parties, the above warranty is rendered null and void. Beyond this, problems arising due to accidents, improper use, incorrect application, damage due to incorrect storage, negligence or modification of the system or its components are excluded from the above warranty.

The above warranty is granted in place of all other warranties, rights or conditions, and the system is delivered "without warranty for defects" except for this limited warranty. Tyromotion GmbH and its third-party suppliers refuse, specifically and without limitation, all other explicit or tacit warranties towards the owner, his personnel and patients, customers and users and absolutely any third party, including without restriction absolutely any guarantees of marketability, suitability for a certain purpose or non-infringement, and absolutely any guarantees arising from performance development, a business development or a commercial custom.

Tyromotion GmbH and its third-party suppliers make no declarations or guarantees regarding whether the system meets the requirements of the owner or runs free from interruption, errors or defects.

Tyromotion GmbH is under no circumstances liable for indirect, incidental, special or consequential losses or for punitive damages including but not limited to the loss or shortfall of profit, income, goodwill or use that arise for the owner or third parties or for damage to related equipment, for the costs of replacement products, setting-up, services, replacement parts or downtime or for claims of patients, customers, visitors or employees of the owner or of other people, irrespective of whether in the framework of a contractual complaint, arising from an unlawful act, strict liability or as imposed by the law or otherwise, even if Tyromotion GmbH has been informed of the possibility of such losses.

The liability of Tyromotion GmbH for losses that occur from or in connection with this contract must under no circumstances exceed the purchasing price of the system.

Some jurisdictions restrict the scope of limitations or ruling-out of legal remedy or compensation or liability, or rule it out, such as liability for gross negligence or deliberate wrong conduct according to the above or to the above described extent or do not permit the exclusion of tacit guarantees. In such jurisdictions, the limitation or the ruling-out of guarantees, legal remedy, compensation or liability as described above, may possibly have no validity for the owner. Even if such limitations or exclusions do not apply in the legally forbidden degree, they apply to the greatest extent legally permitted. The owner can also have other rights which vary according to the state or other jurisdictions.

1.1.11 Service Life of medical device

Service life of this basic equipment is 7 years. Some parts of this like as Table, Stool, Mouse, Keyboard have different service life (see para. 5.2.2 and 5.3.2)

1.1.12 Environmental Protection

For the purposes of the environmental protection it is forbidden to dispose this device and its accessories with domestic waste. Please return the device to the Manufacturer or its authorized representative in the territory of your country. Computer should be disposed according to local regulations.

1.1.13 Utilization

Special waste: The device must not be disposed of with the normal household waste. Please contact with Tyromotion GmbH or its representatives for disposal.

2 Technical Specifications

2.1 Tymo Therapy Plate

Device:

- Sensor plate for the detection of forces and angles of inclination
- TyroS Basic Software V. 2.0 with integrated measuring modules and interactive therapies


Fig. 1. Tymo Therapy Plate

Features:

- Total weight <5kg;
- Dimensions: 51 x 30 cm (LxW);
- Height: 1.7 cm (from floor with legs), 0.8 cm (without legs);
- Maximum load: 150 kg;
- Upper part of plate is made of stainless steel, AISI 201;
- Evaluation and comparison of the two upper extremities;
- Evaluation and comparison of the two lower extremities;
- Interactive movement training with audiovisual feedback;
- Connection to the supplied Software via Bluetooth technology (transmission range of 3 m);
- Display of the current status of the device by built-in RGB-LED;
- Optionally available with All-in-One PC and pre-installed Software.

Measurement Ranges:

- Force measurement range: 0 – 120 kg;
- Inclination measurement range: approx. 10°;
- Resolution of force measurement: 0.1kg / force measurement accuracy: <1%;

TyroS Basic V.2.0 and TyroS Professional 2.0 Software Features:

- Manual switching option for the Software whether tilt mode or static mode is used.
- Storage and selection of individual therapy aims for each of the available measurement methods.
- Issue date TyroS Basic 2.0, TyroS Professional 2.0: 09.07.2014.
- **Measurement programs**
 - Weight Distribution;
 - Force;
 - Posture Test;
 - Active Weight Transfer;
 - Pedometer.
- **Therapy programs**
 - Applehunter;
 - Firefighters;
 - Balloon;

- Motorway;
 - Chicken and worm;
 - Labyrinth;
 - Dinner;
 - Trampoline;
 - Symbols;
 - Math;
 - Words.
- Task-oriented training with motivating aspect.
- Cognitive and coordinative tasks with audiovisual feedback.
- Database with detailed therapy progress and measurement data of each patient.
- Therapy report for the documentation of therapy and measurements.

Minimum System Requirements:

- Windows XP, Windows 7 or Windows 8; Home or Professional;
- Intel or AMD processor with min. 1.6GHz;
- Min. 1GB RAM;
- Free USB port;
- Screen resolution: 1 280x768 pxl.;

- Microsoft .NET framework 1.1 (with Windows installed as standard);
- Adobe Flash Player 10 ActiveX (is installed with the Software).

Approval:

- EU: Medical Device Class I in accordance with 93/42/EEC – CE;
- According to the Rule 12 of the Commission Directive 93/42/EWG, Annex IX, Tymo Therapy Plate is a class I medical device.

Tymo Battery:

- Rated capacity: 2 500 mAh;
- Nominal voltage: 3.7 V;
- Charging current (max.): 2.5 A (1C);
- Charging voltage (max.): 4.2 V (4.2 V +- 0.03 V);
- Weight: approx. 50 gm;
- Dimensions (LWH): 85 x 55 x 5 mm;
- Battery life in full operation: min. 8 hours.

tyromotion

Page intentionally is left blank

2.2 Markings on the Tymo Therapy Plate

Marking	Meaning
TYMO® Designed and produced by tyromotion in Austria     tyromotion Bahnhofgürtel 50 8020 Graz Austria	The type label of Tymo Therapy Plate is attached to the bottom side of the Therapy Plate, in the form of a sticker. The symbols included in it are explained in the following: Special waste: the device must not be disposed of with the normal household waste. For disposal, please contact Tyromotion GmbH.
	
	Read instructions: they are to be read through carefully before using the device for the first time.
	User part: electrical safety class II.
	CE label – the product corresponds to the product-specifically applicable European directives.
Tymo Therapy Plate serial number for example: TR1-2013-001	From this serial number you can derive the year of manufacture.

3 Preparation/Basic Information

3.1 Operation, transportation, storage

- **Operation conditions**

The Tymo Therapy Plate is only intended for use in therapeutic rehabilitation centers and private households with moderate environmental conditions. Training outside the premises intended for the purpose is forbidden.

Operation temperature	+10°C to +30°C
Relative humidity	30% to 75 %
Ambient pressure	700 hPa to 1060 hPa

The Tymo Therapy Plate may not be used in hazardous areas in areas AP and APG in accordance with the requirements of IEC 60601-1. It is not allowed to use flammable or explosive anesthetic gases and mixtures thereof, such as ethyl alcohol and cyclopropane.

- **Storage conditions**

Storage temperature	-20°C to +40°C
Relative humidity	30% to 90%
Ambient pressure	550 hPa to 1060 hPa

- Keep out of the reach of children.
- Store the Tymo Therapy Plate only indoors, only use in premises intended for medical purposes.
- Keep the plate in a place protected from damage.

- **Transportation conditions**

Transported by all modes of transport in covered vehicles in accordance with the current rules of this type of transport services.

Transportation temperature	-20°C to +40°C
Relative humidity	30% to 90%
Ambient pressure	550 hPa to 1060 hPa

tyromotion

Packaging requirements:

TYMO device is supplied with accessories in the original packaging: a wooden box with a nested cardboard boxes.


Each single element is placed in a separate box (except table). Devices are fixed in the cardboard box by foam inserts with the slots.

The table is wrapped with stretch film.

tyromotion

The label is affixed to the each boxes with the designation included in the composition of equipment and its serial number.

All equipment is securely anchored.

3.2 Installation of the Software

If you purchased a PC in Tyromotion GmbH, it is already integrated the latest version of the Software, and you can ignore the points listed below. If you want to install the Software on another computer, follow definite steps.

Before installing the Software, make sure that your system meets the minimum system requirements (see Chapter 2.1).

Attention! To install and use the program you need unlimited access to the records in the installation directory. Make sure that you have access. If you are not sure, contact your system administrator. Then perform the following setting points:

1. Enter your PC to your account by entering your username and password.
2. Insert the Software CD into your PC.
3. If your system supports AutoPlay, the installation will start automatically. If it does not, start the installation by double clicking the left mouse button on the icon "install.bat".
4. Follow the instructions of the program. If you have problems, contact your system administrator.
5. After installation is complete, restart your computer, as required by the system.

Please follow the instructions provided in Chapter 4.12 "Upgrading modules" to activate your purchased module Tymo Therapy Plate. (It is necessary if you install the Software on your own PC).

3.3 Installation of the Bluetooth Connector

The Tymo Therapy Plate is connected to your PC with Bluetooth technology. In order to produce this connection it is first necessary to connect the Bluetooth connector included in the scope of delivery to your PC and to install it. If you have purchased a PC from Tyromotion GmbH, this installation will already have been carried out.

Otherwise, connect the Bluetooth connector to a free USB port of your PC and follow the installation instructions of your PC.

Important points before the installation:

- The Tymo Therapy Plate must be fully charged.
- The Tymo Therapy Plate must not be more than 3m away from the PC.
- The LED on the top of the Tymo Therapy Plate must pulse green after the battery has been charged and the charging unit has been unplugged.

tyromotion

Plug Bluetooth adapter into the USB port:

Step 1



Step 2

3.3.1 Installation on Windows 8


The drivers are now automatically installed. The symbol for the device installation appears for a short time on the task bar:

Open the active edge by moving the cursor through the corner on the bottom right and click on „Settings“:


Select “Devices and Printers” and click on “Add Device”:


After a short while, all nearby Bluetooth devices appear. Select the entry “Serial Adapter” or the device with the Tymo Therapy Plate serial number and click on “Next”. Enter the pairing code “1234” and click on “Next”:

tyromotion

Enter the passcode for your device

The settings are now carried out automatically. Then click on “Close”. Installation is completed and the Tymo Therapy Plate is now ready for use.

Installing Serial Adaptor

Please wait while Setup installs necessary files on your system. This may take several minutes.


3.3.2 Installation on Windows 7


The drivers are now automatically installed:

Open the menu bar by clicking on the Windows symbol on the far left and select "Devices and Printers". Click now on "Add a device":

After a short while, all nearby Bluetooth devices appear. Select the entry "Serial Adapter" or the device with the Tymo Therapy Plate serial number and click on "Next". Click on the 2nd entry "Enter the device's pairing code":

Enter the pairing code "1234" and click on "Next":

tyromotion

The settings are now carried out automatically. Then click on "Close". Installation is completed and the Tymo Therapy Plate is now ready for use:


3.3.3 Installation on Windows XP



The drivers are now automatically installed:


Open the Menu bar by clicking on the Windows symbol on the far left and select "Control Panel" and Click on "Printers and other Hardware":


Click on "Bluetooth Devices". A new window pops up. Now click "Add...":

Again a new window pops up. Now tick "My device is set up and ready to be found" and then click "Next >":

After a short while, all nearby Bluetooth devices appear. Select the entry "Serial Adapter" or the device with the Tymo Therapy Plate serial number and click on "Next". Click on the 2nd entry "Let me choose my own passkey". Then enter the passkey "1234" and click on "Next":

The settings are now carried out automatically. To accomplish the installation click "Finish":

3.4 Preparation of the Tymo Therapy Plate

1. Take the Tymo Therapy Plate and place it on an even surface (depending on use) such e.g. on the floor or a table.
2. Ensure that the Tymo Therapy Plate stands steadily on the respective surface. If you work with the Force Mode, please be sure that all 4 Sensors (this are the 4 feet of the Tymo Therapy Plate) are on a stable surface.
3. When the Tymo Therapy Plate is supplied, the integrated battery is uncharged. Before first use, charge the battery according to the description in Chapter 3.6. Then proceed as described here.
4. Insert the provided Bluetooth connector into a free USB port on your PC.
5. Install the driver for the Bluetooth connector (the driver for this device is included in the scope of supply but is automatically recognized by Windows). For the operating systems Windows 7 and Windows XP an exact description of how to proceed with the settings is included in the scope of supply.
6. Start the Tymo Therapy Plate Software by double-clicking on the Tymo Therapy Plate symbol on your desktop.
7. The Tymo Therapy Plate connects itself to your PC automatically via Bluetooth adapter as soon as you start a therapy session with Tymo Therapy Plate (Chapter 4.3 of the Instruction Manual) or open the hardware configuration (Chapter 4.12 of the Instruction Manual). Please note the LED display on the Tymo Therapy Plate (Chapter 3.5). It gives information on which status the Tymo Therapy Plate is currently in.

3.5 Explanation of the LED Display on the Therapy Plate

The Tymo Therapy Plate has an LED display attached to the top of the plate. By using different colored and flashing signals, a large number of different operational states can be indicated by means of this display. Please see the following table for a description of this.

LED Display	Meaning	Necessary Action
Shining red	Connection with Software okay, battery level under 15%	Charging of the battery (Chapter 3.6)
Shining green	Connection with Software okay, battery level okay	None
Pulsing red (brief flaring up)	Device is in energy-saving mode, battery level is less than 15%	Charging of the battery (Chapter 3.6)
Flashing blue	Device is in energy-saving mode, battery is being charged	None
Shining blue	Device is in energy-saving mode, battery-charging process completed.	Removal of the charging cable.
Flashing green	Battery is charged, device is in energy-saving mode,	None

Alternating
between
green / blue

Connection with Software okay,
battery is being charged

Follow instructions
from Chapter 3.6.

Off

Battery of the device is
completely flat or the device is
switched off

Charging of the battery or
switch on the device
(Chapter 3.6)

Because the Tymo Therapy Plate can be in several operational modes at the same time (e.g. "charging of the battery" (flashing blue) and "connection to PC established" (shining green)), priorities are assigned to the different signals. Higher priorities have precedence over lower priorities. In the example described, the LED display would therefore flash blue.

3.6 Battery and Energy-Saving Mode

The integrated battery of the Tymo Therapy Plate can be charged via the charging unit included in the scope of supply. Because this charging unit has been tested as a medical device, the Tymo Therapy Plate is also electrically safe during a charging process. Follow the instructions, however, about the use of the device during a charging process. The use of charging units that have not been recommended by Tyromotion GmbH is expressly forbidden, because unforeseen risks could arise for people and equipment.

On/Off switch:

It is recommended to switch off the Tymo Therapy Plate when it is not used for a longer period of time as the rechargeable battery service life can be increased considerably in that way.

- When the button on the underside (it is next to the type plate) is pressed, the device is switched off > LED display not illuminated any longer.
- When you press it again, it is switched on > LED display illuminated again.
- When the charger is plugged into the Tymo Therapy Plate, the Tymo Therapy Plate switches on automatically.
- When the charger is unplugged again, the Tymo Therapy Plate remains switched on.
- When the rechargeable battery does not have enough power, it is not possible to switch on the Tymo Therapy Plate.
- After switching on, the Tymo Therapy Plate is ready for operation within 5 seconds.
- Only switch off the Tymo Therapy Plate after having exited the Tymo Therapy Plate Software.
- Do not switch off the Tymo Therapy Plate during the calibration process.

Charging Battery

In order to charge the battery, first connect the charging unit provided to a 230V plug and then insert the other end into the charging socket of the Tymo Therapy Plate. The Tymo Therapy Plate switches into the charging mode automatically (Note the LED display described in Chapter 3.5).

NOTE: The Tymo Therapy Plate is also fully capable of being used while in the charging mode. In this case, however, there is an increased risk of tripping for therapists and patients due to the charging cable. Furthermore, there is an increased risk of a mechanical defect if the Tymo Therapy Plate is used during a charging process. Tyromotion GmbH can accept no liability for harm caused to people or the equipment by this.

Duration of Charging

The duration of charging depends, among other things, upon whether the battery was fully uncharged at the start of the charging process or not. As a general rule, however, the charging process should not last longer than 6 hours. In this connection, see the LED display on the Tymo Therapy Plate. If the LED flashes blue, then the charging process is active. If the LED shines blue, the charging process has been completed and the battery is fully charged again.

Further information about the LED display can be found in Chapter 3.5.

Energy-Saving Mode

The battery runs down at a varying speed, depending on how frequently you use the Tymo Therapy Plate. If Tymo Therapy Plate is not used for several minutes, the device switches automatically into an energy-saving mode. It stays in this mode until the device is automatically awoken by the start of the Tymo Therapy Plate Software.

Charging Cycles and Lifespan:

Due to the use of high-quality lithium-polymer batteries, it is ensured that the battery running times still meet the requirements even after many charging cycles. Nevertheless, no guarantee can be given for the lifespan of batteries going beyond 6 months after the date of purchase.

If you adhere to the following advice, you can certainly save the battery of your device and thereby extend its lifespan:

- Only charge the battery of the device when the battery level is approaching the end (can be recognized by red flashing, 50% on, 50% off, of the LED display).
- Avoid direct exposure of the Tymo Therapy Plate to sunlight in order to prevent the battery inside the device from becoming too hot.
- Use exclusively the charging cable provided by Tyromotion GmbH for the recharging of the battery.
- When the Tymo Therapy Plate is switched off before longer periods of non-use, the rechargeable battery service life is extended to over one week (standby mode).

4 Software

Software may be used only with the Tymo Therapy Plate. Tymo Therapy Plate offers a large variety of therapy variants, which can be independently chosen by a doctor according to therapeutic needs to receive certain results.

However, besides this, training in the Tyromotion GmbH or with a certified specialist of Tymo Therapy Plate is obligatory.

Note that, depending on the purchased Software package, some options may not be available (Chapter 2.1).

Software designed to work in all major operating systems. However, the manufacturer guarantees trouble-free operation only if the Software is installed and entered on a computer system supplied by Tyromotion GmbH. Installing Software may also have a negative impact on the Software and the manufacturer cannot guarantee the functioning and / or perform equipment maintenance.

4.1 Software start

Start by double-clicking on the program icon. The program will start the required program elements before the opening of the interface (Fig. 2). When you first start Tymo Therapy Plate must be calibrated (the procedure should be repeated regularly after the first use). Follow the instructions in Chapter 4.12.

Fig. 2. Loading program components

4.2 Work program

Start screen

After starting the program the start-up screen appears where you can choose the following options:

- **Initiation of therapy with Tymo Therapy Plate - Chapter 4.3.**

There is carried out therapy with Tymo Therapy Plate and stored patient data. In the transition to this Chapter of the program Tymo Therapy Plate automatically switches from power saving mode to the active mode.

- **Changing the settings of the program - Chapter 4.11.**

This Chapter is used to configure the parameters in accordance with the requirements of your establishment and regulation of the system according to your needs.

- **Change the hardware configuration - Chapter 4.12.**

In this Chapter, you can calibrate the platform Tymo Therapy Plate, change the software modules and to test equipment.

Fig. 3. Menu of starting screen

Language selection

Fig. 4. Language selection and right buttons

The standard program includes seven languages. Language selection is performed in the upper left corner of the screen by clicking the left mouse button on the appropriate symbol (see. Fig. 4).

To get information about the current version and the manufacturer, click the question mark.

4.3 Therapy Chapter

When navigating to therapy on the left side of the screen toolbar appears (Fig. 5). This panel is designed to work with the program.

Fig. 5. Toolbar

Important: When you select an item from the toolbar, the current measurement or therapeutic session will be interrupted. To prevent loss of information, save measurement or the current session.

Start Screen. Clicking on this button returns you to the home screen menu (see Chapter 4.2, Fig. 3).

Patient selection. With this button you can go to the patient database to create a new account or select an existing patient account.

Patient data. It displays information about the patient and all the therapy and measurements.

Therapy report. Therapy report displays various indicators that you can use in the patient's chart (see. Chapter 4.10)

Assessment Tools. The button opens a menu of assessment tools (see Chapter 4.8).

Therapeutic programs. The button opens a menu of therapeutic programs (see Chapter 4.9).

LED connection and charging status. Displays the current state of charge of the battery and Bluetooth-connection with Tymo Therapy Plate.

Centering. Press this button before starting therapy or assessment to reset Tymo Therapy Plate sensors. While pressing the plate should be on a flat surface and without load. (This action is not a substitute for calibration as described in Chapter 4.12).

4.4 The main components

To ease operation and recognition, a number of symbols are doubled in the program of Tymo Therapy Plate. The following information will allow you to carry out individual therapy sessions and measurement.

 **Start (therapy).** Load settings, the customization options are disabled and the therapeutic session begins.

 **Pause.** Pause therapy session. Now you can change the settings. Pressing this button to resume therapy session.

 **Start (measurements).** It resets the first measurement and starts a new record values of the sensor.

 **Stop (measurements).** Stop recording of the sensor indications. Although the indications are still displayed on the screen, they are not stored and are not used in the analysis.

 **OK / stop.** If you press this button during the therapy session, it is completed, and you go into the menu.

 **Close / finish.** Measurement is terminated without saving the data, even during the measurement session.

 **Save and Close.** The measurement results are saved and session measurement is completed.

 **Starting position and documentation.** You go to the section of the patient to secure the individual therapeutic goals and results (see Chapter 4.7).

4.5 The patient account

Patients' data should be set before starting the therapeutic session or a session of measurements. Data that had already set in database of patients can call again.

Fig. 6. The patient account selection

The list includes the data of all patients whose parameters have been saved. To facilitate the search of the patient, the program provides filters (name, date of birth, Social Security number and status of backup (yes / no)). Backup Status determines whether the patient receives treatment in the future, or it has already been discharged.

Filter functions allow you to search the patient is not only the first but also the subsequent letters of the name, i.e. letters in the middle of the name. This feature makes it easier to search in case if the order of the first or last name has been broken.

The patient accounts are selected with the mouse. Such database features as "edit data" and "delete account" may also be used. In addition, you can use the main toolbar to generate reports on the treatment of selected patients. If you want to start a therapy session with the selected patient, confirm the account (use the basic elements presented

in Chapter 4.4), or select an account by double-clicking the left mouse button. The following is a detailed description of the patient account's choice.

New patient account

Creating the patient account is carried out in two stages:

1. Enter the original patient data
 - Name;
 - Date of Birth;
 - Number of health insurance;
 - Gender.

Before you go to step 2, save the data by clicking on the diskette symbol. It should be particularly attentive to these input parameters, because they will be used in future applications as filters to select the patient's account.

2. Enter the therapeutic effects of patient
 - Affected parties;
 - The dominant hand;
 - The date of issuance;
 - Assumptions and documentation (Chapter 4.7);
 - Additional information about the patient.

Edit patient account

You can always change the stored patient information. To do this, select the account and click on button [edit data]. The further procedure is similar to creating a new account (see. [New patient account]). If further therapy is not planned for the patient, it is recommended to set the status of a backup on "yes". If the filter of criterion for archiving is set to "no", the patient account is no longer displayed in the account list.

Delete patient account

To delete the patient account, select the account and click [delete patient].

IMPORTANT! This action involves the permanent removal of all patient information.

Existing the patient account selection

Stored the patient account can be selected by double-clicking the left mouse button on an account by pressing the "OK" button or by selecting the item "patient details / patient information" on the toolbar.

4.6 The patient details

The "Patient details" displays information about all the passed therapy sessions and records of the selected patients recorded in the documents. The program saves all passed sessions with date, time and duration of therapy. In addition, it displays the original position and these comments. With double-click on one of the presented therapy sessions or performance (and / or by pressing the [display data] you can redisplay the required values. Unnecessary or incorrect parameters can be removed from the list by clicking [delete therapy result].

When displaying measurement information (by double-click or the [display information] button, you can also navigate to display all previously performed measurements. The measurement results are displayed in chronological order and can be run with the arrow keys (see. Fig. 8).

Fig. 7. The patient details

Fig. 8. Navigation

4.7 Assumptions and documentation

This section is particularly innovative and multi-functional tool and documentation in the program. You can adapt the various measurements made with the platform Tymo Therapy Plate in accordance with your patient needs, appointed these therapeutic targets in the form of tasks (starting position in the program). At the same time, in this section you can leave comments on the therapy and tasks, which facilitates the daily work with documents.

You can always go to this section, if you see the following symbol:

Is not necessary to indicate starting position for the test on hold poses and in the module "Pedometer", because the position is determined in advance.

The comments section presents all ever mentioned comments on the patient. These comments can also be reused. You can add a comment to this current measurement by dragging it with the mouse, using the [+] button or double-click. All new comments can be written directly in the comments box.

Documentation

Fig. 9. Documentation

4.8 Assessment tools

In the menu item "Assessment tools" are a set of measurements developed in collaboration with different doctors. To start treatment, press the corresponding button.

In the assessment tools you can control the sequence using the basic elements described in Chapter 4.4.

Assessment tools can be expanded with individually designed sequences, such that each patient will be optimally worked through its therapeutic problems. Due to the integrated functions of the starting position and the documentation you can create and continue to use an unlimited number of therapy options for their patients. For more information, see. Chapter 4.7 "Starting position and documentation."

The transition into this section by using the symbol:

Assessment Program has been developed for the display of the dynamics of therapy.

Assessment tools

Please select your desired assessment tool


Fig. 10. Assessment tools

4.8.1 Weight distribution

Symmetry assessment occurs within predetermined 5 seconds. During this time the patient stands on Tymo Therapy Plate, and the doctor can give him simple tasks (e.g., sit down, bend your knees, just stand, etc.). After this time the system automatically displays how evenly distribute the load on the plate during the test (for example, the left side of - 25%, right - 75%).


Fig. 11. Weight distribution

4.8.2 Force

Force assessment is designed to calculate the force that the patient can attach to the Tymo Therapy Plate. A method of pressing on the plate can be chosen arbitrarily.

To do this, go to "Starting position and documentation" presented in Chapter 4.7. To calculate the force of patient legs, it should be located in front of the plate to activate display showing the force of the left side.

Now ask the patient as much as possible to put pressure on the plate only left foot. Repeat the whole procedure for the right leg.

Force assessment can be carried out in a similar manner with the upper extremities.

4.8.3 Position holding test

This assessment tool is designed for tests that are widely-distributed in medical centers. With it you can estimate how long a patient can hold position without losing balance.

Estimation takes place in three stages:

- Stand on two legs with opened eyes;
- Stand on two legs with closed eyes;
- Stand on one leg with open eyes.

Put the patient in one of the assumptions on Tymo Therapy Plate and click on the symbol with the appropriate schedule for the start of the estimation. Estimation is carried out within 20 seconds, it may be interrupted only if the evasive compensator motions of the patient too pronounced. Then you can save the results or to delete them by clicking on the appropriate symbols. Estimation may be performed in any order and any number of times.


Fig. 13. Position holder test

4.8.4 Active weight transfer

Balance estimation allows you to calculate the active range of motions of the patient. To do this, place the patient on Tymo Therapy Plate and ask them to transfer the weight of the body as far as possible in all directions. The more the patient transfers the body weight, the more intensive the relevant sectors change color, the higher balance sheet (specified in percentage).

An alternative is to ask the patient to move with the lower amplitude as in testing "Position holder" (Chapter 4.8.3) Thus, it can be evaluated the patient deviation from its position. To do this, enter the settings specified in Chapter 4.7 for the starting position.

Fig. 14. Active weight transfer

4.8.5 Pedometer

Here, the patient is asked to rise on the plate and alternately transfer weight from one foot to the other. The graph displays the direction of movement of the patient. Take the measurement using the basic functions provided in Chapter 4.4. During the measurement program counts the number of performed steps and the average number of steps per second. The assessment shows the highest step frequency.


Fig. 15. Pedometer

4.9 Therapy programs

Eleven specialized modules with audio-visual and tactile feedback contribute to the optimal concentration of the patient and increase motivation. These modules are used in all therapeutic complexes of Tyromotion GmbH, and if you work with different equipment of Tyromotion GmbH, your patients will know them. Due to therapeutic programs you can achieve a variety of purposes when training with patients.



Fig. 16. Therapy programs

Before starting the selected therapeutic module should determine a way to control the tasks. Relevant information is loaded from a database or entered in the control settings, which opens automatically if this information has not been stored for the individual patient.

There are three control modes of therapeutic session, which you can select by clicking on the corresponding symbol. **Power Mode 1** - control of force application EXCLUDING place of application (e.g., pressure with one hand in the center of the plate). **Power Mode 2** – control of force application with the place of application (for example, one hand - the left side of the plate, and the other - on the right). **Dynamic mode** - the use one of the two clips under the plate.

Control Mode

4.9.1.1 Power Mode 1

The patient should press on the Tymo Therapy Plate as strong as possible, regardless of the starting position. Indicators of force are highlighted in orange color. The displayed data is only for the visual presentation of information and do not contain indicators in numbers.

4.9.1.2 Power Mode 2 or dynamic mode for therapy 1D and 2D

It is estimated the transfer of body weight, and parameters are used to control therapy. This mode also takes into position of force application on the plate.

During estimation of motions for therapeutic applications in one plane (1D) the errors curve displays. The longest diagonal of this curve shows the control motions.

In therapeutic applications 2D movement are displayed by a square. This square displays the range of motions, which will be used in therapeutic sessions.

4.9.2 During therapy

During therapy in the game space various symbols are displayed. Symbols in Fig. 19 have the following meanings:

Fig. 19. Display during therapy

The left top angle: the game's level and duration.

The circle shows the current level of the game / level of difficulty. Fig. 19 displays the first level.

Next to the level figure displays the duration of the game. Depending on the type of therapy count may be direct or reversed. In some therapeutic sessions time also affects the score.

The right top angle

It displays information about the game, which is usually indicated by three symbols.

 The number that is next to the exclamation mark means the number of points required to pass the level.

 The number that is next to the "tick" represents the number of earned points.

tyromotion

In most cases, there is a third component in therapy. This number indicates the number of remaining attempts to recruit the required number of points.

GENERALLY, these symbols are also displayed in full-screen mode. To activate this mode, during therapy click in the field of game by the left mouse button. To return to normal mode, press ESC.

4.9.3 After therapy

Points. At the end of therapy the patient sees how many points he earned at the current level. This indicator (1-100) is stored in the database with reference to the appropriate patient and session.

Smiles. Depending on the number of points appropriate smile is displayed. As a rule, the program displays funny emotions.

4.9.4 Balloon

Task: maneuvering of the balloon through the obstacles.

Control: the dynamic sequence of movements and the force for a long period of time.

Training skills: control / coordination of motions, force monitoring, position control, continuous contraction of muscles, concentration, balance.

Control: see. Chapter 4.9.1.

Fig. 20. Balloon

4.9.5 Chicken and worm

Task: control chicken that catches worms on the ground.

Training skills: doing targeted repetitive movements, complex movements performed in daily life, precise control of movements, preventive compensatory movements during training, monitoring and normalization of vitality, orientation and reaction training in space.

Control: see. Chapter 4.9.1.


Fig. 21. Chicken and worm

4.9.6 Applehunter

Task: During the game you need to catch the falling apples in the basket by the active repetition of coordinated motions or force.

Training skills: movement control, force (dosing), vitality control, purposeful movements, coordination, attention, balance and position holding.


Fig. 22. Applehunter

Control: see. Chapter 4.9.1.

4.9.7 Firefighters

Task: During the game you need as close as possible to extinguish the fire by water jets. Achieving and maintaining the required level of force and / or execution of movements.

Training skills: force control, force dosing, coordination, purposeful movements, prolonged contraction of muscles, attention, balance, position holding.

Fig. 23. Firefighters

Control: see. Chapter 4.9.1.

4.9.8 Cars

Task: During the game the patient controls car avoiding obstacles. Targeted achievement end positions and intermediate values.

Training skills: attention, balance, position holding.

Control: see. Chapter 4.9.1.

Fig. 24. Cars

4.9.9 Labyrinth

Task: During the game you need conduct the ball through labyrinth with obstacles. Implementation of effective active motionы without compensatory movements.

Training skills: See “Chicken and worm” module description.

Control: see. Chapter 4.9.1.


Fig. 25. Labyrinth

4.9.10 Dinning

Task: Correct arrangement of dishes on the table (as if the patient is standing at the table). Implementation of effective active motions without compensatory movements.

Training skills: aimed repetitive movements, complicated daily movements, suppression of compensatory movements, orientation in space and cognitive skills. Green point operates as a cursor and disappears after taking an object.

Control: see. Chapter 4.9.1.


Fig. 26. Dinning

4.9.11 Trampoline

Task: Displacement of figures on moving trampoline.

Training skills: focus on the appointed sequence of movements, tilts and / or directional control of movements and force.

Control: see. Chapter 4.9.1.


Fig. 27. Trampoline

4.9.12 Symbols

Task: Recognition and comparison of similar symbols using red cursor. One of the symbols is activated and selected the same symbol. Then, in the center of the playing area the next symbol is activated.

Training skills: Defined movements with hand-eye coordination, posture control and the recognition of simple symbols.

Control: see. Chapter 4.9.1.

Fig. 28. Symbols

4.9.13 Math

Task: doing simple sums in mind, selecting the correct answer with the cursor as in the module with the symbols.

Training skills: Defined movements with hand-eye coordination, posture control and doing simple sums.

Control: see. Chapter 4.9.1.


Fig. 29. Math

4.9.14 Words

Task: reading simple words and their comparison with the corresponding symbols. Control is similar to the game "Symbols."

Training skills: Defined movements with hand-eye coordination, posture control and dyslexia correction.

Control: see. Chapter 4.9.1.


Fig. 30. Words

4.10 Therapy report

The therapy report provides all the information about the patient in the assembled form. You can customize the information displayed in the report, noting the relevant points from the left. In addition to the input of patient data, you can store additional information such as the number of wards, beds, type of insurance, etc. You can also fill in the fields of diagnosis and the results that will allow you to conduct a full-fledged document management.

To change these details, click [Edit record].

After you have changed the report in accordance with your wishes and confirm the data, you can print the report or save it as a file. Pay attention to the following information:

	Print
	Create PDF-file with report (including all diagrams of report).
	Data export in txt-document without diagrams. Then you can import this file in Excel, for example.

Operations in therapy report are carried out using the basic elements presented in Chapter 4.4.

4.11 Settings

You can enter in the section of the program settings by using the appropriate button on the start screen. You can make the following settings:

Change of the force displaye - all measured values are displayed in the selected units. You can select one of the following options:

- kg (1 kg corresponds to approx. 9,81N);
- pounds (1 pound corresponds to approx. 4.45N or approx. 0.45 kg);
- Newtons (1N corresponds to approx. 0.1 kg).

Name of the institution. The appointed name is displayed in the report.


Logo of the institution. You can choose the logo of your institution. This logo will be displayed in the report. Permitted images include BMP, JPEG, JPG, or PNG.

4.12 Hardware configuration

Using the hardware configuration button on the start screen, you can get the section of the program that is responsible for the recognition, selection and calibration of the connected equipment. Under the transition to this section the screen displays a list with the number of Tymo Therapy Plate, detected on your computer (usually only one plate).

To select a sector, choose the appropriate sensor and then click "Configure".

In the next screen, you can calibrate the sensor or change the module of Tymo Therapy Plate - see. below.

The screenshot shows a window titled "Sensor selection" with a table of detected equipment. The table has five columns: "Type", "Serial number", "Firmware Version", "HW Version", and "Details". Two rows are visible, both highlighted in orange. The first row shows "TR1" with serial number "TR1-2012-001", firmware version "V2.1", and hardware version "V2.0". The second row shows "TR1" with serial number "TR1-2012-003", firmware version "V2.1", and hardware version "V2.0".

Type	Serial number	Firmware Version	HW Version	Details
TR1	TR1-2012-001	V2.1	V2.0	
TR1	TR1-2012-003	V2.1	V2.0	

Fig. 31. Sensor selection

4.12.1 Calibrate sensors

Due to high-quality sensors and structures Tymo Therapy Plate is stable concerning to the duration of exposure. Nevertheless, we recommend calibrating the sensors on a regular basis.

During Tymo Therapy Plate calibration do the following points:

1. Make sure that the battery is sufficiently charged (see. Chapter 3.5 and 3.6).
2. Remove the device Tymo-Knobs from the bottom side of Tymo Therapy Plate.
3. Place Tymo Therapy Plate on a flat, stable surface.
4. Make sure that during calibration on the plate there are no secondary objects, nobody stands on it and has no pressure on it.
5. Return to the Start screen and click the [Hardware configuration] button.
6. To begin the calibration, press the [calibrate sensors] button in the configuration.
7. Wait until the screen displays "100%", and there will be a green box with the text "Calibration OK".
8. Save the calibration results by clicking on the corresponding symbol.

If the calibration was incorrect, a message "Calibration failed" appears in the window. In this case, make sure that Tymo Therapy Plate is ready for use and repeat steps 1-8.

4.12.2 Upgrade module

Software is divided into 2 modules: TyroS Basic V. 2.0 and TyroS Professional 2.0. Depending on the purchased module you can be accessed various therapeutic functions. The standard configuration of the program is set to the lowest module at delivery. If the module does not correspond bought module, you can improve it by using the appropriate code.

In order to improve the module, follow these steps:

1. Make sure the battery is charged (see. Chapter 3.5 and 3.6).
2. Return to the Start screen and click the [hardware configuration] button.
3. Select Tymo Therapy Plate.
4. Press the [upgrade Tymo Therapy Plate module] button. The field "The current version of the module" displays your version of the software module.
5. If this module is different from the desired, enter the four-digit code in the free field and confirm this action by pressing the [upgrade] button.
6. You will be notified what module was activated. Go one level back.

4.12.3 Hardware testing

Testing Tymo Therapy Plate is an integrated diagnostic program that automatically checks the validity of indicators of sensors. In normal use, such testing is not necessary to carry out. It is an auxiliary function for the manufacturer to solve problems with the equipment.

5 Software

5.1 Rolling element 1D and 2D

5.1.1 Description

Rolling elements are designed to assist patients during therapy. Rolling element 1D is intended for legs, 2D – for elbow. The product is made of ABS plastic, Terluran GP-22.

5.1.2 Technical specifications

Dimensions and mass:

Rolling element 1D:

- Dimensions: 20x15 cm;
- Mass: 0.25 kg.

Rolling element 2D:

- Dimensions: diameter – 15 cm;
- Mass: 0.30 kg.

5.1.3 Images


Fig. 32. Rolling elements 1D and 2D

5.2 Tyrostation electric height adjustable desk with motor and Tyrostation height adjustable stool

5.2.1 Description

Tyrostation electric height adjustable desk with motor – is a special table with adjustable height and compartments to keep Tymo Therapy Plate and its accessories. The table is supplied with two boxes for accessories in order to store CDs, reports and minor items.

Tyrostation includes an ergonomic height adjustable stool for comfortable use of Tymo Therapy Plate.

Material of desk and stool base is stainless steel, AISI 304.

Surface is made of acrylic artificial stone. Staron.

5.2.2 Technical specifications

Desk specifications:

- Dimensions (LxW): 118x102 cm;
- Height: 114 cm, height adjustable to 195 cm;
- Mass: 80 kg;
- Height adjustment method: electric adjustment;
- Rated load: 200 kg;
- Motor power – 200 W;
- Service life: 5 years.

Stool specifications:

- Complete stool height with castors [floor to the top]: 51 cm from the floor;
- Seat dimensions (LxW): 52.5x31.5 cm;
- Height: 38-51 cm;
- Mass: 9 kg;
- Height adjustment method: mechanic;
- Service life: 5 years;
- Castors' characteristics: 4 castors without breaks made of polyurethane, diameter – 50 mm, color – black.

5.2.3 Images

Fig. 33. Stool

Fig. 34. Desk

5.3 All-in-One PC, wireless mouse and wireless keyboard

5.3.1 Description

Computer is intended for installation of motivating software and provision of feedback during training, and for collection statistical and diagnostic data. Patient has no access to the computer!

5.3.2 Technical specifications

All-in-One PC

- Dimensions (WxHxD): 582 x 459 x 50 cm;
- Mass: 13 kg;
- Power supply: 135 V, power supply adapter;
- Screen diagonal: 23.6" (59.9 cm);
- Resolution: 1920x1080 pxl.;
- Large screen display (16:9), Full HD;
- Processor: Intel® Core™ i5 -2400S, Intel® Core™ i3 -2100;
- Supports MultiTouch;
- Gb – 8 Gb, DDR3, 1333 MHz frequency, 2 x SO-DIMM;
- Data storage: 500 Gb – 1 Tb, 6 Gb/s (7200 rpm);
- WN operation mode: 802.11 b/g/n;
- LAN: 10/100/1000 Mbit/s;
- Sound: SonicMaster;
- Microphone: 2 x 5 V;
- Side I/O ports: 2 x USB 3.0; 1 x 3-in-1 Card Reader; 1 x Headphones; 1 x Microphone;

- Rear panel sockets: 4 x USB 2.0; 1 x HDMI-In; 1 x VGA(D-Sub)-In; 1 x TV Jack (optional); 1 x Kensington lock.
- Card Reader 3-in-1: SD/SDHC/MMC

Wireless mouse

- Dimensions: 95 x 62 x 35.6 mm;
- Mass: 61 gm;
- Connection type: wireless;
- Operating range: 10m;
- Wireless connection type: RF 2.4 GGz;
- Service life: 2 years;
- Number of buttons: 2 + scroll ring;
- Sensor type: optical;
- Response time: 8 m/s;
- OS compatibility: Windows;
- Power supply: batteries;
- Color: black/white.

Wireless keyboard

-
- Dimensions (LxHxW): 25.9 x 0.5 x 16.7 cm;
- Mass: 290 gm;
- Connection type: wireless;
- Wireless communication type: Bluetooth adapter;
- Operating range: up to 10 m;
- OS compatibility: Windows;
- Response time: 8 m/s;
- Service life: 2 years;

tyromotion

- Design type: semi mechanic;
- Number of keys: 86;
- Connection interface: USB;
- Keypad water protection: yes;
- Keypad dust protection: yes;
- Power supply: built-in rechargeable battery (up to 12 hours of autonomous operation);
- Color: black/white.

5.3.3 Images

Fig. 35. All-in-One PC, mouse, keyboard

Fig. 36. Mouse pad

tyromotion

5.3.4 Notes

Also you can purchase with the computer mouse pad with non-slip base (dimensions: 320x270x2 mm) and aluminum holder for the monitor (dimensions: 48x20x14 cm, mass - 4.2 kg). The power cord for the computer (length – 3 m) is supplied if you purchase computer. To connect the cable into the socket it is needed socket type shuko.

5.4 Tymo X-Pad

5.4.1 Description

Patient may use Tymo X-Pad for more comfort during working on Tymo Therapy Plate (if necessary it may be put under elbows, arms, legs, or sit on it during therapy).

5.4.2 Technical specifications

- Dimensions (LxW): 50x30 cm, thickness – 3 cm;
- Mass: 0.35 kg;
- Material: coat is made of 100%-polypropylene; filter – elastic polyurethane foam, ST 2536.

5.4.3 Images

Fig. 37. Tymo X-Pad

Fig. 38. Tymo X-Pad

6 Delivery set

Description	Qty, note
Tymo Therapy Plate for testing and rehabilitation: 51x30 cm	1 pcs., basic completeness
Tymo battery	1 pcs., basic completeness. May be ordered as option.
TyroS Basic Software V. 2.0, on CD	1 pcs., basic completeness
Instruction Manual	1 pcs., basic completeness
Rolling element 1D	Optional. 1 pcs. is enough to work with a patient.
Rolling element 2D	Optional. 1 pcs. is enough to work with a patient.
Tyrostation electric height adjustable desk with motor	Optional. 1 pcs. is enough to work with a patient.
Tyrostation height adjustable stool	Optional. 1 pcs. is enough to work with a patient.
All-in-One PC	1 pcs., optional

tyromotion

Wireless keypad	1 pcs., optional
Wireless mouse	1 pcs., optional
Mouse pad	Optional
Monitor holder	1 pcs., optional
Bluetooth adapter	Optional
Charging Unit "Ansmann AC/DC Adaptor for Medical Applications"	Optional
PC power cable	1 pcs., optional
Tymo X-Pad: 50x30x3 cm	Optional
TyroS Professional Software 2.0, on CD or USB	Optional. 1 pcs. is enough to work with a patient.

7 Tymo Therapy Plate Maintenance

7.1 Cleaning the Tymo Therapy Plate

The following hygiene measures absolutely must be carried out:

- Disinfection of patient's hands before the start of the therapy;
- Desinfection of the surface of the Therapy Plate after the therapy;
- In the event of pronounced build-up of sweat, wiping the Therapy Plate dry after the therapy.

In addition, you should clean all Tymo Therapy Plate components at regular intervals. For this use a dry or slightly damp cloth. In the case of stubborn dirt, the cloth can be soaked in cleaning alcohol.

7.2 Repair Tymo Therapy Plate and its elements

- In case of incorrect operation of the Tymo Therapy Plate or if it does not respond to commands provided in the Instruction Manual, an operator should contact Tyromotion GmbH.
- Tymo Therapy Plate has a limited impact protection. If a device falls to the floor from a table height, sensors or its battery may be damaged. Defects may manifest themselves after the fall. Therefore stop using the Tymo Therapy Plate after its fall and contact Tyromotion GmbH or his distributor.
- If you encounter any fault or any questions, see the Chapter 7 of the Instruction Manual. If you do not find the answer, contact Tyromotion GmbH or its representatives in your country.

8 FAQ - Frequently Asked Questions

- **How do I read the user manual?**

You will need a program to read files in PDF, for example, Adobe Acrobat Reader. The English version of Adobe Acrobat Reader, run the file "Adobe \ AdbeRdr920_en_US.exe" on your disk Tymo Therapy Plate. Versions in other languages can be downloaded from the manufacturer's website.

- **How to change the language of the program?**

The program has a list of different languages in the form of flags of states in the top right corner of the screen. Clicking on one of these symbols, you can select the language you want.

- **How to change the name and logo of the institution?**

In the program, go to Settings. Enter the name of the institution in the respective field. To download a logo, specify the appropriate path to the file in the format of BMP, JPEG, JPG, or PNG.

- **What can I do if the therapeutic program does not start?**

Make sure that you have installed Adobe Flash Player 10 ActiveX controls in your system. Otherwise, run the file "Adobe \ install_flash_player_10_active_x.exe" on the disc.

- **I need a printer?**

Typically, you do not need printer for working, but if you want to print the report - it is indispensable.

- **How to view a report in PDF?**

You will need a program to read files in PDF, for example, Adobe Acrobat Reader. The English version of Adobe Acrobat Reader, run the file "Adobe \ AdbeRdr920_en_US.exe" on your hard drive. Versions in other languages can be downloaded from the manufacturer's website.

- **What does the red flashing of the LED display mean?**

If the LED display starts to flash, this is the instruction to charge the battery because the charge level is only sufficient to last for a short while (but at least an hour from the moment when the LED display starts to flash red).

- **Why will the Tymo Therapy Plate not connect with my PC?**

Make sure that your Bluetooth adapter is correctly connected to your PC and that you have started the therapy area of the Tymo Therapy Plate Software. If you have problems with the installation of your Bluetooth receiver, ask your system administrator for help and read the installation instructions provided for the Bluetooth receiver.

tyromotion

- **In the force measurement function, why does a strength value appear although no active force is being exerted upon the sensor plate?**

The sensors of Tymo Therapy Plate are possibly not correctly calibrated. Open the area "Hardware Configuration" within the Tymo Therapy Plate Software and select your Tymo Therapy Plate from the list. Place Tymo Therapy Plate on an even, stable surface and press on the button "Configure". After you have received the message "Configuration OK" leave the area by using the "Save" symbol.

-
- **In case of any other technical or other problems please contact Tyromotion GmbH or its representative in your country.**

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