02387555

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Zhuhai Vision Medical Technology Co., Ltd on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Sun Danfeng

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(Date: Jun. 11, 2024)

查询网址 Website for verifying the certificate: <http://www.cccpit.com/validate.html>

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
Zhuhai Vision Medical Technology Co., Ltd.

Product Information

- **Product Name:** Flexible Digital Uretero-Renoscope
- **Product Type:** UF25/ UF28/ UF30
- **Product Life:** 5 years or 1000 times peracetic acid disinfection
- **Manufacturer Name:** Zhuhai Vision Medical Technology Co., Ltd.
- **Address:** 301 Building 4, No. 30 Feipeng Road, GMTCM Park, Hengqin New Area, 519031 Zhuhai City, Guangdong Province, PEOPLE' S REPUBLIC OF CHINA
- **Contact information**
 - Tel: +86-756-6910150
 - Fax : +86-756-6910511
 - Zip code : 519031

Version Information

This manual version may be upgraded due to a change of software or technical specifications at any time without notice. The version information of this specification are as follows:

- Document number: TD/CE-UF-A-009
- Version number: V2.5
- Edit date: Jan. 2024

©2024, The copyright of Zhuhai Vision Medical Technology Co., Ltd

CE Mark

EC Representative Name:

Shanghai International Holding Corp. GmbH (Europe)

EC Representative Address:

Eiffestraße 80, 20537 Hamburg, Germany

Statement

The manufacturer owns the copyright of this non-publicly published manual and has the right to treat it as confidential information. This instruction manual is only used as a reference material for operation, maintenance and repair of the product. Others have no right to disclose this instruction manual to others.

This product manual contains proprietary information protected by copyright law; all rights reserved. Any part of this manual may not be photocopied, copied or translated into other languages without the written consent of the manufacturer.

The content contained in this product manual can be changed without notice.

Responsibilities of the Manufacturer

In the following conditions are met, the manufacturer will responsible for the product safety, reliability and performance:

- The installation, maintenance or upgrade of the product shall be carried out by personnel authorized or approved by the manufacturer;
- The storage environment, working environment, and electrical environment of the product conform to product specifications;
- Use this product in accordance with the "Instruction Manual";
- Use this product within the period of use.

Overview

◆ Main Structure

The Flexible Digital Uretero-Renoscope is composed of distal end, insertion part, operation section.

◆ Basic performance of the product

Used in conjunction with video processor product by our company, it can output real-time images without unexpected image orientation

◆ Intended Use

This product is used together with the Video Endoscope Processor produced by Zhuhai Vision Medical Technology Co., Ltd. to provide images for diagnosis and treatment of ureter and renal pelvis for adults.

◆ Contraindication

- 1) Uncontrolled systemic hemorrhagic diseases;
- 2) physically weak and unable to tolerate anesthesia;
- 3) Patients with prostatic hyperplasia affecting ureteroscopic insertion;
- 4) Patients with lower urinary tract tumors;
- 5) Patients with bladder contracture;
- 6) Systemic hemorrhagic diseases that can be controlled;
- 7) uncontrolled urinary tract infection;
- 8) Severe cardiopulmonary insufficiency, unable to tolerate surgery;
- 9) Severe urethral stricture and ureteral stricture;
- 10) Severe hip malformation and difficult lithotomy position;
- 11) Severe gross hematuria.
- 12) The device is not intended for use in pregnant women or children.

◆ Complication

Hematuria; dysuria; lower abdominal pain; flank pain-colic; urethral tear or perforation; damage to the ureter orifice, ureter and-or kidney; ureteral tear; ureteral avulsion or perforation; cystitis or urinary tract infection (UTI) including pyelonephritis or urosepsis; ureteral stricture (the most common); and urethral stricture.

◆ Cautions, Warnings and Notes

- 1) This product does not contain any parts that can be repaired by the user. The internal maintenance of the instrument must be carried out by the technical maintenance personnel authorized by the manufacturer.
- 2) This product shall not be used in the home.
- 3) The connection with the video endoscope processor must be disconnected before cleaning, disinfection .

- 4) Operate in strict accordance with the instructions provided by the manufacturer.
- 5) This product cannot be used with high frequency equipment.
- 6) The equipment can not be used in MRI scan environment.
- 7) The power supply system for this product depends on the Video Endoscope Processor you selected. And the detail operation description please refer to the user manual of the Processor.

About the Manual

This manual introduces the use, function and operation method of the product in detail. Before using this product, please read and understand the contents of this manual carefully to ensure the correct use of this product and the safety of patients and users.

This manual introduces the product according to the most complete configuration, so some of the content may not be applicable to the product you purchased. If you have any questions, please contact our company.

Please place this manual near the product so that you can get it conveniently and in time when you need it.

This manual is suitable for professional clinical medical staff to read. Readers should have knowledge and work experience in medical procedures, practice and terminology necessary for endoscopy.

All illustrations provided in this manual are for reference only, and the settings or data in the illustrations may not be exactly the same as the actual display you see on the product.

In the following description in this manual, the Flexible Digital Uretero-Renoscope will be referred to as "endoscope" for short.

About the product

The main material of endoscope insertion part is TPU. The ingredients and raw materials of ME equipment comply with local laws and regulations. Insertion Section、 distal end、 bending section is an APPLIED PART, see details Fig 2.1 UF series constitute.

Every two years to verify the basic safety of the product, detailed details, please consult the manufacturer.

Catalogue

SYMBOLS	1
KEY POINTS (PLEASE READ BEFORE USE)	2
USER QUALIFICATION	2
SIGNAL WORDS	2
MAINTENANCE	3
INAPPROPRIATE MAINTENANCE AND MODIFICATION PROHIBITED	3
WARNINGS AND CAUTIONS	3
ENDOSCOPIC IMAGES DISAPPEAR OR FREEZE PRECAUTIONS	7
EXAMPLES OF INAPPROPRIATE HANDLING	8
CHAPTER1 CHECKING THE PACKAGE CONTENTS	9
CHAPTER 2 INSTRUMENT NOMENCLATURE AND SPECIFICATIONS	11
2.1 CONSTITUTE	11
2.2 ENDOSCOPE FUNCTIONS	11
2.3 SPECIFICATIONS	12
CHAPTER 3 PREPARATION AND INSPECTION	15
3.1 PREPARATION OF THE EQUIPMENT	15
3.2 INSPECTION OF VIDEO ENDOSCOPE	16
3.3 INSPECTION OF ACCESSORIES	18
3.4 ACCESSORIES INSTALLATION TO THE ENDOSCOPE	19
3.5 INSPECTION OF THE WORKING CHANNEL	20
3.6 EQUIPMENT CONNECTIONS	21
3.7 INSPECTION OF DISTAL TIP	22
3.8 IMAGE CHECK	23
3.9 CHECK THE FREEZING FUNCTION	24
CHAPTER 4 OPERATION	25
4.1 INSERT THE ENDOSCOPE	26
4.3 WITHDRAW THE ENDOSCOPE	29
CHAPTER 5 CLEANING AND DISINFECTION : GENERAL POLICY	30
5.1 KEY POINTS OF CLEANING AND DISINFECTION	30
5.2 PRECAUTION	30
CHAPTER 6 SUITABLE METHODS OF CLEANING, DISINFECTION AND CHEMICAL PREPARATION	32
6.1 COMPATIBILITY SUMMARY	32
6.2 CLEANING REAGENT	33
6.3DISINFECTION SOLUTION	33
6.4 THE INFORMATION OF WASHER	34
6.5 RINSING WATER	34
CHAPTER 7 CLEANING AND DISINFECTION PROCEDURES	35
7.1 MAINTAIN	35
7.2 REQUIRED CLEANING AND DISINFECTION EQUIPMENT	35
7.3 CLEANING AND DISINFECTION PROCEDURES	39
7.4 PRE-CLEANING	41
7.5 TRANSPORTATION OF THE ENDOSCOPE	42
7.6 LEAKAGE TESTING	43
7.7 MANUAL CLEANING	46
7.8 RINSING AND DRYING AFTER MANUAL CLEANING	50
7.9DISINFECTION	51
7.10 RINSING AND DRYING AFTER DISINFECTION	52
CHAPTER 8 STORAGE AND DISPOSAL	56

8.1 STORAGE OF THE ENDOSCOPE	56
8.2 STORAGE OF REUSABLE PARTS AND REPROCESSING EQUIPMENT	56
8.3 TRANSPORTING OUTSIDE THE HOSPITAL	57
8.4 DISPOSAL	57
CHAPTER 9 TROUBLESHOOTING	58
9.1 TROUBLESHOOTING GUIDE	58
9.2 WITHDRAWAL OF THE ENDOSCOPE WITH ANY ABNORMALITY	59
9.3 RETURNING THE ENDOSCOPE FOR REPAIR	59
CHAPTER 10 ELECTROMAGNETIC COMPATIBILITY	61

Symbols

The meanings of the symbols shown on the package with the components, the back cover of this instruction manual and/or this instrument are as follows:

Symbols	The meaning of the symbol	Symbols	The meaning of the symbol
	Follow instructions for use		Medical Device
	Type BF Applied Part		Caution In case of application as a SAFETY SIGN, the rules according to ISO 3864-1 are to be adhered to.
	Serial Number		This way up
	Fragile; handle with care		Keep dry
	Manufacturer		Date of manufacture
IPX7	Protected against the effects of temporary immersion in water		Stacking limit by number
	Symbol for Waste Electrical and Electronic Equipment (Please follow Local Ordinances or Regulations of disposal.)		Authorized Representative in the European Community
	CE mark.		

Key points (Please read before use)

User Qualification

If there are official standards for use established by the medical administration or other organizations, please follow these standards. If there is no official standard, the operator of this product must be a doctor who has passed the audit of the medical safety director of the hospital or a doctor who can safely carry out the established operation according to the official regulations designated by the director of the using department according to the difficulty of using the technology.

Instrument compatibility

Refer to Chapter I "Checking the Package Contents", confirm that this instrument is compatible with the ancillary equipment being used. Using incompatible equipment can result in patient or operator injury and/or equipment damage.

- There is no guarantee that instruments selected solely using maximum insertion portion width and working length will be compatible in combination.
- There is no guarantee that instruments selected solely using this minimum instrument channel width will be compatible in combination.

This product conforms to the EMC standard IEC 60601-1-2.

Signal Words

The following signal words are used throughout this manual:

Indicates an emergency hazardous situation which, if not avoided, could result in death or serious injury.

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to warn against unsafe operations or potential equipment damage

Indicates additional helpful information.

Maintenance

With increasing number of uses and / or total use time, the probability of endoscopic and adjuvant failure increases accordingly. In addition to the inspection before each operation, the person in charge of medical equipment maintenance in each hospital shall regularly check the items specified in this manual. Endoscopes suspected of abnormality cannot be used and should be checked as per Section 9.1 Troubleshooting Guidelines described below. If there is still abnormal after the examination, please contact Zhuhai Vision Medical Technology Co., LTD.

Inappropriate maintenance and modification prohibited

This product does not contain any user-serviceable parts. Do not disassemble, modify or attempt to repair, otherwise it may result in injury to the patient or damage to the equipment. If necessary, please contact the manufacturer.

Zhuhai Vision Medical Technology Co., Ltd. is not responsible for the injury or damage caused by the maintenance equipment of personnel not from Zhuhai Vision Medical Technology Co., Ltd.

Warnings and Cautions

Follow the warnings and cautions given below when handling this instrument. This information is to be supplemented by the warnings and cautions given in each chapter.

- Video Endoscope should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the Video Endoscope should be observed to verify normal operation in the configuration in which it will be used.
- Interconnections conditions require the applied parts of other ME equipment used within the configuration for endoscopic application to be type BF applied parts or type CF applied parts
- Class A equipment is intended to be used in industrial environment. The Video Endoscope may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances.
- Except for the cables sold as spare parts of internal components by the manufacturer of the Endoscope, the use of accessories and cables outside the regulations may result in an increase in the emission of the electronic ureteral pyeloscopy or a decrease in immunity.

- After using this instrument, clean, disinfect and store it according to the instructions given in Chapter 5 to Chapter 8. Using improperly or incompletely reprocessed or stored instruments may cause patient cross-contamination and/or infection.
- When this product is used together with the laser equipment, in order to avoid the damage to the eyes, Take the corresponding protective measures, according to the classification of the laser equipment, such as by wearing suitable anti-filter glasses, or with suitable filters on the eyepiece of the live endoscope.
- Modification of this equipment is not permitted.
- Energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive.
- When this product is used together with the laser equipment, in order to avoid the damage to the eyes, Take the corresponding protective measures, according to the classification of the laser equipment, such as by wearing suitable anti-filter glasses, or with suitable filters on the eyepiece of the live endoscope.
- Be sure to prepare standby equipment or equipment with similar functions to avoid accidents in case of equipment failure or abnormal function during operation.
- When using the endoscope, the working part of the endoscope should be coated with lubricant or the instrument channel should be injected with different liquids. After using the endoscope, the surface of the lens body and the instrument channel should be wiped and scrubbed according to the instructions. Otherwise, rapid aging of the mirror may occur.
- Do not use the endoscope in the presence of inflammable substances or in an oxygen-rich environment. Otherwise, fire or explosion may occur.
- cooperate with the manufacturer's equipment for use, the use of incompatible equipment will lead to patient or operator injury or equipment damage
- Instructions warning: Devices should be placed reliably to avoid falling.
- Heavy objects cannot be placed on the device, otherwise it will cause device damage or injury to the patient
- Cooperate with the manufacturer supporting endoscope image processor

use, using non-manufacturer supporting processor, may lead to patient casualties

- Do not touch the rubber and insert the tube with sharp objects, otherwise the equipment will be damaged
- When this product is used with laser equipment, avoid potential eye damage to the operator and protect it in advance.
- When this product is used with the laser, the laser will reflect back and cause an abnormal temperature rise at the head end. Please note that the head end is not too close to the reflection source when operating.
- Do not use a high frequency electric knife for electric burning surgery and argon plasma coagulation when using this product, because this product cannot be used together with a high frequency electric knife and argon plasma coagulation.
- When this product is used with laser equipment, operators, maintenance personnel or personnel who need to look directly or observe for a long time should wear laser protective goggles with appropriate protection strength according to the laser equipment and power used.
- In the process of examination or surgery, the doctor will match other auxiliary instruments according to the actual situation of the patient. Usually, the instruments used together with this product are stone basket, biopsy forceps and laser fiber. Since the channel of endoscope is 1.2mm, it is recommended that the outer diameter of the stone basket, biopsy forceps and laser fiber used together with this device should be less than or equal to 1mm. And the power of the laser during the operation should be controlled at 15W to 20W.
- Do not hit, bump, or drop the tip, insertion, bending, operation, or display of the endoscope. At the same time, do not forcefully bend, pull or twist the tip, insertion, bending, operation, or display. Otherwise, the endoscope will be damaged. And it may cause the endoscope parts to fall off into the patient's body.
- Do not perform bending operations forcibly or suddenly. Do not forcibly pull, twist or rotate angled bends. Otherwise, the patient may be injured. It will also result in the inability to straighten the bend during the inspection.
- Do not operate the endoscope in the environment with excessive

www.goszdravnadzor.ru

pressure difference. Otherwise, the bending angle may be affected.

- If the Angle fails during the operation, do not be nervous, control the Angle handle in the middle 0 Angle position, slowly pull out the operation department, and never pull it out suddenly, otherwise the patient will be injured.
 - The endoscope uses LED light source lighting, there is the possibility of failure, it is recommended to prepare the same or similar standby equipment before surgery. If there is no light source in the process of use, do not be nervous, control the Angle handle to 0, slowly pull out the operation department, prohibit suddenly pull out, otherwise it will lead to patient injury.
 - If the rubber is found or side leakage, please stop using it, otherwise it will be damaged to the patient.
 - Instructions: This product does not contain any parts that can be repaired by users. Do not remove, modify, or attempt to repair patients that may cause injury or equipment damage
 - When the image is frozen, do not bend the operation part, insert or extract the endoscope insertion part, otherwise it will hurt the patient
 - When the bend is at an angle, do not insert or pull out the insertion part of the endoscope. Otherwise, the patient may be injured.
 - The bend can only bend upwards or downwards. The bending direction of the endoscope should be fully considered when inserting or withdrawing and operating the endoscope. When inserting or withdrawing the endoscope, do not apply excessive force to the right or left. Otherwise, the patient may be injured.
 - When a clear endoscopic image is not obtained, do not manipulate the bending part, insert or pull out the insertion part of the endoscope. Otherwise, the patient may be injured.
 - Do not insert or pull out the insertion part suddenly or forcefully. Otherwise, the patient may be injured.
 - When it is difficult to insert the endoscope, do not insert it forcibly, but stop the operation. Forcible insertion can cause injury to the patient.
 - The following precautions should be followed when lighting.
- Do not keep close observation, or keep the tip of the endoscope in

contact with active tissues for a long time.

The user shall install and use according to the EMC information provided in the attached documents.

- Portable and mobile RF communication equipment may affect the Endoscope. Avoid strong electromagnetic interference when using it, such as close to mobile phones, microwave ovens, etc.;
- Guidance and manufacturer's declaration stated in the appendix.
- Do not make the coiled diameter of the insertion tube less than 12 cm, otherwise, it will cause damage to the equipment.
- Do not try to bend the insertion part of the endoscope forcibly. Otherwise, the insertion part may be damaged.
- Do not bump the tip of the insertion part, especially the surface of the objective lens at the tip, otherwise it will cause framing failure.
- Do not twist or bend the bent part with your hands, otherwise the equipment will be damaged.
- Do not squeeze the bending section forcefully. The covering of the bending section may stretch or break and cause water leaks.
- Electromagnetic interference may occur on this instrument near equipment marked with the following symbol or other portable and mobile RF (Radio Frequency) communications equipment such as cellular phones. If electromagnetic interference occurs, mitigation measures may be necessary, such as reorienting or relocating this instrument, or shielding the location. 
- In order to check the electromagnetic interference from other devices (any devices like X-ray machine, CT, MRT , etc or equipment components of other system), the system should be observed in the operating environment whether it is working properly.

Endoscopic images disappear or freeze precautions

Warning

- If the endoscopic image disappears unexpectedly, or the frozen image cannot be restored to the real-time display working state, you must stop using it immediately and withdraw the endoscope from the patient as described in Section 9.2 "Extracting the Abnormal Endoscope". In this case, continuing to use the operating curve to insert or withdraw the

endoscope may cause injury to the patient.

- Please follow these precautions when operating these instruments. Otherwise, the endoscope image may disappear unexpectedly, or in the inspection process frozen image can't be restored.
- Do not bend, strike or twist the insertion part and the operation section. It may damage the endoscope, causing leakage and / or damage to internal components, such as CMOS, or circuitry.
- If there are continuous air bubbles coming from the endoscope during leak test. Do not use the endoscope. Otherwise, water may enter the instrument and cause a short circuit. This will cause damage to the CMOS and circuit.

Examples of inappropriate handling

Details on clinical endoscopic technique are the responsibility of trained specialists. Patient safety in endoscopic examinations and endoscopic treatment can be ensured through appropriate handling by the physician and the medical facility. Examples of inappropriate handling are given below.

Warning

- Inserting, withdrawing without a clear endoscopic image may cause patient injury.
- Inserting or withdrawing the endoscope, applying suction or operating the bending section without a clear endoscopic image may cause patient injury.

LED light emitting diode

LED is built in the device, which is basically safe for eyes. But looking directly at the light emitted by the LED for a long time may damage the eyes. Please read this manual carefully.

Warning

- If the control or adjustment device cannot be used according to the instructions, or the operation of each step can be carried out, harmful radiation exposure may be caused.




Note


- Avoid direct or scattered light exposure to eyes or skin.
- The endoscope uses LED light source to illuminate, and there is a possibility of failure. It is recommended to prepare the same or similar spare equipment before the operation.

Chapter1 Checking the package contents

The product is only intended to use with the Video Endoscope Processor manufactured by Vision. Match all items in the package with the components shown below. Inspect each item for damage. If the instrument is damaged, a component is missing or you have any questions, do not use the instrument; immediately contact our company. **This instrument was not disinfected before shipment.**

Before using this instrument for the first time, clean, disinfect it according to the instructions given in Chapter 5 “Cleaning and Disinfection : General Policy” to Chapter 7 “Cleaning and Disinfection Procedures”.

Items	Picture	function
Video Endoscope		When used with Video Endoscope Processor, it can be used for clinical examination, diagnosis and adjuvant therapy.
Leakage detector		Connect the operating part to test the tightness of the instrument.
Channel cleaning brush		Brushing the working channel of endoscope
Channel port cleaning brush		Brushing of injection port
Three-Way stopcock		Connect the injection port for instruments and liquid to enter the endoscope
User manual		Follow instructions for use
Suction cleaning adapter		Suitable for cleaning and disinfection connection
Sealing plug		Sealed injection port.

1.waterproof cover 2.balanced valve		1. Connect the waterproof lid to the handle during cleaning, soaking, and gas sterilization. 2. Install the balance valve on the leakage detector before shipment to balance the air pressure inside and outside the mirror during transportation. The balance valve must be removed before the patient examination and immersion into the cleaning fluid.
---	---	--

Form 1.1

Chapter 2 Instrument Nomenclature and Specifications

The constitute of the cooperative-used Video Endoscope Processor please refer to the user manual of the processor.

2.1 Constitute

Distal end, bending section and Insertion Section is an APPLIED PART, see details Fig 2.1 UF series constitute.


Fig 2.1 UF series constitute

2.2 Endoscope Functions

No.	Part Name	Function
1	Dimmer push button	It is used to adjust the brightness of the tip LED, which can be divided into 5 levels of dimming.
2	Photo/video recording button	1. Press this button once in the power on preview mode to record the static picture.

No.	Part Name	Function
		2. Press and hold this button for 2 seconds in the power on preview mode to enter the video mode. 3. In video mode, long press this button for 3 seconds to exit video mode and enter real-time video.
3	Injection port	The injection port is a standard Ruhr connector, from which medical instruments can enter during treatment.
4	Up / Down angulations' control knob	Turn the control knob to 'U' direction to make the distal tip bend up; Turn the control knob to 'D' direction to make the distal tip bend down.
5	Operation section	Handle and operating Insertion tube.
6	Conical sheath	Protect the Insertion tube.
7	Insertion part	It extends into the natural cavity of human body.
8	Bending section	Adjust the position of distal tip via UP/DOWN Angulations Control Knob.
9	Distal end	Endoscope lens end with LED light.

Form 2.1

2.3 Specifications

Note

- This product is equipped with an image sensor chip on the inner lens side to obtain images, convert them into electrical signals and then transmit them to the image and video processing circuit by cables, and display real-time images on the LCD screen. No optical fiber bundles are used.

Environment specifications

Normal working conditions	Environment temperature:	5°C~35°C
	Atmospheric pressure:	700hPa~1060hPa
	Relative humidity:	30%~85% (Non-Condensing)
Transportation and storage	Environment temperature:	-10°C~+55°C
	Atmospheric pressure:	500hPa~1060hPa
	Relative humidity:	≤93% (Non-Condensing)

Form 2.3

Safety Standard


Note

- In accordance with the CE standard classification, this product is class IIa.
This product is classified as follows according to IEC60601-1:

Classification of protection against electric shock	TYPE BF APPLIED PART
Classification of protection against hazards of explosion	Not suitable for use in places with flammable anesthetic gas
Classification of protection against ingress of liquid	IPX7
Classification of running mode	Continuous operation equipment
Mobile level	Handheld

Form 2.4

Product specification

Product specification				
Model		UF30	UF28	UF25
Parameter				
Optical system	Field of view	120°±10%		
	Direction of view	0 ° limit deviation not more than 10 °		
	Depth of field	3mm~150mm		
Insertion tube outer	Diameter	Ø3.0mm, allowed: +10%, lower limit not	Ø2.8mm, allowed: +10%, lower limit not	Ø2.5mm, allowed: +10%, lower limit not
	Distal end enlarge picture	 1. CMOS camera lens 2. LED light 3. Suction Channel exit		
	Working length	680mm ± 10%		
Working channel	Minimum suction channel inner diameter	≥Ø1.2mm		
	Direction from which endo-therapy accessories enter and exit the endoscopic image			

Channel Opening Type		
Bending section	Angulation range	Up275°, Down275°

Form 2.5

System connection diagram


Fig 2.2 System connection diagram

Warning: Additional equipment connected to the medical electrical equipment must comply with the respective IEC or ISO standards (e.g. IEC60950 for data processing equipment). Furthermore all configurations shall comply with the requirements for medical electrical systems (see clause 16 of IEC60601-1 edition 3.2). Anybody connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. When more than one device is connected to use, the accumulation of leakage current may cause the danger of security. If in doubt, consult your local representative or the technical service department.

Chapter 3 Preparation and Inspection


Before each case, prepare and inspect this instrument as instructed below. Inspect other equipment to be used with this instrument as instructed in their respective instruction manuals. If the irregularities are suspected after inspection, follow the instructions given in Chapter 9, "Troubleshooting". If this instrument still has malfunctions, do not use it. Return it to our company for repair as described in Section 9.3, "Returning the endoscope for repair".



- Before endoscopy, it is necessary to check whether the performance of the endoscope is safe to use, such as the display screen can clearly receive real-time images, and can take and record images in real-time.
- The use of a malfunctioning endoscope will endanger the safety of the patient and result in aggravated equipment damage.
- This product has not been cleaned and disinfected before shipment. Before using for the first time, reprocess it according to the instructions given in Chapter 5, "Cleaning and Disinfection : General Policy" to Chapter 7, "Cleaning and Disinfection Procedures" instructions for cleaning and disinfection.

3.1 Preparation of the equipment

Before use, please prepare the equipment and personal protective equipment listed in the table below, such as eye wear, face mask, moisture-resistant clothing and chemical-resistant gloves. Personal protective equipment should be of the right size and long enough to prevent the skin from being exposed. Please refer to the instruction manual of the respective device.

Before the first use, cleaning and disinfection should be carried out in accordance with the instructions in Chapter 5 "Cleaning and Disinfection : General Policy" to Chapter 7 "Cleaning and Disinfection Procedures".

Items	Picture	function
Video Endoscope		When used with Video Endoscope Processor, it can be used for clinical examination, diagnosis and adjuvant therapy.
Leakage detector		Connect the operating part to test the tightness of the instrument.

Channel cleaning brush		Brushing the working channel of endoscope
Channel port cleaning brush		Brushing of injection port
Three-Way stopcock		Connect the injection port for instruments and liquid to enter the endoscope
User manual		Follow the operating instructions
Suction cleaning adapter		Suitable for cleaning and disinfection connection
Sealing plug		Sealed injection port.
1.waterproof cover 2.balanced valve		1. Connect the waterproof lid to the handle during cleaning, soaking, and gas sterilization. 2. Install the balance valve on the leakage detector before shipment to balance the air pressure inside and outside the mirror during transportation. The balance valve must be removed before the patient examination and immersion into the cleaning fluid.

Form 3.1

3.2 Inspection of Video Endoscope

Clean and disinfect as described in Chapter 5 "Cleaning and Disinfection : General Policy" to Chapter 7 "Cleaning and Disinfection Procedures".

Inspection of Video Endoscope

1. Inspect the operation section and the monitor for excessive scratching, deformation, loose parts or other irregularities.
2. Inspect the protective boot and the insertion tube near the protective boot for bends, twists or other irregularities.
3. Inspect the external surface of the entire insertion tube including the bending section and the distal tip for dents, bulges, swelling, scratching, holes, sagging, transformation, bends, and adhesion of foreign bodies, dropout of parts, any protruding objects or other irregularities.
4. Holding the operation section with one hand, the entire surface of the insertion part and the light stroking in two directions with the other hand. (See fig3.4). Confirm no objects or wire from the insert protrusion inserted into the pipe and make sure that no abnormal hard.


Fig 3.4

5. Inspect the covering of the bending section for sagging, swelling, cutting, holes, or other irregularities.
6. Using both hands, bend the insertion tube of the endoscope into a semicircle. Then, moving your hands as shown by the arrows, confirm that the entire insertion tube can be smoothly bent to a semicircle. (Fig 3.5)


Fig 3.5

7. Gently hold the midpoint of the bending section and a point 10 cm from the distal tip. Push and pull

gently to confirm that the border between the bending section and the insertion tube is not loose.

8. Inspect the distal tip of the endoscope's insertion tube for scratching, cracks, stains, gaps around the lens or other irregularities.

Inspection of the bending functions

Perform the following inspections while the bending section is straight.

- If Angle control handle movement is too loose or not smooth, or the bending section does not angulate smoothly, the bending mechanism may be abnormal. In this case, do not use the endoscope because it may be impossible to straighten the bending section during an examination.

1. Slowly rotate the up/down angle control handle to the end points in each direction, and then return to their natural positions. Confirm that the bending operation of the bending part is smooth and normal, can reach the maximum amplitude, and return to the natural position.

2. As shown in Figure 3.6 when moving the UP/DOWN angle control handle to the natural position, confirm that the curved part can be smoothly restored to a nearly straight condition.


Fig 3.6

3.3 Inspection of accessories

Inspection of Three-Way Stopcock



Fig 3.7

- Inspect and make sure the Three-Way Stopcock is without tear, severe deformation, or damage.
- The medical three-way valve with abnormal function will cause leakage and bring the risk of infection.

Note

- If the three-way valve is abnormal, please replace it with a new one.

3.4 Accessories Installation to the endoscope

As shown in the figure below, align the medical three-way valve with the injection port connector and tighten it clockwise.

Insert the endoscope attachment through the working channel port of the three-way lock valve, and push it in slowly to ensure that the attachment extends smoothly from the endoscope tip.

- When inserting the endoscope attachment, always
- select the one of the correct size.
- The endoscope accessories should be checked before use, and if there is any abnormality in the appearance or operation, they must be replaced.
- When inserting the endoscopic attachment, it shall confirm that the endoscopic attachment is in a closed state. Otherwise, it will cause working channel damage and parts loss.
- If it is difficult to insert the endoscopic accessories, do not insert them forcefully. Otherwise, damage to the working channel or shedding of the endoscopic accessory parts will result.

- The attachment should be installed by trained medical personnel. Otherwise, improper installation will affect the sealing effect.


Fig 3.8

3.5 Inspection of the working channel

1. Install the three-way lock valve to the working channel port, connect to the water injection channel interface of the three-way lock valve with the syringe or the water injection pipe joint, and inject the normal saline through the water injection channel.
2. Make sure that water is emitted from the distal tip of the insertion tube.
3. Insert the cleaning brush through the injection port. Make sure the cleaning brush extends smoothly from the distal tip. At the same time, make sure that there is no foreign article at the distal tip.
4. Make sure that the cleaning brush can be pulled out smoothly from the injection port.

Note

- The syringe must be inserted fully and held perpendicular to the valve for proper operation. Angled or incomplete insertion may result in fluid leakage.
- If fluid is not emitted from the endoscope's distal tip, flush air through the channel.
- Keep eyes away from the distal tip when inserting the cleaning brush. A cleaning brush reaches out from the distal tip of endoscope may cause eye injury.
- That before each use, the outer surface of the portions of any energized endotherapy devices which are intended to be inserted into a patient checked to ensure there are no unintended rough surfaces, sharp edges or protrusions which may cause harm.

Caution



Fig 3.9

3.6 Equipment Connections

The connection mode and display setting of this product depends on the Video Endoscope Processor (Type BF Applied Part) you selected. And the detail operation description please refer to the user manual of the Processor. The product can choose the wired connection mode and the wireless connection mode.

Wireless Connection Mode

1. Remove the wireless transmitter from the charging card slot of the electronic endoscope image processor.
2. Press the black mark point on the right side of the wireless transmitter white balance button at the black mark point at the interface of the endoscope connector. Hearing the click, the wireless transmitter connects to the endoscope operator successfully, as shown in the figure below

Wired connection mode

3. Insert one end of the extension cable (Lemo plug) into the video endoscope processor output interface, the Lemo plug (mother plug). Clicks were heard after insertion.

2. Place the extension cable at the other end (endoscopic manipulation handle connection end) from top to bottom above the operant. Align the black dots on the to turn axis to the black dots above the maneuver (endoscopic maneuver handle connection end) rotate from right to left (endoscopic maneuver handle connection end) by gently clasping the (endoscopic maneuver handle connection end) to the maneuver, and after hearing a "click", the illustration follows:


Fig 3.1

3.7 Inspection of distal tip

- Do not stare directly into the light beam. This may result in eye injury.
- This product is of LED lighting, LED inside the equipment, not replaceable.
- If the LED light does not work, do not replace it by yourself, please contact Zhuhai Vision Medical Technology Co., LTD.
- Do not stare directly into the light beam. This may result in eye injury.

1. Connect the processor to the handle according to the guidance in the processor's manual. The power indicator lights, and endoscope image appears on the LCD screen.

2. Confirm that Examination light is emitted from the LCD light on distal tip of the endoscope. (Fig 3.2)


Fig 3.2

3. Point the endoscope's distal tip to a suitable object, and move between 1cm and 6 cm. Confirm that the brightness of the image on the LCD screen remains constant. (Fig 3.3)


Fig 3.3

Note

- If you feel abnormal endoscopic image color, refer to "White Balance Adjustment."

4. Point the distal tip of endoscope to your palm to observe, confirm that endoscopic image is free from noise, blur, fog or other irregularities.

3.8 Image check

1. Entering the preview state by default, that is, the real-time image display interface of the electronic endoscope image processor is implemented. At this time, the real-time image should be observed normally, and the image is clear, without blur and incomplete image abnormalities.

2. Hold the mirror with your left hand and put the thumb on the Angle lever, as shown in the Fig3.6

. At this point, the image is seen in the same direction as the operator . If there is any abnormality, please stop using it immediately and contact the manufacturer.

- Do not use this video endoscope if no real-time endoscopic image is observed. Otherwise, it may lead to patient injury.

3.9 Check the freezing function

- Do not use this video endoscope if no real-time endoscopic image is observed. Otherwise, it may lead to patient injury.

3.9.1 Under the main menu, press the "function switch" button. Enter the function switch setting interface, and select "freeze" below (as shown in Fig 3.13). Then return to the preview interface.


Fig3.13

- Press the light button switch. Confirm the real-time endoscopic image freezing. A frozen icon appears in the display interface, as shown in the figure below.

- Press the light button switch again to confirm that the frozen image returns to the real-time image. And the freeze icon disappears.

Chapter 4 Operation

This instruction manual does not explain or discuss any clinical endoscopic surgery. Only introduce the basic operations and precautions related to the use of this endoscope. Therefore, the operator of this endoscope must be a physician or medical staff under the supervision of a physician, and they must have received sufficient training in clinical endoscopy techniques.

Warning

- That before each use, the compatibility of the endoscopic equipment with any accessories and/or energized endotherapy devices should be checked according to any criteria for safe use defined in the instructions for use.
- In the process of examination or surgery, the doctor will match other auxiliary instruments according to the actual situation of the patient. Usually, the instruments used together with this product are stone basket, biopsy forceps and laser fiber. Since the channel of endoscope is 1.2mm, it is recommended that the outer diameter of the stone basket, biopsy forceps and laser fiber used together with this device should be less than or equal to 1mm. And the power of the laser during the operation should be controlled at 15W to 20W.
- Distal end surface and the light emission part temperature exceeds 41 ° C, temperatures over 41° C (105.8° F) may cause mucosal burns. Always use the minimum time and suitable distance necessarily for adequate viewing. Whenever possible, avoid close stationary viewing and do not leave the distal tip of the endoscope close to the mucous membrane for a long time.
- Never insert or withdraw the endoscope under any of the following conditions. Otherwise, the patient can result in injury.
 - When the accessory protrudes from the tip of the endoscope
 - Insertion or withdrawal while the endoscope is bending.
 - Inserting or pulling out the endoscope with excessive force
- The bending part will only bend up or down, but will not bend to the left or right. Take full consideration of the bending direction when doing insertion or withdrawal. Do not insert or withdraw with an excessive force in the left and right direction. Otherwise, the patient will be injured.
- If any of the following phenomena occurs during an examination, immediately stop the examination and withdraw the endoscope from

the patient as described in Section 9.2“Withdrawal of the Endoscope with any Abnormality”.

——The function of the endoscope is found to be abnormal.

——The endoscopic image on the LCD monitor disappears unexpectedly or freezes.

——The angle control handle is locked.

——The angle control function is abnormal.

Continuous use of endoscope under the situation stated above may cause injury to patients.

- If an abnormal endoscopic image/function occurs and returns to normal condition by itself, the endoscope may have malfunctioned. Continuous use of such an endoscope may cause repetition of the abnormality, which may not return to normal. In this case, stop the procedure immediately and slowly withdraw the endoscope while viewing the endoscopic image. Otherwise, patient can result in injury.
- That before each use or after a change of viewing modes/settings, the operator should check to ensure the view observed through the endoscope provides a live image (rather than a stored one) and has the correct image orientation. Image inspection requires the Video Endoscope Processor, please refer to the Video Endoscope Processor user manual.

4.1 Insert the endoscope

Holding and manipulating the endoscope

The operation section of the endoscope is designed to be held in the left hand. The UP/DOWN angulation control knob can be operated using the left thumb. The right hand is free to operate the insertion tube. (Fig 4.1)


Fig 4.1

Caution

- When inserting the endoscope, first lubricate the insertion tube with medical grade lubricant to ensure the lowest friction when inserting the endoscope into the patient.

1. Enter the insertion section of endoscope into the patient through the adult urethra. When entering the patient, be sure to observe the real-time image displayed on the monitor.
2. By operating the bending control pull rod to control the bending angle of the head end, guide the head end into the anterior and posterior urethra and renal pelvis for diagnosis and treatment.
3. That before each use, the outer surface of the portions of endoscope which are intended to be inserted into a **PATIENT** checked to ensure there are no unintended rough surfaces, sharp edges or protrusions which may cause harm.

Insertion of endoscopic accessories

Warning

- Do not force or suddenly insert the endoscopic accessories. Otherwise, the endoscopic accessories will suddenly protrude from the end of the inner lens, causing injury to the patient.
- When using the endoscope accessory, the curved part of the endoscope should be in a straightened state and the distance between the end of the inner lens and the mucous membrane should be greater than the minimum visible distance of the endoscope.
- Don't operate the accessories if you can't see the accessories in the endoscopic image. Otherwise, the instrument will be damaged.

1. Please refer to the contents of this manual to select the accessories that match the endoscope.
2. First, the endoscope is retracted to the bending section in a natural straight state, and the bending control

knob is not moved.

3. Confirm the status of the endoscopic appendage, and then slowly insert the endoscopic appendage into the endoscopic biopsy valve port or Ruhr connector.


4. While observing the real-time image displayed on the monitor, insert the endoscopic accessory slowly into the working channel and fix the endoscopic accessory.

5. Operate the bending control knob to move the bending section slowly to the observed position for operation.

4.2 Using laser system

This product can be used together with the HO: YAG laser system.

Please refer to the instructions for use of the laser system, and follow the following precautions,

Warning

- Do not use laser systems other than the HO: YAG. Otherwise, the endoscopic image on the monitor cannot be observed, resulting in injury to the patient.
- Do not emit a laser without observing the laser probe tip and the target beam in the endoscopic image. Otherwise, patient injury and forceps pipe damage and water leakage.
- Do not draw the laser probe while emitting the laser; cause patient injury and / or damage to the endoscope. For laser treatment, please wear appropriate protective glasses to protect your eyes.
- Do not insert the laser probe into the clamp pipe when the bend is angled. Otherwise cause damage to the forceps pipe.

Insert the laser probe into the forceps channel of the three-way valve

1. insert the laser probe into the protective sheath.
2. Straighten the bending part and insert the laser probe slowly and vertically into the pliers mouth of the three-way valve
3. Hold the laser probe near the three-way valve, gently tap and slowly insert it into the clamp pipe forward

Operate laser system

Operate the laser system according to the laser system instructions.

Remove the laser probe

Confirm that the laser is not emitted and carefully withdraw the laser probe from the endoscope.

4.3 Withdraw the endoscope

Warning

- If blood unexpectedly adhere to the surface of the insertion tube of the withdraw endoscope, check the conditions of the patient carefully.
- If the endoscope cannot be withdrawn from the patient smoothly, do not attempt to forcibly withdraw it. Forcibly withdrawing the endoscope may cause patient injury.

1. Carefully withdraw the endoscope while observing the endoscopic image.
2. Press the power switch to turn off the whole equipment.

Chapter 5 Cleaning and Disinfection : General Policy

5.1 Key Points of cleaning and disinfection

The medical literature reports incidents of patient cross contamination resulting from improper cleaning or disinfection. It is strongly recommended that all individuals responsible for cleaning should thoroughly understand and follow the instructions in this manual and all manual instructions for auxiliary equipment.

- occupational health and safety regulations in local hospital
- cleaning and disinfection procedures for each individual equipment
- the structure and operation method of endoscopic equipment
- usage of pertinent chemical agents

Please make judgments from a professional point of view on the type and conditions of methods of cleaning and disinfection..

Please refer to the cleaning and disinfection standards of the region/country where the product is used. Contact the local health agency to learn about local standards and regulations.

This product has not been cleaned, disinfected before leaving the factory. It must be cleaned, disinfected before use. Please refer to the method recommended by the manufacturer for cleaning and disinfecting.

Since the product can only withstand 1000 times of peracetic acid disinfection , please record the number of disinfection at each time to avoid accidents caused by product failure after disinfection times exceeding 1000 times.

5.2 Precaution

Warning

- After each case, if the endoscopy equipment is not effectively cleaned, disinfected , the safety of the patient will be endangered. In order to reduce the risk of cross-infection, the endoscope must be thoroughly cleaned according to the contents of Chapter 7 " Cleaning and Disinfection Procedures" after each case inspection, and then disinfected . The outer surface of the endoscope should be cleaned, disinfected.

- If the endoscope is not cleaned meticulously, effective disinfection may not be possible. Before using the endoscope cleaning machine, it should be confirmed that it can clean, disinfect the endoscope. If you are not sure whether the endoscope cleaning machine can clean and disinfect the endoscope, please contact the manufacturer of the endoscope cleaning machine to confirm the specific instructions or connector information. If the endoscope is not cleaned and disinfected properly, it may cause infection risk to the patients who use the endoscope next time.
- The manual specifies fitting/usable and unfitting/unusable detergent, disinfectant, sterilant and another chemical agent and washer/disinfector for the endoscope. As for the detergent, disinfectant, sterilant and other chemical agent and washer/disinfector that are not shown in this operation manual, please contact our company, service center qualified by our company or Sales Representatives of our company. Improper use of chemical agent or washer/disinfector may cause deterioration of endoscope, parts of the endoscope fall off and cause danger to patient's health.
- The chemical agents used for cleaning, disinfection maybe dangerous, so personal protective equipment should be worn to avoid being harmed by dangerous chemical agents and potentially infectious substances. When cleaning, disinfecting , please wear appropriate personal protective equipment, such as goggles, mask, waterproof clothing and protective gloves.
- Be sure to rinse thoroughly with chemicals. Clean and disinfect the external surface of the endoscope thoroughly with (Aq. Dest) water to remove any residual chemicals.
- Be sure to perform the water leakage test before the manual cleaning. If leakage is detected, stop using the endoscope. Just before each procedure, check if the endoscope needs cleaning and disinfection. If necessary, clean, disinfect it by following the instructions given in this manual.

Chapter 6 Suitable methods of cleaning, disinfection and chemical preparation

6.1 Compatibility Summary

Certain materials and structure of our endoscope may not be compatible with some cleaning, disinfection methods.

We have below two standards to decide whether the method is effective:

- Microbiological effectiveness
- Material durability

Microbiological effectiveness

The standard methods described in this manual may be useful for cleaning and disinfection of this equipment.

Material durability

This product can be cleaned, disinfected for up to 1000 times with this method.

Cleaning and disinfection methods

Choose appropriate cleaning and disinfection methods in compliance with the recommendations of your infection control committee and all national and local hospital guidelines and policies.

In summary, adopt the cleaning methods that can achieve appropriate efficacy. Our company would recommend automatic cleaning methods. Please contact your local infection control agency to determine the appropriate cleaning method for the equipment.

Monitoring

In accordance with national and local regulations or requirements, inspect and confirm cleaning and disinfection methods for reusable instruments, record the regular interval and cleaning and disinfection results.

Although there is no available biological indicator to confirm the disinfection process, test strips can be used to monitor the concentration of disinfectants.

Monitoring the concentration of preparations in daily use can ensure that the solution does not fall below

the effective concentration.

In order to monitor the process, appropriate biological or chemical indicators should be used.

List of compatible methods and chemical preparations

The methods and preparations for washing, disinfection, listed in this table, have been previously verified on endoscopes and their accessories. So that it's suitable for endoscopes and fittings and joints of leakage detector

Recommended cleaning and disinfection solution	
Cleaning Reagent	EmPower multienzyme cleaning agent
Disinfection Solution	Peracetic acid disinfection(Dilute the peracetic acid with pure water at the ratio of 1:22.1 to make the final concentration between 2380mg/L-2814 mg/L.)

Form 6.1

- Never use steam sterilization (high-temperature autoclaving) for endoscopes. Steam sterilization (high temperature autoclaving) may damage endoscopy.

6.2 Cleaning Reagent

Use cleaning reagent which manufacturer recommend EmPower multienzyme cleaning agent.

- The cleaning reagent must comply with national and local regulations or requirements and has been certified. Inadequate cleaning agents will bring the risk of infection.

- Do not reuse cleaning reagent.
- The cleaning reagent must be approved by our company for compatibility with the endoscope and certified by authority organization. Uncompatible detergent slution may bring great damage to the endoscope and accessories.

6.3 Disinfection Solution

Peracetic acid disinfection: According to the manufacturer's instruction manual, in general, 2380mg/L-2814mg/L peracetic acid disinfection solution is suitable for the Endoscope produced by our company. If the disinfectant is reused, the test paper should be used to regularly test its effectiveness. Do

not use solutions beyond the expiration date.

Caution

- Do not use endoscope only disinfected by alcohol.
- The disinfection solution must be confirmed by our company for compatibility with the endoscope and certified by authority organization. Uncompatible disinfection solution may pose great damage to the endoscope and accessories.
- Refer to the manufacturer's instructions for the concentration, temperature and immersion time. Do not exceed the level specified by the manufacturer.

6.4 The Information Of Washer

Name	Automatic soft endoscope washer-disinfector
Model	Rider 60A
Manufacturer	SHINVA MEDICAL INSTRUMENT CO., LTD.

6.5 Rinsing Water

Use water that has been processed (e.g. filtered) to improve its microbiological quality. When non-sterile water is used after manual disinfection, wipe the endoscope and flush the channels with 75%~95% ethyl alcohol, then air-dry all internal channels to inhibit the growth of residual bacteria. The selected water quality should be tested by national and local regulations or requirements, and certified by Infection Control Committee in the hospital.

Warning

- Due to possible chlorination, do not use tap water to wash the equipment.
- Do not reuse rinsing water.

Chapter 7 Cleaning and Disinfection Procedures

Warning

- Medical institutions need to have the operation specifications and related equipment for endoscopic cleaning, disinfection and sterilization.
- ALL channels of the endoscope must be cleaned, disinfection during every reprocessing cycle, even if the channels were not used during the previous patient procedure. Otherwise, insufficient cleaning or disinfection of the endoscope may bring an infection risk to the patient and/or operators when using the endoscope.
- Do not bend the endoscope's insertion tube with a diameter of less than 12 cm. The endoscope can be damaged if bending too tightly.

Caution

This chapter covers the cleaning and disinfection steps of endoscopic equipment related to patient safety and health.

7.1 Maintain

Chapter 7 focuses on the endoscope maintenance and maintenance methods from the endoscope preparation, cleaning, and disinfection.

Warning

- Do not perform maintenance while the operator is using the device for the patient, which may harm the patient and the operator.



7.2 Required Cleaning and Disinfection Equipment

Preparation of the Equipment

Prior to cleaning and disinfection, prepare the equipment shown in Fig 7.1.

Warning

- Use basins which are at least 40 cm × 40 cm (16" × 16") in size and deep enough to allow the endoscope to be completely immersed. Incompletely immersed can compromise the effect of cleaning.

Items	Pictures
Channel cleaning brush	
Channel port cleaning brush	
leakage detector	
waterproof cover and balanced valve	
Personal protective equipment	Sterile water (Rinsed in sterile water)
500cm ³ (500ml)	70% ethyl or isopropyl alcohol (Rinsed in non-sterile water)
Soft brush	Clean, lint-free cloths / paper
Clean water	Sterile, lint-free cloths
Detergent solution	Sterile, cotton swabs
Disinfectant solution	30cm ³ (30 ml) syringe
Pressure water gun	Pressure air gun
40cm×40cm(16"×16")Basins with tight fitting lids, 40 cm × 40 cm (16" × 16") minimum size	Basins with tight fitting lids, 25 cm × 25 cm (10" × 10") minimum size
(Aq.Dest) Deionized water (Aq. Dest)	Small basins

Form 7.1

Cleaning and Disinfection Equipment Parts and Functions

For inspection of other equipment than that mentioned below, refer to the instruction manual for the equipment being used.

Waterproof cover

1. When cleaning and Disinfection, the waterproof cover must be attached to the electrical connection part above the endoscope. (Fig 7.1)


Fig 7.1

Channel port cleaning brush

The Channel port cleaning brush is used to brush the endoscopic injection opening. (Fig 7.2)


Fig 7.2

Channel Cleaning Brush

Warning

- Never use cleaning brush for cell tissue sampling or other diagnostic or therapeutic purposes. Otherwise, it can result in injury to the patient, cross infection, or damage to the device.

Channel cleaning brush is used to brush the inside of the working channel. (Fig 7.3)


Fig 7.3

Leakage detector

During endoscopy leak testing, connect the cap of leakage detector with the leak valve connector of the endoscope handle and press the airbag to check for leakage .(Fig7.4)


Fig 7.4

Inspection of the Reusable Equipment

For inspection of other equipment than that mentioned below, refer to the instruction manual for the equipment being used.

Warning

- All equipments mentioned below are consumable items. Whenever the slightest irregularity is suspected, use a spare instead. Using a defective equipment may make it difficult to effectively reprocess the endoscope, and cause endoscope and/or equipment damage.
- Please clean and disinfect the equipment before each use, and check the sealing condition of the equipment to ensure that the endoscope can be used normally. Refer to the following steps.

Inspection of the Pressure Compensation Cap

Make sure that is without scratches, flaws and debris. (See Fig7.4)

Inspection of the Channel Cleaning Brushes


1. Make sure that the brush section and the metal tip of the distal tip are securely in place. Check for loose or missing bristles.
2. Check for bends, scratches and other damage to the shaft.
3. Check for debris on the shaft and/or in the bristles of the brush.

Inspection of the Channel port cleaning brush

1. Check for loose or missing bristles. (See Fig 7.2)
2. Check for debris on the shaft and/or in the bristles of the brush.
3. Check for bends, scratches and other damage to the shaft.

Inspection of the leakage detector

1. Make sure that the leakage detector is without cracks, scratches, flaws, debris and/or other damage. (See Fig 7.4)
2. Press the airbag, make sure that the leak detect or release the air from the connector. (see Fig 7.5)

(leakage detector) Fig 7.5

7.3 Cleaning and Disinfection Procedures

After the procedure, clean, disinfect the endoscope according to the procedure described below.

Endoscope reprocessing summary chart


Fig 7.6

Warning

- All channels of the endoscope must be cleaned, disinfected during every reprocessing cycle, even if the channels were not used during the previous patient procedure. Otherwise, insufficient cleaning of the endoscope may bring an infection risk to the patient and/or operators using the endoscope.
- Do not skip any of the steps described below. Otherwise, it may lead to infection.
- Do not store disinfected instruments. Use it immediately after disinfection, or send it to the sterilization room. If you want to use stored instruments, please refer to local regulations and guidelines.

7.4 Pre-cleaning

Water wash

Warning

- If the endoscope is not immediately cleaned after each procedure, residual organic debris will begin to solidify and it may be difficult to effectively clean, disinfect the endoscope.

The endoscope should be pre-cleaned immediately after each examination.

pretreatment

1. After the endoscope is used, unload the display or remove the extension cable . The display / extension cord is not waterproof.

Fig 7.7

2. Cover the waterproof cover and immediately wipe off all visible dirt on the outer surface of the endoscope with lint free gauze containing cleaning solution. Lint free gauze shall be used once. Then repeatedly wash the working channel with syringe and detergent solution for 10 seconds. Cover the waterproof cover , put it into the conveying container, and then send it to the cleaning and disinfection room.

Warning

- Without a waterproof cover , poor contact with metal contacts.

Fig 7.8

7.5 Transportation of the Endoscope

Warning

- Do not spill out debris or pre-cleaning liquid during transportation.

Note

- Please use the handling equipment to transport endoscope to the cleaning, disinfection area.

Transporting within the hospital

1. Place the endoscope on the transport trolley or use a bracket to secure the device. The specific container will protect the endoscope from damage during transport.

2. Transport the endoscope and reusable components from using site to the cleaning, disinfection area.

Caution

- The equipment should be cleaned, disinfected immediately after use. Do not let the instruments without cleaning, disinfection overnight. If you place the instrument for a long time without cleaning and disinfection processing, equipment of debris will become hard and difficult to remove. If you leave the equipment into the liquid for a long time without cleaning and disinfection processing, internal endoscope will become too humid and result in equipment damage.
- Do not hit the distal tip of endoscope during transportation. Otherwise it will damage the distal tip.
- Do not only hold the endoscope insertion section to take up the whole endoscope. Otherwise it will damage the endoscope.
- Please refer to Section 2.3 of the manual for transportation environment

7.6 Leakage Testing

After pre-cleaning, perform leakage testing on the endoscope to ensure that it is waterproof.

Equipment needed

Please prepare the following equipment:

- Personal protective equipment
- The water basin is at least 40cm x 40cm(16" x 16") and deep enough to completely soak the entire endoscope
- Deionized water (Aq. Dest.)
- Leakage tester

Caution

- The balance valve shall be removed before soaking. Otherwise, the water will enter the endoscope and cause damage to the equipment
- Rotate the Pressure compensation cap until it stops. If it is not fully and properly attached, the endoscope's interior will not be pressurized and accurate leakage testing will be impossible.
- For leak detection, if continuous bubbles emerge from one part of the endoscope, it indicates the leak. This suggests that the water is able to enter the interior of the endoscope. If a leak is found, remove the endoscope from the water and contact the company.
- When connecting the Pressure compensation cap, make sure that the inside of the leakage tester's and the outside of the water-resistant cap's leakage valve are thoroughly wipe and dry. Water remaining on either component may penetrate the inside of the water-resistant cap and could cause endoscope malfunction.
- Be sure to remove the monitor from the operation section before immersion, because the monitor can not be immersed, otherwise it will cause water penetration in the monitor and malfunction.

Note

- When the Pressure compensation cap is in place, the covering of the bending section will expand as the air pressure inside the endoscope increases. This is normal.

Leakage Testing With Gasbag

- Always keep the leakage tester dry. Any water remaining on the

leakage tester may cause damage to the endoscope.

- Press the pressure release valve before taking the Pressure compensation cap off the leakage detector. Discharge air in the endoscope until the indicator on the screen points to 0 mmHg. Otherwise, it may lead to endoscope damage.
- Only the duct of leakage tester and the Pressure compensation cap can be immersed.(as shown in Fig 7.11).Immersion of other components of the leakage tester may lead to its damage.
- The pressure shall not exceed 240 mmHg. If the pressure increases too much, it might damage the endoscope. Press the pressure release valve to discharge the air (as shown in Fig 7.11).

The endoscope is completely waterproof and can be thoroughly soaked and cleaned. However, the tightness of the endoscope should be checked before immersion.

1. **Connect the leakage detector:** Point the cracks in the leak cover to the slot of the leak joint on the video endoscope. Push into the leak cover and rotate clockwise to the head.


Fig 7.10

2. **Pressurization:** use the random leakage detector to pump air into the inside of the endoscope body through the leak detection cover, and maintain it to 220-240 mm Hg. (Note: Pump slowly when pumping)Stabilize the indicator needle for a few seconds after compression. Read the pressure when the indicator needle is stable.


Fig 7.11

3.Observation: immerse the endoscope in clear water, and observe for about 30 seconds while making an angle, and confirm that there is no continuous bubble on the whole endoscope. If leakage is detected, please contact the manufacturer.


Fig 7.12

Warning

- If the indicator needle drops to 0 mmHg continuously, it means that the endoscope may have serious water leakage or the leakage detector may be damaged. The leak test should be stopped immediately. If the endoscope is still immersed in water, water will enter the endoscope without pressure inside. As a result, more serious problems will occur.
- During the leakage test, if bubbles continuously emerge from the leakage test cover, the leakage test cover or joint may be damaged. Please replace the leak detection cover in time or contact the manufacturer for maintenance.

4. Remove the leakage detector: take out the endoscope from the clear water, unscrew the vent valve on the left side of the middle of the leakage detector, release the air pressure inside the endoscope, and remove the leakage detector.


Fig 7.13

7.7 Manual Cleaning

Warning

- Clean the endoscope and accessories thoroughly to remove microorganisms or organic material. If the endoscope is not immediately cleaned after each procedure, residual organic debris will begin to solidify and it may be difficult to effectively clean, disinfect the endoscope.

Caution

- To avoid water leakage, do not use excessive force when cleaning the endoscope.
- Remove balance valve from endoscope. If a balance valve is installed, water may enter the endoscope and damage the endoscope.
- For details of detergent, please refer to section 6.2 "Detergent".

After completing the leakage test, perform manual cleaning according to the procedures described below.

Equipment needed

Prepare the following equipment

- Personal protective equipment.
- Basin which is at least 40 cm × 40 cm (16" × 16") in size and deep enough to allow the endoscope to be completely immersed.
- clean, lint-free cloths
- EmPower multienzyme cleaning agent
- 30cm³ (30 ml) syringe
- cotton swab

- channel cleaning brush
- Channel connector cleaning brush
- soft brush

Caution

- In order to avoid damage to the endoscope, do not soak the endoscope with other instruments except the above instruments.

Preparation

1. Wear appropriate personal protective equipment.
2. Fill a basin with EmPower solution at the room temperature and the concentration is 1:128. Use a basin which is at least 40cm x 40 cm (16" x 16") in size and deep enough to allow the endoscope to be completely immersed.

Cleaning the external surface

1. Soak the endoscope in the washing liquid.
2. Use a clean lint free cloth to wipe the external surface of the whole endoscope, especially the insertion part and the control part. Use and replace the lint free cloth once to ensure that the external surface is completely clean. Do not infiltrate liquid into the monitor. Liquid infiltrating into the monitor may damage the internal circuits.
3. Visually check that there is no residual debris on the surface of the endoscope.

Brushing the channels

Warning

- Be sure to brush the inside of the instrument and suction channel. Otherwise, insufficient cleaning and disinfection of the endoscope may pose an infection control risk to the patient and/or operators performing the next procedure with the endoscope.
- To avoid splattering detergent solution when the channel cleaning brush is pulled out, keep the endoscope submerged.

Caution

- The channel cleaning brush is a consumable item. Disposable brush belongs to the disposable products that cannot be repeated use. Repeated using may cause the brush head being bent or kinked, as a result the brush head may come off. Confirm that the brush is free from any damages or other irregularities before and after each use. If the brush head comes off after brushing, immediately retrieve it and carefully check that it is not left inside the channel of the endoscope by passing a new cleaning brush or other endo-therapy accessories.

If the brush head is left in the channel, it can drop into the patient body during the procedure. Due to the difference in location, the new cleaning brush passing through the pipe may not be able to remove the missing parts. In this case, please contact our company as soon as possible.

- Do not try to pass the cleaning brush through the tip of the insertion part. Otherwise, the cleaning brush will get stuck and cannot be pulled out.

Follow the steps below to scrub all pipes and injection ports while the endoscope is still immersed in the washing liquid (as shown in Fig7.14).


Fig 7.14

Brushing injection forceps tubing (location A)

1. Hold the cleaning brush at a distance of 3 cm from the bristles.
2. Insert the channel cleaning brush into the suction valve Housing as illustrated by A in Fig. Using short strokes, feed the brush through the channel until it emerges from the distal tip of the endoscope.
3. Clean the bristles with your fingertips in the detergent solution. Carefully pull the brush back through the channel.
4. Clean the bristles in detergent solution again.
5. Visual inspection, if there was residual debris on the sample, repeat Step 1 to 4 until all debris is removed.

Brushing injection / forceps tube opening (location B)

1. As shown in Fig B, inserted the tube channel port cleaning brush into the injection tube opening until the brush handle touches the tube opening.
2. Turned the inserted brush once.
3. Pulled the brush out and clean the bristles with your fingertips in the detergent solution.
4. Visual inspection, if there was residual debris on the sample, repeat Step 1 to 3 until all debris is removed.

Soak the endoscope in the detergent solution

1. Completely immerse the endoscope in the washing solution.
2. Attach the 30cm³ (30 ml) syringe on the injection port.
3. Pull up the push rod of the syringe and fill the working channel with detergent. (See Fig 7.15)
4. Completely immersed the endoscope in the detergent solution and detach the syringe.


Fig 7.15

5. Using a lint-free cloth, wipe all debris from the endoscope's external surfaces while the endoscope is immersed in the detergent solution.

6. Soak the endoscope for 5 min in EmPower solution at the room temperature and the concentration is 1:128

7. Remove the endoscope from the detergent solution.

7.8 Rinsing and Drying after Manual Cleaning

Warning

- Rinse the external surface and channels of the endoscope, and wash thoroughly the equipment with rinse water to remove the residual detergent solution. Otherwise, the residual detergent solution may cause mucosal inflammation.

Caution

- do not use tap water to wash the equipment

Equipment needed

Please prepare the following equipment:

- personal protective equipment
- basin, 40 cm × 40 cm (16" × 16") minimum size and deep enough to immerse the endoscope completely
- deionized water (Aq. Dest.)
- clean, lint-free cloth or paper
- Pressure water gun
- Pressure air gun

Rinse the outer surface

1. Wear appropriate personal protective equipment.
2. Fill a basin with the deionized water (Aq. Dest.). Use a basin which is at least 40 cm x 40 cm (16"x 16") in size and deep enough to allow the endoscope to be completely immersed.
3. Soak the endoscope and suction cleaning adapter in the deionized water (Aq. Dest.), and gently agitate to thoroughly rinse them.

Rinsing all channels

1. Aim the pressure water gun at the injection port to flush the working channel until there is no residual cleaning solution.

2. Flush the external surface of the endoscope and all pipe ports with flowing water.
3. Use a pressure air gun to inflate the pipe for 30 s to remove the water in the pipe.

Dry external surfaces

1. Place the endoscope, buttons and valves on a special drying table with sterile towels. Sterile towels should be replaced every 4 hours.
2. Fill all pipes with 75%~95% ethanol or isopropyl alcohol.
3. Use a pressure air gun to inflate all pipes with clean compressed air for at least 30 s until they are completely dry.
4. Dry the outer surface of the endoscope, buttons and valves with sterile wiping cloth and pressure air gun.
5. Install buttons and valves.

7.9 Disinfection

Warning

- During the entire process of disinfection, the endoscope and all equipment should be completely immersed in the disinfectant. If the equipment is not removed from the endoscope during immersion, or all parts are not completely immersed, the disinfectant will not be able to fully contact all surfaces of the equipment. This may reduce the effectiveness of disinfection.
- Refer to "Disinfectant solution" in Section 6.3 for details of disinfectant solution.

After manual cleaning, disinfect the endoscope according to the procedures described below.

Equipment needed

Please prepare the following equipment

- personal protective equipment (gloves, glasses and mask)
- Peracetic Acid (Dilute the peracetic acid with pure water at the ratio of 1:22.1 to make the final concentration between 2380mg/L-2814 mg/L.)
- Automatic soft endoscope washer-disinfector
- Suction cleaning adapter
- Pressure air gun

Note

The information of disinfectant:

Name: Automatic soft endoscope washer-disinfector

Model: Rider 60A

Manufacturer: SHINVA MEDICAL INSTRUMENT CO., LTD.

Process

1. Wear appropriate personal protective equipment
2. Put device into the Automatic soft endoscope washer-disinfector until it is fully loaded.
3. Soak the Flexible Digital Uretero-Renoscope in peracetic acid for 10 min.
4. Disinfection: Time: 5 min.
5. Rise 1: Time: 30 s.
6. Rise 2: Time: 30 s.

Warning

- All disinfection steps should be performed with the endoscope and all equipment completely immersed. If the equipment is disconnected respectively while not immersed, disinfectant solution may not adequately contact all surfaces of the equipment. As a result, the effectiveness of disinfection may be reduced.

Removing the endoscope and all equipment from automatic soft endoscope washer-disinfector

1. Remove the endoscope from the disinfectant solution.
2. Use a pressure air gun to inject air into the pipe through the injection port for 30s to remove the disinfectant in the pipe.

7.10 Rinsing and Drying after Disinfection**Warning**

- After cleaning and disinfection, thoroughly clean and dry the tube of the endoscope. Otherwise, bacteria will grow in the pipeline, posing a risk of infection to the patient or operator who uses the endoscope next time.
- Rinse the external surface and channels of the endoscope, and wash thoroughly the equipment with rinse water to remove the residual

detergent solution. Otherwise, the residual detergent solution may cause mucosal inflammation.

Caution

- Alcohol is flammable. Handle it with care.
- About the quality of the rinsing water, please refer to Section 6.5 "Rinsing Water".

After disinfection, rinse and dry according to the following steps endoscope and all equipment.

Equipment needed

Please prepare the following equipment

Warning

- Never reuse the rinsing water.
- Personal protective equipment
- Sterile basins, 40 cm × 40 cm (16" × 16") minimum size and deep enough to immerse the endoscope completely
- Sterile deionized water (for sterile water rinse)
- Sterile, lint-free cloths
- Pressure water gun
- Pressure air gun

If sterile water is not available, prepare the following equipment.

- Deionized water (aq. Dest.) (For non-sterile water rinse)
- Small basin
- 70% ethyl or isopropyl alcohol
- Sterile cotton swabs
- small basin
- Pressure water gun
- Pressure air gun

Sterile water rinse

1. Wearing appropriate personal protective equipment.
2. Put sterile water in the basin. The water basin should be at least 40 cm × 40 cm (16" × 16") deep enough to soak the entire endoscope.
3. Immerse the endoscope in the sterile water. Using a sterile, lint-free cloth, thoroughly rinse and wipe

all external surfaces.

4. Aim the pressure water gun at the injection port and flush the pipeline with sterile water for at least 2 minutes to flush the working channel until there is no residual disinfectant.
5. Place the endoscope, buttons and valves on a special drying table with sterile towels. Sterile towels should be replaced every 4 hours.
6. Fill all pipes with 75%~95% ethanol or isopropyl alcohol.
7. Use a pressure air gun to inflate all pipes with clean compressed air for at least 30 s until they are completely dry.
8. Dry the outer surface of the endoscope, buttons and valves with sterile wiping cloth and pressure air gun.
9. Install buttons and valves.

Note

- Flushing the channels with 70% ethyl or isopropyl alcohol after rinsing them with sterile water facilitates drying inside the channels.

Non-sterile water rinse and alcohol flush

1. Wearing appropriate personal protective equipment.
2. Fill a basin with deionized water (Aq. Dest.). Use a basin which is at least 40 cm × 40 cm (16" × 16") in size and deep enough to allow the endoscope to be completely immersed.
3. Immerse the endoscope in deionized water. Use a sterile lint-free cloth to thoroughly rinse and wipe all external surfaces.
4. Aim the pressure water gun at the injection port and flush the pipeline with deionized water for at least 2 minutes to flush the working channel until there is no residual disinfectant.
5. Take out the endoscope and put the endoscope and corresponding accessories on the special drying table for laying sterile towel.
6. Put 70% ethanol or isopropyl alcohol in a small bowl.
7. Pull up the syringe push rod and inject 70% ethanol or isopropyl alcohol into the working channel from the injection port. Make sure that the ethanol or isopropyl alcohol injected from the injection port flows out of the inner lens end.
8. Use a pressure air gun to inflate all pipes with clean compression for at least 30s until they are completely dry.
9. Thoroughly wipe and dry the external surfaces of the endoscope and the injection port with a sterile lint

free cloth dipped in alcohol.

10. Using a sterile, lint-free cotton swabs, wipe and dry the inside of the injection port.

11. Store the components as described in Chapter 8 "Storage, Handling and Disposal Outside the Hospital".

Chapter 8 Storage and Disposal

After cleaning, and disinfection the clean endoscope and accessories should be kept away from any contaminated equipment. If the clean endoscope and accessories are contaminated between two inspections, it will pose a risk of infection to the patient or operator who will use it next time.

- To prevent contamination of clean endoscope and accessories, be sure to keep storage cabinet also clean.
- The storage cabinet must be clean, dry, well ventilated and maintained at ambient temperature. Do not store the endoscope in direct sunlight, at high temperature, in high humidity or exposed to X-rays and ultraviolet-rays. Otherwise, it may damage the endoscope or present an infection control risk.
- Do not store the endoscope in the carrying case. Use the carrying case only for shipping the endoscope. Routinely storing the endoscope in a humid, non-ventilated environment such as the carrying case may present an infection control risk.

8.1 Storage of the Endoscope

1. Confirm that the outer surfaces of the endoscope (especially all internal channel, distal tip, lens and electrical contacts) are dry.
2. Carefully wipe the objective lens of the distal tip with cotton swab moistened with 70% alcohol.
3. Hang the endoscope in the storage cabinet with the distal tip hanging freely. Make sure that the insertion tube hangs vertically and as straight as possible.
4. Store the endoscope in a clean and dry place at room temperature (5-35 ° C, humidity of 30-85%).

Note

- According to the differences of the packaging type and storage conditions (refer to DIN 58 953, Section 9, or local regulations) the storage time of equipment is also different.

8.2 Storage of Reusable Parts and Reprocessing Equipment

1. Confirm that all the reusable parts and reprocessing equipment are thoroughly dry.
2. Store all reusable parts loosely in the storage cabinet. Make sure these parts during storage are not in

contact with each other.

3. Store all cleaning equipment in a container, and then store the container in the storage cabinet.

8.3 Transporting Outside the Hospital

Transport the endoscope in the carrying case.

Warning

- The endoscope must be cleaned, disinfected before being transported outside the hospital. Otherwise, there will be a risk of infection

Note

- The carrying case cannot be cleaned, disinfected. Clean and disinfect the endoscope before placing it in the carrying case.

8.4 Disposal

When handling the endoscope and any of its components (such as cleaning brushes), all applicable national and local laws and regulations should be followed.

Chapter 9 Troubleshooting

If it is found that the endoscope is damaged and cannot function normally, or if there are other abnormalities found during the inspection according to the method described in Chapter 3 "Preparation and Inspection", do not use the endoscope. Please contact our company.

Some non-functional faults can be solved by the content of Section 9.1 "Troubleshooting Guide". If the fault still cannot be eliminated after taking the above-mentioned treatment methods, please stop using the endoscope and send it back to our company for repair.

Our company is not responsible for repairing accessories. If the accessory is damaged, please contact our company to purchase a new one.

- If an abnormality of the endoscope is observed, do not use it on the patient. The use of damaged or abnormal endoscopes will endanger the lives of patients and users, as well as cause more serious equipment damage.
- If parts of the endoscope fall into the patient's body due to equipment damage or malfunction, immediately stop using the endoscope and take appropriate measures to retract the parts.

If any abnormality in the function of the endoscope and/or endoscopic image is suspected during use, stop the examination immediately and carefully withdraw the endoscope from the patient as described in Section 9.2, "Withdrawal of the Endoscope with any Abnormality".

9.1 Troubleshooting Guide

The following table shows faults and countermeasures caused by incorrect settings or damaged consumables.

Failures caused by reasons other than the following should be repaired. Repairs performed by maintenance personnel not authorized by the company may cause injury to the patient or user and/or damage to the equipment. Therefore, please contact the company for maintenance in accordance with the provisions in section 9.3 "Repair of Endoscopes".

Irregularity description	Possible cause	Solution
Fluid leakage at the working channel port	Blocking channel	Cleaning the channels with Channel Cleaning Brush.
	The Channel Port is damaged.	Replace it with a new one.
Leakage of fluid between injection port	Blocking channel	Cleaning the channels with Channel Cleaning Brush.

and syringe	The syringe was not inserted firmly	Insert syringe firmly
An image is not clear.	The endoscope lens is dirty.	Wipe with cotton swab moistened with 70% ethyl or isopropyl alcohol; If the image can't be observed, return the endoscope to our company for repair.

Form 9.1

9.2 Withdrawal of the Endoscope with any Abnormality

If an abnormality occurs while the endoscope is in use, take a proper measure as described in either “When the endoscopic image appears on the monitor” or “When the endoscopic image does not appear on the monitor or the frozen image cannot be restored” below. After withdrawal, return the endoscope for repair as described in Section 10.3, “Returning the endoscope for repair”.

- If the endoscope cannot be withdrawn from the patient smoothly, do not attempt to forcibly withdraw it. If any irregularities are suspected, immediately contact us. Forcibly withdrawing the endoscope may cause patient injury.

When the endoscopic image appears on the monitor

1. Turn off the power supply for all devices except the monitor and infinite transmitter.
2. Carefully withdraw the endoscope while observing the endoscopic image.

When the endoscopic image does not appear on the monitor or the frozen image cannot be restored

1. Turn off the power supply for all devices except the monitor and infinite transmitter.
2. Turn the monitor OFF and then ON again. If the endoscopic image appears or the frozen image restored, follow the procedure of step 2 and steps below in “When the endoscopic image appears on the monitor”.
When the endoscopic image still does not appear or the frozen image cannot be restored, perform the following steps.
3. Turn the UP/DOWN angulation control knob to neutral positions.
4. Release the angulation control knobs and carefully withdraw the endoscope.

9.3 Returning the Endoscope for Repair

- Before returning the endoscope for repair, thoroughly clean and disinfect it. Improper cleaning, and disinfection can bring infection

risks to the hospital or the company's personnel who dispose of the equipment.

●

- Our company is not responsible for the equipment damage and/or person injury which are not caused by our company personnel trying to repair the equipment.

Please contact our company before returning the endoscope for repair. When returning for repair, please attach a description of the malfunction or damage of the endoscope, as well as the name and contact number of the person most familiar with the malfunction of the endoscope in your organization. Please also attach the warranty card.

When returning the endoscope for repair, please carry it according to the instructions in section 9.3 "Transporting the Endoscope Outside the Hospital".

Chapter 10 Electromagnetic compatibility

Electromagnetic compatibility

Caution:

- **UF series Video Endoscope meets the requirements of IEC 60601-1-2 standard for electromagnetic compatibility**
- **The user shall install and use the electromagnetic compatibility information provided by the random files**
- **Portable and mobile RF communications devices may affect the performance of the UF series Video Endoscope and avoid strong electromagnetic interference such as near mobile phones and microwave ovens;**
- **The guidelines and manufacturer's declarations are detailed in the annex**

Warning:

- **A device or system should not be used with other devices placed close to or over stacked, if must be close to or stacked, should observe and verify the normal operation in its use of the configuration.**
- **Type A equipment for the industrial environment, since the UF series of Video Endoscope conduction harassment and radiation harassment, in other environments to ensure electromagnetic compatibility may have potential difficulties.**
- **The provisions of the accessories, transducer or cables and equipment and systems used together, which may cause the equipment or system emission increase or decrease of immunity.**

Annex:

Guide and manufacturer's statement - Electromagnetic Emission		
UF series of Video Endoscope used in the electromagnetic environment of the expected performance of the following provisions, the performance of UF series of Video Endoscope purchasers or users should ensure that it is used in the electromagnetic environment:		
Emissions test	Compliance	Electromagnetic environment - Guide
RF emissions CISPR 11	Group 1	The performance of UF series Video Endoscope is only its internal function and the use of RF energy. Therefore, its RF emission is very low, and may not cause any interference to nearby electronic equipment.
RF emissions CISPR 11	Class A	The UF series Video Endoscope is suitable for use in non-household and not directly connected to the public facilities for all residential low-voltage power supply network in the home
Harmonic emissions IEC 61000-3-2	NA	
Voltage Fluctuations / Flicker Emissions IEC 61000-3-3	NA	

Guide and manufacturer's statement - Electromagnetic Immunity			
UF series of Video Endoscope used in the electromagnetic environment of the expected performance of the following provisions, the performance of UF series of Video Endoscope purchasers or users should ensure that it is used in the electromagnetic environment:			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guide
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	The ground should be wood, concrete or tile, if the ground with synthetic material covering, the relative humidity should be at least 30%
Electrical Fast Transient IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	The power network should have the quality of typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 0,5$ kV, ± 1 kV, ± 2 kV	$\pm 0,5$ kV, ± 1 kV, ± 2 kV	The power network should have the quality of typical commercial or hospital environment.
Voltage dips IEC 61000-4-11	0 % UT; 1 cycle And 70 % UT; 25/30 cycles Single phase: at 0°	0 % UT; 1 cycle And 70 % UT; 25/30 cycles Single phase: at 0°	The power network should have the quality of typical commercial or hospital environment. If the UF series of Video Endoscope performance users need to run continuously in the power interruption, it is recommended that the UF series of Video Endoscope performance using the uninterruptible power supply or battery power.

Power frequency magnetic field (50/60Hz) IEC 61000-4-8	3A/m(50Hz/60Hz)	3A/m(50Hz/60Hz)	The power frequency magnetic field shall have the horizontal magnetic field characteristics of a typical place in a typical commercial or hospital environment
---	-----------------	-----------------	--

Note: UT refers to an alternating voltage test voltage before

Guide and manufacturer's statement - Electromagnetic Immunity

UF series of Video Endoscope used in the electromagnetic environment of the expected performance of the following provisions, the performance of UF series of Video Endoscope purchasers or users should ensure that it is used in the electromagnetic environment:

immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guide
RF transmission IEC 61000-4-6	3 V 150 kHz to 80 MHz	3 V	Portable and mobile communication equipment should not be part of any RF closer to the UF series of Video Endoscope performance than the isolation distance recommended use, including cable. The distance should be with the transmitter frequency corresponding to the formula. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Among them, P is based on maximum output power supplied by rating of the transmitter manufacturers, Watt (W) as a unit, D is recommended isolation distance, mas a unit. Fixed field RF transmitter based on the electromagnetic field survey to determine, in each frequency range should be lower than compliance voltage In the following symbols marking equipment, interference may occur. 
RF radiation IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	

Note 1: in 80MHz and 800MHz frequency, adopt the high frequency formula.

Note 2: these guidelines may not be appropriate in all cases, the electromagnetic propagation may be affected by buildings, absorption and reflection of objects and people.

- a Fixed transmitter electric field strength, such as base station of wireless (cellular / cordless) telephone and ground mobile radio, amateur radio, AM (AM) and FM (FM) radio and television broadcasting etc. their emission electric field strength is impossible to be predicted in theory. In order to evaluate the electromagnetic environment of a fixed RF transmitter, the investigation of electromagnetic field should be conducted. If the field strength of the UF series Video Endoscope is measured higher than the RF compliance voltage of the above-mentioned application, the performance of the UF series Video Endoscope should be observed to verify its normal operation. If the observed performance is abnormal, then additional measures may be required, such as the new series of UF Video Endoscope performance orientation or location.
- b In the whole frequency range of 150KHz ~ 80MHz, the emission electric field strength should be less than 3 V/m.

Recommended isolation distance between portable and mobile communication equipment RF and UF series of Video Endoscope performance

The performance of the UF series Video Endoscope is expected to be used in the controlled electromagnetic environment of radiated RF harassment. On the basis of the maximum output power of communication equipment, the purchasers or users can maintain portable and mobile RF communication equipment through the following recommended (transmitter) and UF series Video Endoscope performance minimum distance to prevent electromagnetic interference.

Rated maximum output power of the transmitter/W	Corresponding to different frequency transmitter isolation distance/m		
	150 kHz ~ 80 MHz $d = 1.2\sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2\sqrt{P}$	800 MHz ~ 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For the transmitter rated maximum output power not listed in the above table, the recommended isolation distance d, in meters (m), can be determined by the formula in the corresponding transmitter frequency column, where P is provided by the transmitter manufacturer The maximum output rated power of the transmitter, in watts (W).

Note 1: At 80 MHz and 800 MHz, the higher frequency range formula is used.

Note 2: These guidelines may not be suitable for all situations. Electromagnetic propagation is affected by absorption and reflection from buildings, objects and humans.

Manufacturers: Zhuhai Vision Medical Technology Co., Ltd.

**Address: 301 Building 4, No. 30 Feipeng Road, GMTCM Park, Hengqin New Area,
519031 Zhuhai City, Guangdong Province, PEOPLE' S REPUBLIC OF CHINA**

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