897553

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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
THIS IS TO CERTIFY THAT: the seal of Zhuhai Vision Medical Technology Co., Ltd on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

授权签字: 
Authorized
Signature: Sun Danfeng

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证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Zhuhai Vision Medical Technology Co., Ltd.

Product Information

Product Name: Flexible Digital Bronchoscope
Product Type: BF20/BF25/BF30/BF42/BF45/BF52/BF58/BF62/BF20H/BF25H/BF30H
BF42H/BF45H/BF52H/BF58H/BF62H/BF28P/BF30P/BF42P/BF45P/BF52P
Product Life: 5 years or 1000 times peracetic acid disinfection
Manufacturer Name: Zhuhai Vision Medical Technology Co., Ltd.
Address: 301 Building 4, No. 30 Feipeng Road, GMTCM Park, Hengqin New Area,
519031 Zhuhai City, Guangdong Province, PEOPLE' S REPUBLIC OF CHINA
Contact information
Tel: +86-756-6910150
Fax : +86-756-6910511
Zip code : 519031

Version Information

This manual version may be upgraded due to a change of software or technical specifications at any time without notice. The version information of this specification are as follows:

- Document number: TD/CE-BF-A-009
- Version number: V2.4
- Edit date: Jan. 2024

©2024. The copyright of Zhuhai Vision Medical Technology Co., Ltd

CE Mark

EC Representative Name:

Shanghai International Holding Corp. GmbH (Europe)

EC Representative Address:

Eiffestraße 80, 20537 Hamburg, Germany

Statement

The manufacturer owns the copyright of this non-publicly published manual and has the right to treat it as confidential information. This instruction manual is only used as a reference material for operation, maintenance and repair of the product. Others have no right to disclose this instruction manual to others. This product manual contains proprietary information protected by copyright law; all rights reserved. Any part of this manual may not be photocopied, copied or translated into other languages without the written consent of the manufacturer.

The content contained in this product manual can be changed without notice.

Responsibilities of the Manufacturer

In the following conditions are met, the manufacturer will responsible for the product safety, reliability and performance:

The installation, maintenance or upgrade of the product shall be carried out by personnel authorized or approved by the manufacturer;

The storage environment, working environment, and electrical environment of the product conform to product specifications;

Use this product in accordance with the "Instruction Manual";

Use this product within the period of use.

Overview

◆ Main Structure

The product is composed of distal end, insertion part, operation section.

◆ Basic performance of the product

Used in conjunction with video processor product by our company, it can output real-time images without unexpected image orientation

◆ Intended Use

This product is used together with the Video Endoscope Processor produced by Zhuhai Vision Medical Technology Co., Ltd. for endoscopy and endoscopic surgery within the airways and tracheobronchial tree.

◆ Contraindication

- 1) Cardiovascular instability
- 2) Bleeding diathesis
- 3) Severe bronchospasm
- 4) Hypoxemia

◆ Complication

Throat pain and throat spasms, difficulty breathing, blood oxygen saturation decreased, bronchospasm and cough, laryngeal edema, vocal cord damage, hemorrhage, infect ect.

◆ Cautions, Warnings and Notes

- 1) This product does not contain any parts that can be repaired by the user. The internal maintenance of the instrument must be carried out by the technical maintenance personnel authorized by the manufacturer.
- 2) This product shall not be used in the home.
- 3) The connection with the video endoscope processor must be disconnected before cleaning, disinfection .
- 4) Operate in strict accordance with the instructions provided by the manufacturer.
- 5) This product cannot be used with high frequency equipment.
- 6) The equipment can not be used in MRI scan environment.
- 7) The power supply system for this product depends on the Video Endoscope Processor you selected. And the detail operation description please refer to the user manual of the Processor.

About the Manual

This manual introduces the use, function and operation method of the product in detail. Before using this product, please read and understand the contents of this manual carefully to ensure the correct use of this product and the safety of patients and users.

This manual introduces the product according to the most complete configuration, so some of the content may not be applicable to the product you purchased. If you have any questions, please contact our

company.

Please place this manual near the product so that you can get it conveniently and in time when you need it.

This manual is suitable for professional clinical medical staff to read. Readers should have knowledge and work experience in medical procedures, practice and terminology necessary for endoscopy.

All illustrations provided in this manual are for reference only, and the settings or data in the illustrations may not be exactly the same as the actual display you see on the product.

In the following description in this manual, the Flexible Digital Bronchoscope will be referred to as "Video endoscope" for short.

About the product

The main material of Video endoscope insertion part is TPU. The ingredients and raw materials of ME equipment comply with local laws and regulations. Insertion Section, distal end, bending section is an APPLIED PART, see details Fig 2.1, Fig 2.2 BF series constitute.

Every two years to verify the basic safety of the product, detailed details, please consult the manufacturer.



Catalogue

SYMBOLS	1
KEY POINTS (PLEASE READ BEFORE USE)	2
USER QUALIFICATION	2
INSTRUMENT COMPATIBILITY	2
SIGNAL WORDS	2
MAINTENANCE	2
INAPPROPRIATE MAINTENANCE AND MODIFICATION PROHIBITED	3
WARNINGS AND CAUTIONS	3
ENDOSCOPIC IMAGES DISAPPEAR OR FREEZE PRECAUTIONS	7
EXAMPLES OF INAPPROPRIATE HANDLING	7
LED LIGHT EMITTING DIODE	8
CHAPTER1 CHECKING THE PACKAGE CONTENTS	9
CHAPTER 2 INSTRUMENT NOMENCLATURE AND SPECIFICATIONS	11
2.1 CONSITUTE	11
2.2 SPECIFICATIONS	15
CHAPTER 3 PREPARATION AND INSPECTION	19
3.1 PREPARATION OF THE EQUIPMENT	19
3.2 INSPECTION OF VIDEO ENDOSCOPE	21
3.3 INSPECTION OF ACCESSORIES	23
3.4 ACCESSORIES INSTALLATION TO THE ENDOSCOPE	23
3.5 INSPECTION OF THE ENDOSCOPIC	24
3.6 EQUIPMENT CONNECTIONS 25	
3.7 Inspection of distal end 27	
3.8 IMAGE CHECK	29
3.9 CHECK THE FREEZING FUNCTION	29
CHAPTER 4 OPERATION	31
4.1 INSERT THE ENDOSCOPE	32
4.2 WITHDRAW THE ENDOSCOPE	35
CHAPTER 5 CLEANING AND DISINFECTION: GENERAL POLICY	37
5.1 KEY POINTS OF CLEANING AND DISINFECTION	37
5.2 PRECAUTION	37
CHAPTER 6 SUITABLE METHODS OF CLEANING, DISINFECTION AND CHEMICAL PREPARATION	40
6.1 COMPATIBILITY SUMMARY	40
6.2 CLEANING REAGENT	41
6.3 DISINFECTION SOLUTION	41
6.4 THE INFORMATION OF WASHER	42
6.5 RINSING WATER	42
CHAPTER 7 CLEANING AND DISINFECTION PROCEDURES	43
7.1 REQUIRED CLEANING AND DISINFECTION EQUIPMENT	43
7.2 CLEANING AND DISINFECTION PROCEDURES	48
7.3 PRE-CLEANING	50
7.4 TRANSPORTATION OF THE ENDOSCOPE	51
7.5 LEAKAGE TESTING	51
7.6 MANUAL CLEANING	55
7.7 RINSING AND DRYING AFTER MANUAL CLEANING	58
7.8 DISINFECTION	59
7.9 RINSING AND DRYING AFTER DISINFECTION	60
CHAPTER 8 STORAGE AND DISPOSAL	63
8.1 STORAGE OF THE ENDOSCOPE	63

8.2 STORAGE OF REUSABLE PARTS, REPROCESSING EQUIPMENT AND LEAKAGE TESTER	63
8.3 TRANSPORTING OUTSIDE THE HOSPITAL	64
8.4 DISPOSAL	64
CHAPTER 9 TROUBLESHOOTING	65
9.1 TROUBLESHOOTING GUIDE	65
9.2 WITHDRAWAL OF THE ENDOSCOPE WITH ANY ABNORMALITY	66
9.3 RETURNING THE ENDOSCOPE FOR REPAIR	66
CHAPTER 10 ELECTROMAGNETIC COMPATIBILITY	68

Symbols

The meanings of the symbols shown on the package with the components, the back cover of this instruction manual and/or this instrument are as follows:

Symbols	The meaning of the symbol	Symbols	The meaning of the symbol
	Follow instructions for use		Medical Device
	Type BF Applied Part		Caution In case of application as a SAFETY SIGN, the rules according to ISO 3864-1 are to be adhered to.
	Serial Number		This way up
	Fragile; handle with care		Keep dry
	Manufacturer		Date of manufacture
IPX7	Protected against the effects of temporary immersion in water		Stacking limit by number
	Symbol for Waste Electrical and Electronic Equipment (Please follow Local Ordinances or Regulations of disposal.)		Authorized Representative in the European Community
	CE mark.		

Key points (please read before use)

User Qualification

If there are official standards for use established by the medical administration or other organizations, please follow these standards. If there is no official standard, the operator of this product must be a doctor who has passed the audit of the medical safety director of the hospital or a doctor who can safely carry out the established operation according to the official regulations designated by the director of the using department according to the difficulty of using the technology.

Instrument compatibility

Refer to Chapter 1 "Checking the Package Contents", confirm that this instrument is compatible with the ancillary equipment being used. Using incompatible equipment can result in patient or operator injury and/or equipment damage.

- There is no guarantee that instruments selected solely using maximum insertion portion width and working length will be compatible in combination.

- There is no guarantee that instruments selected solely using this minimum instrument channel width will be compatible in combination.


This product conforms to the EMC standard IEC 60601-1-2.

Signal Words

The following signal words are used throughout this manual:

- Indicates an emergency hazardous situation which, if not avoided, could result in death or serious injury.

- Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

- Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to warn against unsafe operations or potential equipment damage.

- Indicates additional helpful information.

Maintenance

With increasing number of uses and / or total use time, the probability of endoscopic and adjuvant failure increases accordingly. In addition to the inspection before each operation, the person in charge of medical equipment maintenance in each hospital shall regularly check the items specified in this manual.

Endoscopes suspected of abnormality cannot be used and should be checked as per Section 9.1 troubleshooting Guidelines described. If there is still abnormal after the examination, please contact manufacturer .

Inappropriate maintenance and modification prohibited

This product does not contain any user-serviceable parts. Do not disassemble, modify or attempt to repair, otherwise it may result in injury to the patient or damage to the equipment. If necessary, please contact the manufacturer.

Manufacturer is not responsible for the injury or damage caused by the maintenance equipment of personnel not from manufacturer .

Warnings and Cautions

Follow the warnings and cautions given below when handling this instrument. This information is to be supplemented by the warnings and cautions given in each chapter.

Warnings

- Video endoscope should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the Bronchoscope should be observed to verify normal operation in the configuration in which it will be used.
- Interconnections conditions require the applied parts of other ME equipment used within the configuration for endoscopic application to be type BF applied parts or type CF applied parts.
- Class A equipment is intended to be used in industrial environment. The Video Endoscope may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances.
- When this product is used together with the laser equipment, in order to avoid the damage to the eyes, Take the corresponding protective measures, according to the classification of the laser equipment, such as by wearing suitable anti-filter glasses, or with suitable filters on the eyepiece of the live endoscope.
- Modification of this equipment is not permitted.
- Energized endoscopes are used with energized endotherapy devices, patient leakage currents may be increases.
- Be sure to prepare standby equipment or equipment with similar

functions to avoid accidents in case of equipment failure or abnormal function during operation.

- When using the endoscope, the working part of the endoscope should be coated with lubricant or the instrument channel should be injected with different liquids. After using the endoscope, the surface of the lens body and the instrument channel should be wiped and scrubbed according to the instructions. Otherwise, rapid aging of the mirror may occur.
- Do not use the endoscope in the presence of inflammable substances or in an oxygen-rich environment. Otherwise, fire or explosion may occur.
- Cooperate with the manufacturer's equipment for use, the use of incompatible equipment will lead to patient or operator injury or equipment damage.
- Instructions warning: Devices should be placed reliably to avoid falling.
- Heavy objects cannot be placed on the device, otherwise it will cause device damage or injury to the patient.
- Cooperate with the manufacturer supporting endoscope image processor use, using non-manufacturer supporting processor, may lead to patient Casualties.
- Do not touch the rubber and insert the tube with sharp objects, otherwise the equipment will be damaged.
- When this product is used with laser equipment, avoid potential eye damage to the operator and protect it in advance.
- When this product is used with the laser, the laser will reflect back and cause an abnormal temperature rise at the head end. Please note that the head end is not too close to the reflection source when operating.
- Do not use a high frequency electric knife for electric burning surgery and argon plasma coagulation when using this product, because this product cannot be used together with a high frequency electric knife and argon plasma coagulation.

- When this product is used with laser equipment, operators, maintenance personnel or personnel who need to look directly or observe for a long time should wear laser protective goggles with appropriate protection strength according to the laser equipment and power used.
- In the process of examination or surgery, the doctor will match other auxiliary instruments according to the actual situation of the patient. Usually, the instruments used together with this product is biopsy forceps.
- Do not hit, bump, or drop the tip, insertion, bending, operation, or display of the endoscope. At the same time, do not forcefully bend, pull or twist the tip, insertion, bending, operation, or display. Otherwise, the endoscope will be damaged. And it may cause the endoscope parts to fall off into the patient's body.
- Do not perform bending operations forcibly or suddenly. Do not forcibly pull, twist or rotate angled bends. Otherwise, the patient may be injured. It will also result in the inability to straighten the bend during the inspection.
- Video endoscope may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances.
- Except of the cable supplied by manufacturers, the others use of accessories and cable may result in increased product launch interference or low immunity.
- When a clear endoscopic image is not obtained, do not manipulate the bending part, insert or pull out the insertion part of the endoscope. Otherwise, the patient may be injured.
- After using this device, clean, disinfect and store it according to the instructions given in Chapter 5 to Chapter 7. Using improperly or incompletely reprocessed or stored device may cause patient cross-contamination and/or infection.
- Don't insert or withdraw the endoscope's insertion tube while the bending section is angulated. Patient injury can result.
- Bending section can be bent upward or downward only. The operator should

take full account of the bending direction of the endoscope while doing the insertion and withdrawal. Do not insert or withdraw to right or left direction with excessive force. Otherwise, patient injury can result.

- Don't apply suction, operate the bending section, insert or withdraw the endoscope's insertion tube without viewing the clear endoscopic image. Otherwise, patient injury can result.
- Don't insert or withdraw the insertion tube with forcefully or suddenly. Otherwise, patient injury can result.
- When it is difficult to insert the endoscope, do not force inserting the endoscope but stop the endoscopy. Forcible insertion can result in patient injury.
- If liquid has infiltrated into endoscope, do not use the endoscope, because its internal structure may be damaged.
- Do not keep close observation, or keep the tip of the endoscope in contact with active tissues for a long time.
- The user needs to install and use according to electromagnetism compatibility information which is attached with it.
- Portable and mobile RF communication equipment may affect the Endoscope. Avoid strong electromagnetic interference when using it, such as close to mobile phones, microwave ovens, etc.;
- Do not twist or bend the bending section with your hands. Equipment damage may result.
- Guidance and manufacturer's declaration stated in the appendix.
- Do not make the coiled diameter of the insertion tube less than 12 cm, otherwise, it will cause damage to the equipment.
- Do not try to bend the insertion part of the endoscope forcibly. Otherwise, the insertion part may be damaged.
- Do not bump the tip of the insertion part, especially the surface of the objective lens at the tip, otherwise it will cause framing failure.
- Do not squeeze the bending section forcefully. The covering of the bending section may stretch or break and cause water leaks.
- Electromagnetic interference may occur on this instrument near equipment marked with the following symbol or other portable and

mobile RF (Radio Frequency) communications equipment such as cellular phones. If electromagnetic interference occurs, mitigation measures may be necessary, such as reorienting or relocating this instrument, or shielding the location. 

- In order to check the electromagnetic interference from other devices (any devices like X-ray machine, CT, MRT , etc or equipment components of other system), the system should be observed in the operating environment whether it is working properly.

Endoscopic images disappear or freeze precautions

- If endoscopic image disappears unexpectedly, or frozen image cannot recover to real-time display working state, immediately stop using the endoscope, and pull out the endoscope insertion tube as specified in Section 9.2 "Withdrawal of the Endoscope with any Abnormality". Continuously operate the endoscope accessories and bending section, doing suction and/or withdraw or insert the endoscope may lead to patient injury.
- Please follow these precautions when operating these instruments. Otherwise, the endoscope image may disappear unexpectedly, or in the inspection process frozen image can not be restored.
 - Do not bend, strike or twist the insertion part and the operation section. It may damage the endoscope, causing leakage and / or damage to internal components, such as CMOS, or circuitry.
 - Before the endoscope is inserted through the mouth, if necessary, put an oral pad or guide blade into the mouth of the patient in order to prevent the patient from accidentally biting the insertion part. Chewed insertion part may cause the cable damage and broken the insertion tube.
 - Ensure that the monitor is removed before the endoscope is immersed. Otherwise, water will enter into the monitor, and may cause a short circuit. This will cause damage to the monitor.
 - If there are continuous air bubbles coming from the endoscope during leak test. Don't use the endoscope. Otherwise, water may enter the instrument and cause a short circuit. This will cause damage to the CMOS and circuit.

Examples of inappropriate handling

Details on clinical endoscopic technique are the responsibility of trained specialists. Patient safety in endoscopic

examinations and endoscopic treatment can be ensured through appropriate handling by the physician and the medical facility. Examples of inappropriate handling are given below.

- Applying suction with the distal end in prolonged contact with the mucosal surface, with higher suction pressure than required or with prolonged suction time may cause bleeding and/or lesions.
- Inserting, withdrawing without a clear endoscopic image may cause patient injury.
- Inserting or withdrawing the endoscope, applying suction or operating the bending section without a clear endoscopic image may cause patient injury.

LED light emitting diode

LED is built in the device, which is basically safe for eyes. But looking directly at the light emitted by the LED for a long time may damage the eyes. Please read this manual carefully.

WARNING

- If the control or adjustment device cannot be used according to the instructions, or the operation of each step can be carried out, harmful radiation exposure may be caused.





Note

- Avoid direct or scattered light exposure to eyes or skin.
- The endoscope uses LED light source to illuminate, and there is a possibility of failure. It is recommended to prepare the same or similar spare equipment before the operation.

Chapter1 Checking the package contents

The product is only intended to use with the Video Endoscope Processor manufactured by Vision. Match all items in the package with the components shown below. Inspect each item for damage. If the instrument is damaged, a component is missing or you have any questions, do not use the instrument; immediately contact our company. This instrument was not disinfected before shipment.

Before using this instrument for the first time, clean, disinfect it according to the instructions given in Chapter 5 "Cleaning and disinfection: General Policy" to Chapter 7 "Cleaning and disinfection Procedures".

Items	Picture	function
Video Endoscope(with channel)		When used with Video Endoscope Processor, it can be used for clinical examination, diagnosis and adjuvant therapy.
Video Endoscope(without channel)		When used with Video Endoscope Processor, it can be used for clinical examination, diagnosis and adjuvant therapy.
1. Waterproof cover 2. Balanced valve		1. Connect the waterproof lid to the handle during cleaning, soaking, and gas sterilization. 2. Install the balance valve on the leakage detector before shipment to balance the air pressure inside and outside the mirror during transportation. The balance valve must be removed before the patient examination and immersion into the cleaning fluid.
Biopsy valve		Sealing the pliers crossing.
Suction valve		Connect the negative pressure machine and the suction port of the operation department for suction operation.
Leakage detector		Connect the operating part to test the tightness of the instrument.

Suction cleaning adapter		Cleaning and disinfection connection.
Channel cleaning brush		Brushing the working channel of endoscope.
Channel port cleaning brush		Brushing of Clamp Port.
User manual		Follow instructions for use

Fig 1.1

Chapter 2 Instrument Nomenclature and Specifications

The constitute of the cooperative-used Video Endoscope Processor please refer to the user manual of the processor.

2.1 Constitute

Distal end, bending section and Insertion Section is an APPLIED PART, see details Fig 2.1、 Fig 2.2 BF series constitute.


Fig 2.1 BF series constitute (without channel)

No.	Part Name	Function
1	Dimmer push button	It is used to adjust the brightness of the leading led, which can be divided into five levels of dimming.
2	Photo/video recording button	1. Press this button once in the power on preview mode to record the static picture. 2. Press and hold this button for 2 seconds in the power on preview mode to enter the video mode. 3. In video mode, long press this button for 3 seconds to exit video mode and enter real-time video.
3	Bending section	Adjust the position of distal end via UP/DOWN Angulations Control Knob.
4	Distal end	Endoscope lens end with LED light.
5	Insertion part	It extends into the natural cavity of human body.
6	Conical sheath	Protect the Insertion tube.
7	Operation section	Handle and operating Insertion tube.
8	Up / Down angulations' control knob	Turn the control knob to 'U' direction to make the distal end bend up; Turn the control knob to 'D' direction to make the distal end bend down.

Form 2.1


Fig 2.2 BF series constitute (with channel)

No.	Part Name	Function
1	Suction Port	This valve is depressed to activate suction. This valve is used to remove any liquid and debris that block the line of sight.
2	Dimmer push button	It is used to adjust the brightness of the leading led, which can be divided into five levels of dimming.
3	Photo/video recording button	1. Press this button once in the power on preview mode to record the static picture. 2. Press and hold this button for 2 seconds in the power on preview mode to enter the video mode. 3. In video mode, long press this button for 3 seconds to exit video mode and enter real-time video.
4	Clamp Port	The Clamp Port is a standard Ruhr connector, from which medical instruments can enter during treatment.
5	Bending section	Adjust the position of distal end via UP/DOWN Angulations Control Knob.
6	Distal end	Endoscope lens end with LED light.
7	Insertion part	It extends into the natural cavity of human body.
8	Conical sheath	Protect the Insertion tube.
9	Operation section	Handle and operating Insertion tube.
10	Up / Down angulations' control knob	Turn the control knob to 'U' direction to make the distal end bend up; Turn the control knob to 'D' direction to make the distal end bend down.

Form 2.2

2.2 Specifications

Note

- Installed in the distal end of endoscope, the image sensor is to acquire images, and convert to electrical signals by the cable transmission to video image processing circuit. The real-time image is displayed in the LCD screen. This endoscope does not use the optical fiber bundle.

Environment specifications

Normal working conditions	Environment temperature:	5°C~40°C
	Atmospheric pressure:	700hPa~1060hPa
	Relative humidity:	30%~85% (Non-Condensing)
Transportation and storage	Environment temperature:	-10°C~+55°C
	Atmospheric pressure:	500hPa~1060hPa
	Relative humidity:	≤93% (Non-Condensing)

Form 2.3

Safety Standard




Note

- In accordance with the CE standard classification, this product is class IIa. This product is classified as follows according to IEC60601-1:




Classification of protection against electric shock	TYPE BF APPLIED PART
Classification of protection against hazards of explosion	Not suitable for use in places with flammable anesthetic gas
Classification of protection against ingress of liquid	IPX7
Classification of running mode	Continuous operation equipment
Mobile level	Handheld

Form 2.4



Product specification

Model specification									
Model		BF62	BF58	BF52	BF45	BF42	BF30	BF25	BF20
Optical system	Field of view	120°							
	Direction of view	0°							
	Depth of field	3mm~200mm							
	Resolution	400*400							
Insertion tube outer	Diameter	Ø 6.2mm±0.3 mm	Ø 5.8mm±0.3 mm	Ø 5.2mm±0.3 mm	Ø 4.5mm±0.3 mm	Ø 4.2mm±0.3 mm	Ø 3.0mm±0.3 mm	Ø 2.5mm±0.3 mm	Ø 2.0mm±0.3 mm
	Distal end enlarge picture	 1. CMOS camera lens 2. LED light 3. Suction Channel exit							 1. CMOS camera lens 2. LED light
	Working length	600±10%mm							
Working channel	Minimum suction channel inner diameter	3.2±0.2mm	3.0±0.2mm	2.6±0.2mm	2.2±0.2mm	2.0±0.2mm	1.2±0.2mm	1.2±0.2mm	/
	Direction from which endo-therapy accessories enter and exit the endoscopic image								/
	Channel Opening Type								/
Bending section	Angulation range	Up180°,Down130°							

Form 2.5

Model specification									
Model		BF62H	BF58H	BF52H	BF45H	BF42H	BF30H	BF25H	BF20H
Optical system	Field of view	120°							
	Direction of view	0°							
	Depth of field	3mm~200mm							
	Resolution	Not less than 720*720							
Insertion tube outer	Diameter	Ø 6.2mm±0.3 mm	Ø 5.8mm±0.3 mm	Ø 5.2mm±0.3 mm	Ø 4.5mm±0.3 mm	Ø 4.2mm±0.3 mm	Ø 3.0mm±0.3 mm	Ø 2.5mm±0.3 mm	Ø 2.0mm±0.3 mm
	Distal end enlarge picture	 1. CMOS camera lens 2. LED light 3. Suction Channel exit							 1. CMOS camera lens 2. LED light
	Working length	600±10%mm							
Working channel	Minimum suction channel inner diameter	3.2±0.2mm	3.0±0.2mm	2.6±0.2mm	2.2±0.2mm	2.0±0.2mm	1.2±0.2mm	1.2±0.2mm	/
	Direction from which endo-therapy accessories enter and exit the endoscopic image								/
	Channel Opening Type								/
Bending section	Angulation range	Up180°,Down130°							

Form 2.6

Model specification						
Model Parameter		BF52P	BF45P	BF42P	BF30P	BF28P
Optical system	Field of view	120°				
	Direction of view	0°				
	Depth of field	3mm~200mm				
	Resolution	400*400				
Insertion tube outer	Diameter	Ø 5.2mm±0.3mm	Ø 4.5mm±0.3mm	Ø 4.2mm±0.3mm	Ø 3.0mm±0.3mm	Ø 2.8mm±0.3mm
	Distal end enlarge picture	 1. CMOS camera lens 2. LED light 3. Suction Channel exit				
	Working length	600±10%mm				
Working channel	Minimum suction channel inner diameter	2.8±0.2mm	2.4±0.2mm	2.2±0.2mm	1.6±0.2mm	1.4±0.2mm
	Direction from which endo-therapy accessories enter and exit the endoscopic image					
	Channel Opening Type					
Bending section	Angulation range	Up180°,Down130°				

Form 2.7

Chapter 3 Preparation and Inspection

Before each case, prepare and inspect this instrument as instructed below. Inspect other equipment to be used with this instrument as instructed in their respective instruction manuals. If the abnormality are suspected after inspection, follow the instructions given in Chapter 9, "Troubleshooting". If this instrument still malfunctions, do not use it. Return it to our company for repair as described in Section 9.3, "Returning the endoscope for repair".



Warning



- Before video endoscopy, it is necessary to check whether the performance of the video endoscope is safe to use, such as the display screen can clearly receive real-time images, and can take and record images in real-time.
- The use of a malfunctioning endoscope will endanger the safety of the patient and result in aggravated equipment damage.
- This instrument was not cleaned, disinfected before shipment. Before using this instrument for the first time, reprocess it according to the instructions given in Chapter 5, "Cleaning and disinfection: General policy" to Chapter 7, "Cleaning and disinfection procedure" instructions for cleaning and disinfection.

3.1 Preparation of the equipment

Before use, please prepare the equipment and personal protective equipment listed in the table below, such as eye wear, face mask, moisture-resistant clothing and chemical-resistant gloves. Personal protective equipment should be the appropriate size, and long enough to avoid skin exposed. Please refer to the instruction manual of the respective device.

Before using this instrument for the first time, reprocess it according to the instructions given in Chapter 5, "Cleaning and disinfection: General policy" to Chapter 7, "Cleaning and disinfection procedure" instructions for cleaning and disinfection .

Items	Picture	function
Video Endoscope(with channel)		When used with Video Endoscope Processor, it can be used for clinical examination, diagnosis and adjuvant therapy.
Video Endoscope(without channel)		When used with Video Endoscope Processor, it can be used for clinical examination, diagnosis and adjuvant therapy.

1. Waterproof cover 2. Balanced valve		1. Connect the waterproof lid to the handle during cleaning, soaking, and gas sterilization. 2. Install the balance valve on the leakage detector before shipment to balance the air pressure inside and outside the mirror during transportation. The balance valve must be removed before the patient examination and immersion into the cleaning fluid.
Biopsy valve		Sealing the pliers crossing.
Suction valve		Connect the negative pressure machine and the suction port of the operation department for suction operation.
Leakage detector		Connect the operating part to test the tightness of the instrument.
Suction cleaning adapter		Cleaning and disinfection connection.
Channel cleaning brush		Brushing the working channel of endoscope.
Channel port cleaning brush		Brushing of Clamp Port.
User manual		Follow instructions for use

Form 3.1

3.2 Inspection of Video Endoscope

Clean and disinfect as described in Chapter 5 "Cleaning and Disinfection : General Policy" to Chapter 7 " Cleaning and Disinfection Procedures", other steps are as follows:

Inspection of Video Endoscope

1. Inspect the operation section for excessive scratching, deformation, loose parts or other abnormality.
2. Inspect the protective boot and the insertion tube near the protective boot for bends, twists or other abnormality.
3. Inspect the external surface of the entire insertion tube including the bending section and the distal end for dents, bulges, swelling, scratching, holes, sagging, transformation, bends, and adhesion of foreign bodies, dropout of parts, any protruding objects or other abnormality.
4. Holding the operation section with one hand, the entire surface of the insertion part and the light stroking
in two directions with the other hand. (See fig3.1). Confirm no objects or wire from the insert protrusion inserted into the pipe and make sure that no abnormal hardness.


Fig 3.1

5. Inspect the covering of the bending section for sagging, swelling, cutting, holes, or other abnormality.
6. Using both hands, bend the insertion tube of the endoscope into a semicircle. Then, moving your hands as shown by the arrows, confirm that the entire insertion tube can be smoothly bent to a semicircle (Fig 3.2)


Fig 3.2

7. Gently hold the midpoint of the bending section and a point 10 cm from the distal end. Push and pull gently to confirm that the border between the bending section and the insertion tube is not loose.
8. Inspect the distal end of the endoscope's insertion tube for scratching, cracks, stains, gaps around the lens or other abnormality.

Inspection of the bending functions

Perform the following inspections while the bending section is straight.

- If Angle control handle movement is too loose or not smooth, or the bending section does not angulated smoothly, the bending mechanism may be abnormal. In this case, do not use the endoscope because it may be impossible to straighten the bending section during an examination.

1. Slowly rotate the up/down angle control handle to the end points in each direction, and then return to their natural positions. Confirm that the bending operation of the bending part is smooth and normal, can reach the maximum amplitude, and return to the natural position.
2. As shown below in Fig 3.3, turn the UP/DOWN angulation control knob slowly to its neutral position, confirm that the bending section returns smoothly to an approximately straight condition.


Fig 3.3

3.3 Inspection of accessories

Inspection of Biopsy Valve

Inspect the Biopsy Valve for tear, severe deformation, or damage.

Warning

- Dysfunction Biopsy Valve will not only compromise the suction effect, but also result in leakage, or spread of patients' debris, bringing the risk of infection.

Note

- If the Biopsy Valve deterioration, or abnormal, replace the Biopsy Valve.

3.4 Accessories Installation to the endoscope

Installing Biopsy Valve

Warning

- If the Biopsy Valve is not properly attached to the instrument channel port, it will reduce the effect of endoscope suction, but also result in leakage, or spread of patients' debris, bringing the risk of infection.
- If the Biopsy Valve can not achieve sealing effect, it will lead to failed suction without pressure difference.

1. Install the Biopsy Valve to the instrument channel port, and then push to the end, until it is fixed to the body (Fig 3.4)
2. Ensure that the Biopsy Valve is firmly fixed on the instrument channel port.


Fig 3.4

3.5 Inspection of the endoscopic

Inspection of the suction function

- Set the pressure of suction device within -34 kPa to 0 kPa. Excessive pressure may make it difficult to stop suction.
- If the suction valve does not operate smoothly, detach it and reattach it, or replace it with a new one.
- Stop using the endoscope while the suction valve is not working properly, because it may be impossible to stop suction, which could cause patient injury. If the reattached or replaced suction valve still fails to operate smoothly, the endoscope may be malfunctioning; stop using it and contact our company.

1. Adjust the suction pressure to the same level as it will be during the procedure.
2. Immerse the distal end of the insertion tube in sterile water. Press the suction valve and confirm water is continuously aspirated into the suction bottle of the suction device.
3. Release the suction valve. Confirm that suction stops and the valve returns to its original position.
4. Remove the distal end of the endoscope from the water. Depress the suction valve and aspirate air for a few seconds to remove any water from the working channel.

Inspection of the water feeding function

1. Insert a syringe filled with sterile water into the instrument channel port (Biopsy Valve) and depress the plunger.
2. Confirm that water is emitted from the distal end of the insertion tube.
3. Insert the cleaning brush through the suction port and Clamp Port. Make sure the cleaning brush extends smoothly from the distal end. At the same time, make sure that there is no foreign article at the distal end.

4. Make sure that the cleaning brush can be pulled out smoothly from the Clamp Port.

Note

- The syringe must be inserted fully and held perpendicular to the valve for proper operation. Angled or incomplete insertion may result in fluid leakage.
- Do not depress the suction valve during feeding. If the suction valve is depressed during water feeding, water will be aspirated into the suction tube and not emitted from the endoscope's distal end
- If fluid is not emitted from the endoscope's distal end, flush air through the channel.
- Keep eyes away from the distal end when inserting the cleaning brush. A cleaning brush reaches out from the distal end of endoscope may cause eye injury.


Fig 3.5

Inspection of the working channel

- Keep your eyes away from the distal end when inserting endo-therapy accessories. Extending the endo-therapy accessory from the distal end could cause eye injury.

1. Insert the channel cleaning brush through the instrument channel. Confirm that the channel cleaning brush extends smoothly from the distal end and that no foreign objects come out of the distal end.
2. Confirm that the channel cleaning brush is withdrawn smoothly from the working channel.

3.6 Equipment Connections

The connection mode and display setting of this product depends on the Video Endoscope Processor (Type BF Applied Part) you selected. And the detail operation description please refer to the user manual of the Processor. The product can choose the wired connection mode and the wireless connection mode.

Wireless Connection Mode

1. Remove the wireless transmitter from the charging card slot of the electronic endoscope image processor.
2. Press the black mark point on the right side of the wireless transmitter white balance button at the black mark point at the interface of the endoscope connector. Hearing the click, the wireless transmitter connects to the endoscope operator successfully, as shown in the figure below.


Fig 3.6a


Fig 3.6b

Wired connection mode

1. Insert one end of the extension cable (Lemo plug) into the video endoscope processor output interface, the Lemo plug (mother plug). Clicks were heard after insertion.
2. Place the extension cable at the other end (endoscopic manipulation handle connection end) from top to bottom above the operant. Align the black dots on the to turn axis to the black dots above the maneuver (endoscopic maneuver handle connection end) rotate from right to left (endoscopic maneuver handle connection end) by gently clasp the (endoscopic maneuver handle connection end) to the maneuver, and after hearing a "click", the illustration follows:


Fig 3.7

3.7 Inspection of distal end

- Do not stare directly into the light beam. This may result in eye injury.
- This product is of LED lighting, LED inside the equipment, not replaceable.
- If the LED light does not work, do not replace it by yourself, please contact manufacturer.

1. Connect the processor to the handle according to the guidance in the processor's manual. The power indicator lights, and endoscope image appears on the LCD screen.
2. Confirm that Examination light is emitted from the LCD light on distal end of the endoscope. (Fig 3.8)


Fig 3.8

3. Put the endoscope's distal end to a suitable object, and move between 1cm and 6 cm. Confirm that the brightness of the image on the LCD screen remains constant. (Fig 3.9)


Fig 3.9

Note

- If you feel abnormal endoscopic image color, refer to "White Balance Adjustment."

4. Put the distal end of endoscope to your palm to observe, confirm that endoscopic image is free from noise, blur, fog or other abnormality.

Warning

- If it is loose and/ or is not smooth when you rotate the monitor and/ or adjust the tilt angle of the LCD screen, the structure may be abnormal. In this case, do not use the endoscope, because it may not move these structures during the operation.
- Slowly rotate the monitor and tilt the LCD screen to adjust to its end. Confirm that no other abnormalities like difficult to move during the move.

3.8 Image check

1. Entering the preview state by default, that is, the real-time image display interface of the electronic endoscope image processor is implemented. At this time, the real-time image should be observed normally, and the image is clear, without blur and incomplete image abnormalities.

2. Hold the mirror with your left hand and put the thumb on the Angle lever, as shown in the Fig3.3. At this point, the image is seen in the same direction as the operator. If there is any abnormality, please stop using it immediately and contact the manufacturer.

Warning

- Do not use this video endoscope if no real-time endoscopic image is observed. Otherwise, it may lead to patient injury.

3.9 Check the freezing function

Warning

- Do not use this video endoscope if no real-time endoscopic image is observed. Otherwise, it may lead to patient injury.

3.10.1 Under the main menu, press the "function switch" button. Enter the function switch setting interface, and select "freeze" below (as shown in Fig 3.10). Then return to the preview interface.


Fig3.10

- Press the light button switch. Confirm the real-time endoscopic image freezing. A frozen icon appears in the display interface, as shown in the figure below.


Fig3.11

- b. Press the light button switch again to confirm that the frozen image returns to the real-time image. And the freeze icon disappears.

Chapter 4 Operation

This instruction manual does not explain or discuss any clinical endoscopic surgery. Only introduce the basic operations and precautions related to the use of this endoscope. Therefore, the operator of this endoscope must be a physician or medical staff under the supervision of a physician, and they must have received sufficient training in clinical endoscopy techniques.

WARNING

- That before each use, the compatibility of the endoscopic equipment with any accessories and/or energized endotherapy devices should be checked according to any criteria for safe use defined in the instructions for use.
- In the process of examination or surgery, the doctor will match other auxiliary instruments according to the actual situation of the patient. Usually, the instruments used together with this product is biopsy forceps.
- Distal end surface and the light emission part temperature exceeds 41 ° C. temperatures over 41 ° C (105.8 ° F) may cause mucosal burns Always use the minimum time and suitable distance necessarily for adequate viewing. Whenever possible, avoid close stationary viewing and do not leave the distal end of the endoscope close to the mucous membrane for a long time.
- Never insert or withdraw the endoscope under any of the following conditions. Otherwise, it may result patient's injury.
 - Insertion or withdrawal while the endo-therapy accessory extends from the distal end of the endoscope.
 - Insertion or withdrawal while the bending section are in an shape angle.
 - Insertion or withdrawal with an excessive force, or forcible insertion or withdrawal.
- The bending section has up/down deflection only, but not right/left deflection. Take full account of the bending direction when doing insertion or withdrawal .Do not insert or withdraw with an excessive force in the left and right direction. Otherwise, the patient will be injured.
- If any of the following situations occurs during an examination, immediately stop the examination and withdraw the endoscope from the patient as described in Section 9.2“Withdrawal of the Endoscope with any Abnormality”.
 - The function of the endoscope is found to be abnormal.

——The endoscopic image on the LCD monitor disappears unexpectedly or freezes.

——The angle control handle is locked.

——The angle control function is abnormal.

Continuous use of endoscope under the situation stated above may cause the patient's injury.

- If an abnormal endoscopic image/function occurs and returns to normal condition by itself, the endoscope may have malfunctions. Continuous use of such an endoscope may cause repetition of the abnormality, which may not return to normal. In this case, stop the procedure immediately and slowly withdraw the endoscope while viewing the endoscopic image. Otherwise, patient injury can be resulted.
- That before each use or after a change of viewing modes/settings, the operator should check to ensure the view observed through the endoscope provides a live image (rather than a stored one) and has the correct image orientation. Image inspection requires the Video Endoscope Processor, please refer to the Video Endoscope Processor user manual.

4.1 Insert the endoscope

Holding and manipulating the endoscope

The operation section of the endoscope is designed to be held in the left hand. The UP/DOWN angulation control knob can be operated using the left thumb. The right hand is free to operate the insertion tube. (Fig 4.1).


Fig 4.1

Insertion of the endoscopic accessories

Precautions

- When inserting the endoscope, first lubricate the insertion tube with medical grade lubricant to ensure the lowest friction when inserting the endoscope into the patient.
- In order to prevent the patient from biting into the insertion section during the inspection, it is recommended to put the tooth pad in the patient's mouth before inserting the endoscope.
- Do not apply olive oil or products containing petroleum-based lubricants (e.g. aseline). These products may cause stretching and deterioration of the bending section's covering.
- Do not allow the insertion tube to be bent within a distance of 10 cm or less from the junction of the protective cover. Insertion tube damage can occur (Fig 4.2).


Fig 4.2

1. If necessary, apply a medical-grade water-soluble lubricant to the insertion tube.
2. If necessary, place the tooth pad between the patient's teeth and outside the patient's mouth.
3. Always view the endoscopic image when passing the distal end of the endoscope from the mouth or nose to the pharynx.

Observation of the endoscopic image

Keep eyes on endoscopic image all the time when operating the angulation control knob.

Angulation of the distal end

- Avoid forcible or excessive angulation, as this imposes load on the wire controlling the bending section and may cause stretching or tearing of the wire and trouble in the action of the bending section.
- Operate the UP/DOWN angulation control lever as necessary to guide the distal end for the insertion and observation.

Suction

- Aspirating thick fluids may lead to channel or valve clogging. If the suction valve gets clogged and suction cannot be stopped, disconnect the suction tube from the suction connector. Turn the suction pump OFF, detach the suction valve and remove thick fluids.
- When aspirating, maintain the suction pressure at the lowest level necessary to perform the procedure. Excessive suction pressure could cause aspiration of and/or injury to the mucous membrane. In addition, patient fluids could leak or spray from the biopsy valve, posing an infection-control risk.
- When aspirating, attach the Biopsy Valve to the instrument channel port. The uncapped instrument channel port can reduce the efficacy of the endoscope's suction system.

1. Depress the suction valve to aspirate excess fluid or other debris obscuring the endoscopic view. (Fig 4.3)


Fig 4.3

Using a Syringe to flush the working channel

- If the direction is not properly installed, or the installation is not complete, the liquid may leak from the instrument channel port.

1. Securely insert a syringe filled with lavage fluid (e.g. saline) into the instrument channel port and press the plunger to feed lavage fluid. (Fig 4.4(a))
2. With the syringe attached, slowly withdraw the plunger to aspirate lavage fluid. (Fig 4.4(b))


Fig 4.4

4.2 Withdraw the endoscope

- If blood unexpectedly adhere to the surface of the insertion tube

of the withdraw endoscope, check the conditions of the patient carefully.

- If the endoscope cannot be withdrawn from the patient smoothly, do not attempt to forcibly withdraw it. Forcibly withdrawing the endoscope may cause patient injury.

1. Aspirate accumulated blood, mucus or other debris by depressing the suction valve.
2. Carefully withdraw the endoscope while observing the endoscopic image. Remove the guide blade from the patient's mouth.
3. Press the power button to turn off the whole equipment.

Chapter 5 Cleaning And Disinfection: General Policy

5.1 Key Points of cleaning and disinfection

The medical literature reports incidents of patient cross contamination resulting from improper cleaning or disinfection. It is strongly recommended that all individuals responsible for cleaning, disinfection should thoroughly understand and follow the instructions in this manual and all manual instructions for auxiliary equipment.

- occupational health and safety regulations in local hospital
- cleaning, disinfection procedures for each individual equipment
- the structure and operation method of endoscopic equipment
- usage of pertinent chemical agents

Please make judgments from a professional point of view on the type and conditions of methods of cleaning and disinfection.

Please refer to the cleaning and disinfection standards of the region/country where the product is used. Contact the local health agency to learn about local standards and regulations.

This product has not been cleaned, disinfected before leaving the factory. It must be cleaned, disinfected before use. Please refer to the method recommended by the manufacturer for cleaning and disinfecting.

Since the product can only withstand 1000 times of peracetic acid disinfection, please record the number of disinfection at each time to avoid accidents caused by product failure after disinfection times exceeding 1000 times.

5.2 Precaution

Warning

- After each case, if the endoscopy equipment is not effectively cleaned, disinfected, the safety of the patient will be endangered. In order to reduce the risk of cross-infection, the endoscope must be thoroughly cleaned according to the contents of Chapter7 "Cleaning and Disinfection Procedures" after each case inspection, and then disinfected. The outer surface of the endoscope should be cleaned, disinfected.
- If the endoscope is not cleaned meticulously, effective disinfection may not be possible. Before using the endoscope

cleaning machine, it should be confirmed that it can clean, disinfect the endoscope. If you are not sure whether the endoscope cleaning machine can clean and disinfect the endoscope, please contact the manufacturer of the endoscope cleaning machine to confirm the specific instructions or connector information. If the endoscope is not cleaned and disinfected properly, it may cause infection risk to the patients who use the endoscope next time.

- The manual specifies fitting/usable and unfitting/unusable detergent, disinfectant, sterilant and another chemical agent and washer/disinfector for the endoscope. As for the detergent, disinfectant, sterilant and other chemical agent and washer/disinfector that are not shown in this operation manual, please contact our company, service center qualified by our company or Sales Representatives of our company. Improper use of chemical agent or washer/disinfector may cause deterioration of endoscope, parts of the endoscope fall off and cause danger to patient's health.
- The chemical agents used for cleaning, disinfection maybe dangerous, so personal protective equipment should be worn to avoid being harmed by dangerous chemical agents and potentially infectious substances. When cleaning, disinfecting , please wear appropriate personal protective equipment, such as goggles, mask, waterproof clothing and protective gloves.
- Be sure to rinse thoroughly with chemicals. Clean and disinfect the external surface of the endoscope thoroughly with (Aq. Dest) water to remove any residual chemicals.
- Be sure to perform the water leakage test before the manual cleaning. If leakage is detected, stop using the endoscope. Just before each procedure, check if the endoscope needs cleaning and disinfection. If necessary, clean, disinfect it by following the instructions given in this manual.
- With the cleaning and disinfection methods stated in this instruction manual, prions, which are considered to be the pathogenic substance of the Creutzfeldt-Jakob disease (CJD) cannot be destroyed or inactivated. When using this instrument on a patient with CJD or variant Creutzfeldt-Jakob disease (vCJD), be sure to use this product for such patient only and/or

immediately dispose of this product after use in an appropriate manner. For methods to handle CJD, please follow the respective guidelines in your country.

Chapter 6 Suitable methods of cleaning, disinfection and chemical preparation

6.1 Compatibility Summary

Certain materials and structure of our endoscope may not be compatible with some cleaning, disinfection methods.

To judge whether the method is effective, we base on the following two points:

- Microbiological effectiveness
- Material durability

Microbiological effectiveness

The standard methods described in this manual may be useful for cleaning and disinfection of this equipment.

Material durability

This product can be cleaned, disinfected for up to 1000 times with this method.

Cleaning and disinfection methods

Choose appropriate cleaning and disinfection methods in compliance with the recommendations of your infection control committee and all national and local hospital guidelines and policies.

Our company would recommend automatic cleaning methods. Please contact your local infection control agency to determine the appropriate cleaning method for the equipment. Please contact the local infection control institutions to determine suitable cleaning methods specified by the institutions.

Monitoring

In accordance with national and local regulations or requirements, inspect and confirm cleaning and disinfection methods for reusable instruments, record the regular interval and cleaning and disinfection results.

Although there is no available biological indicator to confirm the disinfection process, test strips can be used to monitor the concentration of disinfectants.

Monitoring the concentration of preparations in daily use can ensure that the solution does not fall below the effective concentration.

In order to monitor the process, appropriate biological or chemical indicators should be used.

List of compatible methods and chemical preparations

The methods and preparations for washing, disinfection, listed in this table, have been previously verified on endoscopes and their accessories. So that it's suitable for endoscopes and fittings and joints of leakage detector.

Recommended cleaning and disinfection solution	
Cleaning Reagent	EmPower multienzyme cleaning agent
Disinfection Solution	Peracetic acid disinfection(Dilute the peracetic acid with pure water at the ratio of 1:22.1 to make the final concentration between 2380mg/L-2814 mg/L.)

Form 6.1

Warning

- Do not use steam disinfection (high-temperature autoclaving) for endoscopes. Steam disinfection (high temperature autoclaving) may damage endoscopy

6.2 Cleaning Reagent

Use cleaning reagent which manufacturer recommend EmPower multienzyme cleaning agent.

Warning

- The cleaning reagent must comply with national and local regulations or requirements and has been certified. Inadequate cleaning agents will bring the risk of infection.
- Do not reuse cleaning reagent.

Caution

- The cleaning reagent must be approved by manufacturer for compatibility with the endoscope and certified by authority organization. Incompatible cleaning reagent may bring great damage to the endoscope and accessories.

6.3 Disinfection Solution

Peracetic acid disinfection: According to the manufacturer's instruction manual, in general, 2380mg/L-2814mg/L peracetic acid disinfection solution is suitable for the Endoscope produced by our company. If the disinfectant is reused, the test paper should be used to regularly test its effectiveness. Do not use solutions beyond the expiration date.

Warning

- Do not use endoscope only disinfected by alcohol.

Caution

- The disinfection solution must be confirmed by our company for compatibility with the endoscope and certified by authority organization. Uncompatible disinfection solution may pose great damage to the endoscope and accessories.
- Refer to the manufacturer's instructions for the concentration,

temperature and immersion time. Do not exceed the level specified by the manufacturer.

6.4 The Information Of Washer

Name	Automatic soft endoscope washer-disinfector
Model	Rider 60A
Manufacturer	SHINVA MEDICAL INSTRUMENT CO., LTD.

Form 6.2

6.5 Rinsing Water

Use water that has been processed (e.g. filtered) to improve its microbiological quality. When non-sterile water is used after manual disinfection, wipe the endoscope and flush the channels with 75%–95% ethyl alcohol, then air-dry all internal channels to inhibit the growth of residual bacteria. The selected water quality should be tested by national and local regulations or requirements, and certified by Infection Control Committee in the hospital.

- Due to possible chlorination, do not use tap water to wash the equipment.
- Do not reuse rinsing water.

Chapter 7 Cleaning and Disinfection Procedures

This chapter covers the cleaning and disinfection steps of endoscopic equipment related to patient safety and health.

This chapter focuses on the endoscope maintenance and maintenance methods from the endoscope preparation, cleaning, and disinfection.

WARNING

- Do not perform maintenance while the operator is using the device for the patient, which may harm the patient and the operator.
- Medical institutions need to have the operation specifications and related equipment for endoscopic cleaning, disinfection.
- ALL channels of the endoscope MUST be cleaned and high-level disinfected during every reprocessing cycle, even if the channels were not used during the previous patient procedure. Otherwise, insufficient cleaning and disinfection of the endoscope may pose an infection risk to the patient and/or operators performing the next procedure with the endoscope.

CAUTION

- Do not coil the endoscope's insertion tube with a diameter of less than 12 cm. The endoscope can be damaged if coiled too tightly.

7.1 Required Cleaning and Disinfection Equipment



Preparation of the Equipment

Prior to cleaning and disinfection, prepare the equipment shown in Fig 7.1.

WARNING

- Use basins which are at least 40 cm × 40 cm (16" × 16") in size and deep enough to allow the endoscope to be completely immersed. Incompletely immersed can compromise the effect of cleaning and disinfection.

Items	Pictures
-------	----------

waterproof cover and balanced valve	
Biopsy valve	
Suction cleaning adapter	
Channel cleaning brush	
Channel connector cleaning brush	
Suction valve	
Leakage detector	
Personal protective equipment	Sterile water (Rinsed in sterile water)
500cm ³ (500ml)	70% ethyl or isopropyl alcohol (Rinsed in non-sterile water)
Soft brush	Clean, lint-free cloths / paper
Clean water	Sterile, lint-free cloths
Cleaning Reagent	Sterile, cotton swabs
Disinfectant solution	30cm ³ (30 ml) syringe
Pressure water gun	Pressure air gun

40cm×40cm(16"×16")Basins with tight fitting lids, 40 cm × 40 cm (16" × 16") minimum size	Basins with tight fitting lids, 25 cm × 25 cm (10" × 10") minimum size
(Aq.Dest) Deionized water (Aq. Dest)	Small basins

Form 7.1

Cleaning and disinfection Equipment Parts and Functions

For inspection of other equipment not mentioned below, refer to the instruction manual for the equipment being used.

Waterproof cover

1. When cleaning and disinfection, the Waterproof cover must be attached to the electrical connection part above the endoscope. (Fig 7.1)


Fig 7.1

2. Suction cleaning adapters are used to flush, clean, disinfect and disinfect the liquid through the injection channel and suction channel, and remove the liquid from channels(Fig 7.2)


Fig7.2

Channel Cleaning Brush

- Do not use the channel cleaning brush for cytological tissue sampling.

or other diagnostic purposes. Otherwise, it will cause patient injury, cross-infection or equipment damage.

Channel cleaning brush is used to brush the inside of the working channel. ((Fig 7.3).


Fig 7.3

Channel port cleaning brush

The Channel port cleaning brush is used to brush the endoscopic injection opening (Fig 7.4)


Fig 7.4

Leakage detector

During endoscopy leak testing, connect the cap of leakage detector with the leak valve connector of the endoscope handle and press the airbag to check for leakage (Fig 7.5).


Fig 7.5

Inspection of the Reusable Equipment

For inspection of other equipment not mentioned below, refer to the instruction manual for the equipment being used.

- All equipments mentioned below are consumable items. Whenever the slightest irregularity is suspected, use a spare instead. Using a defective equipment may make it difficult to effectively reprocess the endoscope, and cause endoscope and/or equipment damage.
- Please clean and disinfect the equipment before each use, and check the sealing condition of the equipment to ensure that the endoscope can be used normally. Refer to the following steps.

Inspection of the Pressure Compensation Cap

Confirm that it is free from scratches, flaws and debris. (See Fig 7.5)

Inspection of the Suction Cleaning Adapter

Confirm that the adapter is free from cracks, scratches, flaws, debris and/or other damage. (See Fig 7.2)

Inspection of the Channel Cleaning Brushes


1. Confirm that the brush section and the metal tip of the distal end are securely in place. Check for loosening or missing bristles.
2. Check for bends, scratches and other damage to the shaft.
3. Check for debris on the shaft and/or in the bristles of the brush.

Inspection of the Channel port cleaning brush

1. Check for loose or missing bristles (See Fig 7.4)
2. Check for debris on the shaft and/or in the bristles of the brush
3. Check for bends, scratches and other damage to the shaft.

Inspection of the leakage detector

1. Make sure that the leakage detector is without cracks, scratches, flaws, debris and/or other damage. (See Fig 7.5)
2. Press the air bag, make sure that the leak detect or release the air from the connector. (see Fig 7.6).

(leakage detector) Fig 7.6

7.2 Cleaning and Disinfection Procedures

After the procedure, clean, disinfect the endoscope according to the procedure described below.

Endoscope reprocessing summary chart


Fig 7.7

Warning

- All channels of the endoscope must be cleaned, disinfected during every reprocessing cycle, even if the channels were not used during the previous patient procedure. Otherwise, insufficient cleaning of the endoscope may bring an infection risk to the patient and/or operators using the endoscope.
- Do not skip any of the steps described below. Otherwise, it may lead to infection.
- Do not store disinfected instruments. Use it immediately after disinfection, or send it to the disinfection room. If you want to use

stored instruments, please refer to local regulations and guidelines.

- Never immerse the device in any liquid more than 60 minutes. If immersed for a long time, it will cause excessive moisture inside endoscope, cause dysfunction of endoscope.

7.3 Pre-cleaning

Water wash

WARNING

- If the endoscope is not immediately cleaned after each procedure, residual organic debris will begin to solidify and it may be difficult to effectively clean, disinfect the endoscope.

The endoscope should be pre-cleaned immediately after each examination.

pretreatment

1. After the endoscope is used, unload the display or remove the extension cable . The display / extension cord is not waterproof.

Fig 7.8

2. Cover the waterproof cover and immediately wipe off all visible dirt on the outer surface of the endoscope with lint free gauze containing cleaning solution. Lint free gauze shall be used once. Then repeatedly wash the working channel with syringe and detergent solution for 10 seconds. Cover the waterproof cover, put it into the conveying container, and then send it to the cleaning and disinfection room.


Fig 7.9

- Without a waterproof cover, poor contact with metal contacts.

7.4 Transportation of the Endoscope

Warning

Note

- Do not spill out debris or pre-cleaning liquid during transportation.
- Please use the handling equipment to transport endoscope to the cleaning, disinfection area.

Transporting within the hospital

1. Place the endoscope on the transport trolley or use a bracket to secure the device. The specific container will protect the endoscope from damage during transport.
2. Transport the endoscope and reusable components from using site to the cleaning, disinfection area.

Caution

- The equipment should be cleaned, disinfected immediately after use. Do not let the endoscopes without cleaning and disinfection overnight. If you leave the endoscopes for a long time without cleaning and disinfection, debris will become hard and difficult to remove. If the endoscopes were kept into the liquid for a long time without cleaning and disinfection, internal endoscope will become too humid and result in equipment damage.
- Do not hit the distal end of endoscope during transportation. Otherwise, it will damage the distal end.
- Do not only hold the endoscope insertion section to take up the whole endoscope. Otherwise, it will damage the endoscope.
- Please refer to Section 2.3 of the manual for transportation environment

7.5 Leakage testing

After pre-cleaning, perform leakage testing on the endoscope to ensure that it is waterproof.

Equipment needed

Please prepare the following equipment:

- Personal protective equipment
- The water basin is at least 40cm x 40cm(16" x 16") and deep enough to completely soak the entire endoscope
- Deionized water (Aq. Dest.)
- leakage detector

Caution

- Disconnect the balanced valve before immersed. Otherwise it could allow water to enter the endoscope and may damage the endoscope.
- Rotate the pressure compensation cap until it stops. If it is not fully and properly attached, the endoscope's interior will not be pressurized and accurate leakage testing will be impossible.
- During leakage testing, a continuous series of bubbles emerging from a location on the endoscope indicates a leak at that location. This means that water will be able to penetrate the endoscope. If you locate a leak, remove the endoscope from the water and contact our company.
- When connecting the pressure compensation cap connector, make sure that the inside of the leakage tester's and the outside of the water-resistant cap's leakage valve are thoroughly wipe and dry. Water remaining on either component may penetrate the inside of the water-resistant cap and could cause endoscope malfunction.
- Be sure to remove the monitor from the operation part before immersion, because the monitor can not be immersed, otherwise it will cause water penetration in the monitor and malfunction.

Note

- When the pressure compensation cap is in place, the covering of the bending section will expand as the air pressure inside the endoscope increases. This is normal.

Leakage Testing With Gasbag

- Always keep the leakage tester dry. Any water remaining on the leakage tester may cause damage to the endoscope.
- Press the pressure release valve before taking the pressure compensation cap off the leakage detector. Discharge air in the endoscope until the indicator on the screen points to 0 MMHG. Otherwise, it may lead to endoscope damage.

- Only the duct of leakage tester and the pressure compensation cap can be immersed.(as shown in Fig 7.11).Immersion of other components of the leakage tester may lead to its damage.
- The pressure shall not exceed 240 mmHg. If the pressure increases too much, it might damage the endoscope. Press the pressure release valve to discharge the air (as shown in Fig 7.11).
- The leakage detector is only dedicated for endoscope leak detection and must not be used for other purposes, otherwise it may cause damage to tested object.

The endoscope with Waterproof cover is completely waterproof and can be thoroughly soaked and cleaned. However, the tightness of the endoscope should be checked before immersion.

1. **Connect the leakage detector:** Point the cracks in the leak cover to the slot of the leak joint on the video endoscope. Push into the leak cover and rotate clockwise to the head.

Fig 7.10

2. **Pressurization:** use the random leakage detector to pump air into the inside of the endoscope body through the leak detection cover, and maintain it to 220-240 mm Hg. (Note: Pump slowly when pumping)
Stabilize the indicator needle for a few seconds after compression. Read the pressure when the indicator needle is stable.


Fig 7.11

2. Observation: immerse the endoscope in clear water, and observe for about 30 seconds while making an angle, and confirm that there is no continuous bubble on the whole endoscope. If leakage is detected, please contact the manufacturer.


Fig 7.12

Warning

- If the indicator needle drops to 0 mmHg continuously, it means that the endoscope may have serious water leakage or the leakage detector may be damaged. The leak test should be stopped immediately. If the endoscope is still immersed in water, water will enter the endoscope without pressure inside. As a result, more serious problems will occur.
- During the leakage test, if bubbles continuously emerge from the leakage test cover, the leakage test cover or joint may be damaged. Please replace the leak detection cover in time or contact the manufacturer for maintenance.

4. Remove the leakage detector: take out the endoscope from the clear water, unscrew the vent valve on the left side of the middle of the leakage detector, release the air pressure inside the endoscope, and remove the leakage detector.


Fig 7.13

7.6 Manual Cleaning

Warning

- Clean the endoscope and accessories thoroughly to remove microorganisms or organic material. If the endoscope is not immediately cleaned after each procedure, residual organic debris will begin to solidify and it may be difficult to effectively clean, disinfect the endoscope.

Caution

- To avoid water leakage, do not use excessive force when cleaning the endoscope.
- Remove balance valve from endoscope. If a counterbalance valve is installed, water may enter the endoscope and damage the endoscope.
- For details of detergent, please refer to section 6.2 "Detergent".

After completing the leakage test, perform manual cleaning according to the procedures described below.

Equipment needed

Prepare the following equipment

- Personal protective equipment.
- Basin which is at least 40 cm × 40 cm (16" × 16") in size and deep enough to allow the endoscope to be completely immersed.
- clean, lint-free cloths
- EmPower multienzyme cleaning agent
- 30cm³ (30 ml) syringe
- cotton swab
- channel cleaning brush
- Channel connector cleaning brush
- soft brush

Caution

- In order to avoid damage to the endoscope, do not soak the endoscope with other instruments except the above instruments.

Preparation

1. Wear appropriate personal protective equipment.
2. Fill a basin with EmPower solution at the room temperature and the concentration is 1:128. Use a basin which is at least 40 cm × 40 cm (16" × 16") in size and deep enough to allow the endoscope to be completely immersed.

Cleaning the external surface

1. Soak the endoscope in the washing liquid.

2. Use a clean lint free cloth to wipe the external surface of the whole endoscope, especially the insertion part and the control part. Use and replace the lint free cloth once to ensure that the external surface is completely clean. Do not infiltrate liquid into the monitor. Liquid infiltrating into the monitor may damage the internal circuits.

3. Visually check that there is no residual debris on the surface of the endoscope.

Brushing the channels

- Be sure to brush the inside of the instrument and suction channel. Otherwise, insufficient cleaning and disinfection of the endoscope may pose an infection control risk to the patient and/or operators performing the next procedure with the endoscope.

- To avoid splattering detergent solution when the channel cleaning brush is pulled out, keep the endoscope submerged.
- The channel cleaning brush is a consumable item. Disposable brush belongs to the disposable products that can't be repeated use. Repeated using may cause the brush head being bent or kinked, as a result the brush head may come off. Confirm that the brush is free from any damages or other irregularities before and after each use. If the brush head comes off after brushing, immediately retrieve it and carefully check that it is not left inside the channel of the endoscope by passing a new cleaning brush or other endo-therapy accessories. If the brush head is left in the channel, it can drop into the patient body during the procedure. Due to the difference in location, the new cleaning brush passing through the pipe may not be able to remove the missing parts. In this case, please contact manufacturer as soon as possible.
- Do not try to pass the cleaning brush through the tip of the insertion part. Otherwise, the cleaning brush will get stuck and cannot be pulled out.

Follow the steps below to scrub all pipes and Clamp Ports while the endoscope is still immersed in the washing liquid (as shown in Figure 7.14).


Fig 7.14

Brushing injection forceps tubing (location A)

1. Hold the cleaning brush at a distance of 3 cm from the bristles.
2. Insert the channel cleaning brush into the suction valve Housing as illustrated by A in Fig 7.14. Using short strokes, feed the brush through the channel until it emerges from the distal end of the endoscope.
3. Clean the bristles with your fingertips in the detergent solution. Carefully pull the brush back through the channel.
4. Clean the bristles in detergent solution again.
5. Visual inspection, if there was residual debris on the sample, repeat Step 1 to 4 until all debris is removed.

Brushing injection / forceps tube opening (location B)

1. As shown in Figure 7.14, B, insert the tube channel port cleaning brush into the injection tube opening until the brush handle touches the tube opening.
2. Turn the inserted brush once.
3. Pull the brush out and clean the bristles with your fingertips in the detergent solution.
4. Visual inspection, if there was residual debris on the sample, repeat Steps 1 to 3 until all debris is removed.

Soak the endoscope in the detergent solution

1. Completely immerse the endoscope in the washing solution.
2. Attach the 30cm³ (30ml) syringe on the Clamp Port.
3. Pull up the push rod of the syringe and fill the working channel with detergent. (See Fig 7.15)
4. Completely immersed the endoscope in the detergent solution and detach the syringe.


Fig 7.15

5. Using a lint-free cloth, wipe all debris from the endoscope's external surfaces while the endoscope is immersed in the detergent solution.
6. Soak the endoscope for 5 min in EmPower solution at the room temperature and the concentration is 1:128
7. Remove the endoscope from the detergent solution.

7.7 Rinsing and Drying after Manual Cleaning

- Rinse the external surface and channels of the endoscope, and wash thoroughly the equipment with rinse water to remove the residual detergent solution. Otherwise, the residual detergent solution may cause mucosal inflammation.

- Do not use tap water to wash the equipment

Equipment needed

Please prepare the following equipment:

- personal protective equipment
- basin, 40 cm × 40 cm (16" × 16") minimum size and deep enough to immerse the endoscope completely
- deionized water (Aq. Dest.)
- clean, lint-free cloth or paper
- Pressure water gun
- Pressure air gun

Rinsing the outer surface

1. Wear appropriate personal protective equipment.
2. Fill a basin with the deionized water (Aq. Dest.). Use a basin which is at least 40 cm × 40 cm (16" × 16") in size and deep enough to allow the endoscope to be completely immersed.
3. Soak the endoscope and suction cleaning adapter in the deionized water (Aq. Dest.), and gently agitate to thoroughly rinse them.

Rinsing all channels

1. Aim the pressure water gun at the Clamp Port to flush the working channel until there is no residual cleaning solution.
2. Flush the external surface of the endoscope and all pipe ports with flowing water.
3. Use a pressure air gun to inflate the pipe for at least 30s to remove the water in the pipe.

Dry external surfaces

1. Place the endoscope, buttons and valves on a special drying table with sterile towels. Sterile towels should be replaced every 4 hours.
 2. Fill all pipes with 75%–95% ethanol or isopropyl alcohol.
 3. Use a pressure air gun to inflate all pipes with clean compressed air for at least 30 s until they are completely dry.
 4. Dry the outer surface of the endoscope, buttons and valves with sterile wiping cloth and pressure air gun.
- Install buttons and valves.

7.8 Disinfection

- During the entire process of disinfection, the endoscope and all equipment should be completely immersed in the disinfectant. If the equipment is not removed from the endoscope during immersion, or all parts are not completely immersed, the disinfectant will not be able to fully contact all surfaces of the equipment. This may reduce the effectiveness of disinfection.

- Refer to "Disinfectant solution" in Section 6.3 for details of disinfectant solution.

After manual cleaning, disinfect the endoscope according to the procedures described below.

Equipment needed

Please prepare the following equipment

- personal protective equipment (gloves, glasses and mask)
- Peracetic Acid (Dilute the peracetic acid with pure water at the ratio of 1:22.1 to make the final concentration between 2380mg/L~2814 mg/L.)
- Automatic soft endoscope washer-disinfector
- Suction cleaning adapter
- Pressure air gun

Note

The information of disinfecter:

Name: Automatic soft endoscope washer-disinfector

Model: Rider 60A

Manufacturer: SHINVA MEDICAL INSTRUMENT CO., LTD.

Process

1. Wear appropriate personal protective equipment
2. Put device into the Automatic soft endoscope washer-disinfector until it is fully loaded.
3. Soak the Flexible Digital Bronchoscope in peracetic acid for 10 min.
4. Disinfection: Time: 5 min.
5. Rise 1: Time: 30 s.
6. Rise 2: Time: 30 s.

Warning

- All disinfection steps should be performed with the endoscope and all equipment completely immersed. If the equipment is disconnected respectively while not immersed, disinfectant solution may not adequately contact all surfaces of the equipment. As a result, the effectiveness of disinfection may be reduced.

Removing the endoscope and all equipment from automatic soft endoscope washer-disinfector

1. Remove the endoscope from the disinfectant solution.
2. Use a pressure air gun to inject air into the pipe through the Clamp Port for 30s to remove the disinfectant in the pipe.

7.9 Rinsing and Drying after Disinfection

Caution

- After cleaning and disinfection, thoroughly clean and dry the tube of the endoscope. Otherwise, bacteria will grow in the pipeline, posing a risk of infection to the patient or operator who uses the endoscope next time.
- Rinse the external surface and channels of the endoscope, and wash thoroughly the equipment with rinse water to remove the residual detergent solution. Otherwise, the residual detergent solution may cause mucosal inflammation.

Caution

- Alcohol is flammable. Handle it with care.
- About the quality of the rinsing water, please refer to Section 6.5 "Rinsing Water".

After disinfection, rinse and dry according to the following steps endoscope and all equipment.

Equipment needed

Please prepare the following equipment

Caution

- Never reuse the rinsing water.

- Personal protective equipment
- Sterile basins, 40 cm × 40 cm (16" × 16") minimum size and deep enough to immerse the endoscope completely
- Sterile deionized water (for sterile water rinse)
- Sterile, lint-free cloths
- Pressure water gun
- Pressure air gun

If sterile water is not available, prepare the following equipment.

- Deionized water (aq. Dest.) (For non-sterile water rinse)
- Small basin
- 70% ethyl or isopropyl alcohol
- Sterile cotton swabs
- small basin
- Pressure water gun
- Pressure air gun

Sterile water rinse

1. Wearing appropriate personal protective equipment.
2. Put sterile water in the basin. The water basin should be at least 40 cm×40 cm (16"×16") deep enough to soak the entire endoscope.
3. Immerse the endoscope in the sterile water. Using a sterile, lint-free cloth, thoroughly rinse and wipe all external surfaces.
4. Aim the pressure water gun at the Clamp Port and flush the pipeline with sterile water for at least 2 minutes to flush the working channel until there is no residual disinfectant.
5. Place the endoscope, buttons and valves on a special drying table with sterile towels. Sterile towels should be replaced every 4 hours.
6. Fill all pipes with 75%~95% ethanol or isopropyl alcohol.
7. Use a pressure air gun to inflate all pipes with clean compressed air for at least 30 s until they are completely dry.
8. Dry the outer surface of the endoscope, buttons and valves with sterile wiping cloth and pressure air gun.
9. Install buttons and valves.

Note

- Flushing the channels with 70% ethyl or isopropyl alcohol after rinsing them with sterile water facilitates drying inside the channels.

Non-sterile water rinse and alcohol flush

1. Wearing appropriate personal protective equipment.
2. Fill a basin with deionized water (Aq. Dest.). Use a basin which is at least 40 cm × 40 cm (16" × 16") in size and deep enough to allow the endoscope to be completely immersed.
3. Immerse the endoscope in deionized water. Use a sterile lint-free cloth to thoroughly rinse and wipe all

external surfaces.

4. Aim the pressure water gun at the Clamp Port and flush the pipeline with deionized water for at least 2 minutes to flush the working channel until there is no residual disinfectant.
5. Take out the endoscope and put the endoscope and corresponding accessories on the special drying table for laying sterile towel.
6. Put 70% ethanol or isopropyl alcohol in a small bowl.
7. Pull up the syringe push rod and inject 70% ethanol or isopropyl alcohol into the working channel from the Clamp Port. Make sure that the ethanol or isopropyl alcohol injected from the Clamp Port flows out of the inner lens end.
8. Use a pressure air gun to inflate all pipes with clean compression for at least 30s until they are completely dry.
9. Thoroughly wipe and dry the external surfaces of the endoscope and the Clamp Port with a sterile lint free cloth dipped in alcohol.
10. Using a sterile, lint-free cotton swabs, wipe and dry the inside of the Clamp Port.
11. Store the components following the instructions given in Chapter 8, "Storage and Disposal".

Chapter 8 Storage and Disposal

After cleaning, and disinfection the clean endoscope and accessories should be kept away from any contaminated equipment. If the clean endoscope and accessories are contaminated between two inspections, it will pose a risk of infection to the patient or operator who will use it next time.

- To prevent contamination of clean endoscope and accessories, be sure to keep storage cabinet also clean.
- The storage cabinet must be clean, dry, well ventilated and maintained at ambient temperature. Do not store the endoscope in direct sunlight, at high temperature, in high humidity or exposed to X-rays and ultraviolet-rays. Otherwise, it may damage the endoscope or present an infection risk.
- Do not store the endoscope in the carrying case. Use the carrying case only for shipping the endoscope. Routinely storing the endoscope in a humid, non-ventilated environment such as the carrying case may present an infection risk.

8.1 Storage of the Endoscope

1. Detach all equipment (the suction valve cover, Biopsy Valve, and waterproof cover and balanced valve).
2. Confirm that the outer surfaces of the endoscope (especially all internal channel, distal end, lens and electrical contacts) are dry.
3. Carefully wipe the objective lens of the distal end with cotton swab moistened with 70% alcohol.
4. Hang the endoscope in the storage cabinet with the distal end hanging freely. Make sure that the insertion tube hangs vertically and as straight as possible.
5. Store the endoscope in a clean and dry place at room temperature .

Note

- According to the differences of the packaging type and storage conditions (refer to DIN 58 953, Section 9, or local regulations) the storage time of equipment is also different.

8.2 Storage of Reusable Parts, Reprocessing Equipment and Leakage Tester

1. Confirm that all the reusable parts and reprocessing equipment are thoroughly dry.
2. Store all reusable parts loosely in the storage cabinet. Make sure these parts during storage are not in contact with each other.
3. Store all cleaning equipment in a container, and then store the container in the storage cabinet.

8.3 Transporting Outside the Hospital

Transport the endoscope in the carrying case.

Warning

- Always clean, disinfect the endoscope after removing it from the carrying case. If the endoscope is not cleaned, disinfected, the patient might be infected.

Note

- The carrying case cannot be cleaned, disinfected. Clean and disinfect the endoscope before placing it in the carrying case.
- Attach the waterproof cover and balanced valve when transporting the endoscope, to avoid damage to the endoscope caused by changes in air pressure.

8.4 Disposal

When disposing of the endoscope, or any of its components (such as cleaning brush, valves), follow all applicable national and local laws and guidelines.

Chapter 9 Troubleshooting

If it is found that the endoscope is damaged and cannot function normally, or if there are other abnormalities found during the inspection according to the method described in Chapter 3 "Preparation and Inspection", do not use the endoscope. Please contact our company.

Some non-functional faults can be solved by the content of Section 9.1 "Troubleshooting Guide". If the fault still cannot be eliminated after taking the above-mentioned treatment methods, please stop using the endoscope and send it back to our company for repair.

Our company is not responsible for repairing accessories. If the accessory is damaged, please contact our company to purchase a new one.

- If an abnormality of the endoscope is observed, do not use it on the patient. The use of damaged or abnormal endoscopes will endanger the lives of patients and users, as well as cause more serious equipment damage.
- If parts of the endoscope fall into the patient's body due to equipment damage or malfunction, immediately stop using the endoscope and take appropriate measures to retract the parts.

If any abnormality in the function of the endoscope and/or endoscopic image is suspected during use, stop the examination immediately and carefully withdraw the endoscope from the patient as described in Section 9.2, "Withdrawal of the Endoscope with any Abnormality".

9.1 Troubleshooting Guide

The following table shows faults and countermeasures caused by incorrect settings or damaged consumables.

Failures caused by reasons other than the following should be repaired. Repairs performed by maintenance personnel not authorized by the company may cause injury to the patient or user and/or damage to the equipment. Therefore, please contact the company for maintenance in accordance with the provisions in section 9.3 "Repair of Endoscopes"

Irregularity description	Possible cause	Solution
Fluid leakage at the working channel port	Blocking channel	Cleaning the channels with Channel Cleaning Brush.
	The Channel Port is damaged.	Replace it with a new one.
Leakage of fluid between Clamp Port and syringe	Blocking channel	Cleaning the channels with Channel Cleaning Brush.
	The syringe was not inserted firmly	Insert syringe firmly

An image is not clear.	The endoscope lens is dirty.	Wipe with cotton swab moistened with 70% ethyl or isopropyl alcohol; If the image can't be observed, return the endoscope to our company for repair.
------------------------	------------------------------	---

9.2 Withdrawal of the Endoscope with any Abnormality

If an abnormality occurs while the endoscope is in use, take a proper measure as described in either "When the endoscopic image appears on the monitor" or "When the endoscopic image does not appear on the monitor or the frozen image cannot be restored" below. After withdrawal, return the endoscope for repair as described in Section 10.3, "Returning the endoscope for repair".

- If the endoscope cannot be withdrawn from the patient smoothly, do not attempt to forcibly withdraw it. If any irregularities are suspected, immediately contact us. Forcibly withdrawing the endoscope may cause patient injury.

When the endoscopic image appears on the monitor

1. Aspirate blood, mucus or other debris by depressing the suction valve.
2. Carefully withdraw the endoscope while observing the endoscopic image.

When the endoscopic image does not appear on the monitor or the frozen image cannot be restored

1. Turn the monitor OFF and then ON again. If the endoscopic image appears or the frozen image restored, follow the procedure of step 2 and steps below in "When the endoscopic image appears on the monitor". When the endoscopic image still does not appear or the frozen image cannot be restored, perform the following steps.
2. Turn OFF the monitor and suction pump.
3. Turn the UP/DOWN angulation control knob to neutral positions.
4. Release the angulation control knobs and carefully withdraw the endoscope. Remove the Guide blade from the patient's mouth.

9.3 Returning the Endoscope for Repair

- Before returning the endoscope for repair, thoroughly clean and disinfect it. Improper cleaning, and disinfection can bring infection risks to the hospital or the company's personnel who dispose of the equipment.

- Our company is not responsible for the equipment damage and/or person

injury which are not caused by our company personnel trying to repair the equipment.

Please contact our company before returning the endoscope for repair. When returning for repair, please attach a description of the malfunction or damage of the endoscope, as well as the name and contact number of the person most familiar with the malfunction of the endoscope in your organization. Please also attach the warranty card.

When returning the endoscope for repair, please carry it according to the instructions in section 8.3

"Transporting the Endoscope Outside the Hospital".

Chapter 10 Electromagnetic compatibility

Electromagnetic compatibility

Caution:

- BF series Bronchoscope meets the requirements of IEC60601-1-2 standard for electromagnetic compatibility
- The user shall install and use the electromagnetic compatibility information provided by the random files
- Portable and mobile RF communications devices may affect the performance of the BF series Bronchoscope and avoid strong electromagnetic interference such as near mobile phones and microwave ovens;
- The guidelines and manufacturer's declarations are detailed in the annex

Warning:

- A device or system should not be used with other devices placed close to or over stacked, if must be close to or stacked, should observe and verify the normal operation in its use of the configuration.
- Type A equipment for the industrial environment, since the BF series of Bronchoscope conduction harassment and radiation harassment, in other environments to ensure electromagnetic compatibility may have potential difficulties.
- The provisions of the accessories, transducer or cables and equipment and systems used together, which may cause the equipment or system emission increase or decrease of immunity.

Annex:

Guide and manufacturer's statement - Electromagnetic Emission			
BF series of Video Endoscope used in the electromagnetic environment of the expected performance of the following provisions, the performance of BF series of Video Endoscope purchasers or users should ensure that it is used in the electromagnetic environment:			
Test item	EMC standard or test method	Immunity IEC 60601 test level	Compliance level
Conducted Emissions	CISPR 11	Group 1 Class A	N/A
Radiated RF Emissions	CISPR 11	Group 1 Class A	Group 1 Class A
Harmonic distortion	IEC 61000-3-2	Class A	N/A
Voltage fluctuations and flicker	IEC 61000-3-3	d_{max}<4%, P_s≤ 1, dc<3.3%	N/A

Guide and manufacturer's statement - Electromagnetic Immunity

BF series of Bronchoscope used in the electromagnetic environment of the expected performance of the following provisions, the performance of BF series of Bronchoscope purchasers or users should ensure that it is used in the electromagnetic environment:

Test item	EMC standard or test method	Immunity IEC 60601 test level	Compliance level
Electrostatic discharge (ESD)	IEC 61000-4-2	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields	IEC 61000-4-3	3V/m rms 80MHz-2.7GHz 80% AM at 1 kHz	3V/m rms 80MHz-2.7GHz 80% AM at 1 kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	9V/m, 217 Hz LTE Band 13,17 for 710-780MHz;WLAN 802.11 a/n for 5240-5785 MHz; 27V/m,18Hz for 385MHz; 28V/m,18Hz for 810-930MHz, ± 5 kHz deviation 1 kHz sine for 450MHz, 217Hz for 1720-2450.	9V/m, 217 Hz LTE Band 13,17 for 710-780MHz;WLAN 802.11 a/n for 5240-5785 MHz; 27V/m,18Hz for 385MHz; 28V/m,18Hz for 810-930MHz, ± 5 kHz deviation 1 kHz sine for 450MHz, 217Hz for 1720-2450.
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50Hz or 60Hz	30 A/m 50Hz or 60Hz
Proximity magnetic fields	IEC 61000-4-39	8A/m, CW for 30kHz. 65A/m, 2,1kHz Pulse modulation for 134,2kHz. 7,5A/m, 50kHz Pulse modulation for 13,56MHz	8A/m, CW for 30kHz. 65A/m, 2,1kHz Pulse modulation for 134,2kHz. 7,5A/m, 50kHz Pulse modulation for 13,56MHz
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V rms 150 kHz- 80 MHz 6 V rms in ISM and amateur radio bands between 0,15 MHz and 80 MHz	3 V rms 150 kHz- 80 MHz 6 V rms in ISM and amateur radio bands between 0,15 MHz and 80 MHz

Note: UT refers to an alternating voltage test voltage before

Guide and manufacturer's statement - Electromagnetic Immunity

BF series of Bronchoscope used in the electromagnetic environment of the expected performance of the following provisions, the performance of BF series of Bronchoscope purchasers or users should ensure that it is used in the electromagnetic environment:

Test item	EMC standard or test method	Immunity IEC 60601 test level	Compliance level
RF transmission	IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz
RF radiation	IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz

Note 1: in 80MHz and 800MHz frequency, adopt the high frequency formula.

Note 2: these guidelines may not be appropriate in all cases, the electromagnetic propagation may be affected by buildings, absorption and reflection of objects and people.

a Fixed transmitter electric field strength, such as base station of wireless (cellular / cordless) telephone and ground mobile radio, amateur radio, AM (AM) and FM (FM) radio and television broadcasting etc. their emission electric field strength is impossible to be predicted in theory. In order to evaluate the electromagnetic environment of a fixed RF transmitter, the investigation of electromagnetic field should be conducted. If the field strength of the BF series Bronchoscope is measured higher than the RF compliance voltage of the above-mentioned application, the performance of the BF series Bronchoscope should be observed to verify its normal operation. If the observed performance is abnormal, then additional measures may be required, such as the new series of BF Bronchoscope performance orientation or location.

b In the whole frequency range of 150KHz ~ 80MHz, the emission electric field strength should be less than 3 V/m.

Recommended isolation distance between portable and mobile communication equipment RF and BF series of Bronchoscope performance

The performance of the BF series Bronchoscope is expected to be used in the controlled electromagnetic environment of radiated RF harassment. On the basis of the maximum output power of communication equipment, the purchasers or users can maintain portable and mobile RF communication equipment through the following recommended (transmitter) and BF series Bronchoscope performance minimum distance to prevent electromagnetic interference.

Rated maximum output power of the transmitter/W	Corresponding to different frequency transmitter isolation distance/m		
	150 kHz ~ 80 MHz $d = 1.2\sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2\sqrt{P}$	800 MHz ~ 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

The transmitter rated maximum output power are not listed in the form, the recommended isolation distance D, in meters (m) as a unit, can use the corresponding transmitter frequency bar formula to determine, where P is the maximum rated output power provided by the transmitter of transmitter manufacturers, by Watt (W) as a unit.

Note 1: in 80MHz and 800MHz frequency, adopt the high frequency formula.

Note 2: these guidelines may not be appropriate in all cases, the electromagnetic propagation may be affected by buildings, absorption and reflection of objects and people.

Manufacturers: Zhuhai Vision Medical Technology Co., Ltd.

**Address: 301 Building 4, No. 30 Feipeng Road, GMTCM Park, Hengqin New Area,
Zhuhai City, Guangdong Province 519031 P.R. CHINA**

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