
Dentsply Tulsa Dental Specialties
608 Rolling Hills Drive, Johnson City,
Tennessee, 37604, USA

Approved	Terri Stanton – Senior QA Manager	Signature: <i>Terri Stanton</i>	Date: <i>Aug 30, 2024</i>
Approved	Anthony Skeans – Notary	Signature: <i>Anthony Skeans</i>	Date: <i>Aug 30, 2024</i>

My Commission Expires
08/02/2026

Instructions for Use of the medical device

Version 2.0

ProTaper Ultimate Conform Fit Gutta-Percha Points in variants

2024

1. Name of the medical device

ProTaper Ultimate Conform Fit Gutta-Percha Points in variants (hereinafter referred to as "Gutta-Percha Points", "Points", "medical devices", "devices", "product").

2. Information about the developer, manufacturer, manufacturing site:

- **Developer, manufacturer**

DENTSPLY Tulsa Dental Specialties, 608 Rolling Hills Drive, Johnson City, Tennessee, 37604, USA

- **Manufacturing site**

DENTSPLY Tulsa Dental Specialties, 608 Rolling Hills Drive, Johnson City, Tennessee, 37604, USA

- **Authorized representative**

Name: "Dentsply Sirona" LLC

Address: 115035, Moscow, Ovchinnikovskaya nab., 18/1, bldg 2, room 3n

Telephone: + 7(495)725-10-87

3. Classification of the medical device

Potential risk class: 2a

Medical device type in accordance with the Nomenclature Classification of Medical Devices: 174010

All-Russian Classifier of Products by Types of Economic Activity Code: 21.20.24.180

Sterility: non-sterile

Number of uses: single use

Disinfect before use

4. General information

1. Intended use of the medical device

Dental gutta-percha points are designed to fill the formed, cleaned and rinsed lumen of the root canal in the course of endodontic treatment.

2. Indications for use of the medical device

Dental gutta-percha points are indicated for the treatment of endodontic diseases (also used as a support for the restored dental crown).

Indications for root canal treatment:

1. Teeth with improperly filled root canals with radiologic evidence of developing or persistent apical periodontitis and/or its symptoms.

2. Teeth with improperly filled root canals when replacement of the crown restoration or whitening of the dental crown tissue is required.

Specific indications for root canal treatment:

1. Permanently damaged or necrotic pulp with or without clinical and/or radiologic evidence of apical periodontitis.
2. Selective devitalization, for example, to provide space for a post before making an overdenture, doubtful pulp health prior to restoration procedures, likelihood of pulp exposure during restoration of a displaced tooth, and prior to root resection or hemisection.

3. Contraindications for root canal treatment:

1. Teeth that cannot be restored (impaired integrity and thinning of the root walls; difficult or impossible root canal passage, its curvature, granulomas and cysts formation).
2. Insufficient periodontal support (severe periodontitis).
3. Hypersensitivity to any materials of the product.

4. Adverse reactions

Irritation or allergic reaction due to hypersensitivity to product materials.

5. Intended user

Dental gutta-percha points should only be used in a clinical or hospital setting, in accordance with good dental practice, by qualified dentists, such as general practitioners, and by endodontists and dental assistants.

6. Patient population

Patients who require root canal treatment.

7. Area of application

Dentistry

8. Warnings and precautions

Warning when using the product:

Gutta-percha points are single-use products. Disposable products are not intended to be reprocessed or reused. Reuse of these products increases the risk of cross-contamination and/or breakage.

Use of the products may result in the following residual risks:

- User dissatisfaction and delayed treatment / patient discomfort.
- Retreatment (original infection not cured)
- Serious infection / serious consequences for patient's health (infection, toxic or allergic reaction)
- Tooth fracture due to obturation instruments

Precautions:

- Safety and efficacy of use in children, pregnant and
- breast-feeding women have not been established.
- Inspect the packaging and product before use: if the packaging or product is damaged, do not use the product and discard it.

- Dispose of the product in accordance with local regulations for the safe disposal of sharp and contaminated products.
- Do not sterilize at high temperatures.
- Do not use disinfecting solutions containing phenol or any products incompatible with gutta-percha points.
- Use rubberdam during endodontic treatment.

5. Design variants of the medical device

ProTaper Ultimate Conform Fit Gutta-Percha Points in variants:

1. ProTaper Ultimate Conform Fit, size F1, 60 pcs. per package.
2. ProTaper Ultimate Conform Fit, Size F2, 60 pcs. per package.
3. ProTaper Ultimate Conform Fit, Size F3, 60 pcs. per package.
4. ProTaper Ultimate Conform Fit, Size FX, 60 pcs. per package.
5. ProTaper Ultimate Conform Fit, Size FXL, 60 pcs. per package.
6. ProTaper Ultimate Conform Fit, Sizes: F1 - 12 pcs, F2 - 24 pcs, F3 - 24 pcs, 60 pcs. total per package.

6. Description of the medical device

The product is a one-piece injection molded product in the form of a gutta-percha cone, a color-coding ring on the coronal part (as a visual aid and in accordance with the color of the corresponding instrument size) and a cone head with a molded size. During the root canal obturation procedure, the colored ring and cone head will be removed.



Figure 1 ProTaper Ultimate Conform Fit Gutta-Percha Points

7. Main technical characteristics

Graphic presentation:

Overall point length 29 mm ± 0.5 mm

Note:

Dimensions in mm.
unless otherwise specified

X. ± 0.5 mm
.X ± 0.1 mm
.XX ± 0.01 mm

Size	Text 1	Text 2	Ring color	ØD0 ± 0.03 mm	ØD3 ± 0.03 mm	ØD9 ± 0.03 mm	ØD16 ± 0.03 mm
F1	F1	A-Z a-z	Yellow	0.170	0.390	0.692	0.937
F2	F2	A-Z a-z	Red	0.230	0.455	0.729	0.939
F3	F3	A-Z a-z	Blue	0.280	0.535	0.787	0.940
FX	FX	A-Z a-z	Green	0.330	0.660	1.018	1.134
FXL	FXL	A-Z a-z	Yellow	0.480	0.750	1.023	1.134

Size	Taper, $\pm 10\%$	Weight, g, $\pm 10\%$
F1	5.0%	0.037
F2	4.6%	0.039
F3	4.2%	0.041
FX	6.0%	0.058
FXL	4.6%	0.061

Meting point: $< 100^\circ\text{C}$

8. Description of use of the medical device

Gutta-percha points should only be used in a clinic or hospital setting by qualified dentists. After forming and cleaning the root canal, the specialist must perform obturation technique. This procedure involves the use of gutta-percha points and root canal filling paste. The gutta-percha points selected must match the color code and shaping file size of the last tool used in the root canal preparation. Points must be disinfected before use. The heated tool can be used to soften the material and provide a perfect seal and adaptation of the root canal dimensions (except for ISO standardized and auxiliary gutta-percha posts designed for lateral/cold compaction). When using a heated tool, the instructions for use of the specified tool must be followed.

It is necessary to use the appropriate tool to cut the gutta-percha points, including the colored coding ring, if any, after inserting the gutta-percha point into the root canal. The colored ring

must not remain in the root canal after insertion.

Branded gutta-percha points

The described gutta-percha points are a separate device as well as a part of the complete ProTaper Ultimate system and can be used in combination with shaping files and paper points of this system.

The medical device “ProTaper Ultimate Conform Fit Gutta-Percha Points in variants” is radiopaque.

9. Materials used in the medical device

10. Disinfection

Name	Material	Supplier / Manufacturer
Design variant of the medical device		
- ProTaper Ultimate Conform Fit, Size F2	Polyisoprene	Thermofisher Scientific
	Zinc oxide	American International Chemical
	Barium sulfate	
	Titanium dioxide	
	Red ink 7224 RED 456-C	Markem-Imaje
ProTaper Ultimate Conform Fit, Size F1	Polyisoprene	Thermofisher Scientific
	Zinc oxide	American International Chemical
	Barium sulfate	
	Titanium dioxide	
	Yellow ink 7224 YELLOW 429-E	Markem-Imaje
- ProTaper Ultimate Conform Fit, Size FXL	Polyisoprene	Thermofisher Scientific
	Zinc oxide	American International Chemical
	Barium sulfate	
	Titanium dioxide	
	Blue ink 7224 BRITE BLUE	Markem-Imaje
ProTaper Ultimate Conform Fit, Size F3	Polyisoprene	Thermofisher Scientific
	Zinc oxide	American International Chemical
	Barium sulfate	
	Titanium dioxide	
	Blue ink 7224 BRITE BLUE	Markem-Imaje
- ProTaper Ultimate Conform Fit, Size FX	Polyisoprene	Thermofisher Scientific
	Zinc oxide	American International Chemical
	Barium sulfate	
	Titanium dioxide	
	Green ink 7224 GREEN 575-F	Markem-Imaje

- Immerse the product in 5.25% NaOCl for 1-5 minutes at ambient temperature. Remove all bubbles from the surface of the obturator.
- Gently soak the gutta-percha point with sterile gauze moistened with 70% isopropyl alcohol.
- Dry the products in air.
- Do not use disinfecting solutions containing phenol or any products incompatible with gutta-percha points.

11. Shelf life

The shelf life is 4 years. Do not use after the expiration date indicated on the package.

12. Storage, transportation and operating conditions

Storage conditions	Air temperature: from +10 °C to +30°C Relative humidity: 20%-85%. Atmospheric pressure: 48 kPa - 108 kPa Keep away from sunlight and heat. Store in a dry place.
Transportation conditions	Air temperature: from +5 °C to +25°C Relative humidity: 20%-85% Atmospheric pressure: 48 kPa - 108 kPa
Operating conditions	The devices are properly operated inside the oral cavity, in the presence of body fluids at nominal temperatures: from +32 °C to 42 +°C.

	Relative humidity: no specific requirements. Atmospheric pressure: no specific requirements. Products should only be used by qualified medical personnel in clinical settings.
--	--

13. Information about packaging of the medical device

Consumer packaging: Dental gutta-percha points are packaged in plastic box (bottom and lid) with a shrink wrap.

Item	Description
Dimensions (LxWxH), mm ($\pm$ 5%)	80x48x12
Weight, g, not more	23

Materials of packaging for all variants of the medical device:

	Material	Supplier / Manufacturer
Plastic packaging lid	Polystyrene POLYREX® GPPS, PG-33	Chi Mei Corporation
Plastic packaging bottom	Styrene Butadiene Copolymer POLYREX® HIPS, PH-88	Chi Mei Corporation
Holds gutta-percha points in place to avoid mixing	Thermo-Plastic-Rubber material (Styrene Butadiene copolymer) SBS 230	Ningbo Poly U plastic product Co., LTD
Protection of gutta-percha points	Food grade paper	Zhejiang Wuxing Paper Co., Ltd

14. Labeling

Explanation of symbols on the consumer packaging of the medical device.

	Reference Number
	Do not reuse
	Non-sterile
	Do not use if packaging is damaged
	Caution! Consult instructions for use
	Upper limit of temperature range

	Keep dry
	Keep away from sunlight
	Consult instructions for use
	Manufacturer
	Authorized representative in the European Community
	CE mark
	Use by
	Date of manufacture
	Batch code
	Quantity (Indicates the number of pieces in packaging)
	Do not sterilize
	This product does not contain natural rubber latex.
	Gutta-percha
	Medical device
	Assortment*
	Temperature range
	Atmospheric pressure range
	Humidity range

* Only for ProTaper Ultimate Conform Fit, Sizes: F1 - 12 pcs, F2 - 24 pcs, F3 - 24 pcs, 60 pcs. total per package.

Transport marking:

The cardboard box is labeled with information about the parcel (sender information, recipient information, total number of items, delivery address information including product name, lot

number and catalog number), a label with information about the country of origin and package dimensions, as well as a label, indicating the method of handling goods:

Symbol	Meaning
	Caution, fragile! The symbol must be present on fragile cargoes. Cargoes marked with this symbol must be handled with care and must not be tipped over or thrown.
	Keep dry It is necessary to protect the cargo from moisture. Packages with this marking must always be protected from excessive moisture and, accordingly, must be stored under shelter. If particularly large or bulky packages cannot be stored in warehouses or under shelter, they should be carefully covered with tarpaulins.
	This side up The marking indicates the correct vertical position of the load. The packaging must be transported, stored and used in such a way that the arrows always point upwards. Rotation, rocking, overturning or other similar handling of the load should be avoided.
	Temperature range
	Humidity range
	Atmospheric pressure range

15. Information on disposal or destruction

In the Russian Federation, this product should be disposed of as medical waste of Class B in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures".

16. Warranty

The manufacturer does not bear any responsibility and does not pay compensation for possible damage or accidents caused by:

- **Use of tools unrelated to the system which can affect normal operation.**
- **Failure to follow the instructions for use.**

The responsibility for testing the material for suitability and use for any purpose not explicitly stated in the instructions shall be borne by the user.

For product quality issues, please contact the manufacturer's authorized representative.

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