02402730

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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
THIS IS TO CERTIFY THAT: the seal of Zhuhai Vision Medical Technology Co., Ltd on the annexed DOCUMENT is genuine.

China Council for the Promotion of International Trade

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证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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
Zhuhai Vision Medical Technology Co., Ltd.

Product Information

- Product Name: Video Endoscope Processor
- Product Type: VS610 (With Screen)
- Manufacture Name: Zhuhai VisionMedical Technology Co., Ltd.
- Address: 301 Building 4, No. 30 Feipeng Road, GMTCM Park, Hengqin New Area,
519031 Zhuhai City, Guangdong Province, PEOPLE' S REPUBLIC OF CHINA
- Contact information
Tel: +86-756-6910150
Fax : +86-756-6910511
Zip code : 519031

Version Information

This manual version may be upgraded due to a change of software or technical specifications at any time without notice. The version information of this specification are as follows:

- Document number: TD/CE-V6-A-009
- Document version: V1.2
- Released date: Jan. 2024

©2024, The copyright of Zhuhai Vision Medical Technology Co., Ltd

CE Mark

EC Representative Name:

Shanghai International Holding Corp. GmbH (Europe)

EC Representative Address:

Eiffestraße 80, 20537 Hamburg, Germany

Statement

The manufacturer owns the rights to the unpublished use of the instructions and has the right to treat

them as confidential information. The use of the manual is only a reference for the operation and maintenance of the products, and no other person has the right to disclose the use of the instruction to others.

This product specification contains by copyright law protection of proprietary information, all rights reserved, without manufacturer's written consent shall not be any part of this instruction manual for camera copy, copy, or translated into other languages.

The contents of this manual shall be subject to change without notice.

Responsibilities of the Manufacturer

In the following conditions are met, the manufacturer will responsible for the product safety, reliability and performance:

- Product installation, maintenance or upgrades were carried out by the manufacturer authorized or approved personnel;
- The storage environment, product environment, electrical environment shall meet product specifications;
- Use this product in accordance with < use instruction >
- Use this product within the period of use

Overview

♦ Main Structure

The product mainly consists of the processor, power adapter, extension cable, wireless video transmitter and video signal transmission cables.

Caution

The power adapter is part of the endoscopic image processor.

♦ Basic performance of the product

When used with the video endoscope produced by our company, it can normally accept the image signal output by the video endoscope and ensure that the image signal is not lost.

♦ Intended Use

The product is intended used to work with the Video Endoscope manufactured by Zhuhai Vision Medical Technology Co., Ltd. to provide the images for observation, photography and treatment.

♦ Contraindication

Video Endoscope Processor is intended used in conjunction with Video Endoscope, to provide image display and there is no absolute contraindication.

♦ Cautions, Warnings and Notes

- 1) No user self-service components are contained within this product, and internal repair of instruments must be performed by the manufacturer. Modification of equipment may result in equipment damage and risk to operators and patients.
- 2) This product is not part of the therapeutic armamentarium and cannot be applied to families.
- 3) Operate in strict accordance with the product instructions provided by the manufacturer.

About the Manual

The manual gives a detailed account of the usage, functions and operation methods of the product. Please carefully read and understand the contents of the manual before using the product, so that the product can be correctly used and the safety of both patients and operators can be ensured.

The manual gives an introduction of the product according to the most complete configuration. It is possible that some information herein is not applicable to the product purchased by you. For any queries, please contact our company.

Please place the manual near the product in necessary. So that it can be obtained conveniently and in time when needed.

The manual is suitable for professional clinical medical personnel. The readers shall have the knowledge of medical procedures, practices, terms and working experience that is necessary for endoscope.

All the illustrations provided in the manual are for reference only. The settings or data in the illustrations are possibly different from what actually displayed on the product.

About the product

The product is only in contact with the user, and the material of the contact part is ABS, PVC, PC+ABS, PET. The ingredients and raw materials of ME equipment comply with local laws and regulations.

Every two years to verify the basic safety of the product, detailed details, please consult the manufacturer.

Catalogue

SYMBOLS WORDS.....	1
KEY POINTS (PLEASE READ BEFORE USE)	2
USER QUALIFICATION.....	2
COROLLARY EQUIPMENT	2
STANDBY EQUIPMENT	2
INAPPROPRIATE MAINTENANCE AND MODIFICATION PROHIBITED.....	2
SIGNAL WORDS.....	3
WARNINGS AND CAUTIONS.....	3
LCD DISPLAY (FOR VS610 ONLY)	5
CHAPTER1 CHECKING THE PACKAGE CONTENTS	6
CHAPTER 2 INSTRUMENT NOMENCLATURE AND SPECIFICATIONS.....	9
2.1 COMPONENT AND FUNCTION.....	9
2.1.1 <i>Front panel (for VS600, not applied)</i>	9
2.1.2 <i>Front panel (for VS610)</i>	10
2.1.3 <i>Rear panel</i>	11
2.2 MENU FUNCTIONS	13
2.3 SPECIFICATIONS.....	14
2.3.1 <i>Safety Standard</i>	14
2.3.2 <i>Model specifications</i>	15
CHAPTER 3 CONNECT AND INSPECTION	17
3.1 CONNECT	17
3.1.1 <i>Connect the monitor</i>	18
3.1.2 <i>Connection of the power adapter</i>	22
3.2 INSPECTION	22
3.2.1 <i>Inspection precautions</i>	23
3.2.2 <i>Connected endoscopy</i>	23
3.2.3 <i>Inspection power adapter</i>	25
3.2.4 <i>Check image output</i>	25
3.2.5 <i>Lighting and automatic dimming inspection</i>	25
3.2.6 <i>Check the freezing function</i>	26
3.2.7 <i>Check the zoom function</i>	27
3.2.8 <i>Check the EDI, HBE function</i>	28
3.2.9 <i>Check the photo and video recording function</i>	28
CHAPTER 4 FUNCTION AND OPERATION.....	29
4.1 PREVIEW MODE	29
4.2 PLAYBACK MODE	32
4.3SETUP	35
4.4FILE MANAGER	39
4.5OTHER	42
4.6TAKINGPICTURES.....	45
4.7RECORDING VIDEO	46
4.8 WHITE BALANCE ADJUSTMENT	47
4.9DELETE PHOTO OR VIDEO.....	47
4.10 POWER OFF.....	48
CHAPTER 5 STORAGE AND DISPOSAL	49
5.1 STORAGE OF THE VIDEO ENDOSCOPE PROCESSOR	49
5.2TRANSPORTING OUTSIDE THE HOSPITAL	49
5.3DISPOSALOF EXPIRED OR SCRAPPED.....	50

CHAPTER 6 TROUBLESHOOTING.....	51
6.1 TROUBLESHOOTING GUIDE.....	51
6.2 RETURNING THE VIDEO ENDOSCOPE PROCESSOR FOR REPAIR.....	52
CHAPTER 7 CLEANING AND DISINFECTION.....	53
7.1 MAINTAIN.....	53
7.2 CLEANING	53
7.3 DISINFECTION.....	54
CHAPTER 8 ELECTROMAGNETIC COMPATIBILITY.....	55

Symbols Words

The meaning(s) of the symbol(s) shown on the package with the components, the back cover of this instruction manual and/or this instrument are as follows:

Symbols	The meaning of the symbol	Symbols	The meaning of the symbol
	Follow instructions for use		Caution In case of application as a SAFETY SIGN, the rules according to ISO 3864-1 are to be adhered to.
	Type BF Applied Part		CE Mark
	Serial Number		This way up
	Fragile; handle with care		Keep away from rain
	Stacking limit by number		Symbol for Waste Electrical and Electronic Equipment (Please follow Local Ordinances or Regulations of disposal.)
	Manufacturer		Date of manufacture
	Non-ionizing electromagnetic radiation		Authorized Representative in the European Community
IPX3	Protected against spraying water		Medical Device
	Direct current		Alternating current
	Standby		Use indoors only
RoHs	Comply with RoHs directive		

Key points (please read before use)

User qualification

If there are official standards for use established by the medical administration or other organizations, please follow these standards. If there is no official standard, the operator of this product must be a physician who has passed the audit of the medical safety director of the hospital or a doctor who can safely carry out the established operation according to the official regulations designated by the director of the using department according to the difficulty of using the technology.

The physician followed the legal standards set by the Endoscopy Society and other institutions while considering the difficulties of various endoscopic examinations and treatment, and was able to safely perform the endoscopy and treatment procedures. Therefore, this instruction manual does not explain or discuss any clinical endoscopic technique.

Corollary equipment

Refer to Chapter 1 “Checking the Package Contents”, confirm that this instrument is compatible with the ancillary equipment being used. Using incompatible equipment can result in patient or operator injury and/or equipment damage.

This product conforms to the EMC standard IEC 60601-1-2.

Standby equipment

Be sure to prepare standby equipment or equipment with similar functions to avoid accidents in case of equipment failure or abnormal function during operation.

Please check the equipment according to the instructions in the manual before each use. Furthermore, in addition to the checks before each use, the person responsible for the maintenance of the medical device should also regularly check the parts mentioned in this instruction. If the color of the image is abnormal, please adjust it with the white balance button in conjunction with the endoscope and see Section 4.8 for details. If any abnormalities are observed, do not use them and follow the instructions in Section 6.1, Troubleshooting Guidelines. If any abnormalities are observed after examination, please contact with manufacturer.

Inappropriate maintenance and modification prohibited

This product does not contain any user-serviceable parts. Do not disassemble, modify or attempt to repair, otherwise it may result in injury to the patient or damage to the equipment. If necessary, please contact the manufacturer.

Zhuhai Vision Medical Technology Co., Ltd. is not responsible for the injury or damage caused by the maintenance equipment of personnel not from Zhuhai Vision Medical Technology Co., Ltd.

Signal Words

The following signal words are used throughout this manual:

Danger

Indicates an emergency hazardous situation which, if not avoided, could result in death or serious injury.

Warning

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

Caution

Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to warn against unsafe operations or potential equipment damage.

Note

Indicates additional helpful information.

Warnings and Cautions

Follow the warnings and cautions given below when handling this instrument. This information is to be supplemented by the warnings and cautions given in each chapter.

Warning

- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- Modification of this equipment is not allowed
- Avoid scald operator, continuous contact with the product for less than 10 minutes
- A Class equipment is intended to be used in industrial environment. The Video Endoscope Processor may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances.
- Except of the cable supplied sold as spare parts of internal components by manufacturers of the video endoscope processor, the others use of accessories and cable may result in increased product launch interference or low immunity.
- No non-manufacturer adapter. Please confirm whether it is the manufacturer's accessories before accessing the adapter
- Instructions warning: the waterproof grade of the equipment is IPX1, do not say that the equipment is placed in a liquid dangerous environment,

and the liquid intake will cause equipment damage or patient injury

- Replacement of battery by untrained personnel causes hazard (source) (e. g. overtemperature, fire or explosion)
- Devices should be placed reliably to avoid falling.
- Heavy objects cannot be placed on the device, otherwise it will cause device damage or injury to the patient
- Battery replacement is performed by the manufacturer
- Device output is abnormal, please stop using it immediately
- If the liquid enters the video endoscope processor, please turn off the power immediately to prevent fire or electric shock. If liquid has penetrated into the video endoscope processor, do not use it because its internal circuit may be damaged.
- The wireless transmitter contains rechargeable batteries inside, but the batteries do not require professional maintenance. If it is not in use for a long time, the wireless transmitter needs to be charged every three months, otherwise the lithium battery will cause irreversible capacity loss due to the low voltage caused by self-discharge. Please contact the manufacturer with any questions.
- The video endoscope processor contains a signal input and output part. The signal input port (Lemo interface) of the operation Department can only be connected with the endoscope produced by the company through the video connection line attached to the product. The DVI, CVBS, SDI and S-Video video output interfaces must be connected with monitors meeting the requirements of IEC60601-1-2 and IEC60601-1 standards.
- When it is used with other medical electrical equipment, medical equipment meeting the performance requirements of IEC60601-1 and electromagnetic compatibility shall be used to ensure the safety when used with video endoscope processor.
- The user shall install and use according to the EMC information provided in the attached documents.
- Portable and mobile RF communication devices may influence video endoscope processor performance and avoid strong electromagnetic interference when used, such as close to mobile phones, microwave ovens.

Caution

- Electromagnetic interference may occur on this instrument near equipment marked with the following symbol or other portable and mobile RF (Radio Frequency) communications equipment such as cellular phones. If electromagnetic interference occurs, mitigation measures may be necessary, such as reorienting or relocating this instrument, or shielding the location. 
- In order to check the electromagnetic interference from other devices (any devices like X-ray machine, CT, MRT, etc or equipment components of other system), the system should be observed in the operating environment whether it is working properly.

LCD display




- Do not squeeze the LCD display. Otherwise, the image may become blurred or damage the LCD display.
- If the ambient temperature exceeds the normal operating range, the LCD display may take a long time to open, or its color may temporarily change. LCD displays that exhibit poor performance due to high and low temperatures recover at normal temperatures.
- The LCD used for display is made of high-precision technology. Due to its characteristics or different viewing angles on the LCD display, the color and brightness of spots may be inconsistent. This is not a functional exception.

Chapter1 Checking the package contents

Match all items in the package with the components shown below. Inspect each item for damage. If the instrument is damaged, a component is missing or you have any questions, do not use the instrument; immediately contact our company.

All of the following accessories have an expected service life of 5 years

Items	Picture
Video Endoscope Processor	 VS610

Power Plug		Power Adapter	
USB flashdisk		DVI-D to DVI-D Cable	
USB Male to Male Connector		HDMI To DVI Cable (Optional)	
Equipotential Cable		Video connection Cable	
DVI to DVI Cable (Optional)		S-Video HD Cable (Optional)	

CVBS Video Cable (BNC to BNC) (Optional)		Wireless Transmitter (Optional)	
HD-SDI Cable (Optional)		Certificate Of Inspection, Packing List, Product Acceptance List, Warranty List, Instructions for Use	

Form 1.1

Chapter 2 Instrument Nomenclature and Specifications

2.1 Component and function

2.1.1 Front panel (for VS600, not applied)


Fig 2.1VS600 front panel structure schematic diagram

No.	Name	Function
1	Turn on/off Button	Touch the button to turn on the device. Long touch the button for 2.3s to turn off the device. The blue indicator light is power on, and the green indicator light is on ready for use.
2	Reserved Ports	The reserved port is unavailable.
3	Lemo Connection Port	In wired connection mode, extension cable is used to electrically connect the processor with the endoscope. In wire connection mode, the indicator is blue flashing when the wire connection is unsuccessful. When the wire connection is successful, the blue indicator light is always on
4	Wireless connection status indicator	In wired connection mode, the indicator is off. In wireless connection mode, the indicator is blue flashing when the wireless connection is unsuccessful. When the wireless connection is

No.	Name	Function
		successful,the blue indicator light is always on
5	Up Button	Press this key to go up.
6	Left Button	In the interface of preview, touch the button to decrease the brightness in the usual light source mode. In the interface of menu, touch the button to move the selection box to the left.
7	Down button	Press this key to go down.
8	Right button	In the interface of preview, touch the button to increase the brightness in the usual light source mode. In the interface of menu, touch the button to move the selection box to the right.
9	Back button	Returns to the previous menu
10	Main menu button	The main menu is displayed
11	Wireless video transmitter charging port	Used for charging wireless transmitter, providing charging interface.

Form 2.1VS600 Front panel Components Functions

2.1.2 Front panel


Fig 2.2 VS610 front panel structure schematic diagram

No.	Name	Function
1	Turn on/off Button	Touch the button to turn on the device. Long touch the button for 2.3s to turn off the device. The blue indicator light is on when turning off, and the green indicator light is on after turning on.
2	Reserved Ports	The reserved port is unavailable.
3	LemoConnection Port	In wired connection mode, extension cable is used to electrically connect the processor with the endoscope. In wire connection mode, the indicator is blue flashing when the wire connection is unsuccessful. When the wire connection is successful, the blue indicator light is always on
4	Wireless connection status indicator	In wired connection mode, the indicator is off. In wireless connection mode, the indicator is blue flashing when the wireless connection is unsuccessful. When the wireless connection is successful, the blue indicator light is always on
5	LCD display	Display information such as endoscope image and time.
6	Wireless video transmitter charging port	Used for charging wireless transmitter, providing charging interface. When the wireless transmitter is charged, the blue light flashes, and when the battery is fully charged, the blue light goes out

Form 2.2VS610 Front panel Components Functions

2.1.3 Rear panel



Fig 2.3 Video Endoscope Processor rear panel structure schematic diagram

No.	Name	Function
1	SD card	It is used to store photos and videos taken by endoscopy.
2	Reserved Ports	The reserved port is unavailable.
3	Equipotential grounding terminal	Connect the shell or metal part of the equipment with the ground wire
4	USB interface	The pictures and recorded videos stored in the image processor can be copied by connecting with the computer with USB data cable.
5	Foot switch interface	For connecting foot switch
6	CVBS interface	For CVBS data transmission connection
7	SDI interface	For SDI data transmission connection
8	S-Video interface	For S-Video data transmission connection
9	DVI interface	For DVI data transmission connection
10	Power Supply Socket	Power on and off

Form 2.3 Video Endoscope Processor Rear Panel Components Functions

2.1.4 Wireless transmitter

No.	Name	Function
1	Release button	It is used to separate from the operation department, and can easily separate the wireless transmitter and the operation department after pressing the press
2	power on/off button	Touch power on, press for 2.3s power off, with an indicator light. The green light flashes when the WiFi is not connected, and the green light is on when the WiFi is connected,. At low power when the blue light twinkle.
3	White balance button	White balance after a long press for 2.3s
4	Connection port of operation department	Connect with the operation department

2.2 Menu Functions


Fig2.4

No.	Part Name	Function
1	Date and time	Display real-time wireless status
2	SDcard indicator	SD card status is shown, representing an uninserted SD card when the indicator is black and an inserted indicator card when the indicator is white.
3	Wireless transmitter connection status indicator	Displays the current network state. When the indicator is black, the representative camera system is not connected to

No.	Part Name	Function
		a wireless transmitter, and when the indicator is white, the representative camera system is connected to a wireless transmitter.
4	Battery capacity (Wireless Transmitter)	Displays current power of Wireless Transmitter
5	Main menu	Return to main menu.
6	Setup	Product video source, time and date, language, wireless, storage and other Settings entry.
7	Play back	Replay video playing or image.
8	Preview	Observe the real-time image.
9	File Manager	Product internal document viewing, transmission and other operations.
10	Other	Product version number information, etc.

Form 2.4 Video Endoscope Processor Menu Functions

2.3 Specifications

2.3.1 Safety Standard

Note

- According to classification rules, this product is classified as Class I.
This product is classified according to IEC60601-1 is as follows:

Classification of protection against electric shock	Class I equipment external power supply equipment
Classification of protection against electric shock	Type BF
Classification of protection against hazards of explosion	Not suitable for use in places where flammable anesthetic gas is available
Classification of protection against ingress of liquid	Normal equipment (Wireless Transmitter: IPX3)
Classification of running mode	Continuous operation equipment
Classification of disinfection	Equipment of disinfection method recommended by manufacturer

method	
Methods of disinfection	See relevant chapters in the manual

Form 2.5

2.3.2 Model specifications

Model:	VS610
Screen Size:	5 inches
Resolving Power:	1920×1080

Adapter Power Supply

Input: AC100~240V, 1.4~0.3A, 50/60Hz

Output: DC 12V, 2.5A

Device Power Supply	Input: DC 12V
Shape Size:	280×250×90mm
Software Version:	V1.0
Expected service life	5years

Form 2.6

System connection diagram


Fig 2.2 System connection diagram

Warning: Additional equipment connected to the medical electrical equipment must comply with the respective IEC or ISO standards (e.g. IEC60950 for data processing equipment). Furthermore all configurations shall comply with the requirements for medical electrical systems (see clause 16 of IEC60601-1 edition 3.2). Anybody connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. When more than one device is connected to use, the accumulation of leakage current may cause the danger of security. If in doubt, consult your local representative or the technical service department.

Chapter 3 Connect and Inspection

3.1 Connect

Prepare this endoscope video image processing device and compatible device before each use. connect the equipment according to the procedures described in this chapter and instructions.

Caution

- When the endoscope camera system is connected to other equipment, please confirm the improper connection, which may cause uncertain electrical safety risks.
- Terminal is accessible to operator with the video endoscope processor in any position of normal use.
- Terminal for connection of a potential equalization conductor use the magnetic suction interface can be effectively fixed to avoid accidental disconnection.
- Because of the magnetic suction interface, unplugging does not require any tools.
- Terminal not used for a protective earth connection
- Please confirm that the monitor can be used with the video endoscope processor. Using an incompatible monitor will damage the device or display the image correctly.
- Please ensure that the supporting monitor meets the electromagnetic compatibility standards and safety standards of medical electrical equipment.
- Power supply cord does not incorporate a potential equalization conductor, When the endoscope camera system is connected to other equipment, please confirm the improper connection, which may cause uncertain electrical safety risks.

3.1.1 Connect the monitor

Caution

- Do not place any object on top of the endoscope video image processing device. Otherwise, it may cause equipment deformation, abnormal function and damage.
- Do not install video endoscope processor devices (such as near microwave therapy, SMS, MRI, radio, or portable phone) where exposed to strong electromagnetic field radiation. Otherwise, the video endoscope processor device will be abnormal.
- Read this chapter thoroughly and prepare the equipment correctly before use. Otherwise, it may affect the performance of the device, and cause electric shocks, patient and operator burns, or fire.
- Only the supporting cables may be used. Otherwise, it may cause equipment damage, or abnormal function.
- If a non-medical electrical peripheral device is used, connect the power cord through an insulated converter before connecting to an endoscope video image processing device. Otherwise, electrical shocks, burns, and / or fires can occur.
- Connect all the cables correctly and firmly. If the cable connector has connection screws, screw it down. Otherwise, it may cause equipment damage, or abnormal function.
- Do not overbend, pull, twist, or squeeze the cable. Otherwise, it may cause cable damage.
- Do not overpress the connector. Otherwise, the joint damage may be caused.

Applicable monitor: E190DH and E240D

3.1.1.1 The E190DH signal lines were connected by using one of the three medium methods below.

3.1.1.1.1 Use the DVI-D to DVI-D Cable

Use the DVI-D to DVI-D Cable connect the monitor (E190DH) to the video endoscope processor, as shown in Fig 3.1.


Fig 3.1

3.1.1.1.2 Use the CVBS Video Cable (BNC to BNC)

Use the CVBS Video Cable (BNC to BNC) connect the monitor (E190DH) to the video endoscope processor, as shown in Fig 3.2.


Fig 3.2

3.1.1.1.3 Use the S-Video HD Cable

Use the S-Video HD Cable connect the monitor (E190DH) to the video endoscope processor, as shown in Fig 3.3.


Fig 3.3

3.1.1.2 The E240D signal lines were connected by using one of the two medium methods below.

3.1.1.2.1 Use the DVI-D to DVI-D Cable

Use the DVI-D to DVI-D Cable connect the monitor (E240D) to the video endoscope processor, as shown in Fig 3.4.


Fig 3.4

3.1.1.2.2 Use the HD-SDI Cable

Use the HD-SDI Cable connect the monitor (E240D) to the video endoscope processor, as shown in Fig 3.5.


Fig 3.5

3.1.1.3 Connection of the equipotential lines

3.1.1.3.1 The E240D and E190DH equipotential cable is connected to the Equipotential ground as shown in Fig 3.6.


Fig 3.6

3.1.2 Connection of the power adapter

Connect the secondary connector of the power adapter to the video endoscope processor as shown in Fig 3.8.


Fig3.8

Caution

- Appliance coupler or mains plug of power adapter is used as the isolation means from supply mains, not to position the device so that it is difficult to operate the disconnection device.

3.2 Inspection

Every 6 months, please check the adapter plug for dust every six months. If there is any dust, please clean it up before you use it.

Before each operation, please prepare the video endoscope processor devices and other peripheral devices.

Check the video endoscope processor.

3.2.1 Inspection precautions

Caution

- Read chapter 3.1, "Connect", thoroughly before checking and prepare the equipment correctly. If the equipment is not properly prepared before each use, it may lead to improper performance, electric shock, patient and operator injuries, or a fire.
- The video endoscope processor device should be examined before each procedure as described below. Other devices used with this video endoscope processor device are checked as described in the respective instructions. Do not use the video endoscope processor, even if any abnormality is found, and please consult Chapter 9, "Troubleshooting. If there are still abnormalities after treatment as described in Chapter 9, please contact Zhuhai Vision Medical Technology Co., Ltd. Ltd. Otherwise may cause the device to perform performance, electric shock, patient and operator injury, or fire.
- Other devices used with this video endoscope processor are checked as described in the respective instructions. Do not use the device even if any abnormality is found. If an abnormal equipment is used, it may cause abnormal equipment function, or damage, as well as electric shock, burn, or fire.

3.2.2 Connected endoscopy

Caution

- Before connecting the video connector to the endoscope video image processing device, confirm that the video connector, including the electrical contact, is completely dry and free of foreign objects, such as residual detergent, hard water residues, finger prints, dust, and lint. If the endoscope is operated with wet hands, or if the electrical contact is contaminated, it may cause functional abnormalities of the endoscopic and video endoscope processor.
- Check the extension cable appearance and determine that has no blistering, cracking, breakage, wire has no kinking, the Lemo plug (male) has no mechanical damage, and the plug connected to the endoscope operation department can rebound after pressing the tip.
- Insert one end of the extension cable (Lemo plug) into the video endoscope processor output interface, the Lemo plug (mother plug). Clicks were heard after insertion.

- Place the extension cable at the other end (endoscopic manipulation handle connection end) from top to bottom above the operant. Align the black dots on the to turn axis to the black dots above the maneuver (endoscopic maneuver handle connection end) rotate from right to left (endoscopic maneuver handle connection end) by gently clasping the (endoscopic maneuver handle connection end) to the maneuver, and after hearing a "click", the illustration follows:

3.2.2.1 The Connection for the extension cable

Connect the LEMO connector of the extension cable to the Video Endoscope Processor as shown in Figure 3.9.


Fig3.9

3.2.2.2 Connection of the wireless transmitter and extension cable to the endoscope

Fig 3.10, Place the wireless video transmitter above the operation unit, align the black mark point ● " on the right side of the white balance button of the wireless transmitter with the black mark point ● " at the electrical connection port on the handle of the endoscope operation section, and then press it in opposite directions. After hearing a click, the wireless transmitter is successfully connected with the endoscope operation section.


Fig3.10

3.2.3 Inspection power adapter

1. Check the appearance of the power adapter to ensure that the plastic shell of the power adapter and the outer skin of the power cord are not damaged. The power adapter with abnormality must be stopped immediately, because the faulty equipment may cause electric leakage.
2. Hold both ends of the power adapter by hand and pull gently to ensure that the power adapter can be easily straightened without bending of the wire.
3. The output head of the power adapter is inserted into the power interface of the video endoscope processor, and the other head of the power adapter is inserted into the mains power.
4. The Power adapter output connector acts as a separate plug for network power separation, and do not place the device where it is difficult to operate the unplug connector
5. After connecting the adapter, the switch button is on the red light, press the switch button, the power indicator light is off, and the blue light is on.

Warning

- When installing the power adapter, first connect the electronic endoscope image processor and the adapter, and then connect the adapter to the power grid. Otherwise, it will cause electric shock to the host and damage the host.
- The power adapter with damaged appearance or abnormal structure must be stopped immediately, otherwise it may

3.2.4 Check image output

3.2.4.1. After the logo is displayed for 10s, the preview state enters by default, that is, the real-time image display interface of the electronic endoscope image processor is implemented. At this time, the real-time image should be observed normally, and the image is clear, without blur and incomplete image abnormalities.

Warning

- Do not use this video endoscope processor if the image displays abnormally or the real-time image cannot be observed. Otherwise, it may lead to patient injury.

3.2.5 Lighting and automatic dimming inspection

Warning

- Do not look directly at the column. Otherwise, it can lead to an eye injury.

3.2.5.1 Under the main menu, press the Function Switch " button. Enter the function switching setting interface, and select the dimming key function selection option under the word "dimming" (as shown in

Fig 3.11). Then return to the preview interface.


Fig3.11

3.2.5.2 In Figure 3.12, confirm that the LED will glow at the foreend end, and press the dimming button to confirm that the LED brightness changes with the adjustment.


Fig3.12

3.2.6 Check the freezing function

Warning

- Do not use this video endoscope processor if no real-time endoscopic image is observed. Otherwise, it may lead to patient injury.

3.2.6.1 Under the main menu, press the "function switch" button. Enter the function switch setting interface, and select "freeze" below (as shown in Fig 3.13). Then return to the preview interface.


Fig3.13

- a. Press the light button switch. Confirm the real-time endoscopic image freezing.
- b. Press the light button switch again to confirm that the frozen image returns to the real-time image.

3.2.7 Check the zoom function

3.2.7.1. Under the main menu, press the "function switch" button. Enter the function switching setting interface, and select "Scale" below (as shown in Fig 3.14). Then return to the preview interface.


Fig3.14

- 3.2.7.2. Press the light button switch. Verify that the real-time endoscopic image display screen is full.
- 3.2.7.3. Press and lower the light button switch again to confirm that the real-time endoscopic image display screen is full, up and down.
- 3.2.7.4. Press the light button switch again to confirm that the real-time endoscopic image restores the original size.

3.2.8 Check the EDI, HBE function

3.2.8.1. Under the main menu, press the "function switch" button. Enter the function switch setting interface, and select "EDI, HBE" below it (as shown in Figure 3.15). Then return to the preview interface.


Fig3.15

3.2.8.2. Press the light button switch. Verify that the real-time endoscopic images are switched to the EDI display mode.

3.2.8.3. Press the light button switch again to confirm that the real-time endoscopic image is switched to HBE display mode.

3.2.8.4. Press the light button switch again to confirm that the real-time endoscopic image returns to normal display mode.

3.2.9 Check the photo and video recording function

Warning

- Do not use this video endoscope processor if not properly photographed and recorded. Otherwise, it may lead to patient injury.

3.2.9.1. Touch the camera button switch under the preview screen. Verify that the photo icon appears on the processor.

3.2.9.2. Long press the photo button for 2~3s. Verify that the video icon appears on the processor, enter the video mode, and start the timing.

3.2.9.3. Long press the photo button again for 2~3s, exit the video mode, and enter the file management interface to view the photographed photos and video files.

Chapter 4 Function and Operation

Caution

- The VS610 can be operated by directly selecting the menu touch on the touch screen, but the VS610 needs to select the corresponding menu by touching the upper and down

4.1 Preview mode

4.1.1. After connecting the video endoscope processor with the endoscope operating handle produced by our company as required, gently press the power button on the left side of the display:  the power indicator lights up. Video endoscope processor or external monitor display company Logo for 10 seconds (see Fig4.1 "Boot image").

Fig4.1Boot image

4.1.2. After starting up, the LOGO shows 10 seconds, and the default is in the preview mode, means real-time endoscope image. (See Fig 4.2, Fig 4.3)


Fig4.2 Real time image (VS600) – not applied


Fig4.3 Real time image (VS610)

4.1.3. The display has two display modes to choose. "standard display mode" and "full-screen display mode": Press the power button on the left side of the monitor, the monitor light lights up. After starting up, the LOGO shows 10 seconds, and the default is in the preview mode, the default is the "standard display mode". (Fig 4.4)

Fig4.4 Standard display mode

In the "standard display mode" mode, press on the right display screen full screen button , you can enter the "full screen mode" (Fig 4.5)

Fig4.5 Full-screen Display Mode

Caution

- Monitor is used of multi touch capacitive touch screen, do not use sharp objects to move on the screen. Otherwise, there may be damage to the monitor.

4.1.4. Press the "video settings" button to enter the video interface settings, it displayed the "Brightness", "contrast", "saturation", "clarity", "default ", "save" 6 options. The "Brightness", "contrast", "saturation", "definition" with level of 20adjustable (maximum value is 10, the minimum value -10). The operator can select brightness, contrast, color saturation and resolution according to the actual needs of the clinic, press "+" or "-" on both sides can adjust the brightness, contrast, color saturation and sharpness. "+" represents an increase in value, "-" represents a reduction in value (as shown in the following picture).


Fig4.6 Video Setup interface (VS600) – not applied

Fig4.7 Video Setup interface (VS610)

4.2 Playback mode

4.2.1. In the main menu mode, click the "playback"  icon on the screen. Enter the playback interface(see Fig 4.8, Fig 4.9), review image or video that has been stored.(See Fig 4.10, Fig 4.11, Fig 4.12, Fig 4.13)


Fig4.8 Playback menu interface (VS600) – not applied


Fig4.9 Playback menu interface (VS610)


Fig4.10 Replay A (Image Replay for VS600) – not applied


Fig4.11 Replay A (Image Replay for VS610)


Fig4.12 Replay B (Video Replay for VS600) – not applied


Fig4.13 Replay B (Video Replay for VS610)

4.3 Setup

In the main menu mode, click the "setup"  icon on the screen. there are five function options: "Video Source", "Date & Time", "Language", "Wireless" and "Storage".

4.3.1. In the "Setup" page, click "video source", there are two modes of "cable" and "wireless". (As shown in the figure below)

Caution

- When choosing wired or wireless connection, it needs to be set in the video endoscope processor. When using wireless connection mode, click "wireless" connection at the video source. When using wired connection mode, click "cable" connection at the video source.
- The mismatched video source selection will lead to no endoscopic image in the video endoscope processor.
- When using the wireless transmitter to connect with the video endoscope processor for the first time, the pairing password "12345678" should be input.
- Switching video source signal equipment software will restart.


Fig4.14 Video Source Setting interface(VS600) – not applied


Fig4.15 Video Source Setting interface (VS610)

4.3.2. In the "Setup" page, press "Date & Time" button, enter the time and date setting interface. (As shown in the figure below)


Fig4.16 Date & Time Setting interface (VS600) – not applied


Fig4.17 Date & Time Setting interface (VS610)

4.3.3. In the "Setup" page, press the "Language" button. Enter the language setting interface. (As shown in the figure below)

Fig 4.18 Language Settings interface (VS600) – not applied

Fig 4.19 Language Settings interface (VS610)

4.3.4. In the "Setup" page, press the " wireless " button. Enter the WIFI network selection interface. (As shown in the figure below)

Note

- The WIFI used in the video endoscope processor is a special channel for the product and will not be mixed with the WIFI of mobile communication.
- When wireless is enabled for the first time, the wireless transmitter should be paired with the video endoscope processor. The pairing password is "12345678".


Fig4.20 Wireless setting interface (VS600) – not applied


Fig4.21 Wireless setting interface (VS610)

4.3.5. In the "Setup" page, press the "Storage" button. Enter the storage setting interface. (As shown in the figure below)

Caution

- There are two kinds of storage modes, flash and SD card, in the video endoscope processor. When using SD card for storage, it is necessary to insert a qualified SD card into the SD card slot.

Note

- Don't remove the SD card while operating, otherwise files may be damaged.


Fig4.22Storage setting interface (VS600) – not applied


Fig4.23Storage setting interface (VS610)

4.4 File Manager

In the main menu mode, click "File Manager"  button. Enter the file management setting interface. You can create a new folder, delete a folder, rename a patient folder (remarks), upload photos and videos stored on the device, comment, delete and other operations under the current date.

4.4.1. In the "File Manager" page, the current date can be used for new folders, folder deletion, upload and other operations. (As shown in the figure below)


Fig4.24 File manager interface (VS600) – not applied

Fig4.25 File manager interface (VS610)

4.4.2. In the "File Manager" page, Click the number at the bottom of the folder, the system will automatically pop up the "patient information" box, and the patient information can be noted. (See Fig4.26, Fig4.27)

Fig4.26Folder rename interface (for VS610 only)


Fig4.27 Folder rename interface (for VS610 only)

4.4.3. In the file management menu, find the video or photo that needs to be commented. Click the video or photo to comment on the video or photo. (See Fig4.28, Fig4.29)


Fig4.28Photo remarks interface


Fig4.29Video remarks interface

4.5Other

In the main menu mode, click "Other" . There are six function options under this function menu: "External Device", "Get App", "Upgrade", "About", "Teaching" and "Help".

4.5.1. In the " Other " page, click " External Device", there are two modes of "Connect PC" and " Connect USB Stick". (As shown in the figure below)

Note

- After connecting with the computer through the USB port through the USB data cable, select "  Connect PC " to copy the pictures or videos stored on the video endoscope processor.
- When using the USB disk connection, select"  Connect USB Stick " to copy the pictures or videos stored on the video endoscope processor. The software can also be upgraded.


Fig4.30 External Device interface

4.5.2. In the "Other" page, Click the "Get App" button to have two-dimensional code for "Android" and "IOS" systems (as shown in the figure below)

Note

- After obtaining the app of the mobile platform by scanning the QR code, it can realize the interconnection with the video endoscopeprocessor and obtain the real-time image of the operation.
- After the mobile terminal is interconnected with the video endoscopeprocessor, it is only used for synchronous observation with the video endoscopeprocessor. You can't control any operation.


Fig4.31 Get App interface

4.5.3. In the "Other" page, Click the "Upgrade" button, with "Upgrade Software" and "Upgrade Firmware" two options. (As shown in the figure below)

Note

- Software update or firmware update can be carried out through the USB disk with upgrade package.

- When it is necessary to upgrade the software, insert the USB disk with the software upgrade package into the video endoscopeprocessor and click "Upgrade software".
- When you need to upgrade the firmware, insert the USB disk with the firmware upgrade package into the video endoscopeprocessor, and click "Upgrade firmware".


Fig4.32 Upgrade interface

4.5.4. In the "Other" page, Click the "About" button to display the software version, system version and operation section information. (As shown in the figure below)

- Note**
- The "operation section" on the interface refers to the endoscope used with the video endoscopeprocessor.
 - The video endoscopeprocessor can display the model, serial number, production date and used times of the endoscope.

Fig4.33About interface

4.5.5. In the "Other" page, Click the "Teaching" button to play the video about the use and operation of endoscope (as shown in the figure below)

Caution The teaching video is only used as the basic information of learning and understanding the use and operation of products, not as the diagnosis and treatment procedures.


Fig4.34Teaching interface

4.5.6. In the "Other" page, Click the "Help" button to view the help content. (As shown in the figure below)


Fig4.35Help interface

4.6 Taking pictures

4.6.1. Press the record button  on the operating handle, the Photo icon at the top of the screen is shown in red, about a second after, the red photo icon is gone. This picture has been added on time and the file name in the built-in SD card successfully saved. (Fig 4.4)

Fig 4.36

Note

- The instrument automatically generates the real-time time and file name for recording photos. File names are mainly recorded in numbers. Visible in playback mode.

4.7Recording video

4.7.1. In the real time mode keep pressing video/camera button  for 3 seconds, enter into record mode, the video icon at the top of screen is shown in red  00: 00: 00 . If you need to exit the recording mode, press the video/camera button  for 3 seconds: you can switch to preview mode to display the image of the current endoscope. (See Fig4.37)

Fig4.37

4.8 White balance adjustment

Warning

- White balance adjustment is used to display the image color correctly in the LCD display. Without the correct white balance may cause color distortion.

4.8.1. Press the power button on the left side of the monitor, the LCD display light lights up, enter preview mode. Let the distal tip of the insertion section aim to white objects or white gauze. Keep pressing white balance  button on the extension cable or wireless transmitter for 3 seconds, until the icon  appear, white balance adjustment succeeds.

4.8.2. Repeat step 1, if adjust fault.

4.9 Delete Photo or Video

This function deletes recorded photos / videos. Photos / videos in the product can be deleted one by one, or all at once.

Warning

- Photo or video once deleted that cannot be recovered. Please check it before deleting, in order to avoid accidentally deleting it that you want to keep.
- In order to avoid photo or video is deleted, transfer all data into computer by USB.

Note

- The devices connected with computer through the USB does not need to install any driver.

Delete photo

In the file management mode, select the image to be deleted and press the "delete" key to delete it. (See fig 4.38)


Fig 4.38

Delete video

In the file management mode, use the key to select the video to be deleted, and press the "delete" key to delete it. (See Fig 4.39)


Fig 4.39

4.10 power off

When the equipment is not in use, gently press the power button on the left side of the display:  Then disconnect the power adapter and power plug.

The device should be placed in a suitable position that the operator can disconnect it.

Chapter 5 Storage and Disposal

Keep the clean video endoscope processor and accessories after cleaning and sterilization away from the contaminated equipment. If the clean video endoscope processor and accessories are contaminated, their use may pose an infection control risk to the patient and/or operators performing the next procedure with the video endoscope processor.

Environment specifications

Normal working conditions	Environment temperature:	5°C ~ +40°C
	Atmospheric pressure:	700hPa ~ 1060hPa
	Relative humidity:	30% ~ 85% (Non – Condensing)
Transportation and storage	Environment temperature:	-10°C ~ +55°C
	Atmospheric pressure:	500hPa ~ 1060hPa
	Relative humidity:	≤93% (Non – Condensing)

Form 5.1

5.1 Storage of the video endoscope processor

5.1.1. Remove all attachments from the video endoscope processor.

5.1.2. Verify that all surfaces of the electronic endoscopic image processor are completely dry.

Caution

To prevent contamination clean video endoscope processor and accessories after cleaning are always kept in place.

- The places where the video endoscope processor is stored must be clean, dry, well ventilated and kept at room temperature. Never store it in an environment that has direct sunlight, high temperature, humidity, or is exposed to ozone, X-rays, or UV light, which can cause damage to the video endoscope processor.
- Never store the video endoscope processor inside the porter box. Porter boxes should only be used to transport the video endoscope processor.

5.2 Transporting Outside the Hospital

Warning

- Transport the video endoscope processor in the carrying case.
- Be sure to clean the surface before transporting the video endoscope processor outside the hospital. Otherwise, the patient might be infected.

5.3 Disposal of expired or scrapped

When the video endoscope processor reaches the expiration date or can no longer be used, the disposal of the video endoscope processor and any components (such as power adapter, etc.) shall comply with all applicable national and local laws and regulations.

-  □ The wireless transmitter contains a built-in battery. When scrapping, do not put it into fire or dispose with domestic waste, otherwise it may cause explosion or environmental pollution.
- The circuit board of this product contains trace heavy metals and other harmful substances. Do not discard it with domestic waste, otherwise it will cause environmental pollution.

5.1.1. If possible, please disassemble the wireless transmitter and remove the battery before handling it separately.

5.1.2. Please comply with local laws and regulations for the disposal of abandoned electrical and electronic equipment.

Chapter 6 Troubleshooting

Some problems that appear to have irregularities may be correctable by referring to Section 6.1, "Troubleshooting guide". If the problem cannot be resolved by the described remedial action, stop using the video endoscope processor and send it to us for repair.

Our company does not repair accessory parts. If an accessory part becomes damaged, contact us to purchase a new one.

Warning

- Never use the video endoscope processor on a patient if an abnormality is suspected. Damage or irregularity in the instrument may compromise patient or user safety and may result in more severe equipment damage.

6.1 Troubleshooting Guide

Troubles or failures due to other causes than those listed below should be serviced. As repair performed by persons who are not qualified by our company could cause patient or user injury and/or equipment damage, be sure to contact us for repair according to Section 6.2, "Returning the endoscope for repair".

Image Quality

Irregularity description	Possible cause	Solution
There is no video image.	Not all power switches are ON.	Turn ON all the power switches.
	Improper connection between video endoscope processor and endoscope operation unit	Correctly connect the endoscope operation unit and the video endoscope processor according to the requirements of the manual.
An image is not clear.	The objective lens is dirty.	Wipe with cotton swab moistened with 70% ethyl or isopropyl alcohol. If the image cannot be observed, return the video endoscope processor and the supporting endoscope operation to our company for repair.

Form 6.1

Key

Irregularity description	Possible cause	Solution
Cannot Soft Shutdown	Software system crash	Press the left power switch button to restart
No response when touching virtual key	Touch screen failure	

Form 6.2

Photograph

Irregularity description	Possible cause	Solution
No picture generated	Internal memory or SD card is full.	Delete unnecessary folders or files. Download the important pictures to PC before delete.
The image data are not recorded in accordance with the date.	The date and time are not set.	Set the date and time.
	Low battery power with built-in clock.	If it is still wrong with the date and time after re-setting. Return the endoscope to our company for repair.

Form 6.3

6.2 Returning the video endoscope processor for Repair

Warning

- Before repairing the video endoscope processor, please thoroughly clean its surface with 70% alcohol. Improper cleaning will bring infection risk to the hospital or the company's personnel handling the equipment.

Caution

- Our company is not responsible for the equipment damage and/or person injury which are not caused by our company personnel trying to repair the equipment.

Before returning the instrument for repair, contact our company. With the instrument, include a description of the malfunction or damage and the name and telephone number of the individual at your location who is most familiar with the problem. Also include the warranty card. When returning the instrument for repair, follow the instructions given in Section 5.2 "Transporting outside the hospital".

Chapter 7 Cleaning and Disinfection

7.1 Maintain

The expected service life of the product is 5 years, the entire product life period, the product can work normally, but does not exclude the use of the process or the damage in the transportation process, if the power line is found to be squeezed, the plug is deformed or the appearance of the adapter or the appearance of the processor is damaged, you need to contact the manufacturer for maintenance confirmation.

Caution

- Do not provide service or maintain while it is being used by the patient, which may harm the patient and the operator.
- Disconnect the video endoscope processor from the endoscope, remove any accessories, and completely turn off the video endoscope processor before cleaning and disinfection.
- Do not let water or any liquid drip into the video endoscope processor. When using cleaning or disinfectant, follow the manufacturer's specifications for exposure time and product life.
- After use, the video endoscope processor must be stored according to local guidelines until the next use.

7.2 Cleaning

7.2.1 Recommended cleaning agents

Please use the detergent recommended by manufacturer as below.

Cleaning agent	Metrex Empower 1
Manufacturer	Metrex

Form 7.1

7.2.2. Manual cleaning step at room temperature

- 1) Turn off the electronic endoscopic image processor.
- 2) Prepare cleaning solution using recommended enzyme detergent.
- 3) Permeabilize a piece of sterile gauze in the enzymatic cleaning solution and ensure that the gauze is wet but does not drip.
- 4) Clean the electronic endoscope image processor's buttons, screens, housing thoroughly with moistened gauze. Avoid wetting the unit to prevent damage to the internal electronic components.

- 5) Wait 10 min (or press the recommended time from the detergent manufacturer) to allow the enzyme to work.
- 6) Wipe the electronic endoscopic image processor using sterile gauze that has been moistened with Rodi water. Ensure that all traces of cleaning fluid have been removed.
- 7) Repeat steps 1-6 until the electronic endoscopic image processor surface is clean and visually inspected free of stains and water.

7.3 Disinfection

Caution

- It is ineffective to do disinfection before active cleaning, which easily leads to coagulation of proteins and organic residues on the equipment and prevents correct disinfection.

7.3.1. Disinfectant

Disinfectant	Concentration
Ethyl alcohol	75%

Form 7.2

7.3.2. Manual disinfection procedure at room temperature

- 1) Prepare a piece of gauze dipped with 75% Ethyl Alcohol to ensure that the gauze is wet but not dripping, so as to avoid any liquid dripping into the video endoscope processor.
- 2) Wipe the buttons, screens, housings, slots and gaps on the video endoscope processor.
- 3) Place the product at room temperature for 10 minutes.
- 4) Wipe the surface of the video endoscope processor with dry gauze.

Chapter 8 Electromagnetic compatibility

Electromagnetic compatibility

Caution:

- VS610 series video endoscope processor meets the requirements of IEC 60601-1-2 standard for electromagnetic compatibility.
- The user shall install and use the electromagnetic compatibility information provided by the random files.
- Portable and mobile RF communications devices may affect the performance of the VS610 series video endoscope processor and avoid strong electromagnetic interference such as near mobile phones and microwave ovens.
- The guidelines and manufacturer's declarations are detailed in the annex.

Warning:

- A device or system should not be used with other devices placed close to or over stacked, if must be close to or stacked, should observe and verify the normal operation in its use of the configuration.
- Type A equipment for the industrial environment, since the UF series of Video Endoscope conduction harassment and radiation harassment, in other environments to ensure electromagnetic compatibility may have potential difficulties.
- In addition to the cables sold as spare parts of internal components by the manufacturer of VS610 series video endoscope processor, the use of accessories and cables outside the regulations may increase the emission or reduce the immunity.
- Even if other devices meet the emission requirements of corresponding national standards, VS610 series video endoscope processor may still be disturbed by other devices.

Type	Launch frequency	modulation type	Radiated Power
WIFI	2.4G/5.8G	ISM/UNII	< = 2dB

Sample cable

No.	Item	Length of cable	Shielding(Y/N)	Remark
1	Power cable	1.5	N	
2	Adapter connecting cable	1.4	N	
3	Equipotential Cable	1.8	N	
4	DVI to DVI Cable	1.5	Y	

No.	Item	Length of cable	Shielding(Y/N)	Remark
5	USB-A to USB-A Cable	1.5	Y	
6	Extension Cable	1.5	Y	(EC182)
7	Extension Cable	2.2	Y	(EC181)
8	HDMI To DVI Cable	1.5	Y	
9	HD-SDI Cable (BNC to BNC)	1.5	Y	
10	CVBS Video Cable (BNC to BNC)	1.5	Y	
11	S-Video HD Cable	1.5	Y	

Annex:

Guide and manufacturer's statement - Electromagnetic Emission		
VS610 series of Video Endoscope Processor used in the electromagnetic environment of the expected performance of the following provisions, the performance of VS610 series of Video Endoscope Processor purchasers or users should ensure that it is used in the electromagnetic environment:		
launching test	compliance	Electromagnetic environment - Guide
IEC/CISPR 11 RF	Group 1	The performance of VS610 series video endoscope processor is only its internal function and the use of RF energy. Therefore, its RF emission is very low, and may not cause any interference to nearby electronic equipment. The VS610 series video endoscope processor is suitable for use in non-household and not directly connected to the public facilities for all residential low-voltage power supply network in the home
IEC/CISPR 11 RF	Class A	
IEC6100-3-2	NA	
IEC6100-3-3	NA	

Guide and manufacturer's statement - Electromagnetic Immunity			
VS610 series of Video Endoscope Processor used in the electromagnetic environment of the expected performance of the following provisions, the performance of VS610 series of Video Endoscope Processor purchasers or users should ensure that it is used in the electromagnetic environment:			
immunity test	IEC60601-1 test level	Compliance level	Electromagnetic environment - Guide
Electrostatic discharge (ESD) IEC6100-4-2	±6 kV Contact discharge ±8 kV air discharge	±6 kV Contact discharge ±8 kV air discharge	The ground should be wood, concrete or tile, if the ground with synthetic material covering, the relative humidity should be at least 30%
electrical fast transient IEC6100-4-4	±2kV for the power cord ±1kV the input/output line	±2kV for the power cord	The power network should have the quality of typical commercial or hospital environment.
Surge IEC6100-4-5	±1 kV Differential-Mode Voltage ±2 kV common-mode voltage	±1 kV Differential-Mode Voltage	The power network should have the quality of typical commercial or hospital environment.
Voltage sags, short interruptions, and voltage variations on the power input cable IEC6100-4-11	<5% UT(>95% dip in UT) for 0.5 cycle, 40% UT (60% dip in UT) for 5 cycles, 70% UT (30% dip in UT) for 25 cycles, <5% UT (>95% dip in UT) for 5 sec	<5% UT(>95% dip in UT) for 0.5 cycle, 40% UT (60% dip in UT) for 5 cycles, 70% UT (30% dip in UT) for 25 cycles, <5% UT (>95% dip in UT) for 5 sec	The power network should have the quality of typical commercial or hospital environment. If the VS610 series of video endoscope processor performance users need to run continuously in the power interruption, it is recommended that the VS610 series of videoendoscope processor performance using the uninterruptible power supply or battery power.
(50/60Hz) IEC6100-4-8 Power frequency magnetic field	3A/m	3A/m(50Hz/60Hz)	The power frequency magnetic field shall have the horizontal magnetic field characteristics of a typical place in a typical commercial or hospital environment
Note: UT refers to an alternating voltage test voltage before			

Guide and manufacturer's statement - Electromagnetic Immunity			
VS610 series of Video Endoscope Processor used in the electromagnetic environment of the expected performance of the following provisions, the performance of VS610 series of Video Endoscope Processor purchasers or users should ensure that it is used in the electromagnetic environment:			
immunity test	IEC60601-1 test level	Compliance level	Electromagnetic environment - Guide
RF transmission IEC6100-3	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile communication equipment should not be part of any RF closer to the VS610 series of video endoscope processor performance than the isolation distance recommended use, including cable. The distance should be with the transmitter frequency corresponding to the formula. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Among them, P is based on maximum output power supplied by rating of the transmitter manufacturers, Watt (W) as a unit, D is recommended isolation distance,mas a unit. Fixed field RF transmitter based on
RF radiation IEC6100-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

			the electromagnetic field survey to determine, in each frequency range should be lower than compliance voltage In the following symbols marking equipment, interference π :ur.
--	--	--	--

Note 1: in 80MHz and 800MHz frequency, adopt the high frequency formula.

Note 2: these guidelines may not be appropriate in all cases, the electromagnetic propagation may be affected by buildings, absorption and reflection of objects and people.

a Fixed transmitter electric field strength, such as base station of wireless (cellular / cordless) telephone and ground mobile radio, amateur radio, AM (AM) and FM (FM) radio and television broadcasting etc.their emission electric field strength is impossible to be predicted in theory. In order to evaluate the electromagnetic environment of a fixed RF transmitter, the investigation of electromagnetic field should be conducted. If the field strength of the VS610 series video endoscope processor is measured higher than the RF compliance voltage of the above-mentioned application, the performance of the VS610 series video endoscope processor should be observed to verify its normal operation. If the observed performance is abnormal, then additional measures may be required, such as the new series of VS610video endoscope processor performance orientation or location.

b In the whole frequency range of 150KHz ~ 80MHz, the emission electric field strength should be less than 3 V/m.

Recommended isolation distance between portable and mobile communication equipment RF and VS610 series of Video Endoscope Processor performance			
The performance of the VS610 series video endoscope processor is expected to be used in the controlled electromagnetic environment of radiated RF harassment. On the basis of the maximum output power of communication equipment, the purchasers or users can maintain portable and mobile RF communication equipment through the following recommended (transmitter) and VS610 series video endoscope processor performance minimum distance to prevent electromagnetic interference.			
Rated maximum output power of the transmitter/W	Corresponding to different frequency transmitter isolation distance/m		
	150 kHz ~ 80 MHz $d = 1.2\sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2\sqrt{P}$	800 MHz ~ 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
The transmitter rated maximum output power are not listed in the form, the recommended isolation distance D, in meters (m) as a unit, can use the corresponding transmitter frequency bar formula to determine, where P is the maximum rated output power provided by the transmitter of transmitter manufacturers, by Watt (W) as a unit. Note 1: in 80MHz and 800MHz frequency, adopt the high frequency formula. Note 2: these guidelines may not be appropriate in all cases, the electromagnetic propagation may be affected by buildings, absorption and reflection of objects and people.			

Manufacturers: Zhuhai Vision Medical Technology Co., Ltd.

Address: 301 Building 4, No. 30 Feipeng Road, GMTCM Park, Hengqin New Area,

519031 Zhuhai City, Guangdong Province, PEOPLE' S REPUBLIC OF CHINA

Zip code: 519031

Tel: 4000756365

Fax: +86-756-6910511

E-mail: info@vision-medt.com

service@vision-medt.com

Website: <https://www.vision-medt.com/>

MINDSION

MINDSION

User manual

—Video Endoscope Processor

VS50
VS100

Zhuhai Vision Medical Technology Co., Ltd.

Product Information

- Product Name: Video Endoscope Processor
- Product Type: VS50/VS100
- Manufacture Name: Zhuhai VisionMedical Technology Co., Ltd.
- Address: 301 Building 4, No. 30 Feipeng Road, GMTCM Park, Hengqin New Area, 519031 Zhuhai City, Guangdong Province, PEOPLE' S REPUBLIC OF CHINA
- Contact information
 - Tel: +86-756-6910150
 - Fax : +86-756-6910511
 - Zip code : 519031

Version Information

This manual version may be upgraded due to a change of software or technical specifications at any time without notice. The version information of this specification are as follows:

- Document number: TD/CE-VS-A-009
- Document version: V1.2
- Released date: Jan.2024

©2024, The copyright of Zhuhai Vision Medical Technology Co., Ltd

CE Mark

EC Representative Name:

Shanghai International Holding Corp. GmbH (Europe)

EC Representative Address:

Eiffestraße 80, 20537 Hamburg, Germany

Statement

The manufacturer owns the rights to the unpublished use of the instructions and has the right to treat them as confidential information. The use of the manual is only a reference for the operation, maintenance and repair of the products, and no other person has the right to disclose the use of the instruction to others.

This product specification contains by copyright law protection of proprietary information, all rights

reserved, without manufacturer's written consent shall not be any part of this instruction manual for camera copy, copy, or translated into other languages.

The contents of this manual shall be subject to change without notice.

Responsibilities of the Manufacturer

In the following conditions are met, the manufacturer will responsible for the product safety, reliability and performance:

- Product installation, maintenance or upgrades were carried out by the manufacturer authorized or approved personnel;
- The storage environment, product environment, electrical environment shall meet product specifications;
- Use this product in accordance with < use instruction >
- Use this product within the period of use

Overview

◆ Main Structure

The product mainly consists of the processor, power adapter, extension cable, wireless video transmitter and video signal transmission cables.

◆ Basic performance of the product

When used with the video endoscope produced by our company, it can normally accept the image signal output by the video endoscope and ensure that the image signal is not lost.

◆ Intended Use

The product is intended used to work with the Video Endoscope manufactured by Zhuhai Vision Medical Technology Co., Ltd. to provide the images for observation, photography and treatment.

◆ Contraindication

Video Endoscope Processor is intended used in conjunction with VideoEndoscope, to provide image display and there is no absolute contraindication.

◆ Cautions, Warnings and Notes

- 1) No user self-service components are contained within this product, and internal repair of instruments must be performed by technical service personnel authorized by the manufacturer.
- 2) This product is not part of the therapeutic armamentarium and cannot be applied to families.
- 3) Operate in strict accordance with the product instructions provided by the manufacturer.

About the Manual

The manual gives a detailed account of the usage, functions and operation methods of the product. Please carefully read and understand the contents of the manual before using the product, so that the product can be correctly used and the safety of both patients and operators can be ensured.

The manual gives an introduction of the product according to the most complete configuration. It is possible that some information herein is not applicable to the product purchased by you. For any queries, please contact our company.

Please place the manual near the product in necessary. So that it can be obtained conveniently and in time when needed.

The manual is suitable for professional clinical medical personnel. The readers shall have the knowledge of medical procedures, practices, terms and working experience that is necessary for endoscope.

All the illustrations provided in the manual are for reference only. The settings or data in the illustrations are possibly different from what actually displayed on the product.

About the product

The product is only in contact with the user, and the material of the contact part is ABS, PVC, PC+ABS, PET. The ingredients and raw materials of ME equipment comply with local laws and regulations.



Every two years to verify the basic safety of the product, detailed details, please consult the manufacturer.

Catalogue

SYMBOLS WORDS	1
KEY POINTS (PLEASE READ BEFORE USE)	2
USER QUALIFICATION	2
COROLLARY EQUIPMENT	2
STANDBY EQUIPMENT	2
MAINTENANCE	2
INAPPROPRIATE MAINTENANCE AND MODIFICATION PROHIBITED	2
SIGNAL WORDS	3
WARNINGS AND CAUTIONS	3
BATTERY AND ADAPTER	5
LCD DISPLAY	6
CHAPTER1 CHECKING THE PACKAGE CONTENTS	7
CHAPTER 2 INSTRUMENT NOMENCLATURE AND SPECIFICATIONS	10
2.1 COMPONENTAND FUNCTION	10
2.2 MENU FUNCTIONS	13
2.3 SPECIFICATIONS	15
CHAPTER 3 PREPARATION AND INSPECTION	17
3.1 BATTERY CHARGING	17
3.2 INSPECTION	17
3.2.1 Inspection precautions	18
3.2.2 Connection of VS50 to endoscope	18
3.2.3 Inspection and connectionof extension cable	19
3.2.4 Connection of wireless transmitter	19
3.2.5 Inspection and connectionof power adapter	20
3.2.6 Check image output	20
3.2.7 Check the photo and video recording function	21
CHAPTER 4 FUNCTION AND OPERATION	22
4.1 PREVIEW MODE	22
4.2 PLAYBACK MODE	24
4.3SETUP	26
4.4 FILE MANAGER	30
4.5 OTHER	32
4.6 TAKINGPICTURES	36
4.7RECORDING VIDEO	37
4.8 WHITE BALANCE ADJUSTMENT	37
4.9 DATASTORAGE	37
4.10 POWER OFF	38
CHAPTER 5 STORAGE AND DISPOSAL	39
5.1 STORAGE OF THE VIDEO ENDOSCOPE PROCESSOR	39
5.2TRANSPORTING OUTSIDE THE HOSPITAL	39
5.3DISPOSALOF EXPIRED OR SCRAPPED	40
CHAPTER 6 TROUBLESHOOTING	41
6.1 TROUBLESHOOTING GUIDE	41
6.2RETURNING THE VIDEO ENDOSCOPE PROCESSOR FOR REPAIR	42
CHAPTER 7 CLEANING AND DISINFECTION	43
7.1 MAINTAIN	43
7.2CLEANING	43
7.3DISINFECTION	44
CHAPTER 8 ELECTROMAGNETIC COMPATIBILITY	45

Symbols Words

The meaning(s) of the symbol(s) shown on the package with the components, the back cover of this instruction manual and/or this instrument are as follows:

Symbols	The meaning of the symbol	Symbols	The meaning of the symbol
	Refer to instruction manual/booklet		Caution
	Type BF Applied Part		CE Mark
	Serial Number		This way up
	Fragile; handle with care		Keep away from rain
	Stacking limit by number		Symbol for Individually Disposal Electrical Appliance and Electronic Equipment (Please follow Local Ordinances or Regulations of disposal.)
	Manufacturer		Date of manufacture
	Non-ionizing electromagnetic radiation		Authorized Representative in the European Community
IPX3	The degree of protection against ingress of water is 3 (Wireless Transmitter) .	IPX1	The degree of protection against ingress of water is 1 (equipment) .
	Standby		Medical Device

Key points (please read before use)

User qualification

If there are official standards for use established by the medical administration or other organizations, please follow these standards. If there is no official standard, the operator of this product must be a doctor who has passed the audit of the medical safety director of the hospital or a doctor who can safely carry out the established operation according to the official regulations designated by the director of the using department according to the difficulty of using the technology.

The physician followed the legal standards set by the Endoscopy Society and other institutions while considering the difficulties of various endoscopic examinations and treatment, and was able to safely perform the endoscopy and treatment procedures. Therefore, this instruction manual does not explain or discuss any clinical endoscopic technique.

Corollary equipment

Refer to Chapter 1 "Checking the Package Contents", confirm that this instrument is compatible with the ancillary equipment being used. Using incompatible equipment can result in patient or operator injury and/or equipment damage.

This product conforms to the EMC standard IEC 60601-1-2.

Standby equipment

Be sure to prepare standby equipment or equipment with similar functions to avoid accidents in case of equipment failure or abnormal function during operation.

Maintenance

Please check the equipment according to the instructions in the manual before each use. Furthermore, in addition to the checks before each use, the person responsible for the maintenance of the medical device should also regularly check the parts mentioned in this instruction. If any abnormalities are observed, do not use them and follow the instructions in Section 6.1, Troubleshooting Guidelines. If any abnormalities are observed after examination, please contact with manufacturer.

Inappropriate maintenance and modification prohibited

This product does not contain any user-serviceable parts. Do not disassemble, modify or attempt to repair, otherwise it may result in injury to the patient or damage to the equipment. If necessary, please contact the manufacturer.

Zhuhai Vision Medical Technology Co., Ltd. is not responsible for the injury or damage caused by the maintenance equipment of personnel not from Zhuhai Vision Medical Technology Co., Ltd.

Signal Words

The following signal words are used throughout this manual:

Danger

- Indicates an emergency hazardous situation which, if not avoided, could result in death or serious injury.

Warning

- Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

Caution

- Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to warn against unsafe operations or potential equipment damage.

Note

- Indicates additional helpful information.

Warnings and Cautions

Follow the warnings and cautions given below when handling this instrument. This information is to be supplemented by the warnings and cautions given in each chapter.

Warning

- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- Modification of this equipment is not allowed
- Avoid scald operator, continuous contact with the product for less than 10 minutes.
- A Class equipment is intended to be used in industrial environment. The Video Endoscope Processor may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances.
- Except of the cable supplied sold as spare parts of internal components by manufacturers of the video endoscope processor, the others use of accessories and cable may result in increased product launch interference or low immunity.
- No non-manufacturer adapter. Please confirm whether it is the manufacturer's accessories before accessing the adapter.

Caution

- Instructions warning: the waterproof grade of the equipment is IPX1, do not stay that the equipment is placed in a liquid dangerous environment, and the liquid intake will cause equipment damage or patient injury.
- Replacement of battery by untrained personnel causes hazard (source) (e. g. overtemperature, fire or explosion)
- Devices should be placed reliably to avoid falling.
- Heavy objects cannot be placed on the device, otherwise it will cause device damage or injury to the patient
- Battery replacement is performed by the manufacturer
- Device output is abnormal, please stop using it immediately
- If the liquid enters the video endoscope processor, please turn off the power immediately to prevent fire or electric shock. If liquid has penetrated into the video endoscope processor, do not use it because its internal circuit may be damaged.
- The power adapter is part of the endoscopic image processor.
- The wireless transmitter contains rechargeable batteries inside, but the batteries do not require professional maintenance. If it is not in use for a long time, the wireless transmitter needs to be charged every three months, otherwise the lithium battery will cause irreversible capacity loss due to the low voltage caused by self-discharge. Please contact the manufacturer with any questions.
- The video endoscope processor contains a signal input and output part. The signal input port (Lemo interface) of the operation Department can only be connected with the endoscope produced by the company through the video connection line attached to the product. The HDMI output interfaces must be connected with monitors meeting the requirements of IEC60601-1-2 and IEC60601-1 standards.
- When it is used with other medical electrical equipment, medical equipment meeting the performance requirements of IEC60601-1 and electromagnetic compatibility shall be used to ensure the safety when used with video endoscope processor.
- The user shall install and use according to the EMC information provided in the attached documents.
- Portable and mobile RF communication devices may influence video

endoscope processor performance and avoid strong electromagnetic interference when used, such as close to mobile phones, microwave ovens.

- Electromagnetic interference may occur on this instrument near equipment marked with the following symbol or other portable and mobile RF (Radio Frequency) communications equipment such as cellular phones. If electromagnetic interference occurs, mitigation measures may be necessary, such as reorienting or relocating this instrument, or shielding

the location. 

- In order to check the electromagnetic interference from other devices (any devices or equipment components of other system), the system should be observed in the operating environment whether it is working properly.

Battery and adapter

Operate according to the following regulations to prevent battery leakage, overheating, burns, explosion or electric shock.

Danger

- The battery used in this instrument is a built-in lithium battery. Please use the specified battery charger to charge the battery. Do not use any other charger.

Warning

- Do not disassemble or modify the battery in any way.
- Do not use any charger other than those recommended by the company. Failure to do so may result in fire, explosion, leakage or overheating.
- If the rechargeable batteries are not fully charged within the specified time, stop charging immediately and do not use them again. Because the battery may have been damaged.

Caution

- Before placing the battery in the charger and endoscope camera system, be sure to carefully check the battery for leakage, discoloration, crimping or any other abnormalities.

Note

- Generally speaking, the battery performance will temporarily decline with the decrease of ambient temperature. When the instrument is used in cold environment, the performance of the battery will recover after returning to normal temperature.
- It takes about 120 minutes to charge the battery to the full state.

- The power consumption of video endoscope processor varies with different functions.

LCD display




Note


- Do not squeeze the LCD display. Otherwise, the image may become blurred or damage the LCD display.
- If the ambient temperature exceeds the normal operating range, the LCD display may take a long time to open, or its color may temporarily change. LCD displays that exhibit poor performance due to high and low temperatures recover at normal temperatures.
- The LCD used for display is made of high-precision technology. Due to its characteristics or different viewing angles on the LCD display, the color and brightness of spots may be inconsistent. This is not a functional exception.

Chapter1 Checking the package contents

Match all items in the package with the components shown below. Inspect each item for damage. If the instrument is damaged, a component is missing or you have any questions, do not use the instrument; immediately contact our company.

Items	Picture
Video Endoscope Processor	 A compact, silver-colored video endoscope processor with a black screen and a small control panel below it. The brand name 'MINDSION' is visible on the front. VS50  A larger, silver-colored video endoscope processor with a black screen and a prominent horizontal control panel below it. The brand name 'MINDSION' is visible on the front. VS100

USB Sharp (VS50 only)		HDMI conversion cable (VS50 only)	
Lithium battery (VS50 only)		Android data cable (VS50 only)	
Battery charging stand (VS50 only)		Battery charger (VS50 only)	
HDMI to DVI adapter (VS50 only)		Power Adapter (VS100 only)	
Power cord (VS100 only)		HDMI Converter Cable (VS100 only)	

Wireless video transmitter (VS100 only, optional)		Extension cable (VS100 only)	
U Disk		USB Male To Male Connector	
Certificate Of Inspection, Packing List, Product Acceptance List, Warranty Card, Instructions for Use, Video cables are optional.			

Form 1.1 Package contents

Chapter 2 Instrument Nomenclature and Specifications

2.1 Component and function

2.1.1 Component and function (for VS50)


Fig 2.1 VS50 structure schematic diagram

No.	Name	Function
1	LCD display	Display information such as endoscope image and time.
2	SD card	It is used to store photos and videos taken by endoscopy.
3	Mini USB	The pictures and recorded videos stored in the image processor can be copied by connecting with the computer with USB data cable.
4	Mini HDMI	For Mini HDMI data transmission connection
5	Power Supply Socket	Power on and off. This port is only used by manufacturer.
6	White balance adjustment button	Press the button to adjust the white balance.
7	Release Button	Press the button to release the endoscope.
8	Turn on/off Button	Press the button to turn on the device. Long press the button for about 3s to turn off the device.
9	Battery socket	To insert or remove the battery.

Form 2.1 VS50 Component and function

2.1.2 Component and function (for VS100)



Fig 2.2 VS100 structure schematic diagram

No.	Name	Function
1	LCD display	Display information such as endoscope image and time.
2	SD card	It is used to store photos and videos taken by endoscopy.
3	USB interface	The pictures and recorded videos stored in the image processor can be copied by connecting with the computer with USB data cable.
4	HDMI	For HDMI data transmission connection
5	Audio Port	For audio signal transmission connection
6	IO interface	Connect and control the foot switch
7	Lemo Connection Port	In wired connection mode, extension cable is used to electrically connect the processor with the endoscope.
8	Wireless video transmitter charging port	Used for charging wireless transmitter, providing charging interface.
9	Power Supply Socket	Power on and off
10	Turn on/off Button	Touch the button to turn on the device. Long touch the button for 2.3s to turn off the device.

Form 2.2VS100 Component and function

2.1.3 Wireless transmitter

No.	Name	Function
1	Release button	It is used to separate from the operation department, and can easily separate the wireless transmitter and the operation department after pressing the press
2	power on/off button	Touch power on, press for 2.3s power off, with an indicator light. The green light flashes when the WiFi is not connected, and the green light is on when the WiFi is connected,. At low power when the blue light twinkle.
3	White balance button	White balance after a long press for 2.3s
4	Connection port of operation department	Connect with the operation department

2.2 Menu Functions


Fig2.3Menu interface for VS50

No.	Part Name	Function
1	Date and time	Display real-time wireless status
2	Battery capacity	Displays current power of the device
3	SDcard indicator	SD card status is shown, representing an uninserted SD card when the indicator is black and an inserted indicator card when the indicator is white.
4	Main menu	Return to main menu.
5	Setup	Product video source, time and date, language, wireless, storage and other Settings entry.
6	Play back	Replay video playing or image.
7	Preview	Observe the real-time image.
8	File Manager	Product internal document viewing, transmission and other operations.
9	Other	Product version number information, etc.

Form 2.3 VS50 Video Endoscope ProcessorMenu Functions


Fig2.4Menu interface for VS100

No.	Part Name	Function
1	Date and time	Display real-time wireless status
2	Battery capacity	Displays current power of the device
3	SD card indicator	SD card status is shown, representing an uninserted SD card when the indicator is black and an inserted indicator card when the indicator is white.
4	Battery capacity	Displays current power of Wireless Transmitter
5	Main menu	Return to main menu.
6	Setup	Product video source, time and date, language, wireless, storage and other Settings entry.
7	Play back	Replay video playing or image.
8	Preview	Observe the real-time image.
9	File Manager	Product internal document viewing, transmission and other operations.
10	Other	Product version number information, etc.

Form 2.4 VS100 Video Endoscope ProcessorMenu Functions

2.3 Specifications

Safety Standard

Note

- According to classification rules, this product is classified as Class I. This product is classified according to IEC60601-1 is as follows:

Classification of protection against electric shock	VS50 is internal power supply equipment. VS100 is external power supply equipment, class I.
Classification of protection against electric shock	Type BF
Classification of protection against hazards of explosion	Not suitable for use in places where flammable anesthetic gas is available
Classification of protection against ingress of liquid	Normal equipment (Wireless Transmitter: IPX3)
Classification of running mode	Continuous operation equipment
Classification of disinfection method	Equipment of disinfection method recommended by manufacturer
Methods of disinfection	See relevant chapters in the manual

Form 2.5 Specifications

Model specifications

Model VS50

Model	VS50
Power Supply	DC 3.7V Li-ion Battery
Shape Size:	130×125×44mm
Screen Size:	5 inches
Resolving Power:	1920×1080
Software Version:	V4.0
Product Life:	5years

Form 2.6 VS50 Specifications

Model VS100

Model:	VS100
	Input: 100~240V,50/60Hz
Power Supply:	Output: DC 12V, 2.5A
	DC 3.7V Li-ion Battery
Shape Size:	254×246×60mm
Screen Size:	10 inches
Resolving Power:	1280×800
Software Version:	V4.0
Product Life:	5years

Form 2.7 VS100 Specifications

Chapter 3 Preparation and Inspection

Prepare this endoscope video image processing device and compatible device before each use. connect the equipment according to the procedures described in this chapter and instructions.

Caution

When the endoscope camera system is connected to other equipment, please confirm the proper connection, otherwith may cause uncertain electrical safety risks.

- Terminal is accessible to operator with the video endoscope processor in any position of normal use.
- Terminal for connection of a potential equalization conductor use the magnetic suction interface can be effectively fixed to avoid accidental disconnection,
- Because of the magnetic suction interface, unplugging does not require any tools.
- Terminal not used for a protective earth connection
- Power supply cord does not incorporate a potential equalization conductor, When the endoscope camera system is connected to other equipment, please confirm the proper connection, otherwith may cause uncertain electrical safety risks.

3.1 Battery Charging

Note

- The battery will discharge automatically when not in use. Charge the battery regularly before use.
- When the battery is fully charged, the battery's operating time is still shortening which means its service life may have expired. Please replace the battery with a new one by contact the manufacturer.
- Charging time of lithium battery is approximately 120 minutes.
- It can be charged about 500 times. This varies according to usage.
- When using, battery may be hot. This is not a malfunction.

3.2 Inspection

Every 6 months, please check the adapter plug for dust every six months. If there is any dust, please clean it up before you use it.

Before each operation, please prepare the video endoscope processor devices and other peripheral devices.

3.2.1 Inspection precautions

Caution

- Before use, thoroughly before checking and prepare the equipment correctly. If the equipment is not properly prepared before each use, it may lead to improper performance, electric shock, patient and operator injuries, or a fire.
- The video endoscope processor device should be examined before each procedure as described below. Other devices used with this video endoscope processor device are checked as described in the respective instructions. Do not use the video endoscope processor, even if any abnormality is found, and please consult Chapter 9, "Troubleshooting. If there are still abnormalities after treatment as described in Chapter 9, please contact Zhuhai Vision Medical Technology Co., Ltd. Ltd. Otherwise may cause the device to perform performance, electric shock, patient and operator injury, or fire.
- Other devices used with this video endoscope processor are checked as described in the respective instructions. Do not use the device even if any abnormality is found. If an abnormal equipment is used, it may cause abnormal equipment function, or damage, as well as electric shock, burn, or fire.

3.2.2 Connection of VS50 to endoscope

Place the VS50 processor above the operation unit, align the black dots "●" on the to turn axis to the black dots "●" above the maneuver (endoscopic maneuver handle connection end) rotate from right to left (endoscopic maneuver handle connection end) by gently clasping the (endoscopic maneuver handle connection end) to the maneuver, and after hearing a "click", the illustration follows:


Fig3.1 Connection of VS50 to endoscope

3.2.3 Inspection and connection of extension cable

1. Check the extension cable appearance and determine that has no blistering, cracking, breakage, wire has no kinking, the Lemo plug (male) has no mechanical damage, and the plug connected to the endoscope operation department can rebound after pressing the tip.
2. Insert one end of the extension cable (Lemo plug) into the video endoscope processor output interface, the Lemo plug (mother plug). Clicks were heard after insertion.
3. Place the extension cable at the other end (endoscopic manipulation handle connection end) from top to bottom above the operant. Align the black dots on the to turn axis to the black dots above the maneuver (endoscopic maneuver handle connection end) rotate from right to left (endoscopic maneuver handle connection end) by gently clasping the (endoscopic maneuver handle connection end) to the maneuver, and after hearing a "click", the illustration follows:


Fig3.2 Connection of VS100 to endoscope

3.2.4 Connection of wireless transmitter

Place the wireless video transmitter above the operation unit, align the black mark point "●" on the right side of the white balance button of the wireless transmitter with the black mark point "●" at the electrical connection port on the handle of the endoscope operation section, and then press it in opposite directions. After hearing a click, the wireless transmitter is successfully connected with the endoscope operation section. the illustration follows:


Fig3.3 Connection of wireless transmitter to endoscope

3.2.5 Inspection and connection of power adapter

1. Check the appearance of the power adapter to ensure that the plastic shell of the power adapter and the outer skin of the power cord are not damaged. The power adapter with abnormality must be stopped immediately, because the faulty equipment may cause electric leakage.
2. Hold both ends of the power adapter by hand and pull gently to ensure that the power adapter can be easily straightened without bending of the wire.
3. The output head of the power adapter is inserted into the power interface of the video endoscope processor, and the other head of the power adapter is inserted into the mains power.
4. The Power adapter output connector acts as a separate plug for network power separation, and do not place the device where it is difficult to operate the unplug connector
5. After connecting the adapter, the switch button is on the purple light, press the switch button, the power indicator light , and the blue light is on.

Warning

- When installing the power adapter, first connect the electronic endoscope image processor and the adapter, and then connect the adapter to the power grid. Otherwise, it will cause electric shock to the host and damage the host.
- The power adapter with damaged appearance or abnormal structure must be stopped immediately, otherwise it may cause injury to patients or operators.

3.2.6 Check image output

1. After the logo is displayed for 10s, the preview state enters by default, that is, the real-time image display interface of the electronic endoscope image processor is implemented. At this time, the real-time image

should be observed normally, and the image is clear, without blur and incomplete image abnormalities.

- Do not use this video endoscope processor if the image displays abnormally or the real-time image cannot be observed. Otherwise, it may lead to patient injury.

3.2.7 Check the photo and video recording function

1.Touch the camera button switch under the preview screen. Verify that the photo icon appears on the processor.

2.Long press the photo button for 2.3s. Verify that the video icon appears on the processor, enter the video mode, and start the timing.

3.Long press the photo button again for 2.3s, exit the video mode, and enter the file management interface to view the photographed photos and video files.

- Do not use this video endoscope processor if not properly photographed and recorded. Otherwise, it may lead to patient injury.

Chapter 4 Function and Operation

4.1 Preview mode

1. After connecting the video endoscope processor with the endoscope operating handle produced by our company as required, press the power button on the right side of the display:  the power indicator lights up. Video endoscope processor or external monitor display companyLogo for 10 seconds (see Fig4.1 "Boot image")

Fig4.1 Boot image

2. After starting up, the LOGO shows 10 seconds, and the default is in the preview mode, means real-time endoscope image. (See Fig 4.2)


Fig4.2 Real time image (VS50)


Fig4.3 Real time image (VS100)

3. Press the "video settings" button to enter the video interface settings, it displayed the "Brightness", "contrast", "saturation", "clarity", "default ", "save" 6 options. The "Brightness", "contrast", "saturation", "definition" with level of 20 adjustable (maximum value is 10, the minimum value -10). The operator can select brightness, contrast, color saturation and resolution according to the actual needs of the clinic, press "+" or "-" on both sides can adjust the brightness, contrast, color saturation and sharpness. "+" represents an increase in value, "-" represents a reduction in value (as shown in the following picture).


Fig4.4 Video Setup interface (VS50)

Fig4.5Video Setup interface (VS100)

4.2 Playback mode

1. In the main menu mode, click the "playback"  icon on the screen. Enter the playback interface(see Fig 4.6), review image or video that has been stored.(See Fig 4.7, Fig 4.8, Fig 4.9, Fig 4.10, Fig 4.11)


Fig4.6Playback menu interface (VS50)


Fig4.7 Playback menu interface (VS100)


Fig4.8 Replay A (Image Replay for VS50)


Fig4.9 Replay A (Image Replay for VS100)


Fig4.10 Replay B (Video Replay for VS50)


Fig4.11 Replay B (Video Replay for VS100)

4.3 Setup

In the main menu mode, click the "setup"  icon on the screen. there are five function options: "Video Source", "Date & Time", "Language", "Wireless" and "Storage".

1. In the "Setup" page, click "video source", there are two modes of "cable" and "wireless". (As shown in the figure below)

Caution

- When choosing wired or wireless connection, it needs to be set in the video endoscope processor. When using wireless connection mode, click "wireless" connection at the video source. When using wired connection mode, click "cable" connection at the video source.

- The mismatched video source selection will lead to no endoscopic image in the video endoscopeprocessor.
- When using the wireless transmitter to connect with the video endoscopeprocessor for the first time, the pairing password
- "12345678" should be input.
- Switching video source signal equipment software will restart.


Fig4.12 Video Source Setting interface(for VS100 only)

2. In the "Setup" page, press "Date & Time" button, enter the time and date setting interface. (As shown in the figure below)


Fig4.13 Date & Time Setting interface (VS50)


Fig4.14 Date & Time Setting interface (VS100)

3. In the "Setup" page, press the "Language" button. Enter the language setting interface. (As shown in the figure below)

Fig 4.15 Language Settings interface (VS50)


Fig 4.16 Language Settings interface (VS100)

4. In the "Setup" page, press the " wireless " button. Enter the WIFI network selection interface. (As shown in the figure below)

Note

- The WIFI used in the video endoscope processor is a special channel for the product and will not be mixed with the WIFI of mobile communication.
- When wireless is enabled for the first time, the wireless transmitter should be paired with the video endoscope processor. The pairing password is "12345678".


Fig4.17 Wireless setting interface (for VS100 only)

5. In the "Setup" page, press the "Storage" button. Enter the storage setting interface. (As shown in the figure below)

Note

- There are two kinds of storage modes, flash and SD card, in the

video endoscope processor. When using SD card for storage, it is necessary to insert a qualified SD card into the SD card slot.

Caution

- Don't remove the SD card while operating, otherwise files may be damaged.


Fig4.18Storage setting interface (VS50)


Fig4.19Storage setting interface (VS100)

4.4 File Manager

In the main menu mode, click "File Manager"  button. Enter the file management setting interface. You can create a new folder, delete a folder, rename a patient folder (remarks), upload photos and videos stored on the device, comment, delete and other operations under the current date.

1. In the "File Manager" page, the current date can be used for new folders, folder deletion, upload and other

operations. (As shown in the figure below)


Fig4.20File manager interface (VS50)


Fig4.21File manager interface (VS100)

2.In the "File Manager" page,Click the number at the bottom of the folder, the system will automatically pop up the "patient information" box, and the patient information can be noted. (See Fig4.22, Fig4.23)

Fig4.22 Folder rename interface (VS50 & VS100)


Fig4.23 Folder rename interface (VS50 & VS100)

4.5 Other

In the main menu mode, click "Other" . There are six function options under this function menu: "External Device", "Get App", "Upgrade", "About", "Teaching" and "Help".

1. In the "Other" page, click "External Device", there are two modes of "Connect PC" and "Connect USB Stick". (As shown in the figure below)

Note

- After connecting with the computer through the USB port through the USB data cable, select "  Connect PC " to copy the pictures or videos stored on the video endoscope processor.
- When using the USB disk connection, select "  Connect USB Stick " to copy the pictures or videos stored on the video endoscope processor. The

software can also be upgraded.


Fig4.24External Device interface (VS50 & VS100)

2.In the "Other" page,Click the "Get App" button to have two-dimensional code for "Android" and "IOS" systems (as shown in the figure below)

Note

- After obtaining the app of the mobile platform by scanning the QR code, it can realize the interconnection with the video endoscopeprocessor and obtain the real-time image of the operation.
- After the mobile terminal is interconnected with the video endoscopeprocessor, it is only used for synchronous observation with the video endoscopeprocessor. You can't control any operation.

Fig4.25Get App interface (VS50 & VS100)

3. In the "Other" page, Click the "Upgrade" button, with "Upgrade Software" and "Upgrade Firmware" two options. (As shown in the figure below)

Note

- Software update or firmware update can be carried out through the USB disk with upgrade package.
- When it is necessary to upgrade the software, insert the USB disk with the software upgrade package into the video endoscopy processor and click "Upgrade software".
- When you need to upgrade the firmware, insert the USB disk with the firmware upgrade package into the video endoscopy processor, and click "Upgrade firmware".


Fig4.26 Upgrade interface (VS50 & VS100)

4. In the "Other" page, Click the "About" button to display the software version, system version and operation section information. (As shown in the figure below)

Note

- The "operation section" on the interface refers to the endoscope used with the video endoscopy processor.
- The video endoscopy processor can display the model, serial number, production date and used times of the endoscope.


Fig4.27 About interface (VS50 & VS100)

5. In the "Other" page, Click the "Teaching" button to play the video about the use and operation of endoscope (as shown in the figure below)

Caution

- The teaching video is only used as the basic information of learning and understanding the use and operation of products, not as the diagnosis and treatment procedures.


Fig4.28 Teaching interface (VS50 & VS100)

6. In the "Other" page, Click the "Help" button to view the help content. (As shown in the figure below)


Fig4.29Help interface (VS50 & VS100)

4.6 Taking pictures

4.6.1. Press the record button  on the operating handle, the Photo icon at the top of the screen is shown in red, about a second after, the red photo icon is gone. This picture has been added on time and the file name in the built-in SD card successfully saved.

Fig 4.36

Note

- The instrument automatically generates the real-time time and file name for recording photos. File names are mainly recorded in numbers. Visible in playback mode.

4.7 Recording video

4.7.1. In the real time mode keep pressing video/camera button  for 3 seconds, enter into record mode, the video icon at the top of screen is shown in red  00: 00: 00 . If you need to exit the recording mode, press the video/camera button  for 3 seconds: you can switch to preview mode to display the image of the current endoscope. (See Fig4.37)

Fig4.37

4.8 White balance adjustment

Warning

- White balance adjustment is used to display the image color correctly in the LCD display. Without the correct white balance may cause color distortion.

1. Press the power button on the left side of the monitor, the LCD display light lights up, enter preview mode. Let the distal tip of the insertion section aim to white objects or white gauze. Keep pressing white balance  button on the extension cable or wireless transmitter for 3 seconds, until the icon  appear, white balance adjustment succeeds.
2. Repeat step 1, if adjust fault.

4.9 DataStorage

To avoid videos or images, if the film is deleted by mistake, U disk can be inserted into the USB interface of the endoscope camera system for backup Copies. You do not need to install any driver software to back up data using a USB flash drive.

4.10 power off

When the equipment is not in use, gently press the power button on the left side of the display:  Then disconnect the power adapter and power plug.

The device should be placed in a suitable position that the operator can disconnect it.

Chapter 5 Storage and Disposal

Keep the clean video endoscope processor and accessories after cleaning and sterilization away from the contaminated equipment. If the clean video endoscope processor and accessories are contaminated, their use may pose an infection control risk to the patient and/or operators performing the next procedure with the video endoscope processor.

Environment specifications

Normal working conditions	Environment temperature:	5℃ ~ +40℃
	Atmospheric pressure:	700hPa ~ 1060hPa
	Relative humidity:	30% ~ 85% (Non – Condensing)
Transportation and storage	Environment temperature:	-10℃ ~ +55℃
	Atmospheric pressure:	500hPa ~ 1060hPa
	Relative humidity:	≤93% (Non – Condensing)

Form 5.1 Environment specifications

5.1 Storage of the video endoscope processor

1. Remove all attachments from the video endoscope processor.
2. Verify that all surfaces of the electronic endoscopic image processor are completely dry.

Caution

- To prevent contamination clean video endoscope processor and accessories after cleaning are always kept in place.
- The places where the video endoscope processor is stored must be clean, dry, well ventilated and kept at room temperature. Never store it in an environment that has direct sunlight, high temperature, humidity, or is exposed to ozone, X-rays, or UV light, which can cause damage to the video endoscope processor.
- Never store the video endoscope processor inside the porter box. Porter boxes should only be used to transport the video endoscope processor.

5.2 Transporting Outside the Hospital

Warning

- Transport the video endoscope processor in the carrying case.
- Be sure to clean the surface before transporting the video endoscope

processor outside the hospital. Otherwise, the patient might be infected.

5.3 Disposal of expired or scrapped

When the video endoscope processor reaches the expiration date or can no longer be used, the disposal of the video endoscope processor and any components (such as power adapter, etc.) shall comply with all applicable national and local laws and regulations.

Caution

- The wireless transmitter contains a built-in battery. When scrapping, do not put it into fire or dispose with domestic waste, otherwise it may cause explosion or environmental pollution.
- The circuit board of this product contains trace heavy metals and other harmful substances. Do not discard it with domestic waste, otherwise it will cause environmental pollution.

1. If possible, please disassemble the wireless transmitter and remove the battery before handling it separately.
2. Please comply with local laws and regulations for the disposal of abandoned electrical and electronic equipment.

Chapter 6 Troubleshooting

If the video endoscope processor is visibly damaged, does not function as expected or is found to have irregularities during the inspection described in Chapter 3, "Preparation and Inspection", do not use it. Contact our company.

Some problems that appear to have irregularities may be correctable by referring to Section 6.1, "Troubleshooting guide". If the problem cannot be resolved by the described remedial action, stop using the video endoscope processor and send it to us for repair.

Our company does not repair accessory parts. If an accessory part becomes damaged, contact us to purchase a new one.

Warning

- Never use the video endoscope processor on a patient if an abnormality is suspected. Damage or irregularity in the instrument may compromise patient or user safety and may result in more severe equipment damage.

6.1 Troubleshooting Guide

Troubles or failures due to other causes than those listed below should be serviced. As repair performed by persons who are not qualified by our company could cause patient or user injury and/or equipment damage, be sure to contact us for repair according to Section 6.2, "Returning the endoscope for repair".

Image Quality

Irregularity description	Possible cause	Solution
There is no video image.	Not all power switches are ON.	Turn ON all the power switches.
	Improper connection between video endoscope processor and endoscope operation unit	Correctly connect the endoscope operation unit and the video endoscope processor according to the requirements of the manual.
An image is not clear.	The objective lens is dirty.	Wipe with cotton swab moistened with 70% ethyl or isopropyl alcohol. If the image cannot be observed, return the video endoscope processor and the supporting endoscope operation to our company for repair.

Form 6.1 Image Quality

Key

Irregularity description	Possible cause	Solution
Cannot Soft Shutdown	Software system crash	Press the left power switch button to restart
No response when touching virtual key	Touch screen failure	

Form 6.2 Key

Photograph

Irregularity description	Possible cause	Solution
No picture generated	Internal memory or SD card is full.	Delete unnecessary folders or files. Download the important pictures to PC before delete.
The image data are not recorded in accordance with the date.	The date and time are not set.	Set the date and time.
	Low battery power with built-in clock.	If it is still wrong with the date and time after re-setting. Return the endoscope to our company for repair.

Form 6.3 Photograph

6.2 Returning the video endoscope processor for Repair

Warning

- Before repairing the video endoscope processor, please thoroughly clean its surface with 70% alcohol. Improper cleaning will bring infection risk to the hospital or the company's personnel handling the equipment.

Caution

- Our company is not responsible for the equipment damage and/or person injury which are not caused by our company personnel trying to repair the equipment.

Before returning the instrument for repair, contact our company. With the instrument, include a description of the malfunction or damage and the name and telephone number of the individual at your location who is most familiar with the problem. Also include the warranty card. When returning the instrument for repair, follow the instructions given in Section 5.2 "Transporting outside the hospital".

Chapter 7 Cleaning and Disinfection

7.1 Maintain

The expected service life of the product is 5 years, the entire product life period, the product can work normally, but does not exclude the use of the process or the damage in the transportation process, if the power line is found to be squeezed, the plug is deformed or the appearance of the adapter or the appearance of the processor is damaged, you need to contact the manufacturer for maintenance confirmation.

Caution

- Do not provide service or maintain while it is being used by the patient, which may harm the patient and the operator.
- Disconnect the video endoscope processor from the endoscope, remove any accessories, and completely turn off the video endoscope processor before cleaning and disinfection.
- Do not let water or any liquid drip into the video endoscope processor. When using cleaning or disinfectant, follow the manufacturer's specifications for exposure time and product life.
- After use, the video endoscope processor must be stored according to local guidelines until the next use.

7.2 Cleaning

1. Recommended cleaning agents

Cleaning agent	Metrex Empower
Manufacturer	Metrex

Form 7.1

2. Manual cleaning step at room temperature

- 1) Turn off the electronic endoscopic image processor.
- 2) Prepare cleaning solution using recommended enzyme detergent.
- 3) Permeabilize a piece of sterile gauze in the enzymatic cleaning solution and ensure that the gauze is wet but does not drip.
- 4) Clean the electronic endoscope image processor's buttons, screens, housing thoroughly with moistened gauze. Avoid wetting the unit to prevent damage to the internal electronic components.

- 5) Wait 10 min (or press the recommended time from the detergent manufacturer) to allow the enzyme to work.
- 6) Wipe the electronic endoscopic image processor using sterile gauze that has been moistened with Rodi water. Ensure that all traces of cleaning fluid have been removed.
- 7) Repeat steps 1-6 until the electronic endoscopic image processor surface is clean and visually inspected free of stains and water.

7.3 Disinfection

Caution

- It is ineffective to do disinfection before active cleaning, which easily leads to coagulation of proteins and organic residues on the equipment and prevents correct disinfection.

1. Disinfectant

Disinfectant	Concentration
Ethyl Alcohol	75%

Form 7.2

2. Manual disinfection procedure at room temperature

- 1) Prepare a piece of gauze dipped with 75% Ethyl Alcohol to ensure that the gauze is wet but not dripping, so as to avoid any liquid dripping into the video endoscope processor.
- 2) Wipe the buttons, screens, housings, slots and gaps on the video endoscope processor.
- 3) Place the product at room temperature for 10 minutes.
- 4) Wipe the surface of the video endoscope processor with dry gauze.

Chapter 8 Electromagnetic compatibility

Electromagnetic compatibility

Caution:

- Video endoscope processor meets the requirements of IEC 60601-1-2 standard for electromagnetic compatibility.
- The user shall install and use the electromagnetic compatibility information provided by the random files.
- Portable and mobile RF communications devices may affect the performance of the Video endoscope processor and avoid strong electromagnetic interference such as near mobile phones and microwave ovens.
- The guidelines and manufacturer's declarations are detailed in the annex.

Warning:

- A device or system should not be used with other devices placed close to or over stacked, if must be close to or stacked, should observe and verify the normal operation in its use of the configuration.
- Type A equipment for the industrial environment, since the product conduction harassment and radiation harassment, in other environments to ensure electromagnetic compatibility may have potential difficulties.
- In addition to the cables sold as spare parts of internal components by the manufacturer of Video endoscope processor, the use of accessories and cables outside the regulations may increase the emission or reduce the immunity.
- Even if other devices meet the emission requirements of corresponding national standards, Video endoscope processor may still be disturbed by other devices.

Type	Launch frequency	modulation type	Radiated Power
WIFI (for VS100 only)	2.4G/5.8G	ISM/UNII	< = 2dB

Sample cable (for VS100 only)

No.	Item	Length of cable	Shielding(Y/N)	Remark
1	Power cable	1.5	N	/
2	Adapter connecting cable	1.4	N	/
3	Extension Cable	1.5	Y	(EC182)
4	Extension Cable	2.2	Y	(EC181)

Annex:

Guide and manufacturer's statement - Electromagnetic Emission			
Video Endoscope Processor used in the electromagnetic environment of the expected performance of the following provisions, the performance of Video Endoscope Processor purchasers or users should ensure that it is used in the electromagnetic environment:			
Test item	EMC standard or test method	Immunity IEC 60601 test level	Compliance level
Conducted Emissions	CISPR 11	Group 1 Class A	Group 1 Class A
Radiated RF Emissions	CISPR 11	Group 1 Class A	Group 1 Class A
Harmonic distortion	IEC 61000-3-2	Class A	N/A
Voltage fluctuations and flicker	IEC 61000-3-3	$d_{max} < 4\%$, $P_{st} \leq 1$, $d_c < 3.3\%$	N/A

Guide and manufacturer's statement - Electromagnetic Immunity			
Video Endoscope Processor used in the electromagnetic environment of the expected performance of the following provisions, the performance of Video Endoscope Processor purchasers or users should ensure that it is used in the electromagnetic environment:			
Test item	EMC standard or test method	Immunity IEC 60601 test level	Compliance level
Electrostatic discharge (ESD)	IEC 61000-4-2	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact, ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields	IEC 61000-4-3	3V/m rms 80MHz-2.7GHz 80% AM at 1 kHz	3V/m rms 80MHz-2.7GHz 80% AM at 1 kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	9V/m, 217 Hz LTE Band 13,17 for 710-780MHz; WLAN 802.11 a/n for 5240-5785 MHz; 27V/m, 18Hz for 385MHz; 28V/m, 18Hz for 810-930MHz, ± 5 kHz deviation 1 kHz sine for 450MHz, 217Hz for 1720-2450.	9V/m, 217 Hz LTE Band 13,17 for 710-780MHz; WLAN 802.11 a/n for 5240-5785 MHz; 27V/m, 18Hz for 385MHz; 28V/m, 18Hz for 810-930MHz, ± 5 kHz deviation 1 kHz sine for 450MHz, 217Hz for 1720-2450.
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50Hz or 60Hz	30 A/m 50Hz or 60Hz
Proximity magnetic fields	IEC 61000-4-39	8A/m, CW for 30kHz. 65A/m, 2,1kHz Pulse modulation for 134,2kHz. 7,5A/m, 50kHz Pulse modulation for 13,56MHz	8A/m, CW for 30kHz. 65A/m, 2,1kHz Pulse modulation for 134,2kHz. 7,5A/m, 50kHz Pulse modulation for 13,56MHz
Electrical fast transients / bursts	IEC 61000-4-4	± 2 KV (Input AC power port) ± 1 kV (SIP/SOP port) 100 kHz repetition frequency	± 2 KV (Input AC power port) ± 1 kV (SIP/SOP port) 100 kHz repetition frequency
Surges	IEC 61000-4-5	Line-to-line: $\pm 0,5$ kV, ± 1 kV Line-to-ground: $\pm 0,5$ kV, ± 1 kV, ± 2 kV (Input AC power port) Line-to-ground: ± 2 kV (SIP/SOP port)	Line-to-line: $\pm 0,5$ kV, ± 1 kV Line-to-ground: $\pm 0,5$ kV, ± 1 kV, ± 2 kV (Input AC power port) Line-to-ground: ± 2 kV (SIP/SOP port)
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V rms 150 kHz- 80 MHz 6 V rms in ISM and amateur radio bands between 0,15 MHz and 80 MHz	3 V rms 150 kHz- 80 MHz 6 V rms in ISM and amateur radio bands between 0,15 MHz and 80 MHz
Voltage dips	IEC 61000-4-11	0 % UT; 0,5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles, Single phase: at 0°	0 % UT; 0,5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles, Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0 % UT; 250/300 cycle	0 % UT; 250/300 cycle
Note: UT refers to an alternating voltage test voltage before			

Guide and manufacturer's statement - Electromagnetic Immunity			
Video Endoscope Processor used in the electromagnetic environment of the expected performance of the following provisions, the performance of Video Endoscope Processor purchasers or users should ensure that it is used in the electromagnetic environment:			
Test item	EMC standard or test method	Immunity IEC 60601 test level	Compliance level
RF transmission	IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz
RF radiation	IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz
Note 1: in 80MHz and 800MHz frequency, adopt the high frequency formula. Note 2: these guidelines may not be appropriate in all cases, the electromagnetic propagation may be affected by buildings, absorption and reflection of objects and people.			
a Fixed transmitter electric field strength, such as base station of wireless (cellular / cordless) telephone and ground mobile radio, amateur radio, AM (AM) and FM (FM) radio and television broadcasting etc.their emission electric field strength is impossible to be predicted in theory. In order to evaluate the electromagnetic environment of a fixed RF transmitter, the investigation of electromagnetic field should be conducted. If the field strength of the video endoscope processor is measured higher than the RF compliance voltage of the above-mentioned application, the performance of the video endoscope processor should be observed to verify its normal operation. If the observed performance is abnormal, then additional measures may be required, such as the new series of video endoscope processor performance orientation or location.			
b In the whole frequency range of 150KHz ~ 80MHz, the emission electric field strength should be less than 3 V/m.			

Recommended isolation distance between portable and mobile communication equipment RF and Video Endoscope Processor performance

The performance of the video endoscope processor is expected to be used in the controlled electromagnetic environment of radiated RF harassment. On the basis of the maximum output power of communication equipment, the purchasers or users can maintain portable and mobile RF communication equipment through the following recommended (transmitter) and video endoscope processor performance minimum distance to prevent electromagnetic interference.

Rated maximum output power of the transmitter/W	Corresponding to different frequency transmitter isolation distance/m		
	150 kHz ~ 80 MHz $d = 1.2\sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2\sqrt{P}$	800 MHz ~ 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

The transmitter rated maximum output power are not listed in the form, the recommended isolation distance D, in meters (m) as a unit, can use the corresponding transmitter frequency bar formula to determine, where P is the maximum rated output power provided by the transmitter of transmitter manufacturers, by Watt (W) as a unit.

Note 1: in 80MHz and 800MHz frequency, adopt the high frequency formula.

Note 2: these guidelines may not be appropriate in all cases, the electromagnetic propagation may be affected by buildings, absorption and reflection of objects and people.

Manufacturers: Zhuhai Vision Medical Technology Co., Ltd.

**Address: 301 Building 4, No. 30 Feipeng Road, GMTCM Park, Hengqin New Area,
519031 Zhuhai City, Guangdong Province, PEOPLE' S REPUBLIC OF CHINA**

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




Наименование	Изображение
Видеопроцессор эндоскопический VS100	
Видеопроцессор эндоскопический VS100	 The image shows the VS100 endoscopic video processor, which consists of a monitor and a video processor unit. The monitor is a small, square, black screen with a silver frame, displaying a dark image. Below the monitor is the video processor unit, which is a long, cylindrical, silver-colored device with two large, hexagonal, silver-colored connectors at each end. The brand name 'MINDSIGN' is visible on the monitor's bezel, and 'VS100' is printed below the processor unit.
Кабель питания	 The image shows a black power cable with a standard three-pronged AC power plug on one end and a multi-pin connector on the other. The cable is coiled and lies on a light-colored surface.
Адаптер питания	 The image shows a black power adapter with a multi-pin connector on one end and a USB-A connector on the other. The adapter has a label with technical specifications and safety warnings. It is connected to a black cable that lies on a light-colored surface.

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