Approved by

Vladislav Zaitseva

[Signature] Director of RA

"04" June 2021

[Signature] 04, June 2021

INSTRUCTION FOR USE
Applicators for BioGlue Surgical Adhesive
CryoLife, Inc., USA

Version 04-2021

2021

Index

1. Product Name.....	3
2. Intended use:	3
3. Special properties of the product	3
4. Technical specification.....	3
4.1. Product Variants:	3
4.2. Materials used in the product. Type of contact with the body:.....	4
4.3. Product Image	5
4.4. Design and dimensions of the product	6
5. Labelling and packaging.....	11
6. Compliance with the standards of the Russian Federation.	18
7. General indications, contraindications, potential complications, precautions	19
7.1. Indications	19
7.2. Contraindications	19
7.3. Potential complications	19
7.4. Precautions	19
9. Transportation, storage and operation.....	22
10. Information on the sterility of the product	23
11. Disposal information.....	23
12. Warranty.....	23

This instruction for use has been prepared for the purpose of registration of the medical device in The Russian Federation. Information in the file corresponds to the requirements of the Resolution of the Government of the Russian Federation dated December 27, 2012 No. 1416 "Rules for Registration of Medical Products".

1. Product Name

"Applicator for BioGlue Surgical Adhesive." (Hereinafter: "product", "applicator", "tip").

Manufacturer:	CryoLife, Inc., USA
Manufacturer's Address:	1655 Roberts Blvd., NW, Kennesaw, Georgia, 30144 USA
Telephone:	+1 (770) 419 33 55
E-mail:	usbioglue@cryolife.com
Place of production:	CryoLife, Inc., USA 1655 Roberts Blvd., NW, Kennesaw, Georgia, 30144 USA

Type of contact with the body: short-term contact with the contact with soft tissue.

Method of application: in accordance with instructions for use.

Conditions of use: intended for use in hospitals and clinics, is to be used only by qualified medical personnel.

2. Intended use:

The product is intended for mixing two components of BioGlue Surgical Adhesive and applying to the surgical field.

3. Special properties of the product

These applicators are intended only for use with BioGlue Surgical Adhesive, which is registered on the territory of the Russian Federation (RU No. FCZ 2010/07172 dated June 09, 2010). The range of applicators * and extensions ** is intended for use with a standard applicator. They are worn on standard applicators for applying surgical glue according to the instructions for using BioGlue Surgical Adhesive.

* in this case and hereinafter, the assortment of applicators refers to all options for the performance of applicators, except for the standard applicator.

** common name for:

- *Delivery Tip Extension 10 cm;*
- *Delivery Tip Extension 27 cm;*
- *Delivery Tip Extension 35 cm;*

4. Technical specification

4.1. Product Variants:

Applicator for BioGlue Surgical Adhesive.

I. In the following variants:

1. Syringe Standard Applicator Tip;
2. Syringe Spreader Tip 12 mm; *
3. Syringe Spreader Tip 16 mm; *
4. Syringe Extender Tip 10 cm; **

5. Syringe Extender Tip 27 cm; **

II. Accessories

- Delivery Tip Extension 10 cm; ***

- Delivery Tip Extension 27 cm; ***

- Delivery Tip Extension 35 cm; ***

* hereinafter: "Spreader Tip", "applicators with an Spreader Tip"

** hereinafter: "Extender Tip"

*** hereinafter: "extenders", "Tip Extension"

4.2. Materials used in the product. Type of contact with the body:

Table 4.1 - Materials used in the product

Accessory Part Name	Material Name/ Model /Grade/ Number/ Colourant Name	Material Manufacturer / Vendor	Patient Contacting	BioGlue Surgical Adhesive contacting
Syringe Standard Applicator Tip				
Tip Housing	Polypropylene copolymer ExxonMobil, PP9074	EXXONMOBIL CHEMICAL COMPANY	Short-term soft tissue contact	Yes
Mixing Element	Polyester Vectra MT1345 VF3001 (natural)	Celanese Co.	Not in contact	Yes
Retaining Ring	Polypropylene w/ Green Dye Bormed RF830MO Colourant Name: MB 11962 PP Green	Borealis	Short-term soft tissue contact	No
Syringe Extender Tip 10 cm / Syringe Extender Tip 27 cm				
Tubing (External)	Polyamide Rilsan Besno A61297	Arkema Inc.	Short-term soft tissue contact	Yes
Tubing (Internal)	Polyurethane Thermedics TT-1055D / Pellethane 2363-55DE NAT	THE LUBRIZOL CORPORATION	Not in contact	Yes
Connector	Acrylonitrile-butadiene styrene (ABS) Lustran ABS 633 WT013003 Colourant Name: White 013003	INEOS Styrolution America LLC	Short-term contact with soft tissue is possible	Yes
Syringe Spreader Tip 12mm / Syringe Spreader Tip 16mm				
Tip	Polypropylene Bormed RF830MO	Borealis	Short-term soft tissue contact	Yes
Tip Housing	Polypropylene copolymer Bormed RF830MO	Borealis	Short-term soft tissue contact	Yes
Mixing Element	Polypropylene w/ Green Dye Vectra MT1345 VF3001 (natural)	Celanese Co.	Not in contact	Yes
Retaining Ring	Polypropylene Bormed RF830MO Краситель: MB 11962 PP Зеленый	Borealis	Short-term contact with soft tissue is possible	No
Spreader tip of 12 mm	Polycarbonate (PC) Lexan HPS1 -1125	Sabic Innovative Plastics	Short-term soft tissue contact	Yes
Spreader tip of 16 mm	Polycarbonate (PC) Lexan HPS1 -1125	Sabic Innovative Plastics	Short-term soft tissue contact	Yes
Delivery Tip Extension (10 cm, 27 cm, 35 cm)				
Proximal Flare	Soft polyether block amide Pebax 4033 SA 01 MED	Arkema Inc.	Short-term contact with soft tissue is possible	Yes

Extrusion Tube	Polyurethane Isoplast 2530 NAT	THE LUBRIZOL CORPORATION	Short-term soft tissue contact	Yes
Distal Tip	Soft polyether block amide Pebax 4033 SA 01 MED	Arkema Inc.	Short-term soft tissue contact	Yes
Packages				
Primary packaging	Tyvek 1073B; Poly, TPF-0501A	Oliver-Tolas	Short-term soft tissue contact	No

4.3. Product Image

Fig. 4.1 - Syringe Standard Applicator Tip


Fig. 4.2. - Syringe Extender Tip 10 cm

Fig. 4.2. - Syringe Extender Tip 27 cm

Fig. 4.3. - Syringe Spreader Tip 12 mm


Fig. 4.4. - Syringe Spreader Tip 16 mm


Fig. 4.5.- Delivery Tip Extension 10 cm



Fig. 4.6.- Delivery Tip Extension 27 cm*

* the appearance of the Delivery Tip Extension 27 cm and the Delivery Tip Extension 35 cm is completely identical except for the length.


4.4. Design and dimensions of the product

Drawing 4.1.-Drawing of Syringe Spreader Tip 12 mm


Size A (mm)	Size B (mm)	Size C (mm)	Size D (mm)	Size E (mm)	Size F (mm)	Size G (mm)	Size H (mm)	Size I (mm)	Size J (mm)	Size K (mm)
$\varnothing 1.6 \pm 0.1$	$\varnothing 2.25 \pm 0.1$	$\varnothing 2.55 \pm 0.1$	10.5 ± 0.5	3.8 ± 0.1	0.2 ± 0.05	1.5 ± 0.1	6.6 ± 0.1	18.3 ± 0.5	12 ± 0.5	4.7 ± 0.1

Drawing 4.2.- Drawing of Syringe Spreader Tip of 16 mm


Size A (mm)	Size B (mm)	Size C (mm)	Size D (mm)	Size E (mm)	Size F (mm)	Size G (mm)	Size H (mm)	Size I (mm)	Size J (mm)	Size K (mm)
$\varnothing 1.6 \pm 0.1$	$\varnothing 2.25 \pm 0.1$	$\varnothing 2.55 \pm 0.1$	10.5 ± 0.5	3.8 ± 0.1	0.2 ± 0.05	1.5 ± 0.1	6.6 ± 0.1	18.3 ± 0.5	16 ± 0.5	4.7 ± 0.1

*Drawing 4.3.-Drawing of an applicator tip with an extended tip of 12 mm /
applicator with an extended tip of 16 mm*


Size A (mm)	Size B (mm)	Size C (mm)	Size D (mm)	Size E (mm)	Size F (mm)	Size G (mm)	Size H (mm)	Size I (mm)	Size J (mm)	Size K (mm)
$\varnothing 13.6 \pm 0.1$	5.2 ± 0.1	$\varnothing 4 \pm 0.1$	$\varnothing 2.25 \pm 0.1$	8.5 ± 0.1	67.9 ± 0.1	$\varnothing 11.8 \pm 0.1$	$\varnothing 5.05 \pm 0.1$	$\varnothing 3.8 \pm 0.07$	$\varnothing 2.6 \pm 0.07$	0.7 ± 0.05

Drawing 4.4.-Drawing of Syringe Standard Applicator Tip


Size A (mm)	Size B (mm)	Size C (mm)	Size D (mm)	Size E (mm)	Size F (mm)	Size G (mm)	Size H (mm)	Size I (mm)
$\varnothing 1 \pm 0.1$	74.25 ± 0.7	$\varnothing 3.25 \pm 0.1$	$\varnothing 9.8 \pm 0.5$	13.6 ± 0.5	$\varnothing 2.6 \pm 0.1$	$\varnothing 3.8 \pm 0.1$	$\varnothing 4.5 \pm 0.1$	$\varnothing 1.5 \pm 0.1$

Drawing 4.5.-Drawing of Syringe Extender Tip 10 cm / Syringe Extender Tip 27 cm

Size A (mm)	Size B (mm)	Size C (mm)	Size D (mm)	Size E (mm)	Size F (mm)
4.8 ± 0.1	$\varnothing 1.02 \pm 0.08$	$\varnothing 2.39 \pm 0.08$	$\varnothing 5.2 \pm 0.1$	$\varnothing 2.5 \pm 0.1$	100 ± 0.1

Drawing. 4.6.-Drawing of Delivery Tip Extension 10 cm / Delivery Tip Extension 27 cm / Delivery Tip Extension 35 cm

Size A (mm)	Size B (mm)	Size C (mm)	Size D (mm)	Size E (mm)
$\varnothing 1.5 \pm 0.1$	100 ± 0.5	270 ± 0.5	350 ± 0.5	$\varnothing 4.5 \pm 0.1$

4.5. Primary packaging characteristics

The number of products in the primary packaging, the dimensions of the primary packaging, net weight and gross weight of the primary packaging:

- Syringe Standard Applicator Tip:

Dimensions of the primary packaging: 30.50 x 16.50 cm $\pm$ 5%

Number of products: 4 pcs.

Net weight: 6g $\pm$ 0.5g

Gross weight: 16.5g $\pm$ 1g

- Syringe Spreader Tip 12 mm:

Dimensions of the primary packaging: 30.50 x 16.50 cm $\pm$ 5%

Number of products: 3 pcs.

Net weight: 5.4g $\pm$ 0.5g

Gross weight: 16.1g $\pm$ 1g

- Syringe Spreader Tip 16 mm:

Dimensions of the primary packaging: 30.50 x 16.50 cm $\pm$ 5%

Number of products: 3 pcs.

Net weight: 5.7g $\pm$ 0.5g

Gross weight: 16.4g $\pm$ 1g

- Syringe Extender Tip 10 cm:

Dimensions of the primary packaging: 60.96x17.78cm $\pm$ 5%

Number of products: 4 pcs.

Net weight: 11.6g $\pm$ 1g

Gross weight: 34.9g $\pm$ 1 g

- Syringe Extender Tip 27 cm:

Dimensions of the primary packaging: 60.96x17.78cm $\pm$ 5%

Number of products: 4 pcs.

Net weight: 14g $\pm$ 1g

Gross weight: 40.9g $\pm$ 1g

- Delivery Tip Extension 10 cm:

Dimensions of the primary packaging: 32.08 x 11.43 cm $\pm$ 5%

Number of products: 2 pcs.

Net weight: 7g $\pm$ 0.5g

Gross weight: 14.3g $\pm$ 1g

- Delivery Tip Extension 27 cm:

Dimensions of the primary packaging: 56.8325x11.43cm $\pm$ 5%

Number of products: 2 pcs.

Net weight: 14.6g $\pm$ 1g

Gross weight: 23.2g $\pm$ 1g

- Delivery Tip Extension 35 cm:

Dimensions of the primary packaging: 56.8325x11.43cm $\pm$ 5%

Number of products: 2 pcs.

Net weight: 18 g ± 1 g

Gross weight: 28.2g ± 1g

Characteristics of primary packaging materials Tyvek 1073B; Poly TPF-0501A:

Parameter	Material Tyvek 1073B	Material Poly TPF-0501A
Tensile strength (sealing strength)	≥ 87.56 N / m for the left and right sides of the individual packaging. From 87.56 to 700.51 N / m for the seam.	
Density	74,7 g/m ²	63.5 g/m ²
Tensile strength	205 H (longitudinal) 219 H (transverse)	544,69 kPa (longitudinal) 43 kPa (transverse)
Thickness deviation	199 mkm	63,5 mkm
Breathability	540 ml/min	69,8 (cm ³ /m ²)/day
Punching strength	10324 Dj/m ²	----

5. Labelling and packaging


Standard and extended applicators are packaged in primary packaging of 4 pieces. Applicators with an expanded tip of 12 mm / applicators with an expanded tip of 16 mm are packed in primary packages of 3 pieces. Applicator extenders 10 cm / applicator extenders 27 cm / Applicator extenders 35 cm are packed in primary packages of 2 pieces.

The labeling of the primary packaging contains the following information in English:

-Product name;

- Information about the version of the product;

-Symbol "Manufacturer" and information about the manufacturer ;

-Symbol "Authorized Representative in the European Union" ;

-Symbol "No reuse" ;

-Symbol "Do not use if packaging is damaged" ;

-Symbol "Caution! Refer to the instructions for use" ;

-Symbol "Refer to the instructions for use" ;


- Information on non-toxicity;

-Symbol "Mark of conformity" ;

-Image of the product;



-Number of products in a package;

- The symbol "Radiation sterilization"  or "Sterilization with ethylene oxide" 
(depending on the model of the product);

-Symbol "The purchase of this product is available exclusively for doctors and other qualified medical personnel and in accordance with their requirements"  ;

-Symbol "Does not contain latex"  ;

A sticker with information in Russian containing the following information is also glued to this individual package.

- Product name;
- Information about the version of the product;
- Information about AR;
- Information about RU;
- REF;
- LOT;
- Date of manufacture  ;
- Expiration date  ;
- Symbol "Manufacturer" and information about the manufacturer  ;
- Symbol "No reuse"  ;
- Symbol "Do not use if packaging is damaged"  ;
- Symbol "Caution! Refer to the instructions for use"  ;
- Symbol "Refer to the instructions for use"  ;
- Information on non-toxicity;
- The symbol "Radiation sterilization"  or "Sterilization with ethylene oxide"  (depending on the model of the product);
- Symbol "The purchase of this product is available exclusively for doctors and other qualified medical personnel and in accordance with their requirements"  ;
- Symbol "Does not contain latex"  ;

* REFs vary depending on the variants:

- BGAT-10-SY for Syringe Extender Tip 10 cm;
- BGAT-27-SY for Syringe Extender Tip 27 cm;
- BGAT-SY for Syringe Standard Applicator Tip;
- BGDTE-10 for Delivery Tip Extension 10 cm;
- BGDTE-27 for Delivery Tip Extension 27 cm;
- BGDTE-35 for Delivery Tip Extension 35 cm;
- BGST-12 for Syringe Spreader Tip 12 mm;
- BGST-16 for Syringe Spreader Tip 16 mm;

Individual packaging in the amount of 10 pcs. packed in consumer packaging.
Dimensions and mass of consumer packaging:

- Standard Applicator Tip:

Dimensions: Length: 16.5 cm $\pm$ 0.5 cm Width: 7.9375 cm $\pm$ 0.5 cm Height: 25.4 cm $\pm$ 0.5 cm

Gross weight: 294.8g $\pm$ 20g

- The applicator with an spreader tip 12 mm:

Dimensions: Length: 16.5 cm $\pm$ 0.5 cm Width: 7.9375 cm $\pm$ 0.5 cm Height: 25.4 cm $\pm$ 0.5 cm

Gross weight: 294.8g $\pm$ 20g

- The applicator with an spreader tip 16 mm:

Dimensions: Length: 16.5 cm $\pm$ 0.5 cm Width: 7.9375 cm $\pm$ 0.5 cm Height: 25.4 cm $\pm$ 0.5 cm

Gross weight: 295.1g $\pm$ 20g

- The applicator, extended tip by 10 cm:

Dimensions: Length: 62.23 cm $\pm$ 0.5 cm Width: 18.415 cm $\pm$ 0.5 cm Height: 10.16 cm $\pm$ 0.5 cm

Gross weight: 725.7 g $\pm$ 20 g

- The applicator, extended tip by 27 cm:

Dimensions: Length: 62.23 cm $\pm$ 0.5 cm Width: 18.415 cm $\pm$ 0.5 cm Height: 10.16 cm $\pm$ 0.5 cm

Gross weight: 739.4 g $\pm$ 20 g

- Delivery Tip Extension of applicator 10 cm:

Dimensions: Length: 32.0675 cm $\pm$ 0.5 cm Width: 6.35 cm $\pm$ 0.5 cm Height: 11.43 cm $\pm$ 0.5 cm

Gross weight: 362.9 g $\pm$ 20g

- Delivery Tip Extension of applicator 27 cm:

Dimensions: Length: 56.8325 cm $\pm$ 0.5 cm Width: 6.35 cm $\pm$ 0.5 cm Height: 11.43 cm $\pm$ 0.5 cm

Gross weight: 589.7 g $\pm$ 20 g

- Delivery Tip Extension of applicator 35 cm:

Dimensions: Length: 56.8325 cm $\pm$ 0.5 cm Width: 6.35 cm $\pm$ 0.5 cm Height: 11.43 cm $\pm$ 0.5 cm

Gross weight: 612.3 g $\pm$ 20g





Рис 5.3. Consumer package of the product

Consumer packaging design contains:

- Name of the product in various foreign languages and a brief information about it;
- The name of the manufacturer in English, as well as information about it;
- Date of expiration;
- Symbol "Prohibition of reuse" ;

- Symbol "Refer to the instructions for use" ;
- Symbol "Mark of conformity" ;
- Image of products;
- Equipment;
- The symbol "Radiation sterilization"  or "Sterilization with ethylene oxide" 
(depending on the model of the product);
- Lot number (LOT);
- Product Code (REF);
- A sticker is also glued on consumer packaging with information in Russian, containing the following information:
- Product name;
- Variant of execution;
- Information about AR;
- Information about RU;
- REF; *
- LOT;
- Date of production ;
- Date of expiration ;
- Symbol "Manufacturer" and information about the manufacturer ;
- Symbol "Prohibition of reuse" ;
- Symbol "Do not use if packaging is damaged" ;
- Symbol "Caution! Refer to the instructions for use" ;
- Symbol "Refer to the instructions for use" ;
- Information on non-toxicity;
- The symbol "Radiation sterilization"  or "Sterilization with ethylene oxide" 
(depending on the model of the product);
- Symbol "The purchase of this product is available exclusively to doctors and other qualified medical personnel and in accordance with their requirements" ;
- Symbol "Not made with natural rubber latex" ;

Consumer packaging is packaged in one shipping package.
 Sizes of transport packaging:

Shipping packaging	Information on the number of products in transport packaging
290,51x239,71x311,2(±10 MM) 	3 pcs. consumer packaging: either a Standard Applicator Tip , or an applicator with an spreader tip 16 mm, or an applicator with an spreader tip 12 mm.
425,4x174,6x2635,5(±10 MM) 	5 pcs. consumer packs of either a Standard Applicator Tip , or an applicator with an spreader tip 16 mm, or an applicator with an spreader tip 12 mm.
889x269,9x177,8 (±10 MM) 	8 pcs. consumer packs of either a Standard Applicator Tip , or an applicator with an spreader tip 16 mm, or an applicator with an spreader tip 12 mm.

425,4x425,4x419,1(±10 mm)


16pcs consumer packs of either a Standard Applicator Tip , or an applicator with an spreader tip 16 mm, or an applicator with an spreader tip 12 mm.

568,3x368,3x477,8(±10 mm)

20 pcs. consumer packs of either a Standard Applicator Tip , or an applicator with an spreader tip 16 mm, or an applicator with an spreader tip 12 mm.

685,8x571,5x635 mm (±10 mm)

46pcs consumer packs of either a Standard Applicator Tip , or an applicator with an spreader tip 16 mm, or an applicator with an spreader tip 12 mm.

Or 15 pcs. consumer packs of the applicator, extended tip by 10 cm; or the applicator, extended tip by 27 cm

813,2x203,2x203,2 MM (±10 MM)

2 pcs. consumer packs of the applicator,
extended tip by 10 cm; or the applicator,
extended tip by 27 cm.

Labelling of transport packaging contains:

- Product name;
- Variant of execution;
- Lot number (LOT);
- Number of products in transport packaging;
- Information about the manufacturer (name and address) ;
- Date of expiration ;
- The symbol "Fragile, handle with care"  ;
- Symbol "Avoid exposure to sunlight"  ;
- Symbol "Keep away from moisture"  ;
- Catalog number (REF);
- Symbol "Refer to the instructions for use"  ;
- Information on SCP;
- Information about the registration certificate;

6. Compliance with the standards of the Russian Federation.

GOST R 50444-92 "Medical devices, apparatus and equipment. General specifications. "

GOST ISO 11737-2-2011. "Sterilization of medical devices. Microbiological methods. Part 2. Sterility tests carried out during validation of sterilization processes "

GOST ISO 11607-2011. "Packaging for medical devices subject to final sterilization."

GOST 31814-2012. "Conformity assessment. General rules for sampling for product testing in conformity confirmation."

GOST 31214-2003. "Medical products. Requirements for samples and documentation provided for toxicological, sanitary and chemical tests, sterility and pyrogenicity tests."

GOST R 52770-2007. "Medical products. Safety requirements. Methods of sanitary-chemical and toxicological tests."

GOST ISO 10993-1-2011. "Medical products. Assessment of the biological effects of medical devices. Part 1. Assessment and research."

GOST R ISO 10993-2-2009. "Medical products. Assessment of the biological effects of medical devices. Part 2. Requirements for handling animals."

GOST ISO 10993-3-2011 "Medical devices. Assessment of the biological effects of medical devices. Part 3. Studies of genotoxicity, carcinogenicity and toxic effects on reproductive function."

GOST ISO 10993-4-2011. "Medical products. Assessment of the biological effects of medical devices. Part 4. The study of products that interact with blood."

GOST ISO 10993-5-2011. "Medical products. Assessment of the biological effects of medical devices. Part 5. Cytotoxicity test: in vitro methods."

GOST ISO 10993-6-2011. "Medical products. Assessment of the biological effects of medical devices. Part 6. Investigation of local action after implantation."

GOST ISO 10993-9-2015 "Medical devices. Assessment of the biological effects of medical devices. Part 9. Basic principles of identification and quantification of potential destruction products."

GOST ISO 10993-10-2011. "Medical products. Assessment of the biological effects of medical devices. Part 10. Studies of irritating and sensitizing effects."

GOST ISO 10993-12-2015. "Medical products. Assessment of the biological effects of medical devices. Part 12. Sample preparation and control samples."

GF RF edition XIII. OFS.1.2.4.0003.15. Sterility.

GF RF edition XIII. OFS.1.2.4.0005.15. Pyrogenicity.

7. General indications, contraindications, potential complications, precautions

7.1. Indications

The product is intended for mixing of two components of BioGlue Surgical Adhesive and applying to the surgical field.

7.2. Contraindications

- Not for intravascular use.
- Do not use in patients with known hypersensitivity to materials.

7.3. Potential complications

- Allergic reactions are possible.

7.4. Precautions

- The assortment of applicators and extension cords is intended for use in one patient only. Do not re-sterilize.
- Do not use if packaging has been opened or damaged.
- BioGlue Surgical Adhesive quickly polymerizes. The tip should be filled quickly, and then immediately apply BioGlue Surgical Adhesive to the fabric. The delay in applying glue to the fabric after filling the tip can lead to its polymerization in the tip.

8. Directions for use

(applicable in conjunction with BioGlue Surgical Adhesive, which is registered on the territory of the Russian Federation (RU No. FSZ 2010/07172 dated June 09, 2010))

The BioGlue Surgical Adhesive delivery system is a joint use of the BioGlue Surgical Adhesive applicator and BioGlue Surgical Adhesive. The syringe and piston are part of the BioGlue Surgical Adhesive, which is registered on the territory of the Russian Federation (RU No. ФСЗ 2010/07172 dated June 09, 2010)

8.1. Preparation

NOTE. Read the information in the instructions for use of BioGlue Surgical Adhesive, which is registered in the Russian Federation (RU No. ФСЗ 2010/07172 dated June 09, 2010)

1. Remove the syringe, piston, and applicators from the packaging. Take the syringe with BioGlue Surgical Adhesive in your hand, lifting it with the tip up, and tap it several times so that the air bubbles in the BioGlue Surgical Adhesive collect at the top of the syringe.

Fig. 8.1

NOTE: Hold the syringe in this position until all components of the BioGlue Surgical Adhesive delivery system are fully assembled so that air bubbles remain at the top of the syringe.

2. Remove the applicator from its packaging. Inspect the circlip: the pointer should be directly above the large hole. Otherwise, turn the ring on the tip body until the pointer is positioned opposite the large hole.

NOTE: Continue to hold the applicator upright while assembling the BioGlue Surgical Adhesive Delivery System and removing air from the applicator.

Fig. 8.2

3. Hold the syringe firmly with the tip up, rotate the cap 90 ° counterclockwise and, slightly swaying from side to side, remove it.

Align the applicator with the syringe along the notches on them and slide the applicator onto the syringe.


Fig. 8.3

ATTENTION: Try not to spill the solution from the syringe during assembly of the BioGlue Surgical Adhesive delivery system.

4. Lock the tip against the syringe and lock it by turning the tip retaining ring 90 ° clockwise.

Fig. 8.4

5. Continue to hold the syringe with the applicator up. Match the piston heads to the large and small openings of the syringe body, and insert the piston into the syringe until you feel resistance. The syringe is now ready for use.

Fig. 8.5

WARNING: Do not lay the syringe on its side until all air has been removed from it.

6. When using an elongated applicator, the tip can be given the desired angle by bending it in the right place at the right angle and holding it for 3 to 5 seconds. The created angle will last up to 5 minutes.

7. To remove the clogged tip, turn the tip retaining ring counterclockwise and remove the tip from the syringe, swaying it from side to side.

8.2 Attaching an extension of the applicator

1. While inside the sterile zone, carefully remove the extension cord from the individual packaging.
2. Slide the wide end of the extension onto the standard applicator and press until it is full. The extension cord is considered fully worn when the tapered end of the standard applicator fits into the inside of the tapered extension, as shown in Figure 1.

Fig. 8.1 - Fully attached

Fig. 8.2 - Not fully attached

ATTENTION! Do not press the piston while attaching the extension cord to the applicator. Pressing the piston can be avoided: (1) removing the piston before attaching the extension cord and returning it to its place after attaching, or (2) carefully taking care not to touch the piston and not to press it against any surface.

Fig. 8.3 - Wrong

Fig. 8.4 - Right

8.3 Preparing and filling applicators

1. Before applying BioGlue Surgical Adhesive, the applicator must be filled. The surgeon must press the plunger and release a narrow strip of BioGlue glue (approximately 3 cm long) onto a sterile disposable surface (such as a swab, gauze ball or tissue).

2. The surgeon should inspect the material that came out of the tip during filling, and make sure that the material is uniform yellow or amber in color and does not contain air bubbles. If this material is colorless or contains bubbles, repeat the filling until the product provides a homogeneous liquid without bubbles.

ATTENTION! Avoid direct contact with material extruded from the tip when filling it.

3. After filling the tip correctly, proceed immediately to the application of glue.

ATTENTION! BioGlue Surgical Adhesive cures very quickly. The surgeon must apply BioGlue immediately after filling the tip. The delay in applying glue to the fabric after filling the tip can lead to its polymerization in the tip itself. If this happens, replace the tip / extension with a new one and repeat the filling procedure.

Do not press the piston if the applicator is clogged.

NOTE. The extension cord is designed to simplify the delivery of BioGlue Surgical Adhesive and cannot be used as a surgical tool (i.e. to compress, raise, or hold tissue)

9. Transportation, storage and operation

Storage conditions: This product should be stored in a well-ventilated room with a humidity of not more than 80% at a temperature of +15 to + 25 ° C. Handle with care. Store the product in the manufacturer's packaging in a dry, clean place, away from heat and direct sunlight, at a controlled temperature. Shelf life is indicated on the packaging of each product.

The product is transported by all types of vehicles in accordance with the conditions:

Transportation conditions: This product should be stored in a well-ventilated room with a humidity of not more than 80% at a temperature of +25 to + 35 ° C.

Operating conditions: air temperature from + 25 ° C to + 42 ° C at relative humidity from 35% to 65% (non-condensing) and atmospheric pressure 500 - 1060 hPa.

10. Information on the sterility of the product

Standard applicators and extended-tip applicators are sterilized by gamma radiation and are intended for use in a single procedure and in only one patient. Elongated applicators and applicator extensions are sterilized in primary packaging with ethylene oxide gas and are intended for use in a single procedure and in only one patient. Standard applicators, extended-tip applicators, elongated applicators, and applicator extensions are used for a short time and do not involve permanent implantation into the patient's body.

Sterility level not less than 10^{-6}

The value of the minimum absorbed dose of 25 kGy.

11. Disposal information

Dispose the product in accordance with SanPiN 2.1.7.2790-10 for class B waste (epidemiologically hazardous waste). After hardware-based disinfection methods using physical methods and changing the appearance of the waste, which excludes the possibility of their reuse, class B waste can be accumulated, temporarily stored, transported, disposed of and disposed of together with class A waste (epidemiologically safe waste in composition close to MSW). The packaging of Class B medical disinfected medical waste must be labeled to indicate that the waste has been decontaminated. In case of non-use of the product, damage to the individual packaging and expiration of the storage period, dispose in accordance with San PiN 2.1.7.2790-10 for class A waste.

12. Warranty

The manufacturer makes no warranties, either express or implied, including, but not limited to, the warranties of merchantability or fitness for a particular purpose. The manufacturer assumes no liability for incidental or consequential damages. If this statement is declared invalid or invalid, for any reason: (i) any action related to the violation of this guarantee must be started within one year after the occurrence of the event that constitutes the basis for the claim, and (ii) liability limited to product replacement. Prices, specifications and availability are subject to change without notice.

Guaranteed shelf life is 5 years, the shelf life of sterility products is 5 years.

Manufacturer: CryoLife, Inc., USA

Manufacturer Address: 1655 Roberts Blvd., NW, Kennesaw, Georgia, 30144 USA

Phone: (770) 419 33 55

Email: usbioglue@cryolife.com

Authorized representative of the manufacturer in the Russian Federation:

Firm Finco Limited Liability Company

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