Date: November 2, 2021

To whom it may concern;

DECLARATION


We, Nipro Corporation, hereby declare that the attached document is a true and original of Blood tubing set's IFU for Russia registration.

Sincerely yours,

Y. Nagasawa

YOSHIKI Nagasawa,
General Manager,
Medical Regulatory Affairs Department,
Quality Assurance & Regulatory Compliance Headquarters
NIPRO CORPORATION

APPROVED BY

General Manager

Yoshiki Nagasawa

Y. Nagasawa

Instuction for Use
Blood tubing set for hemodialysis
(Nipro Corporation, Japan)

2021

1. Product name, information about manufacturer:

Blood tubing set for hemodialysis (hereinafter – blood lines, kit, lines).

Design variants of the product:

1) Blood tubing set for hemodialysis A146R-V715R (REF: A146R-V715R-GT) including:

- Arterial blood line A146R with connected plastic “spike” needle with cap;
- Venous blood line V715R with connected recirculation connector;
- Drain bag.

2) Blood tubing set for hemodialysis A165R-V733R (REF: A165R-V733R-GT) including:

- Arterial blood line A165R with connected plastic “spike” needle with cap;
- Venous blood line V733R with connected drain bag;
- Recirculation connector.

3) Blood tubing set for hemodialysis A196R-V764R (REF: A196R-V764R-GT) including:

- Arterial blood line A196R with two connected plastic “spike” needles with caps;
- Venous blood line V764R with connected drain bag;
- Luer lock protection cap;
- Recirculation connector.

4) Blood tubing set for hemodialysis A197R-V765R (REF: A197R-V765R-GT) including:

- Arterial blood line A197R with two connected plastic “spike” needles with caps;
- Venous blood line V765R with connected drain bag;
- Line for substitution solution with connected recirculation connector;
- Luer lock protection cap;
- Recirculation connector.

5) Blood tubing set for hemodialysis A329R-V827R (REF: A329R-V827R-GT) including:

- Arterial blood line A329R with connected plastic “spike” needle with cap;
- Venous blood line V827R with connected recirculation connector;
- Drain bag.

6) Blood tubing set for hemodialysis A316R-V814R (REF: A316R-V814R-GT) including:

- Arterial blood line A316R with connected plastic “spike” needle with cap;
- Venous blood line V814R with connected drain bag;
- Recirculation connector.

7) Blood tubing set for hemodialysis A317R-V815R (REF: A317R-V815R-GT) including:

- Arterial blood line A317R with connected plastic “spike” needle with cap;
- Venous blood line V815R with connected recirculation connector;
- Infusion line with connected recirculation connector and plastic “spike” needle with cap;
- Drain bag.

8) Blood tubing set for hemodialysis A363R-V849R (REF: A363R-V849R-GT) including:

- Arterial blood line A363R with connected recirculation connector;
- Venous blood line V849R with two connected extension lines;
- Plastic “spike” needle with cap;
- Luer lock protection cap.

9) Blood tubing set for hemodialysis A364R-V850R (REF: A364R-V850R-GT) including:

- Arterial blood line A364R;
- Venous blood line V850R with connected extension line;
- Line for substitution solution with connected recirculation connector;
- Recirculation connector.

10) Blood tubing set for hemodialysis A365R-V851R (REF: A365R-V851R-GT) including:

- Arterial blood line A365R with connected recirculation connector;
- Venous blood line V851R with two connected extension lines;
- Plastisc "spike" needle with cap;
- Luer lock protection cap.

11) Blood tubing set for hemodialysis A366R-V852R (REF: A366R-V852R-GT) including:

- Arterial blood line A366R;
- Venous blood line V852R with connected extension line;
- Line for substitution solution with connected recirculation connector;
- Recirculation connector.

12) Blood tubing set for hemodialysis A274R-V785R-UNO (REF: A274R-V785R-UNO-GT) including:

- Arterial blood line A274R with connected plastic "spike" needle with cap;
- Venous blood line V785R-UNO with connected drain bag;
- Recirculation connector.

13) Blood tubing set for hemodialysis A434R-V914R (REF: A434R-V914R-GT) including:

- Arterial blood line A434R;
- Venous blood line V914R with connected drain bag;
- Line for substitution solution;
- Recirculation connector.

Manufacturer: Nipro Corporation, Japan
3-9-3, Honjo-Nishi, Kita-ku, Osaka 531-8510, Japan

Production site:

Nipro (Thailand) Corporation Limited, Thailand.
10/2 Moo 8, Bangnomko, Sena, Phra Nakhon Si Ayutthaya 13110, Thailand

2. Intended use and application of product

Intended Use

This blood tubing set is indicated for use during hemodialysis to provide blood flow.

Application condition

Used in medical institutions by specialists nephrologists

3. Special characteristics

The various bloodline sets are to be used depend on various models of dialysis machines. The blood tubing sets include Arterial blood lines, Venous blood lines and may also include drain bags, infusion lines, lines for substitution solution, extension lines, caps, "spike" needles, recirculation connectors which are used to pump blood, infuse medications or fluids, sample blood and monitor pressure during dialysis treatment. The appearance of products is given in a attachment 1.

4. General indications, contradictions, potential complications,

Indications:

When used in hemodialysis, it is part of an artificial kidney system for treatment of patients with renal failure or toxemic conditions. Compliance with possible configurations of a blood tubing set is determined by the responsible specialist.

Contraindications:

Allergic reaction for PVC and silicone rubber.

Potential complications:

As this blood tubing set is not, in itself, a medical treatment device, there is no absolute potential complications to its use.

5. Warnings & Cautions

1. This blood tubing set is for a single use only. Do not reuse.
*Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity.
2. Do not use if the package is damaged.
3. Aseptic technique is required throughout set up and use.
4. Inspection prior to use and constant monitoring during the dialysis procedure is required in order to insure the product is in proper condition for use.
5. Air entering the extracorporeal blood circuit may cause fatal air embolism. Use of an air detector is recommended, in any case, continually observe the air trap for adequate blood level (at least filled to $\frac{3}{4}$ guide line of the air trap).
6. Leakage at joints, connections, or any point on the extracorporeal circuit may cause blood loss or air embolism. Observe carefully for leaks before and during treatment and take corrective measures as necessary.
7. Avoid kinking or occlusion of tubing during treatment. Excessive pressure may damage the extracorporeal circuit or blood access site, or may cause disconnection and/or blood loss.
8. Fluid in the monitor lines should be avoided.
9. Keep the infusion line clamped except when administering fluid.
10. Serrated metal hemostats can cut or break heparin lines and should not be used with these blood tubing sets.
11. Keep air adjustment lines clamped except when adjusting the air trap level.
12. There may be slight variations in the vascular access devices and connectors provided by different manufacturers. Check to make certain that the vascular access utilized is seated and taped securely. Failure of connectors to fit tightly can result in serious blood loss and/or air embolism.
13. Do not use any needle larger than 20 gauge to puncture blood sample injection site.
14. Make certain that all connections are firm and secure.
15. Immediately prior to puncture, thoroughly disinfect the injection port (Ethyl alcohol 70% or Isopropyl alcohol 70%).
16. Do not store at extreme temperature and humidity.
17. No conclusive or corroborated scientific evidence exists that phthalate exposure via medical treatments has harmful effects in humans. Multiple or prolonged medical treatments may increase the exposure level in patients. Special consideration should be given to potentially sensitive patients, such as children, and nursing or pregnant woman. (For blood tubing set containing phthalate: A146R-V715R, A329R-V827R, A165R-V733R, A274R-V785R-UNO, A316R-V814R, A317R-V815R, A196R-V764R, A197R-V765R, A434R-V914R)

6. Instruction for use

Set up the dialyzer and blood tubing set according to the instructions provided with the dialyzer and blood tubing set.

1. Remove blood tubing set from protective plastic bag and connect arterial (red) and venous (blue) lines to the blood pump, monitoring equipments, dialysis machine and dialyzer. Make sure the pump segment is not kinked or twisted.
2. Connect the monitor lines to the pressure monitors utilizing hydrophobic filters.
3. Make sure all connections are secure and made in the proper manner before proceeding.
4. Insert the air trap or line into an air detector.
5. Prime the system according to dialyzer manufacturer's direction sheet.
6. Adjust air trap levels to $\frac{3}{4}$ guide line of the air trap using a sterile syringe connected to the air adjustment line.

7. Check blood tubing set carefully making sure that all connections are secure and there are no leaks in the blood lines.

8. Initiate dialysis in accordance with dialyzer and equipment manufacturers' instructions.

7. Transportation, storage and operation

Products are transported by all modes of transport in closed vehicles in accordance with the rules for the carriage of goods applicable to this type of transport in accordance with the conditions:

Transportation conditions:	Temperature: +5°C - +40° C
	Humidity: 20-80 %
	Atmosphere pressure: 0-0,1 MPa
Storage conditions and Operating conditions:	Temperature: +5°C - +40° C
	Humidity: 20-80 %
	Atmosphere pressure: 0-0,1 MPa

Store the product in the manufacturer's packaging in a dry, clean place, away from heat and direct sunlight, at a controlled temperature. Shelf life is indicated on the label of each product.

8. Information on the sterility of the product

Bloodlines of the sets are supplied sterile.

Sterilization parameters:

- radiation sterilization method, gamma radiation in accordance with the requirements of ISO 11137-1.

The product is apyrogenic.

Re-sterilization of the product is prohibited. Reuse is prohibited

9. Disposal Information

Dispose of the product in accordance with SanPiN 2.1.3684-21 for waste class B (epidemiologically hazardous waste). After hardware methods of disinfection using physical methods and changing the appearance of waste, eliminating the possibility of their reuse, waste of class B can be accumulated, temporarily stored, transported, destroyed and disposed of together with waste of class A (epidemiologically safe waste similar in composition to MSW). The packaging of disinfected medical waste of class B must be labeled to indicate that the waste has been disinfected. Dispose should be performed in accordance with SanPiN 2.1.3684-21 for wastes class A if product was not used, individual package damaged or product expired.

10. Guarantee

1. Our blood tubing set is manufactured under strict quality control and the quality is assured. If our blood tubing set is defective (broken package, damaged product), however, it shall be replaced with a new one at our cost upon return of the broken package or damaged product. We will not be responsible, however, for the injury on a patient or any person or the damage to any object that is attributed to transport, storage and operation in your institution.

2. If a patient or any person is injured or any object is damaged by use of our blood tubing set, we will not be responsible for the injury or damage unless we are clearly identified as being at fault.

3. If a patient or any person is injured or any object is damaged by reuse of our blood tubing set, we will not be responsible for the injury or damage of any nature.

4. We will not be responsible for any injury or damage caused by use of our blood tubing set after the expiry date mentioned on the packages.

Information about expire storage date see on package. Expire date shown in YYYY-MM-DD format, where DD – last day of the month despite the production date. Maximum shelf life/preservation sterility is 5 years.

Manufacturer: Nipro Corporation, Japan

1-2 Hanjo-Nishi, Kita-ku, Osaka 531-8510, JAPAN

Authorized representative of the manufacturer in the Russian Federation:

Company «BW» LLC

125035, Moscow, Prospect Mira, 101 bld.1, room 17

+7(499) 2816768

1. Marking and packaging

Each product is sealed in a package for final sterilization (individual packaging). Further, together with the instructions for use - in a box of cardboard. The product is sterile. Application in the event of packaging failure is strictly prohibited.

Packaging complies with ISO 11607-1-2018 and ISO 11607-2-2018.

The marking of individual package contains the lot number. The decoding number is shown in the figure 2.

LOT / 

19 A 21/2023-12-31

Expiry date of month e.g.30,31

Expiry month e.g. 01 (January) - 12 (December)

Expiry year e.g. 2023, 2024, ...

Manufacturing date e.g. 01 - 31



Manufacturing month e.g. A (January), ... , L (December)

Manufacturing year e.g. 19 (2019), 20 (2020), ...

Figure 2 – Lot number decoding

The marking of the individual packaging of the product contains the following information (including information on the package itself (original marking) and sticker):

Product name	Blood tubing set for hemodialysis A364R-V850R (REF: A364R-V850R-GT) including*: - Arterial blood line A364R; - Venous blood line V850R with connected extension line; - Line for substitution solution with connected recirculation connector; - Recirculation connector. *The exact design variant indicated as an example
Identification code of the product (Design variant Code)	CODE
Batch code (LOT)	LOT
REF number	REF
Expiration Date	
Pump segment size (internal diameter x length)	PUMP SEGMENT (ID x Length)
Priming volume	PRIMING VOLUME
Indication of compatible equipment	EQUIPMENT

Symbol "Radiation sterilization"	 STERILE R
Symbol "Prohibition of reuse"	
The symbol of "Caution! Refer to the instructions for use"	
The symbol "He use when damaged packaging"	
do not contain latex	
Temperature Range	 On sticker: Transport at temperature from +5 to +40°C. Store at temperature from +5 to +40°C.
Symbol "Keep away from moisture"	
Symbol "Avoid exposure to sunlight"	
Symbol "not pyrogenic"	
only as prescribed	Rx only
Phthalates (For bloodlines: A146R-V715R, A329R-V827R, A165R-V733R, A274R-V785R-UNO, A316R-V814R, A317R-V815R, A196R-V764R, A197R-V765R, A434R-V914R).	
Do not contain phthalates (For bloodlines: A363R-V849R, A364R-V850R, A365R-V851R, A366R-V852R)	
Manufacturer Information	
Information about the place of production	Nipro (Thailand) Corporation Limited 10/2 Moo 8, Bangnomko, Sena, Phra Nakhon Si Ayutthaya 13110, Thailand On sticker: Made in Thailand
Authorized representative at the territory of Russian Federation	Authorized representative at the territory of Russian Federation: Company «BW» LLC 129085, Moscow, Prospect Mira, 101 bld.1, room 17
Not toxic	Not toxic
Product use	Blood tubing set for hemodialysis are intended to provide blood flow during hemodialysis
Information about registration certificate	Registration certificate: No. from

For transportation, products are packed in corrugated cardboard group package. The number of packages in the box is shown in table 10. Product is transported in group package.

12. Compliance with the standards of the Russian Federation

GOST ISO 8638-2012 «Blood tubing sets for haemodialysers, haemodiafilters and haemoconcentrators.

Technical requirements and test methods».

GOST R IEC 62366-2013 «Medical devices. Designing medical devices with usability».

GOST ISO 11607-1-2018. «Packaging for terminally sterilized medical devices. Part 1. Requirement for materials, barrier systems for sterilization and packing systems».

GOST ISO 11607-2-2018. «Packaging for terminally sterilized medical devices. Part 2. Requirements for the validation of the processes of formation, sealing and assembly».

GOST EN 556-1-2011 «Sterilization of medical devices. Requirements for medical devices of the category «sterile». Part 1. Requirements for medical devices subject to final sterilization».

GOST R ISO 15223-1-2020 «Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements»

GOST 31214-2016. «Medical devices. Requirements for samples and documentation presented for toxicological tests, sanitary and chemical analyses, tests for sterility and pyrogenicity».

GOST R 52770-2016. «Medical devices. Safety requirements. Methods of sanitation-chemical and toxicological tests».

GOST 31209-2003 «Containers for blood and its components. Requirements for chemical and biological safety and methods of testing».

GOST ISO 10993-1-2011. «Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing».

GOST R ISO 10993-2-2009. «Medical devices. Biological evaluation of medical devices. Part 2. Animal welfare requirements ».

GOST ISO 10993-4-2011. «Medical devices. Biological evaluation of medical devices. Part 4. Selection of tests for interactions with blood».

GOST ISO 10993-5-2011. «Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity».

GOST ISO 10993-10-2011. «Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and delayed-type hypersensitivity».

GOST ISO 10993-11-2011 «Medical devices. Biological evaluation of medical devices. Part 11. Tests for systemic toxicity».

GOST ISO 10993-12-2015. «Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference materials».

State pharmacopeia. RF OFS.1.2.4.0005.15. Pirogenity.

PND F 14.1: 2: 4.139-98 «Quantitative chemical analysis of water. Methodology for measuring mass concentrations of cobalt, nickel, copper, zinc, chromium, manganese, iron, silver, cadmium and lead in samples of drinking, natural and waste water by atomic absorption spectrometry».

MU 4628-88 «Guidelines for gas chromatographic determination of residual monomers and non-polymerizable impurities released from polystyrene plastics in water, model environments and food products».

MUK 4.1.650-96 «Guidelines for the gas chromatographic determination of acetone, methanol, benzene, toluene, ethylbenzene, pentane, o-, m-, p-xylene, hexane, octane and decane in water».

MUK 4.1.653-96 «Guidelines for the reaction-chromatographic determination of formaldehyde in water».

Attachment 1

The appearance of products

The appearance of the Blood tubing set on example of Blood tubing set for hemodialysis A146R-V715R:


Diagram of the Blood tubing set on example of Blood tubing set for hemodialysis A165R-V733R:


Diagram of the Blood tubing set on example of Blood tubing set for hemodialysis A196R-V764R:





Diagram of the Blood tubing set on example of Blood tubing set for hemodialysis A197R-V765R:


of the Blood tubing set on example of Blood tubing set for hemodialysis A329R-V827R:



The appearance of the Blood tubing set on example of Blood tubing set for hemodialysis A316R-V814R:

Appearance of the Blood tubing set on example of Blood tubing set for hemodialysis A317R-V815R:


appearance of the Blood tubing set on example of Blood tubing set for hemodialysis A363R-V849R:


appearance of the Blood tubing set on example of Blood tubing set for hemodialysis A364R-V850K:


ARTERIAL BLOOD LINE A364R


VENOUS BLOOD LINE V850R

appearance of the Blood tubing set on example of Blood tubing set for hemodialysis A365R-V851R:


appearance of the Blood tubing set on example of Blood tubing set for hemodialysis A366R-V852R:


the appearance of the Blood tubing set on example of Blood tubing set for hemodialysis
A274R-V785R-UNO:


ARTERIAL BLOOD LINE A274R

VENOUS BLOOD LINE V785R-UNO

The appearance of the Blood tubing set on example of Blood tubing set for hemodialysis A434R-V914R:

Attachment 2

Table 1 - Material list for main component of product

Componenet	Material
Air adjusting line	Polyvinylchloride
Saline line	Polyvinylchloride
Monitor line	Polyvinylchloride
Heparine line	Polyvinylchloride
Pump segment	Polyvinylchloride
Extension line	Polyvinylchloride
Line for substitution solution	Polyvinylchloride
Recirculation connector	Polypropylene
Luer lock protection cap	Polypropylene
Dialyzer connector and Dialyzer connector cap	Polyvinylchloride
Air trap	
Injection port	
T-valve	Acrylonitrile-butadiene-styrene
Drain bag	Polyvinylchloride
plastic "spike" needle	Acrylonitrile-butadiene-styrene

*Phtalates are part of PVC material

Table 2– Weight and size parameters of package with product

	Design variants												
	A146R - V715R	A329R - V827R	A165R - V733R	A316R - V814R	A317R - V815R	A363R - V849R	A364R - V849R	A365R - V851R	A366R - V852R	A196R - V764R	A197R - V765R	A274R - V785R -UNO	A434R - V914R
Size: Individual package, (± 5%)	250x320 mm												
SizeCarton, (± 5%)	520x430x250 mm												
Total carton weight (kg), (± 5%)	8.28	8.15	9.34	9.15	8.54	8.27	7.91	9.02	8.88	8.22	8.89	9.22	9.49

Mechanical characteristics

and tubing set:

sealed at an overpressure in the internal cavity that is 1.5 times the maximum pressure set by the manufacturer
 $600 \text{ mmHg} \times 1.5 = 900 \text{ mmHg} = 1.223 \text{ kgf/cm}^2$

sealed at vacuum pressure in the internal cavity minus 200 mmHg

the pump segment of the blood tubing provides continuous operation for 24 hours under the worst operating conditions.

the connecting dimensions of the parts that provide the connection of lines with hemodialyser, hemofilter or hemoconcentrator comply with GOST ISO 8638-2012 (Figure 1-3)


Figure 1 - The main connecting dimensions of the parts that provide the connection of the blood tubings with hemodialyser, hemofilter or hemoconcentrator


Figure 2 - The length of the cone mating parts that provide connection with the blood supplying cavity


Figure 3 - Details of conical devices with a 6% cone to control the connection with the blood-supplying cavity of hemodialyser, hemofilter or hemoconcentrator

All elements of the blood tubings that are designed for use with assistive devices, such as heparin lines and pressure sensor lines, end in a threaded connection with an internal cone of 6% (Luer) injection units that are part of the kit's lines remain airtight after piercing with a needle.

Physical characteristics

Table 3 - Main characteristics

Design variants	Pump Segment Size for 100-300 mm. ± 15 mm for 301-500 mm. ± 25 mm	Priming Volume (ml) $\pm 20\%$.	Quantity of individual packages with product in carton (pcs)
A146R-V715R	8.0x12.0x240	133	24
A165R-V733R	8.0x12.0x326	169	24
A196R-V764R	8.0x12.0x234	191	20
A197R-V765R	8.0x12.0x234	192	18
A329R-V827R	8.0x12.0x350	126	24
A316R-V814R	8.0x12.0x350	156	24
A317R-V815R	8.0x12.0x240	110	24
A363R-V849R	8.0x12.0x233	132	24
A364R-V850R	8.0x12.0x233	132	20
A365R-V851R	8.0x12.0x233	231	20
A366R-V852R	8.0x12.0x233	231	18
A274R-V785R-UNO	8.0x12.0x350	171	24
A434R-V914R	8.0x12.0x310	153	20

Maximum pressure value inside circuit -200 mmHg - +600 mmHg

Compatible equipment indicated on each individual label (table 4).

Table 4- Compatible Equipment (not included into supply set)


№	Design variant	As indicated on blood tubing set package	Manufacturer of the equipment	Description	
				List of equipment models	Treatment
1	A146R-V715R	SURDIAL55PLUS	Nipro Corporation, Japan	Surdial 55, Surdial 55 Plus	Hemodialysis
2	A165R-V733R	FOR B.BRAUN DIALOG	B. Braun Avitum AG, Germany	Dialog, Dialog +	Hemodialysis
3	A196R-V764R	FOR NIKKISO DBB-03/DBB-05/DBB-07 HD DN	NIKKISO Europe GmbH, Germany	DBB-03, DBB-05, DBB-07	Hemodialysis
4	A197R-V765R	FOR NIKKISO DBB-07 ON LINE HDF	NIKKISO CO., LTD, Japan	DBB-07	Hemodialysis and Hemodiafiltration On-line
5	A329R-V827R	UNIVERSAL 1 BLOOD LINE	Nipro Corporation, Japan	Surdial 55, Surdial 55 Plus	Hemodialysis
			OOO Zavod Dizet, Russia	Malahite	Hemodialysis
			Fresenius Medical Care AG & Co. KGaA, Germany	4008S, 4008B	Hemodialysis
			Bellco S.r.l., Italy	Formula, Formula 2000, Formula 2000 plus	Hemodialysis
			B. Braun Avitum AG, Germany	Dialog, Dialog +	Hemodialysis
6	A316R-V814R	UNIVERSAL LINE	Nipro Corporation, Japan	Surdial 55, Surdial 55 Plus	Hemodialysis
			OOO Zavod Dizet, Russia	Malahite	Hemodialysis
			Fresenius Medical Care AG & Co. KGaA, Germany	4008S, 4008B	Hemodialysis
			Bellco S.r.l., Italy	Formula, Formula 2000, Formula 2000 plus	Hemodialysis
			B. Braun Avitum AG, Germany	Dialog, Dialog +	Hemodialysis
7	A317R-V815R	SURDIAL55PLUS	OOO Zavod Dizet, Russia	Malahite	Hemodialysis
			Nipro Corporation, Japan	Surdial 55, Surdial 55 Plus	Hemodialysis
8	A363R-V849R	SURDIAL X HD DIF	Nipro Corporation, Japan	Surdial X	Hemodialysis
9	A364R-V850R	SURDIAL X HDF	Nipro Corporation, Japan	Surdial X	Hemodiafiltration on line
10	A365R-V851R	SURDIAL X HD SN DIF	Nipro Corporation, Japan	Surdial X	Single needle hemodialysis
11	A366R-V852R	SURDIAL X HDF SN	Nipro Corporation, Japan	Surdial X	Single needle hemodiafiltration On-line
12	A434R-V914R-GT	HDF DN FOR B.BRAUN DIALOG	B. Braun Avitum AG, Germany	Dialog, Dialog +	Hemodialysis and Hemodiafiltration On-line
13	A274R-V785R-UNO	UNIVERSAL 3 BLOOD LINE	Nipro Corporation, Japan	Surdial 55, Surdial 55 Plus	Hemodialysis
			OOO Zavod Dizet, Russia	Malahite	Hemodialysis
			Fresenius Medical Care AG & Co. KGaA, Germany	4008S, 4008B	Hemodialysis
			Bellco S.r.l., Italy	Formula, Formula 2000, Formula 2000 plus	Hemodialysis
			B. Braun Avitum AG, Germany	Dialog, Dialog +	Hemodialysis


Registered No. 2175 -2021

NOTARIAL CERTIFICATE

This is to certify that Mr. ITO Hisayuki an agent of Mr. Yoshiki NAGASAWA, General Manager, Medical Regulatory Affairs Department, Quality Assurance & Regulatory Compliance Headquarters of NIPRO CORPORATION has stated in my very presence that said Mr. Yoshiki NAGASAWA acknowledged himself to have signed to the attached document.

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This public document

2. has been signed by TAMAKI Shunji

3. acting in the capacity of Notary of the Osaka Legal Affairs Bureau

4. bears the seal/stamp of TAMAKI Shunji, Notary
Certified

5. at Osaka

6. NOV. 12, 2021

7. by the Ministry of Foreign Affairs

8. 21-No 010159

9. Seal/stamp:

10. Signature

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FUKUTANI Katsumi

For the Minister for Foreign Affairs

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