
Managing Director (General Director)

"Valentine International Ltd."

Chi Yui  Yuhsin Chiu
Signature, Place for seal

"October 07", 2024

Operating Instructions for a medical device

"Barry anti-bedsore system with accessories"

Manufacturer:
"Valentine International Ltd",

8F, No. 149, Sec. 2, Ta Tung Road, Hsichih Dist., New Taipei City 221, Taiwan,
R.O.C.

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北院民認麟字	日期
案號 116102	Date: OCT. 21 2024
Case No.	
本文件之簽名或蓋章，於臺灣臺北地方法院所屬民間公證人南京聯合事務所認證。公證人陳幼麟	
Attested at the NanJing Notary Public Office of Taiwan	
Taipei District Court, R.O.C., that the signature(s)/seal(s) in	
this document is/are authentic.	
Chen, Yu-Lin	
Notary Public 	
3F-1, 285, Sec.3, Nanjing E. Rd., Taipei, Taiwan, R.O.C.	

2024

Table of Contents

1. Name and version of the medical device.....	3
2 Information on the manufacturer and the manufacturer's authorized representative	3
3 Purpose of the medical device	4
4. Potential consumers of the product	4
5. Special properties of the medical device. Indications and contraindications	4
6. Design, principle of operation of a medical device	6
6.1 Design and operating principle.....	錯誤! 尚未定義書籤。
6.2 Design variants, composition and appearance of the product	6
6.3 Description of structural components and composition elements	錯誤! 尚未定義書籤。
7. Technical specification	12
7.1 Basic technical characteristics	錯誤! 尚未定義書籤。
8. Medical device labeling.....	20
9. Packaging, transportation, marking of medical device packaging	22
10. Disposal	24
11. List of normative documentation.....	24
12. Operation	25
11.1 Preparing for operation.....	錯誤! 尚未定義書籤。
11.2 Operating procedures.....	錯誤! 尚未定義書籤。
13. Sterility, cleaning and disinfection	26
14. Maintenance.....	26
15. Warranty obligations	27
16. Contact info	27
Annex 1: Electromagnetic compatibility	27

Valentine International Ltd.	Operating Instructions	Version No. 2.01	Date 07/10/2024.
-------------------------------------	-------------------------------	-------------------------	-------------------------

Dear User!

Thank you for the trust you have placed in us by purchasing the Barry anti-bedsore system with accessories.

Please read this Operation Manual carefully and completely before starting work and strictly follow all procedures described in it.

This User Manual is intended to study the device, principle of operation, technical characteristics of the bed and contains information necessary for its proper operation and maintenance.

1. Name and version of the medical device

Barry anti-bedsore system with accessories (hereinafter referred to as "system").

Barry anti-bedsore system with accessories is available in the following designs:

Barry Primo; Barry Mezzo Basic; Barry Serio; Barry Serio Plus; Barry Serio Forte; Barry Serio Forte 2; Barry Mezzo Plus; Barry Mezzo Forte; Barry Mezzo Forte 2; Barry Mezzo Premium; Barry Mezzo Pro.

The medical device may be operated by persons who have familiarized themselves with this operating manual and who have undergone appropriate safety training.

2. Information about the manufacturer and the manufacturer's authorized representative

Manufacturer:

Valentine International Ltd.

8F, No. 149, Sec. 2, Ta Tung Road, Hsichih Dist., New Taipei City 221, Taiwan, R.O.C.

tel: +886.2.8692.6662

e-mail: info@valentine-dynalife.com

website: <http://www.valentine-dynalife.com/>

Place of production:

Xiamen Senyang Co., Ltd.

Floor 2, No. 189-3 and No. 189-1, Tongji South Road, Jimei District, Xiamen City, Fujian Province, 361000, P.R. China.

Tel: +86 15805935377

e-mail: vip@medicalxm.com

website: www.medicalxm.com

Authorized representative of the manufacturer in the territory of the Russian Federation (APR):
Limited Liability Company SIMS-2 (LLC SIMS-2).

125430, Moscow, Mitinskaya street. 16, est. 10, pom. 1012B, rooms 15 to 18


tel.: +7(495)792-31-90

e-mail: info@sims2.ru

website: www.sims2.ru

The design, manufacture and sale of the product complies with the requirements: ISO 13485:2016.
Certificate number SN23/00006519. Date of issue: 2023-12-05, valid until 2026-12-04, notified body TUV Rheinland LGA Products GmbH.

Class depending on the degree of potential risk of use for medical purposes in accordance with the nomenclature classifier of medical devices approved by the order of the Ministry of Health of the Russian Federation dated June 06, 2012 № 4n: 1.

Classification of the product in accordance with the nomenclature classifier of medical devices approved by the order of the Ministry of Health of the Russian Federation dated June 06, 2012 № 4n: 266750.

OKPD 2 code: 31.03.12.140

3. Purpose of medical device

Pressure relief support system is designed for prevention and treatment of medium and high risk bedsores in patients with injuries and diseases of the musculoskeletal system and central nervous system, severe somatic diseases accompanied by prolonged forced immobility.

Field of application: hospitals of medical and preventive institutions, rehabilitation, traumatology and orthopedics departments, palliative care units, hospices, residential homes for the elderly and disabled, for care of bedridden patients or elderly people at home.

Conditions of use: used as an additional mattress, for patient positioning.

Method of application: according to the operating manual.

4. Potential consumers of the product

Skilled nursing facilities, patients in various health care facilities, nursing homes and comparable facilities, people with musculoskeletal injuries and diseases who need care at home.

5. Special properties of the medical device. Indications and contraindications

5.1 Sterility

The medical device is supplied non-sterile and is not subject to subsequent sterilization.

5.2 Presence of medicinal substances

This medical device does not contain as an ingredient substances which, when used alone, can be evaluated as medicinal substances or substances as defined in Article 1 of Directive 65/65/EEC.

5.3 Content of human blood derivatives

This medical device does not contain human blood or blood components.

5.4 Content of components from animal tissues

This medical device does not contain components derived from animal tissue.

5.5 Phthalate content

This medical device does not contain phthalates.

5.6 Content of mutagenic and carcinogenic components


This medical device does not contain any mutagenic or carcinogenic components.

5.7 Treatment of AP or APG products

This medical device is not an AP or APG product.

5.8 Software

This medical device does not contain any software.

5.9 Indications

It is used for the prevention and treatment of pressure sores in patients on bed rest for long periods of time, allowing to maintain normal tissue trophicity by regularly changing the pressure on the soft tissues adjacent to the mattress.

- Care of bedridden patients;
- is used for patients with severe stages of oncology;
- It is used to improve blood circulation in tissues;
- It is used to neutralize muscle tension at the moment of rest;
- It is used to prevent the formation and progression of pressure sores;
- It is used to prevent the occurrence of defects of the musculoskeletal system;
- It is used to normalize sleep;
- It is used to resist the appearance and development of ulcers.

5.10 Contraindications


- Cervical and skeletal traction;
- Unstable spinal cord injuries.

5.11 Possible side effects

When properly cared for and maintained according to the instruction manual, this medical device has no side effects.

5.12 Precautions

- Always unplug the product when not in use;
- Do not use the product while bathing;
- Do not store the product in the bathroom;
- Protect the device from moisture;
- If the product falls into water, disconnect it from the mains immediately;
- Do not leave the appliance unattended when it is switched on;
- Continuous monitoring is required if the product is used around children or people with disabilities;
- Do not use additional accessories and attachments not specified by the manufacturer;
- Do not install the product near heated surfaces;
- Do not clog or block the openings on the product;
- Do not use the product compressor outdoors or in areas with flammable or chemical substances;

6. Design, composition, appearance of a medical device

6.1 Design variants, composition and appearance of the product

I. Barry anti-bedsore system with accessories, Barry Primo version, as part of:

1. Barry Primo compressor – 1 pc.
2. Barry Primo mattress – 1 pc.
3. Connecting tubes – 1 pc.
4. Instruction manual – 1 pc.

Accessories:

1. Barry Primo repair kit, as part of:

- self-adhesive patch - 2 pcs.

II. Barry anti-bedsore system with accessories, version Barry Mezzo Basic, as part of:

1. Barry Mezzo Basic compressor – 1 pc.
2. Barry Mezzo Basic mattress – 1 pc.
3. Connecting tubes – 1 pc.
4. Instruction manual – 1 pc.

Accessories:

1. Removable case B – 1 pc.

2. Barry Mezzo Basic repair kit, as part of:

- self-adhesive patch - 2 pcs.:
- spare cylinder B - 1 pc.

III. Barry anti-bedsore system with accessories, Barry Serio version, as part of:

1. Barry Serio compressor – 1 pc.
2. Barry Serio mattress – 1 pc.
3. Connecting tubes – 1 pc.
4. Instruction manual – 1 pc.


Accessories:

1. Barry Serio repair kit, as part of:

- self-adhesive patch - 2 pcs.

IV. Barry anti-bedsore system with accessories, Barry Serio Plus version, as part of:

1. Barry Serio Plus compressor – 1 pc.

3. Connecting tubes -- 1 pc.

4. Instruction manual - 1 pc.

Accessories:

1. Barry Serio Plus repair kit, as part of:

- self-adhesive patch - 2 pcs.

V. Barry anti-bed sore system with accessories, Barry Serio Forte version, as part of:

1. Barry Serio Forte compressor – 1 pc.

2. Barry Serio Forte mattress – 1 pc.

3. Connecting tubes – 1 pc.

4. Instruction manual – 1 pc.

Accessories:

1. Barry Serio Forte repair kit, as part of:

- self-adhesive patch - 2 pcs.

VI. Barry anti-bed sore system with accessories, version Barry Serio Forte 2, as part of:

1. Barry Serio Forte compressor 2 -- 1 pc.

2. Barry Serio Forte mattress 2 – 1 pc.

3. Connecting tubes – 1 pc.

4. Instruction manual – 1 pc.

Accessories:

1. Barry Serio Forte 2 repair kit as part of:

- self-adhesive patch - 2 pcs.

VII. Barry anti-bed sore system with accessories, version Barry Mezzo Plus, as part of:

1. Barry Mezzo Plus compressor – 1 pc.

2. Barry Mezzo Plus mattress – 1 pc.

3. Connecting tubes -- 1 pc.

4. Instruction manual – 1 pc.

Accessories:

1. Removable case B – 1 pc.

2. Barry Mezzo Plus repair kit, as part of:

- self-adhesive patch - 2 pcs.:

- spare cylinder B – 1 pc.

VIII. Barry anti-bedsore system with accessories, version Barry Mezzo Forte, as part of:

1. Barry Mezzo Forte compressor – 1 pc.
2. Barry Mezzo Forte mattress – 1 pc.
3. Connecting tubes – 1 pc.
4. Instruction manual – 1 pc.

Accessories:

1. Removable case B – 1 pc.
2. Barry Mezzo Forte repair kit, as part of:
 - self-adhesive patch - 2 pcs.;
 - spare cylinder B - 1 pc.

IX. Barry anti-bedsore system with accessories, version Barry Mezzo Forte 2, as part of:

1. Barry Mezzo Forte compressor 2 – 1 pc.
2. Barry Mezzo Forte mattress 2 – 1 pc.
3. Connecting tubes – 1 pc.
4. Instruction manual – 1 pc.

Accessories:

1. Removable case B – 1 pc.
2. Barry Mezzo Forte 2 repair kit, as part of:
 - self-adhesive patch - 2 pcs.;
 - spare cylinder B - 1 pc.

X. Barry anti-bedsore system with accessories, Barry Mezzo Premium version, as part of:

1. Barry Mezzo Premium compressor – 1 pc.
2. Barry Mezzo Premium mattress – 1 pc.
3. Connecting tubes – 1 pc.
4. Removable 2 m power cable – 1 pc.
5. Instruction manual – 1 pc.

Accessories:

1. Removable MP case – 1 pc.
2. Barry Mezzo Premium repair kit, as part of:

- self-adhesive patch - 2 pcs.;

- spare MP cylinder – 1 pc.

XI. Barry anti-bedsore system with accessories, Barry Mezzo Pro version, as part of:

1. Barry Mezzo Pro compressor – 1 pc.

2. Barry Mezzo Pro mattress – 1 pc.

3. Connecting tubes – 1 pc.

4. Removable 2 m power cable – 1 pc.

5. Instruction manual – 1 pc.

Accessories:

1. Removable MP case – 1 pc.

2. Barry Mezzo Pro repair kit, as part of:

- self-adhesive patch - 2 pcs.;

- spare MP cylinder – 1 pc.

Appearance of the medical product "Barry anti-bedsore system with accessories"	
Barry Primo version	
Barry Mezzo Basic version	

Barry Serio version	
Barry Serio Plus version	
Barry Serio Forte version	
Barry Serio Forte 2 version	

Barry Mezzo Plus version


Barry Mezzo Forte version

Barry Mezzo Forte 2 version

Barry Mezzo Premium version	 
Barry Mezzo Pro version	 

7. Technical specification

The main technical parameters are summarized in the table below:

Name of the parameter	Value
Barry Primo version	
Compressor	
Overall dimensions (L×W×C), mm	240±10×120±10×90±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0,5) Hz
Power consumption, VA, no more	12
The fuse	Fusible, 3.6×10 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. rs. art.	70-130
Weight, kg	0,75±0,1
Cycle duration, min	12±1

Noise level, dBA, no more	33
Air consumption, l/min, no more	6
The Control Panel	The working part type B
Power cable length, mm	2000±40
The mattress	
Mattress pumping time, min., no more	25±1
The shape of the mattress	Cellular
Overall dimensions (L×W×C), mm	2265±20×970±20×70±20
Maximum load, not more than, kg	135
Weight, kg	2,2±0,2
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1
Barry Mezzo Basic version	
Compressor	
Overall dimensions (L×W×C), mm	275±10×125±10×100±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz
Power consumption, VA, no more	12
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. rs. art.	20-120
Weight, kg	0,8±0,1
Cycle duration, min	12±1
Noise level, dBA, no more	35
Air consumption, l/min, no more	6
The Control Panel	The working part of type B
Power cable length, mm	2000±40
The mattress	
Mattress pumping time, min., no more	15±1
The shape of the mattress	Tubular
Overall dimensions (L×W×C), mm	2260±20×820±20×80±20
Maximum load, not more than, kg	145
Weight, kg	4±0,5
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Overall dimensions of the removable cover (L×W), mm	1965±20×935±20
The weight of the removable cover, g	160±5
Overall dimensions of the spare cylinder (L×W), mm	1025±10×170±10
Weight of the spare cylinder, g	135±5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1

Special functions: - Ventilation.	
Barry Serio version	
Compressor	
Overall dimensions (L×W×C), mm	275±10×140±10×90±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz
Power consumption, VA, no more	12
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	50-110
Weight, kg	1,25±0,3
Cycle duration, min	12±1
Noise level, dBA, no more	38
Air consumption, l/min, no more	6
The Control Panel	The working part of type B
Power cable length, mm	1700±40
The mattress	
Mattress pumping time, min., no more	25±1
The shape of the mattress	Cellular
Overall dimensions (L×W×C), mm	2265±20×970±20×70±20
Maximum load, not more than, kg	135
Weight, kg	2±0,2
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1
Special functions: - Ventilation; - Static function.	
Barry Serio Plus version	
Compressor	
Overall dimensions (L×W×C), mm	240±10×120±10×90±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz
Power consumption, VA, no more	12
The fuse	Fusible, 3.6×10 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	70-130
Weight, kg	0,75±0,1
Cycle duration, min	12±1
Noise level, dBA, no more	33
Air consumption, l/min, no more	6
The Control Panel	The working part of type B

Power cable length, mm	2000±40
The mattress	
Mattress pumping time, min., no more	25±1
The shape of the mattress	Cellular
Overall dimensions (L×W×C), mm	2265±20×970±20×70±20
Maximum load, not more than, kg	135
Weight, kg	2±0,2
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1
Special functions: - Ventilation.	
Barry Serio Forte version	
Compressor	
Overall dimensions (L×W×C), mm	220±10×110±10×80±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz
Power consumption, VA, no more	12
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	70-130
Weight, kg	0,85±0,1
Cycle duration, min	12±1
Noise level, dBA, no more	35
Air consumption, l/min, no more	6
The Control Panel	The working part of type B
Length of the power cable, mm	2000±40
The mattress	
Mattress pumping time, min., no more	30±1
The shape of the mattress	Tubular
Overall dimensions (L×W×C), mm	2520±20×820±20×80±20
Maximum load, not more than, kg	130
Weight, kg	3±0,3
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1
Barry Serio Forte 2 version	
Compressor	
Overall dimensions (L×W×C), mm	215±10×110±10×95±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz
Power consumption, VA, no more	12
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA

The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	70-130
Weight, kg	0,75±0,1
Cycle duration, min	12±1
Noise level, dBA, no more	35
Air consumption, l/min, no more	6
The Control Panel	The working part of type B
Length of the power cable, mm	1920±40
The mattress	
Mattress pumping time, min., no more	30±1
The shape of the mattress	Tubular
Overall dimensions (L×W×C), mm	2520±20×820±20×80±20
Maximum load, not more than, kg	130
Weight, kg	3±0,3
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1
Special functions: - Ventilation.	
Barry Mezzo Plus version	
Compressor	
Overall dimensions (L×W×C), mm	215±10×110±10×95±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz
Power consumption, VA, no more	12
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	70-130
Weight, kg	0,75±0,1
Cycle duration, min	12±1
Noise level, dBA, no more	35
Air consumption, l/min, no more	6
The Control Panel	The working part of type B
Length of the power cable, mm	2000±40
The mattress	
Mattress pumping time, min., no more	15±1
The shape of the mattress	Tubular
Overall dimensions (L×W×C), mm	2060±20×820±20×80±20
Maximum load, not more than, kg	135
Weight, kg	4±0,5
Length of the connecting tube, mm	1200±40
The thickness of the connecting tube, mm	2±0,1
The outer diameter of the connecting tube, mm	9±0,5

The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1
Overall dimensions of the removable cover (L×W×C), mm	1965±20×935±20
The weight of the removable cover, g	160±5
Overall dimensions of the spare cylinder (L×W×C), mm	1025±10×170±10
Weight of the spare cylinder, g	135±5
Barry Mezzo Forte version	
Compressor	
Overall dimensions (L×W×C), mm	270±10×130±10×90±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz
Power consumption, VA, no more	12
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3×0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	20-120
Weight, kg	1,25±0,3
Cycle duration, min	12±1
Noise level, dBA, no more	30
Air consumption, l/min, no more	6
The Control Panel	The working part of type B
Length of the power cable, mm	2000±40
The mattress	
Mattress pumping time, min., no more	15±1
The shape of the mattress	Tubular
Overall dimensions (L×W×C), mm	2060±20×820±20×80±20
Maximum load, not more than, kg	145
Weight, kg	4±0,5
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1
Overall dimensions of the removable cover (L×W×C), mm	1965±20×935±20
The weight of the removable cover, g	160±5
Overall dimensions of the spare cylinder (L×W×C), mm	1025±10×170±10
Weight of the spare cylinder, g	135±5
Special functions: - Static mode.	
Barry Mezzo Forte 2 version	
Compressor	
Overall dimensions (L×W×C), mm	275±10×120±10×100±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz

Power consumption, VA, no more	12
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Non-removable LAN-PC-10A-1.5 (Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	20-120
Weight, kg	1,2±0,3
Cycle duration, min	12±1
Noise level, dBA, no more	38
Air consumption, l/min, no more	6
The Control Panel	The working part of type B
Length of the power cable, mm	2000±40
The mattress	
Mattress pumping time, min., no more	15±1
The shape of the mattress	Tubular
Overall dimensions (L×W×C), mm	2060±20×820±20×80±20
Maximum load, not more than, kg	145
Weight, kg	4±0,5
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	160±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1±0,1
Overall dimensions of the removable cover (L×W×C), mm	1965±20×935±20
The weight of the removable cover, g	160±5
Overall dimensions of the spare cylinder (L×W×C), mm	1025±10×170±10
Weight of the spare cylinder, g	135±5
Special functions: - Static mode.	
Barry Mezzo Premium version	
Compressor	
Overall dimensions (L×W×C), mm	270±10×130±10×90±10
Power supply characteristics: power supply voltage, V; mains frequency, Hz	220(±22) V; 50(±0.5) Hz
Power consumption, VA, no more	15
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Removable LAN-PPM-10A-2.0 (C13-Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	14-62
Weight, kg	1,5±0,3
Cycle duration, min	14±1/25±1
Noise level, dBA, no more	30
Air consumption, l/min, no more	6
The Control Panel	The working part of type B
	2000±40
The mattress	

Mattress pumping time, min., no more	30±1
The shape of the mattress	Tubular
Overall dimensions (L×W×C), mm	2080±20×840±20×110±20
Maximum load, not more than, kg	145
Weight, kg	4,58±0,5
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	1,75±0,1
The outer diameter of the connecting tube, mm	8,5±0,5
The inner diameter of the connecting tube, mm	5±0,5
Weight of the connecting tube, g	225±5
The diameter of the self-adhesive patch, mm	69±1
Weight of the self-adhesive patch, g	1,5±0,1
Overall dimensions of the removable cover (L×W×C), mm	2045±20×915±20
The weight of the removable cover, g	650±5
Overall dimensions of the spare cylinder (L×W×C), mm	935±10×345±10
Weight of the spare cylinder, g	145±5
Special functions:	
- Emergency air vent valve;	
- Static function;	
- Notification;	
- Seat mode;	
- Lock the control panel.	
Barry Mezzo Pro version	
Compressor	
Overall dimensions (L×W×C), mm	345±10×115±10×205±10
Power supply characteristics:	220(±22) V; 50(±0.5) Hz
power supply voltage, V; mains frequency, Hz	
Power consumption, VA, no more	30
The fuse	Fusible, 5.2×20 mm; 250 V; 500mA
The degree of protection of the IP shell	IP 21
The cross-sectional area of the power cable wires, mm ²	0,75
Type of power cord	Removable LAN-PPM-10A-2.0 (C13-Schuko, 3x0.75, 220V, 10A, black)
Pressure, mm. Rs. St.	10-56/30-80
Weight, kg	2,5±0,3
Cycle duration, min	6±1/9±1/12±1/25±1
Noise level, dBA, no more	35
Air consumption, l/min, no more	7-14
The Control Panel	The working part of type B
Length of the power cable, mm	2000±40
The mattress	
Mattress pumping time, min., no more	25±1
The shape of the mattress	Трубчатая
Overall dimensions (L×W×C), mm	2000±20×860±20×130±20
Maximum load, not more than, kg	210
Weight, kg	8,35±0,5
Length of the connecting tube, mm	1500±40
The thickness of the connecting tube, mm	2,75±0,1
The outer diameter of the connecting tube, mm	12,5±0,5
The inner diameter of the connecting tube, mm	7±0,5
Weight of the connecting tube, g	225±5
The diameter of the self-adhesive patch, mm	69±1

Weight of the self-adhesive patch, g	1,5±0,1
Overall dimensions of the removable cover (L×W×C), mm	2045±20×915±20
The weight of the removable cover, g	650±5
Overall dimensions of the spare cylinder (L×W×C), mm	935±10×345±10
Weight of the spare cylinder, g	145±5
Special functions: - Static mode; - Seat mode; - Display with weight and pressure adjustment parameters; - Alerts; - Lock the control panel.	

8. Medical device labeling

A brief designation of MI versions is presented in the table below.

Name of MI version	Designation
Barry Primo version	P
Barry Mezzo Basic version	MB
Barry Serio version	S
Barry Serio Plus version	SP
Barry Serio Forte version	SF
Barry Serio Forte 2 version	SF2
Barry Mezzo Plus version	MP
Barry Mezzo Forte version	MF
Barry Mezzo Forte 2 version	MF2
Barry Mezzo Premium version	MPM
Barry Mezzo Pro version	MPR


The marking is typographic or printed on a paper label, indicating:

- product name;
- product model designation;
- name and (or) trademark of the manufacturer;
- address of the manufacturing company;
- serial number or lot;
- date of manufacture;
- symbol "Refer to the instructions for use";
- registration certificate number;
- symbols for electrical safety according to GOST R IEC 60601-1-2022;
- power supply; power consumption; duty cycle; air volume;
- sign "Type B working part";
- protection classification against dangerous ingress of water or solid particles IP21;
- bar code;
- "Special Disposal" sign.

The marking is applied in the form of a plate on the mattress of the anti-decubitus system in the lower left corner. The product labeling is as follows:

Barry anti-bedsore system with accessories

T P O

Version: Barry _____					
	_____		____/20____		
	Valentine International Ltd., 8F, No. 149, Sec. 2, Ta Tung Road, Hsichih Dist., New Taipei City 221, Taiwan, R.O.C.				
Authorized manufacturer's representative: OOO SIMS-2, 125430, Moscow, 16, Mitinskaya St., est. 10, pom. 1012B, room 15 to 18 tel.: +7 (495) 792-31-90 info@sims2.ru 					
Registration certificate № _____				 2 000988 666470	

The compressor marking is applied to the back of the compressor casing. The compressor markings are as follows:

Compressor for Barry anti-bedsore system with accessories, Barry _____			
Power: _____			____/20____
Power consumption: _____		  IP21 	
Operating cycle: _____			
Air Volume: _____			
	Valentine International Ltd., 8F, No. 149, Sec. 2, Ta Tung Road, Hsichih Dist., New Taipei City 221, Taiwan, R.O.C.		

If it is not possible to apply the marking to the product body, the marking is applied to the product packaging or inserted into the packaging together with the product.

Table of symbols used in marking:

Graphic symbol	Value Description
	Manufacturer's trademark
Compressor for Barry anti-bedsore system with accessories	Name or designation of the product type (kind);
Performance Option:	Designation of MI version
	Identification (serial/factory) number
	The date of manufacture of the product (month / year in xx/xxxxx format);

	Manufacturer information: Name of the manufacturer Manufacturer's address
 7 1000988 666470	The barcode is used to identify information about the manufacturing company and information about the product, place of production
	Special utilization
Registration Certificate No. XXXXXX/XXXXXX	Number of registration certificate MI
	Refer to the operating instructions
	Type B working part
IP21	Classification of protection against dangerous penetration of water or solid particles

9. Packaging, transportation, labeling of medical device packaging

The following conditions must be observed during storage or transportation to maintain the performance of the bed:

Transportation Conditions:

Temperature: -40°C to +70°C

Humidity: 10% to 95% (non-condensing)

Pressure: 50 kPa to 106 kPa

Storage Conditions:

Temperature: -40° C to +40° C

Humidity: up to 75% (non-condensing)

Pressure: 50 kPa to 106 kPa

Packaging is performed in accordance with the requirements of GOST R 50444, GOST R 51474 and technical documentation of the Manufacturer. The medical device and all components of the composition are packed in transparent polyethylene film HDPE (LDPE), density 925 kg/m³, thickness 0.5 mm, manufactured by "Nan Ya Plastic Corporation". Taiwan or similar. Then the product is packed in a box made of five-layer cardboard 7 mm thick, density 150 g/m³, color brown, manufactured by "Shenzhen Xing Sheng Li Imports and Exports Co., Ltd." (China).

The dimensions of the packaging containers for the versions of the medical product "Barry anti-bedsore system with accessories" are presented in the table below:

Product version	Dimensions of packing tare, mm	Gross weight, kg	Net weight, kg
Barry Primo version	310±50×160±50×220±50	3,5±0,5	3,8±0,5
Barry Mezzo Basic version	360±50×295±50×275±50	5,5±0,5	5,8±0,5
Barry Serio version	310±50×160±50×220±50	3,8±0,5	4,1±0,5
Barry Serio Plus version	310±50×160±50×220±50	3,4±0,5	3,7±0,5
Barry Serio Forte version	310±50×160±50×220±50	4,5±0,5	4,8±0,5
Barry Serio Forte 2 version	310±50×160±50×220±50	4,4±0,5	4,7±0,5

T P O

Barry Mezzo Plus version	340±50×290±50×250±50	5,4±0,5	5,7±0,5
Barry Mezzo Forte version	340±50×290±50×250±50	5,8±0,5	6,1±0,5
Barry Mezzo Forte 2 version	340±50×290±50×250±50	6±0,5	6,3±0,5
Barry Mezzo Premium version	370±50×295±50×345±50	6,9±0,5	7,2±0,5
Barry Mezzo Pro version	655±50×255±50×440±50	11,7±0,5	12,2±0,5

Labeling of the product packaging:

Transportation markings include:

- product model designation;
- name and (or) trademark of the manufacturer;
- address of the manufacturing company;
- date of manufacture;
- overall dimensions of the transportation package (L×W×H), mm;
- gross weight, kg;
- net weight, kg;
- bar code;
- "Special Disposal" sign;
- number of registration certificate MI;
- information on the manufacturer's authorized representative (name, address);

Labeling of transportation containers should be made with the indication of manipulation signs "Top", "Keep away from moisture", "Fragile. Caution."

The operating instructions must be enclosed in each consumer package.

Other symbols according to ISO 15223-1 may be used.

		Barry anti-bedsore system with accessories	
Version:		Size (DxWxH), mm:	
	___/20___		___x___x___x___ Gross weight, kg: ___ Net weight, kg: ___
	Valentine International Ltd., 8F, No. 149, Sec. 2, Ta Tung Road, Hsichih Dist., New Taipei City 221, Taiwan, R.O.C.		
		Authorized representative of the manufacturer: LLC "SIMS-2", 125430, Moscow, Mitinskaya street, 16, est. 10, room 1012B, rooms 15 to 18 tel.: +7 (495) 792-31-90 info@sims2.ru	
		Registration certificate № _____	

10. Disposal

At the end of the service life the mattress, compressor and accessories should be disposed of. On the territory of the Russian Federation disposal is carried out in accordance with SanPiN 2.1.3684-21 "Sanitary and Epidemiological Requirements for the maintenance of urban and rural settlements, water bodies, drinking water and drinking water supply, air, soil, living quarters, operation of industrial and public premises, organization and conduct of sanitary and anti-epidemic (preventive) measures". This medical device belongs to class A - epidemiologically safe waste, by composition to solid domestic waste.

11. List of regulatory documents

GOST R 50444-2020 "Medical devices, apparatuses and equipment. General technical requirements".

GOST R IEC 60601-1-2022 "Medical electrical products. Part 1. General safety requirements taking into account the basic functional characteristics".

GOST R IEC 60601-1-2-2014 "Medical electrical products. Part 1-2. General safety requirements, taking into account the main functional characteristics. A parallel standard. Electromagnetic compatibility. Requirements and tests".

GOST R IEC 60601-1-6-2014 "Medical electrical products. Part 1-6. General safety requirements, taking into account the main functional characteristics. An additional standard. Operational suitability".

GOST R IEC 62366-1-2023 Medical devices. Part 1. Design of medical devices taking into account their operational suitability.

GOST ISO 14971-2021 "Medical devices. Application of risk management to medical devices".

GOST R 57769-2021 "Anti-bedsores mattresses and pillows. Types and basic parameters".

GOST R ISO 2439-93 "Elastic cellular polymer materials. Determination of indentation hardness".

GOST R ISO 15223-1-2023 "Medical devices. Symbols used for labeling medical devices, on labels and in accompanying documentation. Part 1. Basic requirements".

GOST 31214-2016 "Medical products. Requirements for samples and documentation submitted for toxicological, sanitary and chemical studies, sterility and pyrogenicity tests".

GOST R 52770-2023 "Medical products. Safety requirements. Methods of sanitary, chemical and toxicological tests".

GOST ISO 10993-1-2021 "Medical devices. Assessment of the biological effect of medical devices. Part 1. Assessment and research".

GOST R ISO 10993-2-2009. "Medical devices. Assessment of the biological effect of medical devices. Part 2. Requirements for the treatment of animals".

GOST ISO 10993-5-2023 "Medical products. Assessment of the biological effect of medical devices. Part 5. Cytotoxicity testing: in vitro methods".

GOST ISO 10993-10-2023 "Medical products. Assessment of the biological effect of medical devices. Part 10. Studies of irritating and sensitizing effects".

GOST ISO 10993-12-2023 "Medical products. Assessment of the biological effect of medical devices. Part 12. Sample preparation and control samples".

GOST 31209-2003 "Containers for blood and its components. Chemical and biological safety requirements and test methods".

GOST 4011-72 "Drinking water. Methods for measuring the mass concentration of total iron".

GOST 31870-2012 "Drinking water. Determination of the content of elements by atomic spectrometry methods" item 4.

MUK 4.1.3166-14 "Gas chromatographic determination of hexane, heptane, acetaldehyde, acetone, methyl acetate, ethyl acetate, methanol, isopropanol, acrylonitrile, n-propanol, n-propyl acetate, butyl acetate, isobutanol, n-butanol, benzene, toluene, ethylbenzene, m-, o- and p-xylenes, isopropylbenzene, styrene, α -methylstyrene in water and aqueous extracts from materials of various compositions."

MUK 4.1.3086-13 "Gas chromatographic determination of hexamethylenediamine in aqueous extracts from polymeric materials used in the food industry".

GOST R 55227-2012 "Water. Methods for determining the content of formaldehyde" item 6.

MP 1941-78 "Methodological recommendations for the determination of vinyl chloride in polyvinyl chloride and polymer materials based on it, in model media simulating food products, in food products".

MUC 4.1.3169-14 "Gas chromatographic determination of dimethyl phthalate, dimethyl terephthalate, diethyl phthalate, dibutyl phthalate, butylbenzyl phthalate, bis(2-ethylhexyl) phthalate and dioctyl phthalate in water and aqueous extracts from materials of various compositions".

12. Operation

MI Operating Conditions:

Temperature: 10°C to 40°C

Humidity: 30% to 75%

Pressure: 70 kPa to 106 kPa

Service life of the medical device: 5 years.

Operation mode of the medical device: continuous operation mode.

12.1 Procedure for use of the medical device

- Unpack the product
- Secure the compressor to the foot of the bed with brackets, or place it on a vibration-resistant horizontal surface next to the patient's bed at the foot of the bed.
- Place the system mattress on top of the bed mattress. The area of the bed mattress should not be smaller than the system mattress.
- Tuck the loose edges of the system mattress under the bed mattress
- Put a protective cover on the system mattress
- Connect the system mattress to the compressor by means of the air connection tubes
- Connect the compressor to a serviceable power supply.

Turn on the compressor. Air will start to flow into the mattress. The red indicator on the compressor will light up. When the pressure in the mattress reaches the required level, the color

of the indicator will change to green. The complete filling of the mattress with air depends on the version of the product.

- The patient should be placed on the mattress after the mattress is completely filled with air.
- After the patient is placed on the mattress, the compressor capacity should be reduced to 60-70% of maximum. (discharge capacity depends on the weight of the patient)
- Adjust the compressor power so that two fingers can pass freely between the inflated mattress cylinders and the patient's body. The compressor power should be monitored daily, adjusting if necessary.

Use the dedicated button to select the compressor cycle time when this function is available.

13. Sterility, cleaning and disinfection

The product and accessories are supplied non-sterile and cannot be sterilized afterwards.

13.1 Cleaning recommendations

Compressor:

Do not immerse the compressor in water;

Check the product for external damage and then proceed to cleaning;

Wipe the compressor surface with a cloth dampened with 25% ammonium hydroxide solution.

Do not spray the disinfection solution directly onto the surface of the product. If the solution penetrates into the compressor housing, the compressor will be damaged;

Allow surfaces to dry after cleaning;

Do not use hypochlorite or phenol-based cleaning solutions, as this may damage the product case.

Mattress:

The outer surfaces of mattresses and covers must be resistant to treatment by double wiping with a calico or gauze cloth soaked in a solution of disinfectant (type "Lotus") according to MU 287-113.

Any obvious stains should be thoroughly soaked in a 1:9 hypochlorite solution (1 part bleach to 9 parts water) and allowed to dry for at least 10 minutes. Then wipe with a clean, damp cloth. Allow the mattress to dry completely. Once the mattress is dry, wipe both sides with a disinfectant (Lotus type). After use, store the mattress in a dry place.

The cover can be washed in warm water (temperature up to 40 degrees Celsius) with mild detergents, wrung out without twisting, dried away from heating devices and heating batteries.

14. Maintenance

Maintenance of the bed is aimed at controlling the technical and operational condition of the bed, in order to avoid unplanned breakdowns, including those related to possible harm to human health. Routine maintenance should be carried out annually by qualified technical personnel of the medical institution or by specialists of an authorized service center.

In addition to regular comprehensive inspections by qualified engineering personnel, routine users (medical institutions) are also required to perform daily visual inspections of the products in use and functional checks prior to the use of the medical device by a new patient.

Possible faults and remedies:

Problem	Problem solving
Compressor does not work	1. Check whether the plug is plugged into the socket outlet 2. Press the on/off button. - if the button is on "ON"/"ON" and the compressor still does not work, contact the service center: - If the button is set to OFF, check the outlet. Connect the plug to another outlet. If the button is still OFF, contact an electrical engineer.
Insufficient inflation (low pressure)	1. For a quick check, switch the pressure to the high stiffness position; 2. Check the connection of the tubes to each other and to the compressor: - the tubes must always be straightened; - replacement of the tubes is necessary in case of air leakage; - make sure the valves are closed.
Weak air flow	A clogged compressor can cause poor air flow. Wash the filter with a mild detergent. Check the air filter once a month.

15. Warranty obligations

The manufacturer guarantees the compliance of the medical product "Barry anti-bedsore system with accessories" with the requirements set forth in the technical and operational documentation, provided that the consumer complies with the conditions of use, transportation and storage.

The warranty period is 12 months from the date of sale.

Guaranteed shelf life is 3 months.

Non-observance of the operating instructions and repair work carried out by unqualified persons, especially technical modifications and retrofitting without the manufacturer's consent may invalidate the warranty.

16. Contact information

In case of questions related to the use of the product, as well as in case of claims to the manufacturer, the consumer can contact the following addresses:

Authorized representative of the manufacturer in the Russian Federation on quality issues:

Limited Liability Company "SIMS-2" (LLC "SIMS-2"). Address: 16, Mitinskaya St., Moscow, 125430, Russia, est. 10, room 1012B, comp. 15 to 18.

Tel: +7 (495)-792-31-90

E-mail address: info@sims2.ru

Website: sims2.ru

Taiwan:

Valentine International Ltd.

Address: 8F, No. 149, Sec. 2, Ta Tung Road, Hsichih Dist., New Taipei City 221, Taiwan, R.O.C.

tel: +886.2.8692.6662

e-mail: info@valentine-dynalife.com

website: <http://www.valentine-dynalife.com/>

Annex 1: Electromagnetic compatibility

This product meets the following EMC requirements regarding electromagnetic emission and immunity:

The Barry anti-bedsore system with accessories is intended for use in the electromagnetic environment defined below. The purchaser or user of the product should ensure that it is used in the specified environment.			
Electromagnetic emission test	Compliance	Electromagnetic environment - guidelines	
Radio interference by SISPR 11	Group 1	The products use RF energy only for internal functions. RF emissions are low and are not likely to cause malfunctions, located near electronic equipment.	
Radio interference by SISPR 11	Class B	The product is suitable for use in any location, including residences and buildings directly connected to the electrical distribution network supplying residences.	
Harmonic current components according to IEC 61000-3-2	Class A		
Voltage fluctuations and flicker according to IEC 61000-3-3	Relevant		
Manual and manufacturer's declaration - immunity for non-LIFE MI.			
The product is intended for use in the electromagnetic environment defined below. The purchaser or user of the product should ensure that it is used in the specified environment.			
Immunity test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharges (ESD) according to IEC 61000-4-2	±6 kV - contact discharge ±8 kV - air discharge	±6 kV - contact discharge ±8 kV - air discharge	Floors should be wood concrete or tiled. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
Nanosecond pulse interference according to IEC 61000-4-4	±2 kV - for power supply lines ±1 kV - for I/O lines	±2 kV - for power supply lines ±1 kV - for I/O lines	The quality of electrical energy in the building's electrical grid should be consistent with typical commercial or hospital environments.
High energy microsecond pulse interference to IEC 61000-4-5	±1 kV - when interference is applied by the "wire-to-wire" scheme ± 2kV - when interference is applied by the "wire-to-ground" scheme	±1 kV - when interference is applied by the "wire-to-wire" scheme ±2 kV - when interference is applied by the "wire-to-ground" scheme	The quality of electrical energy in the building's electrical grid should be consistent with typical commercial or hospital environments.
Dips, interruptions and voltage changes in input power lines to IEC 61000-4-11	<5% U_T^* (voltage dip >95% U_T) for 0.5 period; 40% U_T (voltage dip 60% U_T) for five periods; 70% U_T (voltage dip 30% U_T) for 25 periods; <5% U_T (voltage dip >95% U_T) for 5 s	<5% U_T (voltage dip >95% U_T) for 0.5 period; 40% U_T (voltage dip 60% U_T) for five periods; 70% U_T (voltage dip 30% U_T) for 25 periods;	The quality of electrical energy in the building's electrical grid should be consistent with typical commercial or hospital environments. If the user needs to ensure continuous operation under

INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

		<5% U_T (voltage dip >95% U_T) for 5 s	conditions of possible interruptions of the mains voltage, it is recommended to power the MI from an uninterruptible power supply or battery.
Note - U_T is the voltage level of the electrical network before the test impact is applied.			
Magnetic field of industrial frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m.	3 A/m.	Industrial frequency magnetic field level should be provided in accordance with typical commercial or hospital environments.
Conductive disturbances induced by radio frequency electromagnetic fields according to IEC 61000-4-6	3 V (RMS value) in the band from 150 kHz to 80 MHz (RMS value) in the band from 150 kHz to 80 MHz outside the frequency bands allocated for industrial, scientific and medical high-frequency (HF HF) devices ^{a)}	3 V ($[V_I]$, V) (RMS value) in the band from 150 kHz to 80 MHz	The distance between the mobile radiotelephone systems in use and any element of the MI including cables, shall not be less than the recommended spatial separation, which is calculated in accordance with the following expression for the transmitter frequency Recommended spatial separation: calculated according to the formula for the frequency band from 150 kHz to 80 MHz: $d = [3.5/V_I] \sqrt{P}$ $d = 1.2 \sqrt{P}$ for the frequency band from 80 MHz to 800 MHz: $d = [3.5/E_I] \sqrt{P}$ $d = 1.2 \sqrt{P}$ for the frequency band from 800 MHz to 2.5 GHz: $d = \left[\frac{7}{E_I} \right] \sqrt{P}$ $d = 2.3 \sqrt{P}$ where d - recommended spatial spacing, m ^{b)} P - nominal maximum output power of the
Emitted radio frequency electromagnetic field according to IEC 61000-4-3	3 V/m in the band from 80 MHz to 2.5 GHz	3 V/m [E_I], V/m in the band from 80 MHz to 2.5 GHz	

			transmitter, W, set by the manufacturer. The field strength during radio wave propagation from stationary radio transmitters, according to electromagnetic observations^{a)}, shall be lower than the compliance level in each frequency band^{b)}. Interference may occur in the vicinity of equipment marked with the symbol: 
--	--	--	---

^{a)} The field strength for radio wave propagation from fixed radio transmitters such as base stations of radiotelephone networks (cellular/wireless) and land mobile radio stations, amateur radio stations, AM and FM broadcasting transmitters, television transmitters cannot be determined by calculation with sufficient accuracy. For this purpose practical field strength measurements must be carried out. If the measured values at the location of the [ME DEVICE or ME SYSTEM] exceed the applicable compliance levels, observation of the [ME DEVICE or ME SYSTEM] shall be made to verify its normal operation. If the observation process reveals a deviation from normal operation, additional measures such as reorienting or relocating the [ME DEVICE or ME SYSTEM] may be necessary.

^{b)} Outside the 150 kHz to 80 MHz band, the field strength should be less than V_I , V/m.

Notes:

- 1) At 80 and 800 MHz, the higher field strength value is applied.
- 2) Expressions are not applicable in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people.

Recommended values of spatial separation between portable and mobile RF communications and MI.

The medical device is intended for use in an electromagnetic environment where the levels of radiated interference are controlled. The purchaser or user of the MI can avoid the effects of electromagnetic interference by ensuring the minimum spatial separation between portable and mobile radio frequency communication equipment (transmitters) and the MI, as recommended below, taking into account the maximum output power of the communication equipment.

Rated maximum output power of the transmitter P , W	Spatial separation d , m, depending on the transmitter frequency		
	$d=1.2\sqrt{P}$ in the band from 150 kHz to 80 MHz	$d=1.2\sqrt{P}$ in the band from 80 MHz to 800 MHz	$d=2.3\sqrt{P}$ in the band from 800 MHz to 2.5 GHz
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,20	2,30
10	3,80	3,80	7,30
100	12,00	12,00	23,00

When determining the recommended values of spatial separation for transmitters with rated maximum output power not specified in the table, the rated maximum output power in watts specified in the transmitter manufacturer's documentation shall be substituted into the above expressions.

Notes:

- 1) At 80 and 800 MHz, the higher field strength value is applied.
- 2) The above expressions are not applicable in all cases. The propagation of electromagnetic waves is affected by structures, objects and people.

(3) When determining the recommended values of spatial separation for transmitters with a nominal maximum output power not specified in the table, the nominal maximum output power in watts specified in the transmitter manufacturer's documentation shall be substituted into the above expressions.

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