
Laerdal

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
© 2012, Laerdal Medical AS. All rights reserved.

20-04466 Rev E

I. The BAG Disposable Resuscitator configurations and parts.

Instructions for Use

The BAG Disposable Resuscitator is a self-inflating, manual resuscitator that is intended for patients requiring total or intermittent ventilation support. The BAG Disposable Resuscitator provides positive pressure ventilation and also allows for spontaneous breathing either with an artificial airway or with a mask (3)(a), (b) or (c).

Contraindications:

Contraindications for using the BAG Disposable Resuscitator practically are not present. Contraindications for using the BAG Disposable Resuscitator for adult on infants and children, because in this model there is no pressure-limiting device. Also contraindications for using the BAG Disposable Resuscitator are dependent on situations where there are contraindications for using method ALV, for example, obstruction to venous return (where there are contraindications for using method ALV, which more destroy it); lung injury where there are contraindications for using method ALV on principle injection with intermittent positive pressure breathing); the presence of foreign object (solid or liquid) in the upper section of trachea or larynx; and so on.

Precautions and Warnings

The BAG Disposable Resuscitator should only be used by persons who have received specific training in their use.

When using supplemental oxygen, do not allow smoking or use unit near sparking equipment, open flame, oil or other flammable chemicals.

The BAG Disposable Resuscitator should not be used in toxic or hazardous atmospheres.

Infant and Child units are equipped with a pressure-limiting device (2) which opens at a pressure of approx. 35 ± 5 cmH₂O. However, an abrupt, high-volume manual breath may cause the unit to exceed this level. The pressure relief valve can be overridden by depressing the plunger momentarily with light finger pressure or locked in place by using the push-twist locking feature.

The resuscitator is for single patient use only. Do not reuse. Do not sterilize.

Select the appropriate resuscitator size (8)(a), (b) or (c). The wrong size can result in inadequate or excessive air pressure being delivered to the patient. See the Performance and Specifications table for sizing.

When using supplemental oxygen, the flow from its source should be monitored. Delivery of supplemental oxygen greater than 30 LPM (Liters Per Minute) may result in inadvertent positive end expiratory pressure (PEEP).

use of third party products and oxygen delivery devices (e.g. filters and demand valves) with The BAG Disposable Resuscitator may have an affect on product performance. Please consult with the manufacturer of the third party device to verify compatibility with The BAG Disposable Resuscitator and obtain information on the possible performance changes.

Warranty

BAG Disposable Resuscitator is warranted against defects in workmanship and materials only. Please refer to the Global Warranty statement for additional terms and conditions (www.laerdal.com).

Initial Operation

The product is delivered in first-aid protective kits, completely suitable for use, additional sanitary processing isn't required.

Remove the resuscitator from the outer protective poly bag. Expand the Adult or Child resuscitator from collapsed configuration to its operating position.

Inspect the unit to be sure the system is complete.

PRE-OPERATIVE FUNCTIONAL TEST:

Compress the ventilation bag (8) with one hand then release the grip on the bag. Rapid bag re-expansion confirms efficient air intake.

Block the patient valve/mask connector part and try to compress the bag. If the bag cannot be compressed with reasonable force, the valve is efficiently preventing backward escape of air.

Place a reservoir bag, (6) or (7), or test lung (if available) over the patient valve. Compress the bag several times. This should fill the reservoir or test lung and confirm that the patient valve is able to efficiently direct air to the patient.

Compress the filled reservoir bag. Air should vent to the atmosphere as indicated by lifting of the disk membrane at the base of the mask connector and not return to the ventilation bag (8).

PRESSURE RELIEF REGULATOR: The Infant and Child resuscitators feature a patient valve with a special pressure limiting device (2) mounted on the upper valve housing. If inspiration meets with pulmonary resistance of approximately 35 ± 5 cmH₂O, the device opens, reducing the risk of stomach distension and / or barotrauma. A hissing sound can be heard when the device opens. If higher ventilation pressure is required, the pressure limiting device can be overridden with finger pressure on the plunger or disabled by depressing and turning the plunger.

If high concentrations of oxygen are needed, attach oxygen tubing (4) to the bag and an adjustable oxygen source, and attach an oxygen reservoir bag (6) or (7) to the bag.

Adjust oxygen flow to insure the oxygen reservoir bag remains fully or partially inflated during use.

When using a mask (3) attached to the resuscitator, be sure you have a tight seal and secure fit.

When using an endotracheal tube or tracheotomy tube, remove mask (3) and attach the resuscitation unit directly to the tube. The mask elbow provides a 15 mm I.D. port for this purpose.

HAND STRAP: Grab the bag with your hand. Adjust the hand strap (9) by loosening the loop from the hook then pull the strap to fit your hand, and finally reattach the loop to the hook.

Operation Instructions

If supplemental oxygen is needed, connect the oxygen source to the resuscitation bag oxygen connector.

Open patient's airway.

Clear patients mouth of foreign matter.

Apply mask firmly to the face. If the patient is intubated, attach the patient valve connector directly to the endotracheal tube. Squeeze and release the bag allowing enough time between inspirations for the patient to exhale and the bag to re-expand. Follow local protocol.

Observe the rise and fall of the patient's chest and listen for the air flow from the valve as the patient exhales.

IMPORTANT: If the patient's chest does not rise and fall with each breath or no airflow is present, the patient's airway or the patient valve itself may be blocked.

WARNING: TAKE IMMEDIATE ACTION by performing artificial respiration (mouth-to-mouth, mouth-to-barrier, mouth-to-tube resuscitation) or follow local protocol. The airway must be cleared before proceeding.

Performance and Specifications¹

SUSCITATOR MODEL	ADULT (>20 KG)	CHILD (10-20 KG)	INFANT (2.5-12 KG)
MAXIMUM RATE	85 BPM	144 BPM	180 BPM
DELIVERED TIDAL VOLUME PER HAND	830 ml	330 ml	180 ml
PRESSURE RELIEF VALVE	Not provided	(Nominal) 35±5cm H ₂ O	(Nominal) 35±5 cm H ₂ O
OXYGEN CONCENTRATION % FLOWRATE (LPM) FREQUENCY (BPM) % O ₂ W/RESERVOIR % O ₂ W/O RESERVOIR	3 5 10 12 12 12 500 500 500 90% 98% 99% 34 % 47% 66%	10 20 250 100% 70%	4 30 40 100% 85%
MAXIMUM MEASURED VOLUME BAG	1 650 ml 2 900 ml	500 ml 2 900 ml	230 ml 810 ml
SIZES	470g	420g	360g
WEIGHT	60 cm	50 cm	45 cm
Overall dimensions			
Compressible bag	length – 26 cm. height – 13 cm.	length – 15 cm. height – 8,5 cm.	length – 13 cm. height – 7 cm.
Facial Mask	length – 15 cm. width – 12 cm. height – 9 cm.	length – 10 cm. width – 8,5 cm. height – 6 cm.	length – 7 cm. width – 6 cm. height – 4,5 cm.
Timing for fixing the wrist with Velcro	length – 45 cm. width – 2,5 cm.	length – 29 cm. width – 2,5 cm.	length – 24 cm. width – 2,5 cm.
Oxygen Tank	length – 31 cm. width – 24 cm.	length – 22 cm. width – 15 cm.	length – 16 cm. width – 19,5 cm.
operating mode setting time	ready for operation		
Parameters holes for attachment to the patient	outer diameter – 22 mm inner diameter – 15 mm The length of the valve (embedment depth) 15 mm Overall length – 20 mm <i>Connector is unified for all models</i>		
Parameters of connectors facial masks <i>Adults worn Children are inserted)</i>	outer diameter - 26 mm inner diameter - 22 mm the length of the connector - 21 mm	outer diameter - 15 mm inner diameter - 13 mm the length of the connector - 15 mm	

Resuscitator Settings valve bag	outer diameter – 25 mm inner diameter – 22 mm The length of the valve (embedment depth) 10 mm		
Resuscitator Settings exhaust valve bag	diameter 19 mm		
Functioning of the patient after contamination	Double-lobe (bicuspid) non-return (non-reversible), the valve can be easily cleaned, if it is dirty vomit, protects a person which assisting, from vomiting, blood and secretions of the patient. If a non-return valve is clogged breathing vomiting blood or secret during mechanical ventilation, intensive care it is necessary to remove the breathing bag and mask away from the patient and other people and to send him to the ground and sharp and strong compression to remove foreign objects from the non-return valve breathing. It is necessary to repeat this procedure until the complete purification of the valve.		
Dimensions of filling the valve breathing bag	breathable, diameter - 18 mm		
Ability to pass through the hole 600 mm	yes	yes	yes
Features and / or dimensions of input gas nozzle	absent		
Sterility	not sterile		
Multiplicity of application	single		

Tested according to ISO 10651 - 4:2002 and ASTM-F920-93 product standards.

Technical Specifications:

Operating Environmental Temperature Limits: -18°C to +50°C

Storage Environmental Temperature Limits: -40°C to +60°C

Respiratory Resistance: 1.8 cmH₂O @ 50 LPM

Respiratory Resistance: 1.5 cmH₂O @ 50 LPM

Patient Valve Dead Space: 6.8 ml

Materials List:

Parts	Material
Facial Mask	Polyvinylchloride (PVC) SG-3, unpainted
Pressure valve	Silicone Rubber (SI) Mold Max 15T, unpainted Polycarbonate (PC) Lexan 1278R, unpainted
Compressible bag	Polyvinylchloride (PVC) SG-3, Cromofix 01B019 yellow
Pressure limiter valve	Polycarbonate (PC) Lexan 1278R, unpainted
Oxygen Tank	Polyvinylchloride (PVC) SG-3, unpainted
Oxygen tube	Polyvinylchloride (PVC) SG-3, Cromofix 06B017 green
Timing for fixing the wrist with Velcro	Nylon 150D

Period of storage:

The period of storage of the disposable Medical ventilator of manual Bag Disposable Resuscitator, at admissible climatic conditions, and temperature which isn't exceeding + 40° With shouldn't exceed 5 years.

ation

life, or in the case of The BAG Disposable Resuscitator should be disposed of in accordance with applicable national legislation. Waste collection provides distributor (contractual service organization or contractual service technicians).

Transportation

BAG Disposable Resuscitator transported by all modes of transport in covered vehicles in accordance with the rules for transport at temperatures from -40 to +60 C and humidity of 60 to 95%.

Labels:

REF

Unique product reference.

LOT

Manufacturing batch code.

Do not reuse.

Warranties on quality products are accepted to:

Warranty liability Company «REEPL» (Ltd REEPL) 125130, Moscow, Staropetrovsky passage, House 7A, Building 3.

Package Includes:

BAG Disposable Resuscitator, comes flat packed in individual packaging.

The BAG Disposable Resuscitator for adults comprising:

- Oxygen Tank – 1pc.
- Pressure valve – 1pc.
- Facial Mask – 1 pc.
- Timing for fixing the wrist with Velcro – 1pc.
- Compressible bag – 1pc.
- Oxygen tube – 1pc.
- Operating Manual – 1pc.

II. The BAG Disposable Resuscitator for child comprising:

- Oxygen Tank – 1pc.
- Pressure valve – 1pc.
- Pressure limiter valve – 1pc.
- Facial Mask – 1 pc.
- Timing for fixing the wrist with Velcro – 1pc.
- Compressible bag – 1pc.
- Oxygen tube – 1pc.
- Operating Manual – 1pc.

III. The BAG Disposable Resuscitator for infant comprising:

- 1. Oxygen Tank – 1pc.
- 2. Pressure valve – 1pc.
- 3. Pressure limiter valve – 1pc.
- 4. Facial Mask – 1 pc.
- 5. Timing for fixing the wrist with Velcro – 1pc.
- 6. Compressible bag – 1pc.
- 7. Oxygen tube – 1pc.
- 8. Operating Manual – 1pc.

LAERDAL MEDICAL AS

William Clive Patrickson
Chief Executive Officer

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