002202410150202

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WYF2412692

241012 183 00008 1/2

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240411389-001 1/2 ⑦

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2024/10/17

002

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CERTIFICATE

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China Council for the Promotion of International Trade

China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of FUJIAN PUTIAN LICHENG PAPER INDUSTRY CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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

证 查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

APPROVED BY
Legal representative
Fujian Putian Licheng Paper Industry Co., Ltd

L.S. Huang Limei

Date: August 28, 2024

**INSTRUCTIONS FOR USE
for the medical device:**

Disposable absorbent medical products Dr.Skipp®

2024

The instructions for use below are drawn up to provide information for the purpose of state registration of the medical device in the territory of the Russian Federation.

Information is provided to the extent necessary in accordance with the Rules for State Registration of a Medical Device (Decree of the Government of the Russian Federation No. 1416 dated December 27, 2012), as well as in accordance with Order of the Ministry of Health of the Russian Federation No. 11n dated January 19, 2017 On Approval of Requirements for the Content of Technical and Operational Documentation of the Medical Device Manufacturer.

Information about the release of the instructions for use: Version No. 2. The latest version of the instructions for use was approved on August 28, 2024.

1. NAME OF THE MEDICAL DEVICE AND DESIGN VERSIONS

1.1 Description of the medical device:

Disposable absorbent medical products Dr.Skipp*.

1.2 Design versions of the medical device:

1. Dr.Skipp* adult pull-up diapers, size S, containing:
 - Pull-up diapers in consumer packaging, 20 pcs/pack;
 - Instructions for use, 1 pc. per transport packaging unit.
2. Dr.Skipp* adult pull-up diapers, size M, containing:
 - Pull-up diapers in consumer packaging, 20 pcs/pack;
 - Instructions for use, 1 pc. per transport packaging unit.
3. Dr.Skipp* adult pull-up diapers, size L, containing:
 - Pull-up diapers in consumer packaging, 20 pcs/pack;
 - Instructions for use, 1 pc. per transport packaging unit.
4. Dr.Skipp* adult pull-up diapers, size XL, containing:
 - Pull-up diapers in consumer packaging, 20 pcs/pack;
 - Instructions for use, 1 pc. per transport packaging unit.
5. Dr.Skipp* adult diapers, size S, containing:
 - Diapers in consumer packaging, 30 pcs/pack;
 - Instructions for use, 1 pc. per transport packaging unit.
6. Dr.Skipp* adult diapers, size M, containing:
 - Diapers in consumer packaging, 30 pcs/pack;
 - Instructions for use, 1 pc. per transport packaging unit.
7. Dr.Skipp* adult diapers, size L, containing:
 - Diapers in consumer packaging, 30 pcs/pack;
 - Instructions for use, 1 pc. per transport packaging unit.
8. Dr.Skipp* adult diapers, size XL, containing:
 - Diapers in consumer packaging, 30 pcs/pack;
 - Instructions for use, 1 pc. per transport packaging unit.

(hereinafter referred to as the "product medical device").

2. INFORMATION ABOUT THE MANUFACTURER OF THE MEDICAL DEVICE, THE DEVELOPER OF THE MEDICAL DEVICE, THE ADDRESS OF THE MEDICAL DEVICE MANUFACTURING SITE, THE AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE TERRITORY OF THE RUSSIAN FEDERATION

2.1 Information about the medical device manufacturer

Name: Fujian Putian Licheng Paper Industry Co., Ltd

Location address: Jiaoxi Village, Huating, Chengxiang District, Putian City, Fujian Province, China;

Phone: +86 0594 2029539

Email: lichengdiapers@163.com

2.2 Information about the medical device developer

Name: Fujian Putian Licheng Paper Industry Co., Ltd

Location address: Jiaoxi Village, Huating, Chengxiang District, Putian City, Fujian Province, China;

Phone: +86 0594 2029539

Email: lichengdiapers@163.com

2.3 Address of the medical device manufacturing site

Fujian Putian Licheng Paper Industry Co., Ltd; Jiaoxi Village, Huating, Chengxiang District, Putian City, Fujian Province, China.

2.4 Authorized representative of the manufacturer in the territory of the Russian Federation

Name: RM-plus, Limited Liability Company (RM-plus, LLC);

Location address: Russia, 117246, Moscow, Nauchnyi proyezd d.19, floor 8/ room 7;

Phone: +7 (495) 642-82-60;

Email: msk@rmplus.ru.

3. INTENDED PURPOSE OF THE MEDICAL DEVICE:

Disposable absorbent medical products Dr.Skipp[®] are designed to absorb and retain urine when caring for adult patients suffering from severe incontinence and bedridden patients.

4. SCOPE OF USE OF THE MEDICAL DEVICE:

Caring for patients suffering from urinary incontinence and caring for bedridden patients. The products can be used in general medical institutions as well as at home.

5. FREQUENCY OF USE OF THE MEDICAL DEVICE:

Non-sterile and non-sterilizable single-use medical device. The product is disposable and intended for single use by one patient. The recommended time of using the product is no more than 8 hours (in case of polyuria syndrome, no more than 5 hours).

6. POTENTIAL USERS OF THE MEDICAL DEVICE:

The potential users of this medical device are adults suffering from severe urinary incontinence and bedridden patients.

7. INDICATIONS FOR USE OF THE MEDICAL DEVICE:

The indications for use of the products are persistent, pronounced or significantly pronounced impairment of the excretory function as a result of:

- diseases, consequences of injuries and malformations of the central, peripheral nervous system, including conditions resulting from cerebrovascular diseases;
- diseases, consequences of injuries, chronic infections and malformations of the genitourinary system;
- malignant neoplasms (clinical group IV);
- mental disorders with severe or profound mental retardation and dementia;
- other diseases accompanied by severe urinary incontinence.

8. CONTRAINDICATIONS FOR USE OF THE MEDICAL DEVICE:

- individual intolerance to the materials of which the products are made.
- diseases of the external genitalia;
- presence of skin lesions in places where the patient's body comes into contact with the product.

9. POSSIBLE SIDE EFFECTS WHEN USING THE MEDICAL DEVICE:

An allergic reaction may occur due to individual intolerance to the materials of which the products are made. If signs of an allergic reaction appear, stop using the products and seek medical advice. If you select the wrong size or use excessive tension during fastening, deformation or damage to the skin may occur.

10. INTERACTION WITH OTHER MEDICAL DEVICES:

Disposable absorbent medical products Dr.Skipp® are used independently of other medical devices.

11. PRESENCE OF MEDICINAL PRODUCTS:

The medical device does not contain medicinal products.

12. PRESENCE OF MATERIALS OF ANIMAL AND/OR HUMAN ORIGIN:

The medical device contains no materials of animal and/or human origin.

13. USE OF BLOOD AND BLOOD COMPONENTS:

The medical device does not contain derivatives of human blood or its components as an integral part necessary for operation.

14. DESCRIPTION AND PRINCIPLE OF OPERATION OF THE MEDICAL DEVICE

The products are multi-layer disposable incontinence products with an absorbent layer of fibrous pulp containing moisture-absorbing substances (superabsorbents). The products allow patients suffering from incontinence to maintain independent active lifestyle (from providing the opportunity to perform family and household duties to carrying out professional and business activities) despite the presence of pathological changes in the body. The products retain moisture and do not let it out, thereby protecting the skin from irritation upon contact with secretions and preventing contamination of clothing.

The products ensure that skin stays dry and prevent the proliferation of pathogenic flora, i.e. they are used as a means of rehabilitation and care to prevent damage to the skin and create physical and psychological comfort.

The disposable medical incontinence products are manufactured:

- 1) in the form of underpants (pull-up diapers)
- 2) in the form of cut underpants with fixing elements (diapers).

Depending on the amount of fluid absorbed, the products can be used for people with severe incontinence.

The shape of diapers and pull-up diapers corresponds to the shape of the part of the human torso with an additional increase in the area to properly cover the sides. The absorbent layer of the products has a shape that allows the products to be used by both men and women. The diapers and pull-up diapers are a multi-layer product.

The design of the products includes (starting from the layer in contact with human skin):

- **top cover layer:** a layer that is in direct contact with human skin and allows liquid to pass into the product.
- **distribution layer:** a layer that is located behind the top cover layer of the diaper and promotes uniform distribution of liquid inside the diaper. The distribution layer is present only in diapers; it is absent in pull-up diapers.
- **absorbent layer:** the inner main layer of the product, which absorbs and retains the absorbed liquid inside the product.
- **protective layer:** a layer that is located directly behind the absorbent layer of the product and prevents liquid from leaking out.
- **bottom cover layer:** the layer that is located behind the protective layer of the product and is in contact with a person's clothing or bed sheet.
- **barrier elements:** barriers and barrier cuffs in the form of an arc through the groin with elastic bands that tighten them, preventing the leakage of liquid onto human skin as well as ensuring the anatomical shape of the product and the best fit.
- **fixing elements (only for diapers):** details of the diaper on its front and back parts for fixing it in the desired position and securing it to the person's waist (4 reusable hook-and-loop fasteners making it possible to ideally adjust the diaper to the anatomical features of any patient). The diapers have an elastic corrugated tie in the lumbar area, which makes it possible to adapt the diaper to any body shape.

- **wetness indicator:** signals in the form of blue colored stripes on the protective layer, indicating complete filling of the product with absorbed liquid.

Depending on the waist/hip circumference of a person, **diapers** are made in the following sizes:

- 1) size S: person's waist/hip circumference 60-90 cm;
- 2) size M: person's waist/hip circumference 80-120 cm;
- 3) size L: person's waist/hip circumference 100-150 cm;
- 4) size XL: person's waist/hip circumference 130-170 cm;

Depending on the waist/hip circumference of a person, **pull-up diapers** are made in the following sizes:

- 1) size S: person's waist/hip circumference 60-90 cm;
- 2) size M: person's waist/hip circumference 80-120 cm;
- 3) size L: person's waist/hip circumference 100-140 cm;
- 4) size XL: person's waist/hip circumference 130-170 cm;

The appearance of the product must correspond to the image in Figures 1-2.

Figure 1 – Appearance (schematic illustration) of Dr.Skipp® adult pull-up diapers


Figure 2 – Appearance (schematic illustration) of Dr.Skipp® adult diapers

15. SPECIFICATIONS

15.1 The bottom cover layer of the products is finished with a printed pattern.

15.2 External defects are not allowed in products: mechanical damage (torn edges, cuts, damage to fixing elements, etc.), stains of various origins, foreign inclusions visible to the naked eye.

15.3 The printed image on the products must be clear, without distortions or gaps. There should be no plucking of fibers from the surface of the product and no paint smearing.

15.4 The layers of the products are held together with hot melt adhesive, which ensures the strength of the gluing of the layers (seams) of the products. The seams should be continuous.

15.5 The fastening elements of diapers must provide a stable connection to the bottom cover layer after fixing the diapers at least 5 times. The product is disposable, however, it may require repeated fixation to find the desired position and secure it to the person's waist.

15.6 The force that must be applied to separate the fastening elements of the diapers and the bottom cover layer must be at least 6.5 N.

15.7 The main parameters of the product (diapers and pull-up diapers) must correspond to the parameters specified in Table No. 1 and Table No. 2.

Table 1 – Parameters for Dr.Skipp® adult diapers.

Parameter	Dimensions according to size designation			
	S	M	L	XL
Waist/hip circumference, mm	600-900	800-1200	1000-1500	1300-1700
Length when unfolded, mm	700±70	800±80	900±90	1000±100
Width at the edges when unfolded, mm	520±52	630±63	750±75	830±83
Width in the center when unfolded, mm	290±29	290±29	290±29	290±29
Length of absorbent layer, mm	500±50	560±56	640±64	740±74
Width of absorbent layer, mm	170±17	170±17	170±17	170±17
Number of fixing elements, pcs	4	4	4	4
Length of fixing element, mm	Along the long edge – 50±5 Along the short edge – 35±5	Along the long edge – 50±5 Along the short edge – 35±5	Along the long edge – 50±5 Along the short edge – 35±5	Along the long edge – 50±5 Along the short edge – 35±5
Width of fixing element, mm	20±2	20±2	20±2	20±2
Maximum distance between the upper and lower fixing elements (Velcro fasteners), mm, NMT	50	50	50	50
Width of side cuffs, mm	50±5	55±5	55±5	60±6
Width of internal barriers, mm	40±4	40±4	40±4	40±4
Number of elastic bands on the barrier cuffs in the form of an arc through the groin, pcs.	3	3	3	3
Number of elastic bands on barriers, pcs.	2	2	2	2
Length of corrugated tie, mm	130±13	130±13	130±13	130±13
Width of corrugated tie, mm	35±5	60±6	60±6	60±6
Length of corrugated tie in stretched state, mm	215±22	215±22	215±22	215±22
Width of corrugated tie in stretched state, mm	35±5	60±6	60±6	60±6
Weight of dry diaper, g, NLT	70	80	90	100
Availability of wetness indicator	available	available	available	available
Absorbent mass, g, NLT	5	10	12	16
Complete moisture absorption, g, NLT	2000	2300	2400	3100
Sorption capacity after centrifugation, g, NLT	440	480	540	650
Resorption, g, NMT	4.4	4.4	4.4	4.4
Absorption rate, sm^3/s , NLT	2.3	2.3	2.3	2.3

Table 2 – Parameters for Dr.Skipp® adult pull-up diapers.

Parameter	Dimensions according to size designation			
	S	M	L	XL
Waist/hip circumference, mm	600-900	800-1200	1000-1400	1300-1700
Length when unfolded, mm	780±78	820±82	890±89	960±96
Width at the edges (at the waist), mm	530±53	650±65	810±81	880±88
Width at the center, mm	240±24	250±25	250±25	270±27
Length of absorbent layer, mm	430±43	440±44	550±55	550±55
Width of absorbent layer at the	130±15	135±15	135±15	135±15

edges, mm				
Width of absorbent layer in the center, mm	140±14	145±14	145±14	145±14
Width of side cuffs, mm	40±4	40±4	40±4	60±6
Width of internal barriers, mm	60±6	60±6	60±6	60±6
Number of elastic bands on the barrier cuffs in the form of an arc through the groin, pcs.	3	3	3	3
Number of elastic bands on barriers, pcs.	3	3	3	3
Number of elastic bands in the waistband, pcs.	17	18	19	22
Waistband width (not stretched), mm	140±14	160±16	160±16	190±19
Weight of dry product, g, NLT	65	70	75	80
Availability of wetness indicator	available	available	available	available
Absorbent mass, g, NLT	5	10	12	16
Complete moisture absorption, g, NLT	1600	2000	2400	2800
Sorption capacity after centrifugation, g, NLT	400	440	500	600
Resorption, g, NMT	4.4	4.4	4.4	4.4
Absorption rate, sm^3/s , NLT	2.3	2.3	2.3	2.3

16. TYPE OF CONTACT WITH THE HUMAN BODY

The type of contact with the body is short-term contact (less than 24 hours) with intact skin.

The products are resistant to biological fluids and body tissue secretions.

17. DELIVERY PACKAGE

17.1 The delivery package for diapers shall include:

- Diapers in consumer packaging, 30 pcs/pack;
- Instructions for use, 1 pc. per transport packaging unit.

17.2 The delivery package for pull-up diapers shall include:

- Pull-up diapers in consumer packaging, 20 pcs/pack;
- Instructions for use, 1 pc. per transport packaging unit.

18. MARKING

18.1 The marking shall be valid, verifiable and legible. The marking is affixed to the packaging or a label (tag) attached to the packaging. The marking is affixed by any method (printing, embossing, stamping) that ensures its clarity, sharpness and legibility. When the marking is printed, the ink shall not set-off. The marking shall be in Russian.

18.2 Each consumer packaging or label affixed to the consumer packaging must be marked with the following information:

- manufacturer's name;
- manufacturer's legal address;
- address of the manufacturing site;
- name of the authorized representative of the manufacturer in the territory of the Russian Federation;
- legal address of the authorized representative of the manufacturer in the territory of the Russian Federation;
- contact details of the authorized representative of the manufacturer in the territory of the Russian Federation (e-mail, phone number);

- designation of the number and date of issue of the marketing authorization;
- trademark;
- full name of the medical device;
- name of the product design version;
- size of the product;
- product dimensions according to waist/hip circumference;
- type of product depending on the purpose (degree of urinary incontinence);
- instructions for use regarding the method of use of the medical device (in the form of pictures);
- instructions for disposal of the product;
- number of products in consumer packaging, pcs;
- lot code;
- date (day, month, year) of manufacture;
- manufacturer's shelf life;
- inscription "COMPLIES WITH GOST R ISO 55082-2012";
- inscription "DO NOT PUT DOWN THE DRAIN";
- symbols in accordance with Table 3:

Table 3

Symbol	Symbol meaning
	- manufacturer
	- lot code
	- manufacturing date (date, month, year)
	- use before (day, month, year)
	- do not re-use
	- keep away from sunlight
	- avoid exposure to moisture
	- lower limit of the temperature range
	- upper limit of the temperature range

Information about the product (trademark, advertising slogan "ULTRA") must be printed on the lower cover layer of the product.

It is not allowed to affix information about special properties of the diaper, i.e. "protects skin from irritation", "absorbs odor", "breathable", "environmentally friendly" etc., without appropriate documentary confirmation.

18.3 Each transport packaging or label affixed to the transport packaging must be marked with the following information:

- manufacturer's name;
- manufacturer's legal address;
- address of the manufacturing site;
- name of the authorized representative of the manufacturer in the territory of the Russian Federation;

- legal address of the authorized representative of the manufacturer in the territory of the Russian Federation;
- contact details of the authorized representative of the manufacturer in the territory of the Russian Federation (e-mail, phone number);
- designation of the number and date of issue of the marketing authorization;
- trademark;
- full name of the medical device;
- name of the product design version;
- lot code;
- size of the product;
- date (day, month, year) of manufacture;
- manufacturer's shelf life;
- total number of products in transport packaging;
- number of consumer packaging units in the transport packaging;
- inscription "STORAGE CONDITIONS: The products must be stored at an ambient temperature of +5 °C to +40 °C and a relative humidity of no more than 80%."
- inscription "TRANSPORTATION CONDITIONS: The transportation of the products must be carried out at ambient temperatures from -50 °C to +50 °C and a relative humidity of no more than 80%";
- instructions for disposal of the product;
- symbols in accordance with Table 3:
 - manufacturer;
 - lot code;
 - manufacturing date (date, month, year);
 - use before (day, month, year);
 - keep away from sunlight;
 - avoid exposure to moisture.

19. PACKAGING

19.1 The products in quantities determined by the manufacturer are packaged in plastic film bags (consumer packaging), which ensures the safety of the products during transportation and storage. The products in the form of ready-made underpants (pull-up diapers) are packaged in plastic film bags in the amount of 20 pcs. The dimensions of consumer packaging for pull-up diapers (L · W · H) are:

Size S: (425±20) mm x (155±20) mm x (160±20) mm;

Size M: (440±20) mm x (155±20) mm x (165±20) mm;

Size L: (475±20) mm x (155±20) mm x (185±20) mm;

Size XL: (485±20) mm x (155±20) mm x (190±20) mm.

The products in the form of cut underpants with fixing elements (diapers) are packaged in plastic film bags in the amount of 30 pieces. The dimensions of consumer packaging for diapers (L×W×H) are: Size S: (300±20) mm x (170±20) mm x (210±20) mm;

Size M: (320±20) mm x (180±20) mm x (220±20) mm;

Size L: (320±20) mm x (190±20) mm x (240±20) mm;

Size XL: (350±20) mm x (180±20) mm x (250±20) mm.

19.2 The consumer packaging is made in the form of a polyethylene bag manufactured by Jinjiang ZW Printing Co., Ltd, China. Seams in plastic film (polyethylene) bags should be welded.

19.3 One consumer packaging contains products of the same group, type, size, design, including technical and decorative design, made of the same materials, with the same quality parameters and the same manufacture date (day, month, year).

19.4 The consumer packaging must provide:

- protection of the product from external mechanical damage during storage and transportation;
- the impossibility of resealing the packaging if its integrity has been compromised as the fact of opening the packaging shall be obvious;

19.5 The consumer packaging with the products inserted into it should be intact, without tears or other mechanical damage on its surface that would allow access to the surface of the products.

19.6 The medical devices in consumer packaging are placed in transport packaging – a polyethylene pack manufactured by Jinjiang ZW Printing Co., Ltd, China. The number of products in transport packaging is established by the supply agreement. The gross weight of the transport packaging should not exceed 30 kg. One transport packaging contains products of the same group, type, size, design, including technical and decorative design, made of the same materials, with the same quality indicators and the same manufacture date (day, month, year).

19.7 The pack is sealed with a polyethylene pressure-sensitive adhesive tape (duct tape).

20. METHOD OF USE

Preparation for use:

- Inspect the packaging carefully to make sure it is not damaged. Do not use the product with damaged packaging;
- Before using the product, hands must be treated with an antiseptic;
- Before use, the patient should be washed with warm water and wiped dry.

WARNING! Using the products for purposes other than its intended purpose is strictly prohibited. Before using the product, you must check its guaranteed shelf life. It is not allowed to use the product after the expiration of the guaranteed shelf life.

Algorithm for using diapers:

The procedure for putting on a diaper in a vertical position (independent use):


1. Open the consumer packaging, then remove a diaper from the consumer packaging, straighten all its parts, smooth out the folds and stretch slightly.
2. Place the back of the diaper between your legs with the top cover layer toward your body from front to back.
3. Flatten the front of the diaper over your stomach.
4. Flatten out the back and sides of the diaper.
5. Wrap the back of the diaper around the front and secure the fasteners, starting from the bottom.
6. Make sure that the diaper fits snugly to the body but does not squeeze it.

The procedure for putting a diaper on a patient in a lying position:

1. Open the consumer packaging and remove a diaper from the consumer packaging, straighten all its parts, smooth out the folds and stretch slightly.
2. Lay the patient on their side with their knees bent;
3. Fold the diaper and place it between slightly spread legs, covering the crotch (you should start from the front, pushing the diaper backward and flattening out the back);
4. Turn the patient on their back and align the diaper on their stomach so that the creases are minimal;
5. The edges of the product should fit snugly to the body; this will help to avoid leakage; after making sure that the diaper fits well, secure it with the fasteners – first with the lower ones, to create a tighter fit around the legs, then with the upper ones.

The procedure for removing the diaper from the patient:


1. With the patient lying on their back, unfasten the diaper fasteners.
2. Unfold the front of the diaper, place the patient on their side, and carefully roll up the bottom of the diaper.
3. When removing a diaper from the ordinary patient, ask them to stand up straight and make the same steps.
4. Remember to put on and take off the diaper only from the back, as this prevents accidental transfer of bacteria from the intestines.
5. To avoid leakage of the diaper, strictly adhere to all the recommendations of the manufacturer.
6. Dispose of the used diaper in the prescribed manner.

Algorithm for using diapers

Algorithm for using pull-up diapers:

- 1) Open the consumer packaging and remove a pull-up diaper from the consumer packaging;
- 2) Put on like usual underwear and straighten the protective barriers.
- 3) Remove used pull-up diapers like usual underwear or rip the side seams.
- 4) Dispose of used pull-up diapers in the prescribed manner.

Algorithm for using pull-up diapers

21. DIRECTIONS FOR USE

Conditions for use:

The product can be used in general medical institutions and at home at ambient temperatures from +5°C to +40°C and relative air humidity up to 98%.

After transportation at negative temperatures, the products must be kept unpacked at an ambient temperature of +5°C to +40°C for at least 2 hours.

The product is intended for single use only. Any re-use is PROHIBITED.

Precautions and warnings for use.

- 1) The product must not be used if the consumer packaging is damaged or if the guaranteed shelf life has expired.
- 2) The product is intended for external use only.
- 3) The product is intended for single use only.
- 4) Do not use the product if it has visible traces of contamination with biological fluids.
- 5) It is prohibited to cut or change the shape of the product as this may affect its effectiveness;
- 6) Products should not be stored in conditions of high humidity, near heating devices or in direct sunlight.
- 7) Products should not be stored next to washing powders, machine oil and other strong-smelling substances because products even in closed packaging may become saturated with foreign odors.
- 8) Products should not be stored near open fire as this may cause the product to deform and fade.
- 9) If the condition of the skin in the area of contact with the product worsens and irritation occurs, stop using the product and consult a doctor.

- 10) Avoid getting moisture on clean products.
- 11) Keep out of the reach of children and incapacitated patients.

WARNING! Do not use the product for more than 8 hours (in case of polyuria syndrome, no more than 5 hours). It is recommended to change the product every 3-4 hours or more often if the product gets wet faster.

WARNING! The product is intended for single use only.

DO NOT USE THE PRODUCT IN THE FOLLOWING CASES: external defects of the product and mechanical damage are detected: compromised integrity of consumer packaging, foreign inclusions, tears, cracks, folds and/or delamination.

Precautions to be taken by the consumer when disposing of the medical device

After use, dispose of the product according to the instructions given in the DISPOSAL section.

22. REQUIREMENTS FOR CLEANING AND DISINFECTION

The products are classified as non-sterile single-use medical devices. They are not subject to cleaning and disinfection.

23. MAINTENANCE, CARE AND REPAIR

The products are classified as non-sterile single-use medical devices. The products have a completely finished design. They are not subject to maintenance and repair.

24. TRANSPORTATION CONDITIONS

The transportation of packaged products is carried out by any type of transport, provided they are protected from mechanical influences, precipitation and aggressive environments, in accordance with the rules for the transportation of goods in force for this type of transport.

The transportation of products must be carried out at ambient temperatures from -50°C to +50°C and a relative humidity of no more than 80%.

When loading and transporting the products in transport packaging, the requirements of handling signs and warning labels shall be met.

25. STORAGE CONDITIONS

The products in the manufacturer's packaging are stored in closed, clean and dry rooms that are protected from moisture and direct sunlight, at a distance of at least 1 m from heating devices, and exclude exposure to aggressive environments (acidic, alkaline, etc.) as well as flammable and combustible liquids. The products should not be stored near washing powders, machine oil and other strong-smelling substances because the products even in closed packaging may become saturated with foreign odors. Products should not be stored near open fire. Keep out of the reach of children and incapacitated patients.

The products must be stored at an ambient temperature of +5°C to +40°C and a relative humidity of no more than 80%.

26. DISPOSAL REQUIREMENTS

The products used for their intended purpose are disposed of in the manner prescribed by SanPiN 2.1.3684-21 for class B medical waste.

The products not used for their intended purpose, expired products and products with damaged packaging are disposed of in the manner prescribed by SanPiN 2.1.3684-21 for class A medical waste.

The procedure and conditions for disposing of the product at home: when at home, the used disposable product must be placed in a separate bag, hermetically sealed and then put in the trashcan. If there are containers for separate collection of solid municipal waste in the yard of the house, make sure that used products do not end up

in the waste container for recycling.

The disposal of the products must be carried out in accordance with the rules and regulations in force in the Russian Federation at the time of disposal.

27. GUARANTEED SHELF LIFE OF THE PRODUCT

The guaranteed shelf life of the products is 3 years from the date of manufacture.

28. MANUFACTURER'S WARRANTY

The guaranteed shelf life of the products is 3 years from the date of manufacture.

The manufacturing company guarantees the quality, safety and efficacy of the device if it is used in accordance with the instructions for use before the guaranteed shelf life has expired, subject to the conditions of use, transportation and storage.

29. CLAIMS:

Regarding matters concerning the quality of the medical device, please contact the authorized representative of the manufacturer in the territory of the Russian Federation. In the event of undesirable events that have signs of an adverse event (incident), please send the information (message) about these events to the authorized representative of the manufacturer in the territory of the Russian Federation.

Authorized representative of the manufacturer in the Russian Federation:

Name: RM-plus, Limited Liability Company (RM-plus, LLC);

Location address: Russia, 117246, Moscow, Nauchniy proyezd d.19, floor 8/ room 7;

Phone: +7 (495) 642-82-60;

Email: msk@rmplus.ru.

30. LIST OF INDUSTRY STANDARDS THAT THE MEDICAL DEVICE COMPLIES WITH

Order No. 4н of the Ministry of Health of the Russian Federation dated 06.06.2012 On the approval of the nomenclature classification of medical devices.

GOST R 55082-2012 Articles of paper for medical purposes. Diapers for adults. General specifications.

GOST R ISO 15223-1-2023 Medical devices. Symbols to be used with information to be supplied by the manufacturer. Part 1. General requirements.

GOST ISO 14971-2021 Medical devices. Application of risk management to medical devices.

SanPiN 2.1.3684-21 Sanitary and epidemiological requirements for the maintenance of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial and public premises, and organization and conduct of sanitary and anti-epidemic (preventive) measures.

GOST R 52770-2023 Medical products. A system for evaluating biological effects. General safety requirements.

GOST ISO 10993-5-2023 Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity.

GOST ISO 10993-10-2023 Medical products. Assessment of the biological effect of medical devices. Part 10. Studies of the sensitizing effect.

GOST ISO 10993-23-2023 Medical products. Assessment of the biological effect of medical devices. Part 23. Studies of irritating effects.


31. INTERNATIONAL CERTIFICATES AND DECLARATIONS

The enterprise has implemented, operates and applies:

- ISO 13485 quality management system; based on the results of audits, confirmation has been received that the product meets the requirements of the ISO 13485 standard.

The certificate was issued by the certification body Shanghai POSI Certification Co., Ltd, Room 1301-C-7, No. 1500, Century Avenue, Pudong New Area, Shanghai, China.

Certificate ISO 13485:2016 No. 381190006R1M. Initial issue date: 30.07.2019 (valid until 29.07.2025).

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