LEGALISATION

Seen for legalisation by me, Antal Camillus Theodoor De Witte, Civil-Law Notary, practising in Katwijk, the signature of:

Mr Mohammed Aldirawi, residing at 2272 XD Voorburg, Leidschendam-Voorburg, Rozentuin 74, born in Gaza on the fourth day of July nineteen hundred and seventy-six, in this respect acting in his capacity as Holder of a Power of Attorney with representative powers of **Mentor Medical Systems B.V.**, a private company with limited liability (besloten vennootschap met beperkte aansprakelijkheid), incorporated and existing under the Laws of the Netherlands, having its corporate seat in Leiden, with offices at Zernikedreef 2, 2333 CL Leiden, registered with the Traderegister of the Chamber of Commerce under number 28062039, and as such authorised to validly represent Mentor Medical Systems B.V., according to an extract of the commercial register of the Chamber of Commerce, which is attached to this document.

The aforementioned statement only relates to the signature of Mr Mohammed Aldirawi, aforementioned, and his capacity as Holder of a Power of Attorney of the aforementioned company Mentor Medical Systems B.V., and does not relate to the contents of the document on which his signature was endorsed.

Signed in Katwijk, the Netherlands, on the twenty-third day of May two thousand twenty-four.

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Instruction for use CPG Gel Breast Implant

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1. Pregnancy or nursing mothers.
2. Lupus (e.g., SLE and DLE).
3. Scleroderma (e.g., progressive systemic sclerosis).
4. Currently has a condition that could compromise or complicate wound healing (except reconstruction patients).
5. Infection or abscess anywhere in the body.
6. Demonstrates tissue characteristics, which are clinically incompatible with implant (e.g., tissue damage resulting from radiation, inadequate tissue, or compromised vascularity).
7. Possesses any condition, or is under treatment for any condition which, in the opinion of the consulting physician(s), may constitute an unwarranted surgical risk.
8. Anatomic or physiologic abnormality, which could lead to significant postoperative complications.
9. A history of sensitivity to foreign materials or repeated attempts and failures at breast augmentation or reconstruction.
10. An unwillingness to undergo any further surgery for revision.
11. Unrealistic expectations, such as inappropriate attitude or motivation, or a lack of understanding of the risks involved with the surgical procedure and implants.
12. Premalignant breast disease without a subcutaneous mastectomy.
13. Untreated or inappropriately treated breast malignancy, without mastectomy.
14. Allergic reactions to any of the implant materials or filling agents.

NOTE: The satisfactory use of gel-filled prostheses for tissue replacement following mastectomy or trauma may require special reconstructive procedures, particularly in the presence of radiation damage on the chest

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
DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS OF MD AND PRINCIPLE OF OPERATION

General description of functional elements of MD and principle of operation

The CPG Gel Breast Implants Cohesive III consist of silicone gel encased in a silicone elastomer shell and can be used for reconstruction and augmentation of breast tissue. At this time, manufacturer offers 9 styles of CPG Gel Breast Implants Cohesive III, the 300 series.

1. CPG Gel Breast Implant 311 - CPG 311 Gel Breast Implant Cohesive III.
2. CPG Gel Breast Implant 312 - CPG 312 Gel Breast Implant Cohesive III.
3. CPG Gel Breast Implant 313 - CPG 313 Gel Breast Implant Cohesive III.
4. CPG Gel Breast Implant 321 - CPG 321 Gel Breast Implant Cohesive III.
5. CPG Gel Breast Implant 322 - CPG 322 Gel Breast Implant Cohesive III.
6. CPG Gel Breast Implant 323 - CPG 323 Gel Breast Implant Cohesive III.
7. CPG Gel Breast Implant 331 - CPG 331 Gel Breast Implant Cohesive III.
8. CPG Gel Breast Implant 332 - CPG 332 Gel Breast Implant Cohesive III.
9. CPG Gel Breast Implant 333 - CPG 333 Gel Breast Implant Cohesive III.

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CONSTRUCTION OF THE CPG GEL BREAST IMPLANT FAMILY

The anatomical shaped Breast Implant devices consist of a shell, gel filler and a patch assembly that seals the patch hole, and contains a button with identification of the implant. The devices have a layer of roughened elastomer on their exterior surfaces: MED 4735 sheeting. In addition orientation markings are placed on the shell to aid the surgeon in proper placement of the CPG Gel Breast Implant to prevent rotation of the implant. All components are manufactured from silicones. The Cohesive III silicone gel filler used for the CPG Gel Breast Implant has a firmer consistency than the silicone gel in the Breast Implants due to a different base to cross-linker ratio. This difference is regarded as insignificant in terms of finished product performance. The prostheses are available in a variety of sizes from 100cc until 775cc.

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Layout of MD labeling in Russian

Product code:

CPG anatomically shaped breast implant with silicone gel

Name of design variant

Volume: ml

Area of usage: Surgery

Intended use: The intended purpose of medical device is to replace missing breast tissue or to enhance a breast in females undergoing breast augmentation or breast reconstruction.

Registration certificate: RZN XXXX/XXXX or XX.XX.XXXX

Manufacturer: Mentor Medical Systems B.V., Zernikedreef 2, 2333 CL Leiden, The Netherlands

Authorized representative of the manufacturer, address: Johnson & Johnson LLC, 121614, Moscow, Krylatskaya str., 17/2.

Tel. +7 (495) 580-7777, Fax: +7 (495) 580-7878

Model: see on the package

Height: see packaging

Width: see packaging

Projection: see packaging

Shape: see on the package

Shell: polydimethylsiloxane MED1-6609, polydimethylsiloxane MED-6649

Filling: Medical Implant Grade (Cohesive III) MED3-6309-1.

Transport at temperature from -20 °C to +50°C. Relative humidity not more 75%. Barometric pressure 700-1060 hPa.

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Prospective patients should be made aware of this potentiality. (See ADVERSE REACTIONS/EVENTS section of this insert.)

- Only one prosthesis should be implanted per breast. Mentor recommends against the stacking of implants, one upon the other. The devices have not yet been tested for this use and the integrity of the implants cannot be guaranteed as the materials may abrade and wear. Such abnormal stress may result in weakening or rupture of the prostheses.

- Do not insert or attempt to repair a damaged or altered prosthesis.

- The action of drugs (e.g., antibiotics and steroids) in contact with the prosthesis has not been tested by the manufacturer, and their use cannot be recommended. Each physician who chooses to use chemotherapeutic drugs with this prosthesis must assure compatibility of the drug with the silicone elastomer.

- Do not introduce or make injections of drugs or other substances into the implant. Injections through the implant shell will compromise the product's integrity.

- If povidone-iodine is used within the implant space, the site should be carefully rinsed with saline to remove the residual solution.

- Preoperative evaluation of the implant design, size and implant site should include allowances for adequate tissue coverage. Pressure, force, tension and other stresses to which the implant site will be susceptible must be considered.

- Sepsis, hemorrhage or thrombosis may result from the placement of any foreign object in the body.

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- The use of microwave diathermy in patients with breast implants has been reported to cause tissue necrosis, skin erosion and extrusion of the implant. Its use in patients with breast implants is not recommended.
- The patient should be made aware that any abnormal stress or trauma to the breast could result in rupture of the prosthesis.
- Mentor strongly recommends against the treatment of capsule firmness by forceful external stress (such as closed capsulotomy) and is not responsible for the structural integrity of the implant should the surgeon elect to perform such a procedure. If the physician uses this technique, several complications may occur: hematoma, displacement of the implant and/or shell rupture. The physician should inform the patient of these potential complications and of alternatives to the procedure. Such abnormal stress or trauma to the breast and the prosthesis could result in rupture of the prosthesis.
- Mentor has not tested the *in vivo* effects of radiation therapy on tissues of patients who have breast implants; however, the literature suggests that radiation therapy may increase the likelihood of experiencing capsular contracture. In addition, the literature notes the following about the effects of radiation therapy on implanted breasts: “(a) when the implanted breast was free of fibrotic changes, radiotherapy produced acceptable results, (b) whenever feasible, 45 Gy/5 weeks seemed preferable over higher doses, (c) irradiation immediately after the reconstructive surgery appeared to produce poorer cosmetic results.” The decision regarding the use of radiation therapy following breast implantation should be made by the surgeon and the radiation oncologist.
- The use of the periumbilical approach in placement of the implant has not been established and is not recommended.

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- Careful hemostasis is important to prevent postoperative hematoma formation. Should excessive bleeding persist, it is recommended that the device not be implanted until the bleeding is controlled.
- If a physician treats a hematoma or serous fluid accumulation by aspiration, or if a biopsy or lumpectomy is performed, care must be taken to avoid damaging the implant. These procedures present possible risk of implant puncture.
- The incidence of extrusion of the prosthesis has been shown to increase when the prosthesis has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has previously been performed.
- Excessive fibrous capsular formation or contracture may occur around any implant placed in contact with soft tissues. The incidence and severity of this occurrence may increase if postoperative local hematoma or infection occurs.
- The physician should use personal discretion when deciding to use these prostheses regarding patients who exhibit psychological instability.
- Surgical implantation of a breast implant may interfere with the ability to breast-feed. The Institute of Medicine (IOM) concluded that there is limited evidence that implantation, especially through a peri-areolar incision, may interfere with lactation and breast-feeding. However, it should be noted that previous breast reconstruction surgery, such as mastectomy, may be the initial cause of this interference.
- Patients with connective tissue diseases (CTD) may experience an increased risk of wound dehiscence, infection and bleeding (likely due to their ongoing medical management) that may require further treatment.

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- A patient should seek evaluation by a plastic surgeon to evaluate for BIA-ALCL if they experience delayed onset seroma in the breast (delayed means developing more than 12 months after implantation), a new breast mass or lymphadenopathy. Please see the BIA-ALCL section below for additional information.

Peri-operative Infection and Inflammation Source Reduction Practices:

The following practices could potentially reduce peri-operative bacterial infection and particulate based inflammation:

- Perform a layered closure to separate skin/wound bacteria from the implant cavity.
- Avoid, where possible, pocket dissection that communicates with breast ducts and glands
- Avoid unnecessary adjacent tissue trauma and maintain meticulous hemostasis in the creation of the implant pocket.
- Minimize intra-operative device exposure, skin contact, and pocket placement trials.
- Administer preoperative intravenous antibiotic prophylaxis.
- Lint, dust, talc, surgical glove powder, drape and sponge lint, fingerprints, skin oils and other surface contaminants deposited on an implant by improper handling may cause foreign body reactions. Strict adherence to clean, aseptic techniques should be maintained to prevent contamination of the implant and possible complications. Surgical instruments and gloves should be rinsed clean of any impurities before handling the implant

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Contour Profile

- Ensure that the pocket is not too large as this may result in excessive movement or rotation of the implant.
- Mentor recommends the surgeon consider the firmer nature of the gel, as well as the size of the implant and the firmer nature and high profile of the SILTEX shell when choosing optimal incision size and surgical approach. (See also Implant Selection section of this insert.)

Precautions

It is the responsibility of the surgeon to advise prospective patients or their representatives, prior to surgery, of the possible complications associated with the use of this product.

- Preexisting infection should be treated and resolved before implantation of the prosthesis.
- Any surgeon performing augmentation or reconstructive mammoplasty with implants should be familiar with the currently available techniques for measuring the patient, determining the implant size and performing surgery.
- The silicone elastomer shell may be easily cut by a scalpel or ruptured by excessive stress, manipulation with blunt instruments or penetration by a needle. Subsequent rupture will result. All prostheses should be carefully inspected for structural integrity prior to and during implantation. Meticulous care must be exercised in handling and implanting the device.
- Excessive force and over manipulation during the implantation procedure are to be avoided as this may damage the implant.

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



- A well-defined dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.

Contour Profile only:


The CPG Gel Breast Implant Cohesive III contains raised orientation marks on the anterior and posterior of the implant. These marks help the physician orient the implant and ensure proper placement during implantation.

Anterior Orientation Marks

Posterior Orientation Marks

NOTE: It is advisable to have more than one size mammary implant in the operating room at the time of surgery to allow for flexibility in determining the appropriate size implant to be used. A backup implant should also be available. Reports in the medical literature suggest that prophylactic antibiotic treatment may be indicated during subsequent dental or other surgical procedures for patients with some types of silicone implants.

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Calcified capsules may require removal if the patient wishes relief from her contracture but otherwise seem to be harmless. Small foci of calcification are commonly seen anywhere in the breast parenchyma. They can usually be identified as benign by the radiologist but on occasion may require biopsy to rule out a malignancy.

Rupture of the Implant

Breast implants may not last a lifetime. Although rupture can occur at any time following implantation, a number of studies examining rupture of current generation, single-lumen, CE marked, round and shaped silicone gel-filled breast implants of various manufacturers over time (using MRI evaluations) have consistently reported similar results indicating the average expected lifetime exceeds 10 years.

While the silicone material itself has not been shown to biodegrade, the shell may rupture due to wear and tear, or direct injury. If the implant shell is ruptured, the escaping gel is usually contained by the scar envelope in the surgical pocket (intracapsular) and may be undetectable except by MRI (silent rupture), which is about 85% effective in detecting rupture. If the scar envelope is torn, the gel can be driven into the local tissue planes and breast tissue (extracapsular). Most of the escaped gel remains in the immediate environment of the breast but on rare occasions it has been reported to migrate down the arm, into nerve sheaths or into the abdominal wall. Ultrasound, mammography and physical examination may also diagnose these ruptures, which have escaped the scar

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envelope. Most of the reported cases occurred in the more fragile, thinner shell devices implanted in the late 1970s.

Rupture should be suspected if there is a change in character of the device such as a new, persistent unilateral burning sensation or a change in softness, texture or shape of the implant. Because of the silent nature of most ruptures and the difficulty of diagnosis without surgical exploration, the true incidence is unknown. Current products have thicker and stronger shells and more cohesive gel contents. Caution should be used when comparing expected or actual rupture rates of current devices to historical incidences, especially when, as is often the case, the brand, vintage and type of device is unknown. Explantation and/or replacement may be indicated if the implant fails, especially if it is seen in the breast parenchyma as it could be confused with or mask a tumor.

Causes of implant rupture include, but are not limited to: damage from surgical instruments, intraoperative or postoperative trauma, excessive stresses or manipulations as may occur during daily routines such as vigorous exercise, contact athletics, routine manual massage, intimate physical contact and from compression required during mammography.

Changes in Nipple and Breast Sensation/Breast Pain

Any surgery that undermines the skin of the breast can result in hypersensitivity or loss of sensation in the nipple-areolar region. These changes can vary in degree and may be temporary or permanent. Changes in nipple/breast sensation may, on occasion, affect sexual response or comfort while nursing. These

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changes are believed to be a result of nerve damage or stretching of the nerves from the surgery. There is no specific treatment for this condition.

Most women undergoing augmentation or reconstruction with a breast implant will experience some breast and/or chest pain postoperatively. While this pain normally subsides in most women as they heal after surgery, it can become a chronic problem in other women. Chronic pain can be associated with hematoma, migration, infection, and implants that are too large or capsular contracture.

Sudden severe pain may be associated with implant rupture.

Interference with Mammography in Detection of Cancer/Calcium Deposits

As silicone is opaque to x-rays, an implant may theoretically interfere with the early detection of cancer by mammography as it may obscure part of the breast. Newer techniques of breast compression improve the amount of breast that can be visualized. Alternatively, most surgeons feel that the device may improve the detection of tumors by palpation. While of considerable theoretical concern, delayed detection due strictly to the presence of an implant has not been reported. Women at high risk of developing breast cancer should consider getting implants with caution. Since the breast is compressed during mammography, it is possible for the implant to rupture, but this is rare and should not deter a woman from regular, routine mammographic screening. Before the mammography exam, women should inform the technologist that they have implants.

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Calcium deposits are seen occasionally in old scars anywhere in the body and this is true of the implant capsule. This usually does not occur until years after implant surgery. Benign calcifications are also commonly seen on mammography in otherwise normal breast parenchyma even in breasts that have never been operated on. These benign calcium deposits usually have a different x-ray appearance than the calcifications that signal a malignancy. An expert radiologist can usually determine if a calcium spot is benign or malignant but occasionally a biopsy may be necessary to rule out malignancy. There is no evidence that these deposits occur more or less frequently in women with implants than those who have not received implants. After many years, some patients may develop a thin layer of calcium in the scar capsule that surrounds the implant. This is almost always associated with capsular contracture but otherwise causes no known problem.

Extrusion of the Implant/Interruption of Wound Healing

Skin necrosis and/or sloughing resulting in implant exposure may occur from undue tension of the tissues overlying the implant, the overdosage of steroids placed in the implant pocket or surgical or external trauma.

This is quite rare in the augmentation patient but occurs occasionally in the more challenging tissue environment of reconstruction after mastectomy, local tissue injury and/or radiation. The areas of the scar, especially after radiation to that area, seem to be the most vulnerable to disruption. Extruded implants may need to be surgically removed.

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Wrinkling of the Implant/Dissatisfaction with Cosmetic Results/Asymmetry

Visible and/or palpable wrinkles in an implant are related to the thinness of the overlying tissue cover, the degree of capsular contracture and the texturing of the implant shell surface. Traditional smooth-walled gel devices rarely demonstrate wrinkles. Some have described camouflaging of wrinkling with autologous fat transfer superficial to the implant. Surgical error, preexisting asymmetry or deformity, keloid formation of the incisional scar, vagaries of time, weight gain or loss, pregnancy and nursing can all contribute to an immediate or late poor aesthetic result. With time, most breasts, with or without implants, become ptotic to some degree. Asymmetry is usually related to the inability to totally correct for preexisting disparity between the two breasts. It can also be attributed to surgical error, asymmetric contracture, or rupture of the implant.

Excessive buildup of collagen at the incision site during the healing process causes some patients to develop scars of cosmetic concern. Keloid scars, which do not respond well to treatment, often extend beyond the edges of the original scars and can continue to enlarge over time. Hypertrophic scars are generally confined to the original site and respond well to scar revision treatment, which may include steroid injections to break down the collagen or surgery to revise the position, direction or line of the scar.

Possible Reactions to Silicone

This text contains a brief summary of information from the medical literature. Mentor recognizes that the information contained within this text is highly technical. However, medical ethics and practice dictate that the physician must be an intervening party between the manufacturer of prescriptive medical devices and the patient. The issue of the possible relationship between silicone

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(and other implantable materials) and various diseases has been the subject of significant scientific debate. Concerns include immunological and neurological disorders, carcinogenicity and connective tissue disorders.

Despite numerous anecdotal reports of established and newly described connective tissue diseases attributed to the device, multiple epidemiological studies have very consistently shown no significant association between silicone breast implants and any established or new connective tissue or other immunological diseases. Four prestigious multispecialty scientific expert panels have reviewed the published literature on this topic, specifically as it relates to silicone breast implants, and have issued lengthy reports on their findings. These expert panels include the Independent Review Group (commissioned by the Chief Medical Officer of the U.K.), the National Science Panel (appointed by Judge Pointer for MDL 926), the IOM and the Scientific and Technical Options Assessment (STOA) Programme (commissioned by the European Parliament). These four panels have uniformly concluded that there was no discernible evidence of causal association or positive risk ratio between exposure to silicone breast implants and recognized or new autoimmune or connective tissue diseases. Essentially the same finding was reported in a 2011 updated review of the epidemiology literature published in the peer-reviewed medical literature in 2011. Coincident neurological problems such as multiple sclerosis and amyotrophic lateral sclerosis (ALS) have occurred in a small number of breast implant patients. Subsequent studies have not demonstrated a relationship between any neurological disease and breast implants.

Several studies have been conducted to determine the carcinogenic risk of breast implants and no evidence of increased risk of cancer has been demonstrated.

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Continuing assessment of both known and possible risks associated with breast implant surgery is ongoing.

No credible reports on birth defects or other reproductive effects in humans associated with implantation of silicone breast implants of any type have been identified in the literature. Recent studies sponsored by Mentor provide further evidence that silicone materials used in breast implants do not cause adverse reproductive effects in experimental animals.

Although any breast surgery, including breast implantation, could theoretically interfere with the adequacy of a woman's milk supply, many women with breast implants have nursed their babies successfully. It is known that any breast surgery such as breast biopsy can affect the quantity of milk produced. In recent years, the question has been raised regarding the potential transfer of silicone into the breast milk of women with silicone breast prostheses and possible effects on the health of the breast-fed children. However, more recent studies have provided strong evidence of the lack of association between silicone breast implants and adverse effects in breast-fed children. The American Academy of Pediatrics has stated that there is no reason why a woman with implants should refrain from nursing. The European Committee on Quality Assurance and Medical Devices in Plastic Surgery (EQUAM) concluded in 2000 that silicone gel-filled breast implants do not adversely affect pregnancy, fetal development, breast-feeding, nor the health of breast-fed offspring.

The report sponsored by the IOM, "Safety of Silicone Breast Implants," released in July 1999 states that women with silicone breast implants are no more likely than the rest of the population to develop cancer, immunologic diseases, or neurological problems. The committee also concluded that there is no evidence that mothers with implants pass silicone on to infants when breast-feeding.

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Anaplastic Large Cell Lymphoma

Based on information reported to global regulatory agencies and found in medical literature, an association has been identified between breast implants and the development of anaplastic large cell lymphoma (ALCL), a type of non-Hodgkin's lymphoma. Women with breast implants have a very small but increased risk of developing Breast Implant Associated ALCL (BIA-ALCL) in the fluid or scar capsule adjacent to the implant, with documented potential for local, regional, and distant spread of the cancer with mortality reported in rare cases.

BIA-ALCL has been reported globally in patients with an implant history that includes Mentor's and other manufacturers' breast implants with various surface properties, styles, and shapes. Most of the cases in the literature reports describe a history of the use of textured implants. Several journal articles explore the risk

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factors for BIA-ALCL, including the varied methods used to create surface texture of the implant and the role of biofilm in causing disease, among others.

You should consider the possibility of BIA-ALCL when a patient presents with a late onset seroma (at least 1 year and on average 7 to 10 years following implantation), breast mass adjacent to the implant capsule or lymphadenopathy. Work-up of delayed enlargement of a previously implanted breast includes an ultrasound or MRI to look for the presence of seroma or lymphoma. Diagnosis starts at presentation, with needle aspiration of as much of the fluid as possible (minimum 50 mL), using ultrasound for implant displacement and protection. A suspicious mass requires needle or open tissue biopsy. Specimens should be sent for cell morphology by cytology, CD30 immunohistochemistry and ALK (anaplastic lymphoma kinase), and flow cytometry for evaluation, quantification, and characterization of T cells within the specimen. Hematopathology consultation at a tertiary cancer center is strongly encouraged to establish or exclude a diagnosis of BIA-ALCL. Once the diagnosis of BIA-ALCL has been established, physicians are strongly encouraged to consult with a multidisciplinary team including oncologists, pathologists, surgical oncologists, and plastic surgeons. The United States National Comprehensive Cancer Network (NCCN) guidelines for treatment of established BIA ALCL includes the removal of the implant and the entire surrounding implant capsule in addition to associated masses and lymph node involvement. 7

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Please note for latest diagnostic and treatment options refer The National Comprehensive Cancer Network (NCCN) guidelines

Gel Bleed

The gel in an implant consists of large three-dimensional mesh-like molecules that constitute about 20% of the total weight of the gel. The spaces within these molecules are filled with medical grade viscous silicone fluid. This viscous fluid is similar to the materials available in many products including antifatulence medication, available without prescription for infants and adults. A small amount of this material can diffuse or bleed through the shell of the implant. The major portion of this material stays on the implant wall. A small amount moves into the scar capsule where it is gradually picked up by certain "scavenger" cells of the body's immune system, called macrophages. Normally, these cells try to destroy foreign material such as bacteria. However, if the material (such as silicone) cannot be destroyed, it is carried to the lymph glands by the macrophages. As noted in the section on Possible Reactions to Silicone, the STOA report published in 2000 concluded: "Studies do not point to an association between silicone implants and serious health risks, such as cancer and connective tissue diseases."

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Granulomas

It is possible for a granuloma to form around a tiny amount of silicone. Although these lumps are noncancerous, they can be difficult to distinguish from cancerous lumps without being removed (biopsied) and examined. Granulomas, if large or suspicious for malignancy, may need to be biopsied or surgically removed and examined.

Other Possible Reactions*

- Thrombosed veins, resembling large cords, have temporarily developed in the area of the prosthesis and have resolved without surgical or medical therapy.
- Pain from an improperly sized and/or placed implant, such as from compression of nerves or interference with muscle movement, may occur.
- Hypertrophic scarring has been reported.
- The prosthesis may become difficult to explant if the degree of tissue adhesion is significant.

Instructions and Precautions for Implant Removal

Standard surgical approaches and practices should be used should implant removal be necessary. The entirety of the device needs to be removed. In the case of a ruptured gel implant, all of the gel should be removed as far as possible.

**According to Order 1113n of 19.10.2020 corresponds to the term "adverse events and adverse reactions".*

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ADVERSE EVENT REPORTING

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

PRODUCT EVALUATION

Mentor requires that any complications and/or explantation resulting from the use of this device be brought to the immediate attention of your local Mentor representative, who will be responsible for informing the Mentor Complaint Department. If explantation is necessary, analysis will be performed on the explanted device and the patient and the physician must be asked for permission to allow tests to be performed which might alter the condition of the device.

RETURNED GOODS AUTHORIZATION

Authorization for return of merchandise should be obtained from your local Mentor representative prior to the return of merchandise. Merchandise must have all manufacturer's seals intact to be eligible for credit or replacement. Returned products may be subject to a restocking charge.

SAFE DISPOSAL

Devices that have been in contact with material of human origin may be infectious. Dispose with the necessary precautionary measures in accordance with local regulations. Note: devices associated with any complaints or potential adverse event shall be brought to the immediate attention of your local Mentor representative for return (see "PRODUCT EVALUATION" section above).

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INFORMATION A PHYSICIAN SHOULD PROVIDE TO THE PATIENT

Breast implantation is an elective procedure and the patient must be well counseled on the risk-benefit relationship. The surgeon should provide each prospective patient with the following:

1. Device Description

MENTOR Gel Breast Implants is constructed of silicone. Mentor offers Gel Breast Implants with shell surfaces: SILTEX. The SILTEX shell is textured. All CPG Gel Breast Implants are textured.

All MENTOR Gel Breast Implants contain silicone gel which is cohesive. Mentor offers Gel Breast Implants with a variety of cohesivity levels of the filling material. The devices are available in a round shape with different projections and in several contour shapes with different heights and projections. The volume indicated on the product label is the fill volume of the gel.

2. Device Materials

MENTOR Gel-Filled Breast Implants have been tested for safe use. Materials to which patients can come in contact with are 100% medical implant grade silicone. The shells of the MENTOR Gel-Filled Breast Implants are constructed from medical-grade silicone elastomer. The shells are filled with medical-grade, cohesive silicone gel.

3. Information for Safe Use

3.1 General Safety Instructions for patient

Breast implants are for women 18 and older who want breast augmentation or reconstruction.

Breast implants surgery is not for women who have any of the following conditions:

- Pregnancy or nursing mothers.
- Autoimmune disease
- Conditions that would increase surgical risk*
- A history of sensitivity to foreign materials

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- Repeated attempts and failures at breast augmentation or reconstruction.
- An unwillingness to have further surgery
- Unrealistic expectations
 - inappropriate attitude or motivation
 - lack of understanding of the risks involved with the surgery and implants.
- Untreated breast cancer or premalignant breast disease *
- * Discuss with your surgeon any possible conditions.
- **Speak to your breast surgeon and oncologist about your cancer treatment plan.

3.2 Postoperative Warnings, Precautions, or Measures Breast Self-Examination Measures

You should learn to perform breast self-examinations monthly. You should also be shown how to distinguish the implant from their breast tissue.

Warnings

Silicone gel can leak or “bleed” through the undamaged shell into the capsule and adjacent breast tissue. The long-term effects of such “bleed” are unknown. You should be aware of this potentiality.

- The use of microwave diathermy is NOT recommended. It may cause negative side effects to your surgical area.
- You should be aware that any abnormal stress or trauma to the breast could result in damage to the implant.
- Mentor strongly recommends AGAINST the treatment of capsule firmness by forceful external stress (such as closed capsulotomy). It has been reported to cause patient injury and/or implant damage.

Precautions

Advise your Healthcare Provider about your existing implants if

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- you have an autoimmune disease
 - have a weakened immune system
 - currently taking drugs that weaken the body's natural resistance to disease.
- Advise your Healthcare Provider if you are planning any of the following:

- Surgery in the breast area
- Chemotherapy
- Radiation*
- MRI
- Breastfeeding**
- Reoperation***.
- Mammography****

* Radiation treatment may possibly cause damage to the breast implant. Radiation therapy may increase the chance of capsular contracture and/or change how your breasts look

** Breast implants may interfere with breastfeeding

*** Breast implants are not lifetime devices. You will likely undergo implant removal, with or without replacement, in the future. Changes to your breast are irreversible after explant.

**** Breast implants can interfere with the detection of cancer/calcium deposits.

Risk and complications to be aware of are:

- Infection
- Hematoma
- Delayed wound healing
- Capsular contracture III/IV,
- Chest wall deformity
- Enlarged lymph nodes
- Rupture
- Changes in nipple and breast sensation/Breast pain
- Possible reactions to silicone and infection.
- Extrusion
- Wrinkling to the implant

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- Granulomas
- Thrombosed veins.
- Pain
- Hypertrophic scarring
- Gel Breed

If you are experiencing any breast pain or discomfort, contact your breast surgeon. Discuss any questions, concerns, or potential side effects with your physician.

Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL)

Breast implants have been associated with the development of a cancer of the immune system (BIA-ALCL). This cancer occurs more commonly in patients with textured breast implants than smooth implants, although rates are not well defined. If you have breast implants, you have a very small, but increased risk of developing breast implant associated anaplastic large cell lymphoma, or BIA-ALCL. BIA-ALCL is not breast cancer.

It is a rare type of non-Hodgkin's lymphoma (cancer of the immune system). In most cases, BIA-ALCL is found in the scar tissue and fluid near the implant. However, in some cases, it can spread throughout the body. In the cases that have spread beyond the scar tissue and fluid near the implant, rare cases of death have been reported.

You should seek evaluation by a plastic surgeon to evaluate for BIA-ALCL if you experience delayed onset seroma in the breast (delayed means developing more than 12 months after implantation), a new breast mass or lymphadenopathy (disease of the lymph nodes).

If you have breast implants, you should monitor them and follow your routine medical care. It is not necessary to remove your breast implants if you have no symptoms and you have not been diagnosed with BIA-ALCL.

4. Expected Device Lifetime and Follow-up

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Breast implantation may not be a one-time surgery. After the surgery, doctor visits and additional surgery may be needed. Implant removal with or without replacement may be needed as well. Breast implants may not last a lifetime. A rupture can occur at any time following implantation. Studies have indicated that the gel breast implant can have an expected average lifetime of over 10 years.

You should have periodic imaging (e.g., MRI, ultrasound) of your implants to screen for implant rupture. This does not replace additional imaging that may be required, such as screening mammography for breast cancer. Even if you have no symptoms, you should have your first ultrasound or MRI at 5-6 years after your initial implant surgery. Then every 2-3 years thereafter. If you have symptoms at any time or uncertain ultrasound results for breast implant rupture, an MRI is recommended.

Additional doctor visits would be needed to monitor if rupture of the implant has occurred. Speak to your surgeon about regular basis monitoring screening.

5. Implant ID Card

You should have received an International Implant Card from your surgeon. This card will have the catalog and serial number of your implants, along with other information. Carry the card with you and show it to doctors or other healthcare providers when you visit them. It will help them treat you appropriately and protect your breast implants during any medical treatment you need in the future. Your doctor should keep a copy of the Implant ID Card with your medical records.

6. Patient Information Website

You can find any updates about the device on this website <http://ic.jnjmedicaldevices.com>. Your implant card will contain information about your specific device. This can also be found in your medical record

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Magnetic Resonance Imaging (MRI) Compatibility

MRI is one of the imaging methods available to assess the integrity of breast implants and the extent of possible silicone leakage. Owing to its greater sensitivity, MRI is typically the recognized standard method used to evaluate the integrity of breast implants, usually in patients that are symptomatic for rupture (Ahnet al., 1993). MRI offers an operator-independent, global view of the entire implant and surrounding tissues in a variety of planes, providing a modality for implant evaluation and preoperative planning. In a study of 100 symptomatic patients, Ahnet al. (1993) established magnetic resonance imaging (MRI) as a definitive, reliable, and reproducible technique to diagnose both intracapsular and extracapsular ruptures. MRI also enabled the identification of the artifacts that may interfere with the definitive diagnosis of implant rupture.

MR SAFE – Implantable devices Mentor which don't contain metal parts are MR SAFE.

Estimated rupture rates over time for MENTOR Breast Implants, as observed in MENTOR clinical studies.

These rates are calculated by the Kaplan-Meier method of analysis. The Kaplan-Meier method attempts to account for patients not returning for follow-up and adjusts the estimated rupture rate accordingly. These rupture rates include both ruptures that are “suspect” and those that are “confirmed” as defined below.

Suspected ruptures are ruptures that are suspected based on the MRI imaging findings, but which have not been confirmed by surgical explantation and physical examination of the implant, usually because the patient chose not to have the implant removed.

Confirmed ruptures are those verified by examining the implant once removed from the patient.

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Of the overall ruptures (those confirmed as well as those suspected but not confirmed) reported for each device, the proportion of confirmed ruptures was 64% (49 of 77 suspected or confirmed ruptures) for MENTOR MemoryGel Breast Implants, 47% (8 of 17 suspected or confirmed ruptures) for Contour Profile Breast Implants (CPG).

10-year suspected or confirmed implant rupture rates for MENTOR MemoryGel Breast Implants:

Primary Augmentation	Revision Augmentation	Primary Reconstruction	Revision Reconstruction
14.9%	16.5%	24.3%	25.8%

Overall, the 10-year follow-up rate was 53%, thus data accuracy is limited.

10-year suspected or confirmed implant rupture rates for Contour Profile Breast Implants (CPG):

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Patient ID Card

Enclosed with each gel-filled breast implant is a Patient ID Card. To complete the Patient ID Card, stick one Patient Record Label for each implant on the back of the Patient ID Card. Patient Record Labels are located on the internal product packaging attached to the label and contain all the information needed for the product identification: the brand name of the implant, the name and address of the manufacturer, the size of the implant, serial number, batch code, catalog number, expiration date of the implant, patient data, date of implantation, physician's data, indication of which breast the implant is installed (left/right).

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This product has been sterilized by Dry Heat sterilization. Do not resterilize the device.

All MENTOR Gel Breast Implants are indicated for use on one individual patient during a single procedure. If a device is re-used against the manufacturer's instruction, there is a risk of infection (microbial, viral and other transmissible agents) as well as immune responses and the sterility of the device can no longer be guaranteed. Furthermore, the integrity of the device cannot be guaranteed due to the risk of damage to the implant. The established shelf life of the product is compromised and thus null and void if compliance with single-use only indication is not followed. Sterility, safety and efficacy cannot be assured for damaged devices. In the event the product becomes contaminated, contact your local Mentor representative.

REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD

Not applicable.

CONDITIONS OF TRANSPORTATION AND STORAGE

Table 3 – Transportation conditions

Designation	Range
Temperature	-20 °C to +50°C
Relative humidity	not more 75%
Barometric pressure	700-1060 hPa

Table 4 – Storage conditions

Designation	Range
Temperature	+10°C to +25°C

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Gel-Filled Breast Implant Surgery Making an Informed Decision

TABLE OF CONTENTS

1. INTRODUCTION	3
Purpose of the Brochure	3
What Gives the Breast Its Shape?	3
What is a Breast Implant?	3
What is Silicone?	3
Is Silicone Safe?	3
The Independent Review Group (IRG)	3
The Scientific and Technical Options Assessment (STOA) Programme	4
The Institute of Medicine (IOM)	4
2. DECIDING TO HAVE BREAST IMPLANTS	4
Reasons for Considering Breast Implants	4
Reasons Why Breast Implants May Not Be an Option	4
What Types of Silicone Gel-Filled Breast Implants Are Available From Mentor?	5
Factors to Consider in Breast Implantation	6
Potential Benefits	6
Living With Breast Implants	6
Steps to Implantation	7
3. SPECIAL CONSIDERATIONS FOR BREAST AUGMENTATION	9
What Are the Alternatives to Breast Augmentation?	9
Should You Have Breast Augmentation?	9
Surgical Setting and Anesthesia	9
Postoperative Care	9
What Questions Should You Ask Your Surgeon About Breast Augmentation?	9
4. SPECIAL CONSIDERATIONS FOR BREAST RECONSTRUCTION	10
What Are the Alternatives to Breast Reconstruction?	10
Should You Have Breast Reconstruction?	10
What Are the Choices in Reconstructive Procedures?	10
Breast Reconstruction With Breast Implants	10
The Timing of Your Breast Implant Reconstruction	10
Surgical Considerations to Discuss With Your Doctor	11
What Is the Breast Implant Reconstruction Procedure?	11
One-Stage (Immediate or Delayed) Breast Implant Reconstruction	11
Two-Stage (Immediate or Delayed) Breast Implant Reconstruction	11
Stage 1: Tissue Expansion	12
Stage 2: Placing the Breast Implant	12
Breast Reconstruction Without Implants: Tissue Flap Procedures	12
The TRAM Flap (Pedicle or Free)	13
The Latissimus Dorsi Flap With or Without Breast Implants	13
What Questions Should You Ask Your Surgeon About Breast Reconstruction?	14
5. WHAT ARE CONTRAINDICATIONS, WARNINGS AND PRECAUTIONS TO BE CONSIDERED?	14
Surgical practices that are contraindicated in breast implantation	14
Safety has not been established in patients with the following conditions	14
Further considerations	14
Pre-implantation Mammography	14

<i>Interference With Mammography</i>	14
<i>Distinguishing the implant from breast tissue during breast self-examination</i>	15
<i>Long-Term Effects</i>	15
<i>Capsule Procedures</i>	15
<i>Radiation Therapy</i>	15
<i>MRI</i>	15
6. WHAT ARE THE BREAST IMPLANT COMPLICATIONS?	15
Deflation/Rupture	15
Capsular Contracture	15
Pain	16
Additional Surgeries	16
Dissatisfaction With Cosmetic Results	16
Infection	16
Hematoma/Seroma	16
Changes in Nipple and Breast Sensation	16
Breast-Feeding	16
Calcium Deposits in the Tissue Around the Implants	17
Delayed Wound Healing	17
Extrusion	17
Necrosis	17
Breast Tissue Atrophy/Chest Wall Deformity	17
Connective Tissue Disease	17
Cancer	17
Anaplastic Large Cell Lymphoma	17
Second-Generation Effects	18
Gel Bleed and Granulomas	18
7. IF YOU EXPERIENCE A PROBLEM, SHOULD YOU CONTACT YOUR DOCTOR?	18
8. ADDITIONAL RECOMMENDATIONS FOR PATIENTS	18
9. PROBLEM REPORTING FOR AUSTRALIAN PATIENTS	19
10. STATUS/LEGAL POSITION OF BREAST IMPLANTS FOR EUROPEAN PATIENTS	19
11. GLOSSARY	20
PATIENT SIGNATURE FORM	22


1. INTRODUCTION

Purpose of the Brochure

The purpose of this brochure is to assist patients in making informed decisions about breast augmentation and breast reconstruction surgery. This educational brochure is designed to help patients talk with their doctors, as well as provide patients with general information on breast implant surgery and give them specific details about MENTOR® Breast Implants.

What Gives the Breast Its Shape?

The breast consists of milk ducts and glands, surrounded by fatty tissue that provides its shape and feel. Situated beneath the breast is the pectoralis major muscle or chest muscle. Factors such as pregnancy (when milk glands are temporarily enlarged), rapid weight loss, and the effects of gravity as a woman ages combine to stretch the skin, which may cause the breast to droop or sag.

What is a Breast Implant?

Currently, there are three main types of breast implants:

- saline-filled breast implants,
- cohesive silicone gel-filled breast implants and
- combination cohesive gel/saline-filled breast implants.

All of these implants consist of silicone elastomer shells, filled with saline (salt water), cohesive silicone gel or both saline and cohesive silicone gel. The surface of the shells may be either textured or smooth.

Saline-Filled Breast Implant

Gel-Filled Breast Implant

Combination Gel/Saline-Filled Breast Implant

What is Silicone?

Silicones are polymers made up of silicon, oxygen, carbon and hydrogen. Silicone can come in the form of liquid, gel or solid. There are many uses of silicone, including cosmetics, cures for stomach gas and in many medical devices, such as foldable intraocular lenses, heart pacemakers and tissue expanders.

Silicon is a common chemical element occurring in nature. It is found in sand, rocks and glass.

Is Silicone Safe?

There have been several prestigious multispecialty panels of scientific experts who have reviewed and studied the available data on silicone breast implants. The following describes the findings of three such panels.

Continued on Page 2

The Independent Review Group (IRG)

This panel commissioned by the British Minister of Health with experts in many scientific disciplines was charged with reviewing the scientific evidence of the possible health risks associated with silicone breast implants and examining issues related to preoperative patient information. The IRG report released in 1998 made the following conclusions:

- “There is no epidemiological evidence for any link between silicone gel breast implants and any established connective tissue disease. If there is a risk of connective tissue disease, it is too small to be quantified. The IRG cannot justify recommending further epidemiological studies to investigate this hypothesis.” For more information on connective tissue disease, please refer to the section on BREAST IMPLANT COMPLICATIONS.

- The overall biological response to silicone is consistent with conventional forms of response to foreign materials, rather than an unusual toxic reaction.

- There is no evidence that children of women with breast implants are at increased risk of connective tissue disease.

The Scientific and Technical Options Assessment (STOA) Programme

This panel was commissioned by the European Parliament to issue a report on the policy options on silicone gel breast implants based on a comprehensive, unbiased analysis of the scientific literature. The STOA report published in 2000 made the following conclusions:

- “Studies do not point to an association between silicone implants and serious health risks, such as cancer and connective tissue diseases.”

- “... studies have concluded that there does not appear to be evidence of harmful effects to infants breast fed by women with silicone breast implants ...”

The Institute of Medicine (IOM)

This panel commissioned by the US Congress with experts in many scientific disciplines was charged with conducting an independent review of past and ongoing research on silicone and other types of breast implants. The IOM report released in 1999 made the following conclusions:

- There is no increased risk of connective tissue disease in women with breast implants. Women with breast implants are no more likely than the rest of the population to develop cancer, immunologic diseases or neurological problems.

- Breast-feeding is safe and beneficial.

- There are no second-generation effects on children of women with breast implants.

Please note, that Mentor’s evaluation from clinical studies and the latest information to consider are:

- Patients who undergo breast implant surgery with connective tissue diseases (CTD) may experience an increased risk of wound dehiscence, infection and bleeding (likely due to their ongoing medical management) that may require further treatment.

- Some patients with breast implants have reported experiences of neurological and/or rheumatological diseases. Mentor is currently collecting data to further understand these potential risks and their possible association with breast implants.

2. DECIDING TO HAVE BREAST IMPLANTS

Reasons for Considering Breast Implants

Implants are to be used for the following indications:

- **Breast Augmentation** – This procedure is performed to increase the size and proportions of a woman’s breasts.

- **Breast Reconstruction** – This procedure is performed to restore a woman’s breast shape after a mastectomy or injury that resulted in either partial or total loss of the breast(s), or to correct a birth defect.

- **Replacement or Revision** – This procedure is performed as replacement or revision surgery for patients with previous augmentation or reconstruction with silicone gel-filled or saline-filled implants.

Reasons Why Breast Implants May Not Be an Option

The use of this prosthesis is contraindicated in patients who have any of the following conditions:

- Pregnancy or nursing mothers.
- Lupus (e.g., SLE and DLE).
- Scleroderma (e.g., progressive systemic sclerosis).

- Currently has a condition that could compromise or complicate wound healing (except reconstruction patients).
- Infection or abscess anywhere in the body.
- Demonstrates tissue characteristics, which are clinically incompatible with implant (e.g., tissue damage resulting from radiation, inadequate tissue, or compromised vascularity).
- Possesses any condition, or is under treatment for any condition which, in the opinion of the consulting physician(s), may constitute an unwarranted surgical risk.
- Anatomic or physiologic abnormality, which could lead to significant postoperative complications.
- A history of sensitivity to foreign materials or repeated attempts and failures at breast augmentation or reconstruction.
- An unwillingness to undergo any further surgery for revision.
- Unrealistic expectations, such as inappropriate attitude or motivation, or a lack of understanding of the risks involved with the surgical procedure and implants.
- Premalignant breast disease without a subcutaneous mastectomy.
- Untreated or inappropriately treated breast malignancy, without mastectomy.

What Types of Silicone Gel-Filled Breast Implants Are Available from Mentor?

All MENTOR® Gel Breast Implants contain silicone gel which is cohesive. The breast implants are offered with three degrees of cohesive filling material: Cohesive I (standard), Cohesive II (moderate) and Cohesive III (high). The implants come in a variety of shapes, surface textures, and sizes. Breast implants are available with either a smooth or textured surface. Mentor has a textured breast implant surface called SILTEX™. Also, implants are available with a single silicone shell filled with cohesive silicone gel or as expander/breast implants with an inner silicone shell filled with saline and an outer silicone shell filled with cohesive silicone gel.

The following is a description of MENTOR® Breast Implant styles. Be sure to familiarize yourself with the different features of breast implants and to discuss the most appropriate type(s) of implants for you with your surgeon.


Gel-Filled Breast Implant Family (fixed volume)

Round Styles:


Smooth

SILTEX™

CONTOUR PROFILE™

BECKER™ Expander/Breast Implant Family (postoperative adjustment of volume)

Factors to Consider in Breast Implantation

Potential Benefits

Women choose primary breast augmentation surgery to increase the size and proportion of their breast(s). Women choose primary breast reconstruction surgery to replace breast tissue that has been removed because of cancer or injury, or to replace breast tissue that has failed to develop properly because of a severe breast abnormality. In addition, women choose revision surgery (replacement of an existing breast implant) to correct or improve the result of a primary augmentation or primary reconstruction surgery.

Living With Breast Implants

- Whether undergoing augmentation or reconstruction, be aware that breast implantation may not be a one-time surgery. Additional doctor visits are likely needed and additional surgery may be needed over the course of a breast implant patient's life.

- Breast implants are not considered lifetime devices. Implant removal with or without replacement may be needed over the course of a breast implant patient's life.
- The expected duration of MENTOR® Breast Implants is based on the estimated rupture rates of the implants. There are many causes of breast implant rupture, including wear and tear on the device over time, unintentional sharp instrument damage during surgery, and physical stress such as that caused by car accidents or strenuous activity. Rupture rates can vary depending upon the type of breast implant surgery performed. In MENTOR® clinical studies, patients undergo MRI screening on a regular basis to assess whether rupture of the implant has occurred. It is important to note that most patients had no symptoms related to their reported implant rupture.

Below are the estimated rupture rates over time for MENTOR® Breast Implants, as observed in MENTOR® clinical studies. These rates are calculated by the Kaplan-Meier method of analysis. These rupture rates include both ruptures that are "suspect" and those that are "confirmed" as defined below.

Suspected ruptures are ruptures that are suspected due to MRI results but have not been confirmed by physical examination, usually because the patient chose not to have the implant removed.

Confirmed ruptures are confirmed by examining the implant once removed from the patient.

The proportion of confirmed ruptures was 64% (49 of 77 suspected or confirmed ruptures) for MENTOR® MemoryGel® Breast Implants, 47% (8 of 17 suspected or confirmed ruptures) for Contour Profile® Breast Implants (CPG), and 56% (10 of 18 suspected or confirmed ruptures) for BECKERTM Implants.

10-year Kaplan-Meier estimated suspected or confirmed implant rupture rates for MENTOR® MemoryGel® Breast Implants:

Primary Augmentation	Revision Augmentation	Primary Reconstruction	Revision Reconstruction
3.3%	4.7%	11.2%	0%

Overall, the 10-year follow-up rate was 45%, thus data accuracy is limited.

Based on the above information, MENTOR® MemoryGel® and Contour Profile® Breast Implants (CPG) can be expected to have an average lifetime of over ten years.

BECKERTM Expander/Breast Implant devices also contain saline, so deflation can occur if there is a tear in the implant that allows the saline to leak out. Deflation rates and suspected or confirmed rupture rates are shown below.

5-year Kaplan-Meier estimated suspected or confirmed implant rupture rates and deflation rates for BECKERTM Expander/Breast Implants for Reconstruction Patients:

Suspected or confirmed ruptures	7.2%
Deflation	1.3%

Smaller sample sizes in the BECKERTM Implants clinical study limit data accuracy.

BECKERTM Expander/Breast Implant devices have data out to 5 years at this time with less than a 10% rupture/deflation rate for reconstruction patients.

- Many of the changes to a patient's breast following implantation are irreversible (cannot be undone). If a patient later chooses to have her implant(s) removed, unacceptable dimpling, puckering, wrinkling, or other cosmetic changes of the breast may be experienced.
- Breast implants may affect a patient's ability to produce a sufficient supply of milk for breast-feeding. Also, breast implants will not prevent a patient's breast from sagging after pregnancy.
- With breast implants, routine screening mammography may be more difficult, and a patient may need to have additional views, which means more time and radiation.

Steps to Implantation

- **Choosing a Surgeon**

When choosing a surgeon who is experienced with breast implantation, the following questions should be answered:

1. How many breast augmentation or reconstruction implantation procedures does the surgeon perform per year?
2. How many years has the surgeon performed breast implantation procedures?
3. What is the most common complication the surgeon encounters with breast implantation?
4. What is the surgeon's reoperation rate with breast implantation, and what is the most common type of reoperation the surgeon performs?

The following options in breast implant surgery should be understood and discussed with the surgeon:

- **Implant Size**

Generally, the larger the cup size wanted, the larger the breast implant (measured in cubic centimeters, or cc's) the surgeon will consider implanting.

The surgeon will also evaluate the existing tissue to determine if there is enough to cover the breast implant. If a breast implant size too large for the available tissue is desired, the doctor may warn that the breast implant edges may be apparent or visible postoperatively. This may even risk surgical complications. Also, excessively large breast implants may speed up the effects of gravity and result in earlier droop or sag.

- **Surface Texturing**

SILTEX™ is the trade name for the surface texturing on MENTOR® Breast Implants. The textured surface was designed to provide a disruptive or rough surface for the body's collagen tissue to interface with the implant.

Mentor also offers nontextured or smooth-surfaced implants.

- **Implant Shape, Projection and Height**

Mentor offers cohesive gel-filled breast implants with a round or contour profile shape. The round implants are available in various projections (low, moderate, moderate plus, high and ultra high profile) and the contour profile shaped implants are available in various heights (low, medium and tall).

- **Gel Cohesivity Levels**

Mentor's silicone gels are cohesive polymers, not liquids. Although soft and fluid feeling, they act as a unit. Mentor's gels hold together uniformly, while retaining the natural 'give' that resembles actual breast tissue. Mentor's gel breast implants are available in three different levels of cohesivity; however, they are all cohesive.

Cohesive I

The standard cohesive level gel used in Mentor's implants. This is Mentor's softest gel.

Cohesive II

A slightly firmer gel, for a firmer feeling implant.

Cohesive III

Mentor's most cohesive gel, providing shape retention with a pleasing level of firmness.

- **Palpability**

The following may cause implants to be more palpable (more easily felt): textured implants, larger implants, subglandular placement, and the lack of adequate skin/tissue available to cover the implant.

- **Implant Placement**

The breast implant can be placed either partially under the pectoralis major muscle (submuscular) or on top of the muscle and under the breast glands (subglandular). The pros and cons of the implant placement should be discussed with the surgeon.

- **The submuscular placement** may make surgery and recovery last longer, may be more painful, and may make it more difficult to have some reoperation procedures than the subglandular placement. The possible benefits of this placement are that it may result in less palpable implants and easier imaging of the breast with mammography.

The subglandular placement may make surgery and recovery shorter, may be less painful, and may be easier to access for reoperation than the submuscular placement. However, this placement may result in more palpable implants and more difficult imaging of the breast with mammography.

- **Incision Sites**

Subglandular

Submuscular

The pros and cons of the incision site should be discussed with the surgeon. There are three common incision sites: under the arm (axillary), around the nipple (periareolar), or within the breast fold (inframammary). If the incision is made under the arm, the surgeon may use a probe fitted with a miniature camera, along with minimally invasive instruments, to create a "pocket" for the breast implant. A fourth incision site: through the navel (umbilicus) using an endoscopic technique has not been studied and is therefore not recommended.

Axillary


This incision is less concealed than periareolar and associated with less difficulty than the periareolar incision site when breast-feeding.

Periareolar

This incision is most concealed, but is associated with a higher likelihood of inability to successfully breast-feed, as compared to the other incision sites.

Inframammary

This incision is less concealed than periareolar and associated with less difficulty than the periareolar incision site when breast-feeding.

3. SPECIAL CONSIDERATIONS FOR BREAST AUGMENTATION

What Are the Alternatives to Breast Augmentation?

- Accept your breasts as they are.
- Wear a padded bra or external prostheses.

Should You Have Breast Augmentation?

You may consider consulting your family, friends, breast implant support groups, and counselor(s) to help you in making this decision. You are advised to wait a minimum of two to four weeks after reviewing and considering the information in this brochure before deciding whether to have augmentation surgery.

Surgical Setting and Anesthesia

Augmentation surgery is usually performed on an outpatient basis, either in a hospital operating room (theater) or surgery center. General anesthesia is commonly used, although local anesthesia is also an option. The surgery usually lasts one to two hours. Your surgeon will make an incision and create a pocket for the breast implant. Then, the breast implant will be placed in the pocket and positioned. Finally, the incision will be closed, usually with stitches, and possibly taped.

Postoperative Care

You will probably feel somewhat tired and sore for several days following the operation, and your breasts may remain swollen and sensitive to physical contact for a month or longer. You may also experience a feeling of tightness in the breast area as your skin adjusts to the new breast size.

Postoperative care may involve the use of a postoperative bra, compression bandage, or sports bra for extra support and positioning while you heal. At your surgeon's recommendation, you will most likely be able to return to work within a few days, although any strenuous activities that could raise your pulse and blood pressure should be avoided for at least a couple of weeks. Your surgeon may also recommend breast massage exercises.

Note: If you experience fever or noticeable swelling and/or redness in your implanted breast(s), you should contact your surgeon immediately.

What Questions Should You Ask Your Surgeon About Breast Augmentation?

The following list of questions may help to remind you of topics to discuss with your doctor.

1. What are the risks and complications associated with having breast implants?
2. How many additional operations in my implanted breast(s) can I expect over my lifetime?
3. How will my breasts look if I opt to have the implants removed without replacement?
4. What shape, size, surface texturing, incision site, and placement site is recommended for me?
5. How will my ability to breast-feed be affected?
6. How can I expect my implanted breasts to look over time?

7. How can I expect my implanted breasts to look after pregnancy? After breast-feeding?
8. What are my options if I am dissatisfied with the cosmetic outcome of my implanted breasts?
9. What alternate procedures or products are available if I choose not to have breast implants?
10. Do you have before and after photos I can look at for each procedure, and what results are reasonable for me?

4. SPECIAL CONSIDERATIONS FOR BREAST RECONSTRUCTION

What Are the Alternatives to Breast Reconstruction?

You may choose not to undergo breast reconstruction. In this case, you may or may not decide to wear an external breast form (prosthesis) inside your bra. Breast forms are available in a variety of shapes, sizes, and materials such as foam, cotton, and silicone. Custom prostheses are also available to match the size and shape of your breast.

Should You Have Breast Reconstruction?

Whether you decide to have breast reconstruction depends on your own individual case, medical condition, general health, lifestyle, emotional state, and breast size and shape. You may consider consulting your family, friends, breast implant support groups, breast cancer support groups and counselor(s) to help you in making this decision. You are advised to wait a minimum of two to four weeks after reviewing and considering the information in this brochure before deciding whether to have reconstruction surgery.

If you are considering breast reconstruction and do not have a plastic surgeon, ask your general surgeon for a referral for the names of experienced plastic surgeons in your area. Your general surgeon, plastic surgeon, and oncologist should work together to plan your mastectomy and reconstruction procedure to give you the best possible results.

What Are the Choices in Reconstructive Procedures?

The type of breast reconstruction procedure available to you depends on your medical situation, breast shape and size, general health, lifestyle, and goals. Women with small- or medium-sized breasts are the best candidates for breast reconstruction.

Breast reconstruction can be accomplished by the use of a prosthesis (a breast implant, either silicone gel-filled, saline-filled or a combination gel/saline-filled), your own tissues (a tissue flap), or a combination of the two. A tissue flap is a section of skin, fat and/or muscle which is moved from your stomach, back or other area of your body, to the chest area, and shaped into a new breast.

Whether or not you have reconstruction with or without breast implants, you will probably undergo additional surgeries to improve symmetry and appearance. For example, because the nipple is often removed with the breast tissue in mastectomy, the nipple is often reconstructed by using a skin graft from the opposite breast or by tattooing the area. Nipple reconstruction is usually done as a separate outpatient procedure after the initial reconstruction surgery is complete.

Breast Reconstruction With Breast Implants

Your surgeon will decide whether your health and medical condition makes you an appropriate candidate for breast implant reconstruction. Women with larger breasts may require reconstruction with a combination of a tissue flap and an implant. Your surgeon may recommend breast implantation of the opposite, unaffected breast in order to make your breasts more alike (maximize symmetry) or he/she may suggest breast reduction (reduction mammoplasty) or a breast lift (mastopexy) to improve symmetry. Reduction mammoplasty involves removal of breast tissue and skin. Mastopexy involves removing a strip of skin from under the breast or around the nipple and using it to lift and tighten the skin over the breast. If it is important to you not to alter the unaffected breast, you should discuss this with your surgeon, as it may affect the breast reconstruction methods considered for your case.

The Timing of Your Breast Implant Reconstruction

The following description applies to reconstruction following mastectomy but similar considerations apply to reconstruction following breast trauma or for reconstruction for congenital defects. The breast reconstruction process may begin at the time of your mastectomy (immediate reconstruction) or weeks to years afterwards (delayed reconstruction). Immediate reconstruction may involve placement of a breast implant, but typically involves placement of a tissue expander, which is eventually replaced with a breast implant or immediate reconstruction may involve placement of an expander/breast implant (BECKER™ Expander/Implant). A BECKER™ Expander/Implant serves as both a tissue expander and a breast implant. It is expanded over time by adding saline through a filling port located beneath the skin. Once the BECKER™ Expander/Implant has

expanded to the desired fullness, the filling port is removed. It is important to know that any type of surgical breast reconstruction may take several steps to complete. Two potential advantages to immediate reconstruction are that your breast reconstruction starts at the time of your mastectomy and that there may be cost savings in combining the mastectomy procedure with the first stage of the reconstruction. However, there may be a higher risk of complications such as rupture with immediate reconstruction, and your initial operative time and recuperative time may be longer.

A potential advantage to delayed reconstruction is that you can delay your reconstruction decision and surgery until other treatments, such as radiation therapy and chemotherapy, are completed. Delayed reconstruction may be advisable if your surgeon anticipates healing problems with your mastectomy, or if you just need more time to consider your options.

There are medical, financial and emotional considerations to choosing immediate versus delayed reconstruction. You should discuss with your surgeon, plastic surgeon, and/or oncologist, the pros and cons of the options available in your individual case.

Surgical Considerations to Discuss With Your Doctor

Discuss the advantages and disadvantages of the following options with your surgeon and/or your oncologist.

- **Immediate Reconstruction**
One-stage immediate reconstruction with a breast implant (breast implant only) or with an expander/breast implant (BECKER™ Expander/Implant).
Two-stage immediate reconstruction with a temporary tissue expander followed by delayed reconstruction several months later with a breast implant.
- **Delayed Reconstruction**
One-stage delayed reconstruction using an expander/breast implant (BECKER™ Expander/Implant).
Two-stage delayed reconstruction with a temporary tissue expander followed several months later by replacement with a breast implant.

What Is the Breast Implant Reconstruction Procedure?

One-Stage (Immediate or Delayed) Breast Implant Reconstruction

Immediate one-stage breast reconstruction may be done at the time of your mastectomy. After the general surgeon removes your breast tissue, the plastic surgeon will then implant a breast implant or an expander/breast implant (BECKER™ Expander/Implant) that completes the one-stage reconstruction. If a BECKER™ Expander/Implant is used, saline is gradually added over time (up to six months) until the expansion is completed at which time the filling port is removed. The filling port removal is usually done using local anaesthesia at the outpatient or day unit section of the hospital. One-stage delayed breast reconstruction is done several months or years later by placement of a breast implant or a BECKER™ Expander/Implant.


Two-Stage (Immediate or Delayed) Breast Implant Reconstruction

Breast reconstruction with a MENTOR® Gel-Filled Breast Implant usually occurs as a two-stage procedure, starting with the placement of a temporary breast tissue expander, which is replaced several months later with a breast implant.

The tissue expander placement may be done immediately, at the time of your mastectomy, or be delayed until months or years later.

The tissue expander may have a filling port or an injection dome that is integral (attached) to the expander shell or an injection dome that is remote (attached to the expander by a fill tube). If an expander with a remote dome is used, the dome is situated under the arm or axilla.

Stage 1: Tissue Expansion

During a mastectomy, the general surgeon often removes skin as well as breast tissue, leaving the chest tissues flat and tight. To create a breast-shaped space for the breast implant, a tissue expander is placed under the remaining chest tissues.

The tissue expander is a balloon-like device made from elastic silicone rubber. It is inserted unfilled, and over time, sterile saline fluid is added by inserting a small needle through the skin to the filling port of the device. As the tissue expander fills, the tissues over the expander begin to stretch, similar to the gradual expansion of a woman's abdomen during pregnancy. The tissue expander creates a new breast-shaped pocket for a breast implant.

Tissue expander placement usually occurs under general anaesthesia in an operating room (theatre). Operative time is generally one to two hours. The procedure may require a brief hospital stay, or be done on an outpatient basis. Typically, you can resume normal daily activity after two to three weeks.

Because the chest skin is usually numb from the mastectomy surgery, it is possible that you may not experience pain from the placement of the tissue expander. However, you may experience feelings of pressure or discomfort after each filling of the expander, which subsides as the tissue expands. Tissue expansion typically lasts four to six months.

Stage 2: Placing the Breast Implant

After the tissue expander is removed, the breast implant is placed in the pocket. The surgery to replace the tissue expander with a breast implant (implant exchange) is usually done under general anaesthesia in an operating room (theatre). It may require a brief hospital stay or be done on an outpatient basis.

Breast Reconstruction Without Implants: Tissue Flap Procedures

The breast can be reconstructed by surgically moving a section of skin, fat and muscle from one area of your body to another. The section of tissue may be taken from such areas as your abdomen, upper back, upper hip, or buttocks.

The tissue flap may be left attached to the blood supply and moved to the breast area through a tunnel under the skin (a pedicled flap), or it may be removed completely and reattached to the breast area by microsurgical techniques (a free flap). Operating time is generally longer with free flaps, because of the microsurgical requirements.

Flap surgery requires a hospital stay of several days and generally a longer recovery time than implant reconstruction. Flap surgery also creates scars at the site where the flap was taken and possibly on the reconstructed breast. However, flap surgery has the advantage of being able to replace tissue in the chest area. This may be useful when the chest tissues have been damaged and are not suitable for tissue expansion. Another advantage of flap procedures over implantation is that alteration of the unaffected breast is generally not needed to improve symmetry.


The most common types of tissue flaps are the TRAM (Transverse Rectus Abdominus Musculocutaneous flap) (which uses tissue from the abdomen) and the Latissimus Dorsi flap (which uses tissue from the upper back). It is important for you to be aware that flap surgery, particularly the TRAM flap, is a major operation, and is more extensive than your mastectomy operation. It requires good general health and strong emotional motivation. If you are very overweight, smoke cigarettes, have had previous surgery at the flap site, or have any

circulatory problems, you may not be a good candidate for a TRAM flap procedure. Also, if you are very thin, you may not have enough tissue in your abdomen or back to create a breast mound with this method.

The TRAM Flap (Pedicle or Free)

Step 1: Mastectomy is performed and the donor site is marked

Step 2: The flap of rectus muscle and tissue is tunneled to the breast

Step 3: Final Result

During a TRAM flap procedure, the surgeon removes a section of tissue from your abdomen and moves it to your chest to reconstruct the breast. The TRAM flap is sometimes referred to as a “tummy tuck” reconstruction, because it may leave the stomach area flatter.

A pedicle TRAM flap procedure typically takes three to six hours of surgery under general anaesthesia; a free TRAM flap procedure generally takes longer. The TRAM procedure may require a blood transfusion. Typically, the hospital stay is two to five days. You can resume normal daily activity after six to eight weeks. Some women, however, report that it takes up to one year to resume a normal lifestyle. You may have temporary or permanent muscle weakness in the abdominal area. If you are considering pregnancy after your reconstruction, you should discuss this with your surgeon. You will have a large scar on your abdomen and may also have additional scars on your reconstructed breast.

The Latissimus Dorsi Flap With or Without Breast Implants

Post-Mastectomy

View Showing Section of Tissue to Be Removed

Latissimus Dorsi Flap Procedure

During a Latissimus Dorsi flap procedure, the surgeon moves a section of tissue from your back to your chest to reconstruct the breast. Because the Latissimus Dorsi flap is usually thinner and smaller than the TRAM flap, this procedure may be more appropriate for reconstructing a smaller breast.

The Latissimus Dorsi flap procedure typically takes two to four hours of surgery under general anaesthesia. Typically, the hospital stay is two to three days. You can resume daily activity after two to three weeks. You may have some temporary or permanent muscle weakness and difficulty with movement in your back and shoulder. You will have a scar on your back, which can usually be hidden in the bra line. You may also have additional scars on your reconstructed breast.

What Questions Should You Ask Your Surgeon About Breast Reconstruction?

The following list of questions may help to remind you of topics to discuss with your doctor.

1. What are all my options for breast reconstruction?
2. What are the risks and complications of each type of breast reconstruction surgery and how common are they?
3. What if my cancer recurs or occurs in the other breast?
4. Will reconstruction interfere with my cancer treatment?
5. How many steps are there in each procedure, and what are they?
6. How long will it take to complete my reconstruction?
7. How much experience do you have with each procedure?
8. Do you have before and after photos I can look at for each procedure, and what results are reasonable for me?
9. What will my scars look like?
10. What kind of changes in my implanted breast can I expect over time?
11. What kind of changes in my implanted breast can I expect with pregnancy?
12. What are my options if I am dissatisfied with the cosmetic outcome of my implanted breast?
13. Can I talk with other patients about their experiences?
14. For staged reconstruction, what is the estimated total cost of each procedure?
15. How much pain or discomfort will I feel, and for how long?
16. How long will I be in the hospital?
17. Will I need blood transfusions, and can I donate my own blood?
18. When will I be able to resume my normal activity (or sexual activity, or athletic activity)?

5. WHAT ARE CONTRAINDICATIONS, WARNINGS AND PRECAUTIONS TO BE CONSIDERED?

A contraindication is a condition or circumstance that, if present, means a procedure should not be done. MENTOR® Breast Implants are contraindicated in the following circumstances because the risk of undergoing breast augmentation or reconstruction with implants outweighs the benefits:

- Women with active infection anywhere in their bodies,
- Women with existing cancer or precancer of their breast who have not received adequate treatment for those conditions, and in
- Women who are pregnant or nursing.

Surgical practices that are contraindicated in breast implantation:

- Injection through implant shell.
- Stacking of implants: more than one implant per breast.

Safety has not been established in patients with the following conditions:

- Autoimmune diseases such as lupus and scleroderma.
- Conditions that interfere with wound healing and blood clotting.
- A weakened immune system (e.g., currently receiving immunosuppressive therapy).
- Reduced blood supply to breast tissue.

Further considerations

Pre-implantation Mammography

You may wish to undergo a preoperative mammogram and another one 6 months to one year after implantation to establish a baseline.

Interference With Mammography

The implant may interfere with finding breast cancer during mammography and also may make it difficult to perform mammography. Therefore, it is essential that you tell your mammography technologist that you have an implant before the procedure. The technologist can use special techniques to minimize the possibility of rupture and to get the best possible views of the breast tissue. Because the breast is squeezed during mammography, it is possible for an implant to rupture during the procedure. More x-ray views are necessary with these special techniques; therefore, women with breast implants will receive more radiation. However, the benefit of the mammogram in finding cancer outweighs the risk of the additional x-rays.

Distinguishing the implant from breast tissue during breast self-examination

You should perform breast self-examination monthly on your implanted breast. In order to do this effectively, you should ask your surgeon to help you distinguish the implant from your breast tissue. Any new lumps or suspicious lesions (sores) should be evaluated with a biopsy. If a biopsy is performed, care must be taken to avoid puncturing the implant. The presence of breast implants may improve the detection of tumors by self-exam (palpation).

Long-Term Effects

The long-term safety and effectiveness of breast implants are currently being studied; however, Mentor is monitoring the long-term (i.e., 10 years) risk of implant rupture, reoperation, implant removal, and capsular contracture (hardening of the tissues around the implant). As new data becomes available, Mentor will include this information in the product inserts, which are packaged with the products and sent to the doctors. You should contact your doctor about updated information.

Capsule Procedures

You should be aware that closed capsulotomy, the practice of forcible squeezing or pressing on the fibrous capsule around the implant to break the scar capsule, is not recommended as this may result in rupture or breakage of the implant.

Radiation Therapy

You should be aware that testing has not been done to show the effects of radiation therapy on tissues of patients who have breast implants; however, the literature suggests that radiation therapy may increase the likelihood of experiencing capsular contracture. The decision regarding the use of radiation therapy following breast implantation should be made by your surgeon and radiation oncologist.

MRI

The MRI scan has been determined to be the best way to find out if your implant has ruptured without performing surgery. If you undergo an MRI scan, you will have to lie on your stomach with your breast in a special holder. You will then be placed in the machine, which may be open or may be like going into a tunnel. Some patients experience an uneasiness at being in a closed space. While the machine is taking images of your breast, it will make a noise. In order to have an MRI scan, you should not have had metal devices implanted in your body, or a history of claustrophobia.

6. WHAT ARE THE BREAST IMPLANT COMPLICATIONS?

Undergoing any surgical procedure may involve the risk of complications such as the effects of anaesthesia, infection, swelling, redness, bleeding, and pain.

In addition, there are potential complications specific to breast implants. These complications include:

Deflation/Rupture

When the shell of a gel-filled breast implant breaks, the gel can leak from the implant. BECKER™ Expander/Implants deflate when saline solution leaks either through an unsealed or damaged valve, or through a break in the implant shells. This may be noticed by a change in shape of the implant or it might be unnoticed ("silent" rupture) and detected by mammogram or MRI. Deflation or rupture may occur in the first few months after being implanted or after several years. Causes of deflation or rupture include damage by surgical instruments during surgery, overfilling or underfilling of the implant with saline solution, capsular contracture, closed capsulotomy, stresses such as trauma or intense physical manipulation, excessive compression during mammographic imaging, umbilical incision placement, and unknown/unexplained reasons. You should also be aware that the breast implant may wear out over time and deflate/rupture.

Deflated or ruptured implants necessitate additional surgery to remove and to possibly replace the implants.

Capsular Contracture

The scar tissue or capsule that normally forms around the implant may tighten and squeeze the implant and is called capsular contracture. Capsular contracture is more common following infection, hematoma, and seroma. Symptoms range from firmness and mild discomfort, to pain, distortion, palpability of the implant, and/or displacement of the implant.

Additional surgery is needed in cases where pain and/or firmness is severe. This surgery ranges from removal of the implant capsule tissue to removal and possibly replacement of the implant itself. Capsular contracture may recur after these additional surgeries.

Pain

Pain of varying intensity and duration may occur and persist following breast implant surgery. In addition, improper size, placement, surgical technique, or capsular contracture may result in pain associated with nerve entrapment or interference with muscle motion. You should tell your doctor if you experience severe pain.

Additional Surgeries

Women should understand there is a chance they will need to have additional surgery(ies) at some point to replace or remove the implant. Also, problems such as deflation/rupture, capsular contracture, infection, shifting, and calcium deposits can require removal of the implants. Many women decide to have the implants replaced, but some women do not. Those who do not may have cosmetically unacceptable dimpling and/or puckering of the breast following removal of the implant without implant replacement.

Dissatisfaction With Cosmetic Results

Dissatisfying results such as wrinkling, asymmetry, implant displacement (shifting), incorrect size, unanticipated shape, implant palpability, scar deformity, hypertrophic (irregular, raised) scarring, and/or sloshing (with implants containing saline) may occur. Careful surgical planning and technique can minimize but not always prevent such results.

Infection

Infection can occur with any surgery. Most infections resulting from surgery appear within a few days to weeks after the operation. However, infection is possible at any time after surgery. Infections with an implant present are harder to treat than infections in normal body tissues. If an infection does not respond to antibiotics, the implant may have to be removed, and another implant may be placed after the infection is resolved.

In rare instances, Toxic Shock Syndrome has been noted in women after breast implant surgery, and it is a life-threatening condition. Symptoms include sudden fever, vomiting, diarrhea, fainting, dizziness, and/or sunburn-like rash. A doctor should be seen immediately for diagnosis and treatment.

Hematoma/Seroma

A hematoma is a collection of blood inside a body cavity, and a seroma is a collection of the watery portion of the blood (in this case, around the implant or around the incision). Postoperative hematoma and seroma may contribute to infection and/or capsular contracture. Swelling, pain, and bruising may result. If a hematoma occurs, it will usually be soon after surgery, however this can also occur at any time after injury to the breast. While the body absorbs small hematomas and seromas, large ones will require the placement of surgical drains for proper healing. A small scar can result from surgical draining. Implant deflation or rupture can occur from surgical draining if damage to the implant occurs during the draining procedure.

Changes in Nipple and Breast Sensation

Feeling in the nipple and breast can increase or decrease after implant surgery. The range of changes varies from intense sensitivity to no feeling in the nipple or breast following surgery. Changes in feeling can be temporary or permanent and may affect sexual response or the ability to breast-feed.

Breast-Feeding

At this time, it is not known if a small amount of silicone may diffuse (pass through) from the silicone gel-filled breast implant shell and may find its way into breast milk. If this occurs, it is not known what effect it may have on the breast-feeding infant. Although there are no current methods for detecting silicone levels in breast milk, a study measuring silicon (one component in silicone) levels did not indicate higher levels in breast milk from women with silicone-filled gel implants when compared to women without implants. The IOM report concluded that mothers with breast implants should attempt breast-feeding due to the benefits to their children. The periareolar incision site may significantly reduce the ability to successfully breast-feed by decreasing the amount of milk produced.

Calcium Deposits in the Tissue Around the Implants

Deposits of calcium can be seen on mammograms and can be mistaken for possible cancer, resulting in additional surgery to biopsy and/or removal of the implants to distinguish them from cancer.

Delayed Wound Healing

In some cases, the incision site fails to heal normally.

Extrusion

Unstable or compromised tissue covering and/or interruption of wound healing may result in extrusion, which is when the breast implant comes through the skin.

Necrosis

Necrosis is the formation of dead tissue around the implant. This may prevent wound healing and require surgical correction and/or implant removal. Permanent scar deformity may occur following necrosis. Factors associated with increased necrosis include infection, use of steroids in the surgical pocket, smoking, chemotherapy/radiation, and excessive heat or cold therapy.

Breast Tissue Atrophy/Chest Wall Deformity

The pressure of the breast implants may cause the breast tissue to thin and shrink. This can occur while implants are still in place or following implant removal without replacement.

In addition to these common complications, there have been concerns with certain diseases, of which you should be aware:

Connective Tissue Disease

Concern over the association of breast implants to the development of autoimmune or connective tissue diseases, such as lupus, scleroderma, or rheumatoid arthritis, was raised because of cases reported in the literature in small numbers of women with implants. A review of several large epidemiological studies of women with and without implants indicates that these diseases are no more common in women with implants than in women without implants. However, some women with breast implants believe that their implants caused a connective tissue disease.

Cancer

Published studies indicate that breast cancer is no more common in women with implants than in those without implants.

Anaplastic Large Cell Lymphoma

Based on information reported to global regulatory agencies and found in medical literature, an association has been identified between breast implants and the development of anaplastic large cell lymphoma (ALCL), a type of non-Hodgkin's lymphoma. Women with breast implants have a very small but increased risk of developing Breast Implant Associated ALCL (BIA-ALCL) in the fluid or scar capsule adjacent to the implant, with documented potential for local, regional, and distant spread of the cancer with mortality reported in rare cases. BIA-ALCL has been reported globally in patients with an implant history that includes Mentor's and other manufacturers' breast implants with various surface properties, styles, and shapes. Most of the cases in the literature reports describe a history of the use of textured implants. Several journal articles explore the risk factors for BIA-ALCL, including the varied methods used to create surface texture of the implant and the role of biofilm in causing disease, among others.

You should consider the possibility of BIA-ALCL when a patient presents with a late onset seroma (at least 1 year and on average 7 to 10 years following implantation), breast mass adjacent to the implant capsule or lymphadenopathy. Work-up of delayed enlargement of a previously implanted breast includes an ultrasound or MRI to look for the presence of seroma or lymphoma. Diagnosis starts at presentation, with needle aspiration of as much of the fluid as possible (minimum 50 mL), using ultrasound for implant displacement and protection. A suspicious mass requires needle or open tissue biopsy. Specimens should be sent for cell morphology by cytology, CD30 immunohistochemistry and ALK (anaplastic lymphoma kinase), and flow cytometry for evaluation, quantification, and characterization of T cells within the specimen. Hematopathology consultation

at a tertiary cancer center is strongly encouraged to establish or exclude a diagnosis of BIA-ALCL. Once the diagnosis of BIA-ALCL has been established, physicians are strongly encouraged to consult with a multidisciplinary team including oncologists, pathologists, surgical oncologists, and plastic surgeons. The United States National Comprehensive Cancer Network (NCCN) guidelines for treatment of established BIA ALCL includes the removal of the implant and the entire surrounding implant capsule in addition to associated masses and lymph node involvement.

Please report all confirmed cases of BIA-ALCL per the details provided in the ADVERSE EVENT REPORTING.

FDA provides additional information regarding questions and answers for BIA-ALCL, please visit:

<http://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/ImplantsandProsthetics/BreastImplants/ucm239995.htm>

Please note for latest diagnostic and treatment options refer The National Comprehensive Cancer Network (NCCN) guidelines.

Second-Generation Effects

There have been concerns raised regarding potential damaging effects on children born of mothers with implants. The IOM report concluded that there are no second-generation effects on children of women with breast implants.

Gel Bleed and Granulomas

The gel in an implant consists of large three-dimensional mesh-like molecules that constitute about 20% of the total weight of the gel. The spaces within these molecules are filled with a mix of medical grade silicone oils. These oils are similar to the materials available in many products including anti-gas medication, available without prescription for infants and adults. A small amount of these oils can bleed through the shell of the implant. The major portion of this oil stays on the implant wall. A smaller amount moves into the scar capsule where it is gradually picked up by certain body cells, called macrophages, the probable mechanism by which it is transported into regional lymph nodes.

It is possible for a granuloma to form around a tiny amount of silicone. Although these lumps are noncancerous, they can be difficult to distinguish from cancerous lumps without being removed (biopsied) and examined.

7. IF YOU EXPERIENCE A PROBLEM, SHOULD YOU CONTACT YOUR DOCTOR?

Consult a physician if you suspect any complications with your breast implant(s), in particular in the case of trauma or compression caused, for example, by extreme massaging of the breast region, by some sport activities or by using seat belts.

8. ADDITIONAL RECOMMENDATIONS FOR PATIENTS

The following are additional recommendations for patients with breast implants:

- Patients should return to their surgeon for necessary postoperative visits.
- Consult a physician or pharmacist before the use of topical medicines such as steroids in the breast area.
- Continue to consult a physician for routine examinations to detect for breast cancer.
- Inform a physician or surgeon of the presence of an implant if any surgery of the breast area is scheduled.
- Keep the Patient ID Card provided by your surgeon (with the style and lot number of your breast implant[s]) permanently and carry it with you to facilitate medical care in case of emergency (e.g., in case of a road accident).

9. PROBLEM REPORTING FOR AUSTRALIAN PATIENTS

For Australian patients, if you experience difficulties with your implant(s), you are encouraged to report those difficulties to the Therapeutic Goods Administration (TGA) via the Device Incident Reporting Scheme. For information contact: Reply Paid 32, Medical Device Incident Report Investigation Scheme, P.O. Box 100, Woden ACT 2606, phone 1800 809 361, email: iris@health.gov.au.

10. STATUS/LEGAL POSITION OF BREAST IMPLANTS FOR EUROPEAN PATIENTS

All breast implants sold within the European Community are governed under the Medical Device Directive (93/42/EEC) and are Class III devices according to the Commission Directive 2003/12/EC of 3 February 2003.

11. GLOSSARY

ALCL – Anaplastic Large Cell Lymphoma – a rare type of non-Hodgkin's lymphoma (cancer of the immune system).

Areola – The pigmented or darker colored area of skin surrounding the nipple of the breast.

Asymmetry – Lack of proportion of shape, size and position between the two breasts.

Autoimmune disease – A disease in which the body mounts an "attack" response to its own tissues or cell types. Normally, the body's immune mechanism is able to distinguish clearly between what is a normal substance and what is foreign. In autoimmune diseases, this system becomes defective and produces antibodies against normal parts of the body, causing tissue injury. Certain diseases such as rheumatoid arthritis and scleroderma are considered to be autoimmune diseases.

Axillary – Pertaining to the armpit area.

Bilateral – Relating to, or affecting, the right and left breast.

Biopsy – The removal and examination of tissue, cells or fluid from the body.

Breast augmentation – A surgical procedure that increases the size and proportions of a woman's breast.

Breast reconstruction – A surgical procedure that restores the natural breast contour and mass following mastectomy, trauma, or injury.

Capsular contracture – A tightening of the scar tissue surrounding an implant, resulting in firmness or hardening of the breast.

Capsulectomy – Surgical removal of the capsule (scar tissue).

Capsulotomy (closed) – A breakage in the capsule (scar tissue) by massage or compression on the outside of the breast.

Capsulotomy (open) – Incision or opening in the capsule (scar tissue) made by an open surgical approach.

Carcinoma – A malignant (cancerous) tumor.

Congenital anomaly – A deviation from a normal body part, existing at or dating from birth.

Connective tissue disease – A disease or group of diseases affecting connective tissue. The cause of these diseases is unknown. The diseases are grouped together on the basis of clinical signs, symptoms, and laboratory abnormalities.

Deflation/rupture – Leakage of saline solution from the implant, often due to a valve leak or a tear or cut in the implant shell, with partial or complete collapse of the implant.

Delayed reconstruction – Breast reconstruction that takes place weeks, months, or years after a mastectomy.

Displacement – Movement from the usual or proper place.

Epidemiological – Relating to the incidence, distribution and control of disease in a population.

Extrusion – The pressing out of the implant through the surgical wound.

Fibrous tissues – Connective tissues composed mostly of fibers.

Flap – A portion of tissue (which may include muscle, fat and skin) with its blood supply moved from one part of the body to another.

Hematoma – A mass or swelling containing blood.

Immune response – A bodily response to the presence of a foreign substance.

Inframammary – Below the breast.

Inframammary fold/incision – An incision made in the fold below the breast.

Inpatient surgery – A surgical procedure in which the patient is required to stay overnight in the hospital.

Latissimus dorsi – Two triangular muscles running from the spinal column to the shoulder.

Mammoplasty – Plastic surgery of the breast.

Mammary – Pertaining to the breast.

Mammography – X-ray examination of the breasts (as for early detection of cancer).

Mastectomy – The removal of breast tissue due to the presence of a cancerous or precancerous growth.

Mastopexy – Plastic surgery to move sagging breasts into a more elevated position.

Modified radical mastectomy – Surgical removal of the entire breast including the nipple, areola, and overlying skin, as well as the lymphatic-bearing tissue in the axilla.

Necrosis – Death of living tissue.

Oncologist – A doctor specializing in the study of tumors.

Outpatient surgery – A surgical procedure in which the patient is not required to stay in the hospital overnight.

Palpability – Capability of being touched or felt.

Palpate/palpability – To feel with the hand.

Pectoralis – Major muscle of the chest.

Plastic surgery – Surgery intended to repair, restore, or improve the body following trauma, injury, or illness.

Prosthesis – Any artificial device used to replace or represent a body part.

Ptosis – Breast sagging that is usually the result of normal ageing, pregnancy, or weight loss.

Radical mastectomy – Surgical removal of the entire breast including the nipple, areola, and overlying skin, as well as the pectoral muscles, lymphatic-bearing tissue in the axilla, and various other neighboring tissue.

Rectus abdominus – A long, flat muscle extending the whole length of the front of the abdomen (stomach).

Saline – A solution that is made up of water and a small amount of salt. Approximately 70% of an adult's body weight consists of this saltwater solution.

Seroma – Accumulation of fluid in tissue.

Serratus – Muscle located beneath the chest's pectoralis major and minor muscles and the rib cage.

Silicone elastomer – A type of silicone that has elastic properties similar to rubber.

Subcutaneous mastectomy – Surgical removal of the breast tissues, but sparing the skin, nipple, and areola.

Subglandular placement – Placement underneath the mammary gland and on top of the chest muscle.

Submuscular placement – Placement wholly or partially underneath the pectoralis major (chest) muscle.

Surgical incision – A cut or wound of body tissue made during surgery.

Tissue expander – An adjustable implant that can be inflated with saltwater to stretch the tissue at the mastectomy site to create a new tissue flap for implantation of the breast implant.

Total mastectomy – Surgical removal of the breast including the nipple, areola, and most of the overlying skin.

Transaxillary incision – An incision made across the long axis of the armpit.

Umbilical – Relating to the navel.

Unilateral – Affecting only one side.

**SILTEX™ AND SMOOTH GEL BREAST IMPLANTS,
SILTEX™ BECKER™ EXPANDER/BREAST IMPLANTS**

CE
2797

PATIENT SIGNATURE FORM

It is the responsibility of the individual surgeon to decide the best method by which to counsel a patient prior to surgery.

Mentor relies upon the surgeon to advise the patient of all potential complications and risks associated with the use of breast implants.

The surgical procedures associated with the use of gel breast implants are not without potential complications and risks. The use of this product is an elective procedure. I (the patient) have been counselled prior to surgery regarding the benefits and possible risks associated with tissue reconstruction and/or breast augmentation using breast prostheses and alternative procedures. I (the patient) have been advised by my doctor that breast implants should not be considered lifetime implants.

I was provided with the brochure, "Gel-Filled Breast Implant Surgery: Making an Informed Decision" in advance prior to surgery (except in cases of emergency surgery for cancer or other reasons).

I have read this brochure. Yes ☐ No ☐

The information from the brochure was discussed with me by my surgeon and understood by me. All my questions have been answered to my satisfaction and I have been provided a copy of this form.

Patient's Signature

Date

Patient's Printed Name

Physician's Signature

Date

Physician's Printed Name

(Note to Physician: A copy of the signed form is to be kept in the patient's file.)

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