THORSTEN WOLF  
NOTARY PUBLIC

## CERTIFIED COPY

I HEREBY CERTIFY THAT THE FOLLOWING IS A TRUE COPY OF THE ORIGINAL.

D-OFFENBACH AM MAIN, THE 24TH DAY OF OCTOBER 2024

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
  
Signature, L.S.

## Operational Documentation for the medical device:

### Transaroll Patient Transfer Aids

**Manufacturer:**

Transatlantic Handels-Gesellschaft Stolpe & Co. mbH, Germany

Rev. 2.2

2024

## **General information:**

The Operational Documentation consists of three Instructions for use of the medical device "Transaroll Patient Transfer Aids", each of which is intended for certain designs of the medical device. The Instructions for use contain all the necessary and detailed information on the proper and safe use, conditions of use, and specifications of the medical device.

Instructions for use of each design of the medical devices are placed in the packaging of the corresponding devices.

Instructions for use are available in different languages, including Russian.

## **Table of contents**

### **I. Instruction for use - Transaroll Patient Transfer Aids, in designs**

1. Transaroll Silverboard Mini, foldable
2. Transaroll Silverboard Mini, non-foldable
3. Transaroll Silverboard Standard, foldable
4. Transaroll Silverboard Special, foldable
5. Transaroll Silverboard Special, non-foldable

### **II. Instruction for use - Transaroll Patient Transfer Aids, in designs**

1. Transaroll Blueboard Mini, foldable
2. Transaroll Blueboard Mini, non-foldable
3. Transaroll Blueboard Standard, foldable
4. Transaroll Blueboard Long, foldable

### **III. Instruction for use - Transaroll Patient Transfer Aids, in designs**

1. Transaroll Hygieneboard Mini, foldable, with accessories – 1 pc.:

Accessories:

Disposable sliding film, design:

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs

2. Transaroll Hygieneboard Mini, non-foldable, with accessories – 1 pc.:

Accessories:

Disposable sliding film, design:

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs

3. Transaroll Hygieneboard Standard, foldable, with accessories – 1 pc.:

Accessories:

Disposable sliding film, design:

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs
- 53 \* 180 cm, 60 pcs per roll, perforated - from 1 to 10 pcs

## Instruction for use

### Transaroll Patient Transfer Aids,

in various designs


**Transaroll Silverboard Mini, foldable**

**Transaroll Silverboard Mini, non-foldable**

**Transaroll Silverboard Standard, foldable**

**Transaroll Silverboard Special, foldable**

**Transaroll Silverboard Special, non-foldable**

Please read the usage instructions carefully before first use!

**Symbols:**

|                                                                                   |                                    |                                                                                   |                                                                                                                              |                                                                                    |                                                             |
|-----------------------------------------------------------------------------------|------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------|
|  | Bed-to bed transfer                |  | Max. weight load                                                                                                             |  | Observe operating instructions                              |
|  | Warning                            |  | Manufacturer                                                                                                                 |  | Mark of conformity                                          |
|  | Temperature range                  |  | Keep away from sunlight                                                                                                      |  | Medical device                                              |
|  | Catalogue number                   |  | Batch no./Lot no.                                                                                                            |  | The bridging of free transfer levels must not exceed 15 cm. |
|  | Keep dry                           |  | Slip hazard                                                                                                                  |  | Latex-free                                                  |
|  | UDI - Unique Device Identification |  | NFC ( Near Field Communication): devices with NFC support technology can exchange the data at close range without connecting |                                                                                    |                                                             |

**Manufacturer:**

**Transatlantic Handels-Gesellschaft Stolpe & Co. mbH**

Siemensstraße 7  
61267 Neu-Anspach, Germany  
Tel.: +49 (0)6081 9430-50  
Fax: +49 (0)6081 9430-80  
info@transat.de  
[www.transatlantic.de](http://www.transatlantic.de)

Abbreviated name: Transatlantic

**Authorized representative of the manufacturer in the Russian Federation, providing warranty/post-warranty service and support:**

**BioSystems LLC**

Russia, 197022, Saint Petersburg, 23 Professora Popova Street, Let. E  
Tel.: (812) 320-49-49, Fax (812) 320-49-40,  
e-mail: main@bioline.ru  
[www.bioline.ru](http://www.bioline.ru)  
Hotline 8-800-444-33-49

## **Name of the medical device and its designs:**

### **Transaroll Patient Transfer Aids**

#### **I. Various designs:**

1. Transaroll Silverboard Mini, foldable – 1 pc.
2. Transaroll Silverboard Mini, non-foldable – 1 pc.
3. Transaroll Silverboard Standard, foldable – 1 pc.
4. Transaroll Silverboard Special, foldable – 1 pc.
5. Transaroll Silverboard Special, non-foldable – 1 pc.
6. Transaroll Blueboard Mini, foldable – 1 pc.
7. Transaroll Blueboard Mini, non-foldable – 1 pc.
8. Transaroll Blueboard Standard, foldable – 1 pc.
9. Transaroll Blueboard Long, foldable – 1 pc.
10. Transaroll Hygieneboard Mini, foldable, with accessories – 1 pc.:

##### **Accessories:**

##### **Disposable sliding film, in designs:**

- 60 cm \* 100 m, a roll – from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs

11. Transaroll Hygieneboard Mini, non-foldable, with accessories – 1 pc.:

##### **Accessories:**

##### **Disposable sliding film, in designs:**

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs

12. Transaroll Hygieneboard Standard, foldable, with accessories – 1 pc.:

##### **Accessories:**

##### **Disposable sliding film, in designs:**

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs
- 53 \* 180 cm, 60 pcs per roll, perforated - from 1 to 10 pcs

#### **II. Instruction for use**


Hereafter referred to as the Transaroll Patient Transfer Aids, the board, the medical device.

## 1. Technical Specifications

The principle of the Patient Transfer Aids operation is based on the movement of two interacting surfaces relative to each other, one of which is in a static state, and the other rotates around it.

The Transaroll Silverboard consists of 3 components - basic body, undercover and cover. The blue undercover fits tightly onto the basic body, the silver cover fits over the blue undercover and rotates freely around the basic body. The basic body may be foldable depending on the board design.

For the foldable board design, the undercover is divided into two parts along the folding line to prevent the cover from twisting and wrinkling. Both covers are removable and reusable.

### 1.1 Main technical specifications

| REF     | Item designation / designs                   | Dimensions *, L x W x T, mm | Weight*, kg | Max. rolling load, kg | Max. tension load**, kg | UDIDI                           |
|---------|----------------------------------------------|-----------------------------|-------------|-----------------------|-------------------------|---------------------------------|
| 4500028 | Transaroll Silverboard Mini, foldable        | 950 x 500 x 22              | 1.300       | 200                   | 300                     | (01)04050856034557(11)x x(10)xx |
| 4500029 | Transaroll Silverboard Mini, non-foldable    | 950 x 500 x 22              | 1.300       | 200                   | 300                     | (01)04050856034373(11)x x(10)xx |
| 4500030 | Transaroll Silverboard Standard, foldable    | 1800 x 500 x 22             | 2.200       | 200                   | 300                     | (01)04050856034366(11)x x(10)xx |
| 4500033 | Transaroll Silverboard Special, foldable     | 1500 x 600x 22              | 1.750       | 200                   | 300                     | (01)04050856043900(11)x x(10)xx |
| 4500034 | Transaroll Silverboard Special, non-foldable | 1500 x 600x 22              | 1.750       | 200                   | 300                     | (01)04050856043924(11)x x(10)xx |

\* Dimensions and weight are indicated taking into account the acceptable tolerance of  $\pm 10\%$

\*\* Lifting for a short time no more than 30 seconds

The force applied by medical personnel to move the board horizontally with the patient placed on it:

| Design                                    | Maximum pushing force, no more than: |         | Maximum traction force, no more than: |         |
|-------------------------------------------|--------------------------------------|---------|---------------------------------------|---------|
| Load                                      | 150 kg                               | 200 kg  | 150 kg                                | 200 kg  |
| For all designs of Transaroll Silverboard | 295.0 N                              | 325.0 N | 342.0 N                               | 360.5 N |

## 1.2. Material

| Part of the device                | Material                         | Color  | Contact type              | Property    |
|-----------------------------------|----------------------------------|--------|---------------------------|-------------|
| Transaroll Silverboard basic body | Polyethylene, latex-free         | Black  | Direct short-term contact | Radiolucent |
| Cover                             | PES fabric with PVC coating      | Silver |                           |             |
| Undercover                        | PES fabric with silicone coating | Blue   |                           |             |

## 1.3 Packaging

The board is packaged with instructions for use in a clear plastic bag and then placed in a brown cardboard box and sealed with a Transatlantic adhesive tape.

Individual packaging is also a shipping packaging.

Box size varies depending on the board dimensions:

| REF     | Item designation / variants                  | Box size*, mm   | Gross weight*, kg |
|---------|----------------------------------------------|-----------------|-------------------|
| 4500028 | Transaroll Silverboard Mini, foldable        | 975 x 520 x 57  | 2.3               |
| 4500029 | Transaroll Silverboard Mini, non-foldable    | 975 x 520 x 57  | 2.3               |
| 4500030 | Transaroll Silverboard Standard, foldable    | 920 x 535 x 75  | 3.1               |
| 4500033 | Transaroll Silverboard Special, foldable     | 750 x 600 x 70  | 2.7               |
| 4500034 | Transaroll Silverboard Special, non-foldable | 1500 x 600 x 70 | 3.0               |

\* Dimensions and gross weight (in a shipping packaging) are indicated taking into account the acceptable tolerance of  $\pm 10\%$

## 1.4 Use, storage and transportation conditions

Use:

|                         |                          |
|-------------------------|--------------------------|
| Operating temperature*: | +10°C+30°C               |
| Humidity:               | Not determined, keep dry |

Storage and transportation:

|                             |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Storage temperature*:       | Board: +10°C+30°C                                                      |
| Storage conditions:         | Store in a dry place, protected from moisture; keep away from sunlight |
| Humidity:                   | Board: Not determined, keep dry                                        |
| Transportation conditions*: | +10°C+30°C                                                             |

\*Storage and transportation conditions are set by the manufacturer to preserve the structural behaviour of the basic body and the sliding properties of the film.

Transportation is allowed by any type of covered means of transport subject to the recommended temperature conditions.

---

## 2. Indication / areas of application / intended use

---

**Indication:** No specific clinical picture can be given as an indication because of the various possible applications. The Transaroll Silverboard is suitable for all applications on persons requiring assistance with lateral transfers.

**Areas of application:** The Silverboard can be used to bridge gaps when two lying surfaces are positioned next to each other (e.g. nursing bed - nursing bed, nursing bed - shower trolley, nursing bed - operating table, stretcher - examination table).

**Intended use:** The Transaroll Silverboard is suitable for gentle, pain-free and low-friction lateral transfer (by sliding) of persons whose movement is limited.

**Potential users:** No special education or knowledge is required to use this board. Please read the instructions carefully before use. Potential users are junior and middle grade medical staff after reading the instruction for use.

---

## 3. Redeployment

---

The Transaroll Silverboards are intended for reuse. They are not a disposable medical device. It is essential to conduct a thorough functional check and to clean or disinfect the basic structure before using it again.

---

## 4. Contraindications / risks of use

---

There are no known contraindications.

Before using the board, users must carefully read the Instruction for use, which describes in detail the procedure and features of handling the device, including visual drawings. If users do not read the instructions for handling the Transaroll Silverboard, the board may be handled incorrectly, which could result in damage to the health of the patient and/or user. Observe the application instructions under point 5 and the safety instructions under point 7.

The application risk analysis shows that the provided preventive measures do not require additional risk management measures, provided that the handling rules and safety regulations described in this instruction are observed.

It is the sole responsibility of the user to determine the suitability of the board for the specific application. Transatlantic shall not be held accountable for any improper use of Transaroll products.

## 5. Use / placement of the board

Position the board to the side, next to the patient.

Position yourself next to the patient, find a support point, and, resting on the patient's shoulder and hip, push him in the transfer direction.

Ensure that the brake is on!

To remove the board, reposition the patient in a short-term gentle lateral position.

Tilt the patient to the side gently and place the edge of the board under the patient so that the patient's head, shoulder girdle, and pelvic region are on the edge of the board.

Keeping the patient in lateral position, slide the board from under the patient.

The Transaroll board is radiolucent, has no metal parts and can be used during X-ray, MRI, CT procedures without the need to reposition the patient.

## 6. Service life / life cycle

The expected service life of the Transaroll Silverboard is 5 years if it is used in accordance with the safety instructions. There is no limitation to this period.

A longer or shorter service life is possible depending on the load. Usability may also end prematurely in the event of relevant damage to the patient transfer aids.

## 7. Safety instructions / functional check

**WARNING** signifies a possible hazard that could lead to personal injury or material damage.

- Defective products must no longer be used.
- Slip hazard! Ensure that the patient does not slide off the bed, recliner, etc. during transfer or repositioning.
- Risk of accident! Always make sure when using the Transaroll Silverboard that the bed, recliner or any movable base is stable, i.e. that any brakes on it are locked.
- Ensure that the Transaroll Silverboard is stable on both sides at the transfer levels. The bridging of free transfer levels must not exceed 15 cm.

- It is the responsibility of the user/institution to evaluate if multiple users may be needed for the transfer or if additional safety measures must be taken.
- Never leave the patient unsupervised on the Transaroll Silverboard.
- The Transaroll Silverboard must not remain under the patient after transfer.
- Pay attention, particularly when lifting the Transaroll Silverboard for a short time, that the patient is lying securely on the board and cannot slide down.
- Make sure that the transfers are performed by the users in a way that is back-friendly.
- **Weight limits must be observed (200 kg).**

Repositioning products are also subject to a natural ageing process. Inspect the device regularly and carefully. Ensure that the device is in impeccable condition and that the material has not been damaged (no tears or cuts in the cover, joint opening, tears, cuts, or spalling of the basic body).

**Defective products must no longer be used!**

**Precautions:** do not use the board if it is damaged.

Unsuitability criteria: there is damage to the cover and/or undercover (tears, cuts, joint opening, adhesion of covers to each other) or to the basic body of the board (tears, cuts, spalling).

The Transaroll Patient Transfer Aids are subject to the highest quality requirements.

The Transaroll Silverboard Patient Transfer Aids do not have accessories, are not intended for use in combination with other medical or non-medical products and are used independently.

The board does not contain any medicines, materials of animal or human origin.

The board is ready for use immediately after unpacking and does not require installation, setting up or other steps for commissioning.

The board is not sterile and does not require sterilization before use.

The board is not a source of any radiation and therefore does not affect other medical devices. It is also not affected by external factors from other medical devices.

Any serious incidents must be reported to the manufacturer and the competent authority in the Member State where the user and/or patient is based.

---

## 8. Cleaning, care and disinfection

---

The silver cover must be cleaned with medical disinfectants prior to and after use each time a patient is changed.

Disinfection is carried out by wiping with a cloth moistened with a disinfectant or spraying a disinfectant solution onto the surface. Surfaces to be disinfected must be left to dry completely.

### 8.1. Undercover (blue)

When cleaning, do not use abrasive or bleaching agents and objects with sharp edges. The use of bleaching agents or cleaning agents, including those that contain alcohol, can permanently impair the functional properties of the blue undercover. The blue undercover can be easily disinfected using oxygen-based disinfectant solutions or solutions containing alkyldimethylbenzylammonium chloride, polyhexamethylene guanidine chloride.

**For users in the Russian Federation**, it is recommended to use for disinfection - Detergent disinfectant "DeMoS", produced by LiG LLC, Russia (certificate of state registration RU 77.99.01.002.E.035064.09.11) or similar products based on alkyldimethylbenzylammonium chloride and/or polyhexamethylene guanidine chloride, approved for use.

If necessary, the undercover can be removed and disinfected separately.

### 8.2. Cover (silver)

When cleaning, do not use abrasive or bleaching agents and objects with sharp edges. The use of bleaching agents can permanently impair the functional properties of the cover. The silver cover can be easily disinfected

using alcohol-based disinfectant solutions or solutions containing alkyldimethylbenzylammonium chloride, polyhexamethylene guanidine chloride.

**For users in the Russian Federation**, it is recommended to use for disinfection:

- Detergent disinfectant "DeMoS", produced by LiG LLC, Russia (certificate of state registration RU 77.99.01.002.E.035064.09.11) or similar products based on alkyldimethylbenzylammonium chloride and/or polyhexamethylene guanidine chloride, approved for use.
- Disinfectant "Avandez-A", produced by Dezos LLC, Russia (certificate of state registration RU 77.99.32.002.E.002145.06.21) or similar products based on ethyl, isopropyl or n-propyl alcohols, approved for use.

If necessary, the cover can be removed and disinfected separately.

### 8.3. Transaroll basic body

The board itself withstands a lot of wear and tear and can be cleaned using normal household cleaners.

Conventional cleaning: Use warm water and a mild cleaning agent for cleaning. If necessary, wipe the board using a damp cloth and then leave it to dry thoroughly.

Disinfection: The board can be effortlessly disinfected using an alcohol-based disinfectant solution (based on isopropyl, n-propyl, ethyl or other alcohols). Always leave to dry thoroughly.

The board is unsuitable for use in an autoclave due to the thermo-sensitive structure of the material.

**For users in the Russian Federation**, it is recommended to use for disinfection:

- Disinfectant "Avandez-A", produced by Dezos LLC, Russia (certificate of state registration RU 77.99.32.002.E.002145.06.21) or similar products based on ethyl, isopropyl or n-propyl alcohols, approved for use.

---

## 9. Product marking /CE

---

The Transaroll Silverboard is a Class I medical device and is compliant with EU Regulation 2017/745 of the European Parliament and of the Council dated 05.04.2017 regarding medical devices (MDR).

We are certified as a manufacturer and distributor of medical devices in accordance with DIN EN ISO 9001:2015 as well as DIN EN ISO 13485:2016. Our medical products have, of course, passed the relevant conformity assessment procedure and therefore bear the CE mark.

**Individual product marking** consists of two film labels of a special shape corresponding to the board form, which [labels] are placed on the top and bottom sides of the board and contain the following information:

Information on the top label: name and address of the manufacturer, catalog number, batch code, manufacturing date, maximum safe load, symbols of use, storage and transportation conditions, caution symbols, instruction symbol, medical device symbol, risk class, UDI symbol, DataMatrix code, CE marking, (details see in table "Symbols"), label revision (rev.), company logo, country of origin (Made in Germany).

Information on the bottom label: name of the design, manufacturing date, label revision (rev.), QR code, NFC symbol, company logo.


The board marking contains a symbol indicating that there is an NFC tag containing information on the board serial number. To read this information, use the pre-installed NFC application on your smartphone (if available). To do this, you need to open the NFC application and bring your smartphone no more than 10 cm away from the NFC symbol on the board; information on the serial number will be displayed on the smartphone screen.

Marking on the board has a DataMatrix code: when reading the code, the following saved information is displayed: (01) GTIN (Global Trade Item Number GS1), (11) manufacturing date and (10) batch code.

Marking on the board has a QR code: when reading, it contains a link to the manufacturer's website.

#### Sample of board marking

#### Marking on the board packaging

The label contains the following information: name and address of the manufacturer, catalog number, batch code, name of the design, size, foldable/non-foldable, manufacturing date, symbols of storage and transportation conditions, instruction symbol, medical device symbol, risk class, UDI symbol, DataMatrix code, CE marking, (see details in table "Symbols"), label revision (rev.), company logo, country of origin (Made in Germany), web-page.

#### Sample of packaging marking

#### Draft of additional marking on the packaging for the Russian Federation

In order to comply with Russian legislation, the marking in the form of a Russian-language stickers will be applied in addition to the manufacturer's marking.

The marking in the form of a Russian-language sticker is additional to the manufacturer's marking, which contains the catalog number, lot/batch number and warning symbols. Additional Russian-language marking is applied next to the manufacturer's marking, without overlapping it. The data thereon complement each other and must be taken into account as a whole.

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General view of marking using an example of Transaroll Silverboard Standard, foldable

#### 10. Standards and specifications applicable in whole or in part

|                             |                                                                                                                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------|
| DIN EN ISO 13485:2021-12    | Medical devices - Quality management systems - Requirements for regulatory purposes                                        |
| DIN EN ISO 14971:2020-07    | Medical devices - Application of risk management to medical devices                                                        |
| DIN EN ISO 20417:2022-03    | Medical devices - Information to be supplied by the manufacturer                                                           |
| DIN EN ISO 15233-1:2022-02  | Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements    |
| DIN EN ISO 10993-1:2021-05  | Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing within a risk management process |
| DIN EN ISO 10993-10:2014-10 | Medical devices - Biological evaluation of medical devices. Part 10: Tests for irritation and skin sensitization           |
| DIN EN ISO 10993-11:2018-09 | Medical devices. Biological evaluation of medical devices. Part 11. Tests for systemic toxicity                            |

#### For the Russian Federation:

##### Russian standards

|                         |                                                                                                                         |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------|
| GOST R 50444-2020       | Medical instruments, apparatus and equipment. General technical requirements                                            |
| GOST R ISO 15223-1-2023 | Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |

#### 11. Warranty

The warranty period is 3 years. This shall not extend to defects that are attributable to improper use or normal wear and tear. Transatlantic shall not be held liable for accidental or indirect injury or material damage.  
The warranty shelf life is 3 years.

---

## 12. Waste disposal

---

The device must be disinfected before being taken out of operation and disposed of. The Transatlantic Silverboard is not subject to any special disposal regulations. It must be disposed of in accordance with national and local legislation.

**For the Russian Federation**, it is recommended to dispose of the device in accordance with the rules implemented by SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" as class A (epidemically non-hazardous waste).

**Transatlantic Handels-Gesellschaft Stolpe & Co. mbH**

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61267 Neu-Anspach, Germany

Tel.: +49 (0)6081 9430-50

Fax: +49 (0)6081 9430-80

info@transat.de

[www.transatlantic.de](http://www.transatlantic.de)

## Instruction for use

### Transaroll Patient Transfer Aids

In various designs

**Transaroll Blueboard Mini, foldable**

**Transaroll Blueboard Mini, non-foldable**

**Transaroll Blueboard Standard, foldable**

**Transaroll Blueboard Long, foldable**

Please read the usage instructions carefully before first use!

**Symbols:**

|                                                                                   |                                    |                                                                                   |                                                                                                                              |                                                                                    |                                                             |
|-----------------------------------------------------------------------------------|------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------|
|  | Bed-to bed transfer                |  | Max. weight load                                                                                                             |  | Observe operating instructions                              |
|  | Warning                            |  | Manufacturer                                                                                                                 |  | Mark of conformity                                          |
|  | Temperature range                  |  | Keep away from sunlight                                                                                                      |  | Medical device Designation                                  |
|  | Catalogue number                   |  | Batch no./Lot no.                                                                                                            |  | The bridging of free transfer levels must not exceed 15 cm. |
|  | Keep dry                           |  | Slip hazard                                                                                                                  |  | Latex-free                                                  |
|  | UDI - Unique Device Identification |  | NFC ( Near Field Communication): devices with NFC support technology can exchange the data at close range without connecting |                                                                                    |                                                             |

**Manufacturer:**

**Transatlantic Handels-Gesellschaft Stolpe & Co. mbH**  
 Siemensstraße 7  
 61267 Neu-Anspach, Germany  
 Tel.: +49 (0)6081 9430-50  
 Fax: +49 (0)6081 9430-80  
 info@transat.de  
[www.transatlantic.de](http://www.transatlantic.de)

Abbreviated name: Transatlantic

**Authorized representative of the manufacturer in the Russian Federation, providing warranty/post-warranty service and support:**

**BioSystems LLC**

Russia, 197022, Saint Petersburg, 23 Professora Popova Street, Let. E  
 Tel.: (812) 320-49-49, Fax (812) 320-49-40,  
 e-mail: main@bioline.ru  
[www.bioline.ru](http://www.bioline.ru)  
 Hotline 8-800-444-33-49

## of the medical device and its designs:

### ansaroll Patient Transfer Aids

#### I. Various designs:

1. Transaroll Silverboard Mini, foldable – 1 pc.
2. Transaroll Silverboard Mini, non-foldable – 1 pc.
3. Transaroll Silverboard Standard, foldable – 1 pc.
4. Transaroll Silverboard Special, foldable – 1 pc.
5. Transaroll Silverboard Special, non-foldable – 1 pc.
6. Transaroll Blueboard Mini, foldable – 1 pc.
7. Transaroll Blueboard Mini, non-foldable – 1 pc.
8. Transaroll Blueboard Standard, foldable – 1 pc.
9. Transaroll Blueboard Long, foldable – 1 pc.
10. Transaroll Hygieneboard Mini, foldable, with accessories – 1 pc.:

#### Accessories:

##### Disposable sliding film, design:

- 60 cm \* 100 m, a roll – from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs

11. Transaroll Hygieneboard Mini, non-foldable, with accessories – 1 pc.:

#### Accessories:

##### Disposable sliding film, design:

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs

12. Transaroll Hygieneboard Standard, foldable, with accessories – 1 pc.:

#### Accessories:

##### Disposable sliding film, design:

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs
- 53 \* 180 cm, 60 pcs per roll, perforated - from 1 to 10 pcs

#### II. Instruction for Use

Hereafter referred to as the Transaroll Patient Transfer Aids, the board, the medical device.


## 1. Technical Specifications

The principle of the Patient Transfer Aids operation is based on the movement of two interacting surfaces relative to each other, one of which is in a static state, and the other rotates around it.

The Transaroll Blueboard consists of 2 components - basic body and cover. The blue cover fits onto the basic body and rotates freely around it.

The basic body may be foldable depending on the board design. The cover is removable and reusable.

### Transaroll Blueboard

Transaroll Blueboard Standard, foldable  
REF 4501100

#### 1.1 Main technical specifications

| REF     | Item designation / variants             | Dimensions*, L x W x T, mm | Weight *, kg | Max. rolling load, kg | Max. tension load**, kg | UDI-DI                         |
|---------|-----------------------------------------|----------------------------|--------------|-----------------------|-------------------------|--------------------------------|
| 4501100 | Transaroll Blueboard Standard, foldable | 1800 x 500 x 20            | 1.800        | 200                   | 300                     | (01)04050856043054(11)xx(10)xx |
| 4500400 | Transaroll Blueboard Mini, non-foldable | 950 x 500 x 20             | 0.990        | 200                   | 300                     | (01)04050856036384(11)xx(10)xx |
| 4500403 | Transaroll Blueboard Mini, foldable     | 950 x 500 x 20             | 0.990        | 200                   | 300                     | (01)04050856033048(11)xx(10)xx |
| 4500406 | Transaroll Blueboard Long, foldable     | 1950 x 500 x 20            | 1.960        | 200                   | 300                     | (01)04050856033049(11)xx(10)xx |

\* Dimensions and weight are indicated taking into account the acceptable tolerance of  $\pm 10\%$

\*\* Lifting for a short time no more than 30 seconds


The force applied by medical personnel to move the board horizontally with the patient placed on it:

| Design                                  | Maximum pushing force, no more than: |         | Maximum traction force, no more than: |         |
|-----------------------------------------|--------------------------------------|---------|---------------------------------------|---------|
| Load                                    | 150 kg                               | 200 kg  | 150 kg                                | 200 kg  |
| For all designs of Transaroll Blueboard | 123.5 N                              | 165.0 N | 148.0 N                               | 195.0 N |

## 1.2 Material

| Part of the device             | Material                 | Color | Contact type              | Property    |
|--------------------------------|--------------------------|-------|---------------------------|-------------|
| Transaroll basic body          | Polyethylene, latex-free | Black | Direct short-term contact | Radiolucent |
| Cover for Transaroll Blueboard | Nylon                    | Blue  |                           |             |

## 1.3 Packaging

The board is packaged with instructions for use in a clear plastic bag and then placed in a brown cardboard box and sealed with a Transatlantic adhesive tape.

Individual packaging is also a shipping packaging.

Box size varies depending on the board dimensions:

| REF     | Item designation / variants             | Box size*, mm  | Gross weight*, kg |
|---------|-----------------------------------------|----------------|-------------------|
| 4501100 | Transaroll Blueboard Standard, foldable | 920 x 535 x 75 | 2.8               |
| 4500400 | Transaroll Blueboard Mini, non-foldable | 975 x 520 x 57 | 2.0               |
| 4500403 | Transaroll Blueboard Mini, foldable     | 975 x 520 x 57 | 2.0               |
| 4500406 | Transaroll Blueboard Long, foldable     | 975 x 520 x 57 | 3.0               |

\* Dimensions and gross weight (in a shipping packaging) are indicated taking into account the acceptable tolerance of  $\pm 10\%$

## 1.4 Use, storage and transportation conditions

Use:

|                         |                          |
|-------------------------|--------------------------|
| Operating temperature*: | +10°C+30°C               |
| Humidity:               | Not determined, keep dry |

Storage and transportation:

|                             |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Storage temperature*:       | Board: +10°C+30°C                                                      |
| Storage conditions:         | Store in a dry place, protected from moisture; keep away from sunlight |
| Humidity:                   | Board: Not determined, keep dry                                        |
| Transportation conditions*: | +10°C+30°C                                                             |

\*Storage and transportation conditions are set by the manufacturer to preserve the structural behaviour of the basic body and the sliding properties of the film.

Transportation is allowed by any type of covered means of transport subject to the recommended temperature conditions.

---

## 2. Indication / areas of application / intended use

---

**Indication:** No specific clinical picture can be given as an indication because of the various possible applications. The Transaroll Blueboard is suitable for all applications on persons requiring assistance with lateral transfers.

**Areas of application:** The Blueboard can be used to bridge gaps when two lying surfaces are positioned next to each other (e.g. nursing bed - nursing bed, nursing bed - shower trolley, nursing bed - operating table, stretcher - examination table).

**Intended use:** The Transaroll board is suitable for gentle, pain-free and low-friction lateral transfer (by sliding) of persons whose movement is limited.

**Potential users:** No special education or knowledge is required to use this board. Please read the instructions carefully before use. Potential users are junior and middle grade medical staff after reading the instruction for use.

---

## 3. Redeployment

---

The Transaroll Blueboards are intended for reuse. They are not a disposable medical device. It is essential to conduct a thorough functional check and to clean or disinfect the basic structure before using it again.

---

## 4. Contraindications / risks of use

---

There are no known contraindications.

Before using the board, users must carefully read the Instruction for use, which describes in detail the procedure and features of handling the device, including visual drawings. If users do not read the instructions for handling the Transaroll Blueboard, the board may be handled incorrectly, which could result in damage to the health of the patient and/or user. Observe the application instructions under point 5 and the safety instructions under point 7.

The application risk analysis shows that the provided preventive measures do not require additional risk management measures, provided that the handling rules and safety regulations described in this instruction are observed.

It is the sole responsibility of the user to determine the suitability of the Transatlantic Blueboard for the specific application. Transatlantic shall not be held accountable for any improper use of Transaroll products.

## 5. Use / placement of the board

Position the board to the side, next to the patient.

Position yourself next to the patient, find a support point, and, resting on the patient's shoulder and hip, push him in the transfer direction.

Ensure that the brake is on!

To remove the board, reposition the patient in a short-term gentle lateral position.

Tilt the patient to the side gently and place the edge of the board under the patient so that the patient's head, shoulder girdle, and pelvic region are on the edge of the board.

Keeping the patient in lateral position, slide the board from under the patient.

The Transaroll board is radiolucent, has no metal parts and can be used during X-ray, MRI, CT procedures without the need to reposition the patient.

## 6. Service life / life cycle

The expected service life of the Transaroll Blueboard is 5 years if it is used in accordance with the safety instructions. There is no limitation to this period.

A longer or shorter service life is possible depending on the load. Usability may also end prematurely in the event of relevant damage to the patient transfer aids.

## 7. Safety instructions / functional check

**WARNING** signifies a possible hazard that could lead to personal injury or material damage.

- Defective products must no longer be used.
- Slip hazard! Ensure that the patient does not slide off the bed, recliner, etc. during transfer or repositioning.
- Risk of accident! Always make sure when using the Transaroll Blueboard that the bed, recliner or any movable base is stable, i.e. that any brakes on it are locked.
- Ensure that the Transaroll Blueboard is stable on both sides at the transfer levels. The bridging of free transfer levels must not exceed 15 cm.
- It is the responsibility of the user/institution to evaluate if multiple users may be needed for the transfer or if additional safety measures must be taken.
- Never leave the patient unsupervised on the Transaroll Blueboard.

- The Transaroll Blueboard must not remain under the patient after transfer.
- Pay attention, particularly when lifting the Transaroll Blueboard for a short time, that the patient is lying securely on the board and cannot slide down.
- Make sure that the transfers are performed by the users in a way that is back-friendly.
- Weight limits must be observed (200 kg).

Repositioning products are also subject to a natural ageing process. Inspect the device regularly and carefully. Ensure that the device is in impeccable condition and that the material has not been damaged (no tears or cuts in the cover, joint opening, tears, cuts, or spalling of the board basic body).

**Precautions:** do not use the board if it is damaged.

Unsuitability criteria: there is damage to the cover (tears, cuts, joint opening) or to the basic body of the board (tears, cuts, spalling).

The Transaroll Patient Transfer Aids in various designs are subject to the highest quality requirements.

The Transaroll Blueboard Patient Transfer Aids do not have accessories, are not intended for use in combination with other medical or non-medical products and are used independently.

The board does not contain any medicines, materials of animal or human origin.

The board is ready for use immediately after unpacking and does not require installation, setting up or other steps for commissioning.

The board is not sterile and does not require sterilization before use.

The board is not a source of any radiation and therefore does not affect other medical devices. It is also not affected by external factors from other medical devices.

Any serious incidents must be reported to the manufacturer and the competent authority in the Member State where the user and/or patient is based.

---

## 8. Cleaning, care and disinfection

---

The cover must be cleaned with medical disinfectants prior to and after use each time a patient is changed.

Disinfection is carried out by wiping with a cloth moistened with a disinfectant or spraying a disinfectant solution onto the surface. Surfaces to be disinfected must be left to dry completely.

### 8.1. Cover (blue)

When cleaning, do not use abrasive or bleaching agents and objects with sharp edges. The use of bleaching agents or cleaning agents, including those that contain alcohol, can permanently impair the functional properties of the blue undercover. The blue cover can be easily disinfected using oxygen-based disinfectant solutions or solutions containing alkyldimethylbenzylammonium chloride, polyhexamethylene guanidine chloride.

For users in the Russian Federation, it is recommended to use for disinfection - Detergent disinfectant "DeMoS", produced by LiG LLC, Russia (certificate of state registration RU 77.99.01.002.E.035064.09.11) or similar products based on alkyldimethylbenzylammonium chloride and/or polyhexamethylene guanidine chloride, approved for use.

If necessary, the cover can be removed. The Transroll Blueboard cover can be washed at 60°C and dried at low temperatures (40-50°C). Use a mild detergent! Don't use a softener!

|                                                                                   |                             |                                                                                   |                                  |
|-----------------------------------------------------------------------------------|-----------------------------|-----------------------------------------------------------------------------------|----------------------------------|
|  | Washing at 60°              |  | Tumble drying at low temperature |
|  | Do not use bleaching agents |  | Do not iron                      |
|  | Do not dry clean            |  | Latex-free                       |

The cover is unsuitable for use in an autoclave due to the thermo-sensitive structure of the material.

## 8.2. Basic body

The board itself withstands a lot of wear and tear and can be cleaned using normal household cleaners.

Conventional cleaning: Use warm water and a mild cleaning agent for cleaning. If necessary, wipe the board using a damp cloth and then leave it to dry thoroughly.

Disinfection: The board can be effortlessly disinfected using an alcohol-based disinfectant solution (based on isopropyl, n-propyl, ethyl or other alcohols). Always leave the solution to dry thoroughly.

The board is unsuitable for use in an autoclave due to the thermo-sensitive structure of the material.

**For users in the Russian Federation**, it is recommended to use for disinfection:

- Disinfectant "Avandez-A", produced by Dezos LLC, Russia (certificate of state registration RU 77.99.32.002.E.002145.06.21) or similar products based on ethyl, isopropyl or n-propyl alcohols, approved for use.

## 9. Product marking /CE


The Transroll Blueboard is a Class I medical device and is compliant with EU Regulation 2017/745 of the European Parliament and of the Council dated 05.04.2017 regarding medical devices (MDR).

We are certified as a manufacturer and distributor of medical devices in accordance with DIN EN ISO 9001:2015 as well as DIN EN ISO 13485:2016. Our medical products have, of course, passed the relevant conformity assessment procedure and therefore bear the CE mark.

**Individual product marking** consists of two film labels of a special shape corresponding to the board form, which [labels] are placed on the top and bottom sides of the board and contain the following information:

Information on the top label: name and address of the manufacturer, catalog number, batch code, manufacturing date, maximum safe load, symbols of use, storage and transportation conditions, caution symbols, instruction symbol, medical device symbol, risk class, UDI symbol, DataMatrix code, CE marking, (details see in table "Symbols"), label revision (rev.), company logo, country of origin (Made in Germany).

Information on the bottom label: name of the design, manufacturing date, label revision (rev.), QR code, NFC symbol, company logo.



The board marking contains a symbol indicating that there is an NFC tag containing information on the board serial number. To read this information, use the pre-installed NFC application on your smartphone (if available). To do this, you need to open the NFC application and bring your smartphone no more than 10 cm away from the NFC symbol on the board; information on the serial number will be displayed on the smartphone screen.

Marking on the board has a DataMatrix code: when reading the code, the following saved information is displayed: (01) GTIN (Global Trade Item Number GS1), (11) manufacturing date and (10) batch code.

### Sample of board marking

QR code

Sample of packaging marking

The marking in the form of a Russian-language sticker is additional to the manufacturer's marking, which contains the catalog number, lot/batch number and warning symbols. Additional Russian-language marking is applied next to the manufacturer's marking, without overlapping it. The data thereon complement each other and must be taken into account as a whole.

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General view of marking using an example of Transaroll Blueboard Mini, non-foldable

#### 10. Standards and specifications applicable in whole or in part

|                             |                                                                                                                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------|
| DIN EN ISO 13485:2021-12    | Medical devices - Quality management systems - Requirements for regulatory purposes                                        |
| DIN EN ISO 14971:2020-07    | Medical devices - Application of risk management to medical devices                                                        |
| DIN EN ISO 20417:2022-03    | Medical devices - Information to be supplied by the manufacturer                                                           |
| DIN EN ISO 15233-1:2022-02  | Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements    |
| DIN EN ISO 10993-1:2021-05  | Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing within a risk management process |
| DIN EN ISO 10993-10:2014-10 | Medical devices - Biological evaluation of medical devices. Part 10: Tests for irritation and skin sensitization           |
| DIN EN ISO 10993-11:2018-09 | Medical devices. Biological evaluation of medical devices. Part 11. Tests for systemic toxicity                            |

#### For the Russian Federation:

##### Russian standards

|                         |                                                                                                                         |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------|
| GOST R 50444-2020       | Medical instruments, apparatus and equipment. General technical requirements                                            |
| GOST R ISO 15223-1-2023 | Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |

#### 11. Warranty

The warranty period is 3 years. This shall not extend to defects that are attributable to improper use or normal wear and tear. Transatlantic shall not be held liable for accidental or indirect injury or material damage. The warranty shelf life is 3 years.

#### 12. Waste disposal

The device must be disinfected before being taken out of operation and disposed of. The Transatlantic Blueboard is not subject to any special disposal regulations. It must be disposed of in accordance with national and local legislation.

**For the Russian Federation**, it is recommended to dispose of the device in accordance with the rules implemented by SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" as class A (epidemically non-hazardous waste).

**Transatlantic Handels-Gesellschaft Stolpe & Co. mbH**

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61267 Neu-Anspach, Germany

Tel.: +49 (0)6081 9430-50

Fax: +49 (0)6081 9430-80

info@transat.de

[www.transatlantic.de](http://www.transatlantic.de)

## Instruction for use

### Transaroll Patient Transfer Aids

In various designs

**Transaroll Hygieneboard Mini, foldable, with accessories**

**Transaroll Hygieneboard Mini, non-foldable, with accessories**

**Transaroll Hygieneboard Standard, foldable, with accessories**

Please read the usage instructions carefully before first use!

## Symbols:

|  |                                    |  |                                                                                                                              |  |                                                             |
|--|------------------------------------|--|------------------------------------------------------------------------------------------------------------------------------|--|-------------------------------------------------------------|
|  | Bed-to bed transfer                |  | Max. weight load                                                                                                             |  | Observe operating instructions                              |
|  | Warning                            |  | Manufacturer                                                                                                                 |  | Mark of conformity                                          |
|  | Temperature range                  |  | Keep away from sunlight                                                                                                      |  | Medical device Designation                                  |
|  | Catalogue number                   |  | Batch no./Lot no.                                                                                                            |  | The bridging of free transfer levels must not exceed 15 cm. |
|  | Keep dry                           |  | Slip hazard                                                                                                                  |  | Latex-free                                                  |
|  | UDI - Unique Device Identification |  | NFC ( Near Field Communication): devices with NFC support technology can exchange the data at close range without connecting |  | No reusable (only for disposable sliding film)              |
|  |                                    |  |                                                                                                                              |  | Distributor (only for disposable sliding film)              |

## Manufacturer:

**Transatlantic Handels-Gesellschaft Stolpe & Co. mbH**

Siemensstraße 7  
61267 Neu-Anspach, Germany  
Tel.: +49 (0)6081 9430-50  
Fax: +49 (0)6081 9430-80  
info@transat.de  
[www.transatlantic.de](http://www.transatlantic.de)

Abbreviated name: Transatlantic

**Authorized representative of the manufacturer in the Russian Federation, providing warranty/post-warranty service and support:**

### BioSystems LLC

Russia, 197022, Saint Petersburg, 23 Professora Popova Street, Let. E  
Tel.: (812) 320-49-49, Fax (812) 320-49-40,  
e-mail: main@bioline.ru  
[www.bioline.ru](http://www.bioline.ru)  
Hotline 8-800-444-33-49

## **Name of the medical device and its designs:**

### **Transaroll Patient Transfer Aids**

#### **I. Various designs:**

1. Transaroll Silverboard Mini, foldable – 1 pc.
2. Transaroll Silverboard Mini, non-foldable – 1 pc.
3. Transaroll Silverboard Standard, foldable – 1 pc.
4. Transaroll Silverboard Special, foldable – 1 pc.
5. Transaroll Silverboard Special, non-foldable – 1 pc.
6. Transaroll Blueboard Mini, foldable – 1 pc.
7. Transaroll Blueboard Mini, non-foldable – 1 pc.
8. Transaroll Blueboard Standard, foldable – 1 pc.
9. Transaroll Blueboard Long, foldable – 1 pc.
10. Transaroll Hygieneboard Mini, foldable, with accessories – 1 pc.:

#### **Accessories:**

Disposable sliding film, design:

- 60 cm \* 100 m, a roll – from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs

11. Transaroll Hygieneboard Mini, non-foldable, with accessories – 1 pc.:

#### **Accessories:**

Disposable sliding film, design:

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs

12. Transaroll Hygieneboard Standard, foldable, with accessories – 1 pc.:

#### **Accessories:**

Disposable sliding film, design:

- 60 cm \* 100 m, a roll - from 1 to 10 pcs
- 60 cm \* 100 m, a roll in a dispenser box - from 1 to 10 pcs
- 53 \* 180 cm, 60 pcs per roll, perforated - from 1 to 10 pcs

#### **II. Instruction for Use**

Hereafter referred to as the Transaroll Patient Transfer Aids, the board, the medical device.

## 1. Technical Specifications


### 1.1 Main technical specifications

The principle of the Patient Transfer Aids operation is based on the movement of two interacting surfaces relative to each other, one of which is in a static state, and the other rotates around it.

The Transaroll Hygieneboard consists of 2 components - basic body and disposable sliding film which is an accessory to the board.

The basic body may be foldable depending on the board design. The green disposable sliding film fits onto the basic body and rotates freely around it. The sliding film is removable and disposable.

The user can select one or more designs of the disposable sliding film: in a roll or in a dispenser box with the ability to cut off any required length of the film, or in a roll with perforations for tearing off the film 180 cm long.

Standard, foldable

REF 4500026

The Transaroll Hygieneboard should be used together with a disposable sliding film on a roll (accessory). The disposable sliding film is a polyethylene sleeve that is tearproof, latex-free, nonirritant, biocompatible and food-safe.

For convenience of the user, the disposable sliding film is available in several designs: in a roll for independently cutting the film to the required length, or in a roll with perforations for tearing off the film of a specified size. Depending on the needs, the medical institution can select one or more of the possible suitable sliding film designs, in quantities from 1 to 10 rolls.

| No. | REF     | Item designation / variants                                         | Dimensions*, L x W x T | Weight* , kg | Max. rolling load, kg | Short-term max. weight load**, kg | UDI-DI                   |
|-----|---------|---------------------------------------------------------------------|------------------------|--------------|-----------------------|-----------------------------------|--------------------------|
| 1   | 4500036 | <b>Hygieneboard Mini</b> , foldable:<br>Accessories:                | 950 x 500 x 20 mm      | <b>0.700</b> | 200                   | 300                               | (01)04050856042156(11)xx |
|     | 4400200 | Disposable sliding film, - 60 cm * 100 m, a roll                    | 100 m x 600 mm         | <b>6.180</b> | -                     | -                                 | (01)4250355360601(10)x   |
|     | 4400203 | Disposable sliding film, - 60 cm * 100 m, a roll in a dispenser box | 100 m x 600 mm         | <b>6.360</b> | -                     | -                                 | (01)4250355360656(10)x   |
| 2   | 4500035 | <b>Hygieneboard Mini</b> non-foldable<br>Accessories:               | 950 x 500 x 20 mm      | <b>0.700</b> | 200                   | 300                               | 01)04050856042149(11)xx  |

|   |         |                                                                           |                          |       |     |     |                               |
|---|---------|---------------------------------------------------------------------------|--------------------------|-------|-----|-----|-------------------------------|
|   | 4400200 | Disposable sliding film,<br>- 60 cm * 100 m, a roll - max.<br>10 pcs      | 100 m x 600 mm           | 6.180 | -   | -   | (01)42503553606<br>01(10)x    |
|   | 4400203 | Disposable sliding film,<br>- 60 cm * 100 m, a roll in a<br>dispenser box | 100 m x 600 mm           | 6.360 | -   | -   | (01)42503553606<br>56(10)x    |
| 3 | 4500026 | <b>Hygieneboard Standard,<br/>foldable</b><br>Accessories:                | 1800 x 500 x 20<br>mm    | 1.750 | 200 | 300 | (01)04050856035<br>493 (11)xx |
|   | 4400200 | Disposable sliding film,<br>- 60 cm * 100 m, a roll - max.<br>10 pcs      | 100 m x 600 mm           | 6.180 | -   | -   | (01)42503553606<br>01(10)x    |
|   | 4400203 | Disposable sliding film,<br>- 60 cm * 100 m, a roll in a<br>dispenser box | 100 m x 600 mm           | 6.360 | -   | -   | (01)42503553606<br>56(10)x    |
|   | 4400204 | Disposable sliding film,<br>- 53 * 180 cm, 60 pcs per<br>roll, perforated | 1800 x 530 mm,<br>60 pcs | 6.300 | -   | -   | (01)42503553947<br>29 (10)x   |

\* Dimensions and weight are indicated taking into account the acceptable tolerance of  $\pm 10\%$

\*\* Lifting for a short time no more than 30 seconds

The force applied by medical personnel to move the board horizontally with the patient placed on it:

| Design                                        | Maximum pushing force, no<br>more than: |         | Maximum traction force, no<br>more than: |         |
|-----------------------------------------------|-----------------------------------------|---------|------------------------------------------|---------|
| Load                                          | 150 kg                                  | 200 kg  | 150 kg                                   | 200 kg  |
| For all designs of<br>Transaroll Hygieneboard | 150.0 N                                 | 200.0 N | 160.0 N                                  | 205.0 N |

### 1.3. Material

| Part of the device                          | Material                 | Color | Contact<br>type               | Property    |
|---------------------------------------------|--------------------------|-------|-------------------------------|-------------|
| Transaroll basic body                       | Polyethylene, latex-free | Black | Direct short-<br>term contact | Radiolucent |
| Sliding film for the board (all<br>designs) | LDPE                     | Green |                               |             |

### 1.3 Packaging

The board is packaged with instructions for use in a clear plastic bag and then placed in a brown cardboard box and sealed with a Transatlantic adhesive tape.

The disposable sliding film for the Transaroll Hygienboard on a roll and is placed in a cardboard box.

Individual packaging is also a shipping packaging.

Box size varies depending on the board dimensions:

| REF     | Variants and accessories                                               | Box size*, mm   | Gross weight*, kg |
|---------|------------------------------------------------------------------------|-----------------|-------------------|
| 4500036 | Transaroll Hygieneboard Mini, foldable                                 | 975 x 520 x 57  | 1.7               |
| 4530035 | Transaroll Hygieneboard Mini, non-foldable                             | 975 x 520 x 57  | 1.7               |
| 4500026 | Transaroll Hygieneboard Standard foldable                              | 920 x 535 x 75  | 2.7               |
| 4400200 | Disposable sliding film,<br>- 60 cm * 100 m, a roll                    | 660 x 155 x 155 | 6.6               |
| 4400203 | Disposable sliding film,<br>- 60 cm * 100 m, a roll in a dispenser box | 660 x 155 x 155 | 6.8               |
| 4400204 | Disposable sliding film,<br>- 53 * 180 cm, 60 pcs per roll, perforated | 660 x 155 x 155 | 6.8               |

\* Dimensions and gross weight (in a shipping packaging) are indicated taking into account the acceptable tolerance of  $\pm 10\%$

#### 1.4 Use, storage and transportation conditions

Use:

|                         |                          |
|-------------------------|--------------------------|
| Operating temperature*: | +10°C+30°C               |
| Humidity:               | Not determined, keep dry |

Storage and transportation:

|                             |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Storage temperature*:       | Board: +10°C+30°C<br>Disposable sliding film: +10°C+25°C               |
| Storage conditions:         | Store in a dry place, protected from moisture; keep away from sunlight |
| Humidity:                   | Board: Not determined, keep dry<br>Disposable sliding film: 50-75%     |
| Transportation conditions*: | +10°C+30°C                                                             |

\*Storage and transportation conditions are set by the manufacturer to preserve the structural behaviour of the basic body and the sliding properties of the film.

Transportation is allowed by any type of covered means of transport subject to the recommended temperature conditions.

**Disposable sliding film:** Protecting the sliding film from dirt, heat, moisture and sunlight is essential, otherwise the sliding properties may be significantly impaired, so it is essential to comply with storage conditions specified in this section.

Opened rolls of disposable sliding film should be used within 6 months.

## 2. Indication / areas of application / intended use

**Indication:** No specific clinical picture can be given as an indication because of the various possible applications. The Transaroll Hygieneboard is suitable for all applications on persons requiring assistance with lateral transfers.

**Areas of application:** The Transaroll Hygieneboard can be used to bridge gaps when two lying surfaces are positioned next to each other (e.g. nursing bed - nursing bed, nursing bed - shower trolley, nursing bed -

operating table, stretcher - examination table). Through the use of disposable sliding film, the Hygieneboard is recommended to use in cases where hygiene is critical (surgical department, X-ray department, intensive care unit, general medicine ward). By virtue of its exceptional sliding properties, the sliding film also simplifies complex transfers.

**Intended use:** The Transaroll board is suitable for gentle, pain-free and low-friction lateral transfer (by sliding) of persons whose movement is limited.

**Potential users:** No special education or knowledge is required to use this board. Please read the instructions carefully before use. Potential users are junior and middle grade medical staff after reading the instruction for use.

### 3. Redeployment

The Transaroll Hygieneboard basic body is intended for reuse.

It is essential to conduct a thorough functional check and to clean or disinfect the basic structure before using it again.

The disposable sliding film can be used several times to transfer the same patient. Then it must be removed from the basic body and disposed of.

### 4. Contraindications / risks of use

There are no known contraindications.

Before using the board, users must carefully read the Instruction for use, which describes in detail the procedure and features of handling the device, including visual drawings. If users do not read the instructions for handling the Transaroll Hygieneboard, the board may be handled incorrectly, which could result in damage to the health of the patient and/or user. Observe the application instructions under point 5 and the safety instructions under point 7.

The application risk analysis shows that the provided preventive measures do not require additional risk management measures, provided that the handling rules and safety regulations described in this instruction are observed.

It is the sole responsibility of the user to determine the suitability of the Transaroll Hygieneboard for the specific application. Transatlantic shall not be held accountable for any improper use of Transaroll products.

### 5. Use / placement of the board

Before transferring the patient, a disposable sliding film must be placed on the Hygieneboard basic body.

Depending on design of the sliding film, perform the following actions:

|         |                                                                   |
|---------|-------------------------------------------------------------------|
| 4400200 | Disposable sliding film,<br>- 60 cm * 100 m, a roll - max. 10 pcs |
|---------|-------------------------------------------------------------------|

- > Open the cover of the box.
- > Turn the roll in the box so that the end of the sliding film is on top.
- > Remove the protective tape from the roll
- > Pull the end of the sliding film with hands out of the box and draw the required length out of the box.
- > Hold the sliding film with one hand.
- > With your other hand, put the cutter knife to the side of the box.
- > Cut the disposable sliding film off in one continuous motion.

|         |                                                                        |
|---------|------------------------------------------------------------------------|
| 4400203 | Disposable sliding film,<br>- 60 cm * 100 m, a roll in a dispenser box |
|---------|------------------------------------------------------------------------|

- > Open the cover of the dispenser box.
- > Remove the sliding film cutter and cord.

- Turn the roll in the box so that the end of the sliding film is on top.
- Remove the protective tape from the roll.
- Using both hands, pull the end of the sliding cover about 50cm out of the box.
- Pull the end of the sliding film out of the dispenser box.
- Close the cover of the box.

#### Cutting the disposable sliding film to length

- Pull the required length of the sliding film with hands out of the box.
- Hold the sliding film with one hand.
- With your other hand, put the cutter knife to the cover of the box.
- Cut the disposable sliding film off in one continuous motion.

Recommendation: To get a precise straight cut, you can also run the cutter along the edge of the dispenser box cover.

The text on the surface of the disposable film is printed at intervals of approximately 50 cm, which helps visually determine the required length.

The cord can be used to secure the cutter to the dispenser box or other place using an adhesive tape.

4400204

Disposable sliding film,  
- 53 \* 180 cm, 60 pcs per roll, perforated

- Open the cover of the box.
- Turn the roll in the box so that the end of the sliding film is on top.
- Remove the protective tape from the roll
- Pull the end of the sliding film with hands out of the box and draw the required length (before perforation holes) out of the box.
- Hold the sliding film with one hand.
- With your other hand, tear off the film along the perforation holes by pulling it from one edge.

Put the disposable sliding film over the basic body of the Hygieneboard. Now it can be used.

Position the board to the side, next to the patient.

Position yourself next to the patient, find a support point, and, resting on the patient's shoulder and hip, push him in the transfer direction.

Ensure that the brake is on!

To remove the board, reposition the patient in a short-term gentle lateral position.

Tilt the patient to the side gently and place the edge of the board under the patient so that the patient's head, shoulder girdle, and pelvic region are on the edge of the board.

Keeping the patient in lateral position, slide the board from under the patient.

The Transaroll Hygieneboard is radiolucent, has no metal parts and can be used during X-ray, MRI, CT procedures without the need to reposition the patient.

---

## 6. Service life / life cycle

---

The expected service life of the Transaroll Hygieneboard is 5 years if it is used in accordance with the safety instructions. There is no limitation to this period.

A longer or shorter service life is possible depending on the load. Usability may also end prematurely in the event of relevant damage to the patient transfer aids.

Accessories (disposable sliding film of all designs) have the following shelf life: shelf life of the packaged sliding film is 4 years from the manufacturing date; upon opening the packaging, the disposable sliding film should be used within 6 months to ensure the required sliding properties.

---

## 7. Safety instructions / functional check

---

**WARNING signifies a possible hazard that could lead to personal injury or material damage.**

- Defective products must no longer be used.
- Slip hazard! Ensure that the patient does not slide off the bed, recliner, etc. during transfer or repositioning.
- Risk of accident! Always make sure when using the Transaroll Hygieneboard that the bed, recliner or any movable base is stable, i.e. that any brakes on it are locked.
- Ensure that the Transaroll Hygieneboard is stable on both sides at the transfer levels. The bridging of free transfer levels must not exceed 15 cm.
- It is the responsibility of the user/institution to evaluate if multiple users may be needed for the transfer or if additional safety measures must be taken.
- Never leave the patient unsupervised on the Transaroll Hygieneboard.
- The Transaroll Hygieneboard must not remain under the patient after transfer.
- Pay attention, particularly when lifting the Transaroll Hygieneboard for a short time, that the patient is lying securely on the board and cannot slide down.
- Make sure that the transfers are performed by the users in a way that is back-friendly.
- Weight limits must be observed (200 kg).
- Risk of injury due to cutter/cutter knife! Follow the safety instructions when using a dispenser box with a disposable sliding film.
- The disposable sliding film can be used several times to transfer the same patient. Then it must be removed from the board basic body and disposed of.
- The disposable sliding film cannot be reused for another patient.
- Remove the disposable sliding film after each patient is transferred.
- Opened rolls of disposable sliding film should be used within 6 months. It is necessary to protect opened rolls from contamination, heat, moisture and sunlight, otherwise the sliding properties of the film may deteriorate significantly.
- A damaged disposable sliding film or cover that has lost its sliding properties due to improper use or storage should not be used.

Repositioning products are also subject to a natural ageing process. Inspect the device regularly and carefully. Ensure that the device is in impeccable condition and that the material has not been damaged (no tears or cuts in the sliding film, or tears, cuts, spalling of the board basic body).

The disposable sliding film is also subject to a natural ageing process. Inspect it regularly and carefully. Ensure that the disposable film is in good condition and that the material has not been damaged (no tears or cuts). Check that the smooth surface of the disposable sliding film allows it to slide and move easily over the surface. Check these properties of the film - in comparison with a new film at its best.

**Defective products must no longer be used!**

**Precautions:** do not use the board if it is damaged.

Unsuitability criteria: damage to the basic body of the board (tears, cuts, spalling), damage to the disposable sliding film (tears, cuts) or sliding property loss.

The Transaroll Patient Transfer Aids are subject to the highest quality requirements.

The Transaroll Hygieneboard Patient Transfer Aids have accessories for its combined application, are not intended for use in combination with other medical or non-medical products and are used independently.

The board does not contain any medicines, materials of animal or human origin.

The board is ready for use immediately after unpacking and does not require installation, setting up or other steps for commissioning.

The board is not sterile and does not require sterilization before use.

The board is not a source of any radiation and therefore does not affect other medical devices. It is also not affected by external factors from other medical devices.

Any serious incidents must be reported to the manufacturer and the competent authority in the Member State where the user and/or patient is based.

---

## 8. Cleaning, care and disinfection

---

### 8.1. Disposable sliding film (green) (accessory to the Transaroll Hygieneboard)

The disposable sliding film can be used several times to transfer the same patient. Then it must be removed from the board basic body and disposed of. Cleaning and disinfection are not required.

### 8.2. Basic body

The board itself withstands a lot of wear and tear and can be cleaned using normal household cleaners.

Conventional cleaning: Use warm water and a mild cleaning agent for cleaning. If necessary, wipe the board using a damp cloth and then leave it to dry thoroughly.

Disinfection: The board can be effortlessly disinfected using an alcohol-based disinfectant solution (based on isopropyl, n-propyl, ethyl or other alcohols). Always leave the solution to dry thoroughly.

The board is unsuitable for use in an autoclave due to the thermo-sensitive structure of the material.

For users in the Russian Federation, it is recommended to use for disinfection: - Disinfectant "Avandez-A", produced by Dezos LLC, Russia (certificate of state registration RU 77.99.32.002.E.002145.06.21) or similar products based on ethyl, isopropyl or n-propyl alcohols, approved for use.

---

## 9. Product marking /CE

---

The Transaroll Hygieneboard is a Class I medical device and is compliant with EU Regulation 2017/745 of the European Parliament and of the Council dated 05.04.2017 regarding medical devices (MDR).

We are certified as a manufacturer and distributor of medical devices in accordance with DIN EN ISO 9001:2015 as well as DIN EN ISO 13485:2016. Our medical products have, of course, passed the relevant conformity assessment procedure and therefore bear the CE mark.

**Individual product marking** consists of two film labels of a special shape corresponding to the board form, which [labels] are placed on the top and bottom sides of the board and contain the following information:

Information on the top label: name and address of the manufacturer, catalog number, batch code, manufacturing date, maximum safe load, symbols of use, storage and transportation conditions, caution symbols, instruction symbol, medical device symbol, risk class, UDI symbol, DataMatrix code, CE marking, (details see in table "Symbols"), label revision (rev.), company logo, country of origin (Made in Germany).

Information on the bottom label: name of the design, manufacturing date, label revision (rev.), QR code, NFC symbol, company logo.


The board marking contains a symbol indicating that there is an NFC tag containing information on the board serial number. To read this information, use the pre-installed NFC application on your smartphone (if

available). To do this, you need to open the NFC application and bring your smartphone no more than 10 cm away from the NFC symbol on the board; information on the serial number will be displayed on the smartphone screen.

Marking on the board has a DataMatrix code: when reading the code, the following saved information is displayed: (01) GTIN (Global Trade Item Number GS1), (11) manufacturing date and (10) batch code.

Marking on the board has a QR code: when reading, it contains a link to the manufacturer's website.


#### Sample of board marking

#### Marking on the board packaging

The label contains the following information: name and address of the manufacturer, catalog number, batch code, name of the design, size, foldable/non-foldable, manufacturing date, symbols of storage and transportation conditions, instruction symbol, medical device symbol, risk class, UDI symbol, DataMatrix code, CE marking, (see details in table "Symbols"), label revision (rev.), company logo, country of origin (Made in Germany), web-page.

### Sample of packaging marking

### Draft of additional marking on the packaging for the Russian Federation

In order to comply with Russian legislation, the marking in the form of a Russian-language stickers will be applied in addition to the manufacturer's marking.

The marking in the form of a Russian-language sticker is additional to the manufacturer's marking, which contains the catalog number, lot/batch number and warning symbols. Additional Russian-language marking is applied next to the manufacturer's marking, without overlapping it. The data thereon complement each other and must be taken into account as a whole.

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General view of marking using an example of **Transaroll Hygieneboard Standard, foldable**

### Marking for packaging of the disposable sliding film

The label contains the following information: name and address of the manufacturer, name and address of the distributor, catalog number, name of the design, size, quantity, batch code, manufacturing date, symbols of use, storage and transportation conditions, non-reusable symbol, instruction symbol, medical device symbol, risk class, UDI symbol, DataMatrix code, CE marking, (see details in table "Symbols"), label revision (rev.), company logo, country of origin (Made in Germany), web-page.

### Example of the marking for packaging of the disposable sliding film

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General view of marking using an example of the Disposable sliding film - 60 cm \* 100 m, a roll

### 10. Standards and specifications applicable in whole or in part

|                             |                                                                                                                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------|
| DIN EN ISO 13485:2021-12    | Medical devices - Quality management systems - Requirements for regulatory purposes                                        |
| DIN EN ISO 14971:2020-07    | Medical devices - Application of risk management to medical devices                                                        |
| DIN EN ISO 20417:2022-03    | Medical devices - Information to be supplied by the manufacturer                                                           |
| DIN EN ISO 15233-1:2022-02  | Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements    |
| DIN EN ISO 10993-1:2021-05  | Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing within a risk management process |
| DIN EN ISO 10993-10:2014-10 | Medical devices - Biological evaluation of medical devices. Part 10: Tests for irritation and skin sensitization           |
| DIN EN ISO 10993-11:2018-09 | Medical devices. Biological evaluation of medical devices. Part 11. Tests for systemic toxicity                            |

## For the Russian Federation:

### Russian standards

|                         |                                                                                                                         |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------|
| GOST R 50444-2020       | Medical instruments, apparatus and equipment. General technical requirements                                            |
| GOST R ISO 15223-1-2023 | Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements |

## 11. Warranty

The warranty period is 3 years. This shall not extend to defects that are attributable to improper use or normal wear and tear. Transatlantic shall not be held liable for accidental or indirect injury or material damage.

The manufacturer guarantees the properties of the Disposable sliding film throughout the entire shelf life (4 years).

The warranty shelf life is 3 years.

## 12. Waste disposal

The disposable sliding film should be removed from the board after its using and disposed of separately in accordance with state and local regulations.

The device must be disinfected before being taken out of operation and disposed of. There are no special recommendations of disposal for Transaroll Hygieneboard. It must be disposed of in accordance with national and local legislation.

**For the Russian Federation**, it is recommended to dispose of the device in accordance with the rules implemented by SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures":

Disposable sliding film – as class B (epidemically hazardous waste).

Basic body of the board (after preliminary disinfection) – as class A (epidemically non-hazardous waste).

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