
Synthes GmbH
Eimattstrasse 3
4436 Oberdorf
Switzerland

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Instruction for use
Surgical instruments for the installation, measurement and extraction of intramedullary implants, with accessories

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Official Certification:

The undersigned Notary certifies, that the preceding signature of


Robertson Philippe, 09.04.1977, from Belgium, in Neuchâtel

is authentic and that the signature currently act in the position of Senior Manager of Regulatory Affairs, Synthes GmbH, a company currently existing and organized under the laws of Switzerland, that the address is Eimattstrasse 3, 4436 Oberdorf.

The Notary Public is aware and attests to the above facts upon examination of the corresponding documents.

Given and signed in 4500 Solothurn, Switzerland

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10.1. Protection sleeve, 03.010.075- 1 piece

10.2. Instruction for use - 1 piece

11. Protection sleeve for Medullary Reamers, 03.033.009, contain:

11.1. Protection sleeve for Medullary Reamers, 03.033.009- 1 piece

11.2. Instruction for use - 1 piece

12. Sleeve of drill, 03.010.076, contain:

12.1. Sleeve of drill, 03.010.076- 1 piece

12.2. Instruction for use - 1 piece

13. Sleeve of drill, multihole, 03.033.008, contain:

13.1. Sleeve of drill, multihole, 03.033.008- 1 piece

13.2. Instruction for use - 1 piece

14. Drill Bit cannulated, flexible, for Quick Coupling DHS/DCS, 03.033.004, contain:

14.1. Drill Bit cannulated, flexible, for Quick Coupling DHS/DCS, 03.033.004- 1 piece

14.2. Instruction for use - 1 piece

15. Trocar for protection tube, 03.033.010, contain:

15.1. Trocar for protection tube, 03.033.010- 1 piece

15.2. Instruction for use - 1 piece

16. Trocar for protection sleeve, 03.010.077, contain:

16.1. Trocar for protection sleeve, 03.010.077- 1 piece

16.2. Instruction for use - 1 piece

17. Awl cannulated, 03.010.041, contain:

17.1. Awl cannulated, 03.010.041- 1 piece

17.2. Instruction for use - 1 piece

18. Screw connecting, cannulated, 03.033.006, contain:

18.1. Screw connecting, cannulated, 03.033.006- 1 piece

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- 18.2. Instruction. ... use – 1 piece
19. Screw for extraction Tibial and Femoral Nails, 03.010.00, contain:
 - 19.1. Screw for extraction Tibial and Femoral Nails, 03.010.00– 1 piece
 - 19.2. Instruction for use – 1 piece
20. Direct Measuring Device for Guide Wires, 03.010.493, contain:
 - 20.1. Direct Measuring Device for Guide Wires, 03.010.493– 1 piece
 - 20.2. Instruction for use – 1 piece
21. Direct Measuring Device for Guide Wires, 03.010.085, contain:
 - 21.1. Direct Measuring Device for Guide Wires, 03.010.085– 1 piece
 - 21.2. Instruction for use – 1 piece
22. Reamer for Hip Screws, 03.010.078, contain:
 - 22.1. Reamer for Hip Screws, 03.010.078– 1 piece
 - 22.2. Instruction for use – 1 piece
23. Fixation Sleeve for reamer for Hip Screws, 03.010.079, contain:
 - 23.1. Fixation Sleeve for reamer for Hip Screws, 03.010.079– 1 piece
 - 23.2. Instruction for use – 1 piece
24. Screwdriver Stardrive T25, self-holding, 03.010.519, contain:
 - 24.1. Screwdriver Stardrive T25, self-holding, 03.010.519– 1 piece
 - 24.2. Instruction for use – 1 piece
25. Combination wrench, 321.160, contain:
 - 25.1. Combination wrench, 321.160– 1 piece
 - 25.2. Instruction for use – 1 piece

Accessories:

1. Tray for Opening Instruments FRN, 60.033.006 (if needed) – 1 piece.
2. Tray for Insertion Instruments FRN, 60.033.007 (if needed) – 1 piece.
3. Tray for Locking Instruments FRN, 60.033.010 (if needed) – 1 piece.
4. Lid for Graphic Case, 60.037.010 (if needed) – 1 piece.
5. Label Sheet for Opening Instruments and Insertion Instruments FRN, 60.033.011 (if needed) – 1 piece.
6. Label Sheet for Locking Instruments FRN, 60.033.014 (if needed) – 1 piece.
7. Spacer for Graphic Case Tray, 60.033.015 (if needed) – 1 piece.

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Handle with Quick Coupling , 03.010.500

Handles with Quick Coupling are hand-held manual surgical instruments designed to attach to the proximal end of a surgical instrument to allow the surgeon to perform manipulations with the instrument. The device has a straight handle grip.

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Handle radiolucent, FRN , 03.033.001

Radiolucent handle with radiographic indicators supports x-ray visualization and assists with guide wire placement. Markers on the handle help visualize the nail height and determination of end cap size.

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Aiming arm for Greater Trochanter, radiolucent, FRN, 03.033.003;

Aiming arm for Piriformis Fossa, radiolucent, FRN, 03.033.002
Aiming instruments are intended for guiding of instruments for alignment with corresponding implants during orthopaedic surgery.

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Protection sleeve, 03.033.007;
 Protection sleeve, 03.025.040;
 Protection sleeve, 03.010.075;
 Protection sleeve, for Medullary Reamers, 03.033.009;
 Sleeve of drill, 03.010.076;
 Sleeve of drill, multihole, 03.033.008;
 Fixation Sleeve for reamer for Hip Screws, 03.010.079
 Sleeve and fixation sleeve - collection of hand-held surgical instruments that are able to protect the instrument and/or its operator and/or the patient from injury, guide a rotating surgical drill, allow the surgeon to perform manipulations with the instrument.

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Protection sleeve 03.033.007

Protection sleeve, 03.010.075

Protection sleeve, for Medullary Reamers, 03.033.009

Protection sleeve, 03.025.040


Sleeve of drill, multihole , 03.033.008

Sleeve of drill, 03.010.076

Fixation Sleeve for reamer for Hip Screws , 03.010.079

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Drill Bit cannulated, flexible for Quick Coupling DHS/DCS, 03.033.004

Drill bits are cylindrical constructs with a cutting tip and spiral flutes along their length. Pressure applied to the drill as it is rotating allows it to cut through the bone and advance, creating a cylindrical hole. The spiral flutes along the length of the drill carry bone debris away as the drill advances.

**Trocar for protection tube, 03.033.010;
Trocar for protection sleeve, 03.010.077**

Trocars are rigid hand-held surgical instruments with a sharp conical point used to manually puncture a bodily entry point to assist in the location and positioning of skeletal fractures and/or, upon entry, gently part the surrounding soft tissue.

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**Awl cannulated .010.041**

Awl is a metal surgical instrument and manually used for opening of the medullary canal before the insertion of the dedicated implant.

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**Screw connecting, cannulated, 03.033.006**

Connecting screw is designed for facilitated attachment of nail to insertion handle. Instrument is intended for positioning an implant into the body.

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**Screw for Extraction Tibial and Femoral Nails, 03.010.000**

Surgical instrument that is intended to be used to grasp and manipulate an implantable orthopedic device during extraction during a surgical intervention.

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
LIST AND DESCRIPTION OF MD MATERIALS

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	Number of Units
MATERIAL Stainless Steel	Material
 1	Contains hazardous substance
	Contains EU regulated chemical substance

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
DEVICE SERVICE LIFE

These products are provided in non-sterile packaging systems do not create a sterile barrier which could degrade over time. End of life of a device is normally determined by wear and damage due to use. Evidence of damage and wear on a device may include but is not limited to corrosion (i.e. rust, pitting), discoloration, excessive scratches, flaking, wear and cracks. Improperly functioning devices, devices with unrecognizable markings, missing or removed (buffed off) part numbers, damaged and excessively worn devices should not be used.


Examples of evidence of damage and wear on a device:

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

damaged device


right device

- Corrosion, Rusting, Pitting.
- Surfaces discolored, flaky or scratched.

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




Devices With Unrecognizable Markings, Etches/Buffered off or Missing Part Numbers.


Cracks at welded joints, cracks at assembly joints (e.g., pins), cracks in cannulated devices (e.g., threads, screwdriver tips).

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

Broken at driver/tip breakage, broken at assembly joints (e.g., plastic handles).

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
Bent shaft or body, tips/drivers deformed or twisted.

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- Blunt, bent, broken, rounding of sharp edges.
- Chipped, gouged, deep scratches.
- Does not rotate concentrically; circular scratching due to intended rotation of the device.
- Knocks and can thus break; damage due to impact.

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


- Threads damaged/peeling.

- Parts cracked or broken off.
- Surface is brittle, soft or has a burnt appearance.
- Discoloration or delamination of handle.

- Will not hold/retain (A)
- Will not release (B)
- Instrument cannot be fitted

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Cautions

All devices must be thoroughly cleaned and inspected prior to sterilization. Long, narrow lumens, blind holes, moving and intricate parts require particular attention during cleaning and inspection. During cleaning, only use detergents that are labelled for use on medical devices and in accordance with the manufacturer's instructions (e.g. temperature, contact time, and rinse time). Cleaning agents with a used dilution pH of within 7–9.5 are recommended. Highly alkaline conditions ($\text{pH} > 11$) can damage components/ devices, such as aluminum materials. Do not use saline, environmental disinfection (including chlorine solutions) or surgical antiseptics (such as iodine- or chlorhexidine-containing products). Do not use a cleaning aid that can damage the surface of instruments such as steel wool, abrasive cleaners or wire brushes.

– Only place devices with items of similar metallic composition together in an ultrasonic cleaner.

– Soiled or used devices should not be loaded into a case for cleaning in a mechanical washer. Soiled devices must be processed separate from trays and cases. Cases are designed to be an organizational tool for the steam sterilization process, a storage tool for all medical devices and an organizational tool for surgery.

– The sterilization parameters are only valid for devices that are adequately cleaned.

– The parameters listed are only valid for properly installed, maintained, calibrated and compliant reprocessing equipment in accordance with standards such as the ISO 15883 and ISO 17665 series.

– Surgical patients identified as at-risk for Creutzfeldt-Jakob disease (CJD) and related infections should be treated with single-use instruments. Dispose of instruments used or suspected of use on a patient with CJD after surgery and/or follow current national recommendations.

– Consult national regulations and guidelines for additional information. Compliance is additionally required with internal hospital policies and procedures and recommendations of manufacturers of detergents, disinfectants, and any clinical processing equipment.

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Limits on reprocessing

- Repeated processing cycles that include ultrasonic, mechanical washing and sterilization have minimal effects on surgical instrumentation.
- End of life of a device is normally determined by wear and damage due to use. Evidence of damage and wear on a device may include but is not limited to corrosion (e.g. rust, pitting), discoloration, excessive scratches, flaking, wear and cracks. Improperly functioning devices, devices with unrecognizable markings, missing or removed (buffed off) part numbers, damaged and excessively worn devices should not be used.
- Further details regarding end of life indicators are available from your sales representative or for download from the website: <http://www.depuysynthes.com/ifu>.

Point of Use Care

Wipe blood and/or debris from device throughout surgical procedure to prevent it from drying onto the surface.

- Flush cannulated devices with sterile or purified water to prevent the drying of soil and/or debris to the inside.
- Soiled devices should be separated from non-contaminated devices to avoid contamination of personnel or surroundings.
- Devices should be covered with a towel dampened with sterile or purified water to prevent blood and/or debris from drying.

Containment and Transportation

- Soiled devices should be transported separate from non-contaminated devices to avoid contamination.

Preparation for Cleaning (for all cleaning methods)

- It is recommended that devices should be reprocessed as soon as reasonably practical following use.
- Disassemble device, if device is able to be disassembled, prior to reprocessing. Further detailed instrument dismantling instructions are available from your local sales representative or for download from the website: <http://www.depuysynthes.com/ifu>.
- Open devices with ratchets, box locks or hinges.
- Care should be taken in the handling and cleaning of sharp devices. These are recommended to be cleaned separately to reduce risks of injury.

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– Lumens/cannula if devices should be manually processed prior to cleaning. Lumens/cannula should first be cleared of debris. Lumens/cannula should be brushed thoroughly using appropriately sized soft-bristled brushes and twisting action. Brushes should be tight-fitting. Brush size should be approximately the same diameter of the lumen/cannulation to be cleaned. Using a brush that is too big or too small for the diameter of the lumen/cannulation may not effectively clean the surface of a lumen/cannulation.

– Soak and/or rinse heavily soiled devices or cannulated devices prior to cleaning to loosen any dried soil or debris. Use a neutral or mild alkaline detergent solution. Follow the detergent manufacturer's instructions for use for correct exposure time, temperature, water quality and concentration. Use cold tap water to rinse devices.

– Devices must be cleaned separately from instrument trays and cases. Lids should be removed from cases for the cleaning process, if applicable.

Cleaning-Manual Method

1. Rinse soiled device under running cold tap water for a minimum of two minutes. Use a soft-bristled brush to assist in the removal of gross soil and debris.

2. Soak device in a neutral or mild alkaline detergent solution for a minimum of ten minutes. Follow the detergent manufacturer's instructions for use for correct exposure time, temperature, water quality and concentration.

3. Rinse device with cold water for a minimum of two minutes. Use a syringe, pipette, or water jet to flush lumens, channels and other hard to reach areas.

4. Manually clean device for a minimum of five minutes in a freshly prepared neutral or mild alkaline detergent solution. Use a soft-bristled brush to remove soil and debris. Actuate joints, handles and other moveable device features to expose all areas to the detergent solution, if applicable. Clean device under water to prevent aerosolization of contaminants. Note: fresh solution is a newly-made, clean solution.

5. Rinse device thoroughly using cold or warm tap water for a minimum of two minutes. Use a syringe, pipette or water jet to flush lumens and

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channels. Actuate joints, handles and other moveable device features in order to rinse thoroughly under running water, if applicable.

6. Visually inspect device. Repeat steps 2–6 until no visible soil remains on device. Ultrasonic process: Pre-cleaning steps 1–6 should occur prior to this step

7. Prepare a fresh detergent solution for the ultrasonic bath using a neutral or mild alkaline detergent solution. Follow the detergent manufacturer's instructions for use for correct exposure time, temperature, water quality and concentration. Note: fresh solution is a newly-made, clean solution.

8. Clean device ultrasonically for a minimum of 15 minutes, using a minimum frequency of 38 kHz.

9. Rinse device thoroughly with deionized (DI) or purified (PURW) water for a minimum of two minutes. Use a syringe, pipette or water jet to flush lumens and channels. Actuate joints, handles and other moveable device features in order to rinse thoroughly under running water, if applicable.

10. Visually inspect device. Repeat steps 2–10 until no visible soil remains on device.

11. Perform a final rinse on device using DI or PURW water for a minimum of 15 seconds.

12. Dry devices using a clean, soft, lint-free single-use cloth or medical grade compressed air. Ensure that all lumens and articulated areas are dried using compressed air.

Cleaning-Disinfection Automated Method

1. Rinse soiled device under running cold tap water for a minimum of one minute. Remove gross soil using a soft-bristled brush or soft, lint-free cloth.

2. Manually clean device for a minimum of two minutes in a freshly prepared neutral or mild alkaline detergent solution. Follow the detergent manufacturer's instructions for the correct dilution, temperature, water quality and exposure time. Use a soft-bristled brush to remove soil and debris. Actuate joints, handles and other moveable device features to expose all areas to detergent solution, if applicable. Clean device under water to prevent aerosolization of contaminants, Note: fresh solution is a newly-made, clean solution.

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3. Rinse device in cold to lukewarm running tap water for a minimum of one minute. Use a syringe, pipette or water jet to flush lumens and channels. Actuate joints, handles and other moveable device features in order to rinse thoroughly under running water, if applicable.

4. Prepare a fresh detergent solution for the ultrasonic bath using a neutral or mild alkaline detergent. Follow the detergent manufacturer's instructions for the correct dilution, temperature, water quality and exposure time. Note: fresh solution is a newly-made, clean solution.

5. Clean devices ultrasonically for a minimum of 15 minutes, using a minimum frequency of 38 kHz.

6. Rinse device using DI or PURW water for a minimum of two minutes. Use a syringe, pipette or water jet to flush lumens and channels. DI or PURW water must be used for final rinse.

7. Visually inspect device. Repeat steps 2–7 until no visible soil remains on device.

8. Automated washing shall be conducted in a validated washer-disinfector in compliance to ISO 15883-1 and -2, or to an equivalent standard. Load the device components in the washer-disinfector in accordance with manufacturer's instructions, ensuring that the devices and lumens can drain freely. Automated washing can be included as part of a validated washing, disinfection, and/or drying cycle in accordance to manufacturer's instructions. An example of a validated cycle used for cleaning validation included:

Cycle		Minimum Time (minutes)	Minimum	Type of Detergent
Type of Detergent	2		Cold tap water	N/A
Wash I	2		Cold tap water (<40°C)	Cleaning agent*
Wash II	5		Warm tap water (>40°C)	Cleaning agent*
Rinse	2		Warm DI or PURW (>40°C)	N/A
Rinse	5		>93°C	N/A
Dry	40		>90°C	N/A

* see section Additional Information

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Thermal disinfection

For automated cleaning-disinfection, thermally disinfect at a minimum of 93 °C for a minimum of 5 minutes (see Cleaning-Disinfection Automated Method). For devices with cannulations or lumens, orient the parts such that the lumen or cannulation are in a vertical position. If this is not possible due to space limitations within the automated/mechanical washer, use an irrigating rack/load carrier with connections designed to ensure an adequate flow of process fluids to the lumen or cannulation of the device if necessary.

Drying

If a dry cycle is not included in the mechanical washer:

- Dry each device thoroughly inside and out to prevent rusting and malfunction.
- Use a clean, soft, lint-free single-use cloth to avoid damage to the surface.

Pay special attention to threads, ratchets and hinges or areas where fluid can accumulate. Open and close devices so that all areas are reached. Dry hollow parts (lumens, cannulations) using the air jet with medical grade compressed air.

Inspection

Instruments should be inspected after processing, prior to sterilization, for end of life indicators such as:

- Cleanliness, if any residual soil is discovered during inspection, repeat the cleaning steps on those devices until all visible soil is removed from the device.
- Lack of moisture, carefully inspect device lumens and moving parts. If moisture is detected, manual drying should be performed.
- Damage, including but not limited to corrosion (e.g. rust, pitting), discoloration, excessive scratches, flaking, cracks and wear.
- Proper function, including but not limited to sharpness of cutting tools, bending of flexible devices, movement of hinges/joints/box locks and moveable features such as handles, ratcheting and couplings. Damaged or worn devices should not be used.

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Further detailed : tion control instructions and end of life indicators are available from your local sales representative or for download from the website: <http://www.depuysynthes.com/ifu>.

Lubricate instruments with moving parts, such as hinges and joints, spring-loaded ball bearings, and threaded parts. It is recommended to lubricate and maintain instruments with Special Oil only.

Disassembled devices should be reassembled prior to sterilization unless otherwise noted or the case is not configured for the assembled device. Further detailed instrument dismantling instructions are available from your local sales representative or for download from the website at <http://www.depuysynthes.com/ifu>.

Packaging

Put cleaned, dry devices into the proper location in the case. Additionally, use an appropriate sterilization wrap or re-usable rigid container system for sterilization, such as a sterile barrier system according to ISO 11607. Care should be taken to protect implants and pointed and sharp instruments from contact with other objects that may damage the surface.

Sterilization

Steam (moist heat) sterilization shall be performed in a locally approved, pre-vacuum (forced air removal) cycle. The steam sterilizer should be validated to the requirements of any local standards and guidance such as EN 285 or AAMI/ANSI ST8. The steam sterilizer should be installed and maintained in compliance to manufacturer's instructions and local requirements. Ensure that a steam sterilizer cycle is chosen that is designed to remove air from porous or lumened device loads in accordance to manufacturer's instructions and does not exceed the maximal sterilizer load

Cycle Type	Minimum Sterilization Exposure Time (minutes)	Minimum Sterilization (minutes) Exposure Temperature	Minimum Dry Time*
Prevacuum Saturated	4	132°C	20 minutes
	3	134°C	20 minutes

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steamforced air removal (pre-vacuum, minimum three pulses)			
--	--	--	--

* When applying dry times to cases and their accessories, dry times outside the standard healthcare prevacuum parameters may be required. This is especially important for polymer-based (plastic) cases/trays used in conjunction with heavy duty nonwoven sterilization wraps. The current recommended dry times for cases can range from a standard 20 minutes to an extended time of 60 minutes. The dry time is most often influenced by the presence of polymer based (plastic) materials; therefore changes such as elimination of silicone mats and/or change in sterile barrier system (i.e. heavy grade to light grade wrap) can reduce the necessary dry time. The user should employ verifiable methods (e.g. visual inspections) to confirm adequate drying. Dry times generally range from 20 to 60 minutes due to differences in packaging materials (Sterile Barrier System, e.g. wraps or re-usable rigid container systems), steam quality, device materials, total mass, sterilizer performance, and varying cool-down time.

Storage

Packaged products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity. Refer to sterilization wrap or rigid container manufacturers IFU for limits on sterile product storage time and storage requirements for temperature and humidity.

Additional Information

Further information regarding the use of specific cleaning agents, ultrasonic washers, washer-disinfector, packaging materials or sterilizers during validation studies are available on request. The manufacturer used the following during validation of these reprocessing recommendations:

– Manual cleaning: Manual Pre-Cleaning with Enzol Enzymatic Detergent 8 mL/L at 16–17 °C and Ultrasonic Cleaning with Enzol Enzymatic Detergent 8 mL/L at 18-25 °C.

– Automated Cleaning: Manual Pre-Cleaning with Enzol Enzymatic Detergent 8 mL/L at 16–17 °C and Ultrasonic Cleaning with Enzol

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Enzymatic Deter 8 mL/L at 18–25 °C. Water-disinfectant cleaning with (Wash 1) Prolystica 2X Concentrate Enzymatic Cleaner 1 mL/L at 23–26 °C and (Wash 2) Prolystica 2X Neutral Detergent 1 mL/L at 44–46 °C.

– Lint-Free Cloth: Berkshire Durx 670. The cleaning and sterilization information is provided in accordance with ISO 17664. The recommendations provided above have been validated by the medical device manufacturer as being capable of preparing a non-sterile medical device. It remains the responsibility of the processor to ensure that the processing is performed using equipment, materials and personnel in the reprocessing facility, and achieves the desired result. This requires validation and routine onitoring of the process. Likewise, any deviation by the processor from the recommendations provided should be properly evaluated for effectiveness and potential adverse consequences.

The chemical quality of the water used during reprocessing can impact device safety. Facilities should use the recommended water quality requirements for device reprocessing in accordance with local guidance (such as AAMI TIR 34, Water for the reprocessing of medical devices) and these instructions for use.

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Material Compatibility of Instruments and Implants in Clinical Processing

Knowledge of the materials used and their properties is essential for ensuring that instruments are proficiently processed and maintained.

Stainless steels

Instruments are made predominantly from corrosion-resistant steels, recognizable by their shiny or dull metallic color. As a result of their high chromium and nickel content, corrosion-resistant steels form a protective chromium oxide layer, known as a passive layer, on the metal surface. This passive layer protects the instrument against corrosion and rust. Incorrect or careless handling (e.g. damage to the surface) and attacks of a chemical, electrochemical or physical nature, can adversely affect the corrosion resistance. Two types of stainless steels are used, differentiated based on their composition and properties:

– Martensitic steels, which are corrosion resistant and whose high hardness can be influenced and adjusted by heat treatment, possess high wear resistance and high cutting edge retention. These steels are used for cutting and sharp-pointed instruments, e.g. drill bits, reamer heads, awl, burrs or cutting edges of pliers.

– Austenitic steels, which cannot be hardened by heat treatment, possess high corrosion resistance, elasticity and toughness, and are generally non-magnetic. These steels are used for non-cutting instruments, e.g. drill guides, gauges and aiming devices.

– the manufacturer recommends disinfectants, cleaners or detergents with pH 7–11 for all stainless steels.

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Aluminum, titanium and its alloys

Since aluminum is a lightweight material it is used, for example, for the graphic cases, instrument handles and certain other instrument parts. An electrochemical surface treatment (anodizing, “Ematal” or hard anodizing) produces a resistant oxide layer on the aluminum, which can be dyed.

Titanium and titanium alloys are widely used as implant materials. On instruments titanium is used for only a few applications, mainly color coding of instruments. The surface of titanium alloys is also treated electrochemically (anodizing), producing a resistant oxide layer. Various color shades can be applied using this layer. Although anodized aluminum, titanium and its alloys have good corrosion resistance, contact with strong alkaline detergents or disinfectants and solutions containing iodine or certain metal salts can lead to chemical attack and dissolution of the surface depending on the specific composition of the detergent.

Consequently, the manufacturer recommends disinfectants, cleaners or detergents with pH 6–9.5. Products with a higher pH value, especially higher than pH 11, should only be used subject to the material compatibility requirements stated on the datasheet and other information from the manufacturer of the detergent.

Plastics

Various plastics are used for certain instrument parts, e.g. handles, radiolucent parts. In addition to pure plastics composite materials are also used in some cases, e.g. wood-looking phenolic resin reinforced with fabric for handles of screwdrivers, raspatories, chisels, etc, or carbon-fiber reinforced plastics for aiming arms.

All used plastics are able to withstand correct processing. Some of the plastics can become soft during steam sterilization, but do not undergo permanent deformation at normal sterilization temperatures below 140 °C. The material can, however, be damaged, for example by repeated immersion in disinfectants outside the pH range of 4–9.5 and by oversteering. Also, some rinsing aids can lead to discoloration or embrittlement of plastics and composites by repeated use.

Recommended Temperatures and pH Levels

Material	Temperature *	pH
Stainless steel	Up to 149°C	7–11

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* The recommended processing temperatures take into account material properties and internally validated parameters for processing.

Causes of Corrosion and Surface Change or Damage

The surface of the instruments can be attacked and damaged by incorrect handling or contact with various substances. Awareness of the following possible causes of corrosion and material damage can help to avoid their occurrence.

Blood, pus, secretions etc.

Most human body fluids and residues contain chlorine ions, which can lead to corrosion if left to adhere to, or dry on, the instrument for prolonged periods.

Instruments should therefore be cleaned and dried immediately after every use.

Saline solutions, iodine tinctures, water

The chlorine and iodine ions in these solutions cause pitting corrosion. Keep any contact with these ions to a minimum. Rinse instruments thoroughly with distilled water* to remove all residues.

Normal tap water often also contains chlorides, as well as high concentrations of other minerals, which can form marks with sharply defined edges on the instrument surface. These can usually be removed with distilled water* and non-abrasive stainless steel cleaners. Never leave instruments wet for long periods; always dry them immediately. Condensation moisture produced during sterilization can be avoided by prolonging the drying phase.

* A conductivity of < 0.5 µS is recommended for distilled water.

Detergents, disinfectants, rinsing aids and other additives

Excessive concentrations of these products or strongly acidic or alkaline detergents can attack the protective oxide layer of stainless steel, titanium and aluminum and lead to corrosion, discoloration or other changes of the materials, properties and surface conditions. When using such products, always follow the manufacturer's recommendations in respect of

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concentrations, contact times, temperatures and material compatibility. Products with pH levels between 7 and 9.5 are recommended. During repeated and prolonged use some rinsing aids can attack certain plastics and lead to discoloration or embrittlement. If the instruments are cleaned in an automated washer-disinfector, follow the directions of the manufacturers of the washer-disinfector, detergents, rinsing aids and other additives.

Steel wool, steel brushes, files and other abrasive cleaning tools

Never use extra fine or normal steel wool, steel brushes, files or other cleaning tools with abrasive effect on metals to clean surgical instruments, as this will result in mechanical damage to the passive layer, leading to corrosion and malfunction.

Contact between instruments made from different metals

If stainless steel instruments are left in contact for long periods with surface-damaged instruments and are simultaneously moistened with an electrolyte, rust can form at the points of contact. Steam, water, ultrasonic cleaning solutions or other liquids and solutions can act as electrolytes. Such phenomena are occasionally observed during automated cleaning. Corrosion products that have already formed can also be transferred to other instruments by electrolytes, thereby producing surface rust. If possible, instruments made from differing materials should be cleaned and sterilized separately. Consequently, instruments with corrosion or rust spots must always be excluded and exchanged for unblemished ones. Instruments should be cleaned in their opened and dismantled state in order to avoid not only insufficient cleaning but also crevice and fretting corrosion. The passive layer in crevices or joint gaps can be damaged by chemical or mechanical action, leading to corrosion.

Inadequate lubrication

Moving instrument parts, e.g. joints, sliding parts, dismantlable threaded connections etc. must be regularly lubricated. Constant metallic abrasion increases the damage to the passive layer and thus greatly increases the risk of corrosion.

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Detergent residues on packing cloths

Cloths used to pack the devices for sterilization must be free of detergent or other residues. Such residues can be transferred to the device surface via steam and can interact with the surface.

Overstressing of instruments

Instruments are designed only for a specific purpose and must be used accordingly. Inappropriate use can lead to mechanical overstressing, malfunction and permanent instrument damage, and this in turn increases their susceptibility to corrosion.

Note on latex

Because instruments do not contain any latex, they can safely be used for patients with a latex allergy.

Note on Special Oil

Synthes Special Oil is a synthetic oil and non-toxic. It is recommended to lubricate and maintain instruments only with Synthes Special Oil.

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REQUIREMENTS FOR TECHNICAL SUPPORT AND REPAIR OF MD

Defective instruments can be sent to your local Synthes GmbH customer service for repair.

Customer service will assess whether the instrument can be repaired. Make sure that you enclose a delivery note with the defective instrument containing the following information:

- Hospital address, contact person and telephone number
- Article number of the defective instrument being returned
- Description of the problem.

Your local customer service can deliver spare parts for defective or missing components of simple, multi-part instruments (e.g. depth gauges, drill sleeves). Consult your local customer service for information on the availability of spare parts.

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CONDITIONS OF TRANSPORTATION AND STORAGE**Transportation conditions**

Designation	Range
Temperature	-29°C to +60°C
Humidity	30% to 85%
Barometric pressure	No special conditions

Storage conditions

Designation	Range
Temperature	-29°C to +60°C

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