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Demo samples fully correspond to similar implantable samples in terms of technical characteristics.

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
Table 1 - Compliance of demo samples with implantable ones.

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








Phone	Surgeon's Telephone Number
REF	Catalogue number
LOT	Batch code
QTY 1	Quantity One
STERILE 	Sterilized using steam or dry heat
SN	Serial number
	Use-by date
	Caution
	Do not reuse
	Width, Height, Projection
	Diameter, Projection

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- This product is intended for use on one individual patient during a single procedure. The possibility of damage to the implant and infection exists if a subsequent procedure is performed (such as open capsulotomy, breast pocket revision, etc.). It is the responsibility of the attending physician to determine if a new implant should be inserted. If the implant is damaged, it must be removed.
- Silicone gel can leak or “bleed” through the semipermeable outer silicone envelope into the capsule and adjacent breast tissue. Migration into capillaries has also been reported. The longterm effects of such “bleed” are unknown. Prospective patients should be made aware of this potentiality. (See ADVERSE REACTIONS/EVENTS section of this insert.)
- Only one prosthesis should be implanted per breast. Mentor recommends against the stacking of implants, one upon the other. The devices have not yet been tested for this use and the integrity of the implants cannot be guaranteed as the materials may abrade and wear. Such abnormal stress may result in weakening or rupture of the prostheses.
- Do not insert or attempt to repair a damaged or altered prosthesis.
- The action of drugs (e.g., antibiotics and steroids) in contact with the prosthesis has not been tested by the manufacturer, and their use cannot be recommended. Each physician who chooses to use chemotherapeutic drugs with this prosthesis must assure compatibility of the drug with the silicone elastomer.
- Do not introduce or make injections of drugs or other substances into the implant. Injections through the implant shell will compromise the product’s integrity.

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produce poorer cosmetic results.” The decision regarding the use of radiation therapy following breast implantation should be made by the surgeon and the radiation oncologist.

- The use of the periumbilical approach in placement of the implant has not been established and is not recommended.
- Careful hemostasis is important to prevent postoperative hematoma formation. Should excessive bleeding persist, it is recommended that the device not be implanted until the bleeding is controlled.
- If a physician treats a hematoma or serous fluid accumulation by aspiration, or if a biopsy or lumpectomy is performed, care must be taken to avoid damaging the implant. These procedures present possible risk of implant puncture.
- The incidence of extrusion of the prosthesis has been shown to increase when the prosthesis has been placed in injured areas: scarred, heavily irradiated or burned tissue, crushed bone areas or where severe surgical reduction of the area has previously been performed.
- Excessive fibrous capsular formation or contracture may occur around any implant placed in contact with soft tissues. The incidence and severity of this occurrence may increase if postoperative local hematoma or infection occurs.
- The physician should use personal discretion when deciding to use these prostheses regarding patients who exhibit psychological instability.
- Surgical implantation of a breast implant may interfere with the ability to breast-feed. The Institute of Medicine (IOM) concluded that there is limited evidence that implantation, especially through a peri-areolar incision, may interfere with lactation and breast-feeding. However, it

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should be noted that previous breast reconstruction surgery, such as mastectomy, may be the initial cause of this interference.

- Patients with connective tissue diseases (CTD) may experience an increased risk of wound dehiscence, infection and bleeding (likely due to their ongoing medical management) that may require further treatment.
- Some patients with breast implants have reported experiences of neurological and/or rheumatological diseases. Mentor is currently collecting data to further understand these potential risks and their possible association with breast implants.

Peri-operative Infection and Inflammation Source Reduction Practices:

The following practices could potentially reduce peri-operative bacterial infection and particulate based inflammation:

- Perform a layered closure to separate skin/wound bacteria from the implant cavity.
- Avoid, where possible, pocket dissection that communicates with breast ducts and glands
- Avoid unnecessary adjacent tissue trauma and maintain meticulous hemostasis in the creation of the implant pocket.

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- The silicone elastomer shell may be easily cut by a scalpel or ruptured by excessive stress, manipulation with blunt instruments or penetration by a needle. Subsequent rupture will result. All prostheses should be carefully inspected for structural integrity prior to and during implantation. Meticulous care must be exercised in handling and implanting the device.
- Excessive force and overmanipulation during the implantation procedure are to be avoided as this may damage the implant.
- Any subsequent surgical procedures in the area of the implant should be undertaken with extreme caution as damage to the implant could occur. In the event that an implant is damaged, it must be removed.
- Each device should be checked for patency prior to surgery and continuously monitored throughout the surgical procedure to ensure the structural integrity of the device is not compromised in any way. This prosthesis should not be implanted following any modifications to its original design. A prosthesis which has been damaged, or on which repairs or modifications have been attempted, should not be implanted. A standby prosthesis should be available at the time of surgery.
- Do not contact the device with disposable, capacitor-type cautery instruments as damage to the outer shell of the prosthesis may result.

Patient education and informed consent

The surgical procedures associated with the use of gel breast implants are not without potential complications and risks. The use of this product is an elective procedure. The patient should be counseled prior to surgery regarding the

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benefits and possible risks associated with tissue reconstruction and/or breast augmentation using breast prostheses and alternative procedures. Patients should be advised that breast implants should not be considered lifetime implants and that revision surgery may be necessary including implant removal or replacement.

It is the responsibility of the individual surgeon to decide the best method by which to counsel a patient prior to surgery. Mentor relies upon the surgeon to advise the patient of all potential complications and risks associated with the use of breast implants.

After surgery, surgeons are to provide the patient with the Patient ID Card with information on the implant(s) used. Labels are supplied with the implant to be applied to the Patient ID Card. Patients should be advised to carry the patient card to facilitate medical care in case of emergency.

Following the surgery, the surgeon should inform the patient of necessary postoperative office visits and the need to continue to consult a physician for routine examinations to detect for breast cancer.

Patients should be advised to inform a physician or surgeon of the presence of an implant if any surgery of the breast area is scheduled, and to consult a physician or pharmacist before the use of topical medicines such as steroids in the breast area. In addition, the patient should be instructed to contact her surgeon should she experience any problem related to her breast implants.

INSTRUCTIONS FOR USE

The implantation of gel-filled prostheses for breast reconstruction or augmentation involves a variety of surgical techniques; therefore, the surgeon is best advised to use the method, which his/her own practice and discretion

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around the implant, which, on occasion, will shrink sufficiently to squeeze the implant, producing varying degrees of firmness. At its worst, the implant can feel hard, be painful and/or distorted. This can occur soon after surgery or years later and may be unilateral, bilateral or asymmetric. Surgical release or excision of the scar is often successful, but recurrence is not uncommon. The cause of the contracture phenomenon is poorly understood. In the past, closed disruption of the scar by squeezing the breast was common, but this is rarely practiced today.

Capsular contracture is graded in severity on a scale of I to IV by Baker classification. Calcification of the capsule can also occur. Calcification is a phenomenon that is occasionally seen with long-term scarring especially if there is irritation such as tight burn scars that cross-joints.

Calcified capsules may require removal if the patient wishes relief from her contracture but otherwise seem to be harmless. Small foci of calcification are commonly seen anywhere in the breast parenchyma. They can usually be identified as benign by the radiologist but on occasion may require biopsy to rule out a malignancy.

Rupture of the Implant

Breast implants may not last a lifetime. Although rupture can occur at any time following implantation, a number of studies examining rupture of current generation, single-lumen, CE marked, round and shaped silicone gel-filled breast implants of various manufacturers over time (using MRI evaluations)

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have consistently reported similar results indicating the average expected lifetime

exceeds 10 years.

While the silicone material itself has not been shown to biodegrade, the shell may rupture due to wear and tear, or direct injury. If the implant shell is ruptured, the escaping gel is usually contained by the scar envelope in the surgical pocket (intracapsular) and may be undetectable except by MRI (silent rupture), which is about 85% effective in detecting rupture. If the scar envelope is torn, the gel can be driven into the local tissue planes and breast tissue (extracapsular). Most of the escaped gel remains in the immediate environment of the breast but on rare occasions it has been reported to migrate down the arm, into nerve sheaths or into the abdominal wall. Ultrasound, mammography and physical examination may also diagnose these ruptures, which have escaped the scar envelope. Most of the reported cases occurred in the more fragile, thinner shell devices implanted in the late 1970s.

Rupture should be suspected if there is a change in character of the device such as a new, persistent unilateral burning sensation or a change in softness, texture or shape of the implant. Because of the silent nature of most ruptures and the difficulty of diagnosis without surgical exploration, the true incidence is unknown. Current products have thicker and stronger shells and more cohesive gel contents. Caution should be used when comparing expected or actual rupture rates of current devices to historical incidences, especially when, as is often the case, the brand, vintage and type of device is unknown. Explantation and/or replacement may be indicated if the implant fails, especially if it is seen in the breast parenchyma as it could be confused with or mask a tumor.

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Causes of implant rupture include, but are not limited to: damage from surgical instruments, intraoperative or postoperative trauma, excessive stresses or manipulations as may occur during daily routines such as vigorous exercise, contact athletics, routine manual massage, intimate physical contact and from compression required during mammography.

Changes in Nipple and Breast Sensation/Breast Pain

Any surgery that undermines the skin of the breast can result in hypersensitivity or loss of sensation in the nipple-areolar region. These changes can vary in degree and may be temporary or permanent. Changes in nipple/breast sensation may, on occasion, affect sexual response or comfort while nursing. These changes are believed to be a result of nerve damage or stretching of the nerves from the surgery. There is no specific treatment for this condition.

Most women undergoing augmentation or reconstruction with a breast implant will experience some breast and/or chest pain postoperatively. While this pain normally subsides in most women as they heal after surgery, it can become a chronic problem in other women. Chronic pain can be associated with hematoma, migration, infection, and implants that are too large or capsular contracture. Sudden severe pain may be associated with implant rupture.

Interference with Mammography in Detection of Cancer/Calcium Deposits

As silicone is opaque to x-rays, an implant may theoretically interfere with the early detection of cancer by mammography as it may obscure part of the

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breast. Newer techniques of breast compression improve the amount of breast that can be visualized. Alternatively, most surgeons feel that the device may improve the detection of tumors by palpation. While of considerable theoretical concern, delayed detection due strictly to the presence of an implant has not been reported. Women at high risk of developing breast cancer should consider getting implants with caution. Since the breast is compressed during mammography, it is possible for the implant to rupture, but this is rare and should not deter a woman from regular, routine mammographic screening. Before the mammography exam, women should inform the technologist that they have implants.

Calcium deposits are seen occasionally in old scars anywhere in the body and this is true of the implant capsule. This usually does not occur until years after implant surgery. Benign calcifications are also commonly seen on mammography in otherwise normal breast parenchyma even in breasts that have never been operated on. These benign calcium deposits usually have a different x-ray appearance than the calcifications that signal a malignancy. An expert radiologist can usually determine if a calcium spot is benign or malignant but occasionally a biopsy may be necessary to rule out malignancy. There is no evidence that these deposits occur more or less frequently in women with implants than those who have not received implants. After many years, some patients may develop a thin layer of calcium in the scar capsule that surrounds the implant. This is almost always associated with capsular contracture but otherwise causes no known problem.

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Extrusion of the Implant/Interruption of Wound Healing

Skin necrosis and/or sloughing resulting in implant exposure may occur from undue tension of the tissues overlying the implant, the overdosage of steroids placed in the implant pocket or surgical or external trauma.

This is quite rare in the augmentation patient but occurs occasionally in the more challenging tissue environment of reconstruction after mastectomy, local tissue injury and/or radiation. The areas of the scar, especially after radiation to that area, seem to be the most vulnerable to disruption. Extruded implants may need to be surgically removed.

Wrinkling of the Implant/Dissatisfaction with Cosmetic Results/Asymmetry

Visible and/or palpable wrinkles in an implant are related to the thinness of the overlying tissue cover, the degree of capsular contracture and the texturing of the implant shell surface. Traditional smooth-walled gel devices rarely demonstrate wrinkles. Some have described camouflaging of wrinkling with autologous fat transfer superficial to the implant. Surgical error, preexisting asymmetry or deformity, keloid formation of the incisional scar, vagaries of time, weight gain or loss, pregnancy and nursing can all contribute to an immediate or late poor aesthetic result. With time, most breasts, with or without implants, become ptotic to some degree. Asymmetry is usually related to the inability to totally correct for preexisting disparity between the two breasts. It can also be attributed to surgical error, asymmetric contracture, or rupture of the implant.

Excessive buildup of collagen at the incision site during the healing process causes some patients to develop scars of cosmetic concern. Keloid scars, which do not respond well to treatment, often extend beyond the edges of the

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original scars and can continue to enlarge over time. Hypertrophic scars are generally confined to the original site and respond well to scar revision treatment, which may include steroid injections to break down the collagen or surgery to revise the position, direction or line of the scar.

Possible Reactions to Silicone

This text contains a brief summary of information from the medical literature. Mentor recognizes that the information contained within this text is highly technical. However, medical ethics and practice dictate that the physician must be an intervening party between the manufacturer of prescriptive medical devices and the patient. The issue of the possible relationship between silicone (and other implantable materials) and various diseases has been the subject of significant scientific debate. Concerns include immunological and neurological disorders, carcinogenicity and connective tissue disorders.

Despite numerous anecdotal reports of established and newly described connective tissue diseases attributed to the device, multiple epidemiological studies have very consistently shown no significant association between silicone breast implants and any established or new connective tissue or other immunological diseases. Four prestigious multispecialty scientific expert panels have

reviewed the published literature on this topic, specifically as it relates to silicone breast implants, and have issued lengthy reports on their findings. These expert panels include the Independent Review Group (commissioned by the Chief Medical Officer of the U.K.), the National Science Panel (appointed by Judge Pointer for MDL 926), the IOM and the Scientific and Technical Options Assessment (STOA) Programme (commissioned by the European Parliament). These four panels have uniformly concluded that there was no discernible evidence of causal association or positive risk ratio between exposure to silicone breast implants and recognized or new autoimmune or connective tissue diseases. Essentially the same finding was reported in a 2011 updated review of the epidemiology literature published in the peer-reviewed medical literature in 2011.

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Coincident neurological problems such as multiple sclerosis and amyotrophic lateral sclerosis (ALS) have occurred in a small number of breast implant patients. Subsequent studies have not demonstrated a relationship between any neurological disease and breast implants.

Several studies have been conducted to determine the carcinogenic risk of breast implants and no evidence of increased risk of cancer has been demonstrated. Continuing assessment of both known and possible risks associated with breast implant surgery is ongoing.

No credible reports on birth defects or other reproductive effects in humans associated with implantation of silicone breast implants of any type have been identified in the literature. Recent studies sponsored by Mentor provide further evidence that silicone materials used in breast implants do not cause adverse reproductive effects in experimental animals.

Although any breast surgery, including breast implantation, could theoretically interfere with the adequacy of a woman's milk supply, many women with breast implants have nursed their babies successfully. It is known that any breast surgery such as breast biopsy can affect the quantity of milk produced. In recent years, the question has been raised regarding the potential transfer of silicone into the breast milk of women with silicone breast

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prostheses and possible effects on the health of the breast-fed children. However, more recent studies have provided strong evidence of the lack of association between silicone breast implants and adverse effects in breast-fed children. The American Academy of Pediatrics has stated that there is no reason why a woman with implants should refrain from nursing. The European Committee on Quality Assurance and Medical Devices in Plastic Surgery (EQUAM) concluded in 2000 that silicone gel-filled breast implants do not adversely affect pregnancy, fetal development, breast-feeding, nor the health of breast-fed offspring.

The report sponsored by the IOM, "Safety of Silicone Breast Implants," released in July 1999 states that women with silicone breast implants are no more likely than the rest of the population to develop cancer, immunologic diseases, or neurological problems. The committee also concluded that there is no evidence that mothers with implants pass silicone on to infants when breast-feeding.

Anaplastic Large Cell Lymphoma

Based on information reported to global regulatory agencies and found in medical literature, an association has been identified between breast implants and the development of anaplastic large cell lymphoma (ALCL), a type of non-Hodgkin's lymphoma. Women with breast implants have a very small but increased risk of developing Breast Implant Associated ALCL (BIA-ALCL) in the fluid or scar capsule adjacent to the implant, with documented potential

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for local, regional, and distant spread of the cancer with mortality reported in rare cases.

BIA-ALCL has been reported globally in patients with an implant history that includes Mentor's and other manufacturers' breast implants with various surface properties, styles, and shapes. Most of the cases in the literature reports describe a history of the use of textured implants. Several journal articles explore the risk factors for BIA-ALCL, including the varied methods used to create surface texture of the implant and the role of biofilm in causing disease, among others.

You should consider the possibility of BIA-ALCL when a patient presents with a late onset seroma (at least 1 year and on average 7 to 10 years following implantation), breast mass adjacent to the implant capsule or lymphadenopathy. Work-up of delayed enlargement of a previously implanted breast includes an ultrasound or MRI to look for the presence of seroma or lymphoma. Diagnosis starts at presentation, with needle aspiration of as much of the fluid as possible (minimum 50 mL), using ultrasound for implant displacement and protection. A suspicious mass requires needle or open tissue biopsy. Specimens should be sent for cell morphology by cytology, CD30 immunohistochemistry and ALK (anaplastic lymphoma kinase), and flow cytometry for evaluation, quantification, and characterization of T cells within the specimen. Hematopathology consultation at a tertiary cancer center is strongly encouraged to establish or exclude a diagnosis of BIA-ALCL. Once the diagnosis of BIA-ALCL has been

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established, physicians are strongly encouraged to consult with a multi-disciplinary team including oncologists, pathologists, surgical oncologists, and plastic surgeons. The United States National Comprehensive Cancer Network (NCCN) guidelines for treatment of established BIA ALCL includes the removal of the implant and the entire surrounding implant capsule in addition to associated masses and lymph node involvement.

Please report all confirmed cases of BIA-ALCL per the details provided in the ADVERSE EVENT REPORTING.

Please note for latest diagnostic and treatment options refer The National Comprehensive Cancer Network (NCCN) guidelines

Gel Bleed

The gel in an implant consists of large three-dimensional mesh-like molecules that constitute about 20% of the total weight of the gel. The spaces within these molecules are filled with medical grade viscous silicone fluid. This viscous fluid is similar to the materials available in many products including antiflatulence medication, available without prescription for infants and

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adults. A small amount of this material can diffuse or bleed through the shell of the implant. The major portion of this material stays on the implant wall. A small amount moves into the scar capsule where it is gradually picked up by certain “scavenger” cells of the body’s immune system, called macrophages. Normally, these cells try to destroy foreign material such as bacteria. However, if the material (such as silicone) cannot be destroyed, it is carried to the lymph glands by the macrophages. As noted in the section on Possible Reactions to Silicone, the STOA report published in 2000 concluded: “Studies do not point to an association between silicone implants and serious health risks, such as cancer and connective tissue diseases.”

Granulomas

It is possible for a granuloma to form around a tiny amount of silicone. Although these lumps are noncancerous, they can be difficult to distinguish from cancerous lumps without being removed (biopsied) and examined. Granulomas, if large or suspicious for malignancy, may need to be biopsied or surgically removed and examined.

Other Possible Reactions

- Thrombosed veins, resembling large cords, have temporarily developed in the area of the prosthesis and have resolved without surgical or medical therapy.
- Pain from an improperly sized and/or placed implant, such as from compression of nerves or interference with muscle movement, may occur.
- Hypertrophic scarring has been reported.
- The prosthesis may become difficult to explant if the degree of tissue adhesion is significant.

Instructions and Precautions for Implant Removal

Standard surgical approaches and practices should be used should implant removal be necessary. The entirety of the device needs to be removed. In the

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case of a ruptured gel implant, all of the gel should be removed as far as possible.

INFORMATION A PHYSICIAN SHOULD PROVIDE TO THE PATIENT

Breast implantation is an elective procedure and the patient must be well counseled on the risk-benefit relationship. The surgeon should provide each prospective patient with the following:

- **Gel-Filled Breast Implant Surgery: Making an Informed Decision**

This brochure can be used to facilitate patient education in the risks and benefits of gel-filled breast implant surgery. The patient should be advised to wait a week after reviewing and considering this information before deciding whether to have augmentation surgery.

- **Device Materials**

MENTOR Gel-Filled Breast Implants have met biocompatibility requirements per ISO 10993. Materials to which patients can come in contact with are 100% medical implant grade silicone. The shells of the MENTOR Gel-Filled Breast Implants are constructed from medical-grade silicone elastomer. The shells are filled with medical-grade, cohesive silicone gel.

- **Expected Device Lifetime and Follow-up**

Whether undergoing augmentation or reconstruction, be aware that breast implantation may not be a one-time surgery. Additional doctor visits are likely needed and additional surgery may be needed over the course of a breast implant patient's life. The longer people have breast implants, the greater the chances are that they will develop complications, some of which may require additional surgery.

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Breast implants are not considered lifetime devices. Implant removal with or without replacement may be needed over the course of a breast implant patient's life.

Based on the below information, MENTOR MemoryGel and Contour Profile® Breast Implants (CPG) can be expected to have an average lifetime of over ten years.

The expected lifetime duration of MENTOR Breast Implants is based on the estimated rupture rates of the implants. There are many causes of breast implant rupture, including wear and tear on the indwelling device over time due to normal body motion, unintentional sharp instrument damage during surgery, and physical stress such as that caused by car accidents or strenuous activity. Rupture rates can vary depending upon the type of breast implant surgery performed. In MENTOR clinical studies, patients undergo MRI screening on a regular basis to assess whether rupture of the implant has occurred. It is important to note that most patients had no symptoms related to their reported implant rupture.

Magnetic Resonance Imaging (MRI) Compatibility

MRI is one of the imaging methods available to assess the integrity of breast implants and the extent of possible silicone leakage. Owing to its greater sensitivity, MRI is typically the recognized standard method used to evaluate the integrity of breast implants, usually in patients that are symptomatic for rupture (Ahn et al., 1993). MRI offers an operator-independent, global view of the entire implant and surrounding tissues in a variety of planes, providing a modality for implant evaluation and preoperative planning. In a study of 100 symptomatic patients, Ahn et al. (1993) established magnetic resonance imaging (MRI) as a definitive, reliable, and reproducible technique to

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diagnose both intracapsular and extracapsular ruptures. MRI also enabled the identification of the artifacts that may interfere with the definitive diagnosis of implant rupture.

MR SAFE – Implantable devices Mentor which don't contain metal parts are MR SAFE.

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Estimated rupture rates over time for MENTOR Breast Implants, as observed in MENTOR clinical studies.

These rates are calculated by the Kaplan-Meier method of analysis. The Kaplan-Meier method attempts to account for patients not returning for follow-up and adjusts the estimated rupture rate accordingly. These rupture rates include both ruptures that are “suspect” and those that are “confirmed” as defined below.

Suspected ruptures are ruptures that are suspected based on the MRI imaging findings, but which have not been confirmed by surgical explantation and physical examination of the implant, usually because the patient chose not to have the implant removed.

Confirmed ruptures are those verified by examining the implant once removed from the patient.

Of the overall ruptures (those confirmed as well as those suspected but not confirmed) reported for each device, the proportion of confirmed ruptures was 64% (49 of 77 suspected or confirmed ruptures) for MENTOR MemoryGel Breast Implants, 47% (8 of 17 suspected or confirmed ruptures) for Contour Profile Breast Implants (CPG).

10-year suspected or confirmed implant rupture rates for MENTOR MemoryGel Breast Implants:

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Overall, the 10-year follow-up rate was 53%, thus data accuracy is limited. 10-year suspected or confirmed implant rupture rates for Contour Profile Breast Implants (CPG):

Primary Augmentation	Revision Augmentation	Primary Reconstruction	Revision Reconstruction
3,3%	4,7%	11,2%	0%

Overall, the 10-year follow-up rate was 45%, thus data accuracy is limited.

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Patient ID Card

Enclosed with each gel-filled breast implant is a Patient ID Card. To complete the Patient ID Card, stick one Patient Record Label for each implant on the back of the Patient ID Card. Patient Record Labels are located on the internal product packaging attached to the label. If a Patient Record Label is unavailable, the lot number, catalog number and description of the device may be copied by hand from the device label. The patient should be provided with the Patient ID Card for personal reference.

Information specific to the implant, are included on the implant card as well as within patient records as kept by the health care provider. The physician shall provide the patient card containing information to identify the product-specific information related to your device.

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CONDITIONS OF TRANSPORTATION AND STORAGE

Table 3 – Transportation conditions

Designation	Range
Temperature	-20°C to +50°C
Relative humidity	not more 75%
Barometric pressure	700-1060 hPa

Table 4 – Storage conditions

Designation	Range
Temperature	+10°C to +25°C
Relative humidity	not more 75%
Barometric pressure	700-1060 hPa

MD DISPOSAL RULES AND CONDITIONS

The product belongs to class B of the classification of medical waste. Uncontaminated packaging belongs to class A of the classification of medical waste. After last use the product and its package must be disposed of in accordance with local regulations as specified in the requirements SanPiN 2.1.7.2790-10 and state requirements.

WARRANTY

Mentor disclaims any warranties, whether expressed in writing or orally, whether express or implied by the circumstances, including, but not limited to, any implied warranties of merchantability, fitness for purpose or design, except for those set forth herein and to the extent permitted by law. Mentor is not responsible for any direct, incidental or indirect loss, damage or expenses, directly or indirectly, as a result of using this device to the extent permitted by

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