B. BRAUN
SHARING EXPERTISE

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Etablissement de CHASSENEUIL-DU-POITOU

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86360 CHASSENEUIL-DU-POITOU

TO THE
FEDERAL SERVICE ON SURVEILLANCE IN HEALTHCARE
(ROSZDRAVNADZOR)

Covering letter

Poitiers

March 13th, 2024

With this letter we

20 MARS 2024

visé n° 025

B. BRAUN MEDICAL, 26, rue Armengaud, 92210 Saint-Cloud, France, as the responsible/ legal manufacturer of the medical device

Celsite® Access Ports with installation sets

declare that we have considered the remarks stated in notification no. 10-10615/24 from 15.02.2024 issued by Federal Service on Surveillance in Healthcare (Roszdravnadzor). These remarks do not affect the Instructions for Use and content of Instructions for Use from 2020 (signed 18/02/2021) is still valid.

Based on abovementioned we provide the IFU in English language with Annex 1 to the IFU:

1. IFU for Celsite® Access Ports with installation sets

The Annex 1 to the IFU contains the following items*:

- 1.1. Name of the medical device
- 1.2. Name of manufacturer and manufacturing site
- 1.3. Operating conditions of the medical device in the internal environment of the human body
- 1.4. Conditions for transportation and storage
- 1.5. Delivery set and packaging data
- 1.6. The information about the symbol the conformity Declaration of conformity in the Russian Federation
- 1.7. The information about the symbol of non-pyrogenic
- 1.8. Shelf life for the medical device
- 1.9. Warranty/Complaints
- 1.10. Procedure for disposal or destruction of a medical device
- 1.11. Intended Use
- 1.12. Target category of patients / Target group of users
- 1.13. Model range of the medical device

*- the description of data, see at Annex 1 to the present IFU.

for updating the documents of Registration Dossier № ФС3 2008/03191 issued on January, 24th, 2022 for the medical device Celsite® Access Ports with installation sets and submission of abovementioned document to the Federal Service on Surveillance in Healthcare (Roszdravnadzor).

For and on behalf of B. BRAUN MEDICAL

Sylvie PLUMET

Competence Center Chasseneuil
Director
B. Braun Medical

Catherine BOISMENU

Competence Center Chasseneuil
Deputy Director in charge of Quality and
delegated Regulatory Affairs
B. Braun Medical

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

1/28

Contents


(en) Contents: Refer to unit package label for actual content list.

(en) Access Ports

Celsite® Discreet

Celsite® IMPLANTOFIX




Celsite® Concept


Celsite®




1 F = 0,33 mm

 Connection rings (x2)	 J guide wire	 Syringe
	 Catheter	  
 Thin wall needle 18G/20G	 Vein lifter	 Winged Surecan®

 Surecan® needle 22G (x2)	 Dilator 6F	 Tunnelling rod
 Peelable introducer	 ECG cable	 Introcan® i.v. cannula 20G
 Screw connector, connection spanner		

Celsite® ST304

 Access port		
 Screw connector, connection spanner	 Catheter PUR, anti-kink device	 Catheter PA, anti-kink device
 Syringe	 Syringe LOR	 Sterican® needle 20G 70 mm
 Surecan® needle 22G (×2)	 Filter 0,2 μ	 Tunnelling rod
 Tuohy needle 16/18G	 Scalpel E10, E11	 Winged Surecan®

**INSTRUCTIONS FOR USE
Celsite® ACCESS PORTS**

**PLEASE READ CAREFULLY
SUMMARY**

- I – INDICATIONS
- II - CONTRAINDICATIONS
- III - POTENTIAL COMPLICATIONS (Immediate and delayed)
- IV - WARNINGS V - SPECIAL ATTENTION V-1 PATIENT
- V-2 CATHETER POSITION
- V-2-1 Venous catheters
- V-2-2 Percutaneously placed Arterial catheters (Celsite® Interventional with Anthron® catheter)
- V-3 DRUG INCOMPATIBILITIES
- V-4 NEEDLES
- V-4-1 Puncture of the Celsite® Access Port
- V-4-2 Puncture of the vein
- V-5 OCCLUSION
- V-6 INFECTION
- VI - IMPLANTATION TECHNIQUES
- VI-1 CENTRAL VENOUS IMPLANTATION TECHNIQUES
- VI-1-1 Percutaneous Technique
- VI-1-2 Surgical Cut-down Technique
- VI-1-3 References requiring modifications to the above techniques
- VI-1-3-1 Over-the-Wire Technique
- VI-1-3-2 Brachial Implantation Technique
- VI-1-3-3 ECG Technique
- VI-1-3-4 Braunule Technique
- VI-1-3-5 Valved Catheter
- VI-1-3-6 Pre-connected catheter and port
- VI-2 ARTERIAL IMPLANTATION TECHNIQUES (FOR H.A.I.T.)
- VI-2-1 Surgical Implantation Technique
- VI-2-2 Percutaneous Implantation Technique
- VI-2-2-1 Catheter tip in gastro-duodenal artery
- VI-2-2-2 Catheter tip in hepatic artery
- VI-3 PERITONEAL AND PLEURAL IMPLANTATION TECHNIQUES
- VI-3-1 Peritoneal – Surgical Implantation Technique
- VI-3-2 Peritoneal – Percutaneous Implantation Technique
- VI-3-3 Pleural – Percutaneous Implantation Technique
- VI-4 SPINAL IMPLANTATION TECHNIQUES
- VII - USE OF THE PORT AND CATHETER
- VII-1 VENOUS PORTS
- VII-2 ARTERIAL PORTS
- VIII - MAINTENANCE OF THE PORT AND CATHETER
- VIII-1 RINSING AND HEPARINISATION
- VIII-1-1 Venous ports
- VIII-1-2 Arterial ports
- VIII-1-3 Epidural or Intra-theal ports
- VIII-1-4 Peritoneal ports
- VIII-1-5 Pleural ports
- IX - DURATION OF IMPLANTATION AND REMOVAL OF THE PORT
- IX-1 DURATION OF IMPLANTATION OF THE SYSTEM
- IX-2 REMOVAL OF THE SYSTEM
- IX-3 SPECIAL ATTENTION: BABYPORT®, BABYPORT® S, BABYPORT® PC AND SMALL PORTS
- X - HIGH PRESSURE INJECTION

I - INDICATIONS

- Venous Access Ports are intended to be used for repeated intra-venous administration of, for example, chemotherapy, antibiotic and anti-viral drugs, parenteral nutrition, blood sampling or transfusion. Certain Celsite® Access Ports may be used for Contrast Enhanced Computerised Tomography (CECT) using high pressure injection (see "X - High pressure injection"). For full list of validated references please refer to the tables.
- Arterial Access Ports are intended to be used for intra-arterial administration of chemotherapy.
- Epidural or intra-theal ports are intended to be used for spinal administration of pain relieving drugs.
- Peritoneal ports are intended to be used for intra-peritoneal administration of chemotherapy.
- Peritoneal/Pleural ports are intended to be used for intra-peritoneal administration of chemotherapy, hydration, drainage of malignant ascites, or drainage of malignant pleural effusion.

II – CONTRAINDICATIONS

The following are contraindications to port placement:

- Known infection, bacteraemia or septicaemia.
- Known allergy to any of the materials contained in the access port or catheter.
- If the medications to be used in the access port are incompatible with any of the materials contained in the access port or catheter.
- If the patient's anatomy does not allow the insertion of the catheter into the chosen access site or if the patient has had previous radiotherapy in the chosen area.
- Previous venous thrombosis.
- Heparin induced thrombocytopenia (Celsite® Interventional with Anthron® catheter only).

III – POTENTIAL COMPLICATIONS (Immediate and delayed)

Air embolism Angiospasm (arterial) Angitis (arterial) Aneurysm/false aneurysm (arterial) Brachial plexus injury Cardiac arrhythmia, tamponade, perforation Catheter disconnection, rupture or fragmentation Catheter occlusion	Cerebral infarction or thrombosis (arterial) Embolism Explantation of the port secondary to inflammatory reaction Fibrin sheath formation Haematoma Haemothorax Heparin Induced Thrombocytopenia (HIT) Infection/Sepsis	Liver failure (arterial) Pericarditis (arterial) Pneumothorax Port/Catheter migration Thromboembolism/Thrombophlebitis Thrombosis Ulcer related to cytotoxic drugs (arterial)
---	--	--

Before implantation of the access port system these, and other, well documented potential complications should be considered.

IV - WARNINGS

During storage avoid freezing or excessively high temperatures, storage at room temperature is recommended. Anthron® catheters should be protected from light.

Celsite® Access Ports are for single use only, do not re-sterilise the product, or any of the components, and destroy after use. The device and its accessories are not re-usable nor designed to be re-used. Any re-use would definitely compromise the performance and safety of the device. The product is sterilised using Ethylene Oxide and is sterile and non-pyrogenic in unopened, undamaged, individual packages.

Refer to unit package label for actual content list.

Check the expiry date before use and never implant a device past the expiry date.

Do not remove from the unit carton until ready for use and ensure that all handling and use is under strict aseptic conditions.

This product should not be implanted or accessed except by personnel experienced in the technical and clinical aspects of access ports.

Special Warning for Anthron® catheter:

- a) The Anthron® catheter is made of PUR with a heparin coating.
- b) The Anthron® catheter should not be used for blood sampling for coagulation studies as the results may be affected by the release, however small, of heparin.
- c) As for all heparin, regular supervision of platelet levels is required during the use of the catheter due to the risk of heparin induced thrombocytopenia (HIT).

V - SPECIAL ATTENTION

V-1 PATIENT

Ensure that use of the access port is indicated taking into consideration the patients' anatomy, needs and proposed use of the port.

Check the condition of the skin over the proposed implantation and puncture sites.

For venous catheter placement, place the patient slightly inclined head down.

V-2 CATHETER POSITION

To ensure that the system works correctly, there should be no kinking of the catheter.

It is recommended that the catheter is tunnelled sub-cutaneously to the port.

The position of the tip of the catheter may change when the patient moves.

V-2-1 Venous catheters

a) The catheter tip should always be positioned in the Superior Vena Cava at the entrance to the right atrium taking into account possible catheter tip movement (see Fig. 1). A radiograph should be taken to verify catheter position and exclude haemothorax or pneumothorax.

b) Particular attention should be paid when the catheter is to be implanted via the sub-clavian route. It is recommended that the catheter be inserted outside the costo-clavicular space (see Fig. 2). Catheter ruptures have been observed with the sub-clavian route, the associated risk of extravasation of the infused drugs and embolisation of the distal extremity may have serious consequences. This is due to rupture of the catheter secondary to the catheter being pinched in the costo-clavicular space (Pinch-Off syndrome). It is necessary to pay particular attention to catheters used for long periods or for ambulatory patients rather than short term catheters. Catheter rupture can occur with silicone and polyurethane catheters.

The following clinical signs may suggest catheter pinching:

- The patient needs to lift the arm to permit infusion.
- Intermittent malfunction of the catheter, such as difficulty with aspiration or infusion.
- Sub-clavicular pain or swelling during infusion.
- Palpitations or chest discomfort may indicate catheter fracture.

Remove any sub-clavian catheter which presents any of the signs of Pinch-Off.

A radiograph with contrast injection may be useful to detect catheter abnormalities at the costoclavicular level.

V-2-2 Percutaneously placed Arterial catheters (Celsite® InterVentional with Anthron® catheter)

a) It is recommended that the catheter be placed via the femoral artery (see Fig. A). The entry point into the artery should be below the inguinal ligament. The access port should be positioned on the patient's thigh, medial to the iliac crest, to avoid kinking of the catheter (see Fig. B). The catheter should not traverse the inguinal ligament sub-cutaneously.

b) If needed, access may be via the axillary (see Fig. B, 2), acromiothoracic (see Fig. B, 3) or

9

brachial artery (see Fig. B, 4). Care should be taken with these routes in patients with diabetes, hyper-lipemia and hyper-tension. The risk of cerebral embolism is higher in these patients.

c) Placement in the sub-clavian artery is not recommended as this may lead to catheter displacement and subsequent fracture, and carries an increased risk of cerebral thrombosis or infarct in case of catheter thrombosis.

V-3 DRUG INCOMPATIBILITIES

Never use incompatible drugs together or in sequence: the port should be rinsed of all residues with 10 mL sodium chloride (NaCl) 0.9% between each medication.

V-4 NEEDLES

V-4-1 Puncture of the Celsite® Access Port

Always use special bevelled needles, such as Surecan®, Winged Surecan® or Cytocan, which are adapted for use with implantable access port systems.

Normal hypodermic needles will damage the septum and may cause leakage of the system or blockage due to small silicone particles.

V-4-2 Puncture of the Vein

Certain references contain a Safety Seldinger needle. The Safety Seldinger needle is designed to reduce the risk of accidental needle stick injury. Care must be taken to avoid needle sticks in accordance with local Infection Control standards (and avoid the risk of exposure to contaminated blood).

Access the vein with the needle in the normal manner. After withdrawal of the Safety Seldinger needle activate the safety clip. To do this grasp the clear clip housing, which houses the metal safety clip, with your free hand.

Advance the clear clip housing forward along the Safety Seldinger needle towards the needle bevel in one continuous motion. Advance the clip housing over the needle bevel activating the safety mechanism. The clip will remain on the needle bevel. The clip housing will separate from the clip and needle. Dispose of the needle and the clear clip housing immediately into sharps container.

Do not attempt to re-insert the needle into the patient or to re-cap the needle with the protective guard.

V-5 OCCLUSION

Always use a syringe equal to or greater than 10 mL volume when rinsing the port to avoid generating excessive pressures in case of occlusion of the system.

In case of obstruction of the system never try to clear the blockage using a fluid under high pressure which carries the risk of catheter fracture and migration. According to validated protocols, and under medical supervision, 2 mL of 70% alcohol may be used to aid the unblocking of silicone catheters when the blockage is due to lipid deposits. The use of alcohol with PUR catheters is not recommended.

According to validated protocols, and under medical supervision, 2 mL of Hydrochloric acid (HCl) 0.1 mol/L may be used to aid the unblocking of both silicone and PUR catheters when the blockage is due to mineral deposits.

V-6 INFECTION

In case of infection, if appropriate treat with antibiotic drugs. If this fails, or is not appropriate, the catheter and port should be removed.

VI – IMPLANTATION TECHNIQUES

General considerations:

- To prevent shearing of the guide-wire never withdraw the J guide-wire through the Seldinger needle.
- Remove the guide-wire and the dilator together. Do not remove the guide-wire through the dilator as this may result in the guide-wire unravelling.
- To prevent shearing of the catheter, never withdraw the epidural/spinal catheters through the Tuohy needle.
- It is recommended that NaCl 0.9% (with or without heparin according to local protocols) be used to flush, and lock the port and catheter during implantation (see Fig. 3).
- During implantation ensure that the catheter is not damaged by unguarded forceps, suture needle or other sharp instruments.
- If using the Celsite® Concept ports (with silicone inserts) it is recommended that the suture be placed within the square indicated on the silicone.
- The catheter should be mounted on the exit cannula along its axis (see Fig. 11 a-b-c) and not across the axis and should be completely mounted on the exit cannula before the connection ring is slid over the catheter.

VI-1 CENTRAL VENOUS IMPLANTATION TECHNIQUES

VI-1-1 Percutaneous Technique

- Insert the Seldinger needle into the chosen vein; verify the position by observing blood reflux (see Fig. 4).
- Insert the J guide-wire into the vein, when the correct position is reached, withdraw the needle (see Fig. 5).
- Thread the assembled dilator and peelable sheath over the guide-wire using a twisting motion to pass through the skin planes (see Fig. 6).
- Remove the dilator and guide-wire from the vein and insert the catheter through the peelable sheath (see Fig. 7) to the desired position at the junction of the Superior Vena Cava and the right atrium (see Fig. 1). Check catheter tip position using fluoroscopy.
- Remove the peelable sheath from the vessel by pushing down on the 2 handles of the sheath, then peel the sheath from the catheter ensuring that the catheter is not dislodged from the vessel (see Fig. 8).
- Prepare the port pocket at the chosen site, the port should lie approximately 1/2 - 1 cm below the skin surface away from the injection site. Cut the luer connector from the catheter and attach to the tunnelling rod (see Fig. 9).
- Tunnel the catheter from the puncture point to the port pocket (see Fig. 10); ensure there is no kinking of the catheter. Cut the excess catheter (at right-angles) prior to connecting to the port.
- Catheter port/connection:

- Celsite®, Celsite® Concept, Celsite® Discreet:

Slide the connection ring over the catheter, firmly push the catheter onto the exit cannula ensuring the catheter covers the length of the exit cannula, slide the connection ring over the catheter and exit cannula. The connection ring should be in contact with the port (see Fig. 11c).

- Celsite® IMPLANTOFIX:

Slide the screw connector over the catheter, firmly push the catheter onto the exit cannula ensuring the catheter covers the length of the exit cannula, slide the screw connector over the catheter and exit cannula. The screw connector should be tightened by hand and the spanner used to give the final tightening to the securement device, until an audible click is heard, to ensure correct attachment between the port and catheter.

- Celsite® Interventional with Anthron® catheter:

If necessary, cut the anti-kink portion of the connection ring to the desired length (a minimum of 1.5 cm should remain), then slide the connection ring over the catheter, firmly push the catheter onto the exit cannula ensuring the catheter covers the length of the exit cannula, slide the connection ring over the catheter and exit cannula. The connection ring should touch the port.

- i) Confirm catheter connection by gently pulling on the catheter.
- j) Insert the port into the port pocket, paying attention not to kink the catheter.
- k) Confirm catheter patency (by ensuring that both aspiration and injection are possible), ensure that the skin incision is 1 cm from the injection site.
- l) It is possible to fix the access port to the fascia with sutures.

- Celsite® DRAINAPORT:

Both the connection ring and the catheter must be threaded together on the exit cannula to facilitate connection.

VI-1-2 Surgical Cut-down Technique

- a) The skin is incised and the tissue dissected (see Fig. 12).
- b) The selected vessel is identified and tapes passed under it. The vessel is incised with a scalpel or with scissors (see Fig. 13).
- c) The previously flushed catheter is inserted into the vessel, using either a small vessel dilator or vein lifter (see Fig. 14).
- d) Follow from point e) - l) of the preceding paragraph: "Percutaneous Technique".

VI-1-3 References requiring modifications to the above techniques

VI-1-3-1 Over-the Wire Technique

This technique can only be used with the references ST301OTW and Celsite® IMPLANTOFIX references 04438604, 04438647.

Follow steps a) and b) above: "Percutaneous Technique".

c) Thread the dilator over the guide-wire using a twisting motion to pass through the skin planes. Remove the dilator (see Fig. 15).

d) Flush the catheter with sodium chloride (NaCl) 0.9 %, and pass the catheter over the guide wire to the desired position at the junction of the Superior Vena Cava and the right atrium (see Fig. 1).

e) Follow from point f) - l) above: "Percutaneous Technique".

VI-1-3-2 Brachial Implantation Technique

Brachial ports and catheters can be placed using either the OTW (Over-the Wire), Percutaneous or Surgical Cut-down Technique, (see Fig. E).

Select the appropriate vein on the anterior aspect of the elbow, normally the Basilic or Cephalic vein.

Ensure that the access port is not placed over the mid-line of the selected vein. The port and catheter should be positioned above the elbow on the internal surface of the arm.

VI-1-3-3 ECG Technique

This technique can only be used with the ECG references which have specially adapted accessories.

Percutaneous Technique: (see Fig. E)

- a) Place ECG pads on the patient's chest.
- b) Follow steps a) to c) above: "Percutaneous Technique".

Note: A separate J-guide wire (0.035") is needed for access to the vein.

c) Prepare the catheter by inserting the specially marked J-guide wire, which is contained in the kit, into the catheter.

12

d) Remove the J-guide wire and dilator from the sheath and insert the catheter and guide-wire into the peelable sheath. The guide-wire should be pulled back into the catheter during placement into the vein.

e) When approximately 10-13 cm of the catheter has been placed into the vein the guide-wire should be advanced into the catheter until the mark on the guide wire reaches the level of the luer connector on the catheter.

f) When the mark on the guide-wire is at the level of the luer connector the tip of the guide-wire is in the desired position at the tip of the catheter.

g) The sterile ECG cable is attached to the guide-wire, connected to the monitor via a universal adapter and switched to internal monitoring mode. The catheter and guide-wire are advanced using the ECG image to monitor progress.

h) As the catheter approaches the right atrium the P wave starts to elevate ❶. Maximal P wave height is reached and maintained when the catheter enters into the right atrium. Note that the amplitude of the P wave will decrease if the catheter is advanced down into the inferior vena cava. After identifying the area where the P wave begins to develop its maximal amplitude (which corresponds anatomically to the junction between superior vena cava and the right atrium) advance the catheter a further 2 cm ❷. This is the final position of the catheter tip. As the procedure is performed with the patient in the supine position this catheter position allows for the 2-3 cm cranial movement of the catheter which occurs when the patient is upright.

i) Remove the guide-wire while maintaining the catheter in place, flush and clamp the catheter. Follow from point e) - l) above: "Percutaneous Technique".

A post-operative radiograph should be taken to confirm catheter position and exclude haemothorax or pneumothorax.

Surgical Cut-down Technique: (see Fig. F)

This technique combines the Surgical Cut-down Technique with the ECG technique described above.

If no ECG change is noted catheter tip position has to be confirmed with fluoroscopy as catheter tip malposition is to be suspected.

The Celsite® Access Port with ECG technique should not be placed if any of the following is applicable to the patient:

- The presence of a cardiac pacemaker.
- Known cardiac abnormalities, or heart rhythm disturbances.
- Access through the Inferior Vena Cava.

With regard to the ECG lead:

Safety regulations according to IEC 601-1 and NF EN 60601-2-25.

Protection against the risks of current leakage:

Only use Class 1 ECG monitor CE marked, with floating input (Type CF).

Connect additional potential equalisation (earth) on the monitor.

Protection against the risks of electrostatic discharge:

- Use only in rooms protected from electrostatic discharge (conductive flooring).
- Theatre shoes should be made of a recommended polymer material.
- Do not allow the J guide-wire to come into contact with external conductive parts.
- Before use, check that all connections are secure.
- Do not autoclave.
- Do not immerse in liquid.
- Repairs must be carried out only by the manufacturer.

VI-1-3-4 Braunule Technique

- a) Puncture the vein with the Splittocan® needle percutaneously or through a skin incision.
 - b) Flush the catheter with sodium chloride (NaCl) 0.9 %. Remove the stilette from the cannula, and insert the catheter through the cannula to the desired position at the junction of the Superior Vena Cava and the right atrium.
 - c) Gently withdraw the cannula over the catheter.
 - d) Follow from point e) - l) above: "Percutaneous Technique".
- For the ST301M and ST305M, remove the mandrin from the catheter prior to cutting the Y connector off the catheter.

VI-1-3-5 Valved catheter

Do not place the catheter over the guide-wire, as this may cause damage to the valve. The catheter must be tunnelled from the puncture point to the port pocket; the tunnelling rod must not be placed through the valve.

VI-1-3-6 Pre-connected catheter and port

The catheter should be measured and trimmed to the correct length before insertion in the vein. The catheter must be tunnelled from the port pocket to the puncture point.

VI-2 ARTERIAL IMPLANTATION TECHNIQUES (FOR H.A.I.T.)

H.A.I.T. = Hepatic Arterial Infusion Therapy

VI-2-1 Surgical Implantation Technique

- a) Perform a pre-operative arteriogram.
- b) Any accessory arteries to adjacent organs must be embolised.
- c) The catheter is sutured to the chosen artery using the suture beads at the tip of the catheter.
- d) Perfusion of the organ should be verified after placement of the catheter.
- e) Tunnel the catheter to the port pocket, the port should lie approximately 1/2 - 1 cm below the skin surface. Flush the port with sodium chloride (NaCl) 0.9 %.

VI-2-2 Percutaneous implantation Technique

Percutaneous placement of an arterial catheter is generally via the femoral artery. Access may also be via the axillary, acromiothoracic or brachial artery, and can be made either percutaneously or by surgical cut-down, depending on physician's choice.

An introducer should not be used during arterial procedures.

System implantation and catheter placement may be accomplished through a variety of techniques and the physician must choose the appropriate method for the patient and the intended result. Wider application, addition or alteration will be based on the experience of the physician. In the following chapter we present two of the most commonly used techniques for hepatic artery infusion.

VI-2-2-1 Catheter tip in gastro-duodenal artery (see Fig. C)

Note: The port/catheter reference R305-A5ST has a side hole in the catheter 19 cm from the catheter tip. As the tip of the catheter will be embolised in place no flow occurs from the catheter tip. The taper of the catheter should be trimmed to the correct length (distal to the side hole) before placement in the gastro-duodenal artery so that the flow of medication is not obstructed and the side hole is in the proper hepatic artery.

A rigid 0.018" angiography guide-wire is required for placement of the tapered Anthron® catheter. As the catheter is embolised in place, removal of the catheter percutaneously may not be possible.

- a) An arteriogram should be performed to confirm the location of the tumour, the anatomy of the blood vessels and blood supply to adjacent organs.
- b) Any accessory arteries to adjacent organs which arise distally to the proposed location of the

side hole should be embolised.

c) Place an exchange guide-wire (minimum length 2 m, diameter 0.018") through the angiography catheter to the desired position in the gastro-duodenal artery.

d) Remove the angiography catheter and place the Anthron® catheter over the guide-wire (do not use an introducer).

e) Remove the guide-wire.

f) A check arteriogram should be performed to confirm catheter position.

g) The tip of the catheter should be fixed in place using embolisation coils, a gluing agent or other means of immobilising the catheter tip, and of embolising the gastro-duodenal artery. Access should be via a contra-lateral approach or through a micro-catheter which permits embolisation. Flow of medication is via the catheter side-hole.

h) In order to prevent displacement of the catheter tip during movement of the patient allow the catheter to sag or loop in the aorta.

Follow from point e) – l): "Percutaneous Technique".

The access port should be placed in the pelvic area ensuring that there is a solid base behind the port.

The port should not be placed over the abdominal cavity.

VI-2-2-2 Catheter tip in hepatic artery (see Fig. D)

a) An arteriogram should be performed to confirm the location of the tumour, the anatomy of the blood vessels and blood supply to adjacent organs.

b) Any accessory arteries to adjacent organs which arise distally to the proposed location of the catheter tip should be embolised.

c) Place a guide-wire through the angiography catheter to the desired position in the hepatic artery.

d) Remove the angiography catheter and place the Anthron® catheter over the guide-wire (do not use an introducer).

e) Remove the guide-wire.

f) A check arteriogram should be performed to confirm catheter position, function and adequate perfusion of the liver.

g) Ensure that a sufficient length of catheter is left in the aorta to prevent catheter tip displacement during patient movement.

Follow from point e) - l): "Percutaneous Technique".

The access port should be placed in the pelvic area ensuring that there is a solid base behind the port.

The port should not be placed over the abdominal cavity.

VI-3 PERITONEAL AND PLEURAL IMPLANTATION TECHNIQUES

The catheter and port may be implanted either surgically or percutaneously, and may be used for locoregional chemotherapy or for drainage of symptomatic malignant ascites.

VI-3-1 Peritoneal - Surgical implantation Technique

a) During laparotomy the tip of the catheter is placed in the area of interest depending on the indication.

b) The catheter may be sutured to the fascia if the cuffed catheter (Drainaport) is used, and then tunnelled to the port pocket at the base of the ribs.

c) Catheter function and position should be confirmed by flushing the port.

d) The catheter is cut to the desired length and attached to the access port, taking care that the skin incision does not lie over the septum and that the catheter is not kinked anywhere along its length.

e) Confirm catheter and port function before closing the skin.

f) It is possible to fix the access port to the fascia with sutures.

VI-3-2 Peritoneal - Percutaneous implantation Technique

(Implantation accessories are available separately).

- a) Identify the desired final position of the catheter tip using ultra-sound guidance.
- b) With the 18G needle access the peritoneal cavity sub-cutaneously, taking care not to puncture the intestine. A curved needle may be used to aid entry into the peritoneal space.
- c) Pass the 0.035" (0.89mm) guide-wire through the needle into the peritoneal cavity, remove the needle.
- d) Pass the 12F dilator over the guide-wire to dilate the catheter route, followed by the 16F peelable sheath. Place the catheter into the desired position in the peritoneal space. The use of a rigid catheter placed telescopically inside the peritoneal catheter and a super-stiff guidewire may facilitate its placement. Drain any fluid into a sterile bowl.
- e) Create a pocket for the access port at the base of the ribs or over the iliac crest. If necessary further tunnel the catheter to the port pocket. It is advisable to have the cuff, if used, in the sub-cutaneous tissue so that tissue in-growth helps maintain catheter position.
- f) Trim the catheter to the desired length and attach to the access port, take care that the skin incision does not lie over the septum and that the catheter is not kinked anywhere along its length.
- g) Confirm catheter and port function before closing the skin.
- h) It is possible to fix the access port to the fascia with sutures.

VI-3-3 Pleural - Percutaneous implantation technique

(Implantation accessories are available separately).

- a) Identify the desired puncture point using ultra-sound (US) guidance. If the procedure is performed without US guidance, the puncture point will be marked by aspiration of pleural fluid through the needle used for injecting local anaesthetic.
- b) With the 18G needle access the pleural space, taking care not to puncture the lung or liver.
- c) Pass the 0.035" (0.89mm) guide-wire through the needle into the pleural cavity, remove the needle.
- d) Pass the 12F dilator over the guide-wire to dilate the catheter route, followed by the 16F peelable sheath. Remove the guide-wire and dilator. Place the catheter into the desired position in the pleural space. Drain any fluid into a sterile bowl.
- e) Create a subcutaneous pocket for the access port at the base of the ribs. Tunnel the catheter to the port pocket. Place the cuff in the sub-cutaneous tissue so that tissue in-growth helps maintain catheter position.
- f) Trim the catheter to the desired length and attach to the access port. Take care that the skin incision does not lie over the septum and that the catheter is not kinked anywhere along its length.
- g) Confirm catheter and port function before closing the skin.
- h) It is possible to fix the access port to the fascia with sutures.

VI-4 SPINAL IMPLANTATION TECHNIQUES

- a) Select the level for introduction of the catheter (L2-L3 or lower), and infiltrate the chosen area with local anaesthetic.
- b) Select the position for the port pocket, this should be over a bony structure (the base of the ribs is usually chosen), and infiltrate local anaesthetic.
- c) Make a small incision at the chosen entry site and enter the epidural or intra-theal space using the Tuohy needle.
- d) Introduce the catheter into the Tuohy needle using the insertion aid on the catheter, and advance into the chosen position, the catheter is marked at 5 cm intervals to aid placement. It is suggested that a minimum of 5 cm should be placed in the epidural space. The PU catheter should be advanced until the mandrin reaches the end of the needle, the mandrin should then be withdrawn while advancing the catheter into the intra-theal or epidural space.

- e) Remove the Tuohy needle from the puncture site while holding the catheter in a fixed position.
 - f) To verify the patency of the catheter attach the catheter to the catheter coupling device* (see Fig. 17). Attempt aspiration and inject 3-5 mL of dilute local anaesthetic for epidural placement. For spinal placement confirm the catheter position by ensuring aspiration of CSF (Cerebro-spinal-fluid) is possible.
 - g) Prepare the port pocket at the chosen site, the port should lie approximately 1/2 - 1 cm below the skin surface.
 - h) Thread the anti-kink device onto the catheter and attach the anti-kink device to the fascia with a suture.
 - i) Attach the catheter to the tunnelling rod, the route of the tunnelling rod may be anaesthetised with local anaesthetic prior to tunnelling. The catheter is then tunnelled to the port pocket. It may be necessary to do this in two or more steps.
 - j) Cut the catheter, at right-angles, to the desired length allowing sufficient length for patient movement.
 - k) Slide the connection ring over the catheter; push the catheter onto the exit cannula of the chamber.
- The connection ring should be tightened by hand and the spanner used to give the final tightening to the securement device.
- l) The assembled port and catheter should again be checked for patency.
 - m) The chamber may then be sutured into the port pocket, taking care that the skin incision is not at the injection site.
- * Attention: confirm that the catheter is pushed to the end of the catheter coupling device before closing to ensure correct function (see Fig. 17).

VII - USE OF THE PORT AND CATHETER

VII-1 VENOUS PORTS

Always verify that the port and catheter are functional by aspirating 2 mL of blood into a syringe and injecting 5 mL of sodium chloride (NaCl) 0.9% before attempting to start an infusion.

If aspiration of blood is not possible, attempt to slowly inject 2 mL of sodium chloride (NaCl) 0.9% into the port. If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected. In this case, device integrity should be verified using X-ray and contrast media studies.

During treatment, the infusion should be stopped immediately if any pain or swelling is noted or if blood return is absent.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VII-2 ARTERIAL PORTS

Before attempting to start an infusion, verify that the port and catheter are functional. Using a syringe containing sodium chloride (NaCl) 0.9 %, attempt to inject into the port. Current practice is not to verify function by aspiration of blood into the syringe.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII - MAINTENANCE OF THE PORT AND CATHETER

General recommendations:

- If heparinised sodium chloride (NaCl) 0.9% is used the system should be rinsed first with sodium chloride (NaCl) 0.9% alone, as some drugs may react with heparin and result in blockage of the port/catheter due to the for-mation of precipitates.
- Always rinse the system between administration of different drugs. Special care should be taken with drugs which have a higher risk of precipitation, with anti-coagulation agents, after blood sampling or after transfusion to reduce the risk of catheter occlusion.

- Failure to maintain the system may result in occlusion of the catheter.

VIII-1 RINSING AND HEPARINISATION

VIII-1-1 Venous ports

Rinsing of the access port is essential.

Rinse the port with 10 mL of sodium chloride (NaCl) 0.9% after each use and every 4 weeks when no treatment is being given.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII-1-2 Arterial ports

Arterial ports and catheters should be rinsed with 10 mL of sodium chloride (NaCl) 0.9% followed by heparinised sodium chloride (NaCl) 0.9% after each treatment and regularly, according to clinical needs (eg every 4 weeks) when no treatment is being given.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII-1-3 Epidural or Intra-thecal ports

Sodium chloride (NaCl) 0.9%, 0,5 mL may be used to flush the system after use. Heparin should not be used in an epidural or spinal port.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII-1-4 Peritoneal ports

Peritoneal ports may be rinsed with 20 mL sodium chloride (NaCl) 0.9% after each treatment.

If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

VIII-1-5 Pleural ports

Pleural ports should be rinsed with 20 mL of heparinised sodium chloride (NaCl) 0.9%, initially this should occur weekly. The interval between rinsing may be extended according to clinical needs. If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.

IX – DURATION OF IMPLANTATION AND REMOVAL OF THE PORT

IX-1 DURATION OF IMPLANTATION OF THE SYSTEM

The port and catheter system should be removed at the end of treatment.

The duration of the system is mostly dependant on the catheter. Fracture is the most common complication of ageing. Risk of fracture may increase over time depending on the material of the catheter, the pathway, and the entry point of the catheter into the vessel.

IX-2 REMOVAL OF THE SYSTEM

When removing the system, care must be taken not to fracture the catheter. If the catheter is sutured into the vessel, the sutures should first be removed. Control the catheter while removing the system from the port pocket. Be vigilant if there is excessive resistance to the removal of the catheter. Catheters may become encapsulated and attached to the vein wall. Should this occur and it is not possible to remove the catheter without risking catheter fracture, or if the catheter is fractured, the advice of an Interventional Radiologist, Surgeon or other Physician with endo-luminal experience should be sought.

Arterial catheters - If the catheter has been immobilised in an artery, using either coils or gluing techniques, removal of the catheter may not be possible. Attempts at removal may result in fracture

of the catheter within the artery and/or vessel injury.

IX-3 SPECIAL ATTENTION: BABYPORT®, BABYPORT® S, BABYPORT® PC AND SMALL PORTS

Following implantation of a port in a child particular attention should be given to the position of the catheter tip. Over time, as the child grows, the position of the catheter tip will change and move higher in the Superior

Vena Cava. As for any catheter, a high tip position in the SVC predisposes the catheter to fibrin sheath formation and other mechanical difficulties.

The position of the catheter tip should be checked radiographically a minimum of every 12 months or more frequently if the child is growing quickly. The system should be changed when the tip reaches the level of T4.

X - HIGH PRESSURE INJECTION

SPECIAL ATTENTION

- Always verify that the Celsite® reference implanted is included in the table, only these references may be used for high flow rate/high pressure injection.
- Always verify that the port and catheter are functional by aspirating 2 mL of blood into a syringe and injecting 5 mL of sodium chloride (NaCl) 0.9% into the port/catheter before attempting to start an infusion of medication.
- Do not exceed the recommended pressure (325 psi- 22.4 bars) and flow rate as Access Port system failure may occur.
- Contrast media should be warmed to 37°C (98.6°F) before use. According to drug manufacturers recommendations. Failure to follow this recommendation will result in up to 50% lower flow rates and/or Access Port or injection system failure.
- Use only Winged Surecan® (without Y-site) or Angled Surecan® needles.
Do not use needles which may not withstand high pressure.
- Use only 20G or 22G needles for Babyport® and Brachial ports. Use of 19G needles may lead to leakage of contrast media.
- Ensure that the needle is correctly placed in the port, securely taped to the skin and covered with an adhesive dressing before commencing high pressure injection.
- Catheter lengths more than 20 cm will result in reduction in flow rates.
- Depending on the technical characteristics of the injector system, the target flow rate might not be attained.
- The Access Port system should be flushed with 10 mL of sodium chloride (NaCl) 0.9% before, and after, using the port for CECT, followed by usual rinsing procedures.

		Recommended maximum flow rates (mL/s)*						Recommended maximum pressure setting (CT function)
		Contrast media at 37°C (98,6°F)						
		Viscosity 5,8 mPa.s (cP)**			Viscosity 11,4 mPa.s (cP)**			
		Winged / Angled Surecan® needle			Winged / Angled Surecan® needle			
		22G	20G	19G	22G	20G	19G	
	Double port	Babyport® (1,5 x 0,8 mm)	2	4		1	3	
		Brachial (1,6 x 1,1 mm)	2	4		1	3	
		Babyport® S (2,0 x 1,2 mm)	2	4		2	4	
Double port		ST405L (3,2 x 1,3mm)	2	5	6	2	4	6
Small Size Port		STL205P - STR205P (2,1 x 1,4 mm)	2	4	6	2	3	5
		ST305P - ST205P - BT305P (2,1 x 1,4 mm)	2	4	6	2	3	4
		ST305C (1,6 x 1,1mm)	2	4	5	1	3	4
		T/ST305 - T/ST205 - ST505 - SNT305F ST315 - ST215 - T/ST205F ECG - SNT215F (2,2 x 1,1 mm)	2	4	5	2	3	4
		STL205F - STR205F (2,2 x 1 mm)	2	4	5	2	3	4
		ST305L - ST505L - T/ST205ECG ST205L - ST315L (2,8 x 1,1 mm)	2	4	5	2	3	5
		ST305H - ST505H - ST205H (2,8 x 1,6 mm)	2	5	7	2	4	6
	Double port		ST401L (3,2 x 1,3 mm)	2	5	7	2	4
Standard Size Port		ST301C - ST201C - ST301OTW (1,6 x 1,1 mm)	2	5	6	2	4	5
		T/ST301F - ST311F - SNT301F T/ST201F - T/ST501F - T/ST201F ECG - SNT201F (2,2 x 1,1 mm)	2	5	6	2	4	6
		T/ST301P - ST201P - BT301P (2,1 x 1,4 mm)	2	5	6	2	4	6
		T/ST301 - ST311 T/ST201 - T/ST501 - T/ST201ECG (2,8 x 1,1 mm)	2	5	6	2	4	6
		STL201L - STR201L (2,8 x 1,1 mm)	2	5	6	2	4	6
		ST201H - T/ST301H - ST311H (2,8 x 1,6 mm)	2	5	7	2	5	7
		STL201H - STR201H (2,8 x 1,6 mm)	2	5	7	2	5	7
		ST301G - ST201G (3,2 x 1,6 mm)	2	5	8	2	5	7
* With a catheter of 20 cm								
** Or any contrast media with a similar viscosity at 37°C								

		Recommended maximum flow rates (mL/s)*						Recommended maximum pressure setting (CT function)
		Contrast media at 37°C (98,6°F)						
		Viscosity 5,8 mPa.s (cP)**			Viscosity 11,4 mPa.s (cP)**			
		Surecan® Safety II			Surecan® Safety II			
		22G	20G	19G	22G	20G	19G	
	Babyport® (1,5 x 0,8 mm)	2	4		1	3		325psi
	Brachial (1,6 x 1,1 mm)	2	5		1	4		
	Babyport® S (2,0 x 1,2 mm)	2	5		1	4		
Double port	ST405L (3,2 x 1,3mm)	2	5	8	1	4	6	
Small Size Port	STL205P - STR205P (2,1 x 1,4 mm)	2	5	8	1	4	5	
	ST305P - ST205P - BT305P (2,1 x 1,4 mm)	2	5	8	1	4	5	
	ST305C (1,6 x 1,1mm)	2	5	7	1	3	5	
	T/ST305 - T/ST205 - ST505 - SNT305F ST315 - ST215 - T/ST205F ECG - SNT215F (2,2 x 1,1 mm)	2	5	8	1	4	6	
	STL205F - STR205F (2,2 x 1 mm)	2	5	8	1	4	6	
	ST305L - ST505L - T/ST205ECG ST205L - ST315L (2,8 x 1,1 mm)	2	5	8	1	3	6	
	ST305H - ST505H - ST205H (2,8 x 1,6 mm)	2	6	9	1	4	6	
Double port	ST401L (3,2 x 1,3 mm)	2	5	8	1	4	6	
Standard Size Port	ST301C - ST201C - ST301OTW (1,6 x 1,1 mm)	2	5	6	1	3	5	
	T/ST301F - ST311F - SNT301F T/ST201F - T/ST501F - T/ST201F ECG - SNT201F (2,2 x 1,1 mm)	2	6	7	1	4	5	
	T/ST301P - ST201P - BT301P (2,1 x 1,4 mm)	2	5	7	1	4	6	
	T/ST301 - ST311 T/ST201 - T/ST501 - T/ST201ECG (2,8 x 1,1 mm)	2	6	7	1	4	6	
	STL201L - STR201L (2,8 x 1,1 mm)	2	6	7	1	4	6	
	ST201H - T/ST301H - ST311H (2,8 x 1,6 mm)	2	6	8	1	4	7	
	STL201H - STR201H (2,8 x 1,6 mm)	2	6	8	1	4	7	
	ST301G - ST201G (3,2 x 1,6 mm)	2	6	9	1	4	6	

* With a catheter of 20 cm

** Or any contrast media with a similar viscosity at 37°C

* With a catheter of 20 cm

** Or any contrast media with a similar viscosity at 37°C

en: Reference (A)	Port/ Catheter material (B)	Dead volume of port	Dead volume of catheter/ 10 cm (D)*	Implantatio n routes	Possible implantation techniques
Celsite®					
ST301V ST305V	Polysulphone /Silicone	0,50 mL 0,25 mL	0,18 mL	Venous	Seldinger
ST401L (double lumen) ST405L (double lumen)	Epoxy Resin /Silicone	0,50 mL x 2 0,25 mL x 2	0,13 mL	Venous	Seldinger
T201F, ST201F, T/ST201F ECG, SNT201F T301F, SNT301F ST301F, ST311F T305, ST305, ST315, SNT305F T205, ST205, ST215, SNT215F T/ST205F ECG	Epoxy Resin /Silicone Polysulphone /Silicone Polysulphone /Silicone Polysulphone /Silicone Epoxy Resin /Silicone Epoxy Resin /Silicone	0,50 mL 0,50 mL 0,50 mL 0,25 mL 0,25 mL 0,25 mL	0,08 mL	Venous	Surgical, Seldinger, ECG Surgical, Seldinger, Seldinger Surgical, Seldinger, Surgical, Seldinger, ECG
T202F (arterial) T302 (arterial)	Epoxy Resin /Silicone Polysulphone /Silicone	0,50 mL	0,08 mL	Arterial	Surgical
ST201G ST301G	Epoxy Resin /Silicone Polysulphone /Silicone	0,50 mL	0,2 mL	Venous	Seldinger
T301H, ST301H, ST311H ST305H ST201H ST205H	Polysulphone /PUR Polysulphone /PUR Epoxy Resin /PUR Epoxy Resin /PUR	0,50 mL 0,25 mL 0,50 mL 0,25 mL	0,2 mL	Venous	Surgical, Seldinger
T201, ST201 T301, ST301, ST311 T/ST201ECG T/ST205ECG ST305L, ST315L ST205L	Epoxy Resin /Silicone Polysulphone /Silicone Epoxy Resin /Silicone Epoxy Resin /Silicone Polysulphone /Silicone Epoxy Resin /Silicone	0,50 mL 0,50 mL 0,50 mL 0,25 mL 0,25 mL 0,25 mL	0,1 mL	Venous	Surgical, Seldinger Surgical, Seldinger ECG ECG Surgical, Seldinger Surgical, Seldinger
ST201C ST301C ST301OTW	Epoxy Resin /PUR Polysulphone /PUR Polysulphone /PUR	0,50 mL	0,1 mL	Venous	Seldinger, Braunule Over-the-Wire, Over-the-Wire
ST305C Brachial	Polysulphone /PUR Epoxy Resin /PUR	0,25 mL 0,15 mL	0,1 mL	Venous	Seldinger, Braunule Seldinger, Over-the-Wire
T301P, ST301P, BT301P ST201P ST305P, BT305P ST205P	Polysulphone /PUR Epoxy Resin /PUR Polysulphone /PUR Epoxy Resin /PUR	0,50 mL 0,50 mL 0,25 mL 0,25 mL	0,15 mL	Venous	Surgical, Seldinger Seldinger Seldinger, Surgical Surgical, Seldinger
T203J (peritoneal) DRAINAPORT T203J-1 (pleural/peritoneal)	Epoxy Resin /Silicone	0,50 mL	0,53 mL	Peritoneal Pleural/peritoneal	Seldinger, Surgical
Babyport® Babyport® S	Epoxy Resin /PUR Epoxy Resin /Silicone	0,15 mL	0,05 mL 0,11 mL	Venous	Seldinger
ST304-19G	Polysulphone PUR /Polyamide	0,33 mL	PUR 0,03 mL PA 0,03 mL	Intrathecal Epidural	Intrathecal Epidural
ST304-20G	Polysulphone PUR /Polyamide	0,33 mL	PUR 0,02 mL PA 0,02 mL	Intrathecal Epidural	Intrathecal Epidural
Celsite® Discreet					
STL201L, STR201L	Epoxy Resin /Silicone	0,50 mL	0,1 mL	Venous	Seldinger
STL205F, STR205F	Epoxy Resin /Silicone	0,25 mL	0,08 mL	Venous	Seldinger
STL201H, STR201H	Epoxy Resin /PUR	0,50 mL	0,2 mL	Venous	Seldinger
STL205P, STR205P	Epoxy Resin /PUR	0,25 mL	0,15 mL	Venous	Seldinger
Celsite® Concept					
T501F, ST501F ST505	Polysulphone + silicone /Silicone	0,50 mL 0,25 mL	0,08 mL	Venous	Surgical, Seldinger Seldinger
ST505H	Polysulphone + silicone /PUR	0,50 mL 0,25 mL	0,2 mL	Venous	Seldinger
T501, ST501 ST505L	Polysulphone)+ silicone /Silicone	0,50 mL 0,25 mL	0,1 mL	Venous	Surgical, Seldinger Seldinger
Celsite® IMPLANTOFIX					
04438604, 04438620, 04430263	Polysulphone /PUR	0,33 mL	0,1 mL	Venous	Over-the-Wire Braunule Surgical
04438817	Polysulphone /PUR arterial, open tip with two suture beads	0,33 mL	0,1 mL	Arterial	Surgical
04438647 04438663, 04433521	Polysulphone /PUR	0,15 mL	0,1 mL	Venous	Over-the-Wire Braunule Surgical
04438704 04438747	Polysulphone /Silicone	0,33 mL 0,15 mL	0,11 mL	Venous	Seldinger
Celsite® InterVentional with Anthron® catheter					
R305-A5ST	Polysulphone /PUR	0,25 mL	0,1 mL	Arterial	Percutaneous

*en: to determine the priming volume of the port / catheter, multiply the catheter length by the dead volume of the catheter and add the dead volume of the port.

XIII. Implantation images


Fig. A


Fig. B


Fig. C


Fig. D


Fig. E



Fig. F

en: Special Attention

CT Marking of Epoxy Access Ports

Contrast Enhanced Computerized Tomography (CECT) - refer to I - Indications

The following venous access ports are marked with "CT" to indicate compatibility with power injection. This marking is visible on X-ray. Please refer to for recommended maximum flow rates (mL/sec).

Paediatric	Small / Double	Standard / Double
Babyport®	STL205P, STR205P	ST201C
Brachial	T/ST205, STL205F, STR205F ST205L, ST205P, ST205H	T/ST201F, SNT201F
Babyport® S	ST215, SNT215	ST201F ECG
	T/ST205F ECG	T/ST201, STL201L, STR201L
	T/ST205ECG	ST201ECG
	ST405L	ST201H, STL201H, STR201H
		ST201G
		ST201P
		ST401L

en: Gravity flow rates established by gravity infusion of NaCl 0,9%, height 1 m. Catheter length 40 cm. According ISO 10555-1.

		Technical characteristics of the Silicone Catheters					
		Needles size	6,5F	8F	8,5F	10F	Double Lumen 10F
Gravity flow rate (NaCl 0,9%) mL/min Standard port	19G	26"	30	28	47	32	55
	20G	19"	21	20	29	24	31
	22G	10"	11	13	13	11	13
Gravity flow rate (NaCl 0,9%) mL/min Small port	19G	24	27	25		29	
	20G	18	19	18		22	
	22G	10	10	10		11	

* includes T302

		Technical characteristics of the PUR Catheters		
		Needles size	5F	6,5F
Gravity flow rate (NaCl 0,9%) mL/min Standard port	19G	24	34	45
	20G	19	24	28
	22G	10	11	12
Gravity flow rate (NaCl 0,9%) mL/min Small port	19G	22	30	37
	20G	18	22	25
	22G	10	11	12

		Technical characteristics of the Catheters for		
		Babyport®	Brachial	Babyport® S
Needles size		4,5F PUR	5F PUR	6F Silicone
Gravity flow rate (NaCl 0,9%) mL/min	19G	12	22	24
	20G	10	18	18
	22G	7	10	10

		Technical characteristics of the Catheters for	
		Celsite® Implantofix	
		Needles size	
Gravity flow rate (NaCl 0,9%) mL/min Standard & Small port	19G	22	26
	20G	18	21
	22G	10	11

		Technical characteristics of the PA/PUR 19G/20G Catheters	
	Needles size	T304-19	T304-20
Gravity flow rate (NaCl 0,9%) mL/min Standard port	19G / 20G / 22G	3	1

en: Description of the Device: Celsite® Access Ports are totally implantable access devices. The system consists of an access port, with a silicone septum, which is connected to a catheter. The catheter may be made of Silicone (Si), Polyurethane (PUR) or Polyamide (PA) depending on the reference. Access to the system is via percutaneous injection using a special non-coring needle. Both the port and the catheter are radio-opaque. For non-connected catheters, the marking is indicated from the distal end of the catheter. For pre-connected catheters, the marking is indicated from the proximal end of the catheter. The catheter is marked every cm starting at 5 cm.

MR Information

General Information

According to IEC Standard 60601-2-33(2008), the scanner must be operated in Normal Operating Mode (defined as the mode of operation of the MR system in which none of the outputs have a value that causes physiological stress to the patient):

The whole body averaged specific absorption rate (SAR) must be ≤ 2.0 W/kg

The head SAR must be < 3.2 W/kg

Device Information

MR Conditional

The Celsite® Access Ports were determined to be MR-Conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation:

24/28

F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

Non-clinical testing demonstrated that the Celsite® Access Ports are MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

Static magnetic field of 3-Tesla and 1.5-Tesla

Maximum spatial gradient magnetic field of 710 Gauss/cm or less

Maximum whole body averaged specific absorption rate (SAR) of 2.9 W/kg for 15 minutes of scanning

Note: The MR information is related to the access port only. The implantation accessories were not tested.

MRI-Related Heating





In non-clinical testing, the Celsite® Access Ports produced the maximum temperature rise during MRI performed for 15-min (i.e., per pulse sequence) in 3-Tesla (Excite, Software G3.0-052B, General Electric Healthcare, Milwaukee, WI) MR systems, as:

MRI condition	MR System Reported, Maximum Whole Body averaged SAR (W/Kg)	Calorimetry Value (W/Kg)	Highest Temperature Change	Time for MRI (per pulse sequence)
3-T / 128-MHZ	2.9	2.7	2.2°C	15 min

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Celsite® Access Ports. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Symbols

			
Do not re-use	Do not re-sterilize	Keep dry	MR Conditional
			
Keep away from sunlight	Do not use if package is damaged	DEHP-free	
			
Caution, consult accompanying documents	Latex free	Sterilized using ethylene oxide	Batch code
			
Date of manufacture	Use-by date	Catalogue number	Manufacturer
			
Standard size low profile port	Small size low profile port	Mini size low profile port	
			
Standard size low profile double lumen port		Small size low profile double lumen port	
			
Power injectable access port			

Annex 1

№	Name of additional item	Description
1.1	Name of the medical device	Name of the medical device: "Celsite® Access Ports for implantation with kit for their installation"
1.2	Name of manufacturer and manufacturing site	Manufacturer: B. BRAUN MEDICAL, 26 rue Armengaud, 92210 Saint-Cloud, France. Manufacturing site: B. BRAUN MEDICAL, 30 avenue des Temps Modernes, 86360 CHASSENEUIL-DU-POITOU, France.
1.3	Operating conditions of the medical device in the internal environment of the human body	Temperature, (C°): From 32 C° to 42 C°; Humidity, (%): 100 %.
1.4	Conditions of transportation and storage	Finished products are stored and transported at the temperature from -20 ° C to + 50 ° C and relative humidity up to 80%. Individual and group packages of the medical device must be stored in a dry place protected from direct sunlight.
1.5	Delivery set and packaging data	The delivery set is includes: 1. Consumer packaging: Cardboard corrugated box - 1pc. Cardboard corrugated box contains: 2. Instruction for use - 1 pc. 3. Patient Card Celsite® - 1 pc. 4. Patient Label Celsite® - 1 × 6 pcs. 5. <u>Sterile Individual packaging:</u> 5.1. Outer packaging: "Breathable" bag - 1 pc. 5.2. Inner packaging: Covering Tyvek + Tray - 1 pc. Sterile Inner packaging contains: 6. Medical device: "Celsite® Access Ports" - 1 pc.
1.6	The information about the symbol the conformity Declaration of conformity in the Russian Federation (if necessary)	
1.7	The information about the symbol of non-pyrogenic (if necessary)	
1.8	Shelf life for the medical device	Shelf life for the medical device is 5 years.
1.9	Warranty/Complaints	The company B. BRAUN MEDICAL, France guarantees that the medical device: "Celsite® Access Ports for implantation with kit for their installation" meets the requirements of the technical documentation, subject to conditions, transporting and storage. Company B. BRAUN MEDICAL, France is not responsible for the use of a medical device for other purposes. Claims for civil liability resulting from improper handling or use in combination with other medical devices or accessories will be rejected. On the circulation issues In the territory of the Russian Federation of the medical device "Celsite® Access Ports for implantation with kit for their installation", contact to the authorized organization (importer) in the Russian Federation: Limited Liability Company "B. Braun Medical", Russia, 191040, St. Petersburg, ul. Pushkinskaya, 10, Tel / fax.: 8 (812) 320-40-04, e-mail: office.spb.ru@bbraun.com.
1.10	Procedure for disposal or destruction of a medical device	Disposal of the medical device "Celsite® Access Ports for implantation with kit for their installation" after use is carried out in accordance with the requirements of for the disposal of medical waste according to class B and any other regulatory legal acts of the Russian Federation. Unused medical devices, unused packages must be disposed of in accordance with the procedure for the disposal of domestic solid waste, established by the legislation of the Russian Federation according to class A. If you have any questions regarding disposal of a medical device, you should contact the authorized organization (importer) in the Russian Federation.

nnex 10

1.11	Intended use	• Implanted Venous Access Ports are intended to be used for catheterization and repeated intra-venous administration of, for example, chemotherapy, antibiotic and anti-viral drugs, parenteral nutrition, blood sampling or transfusion. Certain Celsite® Access Ports may be used for Contrast Enhanced Computerised Tomography (CECT) using high pressure injection. • Implanted Arterial Access Ports are intended to be used for intra-arterial administration of chemotherapy. • Implanted Epidural or intra-theal ports are intended to be used for spinal administration of pain relieving drugs. • Implanted Peritoneal ports are intended to be used for intra-peritoneal administration of chemotherapy. • Implanted Peritoneal/Pleural ports are intended to be used for intra-peritoneal administration of chemotherapy, hydration, drainage of malignant ascites, or drainage of malignant pleural effusion.
1.12	Target category of patients / Target group of users	<i>Target category of patients:</i> Adults and children (taking into account anatomical features and individual needs). <i>Target group of users:</i> Only medical personnel trained in the installation, access to and removal of the port system, as well as the clinical aspects of the use of port systems, can install, access and remove port systems.
1.13	Model range of the medical device	Celsite® Anthron is not included in the current registration process (registration dossier) in the Russian Federation.

18/02/2021

Sylvie PLUMET

CoE Chasseneuil, Director

B. Braun Medical

Catherine BOISMENU

CoE Chasseneuil, Quality & Environmental

Management System and Regulatory Affairs, Manager

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Sylvie PLUMET
Catherine BOISMENU
(seen exclusively to certify the above signature)

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- 1.1. Name of the medical device
- 1.2. Name of manufacturer and manufacturing site
- 1.3. Operating conditions of the medical device in the internal environment of the human body
- 1.4. Conditions for transportation and storage
- 1.5. Delivery set and packaging data
- 1.6. The information about the symbol the conformity Declaration of conformity in the Russian Federation
- 1.7. The information about the symbol of non-pyrogenic
- 1.8. Shelf life for the medical device
- 1.9. Warranty/Complaints
- 1.10. Procedure for disposal or destruction of a medical device
- 1.11. Intended Use
- 1.12. Target category of patients / Target group of users

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for updating the documents of Registration Dossier № ФСЗ 2008/03191 issued on January, 24th, 2022 for the medical device Celsite® Access Ports with installation sets and submission of abovementioned document to the Federal Service on Surveillance in Healthcare (Roszdravnadzor).

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
Tél. +33(0)1 41 10 53 00
Fax +33(0)1 41 10 53 99

2


Instruction for use
Celsite® Safety
en: Venous Access Port Systems

Manufacturer: B. BRAUN MEDICAL, 26 rue Armengaud, 92210 Saint-Cloud, France.
Manufacturing site: B. BRAUN MEDICAL, 30 avenue des Temps Modernes, 86360 CHASSENEUIL-DU-POITOU, France.

2020

Contents






(en) System contents:

- (1) en - Housing
- (2) en - Septum
- (3) en - Exit cannula
- (4) en - Catheter
- (5) en - Titanium plate
- (6) en - Connection ring
- (7) en - CT mark

(en) Measurement Equivalence Table:

1 F = 0,33 mm	1 inch = 25,4 mm	18G/1.3 mm	19G/1.1 mm	20G/0.9 mm	22G/0.7 mm
---------------	------------------	------------	------------	------------	------------

 Access Ports	 Catheter	 Vein puncture needle
 Connection ring	 Rinsing hub	 J-guidewire

 Vein lifter	 Omnifix® syringe 10mL Luer slip	 ECG adaptor
 Surecan® Straight needle	 Surecan® Safety II needle	 Tunneling rod
 Peelable introducer	 Dilator 6F	

en: Refer to unit package label for actual content list

en: *18° for safety needle - **16° for non safety needle

I – DESCRIPTION OF THE DEVICE
II – INDICATIONS
III – CONTRAINDICATIONS
IV – POTENTIAL COMPLICATIONS
V – MR INFORMATION
V-1 General Information
V-2 MR-Related Heating
V-3 Artifact Information
VI – WARNINGS
VII – IMPLANTATION TECHNIQUES
VII-1 General considerations
VII-2 Percutaneous Technique
VII-3 Surgical Cut-down Technique
VII-4 Over-the-Wire (OTW) Technique
VIII – ECG TECHNIQUES
VIII-1 Catheter tip position as determined by ECG assisted catheter placement
VIII-2 Warnings
VIII-3 ECG Percutaneous Technique using J guidewire
VIII-4 ECG Percutaneous Technique using Saline (NaCl 0.9%)
VIII-5 ECG Surgical Cut-down Technique using J guidewire
VIII-6 ECG Surgical Cut-down Technique using Saline
IX – USE AND MAINTENANCE OF THE PORT AND CATHETER
IX-1 General recommendations
IX-2 Special attention for young patients
X – HIGH PRESSURE INJECTION
XI – SAFETY REFERENCES
XI-1 Puncture of the Vein
XI-2 Surecan® Safety II removal
XII – DURATION OF IMPLANTATION OF THE SYSTEM
XIII – REMOVAL OF THE SYSTEM
XIV – TECHNICAL CHARACTERISTIC OF THE ACCESS PORTS
XV – IMPLANTATION IMAGES
XVI – SYMBOLS GLOSSARY

I – DESCRIPTION OF THE DEVICE

The Celsite® Safety Venous Access Port System is a totally implantable device for vascular access. The system consists of an access port, with a silicone septum, which is connected to a catheter. The catheter may be made of Silicone (Si) or Polyurethane (PUR) depending on the reference. Access to the system is via percutaneous injection using a special non-coring needle (see § XI – SAFETY REFERENCES). Both the port and the catheter are radiopaque. The catheter is marked every cm starting at 5 cm.

II – INDICATIONS

• The Celsite® Safety Venous Access Port System is intended to be used for repeated intra-venous administration of, for example, chemotherapy, antibiotic and anti-viral drugs, total parenteral nutrition, blood sampling or transfusion. CT marked access ports may be used for contrast enhanced CT using a high pressure injection (see § X – HIGH PRESSURE INJECTION). This marking is visible on X-ray.

III – CONTRAINDICATIONS

The Celsite® Safety Venous Access Port System is contraindicated for patient therapy whenever:

- The presence of infection, bacteraemia, or septicaemia is known or suspected.
- Known allergy to the materials contained in the kit (main raw materials: Cyanoacrylate, Polycarbonate, Polyetheretherketone (PEEK), Polyethylene, Polypropylene, Polyurethane, Polytetrafluoroethylene (PTFE), Silicone, Stainless Steel, Tantalum, Titanium.
- Medications to be used in the access port are incompatible with any of the materials contained in the access port or catheter.
- The patient's anatomy does not allow the insertion of the catheter at the planned access site (body size, e.g. obesity).
- The patient has severe chronic obstructive pulmonary disease (chest placement only).
- The patient has had previous radiotherapy in the chosen area (condition of the skin, e.g. burns).
- Previous episodes of venous thrombosis.

IV – POTENTIAL COMPLICATIONS

Implantation or access procedures should not be attempted by personnel unfamiliar with the well documented potential complications (qualified healthcare practitioners). Complications may occur at any time during or after the procedure.

Possible complications include but are not limited to the following:

Air embolism	Drug extravasation	Perforation/Laceration/Erosion of vessels or Viscus
Artery or vein damage/Injury	Endocarditis	Pneumothorax
Bleeding	Fibrin sheath formation	Port/Catheter migration
Brachial plexus injury	Formation of deposits	Port/Catheter occlusion
Cardiac arrhythmia, tamponade,	Haematoma	Retraction or malposition of the catheter
Catheter disconnection, rupture or fragmentation	Haemothorax	Risks associated with anaesthesia/Surgery
Catheter damage or breakage due to	Haemorrhage	Thoracic duct injury
Pinch-Off syndrome (see below)	Implant rejection/Intolerance	Thromboembolism/Thrombophlebitis
Device rotation or erosion through the skin and/ or blood vessel	reaction to implanted device	Thrombosis
	Infection/Sepsis/Bacteremia	Twiddler Syndrome
	Nerve damage	

Pinch-Off Syndrome (see Fig. B):

Pinch-Off syndrome could occur if subclavian vein puncture is performed too medially. The following clinical signs may suggest catheter pinching:

- The patient needs to lift the arm to permit infusion.

- Intermittent malfunction of the catheter, such as difficulty with aspiration or infusion.
- Sub-clavicular pain or swelling during infusion.
- Palpitations or chest discomfort may indicate catheter fracture.

Remove any subclavian catheter which presents any of the signs of Pinch-Off syndrome.

A radiography with contrast injection should be done to detect catheter abnormalities at the costoclavicular level.

V – MR INFORMATION

V-1 General Information

MR Conditional

The Celsite® Safety Venous Access Port System was determined to be MR-Conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-13. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

Non-clinical testing demonstrated that the Celsite® Safety Venous Access Port System is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of a 1.5 Tesla and 3-Tesla.
- Maximum spatial gradient magnetic field of 10 000 Gauss/cm (extrapolated) or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode of operation for the MR system.

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V-2 MR-Related Heating

In non-clinical testing, the Celsite® Safety Venous Access Port System produced the maximum temperature rises during MRI performed for 15-min (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner), in 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems, as:

MRI condition	MR System Reported, Maximum Whole Body averaged SAR (W/Kg)	Calorimetry Value (W/Kg)	Highest Temperature Change	Time for MRI (per pulse sequence)
1.5-T / 64-MHZ	2.9	2.1	1.7°C	15 min
3-T / 128-MHZ	2.9	2.7	2.0°C	15 min

V-3 Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or within 5 mm of the position of the Celsite® Safety Venous Access Port System. Therefore, optimisation of MR imaging parameters to compensate for the presence of this device may be necessary.

VI - WARNINGS

Please note that the following warnings contain critical information for the patient safety. Failure to follow this Instruction For Use could result in serious injury to the patient:

- The product is sterilised using Ethylene Oxide and is sterile and non-pyrogenic in unopened, undamaged, individual packaging. For single use only.
- The device and its accessories are not reusable nor designed to be re-used. Any re-use would

definitely compromise the performance and safety of the device. Components may not be cleaned or re-sterilized.

- Destroy system components after access port implantation and explantation.
 - Store at room temperature. Avoid freezing and excessive heat. Do not remove from unit carton until ready for use.
 - Check that the sterile packaging is undamaged before opening. Confirm the date of sterile expiration on the labeling. Never implant a device which has exceeded its sterile expiration date. To maintain sterility, all normal aseptic room sterile procedures should be followed.
 - Always verify the integrity of the device before placement.
 - The Celsite® Safety Venous Access Port System should only be used by physicians who have experience with the technical and clinical requirements of access port implantation.
 - Check haemostasis before implantation: Partial thrombin time, Bleeding time, Platelet count.
 - Do not force insertion of the introducer, dilator or J-guidewire which could result in vessel injury (irritation, puncture or trauma).
 - When using a peelable introducer sheath, take care not to tear the vein when peeling the sheath.
 - Always place the catheter tip in the Superior Vena Cava positioned at the entrance to the right atrium taking into account possible catheter tip movement when the patient moves (see Fig. A). Failure to position the catheter tip at the entrance to the right atrium may lead to tip displacement and favours fibrin sheath formation. A radiograph should be taken to verify catheter position and exclude haemothorax or pneumothorax.
 - Subclavian placement should be considered carefully. It is recommended that the catheter be inserted outside the costo-clavicular space (see Fig. B) to avoid a Pinch-Off syndrome (see Pinch-Off Syndrome). Catheter ruptures have been observed with this route which could have serious consequences. It is necessary to pay particular attention to catheters used for long periods or for ambulatory patients. Catheter rupture can occur with silicone and polyurethane catheters.
 - It is recommended that the access port is fixed to the subcutaneous tissue to avoid any catheter problems or port rotation. Do not cut or occlude the catheter when using sutures.
 - If, for any reason, the access port/catheter assembly is disconnected during the implantation procedure it is recommended that at least 1 cm of the catheter should be cut off before re-connecting. Always trim any part of the catheter which has been connected to either the tunneling rod or the rinsing hub before connecting to the exit cannula. Ensure that the catheter is still in the right position after trimming.
 - Do not use standard hypodermic needles as they will damage the septum and may cause leakage of the system or blockage due to small silicone particles. Always use special bevelled needles (non coring needles), such as Angled Surecan®, Straight Surecan®, Winged Surecan®, Cytocan, Surecan® Safety II to penetrate the septum.
 - Never use incompatible drugs together or in sequence: the port should be rinsed of all residues with at least 10 ml sodium chloride (NaCl) 0.9% between each medication.
 - Always use a syringe equal to or greater than 10 ml volume when rinsing the port to avoid generating excessive pressures in case of occlusion of the system.
- In case of obstruction of the system never try to clear the blockage using fluid under high pressure which carries the risk of catheter fracture and migration. In case of blockage due to mineral or lipid deposits follow local validated protocols. The use of alcohol with PUR catheters is not recommended.
- In case of infection, if appropriate, treat with antibiotic drugs. If this fails, or is not appropriate, the catheter and port should be removed.

VI – IMPLANTATION TECHNIQUES

VII-1 General considerations:

- The use of ultra-sound guidance to identify a suitable vein, and during puncture of the vein is strongly recommended.
- The IFU's should be respected to minimise the risk of infection and thrombosis.
- Ensure that the catheter size is adapted to the diameter of the vein.
- Confirm catheter patency before implantation by flushing the catheter using the rinsing hub. Flush the port before connecting the catheter (see Fig. C). Use a non-coring needle to flush the preconnected systems.
- Remove the J-guidewire and the dilator together. Do not remove the J-guidewire through the dilator as this may result in the J-guidewire unravelling.
- To prevent shearing of the J-guidewire never withdraw the J-guidewire through the vein puncture needle.
- It is recommended that NaCl 0.9% (with or without heparin according to local protocols) be used to flush and lock the port and catheter during implantation.
- During implantation ensure that the catheter is not damaged by unguarded forceps, suture needles or other sharp instruments.
- Avoid kinking or sharp angles in the catheter pathway during implantation so that catheter patency is not affected.
- Ensure that the port pocket size is adapted to the size of the port.
- The catheter should be mounted on the exit cannula along its axis (see Fig. D) and not across the axis and should be completely mounted on the exit cannula before the connection ring is slid over the catheter.
- It is recommended that the catheter is tunnelled sub-cutaneously to the port.
- Immediately after port placement, the patient should be observed for potential complications.
- A post-operative radiograph should be taken to confirm catheter position and to exclude complications e.g. pneumothorax, haemothorax.

Notes:

Patient Identification Card

A Patient Identification Card is supplied with each Celsite® Safety Venous Access Port System. The completed card should be given to the patient. The patient should be instructed to carry it at all times.

Pre-connected catheter and port

The catheter must be tunneled from the port pocket to the puncture point in the vein.
The catheter should be measured and trimmed to the correct length before insertion into the vein.
Pre-connected ports cannot be used for ECG techniques.
Do not attempt to disconnect and/or re-connect pre-connected ports.

Valved catheter

Do not place the catheter over the J-guidewire as this may cause damage to the valve.
The tunneling rod must not be placed through the valve.

Catheter port/connection

Slide the connection ring over the catheter, firmly push the catheter onto the exit cannula ensuring the catheter covers the exit cannula up to the stop, slide the connection ring over the catheter and exit cannula. The connection ring should be in contact with the port (see Fig. D) to avoid the risk of premature disconnection.

Confirm catheter connection by gently pulling on the catheter.

VII-2 Percutaneous Technique

- a) Select the appropriate site for placement.
- b) Prepare the placement site following your current aseptic process.
- c) Insert the vein puncture needle into the chosen vein, verify the position by observing blood reflux (see Fig. E).
- d) Insert the J-guidewire into the vein through the needle, when the correct position is reached, withdraw the needle (see Fig. F).
- e) Thread the assembled dilator and peelable sheath over the J-guidewire using a twisting motion to pass through the skin planes (see Fig. G), a small incision can be made in the skin to facilitate passage of the sheath.
- f) Remove the dilator and J-guidewire together from the vein and insert the catheter through the peelable sheath (see Fig. H) to the desired position at the junction of the Superior Vena Cava and the right atrium (see Fig. A).
- g) Remove the peelable sheath from the vessel by pushing down on the 2 handles of the sheath, then peel the sheath from the catheter ensuring that the catheter is not dislodged from the vessel.
- h) Check catheter tip position using fluoroscopy.
- i) Prepare the port pocket at the chosen site, the port should lie approximately 1/2 - 1 cm below the skin surface. Remove the rinsing hub from the catheter and attach to the tunneling rod (see Fig. I).
- j) Tunnel the catheter from the puncture site to the port pocket (see Fig. J); ensure there is no kinking of the catheter.
- k) Cut the excess catheter length (at right-angles) prior to connecting to the port. Connect the catheter to the port (see Catheter/port connection)
- l) Insert the port into the port pocket, paying attention not to kink or loop the catheter.
- m) Confirm catheter patency (by ensuring that both aspiration and injection are possible), ensure that the skin incision is at least 1 cm from the injection site.
- n) It is recommended that the access port is fixed to the subcutaneous tissue with sutures.

VII-3 Surgical Cut-down Technique

- a) Select the appropriate site for placement.
- b) Prepare the placement site following your current aseptic process.
- c) The skin is incised and the tissue dissected (see Fig. K).
- d) The selected vessel is identified and tapes passed under it. The vessel is incised with a scalpel or with scissors (see Fig. L).
- e) The previously flushed catheter is inserted into the vessel, using either a small vessel dilator or vein lifter (see Fig. M) to the desired position at the junction of the Superior Vena Cava and the right atrium (see Fig. A).
- f) Follow from point h) - n) of the preceding paragraph: VII-2 Percutaneous Technique.

VII-4 Over-the Wire Technique

Follow steps a) and d): VII-2 Percutaneous Technique.

This technique is only possible with the references SNT601C OTW, SNT605C OTW and SST601C OTW.

- a) Thread the dilator over the J-guidewire using a twisting motion to pass through the skin planes.
- b) Remove the dilator (see Fig. N). Flush the catheter with sodium chloride (NaCl) 0.9 %, and pass the catheter over the J-guidewire to the desired position at the junction of the Superior Vena Cava and the right atrium (see Fig. O).
- c) Follow from point h) - n) above VII-2 Percutaneous Technique.

VIII - ECG Techniques

These techniques can only be used with the Celsite® Safety ECG references which have specially adapted accessories.

VIII-1 Catheter tip position as determined by ECG assisted catheter placement

It is recommended that the tip of the catheter is placed at the entrance to the right atrium (see VI – WARNINGS).

As the catheter approaches the right atrium the P wave (1) starts to elevate (see Fig. W). Maximal P wave height is reached and maintained when the catheter enters into the right atrium.

Note that the amplitude of the P wave will decrease if the catheter is advanced down into the inferior vena cava.

After identifying the area where the P wave begins to develop its maximal amplitude (which corresponds anatomically to the junction between Superior Vena Cava and the right atrium), advance the catheter a further 2 cm (2). This is the final position of the catheter tip.

As the procedure is performed with the patient in the supine position this catheter position allows for the 2-3 cm cranial movement of the catheter which occurs when the patient is upright.

VIII-2 Warnings

- To protect against the risks of current leakage, only use Class 1 ECG monitor CE marked, with floating input (Type CF) and connect additional potential equalization (earth) on the monitor.
- This technique can only be used with the ECG references which have specially adapted accessories.
- The Celsite® Safety Venous Access Port System with ECG technique should not be placed if any of the following is applicable to the patient:
 - The presence of a cardiac pacemaker.
 - Known cardiac abnormalities, or heart rhythm disturbances.
 - Access through the Inferior Vena Cava.
- If no ECG change is noted, catheter tip position has to be confirmed with fluoroscopy as catheter tip malposition is suspected. A post-operative radiograph should be taken to confirm catheter position.
- Protection against the risks of electrostatic discharge:
 - Use only in rooms protected from electrostatic discharge (conductive flooring).
 - Theatre shoes should be made of a recommended polymer material.
- Do not allow the J-guidewire to come into contact with external conductive parts.
- Before use, check that all connections are secure.
- Do not autoclave.
- Do not immerse in liquid.
- Repairs must be carried out only by the manufacturer.

VIII-3 ECG Percutaneous Technique using J guidewire

Preparation for ECG technique

- a) Place ECG pads on the patient's chest (see Fig. Q).

Note: A separate J-guidewire (0.035") is needed for percutaneous catheter placement.

b) Prepare the catheter by inserting the marked J-guidewire, which is contained in the kit, into the catheter.

Note: Do not attempt to insert the J-guidewire through the rinsing hub. The J-guidewire should be inserted directly into the catheter.

c) Slide the ECG adapter over the end of the marked J-guidewire ensuring that the J-guidewire has completely passed through the device (see Fig. R).

d) Turn the green handle through 90° to establish contact between the J-guidewire and the conducting pin in the ECG device (see Fig. S).

e) Follow steps a) – e): VII-2 Percutaneous Technique).

f) Remove the J-guidewire and dilator together from the sheath and insert the catheter (with the marked J-guidewire) into the peelable sheath. The J-guidewire should be pulled back into the catheter during placement into the vein.

g) Connect the plug end of the cable to the central socket of the Certodyn® universal adapter (see Fig. U) and switch to internal monitoring mode (see Fig. V).

h) When approximately 10-13 cm of the catheter has been placed into the vein the J-guidewire should be advanced into the catheter until the mark on the J-guidewire reaches the end of the catheter.

i) When the mark on the J-guidewire is at the level of catheter the J of the guidewire is in the desired position at the tip of the catheter.

j) See VIII-1 Catheter tip position as determined by ECG assisted catheter placement for catheter tip placement.

k) Remove the J-guidewire while maintaining the catheter in place, flush and clamp the catheter.

l) Remove the peelable sheath from the vessel by pushing down on the 2 handles of the sheath, then peel the sheath from the catheter ensuring that the catheter is not dislodged from the vessel.

m) Follow steps i) – n): VII-2 Percutaneous Technique.

VIII-4 ECG Percutaneous Technique using Saline (NaCl 0.9%)

Preparation for ECG technique

a) Place ECG pads on the patient's chest (see Fig. Q).

b) Introduce the rinsing hub into the proximal end of the catheter and rinse the catheter with saline.

c) Screw the ECG adapter on the rinsing hub.

d) The syringe filled with saline solution (NaCl) 0.9%, is screwed onto the ECG adapter, ensure there is no air in the system. (see Fig. T)

e) Follow steps a) – e): VII-2 Percutaneous Technique.

f) Remove the J-guidewire and dilator together from the sheath and insert the catheter (into the vein).

g) Connect the plug end of the cable to the central socket of the Certodyn® universal adapter (see Fig. U) and switch to internal monitoring mode (see Fig. V).

h) When approximately 10-13 cm of the catheter has been placed into the vein start to slowly inject the saline into the catheter through the adapter.

i) See VIII-1 Catheter tip position as determined by ECG assisted catheter placement for catheter tip placement.

j) Remove the peelable sheath from the vessel by pushing down on the 2 handles of the sheath, then peel the sheath from the catheter ensuring that the catheter is not dislodged from the vessel.

k) Follow steps i) – n): VII-2 Percutaneous Technique.

VIII-5 ECG Surgical Cut-down Technique using J guidewire

a) Preparation for ECG technique: follow steps a) – d) VIII-3 ECG Percutaneous Technique with J-guidewire.

b) Preparation for vein access: follow steps a) - d) VII-3 Surgical Cut-down Technique.

c) Insert the catheter (with the marked J-guidewire) into the vein. The J-guidewire should be pulled back into the catheter during placement.

d) Catheter placement: follow steps e) – i) VIII-3 ECG Percutaneous Technique with J-guidewire.

e) Follow steps k) – n): VII-2 Percutaneous Technique.

VIII-6 ECG Surgical Cut-down Technique using saline

- a) Preparation for ECG technique: follow steps a) – d) VIII-4 ECG Technique with Saline.
- b) Preparation for vein access: follow steps a) – d): VII-3 Surgical Cut-down Technique.
- c) The catheter is inserted into the vessel, using either a small vessel dilator or vein lifter (see Fig. M).
- d) Catheter placement: follow steps g) – i) VIII-4 ECG Technique with Saline.
- e) Follow steps k) – n): VII-2 Percutaneous Technique.

IX – USE AND MAINTENANCE OF THE PORT AND CATHETER

IX-1 General recommendations

- Always follow strict aseptic techniques when accessing or de-accessing the access port.
- Follow all contraindications and instructions for all infused products as specified by the manufacturer.
- If heparinised sodium chloride (NaCl) 0.9% is used, the system should be rinsed first with sodium chloride (NaCl) 0.9% alone, as some drugs may react with heparin and result in blockage of the port/
catheter due to the formation of precipitates. The use of alcohol with PUR catheters is not recommended.
- Always rinse the system between administration of different drugs. Special care should be taken with drugs which have a higher risk of precipitation, with anti-coagulation agents, after blood sampling or after transfusion to reduce the risk of catheter occlusion.
- At any time during implantation, use, maintenance or explantation of the device, a radiograph should be taken to confirm port/catheter position and to eliminate any doubts or complications if needed.
- Always verify that the port and catheter are functional by aspirating 2 mL of blood into a syringe and injecting 5 mL of sodium chloride (NaCl) 0.9% before attempting to start an infusion. If aspiration of blood is not possible, attempt to slowly inject 2 mL of sodium chloride (NaCl) 0.9% into the port. If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected. In this case, device integrity should be verified using X-ray and contrast media studies.
During treatment, the infusion should be stopped immediately if any pain or swelling is noted or if blood return is absent.
- Rinsing of the access port is essential and should be done with at least 10 mL of sodium chloride (NaCl) 0.9% after each use (injection and/or blood sampling) and every 4 weeks when no treatment is being given. If resistance to injection is noted, or if swelling occurs around the port or along the catheter, device malfunction should be suspected.
- A long term access port needle (Winged Surecan®, Cytocan®, or Surecan® Safety II) may be used from 24 hours up to 7 days, in the absence of infection, redness, swelling or pain, and in accordance with local, country or professional guidelines.
- The length and gauge of the needle should be chosen according to the morphology of the patient and the proposed treatment.
- To minimise blood loss and reduce the risk of air embolism during insertion of the guidewire into the needle or of the catheter into the peelable sheath, the patient should perform the Valsalva manoeuvre or a syringe or end cap may be used on the hub, or a thumb placed over the hub opening.
- The needle should be removed at the end of the treatment cycle.
- Do not tilt the needle once the septum is punctured.
- An occlusive dressing should cover the needle site at all times.
- Check the dressing regularly for cleanliness, change the dressing according to local protocols or if the dressing becomes soiled, wet or non-occlusive.
- It is recommended to use positive pressure techniques when removing the needle from the port (see Fig. P).

IX-2 Special attention for young patients

Particular attention should be given to the position of the catheter tip following implantation of a port in a child. Over time, as the child grows, the position of the catheter tip will change and move higher in the Superior Vena Cava. As for any catheter, a tip position high in the SVC predisposes the catheter to fibrin sheath formation and other mechanical difficulties.

The position of the catheter tip should be checked radiographically a minimum of every 12 months or more frequently if the child is growing quickly. The system should be changed when the tip reaches the level of the 4th thoracic vertebrae (T4).

X – HIGH PRESSURE INJECTION

CT marked access ports may be used for contrast enhanced CT using a high pressure injection. Ensure that the patient is not allergic to the contrast medium used.

Note: the references with valved catheters are not indicated for power injection.

Table for Standard and Small ports with silicone and PUR catheters (excluding Valved catheter)

Maximum recommended flow rates (mL/s) with a catheter of 20 cm			
Contrast media at 37°C (98,6°F)			
Viscosity 11.4 mPa.s (Cp)			
Winged/Angled Surecan® Surecan® Safety II Max. pressure 325 psi	19G x 20 mm	20G x 20 mm	22G x 20 mm
	5 mL/sec	5 mL/sec	2 mL/sec

Maximum recommended flow rates have been determined according to ISO 10555-6: Subcutaneous implanted ports with the following devices: bonded braided high pressure extension line (1200 psi)

(length 1219 mm; ID 1.8 mm x OD 3.5 mm) and Surecan® Safety II non coring needles.

Special attention:

- Always verify that the port and catheter are functional by aspirating 2 mL of blood into a syringe and injecting 5 mL of sodium chloride (NaCl) 0.9% into the port/catheter before attempting to start a high pressure injection.
- Do not exceed the recommended pressure (325 psi- 22.4 bars) and flow rate as Access Port system failure may occur.
- Contrast media should be warmed to 37°C (98.6°F) before use according to drug manufacturers recommendations. Failure to follow this recommendation will result in up to 50% lower flow rates and/or Access Port or injection system failure.
- Use only Winged Surecan® (without Y-site) , Angled Surecan® or Surecan® Safety II needles. Do not use needles which may not withstand high pressure.
- Ensure that the needle is correctly placed in the port, securely taped to the skin and covered with an adhesive dressing before commencing high pressure injection.
- Catheter lengths more than 20 cm will result in reduction of flow rates.
- Depending on the technical characteristics of the injector system, the target flow rate might not be attained.
- A pulsative flush should be performed through the access port system with at least 10 mL of sodium chloride (NaCl) 0.9% before and after using the port for CECT (Contrast Enhanced Computerised Tomography), followed by usual rinsing procedures.

XI – SAFETY REFERENCES

Certain references contain a Safety and echogenic vein puncture needle and a Surecan® Safety II non-coring needle which are designed to reduce the risk of accidental needle stick injury. Care must be taken to avoid needle sticks in accordance with local Infection Control standards and avoid the risk of exposure to contaminated blood.

Handle and dispose of all needles according to local protocols.

XI-1 Puncture of the Vein

Access the vein with the needle in the normal manner. After withdrawal of the Safety vein puncture needle activate the safety clip. To do this, grasp the clear clip housing, which houses the metal safety clip, with your free hand.

Advance the clear clip housing forward along the Safety vein puncture needle towards the needle bevel in one continuous motion. Advance the clip housing over the needle bevel activating the safety mechanism.

The clip will remain on the needle bevel. The clip housing will separate from the clip and needle. Dispose of the needle and the clear clip housing immediately into sharps container.

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XI-2 Surecan® Safety II removal

After completion of treatment, flush the port per local protocols. Stabilise the port by securely holding the base down. Firmly pull the wings up until you feel a firm stop and the needle is locked in the safety position. A green dot on the base of the needle confirms safety device deployment.

XII – DURATION OF IMPLANTATION OF THE SYSTEM

At the end of treatment and after medical approval, the port and catheter system should be removed. The duration of the system is mostly dependant on the catheter (e.g. the maintenance). Fracture

is the most common complication of ageing. Risk of fracture may increase over time depending on the material of the catheter, the pathway, and the entry point of the catheter into the vessel.

XIII – REMOVAL OF THE SYSTEM

When removing the system, care must be taken not to fracture the catheter.

Open the fibrous capsule around the port, cut any sutures which may have been placed to stabilise the port and remove the port.

Localise and stabilise the catheter connection ring with forceps and slowly remove the catheter from the vein. If the catheter is sutured into the vessel, the sutures should first be removed.

Be vigilant if there is excessive resistance to the removal of the catheter. Catheters may become encapsulated and attached to the vein wall. Should this occur and it is not possible to remove the catheter without risking catheter fracture, or if the catheter is fractured, the advice of an Interventional Radiologist, Surgeon or other Physician with experience in endovascular treatment should be sought.

Ensure that all the elements are removed from the patient (access port, connection ring and catheter).

Verify that the length of catheter removed corresponds to the length inserted (noted in patient's file).

If the explantation is due to suspected malfunction, the port and catheter should be returned to the local B. Braun representative for analysis.

XIV – TECHNICAL CHARACTERISTICS OF THE ACCESS PORTS

Technical characteristics of the Access Ports					
	Height (mm)	Length (mm)	Width (mm)	Weight(g)	Priming Volume (mL)
Standard size	13,7	33	26	8	0,5
Small size	11,3	29	22	5	0,3

Technical characteristics of the Catheters							
	Needle size	Silicone	Silicone	Silicone	PUR	PUR	PUR
		6,5F	8,5F	10F	5F	6,5F	8,5F
Priming volume of catheter/ 10 cm (mL)							
Gravity flow rate (NaCl 0,9%) mL/min Standard port	19G	26	34	48	26	37	48
	20G	21	25	31	21	27	31
	22G	10	11	12	10	12	12
Gravity flow rate (NaCl 0,9%) mL/min Small port	19G	26	34	48	26	37	48
	20G	21	25	31	21	27	31
	22G	10	11	12	10	12	12

Catheter diameter	Catheter length (mm)	Peelable sheath size (F)	Peelable sheath ID/length (mm)	Dilator ID/OD/length (mm)
5F	500	6F	2,11 / 140	O.98/2.08/178
6.5F	500	7F	2,44 / 140	O.98/2.41/178
8.5F	500	9F	3,10 / 140	0.98/3.07/178
10F	500	10F	3,43 / 140	0.98/3.40/178

XV – IMPLANTATION IMAGES


Fig A


Fig B


Fig C


Fig D


Fig E


Fig F


Fig G


Fig H


Fig I


Fig J


Fig K


Fig L


Fig M


Fig N


Fig O


Fig. P


Fig. Q


Fig. R


Fig. S


Fig. T


Fig. U


Fig. V


Fig. W

Annex 1

№	Name of additional item	Description
1.1	Name of the medical device	Name of the medical device: "Celsite® Access Ports for implantation with kit for their installation"
1.2	Name of manufacturer and manufacturing site	Manufacturer: B. BRAUN MEDICAL, 26 rue Armengaud, 92210 Saint-Cloud, France. Manufacturing site: B. BRAUN MEDICAL, 30 avenue des Temps Modernes, 86360 CHASSENEUIL-DU-POITOU, France.
1.3	Operating conditions of the medical device in the internal environment of the human body	Temperature, (C°): From 32 C° to 42 C°; Humidity, (%): 100 %.
1.4	Conditions of transportation and storage	Finished products are stored and transported at the temperature from -20 ° C to + 50 ° C and relative humidity up to 80%. Individual and group packages of the medical device must be stored in a dry place protected from direct sunlight.
1.5	Delivery set and packaging data	The delivery set is includes: 1. Consumer packaging: Cardboard corrugated box - 1pc. Cardboard corrugated box contains: 2. Instruction for use - 1 pc. 3. Patient Card Celsite® - 1 pc. 4. Patient Label Celsite® - 1 × 6 pcs. 5. <u>Sterile Individual packaging:</u> 5.1. Outer packaging: "Breathable" bag - 1 pc. 5.2. Inner packaging: Covering Tyvek + Tray - 1 pc. Sterile Inner packaging contains: 6. Medical device: "Celsite® Access Ports" - 1 pc.
1.6	The information about the symbol the conformity Declaration of conformity in the Russian Federation (if necessary)	
1.7	The information about the symbol of non-pyrogenic (if necessary)	
1.8	Shelf life for the medical device	Shelflife for the medical device is 5 years.
1.9	Warranty/Complaints	The company B. BRAUN MEDICAL, France guarantees that the medical device: "Celsite® Access Ports for implantation with kit for their installation" meets the requirements of the technical documentation, subject to conditions, transporting and storage. Company B. BRAUN MEDICAL, France is not responsible for the use of a medical device for other purposes. Claims for civil liability resulting from improper handling or use in combination with other medical devices or accessories will be rejected. On the circulation issues In the territory of the Russian Federation of the medical device "Celsite® Access Ports for implantation with kit for their installation", contact to the authorized organization (importer) in the Russian Federation: Limited Liability Company "B. Braun Medical", Russia, 191040, St. Petersburg, ul. Pushkinskaya, 10, Tel / fax.: 8 (812) 320-40-04, e-mail: office.spb.ru@bbraun.com.
1.10	Procedure for disposal or destruction of a medical device	Disposal of the medical device "Celsite® Access Ports for implantation with kit for their installation" after use is carried out in accordance with the requirements of for the disposal of medical waste according to class B and any other regulatory legal acts of the Russian Federation. Unused medical devices, unused packages must be disposed of in accordance with the procedure for the disposal of domestic solid waste, established by the legislation of the Russian Federation according to class A. If you have any questions regarding disposal of a medical device, you should contact the authorized organization (importer) in the Russian Federation.

1.11	Intended use	Implanted Venous Access Ports are intended to be used for catheterization and repeated intra-venous administration of, for example, chemotherapy, antibiotic and anti-viral drugs, parenteral nutrition, blood sampling or transfusion. Certain Celsite® Access Ports may be used for Contrast Enhanced Computerised Tomography (CECT) using high pressure injection.
1.12	Target category of patients / Target group of users	<i>Target category of patients:</i> Adults and children (taking into account anatomical features and individual needs). <i>Target group of users:</i> Only medical personnel trained in the installation, access to and removal of the port system, as well as the clinical aspects of the use of port systems, can install, access and remove port systems.

18/02/2024

Sylvie PLUMET

CoE Chasseneuil, Director

B. Braun Medical

Catherine BOISMENUCoE Chasseneuil, Quality & Environmental
Management System and Regulatory Affairs, Manager
B. Braun Medical

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