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APPROV
«I CERTIFY ACCURACY, CORRECTNESS AND
RELIABILITY OF THIS DOCUMENT TEXT»


Dr. Wolfgang G.K. Müller-Lierhe

Acting

Chief Executive Officer (CEO)

i.com medical GmbH, Germany

31.10.20

Sign, seal/stamp

i.com medical GmbH

Kaflerstrasse 15

D - 81241 München

INSTRUCTIONS ON USE

“Moisturizing ophthalmic solution
Ocutears®”

1. INDICATION

Moisturizing ophthalmic solution Ocutears® imitates the effect of the natural tears, protects and moistens the eye surface and lubricates it. This provides long-lasting relief for eye dryness as a result of:

- Keratoconjunctivitis sicca (Sjögren's syndrome)
- Injuries to the cornea, erosion, corneal surgery
- Wearing contact lenses
- Hormonal disturbances (menopause)
- Air conditioning, wind, cold, dust, smoke, dry heat, air pollution
- Work at a computer screen, long car journeys, long-haul flights
- Swimming in chlorinated water

2. MAIN CHARACTERISTICS

2.1. Principle of action

Moisturizing ophthalmic solution Ocutears® imitates the effect of the natural tears, protects and moistens the eye surface and lubricates it.

2.2. Composition

Moisturizing ophthalmic solution Ocutears® is a sterile, unpreserved, viscoelastic, protective shield for the cornea.

Moisturizing ophthalmic solution Ocutears® contains: Sodium hyaluronate 0.15%, sodium chloride, sodium hydrogen phosphate, sodium dihydrogen phosphate, manganese chloride tetrahydrate, HCl/NaOH to adjust pH close to physiological human tear pH range (6.8 - 7.4) purified water.

Each bottle of Moisturizing ophthalmic solution Ocutears® contains 10 ml of sterile fluid.

3. SPECIAL WARNINGS

- Do not use Moisturizing ophthalmic solution Ocutears® if you are allergic to any of the components of the product.
- Do not use if the tamper-evident ring of the cap or the bottle are damaged.
- Do not allow the bottle to touch the eye.
- Because of the risk of contamination, avoid bringing the tip of the bottle into contact with ocular surfaces or your hands.
- Do not use simultaneously with other medications for the eye. If you currently use other eye drops, you should wait at least five minutes before using Moisturizing ophthalmic solution Ocutears®.
- Do not apply to infected eyes.
- Write down the date of first opening the bottle on the bottle itself.

Note:



Moisturizing ophthalmic solution Ocutears has no influence on the ability to drive and use machines. As with any ocular treatment, if transient blurred vision occurs at instillation, the patient should wait until the vision clears before driving or using machinery.

Patients with chronic disease should consult doctor before starting to use moisturizing ophthalmic solution Ocutears®.

Pregnant and breast-feeding women should consult with doctor before starting to use moisturizing ophthalmic solution Ocutears®.

No clinical trials in children have been conducted with moisturizing ophthalmic solution Ocutears®.

4. INSTRUCTIONS FOR USE

When you are using the eye drops bottle for the first time, remove the tamper-evident ring of the cap by pulling the tab; Ocutears® can then be used for up to 6 months.	Remove the cap.	By pressing on the sides of the bottle squeeze one drop directly onto the eye or contact lens.	Replace the bottle's cap after use.
			

By gently massaging the eyelids with the eyes closed, the product will spread more quickly across the surface of the eye after application.

5. SIDE EFFECTS

As with all eye drops, it is possible that you have some mild discomfort in the eye for a short while after using the drops. IF YOU NOTICE ANY SIDE EFFECTS THAT ARE NOT LISTED IN THIS LEAFLET, PLEASE TELL YOUR DOCTOR.

6. STORAGE CONDITIONS

Store at 2°C to 25°C. Protect from frost.

Keep out of the reach and sight of children.

Do not use more than 6 months after the first opening.

7. SHELF LIFE

2 years.

Do not use after the expiry date printed on the box.

8. DISPOSAL

An used bottle, a bottle 6 months after its first opening, or expired bottle should be disposed via household waste.

9. GUARANTEE

This product was developed, manufactured, examined and packed according to all relevant requirements.

Manufacturer's guarantee is valid for unexpired, unopened, undamaged and properly stored product.

Developer and CE-mark holder in EU:

i.com medical GmbH

Kaflerstrasse 15

81241 München, Germany

Manufacturer:

i.com medical GmbH
Kaflerstrasse 15
81241 München, Germany
in factory:
Laboratoires Pharmaster
Z.I. de Krafft
67150 Erstein, France

Name and address of legal entity-holder of Registration Certificate:

Santen Oy
Niittyhaankatu 20
33720 Tampere, Finland

Additional information regarding the medication could be received and claims should be sent to the following address in Russia:

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