002202410160224

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SW2404601

241011 011 00002 7/12

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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Weifang Hengcai Digital Photo Materials Co., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明 网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Approved by
Weifang Hengcai Digital Photo
Materials Co., Ltd.

Director

 Xinjian Gong
signature stamp

"10" 10/10/2024

INSTRUCTION FOR USE
NO. CE/MDR-Z130106-01-02
FOR THE MEDICAL DEVICE

FILM FOR MEDICAL PRINTER RM-HL

2024

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

CONTENT	
1 GENERAL INFORMATION	4
2 DESIGNATION AND DESCRIPTION	4
2.1 Purpose and conditions of use	4
2.2 Potential customers	4
2.3 Indications	4
2.4 Contraindications and side effects	4
2.5 Precautions	5
2.6 Decoding of symbols	5
3 DESCRIPTION OF THE PRODUCT	6
4 TECHNICAL REQUIREMENTS	7
4.1 Main technical specifications	7
4.2 Materials used in the manufacture of the product	8
4.3 Set	8
4.4 Information on product labeling and packaging	9
5 OPERATING INSTRUCTIONS	11
6 INFORMATION ON MATERIALS OF HUMAN AND ANIMAL ORIGIN CONTAINED IN THE MEDICAL DEVICE	11
7 INFORMATION ON THE CONTENT OF MEDICINAL PRODUCTS IN THE MEDICAL DEVICE	11
8 ANIMAL STUDIES	11
9 STERILIZATION INFORMATION	11
10 DESIGNATED USE	12
10.1 General provisions	12
10.2 Operational limitations	12
10.3 Preparation before use	12
10.4 How to use the product	12
11 TRANSPORTATION AND STORAGE	13
12 REQUIREMENTS FOR MAINTENANCE AND REPAIR OF THE MEDICAL DEVICE	13
13 REQUIREMENTS FOR DISPOSAL AND ENVIRONMENTAL PROTECTION	13
14 MANUFACTURER'S WARRANTIES	13
15 ACCEPTANCE CERTIFICATE	14
16 PACKING CERTIFICATE	14
17 INFORMATION ABOUT CLAIMS	15
ANNEX A	15

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

This instruction for use applies to the Film for medical printer RM-HL (hereinafter - film, product) and is a document containing information on the purpose, principle of operation and characteristics of the product necessary for its proper use, transportation and storage, as well as information certifying the manufacturer's warranty.

It is recommended to read and follow the instructions for use before use.

Any questions regarding the use of the product should be addressed to the manufacturer or its authorized representative:

Manufacturer:

Weifang Hengcai Digital Photo Materials Co., Ltd.

Address: (Yingzi Village South 500 Meters) Yishan Road East, Dongcheng Subdistrict Office, Linqu County, Weifang City, Shandong, P. R. China

Tel.: +86 0536-2190926

E-mail: david@hcsminkjet.com

Production Site:

Weifang Hengcai Digital Photo Materials Co., Ltd.

Address: (Yingzi Village South 500 Meters) Yishan Road East, Dongcheng Subdistrict Office, Linqu County, Weifang City, Shandong, P. R. China

Authorized representative of the manufacturer in the Russian Federation:

AURA Limited Liability Company (AURA LLC)

Address: 355008, Stavropol Territory, Stavropol urban district, Stavropol, Golenyev Str. 47, apt. 40

Tel.: +7 (905) 338-89-63

E-mail: auramed99@mail.ru

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0
1 GENERAL INFORMATION 1.1 Introduction These instructions for use have been developed in accordance with the requirements of national standards and international normative documents (standards). 1.2 Name of the device FILM FOR MEDICAL PRINTER RM-HL I. Versions: 1. Film for medical printer RM-HL 8*10 inches (20×25 cm) (100 pcs/unit) 2. Film for medical printer RM-HL 10*12 inches (25×30 cm) (100 pcs/unit) 3. Film for medical printer RM-HL 11*14 inches (28×35 cm) (100 pcs/unit) 4. Film for medical printer RM-HL 14*17 inches (35×43 cm) (100 pcs/unit) II. Instructions for use - 1 pc. 2 DESIGNATION AND DESCRIPTION 2.1 Purpose and conditions of use The film is used in medical imaging systems in printers that utilize the principle of imaging by changing the temperature through dry printing. Scope of application: in conditions of therapeutic and preventive medical institutions in special rooms, designed including for X-ray studies, mobile medical complexes. The class of the product depending on the potential risk of use is 2a. No contact with the patient's body, X-ray personnel must wear gloves. 2.2 Potential customers The device is intended for use by qualified personnel of X-ray departments, who directly work with the equipment, as well as supervise its use, and who have familiarized themselves with the instructions for use. 2.3 Indications To use in medical imaging systems in printers that utilize the principle of imaging by changing temperature. The film is compatible with the HYC-4600 medical thermal printer, manufactured by Weifang Hengcai Digital Photo Materials Co., Ltd. 2.4 Contraindications and side effects No contraindications and side effects have been detected.		

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

2.5 Precautions

Attention!

- Handle the film with care, do not crumple, fold, deform or rub, do not touch the film with wet hands.
- The film must be stored in an upright position to avoid damage caused by pressure.
- Do not leave the film near heat or radiating sources.

2.6 Decoding of symbols

An explanation of the warning symbols used on the product packaging and in these instructions is given in Table 1.

Table 1 - Decoding of warning symbols

Symbol	Description	Location
	Date of manufacture	Consumer packaging
	Manufacturer name and address	Consumer packaging / transport packaging
	Batch code	Consumer packaging
	Non-sterile product	Consumer packaging
	Prohibition of repeated use	Consumer packaging
	Do not use if the packaging is damaged	Consumer packaging
	Model number (version)	Consumer packaging / transport packaging
	Medical device	Consumer packaging / transport packaging
	European conformity mark	Consumer packaging
	Authorized representative in the European Union	Consumer packaging
	Shelf life	Consumer packaging
	Refer to the instructions for use	Consumer packaging
	Protect from exposure to heat and radiating sources	Consumer packaging / transport packaging
	Keep away from moisture	Consumer packaging / transport packaging
	Fragile	Consumer packaging / transport packaging
	Temperature range	Transport packaging

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

3 DESCRIPTION OF THE PRODUCT


The film (Figure 1) is a thermographic medical dry-print film made of blue polyester-based film with an image and protection layer. The image layer is made with thermal dyes, polyvinyl alcohol and latex, and the protective layer includes zinc stearate, silicone oil and polyvinyl alcohol. Suitable for printing and outputting medical digital images such as CR, DR, CT, MRI. Can be used for medical printers with direct to lightroom imaging. The film is not photosensitive. Designed for high quality diagnostics and used as an alternative to paper copy. This product is a single-use consumable and is not sterile.

To determine the sensitive side of the film, a prong is provided on the film, which must be located at the bottom of the left corner of the short side of the film, when the sheet is placed with the long edges vertical to the observer.

The film is produced in four versions depending on the size. Each package contains 100 sheets of one size according to the version, placed in a tray.

a)

b)

Figure 1 - Appearance of film for medical printer,
where a) - appearance of the film sheet; b) - appearance of the film in the tray

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

4 TECHNICAL REQUIREMENTS

4.1 Main technical specifications

4.1.1 The film complies with the requirements of this document, the set of design documentation (DD) of the Manufacturer, approved in accordance with the established procedure.

4.1.2 The main technical specifications of the film are given in Table 2.

4.1.3 Dimensions and tolerances of the film are in accordance with ISO 4090:2001, ISO 18939:2013.

Table 2 - Main technical specifications of the film

No.	Name of characteristic	FILM FOR MEDICAL PRINTER RM-HL			
		8*10 inches	10*12 inches	11*14 inches	14*17 inches
1.	Nominal size, cm	20 × 25	25 × 30	28 × 35	35 × 43
	Overall dimensions of one sheet, mm	200,8 × 250,0	251,9 × 302,5	278.4 × 353.2	353,0 × 429,0
	Permissible deviation, mm	± 0.8	± 0.8	± 1.0	± 1.0
2.	Thickness of one sheet, μm	193	196	195	200
	Permissible deviation, μm	± 5	± 5	± 5	± 5
3.	Weight of one sheet, g	13	20	26	40
	Permissible deviation, g	± 3		± 5	
4.	Optical density of the veil	≤ 0,25			
5.	Maximum optical density	≥ 2.8			
6.	Sensitivity to light	Not sensitive			
7.	Film color	Blue transparent			

4.1.4 The shape of the film is rectangular with rounded corners designed within a quarter of a circle radius of 4 ± 1.2 mm.

4.1.5 The film has an image layer and a protective layer. The image layer reacts to high temperatures, ensuring the application of the image. The protective layer makes the film resistant to scratches and moisture.

4.1.6 There is a notch on the short edge of the film near the lower left corner, not more than 3.5 mm deep, designed to define the sensitive side of the film.

4.1.7 The film is placed in a tray, tray dimensions according to Table 3. Permissible deviation of tray dimensions ± 5 mm. On the rear surface of the tray there must be a radio frequency tag, carrying a code for each film format, for identification of the film by a compatible thermal printer.

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

Table 3 - Tray dimensions

No.	Film size	Tray size		
		Length, mm	Width, mm	Height, mm
1.	8*10 inches (20×25 cm)	401	378	15
2.	10*12 inches (25×30 cm)	401	378	15
3.	11*14 inches (28×35 cm)	401	378	15
4.	14*17 inches (35×43 cm)	475	378	15

4.1.8 The film surfaces are free from peeling, scratches, contamination, foreign inclusions and other defects that impair the appearance of the product.

4.1.9 Shelf life of the film is 2 years from the date of its manufacture under observance of storage conditions.

4.2 Materials used in the manufacture of the product

4.2.1 The materials used in the production of the film comply with the requirements of the current standards and normative and technical documentation for these materials, stipulated by the current legislation, and have accompanying documents on the quality of the supplier company.

4.2.2 Input control of raw materials and materials used in film production is carried out in accordance with the requirements of a certified and effectively used quality management system implemented at the Manufacturer.

4.2.3 There is no contact with the patient's body, the personnel must work in gloves.

4.3 Set

4.3.1 Supply set of film for medical printer RM-HL according to Table 4.

Table 4 - Set

No.	Name	Quantity
1.	Film for medical printer RM-HL in variants:	-
1.	Film for medical printer RM-HL 8*10 inches (20×25 cm) (100 pcs/unit)	1 package
2.	Film for medical printer RM-HL 10*12 inches (25×30 cm) (100 pcs/unit)	1 package
3.	Film for medical printer RM-HL 11*14 inches (28×35 cm) (100 pcs/unit)	1 package
4.	Film for medical printer RM-HL 14*17 inches (35×43 cm) (100 pcs/unit)	1 package
2.	Instructions for use	1 Pcs

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0
4.4 Information on product labeling and packaging 4.4.1 The film labeling is designed in accordance with ISO 15223-1:2021 and the requirements of Annex I of Regulation (EU) 2017/745. The labeling is authentic, clear and easy to read. 4.4.2 The labeling is fully applied on one side of the individual package in the form of a label or in the form of a paper label attached to it and contains: – name and address of the manufacturing company; – full name of the device; – date of manufacture (month, year); – batch No.; – identification number of the version (model); – film size in inches; – nominal film size in centimeters; – number of film sheets in a package; – shelf life; – conformity mark of the European quality certificate – “Medical device” symbol; – “Not sterile” symbol; – “Prohibition of repeated use” symbol; – “Do not use if the packaging is damaged” symbol; – “Refer to the instructions for use” symbol; – “Authorized representative in the European Union” symbol; – storage conditions recommended by the manufacturer: “Protect from heat and radiation sources” symbol, “Fragile. Warning” symbol, “Keep away from moisture” symbol, “Temperature limits” symbol. 4.4.3 An explanation of the warning symbols used on the product packaging is given in Table 1. 4.4.4 The details are indicated and filled in full or abbreviated in accordance with the rules of abbreviation in writing. 4.4.5 Transport labeling is applied on paper or cardboard labels, or directly on the container, labels are attached to the package with glue or other materials that ensure the safety of the cargo and labeling. 4.4.6 Labeling of transport packaging with application of manipulation signs: – “Fragile. Warning”; – “Keep away from moisture”; – “Temperature limits”; – “Up”; – “Protect from exposure to heat and radiating sources”; 4.4.7 The film in the amount of 100 sheets with an enclosed white sheet of the corresponding size is placed in a tray and then in an individual package made of aluminum foil under vacuum (foil bag). 4.4.8 The individual packages have the following characteristics:		

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

– for RM-HL 8*10 inch (20×25 cm) film: size (L×W×H) 500×470×24 mm ± 5%, gross weight 1.6 kg ± 5%;
 – for RM-HL 10*12 inch (25×30 cm) film: size (L×W×H) 500×470×24 mm ± 5%, gross weight 2.3 kg ± 5%;
 – for RM-HL 11*14 inch (28×35 cm) film: size (L×W×H) 500×470×24 mm ± 5%, gross weight 2.9 kg ± 5%;
 – for RM-HL 14*17 inch (35×43 cm) film: size (L×W×H) 575×470×24 mm ± 5%, gross weight 4.4 kg ± 5%;

4.4.9 Operation documentation is placed in a polyethylene film bag and placed with the product in a foam mold used for transportation.

4.4.10 The transport package contains 5 individual packages. Each individual package is placed in a foam mold used for transportation.

4.4.11 Gross weight of the transport package:

- for film RM-HL 8*10 inch (20×25 cm): 9.0 kg ± 5%;
- for film RM-HL 10*12 inches (25×30 cm): 12.5 kg ± 5%;
- for film RM-HL 11*14 inches (28×35 cm): 15.5 kg ± 5%;
- for film RM-HL 14*17 inches (35×43 cm): 23.5 kg ± 5%.

4.4.12 The packaging ensures protection of the contents from mechanical and climatic factors during transportation and storage, as well as the fullest use of the capacity of vehicles and convenience of unloading and loading operations.

10

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0
5 OPERATING INSTRUCTIONS 5.1 The film must be used in accordance with the instructions for use. 5.2 The product may be used by qualified personnel of X-ray departments, who work directly with the equipment and supervise its use, and who have familiarized themselves with the instructions for use. 5.3 Handle the film with care, do not crumple, fold, deform or rub, do not touch the film with wet hands. 6 INFORMATION ON MATERIALS OF HUMAN AND ANIMAL ORIGIN CONTAINED IN THE MEDICAL DEVICE The product does not contain any materials of animal or human origin. 7 INFORMATION ON THE CONTENT OF MEDICINAL PRODUCTS IN THE MEDICAL DEVICE The product does not contain medicinal products for medical use or pharmaceutical substances. 8 ANIMAL STUDIES No animal studies were carried out. 9 STERILIZATION INFORMATION Not applicable. The film is not sterile and must not be sterilized.		
11		

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

10 DESIGNATED USE

10.1 General provisions

Each product undergoes a special inspection by the manufacturer before sale. However, unforeseen damage not related to the production of the product may occur during transportation.

Please make a visual assessment of the packaging, the condition of the film tray to make sure that there are no various damages.

10.2 Operational limitations

⚠ Attention! Use the film in accordance with its intended use.

Operating ambient conditions:

Temperature: from 18 °C to 25 °C.

Relative humidity: max. 70 %.

Atmospheric pressure: 86.6 to 106.7 kPa (650 to 800 mmHg).

10.3 Preparation before use

Open the package, carefully remove the film tray. The tray has a recess, simulating a handle, at the bottom of the tray for easy handling.

⚠ Attention! Do not open the package with a knife (such as a paper cutter). This can damage the film and it cannot be used for printing.

10.4 How to use the product

1. Place the film tray and white sheet of paper in the appropriate tray of the medical printer.
2. Remove the white paper and close the medical printer tray.

⚠ Attention! Do not touch the printed side of the film.

3. Print the calibration page of the film.

⚠ Attention! When in doubt whether the film is correctly placed in the tray, the user should carefully remove the top sheet of film and make sure that the prong that serves as a reference for determining the sensitive side of the film is at the bottom of the left corner of the short side of the film. Be careful not to bend or crumple the film.

4. The film is now ready for use.

⚠ Attention! When the film is used with a medical printer, image problems may occur due to equipment malfunction. You should contact an engineer to solve this problem and follow the instructions for the medical printer.

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

11 TRANSPORTATION AND STORAGE

11.1 The film can be transported by all types of covered vehicles in accordance with the requirements of storage conditions. It is forbidden to cause mechanical damage.

11.2 Conditions of storage and transportation of the product: at ambient air temperature from 5°C to 25°C and relative humidity from 30% to 70%.

⚠ Attention! Excessive temperature and humidity during storage will cause film sticking, fog rising and quality defects.

11.3 The film must be stored in a cool and dry place without dust, protected from direct sunlight, acid and alkaline gases (such as hydrogen sulfide, ammonia, sulfur dioxide and formaldehyde).

11.4 The film must be stored in an upright position to avoid damage caused by pressure.

12 REQUIREMENTS FOR MAINTENANCE AND REPAIR OF THE MEDICAL DEVICE

The film is a consumable material and cannot be serviced or repaired.

13 REQUIREMENTS FOR DISPOSAL AND ENVIRONMENTAL PROTECTION

13.1 Disposal must be in accordance with accepted medical practice standards and relevant local, state and federal laws and regulations. In the Russian Federation, SanPiN 2.1.3684-21 is stipulated.

13.2 According to SanPiN 2.1.3684-21 the film belongs to class A - epidemiologically safe waste.

13.3 The film is subject to disposal in case of confirmation of facts and circumstances indicating failure to fulfill the intended purpose.

13.4 All packaging, including transport packaging, must be disposed of. Paper and polyethylene must be disposed of separately.

14 MANUFACTURER'S WARRANTIES

The Manufacturer guarantees conformity of the product to the requirements specified in this document, provided that the consumer complies with the conditions and rules of transportation, storage and operation during the entire shelf life of the product.

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

15 ACCEPTANCE CERTIFICATE

FILM FOR MEDICAL PRINTER RM-HL

Version _____

Batch No. _____

accepted in accordance with the mandatory requirements of state standards, current technical documentation and recognized as fit for operation.

Representative of the quality control department

(signature) (signature in words) (date)

16 PACKING CERTIFICATE

FILM FOR MEDICAL PRINTER RM-HL

Version _____

Batch No. _____

packed by Weifang Hengcai Digital Photo Materials Co., Ltd. in accordance with the requirements stipulated in the applicable technical documentation.

Packed by _____

(signature) (signature in words) (date)

The product after packing was accepted by

(signature) (signature in words) (date)

Instructions for use	Document No.	CE/MDR-Z130106-01-02
	Version	A/0

17 INFORMATION ABOUT CLAIMS

All submitted claims must be submitted to the manufacturer's authorized representative in the Russian Federation, registered by the manufacturer and contain information about the measures taken.

The period of claim consideration by the manufacturer is 1 (one) month from the date of receipt of the claim. The consumer is informed in writing about the measures taken.

To determine the reasons for non-compliance of the delivered product, its packaging, labeling and completeness with the requirements of the accompanying documentation, it is necessary to draw up an act, which should specify:

- name of the product;
- batch No.;
- detailed description of the discrepancy, etc.

The manufacturer accepts the claim if the fault of the recipient for the defect in the product has not been established.


ANNEX A

LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS

Scope of application	Designation of the national (or international) standard (document) in the country of manufacture	Standard (document) name
-	Regulation (EU) 2017/745	Directive 2017/745 of the European Parliament and of the Council of April 5, 2017 on medical devices
Label / Operating instructions	ISO 15223-1:2021	Medical devices. Symbols used to convey information provided by the manufacturer. Part 1. General requirements.
General standards	ISO 18939:2013	Materials for visualization. Digital hard copies for medical imaging. Measurement methods
	ISO 4090:2001	Photography. Medical X-ray cassettes/screens/films and films for printing images. Dimensions and specifications
Disposal of medical devices	SanPiN 2.1.3684-21	Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and potable water supply, atmospheric air, soils, living quarters, operation of industrial and public premises, organization and conduct of sanitary and anti-epidemic (preventive) measures

Approved by
Weifang Hengcai Digital Photo
Materials Co., Ltd.

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