
Celtra® Press Starter Kit dental material for the fabrication of glass-ceramic restorations

Directions for use

For dental use only

Version 2. November 2019.

JS

Meister, 02. Dec. 2019

 Dentsply
Sirona

DequDent GmbH
Rochanbacher Chaussee 4
D 63457 Hanau-Wolfgang

K. Neeb

K. Neeb/HR

Contents

1. Purpose	1	1.	Celtra® F
1.1 Technical Data	1	properties	labs. The
1.2 Indications for Use	1	Celtra® F	The Inge
1.3 Contraindications	1	layer, to	Sirona U
1.5 Compatible stains/glazes & veneering porcelain	1	restorati	Celtra® F
1.6 Compatible luting cements	1	3-unit br	splinted
1.4 Adverse reactions	1		
1.7 Information about disinfection, sterilization and frequency of application	1		
1.8 Conditions of storage and transportation	1		
1.9 Service conditions	1		
2. General Safety Notes	6	Intende	
2.1 Warnings	6	This pro	
2.2 Precautions	6	ceramic	
2.3 Adverse reactions	6	applicat	
3. Preparation	7	Deliver	
3.1 Preparation Guidelines	7	Celtra®	
3.1.1 Ingot types	8	-	
3.1.2 Inlays and Onlays	9	-	
3.1.3 Crowns and bridges	9	-	
3.1.4 Veneers	9	-	
4. Celtra® Press processing	10	-	
4.1 Model preparation (analogue)	10	-	
4.2 Wax-up	11	-	
4.2.1 Waxing & CAD Considerations	12	-	
4.2.2 Design	14	-	
4.2.3 Use wax weight to determine the size of ingot use	14	-	
4.2.4 Sprueing	15	-	
4.3 Investing procedures	16	-	
4.4 Pre-heating	16	-	
4.5 Pressing procedures	17	-	
4.6 Divesting procedures, Divesting procedures (Investment plunger)	18	-	
4.7 Removing the reaction layer	19	-	
4.8 Finishing	20	-	
4.9 Composite die procedures	21	-	
4.10 Build-up, stain and glaze	22	-	
4.10.1 PowerFire	22	-	
4.10.2 Porcelain application	22	-	
4.10.3 Dentin and/or Enamel cut-back/ Layered Technique	23	-	
4.10.4 Staining/full contour Technique	24	-	
5. General firing recommendation for Celtra® Ceram	25	Manufa	DeguD
6. Cementation	27	Manufa	
6.1 Preparation of the Celtra® restoration	27	1.	
6.2 Cementing	27	2.	
7. Technical and physical data	28		
8. Chemical composition	29		
9. Troubleshooting Guide	31		
10. Labelling	32		
11. Environmental Requirements	32		
12. Disposal	32		
13. Requirements to environment protection	32		
14. Shelf Life	33		
15. List of applied national and international regulatory documents/standards the medical product complies with	33		
16. Manufacturer Warranty	33		
17. Authorized representative in the Russian Federation	33		

To download this product information in other languages and to obtain additional firing recommendations, please visit www.celtra-dentsplysirona.com

1. Purpose

Celtra® Press is a high-strength zirconia-reinforced lithium silicate glass ceramic material that thanks to its translucent and opalescent properties can be used for the fabrication of highly aesthetic all-ceramic restorations by using the hot-pressing technique in dental labs. The material is supplied in the form of homogeneous ingots, which are produced in industrial conditions. The material in the Celtra® Press Starter Kit is available in two types of translucency: medium (MT) and low (LT).

The ingots are pressed in pressing furnaces, using ideally Celtra® Press investment material which does not form surface reaction layer, to obtain tooth-colored, highly aesthetic restorations. Subsequently, the pressed substructures can be stained with Dentsply Sirona Universal Stains & Glaze (for full-contour restorations) and/or veneered with Celtra® Ceram veneering porcelain (for cut-back restorations).

Celtra® Press is used to create occlusal veneers, thin veneers, veneers, inlays, onlays, crowns in the anterior and posterior regions, 3-unit bridges in the anterior region, 3-unit bridges in the premolar region up to the second premolar as the terminal abutment, crowns, splinted crowns or 3-unit bridges up to the second premolar placed on top of an implant abutment.

Intended Use

This product is designed for the fabrication of all-ceramic dental prostheses using pressure technique, e.g., partial ceramic and all ceramic crowns, inlays, onlays, veneers as well as bridges. It can also be used as a material for the fabrication of fully contoured applications.

Delivery form

Celtra® Press dental material for fabrication of glass-ceramic restorations in kit, including:

- Celtra® Press glass-ceramic material in ingots, medium translucency (MT), shade A2 or low translucency (MT/LT), shade BL2, 3g (5 pieces/package) – 1 package.
- Celtra® Press glass-ceramic material in ingots, medium translucency (MT), shade A3, 3g (5 pieces per package) – 1 package.
- Dentsply Prosthetics Die Material composite material, 4 g, shades F1 and/or F2 and/or F3 and/or F4 and/or F5 and/or F6 and/or F7 and/or F8 and/or F9 and/or F10 and/or F11 and/or F12 (2 pieces of different colors)
- Celtra® Ceram Enamel, shade E1, Extra-light, 15 g – 1 piece.
- Celtra® Ceram Enamel, shade E2, Light, 15 g – 1 piece.
- DS (Dentsply Sirona) Universal Overglaze - High Flu, 5r – 1 piece.
- Universal Stain & Glaze Liquid, 15 ml – 1 piece.
- Universal modeling liquid, 15 ml – 1 piece.
- Celtra® Undercoating Wax, 60 g – 1 piece.
- Celtra® Modeling Wax, 60 g – 1 piece.
- Celtra® Press investment, 100 g – 2 piece.
- Celtra® Press investment Liquid, 100 ml – 1 piece.
- Celtra® Press Muffle former, silicon based – 1 piece.
- Celtra® Press Muffle ring, 2 parts (base and lid) – 1 piece.
- Celtra® Press investment Plunger – 2 pieces.
- Firing Pad – 3 pieces/package.
- Illustrated Technical Guide¹ – 3 pieces.
- Order form for Celtra® Press investment² – 1 piece.

Manufacturer:

DeguDent GmbH, Rodenbacher Chaussee 4, 63457 Hanau-Wolfgang, Germany

Manufacturing sites:

1. DeguDent GmbH, Rodenbacher Chaussee 4, 63457 Hanau-Wolfgang, Germany
2. Dentsply International, Inc. Prosthetics Division – Also trading as: Dentsply Prosthetics, Dentsply Trubyte, Dentsply Austenal and Dentsply Ceramco and Dentsply Sirona, 570 West College Avenue, York, Pennsylvania 17401 USA

¹ Illustrated technical guide is not intended for users located in the Russian Federation.

² Order form for Celtra® Press investment is not intended for users located in the Russian Federation.

1.1 Technical Data	1.2 Indications for Use
Based on the classification of ISO 6872*, Celtra® Press is a Type II, Class 1-3, zirconia-reinforced lithium silicate (ZLS), dental ceramic substructure material with a CTE of $9.7 \times 10^{-6} \text{ K}^{-1}$ @ 25-500°C; and Celtra® Ceram is a Type I, Class I, leucite-reinforced porcelain for coverage of a ceramic substructure with a CTE of $9.0 \times 10^{-6} \text{ K}^{-1}$ @ 25-500°C. *2015-06	Celtra® Press is designed for use as a pressable material for the fabrication of dental restorations both for substructures and fully contoured applications for indications as follows: > Occlusal veneers > Thin veneers > Veneers > Inlays > Onlays > Crowns in the anterior and posterior region > 3-unit bridges in the anterior region > 3-unit bridges in the premolar region up to the second premolar as the terminal abutment > Crown, splinted crown or 3 unit bridge up to the second premolar placed on top of an implant abutment.
Risk class: 2a	
Scope of supply: The product is supplied together with the Instructions for Use.	

Celtra® Ceram porcelain is indicated for use as veneering porcelain for Celtra® Press substructure with the same indications. The restorations can be cemented with conventional adhesives.	1.3 The Celtra® > three > Sho tem frac > > > duc >
	1.4
	No for ma stru Dat to res

1.3 Contraindications The following are not indicated for Celtra® Press: > Bridges spanning more than three units > Temporary restorations. Should not be placed due to unstable temporary cements which could cause fracture of restorations. > Parafunction (bruxism) > Cantilever bridges > Patients with a substantially reduced residual dentition > Inlay bridges/Maryland bridges	1.5 Compatible stains/glazes & veneering porcelain Staining and glaze firing is recommended for completing Celtra® Press full-contour restorations. Dentsply Sirona Universal Stains & Glaze can also be used for restorations veneered with Celtra® Ceram porcelain. Celtra® Press restorations are compatible with Dentsply Sirona Universal Stains & Glaze for customizing and glazing (available separately). The use of other stain and glaze systems is not recommended (see Warnings). Use of other stain/glaze porcelain with Celtra® Press is at the discretion and sole responsibility of the dental technician. Likewise, Celtra® Press substructures are compatible with Celtra® Ceram veneering porcelain (available separately). The use of other veneering porcelain systems which may not have the correct coefficient of thermal expansion is not recommended (see Warnings). While initial results with some materials may appear acceptable, internal stress can compromise long term success. Use of other veneering porcelain with Celtra® Press is at the discretion and sole responsibility of the dental technician.	1.6 Compatible luting cements Full coverage crown restorations and bridges are compatible with self-adhesive resin cements, including all Dentsply Sirona self-adhesive resin cements (available separately). Inlays and onlays may be cemented with self-adhesive resin cements, however, adhesive bonding is recommended. Alternatively, full crowns and bridges can be cemented with glass- ionomer cement. Use of other cements or cement systems with Celtra® Press is at the discretion and sole responsibility of the dental technician.
1.4 Adverse reactions No adverse reactions have been reported for Celtra® Press. When working with these materials, make sure to comply with the Instructions for Use and the pertinent Safety Data Sheets (SDS). If the patient is allergic to any of the ingredients, Celtra® Press restorations should not be used.		1.7 Information about disinfection, sterilization and frequency of application Celtra® Press is a non-sterile product and is designed for single use, and re-sterilization or re-disinfection are not allowed.
		1.8 Conditions of storage and transportation - The storage and transportation temperature shall be 2-28 ° C³ - Do not use the product after the specified shelf life date expiry. The product shall be conveyed using all types of enclosed carriers in line with the rules of carriage applicable for the given transportation vehicle 1.9 Service conditions - The service temperature may vary from 18 to 25°C (the environment) - For use in dental laboratories of dental clinics - The materials may be also properly used inside the mouth cavity when there are biological fluids in it at the reference temperature of 32 to 42° C

2. General Safety Notes

Be aware of the following general safety notes and the special safety notes in other sections of these directions for use. This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

If properly processed and used, adverse effects of these medical products are highly unlikely. However, reactions of the immune system (such as allergies) or localized paresthesia (such as an irritating taste or irritation of the oral mucosa) cannot be completely ruled out as a matter of principle. In case of skin sensitization or rash, discontinue use and seek medical attention. Celtra® Press restorations are not suitable for patients with clinical symptoms of parafunctional habits or bruxism (see Contraindications). Do not inhale dust particles during grinding. Wear suitable protective mask. Avoid use of thick ceramic pins or third-party auxiliary firing pastes, correction porcelains or stains and glazes. Such use may result in cracking of Celtra® Press restorations. For patients with hypersensitivity to any of the ingredients, this medical device may not be used at all or only under the particular scrutiny of the dentist or physician in charge.

2.2 Precautions

This product is intended to be used only as specifically outlined in these "Directions for Use". Any use of this product that is inconsistent with the "Directions for Use" is at the discretion and is the sole responsibility of the practitioner. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients. Contact with saliva, blood and/or some astringent solutions during adhesive procedures may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended. Devices marked "single use" on the labeling is intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination. Celtra® Press restorations require adequate preparation reduction (see 3.1 Preparation Guidelines). Insufficient wall thickness may lead to premature failure. Celtra® Press restorations are designed to be polished and/or glazed before insertion. Direct insertion without polishing or glazing may lead to excessive wear on opposing dentition and compromise its esthetics.

2.3 Adverse reactions


No adverse reactions have been reported for Celtra® Press. When working with these materials, make sure to comply with the Instructions for Use and the pertinent Safety Data Sheet (SDS). If the patient is allergic to any of the ingredients, Celtra® Press restorations should not be used.

3. Preparation

3.1 Preparation Guidelines


Proper reduction of the hard tissue of the tooth during preparation is essential for maximizing the strength, shade and retention of the finished restoration. When preparing anterior or posterior teeth, the anatomical form has to be reduced as shown below.

Minimum wall thickness: The following diagram shows the specified minimum wall thickness for each indication. The minimum wall thickness must be retained after all manual adjustments have been made.

Important application precautions

All internal line angles of a preparation should be rounded. Sharp internally prepared angles have to be softened. Sharp line angles have to be rounded to prevent stresses in restoration. Celtra® Press substructure cusp tips and incisal edges must be designed to support the porcelain when veneered.

3 unit bridges including 2nd premolar:

Maximum pontic width:

anterior 11 mm

premolar 9 mm

Connectors for 3 unit bridges need a cross section of 16 mm².

> Principle Height > Width

3.1.2 In
con
rende
rat th
egree
rat al
cen
sally
Celtra
verec
Cely,
ons
sin o

3.1.1 Ingot types

The following information is for the reference purposes, the Celtra® Press Starter Kit only supplies medium and/or low trasluc ingots.

Possible shades and degrees of translucency of Celtra® Press ingots.

1. **HT ingots** are high translucency ingots and are available in values of I1, I2 and I3. HT ingot can be used for restorations in incisal region like inlays, onlays and veneers.

2. **MT ingots** are ingots with a medium translucency and are available in A-D shades. MT ingots can be used for crowns and bridges designed as a full contour restoration. MT ingots can also been used for restoration with a cut-back of the incisal region, pressing and veneering with Celtra® Ceram to complete the esthetic restoration.

3. **LT ingots** are shaded to Vita® A-D shades in dentin.

These are used for veneers or crowns and bridges with a cut-back of the incisal region, pressing and veneering with Celtra® Ceram to complete the esthetic restoration.

CELTRA® PRESS ingot portfolio

Restoration type	Translucency	Shade											Individualization technique
Incisal (inlay, onlay, veneer)	HT	I1 I2 I3											Glaze
Full Contour (posterior)	MT	BL2*	A1	A2	A3	B1	B3	C1	C3	D2	D3		Stain & Glaze
Cut-back (anterior)	LT	BL2*	A1	A2	A3	B1	B3	C1	C3	D2	D3		Build-up Stain & Glaze

• covered by one universal ingots MT/LT BL2

HT ingots are value-based and can be used for restorations in the incisal/enamel region.

Shade mapping:

I1: A1, B1, C1

I2: A2, A3, B2, C2, D2

I3: B3, C3, D3, A3.5, A4

3.1.2 Inlays and Onlays

A conventional inlay/onlay design is recommended. Do not prepare undercuts. Ensure that the cavity walls form an angle of 5 to 6 degrees with the long axis of the tooth. Ensure that all sharp edges and angles are rounded. In centric and dynamic occlusion, reduce incisally/occlusally by 1.5 to 2 mm. Celtra® Press inlays and onlays are ideally delivered by fully adhesive cementing. Alternatively, highly retentive inlay or onlay restorations may be cemented with self-adhesive resin cements.

3.1.3 Crowns and bridges

Ensure that there is an axial reduction of 1.0 to 1.5 mm with the walls form an angle of 5 - 6 degrees with the long axis of the tooth. In centric and dynamic occlusion, reduce incisally/occlusally by 1.5mm. The lingual shoulders must be extended at least 1.0 mm into the proximal contacts surfaces. It is recommended to use a shoulder preparation without a bevel: All angles must be rounded, and the preparation surfaces must be smooth.

Given the different masticatory forces, the maximum acceptable pontic width is different in the anterior and posterior region.

The pontic width is determined on the unprepared tooth.

> **In the anterior region (up to the canine), the pontic width should not exceed 11 mm.**

> **In the premolar region (canine up to the second premolar), the pontic width should not exceed 9 mm.**

Always observe the relation between width and height as well as the suitable dimensions (min. 16 mm²) when designing the connectors. Basically, the following applies: **Height > Width.**

Celtra® Press crowns and bridges can be delivered by either fully adhesive or self-adhesive cementing.

3.1.4 Veneers

The standard reduction is 0.6 mm for the labial surface and 0.4 mm in the gingival area (since the enamel is thinner in this region). Reduce the labiolingual incisal angle by 0.6 to 1.5 mm. The preparation margins should be located in enamel. A chamfer or rounded-shoulder preparation is recommended for all veneer margins. Proximal extensions must be located far enough proximally to conceal preparation margins from site and to avoid proximal gingival undercuts. Celtra® Press veneers are delivered by fully adhesive cementing. Self-adhesive cementing is not recommended for veneer restorations.

4. Celtra® Press processing

4.1 Model preparation (analogue)

Producing a die (using a sealer coating (e.g. Cergo® Sienna) for surface hardening). Apply a die spacer to within 1 mm of the preparation margin line, in two layers (for inlays and onlays: three layers). The procedure is analogous to that for natural dies. Prepare the master casts as for the fabrication of crowns and crowns made of precious dental alloys.

Procedure

- >Establish the preparation margin.
- >Use die hardener on the die.
- >Apply Cergo die spacer Sienna to the die.

- >Apply die spacer in one layer or in two layers as a placeholder for the cementing gap.
- >In the case of crowns, apply die spacer to within 1 mm of the preparation margin on the die.
- >For inlays, apply the die spacer over the entire prepared surface to just short of the preparation margin.

- >The die spacer was designed to simplify individual characterization of the working die optimizing the final restoration.

4.2 Wax-up

- > Use only designated organic pressing waxes that burn out without residue (Celtra® Cervical- and Undercoating Wax, 60 g; Celtra® Modeling Wax, 60 g).
- > Avoid sharp internal edges and (line) angles to reduce tension.
- > Follow the requirements for minimum connector size (cross-sections) and layering thicknesses.
- > Do not over contour preparation margins; create precise margins.
- > The model for the cut-back technique is first waxed up to full contour, as for the staining technique and then cut-back in the incisal third. Maintain the minimum wall thickness (a check with a silicone index is recommended).
- > Do not reduce the occlusal surfaces.

4.2.1 Waxing & CAD Considerations

> Enamel cut-back and/or Dentin/Enamel Layered Technique

This technique offers excellent esthetics. For cut-back, wax a crown/bridge or design a (CAD file) to full contour, and cutback of the enamel area. This crown/bridge is pressed using one of the dentin ingots (LT) and Celtra® Ceram is then used to complete the build-up. If too much dentin (core material) has been taken away, dentin porcelain can be applied, and then completed with enamel porcelain (Celtra® Ceram Enamel, shade E1, Extra-light; Celtra® Ceram Enamel, shade E2, Light). (If the core in a specific area is below recommended minimum thickness for the framework, dentin porcelain shall not be applied). For layered technique, the crown/bridge wax-up (or CAD file) are designed based on the preparation geometry and using the working model.

3 unit bridges including 2nd premolar:

Maximum pontic width:


Anterior 11 mm

Premolar 9 mm

Connectors for 3 unit bridges need a cross section of 16 mm².

> Principle Height > Width

The all-ceramic substrate must be 50 % larger than the porcelain being applied to it.

> Full Co
This tech
crowns to
This tech
applicati
Celtra®
Techno
Staining
techniq
Cut-ba
Layer
techniq
full-a

>Full Contour Technique / Staining technique

This technique is used primarily for the posterior application. The dentin/enamel ingots can be used for waxed or CAD/ CAM-designed crowns to full contour and stained and finalized with staining and glazing.

This technique is highly productive and due to the thicker core material, the result is a restoration with higher strength. For this application wax the unit or design a (CAD file) to full contour. MT ingots would be used for creating crowns and bridges.

Celtra® Press - minimum framework wall thickness/veneer thickness (mm)

Technology Area		Inlays	Onlays/ Tabletop	Veneers	Anterior crowns	Posterior crowns	Anterior bridges	Posterior bridges
							Connector 16 mm ²	cross-section
Staining technique	Framework wall thickness (fully contoured)	full-arch* 1,0 > Isthmus width 1,5	1,5	0,6	1,2	1,5	1,2	1,5
		incisal/occlusal	1,5	0,6	1,5	1,5	1,5	1,5
Cut-back	Framework wall thickness	full-arch* -	-	0,6	1,2	1,5	1,2	1,5
		incisal/occlusal	-	0,4	0,8	0,8	0,8	0,8
	Veneering (thickness)	-	-	0,4	0,7	0,7	0,7	0,7
Layering technique	Framework wall thickness	full-arch*/ incisal/occlusal -	-	-	0,8	0,8	0,8	0,8
	Veneering (thickness)	-	-	-	0,4 - 0,7	0,7	0,7	0,7

* „full-arch" refers to the region of the tooth equator.

4.2.2 Design

4.2.2.1 Wax-up

- 1 Apply a thin coat of a conventional die lubricant for a wax-up.
- 2 Wax the crown using a conventional inlay wax. Note:

Note

Use only non-contaminated inlay wax. After burning out, some inlay waxes may leave ash or carbon residue. This will result in dark areas on the pressed crown.

3 Recreate all necessary anatomical features in wax and completely seal the margins.

4 For full coverage restorations, be sure the wax or (CAD file) has a minimum thickness of 0.8 mm. Laminate veneers should have a minimum thickness of 0.4 mm.

3 unit bridges including 2nd premolar:

Maximum pontic width:
anterior 11 mm
premolar 9 mm

Connectors for 3 unit bridges need a cross section of 16 mm².
> Principle Height > Width

4.2.2.2 CAD design

- 1 Establish model with CAD/CAM gypsum or alternatively prepare the model with CAD spray.
- 2 When constructing digital restoration to specifications for the cement gap through the menu system while scanning result.
- 4 For full coverage restorations, be sure the wax or (CAD file) has a minimum thickness of 0.8 mm. Laminate veneers should have a minimum thickness of 0.4 mm

4.2.3 Use wax weight to determine size of ingot use

- 1 Weigh the sprue base.
- 2 Use the following table for ingot selection according table 1.

4.2.4 Sprue
Weigh the
nm [6 ga
Table 1
Max Wei
0.0 - 0.70
0.71 - 1.7
Always at
> S
> A
> -
> -
> F
For 3 uni

4.2.4 Sprueing

Weigh the wax-up including the sprues and base as recommended above. Sprues should be (0 2.5mm - 3.0mm single units or 0 4.0 mm [6 gauge] for bridges).

Table 1

Wax Weight	# Ingots	Ingot Size
0.0 - 0.70 g	1	3 g
0.71 - 1.70 g	1	6 g

Always attach sprues in the direction of flow and at the thickest point of the object in order to ensure complete pressing.

- > Select the 100 g or 200 g Celtra® Press muffle former (use 200 g rings only for bridges).
- > Attach only a single sprue (0 4.0 mm) to bridges (increase pressing temperature by 10 °C to 870 °C)
- > The distance from the Celtra® Press Muffle ring should be at least 10 mm.
- > The wax-up plus sprue should not exceed 16 mm in length.
- > No short "blind" pressing sprue is required for single objects when using Celtra® Press.

For 3 unit bridges one sprue (0 4.0 mm) is sufficient. Pressing temperature must be increased by 10 °C to 870 °C.

4.3 Investing procedures

Note

You must use a phosphate-bonded investment for these procedures. It is recommended to use the Celtra® Press investment. This investment was specially developed for Celtra® Press with the aim to prevent reaction layers from forming during pressing. As a result, it won't be necessary to acid-etch the restoration.

The following procedure is recommended:
Place the investment ring on the crucible for channel (while cold, use a vacuum plunger and mix the Celtra® Press Investment material (following the Instructions for Use) until all objects are completely covered while vibrating slightly to remove all bubbles. Then continue filling without vibration on and check the height and orthogonal position of the investment ring with the help of the ring gauge. Remove excess investment with the ring gauge.

4.4 Pre-heating

Keep the pre-heating oven clean to avoid residual investment in the channel (while cold, use a vacuum small broom to clean). After 20 minutes of setting, remove the investment ring and place it with the opening facing downwards in the heated furnace at 850 °C for 1 hour (100-g ring) or 45 minutes (100-g ring). **possible type investment plungers not require preheating. We recommend the use of Celtra® Press plungers.** We recommend placing the investment ring on a firing support or similar with the opening down during heating to avoid contact with the bottom plate of the chamber. (This ensures that the investment ring is heated uniformly and the wax can run out freely.)

Caution
CAUTION
During the burnout process, opening the furnace door will generate flames.

Note
Do not pre-heat Celtra® Press investment and/or plungers.

4.5 Pressing procedures

Avoid long waiting time between removing investment ring from the burnout furnace and placing into the pressing furnace, to prevent too much cooling of the ring.

Note

Pay attention to a precise pressing temperature of your press furnace by calibrating these at regular intervals. (Calibration kit Dentsply Sirona)

Up to max 0.7 g wax weight	1 press ingot, 3 g	100 g investment ring
Up to max 1.7 g wax weight	1 press ingot, 6 g	200 g investment ring

- 1 Remove a ring from the burnout furnace and place it on a heat-resistant surface with the sprue hole facing up.
- 2 Carefully place the correct size and shade of the ingot into the sprue hole.
- 3 Be aware that Celtra® Press ingots are exactly keyed to the A-D shade guide.
- 4 Do not stack ingots. Use 6 g ingots for a 200 g ring.
- 5 Place plunger it in the sprue hole on top of the ingot.
- 6 Place the ring in the center of the firing platform of the pressing furnace and start pressing procedure according to pressing DFU.
- 7 After the pressing cycle has been completed, remove the ring from the furnace and allow to bench cool. The ring may be fast cooled by placing it in front of a circulating fan.

General pressing recommendations

Low temp	Heating rate	Vacuum level	High temp	Hold time	Pressing time	Cool Time	Pressure
700 °C	40 °C/min	45 hPa	860 °C (100 g ring) 865 °C (200 g ring) (bridge, 200 g ring)	870 °C 30 min	3 min	0:00 min	2.7 or 4.5 bar depending on furnace design

Pressing recommendations Programat EP5000

Stand-by	Heating rate t	High temp T	Hold time H	Stop speed E
700 °C	40 °C/min	860 °C (100 g ring) 865 °C (200 g ring) 870 °C (bridge, 200 g ring)	30 min	250 pm/min

4.6 Divesting procedures, Divesting procedures (Investment plunger)

Note

Do not use pliers or chippers for divesting.

1 Indicate the position of the pressed objects using an investment plunger.

2 Make a deep cut into the investment compound, preferably using a diamond-covered and sintered large carbide disc or a carbide disc for metal castings.

3 Separate the part of the investment ring containing the investment plunger from the rest of the investment ring by turning in opposite directions.

4 Use a sand blaster (aluminum oxide 110 µm, max. 2.5 bars or 50 µm glass beads with 3 bars) to remove the investment. Make sure not to touch the objects (abrade the "internal cylinder" and remove the external ring of investment).

5 Once the objects have become visible, continue abrading across the area using reduced pressure (1.5 bars).

6 Use 50 µm glass beads with 1.5 bars for inlays for ensure the fit.

> *Using the investment press plunger (Celtra® Press plunger)*

> *To save time: Trim the investment including the plunger on the wet or dry trimmer.*

4.7 Rem
If using C

Bridge v

Note

If a rea
reduce

4.7 Removing the reaction layer

If using Celtra® Press investment material, no hydrofluoric acid gel or liquid is required.

Bridge with one 4-mm sprue.

Note

If a reaction layer has been formed, it's an indication of too high pressing temperature. Either calibrate your pressing furnace or reduce the pressing temperature.

4.8 Finishing

1 Conventional abrasives used for dental porcelains may be used for the Celtra® Press all-ceramic. Use a diamond disk to remove the sprues from the object.

2 Do not use carbide burs or coarse grit stones. This will cause a thermal crack and the ceramic coping/crown will have to be re-made.



3 Use a diamond bur, grinder for glass ceramic material or suitable aluminum oxide stones to re-contour the sprue attached area. Prevent overheating off the pressed unit.

4 Closely examine the inside of the restoration for bubbles or irregularities. These can be removed with a thin diamond bur or stone. A cracked restoration can't be repaired by firing porcelain on it and should be discarded.

5 Carefully seat the object on the die. Indicating sprays and materials may be used to aid in seating the object. Be sure indicating materials have been completely removed prior to placing objects in the porcelain furnace.

6 Use a porcelain finishing silicone wheel or point to adjust margins. Other abrasives may chip the pressed ceramic object and decrease the marginal integrity.

7 Grinding with diamond instruments (conventional dental technician hand pieces). Prevent overheating of the restoration. When using turbines, water cooling is recommended.

4.9 Con

The De
this mat

The den
If the d
shade.

1 Apply
2 Place
any voi
area.
3 Light
4 Remo
sonic cl

Due to
account
dental t
the oral

A1
F1

4.9 Composite die procedures

The **Dentsply Sirona Die Material** composite material is designed to mimic the actual shade of the patients prepared tooth. When this material is placed inside the pressed Celtra® crown, it will assist in accurate shade reproductions.

The dentist should take a prep shade of the tooth being restored for the laboratories reference with the Dentsply Sirona Prep Guide. If the dentist did not take the shade of the prepared tooth, the composite die material guide below may be used to verify the final shade. Select the proper composite die material from the appropriate table.

- 1 Apply the Dentsply Sirona Prosthetics Die Release to the inside of the ceramic restoration and allow it to dry.
- 2 Place a small amount of the Dentsply Sirona composite die material on the inside of the restoration. Pack the material to remove any voids. Immediately push a dowel pin into the uncured composite die material. Remove any excess composite from the margin area.
- 3 Lightcure the composite for 1-2 minutes using a hand-held light curing unit or the Triad 2000 curing unit from Dentsply Sirona.
- 4 Remove the composite die material from the restoration and carefully clean using a steam cleaner or in distilled water in an ultra-sonic cleaner for 10 minutes.

Due to the high translucency of Celtra® Press, the influence of the die shade on the shade of the restoration must be taken into account. The aesthetic result is also influenced by the color of the adhesive material. Using the supplied light-curing die material, the dental technician has the ability to map the shade information supplied by the dentist to a control die to replicate information about the oral situation in the shade reproduction. The aim is to simulate the shade of the prepared tooth (follow the working instructions).

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4.10 Build-up, stain and glaze

4.10.1 PowerFire

PowerFire is a firing program that is carried out before the first ceramic firing of the veneering porcelain. PowerFire increases flexural strength of the Celtra® Press restoration to >500 MPa. **After PowerFire, blasting must be avoided, since otherwise reduce the strength of the restoration.**

4.10.2 Porcelain application

Note

When firing a Celtra® Press All-Ceramic restoration it is important to **use only ceramic/ porcelain** type pins/pegs or to place object directly on the firing pad, to prevent issues with Celtra® Press during porcelain and glaze firings. When other than recommended type of pins/pegs were used, while initial results with some firing pins may appear acceptable, internal stress can compromise long term success. Do not fill the entire restoration with refractory putty. Doing so could cause cracking of restoration.

Framework- porcelain ratio standards*

	Veneer			Crowns & Bridges including 2nd premolar						
Overall thickness of restoration (mm)	0.8	1.0	1.1	1.2	1.5	1.7	2.0	2.2	2.4	2.8
Minimum framework thickness (mm)	0.4	0.5	0.6	0.8	0.8	0.9	1.1	1.2	1.3	1.5
Maximum layer porcelain thickness (mm)	0.4	0.5	0.5	0.4	0.7	0.8	0.9	1.0	1.1	1.3

*the strength of the veneering must not exceed 2.0 mm at any point.

- > **The minimum framework wall thickness will always be based on the total thickness of the restoration.**
- > **The thickness ratio of the framework wall to the ceramic layer must be at least 1:1 to ensure framework stability and esthetics.**


4.10.3

Layered

1. Use at 2 ext ful
2. Use fac wa cle per
3. Ap por tion Ce enl the
4. Alw tray Do
5. The the teri
6. To De De des app the tim


4.10.3 Dentin and/or Enamel cut-back/ Layered Technique

1. Use 50-micron sized aluminum oxide at 20 psi pressure and lightly blast the exterior surface of the crown. Be careful not to harm the margins.
2. Use a steam cleaner to clean the surfaces or put the restoration in distilled water and place in an ultrasonic cleaner for 10 minutes. Conclude by performing PowerFire.
3. Apply and fire the dentin and/or enamel porcelains to complete the restorations. All the components of the Celtra® Ceram system may be used to enhance the esthetics and contours of the restoration.
4. Always use the honeycomb sagger trays with the recommended firing pins. Do not use other types of firing pins.
5. The final shade may be checked using the Dentsply Sirona composite die material.
6. To glaze the restoration, mix the Dentsply Sirona Overglaze with Dentsply Sirona Stain & Glaze Liquid to desired consistency, if required, and apply onto the porcelain surface. Note the glazing cycle has 2:00 minutes hold time.

4.10.4 Staining/full contour Technique

1. Use 50-micron sized aluminum oxide at 20 psi pressure and lightly blast the exterior surface of the crown. Be careful not to harm the margins.
2. Use a steam cleaner to clean the surfaces or put the restoration in distilled water and place in an ultrasonic cleaner for 10 minutes.
3. Place the pressed crown on the prepared Dentsply Sirona Prosthetics die material.
4. Place a small amount of the enamel stain or glaze on the palette. Mix in the Dentsply Sirona Universal Stain and Glaze with the liquid to achieve a creamy viscosity and apply the mixture to the porcelain surface.
5. The Celtra® Universal Stains may be used for individual characterizations.
6. Remove the die material from the restoration. To ensure proper fit, remove excessive glaze from the interior of the crown as well as the internal margin areas.
7. Always use the honeycomb sagger trays with ceramic/porcelain type pins/pegs or place the object directly on the firing pad.
8. If necessary, corrections may be made to the crown using Celtra® Ceram Add-On/Correction porcelains.
9. Fire the crown using the recommended PowerFire and Glaze firing cycles.

Note

If a higher sheen is desirable, either raise the high firing temperature 10°C or use an additional 30 second hold time at the high temperature.

5. General firing recommendation for Celtra® Ceram

Cut-back: Layered technique

PowerFire - is a heating program to be done to the pressed units prior to layering ceramic. PowerFire increases flexural strength of the Celtra® Press restoration to its maximum >500 MPa.

PowerFire

Drying	Closing	Start temp	Pre heating	Heating rate	Final temp	Vacuum start	Vacuum stop	Vacuum time	Hold time*	Cooling**
min	min	°C	min	°C/min	°C			min	min	min
1:00	1:00	400	1:00	55	760	Off	Off	0:00	2:00	0:00

2. First Dentin/Incisal Firing

Drying	Closing	Start temp	Pre heating	Heating rate	Final temp	Vacuum start	Vacuum stop	Vacuum time	Hold time*	Cooling**
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	770	400	770	1:00	1:00	5:00

3. Second Dentin/Incisal Firing

Drying	Closing	Start temp	Pre heating	Heating rate	Final temp	Vacuum start	Vacuum stop	Vacuum time	Hold time*	Cooling**
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	760	400	760	1:00	1:00	5:00

Gaze firing

Drying	Closing	Start temp	Pre heating	Heating rate	Final temp	Vacuum start	Vacuum stop	Vacuum time	Hold time*	Cooling**
min	min	°C	min	°C/min	°C			min	min	min
2:00	2:00	400	2:00	55	750	Off	Off	0:00	2:00	0:00

Note

After PowerFire blasting must be avoided, otherwise the strength of the restoration will decline.

Add-on (correction) with 1st glaze firing

Drying	Closing	Start temp	Pre heating	Heating rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Time	Hold Time*	Cooling**
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	760	400	760	1:00	1:00	5:00

Add-on (correction) after glaze firing

Drying	Closing	Start temp	Pre heating	Heating rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Time	Hold Time*	Cooling**
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	750	400	750	1:00	1:00	5:00

Full contour: Staining technique**PowerFire & Glaze**

Drying	Closing	Pre heating	Start temp	Heating rate	Final temp	Vacuum Start	Vacuum Stop	Vacuum Time	Hold time*	Cooling
min	min	min	°C	°C/min	°C				min	min
2:00	2:00	2:00	400	55	1 st : 760 °C 2 nd : 750 °C	Off	Off	Off	2:00	5:00

> For this application, the glaze firing may be conducted together with the PowerFire.

> Firing recommendations for other ceramic furnaces, please visit: Celtra®-dentsplysirona.com

Add-on (correction) after glaze firing

Drying	Closing	Start temp	Pre heating	Heating rate	Final Temp	Vacuum Start	Vacuum Stop	Vacuum Time	Hold Time*	Cooling**
min	min	°C	min	°C/min	°C	°C	°C	min	min	min
2:00	2:00	400	2:00	55	750	400	750	1:00	1:00	5:00

* Hold time w/o vacuum

** In furnaces that cannot constitute a cooling phase, it is recommended to cool down to 600 °C until removal of the object.

Cementation

1.1 Preparation of the Celtra® restoration

Clean the restoration with a steam cleaner, in ultrasonic bath or with alcohol.

Apply 5% - 9% hydrofluoric acid etching gel (Available separately, see manufacturer's complete Directions for Use) to the interior of the restoration only and allow to soak for 30 seconds.

CAUTION: Follow manufacturer's precautions. Do not allow tissue or eyes to come into contact with the acid!

Remove the hydrofluoric acid as per the manufacturer's instructions.

Dry the restoration in an air stream. It is recommended to silanize the etched surfaces immediately.

At chairside, apply silane only to those surfaces required for adhesive cementing.

Allow to soak for 60 seconds. If the silane layer is no longer liquid, add more silane. Blow-dry in a powerful air stream.

Recommended material: Calibra® Silane Coupling Agent, available separately, see complete Directions for Use).

1.2 Cementing

Depending on the indication for Celtra® Press restorations a self-adhesive or full-adhesive cementation can be chosen. Compatible time-proven adhesive cementing materials are available as part of the Dentsply Sirona range of products. Alternatively, full crowns and bridges can also be fixed with glass ionomer cement. Cements are available separately.

	Self-adhesive	Fully adhesive	Glass-ionomer
Inlays	R	HR	-
Onlays	R	HR	-
Veneers	-	HR	-
Crowns	HR	HR	R
Bridges	R	HR	R

R = recommended HR = highly recommended


Directions for use in other languages: Celtra®-dentsplysirona.com

7. Technical and physical data

CELTRA® Press

Sizes

a	b
∅12.6mm	9.8-11.5mm

Storage conditions of kit components:

- Celtra® Press glass-ceramic material in ingots: bc 2-28 °C
- Celtra® Press investment: 5-30 °C
- Celtra® Press investment Liquid: 5-30 °C
- Dentsply Prosthetics Die Material composite material: 2-28 °C
- Universal modeling liquid: 2-28 °C
- DS (Dentsply Sirona) Universal Overglaze
- High Flu: 2-28 °C
- Universal Stain & Glaze Liquid: 2-28 °C
- Celtra® Ceram Enamel: 2-28 °C
- Celtra® Undercoating Wax: 2-28 °C
- Celtra® Modeling Wax: 2-28 °C
- Celtra® Press Muffle former: 2-28 °C
- Celtra® Press Muffle ring: 2-28 °C
- Celtra® Press investment Plunger: 2-28 °C
- Firing Pad: 2-28 °C

CELTRA® Press physical properties

Physical property	Specification limits
ISO 6872 type and classification	Type II, Class 4
Flexural strength	555 ±66.3 MPa
Chemical solubility	≤0.05%
Coefficient of thermal expansion	(9.8 ± 0.06) × 10 ⁻⁶ K ⁻¹
Radioactivity	<0.03 Bq/g
Resistance to coloration, the quantity of spots on samples after holding in a dyeing solution	0
Porosity	The number of pores with a diameter of 30 µm shall not exceed 16 and the number of pores with a diameter of 40 to 150 µm shall not be more than 6 on the surface with an area of 1mm ² . Pores with a diameter over 150 µm shall not present
Linear shrinkage on firing	≤16%

Ingot weight:

6 g ingots	6 g ± 5%
3 g ingots	3 g ± 5%

The consumer package sizes 7.7×1.7×1.7 cm ±5%
For physical properties of all components please see annex

3.	Comp
	silicon
	phosph
	alumin
	lithium
	potass
	zircon
	ceric o
	terbiu
	vanad
	iridium
	yttrium
	manga
	neody
	ferrou
	boron
	sodium
	Dentsply I
	Compo
	Silanate
	Methacr
	Silanate
	Urethane
	Ethyl 4-
	Campho
	Butylate
	2,5-dihy
	hyl ester
	Pigment
	Celtra® C
	Compo
	silica
	lithium
	potassiu
	aluminium
	sodium
	zirconiu
	cerium o
	boron tr
	calcium
	barium
	magnesi
	antimon
	terbium
	titanium
	tin oxide
	fluorine
	Pigment
	DS (Dent
	Compo
	silica
	aluminium
	sodium
	potassiu
	lithium
	barium
	calcium
	boron tr
	cerium
	terbium
	titanium

Chemical composition

Composition	Formula
silicon oxide	SiO ₂
phosphorus oxide	P ₂ O ₅
aluminum oxide	Al ₂ O ₃
lithium oxide	Li ₂ O
potassium oxide	K ₂ O
zirconium oxide	ZrO ₂
ceric oxide	CeO ₂
terbium oxide	Tb ₄ O ₇
vanadium oxide	V ₂ O ₅
iridium oxide	Er ₂ O ₃
yttrium oxide	Y ₂ O ₃
manganese oxide	MnO ₂
neodymium oxide	Nd ₂ O ₃
ferrous oxide	Fe ₂ O ₃
boron trioxide	B ₂ O ₃
sodium oxide	Na ₂ O

Dentsply Prosthetics Die Material composite material, 4 g

Component
Silanated Quartz
Methacrylate resin
Silanated Barium Aluminoborosilicate Glass
Urethane modified Bis-GMA dimethacrylate
Ethyl 4--dimethylamino Benzoate
Camphorquinone
Butylated Hydroxy Tlouene
2,5-dihydroxybenzene-1, 4- dicarboxylic acid diet- hyl ester
Pigments

Cera® Ceram Enamel

Composition	Chemical formula	E1, Extra-light	E2, Light
silica	SiO ₂	+	+
lithium oxide	Li ₂ O	+	+
potassium oxide	K ₂ O	+	+
aluminium oxide	Al ₂ O ₃	+	+
sodium oxide	Na ₂ O	+	+
zirconium oxide	ZrO ₂	+	+
cerium oxide	CeO ₂	+	+
boron trioxide	B ₂ O ₃	+	+
calcium oxide	CaO	+	+
barium oxide	BaO	+	+
magnesium oxide	MgO	+	+
antimony oxide	Sb ₂ O ₃	+	+
terbium oxide	Tb ₂ O ₃	+	+
titanium oxide	TiO ₂	+	+
tin oxide	SnO ₂	+	+
fluorine	F ⁻	+	+
Pigments	-	+	+

DS (Dentsply Sirona) Universal Overglaze - High Flu

Composition	Chemical formula
silica	SiO ₂
aluminium oxide	Al ₂ O ₃
sodium oxide	Na ₂ O
potassium oxide	K ₂ O
lithium oxide	Li ₂ O
barium oxide	BaO
calcium oxide	CaO
boron trioxide	B ₂ O ₃
cerium oxide	CeO ₂
terbium oxide	Tb ₂ O ₃
titanium oxide	TiO ₂

antimony oxide	Sb ₂ O ₃
----------------	--------------------------------

Universal Stain & Glaze Liquid

Composition	CAS#
Butylene glycol	107-88-0

Universal modeling liquid

Composition
Water, deionized
Distilled water
Other components

Celtra® Press investment

Composition
Quartz (SiO ₂)
Cristobalite
Aluminium oxide
Magnesium oxide
Ammonium dihydrogen orthophosphate
Boron Nitride

Celtra® Press investment liquid

Composition
Aqueous solution of alkali-stabilized colloidal silicic acid

Firing pad -SiO₂

Celtra® Undercoating Wax, Celtra® Modeling Wax:
Composition: Mineral wax, synthetic hydrocarbon wax, polymer resin

Celtra® Press Muffle former, silicon based
Material: Silicon

Celtra® Press Muffle ring, 2 parts (base and lid)
Material: POM black

Celtra® Press investment Plunger

Component
Quartz powder (SiO ₂)
Cristobalite
Magnesium oxide

9.

Issue
1. Crack
stacked
2. Fins o
Reason
der-liqui
3. Incom
vestmen
4. Failur

5. Resto
cially in t
apparen
Reason
6. Press
ing and/

7. Crack

8. Bullse

Troubleshooting Guide

Issue	Recommendation for resolution
Cracking of the investment rings due to misalignment of cracked ingots.	> Ingots must not be stacked. If insufficient material is a concern, use 6 g instead of 3 g ingot size.
Fins on restorations, while no ring cracking was noted. Reason: incomplete drying of the ring before pressing, powder-liquid ratio for investment is not perfect.	> Adhere to drying time of the last ring poured (20 min), if multiple rings are being placed in the burn-out furnace. > Check investment powder-liquid ratio to assure it's correct.
Incomplete investment ring burnout causing cracking of investment ring.	> Follow exactly the recommended burn out temperature as stated in the Investment DFUs.
Failure during divesting by using pliers or chipper.	> Avoid using chippers or other similar tools during divesting. > Sandblasting only!
Restoration appears too purple and/or transparent, especially in the case of bleach shades. Substantial reaction layer apparent after divesting Reason: The pressing temperature is too high	> Calibrating the pressing furnace > The pressing temperature was too high. > Repeat the pressing at the correct temperature (860 °C for 100 g ring, 865 °C for 200 g ring; 3-unit bridges: 870 °C for 200 g ring).
Pressed restorations are milky-white, opaque, dead-looking and/or pressing is incomplete.	> Check calibration of the furnace. > Pressing temperature is too low. > Either calibrate pressing furnace or validate that the pressing temp. is correct (860 °C for 100 g ring, 865 °C for 200 g ring; 3-unit bridges: 870 °C for 200 g ring).
Cracked restorations after glazing	> Do not completely fill restorations with Peg putty > Use only ceramic pins > Make sure restoration have recommended thickness according to DFU
Bullseye swelling at the sprue attachment site	> Use a sprue that is wider in one direction while reducing the height of the wax sprue. > See www.celtra-dentspolysirona.com

10. Labelling

Product (ingot) labelling

- Celtra® batch name
- translucency
- shade

- Manufacturer


- The product is in conformity with the requirements of Directive 93/42/EEC on medical devices

- For professional use only


- Manufacturer's trademark

- Consult instruction for use

- The medical device unique identifier

- Catalogue number

- Lot

- Use by

- Date of manufacture

- Temperature range (Celtra® Press investment Liquid)

- Hazard to human health (Celtra® Press investment)

11. Environmental Requirements

The device does not have an adverse impact on the environment during its life cycle.

12. Disposal

Disposal according to local authority regulations.

When disposed, the product itself is not a hazardous material. Please, observe the established national or local rules for waste disposal during disposal. No special requirements to waste disposal are applied.

13. Requirements to environment protection

Class of waste hazard Class B (hazardous wastes with respect to epidemiology)

Labelling layout for supplies in the Russian Federation

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14.

No stability as a glass assigned b

The life-time material in will have a Packaging transporta The shelf Storage co

15.

the

Standard

EN ISO

EN ISO

EEC Dir

cluding S

DIN EN

DIN EN

DIN EN

DIN EN

DIN EN

DIN EN

DIN EN

DIN EN

16

The mar

- Us

- No

The user

For prod

17

Dentsply

Phone: +

info@de

14. Shelf Life

A stability study was conducted for Celtra® Press. In principle, due to its inherent stability, Celtra® Press has an unlimited shelf life as a glass ceramic. However, for operational purposes, a 15-year shelf life will be applied to Celtra® Press. This shelf life was assigned based on experience with similar porcelain products.

The life-time of glasses and glass-ceramics is mainly limited by corrosion attack on the surface, mainly in humid atmosphere. Ceramic material in powder (with large specific surface) form has proven to be stable over 15 years. Therefore, the shelf-life of Celtra® Press will have a life-time exceeding 15 years. So the shelf life can be set to 15 years without further testing.

Packaging must be suitable to provide protection against mechanical deterioration or minimize the possibility of contamination during transportation and storage. Therefore, a detailed shelf-life testing of the packaging for the 15-year period is not required.

The shelf life is 15 years.

Storage conditions of kit components:

- Celtra® Press glass-ceramic material in ingots: 15 years
- Celtra® Press investment: 24 months
- Celtra® Press investment Liquid: 30 months
- Dentsply Prosthetics Die Material composite material: 5 years
- Universal modeling liquid: 5 years
- DS (Dentsply Sirona) Universal Overglaze - High Flu: 4 years
- Universal Stain & Glaze Liquid: 5 years
- Celtra® Ceram Enamel: 10 years
- Celtra® Undercoating Wax: unlimited shelf-life
- Celtra® Modeling Wax: unlimited shelf-life
- Celtra® Press Muffle former: unlimited shelf-life
- Celtra® Press Muffle ring: unlimited shelf-life
- Celtra® Press investment Plunger: unlimited shelf-life
- Firing Pad: unlimited shelf-life

15. List of applied national and international regulatory documents/standards the medical product complies with

Standard #	Standard name
EN ISO 9001:2008	Quality management systems. Requirements.
EN ISO 13485:2012 + AC:2012	Medical devices. Quality management systems. Requirements for regulatory purposes
EEC Directive 93/42/EEC; Annex II, excluding Section 4	Directive 93/42/EEC Medical devices
DIN EN 1041:2013	Information supplied by the manufacturer of medical devices
DIN EN 1641:2010	Dentistry. Medical devices for dentistry. Materials
DIN EN ISO 7405:2013-12	Dentistry. Evaluation of biocompatibility of medical devices used in dentistry
DIN EN ISO 10993-1:2010-04	Biological evaluation of medical devices. Part 1: Evaluation and testing within a risk management process.
DIN EN ISO 13485:2012	Medical devices. Quality management systems. Requirements for regulatory purposes.
DIN EN ISO 14971:2013	Medical devices - Application of risk management to medical devices
DIN EN ISO 15223-1:2013	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1. General requirements
DIN EN ISO/IEC 17050-1:2010	Conformity assessment. Supplier's declaration of conformity. Part 1. General requirements

16. Manufacturer Warranty

The manufacturer is not responsible and shall not indemnify for any possible damage or accidents caused by:

- Use of components unrelated to the product which can affect normal operation.
- Noncompliance with the Directions for Use.

The user is responsible for testing of materials for fitness and use for any purpose not expressly specified in the instructions!

For product quality issues, please contact the Authorized Manufacturer's Representative.

17. Authorized representative in the Russian Federation

Dentsply Sirona LLC, 115432, Moscow, Andropov prospect, 18, building 6, offices 1 39 40

Phone: +7 (495) 725-10-87

info@dentsplysirona.com

Appendix

Physical properties of the components of the Celtra® Press Starter Kit for the manufacture of glass ceramic restorations

Dentsply Prosthetics Die Material composite material

Physical property	Values
Appearance	Homogeneous paste, free of foreign matter Main color: yellowish Catalyst color: white
Use	The material is in working condition for at least 5 minutes from the moment it is extruded from the syringe
Curing	Easily removed from the base of the crown
Syringe dimensions	Length: 142 mm ± 10% Diameter: 10 mm ± 10%
Mass of material	Not less than 4 g

Celtra® Ceram Enamel

Physical property	Values
Coefficient of thermal expansion (CTE)	(CTE): $9,0 \pm 0,5 \times 10^{-6} \text{K}^{-1}$ (25-500°C)
Glass transition temperature – Tg	Tg = 520°C
Flexural strength	112±19 MPa
Radiactivity	< 0,405 Bq/g
Chemical solubility	≤0,05%
Linear shrinkage during firing	Not more than 16 %
Resistance to coloration, the quantity of spots on samples after holding in a dyeing solution	0
Porosity	The number of pores with a diameter of 30 µm shall not exceed 16, the number of pores with a diameter of 40 to 150 µm shall not be more than 6 on the surface with an area of 1 mm². Pores with a diameter of 150 µm shall not be in the place.
Bottle dimensions	Height: 58 mm ± 10% Diameter: 35 mm ± 10%
Mass of enamel	Not less than 15 g

DS (Dentsply Sirona) Universal Overglaze - High Flu

Physical property	Values
Appearance	Opaque paste
Flexural strength	≥100 MPa
Chemical solubility	<0.05%
Coefficient of thermal expansion (CTE)**	$(7.9 \pm 0.1) \times 10^{-6} \text{K}^{-1}$
Glass transition temperature, °C**	517±0.7
Softening temperature, °C**	576±0.1
Porosity	The number of pores with a diameter of 30 µm shall not exceed 16 and the number of pores with a diameter of 40 to 150 µm shall not be more than 6 on the surface with an area of 1 mm². Pores with a diameter over 150 µm shall not be in the place.
Linear shrinkage	Not more than 16%
Resistance to coloration, the quantity of spots on samples after holding in a dyeing solution	0
Bottle dimensions	Height: 25 mm ± 10% Diameter: 30 mm ± 10%
Mass of overglaze	Not less than 5 g

Universal Stain & Glaze Liquid

Physical property	Values
Appearance	Clear colorless liquid
Odour	Odourless
Melting / Freezing Temperature	-77°C
Normal boiling point	207.5°C
Vapour density (air = 1)	3.1
Solubility	Completely soluble in water
Viscosity	104 mPa*s
Bottle dimensions	Height: 70 mm ± 10% Diameter: 25 mm ± 10%
Volume of liquid	Not less than 15 ml

Physical property	Values
Appearance	Clear liquid
Odour	Odourless
Melting / Freezing Temperature	0°C
Normal boiling point	100°C
Vapour density (air = 1)	Completely soluble in water
Solubility	Similar to water
Viscosity	Height: 70 mm ± 10% Diameter: 25 mm ± 10%
Bottle dimensions	Not less than 15 ml

Physical property	Values
Appearance	Hard wax
Odour	Apricot
Colour	Moderate
Spontaneous ignition	> 300 °C
Vapour pressure	<0.01 mbar at 25 ° C
Density	~ 0.94 g/cm³
Solubility	Insoluble in water
Bottle dimensions	Height: 26 mm ± 10% Diameter: 76 mm ± 10%
Mass of wax	Not less than 60 g

Celtra® Modeling Wax

Physical property	Values
Appearance	Hard wax
Color	Light blue
Odour	Moderate
Spontaneous ignition	> 300 °C
Vapor pressure	<0.01 mbar at 25 ° C
Density	~ 0.94 g/cm ³
Solubility	Insoluble in water
Bottle dimensions	Height: 26 mm ± 10% Diameter: 76 mm ± 10%
Mass of wax	Not less than 60 g

Celtra® Press investment

Physical property	Values
Delivery form	Powder
Color	White to grey
Odour	Fruity
pH	Approx. 5
Melting temperature	>1500°C
Solubility in water	Low solubility
Bulk density	1100-1200 kg/m ³
Packaging dimensions	Length: 122 mm ± 10% Width: 130 mm ± 10%
Mass of the investment material	Not less than 100 g

Celtra® Press investment Liquid

Physical property	Values
Composition	Aqueous solution of alkali-stabilized colloidal silicic acid
Appearance	Liquid
Color	Milky white
Smell	Odourless
pH	8-10
Freezing point (t)	Approx. 0°C
Boiling point (t)	Approx. 100°C
Density (at 20°C)	1.1-1.4 g/cm ³
Vapor pressure	Approx. 20 GPa
Bottle dimensions	Height: 118 mm ± 10% Length: 45 mm ± 10% Width: 30 mm ± 10%
Liquid volume	Not less than 100 ml

Celtra® Press Muffle former, silicon-based

Dimensions:

Mass: 105,96 g ± 10%, length - 74 mm ± 10%, diameter - 53 mm ± 10%

Material: Silicone

Celtra® Press Muffle ring, 2 parts (base and lid)

Dimensions:

Mass: (each separate part)

Upper part (lid): weight - 35,45 g ± 10%, height - 43 mm ± 10%, diameter - 53 mm ± 10%

Lower part (base): weight - 48,63 g ± 10%, height - 82 mm ± 10%, diameter - 53 mm ± 10%

Material: POM black

Celtra® Press investment Plunger

Dimensions:

Mass: 8,43 g ± 10%, length - 35 mm ± 10%, diameter - 13 mm ± 10%

Material:

Component
Quartz powder (SiO ₂)
Cristobalite
Magnesium oxide

Firing Pad

Dimensions of the package: 57 mm x 58 mm x 5 mm ± 10%

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Eberhard Uhlig, Notar

Document Register No.: 1369 / 2019 U

I herewith certify the signatures enforced in front of me of

1. Mr. **Jörg Meister**, born on September 23th 1966,
with registered office at Rodenbacher Chaussee 4, 63457 Hanau,
2. Mrs. **Kristina Neeb**, born on June 12th 1971,
with registered office at Rodenbacher Chaussee 4, 63457 Hanau.

Mr. **Jörg Meister** acts as Managing Director and Mrs. **Kristina Neeb** acts as authorized signatory for DeguDent GmbH with registered office at Hanau.

I hereby certify on the basis of today's inspection of the electronic commercial register at the local court in Hanau under HRB 7182, that Mr. **Jörg Meister** and Mrs. **Kristina Neeb** are entitled to represent the company DeguDent GmbH at Hanau jointly.

Mr. **Jörg Meister** and Mrs. **Kristina Neeb** are personally known by the notary.

I do establish at the same time, that the interested parties denied the question of a prior involvement according to § 3 Sec. No. 7 BeurkG.

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