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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of URIT MEDICAL
ELECTRONIC CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized

Signature:

HAN QI

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(Date: Jul. 15, 2024)

验证网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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ADD: No. D-07 Information Industry District, High-Tech Zone, Guilin, Guangxi 541004, P.R. China

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E-MAIL: export@uritest.com

TEL: +86 773 2288583 FAX: +86 773 2288560

HTTP://www.urit.com

Instruction For Use


URIT D43 Probe Cleaner

APPROVED BY:

URIT Medical Electronic Co., Ltd., No. D-07, Information Industry District, High-Tech Zone, Guilin City, Guangxi, China

Name, Surname: Antique Li

Position held: Sales manager of URIT Medical Electronic Co., Ltd.

Signature: 

Place and Date of issue:

Guilin, China

July.1,2024

Instruction For Use

URIT D43 Probe Cleaner

manufacturer
URIT Medical Electronic Co., Ltd.,
P.R. China

for use applies to the medical device "URIT D43 Probe Cleaner" (hereinafter - device,

Use: probe cleaning solution is designed to clean the sample probe and instrument fluid lines. The agent can be used for daily and periodic maintenance. The solution emulsifies fats and dissolves hardened dirt to keep the measuring system clean.

Description of target analyte:

not applicable.

Functional purpose: in vitro diagnostics.

A specific pathology, condition or risk factor that an in vitro diagnostic medical device is designed to detect, identify or differentiate:

not applicable.

Type of sample analysed: not applicable.

Field of application: clinical laboratory diagnostics.

Potential user: the medical device is intended for use in various clinical and laboratory settings, including, without limitation, hospitals and laboratories, for professional use. Personnel must have the necessary qualifications to perform in vitro diagnostic tests.

Indications for use: the medical device must be used in accordance with its intended purpose.

Contraindications: not applicable.

Possible side effects: the medical device has no side effects if the operating instructions are followed.

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