More than
75
Years

Shaping the Future of Endoscopy with you

77C 6771
STORZ
KARL STORZ — ENDOSKOPE

KARL STORZ SE & Co. KG • PO Box 230 • 78503 Tuttlingen/Germany

APPROVED


**Senior Director Regulatory Affairs
Global Patient Health & Regulatory Compliance**

(position)

i. A. Dr. Bettina Hafen

(name)

**Karl Storz SE & Co. KG
Dr.-Karl-Storz-Straße
34, 78532, Tuttlingen,
Germany**

(signature)

18.12.2023
(«day» month (numerals))

(Stamp)

INSTRUCTION MANUAL

Video-esophagoscope flexible CMOS, with accessories

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

Flexibles Video-Oesophagoskop –
Serie 13303 E, Monitor 8403 ZX

INSTRUCTION MANUAL

Flexible Video Esophagoscope –
Series 13303 E, Monitor 8403 ZX

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It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and, where required, sterilization of the Video-Oesophagoscope. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

WARNING: KARL STORZ instruments are provided unsterilized, and must be cleaned, disinfected and/or sterilized prior to the initial use and before each subsequent use.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be deduced from the information, figures and descriptions in this manual.

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1	Pictures of the instruments	4
2	Symbol description	6
2.1	Pictograms	6
2.2	Symbols for safety instructions	8
3	Intended use	9
3.1	Indications	9
3.2	Contraindications	9
4	User qualification	10
5	Safety information	11
6	Unpacking the device	15
6.1	Basic equipment	15
6.2	Optional accessories	15
6.3	Operating conditions	15
6.4	Electromagnetic compatibility	15
7	External system components	16
7.1	Charger	16
7.2	Country specific configuration	16
8	Handling	17
8.1	Operating	17
8.2	Holding correctly	17
8.3	Handling the shaft	18
8.4	Handling the deflection mechanism	19
8.5	Connection to the monitor 8403 ZX	20
8.6	Use of the HDMI output	21
8.7	Functions of the C-MAC® Monitor 8403 ZX	21
8.8	Focus	21
8.9	Function buttons of the endoscope	22
8.10	Monitor function keys	22
8.11	Energy management	35
8.12	Fixing the monitor to a stand	37
9	Preparation prior to use	38
9.1	Inspecting the surfaces	38
9.2	Checking the deflection	39
9.3	Switching the Video-Esophagoscope on/off	40

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9.4	Inspecting the image quality	41
9.5	Switching off the C-MAC® Monitor	41
10	System reset.....	42
11	Monitor bag 8403 YD.....	42
12	Disassembly	42
13	Reprocessing	43
13.1	Overview of reprocessing procedures.....	43
13.2	Accessories	45
13.3	Preliminary cleaning directly after use.....	45
13.4	Leakage test.....	46
13.4.1	Performing the leakage test – dry.....	46
13.4.2	Performing the leakage test – in liquid.....	46
13.5	Manual reprocessing.....	47
13.5.1	Manual cleaning	47
13.5.2	Manual disinfection	48
13.5.3	Manual High Level Disinfection (HLD)	48
13.5.4	Drying	48
13.6	Automated reprocessing.....	49
13.6.1	Manual preparation	49
13.6.2	Machine cleaning/chemical-thermal disinfection	49
13.7	Assembly, inspection and care.....	49
13.8	Packaging.....	50
13.9	Sterilization	51
13.9.1	Hydrogen peroxide (H ₂ O ₂) sterilization - ASP STERRAD®	51
13.9.2	Hydrogen peroxide (H ₂ O ₂) sterilization - STERIS® AMSCO® V-PRO™	52
13.9.3	Ethylene oxide sterilization (EO).....	52
13.9.4	Low-temperature steam and formaldehyde process (LTSP)	53
13.10	Limits of reprocessing.....	53

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


Controls, displays, connectors and their uses

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- ① Connecting cord
- ② Connection for leakage tester or pressure compensation cap (vent port)
- ③ Distal tip
- ④ Shaft
- ⑤ Handle
- ⑥ Deflection lever
- ⑦ Leakage tester 13242 XL
- ⑧ Pressure compensation cap 11025 E
- ⑨ Monitor
- ⑩ Connection 2 for video endoscope
- ⑪ Connection 1 for video endoscope
- ⑫ Mains connection
- ⑬ Reset button
- ⑭ Holder for SD memory card
- ⑮ Status bar
- ⑯ Pilot lamp
- ⑰ On/Off switch
- ⑱ Pilot lamp charging (orange)
- ⑲ Softkeys
- ⑳ Function symbols
- ㉑ USB port
- ㉒ HDMI port
- ㉓ Video recording Start/Stop
- ㉔ White balance
- ㉕ Still image recording

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2 Symbol description

2.1 Pictograms

Refer to instruction manual
White balance
Single image capture
Video recording
This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
Device type BF
Electronic information product pollution control (China RoHS)
Manufacturer
UL approval mark, valid in the USA and Canada
Degree of protection in accordance with International Protection Codes: Protection against contact with tools or similar and small foreign bodies, protection against vertically dripping water
Rated ambient temperature/transducer temperature class
DC voltage connection
Consult instructions for use

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2.2 Symbols for safety
instructions

WARNING: A Warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a Warning could result in injury to the patient or personnel.

CAUTION: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.

NOTE: A Note indicates special information regarding operation or clarifies important information.

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3 Intended use

The video esophagoscope is used for endoscopic imaging and image transmission, thus enabling an endoscopic diagnosis of body cavities and hollow organs in esophagoscopy. It is also used for inspection of the upper and lower respiratory tracts, endotracheal intubation as well as inspection of the tube position during PCT. The video esophagoscope is designed for short-term use in invasive procedures through natural orifices and surgical invasive interventions.

3.1 Indications

The use of the instrument set for esophagoscopy is indicated if esophagoscopy is indicated in the opinion of the attending physician. This includes, amongst other things, the diagnosis and treatment of:

- Removal of foreign bodies from the esophagus
- Eosinophilic esophagitis
- Stenoses
- Diagnosis of disorders of the esophagus, e.g., esophagitis, tumors, bleeding
- Biopsy sample
- Endoscopic ablation with esophagus tumors/squamous cell carcinoma of the hypopharynx and the upper third of the esophagus

3.2 Contraindications

The use of the instrument set for esophagoscopy and diverticuloscopy is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered through such use, due to the general condition of the patient or if the endoscopic method itself is contraindicated.

Esophagoscopes and instrument sets must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

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The following contraindications also apply:

- General anesthesia
- Coagulation disorders
- ASA > 3
- Multiple diverticula
- Esophageal stricture

4 User qualification

The KARL STORZ Video-Esophagoscope may only be used by persons with an appropriate medical qualification and who are acquainted with the esophagoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the Video-Esophagoscope. It is not intended as an introduction to the technique of esophagoscopy.

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5 Safety information

Safety precautions when using the Video-Esophagосcope

The Video-Esophagосcope must be used according to the recognized medical rules and procedures of endoscopy.

NOTE: The patient leakage current limit of the CF application part may be exceeded when CF application parts (application on heart) are used in combination with other BF application parts. In this case, the reduction of the patient leakage current should be effected using only CF application parts as far as this is possible.

NOTE: The flexible endoscope does not offer any electric insulation against high frequency (HF) currents. Protection against high frequency currents should be ensured through the use of correct HF instruments. Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

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Risks of infection

- Warning: The Video-Esophagoscope is delivered unsterilized and must be cleaned and disinfected/sterilized before initial and subsequent use.
- Warning: In order to prevent infection and protect personnel, it is strictly forbidden to ship contaminated medical devices. To avoid contact and aerogen infections, medical devices must be decontaminated on-site.

Risks due to electric energy

- Caution: Do not allow any live components or instruments to come into contact with the Video-Esophagoscope.
- Warning: Danger of glare and burns! Never look directly into the light of the Video-Esophagoscope. The LED light radiated through the Video-Esophagoscope may cause glare or high temperatures at the distal light outlet and therefore lead to burns.
- Warning: Danger of explosion if flammable gases are used in the immediate proximity of the Video-Esophagoscope. The Video-Esophagoscope must not be used in the presence of flammable or combustible media.
- Warning: Always unplug the Video-Esophagoscope before carrying out any maintenance and cleaning work on it.
- Warning: Danger of electric shock! Refer servicing for this device to authorized personnel only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims. Exceptions are those that are described in the instruction manual.

Risks due to damaged parts

- Warning: Every time the Video-Esophagoscope is used, you must first check that it, and any accessories used in combination with it, is in perfect condition. Damaged Video-Esophagoscopes or damaged accessories must not be used.

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- **Warning:** The Video-Esophagосcope may only be repaired by personnel authorized by KARL STORZ. Unauthorized repairs may impair functioning and safety, and are therefore strictly forbidden. KARL STORZ gives no warranty for Video-Esophagосcopes that are repaired by unauthorized personnel.

Risks due to combination of system components

- **Warning:** Combinations of medical devices are only then assured to be safe if
 - they are identified as such in the respective instruction manual or
 - the intended purpose and the interface specifications of the devices used in combination permit this.
- **Warning:** Only use accessories that are specifically suited for use with Video-Esophagосcopes and can be employed in line with generally accepted technical standards and safety guidelines.
- **Caution:** Only use accessories from KARL STORZ. Accessories from other manufacturers may damage fiberscopes and the flexible videoscopes.
- **Warning:** Only KARL STORZ devices (e.g. monitor 8403 ZX) may be used. Disregarding this warning can endanger the patient and user.
- **Warning:** The Video-Esophagосcope must not be used during the discharge of a cardiac defibrillator. Remove the Video-Esophagосcope from the operation area before discharge takes place.
- **Warning:** When in operation, the integrated LED light source heats up the Video-Esophagосcope. Switch off the Video-Esophagосcope before putting it down and do not place it in the proximity of the patient (as per IEC 60601-1).
- **Warning:** Damage to the Video-Esophagосcope and injury to the patient or user resulting from incorrect operation are not covered by the manufacturer's warranty.

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Damage to the instrument

- Warning: If the Video-Esophagoscope should malfunction during application on a patient, stop the application immediately. Place the distal tip of the Video-Esophagoscope in the neutral, nondeflected position and remove the Video-Esophagoscope slowly and carefully from the patient.
- Caution: Do not steam sterilize (autoclave) the Video-Esophagoscope under any conditions.
- Caution: The maximum permissible temperature for decontamination, sterilization, transport and storage is 65 °C.
- Caution: Do not place the Video-Esophagoscope on devices that become very hot. Do not expose it to direct sunlight or excessive heat, to avoid damage to the device.
- Caution: Do not hit the distal tip of the Video-Esophagoscope against hard objects, and never use a sharp object to remove dirt.
- Caution: Never lay heavy objects on the Video-Esophagoscope.
- Caution: To avoid damaging the monitor, do not remove the SD card or connecting cord while a video is recording.
- Caution: To avoid damaging the monitor, always switch it off before replacing the SD card.

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6 Unpacking the device

Check for missing items and evidence of shipping damage.

File any complaints immediately with KARL STORZ or the supplier.

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6.1 Basic equipment

Video-Esophagoscope 13303 E with the following accessories:

- 1x Carrying case 27677 SM
- 1x Instruction manual
- 1x Leakage tester 13242 XL
- 1x Pressure compensation cap 11025 E

6.2 Optional accessories:

- Transportation protection sheet 13990 SFN
- Protective cap 8401 YZ
- Adaptor for leakage test in the washer and disinfectant (green) 13242 XA

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7 External system components

7.1 Charger

Using the medically approved charger (Art. No. ET27-30-0004515), the Li-Ion battery of the C-MAC® Monitor can be recharged extremely quickly.

Moreover, this charger enables the C-MAC® Monitor to be operated with a discharged Li-Ion battery.

The charger has a green LED display (LED), which indicates that the output voltage is correct.

7.2 Country specific configuration

The charger features interchangeable primary adapters, which allow the system to be adapted to the country-specific power supply outlets.

- a Primary adaptor for EUROPE:
ET27-30-0004 516
- b Primary adaptor for UK:
ET27-30-0004 517
- c Primary adaptor for US/Japan:
ET27-30-0004 518
- d Primary adaptor for Australia:
ET27-30-0004 519
- e Primary adaptor for USA:
ET27-30-0004 520
- f Adaptor cable-ET27-30-0004 370

Press the button on the power supply unit to change the adaptor (fig. 1). This releases the adaptor which can then be lifted out.

When inserting a new adaptor, this must be placed in the guide (fig. 2) and pushed toward the cable until there is an audible click (fig. 2, see arrow).

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8 Handling

NOTE: Damage to the Video-Esophagoscope resulting from incorrect operation is not covered by the manufacturer's warranty.

WARNING: The Video-Esophagoscope is delivered unsterilized, so it has to be reprocessed prior to initial use and before each subsequent reuse.

CAUTION: Do not store or disinfect the flexible Video-Esophagoscope together with sharp objects in order to avoid causing damage to the outer casing.

NOTE: The C-MAC® monitor has passed a fall and crash test successfully.

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8.1 Operating

The following instructions outline the correct operation of the Video-Esophagoscope and the corresponding accessories. This is not an introduction to the techniques of esophagoscopy. Instructions on the endoscopic technique can be found in the relevant medical literature.

The Video-Esophagoscope is intended for visualization and diagnostic purposes.

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8.2 Holding correctly

Grip the handle of the Video-Esophagoscope in your left or right hand, as shown in the picture. With this kind of control grip, the Video-Esophagoscope can be operated with one hand.

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8.3 Handling the shaft

CAUTION: Do not carry the Video-Esophagoscope by the shaft. Do not pull, clamp or twist the shaft of the Video-Esophagoscope.

CAUTION: The shaft of the Video-Esophagoscope must never be kinked or coiled up too tightly. It is best to store the Video-Esophagoscope in a suspended position or to keep it in a tray specifically designed for sterilizing and storage (39405 AS).

WARNING: Do not apply lubricant to the lens, to avoid clouding the image and thus posing a risk to the patient!

Use a suitable water-soluble lubricant for the distal section of the shaft (For materials recommended and approved insertion by KARL STORZ, please see www.karlstorz.com).

If, during use, the image fails, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the shaft of the Video-Esophagoscope out of the patient's body.

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
8.4 Handling the deflection mechanism

- !** **NOTE:** The neutral position of the distal tip is straight, without any deflection whatsoever.
- !** **WARNING:** The deflection lever must never be operated with force or abruptly to avoid perforation and injury to the patient or damaging the device irreparably.
- !** **WARNING:** The distal tip of the Video-Esophagoscope must be straight, not angled, when the videoscope is inserted into the patient.
- !** **WARNING:** Before removing the Video-Esophagoscope from the patient, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.
- !** **WARNING:** If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the Video-Esophagoscope must not be used.
- !** **WARNING:** If a functional defect occurs with the Video-Esophagoscope during use on the patient, discontinue use immediately. Put the distal tip of the Video-Esophagoscope into the neutral, non-deflected position and remove the Video-Esophagoscope slowly and carefully from the patient.
- !** **CAUTION:** The resistance that can be felt when deflecting the tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.

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WARNING: The lens at the distal end of the Video-Esophagoscope may become detached, e.g. due to mechanical forces. For this reason, it is essential to check that the lens is undamaged and present before and after decontamination.

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8.5 Connection to the monitor 8403 ZX

CAUTION: The Video-Esophagoscope can be used together with the KARL STORZ monitor 8403 ZX. The instruction manual for the monitor must be observed.

WARNING: It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces, or voltage-carrying elements of other units, at any time.

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

Insert the connecting cord ① of the flexible Video-Esophagoscope into either connection socket ⑩ or ⑪ on the monitor.


The connecting cord ① can be connected or disconnected during operation. Insert the cable from the power supply unit into the connection socket ⑫ of the C-MAC® Monitor and connect the power supply unit with the mains socket. The power cord can be fixed to the connection socket using the union nut.

The battery is charged automatically as soon as the C-MAC® Monitor is connected to the power supply. This will be indicated by the lighting up of the orange colored LED ⑬.

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8.6 Use of the HDMI output

The device features a port for an external monitor. An HDMI cable can be used to send the video signal to an external monitor, if necessary.

! CAUTION: The C-MAC® Monitor must be switched off when an external monitor is connected.

i NOTE: The connected monitor must support a resolution of 1024 x 768 or 1280 x 720.

i NOTE: To connect the C-MAC® Monitor to a KARL STORZ monitor, the HDMI-DVI adaptor with article number 209190 10 is required.

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
8.7 Functions of the C-MAC® Monitor 8403 ZX

A more detailed description of the operation and functions of the monitor 8403 ZX can be found in the corresponding instruction manual supplied with every monitor.

In case of loss of the instruction manual, you can also download it from our homepage at www.karlstorz.com.

! WARNING: The KARL STORZ flexible Video-Oesophagoscope can not be used in combination with camera control units and light sources. The use of non-compatible units can lead to patient injury or unit damage.

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8.8 Focus

The image focus is adjusted automatically. Manual adjustment is not necessary.

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8.9 Function buttons of the endoscope

The endoscope is equipped with 3 keys, which are used to activate the functions displayed as symbols on the endoscope.

The following functions are available:

- Video recording : Starts or stops the recording of a video stream on the SD memory card inserted in the receptacle.
-  White balance : Enables the user to carry out a manual white balance. Point the tip of the videoscope at a white surface and press the button for 1 sec.
-  Single image capture : Saves the live image on the screen on the SD memory card.

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
8.10 Monitor function keys

The monitor is equipped with 4 softkeys. These can be used to activate the functions displayed as symbols on the screen.

NOTE: The **video recording** and **image storage** symbols are only visible if an SD memory card is inserted. The **change the camera unit** symbol is only visible if two video units are connected.

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8.11 Energy management

The charge state of the battery is displayed on the monitor by the charge indicator in the menu bar.

The symbols mean:

Battery empty

25 %: Approx. 30 minutes of operating time

50 %: Approx. 1 hour of operating time

75 %: Approx. 1.5 hours of operating time

100 %: Approx. 2 hours of operating time

Battery is being charged

Battery must be replaced
(Maximum battery capacity is less than 50 %
of the original maximum capacity)

NOTE: The C-MAC® Monitor can be powered using either mains electricity or the integrated battery.

The battery is charged automatically as soon as the monitor is connected to the power supply. This is displayed by the lighting up of the orange-colored LED.

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NOTE: The monitor should be recharged after each use or, at the latest, once the charge state has dropped below 50 %.

CAUTION: Damage to the power cord: To remove the power cord from the socket, pull only on the plug housing, never on the cord itself!

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The lithium-ion battery can be used for approx. 2 hours when charged.

The discharged battery is fully charged after 5 hours (when the device is switched off). During the charging time, the orange-colored LED lights up.

If the system is not used for 12 weeks or longer, you should check the charge state using the progress bar described on page 35 and, if necessary, charge.

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8.12 Fixing the monitor to a stand


Attach the stand holder to the rear of the monitor using 4 screws.

The VESA 75 holder (8401 YC), which is included in the scope of supply, is suitable for securing devices to tubes with a diameter of 25 mm.

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
Press the holder together to allow the stand bar to be pushed through the openings.

Release the holder in the desired position. The selected position is fixed by the holder's basic tension.

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The optionally available clamp 8401 YB (see accessories) can also be secured to the VESA 75 holder and used to attach the monitor to tubes with a diameter of less than 25 mm and/or to rails with a cross section of less than 25 mm.

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
- 9 Preparation prior to use**
- NOTE:** Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding preliminary cleaning.
- WARNING:** Every time the Video-Esophagoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged Video-Esophagoscopes or damaged accessories must not be used.

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- 9.1 Inspecting the surfaces**
- WARNING:** The surface of the flexible shaft must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.
- CAUTION:** Do not squeeze the distal tip of the Video-Esophagoscope too tightly.

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Every time the Video-Esophagoscope is used, you must first check that its surface, and above all the surface of the shaft and distal tip, are in perfect condition and that there are no sharp edges. Video-Esophagoscope with sharp edges must not be used. Otherwise serious injuries may be caused to the patient or the user.

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Inspect the connecting cord for damage (such as breaks, cracks, twisted or squashed sections...) and check the supply connector for loose connections.

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9.2 Checking the deflection

WARNING: Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated with force or abruptly, to avoid injuring the patient or damaging the device! Move the deflection lever smoothly and slowly!

WARNING: Before use, always check that the deflection lever and the distal tip connected to it function correctly. Defective Video-Esophagoscopes must not be used.

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
The distal tip is deflected by actuating the deflection lever slowly. The distal end will deflect in accordance with the direction of actuation (forward or backward).

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9.3 Switching the Video-Esophagoscope on/off

Switching on/off: The Video-Esophagoscope with the integrated LED light source is switched on by pressing the power button 17 on the monitor. Pressing the key a second time switches the unit off.

WARNING: Danger of glare: Do not look directly into the light beam.

WARNING: When in operation, the integrated LED light source heats up the Video-Esophagoscope. Switch off the Video-Esophagoscope before putting it down and do not place it in the proximity of the patient (as per IEC 60601-1).

NOTE: The function keys 17 on the monitor are used for switching the unit on/off. The Video-Esophagoscope itself does not have any function keys.

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
Video recording: When the control key "Video recording" on the C-MAC® Monitor or the flexible Video-Esophagoscope is pressed, the flexible Video-Esophagoscope records a video stream on the SD memory card. Pressing the key a second time stops the recording.


NOTE: The video stream can only be viewed on a PC if an MPEG 4-codec is installed on the PC.

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Single image capture: When the control key "Single image capture" on the C-MAC® Monitor or the flexible Video-Esophagoscope is pressed, the current image on the monitor is saved on the SD memory card.

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9.4 Inspecting the image quality


WARNING: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the Video-Endoscope from the patient's body gently and carefully. If the Video-Endoscope is used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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9.5 Switching off the C-MAC® Monitor

Standby: Briefly pressing the power button (17) will set the C-MAC® Monitor to standby operation mode. In the case of brief interruptions, this enables the C-MAC® Monitor to be ready for operation again quickly. The standby mode is indicated by a flashing green pilot lamp (16) above the power button.

After 5 minutes, the system switches itself from standby mode to the off mode.

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Power off: Pressing the power button (17) for more than approx. 3 seconds will switch off the C-MAC® Monitor. The green pilot lamp (16) will go out.

NOTE: When the C-MAC® Monitor is switched off, the integrated light source in the flexible Video-Endoscope also switches off.

NOTE: Information on the operation and functions of the C-MAC® Monitor 8403 ZX can be found in the corresponding instruction manual supplied with every C-MAC® Monitor (Art. No. 96076008DEF). In case of loss of the instruction manual, you can also download it from our homepage at www.karlstorz.com.

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
10 System reset

If a system error occurs (unit no longer reacts, although connections and power supply are in order), a system reset has to be implemented. This is the equivalent of rebooting a computer. For this, a thin screwdriver or similar is inserted into the round opening next to the SD card holder ⑭ and carefully pushed downward in order to activate the reset button ⑬.

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
11 Monitor bag 8403 YD

The optionally available protective bag 8403 YD facilitates the mobile use of the monitor.

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12 Disassembly

Loosen the union nut and remove the connecting cord from the power supply socket ⑫.

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13 Reprocessing

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13.1 Übersicht

13.1 Overview of reprocessing procedures

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This reprocessing step is part of the validated procedure

- This reprocessing step may be carried out additionally in order to aid the cleaning. The sterilization procedure is approved in respect of material compatibility; however, it must be validated by the operator on site in respect of effectiveness.

- ✗ On the grounds of potential damage, this reprocessing step must NEVER be performed

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WARNING: Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

WARNING: Risk of infection from non-sterile delivery. Reprocess the device before each use. Check the medical device for visible soiling before use.

WARNING: Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

ATTENTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

CAUTION: Never clean flexible endoscopes with metal brushes.

CAUTION: Constant changing between various reprocessing methods is hard on the materials and should be avoided. However, a one-off change of method, due to procurement of new cleaning equipment, for example, is possible.

CAUTION: Flexible endoscopes must never be steam sterilized (autoclaved).

CAUTION: To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.

NOTE: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

NOTE: Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.

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- ❗ **NOTE:** Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be found at www.karlstorz.com.
- ❗ **NOTE:** National laws and regulations must be observed.
- ❗ **NOTE:** We recommend using a plastic basin for disinfection to avoid damaging the endoscope. Furthermore, use of a plastic basin prevents electrolytic corrosion which can occur if there are different metals in a solution.
- ❗ **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

13.2 Accessories

Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/adaptor
Clean surface	27652
Ventilation connector	13242 XA, 11025 E

13.3 Preliminary cleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

Wipe-clean the exterior surfaces of the endoscope with a moistened cloth.

- ❗ **NOTE:** To avoid damaging the shaft, do not use pushing and pulling movements.
- ❗ **NOTE:** Do not use disinfectants for preliminary cleaning, these cause fixation of proteins and organic residues.

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13.4 Leakage test

The endoscope should be tested to ensure it is leakproof after use and before reprocessing.

13.4.1 Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch 25° clockwise.

NOTE: The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

Inflate the leakage tester until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on manometer wait at least 2 minutes to see whether the pressure remains constant or drops.

13.4.2 Performing the leakage test – in liquid

NOTE: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

NOTE: Any lumina must be covered with no air bubbles by filling in a targeted manner.

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in 'Leakage test – dry' and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously)

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from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

NOTE: Never release the pressure in the liquid.

Cleaning can then begin.

If a pressure loss occurs, this means that the endoscope leaks and may no longer be used, immersed in liquid or reprocessed by machine.

CAUTION: In this case, fill out the accompanying repairs form (see Appendix or Art. no. 96096004DER), mark the endoscope, which has been cleaned by wiping to remove organic soiling/residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.

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CAUTION: Never remove dirt with a needle or other sharp objects, since these can cause severe damage.

CAUTION: The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

13.5.2 Manual disinfection

The endoscope must be completely immersed in a cleaning solution. At the end of the exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. For this, water of the best possible quality must be used, taking into account the national regulations. Sterile filtered water must be used for the final rinsing.

13.5.3 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Dialdehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

13.5.4 Drying

Finally, the surfaces are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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CAUTION: The pressure must not exceed 100 kPa (1 bar / 14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

13.6 Automated reprocessing

The following procedure for machine decontamination has been validated and approved subject to compliance with the process parameters described in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments":

13.6.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush.

Finally, it must be rinsed with cold water to ensure neutralization.

13.6.2 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C and a maximum pressure of 100 kPa (1 bar / 14.5 psi) must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off manually afterwards.

13.7 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

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- NOTE:** Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.
- Damaged or corroded medical devices must be withdrawn from use.
- NOTE:** Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.
- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surface is to be dried completely (e.g., plug connections) with (preferably sterile) compressed air and according to national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.
- CAUTION:** The pressure must not exceed 100 kPa (1 bar / 14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.
- Afterwards, a functional check must be carried out.
- NOTE:** Check for a clear view and steerability.

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13.9 Sterilization

WARNING: Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

WARNING: Risk of infection: Sterilize the medical device in a disassembled state.

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CAUTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

he procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

13.9.1 Hydrogen peroxide (H₂O₂) sterilization - ASP STERRAD®

WARNING: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

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i NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S Long Cycle with Boosters
- STERRAD® 100S Short Cycle with Boosters
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

13.9.2 Hydrogen peroxide (H₂O₂) sterilization - STERIS® AMSCO® V-PRO™

i NOTE: The 'STERIS Low Temperature Reprocessing Device Compatibility Matrix' can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® maX Flexible Cycle
- V-PRO® maX Non-Lumen Cycle
- V-PRO® 60 Flexible Cycle
- V-PRO® 60 Non-Lumen Cycle
- V-PRO® 1 PLUS Non-Lumen Cycle

13.9.3 Ethylene oxide sterilization (EO)

The ethylene oxide procedure is validated with 100% ethylene oxide at 55°C and a hold time of 60 minutes.

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The following procedures are approved in respect of material compatibility, but must be validated by the on-site operator with regard to efficacy:

13.9.4 Low-temperature steam and formaldehyde process (LTSF)

KARL STORZ had only the specified sterilization method validated with regard to a reduction in microorganisms on the order of 10^{-6} SAL. However, due to the large number of processes used in the global market, validation in respect of efficacy was not carried out. Therefore, the user must validate his specific processes individually in collaboration with the device manufacturer.

WARNING: Risk of infection: if FO sterilization is performed, the method must be validated on site.

13.10 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

14 Storage

CAUTION: Do not expose the Video-Esophagoscope to direct sunlight or excessive heat.

CAUTION: The distal tip of the Video-Esophagoscope must be in the neutral position during storage.

CAUTION: Do not keep the Video-Esophagoscope in the carrying case.

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CAUTION: The shaft of the Video-Esophagoscope must never be kinked or coiled up too tightly. It is best to suspend the Video-Esophagoscope for storage or to keep it in a large container. Video-Esophagoscope should be packaged and stored in sterile cloths or suitable containers in accordance with EN 868, and must be absolutely dry before being put into storage. It is most suitable to store the Video-Esophagoscope suspended vertically with the distal end pointing down. The room where they are stored should be at a normal temperature and humidity. Store dismantled valves and caps separately.

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15 System description

The flexible Video-Esophagoscope from KARL STORZ has a 750 mm pellethane working shaft, 2.90 mm in diameter. The distal tip made of stainless steel is 2.9 mm in diameter and houses the video module. The distal tip of the flexible Video-Esophagoscope can be deflected 210° up or 130° down. The handle houses the deflection lever, the function keys and the LED light source.

The flexible Video Esophagoscope from KARL STORZ is connected to the recommended monitor 8403 ZX using the integrated 2 m supply cable. Images and videos can be saved using a SD memory card and the control keys on the endoscope or on the monitor, which are assigned, for example, the functions of video recording and image capture.

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Length: 230 mm
Width: 155 mm
Depth: 54 mm
Weight: 1250 g

Drip-water protection: IP54

! **CAUTION:** Drip-water protection can only be ensured if the connection sockets are sealed using the designated rubber seals.

Classification: IEC Class II
Input voltage: 12 V
Current consumption: 1.5 A max.
Battery: Lithium-Ion, 3.6 V 7800 mAh
Resolution: 1280 x 800
Screen width: 150 mm
Screen height: 94 mm
Memory card: SD card
Video out signal: The connected monitor must support the 720p mode.

Operating/storage conditions

Operating conditions: 0 °C to +40 °C
30 % to 70 % rel. humidity
Storage conditions: -10 °C to +60 °C
5 % to 95 % rel. humidity
Atmospheric pressure: 700 to 1080 hPa

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SD card

Throughput: 10 MB/s
Compatibility: Compatible with the Sandisk
Ultra II HD card and SDHC
cards class 10

Color settings

Basic Setting

Brightness 50 %
Contrast 35 %
Saturation 31 %
Hue 50 %

Endoscope 13303 E

Direction of view: 0°
Angle of view: 100°
Deflection: up 210°, down 130°
Working length: 75 cm
Outside diameter: 2.9 mm
Illumination: Integrated LED light source
Electrical power
consumption: 1W
Sensor: CMOS technology
Weight 326 g
Drip water
protection: IPX8
Classification: IIa
Applied part: Type BF

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17 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

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18 Maintenance

Preventive maintenance is not essential. Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

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19 Limitation of liability

KARL STORZ SE & Co. KG shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

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20 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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21 Standard compliance

According to IEC 60601-1, 60601-2-18, UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electric shock:
 - In combination with the monitor 8403 ZX: Internally powered during battery operation, Protection class II during line operation
 - In combination with the C-HUB®: Protection class II

- Degree of protection against electric shock: Applied part type BF 
- Type of protection against moisture: Protection against continuous immersion: IPX 8

According to IEC 60601-1-2:

Please observe the Electromagnetic Compatibility information in the appendix (pp. 64–80).

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22 Directive compliance

According to Medical Device Directive (MDD):

13303 E: This medical product belongs to Class IIa
8403 ZX: This medical product belongs to Class I

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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23 Infection prevention

To avoid the spread of infectious illnesses among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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24 Disposal

This unit has been marked in accordance with the European Directive 2012/19/EC on waste electrical and electronic equipment (WEEE). At the end of its service life, dispose of the device as electronic waste. Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information about your local collection point. Within the scope of application of this directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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Stand, height 120 cm, star-foot roller stand with antistatic castors, crossbar 25 cm x diameter 25 mm for positioning the monitor, including storage tray for laryngoscopes, measurements (w x d x h): 30 x 20 x 10 cm

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Crossbar, for stands 8401 YA, 50 cm x diameter 25 mm, for positioning the monitor 8403 ZX

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Crossbar, for stands 8401 YA, 70 cm x diameter 25 mm, for positioning the monitor 8401 ZX

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Multifunctional clamp with VESA 75 standard holder. Suitable for securing rods with a max. diameter of 24 mm and rails with a max. height of 25 mm

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Multifunctional cap to protect contacts during reprocessing. The cap is reusable.

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Holder for flexible endoscopes, for mounting to standard tubes, incl. installation accessories

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Electromagnetic Compatibility (EMC) Information

25.3 Electromagnetic Compatibility (EMC) Information

The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: The use of an accessory or cable with the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX other than those specified in this manual may result in increased emissions or decreased immunity of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX. When using an accessory or cable with the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX other than those specified in this manual, it becomes the responsibility of the user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX to determine compliance with IEC 60601-1-2 when using this item.

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Appendix

Electromagnetic Compatibility (EMC) Information

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 Guidelines and manufacturer's declaration – electromagnetic immunity			
The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is intended for use in an electromagnetic environment specified below. The user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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Table 3
Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting

The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{ms} 150 kHz to 80 MHz	3 V _{ms}	Portable and mobile RF communications equipment should be used no closer to any part of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ compliance level in each frequency range ^a . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is used exceeds the above compliance levels, the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^a Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 Guidelines and manufacturer's declaration – electromagnetic emissions		
The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is intended for use in the electromagnetic environment specified below. The customer or user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX uses RF energy only for its internal function. The customer or user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5 Recommended separation distances between portable and mobile RF telecommunication devices and the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX			
The flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the flexible Video-Esophagoscope 13303 E in conjunction with monitor 8403 ZX as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies. NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.			

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kse.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ Tra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzecuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6° Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skarholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi Oy
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio st. 59 / Lenkioji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: ausra.kukeniene@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlystøften 33
2840 Høtte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopy Nederland B. V.
Displayweg 2
3821 BT Amerstoft, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

KARL STORZ Endoskop Magyarországi Kft.
Toberek utca 2. fsz. 17/b
HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
E-Mail: info-hu@karlstorz.com

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colanin, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: info-ro@karlstorz.com

KARL STORZ Endoscopy Greece M.E.P.E.*
Patriarhou Grigorou E' 34
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
E-Mail: info-gr@karlstorz.com

*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No. 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy – WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: info-ru@karlstorz.com

TOV KARL STORZ Ukraine
Oholonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20
E-Mail: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG
Representation Office
Sabit Onudschow 1184, apt. 23
1025 Baku, Azerbaijan
Phone: +99 450 613 30 60
E-Mail: info-az@karlstorz.com

KARL STORZ ENDOSKOPE –
East Mediterranean and Gulf (Offshore) S.A.L.
Spark Tower 1st floor
Charles Helou St., Harch Tabet – Sin El Fil
Beirut, Lebanon
Phone: +961 1 501105, Fax: +961 1 501950
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za
TDO KARL STORZ Endoscopy Kazakhstan
Saryarka, 6, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing iFC
No. 8, Jiangmenmenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F, Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82 70 4350-7474, Fax: +82 70 8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

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**Addendum to Instruction manual for medical device:
“Video-esophagoscope flexible CMOS, with accessories”**

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The information contained in this supplement takes precedence in the Russian Federation over the information specified in the corresponding instruction manual for the medical device.

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Videosophagoscopes are endoscopes that conduct light with a built— in LED inside the body of patient. The chip on the tip of the video esophagoscopes transmits digital endoscopic images. These video endoscopes have a flexible body, and the tip can be actively tilted at an angle using a lever located on the handle.

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
The Video-Esophagoscope is intended for visualization and diagnostic purposes.
The video esophagoscope is designed for short-term use in invasive procedures through natural orifices and surgical invasive interventions.

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

The documentation contains a description of the technical characteristics of the device, visualization programs, as well as data on the use and maintenance of the product.

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The cap is used to protect the connectors of electrical contacts during cleaning and disinfection.

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
The adaptor is used to check the tightness

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The following measures must be taken for disposal:

1. After use, transfer the products to a safe, controlled and documented disposal process.
2. It is necessary to comply with national laws and regulations in force in a particular country.

When disposing of products, domestic regulations/laws must be observed. When disposing of products on the territory of the Russian Federation, the norms set out in SanPiN 2.1.3684-21 should be followed. Unused video endoscopes, as well as packaging materials, should be disposed of as epidemiologically safe waste of Class A. Used video endoscopes should be disposed of as epidemiologically hazardous waste of Class B.

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Identification	Designation
ISO 13485	Medical devices Quality management systems – Requirements for regulatory purposes
ISO 14971	Medical devices Application of risk management to medical devices
ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
ISO 20417	Medical devices - Information to be supplied by the manufacturer
ISO 17664-1	Processing of health care products - Information to be provided by the medical device manufacturer for the processing of medical devices - Part 1: Critical and semi-critical medical devices
IEC 60601-1 3.1 Edition	Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance

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Symbol	Description
	Manufacturer
	Date of manufacture
	The catalog number.
	Batch code
	Refer to the instructions for use or to the instructions for use in electronic form
	In accordance with US federal law (according to 21 CFR 801.109), the product is allowed to be sold only to doctors (licensed) or by prescription.
	In accordance with US federal law (according to 21 CFR 801.109), the product is allowed to be sold only to doctors (licensed) or by prescription.
	Warning
	Medical Device
	Serial number
	The number of products in the package
	Unique Device Identifier
	Unsterile

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