IFU
Container + Spoon

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of XIAMEN SCIENDOX BIOLOGICAL TECHNOLOGY CO., LTD. on the annexed DOCUMENT is genuine.

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Authorized
Signature: Zhu Chaoyang

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Xiamen Sciendox Biological Technology Co., Ltd			
Technical Document (Instruction for use)- Single use materials Sciendox for feces analysis: sample container, sampling spoon, applied solution tube			Page 1 of 10
Prepared	Fengmei Zhang	R & D personnel	张凤梅
Verified	Lingling Zheng	Assistant General Manager	郑玲玲
Approved	Lingling Zheng	Assistant General Manager	郑玲玲

INSTRUCTION FOR USE

Sample container and sampling spoon

Contents

1. DESCRIPTION	2
2. DESIGNATION and INDICATION	3
3. CONTRAINDICATIONS AND SIDE EFFECTS	3
4. PRECAUTIONS.....	3
5. INSTRUCTION FOR USE.....	4
6. SYMBOLS ON PACKAGE	4
7. GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.....	7
8. CONDITIONS FOR TRANSPORTATION, STORAGE AND OPERATION. SHELF LIFE.	7
9. WARRANTY OBLIGATION.....	8
10. POTENTIAL USERS	9
11. REQUIREMENTS FOR THE MAINTENANCE AND REPAIR OF MEDICAL DEVICES.....	9
12. FOR THE DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES.	9
13. THE NAME AND LEGAL ADDRESS OF THE MANUFACTURER OF MEDICAL DEVICES AND AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RF:	9
14. AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:	10

1. DESCRIPTION

Sample container

- Made of transparent plastic, the filter is installed at the bottom.
- Two-layer filter, in the first layer a filtration hole with a diameter of 2 mm, in the second layer a hole with a diameter of 250 microns.
- Cover color: white-green.


Figure 1. Sample container

Sampling spoon for feces specimen collection

- Used for sampling.
- If you press the button on the end of the spoon, the front part is automatically reset.
- Convenient for intake, a certain (fixed) amount of sample.
- Color spoons: white and green

2. DESIGNATION and INDICATION

The product is intended for sampling and preparation of stool samples for further research in laboratories and clinical diagnostic departments of medical institutions.

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4. PRECAUTIONS

Take one spoon for forming stool and thin pasty stool. Take 4 spoons for liquid stool.

Preferentially pick the part with blood or mucus, etc. Don't mix with water and foreign matter.

Keep the bottle exterior clean.

Please consult healthcare provider before use.



Biological safety

Feces sample maybe a potential pollution source.

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
5. INSTRUCTION FOR USE


1. Tear the bag
2. Sampling
3. Put into the container
4. Keep the container exterior clean
5. Push the white thin rod at the tail, the spoon head and handle would separate


1. Take appropriate amount of specimen.


2. Especially take specimen with blood or mucous.


3. Insert the spoon into the container.

4. Active the separation.

5. The spoon will separate automatically.

6. Close the cap and send for test.

6. SYMBOLS ON PACKAGE

on individual package

	Manufacturer
	Authorized representative in EC
	MD for in vitro diagnostics
	CE marking
	For single use

Consumer and group package

	Ref number
LOT :	Lot number
EXP :	Expiry date
QTY :	Quantity
	manufacturer
	Authorized representative in EC
	MD for in vitro diagnostics

	CE marking
	For single use
	Protect from moisture
	Do not heat

Transport package

	Fragile, handle with care
	Up
	Protect from moisture
	Do not heat
	Manufacturer
	Authorized representative in EC
	MD for in vitro diagnostics

CE	CE marking
LOT :	Lot
EXP :	Expiry date
QTY :	Quantity

7. GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.

Parameters of Feces specimen collection bottle

Overall dimensions, mm	Collar diameter-25±2, Height-79±2
Volume, ml, not less than	12
Weight, g,	20.7±2

Parameters of sampling spoon

Color	White and green
Overall dimensions, mm	Length-156±2, Width-12 ±2
Cup volume, ml, (for water),	0.4-0.5
Weight, g,	5.0±0.5

8. CONDITIONS FOR TRANSPORTATION, STORAGE AND OPERATION. SHELF LIFE.

The product is transported by all types of closed vehicles, in accordance with the rules of transportation of goods, acting on this type of transport, during transportation should be protected from strong shocks, rain and direct sunlight.

Transportation and storage

Temperature	10-55° C
Atmospheric pressure	50,0kPa -106,0kPa
Relative humidity	Must not exceed 85%

Operation

Temperature	15-30o C
Atmospheric pressure	86,0kPa -106,0kPa
Relative humidity	Must not exceed 70%

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Telephones: +7(495)234-57-31

10. POTENTIAL USERS

patient (during samples collection).

Medical laboratory technician (laboratory assistant) or other specialist. Operators must be trained (during sample processing together with Feces sample liquefaction processor Sciendox 1A).

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Not applicable.

12. FOR THE DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES.

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13. THE NAME AND LEGAL ADDRESS OF THE MANUFACTURER OF MEDICAL DEVICES AND

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RF:

MANUFACTURER:

Xiamen Sciendox Biological Technology Co., Ltd, China

Room 402, Building 3, No.2054 of Wengjiao West Road, Haicang District 361006 Xiamen, Fujian,
People's Republic of China

Manufacturing site:

Xiamen Sciendox Biological Technology Co., Ltd, China

Room 402, Building 3, No.2054 of Wengjiao West Road, Haicang District 361006 Xiamen,
Fujian, People's Republic of China

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Telephones: +7(495)234-57-31

**TOTAL NUMBERED,
BOUND AND SEALED
10 PAGES**

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THIS IS TO CERTIFY THAT: the seal of XIAMEN SCIENDO
BIOLOGICAL TECHNOLOGY CO., LTD. on the annexed DOCUMENT
is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Zhu Chaoyang

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(Date: Jul. 06, 2020)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Xiamen Sciendox Biological Technology Co., Ltd		
Technical Document – (Instruction for use) Single use materials Sciendox for feces analysis: sample container, sampling spoon, applied solution tube		
Prepared	Fengmei Zhang	R & D personnel
Verified	Lingling Zheng	Assistant General Manager
Approved	Lingling Zheng	Assistant General Manager

Page 1 of 2

INSTRUCTION FOR USE

Applied solution tube

Contents


1. DESCRIPTION	2
2. DESIGNATION and INDICATION	2
3. CONTRAINDICATIONS AND SIDE EFFECTS.....	2
4. PRECAUTIONS.....	2
5. INSTRUCTION FOR USE	3
MIXING.....	3
FILTERING	4
6. SYMBOLS ON PACKAGE	4
7. GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.....	6
8. CONDITIONS FOR TRANSPORTATION, STORAGE AND OPERATION. SHELF LIFE.....	7
9. WARRANTY OBLIGATION.....	8
10. POTENTIAL USERS	8
11. REQUIREMENTS FOR THE MAINTENANCE AND REPAIR OF MEDICAL DEVICES	8
12. REQUIREMENTS FOR THE THE DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES.	8
13. THE NAME AND LEGAL ADDRESS OF THE MANUFACTURER OF MEDICAL DEVICES AND AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RF:.....	9
14. AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:	9

1. DESCRIPTION

It is a plastic tube for one-time use, which is used for the treatment of feces with Sciendox equipment and analytical instruments.

The product is intended for preparation of stool samples for further research in laboratories and clinical diagnostic departments of medical institutions.

- Exclusive item for use with sample liquefaction processor and feces analysis system of Sciendox,
- Cannot use when broken or contaminated

Picture 1. Applied solution tube

2. DESIGNATION and INDICATION

The product is intended for preparation of stool samples for further research in laboratories and clinical diagnostic departments of medical institutions.

3. CONTRAINDICATIONS AND SIDE EFFECTS

There are no contraindications and side effects at this time.

4. PRECAUTIONS

biological safety

Feces sample maybe a potential pollution source,

Consumables and waste fluid which contact feces sample are should be disposed according to local law.

5. INSTRUCTION FOR USE

The device is used during sample preparation together with sample liquefaction processor and feces analysis system of Sciendox


It is used during the Mixing & Filtration procedure.

At first the sample container with the biomaterial inside is placed at the dilution position in the sample liquefaction processor Sciendox.


The container is punctured by needle and the saline is added to the bottle.

MIXING

When the device finish adding fluid put the container to mixing position (pic.2). After that, it is necessary to pull out the tongue for the applied solution tube, insert the tube into the special opening and close the tongue to the stop (pic. 3). Start the mixing and filtering process from the control panel.

Picture 2. Sample container in the position of mixing and filtration.

Picture 3. The tube in the position of filtering

The device mixes the sample according to the set parameters. The sample is mixed at high speed for transformation into a homogeneous suspension. The standard number of mixing cycles is 4.

FILTERING

The finished suspension is fed into the applied solution tube under pressure, through a piston mechanism (pic. 4), the standard filtration time is 4 seconds.

Picture 4. Sample filtering

6. SYMBOLS ON PACKAGE

roup package

	Ref number
LOT :	Lot number
EXP :	Expiry date
QTY :	Quantity
	manufacturer
	Authorized representative in EC
	MD for in vitro diagnostics
	CE marking
	For single use
	Protect from moisture
	Do not heat

transport package

	Fragile, handle with care
	Up
	Protect from moisture
	Do not heat
	Manufacturer
	Authorized representative in EC
	MD for in vitro diagnostics
	CE marking
LOT :	Lot
EXP :	Expity date
QTY :	Quantity

**7. GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS,
CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.**

Parameters of Disposable plastic tube for filtrate

Overall dimensions, mm	Diameter – 14 ± 1 Flange (collar) diameter 19 ± 1 Height - 60 ± 1.4
Nominal volume (until flange), ml, not more than	5
Weight, g,	2.5 ± 0.5

8. CONDITIONS FOR TRANSPORTATION, STORAGE AND OPERATION. SHELF LIFE

The product is transported by all types of closed vehicles, in accordance with the rules of transportation of goods; acting on this type of transport, during transportation should be protected from strong shocks, rain and direct sunlight.

Transportation and storage

Temperature	$10-55^{\circ}\text{C}$
Atmospheric pressure	$50,0\text{ kPa} - 106,0\text{kPa}$
Relative humidity	Must not exceed 85%

Operation

Temperature	$15-30^{\circ}\text{C}$
Atmospheric pressure	$86,0\text{ kPa} - 106,0\text{kPa}$
Relative humidity	Must not exceed 70%

shelf life 24 months.

9. WARRANTY OBLIGATION

All questions related to the operation and quality of the product can be addressed to the Authorized representative of the manufacturer in the territory of the Russian Federation.

The manufacturer guarantees that the product meets the requirements of the standards and documentation of the manufacturer under the conditions of transportation, storage, operation.

The manufacturer shall not be liable for any damage or damage caused by the use of the product until proven guilty.

The manufacturer shall not be liable for any damage or damage caused by the use of the product after the expiry date indicated on the packaging.

For questions about the quality of a medical product, manufactured by Xiamen Sciendox Biological Technology Co., Ltd, China, please contact:

REAL CARE LLC

111524, Moscow, Perovskaya street, house 1, building 10, floor 2, room VII room 3

Telephones: +7(495)234-57-31

10. POTENTIAL USERS

Medical laboratory technician (laboratory assistant) or other specialist.

11. REQUIREMENTS FOR THE MAINTENANCE AND REPAIR OF MEDICAL DEVICES

Not applicable.

12. REQUIREMENTS FOR THE THE DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES.

Unused expired products can be disposed of in accordance with medical waste disposal requirements.

In case of clinic, warehouse, etc. medical device is utilized/destroyed in accordance with medical waste circulation rules of SanPin 2.1.7.2790-10 dated 12.12.2010 «Sanitation and epidemiological requirements for medical waste circulation»

Medical waste class B.

Adhere to national regulations when disposing or recycling the product, its components and its packaging.

**13. THE NAME AND LEGAL ADDRESS OF THE MANUFACTURER OF
MEDICAL DEVICES AND AUTHORIZED REPRESENTATIVE OF THE
MANUFACTURER IN RF:**

MANUFACTURER:

Xiamen Sciendox Biological Technology Co., Ltd, China

Room 402, Building 3, No. 2054 of Wengjiao West Road, Haicang District 361006 Xiamen,
Fujian, People's Republic of China

Manufacturing site:

Xiamen Sciendox Biological Technology Co., Ltd, China

Room 402, Building 3, №2054 of Wengjiao West Road, Haicang District 361006 Xiamen,
Fujian, People's Republic of China

14. AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:

REAL CARE LLC

111524, Moscow, Perovskaya street, house 1, building 10, floor 2, room VII room 3

Telephones: +7(495)234-57-31

**TOTAL NUMBERED,
BOUND AND SEALED
9PAGES**

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