Ovum Pickup Needle

Instructions for use for the medical device: Ovum Pickup Needle

APPROVED


Allwin Medical Devices Inc., USA

Business Development Manager

(position)

Digish Mehta

(name)

(signature)

Stamp

July, 30th, 2024

Allwin Medical Devices, Inc.

3305 East Miraloma Avenue, Suite 176, Anaheim, CA 92806, USA

Tel.: +1 714-572-1709 • info@allwinmedical.com • www.allwinmedical.com

Ovum Pickup Needle

Ovum Pickup Needle

Updated: 26.06.2024

Intended Use

Ovum Pickup Needle is used for retrieval and flushing of oocytes in ovarian follicles in the assisted reproductive technologies programs.

Intended user

The product should be used by professionals in medical laboratories specializing in the field of assisted reproductive technologies

Application area

The medical device is used in the assisted reproductive technologies programs in specialized laboratories.

Configuration and characteristics :

1) ACE Single Lumen Ovum Pickup Needle

REF	Needle Outer Diameter G/mm (± 0.05 mm)	Needle Length, mm (± 2 mm)	Aspiration Line, mm (± 10 mm)	Vacuum Line, mm (± 6 mm)
A-OVM-SL16/35	16 / 1,60	350	990	395
A-OVM-SL17/35	17G / 1,40	350	990	395
A-OVM-SL18/35	18G / 1,30	350	990	395

2) ACE-FL Single Lumen Ovum Pickup Needle with Flushing Line

REF	Needle Outer Diameter G/mm (± 0.05 mm)	Needle Length, mm (± 2 mm)	Aspiration Line, mm (± 10 mm)	Vacuum Line, mm (± 6 mm)	Flushing Line, mm (± 5 mm)
A-OVM-SL16/35-FL	16/1,60	350	990	395	130
A-OVM-SL17/35-FL	17/1,40	350	990	395	130
A-OVM-SL18/35-FL	18 / 1,30	350	990	395	130

3) OVUMPICK Double Lumen Ovum Pickup Needle

REF	Needle Outer Diameter G/mm (± 0.05 mm)	Needle Length, mm (± 2 mm)	Aspiration Line, mm (± 10 mm)	Vacuum Line, mm (± 6 mm)	Flushing Line, mm (± 6 mm)
OVM-DL16/35	16 / 1,60	350	1000	395	750
OVM-DL17/35	17 / 1,50	350	1000	395	750

Package contents: 10 individually packaged sterile disposable products of one configuration in a box, instructions - 1 pc.

Quality control

One Cell MEA $\geq 80\%$ embryos expanded to blastocyst at 96 hours
Bacterial Endotoxins (LAL) ≤ 20 EU/device.

Materials not included with the device

Vacuum pump, sterile syringe, oocyte collection tube, transvaginal ultrasound transducer, laparoscopic equipment (if needed)

Indications for Use

Retrieval and flushing of oocytes in ovarian follicles in the assisted reproductive technologies programs

Contraindications:

Ovum pick-up needle should not be used:

- In the presence of or after recent pelvic inflammatory disease
- For intra-fallopian procedures
- In the presence of active vaginal or intrauterine infection
- In the presence of sexually transmitted disease
- In case of recent uterine perforation

Possible complications:

- Hematuria
- Extravasation of urine
- Infection
- Vaginal bleeding

Note:

Ensure that the ultrasound unit has good resolution for better visibility.

Warning/ Precaution:

- Do not use if package is opened or damaged.
- Test the system to ensure all connections are sound and flow can be achieved. Failure to test the system may result in failure to obtain the required sample and/or tissue damage.
- Failure to use ultrasound to locate the follicles may result in failure to obtain the required sample and/or tissue damage.
- Dispose in accordance with recognized medical practice and under observance of applicable law and regulation.
- This device is supplied sterile and intended for single use only.
- Attempt to sterilize a single use device may compromise its structural integrity and/or lead to device failure which may result in patient injury or illness.
- Only experienced qualified specialist should use this device.
- The procedure should be carried out under aseptic conditions. Use sterile surgical clothing (gloves, robe, mask).

Method and operation procedure:

Ultrasound Guided Procedure

1. The procedure should be performed under either general anesthesia or analgesia.
2. Position the patient in the lithotomy position. Introduce the ultrasound transducer into the vagina and scan the pelvis to establish the position of the uterus, the development of the endometrium and the number and position of follicles to be aspirated.
3. Remove the Ovum pickup needle from its packaging and insert the needle into the needle guide of the ultrasound transducer, ensuring that there is smooth movement along its length. Ensure no kinks in the tubing.
4. Connect vacuum pump to the connector on the vacuum line and a sterile collection tube to the silicone bung designed to fit 15ml Falcon type tube. The vacuum pressure required varies with the gauge and the type of needle. The vacuum pressure selection is at the discretion of the performing clinician but should not exceed - 500mm HG.
5. If Flushing is required :
 - a) For configuration 2 – Connect 10ml syringe containing flushing media to the connector of flushing line, and flush the follicle. Then conduct the aspiration of the follicle.
 - b) For configuration 3 – To activate flushing line, turn the 3 way stopcock valve and attach a syringe containing flushing media to 3rd face of stopcock and perform flushing.
6. Insert the ultrasound transducer with the needle into vagina and direct the probe towards the posterolateral aspect of the vaginal fornix. Manipulate the probe until the follicles appear in close proximity to the needle.
7. Advance the needle into the nearest follicle. If the ovaries are mobile apply pressure suprapubically and enter the follicles with a stabbing movement. When the tip of the needle can be seen in the follicle, aspirate the contents until empty.
8. Repeat for all follicles within the ovary. To minimize trauma and possible bleeding keep the needle tip within the ovary between aspirations.
9. Remove the collection tube from the silicon bung when full. Replace it with a new sterile collection tube.
10. When all the follicles in one ovary have been aspirated, repeat the procedure with the other ovary.
11. At the end of the procedure ensure to rinse the needle with flushing solution and discard as per standard procedure.

Laparoscopic Procedure

1. Position the patient in the lithotomy position. Follow steps 3 as previously described.
2. Connect a vacuum unit or a syringe as described in steps 4-5 above.
3. Using standard laparoscopic procedure and the necessary laparoscopic equipment advance the needle into the follicle.
4. Aspirate and/or flush the follicle to obtain the oocyte (step 5 above)
5. When the desired number of oocytes has been aspirated, remove, flush and discard the needle.

Do not reuse, reprocess & resterilize:

Reuse, Reprocessing & resterilization may compromise the structural integrity, can also create risk of contamination & may not give desired

Ovum Pickup Needle

result or create complications, infections which may result in injury, illness or death.

Transportation and storage

- Store between 5-30°C in dry place
- Keep away from sunlight
- Shelf life: 5 years from date of manufacture
- Do not use after expiry date
- Transport stability: maximum 24 hours at temperature -30°C at uncontrolled RH to +60°C at 15% RH (±5% RH)

Procedure of utilization

When a waste classified as class A and B is formed during the use of the product, it should be disposed of according with the requirements of SanPin 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures. Use, expired and unused product is classified as Class B and disposal should happen

according to the requirements of SanPin 2.1.3684-21. Packaging of the product after use of the product belongs to waste of class A and is disposed of with household waste.

Warranty liability

The manufacturer guarantees that the products are in accordance with the requirements declared by the manufacturer, if abide by the rules of use, transportation and storage during the shelf life.

If there are questions regarding of the quality of a medical product, you must contact to the authorized representative: Gera Ltd, Russia, 620026 Ekaterinburg, Karla Marks str., d. 8 of. 416, e-mail: gera.cat96@mail.

Manufacturer:

Allwin Medical Devices Inc. 3305 E Miraloma Ave, Suite 176, Anaheim, CA 92806 (USA)

Tel.: +1 714-572-1709

E-mail: info@allwinmedical.com | website: www.allwinmedical.com

REF	Catalogue Number	LOT	Batch Code	Date of Manufacture	Use By	Do not re-use	Do not Re-Sterilize	Non-toxic	STERILE	EO	Sterilized using ethylene oxide
Consult Instruction for use	Do not use if Packing is damaged	Caution Consult Instruction for use	Keep out of sunlight								
Apyrogenic	Keep Dry	Does not contain natural rubber latex	Temperature Limit								

allwin®
Medical Devices

Allwin Medical Devices, Inc.
3305 E. Miraloma Ave., Suite 176 Anaheim, CA 92806 (USA)
Tel.: +1 714-572-1709
E-mail: info@allwinmedical.com | www.allwinmedical.com

Stickers for the patient's card, warehouse accounting, financial accounting are used when necessary in accordance with the laboratory regulations (they are placed in the patient's card, documents of warehouse or financial accounting, respectively)

ACKNOWLEDGMENT

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California
County of ORANGE

On 9/18/2024 before me, ALEXANDRA S. ZAMBRI, NOTARY PUBLIC
(insert name and title of the officer)

personally appeared Digish Mehta
who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.


I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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