HumaScope Advanced^{LED}

| User Manual

CE

Cat.No. 14600/1

2 2. SEP. 2015

Human

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Dr. Torsten Borchers
Vice President Quality Assurance
and Regulatory Affairs

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Human

REVISION LIST OF THE MANUAL

No.	DATE / Rev.	REVISION DESCRIPTION
01	01/2011-04	First edition
02	02/2015-05	Transfer of DOV data, storage, transport and operation conditions
03	03/2015-09	Correction of typographical errors

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human
Diagnostics Worldwide

1 INTRODUCTION

This manual is considered as a part of the instrument; it has to be at the operator's hand as well as at the maintenance operator's availability. For accurate installation, use and maintenance, please read the following instructions carefully. In order to avoid instrument or personal damages, carefully read the "GENERAL SAFETY WARNINGS", describing the suitable operating procedures. In case of breakdowns or any troubles with the instrument, apply to the local Technical Service.

2 USER WARRANTY

HUMAN warrants that instruments sold by one of its authorised representatives shall be free of any defect in material or workmanship, provided that this warranty shall apply only to defects which become apparent within one year from the date of delivery of the new instrument to the purchaser.

The HUMAN representative shall replace or repair any defective item at no charge, except for transportation expenses to the point of repair.

This warranty excludes the HUMAN representative from liability to replace any item considered as expendable in the course of normal usage, e.g.: lamps, valves, syringes, glassware, fuses, diskettes, tubing etc.

The HUMAN representative shall be relieved of any liability under this warranty if the product is not used in accordance with the manufacturer's instructions, altered in any way not specified by HUMAN, not regularly maintained, used with equipment not approved by HUMAN or used for purposes for which it was not designed.

HUMAN shall be relieved of any obligation under this warranty, unless a completed installation / warranty registration form is received by HUMAN within 15 days of installation of this product.

This warranty does not apply to damages incurred in shipment of goods. Any damage so incurred shall be reported to the freight carrier for settlement or claim.

3 INTENDED USE OF THE INSTRUMENT

The instrument has to be used for the expected purposes and in perfect technical conditions, by qualified personnel, in working conditions and maintenance operations as described in this manual, according to the GENERAL SAFETY WARNINGS. This manual contains instructions for professional qualified operators.

4 GENERAL SAFETY WARNINGS

Use only chemical reagents and accessories specified and supplied by HUMAN and/or mentioned in this manual.

Place the product so that it has proper ventilation.

The instrument should be installed on a stationary flat working surface, free from vibrations.

Do not operate in area with excessive dust.

Work at room temperature and humidity according to the specifications listed in this manual.

Do not operate this instrument with covers and panels removed.

Only use the power cord specified for this product, with the grounding conductor of the power cord connected to earth ground.

Use only the fuse type and rating specified by the manufacturer for this instrument, use of fuses with improper ratings may pose electrical and fire hazards.

To avoid fire or shock hazard, observe all ratings and markings on the instrument.

Do not power the instrument in potentially explosive environment or at risk of fire.

Prior to cleaning and/or maintaining the instrument, switch off the instrument and remove the power cord.

For cleaning use only materials specified in this manual, otherwise parts may become damaged.

It is recommended always to wear protective apparel and eye protection while using this instrument.

Respective warning symbols, if appearing in this manual, should be carefully considered.

5 DISPOSAL

The current
arrange pr
All parts
procedure
carefully c
The instr
the regul
Batteries
off in acc

6 INSTRUMENT

Instrum
be consi
Before c
parts.
decont
person
disinfe
supplie
by the

7 NO

Every
respo
prod
to ch

8 ST

The
The
Dire

Human
Diagnostics Worldwide

5 DISPOSAL MANAGEMENT CONCEPT

The currently valid local regulations governing disposal must be observed. It is in the responsibility of the user to arrange proper disposal of the individual components.

All parts which may comprise potentially infectious materials have to be disinfected by suitable validated procedures (autoclaving, chemical treatment) prior to disposal. Applicable local regulations for disposal have to be carefully observed.

The instruments and electronic accessories (without batteries, power packs etc.) must be disposed off according to the regulations for the disposal of electronic components.

Batteries, power packs and similar power source have to be dismantled from electric/electronic parts and disposed off in accordance with applicable local regulations.

6 INSTRUMENT DISINFECTION

Instruments or parts which may come in contact with biological samples (patient specimens, controls etc.) should be considered at least potentially infectious.

Before doing any servicing on the instrument it is very important to thoroughly disinfect all possibly contaminated parts. Before the instrument is removed from the laboratory for disposal or servicing, it must be decontaminated/disinfected. Decontamination/disinfection should be performed by authorised well trained personnel, observing all necessary safety precautions. Instruments to be returned have to be accompanied by a disinfection certificate completed by the responsible laboratory manager. If a disinfection certificate is not supplied, the returning laboratory will be responsible for charges resulting from non-acceptance of the instrument by the servicing centre, or from authority's interventions.

7 NOTICE

Every effort has been made to avoid errors in text and diagrams, however, HUMAN GmbH assumes no responsibility for any errors which may appear in this publication. It is the policy of HUMAN GmbH to improve products as new techniques and components become available. HUMAN GmbH therefore has to reserve the right to change specifications if necessary in the course of such improvements.

8 STORAGE/TRANSPORTATION

The HumaScope Advanced^{LED} can be stored at -20...+70°C. Intense direct sunlight and humidity should be avoided.

The HumaScope Advanced^{LED} contains no dangerous parts. No restrictions have to be observed for transportation.

Directive on restriction of hazardous substances (RoHS) is met by HumaScope Advanced^{LED}.

CONTENTS

CONTENTS	1
1 NOTICE	3
2 APPLICATIONS	3
3 SPECIFICATIONS	3
4 COMPONENTS	4
5 STRUCTURE	5
5.1 Base	5
5.2 Head	5
5.3 Illumination System	5
5.3.1 Koehler Illumination	6
5.4 Mechanical Stage (Slide Holder)	8
5.5 Ocular, Objective and Nosepiece	8
5.6 Coarse and Fine Focussing Knobs	9
5.7 Condenser	9
5.8 Immersion Oil	9
6 OPERATION	10
6.1 Operating Conditions	10
7 EXCHANGING THE LAMP & THE BATTERY	10
8 CARE AND MAINTENANCE:	11
9 DISPOSAL MANAGEMENT CONCEPT	12
10 MICROSCOPE GLOSSARY	12
10.1 Darkfield Microscopy	12
10.2 Brightfield	12
10.3 Phase Contrast	12
10.4 Fluorescence Microscopy	13
10.5 Lenses	13
10.6 Focus	13

1 NOTICE

Every effort has been made to avoid errors in this text and the accompanying diagrams; however, HUMAN GmbH assumes no responsibility for any errors which may appear in this publication. It is the policy of HUMAN GmbH to improve products as new technologies and components become available. Human GmbH therefore reserves the right to change specifications if necessary in the course of such improvements.

2 APPLICATIONS

The HumaScope Advanced^{LED} biological microscope is suited for use in medical laboratories (examples: biology, bacteriology, cytology, histology, pharmaceutical chemistry), agricultural research institutes, colleges, universities and many more.

3 SPECIFICATIONS

Observation Tube	Compensation-free Siedentopf-type binocular 360° rotatable head, 30° inclined
Magnification	40x-1600x
Eyepiece	WF-10x/18 & P-16x/11 (optional)
Interpupillary Distance	53-75 mm
Objective Revolver	5-position reverse
Objectives	Plan-achromatic (anti-fungal) 4x/0.1, 10x/0.25, 40x/0.65, 100x/1.25 (oil immersion) 40, 100 spring loaded
Condenser	Movable N.A. 1.25 Abbe condenser with iris diaphragm
Filtering	Filter holder and daylight (blue) filter
Stage	Mechanical; size 125 x 135 mm, range 35x75 mm
Focusing	Coaxial coarse/fine focus, with tension adjustment and up stop
Focus Stroke	Gradation 2µm
Collector	Koehler illumination
Illumination	LED lamp 3.6V, 0.3W. adjustable light intensity

4 COMPONENTS


Figure 1

1 = Power Switch	7 = Ocular Adjustment Angle	13 = Condenser
2 = Stage Course Focussing	8 = Objective Revolver	14 = Diaphragm (Koehler)
3 = Fine Focussing	9 = Objective	
4 = Condenser Up/Down	10 = Specimen Holder	
5 = Binocular Head	11 = Stage	
6 = Ocular (Eyepiece)	12 = Longitudinal Stage Adjust.	

5 STRUCTURE

5.1 Base

The base (figure 2) supports the entire weight of the microscope. It has four rubber feet to stabilise the instrument. The base contains the lamp source and the complete electronic circuitry.

Figure 2

5.2 Head

The binocular, compensation-free, Siedentopf-type ergonomic head (figure 3) can be inclined at angles from 10° to 30° and the eyepieces can be extended up to 40 mm. The Siedentopf head does not require compensating adjustment for different interpupillary distances. This ensures an optimum eye level and comfortable viewing position, regardless of the operator's physique. In addition, a digital eyepiece camera can be attached to the tubes.

Figure 3

5.3 Illumination System

In order to provide sufficient brightness to observe the specimen and make full use of the resolution of the objectives, this instrument uses a 3.6V, 0.3W, built-in variable-brightness lamp. For optimal results, the condenser must be adjusted properly. The condenser is composed of two parts. The iris diaphragm is mounted under the stage; the illumination source is in the base. Adjust the condenser iris (mounted under the stage) to match the aperture of the objective. The stage can be raised and lowered using the focussing knobs. The light axis of the condenser must coincide with the light axis of the illumination source. To correct the alignment, adjust the 3 corresponding screws in the condenser.


Figure 4

5.3.1 Koehler Illumination

The HumaScope Advanced^{LED} is equipped with Koehler illumination (see figure 5). In order to get the best image possible from the optical setup of the microscope for brightfield, phase contrast and differential interference contrast, it is crucial that the light path be set up properly.

How to set up HumaScope Advanced^{LED} properly for Koehler illumination:

If the light path is set up properly, the user has the advantage of an evenly illuminated field, a bright image without glare and minimal heating of the specimen.

The following instructions apply to a microscope equipped for transmitted-light brightfield illumination:

1. Switch on the light source and make sure that light is passing through the field diaphragm at the base of the microscope stand.
2. Place your specimen on the stage and turn the nosepiece (which holds the objective lenses) to the 10x or 20x lens. Open the field diaphragm as far as possible.
3. Check if the specimen is illuminated. It will help to place a piece of paper over the top of the specimen to see if light is getting through. If you are using the brightfield condenser stop, open the iris diaphragm (or aperture diaphragm) on the condenser turret (which holds the stops for brightfield and phase, etc.) to the maximum. The front lens should be about 1-3 mm above the specimen. Use the condenser focussing knobs to do this.
4. Bring the specimen into focus with the coarse and fine focussing knobs. If the light is too bright, reduce the light intensity using the light adjustment wheel.
5. When the specimen is in focus, start to close the field diaphragm and also begin to carefully move the condenser up and down with the condenser focussing knobs. Look for a sharp image of the edge of the field diaphragm.
6. When the edge of the field diaphragm silhouette is sharply defined, centre it with the two knurled knobs that extend diagonally from the condenser. Close down the field diaphragm most or all the way to get it centred properly. When it is centred, open the field diaphragm until its edge is outside the field.
7. Your specimen should now be properly illuminated and you should have a good image.

Figure 5

Short Summary on Koehler adjustment:

1. Focus on the specimen.
2. Close the field diaphragm to its most closed state so that you can see the edges of the diaphragm (may be blurry) in the field of view.
3. Use the condenser focus knobs to bring the edges of the field diaphragm into the best focus possible.
4. Use the condenser-centering screws to center the image of the closed field diaphragm in the field of view.
5. Open the field diaphragm just enough so that its edges are just beyond the field of view.
6. Adjust the condenser diaphragm to introduce the proper amount of contrast into your sample.
7. Adjust the light intensity as necessary.

5.4 Mechanical Stage (Slide Holder)

The HumaScope Advanced^{LED} stage mechanism (See figure 6) allows the operator to move the specimen slide along the X or Y axis by turning the corresponding knob. This is very useful at higher magnifications, as it can be difficult to move the slide in small enough increments by hand. Also, moving the slide by hand can be difficult since it must be moved in the opposite direction of the observed movement. The mechanical stage also has a graduated scale so you can see how far the slide has been moved or keep track of the position of various objects on the slide.

Figure 6

5.5 Ocular, Objective and Nosepiece

The microscope imaging system is composed of two oculars and the inward 5-hole objective revolver (see figure 7). The objectives are screwed into the threaded holes in the nosepiece, which can be rotated to quickly achieve the desired magnification. The advanced, high-precision nosepiece design ensures that the field of view is always centred in the observation area for precise focussing (parfocal). The ocular tube is inclined 45° for convenient, fatigue-free observation. Ocular and objectives have an anti-fungal coating.

Figure 7

5.6 Coarse and Fine Focussing Knobs

The HumaScope Advanced^{LED} features coaxial coarse and fine focussing knobs (Figure 8) on the left and right sides of the instrument. The knobs are located below the stage and are easy to operate. The coarse and fine focussing range is 30mm. The fine focussing knobs are the smaller knobs centred in the coarse focussing knobs. The knob on the left, inside the left-side coarse focussing knob, is the coarse focussing stopper. When the specimen is in clear focus, the stopper may be used to lock in the focus. This also allows the stage to be lowered and raised again to exactly the same position. When the position is no longer needed, simply loosen the stopper.


Figure 8

5.7 Condenser

The Abbe condenser is a specially designed lens mounted under the stage (figure 4). It can be raised or lowered. It has an iris-type aperture to control the diameter of the beam of light entering the lens system. By changing the size of the iris and moving the lens toward or away from the stage, the diameter and focal point of the cone of light that goes through the specimen can be controlled. Abbe condensers become most useful at magnifications above 400x.

5.8 Immersion Oil

Immersion oil is a special oil used in microscopy only with the 100x objective lens (usually at 1000x total power). A drop is placed on the cover slip and the objective is lowered until it just touches the drop. Once brought into focus, the oil acts as a bridge between the glass slide and the glass in the lens. This concentrates the light path and increases the resolution of the image.

6 OPERATION

6.1 Operating Conditions

Working: +15...35°C

Humidity: 20...85% rel. Non condensing

During operations HumaScope Advanced^{LED} must be protected against any vibrations and movements. Elevated temperatures, any direct sunlight and humidity should be avoided.

1. Screw the objectives into the threaded holes in the nosepiece in order of their magnification. Insert the oculars into the inclined ocular tube.
2. Place a specimen slide on the stage. Adjust the stage so that the specimen is centred over the hole in the stage.
3. Using first the 10X objective, turn the coarse focussing knob to move the objective close to the specimen. Now observe the image through ocular and turn the coarse focussing knob to lower the objective until the image comes into approximate focus. Next, turn the fine focussing knob until the image is sharp. Thanks to the instrument's precise and reliable nosepiece, a clear focus is maintained even when switching objectives.


Abbildung 9

Adjustment of interpupillary distance

Place the specimen on the stage and bring the specimen into exact focus. Adjust the interpupillary distance of binocular until the right and left fields of view merge to form a single circle. (figure7)

Adjustment of diopter

Put the specimen on the stage. Turn the 40x objective to working position. First, observe the specimen with only the right eye, using the right ocular. Use the coarse and fine focusing knobs to bring the specimen into exact focus. Next, using only the left eye to observe the specimen through the left ocular, adjust the diopter control (1) until the specimen is in exact focus.

7 EXCHANGING THE LAMP & THE BATTERY


Abbildung 10

- A. Lamp exchange is not necessary as the lifetime of a LED lamp is 10.000 hours.
- B. To change the rechargeable batteries First switch off the power and unplug the instrument.
- C. Turn the microscope on its back and loosen the 4 screws beneath the feet.
- D. Carefully pull the case circuit out and turn around.
- E. Loosen the screw of the plastic box. Pull off the lid and pull out the old batteries
- F. Insert the new batteries (type AA, NiMH 1.2V 1800mAh) into the box.
- G. Close the panel and retighten the screw.
- H. Check batteries for correct function and fix the base circuit and retighen the screws.

8 CARE AND MAINTENANCE:

1. As with any other optical instrument, the microscope should be kept in a cool, dry, dust- and acid-free location. The instrument should be covered with a dust cover after use.
2. Never disassemble the lenses; these are precisely adjusted and glued in place at the factory. Stains can be removed from the lens by wiping it carefully with a soft, clean cloth moistened with alcohol. Take care, however, that alcohol does not seep it into the objectives to avoid dissolving the glue. Dust on the lens can be removed using a clean lens brush.
3. The coarse and fine focussing assembly and nosepiece are precisely adjusted at the factory and should not be dismantled.
4. When not in use, keep the objectives in the objective box and cover the ocular tube with the ocular-tube cover.

9 DISPOSAL MANAGEMENT CONCEPT

The currently valid local regulations governing disposal must be observed. It is in the responsibility of the user to arrange for proper disposal of the individual components.

All parts which have come into contact with potentially infectious materials must be disinfected using a valid procedure (autoclaving, chemical treatment) prior to disposal. The instrument and electronic accessories (without batteries, power packs etc.) must be disposed of in accordance with regulations for the disposal of electronic components.

Batteries, power packs etc. must be separated from electronic parts and disposed of in accordance with applicable local regulations.

10 MICROSCOPE GLOSSARY

10.1 Darkfield Microscopy

Darkfield microscopy is a method of acquiring a slide image. It is specific to the way the light is projected onto a slide in the microscope equipped with a darkfield condenser. The method allows elements in the field of view to be distinguished that would be invisible or poorly visible using other methods of lighting (such as brightfield or phase contrast light). This method is commonly used to evaluate slides on which the object and the background are the same colour (such as red blood cells in vital blood or clear bacterial organisms in saline). In darkfield microscopy, objects appear in the form of glowing outlines on a black/dark background.

The advantage of darkfield illumination is that details can be seen that are normally not resolved by the microscope's objective. While the actual details cannot be seen, the light reflected by them can. A nice analogy here is that of dust in a room. Normally, in a well-lit room very small dust particles cannot be seen. However, if the lights are turned out, a beam of light from an acute angle makes these same particles visible. Besides the optical advantages, darkfield illumination is very beautiful and gives an almost science-fiction-like image.

What is the difference between darkfield and brightfield?

The difference between the two is the way light is projected onto the slide. As a result, elements on the slide not visible (or poorly visible) in bright light are visible with darkfield illumination.

10.2 Brightfield

This is the basic microscope configuration.

In brightfield microscopy, the transparent or translucent specimen is either naturally coloured or stained and appears dark against a bright, white background or field.

10.3 Phase Contrast

Phase contrast is a technique for revealing the structural features of microscopic, transparent objects that cannot otherwise be accomplished with brightfield microscopy. Phase contrast achieves the same effect as staining a specimen (which can kill a live specimen). In phase-contrast microscopy, the "staining" is in the optics.

In a phase-contrast microscope, the annular rings in the objective lens and the condenser split the light. The light that passes through the central part of the light path is recombined with the light that travels around the periphery of the specimen. The interference produced by these two paths produces images in which the dense structures appear darker than the background.


Figure 11

10.4 Fluorescence Microscopy

The fluorescence of a substance is seen when a molecule is exposed to a specific wavelength of light (excitation wavelength or spectrum) and the light it emits (the emission wavelength or spectrum) is of a higher wavelength. To view this fluorescence in the microscope, several light-filtering components are needed. Specific filters are needed to isolate the excitation and emission wavelengths of a fluorochrome. A bright light source with proper wavelengths for excitation is also needed. For normal fluorescence applications, this is a mercury vapour arc burner.


Figure 12

10.5 Lenses

objective lens - gathers light from the specimen

ocular - transmits and magnifies the image from the objective lens to your eye

nosepiece - rotating mount that holds the interchangeable objective lenses

tube - holds the ocular at the proper distance from the objective lens and blocks out stray light

10.6 Focus

to position the objective lens at the proper distance from the specimen

coarse-focus knob - used to bring the object into the focal plane of the objective lens

fine-focus knob - used to make fine adjustments to focus the image

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HumaScope Advanced^{LED}

CONTENTS

1	Title	2
2	Purpose	2
3	Area of Application	2
3.1	Description of the function	2
4	Accessories.....	2
5	Special Notes	2
6	Materials and Tools.....	2
7	Referenced Documents	2
8	Procedure.....	2
9	Specifications	3
9.1	Power Supply	3
10	Primary Packaging Instrument	3
10.1	Secondary Packaging Instrument	3
10.2	Tertiary Packing Instrument	3
11	Storage /Transportation	3
12	Operating Conditions	3
13	Maintenance	3
14	CE Mark	3
15	Quality Control	4
16	Documentation	4
17	Disposal.....	4
18	Warranty.....	4

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Dr. Forsten Borchers
Vice President Quality Assurance
and Regulatory Affairs

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1/5

1 Title

Design output und product specification for the HumaScope Advanced^{LED}.

2 Purpose

To specify the HumaScope Advanced in it's design and function.

3 Area of Application

Description	Cat. No.	Content
HumaScope Advanced ^{LED}	14600	• 1 X Rotary head (type Siedentopf) • 4 X Plan achromatic Objectives 4x, 10x, 40x, 100(oil)x • 2 X Ocular 10x /18mm • 1 X Immersion oil (1 Bottle) • 1 X Power supply • 1 X Dust cover Human • 1 X User Manual

3.1 Description of the function

HumaScope Advanced^{LED} is a binocular universal microscope designed for general use in laboratories, schools, hematology, dermatology, urinalysis and histology.

This microscope can be optionally equipped with phase contrast, a dark field condenser, simple polarisation kit, and other accessories as listed in the HumaScope accessory catalogue.

4 Accessories

According to the accessory list, separately available from HUMAN.

5 Special Notes

Not applicable

6 Materials and Tools

7 Referenced Documents

Operating instructions (User Manual), conformity declaration.

8 Procedure

HumaScope Advanced^{LED} has to be unpacked, installed and operated strictly following the instructions in the HUMAN User Manual/Operating Instruction.

9 Specifications

HumaScope Advanced^{LED}

Observation Tube	Compensation-Free Siedentopf-type Binocular 360° rotatable Head, 30° inclined
Magnification	40x-1600x
Eyepiece	WF-10x/18 & P-16x/11 (Optional)
Interpupillary Distance	53-75 mm
Objective Revolver	5 Position reverse
Objectives	Plan-achromatic (anti-fungal) 4x/0.1, 10x/0.25, 40x/0.65, 100x/1.25 (immersion oil) 40,100 spring loaded
Condenser	Movable N.A. 1.25 Abbe condenser with iris diaphragm
Filtering	Filter holder, no filters included
Stage	Mechanical size 125 x 135 mm, range, 35 x 75 mm
Focusing	Coaxial coarse/fine focus, with tension adjustable and up stop
Focus Stroke	Gradation 2µm
Collector	Koehler illumination
Illumination	LED lamp 3.6V 0.3W, adjustable light intensity, 3 rechargeable batteries Ni-Mh AA 1800mAh 1.2V

9.1 Power Supply

The instrument is equipped with a power supply, suitable for 110 to 230 VAC.

10 Primary Packaging Instrument

The instrument including its accessories is packed in a neutral styrofoam box.

10.1 Secondary Packaging Instrument

The neutral styrofoam box is packed into a neutral carton box (inner box).

10.2 Tertiary Packing Instrument

The inner box is packed into an outer carton, appropriately labelled.

11 Storage /Transportation

The HumaScope Advanced^{LED} can be stored at -20 ... +70°C. Intense direct sunlight and humidity should be avoided. The HumaScope Advanced^{LED} contain no dangerous parts. No restrictions have to be observed for transportation. Directive on Restriction of Hazardous Substances (RoHS) is met by HumaScope Advanced^{LED}.

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13 Maintenance

For cleaning the body and movable parts common non-corrosive and non-abrasive laboratory cleaning liquid or disinfecting solvent can be used.

For the lenses and optics the use of lens cleaning solutions is recommended.

14 CE Mark

HumaScope Advanced^{LED} bear the CE mark in compliance with directives 2004/108/EC, 2006/95/EC and RoHS 2002/95/EC and applicable standards.

15 Quality Control

Quality control of each instrument guarantees completeness, correct labelling and function.

16 Documentation

Release documents are prepared and archived in accordance with applicable regulations and internal quality management system requirements.

17 Disposal

For disposal the currently valid local regulations governing disposal of electric and electronic equipment must be observed. The disposal should be done in accordance with regulations of electric and electronic equipment disposal 2002/96/EC.

Please contact our local representative about disposal of the equipment.

18 Warranty

The equipment has 1 year warranty.

CERTIFICATE OF ANALYSIS

CERTIFICATE OF QUALITY AND WARRANTY

HumaScope Advanced^{LED}

Binocular laboratory LED illuminated microscope

Cat. No.:	14600
Serial No.:	LOTNUMMER
Invent. date:	PRODDATUM

Storage Conditions:

Store at ambient conditions!

QUALITY SPECIFICATIONS		
1. Contents:	Complete acc. pack list	Within Specification

This product passed the quality control release criteria and has been released by the responsible Quality Control Department of HUMAN GmbH.

Instruments have 1 year warranty.

Following release of the product, this certificate of analysis is electronically generated and valid without signature.

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