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User's Manual

“Thermal Paper for Video Printers in rolls, model Super ULSTAR-1100HD, ULSTAR-1100S, ULSTAR-1100HG, ULSTAR-2100HD”

Manual No.: MF-1-01 (Rev.0)

Prepared by: H.S. KIM

on April 20, 2013

Approved by: J.W. Kim

on April 25, 2013

Durico C&T, Inc.

1. Product Information:

Thermal Paper for Video Printers in rolls, model Super ULSTAR-1100HD, ULSTAR-1100S, ULSTAR-1100HG, ULSTAR-2100HD (further in the text Thermal Papers) are used for printing Video Images on the Diagnostic Ultrasound Device Monitor with Thermal Type Video Printers. Medical Image Printing requires high precision available from thermal printing method. Thermal printer uses heat to transfer an impression or image onto Thermal papers that are specially coated to make the paper surface sensitive to thermal impulses.

As a kind of Medical Recording Papers, Printed Images are kept in patient's record file before or after being reviewed by doctors.

2. Directions for Use:

a. Storing Unprinted Paper

- 1) Store unused paper in a cool and dark place (below 30°C or 86°F).
- 2) Do not store in a hot or humid place.
- 3) Do not leave unused paper in direct sunlight or other bright places.
- 4) Do not allow any volatile organic solvent or vinyl chloride.
- 5) Do not touch unused paper with finger prints.

b. Using Paper for printing in Thermal Printer

- 1) Check and confirm the width of paper size as 110mm for use in your printer.
- 2) Select the Paper Type Switch in your Printer to the type of Durico Paper.
■ ex.) For Sony Printers, Select TYPE I Switch for ULSTAR-1100 S, and TYPE II Switch for ULSTAR-1100 HD
- 3) Tear Plastic Bag to get Paper Roll, and Remove the label attached to the top of the paper roll.

- 4) Open Paper Loading Bin to load the Paper roll into Printer
- 5) Pull out the first 15 to 20 cm (6 to 8 inches). These first parts may have some dust, and cause poor print quality.
- 6) Do not fold the paper or touch the printing surface. Dust or finger prints may cause poor print quality.

c. Storing Printed Papers

- 1) Keep the paper in a dark, cool and dry place Otherwise the print may fade out gradually.
- 2) Do not apply the followings to the surface of the printed paper:
 - Plastic tapes
 - Plastic erasers
 - Copy toner
 - Vinyl
 - Marker ink
 - Plasticizer such as desk mat.
- 3) Do not put volatile organic solvent, such as alcohol, on the surface of the printed paper.
- 4) Do not put a fingerprint on the surface of the printed paper.
- 5) Do not rub the printed surface hard.
- 6) Do not cover the printed paper with a commercially available laminate. The lamination may cause the print dyes to smudge.
- 7) It is recommended to store each printed sheet face up or to stack with a piece of paper between each sheet.

3. Shelf Life of the product:

- a. Condition:
 - Under Normal Room Condition (Below 30 °C (86F) /80 % RH)
- b. Recommended Maximum Shelf Life:
 - Sealed product - Approx. 3 years
 - Printed paper -Approx. 5 years

4. Cautions:

- a. This Product is for use in Thermal Printers only.
- b. This Product is not, in itself, medical treatment device.
- c. Do not reuse any part of this Product after printing.
- d. Do not use this Product if the plastic bag is torn, or any unused paper is blackened.

5. Maintenance and repair

- The device is for single use only. Maintenance and repair are not applicable.

6. Warranty

- a. Durico warrants that the Products shall be free from defects in material and workmanship.
- b. If any of the Products or any part of a Product is found not to be in compliance with the quality standards, Durico shall supply Distributor free of charge replacement for the products or the part of a Product not meeting the quality standards.

- 7. In case of any problems, please contact the company LLC "Medonika" 115230, Moscow, Hlebozavodskaya passage, 7, bldg 10, office 406 Tel: +7 (499) 800-30-03

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Joowan Kim President of
Durico C&T Inc.

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appeared before me and
admitted said principal's
subscription to the attached

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This is hereby attested on
this 08th day of May, 2015
at this office

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«I certify accuracy, correctness and
reliability of this document»

Durico C&T, Inc., South Korea

Name Company, Country

Joowan Kim

President of the Durico C&T, Inc.

Last name, first name and title of the person authorized to sign such documents

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1. Location of the manufacturing site

67-17, Oedab-dong, Sangju-si, Gyeongsangbuk-do, 742-320, Korea

2. Utilization

When disposing of the thermal paper Super Ulstar follow the instructions and the provisions of national and local laws and regulations.

3. Indications for use

Thermal paper Super Ulstar is used together with the video printer of the ultrasound center, and does not affect the examination, diagnosis and treatment.

4. Contraindications

Contraindications to the use of thermal paper are not revealed.

5. Transportation and operation conditions

Transportation and operation conditions	
Temperature	- 30 °C to + 30 °C
Relative Humidity	10% ~ 80%
Pressure	700 – 1060 hPa

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**«I certify accuracy, correctness and
reliability of this document»**

Durico C&T, Inc., South Korea

Joowan Kim

President of the Durico C&T, Inc.

Additions to the User's Manual

"Thermal Paper for video printers in rolls, model Super ULSTAR-1100HD, ULSTAR-1100S, ULSTAR-1100HG, ULSTAR-2100HD"

☒ Characteristics

1. Size

Models	Size
Super ULSTAR-1100HD	110mm(Width) x 20meter (Length)
Super ULSTAR -1100S	110mm(Width) x 20meter (Length)
Super ULSTAR -1100HG	110mm(Width) x 18meter (Length)
Super ULSTAR -2100HD	210mm(Width) x 25meter (Length)

2. Properties

Below data is the specification (technical standard) of our product's properties for each grade. "Min" stands for minimum, which means our product's quality will at least meet this standard. The number of "Additions to the Technical Data" is the representative figure of each grade, which is the average value of test result by random sampling.

1) Super ULSTAR-1100HD

Item	Specification		Evaluation Method
Physical Properties	Material	Polypropylene	
	Thickness	80 ± 5 μm	Micro-gauge
	Product Size	Width 110±0.1 mm Length 20±0.1 m	Roll unit
	Whiteness	80.0% Min	ASTM E 313 Spectro eye
	Maximum Optical Density (Dmax)	1.60 Min	DIN 16536 Spectro eye
	Gloss	80.0% Min	ASTM D 523 Incident angle ; 85°
Preservation Abilities	Thermal resistance	90% Min	At 50 °C, No Color Change during 1 day (Preservability of Max. Optical Density)

	Humidity Resistance	90% Min	At 40 °C, 90%RH, 26 hours (Preservability of Max. Optical Density)
	Sunlight Resistance	90% Min	Below sunlight, 28 °C, 50%RH, 1 day, (Preservability of Max. Optical Density)

2) Super ULSTAR-1100S

Item	Specification		Evaluation Method
Physical Properties	Material	Polypropylene	
	Thickness	80 ±5 μm	Micro-gauge
	Product Size	Width 110±0.1 mm Length 20±0.1 m	Roll unit
	Whiteness	85.0% Min	ASTM E 313 Spectro eye
	Maximum Optical Density (Dmax)	1.30 Min	DIN 16536 Spectro eye
	Gloss	60.0% Min	ASTM D 523 Incident angle ; 85°
Preservation Abilities	Thermal Resistance	90% Min	At 50 °C, No Color Change during 1 day (Preservability of Max. Optical Density)
	Humidity Resistance	90% Min	At 40 °C, 90%RH, 26 hours (Preservability of Max. Optical Density)
	Sunlight Resistance	90% Min	Below sunlight, 28 °C, 50%RH, 1 day (Preservability of Max. Optical Density)

3) Super ULSTAR-1100HG

Item	Specification		Evaluation Method
Physical Properties	Material	Polypropylene	

	Thickness	$80 \pm 5 \mu\text{m}$	Micro-gauge
	Product Size	Width $110 \pm 0.1 \text{ mm}$ Length $18 \pm 0.1 \text{ m}$	Roll unit
	Whiteness	85.0% Min	ASTM E 313 Spectro eye
	Maximum Optical Density (Dmax)	1.70 Min	DIN 16536 Spectro eye
	Gloss	88.0 % Min	ASTM D 523 Incident angle ; 85°
Preservation Abilities	Thermal resistance	90% Min	At 50°C , No Color Change during 1 day (Preservability of Max. Optical Density)
	Humidity Resistance	90% Min	At 40°C , 90%RH, 26 hours (Preservability of Max. Optical Density)
	Sunlight Resistance	90% Min	Below sunlight, 28°C , 50%RH, 1day, (Preservability of Max. Optical Density)

4) Super ULSTAR-2100HD

Item	Specification		Evaluation Method
Physical Properties	Material	Polypropylene	
	Thickness	$80 \pm 5 \mu\text{m}$	Micro-gauge
	Product Size	Width $210 \pm 0.1 \text{ mm}$ Length $25 \pm 0.1 \text{ m}$	Roll unit
	Whiteness	80.0% Min	ASTM E 313 Spectro eye
	Maximum Optical Density (Dmax)	1.60 Min	DIN 16536 Spectro eye
	Gloss	80.0% Min	ASTM D 523 Incident angle ; 85°
Preservation Abilities	Thermal resistance	90% Min	At 50°C , No Color Change during 1 day (Preservability of Max. Optical Density)

	Humidity Resistance	90% Min	At 40 °C, 90%RH, 26 hours (Preservability of Max. Optical Density)
	Sunlight Resistance	90% Min	Below sunlight, 28 °C, 50%RH, 1day, (Preservability of Max. Optical Density)


Construction of Medical Device

1. Design Drawing

The structure of thermal paper for video printer

- The thermal paper composes of two layers on plastic substrate (polypropylene).

- 1) Thermal ink layer
- 2) Protect Layer



-  : Protect Layer
-  : Thermal Ink Layer
-  : Substrate (YUPO)

The mechanism of coloration of thermal paper

- The main materials of thermal ink are dye called color former, developer and binder. And some additives are added for coating process and printability. The color forming mechanism is as follows:

Materials don't react at room temperature, but with the heat of printer head, among the components of dye, developer and heat sensitive agent (assistance of melting), the lowest melting point material of heat sensitive agent first melts. And in company with melted liquid, dye and developer melt each. The melting point and compatibility of components are important properties.

The reaction of melted dye and developer is acid-base reaction. The color former (dye) is electron donor and the latter is acceptor. With the transfer of electron the ring of dye opens, and it occurs color change. The reaction is irreversible.

Color former (Dye, colorless) Developer Color formation (ring opening)

Because almost heat energy use to melt of materials, with the high efficiency energy it is able to get the high speed and sensitivity.

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