11

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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of MIDEA BIOMEDICAL CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

«APPROVE» / “APPROVE”

Midea Biomedical Co., Ltd.

General Director

(position/position)

Tony Lu

(name/name)

(signature/signature)

M.P. / Stamp

User Manual for medical device

Medical Freezer

1. Description of the medical device

1.1. Name of the medical device

Medical Freezer

Execution options:

1. Medical Freezer BD-86L108;
 2. Medical Freezer BD-86L358;
 3. Medical Freezer BD-86L458;
 4. Medical freezer BD-86L566;
 5. Medical freezer BD-86L568;
 6. Medical freezer BD-86L718;
 7. Medical freezer BD-86L838.
- Operating manual.

Accessories:

1. Stainless steel shelves – from 0 to 100 pcs. (as required);
2. Cryoboxes 5x5 – from 0 to 100 pcs. (as required);
3. Cryoboxes 9x9 – from 0 to 100 pcs. (as required);
4. Cryoboxes 10x10– from 0 to 100 units (as required);
5. Racks for cryoboxes 5x5– from 0 to 100 units (as required);
6. Racks for cryoboxes 9x9– from 0 to 100 units (as required);
7. Racks for cryoboxes 10x10– from 0 to 100 units (as required);

(hereinafter referred to as the device, chamber, or freezer)

1.2. Purpose of the medical device

Intended for use in clinical diagnostic laboratories for the preservation and storage of materials at temperatures below the freezing point of water (biological material, blood, medicines, live viral vaccines, laboratory products, samples, cultures, reagents, and other laboratory materials)

1.3. Scope of Application of the Medical Device

For use in medical facilities, vaccination centers, medical universities, electronic and chemical industry enterprises, as well as other enterprises for the preservation and storage of materials at temperatures below the freezing point of water (biological material, blood, medicines, live viral vaccines, laboratory products, samples, cultures, reagents, and other laboratory materials)

1.4. Operating Conditions for the Medical Device

Ambient air temperature: 10~32 °C; ideal temperature: 18~25 °C; use an air conditioning system if necessary.

Relative humidity: less than 80%. At a maximum operating temperature of 32 °C, the relative humidity must not exceed 60%.

It is not recommended to install the freezer in a room with high dust concentration in the air.

Do not allow mechanical swinging or excessive vibration to affect the freezer.

Installation height above sea level: no more than 2000 m.

Input voltage: 220V~240 V, $\pm 10\%$.

Overvoltage category: II, Capacity: 5 kVA, IPX0.

Pollution degree: 2.

The low-temperature freezer is sensitive to the effects of ambient temperature. If the environmental conditions at the installation site differ, the low-temperature freezer may operate incorrectly. Before using the freezer, ensure that the installation site conditions meet the required specifications.

Do not install the low-temperature freezer outside of a room. Water entering the freezer can cause electrical leakage or electric shock.

1.5. Indications for Use

Preservation and storage of materials at temperatures below the freezing point of water (biological material, blood, medicines, live viral vaccines, laboratory products, samples, cultures, reagents, and other laboratory materials)

1.6. Contraindications

-

1.7. Warnings

- If the low-temperature freezer is not operating correctly, disconnect it from the power supply. Failure to comply with this requirement may result in electric shock or fire.

WARNING: If the low-temperature freezer is to be stored for a long period, take reasonable steps to restrict access by unauthorized persons.

Disposal of the freezer should be carried out by qualified professionals. When disposing of the door, it should be removed to prevent accidents related to suffocation.

WARNING: If a CO₂/LN₂ backup cooling system is used, proper ventilation must be ensured at the installation site.

- Power supply voltage: The units operate from an alternating current power supply of 220–240 V at a frequency of 50/60 Hz. If the line voltage is below 198 V or above 264 V, a suitable voltage stabilizer (5 kVA) must be used.

- The power supply circuit to which the freezer is connected must include a low-voltage automatic circuit breaker and a ground fault protection device.

- To connect the freezer to the power supply, use a separate socket with a reliable ground connection without using an extension cord. If it is necessary to use an extension cord, the extension cord must have copper conductors with a cross-sectional area of at least 2.5 mm². In this case, the wires inside the walls must have copper conductors with a cross-sectional area of at least 4 mm².

- It is strictly prohibited to place flammable, explosive substances, strong corrosive acids, alkalis, and other hazardous materials in the freezer.

- To prevent accidents, store the freezer lock key in a place inaccessible to unauthorized persons.

- It is forbidden to connect the neutral wire (terminal N) and the grounding wire (terminal E) of the power socket, as this can result in voltage being applied to the freezer's housing and electric shock.

- Do not bend the power cord, press it with heavy objects, or place it near heat sources, such as compressors.

- If it is necessary to restart the freezer after disconnecting it from the power supply or a power failure, first check the freezer settings. Changes in settings can lead to the spoilage of stored materials.

- If the freezer is disconnected from the power supply or a power failure occurs, wait at least 5 minutes before restarting to avoid damage to the compressor or the entire system.

- Timely check and clean the mesh filter. Dusty mesh filters can lead to increased temperature or malfunction of the freezer.

- Timely clean the power cord plug from dust. Incorrectly connected or dusty power cord plug can lead to overheating or fire.
- After a malfunction or power outage of the medical freezer, check the settings when turning the power switch back on; otherwise, stored materials may be damaged due to changes in settings.
- During technical maintenance, wear gloves to avoid injuries from sharp edges or corners.
- Do not touch stored materials with bare hands. Direct contact with stored materials or the inner walls of the freezer can cause frostbite.
- When closing the door, exercise caution and hold the handle to avoid pinching your fingers.
- When moving the low-temperature freezer, do not tilt it at an angle greater than 45°.
- Exercise caution when moving the freezer to avoid damaging property or causing injuries.
- Do not use the door handle to lift or move the freezer, as this can result in damage to the freezer or injuries. Do not use electrical freezers inside the freezer's chamber, except those recommended by the manufacturer.

Precautions for using the product:

- During the operation of the low-temperature freezer, the rear walls of the freezer may heat up, which is a normal phenomenon. To prevent condensation around the freezer's body, a heated condensate drain must be installed.
- Before placing materials in the chamber, ensure that the temperature in the chamber has reached the set point, then place the materials in batches.
- To avoid excessive temperature rise in the working chamber, do not place materials in a volume exceeding 1/3 of the working volume of the chamber.
- The temperature displayed on the chamber's display corresponds to the temperature read by the sensor in the working chamber. At the beginning of the freezer's operation, the set and actual temperatures in the working chamber may differ significantly; however, as the temperature inside the chamber stabilizes, the displayed temperature value will gradually approach the actual temperature.
- The working chamber of the unit has an opening for inserting the tester into the chamber for measurements. After removing the tester, it is necessary to reseal the control hole with insulating material; otherwise, the temperature in the freezer may not reach the set point and may cause condensation outside the hole.
- For cleaning the item, use diluted neutral detergent. Do not use brushes, acidic agents, gasoline, soap powder, polishing agents, or hot water to clean the item, as this can lead to damage to the painted surface and plastic/rubber parts. Do not use volatile solvents, such as gasoline, to wipe plastic and rubber parts.
- If the freezer is not used for a long time, disconnect it from the power supply and turn off the battery switch.
- If the freezer does not turn on after a long period of inactivity, the battery may be discharged. After switching the battery switch to the 'ON' position, it will begin to charge. Leave the freezer running for a week to fully charge the battery (applicable to models with a battery).
- Open the door for minimal intervals to avoid significant temperature and humidity fluctuations in the chamber. After opening the door, the temperature inside the chamber temporarily rises sharply, which is a normal condition, and it will return to normal within an hour after the door is closed.
- After operating the low-temperature freezer for some time, a layer of frost will form on the inner walls and the evaporator. If the frost layer is too thick, it negatively affects the insulation of the low-temperature freezer and increases energy consumption. If the frost layer thickness exceeds 5 mm, remove it using the included scraper.
- Before defrosting the freezer, remove stored materials from the chamber and place them in a suitable storage environment; otherwise, this may lead to material damage due to increased temperature in the chamber.
- Behind the side and rear walls of the chamber, several cooling coils are located. Do not remove ice from the inner walls of the chamber with sharp tools such as a knife, chisel, or screwdriver.

Exercise caution to avoid scratching the inner walls during defrosting; otherwise, the low-temperature freezer may malfunction.

1.8. Safety Instructions

CAUTION! Before using the product, remove the lower part of the packaging from the chamber.
CAUTION! When moving the unit, ensure that the feet of the chamber do not touch the power cord.
CAUTION! After moving the unit to the desired location, lock the brakes on the front lower wheels.

CAUTION! Before the first use, leave the chamber in the installation location for 24 hours.
CAUTION! After the first activation of the chamber, do not place materials inside immediately.

Wait until the temperature inside the chamber reaches the set point, then load materials in batches.
Note: the volume of materials placed in the low-temperature freezer at one time should not exceed 1/3 of the chamber volume.

CAUTION! Ensure that all ventilation holes on the cabinet or in the built-in structure are open.
CAUTION! Do not use mechanical devices or other means not recommended by the manufacturer to speed up the defrosting process.

CAUTION! Be careful not to damage the refrigerant circuit.

CAUTION! connect the freezer only to a grounded outlet. If the installation site has only outlets without a grounding contact, consult a technical specialist to replace them with grounded outlets.

CAUTION! the low-temperature freezer must be installed on a solid and level surface with appropriate measures to prevent tipping.

If the floor is uneven or the freezer is installed in an unsuitable location, it may fall and cause injuries.

CAUTION! In case of a gas leak, including flammable gas, close the valve from which the leak is occurring, open doors and windows to ventilate the room. Do not plug in or unplug the freezer's power cord, as this can cause an explosion or fire.

CAUTION! Disassembly and assembly should only be performed by qualified personnel or authorized manufacturer representatives. Otherwise, this can lead to a fire or electric shock.

CAUTION! Connect the freezer to a power supply with characteristics that correspond to the values indicated on the factory label of the freezer. Otherwise, this may cause a fire or electric shock.

CAUTION! The freezer's power cord is equipped with a three-prong (grounded) plug for connection to a standard three-prong (grounded) outlet rated at 16 A. Under no circumstances should you cut or remove the grounding contact of the power cord. Make sure the plug is securely inserted into the outlet; otherwise, this could cause a fire.

CAUTION! When removing the plug from the outlet, hold the plug, not the cord. Otherwise, there is a risk of electric shock or fire due to a short circuit.

CAUTION! Before performing any repairs or maintenance on the freezer, ensure it is disconnected from power; otherwise, this could result in electric shock or injury.

CAUTION! Ensure that medicines or vials with suspended particles inside or outside the freezer are tightly sealed to avoid inhaling them during repair or maintenance, as this can cause health damage.

CAUTION! When storing toxic, harmful, or radioactive substances, place the low-temperature freezer in a suitable safe location. Improper use can cause harm to human health or the environment.

CAUTION! If the low-temperature freezer is not used for a long period of time, unplug the power cord from the outlet to prevent electric shock, current leakage, or fire due to power cord wear.

CAUTION! Increased concentration of carbon dioxide in the air can lead to suffocation or death. In case of inadequate ventilation, appropriate measures must be taken to ensure the safe operation of the equipment.

CAUTION! Do not stand on the freezer and do not place any objects on it, as this may cause it to tip over, result in personnel injuries, or damage the equipment.

CAUTION! Do not operate the low-temperature freezer outside of a room. Water entering the freezer can cause electrical leakage or electric shock.

CAUTION! Do not install the freezer in rooms with high humidity or in places where water can get onto the freezer. Otherwise, a short circuit or electric shock may occur due to insulation damage.

CAUTION! Do not allow water to spill on the low-temperature freezer, as this can cause a short circuit or electric shock.

CAUTION! Do not disassemble, repair, or modify the low-temperature freezer yourself. Otherwise, this can lead to a fire or injuries.

CAUTION! Do not ground the freezer to gas or water pipes, telephone lines, or lightning rods. Otherwise, this can lead to electric shock or other hazards.

CAUTION! Do not touch the electrical parts of the freezer, including the power cord plug, or press switches with wet hands, as this can lead to electric shock.

CAUTION! Do not place containers with water or heavy objects on the freezer. Falling objects can cause injuries, and spilled water can damage insulation, leading to a short circuit or electric shock.

CAUTION! Do not insert any metal objects, such as nails or wires, into the ventilation or exhaust openings of the freezer. This can result in electric shock or injury from accidental contact with moving parts.

CAUTION! Do not tie, pull, wrap, or twist the power cord. Also avoid damaging the power cord plug. Damaged power cords or plugs can lead to fire or electric shock.

CAUTION! Do not use a power cord with a poorly secured plug, as this can lead to fire or electric shock.

CAUTION! Do not place glass bottles or canned materials in the freezer, as they may freeze and crack, causing injury to personnel.

The refrigerant and polyurethane foam with cyclopentane used in the freezer are highly flammable. Therefore, when disposing of the freezer, it should be stored away from sources of fire and handled by a disposal company with the appropriate qualifications (except in cases of disposal by incineration) to prevent environmental damage or any other damage.

IMPORTANT INFORMATION:

The freezer uses a flammable gas and refrigerant. The user of the product is responsible for complying with federal and local regulations when disposing of this freezer.

It is recommended to cut the refrigerant circuit pipeline and drain it in a place where there are no sources of open flame.

Symbol "Caution!".

Symbol "Caution! Do not place near sources of open flame".

CYCLOPENTANE

Symbol "Cyclopentane in the cooling circuit".

1.9. Risk class and applicable classification rules

Type of medical device according to the nomenclature classification of medical devices	122990
Risk class of the medical device according to the nomenclature classification of medical devices	2a
Code of the All-Russian Classifier of Products by Types of Economic Activity	28.25.13.116

1.10. Developer of the medical device

Midea Biomedical Co., Ltd.

No.5 workshop 88 Yulan Main Road, Hi-Tech Development Zone Hefei, Anhui, China

1.11. Manufacturer of the medical device

Midea Biomedical Co., Ltd.


No.5 workshop 88 Yulan Main Road, Hi-Tech Development Zone Hefei, Anhui, China


production site address:

Midea Biomedical Co., Ltd.



Plot A on the south of Mingchuan Road, Baiyan Science and Technology Park, High-tech Development Zone, Hefei, China

1.12. Composition of the medical device, list of versions

General view of the device (for the following versions: Medical freezer BD-86L458, Medical freezer BD-86L568, Medical freezer BD-86L718, Medical freezer BD-86L838, Medical freezer BD-86L358, Medical freezer BD-86L566)

General view (for the version Medical freezer BD-86L108)

The product supply package includes a Medical freezer of one version as per item 1.1. Operating manual – 1 pc. Accessories as agreed with the customer.

1.13. List of accessories

1. Stainless steel shelves – from 0 to 100 pcs.
2. Cryoboxes 5x5 – from 0 to 100 pcs.
3. Cryoboxes 5x4 – from 0 to 100 pcs.
4. Cryoboxes 10x10 – from 0 to 100 pcs.

5. Racks for cryoboxes 5x5 – from 0 to 100 pcs.

6. Racks for cryoboxes 5x4 – from 0 to 100 pcs.

1.14. General description of main functional elements, parts, composition, and functional capabilities

The chamber is designed for the preservation and storage of materials at temperatures below the freezing point of water (biological material, blood, medicines, live viral vaccines, laboratory products, samples, cultures, reagents, and other laboratory materials, etc.)

Freezer functions

Precise temperature control/ Data traceability

- Computer temperature control with an indication accuracy of 0.1°C
- The low-temperature freezer maintains a uniform temperature in the working chamber and allows for temperature adjustment within the range of -86°C to -40°C with high precision.
- The freezer is equipped with a standard USB port and allows for tracking, viewing, and saving temperature data at any time.
- The compressor is equipped with a dual cooling system.

Distinguishing features of the compressor: power, safety, and reliability

- Highly efficient insulation layer
- Multi-layer sealing construction effectively prevents a decrease in cooling performance and frost formation

Extended alarm configuration capabilities, enhanced safety, and stable operation

- Emergency signal system for exceeding upper and lower temperature limits, power failure, and sensor malfunction; extreme temperature values are set by the user
- The door safety lock design enhances safety when working with samples

Ease of operation


- The adjustable shelf structure provides ample space
- The low-temperature freezer display shows the temperature inside the chamber and the ambient temperature, allows setting high and low temperature alarm trigger values, and has a malfunction warning and early warning function.

Electrical Schematic Diagram


BD-86L108


BD-86L358


BD-86L458

BD-86L566

BD-86L568/718

BD-86L838

Indicator Panel

Button and Icon Descriptions

Button 1 — up; button 2 — down; button 3 — sound on/off.

Icon 4 — RS485 communication (interface for communication with PC); icon 5 — button lock; icon 6 — start compressor; icon 7 — power on; icon 8 — door open; icon 9 — emergency alarm sound off; icon 10 — power failure; icon 11 — low battery charge (all the above icons are displayed).

Operating Instructions

1. Check the ambient air temperature

In keyboard lock mode, press the multifunctional 'Sound On/Off' button to display the ambient temperature on the segment LED display. Wait 5 seconds or press the 'Up' or 'Down' buttons to return to the initial screen. In keyboard unlock mode, press the multifunctional 'Sound On/Off' button to display the ambient temperature on the segment LED display. Wait 5 seconds to return to the initial screen.

2. Disable sound signal

When the alarm sounds, you can press the multifunction button 'Sound On/Off' to disable the alarm; after this, the freezer will switch to the ambient temperature measurement state. If the alarm sound is disabled, you can press the multifunction button 'Sound On/Off' to enable the sound; after this, the freezer

will transition to the ambient temperature measurement state. While in the keyboard locked mode and in the ambient temperature measurement state, you can press the multifunctional 'sound on/off' button to turn off or turn on the sound of the emergency alarm.

3. Keyboard Lock

In the unlocked keyboard mode, if the user is inactive for 60 seconds or simultaneously presses and holds the 'Up/Down' buttons for 3 seconds, the keyboard will be locked (a double audible signal will sound). In keyboard lock mode, you can view the ambient temperature and turn the emergency alarm sound on/off. All other functional buttons are locked, and pressing them will produce a double audible signal.

To unlock, press and hold the 'up/down' buttons for 3 seconds — the display will show the code '0000'. Use the 'up/down' buttons to change all zeros (0) to the value '5'. If the password is entered correctly, the display will return to the initial screen, and a single audible signal will sound.

4. Setting Parameters

If the keyboard is unlocked, press and hold the 'Sound On/Off' button for 3 seconds — the display will show the code 'PS1'. Press the 'Sound On/Off' button to switch between menus and corresponding settings, and use the 'Up/Down' buttons to select a menu or adjust settings. To save changes, wait 60 seconds or press and hold the 'Sound On/Off' button for 3 seconds — the display will return to the initial screen and a single audible signal will sound.

Each time the freezer is reconnected to the power supply after a power failure, the symbols «1», «2», «3», «4», «5», «6», «7», «8», «9», n, are displayed sequentially on the display, and users can configure the parameters listed in the table below while the symbol «n» is displayed.

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Complete the setup within 2 minutes. If no action is taken, the system will revert to the standard mode.

5. Exporting temperature data to a USB drive

When a USB drive is connected to the USB port, a single beep sounds, and the display shows the message «on» (Data loading in progress). PDF files with data from the current and previous months will be generated and imported to the USB drive, at which point a single audible signal will be emitted and the display will show the message «End» (Completed), and after 6 seconds, it will return to the standard display mode.

Technical specifications*

Basic information	Appearance	Medical freezer BD-86L108	Medical freezer BD-86L458	Medical freezer BD-86L568	Medical freezer BD-86L718	Medical freezer BD-86L838
						
	Capacity (L)	42	24	24	20	18
	Cabinet type	Vertical	Vertical	Vertical	Vertical	Vertical
	Controller	Microprocessor	Microprocessor	Microprocessor	Microprocessor	Microprocessor
	Display	LED	LED	LED	LED	LED
	Temperature adjustment accuracy (°C)	0.1	0.1	0.1	0.1	0.1
	Temperature range	-40 ~ -86°C	-40 ~ -86°C	-40 ~ -86°C	-40 ~ -86°C	-40 ~ -86°C
	Refrigerator/Freezer	Freezer	Freezer	Freezer	Freezer	Freezer
	Cooling Type	Direct Cooling	Direct Cooling	Direct Cooling	Direct Cooling	Direct Cooling
	Effective Volume	108 L	458 L	568 L	718 L	838 L
	Climate Class	N	N	N	N	N
	Battery Type	sealed lead-acid battery, and it complies with the CE certificate, 12V, 4Ah, operating time 24 hours	sealed lead-acid battery, and it complies with the CE certificate, 12V, 4Ah, operating time 24 hours	sealed lead-acid battery, and it complies with the CE certificate, 12V, 4Ah, operating time 24 hours	sealed lead-acid battery, and it complies with the CE certificate, 12V, 4Ah, operating time 24 hours	sealed lead-acid battery, and it complies with the CE certificate, 12V, 4Ah, operating time 24 hours
	Fuse Type	NB1-63C20	NB1-63C20	NB1-63C20	NB1-63C20	NB1-63C20

	Packaging Type	Brown Box	Brown Box	Brown Box	Brown Box	Brown Box
	Dimensions (mm) (wxdxh)	955*675*815	885*855*1980	885*995*1980	1030*995*1980	1175*995*1980
	Internal dimensions (mm) (wxdxh)	440*435*590	595*580*1310	595*720*1310	740*720*1310	885*720*1310
	Box dimensions (mm) (wxdxh)	1040*775*870	950*970*2030	950*1115*2030	1095*1115*2030	1240*1115*2030
	Package dimensions (mm) (wxdxh)	1040*775*995	980*1000*2155	980*1145*2155	1125*1145*2155	1270*1145*2155
	Net weight, kg	113	237	286	305	332
	Gross weight, kg	141	292	342	368	401
Function of the device	Audible alarm	yes	yes	yes	yes	yes
	Visual alarm	yes	yes	yes	yes	yes
	Remote alarm	yes	yes	yes	yes	yes
	High temperature alarm	yes	yes	yes	yes	yes
	Low temperature alarm signal	yes	yes	yes	yes	yes
	Door ajar	yes	yes	yes	yes	yes
	Power failure alarm signal	yes	yes	yes	yes	yes
	High condenser temperature alarm signal	yes	yes	yes	yes	yes
	High ambient temperature alarm signal	yes	yes	yes	yes	yes
	Low battery charge alarm signal	yes	yes	yes	yes	yes
	Communication error alarm signal	yes	yes	yes	yes	yes
	Sensor error alarm	yes	yes	yes	yes	yes
	Sensor type	PT1000	PT1000	PT1000	PT1000	PT1000
	Number of sensors	1+2	1+2	1+2	1+2	1+2

	USB port	yes	yes	yes	yes	yes
	Door heating	yes	yes	yes	yes	yes
	LN2 backup system	Optional	Optional	Optional	Optional	Optional
	CO2 backup system	Optional	Optional	Optional	Optional	Optional
	Protection against electric shock class/type	I/A	I/A	I/A	I/A	I/A
Product performance	Nominal voltage/Frequency (V/Hz)	220-240/50	220-240/50	220-240/50	220-240/50	220-240/50
	Nominal power (W)	480	780	1520	1300	720
	Nominal current (A)	4.9	5.9	12	12	9.5
	Cooling system	Unified cooling system	Unified cooling system	Dual cooling system	Dual cooling system	Dual cooling system
	Energy consumption (kWh/24 hours)	6.1	11.3	10.95	13	14.3
	Noise level, dB (A)	40	45	56	49	53
	Compressor brand	SCOPE	SCOPE	SCOPE	SCOPE	SCOPE
	Number of compressors	1	1	2	2	2
	Compressor type	Compressor SC21 CL 404a (LBM)	Compressor SC21 CL 404a (LBM)	Compressor SC21 CL 404a (LBM)	Compressor SC21 CL 404a (LBM)	Compressor SC21 CL 404a (LBM)
	Refrigerant/Amount (g)	Mixed refrigerant (R601a/100g R290a/10g R1150/25g)	Mixed refrigerant (R601a/100g R290a/10g R1150/25g)	Mixed refrigerant (R601a/70g R290/7g R1150/18g)	Mixed refrigerant (R601a/70g R290/7g R1150/18g)	Mixed refrigerant (R601a/90g R290/9g R1150/23g)
	Time to reduce temperature	5h	4.45h	2.9h	3.6h	3.5 h
	Time to increase temperature	11.2h(-81°C--10°C) 4 hours(-81°C--40°C) 2.5 hours(-81°C--50°C)	16.4h(-81°C--10°C) 6.5h(-81°C--40°C) 4.3h(-81°C--50°C)	16.4h(-81°C--10°C) 6.3h(-81°C--40°C) 4.5h(-81°C--50°C)	18 hours(-81°C--10°C) 7.3h(-81°C--40°C) 5.1h(-81°C--50°C)	18.4h(-81°C--10°C) 7.6h(-81°C--40°C) 5.1h(-81°C--50°C)

	Temperature uniformity (°C)	≤2°C	≤5°C	≤5°C	≤5°C	≤5°C
	Foaming agent	cyclopentane +VIP(Vacuum insulated panels)	cyclopentane +VIP(Vacuum insulated panels)	cyclopentane +VIP(Vacuum insulated panels)	cyclopentane +VIP(Vacuum insulated panels)	cyclopentane +VIP(Vacuum insulated panels)
	Thickness of foaming layer (insulation thickness) (mm)	90	90	90	90	90
Product composition data	Handle	Grip	Grip	Grip	Grip	Grip
	Spatula	Material – plastic, scraper width 65 mm, handle width 25 mm, length 140 mm	Material – plastic, scraper width 65 mm, handle width 25 mm, length 140 mm	Material – plastic, scraper width 65 mm, handle width 25 mm, length 140 mm	Material – plastic, scraper width 65 mm, handle width 25 mm, length 140 mm	Material – plastic, scraper width 65 mm, handle width 25 mm, length 140 mm
	Number of door seals	4	4	4	4	4
	Door seal material	Silicone	Silicone	Silicone	Silicone	Silicone
	Shelf size / mm	140*413*293	140*549*293	140*685*293	140*685*293	140*685*293
	Shelves/ Boxes	1	3	4	3	3
	Number of internal doors	2	4	4	4	4
	Freezer shelves(number of racks for cryostats)	6	16	16	20	24
	Number of cryotubes 2 ml	9000	32000	40000	50000	60000
	Additional lock	yes	yes	yes	yes	yes
	Swivel casters	yes	yes	yes	yes	yes

Basic information	Appearance	Medical freezer BD-86L358	Medical freezer BD-86L566
--------------------------	------------	---------------------------	---------------------------

			
Capacity (L)	24	24	
Cabinet type	Vertical	Vertical	
Controller	Microprocessor	Microprocessor	
Display	LED	LED	
Temperature adjustment accuracy (°C)	0.1	0.1	
Temperature range	-40 ~ -86°C	-40 ~ -86°C	
Refrigerator/Freezer	Freezer	Freezer	
Cooling Type	Direct Cooling	Direct Cooling	
Effective Volume	358 L	566 L	
Climate Class	N	N	
Battery Type	sealed lead-acid battery, and it complies with the CE certificate, 12V, 4Ah, operating time 24 hours	sealed lead-acid battery, and it complies with the CE certificate, 12V, 4Ah, operating time 24 hours	
Fuse Type	NB1-63C20	Fuse Type	
Packaging Type	Brown Box	Brown Box	
Dimensions (mm) (wxdxh)	866x811x1920	885x995x1980	
Internal dimensions (mm) (wxdxh)	580x691x1280	595x710x1310	

	Box dimensions (mm) (wxdxh)	930x925x1970	950x1115x2030
	Package dimensions (mm) (wxdxh)	960x955x2000	980x1145x2155
	Net weight, kg	185	280
	Gross weight, kg	240	340
Function of the device	Audible alarm	yes	yes
	Visual alarm	yes	yes
	Remote alarm	yes	yes
	High temperature alarm	yes	yes
	Low temperature alarm signal	yes	yes
	Door ajar	yes	yes
	Power failure alarm signal	yes	yes
	High condenser temperature alarm signal	yes	yes
	High ambient temperature alarm signal	yes	yes
	Low battery charge alarm signal	yes	yes
	Communication error alarm signal	yes	yes
	Sensor error alarm	yes	yes
	Sensor type	PT1000	PT1000
	Number of sensors	1+2	1+2
	USB port	yes	yes
	Door heating	yes	yes
	LN2 backup system	Optional	Optional
	CO2 backup system	Optional	Optional
	Protection against electric shock	I/A	I/A

	class/type		
Product performance	Nominal voltage/Frequency (V/Hz)	220-240/50	220-240/50
	Nominal power (W)	710	1520
	Nominal current (A)	5.5	12
	Cooling system	Unified cooling system	Unified cooling system
	Energy consumption (kWh/24 hours)	10.3	10.95
	Noise level, dB (A)	40	49
	Compressor brand	SCOPE	SCOPE
	Number of compressors	1	2
	Compressor type	Compressor SC21 CL 404a (LBM)	Compressor SC21 CL 404a (LBM)
	Refrigerant/Amount (g)	Mixed refrigerant (R601a/60g R290a/6g R1150/15g)	Mixed refrigerant (R601a/70g R290a/10g R1150/25g)
	Time to reduce temperature	4.5 h	2.9h
	Time to increase temperature	11.2h(-81°C--10°C) 4 hours(-81°C--40°C) 2.5 hours(-81°C--50°C)	16.4h(-81°C--10°C) 6.3h(-81°C--40°C) 4.5h(-81°C--50°C)
	Temperature uniformity (°C)	≤2°C	≤5°C
	Foaming agent	cyclopentane +VIP(Vacuum insulated panels)	cyclopentane +VIP(Vacuum insulated panels)
	Thickness of foaming layer (insulation thickness) (mm)	90	90
Product composition data	Handle	Grip	Grip

Spatula	Material – plastic, scraper width 65 mm, handle width 25 mm, length 140 mm	Material – plastic, scraper width 65 mm, handle width 25 mm, length 140 mm
Number of door seals	4	4
Door seal material	Silicone	Silicone
Shelf size / mm	140*330*293	140*685*293
Shelves/ Boxes	1	3
Number of internal doors	2	4
Freezer shelves(number of racks for cryostats)	4	16
Number of cryotubes 2 ml	7000	40000
Additional lock	yes	yes
Swivel wheels	yes	yes

***parameter table deviation $\pm 10\%$**

Maximum permissible time to reach operating mode

BD-86L108	5 h
BD-86L458	4.45 h
BD-86L568	2.9 h
BD-86L718	3.6 h
BD-86L838	3.5 h
BD-86L358	4 h
BD-86L566	2.9 h

Software version:

Software model: EMC-3100-05M

Software version and date: V2.0, 2.2, 08.2018

Software manufacturer: Elitech Technology, Inc., No. 1, Huangshan Road, Tongshan Economic Development Zone, Xuzhou, Jiangsu Province

Class A

The product must comply with GOST R IEC 62366-1. In terms of electromagnetic compatibility, the product must comply with GOST R IEC 61326-1.

1.15. General description of accessories, composition, functional capabilities


1. Stainless steel shelves – used for additional placement of biological materials.

Dimensions (DxWxH), mm $\pm 10\%$

Weight, g $\pm 10\%$

2. Cryoboxes 5x5; Cryoboxes 9x9; Cryoboxes 10x10

Cryoboxes - used for storing 2ml cryovials.


Cryobox example 9x9

Name	Dimensions (DxWxH), mm $\pm 10\%$	Weight, g $\pm 10\%$
------	-----------------------------------	----------------------

Cryoboxes 5x5	75x75x52	65
Cryoboxes 9x9	133x133x53	100
Cryoboxes 10x10	133x133x53	100

3. Racks for cryoboxes 5x5; Racks for cryoboxes 9x9; Racks for cryoboxes 10x10
Racks for cryoboxes - used for storing cryoboxes.

Name	Contents	Dimensions (DxWxH), mm ±10%	Weight, kg ±10%
Racks for cryoboxes 5x5	Contains 25 cryoboxes (cryobox for 2 ml cryovial)	140x675x285	6
Racks for cryoboxes 5x4	Contains 20 cryoboxes cell size (cryobox for 2 ml cryovial)	140x554x285	5

1.16. Material Description

No direct contact with the product, indirect contact through gloves (registered in the established order in the Russian Federation)

Part of the product	Material	Manufacturer, country	Contact
Medical freezer	External material- PSM Internal material - stainless steel	China	No human contact
Dool seal	Silicone	China	No human contact
Compressor	Stainless steel	China	No human contact
Condensate fan	Stainless steel	China	No human contact
Refrigerant	R601a R290a R1150	China	No human contact
Foaming agent	Cyclopentane +VIP	China	No human contact
Swivel wheel	Stainless steel with rubber coating	China	No human contact
Shelf	Stainless steel	China	No human contact
Cryoboxes	PP white dye, blue,	China	No human contact

	yellow, red		
Racks for cryoboxes	Stainless steel, PP blue dye, yellow, red	China	No human contact
Spatula	Plastic	China	No human contact
Packaging	Cardboard box	China	No human contact


2. Technical description of medical device

2.1. Conditions for safe use

Do not install the freezer in a narrow and confined space. The dimensions of the door opening of the room where the freezer is to be installed must allow for its unobstructed movement, or at least provide easy access to the freezer to facilitate maintenance in the event of a malfunction and to prevent damage to stored materials due to untimely repair of the freezer.

The floor surface on which the freezer is installed must be strong and level, non-flammable and have the appropriate load-bearing capacity. The freezer must be installed in a well-ventilated area, protected from direct sunlight.

The freezer must be connected to a separate power outlet. Ensure that the power cord plug is securely inserted into the outlet, and the rated current of the outlet is >16 A. Do not allow the power cord to be twisted or compressed. Before connecting, check the operating voltage. When connecting to a network with unstable voltage, use a voltage stabilizer with the appropriate load, the power of the voltage stabilizer should be more than 4000 W to ensure compliance with input

voltage requirements.


The freezer must be properly grounded. Before using the freezer, ensure that the outlet is properly grounded. If the outlet is not grounded, contact a qualified electrician to perform proper grounding of the outlet. Do not ground the freezer to gas or water pipes, telephone lines, or lightning rods, as this can lead to electric shock. The outlet into which the freezer is plugged should be located in an easily accessible place for quick disconnection of the power cord plug in case of emergencies. Do not cover the ventilation openings of the low-temperature freezer.

Since the ambient temperature has a significant impact on the operation of the freezer, failure to comply with the above installation site conditions can result in improper operation of the freezer. Before using the freezer, ensure that the installation site conditions meet the required specifications.

Preparation for Operation

1. Remove all packaging materials from the freezer.
2. Check for all components. Compare the contents of the package with the packing list. In case of discrepancies, immediately contact the authorized representative of the manufacturer.
3. Install the low-temperature freezer so that the distance from the walls of the room to the freezer is at least 30 cm to ensure proper ventilation and heat dissipation.
4. Turn the horizontal support legs clockwise using a wrench or by hand to extend them downward so that they rest on the floor, to prevent the low-temperature freezer from moving during operation.
5. Open the outer door so that you can remove the shelves from the bottom of the freezer compartment, remove the protective film, and install the shelves one by one on the support strip.

Do not leave plastic packaging bags in places accessible to children.


First startup

Upon the first activation of the freezer, perform the following steps:

1. After installing the freezer in the desired location, leveling, and cleaning, before the first power-up, leave the freezer in the installation location for at least 24 hours. This is necessary to ensure its proper operation in the future;
2. Connect the power cord to a separate outlet with the appropriate rating under no-load conditions.
3. After turning on the power supply, first turn on the freezer's power switch; if the freezer is equipped with a battery, re-engage the battery switch (located on the back side at the bottom of the freezer compartment).
4. After powering on, set the time.
5. Turn on the power. Do not place anything in the low-temperature freezer, start the operation and check if the set temperature is reached in the working chamber. Leave the freezer running for 24 hours and ensure it is functioning properly.
6. To avoid excessive temperature rise in the working chamber, do not place materials exceeding 1/3 of the working volume of the chamber. Ensure the freezer operates normally after power off and power on. The next batch of materials can be placed inside no sooner than 8 hours after the previous loading.
7. Avoid opening the freezer door during the cooling process, as this will lead to an increase in temperature.

Actions after a power failure

The set value is stored in the freezer's memory. In the event of a power failure, after power is restored, the freezer will resume operation according to the settings established before the power outage.

If the low-temperature freezer is disconnected from the power supply or a power failure occurs, wait at least 5 minutes before restarting to avoid damage to the compressor or the entire system.

If the freezer is not used for a long period of time, unplug the power cord from the outlet and turn off the battery switch (for freezer models with a battery) to prevent electric shock, current leakage, or fire due to power cord wear. If you plan to store the low-temperature freezer for a long period of time, take reasonable measures to restrict access to it by unauthorized persons, including children.

Operating rules of the product

This product is a low-temperature freezer. Make sure in advance that the temperature required for storing your materials corresponds to the temperature range of the freezer; otherwise, it may lead to damage of the stored materials due to a mismatch in the temperature regime.

Due to the thermal inertia of the sensors, the temperature and humidity displayed on the screen may differ from the actual temperature and humidity inside the chamber, which is a normal phenomenon. Since all low-temperature freezers are storage equipment, it is strictly prohibited to place several relatively hot materials simultaneously. This will cause the compressor to operate continuously for a long time, reducing its lifespan, and the temperature in the working chamber will decrease slowly. Materials should be loaded in batches to gradually reduce the temperature to the required level!

Exercise caution to avoid damaging the cooling circuit.

Do not use electrical appliances inside the freezer without the manufacturer's permission.

- Do not use mechanical tools or other means to accelerate defrosting without the manufacturer's permission.

Data transfer via USB

The USB drive is connected to the freezer, and temperature data from the device (freezer) is transferred to it. The drive is then removed from the freezer and connected to a PC. The data is provided in an Excel file.

Shutting down the device

Shutting down: If you plan to store the low-temperature freezer for a long period, take reasonable measures to restrict access by unauthorized persons. Also, block the door in such a way that it cannot be fully closed.

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2.2. Maintenance


Operation of the freezer must be carried out by an authorized person with mandatory daily record-keeping (it is recommended to conduct and record inspection results every 2-4 hours). In the event of a power failure, the temperature in the freezer's working chamber will begin to rise. If the malfunction cannot be resolved in a short period, stored materials should be transferred to another suitable environment to prevent their deterioration.

This product is a low-temperature freezer. Make sure in advance that the temperature required for storing your materials corresponds to the temperature range of the freezer; otherwise, it may lead to damage of the stored materials due to a mismatch in the temperature regime.

Due to the thermal inertia of the sensors, the temperature and humidity displayed on the screen may differ from the actual temperature and humidity inside the chamber, which is a normal phenomenon. Since all low-temperature freezers are storage equipment, it is strictly prohibited to place several relatively hot materials simultaneously. This will lead to continuous operation of the compressor for a long time and

a reduction in its service life, while the temperature in the working chamber will decrease slowly. Materials should be loaded in batches to gradually reduce the temperature to the required level!

- Exercise caution to avoid damaging the cooling circuit.
- It is prohibited to use electrical appliances inside the freezer without the manufacturer's permission.
- Do not use mechanical tools or other means to accelerate defrosting without the manufacturer's permission.

To avoid electric shock or injury, it is necessary to disconnect the power supply before any repair and maintenance of the freezer.

Ensure that medicines or vials with suspended particles inside or outside the freezer are tightly sealed to avoid inhaling them during repair or maintenance. Otherwise, it may cause harm to health.

Do not spray water on the freezer, as this can lead to a decrease in the insulation properties of electrical components and the appearance of rust on metal parts.

Do not use hot water and aggressive cleaning agents or organic solvents to clean the internal surfaces of the chamber. Do not place heavy objects on the freezer, as this can lead to deformation of the external panels of the freezer under pressure.

Defrosting of internal walls

A thick layer of ice can lead to a gap between the medical freezer and the door seal, resulting in a decrease in cooling capacity. Recommended defrosting period: once a month.

To clean the inner door, use the scraper attached to the freezer.

To defrost the freezer naturally, perform the following actions.

1. If auxiliary cooling devices are present, turn them off;
2. Remove all stored materials from the freezer and move them to a suitable storage environment;
3. Turn off the power supply.
4. Open the outer and inner doors to allow the freezer to defrost naturally;
5. Place a dry towel at the bottom of the chamber to collect moisture;
6. After cleaning the low-temperature freezer and the inner door, turn the power back on;
7. Wait until the freezer is fully cooled, and then return the products to the working chamber;
8. If auxiliary cooling devices are present, turn them on.

Troubleshooting

In the event of malfunctions, refer to the information provided in the following table. If the malfunction cannot be resolved, contact an authorized manufacturer's representative.

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Note: If the problem persists after all checks, or if the malfunction is not listed in the table above, contact an authorized manufacturer's representative.

Joint Stock Company 'Labtech' (JSC 'Labtech')

RF, 105264, Russia, Moscow, Vn. ter. g. Municipal District Izmailovo, Izmailovsky Blvd, bld. 1/28, unit 1A/1

+7(495) 276-77-00

order@labteh.com

2.3. Cleaning and Disinfection

Clean the low-temperature freezer no less than once a month.

Wipe the freezer body, working chamber, and all accessories with a dry cloth to remove dust. In case of contamination of the freezer, remove dirt with a cloth moistened in a neutral cleaning agent, then remove any residue with a damp cloth and dry thoroughly. Then wipe the surfaces again with a dry cloth.

Do not allow water to enter the casing or the interior of the low-temperature freezer. Otherwise, this can lead to damage to the electrical insulation and, consequently, a fire.

The compressor and other mechanical parts are completely sealed and do not require lubrication.

Clean the inner wall and the condenser's mesh filter of frost or ice at least once a month.

Clean the mesh filter at least once a month. If the mesh filter is clogged, it will lead to a reduction in cooling capacity and the lifespan of the freezer.

To clean the mesh filter, do the following: Do not move the sensor wire.

1. Remove the front lower cover and the mesh filter.
2. Rinse the mesh filter with water.
3. Blow the condenser with compressed air to remove dust;
4. Wait for the mesh filter to dry completely and install it back in place along with the cover.

Sterilization

The medical device is not sterile.

2.4. Environmental protection requirements, including requirements for safe disposal and recycling. Environmental protection.

The products do not release toxic substances into the environment during storage, transportation, and operation, and do not have (upon direct contact) a harmful effect on the human body and the environment. No special environmental protection measures are required when handling the product. Disposal is carried out in accordance with SanPiN 2.1.3684-21 for Class A waste. The electronic and electrical components of the product (excluding batteries, control panel, etc.) are disposed of in accordance with the rules for the disposal of electronic components. Batteries, power supplies, and similar power sources are removed from electrical/electronic components and disposed of in accordance with current local regulations.

Unused items, those past their expiration date, and those that have deteriorated must be disposed of in accordance with SanPiN 2.1.3684-21, Class A – epidemiologically safe waste. Electronic and electrical components of the item (excluding batteries, control panel, etc.) should be disposed of according to the rules for the disposal of electronic components.

For additional information on disposal or recycling, contact local authorities.

Waste and material disposal can be carried out on a contractual basis with a company that has the appropriate license.

2.5. Service life data

Average service life - 10 years.

2.6. Transportation and storage requirements

The product can be transported by all types of covered transport, except unheated compartments of aircraft, and in accordance with the rules approved in the established manner for each type of transport.

Transportation and storage conditions: Ambient temperature $-40^{\circ}\text{C} \sim +55^{\circ}\text{C}$, relative humidity 10~90%.

2.7. Manufacturer's Warranties

From the date of purchase, manufacturing defects of the product, except for unauthorized disassembly and repair during the 'Three Guarantees' period, can be rectified based on the warranty card and sales invoice after evaluation by the customer service department or repair shop in accordance with the Regulations on Responsibility for Replacement and Return of Certain Goods (GJM) ('New Three Guarantees').

The warranty does not cover damages caused by improper transportation, operation, maintenance, storage, lack of a warranty card and valid sales invoice, mismatch between the model on the warranty card and the repaired item or any alteration thereof, as well as damages caused by force majeure.

Warranty Period Definition

This product is covered by a free warranty for the entire warranty period from the date of purchase.

Conditions for Providing a Free Warranty

1. The user purchased a medical freezer and completed the warranty registration procedure.

2. The user provided a valid warranty card (Warranty Card), with the name, model, and barcode of the freezer to be repaired matching those specified in the warranty card. If a valid warranty card cannot be provided, the warranty is extended by 6 months based on the warranty period calculated from the delivery date of the freezer as indicated by the barcode.

3. The freezer must be under warranty.

4. Presence of malfunctions not caused by human factors. Date of manufacture: see. barcode of the freezer.

Warranty on the product.

Complete product – 1 year, Main components – 3 years (Main components include the compressor, main control panel, display control panel, temperature regulator, evaporator, condenser, filter, capillary tube, fan motor).

Warranty on the product:

Warranty period of operation and warranty storage period: 10 years

3. Marking and packaging of the product

3.1 The product must be supplied with a label (located on the back of the product) containing the following data in accordance with GOST IEC 61010-1:


- country of manufacture;
- name and (or) trade mark of the manufacturer (if any), address;
- name of the product;
- date of manufacture (month and year) and serial number (indicated by text or the appropriate symbol);
- electrical characteristics (type of power supply current: 1) alternating current: nominal network frequency or frequency range; nominal voltage or range of nominal voltage values of the power supply (V), nominal current (A));
- maximum nominal power in watts;
- refrigerant/quantity (g);
- effective volume;
- symbols of electrical safety classification (Class I);
- pressure sensor (high/low MPa);
- net weight;
- registration certificate number (upon receipt);
- conformity mark (if required);
- authorized representative of the manufacturer;
- product barcode (if required).

Note – It is permissible to apply to the labeling graphic symbols and diagrams explaining the conditions of storage of items, their disposal, etc., according to GOST R ISO 15223-1 and GOST IEC 61010-1.

Product Labeling

Next to the fuse holder, the type of fuse is indicated, allowing the operator to correctly identify the replacement fuse.

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Midea Biomedical Co., Ltd.
No.5 workshop 88 Yulan Main Road, Hi-Tech
Development Zone Hefei, Anhui, China

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Labeling Layout for the Russian Federation


3.2 Transport Labeling Contains:

- country of manufacture;
- name and (or) trade mark of the manufacturer (if any), address;
- name of the product;
- manufacturing date (month and year) and serial number.

Transport labeling with the handling symbol “Protect from Moisture,” “Fragile. Handle with Care,” “Do Not Expose to Sunlight.” It is permissible to affix handling symbols printed by a printing method.

3.2. Description of symbols marked on the MI

	Caution! Refer to the accompanying documentation!
	Manufacturer
	Date of manufacture

	Keep away from moisture
	Fragile, handle with care
	Do not expose to sunlight
	Flammability hazard
	Fusible link
	Do not dispose of! Requires special disposal!
	CE Marking

3.3. Packaging process of the product

The product is packaged in a cardboard box of the following dimensions:

Product variant	BD-86L108	BD-86L458	BD-86L568	BD-86L718	BD-86L838	BD-86L358	BD-86L566
Box dimensions (mm) (wxdxh) ±10%	1040*775*870	950*970*2030	950*1115*2030	1095*1115*2030	1240*1115*2030	960*955*2000	980*1145*2155
Mass, kg ±10%	141	292	342	368	401	240	340

4. Table of Correspondence of Essential Principles

List of Applied National and International Standards for the Purpose of Ensuring Safety, Effectiveness, and Quality of the Medical Device

GOST ISO 13485-2017 (ISO 13485:2016) Quality Management System/ Medical Devices. Quality Management Systems. Requirements for Regulatory Purposes

GOST ISO 14971-2021 (ISO 14971:2019) Application of risk management to medical devices / Medical devices. Application of risk management to medical devices

GOST R ISO 15223-1-2023 (ISO 15223-1:2021) Medical devices. Symbols to be used with medical device labels / Medical devices. Symbols used to convey information provided by the manufacturer. Part 1. Basic requirements

GOST IEC 61010-1-2014 (IEC 61010-1:2010+A1:2016) Safety of electrical instrumentation and laboratory equipment. Part 1. General requirements / Safety of electrical measuring and laboratory equipment. Part 1. General requirements

GOST IEC 62304-2022 (IEC 62304:2006+A1:2015) Medical device software - Software life-cycle processes. / Medical devices. Software. Life cycle processes

GOST R IEC 62366-1-2023 (IEC 62366:2015+A1:2020) Medical devices - Application of usability engineering to medical devices / Medical devices. Software. Life cycle processes

GOST R IEC 61326-1-2014 (IEC 61326-1:2020) Electrical equipment for measurement, control, and laboratory use. EMC requirements. Part 1. General requirements/ Electrical equipment for measurement, control, and laboratory use. Electromagnetic compatibility requirements. Part 1. General requirements

5. Authorized representative of the manufacturer in the Russian Federation (name, address, phone, email)

Legal form and full name of the legal entity with transliteration, or surname, name, and patronymic (if applicable), and the details of the document verifying the identity of the individual entrepreneur	Joint Stock Company «Labtech»
Short name of the legal entity (if applicable)	JSC «Labtech»
Trade name of the legal entity (if applicable)	-
Taxpayer identification number	9719053826
Legal address of the legal entity	RF, 105264, Russia, Moscow, Vn. ter. g. Municipal District Izmailovo, Izmailovsky Blvd, bld. 1/28, unit 1A/1
Phone numbers	+7(495) 276-77-00
Email address of the legal entity, individual entrepreneur (if applicable)	order@labteh.com

Appendix 1.

Freezers are intended for use in a basic electromagnetic environment. Freezers are intended for use in the electromagnetic environment defined according to GOST R IEC 61326-1-2014 and presented in the table below. Users of the products must ensure their use in the specified electromagnetic environment.

Freezers are not intended for use in conjunction with other medical devices and general-purpose devices.

When connecting equipment not included in the supply package, electromagnetic emissions may occur that exceed the level specified in the GOST R IEC 61326-1-2014 standard.

Manufacturer's recommendations and statement - electromagnetic radiation		
The medical freezer is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in such an environment.		
Radiation testing	Compliance	Electromagnetic environment - recommendations
Radio frequency radiation CISPR 11	Group 1	Radio frequency energy is used in the medical freezer only for its internal functions. Therefore, its radio frequency emission is very low and is unlikely to cause any interference with nearby electronic equipment.
Radio frequency radiation CISPR 11	Class A	The medical freezer is suitable for use in all institutions except for home use, and can be used in residential premises and those directly connected to the public low-voltage power supply network supplying buildings used for residential purposes, provided the following warning is observed: Warning: This equipment/item is intended for use only by healthcare professionals. This equipment/item may cause radio interference or affect the operation of nearby equipment. Measures may need to be taken to mitigate the effects, such as reorientation or relocation of the medical freezer or shielding of the installation site.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies with	

Manufacturer's Recommendations and Statement - Electromagnetic Interference Immunity			
The medical freezer is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in such an environment.			
Immunity Testing	Control level according to IEC 61326-1	Compliance Level	Electromagnetic Environment - recommendations
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	The floor should be wooden, concrete, or ceramic. If the floor is covered with synthetic material, the relative humidity should be at least 30%.
Fast electrical transients/surge immunity according to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input/output lines	± 2 kV for power lines ± 2 kV for input/output lines	The quality of the power supply must comply with the standard for commercial or hospital conditions.
Surge IEC 61000-4-5	± 1 kV in differential	± 1 kV in differential	The quality of the power supply

	mode ± 2 kV in common mode	mode ± 2 kV in common mode	must comply with the standard for commercial or hospital conditions.
Voltage dips, short interruptions, and voltage variations on power supply input lines according to IEC 61000-4-11	$<5\% U_T$ (drop $> 95\% U_T$) for 0.5 cycles $40\% U_T$ (drop $60\% U_T$) for 5 cycles $70\% U_T$ (drop $30\% U_T$) for 25 cycles $<5\% U_T$ (drop $> 95\% U_T$) for 5 sec	Drop $> 95\%$ in 10 ms Drop 60% in 100 ms Drop 30% in 500 ms Drop 95% in 5000 ms	The quality of the power supply must meet the standard for commercial or hospital conditions. If continuous operation of the device is required during power outages, it is recommended to connect the device to an uninterruptible power supply.
Magnetic field of industrial frequency (50/60 Hz) according to IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields of industrial frequency should be at levels characteristic of a standard location in a standard commercial or hospital environment.
NOTE. - U_T - is the AC supply voltage before application of the test level.			

Manufacturer's Recommendations and Statement - Electromagnetic Interference Immunity			
The medical freezer is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in such an environment.			
EMI TESTING	Control level According to IEC 61326-1	Level compliance	Electromagnetic environment - recommendations
Radio frequency testing was performed according to IEC 61000-4-6	3 V RMS From 150 kHz to 80 MHz	[3] V	Portable and mobile radio frequency telecommunications equipment should not be used closer to any part of the product, including cables, than the recommended distance, calculated using the equation applicable to the transmitter frequency. Recommended safe distance From 230 MHz to 1.0 GHz, where P - maximum transmitter output power in watts (W) according to the manufacturer's data, and d - recommended distance in meters (m). Field strength from stationary radio transmitters, determined by electromagnetic site survey, ^a should be below the compliance level in each frequency band. ^B Interference may occur near equipment marked with the following symbol:
Radio frequency radiation level according to IEC 61000-4-3	3 V/m From 80 MHz to 1.0 GHz	[3] V/m	

			
NOTE 1. At 30 MHz and 230 MHz, a higher frequency range is used.			
NOTE 2: These recommendations may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects, and people.			
^a The field strength from fixed transmitters, such as base stations for radiotelephones (cellular/wireless) and terrestrial mobile radio stations, amateur radio stations, AM and FM broadcasting, and television, cannot be predicted with high accuracy theoretically. To assess the electromagnetic environment created by stationary radio transmitters, consider conducting an electromagnetic survey of the site. If the measured field strength in the location where the product is used exceeds the applicable compliance level for radio frequencies specified above, it is necessary to check the [MEDICAL EQUIPMENT or MEDICAL DEVICE] to ensure proper operation. If operational malfunctions are observed, additional measures may be required, such as reorientation or relocation of the device.			
^b In the frequency range from 80 MHz to 230 MHz, the field strength should be less than [3] V/m.			

Recommended safe distances between portable and mobile radio frequency communication equipment and Medical freezer			
The Medical freezer is designed for use in an electromagnetic environment where radiated radio frequency interference is controlled. The buyer or user can help prevent electromagnetic interference by maintaining the minimum distance between portable and mobile radio frequency communication equipment (transmitters) and the medical freezer, as recommended below, in accordance with the maximum output power of the communication equipment.			
Nominal maximum output power transmitter [W]	Safe distance depending on transmitter frequency [m]		
	From 150 kHz to 30 MHz	From 30 MHz to 230 MHz	From 230 MHz to 1.0 GHz
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
For transmitters designed for maximum output power not specified above, the recommended safe distance d in meters (m) can be calculated using the equation applicable to the transmitter frequency, where P is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer's specifications.			
NOTE 1: At 30 MHz and 230 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These recommendations may not apply in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects, and people.			

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